DENTSPLY
TULSA DENTAL
SPECIALTIES

DENTSPLY Tulsa Dental Specialties
608 Rolling Hills Drive
Johnson City, TN 37604
(423) 282-1453
(800) 669-0822
Fax: (423) 262-1000

Directions for Use

For dental use only

GuttaCore gutta-percha endodontic obturators on crosslinked gutta-percha carrier

Manufactured by:

**DENTSPLY Tulsa Dental Specialties, 608 Rolling Hills Drive, Johnson City,
Tennessee 37604, USA**

Contents

1. Description	2
2. Delivery forms	2
3. Photos.....	3
4. Indication for use.....	3
5. Contraindications	3
6. Warnings	4
7. Adverse reactions	4
8. Precautions.....	4
9. Step-by-step instructions	4
10. Shelf life and storage	8
11. Transportation.....	8
12. Environmental requirements	8
13. Conditions of use.....	8
14. Disposal	8
15. Manufacturer's warranties	9
16. Authorized Representative.....	9

Signature: Terri Stanton

Terri Stanton

Site Quality Manager

DENTSPLY Tulsa Dental Specialties

Date: April 10, 2019

Signature: Anthony Skeans

Anthony Skeans – Notary

Date: 4/10/19

My Commission Expires
9-25-2022

1. Description

Name: GuttaCore gutta-percha endodontic obturators on crosslinked gutta-percha carrier.

Description: GuttaCore obturators are gutta-percha endodontic obturators on crosslinked gutta-percha carrier, non-sterile, for single use. The distinctive feature of latest-generation GuttaCore obturators is that they have a carrier of the same natural pink color as the obturator itself.

Crosslinking is a well-known chemical process of linking polymer chains together. This technique allows to achieve the required strength for a thin gutta-percha carrier. Therefore, it became possible to create an obturator with a carrier which is gutta-percha, inside and out. It is important to note that crosslinked gutta-percha is not sticky, does not melt and will not dissolve with solvents.

The main advantage of gutta-percha endodontic obturators on crosslinked gutta-percha carrier is that they move warm gutta-percha 3-dimensionally into the root canal system. Many obturation techniques use lateral or vertical compaction, however the hydraulic force from these techniques generally moves gutta-percha in one or two unequal directions (laterally or apically). When using GuttaCore, due to the presence of a crosslinked gutta-percha carrier, the hydraulic force compacts warm gutta-percha flowing equally in three dimensions.

Information about medical device manufacturer/designer:

DENTSPLY Tulsa Dental Specialties

608 Rolling Hills Drive, Johnson City, Tennessee 37604, USA

Tel.: 1-800-662-1202

Addresses of medical device manufacturing sites:

DENTSPLY Tulsa Dental Specialties

608 Rolling Hills Drive, Johnson City, Tennessee 37604, USA

Tel.: 1-800-662-1202


2. Delivery forms

1. GuttaCore gutta-percha endodontic obturators on crosslinked gutta-percha carrier, sizes 20, 25, 30, 35, 40, 45, 50, 55, 60, 70, 80, 90; 6 pieces in the blister.
2. GuttaCore gutta-percha endodontic obturators on crosslinked gutta-percha carrier for ProTaper NEXT instruments; sizes X2, X3, X4, X5; 6 pieces in the blister.

3. Photos

GuttaCore gutta-percha endodontic obturators on crosslinked gutta-percha carrier

GuttaCore gutta-percha endodontic obturators on crosslinked gutta-percha carrier for Protaper Next instruments

4. Indication for use

GuttaCore gutta-percha endodontic obturators on crosslinked gutta-percha carrier are used for filling root canal systems.

Obturators should be used only in clinical environments by dental practitioners.

5. Contraindications

Not known

6. Warnings

This product contains dry natural rubber.

7. Adverse reactions

Patients with latex sensitivity may experience allergic reactions to gutta-percha, which contains dry natural rubber.

8. Precautions

- Select the GuttaCore obturator according to the size verifier that fits the shaped canal passively.
- GuttaCore obturators should not be used as size verifiers.
- Do not remove gutta-percha from the apical third of the obturator prior to filling the root canal space. Removing the gutta-percha in the apical third damages the obturator.

9. Step-by-step instructions

Prior to using any new product, read all the manufacturer's instructions. Make sure that you are competent when using the product and technique prior to clinical use.

General recommendations:

1. Create straight-line access to all root canal systems for application of both rotary and hand files. Remember, the use of a rubber dam is the standard of care for endodontic procedures.
2. Use only sterile files when entering the canal.
3. Determine the ideal working length position with an apex locator, radiograph or both.
4. Use all files with minimal pressure especially in the apical third and when resistance is felt. Do not force any root canal instruments when negotiating or shaping the canal. It is recommended to use a crown-down technique with rotary instruments.
5. Shape and clean canals to sizes compatible with the root shape, size and curvature. It is necessary to avoid over-enlargement of root canal systems to preserve sound dentin and to prevent possible instrument damage in the canal space, e.g., ledging, creation of false canals, separation of instruments.
6. Flare the canal coronally and ensure the GuttaCore obturator has sufficient space to enter the root canal along the straight-line access by using an orifice shaping instrument as necessary such as EndoAccess Kit, Gates Glidden Bur, GT Accessory File, ProFile, O.S. Orifice Shaper, or a Protaper SX or S1 file.
7. Proceed slowly with all instruments in areas of canal irregularity and curvature.
8. Consider all root canal instruments for "single use/single case" applications, as wear and fatigue may impact the integrity of the instrument and result in separation with excessive use.

9. Verify the canal size prior to obturating.
10. Irrigate frequently throughout the cleaning and shaping procedure.

Shaping

- Proper shaping permits the GuttaCore obturator to move passively through the canal. Shaping begins by creating straight-line access and shaping the canal – crown down.
- Shape canal – crown down.

Description of crown-down technique. The technique includes crown-down canal shaping with changing instruments sequentially from larger to smaller size. Consult the file system's directions for use to determine the starting file in small canals (medial and buccal canals of molars, small premolars and mandibular anteriors). Perform root canal processing using the file to resistance or working length (whichever occurs first). If resistance is encountered before working length is obtained, go to next smaller instrument following the same protocol until working length is achieved. Between each rotary file use a #10 or #15 tip hand file to maintain the glide path to the canal terminus. Frequent irrigation with sodium hypochlorite (NaOCl) is used throughout the shaping procedure.

Consult the file system's directions for use to determine the file size for larger canals (palatal and distal canals of molars, larger premolars, maxillary anteriors). Perform root canal processing using the file to resistance or working length (whichever occurs first). If resistance is encountered before working length is obtained, go to next smaller instrument following the same protocol until working length is achieved. Between each rotary file use hand instrument to maintain a glide path to working length.

Follow the file system's directions for use when creating access, shaping the canal and flaring the coronal portion of the canal.

A properly shaped canal preserves the natural shape of the anatomy and allows for a more effective volume of disinfecting irrigants to penetrate, circulate and clean all areas of the root canal system.

Verify the shape at working length and passive fit.

Following proper cleaning and shaping, confirm the final working length measurement using a size verifier or a small hand file that reaches the apical constriction passively, with no significant resistance or twisting. It is necessary to use a size verifier or a small hand file to confirm the shape at working length and passive fit. Verifying the shape at working length and passive fit ensures that the obturator being used matches the correct apical canal diameter.

- a. Remove the size verifier from its package to confirm a passive fit at working length; set the silicone stopper at your working length using the size verifier's calibration mark. Check passivity by taking the size verifier to working length and confirm by rotating in the canal 180°.

- b. If the fit is passive, remove the obturator from the package and proceed to the next step.
- c. If the fit is not passive, use the size verifier as a finishing file to gently enlarge the canal terminus or select a smaller size verifier. If the smaller size verifier fits passively, remove the obturator of the same size as the verifier and proceed to the next step.

Irrigate the canal

The cleaning of the shaped root canal system is accomplished primarily through the use of copious irrigation with NaOCl 5.25%. The activated irrigants improve debridement and disrupt the smear layer and biofilm inside the canal, promoting deep cleaning and disinfection of lateral canals, fins, webs, isthmuses, anastomoses, and dentinal tubules. Dry the canal thoroughly using absorbent points.

Selecting the obturator and setting working length

Select a GuttaCore obturator from the same package as the size verifier used to verify the shape of the working length and passive fit. Available GuttaCore obturator sizes include:

Verifier size/taper	GuttaCore obturator size
20/0,04	-
25/0,04	20
30/0,04	25/X2
35/0,04	30/X3
40/0,04	35
45/0,04	40/X4
50/0,04	45
55/0,04	50/X5
60/0,04	55
70/0,04	60
80/0,04	70
90/0,04	80

Verifier size/taper	GuttaCore obturator size
20/0,06	20
25/0,06	25/X2
30/0,06	30/X3
35/0,06	35
40/0,06	40/X4
45/0,06	45
50/0,06	50/X5
55/0,06	55
60/0,06	60
70/0,06	70
80/0,06	80
90/0,06	90

When the final shape is a 0.04 taper, select an obturator one size smaller than the last file taken to working length. When the final shape is a 0.06 taper or larger, select an obturator of the same size as the last file taken to working length.

Set the top of the stopper on the obturator to coincide with the working length established. Fix the working length with a stop ring following calibration marks on the obturator shaft. Calibration marks are set at the following working lengths (in millimeters): 18, 19, 20, 22, 24, 27, 29. The last two are on the obturator handle.

Asepsis

Disinfect the obturator in a 5.25% sodium hypochlorite solution (bleach) for one minute. Gently wipe the gutta-percha on the obturator with cotton moistened with 70% alcohol.

Drying the canal and applying sealer

Use a sterile paper point of the same tip and taper as the last file taken to working length. Brush a very light coating of sealer with the paper point. To help ensure the obturator reaches working length, use an additional absorbent point to remove any excess sealer that may have accumulated on the canal walls or pooled in the apical third. If obturating more than one canal in a single tooth, place sealer in all the canals at the same time. This facilitates removal of excess gutta-percha by preventing the gutta-percha from sticking to other orifices or dentin.

Heating obturators in the obturator oven

GuttaCore is heated to 160°C then maintain 140°C until oven turning off. The obturator cools to room temperature (25°C) in about 5 seconds but the material is still liquid. It takes about 25 seconds for the material to become a solid. Please refer to the Obturator Oven Instruction Manual for details.

Placing obturator in the canal

After the first signal “beep” of the oven, the obturator is ready for use. Push the obturator holder, wait until it goes up smoothly and take the holder with your finger. Take the obturator carefully out of the holder by first lifting it up a little and then pulling it toward you. If the obturator is not taken out, the oven will “beep” every 15 seconds to remind you that the obturator is still in the oven. After 90 seconds, the oven will switch off automatically. Insert the obturator directly into orifice of canal. Avoid touching the walls of the occlusal opening.

Note: If the canal has been properly shaped and the gutta-percha properly heated, the obturator should seat to place without twisting or forcing. To avoid overextension, do not force the obturator beyond the working length. You will note a backflow of sealer and gutta-percha resulting in an accumulation at the orifice. This is to be expected, especially in multi-rooted teeth, because the obturator is designed with excess gutta-percha to accommodate even the most widely flared or anatomically irregular canals.

Removing the GuttaCore obturator shaft and handle

Remove the shaft and handle at the orifice by bending to either side of the canal wall. Alternatively, while stabilizing the obturator with your index finger, use a round bur, or an inverted cone bur in a high-speed handpiece. A sharp spoon excavator can be also used. Dispose of obturators in an appropriate biohazard container.

Compacting gutta-percha in large, oval or irregular shaped canals

Insert a small segment of conventional gutta-percha, 4-6 mm long, alongside the shaft of the GuttaCore obturator with a lubricated plugger and compact. The gutta-percha segment compacts the heat softened gutta-percha and becomes an integral part of the filling. This will also prevent inadvertent removal of gutta-percha from the orifice with an inadequately lubricated plugger, as well as compensate for the lack of gutta-percha in

uncommonly large canals or those with internal resorption defects. Due to the tacky nature of the plasticized gutta-percha, you must use a lubricant such as ProLube root canal conditioner, sealer, topical anesthetic, etc., on the plugger.

Removing excess gutta-percha

Using an excavator, remove any excess gutta-percha that may block access to the chamber and other canals. Repeat all of the above steps on each canal of a multi-rooted tooth. If necessary, protect adjacent canal orifices from debris and gutta-percha with tiny cotton pellets or a paper point until time for obturation. Obturation is complete once excess gutta-percha is removed.

Removing GuttaCore obturators

Creating post space and retreatting are accomplished using traditional methods for removing obturation material. For creating post space, remove the GuttaCore obturator by selecting either an appropriately sized Pro-Post, GT Post drill, Gates Glidden or Peeso Reamer. When removing the obturation material for retreatment purposes, use the same sized rotary file used for shaping and engage the obturation material. Take a radiograph to confirm the removal of all obturation material.

Use a hand file with solvent to soften the gutta-percha on the GuttaCore obturator if retreatting and encountering intricate anatomy in the apical 1/3 of the canal.

10. Shelf life and storage

The shelf life is 5 years if the package (blister) remains intact. Store at temperatures below 30°C. Protect from sunlight.

11. Transportation

No special transportation conditions are required. The devices should be transported according to the rules applicable to each transport mode.

12. Environmental requirements

The product does not have an adverse impact on the environment during its life cycle.

13. Conditions of use

Use at: +10 - +45°C. The devices are not subject to repair. No special maintenance is required.

14. Disposal

Disposal is carried out in accordance with local regulations and policies and procedures of the particular medical institution.

15. Manufacturer`s warranties

The manufacturer is not responsible and shall not indemnify for any possible damage or accidents caused by:

- Use of instruments unsuitable for the product and instruments which can affect normal operation.
- Failure to follow the instructions for use.

The user is responsible for testing of materials for fitness and use for any purpose not expressly specified in the instructions!

For product quality issues, please contact the Authorized Manufacturer's Representative.

16. Authorized Representative

Dentsply Sirona LLC, Russia, 115432, Moscow, Andropov prospect 18, building 6, offices 1 39 40, Tel.: +7 (495) 725-10-87.

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