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**«06» June 2016**

## **USER MANUAL**

**Thermal processing system for allogeneic bone grafts  
models Lobator sd-2 and Lobator sd-3 with accessories**

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## 1. General

### Medical device description:

Thermal processing system for allogeneic bone grafts models Lobator sd-2 and Lobator sd-3 with accessories: 1. Transportation and storage Set, 2. Disinfection Set

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Internet address: [www.telos-marburg.de](http://www.telos-marburg.de)

## 2. Product designation

Thermal processing system for allogeneic bone grafts models Lobator sd-2 and Lobator sd-3 with accessories (hereinafter referred to as Unit) is designed for thermal disinfection of allogenic femoral head allografts from surgical (living) donors.

## 3. Field of application

Transplantology; surgery (including faciomaxillary, dental and aesthetic surgery); traumatology; orthopaedic surgery; dentistry; orthodontology.

The Unit is an electronic control unit; it is implemented in the form of the following models:

- Thermal processing system for allogeneic bone grafts model Lobator sd-2 (with an indicator panel);
- Thermal processing system for allogeneic bone grafts model Lobator sd-3 (with a touch screen).

The Unit accessories are: a transportation and storage container; a disinfection set.

The Unit provides inactivation of such viruses as HIV 1, HIV 1/0, HIV 2, human T-cell lymphotropic virus HTLV, Cytomegalovirus CMV, Hepatitis B and C viruses and also vegetative aerobic and anaerobic bacteria and fungi (S. aureus, Staphylococcus spp.(CNS) (S. epidermidis), E. faecium, Non fermenters (Pseudomonas aeruginosa), B. subtilis, M.terrae, Propionibacterium spp. (P. acnes), other anaerobic agent (without Clostridia), Aerobic diphtheroid rods, Streptococcus spp., C. albicans, A. niger).

In disinfection, the Unit maintains the transplants biomechanical and biological value.

The Unit advantages are as follows:

- ability to efficiently disinfect the allografts with maintaining the biological properties;
- efficient obtaining of the grafts from waste tissue harvested during primary hip arthroplasty procedures;
- possibility of extraction, processing and storage of the transplants in the clinics' in-house

bonebanks;

- no quarantine storage required;
- no reexamination of the donor required.

## Advices on the explant

### Maximum size of the femoral head

The Unit is designed for thermal disinfection of femoral heads up to 56 mm in diameter. Therefore, inner dimensions of the disinfection containers are not able to accept any larger femoral heads. It is recommend that larger femoral heads be divided into two parts and disinfected in 2 containers. In this case it is necessary to thoroughly document the fact that 2 different transplant identification numbers belong to one and the same donor.

### Thermal disinfection

Put the femoral head into the disinfection container free of fluid. Tightly close the container, mark its package with the transplant identification number, and put in the deep freezer. Never fill the container with Ringer's solution. In addition, take care not to interconfuse the bones to be disinfected; for this purpose, the deep freezer inner space should be divided into sections. The bone may remain in the deep freezer for as long as a year until thermal disinfection becomes possible. Prior to starting the thermal disinfection, keep in mind that the Ringer's solution will be filled in under operation room conditions according to the actual regulations and that the closed container with the femoral head will be exposed to room temperature for either 3 hours if it has been stored at  $-40 \pm 5^{\circ}\text{C}$  or 5 hours after storage at the temperature of  $-70^{\circ}\text{C}$  or lower.

### Documentation

Each individual explant should be provided with a complete set of documents. The complete documentation set shall comprise documents for both the donor and the recipient, including the informed consent, donor's history and examination report, and also the results of the serological and microbiological laboratory tests.

The requirements for the documentation set are dictated by the legislation of the country where the Unit is exploited.

## 4. Safety precautions

### 4.1. General safety rules

Read attentively the entire **Manual**. It contains important information on the safety of operation and also recommendations for properly using and maintaining the Unit.


Never use the Unit in explosion-prone areas.

Keep the Unit not closer than 10 cm to other devices or walls.


Never cover the Unit while it is switched-on.

Avoid handling liquids above the Unit.


In the process of thermal disinfection, beware of high temperatures of the heating/cooling cavity.

Use for thermal disinfection and defrosting only the special disinfection containers.

Do not clean the Unit with aggressive detergents.

**Telos will not assume any liability for direct or consequential damages arising from the inappropriate installation, transportation or maintenance of the Unit, especially in the event of violation of the User's Manual instructions or improper operation.**

4.2. Electromagnetic compatibility

The Unit user should ensure its operation under the conditions of electromagnetic environment described below.

The Unit may be used in any rooms (including living accommodations), as well as in spaces directly connected to utility low-voltage electric mains supplying electric power to dwelling houses.

is a sign warning about sensitivity to ESD (electrostatic discharge). Never touch the connector contacts marked with the ECD-sensitivity warning sign; never try to perform connection through these contacts without taking measures against ECD effects.

Inform the personnel involved in the Unit operation on the meaning of the ECD warning sign; familiarize the personnel with the measures preventing the ECD action and train them in taking these measures.

Installation warnings

Do not operate or store the Unit near other equipment. If the Unit is still operated or stored near other equipment, its normal operation cannot be guaranteed.

To connect the Unit, use only cables supplied by the Manufacturer. The use of other cables can result in the decrease in the Unit noise immunity and enhancement of electromagnetic emission.

List of cables to be supplied:

| Model        | Cable type                 | Maximum length, m |
|--------------|----------------------------|-------------------|
| Lobator sd-2 | Power-supply cable         | 2,5               |
|              | Computer-connection cables | 3                 |
| Lobator sd-3 | Power-supply cable         | 2,5               |

### 4.3. Noise immunity

The Unit user should ensure its operation under the conditions of electromagnetic environment described below.

| Noise immunity test                                                                                                                             | IEC 6 0601-<br>test load                                                                                                                                                                                                                                                                           | Conformity level                                                                                                                                                                                                                                                                                | Recommendations for electromagnetic<br>compatibility                                                                                                                                                                                                                                                                     |
|-------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Electrostatic discharge<br>(ECD) IEC 61000-4-2                                                                                                  | ± 6 kV, contact<br><br>± 8 kV, air                                                                                                                                                                                                                                                                 | ± 6 kV, contact<br><br>± 8 kV, air                                                                                                                                                                                                                                                              | The floors should be wooden, concrete or lined with ceramic tile If the floor deck is made from a synthetic material, the ambient air relative humidity should be not lower than 30%.                                                                                                                                    |
| Short electric surges/pulses<br><br>IEC 61000-4-4                                                                                               | ± 2 kV<br>for the feeder<br>lines<br><br>±1 kV for<br>incoming and<br>outgoing lines                                                                                                                                                                                                               | ± 2 kV<br>for the feeder lines<br><br>±1 kV for incoming/<br>outgoing lines                                                                                                                                                                                                                     | Quality of electric power-supply lines should meet the requirements for the commercial or hospital environment.                                                                                                                                                                                                          |
| Voltage surges<br>IEC 61000-4-5                                                                                                                 | ± 1 kV<br>differential mode<br>± 2 kV<br>common mode                                                                                                                                                                                                                                               | ± 1 kV<br>differential mode<br>± 0.5 kV<br>common mode                                                                                                                                                                                                                                          | Quality of electric power-supply lines should meet the requirements for the commercial or hospital environment.                                                                                                                                                                                                          |
| Voltage drops, short-term<br>electric power supply<br>interruptions, and voltage<br>fluctuations at the input<br>power trunks<br>IEC 61000-4-11 | <5% U <sub>T</sub><br>during 0.5 cycle<br>(> 95 % U <sub>T</sub> drop)<br><br>40% U <sub>T</sub><br>during 5 cycles<br>(60% U <sub>T</sub> drop)<br><br>70% U <sub>T</sub><br>during 25 cycles<br>(30% U <sub>T</sub> drop)<br><br>5% U <sub>T</sub> during 5<br>sec (95 % U <sub>T</sub><br>drop) | <5% U <sub>T</sub><br>during 0.5 cycle<br>(> 95 % U <sub>T</sub> drop)<br><br>40% U <sub>T</sub><br>during 5 cycles<br>(60% U <sub>T</sub> drop)<br><br>70% U <sub>T</sub><br>during 25 cycles<br>(30% U <sub>T</sub> drop)<br><br>5% U <sub>T</sub> during 5 sec<br>(95 % U <sub>T</sub> drop) | Quality of electric power-supply lines should meet the requirements for the commercial or hospital environment.<br><br>If the user has to continue operating the Unit when the electric power supply is interrupted, it is recommended that the Unit be fed from an Uninterruptible Power Supply or accumulator battery. |
| Industrial-frequency<br>magnetic field<br>IEC 61000-4-8                                                                                         | 3 A/m                                                                                                                                                                                                                                                                                              | 3 A/m                                                                                                                                                                                                                                                                                           | At the assumed place of operation, the industrial-frequency magnetic field should be measured in order to make sure that it is sufficiently low.                                                                                                                                                                         |
| NOTE U <sub>T</sub> is the ac power-supply voltage prior to applying the test load.                                                             |                                                                                                                                                                                                                                                                                                    |                                                                                                                                                                                                                                                                                                 |                                                                                                                                                                                                                                                                                                                          |

The interference can occur in the vicinity of equipment marked with the following symbol:

Never use portable and mobile radio-frequency communication devices closer to any Unit part (including cables) than the recommended separation distance.

The electromagnetic field of the frequencies higher than 150 kHz to 80 MHz should be lower than 3 V/m.

The Unit user can prevent electromagnetic interference by observing the minimal distance between the portable and mobile communication devices (transmitters) depending on the communication device output power according to the following recommendations.

| Recommended separation distances between the Unit and portable/mobile communication devices                                                                                                                                                                                                                        |                                |                               |                                  |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------|-------------------------------|----------------------------------|
| The transmitter nominal output power /Watt                                                                                                                                                                                                                                                                         | Separation distance /mm        |                               |                                  |
|                                                                                                                                                                                                                                                                                                                    | 150 kHz - 80 MHz<br>d= 1.17 VP | 80 MHz - 800 MHz<br>d=1.17 VP | 800 MHz - 0<br>2.5 GHz d=2.33 VP |
| 0.01                                                                                                                                                                                                                                                                                                               | 0.12                           | 0.12                          | 0.23                             |
| 0.1                                                                                                                                                                                                                                                                                                                | 0.37                           | 0.37                          | 0.73                             |
| 1                                                                                                                                                                                                                                                                                                                  | 1.17                           | 1.17                          | 2.30                             |
| 10                                                                                                                                                                                                                                                                                                                 | 3.70                           | 3.70                          | 7.30                             |
| 100                                                                                                                                                                                                                                                                                                                | 11.70                          | 11.70                         | 23                               |
| For transmitters whose maximum output power is beyond the above-specified ranges, the separation distance can be calculated using the formula presented in the respective column, where P is the maximum output power of the transmitter (in W) according to the data presented by the transmitter's manufacturer. |                                |                               |                                  |
| NOTE 1: In case of 80 MHz and 800 MHz, a higher frequency range shall be used.                                                                                                                                                                                                                                     |                                |                               |                                  |
| NOTE 2: These recommendations are not always applicable since the electromagnetic waves propagation depends on their absorption by and reflection from structures, things and people.                                                                                                                              |                                |                               |                                  |



## 5. Main performance characteristics



| Product                                                               | Thermal processing system for allogeneic bone grafts                                                                         |                                                                                                |
|-----------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------|
| Model                                                                 | Lobator sd-2                                                                                                                 | Lobator sd-3                                                                                   |
| Overall dimensions (height x width x depth), mm                       | 140x345 x325                                                                                                                 | 150x300x340                                                                                    |
| Weight, kg                                                            | 7,5                                                                                                                          | 5,8                                                                                            |
| Article No.                                                           | 8700000                                                                                                                      | 5000000                                                                                        |
| Power-supply voltage                                                  | 230 V, 50-60 Hz                                                                                                              |                                                                                                |
| Power consumption                                                     | 230 VA                                                                                                                       |                                                                                                |
| PC Interface                                                          | An IBM-compatible PC with a free serial port                                                                                 |                                                                                                |
| * Unit Software                                                       | EURO-230V                                                                                                                    | EURO-230V                                                                                      |
| Required PC Software                                                  | version 92 (T092), 05.08.2005<br>Windows, version 4.2 and higher                                                             | version 92 (T092), 05.08.2005<br>Windows Embedded for Point of Service, version 1.0 and higher |
| Control panel type                                                    | Indicator panel                                                                                                              | Touch screen                                                                                   |
| Type of PC connection                                                 | Serial port<br>RS232 protocol<br>max.data rate -19200 bit/s                                                                  | USB<br>version 2.0<br>max. data rate - 480 Mbit/s                                              |
| Dust and moisture protection degree                                   | IP 00                                                                                                                        | IP 00                                                                                          |
| ** Fuse type                                                          | 1A, 250 V                                                                                                                    | 1A, 250 V                                                                                      |
| Battery type                                                          | 2 x AA                                                                                                                       | -                                                                                              |
| Noise level                                                           | 35 db                                                                                                                        | 35 db                                                                                          |
| Audible alarm                                                         | The duration of sound signals - 0.3 sec<br>The frequency of the sound signal - 3.5 + 0.5 kHz<br>Sound pressure level - 43 dB |                                                                                                |
| Storage and transportation conditions                                 |                                                                                                                              |                                                                                                |
| Temperature, °C                                                       | -40 - +70                                                                                                                    |                                                                                                |
| Air humidity, %                                                       | 10 - 95 without condensation                                                                                                 |                                                                                                |
| Atmospheric pressure, GPa                                             | 500-1060                                                                                                                     |                                                                                                |
| Operating conditions                                                  |                                                                                                                              |                                                                                                |
| Temperature, °C                                                       | +18 - +30                                                                                                                    | +10 - +40                                                                                      |
| Air humidity, %                                                       | 30-75                                                                                                                        |                                                                                                |
| Atmospheric pressure, GPa                                             | 700 -1069                                                                                                                    |                                                                                                |
| Classification                                                        | Shock Protection Class I, Type B                                                                                             |                                                                                                |
| as per Appendix IV of Guideline 93/42/EEC                             | Active, Class 1, accessories – Class 2a                                                                                      |                                                                                                |
| Approvals as per recommended schedule of labor protection inspections | DIN EN 60601-1 / DIN EN 60601-1-2<br>every 2 years                                                                           |                                                                                                |
| Warranty Period                                                       | 2 years                                                                                                                      |                                                                                                |
| Availability of Spare Parts                                           | 10 years                                                                                                                     |                                                                                                |
| Certification body                                                    | DEKRA Certification GmbH                                                                                                     |                                                                                                |
| Identification Number                                                 | CE 0124                                                                                                                      |                                                                                                |

\* - The manufacturer does not provide an update revision or reprogramming the software settings by the user. The manufacturer or the Authorized Customer Service Representative carry out the software installation in accordance with the manufacturer's service instructions.


\*\* - Maintenance service provides the replace of the fuse in accordance with the service instructions

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


fel s

6. Delivery set

6.1 Thermal processing system for allogeneic bone grafts model Lobator sd-2 Delivery set

| No. | Component parts                                                                 | Manufacturer                                                | Country |
|-----|---------------------------------------------------------------------------------|-------------------------------------------------------------|---------|
| 1   | Thermal processing system for allogeneic bone grafts model Lobator sd-2 – 1 pc. | TELOS Herstellung und Vertrieb med. technischer Geräte GmbH | Germany |
| 2   | Disinfection Set – maximum 1000 pcs.                                            |                                                             |         |
| 3   | Transportation and storage Set – maximum 1000 pcs.                              |                                                             |         |
| 4   | User's Manual – 1 pc.                                                           |                                                             |         |

6.2 Thermal processing system for allogeneic bone grafts model Lobator sd-3 Delivery set

| No. | Component parts                                                                | Manufacturer                                                | Country |
|-----|--------------------------------------------------------------------------------|-------------------------------------------------------------|---------|
| 1   | Thermal processing system for allogeneic bone grafts model Lobator sd-3– 1 pc. | TELOS Herstellung und Vertrieb med. technischer Geräte GmbH | Germany |
| 2   | Disinfection Set – maximum 1000 pcs.                                           |                                                             |         |
| 3   | Transportation and storage Set – maximum 1000 pcs.                             |                                                             |         |
| 4   | User's Manual – 1 pc.                                                          |                                                             |         |

7. Unit connection


Check the power-supply voltage and connect the Unit by using the cable taken from the Delivery Set.


The power-supply voltage will be directly proportional to the heating output. To avoid any risk of undervoltage, please make sure that the Unit has not been connected through a distribution outlet. If the voltage drops below 207V, the process will be aborted (error 85 for Lobator sd-2 and Error xx for Lobator sd-3).

15

## 8. Thermal processing system for allogeneic bone grafts model Lobator sd-2

### 8.1. Controls

## 8.2. Thermal disinfection program

1. The Unit ON/OFF button is on the rear side (Fig. 1)
2. A blinking green arrow prompts you to properly insert the disinfection container into the heating cavity (Fig. 2)
3. After insertion of the disinfection container, the visual display ceases. Push the green START button (Fig. 3) to start the disinfection process (duration: 94 minutes) (Fig. 4)
4. The red LED now displays the three-digit transplant identification number (Fig. 5), which is to be indicated on the disinfection set with a waterproof marker. The timer counts down thus indicating the time remaining till the completion of the disinfection process.
5. A continuous sound signal and red cross in the Error display indicate the error and abortion of the disinfection process (Fig. 6).
6. A continuous sound signal and red blinking arrow indicate the end of the disinfection process. The disinfection container can now be removed.
7. To complete the disinfection process properly, push the blinking "Print" button at the Unit front panel (Fig. 7). 7). Being connected, the printer will print the protocol. The Unit will be available for a new process only after pushing the "Print" button.

Fig. 1

Fig. 2

Fig. 3

Fig. 4

Fig. 5

Fig. 6

Fig. 7

## 8.3. Errors / Malfunction Messages

During the thermal disinfection process, microprocessors are constantly controlling the power supply voltage, heating curve, the process running time, and the magnetic agitator drive.

If any failure takes place in the course of the thermal disinfection, a warning signal sounds (and a cross symbol lights up on the error display thus indicating abortion of the process. In this case the device is to be turned off for a few seconds. After restarting the Unit, the Date and Time Display will show "Error" as well as the corresponding number of the error that has caused the process abortion.

When the "PRINT" button is pushed, the process is saved, and the Unit returns to its initial settings. The Unit will become available again after the heating cavity temperature decreases below 35°C. This shall prevent the Ringer's solution overheating due to the preheated heating/cooling cavity).

## 8.3.1. Potential sources of errors

Potential error sources ~~will be checked~~ automatically in switching on the Unit.

| Error No. | Error                               | Error No. | Error                                                                           |
|-----------|-------------------------------------|-----------|---------------------------------------------------------------------------------|
| 0         | Too low glass temperature           | 18        | Weak signal of the LED indicator                                                |
| 1         | Too low LED voltage                 | 19        | Too bright LED indicator signal                                                 |
| 2         | Too low battery charging voltage    | 20        | Too high voltage of LED display X                                               |
| 3         | Poor software version               | 21        | 160° - control system response                                                  |
| 4         | Too small heating pot               | 22        | IIC bus error - IC32 (the detector port element) is not adequate to the iic bus |
| 5         | 5 V, too low voltage                | 23        | Insufficiently high level of the IC32 port element (detector LED)               |
| 6         | 15 V, too high voltage              | 24        | Unremovable IC32 port element                                                   |
| 7         | Power is too low for LED display X  | 25        | Breaker RTCC is out of operation                                                |
| 8         | Too high glass temperature          | 26        | Interrupting zero crossing circuit is out of operation                          |
| 9         | Too high LED voltage                | 27        | Hardware failure of IIC bus, FRAM or IC32                                       |
| 10        | Too high battery charging voltage   | 28        | The timer does not transfer data to IIC bus                                     |
| 11        | Too high software version           | 29        | The LPT IC1 port element does not transfer data to IIC bus                      |
| 12        | Too large heating pot               | 30        | The LPT IC2 port element does not transfer data to IIC bus                      |
| 13        | 5 V, too high voltage               | 31        | The LPT port element IC3 does not transfer data to IIC bus                      |
| 14        | 15 V, too low voltage               | 32        | The timer detector is not open                                                  |
| 15        | Power is too high for LED display X | 33        | FRAM IC26 does not transfer data to IIC bus                                     |
| 16        | Too loud audible signal             | 34        | The FRAM checknumber error, FRAM variables                                      |
| 17        | Too low audible signal              | 36        | FRAM constants, incorrect checksum                                              |

## 8.3.2. Potential thermal disinfection errors able to cause abortion

| Error No. | Error                                                                                                                                                                                                                                               |
|-----------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 80        | Agitator failure                                                                                                                                                                                                                                    |
| 81        | Control processor failure                                                                                                                                                                                                                           |
| 82        | Actuation of the main processor safety fuse                                                                                                                                                                                                         |
| 83        | Power-supply failure for more than 30 sec<br>Please check whether the plug is properly inserted into the power supply.                                                                                                                              |
| 84        | The container was removed from the heating cavity for more than 30 sec<br>The disinfection process will be continued after the container is placed back into the cavity within 30 seconds. The time of removal is controlled by the microprocessor. |
| 85        | The ADC parameters are beyond the permitted range (207 V or below)<br>Make sure that the Unit has not been plugged into the distribution outlet.                                                                                                    |
| 86        | Control error, no deactivation for more than 5 min. Switch off the Unit and restart it as described above.                                                                                                                                          |
| 87        | The glass temperature has not exceeded 70°C in 19 min                                                                                                                                                                                               |
| 88        | The agitator failure due to wrong ADC parameters                                                                                                                                                                                                    |

In case the display indicates the error No., and You fail to eliminate the problem by yourself, please contact an authorized representative of TELOS-rus LTD.

## 8.4. Disinfection, printing the protocols

The Unit serial No.  
Transplant identification No.

---

Date and time of thermal disinfection

---

Hex - values

---

Temperatures of the heating/cooling cavity  
Measurements every 3 min

---

Process completed

---

Process aborted

---

Signature

**Lobator sd-2**  
 ● Identnr.: 2240\_063  
 ● Date (DMY): 27.09.05  
 ● Time : 11:42

**IntP. :**  
 42 11 67 C9 05 02 00 15  
 2A 3C 4E 5F 6D 78 81 86  
 8B 87 7D 76 73 6F 6D 66  
 61 5D 5E 5D 5D 59 47 3D  
 37 33 2F 2D 2A 28 27 28

**Temp. :**  
 28 45 58 72 87 101  
 116 130 140 136 126 121  
 115 112 111 102 94 92  
 92 93 93 93 86 64  
 52 45 41 39 37 35  
 34 33

**PROCESS COMPLETED**  
 PROCESS ABORTED

SIGNATURE

- The automatic disinfection process consists of three phases: the heat-up phase, the steady state temperature phase, and the cool-down phase. The electronic control system ensures a temperature of 82.5°C in the center of a femoral head up to 56 mm in

diameter for **15 minutes**. The process duration is 94 minutes.

- In order to **expose the femoral head center** to minimum 82.5°C for minimum 15 minutes, the heating/cooling cavity should be heated up to a much higher temperature (up to 140°C).
- For the reasons of **sterility**, the Ringer's solution temperature in the closed disinfection container **cannot be measured** directly. This is why the temperature detector is mounted on the **container outer surface**.
- Therefore, **the temperature** displayed during the process is the temperature measured **between the disinfection container and the heating/cooling cavity** in °C. Each segment corresponds to 5°C. The display will be actuated at 35°C.
- The microprocessor continuously compares the temperatures measured at the bottom of the heating/cooling cavity with those measured between the heating/cooling cavity and disinfection container. Both values are compared with the preprogrammed nominal temperature curve.
- If the disinfection process has been performed without failures, annotation "Process Completed" appears at the end of the protocol. A failure or interruption of the process will be documented as "Process Aborted".
- Using the cables included in the Delivery Set, the Unit can be connected to a PC; the provided software ensures simple Unit-PC communication. The last 200 data blocks will be saved in the Unit memory and can be selected and displayed as a chart. Using a respective button, the desired protocol may be accessed and reprinted for filing.

## 8.5. The DEFROST program

This program will be used to defrost the frozen transplant.

The DEFROST period is 8 minutes; the process may be repeated several times or, if necessary, terminated prematurely by removing the container.

1) Turn on the power switch at the rear side of the Unit to set the device to standby mode.

2) The blinking green arrow prompts the operator to properly and completely insert the transplant container into the heating cavity.

3) After inserting the container, the blinking display ceases, and the defrost process may be activated by pushing the yellow DEFROST button (Fig. 1).

4) The termination of the process is indicated by an audible signal and red blinking arrow. Now the Operator may remove the transplant container from the Unit.

5) To switch off the Unit, push the power switch mounted on the Unit rear side.


Fig 1

**Prior to inserting the container into the heating/cooling cavity for defrosting, remove its outer package.**

## 8.6. Battery change

- The Date and Time display needs accumulator batteries. Before they get empty, the Date and Time display indicates "LOW BATT".
- To prevent the date and time reset while changing the batteries, keep the Unit plugged into the electrical outlet and turned-on.
- Tilt the device laterally and open the battery compartment at the bottom of the Unit by using an appropriate tool. Exchange the batteries.
- 2 Mignon batteries AA are required. Make sure that their poles are positioned correctly.
- Close the battery compartment and return the Unit to the home position.


## 8.7. Date and Time Settings

- To set the date and time, switch on the Unit. However, the thermal disinfection or defrosting process should not be active at this time.
- Use the "select" button at the rear side of the device to select the display for changing.
- Now the blinking display indications may be changed by using the "+" or "-" buttons.
- Having run through all the front displays, the device is again available for operation.


## 9. Thermal processing system for allogeneic bone grafts model Lobator sd-3

### 9.1. Controls and Displays


#### 9.1.1. Main screen


## 9.1.2. Date and Time Input


## 9.1.3. Desktop

## 9.1.4. Administrator


## 9.1.5. Language


## 9.1.6. Volume



## 9.1.7. Archive






## 9.1.8. DEFROST program (1)

## 9.1.9. DEFROST program (2)

## 9.2. Thermal Disinfection Procedure

1. To turn the Unit on, use the special switch mounted on the right from the touch screen. Set the disinfection container into the process panel/heating cavity.


(If the Unit has been already turned on, set the disinfection container into the process panel and then continue the process from the Main Menu, item 4).

2. In 45 seconds the screen displays the date and time. Confirm these data or correct them.
3. The Main Menu appears on the screen in 5 s. The icon of the thermal disinfection start will be in the bottom right corner. The process duration: 94 minutes.


4. An audible signal informs on the process startup. In the process, two temperature curves will be plotted on-line. The blue curve presents the heating cavity temperature, the red one shows that of the container walls. In the top of the screen, the Unit name, process No., data and time are displayed.


5. If the disinfection container is removed during the process, first the verbal message sounds, and then acoustic and optical warnings are generated.

If the disinfection container stays out of the heating cavity for more than 30 seconds, the thermal disinfection will be aborted. If this takes place after the first 64 minutes have passed, only the Scald Warning is generated since in this case the sample will be already disinfected and the only thing it needs is cooling.


The moment and duration of the container removal will be shown graphically on the temperature curve.

6. In 94 minutes, a continuous audible signal sounds, which indicates that the thermal disinfection is completed, and the container may be removed. In addition, the Unit informs on creating a backup PDF-file.

7. Each process that has been started can be found in the archive regardless of whether it has been completed or not. This is true for all the previously accomplished processes. To enter the archive, push button "Verwaltung" (Administrator) and then button "Archiv" (Archive). After that you may look through the process protocol, print it or save on the USB stick.

8. To turn on the Unit, push button "Aus" (OFF) on the main desktop.

## 9.3. Case disinfection process curve


### Measurements processing

telos Lobator sd-3 Measurements estimation

For the PC as at: 22.06.2011 14 10:57 / Interface: 29.10.2010 \*8 33 36

Disinfection startup: 01.07.2011 11:21:11 hr  
 Device No.: 30003  
 Ident. No.: 00347  
 Archive catalogue: sc3-30003  
 Archive file: sc3-30003-00347.PDF  
 Comments: Error-free disinfection.

Signature:

| Time (min) | Cavity t [°C] | Container t [°C] | Error time [s] | Time (min) | Cavity t [°C] | Container t [°C] | Error time [s] | Time (min) | Container t [°C] | Container t [°C] | Error time [s] |
|------------|---------------|------------------|----------------|------------|---------------|------------------|----------------|------------|------------------|------------------|----------------|
| 0          | 29,06         | 28,22            | 0 / 0 / 0      | 32         | 125,19        | 97,51            | 0 / 0 / 0      | 64         | 93,02            | 79,25            | 0 / 0 / 0      |
| 1          | 37,29         | 28,30            | 0 / 0 / 0      | 33         | 121,86        | 97,17            | 0 / 0 / 0      | 65         | 87,06            | 78,22            | 0 / 0 / 0      |
| 2          | 42,01         | 29,40            | 0 / 0 / 0      | 34         | 120,64        | 96,47            | 0 / 0 / 0      | 66         | 79,75            | 75,83            | 0 / 0 / 0      |
| 3          | 46,57         | 30,91            | 0 / 0 / 0      | 35         | 118,54        | 95,92            | 0 / 0 / 0      | 67         | 73,14            | 72,23            | 0 / 0 / 0      |
| 4          | 52,49         | 32,59            | 0 / 0 / 0      | 36         | 116,45        | 95,26            | 0 / 0 / 0      | 68         | 67,90            | 68,30            | 0 / 0 / 0      |
| 5          | 59,29         | 34,74            | 0 / 0 / 0      | 37         | 114,70        | 94,57            | 0 / 0 / 0      | 69         | 63,36            | 64,38            | 0 / 0 / 0      |
| 6          | 66,26         | 37,13            | 0 / 0 / 0      | 38         | 114,01        | 94,05            | 0 / 0 / 0      | 70         | 59,79            | 60,64            | 0 / 0 / 0      |
| 7          | 71,89         | 39,82            | 0 / 0 / 0      | 39         | 113,49        | 93,54            | 0 / 0 / 0      | 71         | 56,73            | 57,33            | 0 / 0 / 0      |
| 8          | 77,71         | 42,68            | 0 / 0 / 0      | 40         | 112,44        | 93,02            | 0 / 0 / 0      | 72         | 54,01            | 54,35            | 0 / 0 / 0      |
| 9          | 83,72         | 45,55            | 0 / 0 / 0      | 41         | 110,87        | 92,50            | 0 / 0 / 0      | 73         | 51,81            | 51,64            | 0 / 0 / 0      |
| 10         | 89,74         | 48,60            | 0 / 0 / 0      | 42         | 109,83        | 91,98            | 0 / 0 / 0      | 74         | 49,95            | 49,44            | 0 / 0 / 0      |
| 11         | 94,23         | 51,64            | 0 / 0 / 0      | 43         | 110,00        | 91,47            | 0 / 0 / 0      | 75         | 48,26            | 47,41            | 0 / 0 / 0      |
| 12         | 99,42         | 54,79            | 0 / 0 / 0      | 44         | 106,99        | 91,12            | 0 / 0 / 0      | 76         | 46,91            | 45,72            | 0 / 0 / 0      |
| 13         | 104,44        | 57,92            | 0 / 0 / 0      | 45         | 105,14        | 90,43            | 0 / 0 / 0      | 77         | 45,72            | 44,37            | 0 / 0 / 0      |
| 14         | 107,74        | 60,98            | 0 / 0 / 0      | 46         | 101,84        | 89,60            | 0 / 0 / 0      | 78         | 44,71            | 43,02            | 0 / 0 / 0      |
| 15         | 112,26        | 64,04            | 0 / 0 / 0      | 47         | 100,45        | 88,71            | 0 / 0 / 0      | 79         | 43,71            | 41,84            | 0 / 0 / 0      |
| 16         | 116,10        | 66,94            | 0 / 0 / 0      | 48         | 98,38         | 87,67            | 0 / 0 / 0      | 80         | 43,02            | 40,88            | 0 / 0 / 0      |
| 17         | 118,54        | 70,01            | 0 / 0 / 0      | 49         | 95,96         | 86,64            | 0 / 0 / 0      | 81         | 42,18            | 39,99            | 0 / 0 / 0      |
| 18         | 121,51        | 72,74            | 0 / 0 / 0      | 50         | 94,92         | 85,61            | 0 / 0 / 0      | 82         | 41,53            | 39,31            | 0 / 0 / 0      |
| 19         | 124,67        | 75,48            | 0 / 0 / 0      | 51         | 93,89         | 84,75            | 0 / 0 / 0      | 83         | 41,00            | 38,54            | 0 / 0 / 0      |
| 20         | 127,64        | 78,22            | 0 / 0 / 0      | 52         | 93,36         | 83,89            | 0 / 0 / 0      | 84         | 40,32            | 37,97            | 0 / 0 / 0      |
| 21         | 130,80        | 80,83            | 0 / 0 / 0      | 53         | 92,85         | 83,03            | 0 / 0 / 0      | 85         | 39,82            | 37,46            | 0 / 0 / 0      |
| 22         | 133,44        | 83,37            | 0 / 0 / 0      | 54         | 93,71         | 82,34            | 0 / 0 / 0      | 86         | 39,37            | 36,96            | 0 / 0 / 0      |
| 23         | 135,02        | 85,79            | 0 / 0 / 0      | 55         | 93,36         | 81,83            | 0 / 0 / 0      | 87         | 38,98            | 36,45            | 0 / 0 / 0      |
| 24         | 137,49        | 88,19            | 0 / 0 / 0      | 56         | 92,85         | 81,31            | 0 / 0 / 0      | 88         | 38,48            | 35,95            | 0 / 0 / 0      |
| 25         | 139,08        | 90,26            | 0 / 0 / 0      | 57         | 93,71         | 80,97            | 0 / 0 / 0      | 89         | 38,14            | 35,61            | 0 / 0 / 0      |
| 26         | 139,78        | 92,33            | 0 / 0 / 0      | 58         | 93,36         | 80,63            | 0 / 0 / 0      | 90         | 37,80            | 35,27            | 0 / 0 / 0      |
| 27         | 139,61        | 94,23            | 0 / 0 / 0      | 59         | 92,92         | 80,45            | 0 / 0 / 0      | 91         | 37,46            | 34,91            | 0 / 0 / 0      |
| 28         | 139,55        | 95,78            | 0 / 0 / 0      | 60         | 93,54         | 80,11            | 0 / 0 / 0      | 92         | 37,13            | 34,52            | 0 / 0 / 0      |
| 29         | 136,18        | 96,99            | 0 / 0 / 0      | 61         | 93,02         | 79,77            | 0 / 0 / 0      | 93         | 36,79            | 34,27            | 0 / 0 / 0      |
| 30         | 131,86        | 97,68            | 0 / 0 / 0      | 62         | 93,71         | 79,60            | 0 / 0 / 0      |            |                  |                  |                |
| 31         | 127,82        | 97,85            | 0 / 0 / 0      | 63         | 93,21         | 79,42            | 0 / 0 / 0      |            |                  |                  |                |

The error time = the moment of the container extraction / agitator current below/above the rating

Archive file saved: 01.07.2011 12:55:15

## 9.4. The DEFROST procedure

This program will be used to defrost the frozen transplant.  
The DEFROST period is 8 minutes; the process may be repeated several times or, if necessary, terminated prematurely by removing the container.

In this case, the transplant identification No. will not be specified, and the protocol will not be drawn.

1. To turn on the Unit, use the special switch mounted on the right from the touch screen.

Insert the disinfection container into the heating cavity.

2. In 45 seconds the screen will display the date and time. Confirm or correct the data.


3. The Defrost icon shall be on the right in the screen. The process duration: 8 minutes.

4. During defrosting, the screen will display the heating cavity temperature, container temperature, elapsed time and remaining time.

In addition, the status audible message will be generated at the beginning and at the end of the process.

5. After an arbitrary number of process repetitions, you may either start the thermal disinfection or turn on the Unit.

## 9.5. Failure/error messages

During the thermal disinfection process, microprocessors are constantly controlling the power supply voltage, heating curve, the process running time, and the magnetic agitator drive.

If any failure takes place in the course of thermal disinfection, a warning audible signal sounds (and a cross symbol lights up on the Error Display indicating abortion of the process. In this case the device is to be turned off for a few seconds. After restarting the Unit, the Date and Time Display will show "Error" as well as the corresponding number of the error that has caused the process abortion.

When the "Print" button is pushed, the process will be saved, and the Unit will return to its initial settings. The Unit will become available again after the heating/cooling cavity temperature has decreased below 35°C. This shall prevent the Ringer's solution overheating due to the preheated heating/cooling cavity).

1. Possible errors to be automatically checked by the Unit after activation.
2. The thermal disinfection errors giving rise to failures that cause the process abortion.

## 10. Accessories

### 10.1. Disinfection Set

#### Description

The special disposable container consists of the outer package and disinfection container. The two-piece outer package is made of polystyrene and can be tightly and sterilely closed due to the thread sealed by Evoprene.

The disinfection container is made from Makrolon and consists of two chambers. The bottom part (container) will accommodate the femoral head and be filled with the disinfectant up to the mark. The top part (lid) that is also referred to as a transfer container will be used as a sterile closure of the bottom chamber and a disinfection fluid receiver after the process completion.

#### Function:


Function: the container is designed to accommodate the transplant for further thermal disinfection and following storage of the allogenic femoral heads under the deep frost conditions.

Disinfection Set is not an independent medical device and is used only in common with the Unit.

## Warnings / Precautions

The red safety thread (seal) on the outer package should be absolutely intact.

Close the valve of the transfer container by turning it by 360° only after the fluid is completely drained. Otherwise, proper functioning of the disinfection container cannot be guaranteed.

## Handling instructions

Prior to using the disinfection set for the first time, in any case it is necessary to be trained and instructed by a consuler competent in medical products or an experienced colleague.

1. The disinfection set will be delivered by an assistant or surgical nurse being non-sterile.

2. Cut the red protective thread and open the outer package.  
Observing the sterility, take the disinfection container out of the package.


3. Open the disinfection container.

4. The proper-size femoral head will be put into the container by the operating surgeon.

Then fill the container with the Ringer's or table-salt solution up to the top mark.

5. Close the container by tightly screwing on it the transfer container.

6. Then install the disinfection container into the Unit process cavity (heating/cooling element).

## On the completion of thermal disinfection:

1. When thermal disinfection is completed, turn over the disinfection container.


The fluid will drain through the open valve into the transfer container.

2. To prevent the backflow and ensure in the container a Class A clean space, thoroughly close the transfer container in its position after draining by turning the shutoff valve up to the stop (by 360°).

Then again turn over the disinfection container.

2


3. Remove the small blue protective cover from the transfer container top and covers from the hemoculture containers; spray twice the disinfection spray.

4. Take 20 ml from the container, change the needles, and transfer 10 ml to each of the two hemoculture bottles (1 aerobic and 1 anaerobic).

Hand over the hemoculture systems to the clinic's laboratory.

4


5. Press the rubber stopper into the transfer container and drain the fluid through this opening into a kidney dish.

Close the opening with the small blue protective cover.

6. Using a waterproof marker, write on the transfer container identification No. of the transplant that have passed the processing and return the disinfection container into the outer package.

When working with several Units, it is necessary to mark the transfer container with the Unit serial No.; otherwise, the process clarity and accountability cannot be guaranteed.

6

7. Seal the outer package with a yellow protective thread.

Put the disinfection set into the deep freezer and store it at the temperature of  $-70^{\circ}\text{C}$  or lower for 5 years or at  $-40 \pm 5^{\circ}\text{C}$  or lower for 2 years.

To provide for a consistent monitoring and documentation of the cooling process, the refrigerator should be equipped with an alarm and a digital temperature detector.

After analysis of the selected samples and obtaining a release from the "responsible person", seal the container with the green protective thread (seal).

In storing the samples in the deep freezer, the already approved and yet unapproved transplants should be thoroughly spatially separated from each other.

## The disinfection set marking:

Certification body: CE 0124

Signs: "Caution! Refer to accompanying documents",

"Never reuse", "Sterilized by radiation"

For Lobator only

Article No.: (8300001)

Lot Number: (30 019-02)

Date of production: (2014-12-12)

Expiry date: (2017-06-30)

Unit identification No. \_\_\_\_\_ (to be filled)

Manufacturer's name and address

Note: "The red seal should be intact"

## 10.2. Transportation and storage Set

### Description

Transportation and storage Set is the special disposable container, consists of the outer package and inner disinfection container (vessel). The outer package is made from polystyrene and can be sterilely closed due to the thread sealed by Evoprene. The inner container is made from Makrolon.

### Function:

transportation and storage of transplants from the moment and place of extraction to the moment and place of thermal disinfection.

Transportation and storage Set is not an independent medical device and is used only for purposes of subsequent thermal disinfection bone grafts together with Unit and Disinfection Set. Transportation and storage Set is not directly intended for thermal disinfection transplant on the Unit, the Disinfection Set is used for this purposes.

## Warnings / Precautions

The red safety thread (~~seal~~) on the outer package should be absolutely intact.

### Transportation and storage Set marking:

Certification body: **CE 0124**

Signs: "Caution! Refer to accompanying documents",  
"Never reuse", "Sterilized by radiation"

Not for Lobator

Article No.: (8250000)

Lot Number: (60 003-03)

Date of production: (2014-11-05)

Expiry date: (2017-06-30)

Product name

Manufacturer's name and address

Note: "The red seal should be intact"

## 11. Documents management and storage

The Unit is provided with two mutually independent data storage devices. Each process will be documented on both devices since the moment of putting the Unit into operation. The memory size is designed for 20 years of intense operation.

The USB-stick automatically saves duplicates of all the process data registered by the Unit. Any other data storage device may be also used.

In addition to simultaneous saving of all the process data on the USB-stick, it is possible to intentionally duplicate separate processes. For this purpose, it is necessary to open the Archive prior to inserting the USB-stick into the USB socket. After the USB-stick has been inserted, an icon appears on the left of the Archive; press this icon after selecting the process. Only the selected process will be copied.

In the case of remote maintenance, as well as in the case of controlling observation of the Unit labor protection rules, all the processes will be stored at an individual server. As the Unit does not fix the patients' data, requirements of the data protection law will be fulfilled. Each Unit user may free request the data on his/her Unit.

## 12. Warranty period

2 years.

## 13. Transportation conditions


The Unit and accessories may be transported in the Manufacturer's package by covered vehicle of any kind according to the cargo transportations rules actual for every transport mode at a temperatures of  $-40^{\circ}\text{C}$  to  $+70^{\circ}\text{C}$  and relative humidity of 10 to 95% w/o condensation.

## 14. Storage conditions

Store the Unit and accessories at indoor storages at air temperatures of  $-40^{\circ}\text{C}$  to  $+70^{\circ}\text{C}$  and relative humidity of 10 to 95% w/o condensation.

## 15. Unit marking

Marking of Thermal processing system for allogeneic bone grafts model Lobator sd-2

### Decoding of the Lobator sd-2 marking:

Product name: Lobator sd-2

Article No.: 8700000

Serial No.: SN 2340

Manufacturer: TELOS

Operating parameters: 230 V, 50/60 Hz, 230 VA

Safety fuse: 1A, 250 V

Made in Germany

## Marking of Thermal processing system for allogeneic bone grafts model Lobator sd-3

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Operating parameters: 230 V, 50/60 Hz, 230 VA

Safety fuse: 1A, 250V

Prior to opening, de-energize the Unit!

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Made in Germany

EC marking

### Warning signs for models Lobator sd-2, Lobator sd-3:

|  |                                                                                                                                             |
|--|---------------------------------------------------------------------------------------------------------------------------------------------|
|  | CE Mark:<br>The product is certified in accordance with the EC Directive 93/42/EEC                                                          |
|  | Waste parts of electric and electronic equipment:<br>At the end of its service life, the device should be disposed via the recycling system |
|  | Electric shock protection (type B) symbol                                                                                                   |
|  | ESD sensitivity sign                                                                                                                        |

## 16. Disposal

Dispose the Unit and accessories in the way adopted by the Medical Treatment and Prevention Institution.

At the end of its service life, the Unit will be disposed via the recycling system as waste electronic and electric equipment.

Accessories that have been used at the M&P Institution should be disposed as Class B medical wastes according to the State regulations and standards actual by the moment of the disposal.

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