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
Medtronic

Medtronic Inc.
710 Medtronic Parkway
Minneapolis, MN 55432-5604
www.medtronic.com
763 505 2900

To whom it may concern:


It is hereby certified that the instructions for use, **INVISx Cranial Fixation System, number 14578-1F,** and **Symbol Definitions 15782-1N,** and to which this certificate are attached, are the original copies.

Sincerely,

Carla P. Carrillo
Regulatory Affairs Specialist


Subscribed and sworn to me before this 16th day of January
month, 2015.

Katherine Perez

Jan 16th, 2015

Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432 USA

*INSTRUCTIONS FOR Use / BRUGSANVISNING / GEBRUIKSAANWIJZING /
MODE D'EMPLOI / GEBRAUCHSANWEISUNG / ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ /
HASZNÁLATI UTASÍTÁS / ISTRUZIONI PER L'USO / BRUGSANVISNING /
INSTRUKCJA UŻYTKOWANIA / INSTRUÇÕES PARA USO /
INSTRUCCIONES DE USO / BRUGSANVISNING / KULLANMA TALIMATI*

INVISx™ CRANIAL FIXATION SYSTEM
INVISx KRANIEFIKSERINGSSYSTEM
INVISx CRANIAAL FIXATIESYSTEEM
SYSTÈME DE FIXATION CRÂNIENNE INVISx
KRANIALES KNOCHENDECKEL-
FIXATIONSSYSTEM INVISx
ΣΥΣΤΗΜΑ ΚΡΑΝΙΑΚΗΣ ΚΑΘΗΛΩΣΗΣ INVISx
INVISx CRANIALIS RÖGZÍTŐRENDSZER
SISTEMA DI FISSAGGIO CRANICO INVISx
INVISx KRANIEFIKSERINGSSYSTEM
ZESTAW DO OSTEOSYNTAZY KOŚCI
POKRYWY CZASZKI INVISx
SISTEMA DE FIXAÇÃO CRANIANA INVISx
SISTEMA DE FIJACIÓN CRANEAL INVISx
INVISx KRANIELLT FIXATIONSSYSTEM
INVISx KRANIYAL FIKSASYON SISTEMI

Medtronic

Conformité européenne (European Conformity). This symbol means that the device fully complies with European Directive 93/42/EEC.

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EN

Description

The INVISx™ Cranial Fixation System is designed to assure the indwelling structural stability of the replaced skull flap following craniotomy. The complete System consists of non-ferromagnetic polymer Locks and two reusable Instruments (Fig. 1). The nonmetallic design of the Locks will not interfere with CT or MR imaging. The INVISx Locks are packaged **sterile** and **non-pyrogenic**. The reusable Instruments are packaged **non-sterile**.

Indications

The INVISx Cranial Fixation System is intended for use in re-fixation of cranial bone flaps following a craniotomy.

Instructions For Use

WARNING: The INVISx Locks are for single use only. DO NOT RESTERILIZE OR REUSE.

Prior to using the INVISx Cranial Fixation System, it is necessary for the attending physician and other responsible personnel to familiarize themselves with the use and function of the various components.

Instruments

The Instruments should be carefully examined prior to use for damage and functionality. Corrosion (freezing-up), scratching, loosening, bending and/or fracture of any or all sections of the Instruments may inhibit or prevent proper function. **Do not use if Instrument is damaged or non-functioning.**

Clean and sterilize the INVISx Tensioning Instrument and the INVISx Trimming Instrument as indicated below or in accordance with the standard procedure of the institution.

Cleaning

Medtronic Neurosurgery INVISx Instruments are provided **NON-STERILE** and must be cleaned and sterilized before each use. Manual cleaning is preferred for all INVISx Instruments. Remove all packaging material prior to sterilization.

Immediately following surgery, rinse the Instruments under warm running water to remove all blood, body fluids, and tissue. To clean, use enzymatic detergent of neutral pH diluted according to the manufacturer's instructions, and soak for at least five (5) minutes. Gently scrub the Instruments with a stiff plastic brush. Do not use wire brushes or steel wool. Rinse device in warm water for at least two minutes. Make sure all Instrument surfaces are clean and free of tissue and stains. Package Instruments in material suitable for steam autoclave processing, or in the original storage/sterilization tray.

Sterilization

Steam Autoclaving

Steam autoclaving is the recommended method of sterilization for the INVISx Instruments.

Follow the autoclave manufacturer's instructions. Do not overload the autoclave chamber, as pockets may form, preventing full steam penetration. Place a towel on the bottom of the autoclave pan to absorb excess water during autoclaving.

CAUTION: ONLY USE TOWELS THAT HAVE NO RESIDUAL DETERGENT AND THAT HAVE BEEN LAUNDERED IN NEUTRAL WATER.

The following recommended cycle parameters are valid for processing the INVISx Instruments.

METHOD	CYCLE	TEMPERATURE	EXPOSURE TIME
Steam	Gravity	121 °C to 123 °C (250 °F to 254 °F)	30 minutes
Steam	Gravity	134 °C to 137 °C (273 °F to 278 °F)	10 minutes
Steam	Pre-Vacuum	134 °C to 137 °C (273 °F to 278 °F)	4 minutes

For instruments contaminated with Transmissible Spongiform Encephalopathies (TSE), the following cycle is referenced in NHS Estates HTM 2010 parts 4 & 6: Appendix 2, "Items Contaminated With TSE Agents and WHO Infection Control Guidelines for Transmissible Spongiform Encephalopathies."

METHOD	CYCLE	TEMPERATURE	EXPOSURE TIME
Steam	Pre-Vacuum	134 °C to 137 °C (273 °F to 278 °F)	Single cycle of 18 minutes or six repeated 3-minute cycles

Store the INVISx Instruments in their original packaging. Store in a clean and dry environment.

Surgical Technique

NOTE: For best results, use a high-speed craniotome (such as MIDAS REX) to perform the craniotomy.

NOTE: Use of a dissecting tool smaller than Ø1.8 mm may prevent use of the Locks.

NOTE: Creating a kerf perpendicular to the surface of the cranium will aid in the ease of placement and stability of the Locks. Locks are not suitable for skull thicknesses less than 2.5 mm.

Lock Placement

INVISx Locks are for single use only and cannot be reused.

Correct handling of the INVISx Locks before and during the operation is important for the success of the cranial flap fixation.

INVISx Locks should be stored in their original packaging until use. Check the expiration date and the integrity of the sterile packaging prior to use. Do not use if package has been previously opened or damaged.

CAUTION: Check Locks for damage. Do not use if bent or damaged.

Proper placement of INVISx Locks (each composed of a Cap and Stem) is the responsibility of the surgeon.

1. Using sterile technique, remove the INVISx Locks from the opened trays. Check to ensure that all components are assembled.
2. It is recommended that dural tack-up sutures be applied every two to three centimeters. Care should be taken to avoid placing a suture where an INVISx Lock will be positioned for closing.

NOTE: Do not place INVISx Locks in the burr hole.

3. A minimum of four INVISx Locks is recommended for secure closure of a craniotomy.
4. With the flap removed, position the INVISx Locks equidistantly around the opening (Fig. 2), taking care to avoid any tack-up sutures. The Stem base should lie between dura and skull (Fig. 3).
5. Replace the flap, ensuring that the Stems remain vertical within the kerf (Fig. 4).
6. Place Tensioning Instrument on one of the Stems as shown in Fig. 5, and tighten just until the Cap touches the skull. Repeat tensioning procedure on each Stem in a diagonal sequence (Fig. 6) until all caps are just touching the skull.

NOTE: Do not proceed with the tensioning procedure until all of the caps are in place and touching the skull.

7. Complete the tensioning procedure by tightening each cap until the Stem safety tab breaks off. Tighten caps in the same diagonal sequence.

NOTE: Each Stem is equipped with a safety tab, which is designed to break off when the Lock has reached the optimum degree of tension needed for cranial flap fixation.

8. Using the Trimming Instrument, cut each Stem flush with the top of the Cap (Fig. 7).

NOTE: After cutting the Stem, a scalpel may be used to trim off any excess or sharp edges that remain. Use caution not to trim Stem below level of Cap.

Lock Removal

Supplies needed:

- Trimming Instrument or Surgical Clamp
- Safety Glasses or Face Shield

1. Using Trimming Instrument or Surgical Clamp, pry Cap off of Stem (Fig. 8).

CAUTION: Removing a Cap under tension may result in its ejection from the surgical site.

2. Repeat on each Lock until flap may be removed.

How Supplied

INVISx Locks are **sterile** and **non-pyrogenic** in a double sterile barrier package. They are intended for **single** (one-time) **use only**. **Do not Resterilize**. Do not use if package has been previously opened or damaged. Resterilization can damage the product, potentially leading to patient injury.

Medtronic Neurosurgery is not responsible for the performance of any product that has been resterilized.

INVISx Instruments are packaged **NON-Sterile** and must be cleaned and sterilized prior to use.

Patient Education

It is the physician's responsibility to educate the patients and/or their representative(s) regarding fixation of a cranial flap following a craniotomy. This should include a description of associated warnings, precautions, complications, and an explanation of potential alternative products and treatments.

Contraindications

Contraindications include, but are not limited to patients with: systemic diseases, metabolic disturbances, sensitivity or allergies to implant materials, active or latent intracranial infections, bone conditions which exclude the application of INVISx Locks, artificial cranial bone flaps, bone tumors in the region of the implant, degenerative bone disease, absence of dura mater, decompressed flap, some orbital or cranial base regions of the skull, skull thicknesses less than 2.5 mm, and children whose skulls are immature.

Warnings and Precautions

Improper fixation of INVISx Locks may result in formation of cranial ridges; changes in the position of the cranium and loosening and/or fracture of implant components. The use of a minimum of four INVISx Locks placed equidistantly is recommended for maximum stability. Refer to the section on Stem Placement. **Do not place INVISx Locks over burr holes. Place only in kerf.**

It is advised that the surgeon be thoroughly familiar with the surgical techniques relative to the use of this device. The surgeon is responsible for proper application of the INVISx Cranial Fixation System. Complications due to incorrect diagnosis and operational technique are the responsibility of the surgeon.

Complications

General risks may include infections, nerve damage, hematomas and impairment of wound healing. The patient should be made aware of these risks and others associated with neurosurgery, general surgery and the use of anesthesia.

Returned Goods Policy

Products must be returned in unopened packages, with manufacturer's seals intact, to be accepted for replacement or credit, unless returned due to a complaint of product defect or mislabeling. Determination of a product defect or mislabeling will be made by Medtronic Neurosurgery, which determination will be final. Products will not be accepted for replacement or credit if they have been in possession of the customer for more than 90 days.


Warranty

A. Standard Limited Warranty. Medtronic Neurosurgery warrants to the original end user purchaser ("Purchaser") that the enclosed single use or reusable product ("Product") purchased by Purchaser, at the time of delivery to Purchaser, shall be substantially free from defects in material and workmanship. Medtronic Neurosurgery makes no warranty (express, implied or statutory) for Products that are modified (except as expressly contemplated herein) or subjected to unusual physical stress, misuse, improper operation, neglect, improper testing, use in combination with other products or components other than those for which the Products were designed, or use in any manner or medical procedure for which the Products are not indicated.







B. Remedy. Purchaser's exclusive remedy and Medtronic Neurosurgery's sole liability for breach of the foregoing warranty shall be, at Medtronic Neurosurgery's sole option and election, to replace the Product or credit Purchaser for the net amount actually paid for any such Product; provided that (i) Medtronic Neurosurgery is notified in writing within ninety (90) days after Purchaser's receipt of the Product that such Product failed to conform, including a detailed

explanation in English of any alleged nonconformity; (ii) such Product is returned to Medtronic Neurosurgery within ninety (90) days after Purchaser's receipt of the Product F.O.B. 125 Cremona Drive, Goleta, California 93117, U.S.A. or as otherwise designated by Medtronic Neurosurgery; and (iii) Medtronic Neurosurgery is reasonably satisfied that the claimed nonconformities actually exist. Except as expressly provided in this paragraph, Purchaser shall not have the right to return Products to Medtronic Neurosurgery without Medtronic Neurosurgery's prior written consent.

C. Exclusion of Other Warranties. EXCEPT FOR THE LIMITED WARRANTY PROVIDED IN (A) ABOVE, MEDTRONIC NEUROSURGERY GRANTS NO OTHER WARRANTIES OR CONDITIONS, EXPRESS OR IMPLIED AND MANUFACTURER SPECIFICALLY DISCLAIMS THE IMPLIED WARRANTIES AND CONDITIONS OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC NEUROSURGERY NEITHER ASSUMES NOR AUTHORIZES ANY OTHER PERSON TO ASSUME ANY OTHER LIABILITIES ARISING OUT OF OR IN CONNECTION WITH THE SALE OR USE OF ANY PRODUCT.

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8

Appendix to Instruction for Use INVISx Cranial Fixation System


INVISx Cranial Fixation System products:

Product	Description	Product Number
INVISx Lock	2 INVISx Locks per box	60100
	1 INVISx Lock per box	60101
	Shelf Pack: 10 boxes of Product Number 60100	60110
Tensioning Instrument	Stainless steel instrument designed to advance cap down onto skull	60710
Trimming Instrument	Stainless steel trimming instrument designed to trim stem flush with the top of the superior cap	60720
Instrument Tray	Instrument Tray	60730
INVISx Cranial Fixation System Instrument Set	Set includes Product Numbers: 60710, 60720, 60730	60700

1. Packaging	Package description for INVISx Locks and Tools: INVISx Locks The INVISx Locks system (CFNs 60100, 60101, and 60110) are packaged in a double cover Tyvek tray and then packaged in a labeled corrugated carton with label and brand name. Each device for the INVISx System is packaged in packaging either a 1 or 2 pieces in pack or in packaging for storage, which contains 10 boxes. IFUs are included inside each carton. The system of fixing elements INVISx Locks are all sterilized via EO; shelf life - 5 years. INVISx Instruments The INVISx Instruments (CFNs 60700, 60710, 60720) are packaged in bubble wrap and then packaged in a labeled corrugated carton along with the IFU. INVISx instruments are non-sterile reusable device that do not have a shelf life.																														
2. Type of contact with the body	The INVISx Locks are implantable devices which have tissue and bone contact for over 30 days (permanent contact).																														
3. Components and materials of devices	<table border="1"> <thead> <tr> <th colspan="2">INVISx Tensioning Tool (CFN 60710)</th> </tr> <tr> <th>Product Part</th> <th>Material</th> </tr> </thead> <tbody> <tr> <td>Entire Tool</td> <td>17-4 Stainless Steel: 17-4 PH Cond A ASTM A564, Heat Treat Cond H900</td> </tr> <tr> <th colspan="2">INVISx Trimmer (CFN 60720)</th> </tr> <tr> <th>Product Part</th> <th>Material</th> </tr> <tr> <td>Cutting Edge</td> <td>Stainless Steel WN-1.4125</td> </tr> <tr> <td>Handle of Trimmer</td> <td>Stainless Steel WN-1.4021</td> </tr> <tr> <th colspan="2">INVISx Locks (CFN 60100, 60101, 60110)</th> </tr> <tr> <th>Product Part</th> <th>Material</th> </tr> <tr> <td>Stem</td> <td>6331 Polypropylene Homopolymer Resin, Virgin Grade</td> </tr> <tr> <td>Round Cap</td> <td>Delrin Resin 500P Natural</td> </tr> <tr> <th colspan="2">Instrument Tray INVISx (CFN 60730)</th> </tr> <tr> <td>Bottom Tray</td> <td>Radel® R-5100 Series Resin Polyphenylsulfone</td> </tr> <tr> <td>Headblock (cover of the tray)</td> <td>Radel® R-5000 Series Resin polyphenylsulfone</td> </tr> <tr> <td>Fasteners of tray body</td> <td>304 Stainless Steel</td> </tr> </tbody> </table>	INVISx Tensioning Tool (CFN 60710)		Product Part	Material	Entire Tool	17-4 Stainless Steel: 17-4 PH Cond A ASTM A564, Heat Treat Cond H900	INVISx Trimmer (CFN 60720)		Product Part	Material	Cutting Edge	Stainless Steel WN-1.4125	Handle of Trimmer	Stainless Steel WN-1.4021	INVISx Locks (CFN 60100, 60101, 60110)		Product Part	Material	Stem	6331 Polypropylene Homopolymer Resin, Virgin Grade	Round Cap	Delrin Resin 500P Natural	Instrument Tray INVISx (CFN 60730)		Bottom Tray	Radel® R-5100 Series Resin Polyphenylsulfone	Headblock (cover of the tray)	Radel® R-5000 Series Resin polyphenylsulfone	Fasteners of tray body	304 Stainless Steel
INVISx Tensioning Tool (CFN 60710)																															
Product Part	Material																														
Entire Tool	17-4 Stainless Steel: 17-4 PH Cond A ASTM A564, Heat Treat Cond H900																														
INVISx Trimmer (CFN 60720)																															
Product Part	Material																														
Cutting Edge	Stainless Steel WN-1.4125																														
Handle of Trimmer	Stainless Steel WN-1.4021																														
INVISx Locks (CFN 60100, 60101, 60110)																															
Product Part	Material																														
Stem	6331 Polypropylene Homopolymer Resin, Virgin Grade																														
Round Cap	Delrin Resin 500P Natural																														
Instrument Tray INVISx (CFN 60730)																															
Bottom Tray	Radel® R-5100 Series Resin Polyphenylsulfone																														
Headblock (cover of the tray)	Radel® R-5000 Series Resin polyphenylsulfone																														
Fasteners of tray body	304 Stainless Steel																														

For Russia

4. Specifications

5. Shelf life	For sterile products: 5 years. For non-sterile: unlimited.
6. Utilization	Utilize (dispose) of all single-use products in accordance to local environmental requirements.
7. Operating conditions	The product can be used under the standard conditions. Operation room temperature range from 15 °C to 30 °C. During operation at nominal temperature values from 32°C to 42°C, humidity range from 80% to 100% (+2%).

For Russia

8. Transportation conditions	Recommended transit temperature -35 °C to 58 °C (-31 °F to 136 °F); 85% relative humidity (non-condensing)
9. Storage conditions	Recommended storage temperature Average of 22 °C (71 °F) with a standard deviation of 2.7 °C (4.86 °F) and actionable control limits of 15 °C to 30 °C (59 °F to 86 °F)
10. Customer Service	If you have any INVISx devices returns or questions please contact with Russia Medtronic Customer Service at tel.: +7 (495) 580-7377.
	Authorized Representative in Russia: Medtronic LLC Russia, 123317 Moscow, Presnenskaya Naberezhnaya, 10, Tel.: +7 495 580 7377; +7 495 580 7377 FREE Fax: + 7 495 580 7378 E-mail: info.russia@medtronic.ru Internet: www.medtronic.ru

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INVISIX® is a registered trademark of Medtronic, Inc.