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
INSTRUCTIONS FOR USE / BRUGSANVISNING / GEBRUIKSAANWIJZING / MODE D'EMPLOI / GEBRAUCHSANWEISUNG /
ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ / HASZNÁLATI UTASÍTÁS / ISTRUZIONI PER L'USO / BRUGSANVISNING / INSTRUKCJA UŻYTKOWANIA /
INSTRUÇÕES PARA USO / INSTRUCCIONES DE USO / BRUGSANVISNING / KULLANMA TALIMATI

CLIP GUN
CLIPSEPISTOL
CLIP-PISTOOL
PISTOLET AGRAFEUR
KLAMMERPISTOLE
ΠΙΣΤΟΛΙ ΑΓΚΤΗΡΩΝ
КАРОЦСПИЗТОЛ
PISTOLA PER CLIP
KLIPSPISTOL
APLIKATOR DO ZAKŁADANIA KLIPSÓW
PISTOLA DE CLIPES
PISTOLA DE GRAPAS
CLIPSPISTOL
KLIPS TABANCASI

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
Clip Gun

Description

The single-use disposable Clip Gun Kit contains a Clip Gun, pre-loaded magazines with scalp clips, and removal forceps. The Raney-type scalp clips are used to provide hemostasis by preventing arterial bleeding from small vessels of the skin and deeper vessels under the scalp. Scalp clips apply force to the scalp through the spring action built into the clip explicitly as a spring and implicitly as a result of the pliable plastic clip composition. Each clip controls bleeding locally along the segment to which it is applied; therefore, the number of clips required for hemostasis is dictated by incision length. Clips are placed in tandem along the incision edge, with no more than a 1 cm gap between clips. Scalp clips are green to enhance visibility and are designed with a serrated edge and a round "hinge" area.

The Clip Gun, magazine with clips, and removal forceps are packaged **sterile** and **non-pyrogenic** and are for **single (one-time) use only**.

Indications

The Clip Gun Raney-type clips provide hemostasis during cranial surgery by preventing arterial bleeding from small vessels of the skin and deeper vessels under the scalp.






Instructions for Use


WARNING: The Clip Gun Kit is designed for single use only. DO NOT RESTERILIZE OR REUSE.

Note: Do not remove clip retaining pin until Step 3.

To use Clip Gun:

- Align magazine rails to gun slots (Fig. 1). Slide the magazine toward the back of the Clip Gun between the gun ears.
- Be sure rear catch is secured in notch on back of gun (Fig. 2).
Note: Do not release rear catch until magazine is completely empty.
- Remove clip retaining pin (Fig. 3).
Note: Once pin is removed, do not remove magazine until completely empty.
- Squeeze trigger firmly until stop occurs, without releasing trigger (Fig. 4).
- Push clip **FULLY** onto the scalp until the "hinge" edge of the clip is touching the scalp edge (Fig. 5).
- Squeeze trigger harder to overcome stop and release the clip (Fig. 6).
- Release trigger to allow jaws to completely retract (Fig. 7). Another clip will self-load into the jaws.
- To use remaining clips, repeat Steps 4-7 until magazine is empty. Magazine may now be released and removed from the gun.

Note: Do not release the magazine until it is completely empty. Failure to completely empty the magazine before removal will cause clips to eject from the gun.

9. To remove empty magazine, lift rear catch and slide magazine forward (Fig. 8).

How Supplied


The Clip Gun, Clip Gun magazine with scalp clips, and removal forceps are **sterile** and **non-pyrogenic**. They are intended for **single (one-time) use only**. Do not Resterilize.

Do not use if package has been previously opened

or damaged. Resterilization can damage the product, potentially leading to patient injury. Medtronic Neurosurgery is not responsible for the performance of any product that has been resterilized.

Patient Education

It is the physician's responsibility to educate the patients and/or their representative(s) regarding the use of these products during surgical procedures. This should include a description of associated warnings, precautions, complications, and an explanation of potential alternative products and treatments.

Contraindications

Contraindications include, but are not limited to patients with: systemic diseases, metabolic disturbances, sensitivity or allergies to scalp clip materials, active or latent infections, and scalp conditions which exclude the application of scalp clips.

Warnings and Precautions

It is advised that the surgeon be thoroughly familiar with the surgical techniques relative to the use of this device. The surgeon is responsible for proper application of scalp clips. Complications due to incorrect diagnosis and operational technique are the responsibility of the surgeon.

Dispose of Clip Gun system components in accordance with local regulations.

Complications



General risks may include infections and necrosis. The patient should be made aware of these risks and others associated with neurosurgery, general surgery and the use of anesthesia.

Returned Goods Policy

Products must be returned in unopened packages, with manufacturer's seals intact, to be accepted for replacement or credit, unless returned due to a complaint of product defect or mislabeling. Determination of a product defect or mislabeling will be made by Medtronic Neurosurgery, which determination will be final. Products will not be accepted for replacement or credit if they have been in possession of the customer for more than 90 days.

Warranty

A. Standard Limited Warranty. Medtronic Neurosurgery warrants to the original end user purchaser ("Purchaser") that the enclosed single use or reusable product ("Product") purchased by Purchaser, at the time of delivery to Purchaser, shall be substantially free from defects in material and workmanship. Medtronic Neurosurgery makes no warranty (express, implied or

statutory) for Products that are modified (except as expressly contemplated herein) or subjected to unusual physical stress, misuse, improper operation, neglect, improper testing, use in combination with other products or components other than those for which the Products were designed, or use in any manner or medical procedure for which the Products are not indicated.

B. Remedy. Purchaser's exclusive remedy and Medtronic Neurosurgery's sole liability for breach of the foregoing warranty shall be, at Medtronic Neurosurgery's sole option and election, to replace the Product or credit Purchaser for the net amount actually paid for any such Product; provided that (i) Medtronic Neurosurgery is notified in writing within ninety (90) days after Purchaser's receipt of the Product that such Product failed to conform, including a detailed explanation in English of any alleged nonconformity; (ii) such Product is returned to Medtronic Neurosurgery within ninety (90) days after Purchaser's receipt of the Product F.O.B. 125 Cremona Drive, Goleta, California 93117, U.S.A. or as otherwise designated by Medtronic Neurosurgery; and (iii) Medtronic Neurosurgery is reasonably satisfied that the claimed nonconformities actually exist. Except as expressly provided in this paragraph, Purchaser shall not have the right to return Products to Medtronic Neurosurgery without Medtronic Neurosurgery's prior written consent.

C. Exclusion of Other Warranties. EXCEPT FOR THE LIMITED WARRANTY PROVIDED IN (A) ABOVE, MEDTRONIC NEUROSURGERY GRANTS NO OTHER WARRANTIES OR CONDITIONS, EXPRESS OR IMPLIED AND MANUFACTURER SPECIFICALLY DISCLAIMS THE IMPLIED WARRANTIES AND CONDITIONS OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC NEUROSURGERY NEITHER ASSUMES NOR AUTHORIZES ANY OTHER PERSON TO ASSUME ANY OTHER LIABILITIES ARISING OUT OF OR IN CONNECTION WITH THE SALE OR USE OF ANY PRODUCT.

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Authorized EC-representative:

Medtronic B.V.
Earl Bakkenstraat 10
6422 PJ Heerlen
The Netherlands

Australian Sponsor:
Medtronic Australiasia Pty Ltd
97 Waterloo Rd
North Ryde, NSW 2113 Australia

Medtronic, Inc.
710 Medtronic Parkway NE
Minneapolis, MN 55432-5604 USA

Medtronic

Design & Manufacturing Facility
Medtronic Neurosurgery
125 Cremona Drive
Goleta, California 93117-5500 USA
(800) 468-9710 USA/Canada
(901) 344-0645 International
(800) 468-9713 FAX USA/Canada
(901) 396-2698 FAX International

15802-1F 26026
© Medtronic, Inc. 2009
All Rights Reserved
Printed in USA

Appendix to Instruction for Use
Clip Gun

Table 1: Clip Gun and Associated Component Devices


Device	Catalog Number
Clip Gun (10 pack)	CG8901A
Clip Gun Kit (1 pack)	CG8900
Clip Gun Kit (10 pack)	CG8901
Magazine, 5 pack	CM8902
Clips (10 each per pack), 20 packs	CL8903

1. Packaging	Package description The Clip Gun products are packaged on a pouch for the gun kit made from a bundle in double polyurethane bags within cardboard box to prevent contamination or damage.																																																																																														
2. Type of contact with the body	The green clips are considered a surface device with contact to breached or compromised surfaces. The contact duration for the green clips is considered limited (<24hrs.) The green clips are used with Clip Guns during cranial surgical procedures. During the surgical operation the green clips are applied to the edge of the incision scalp to help control the bleeding. Once the operation is complete the clips are removed and discarded. The green clips are the only part of the device that has any patient contact.																																																																																														
3. Components and materials of devices	<table><tr><th>Item</th><th>Description</th><th>Mold / Purchased</th><th>Material</th></tr><tr><td>CG8900</td><td>SINGLE CG PACKAGED KIT</td><td>Purchase</td><td>N/A</td></tr><tr><td>CG8901</td><td>10PK CG PACKAGED KIT</td><td>Purchase</td><td>N/A</td></tr><tr><td>CG8901A</td><td>10PK CG GUN ONLY</td><td>Purchase</td><td>N/A</td></tr><tr><td>CM8902</td><td>5PK MAGAZINE W/CLIPS</td><td>Purchase</td><td>N/A</td></tr><tr><td>CL8903</td><td>20-10PK CLIPS ONLY</td><td>Purchase</td><td>N/A</td></tr><tr><td colspan="4">CLIP GUN ASSEMBLY</td></tr><tr><td>17235-1</td><td>LOCATION ARM -RIGHT-CLIP-GUN</td><td>Mold</td><td rowspan="5">Material: - Lexan 10% Glass Filled</td></tr><tr><td>17234-1</td><td>LOCATION ARM -LEFT-CLIP-GUN</td><td>Mold</td></tr><tr><td>17217-1</td><td>TRIGGER SEAR, RIGHT SIDE</td><td>Mold</td></tr><tr><td>17218-1</td><td>TRIGGER SEAR, LEFT SIDE</td><td>Mold</td></tr><tr><td>17175-1</td><td>RAMP CLIP GUN</td><td>Mold</td></tr><tr><td>14929-1</td><td>DOOR, GUN</td><td>Mold</td><td rowspan="4">Material: -Plastic CYCOLAC Resin ABS MG47</td></tr><tr><td>14932-1</td><td>TRIGGER</td><td>Mold</td></tr><tr><td>14928-1</td><td>BODY, HALF, LEFT</td><td>Mold</td></tr><tr><td>15875-1</td><td>BODY, HALF, RIGHT, CE</td><td>Mold</td></tr><tr><td>14944-1</td><td>TRIGGER STOP</td><td>Mold</td><td>Material: - Santoprene Rubber</td></tr><tr><td>14933-1</td><td>TRIGGER SPRING</td><td>Purchase</td><td>302 SS Spring Wire</td></tr><tr><td>14930-1</td><td>DOOR SPRING</td><td>Purchase</td><td>A-228 Music Wire</td></tr><tr><td colspan="4">INNER JAW ASSEMBLY</td></tr><tr><td>15544-1</td><td colspan="2">INNER JAW ASSEMBLY</td><td></td></tr><tr><td>17214-1</td><td>LOWER JAW</td><td>Mold</td><td rowspan="2">Material: - PPS SF-40CF</td></tr><tr><td>17176-1</td><td>UPPER JAW</td><td>Mold</td></tr><tr><td>N/A</td><td>DOWEL PIN, Ø.125 X .625"</td><td>Purchase</td><td>18-8 SS</td></tr><tr><td>17683-1</td><td>CLIPS, GREEN</td><td>Mold</td><td>Material(s): - Nylon #6</td></tr></table>			Item	Description	Mold / Purchased	Material	CG8900	SINGLE CG PACKAGED KIT	Purchase	N/A	CG8901	10PK CG PACKAGED KIT	Purchase	N/A	CG8901A	10PK CG GUN ONLY	Purchase	N/A	CM8902	5PK MAGAZINE W/CLIPS	Purchase	N/A	CL8903	20-10PK CLIPS ONLY	Purchase	N/A	CLIP GUN ASSEMBLY				17235-1	LOCATION ARM -RIGHT-CLIP-GUN	Mold	Material: - Lexan 10% Glass Filled	17234-1	LOCATION ARM -LEFT-CLIP-GUN	Mold	17217-1	TRIGGER SEAR, RIGHT SIDE	Mold	17218-1	TRIGGER SEAR, LEFT SIDE	Mold	17175-1	RAMP CLIP GUN	Mold	14929-1	DOOR, GUN	Mold	Material: -Plastic CYCOLAC Resin ABS MG47	14932-1	TRIGGER	Mold	14928-1	BODY, HALF, LEFT	Mold	15875-1	BODY, HALF, RIGHT, CE	Mold	14944-1	TRIGGER STOP	Mold	Material: - Santoprene Rubber	14933-1	TRIGGER SPRING	Purchase	302 SS Spring Wire	14930-1	DOOR SPRING	Purchase	A-228 Music Wire	INNER JAW ASSEMBLY				15544-1	INNER JAW ASSEMBLY			17214-1	LOWER JAW	Mold	Material: - PPS SF-40CF	17176-1	UPPER JAW	Mold	N/A	DOWEL PIN, Ø.125 X .625"	Purchase	18-8 SS	17683-1	CLIPS, GREEN	Mold	Material(s): - Nylon #6
Item	Description	Mold / Purchased	Material																																																																																												
CG8900	SINGLE CG PACKAGED KIT	Purchase	N/A																																																																																												
CG8901	10PK CG PACKAGED KIT	Purchase	N/A																																																																																												
CG8901A	10PK CG GUN ONLY	Purchase	N/A																																																																																												
CM8902	5PK MAGAZINE W/CLIPS	Purchase	N/A																																																																																												
CL8903	20-10PK CLIPS ONLY	Purchase	N/A																																																																																												
CLIP GUN ASSEMBLY																																																																																															
17235-1	LOCATION ARM -RIGHT-CLIP-GUN	Mold	Material: - Lexan 10% Glass Filled																																																																																												
17234-1	LOCATION ARM -LEFT-CLIP-GUN	Mold																																																																																													
17217-1	TRIGGER SEAR, RIGHT SIDE	Mold																																																																																													
17218-1	TRIGGER SEAR, LEFT SIDE	Mold																																																																																													
17175-1	RAMP CLIP GUN	Mold																																																																																													
14929-1	DOOR, GUN	Mold	Material: -Plastic CYCOLAC Resin ABS MG47																																																																																												
14932-1	TRIGGER	Mold																																																																																													
14928-1	BODY, HALF, LEFT	Mold																																																																																													
15875-1	BODY, HALF, RIGHT, CE	Mold																																																																																													
14944-1	TRIGGER STOP	Mold	Material: - Santoprene Rubber																																																																																												
14933-1	TRIGGER SPRING	Purchase	302 SS Spring Wire																																																																																												
14930-1	DOOR SPRING	Purchase	A-228 Music Wire																																																																																												
INNER JAW ASSEMBLY																																																																																															
15544-1	INNER JAW ASSEMBLY																																																																																														
17214-1	LOWER JAW	Mold	Material: - PPS SF-40CF																																																																																												
17176-1	UPPER JAW	Mold																																																																																													
N/A	DOWEL PIN, Ø.125 X .625"	Purchase	18-8 SS																																																																																												
17683-1	CLIPS, GREEN	Mold	Material(s): - Nylon #6																																																																																												



FORCEP ASSEMBLY			
14948-1	FORCEP ASSEMBLY		
14946-1	TOP ARM, FORCEPS	Mold	Material: - Zelux 30% Glass Filled Polycarbonate (PC)
14945-1	BOTTOM ARM, FORCEPS	Mold	
14947-1	DOWEL PIN	Purchase	304 SS
MAGAZINE ASSEMBLY			
14942-1	MAGAZINE ASSEMBLY		
N/A	LEFT SIDE, MAGAZINE	Mold	Material: - Terluc transparent ABS polymers SDM-050
N/A	RIGHT SIDE, MAGAZINE	Mold	
14940-1	LOCK, CLIP	Mold	Material: - Polycarbonate, Lexan
N/A	CATCH LOCK	Mold	
N/A	SPRING GUIDE	Mold	
N/A	MAGAZINE SPRING	Purchase	
			301 SS





4. Specifications

Specification (nominal) for Clip Gun.		
The clips are used to provide hemostasis by preventing arterial bleeding from small vessels of the skin and deeper vessels under the scalp. The clips apply force through the spring action built into the clip explicitly as a spring and implicitly as a result of the pliable plastic clip composition.		
Ergonomics/Human Factors • Clip Gun is intuitive and ambidextrous • Designed to resemble other currently available surgical instruments • Device does not require repositioning of the surgeon or patient for use • Moderate hand activation is required to use device		
Physical Dimensional • Opening width – minimum 0.25" (6,35 mm) • Clip length – maximum 0.75" (19,05 mm)		
Catalog Number	Description	Mass, no more then
CG8900	Clip Gun + 3 Magazines + Forceps	528.4 g
CG8901	CG8900 x10	528.4 g
CG8901A	Single Clip Gun	382.1 g
CM8902	3 Magazines	272.5 g
CL8903	20 pouches of 10 clips each	282.1 g
Specifications drawing (nominal): Clip gun Body Half Left		

Body Half Right

	<table><tr><th colspan="3">COLORANT CHART</th></tr><tr><th>PIGMENT</th><th>% USED</th><th>MANUFACTURER</th></tr><tr><td>WHITE 6</td><td>0.3084</td><td>CRISTAL</td></tr><tr><td>RED 7</td><td>0.0705</td><td>JECO</td></tr><tr><td>YELLOW 53</td><td>1.5419</td><td>HEMOTEC</td></tr><tr><td>BARIUM SULFATE</td><td>3.00</td><td>SACHTLESEN CHROME BLANC</td></tr><tr><td>IR-6-441</td><td>0.2202</td><td>ACRAWAX C</td></tr></table> 2X EJECTOR PIN PAD	COLORANT CHART			PIGMENT	% USED	MANUFACTURER	WHITE 6	0.3084	CRISTAL	RED 7	0.0705	JECO	YELLOW 53	1.5419	HEMOTEC	BARIUM SULFATE	3.00	SACHTLESEN CHROME BLANC	IR-6-441	0.2202	ACRAWAX C
COLORANT CHART																						
PIGMENT	% USED	MANUFACTURER																				
WHITE 6	0.3084	CRISTAL																				
RED 7	0.0705	JECO																				
YELLOW 53	1.5419	HEMOTEC																				
BARIUM SULFATE	3.00	SACHTLESEN CHROME BLANC																				
IR-6-441	0.2202	ACRAWAX C																				
5. Shelf life	The shelf life is 2 years and 6 months																					
6. Utilization	Dispose of all single-use accessories according to local environmental requirements.																					
7. Operating conditions	15 °C to 30 °C (59 °F to 86 °F) at altitudes less than 2400 meters (8000 feet) above sea level																					
8. Transportation conditions	Recommended transit temperature -35 °C to 58 °C (-31 °F to 136 F); 85% relative humidity (non-condensing)																					
9. Storage conditions	Recommended storage temperature Average of 22 °C (71 °F) with a standard deviation of 2.7 °C (4.86 °F) and actionable control limits of 15 °C to 30 °C (59 °F to 86 °F)																					
10. Customer Service	Customer Service Medtronic Neurosurgery 125 Cremona Drive Goleta, California 93117-5500 USA (800) 468-9710 USA/Canada (901) 344-0645 International (800) 468-9713 FAX USA/Canada (901) 396-2698 FAX International Representative company of Medtronic Inc. in Moscow																					
	Authorized Representative in Russia: Medtronic LLC Russia, 123317 Moscow, Presnenskaya Naberezhnaya, 10, Tel.: +7 495 580 7377; +7 495 580 7377 FREE Fax: + 7 495 580 7378 E-mail: info.russia@medtronic.ru Internet: www.medtronic.ru																					

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PS Medical® is a registered trademark of Medtronic, Inc.