INSTRUCTION FOR USE

EXTENTION LINE - HIGH PRESSURE

MATERIALS USED:

- Poly Vinyl Chloride (PVC), Polypropylene, Polycarbonate

INDICATIONS:

- Infusion of fluid such as Parenteral nutrition and administration of other drugs and medicines.
- Maintaining hydration and/or correct dehydration in patients who are unable to take sufficient volume of oral fluids.
- Suitable for Pressure up to 54 bar (800 PSI)

CONTRA INDICATIONS:

- Administration of highly viscous fluid
- Not to be used in patients with known hypersensitivity to any of the materials used

DIRECTION FOR USE:

- Inspect product package for integrity and expiry and then remove product from package
- Connect the Female luer connector of Extension Tube to the male fitting of Infusion Set.
- Remove the protective cap of Male luer connector and open roller clamp of I. V. Sets.
- Expel the air from PVC tube and close roller clamp.
- Connect the male luer connector to the Injection site (I.V. Cannula etc.) And regulate the flow by adjusting roller clamp to achieve desired flow rate

WARNINGS:

- Read Instructions before use
- This product should be used only by a qualified Doctor or a Paramedic
- Check the Product and Package integrity prior to use and also the expiry date
- The product is sterile, non-pyrogenic and non-toxic if package is not opened or damaged
- For single use only. Do not re-cycle, clean or re-sterile. Discard the product after use as per Medical disposal regulation in your country.
- Use the product immediately after opening the individual pack
- Store under dry and cool condition. Do not store in direct sunlight, at extreme temperature or in high humidity
- The product should be used according to the instructions for use
- This product have been evaluated against Potential Risks & after proposed Risk Control Measures, quality for design & manufacturing. Any residual risks are deemed to be outweighed by the clinical benefits.
- This product is for short term use only.

CE
0120

: See Instruction for use

Ref./Art-No.

: Item Code / Catalogue No.

: For Single Use Only

Lot / Ch-B

: Lot / Batch Number

MFG. DT.: Date of Manufacturing

EXP. DT.: Use by / Expiry date

STERILE EO

: Sterilized by Ethylene Oxide

: Non Pyrogenic

: DEHP

: Latex free

: Do not re-sterile

SURU™ is the Trade Mark of SURU INTERNATIONAL PVT. LTD., INDIA

Manufactured by:
SURU INTERNATIONAL PVT. LTD.

Plot # 3, Panchsheel Centre, Gurgaon, Haryana 122 002, INDIA
www.suru.com • Tel: ISO 9001 2008 ISO 13485 2003 Certified Company
Customer Care No. • 91 22 2820 8707 Email: care@suru.com

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WARRANTY:

SURU warrants its products free from defects in workmanship and material under proper handling and use. This warranty voids if failure or defect is due to improper handling or storage or factors related to the patient, his / her Diagnosis, Treatment, Surgical Procedures or any other matter beyond SURU's control, directly affecting SURU's products and the results obtained from their use. Liability of SURU limited upto replacement of the Product only. SURU disclaims any responsibility for possible consequences resulting from improper use.

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