
Medtronic

PV-TRACKER™ 990045

Guide wire
Fil guide
Führungsdraht
Guia
Voerdraad
Filo guida
Guidewire
Ledare
Οδηγό σύρμα
Guidewire
Fio-guia
Vezetődrót
Prowadnik
Vodici drát
Проводник
Vodiaci drôt
Kilavuz tel
Fio-guia

As a registered Notary Public for the State of Minnesota,
I hereby certify that this is a true and accurate copy of the
original document. The original is filed at Medtronic, Inc.

Donna J. Bergin
May 08, 2015

Technical manual • Manuel technique • Technisches Handbuch • Manual técnico •
Technische handleiding • Manuale tecnico • Teknisk håndbok • Teknisk manual •
Τεχνικό εγχειρίδιο • Teknisk håndbog • Manual técnico • Műszaki leírás • Instrukcja techniczna •
Technická příručka • Техническое руководство • Technická příručka • Teknik el kitabı •
Manual técnico

CE
0050

The following list includes trademarks or registered trademarks of Medtronic in the United States and possibly in other countries. All other trademarks are the property of their respective owners. /

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1 General description

The PV-Tracker guide wire is intended to facilitate the placement of devices during diagnostic and interventional procedures.

The guide wire is manufactured from stainless steel with a polytetrafluoroethylene (PTFE) coating. The dispenser tube is fitted with a luer hub to facilitate flushing of the guide wire prior to use.

The "J" tip configuration on the distal end of the guide wire is packaged with a "J" straightener to aid in the insertion of the guide wire.

1.1 Contents of package

The package contains the following:

- 1 Guide wire
- 1 "J" straightener
- product documentation

2 Indications for use

The PV-Tracker guide wire is intended to facilitate the placement of devices during diagnostic and interventional procedures

3 Contraindications

Use of the PV-Tracker guide wire is contraindicated as follows:

- Atrial septal patch or atrial septal defect repair with percutaneous device
- Venous filtering device (i.e. Greenfield Filter)
- Active sepsis
- Left atrial thrombus or myxoma
- Blood clotting abnormalities
- Known sensitivity to heparin

The guide wire is not recommended for use in patients who cannot undergo standard anticoagulation protocol for a cardiac procedure, or who have had a recent coagulatory or embolic event.

4 Warnings and precautions

For single use only – This device is intended only to be used once for a single patient. Do not reuse, reprocess, or resterilize this device for purpose of reuse. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or create a risk of contamination of the device that could result in patient injury, illness, or death.

Do not resterilize – Do not resterilize this device for purpose of reuse. Resterilization may compromise the structural integrity of the device or create a risk of contamination from the device that could result in patient injury, illness, or death.

Expiration date – Check to verify the guide wire is within its expiration date. Do not use if the product date has expired.

Qualified users – This device should be used only by physicians thoroughly trained in percutaneous procedures.

Disposal of guide wire – Dispose of the guide wire according to hospital biohazard requirements. If returning the guide wire contact your local Medtronic representative.

Sterile package inspection – Inspect the sterile packaging and guide wire prior to use. If the sterile packaging or the guide wire exhibits damage, do not use the guide wire. Contact your local Medtronic representative.

Fluoroscopy required for guide wire placement – Use of fluoroscopy during guide wire manipulation and placement is strongly advised. Manipulating the guide wire without fluoroscopy may result in damage to cardiac and vascular structures.

X-ray and fluoroscopic exposure – Minimize x-ray and fluoroscopic exposure. Due to the intensity of the x-ray beam and the duration of fluoroscopic imaging, patients and laboratory staff may be subjected to acute radiation injury and increased risk for somatic and genetic effects. Take all appropriate measures to minimize x-ray exposure to both patients and clinical staff. The long-term effects of protracted fluoroscopy have not been established.

Storage conditions – store the guide wire in a cool, dry place.

Withdrawing a guide wire back through a metal needle – Avoid withdrawing a guide wire back through a metal needle. The sharp edge of the needle may scrape the PTFE coating or shear the guide wire. A catheter, introducer sheath or vessel dilator should replace the needle as soon as the guide wire has been inserted into the vasculature.

Tissue damage – Care should be taken when manipulating a catheter over the guide wire during the catheter's placement and withdrawal. If resistance is felt during catheter placement, discontinue the procedure and determine the cause of resistance before proceeding. If the cause of resistance cannot be determined, remove the guide wire and catheter as a unit to prevent possible tissue damage and complications.

Thrombus or emboli formation – When using a guide wire potential exists for thrombus formation, plaque dislodgement, and emboli. The physician should be familiar with the complications of interventional procedures.

5 Adverse events

Potential adverse events associated with cardiac guide wire procedures include, but are not limited to, the following conditions:

- Arrhythmias, proarrhythmias
- Atrial septal defect
- Allergic reaction to x-ray contrast media
- Catheter entrapment
- Cardiac tamponade
- Cerebrovascular accident (CVA) or transient ischemic attack (TIA)
- Death
- Heart Failure
- Hematoma
- Hemorrhage
- Hypo/hypertension
- Infection
- Myocardial infarction or ischemia
- Nerve damage
- Obstruction, perforation, damage or spasm of the vascular system including the coronary circulation system
- Pacemaker/defibrillator lead displacement
- Perforation or dissection
- Pericardial/pleural effusion
- Pericarditis or endocarditis
- Pseudoaneurysm in the groin
- Pneumothorax
- Pulmonary edema
- Pulmonary embolism
- Respiratory depression
- Retroperitoneal bleed
- Thromboembolic event
- Valvular insufficiency or damage
- Vasovagal reaction

6 Instructions for use

Caution: Carefully read all instructions prior to use. Observe all contraindications, warnings, and precautions noted in these instructions. Failure to do so may result in suboptimal outcomes.

Note: Administer anticoagulation therapy consistent with procedural indications and the institution's standard of care.

1. Prepare and utilize the guide wire under aseptic conditions.
2. Carefully remove the dispenser tube with loaded guide wire from the packaging.
3. Flush the guide wire by attaching a filled syringe to the luer hub on the guide wire dispenser tube and flush several times.
4. After flushing remove the guide wire from the tubing and inspect the guide wire prior to use to verify it has not been compromised during shipping or handling.
5. Insert the flexible "J" tip end of the guide wire into the percutaneous entry needle.

Note: For a "J" tip guide wire insertion, advance the "J" straightener forward to the distal tip of the wire so the "J" curve is straightened and easily inserted in to the hub of the introducer sheath or needle. After advancing the wire, slide the "J" straightener off the back end of the wire. The straightener may be discarded or flushed and saved for a future wire insertion during the procedure.

6. Advance the guide wire through the needle into the vasculature.
7. Hold the guide wire in place and withdraw the needle.
8. Position the guide wire as needed.
9. Advance and position the catheter over the guide wire.
10. Carefully withdraw the guide wire from the catheter.

11. Clean the guide wire and store in fluid filled basin until completion of the procedure.

7 Specifications

Length	200 cm (79")
Diameter	.8 mm (0.032")
Material	Stainless steel
Coatings	Polytetrafluoroethylene (PTFE)

8 Medtronic disclaimer of warranty

For complete warranty information, see the accompanying disclaimer of warranty document.

9 Service

Medtronic employs highly trained representatives and engineers located throughout the world to serve you and, upon request, to provide training to qualified hospital personnel in the use of Medtronic products. Medtronic also maintains a professional staff to provide technical consultation to product users. For more information, contact your local Medtronic representative, or call or write Medtronic at the appropriate telephone number or address listed on the back cover.

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© Medtronic, Inc. 2014