### Instructions for Use Statement

As a responsible official of Beckman Coulter, Inc., I hereby certify that the attached copies of Product Instruction for Use with this registration dossier for product registration in Russia are the true copies of the original documents on file at Beckman Coulter, Inc.

MODEL  
SX4700

NAME OF DEVICE  
Swinging-Bucket Rotor

Authorized Signatory:

NAME: Jun Maggie Zhang  
TITLE: Sr. Regulatory Affairs Analyst  
Beckman Coulter, Inc.  
State of Indiana

County of Marion

On 02 Feb 2016 before me, the undersigned, a Notary Public, in and for said County and State, Jun Maggie Zhang M personally appeared and proved to me on the basis of satisfactory evidence to be the person whose name is subscribed to the within instrument and acknowledged to me that she executed the same in her authorized capacity, and that by her signature on the instrument, the person executed the instrument.

I certify under PENALTY OF PERJURY under the laws of the State of Indiana that the foregoing paragraph is true and correct.

WITNESS my hand and official seal.

NAME: Janice C. Rice  
Notary Public  
Hamilton County  
State of Indiana

Beckman Coulter, Inc.  
250 S. Kraemer Boulevard  
Brea, CA 92821

Mailing Address:  
250 S. Kraemer Boulevard  
P.O. Box 8000  
Brea, CA 92822-8000

Telephone: (714) 993-5321  
Facsimile: (714) 961-4234  
Internet: [www.beckmancoulter.com](http://www.beckmancoulter.com)

**BECKMAN  
COULTER**

## **SX4700ROTOR IFU**

PROJECT NAME

**ALLEGRA X-5**

PROJECT NUMBER

**CI6515**

DOCUMENT NUMBER

**6515-033**

### **Totally Committed to Quality**

My electronic signature confirms my review and approval of this document for its intended use.

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#### **Approval Authorization:**

The following approvals are required for this documentation:

Bob Indig – Core Team Leader/ISDC

David Rolwing – Marketing & Sales

Kris Klink – Supply Chain

Maggie Zhang – Quality

Guenther Krischer – Service & Support

Lisa Brown – Regulatory

Larry Doriott – PCE

Vinay Goyal – EHS/PE

«version\_label

**See document attributes for additional data control elements.**

#### **CHANGE ABSTRACT:**


Removed incorrect statement about magnetic speed sensor system from the rotor description.

THE EFFECTIVITY OF THIS DOCUMENT IS GUARANTEED ONLY BY VERIFICATION AGAINST THE ON-LINE EDMS DATA. THE USER MUST VERIFY ITS CURRENCY ON THE DAY OF USE.

# Instructions For Use

## SX4700 Swinging-Bucket Rotor

For Use in the Beckman Coulter  
Allegra X-5 Centrifuge

PN B30343AA  
June 2013

Beckman Coulter, Inc.  
250 S. Kraemer Blvd.  
Brea, CA 92821 U.S.A.

## **JX4700 Swinging-Bucket Rotor**

PN B30343AA (June 2013)

Copyright 2013 Beckman Coulter, Inc.

All rights reserved. No part of this document may be reproduced or transmitted in any form or by any means, electronic, mechanical, photocopying, recording, or otherwise, without prior written permission from Beckman Coulter, Inc.

Beckman Coulter and the stylized logo are trademarks of Beckman Coulter, Inc. and are registered in the USPTO.

All other trademarks, service marks, products, or services are trademarks or registered trademarks of their respective holders.

Find us on the World Wide Web at:  
**[www.beckmancoulter.com](http://www.beckmancoulter.com)**

Made in Germany  
Printed in U.S.A.

# Safety Notice


Read all product manuals and consult with Beckman Coulter-trained personnel before attempting to use this equipment. Do not attempt to perform any procedure before carefully reading all instructions. Always follow product labeling and manufacturer's recommendations. If in doubt as to how to proceed in any situation, contact your Beckman Coulter Representative.

## Alerts for Warning, Caution, and Note

---

**WARNING** indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

**CAUTION** indicates a potentially hazardous situation, which, if not avoided, may result in minor or moderate injury. It may also be used to alert against unsafe practices.

**NOTE** NOTE is used to call attention to notable information that should be followed during installation, use, or servicing of this equipment.

## Safety Information for the SX4700 Rotor

---

Handle body fluids with care because they can transmit disease. No known test offers complete assurance that such fluids are free of micro-organisms. Some of the most virulent — Hepatitis (B and C) viruses, HIV (I–V), atypical mycobacteria, and certain systemic fungi — further emphasize the need for aerosol protection. Handle other infectious samples according to good laboratory procedures and methods to prevent spread of disease. Because spills may generate aerosols, observe proper safety precautions for aerosol containment. Do not run toxic, pathogenic, or radioactive materials in this rotor without taking appropriate safety precautions. Biosafe containment should be used when Risk Group II materials (as identified in the World Health Organization *Laboratory Biosafety Manual*) are handled; materials of a higher group require more than one level of protection.

The rotor and accessories are not designed for use with materials capable of developing flammable or explosive vapors, or a hazardous chemical reaction. Do not centrifuge such materials in nor handle or store them near the centrifuge.

This rotor was developed, manufactured, and tested for safety and reliability as part of a Beckman Coulter centrifuge/rotor system. Its safety or reliability cannot be assured if used in a non-Beckman Coulter centrifuge, in a Beckman Coulter centrifuge not specified for use with the rotor, or in a Beckman Coulter centrifuge that has been modified without Beckman Coulter's approval.

## Safety Notice

Safety Information for the SX4700 Rotor

Although rotor components and accessories made by other manufacturers may fit in the SX4700 rotor, their safety in this rotor cannot be ascertained by Beckman Coulter. Use of other manufacturers' components or accessories in the SX4700 rotor may void the rotor warranty and should be prohibited by your laboratory safety officer. Only the components and accessories listed in this publication should be used in this rotor.

Hook all four buckets, loaded or empty, to the rotor for every run. Make sure that filled containers are loaded symmetrically into the rotor and that opposing tubes are filled to the same level with liquid of the same density.

If disassembly reveals evidence of leakage, you should assume that some fluid escaped the rotor. Apply appropriate decontamination procedures to the centrifuge and accessories.

Never exceed the maximum rated speed of the rotor and labware in use. Refer to the section on *Run Speeds*, and derate the run speed as appropriate.

Do not use sharp tools on the rotor that could cause scratches in the rotor surface. Corrosion begins in scratches and may open fissures in the rotor with continued use.

# Contents

## Safety Notice, iii

Alerts for Warning, Caution, and Note, iii

Safety Information for the SX4700 Rotor, iii

## SX4700 Swinging-Bucket Rotor, 1

Specifications, 1

Description, 2

The Rotor, 2

Buckets and Accessories, 2

Tube Breakage, 3

Preparation and Use, 4

Tubes, 4

Temperature Limits, 4

Prerun Safety Checks, 4

Installing the Rotor Yoke, 5

Loading the Buckets, 6

Symmetrical and Balanced Loading, 6

Using Buckets, 6

Operation, 8

Removal and Sample Recovery, 8

Run Speeds, 9

Speed Derating, 9

Care and Maintenance, 10

Maintenance, 10

Cleaning, 10

Rotor Components, 11

Buckets, 11

Adapters, 12

Decontamination, 12

Sterilization and Disinfection, 12

Tube Breakage, 13

Storage, 13

Returning a Rotor, 13

Supply List, 14

Replacement Rotor Parts, 14  
Other, 14

Beckman Coulter, Inc.  
Benchtop Rotor Warranty

# Illustrations


- 1 The SX4700 Rotor and Accessories, 3
- 2 Installing the Rotor Yoke on the Centrifuge Drive Shaft, 5
- 3 Arranging tubes in the Rotor, 7
- 4 Seating Buckets on Yoke Pivot Pins, 8

# Tables

|   |                                                 |
|---|-------------------------------------------------|
| 1 | Available Tube Adapters for the SX4700 Rotor, 3 |
|---|-------------------------------------------------|

# SX4700 Swinging-Bucket Rotor

## Specifications

1. Axis of Rotation

Maximum speed 4700 RPM

Density rating at maximum speed 1.2 g/mL

Critical speed range<sup>a</sup> 500 to 800 RPM

Relative Centrifugal Field<sup>b</sup> at maximum speed at  $r_{\text{max}}$  (181 mm)  $4470 \times g$

Conditions requiring speed reductions see *Run Speeds*

Maximum allowable imbalance of opposing loads 20 grams

Number of buckets 4

Available labware Table 1

Nominal tube capacity (largest tube) 10 mL

Nominal rotor capacity 1000 mL

Approximate acceleration time to maximum speed  
(fully loaded) 33 sec

Approximate deceleration time from maximum speed  
(fully loaded) 27 sec


Weight of fully loaded rotor 10.4 kg (22.8 lbs)

Rotor and bucket material anodized aluminum

- a. The critical speed range is the range of speeds over which the rotor shifts so as to rotate about its center of mass. Passing through the critical speed range is characterized by some vibration.
- b. Relative Centrifugal Field (RCF) is the ratio of the centrifugal acceleration at a specified radius and speed ( $r\omega^2$ ) to the standard acceleration of gravity ( $g$ ) according to the following formula:  $RCF = r\omega^2/g$  — where  $r$  is the radius in millimeters,  $\omega$  is the angular velocity in radians per second ( $2\pi \text{ RPM}/60$ ), and  $g$  is the standard acceleration of gravity ( $9807 \text{ mm/s}^2$ ). After substitution:  $RCF = 1.12r (\text{RPM}/1000)^2$

## Description

---

Beckman Coulter SX4700 rotors are manufactured in a facility that maintains certifications to both ISO 9001:2008 or ISO 13485:2003. They are for use with the Beckman Coulter Allegra X-5 Centrifuge.

## The Rotor

The SX4700 is a four-place swinging bucket rotor. This rotor carries adapters that allow centrifugation of standard 5/7 or 10 mL blood collection tubes.


The rotor yoke and buckets are made of aluminum and are black-anodized for corrosion protection. The buckets are run by placing them over pivot pins on the arms of the yoke. Buckets swing out to a horizontal position during centrifugation. A tie-down screw is used to secure the rotor to the drive shaft during centrifugation. Figure 1 shows the rotor and accessories.

## Buckets and Accessories

Buckets hold several sizes of labware in one-piece adapters.

A variety of tube sizes can be centrifuged in the rotor, supported in tube adapters listed in Table 1. The adapters can also serve as tube racks in the laboratory.

**Figure 1** The SX4700 Rotor and Accessories

1. Bucket                      2. T-handle Rotor Wrench                      3. Rotor Yoke

**NOTE** The latest information regarding labware and accessories used in Beckman Coulter rotors can be obtained at [www.beckmancoulter.com](http://www.beckmancoulter.com).

**Table 1** Available Tube Adapters for the SX4700 Rotor

| Adapters            |                     | Tubes                                                             |                | Maximum Number Tubes per Rotor |
|---------------------|---------------------|-------------------------------------------------------------------|----------------|--------------------------------|
| Part Number (pkg 2) | Tubes per Adapter   | Description                                                       | Nominal Volume |                                |
| B30602              | 35                  | 5/7 mL blood collection tubes                                     | 5/7 mL         | 140                            |
| B30603              | 25                  | 10 mL blood collection tubes                                      | 10 mL          | 100                            |
| B30604              | 4 racks<br>16 tubes | Beckman DxC/DxI 4-place 13/16 mm tube sample racks                | varies         | 16 racks 64 tubes              |
| B32766              | 25                  | Adapter for Beckman Coulter AutoMate 2500 Family, 25 x 13mm tubes | 5/7 mL         | 100                            |
| B32767              | 16                  | Adapter for Beckman Coulter AutoMate 2500 Family, 16 x 17mm tubes | 10 mL          | 64                             |

See the Warranty at the back of this manual for warranty information.

## Tube Breakage

**To reduce the potential for corrosion, clean the buckets thoroughly immediately following a tube breakage.**

If a glass tube breaks, remove the glass very carefully from the adapter and bucket. Imbedded glass particles that remain in the bucket or adapters can cause tube failure during subsequent runs.

## Preparation and Use

*Specific information about the SX4700 rotor is given here. Information about the centrifuge is contained in the centrifuge manual, which should be used together with this manual for complete centrifuge, rotor, and accessory operation.*

**NOTE** Although rotor components and accessories made by other manufacturers may fit in the SX4700 rotor, their safety in this rotor cannot be ascertained by Beckman Coulter. Use of other manufacturers' components or accessories in the SX4700 rotor may void the rotor warranty and should be prohibited by your laboratory safety officer. Only the components and accessories listed in this publication should be used in this rotor.

### Tubes

Tubes should be pretested under anticipated run conditions (using water instead of valuable samples) if operating below 4°C or above room temperature. Refer to *Chemical Resistances* (publication IN-175) for information on the chemical resistances of tube and adapter materials.

### Temperature Limits

- Plastic containers have been centrifuge tested for use at temperatures between 2 and 25°C. For centrifugation at other temperatures, pretest tubes under anticipated run conditions.
- If plastic containers are frozen before use, make sure that they are thawed to at least 2°C prior to centrifugation.

### Prerun Safety Checks

*Read the Safety Notice section at the front of this manual before using the rotor.*

- 1 Make sure that the rotor and buckets are clean and show no signs of corrosion or cracking. If any evidence of damage is present, do not centrifuge the rotor.
- 2 Check the chemical compatibilities of all materials used (refer to *Chemical Resistances*, publication IN-175).
- 3 Verify that the tubes and accessories being used are listed in Table 1.

## Installing the Rotor Yoke

**NOTE** Before the first use of a new rotor, lubricate the pivot pins with Spinkote (306812). After the first use, apply Spinkote to the pivot pins approximately once a week and after cleaning or autoclaving.

**NOTE** Refer to the centrifuge manual for proper installation and removal instructions. Failure to follow printed instructions may result in damage to the rotor and the centrifuge.


- 1 Before installing the yoke in the centrifuge, lightly coat the drive hole with lubricant such as Spinkote (306812) (see instructions under *Maintenance*).
- 2 Carefully lower the yoke straight down onto the centrifuge drive shaft (see Figure 2). Be sure the yoke is properly seated on the shaft.
- 3 Fasten the tie-down screw (368245) onto the drive shaft. Use the T-handle rotor wrench (368246) to tighten the tie-down screw firmly on the shaft.

### CAUTION

**Never drop the rotor yoke onto the centrifuge drive shaft. The drive shaft can be damaged if the rotor is forced sideways or dropped onto it.**

**NOTE** If the yoke is left in the centrifuge between runs, before each run make sure it is properly seated on the drive shaft and that the tie-down screw is tight.

Figure 2 Installing the Rotor Yoke on the Centrifuge Drive Shaft

1. Lower the yoke onto the drive shaft
2. Tighten the tie-down screw with the T-handle wrench

## Loading the Buckets

*For runs at other than room temperature, refrigerate or warm the rotor and precool the centrifuge beforehand for faster equilibration.*

### **WARNING**

Handle body fluids with care because they can transmit disease. No known test offers complete assurance that such fluids are free of micro-organisms. Some of the most virulent—Hepatitis (B and C) viruses, HIV (I–V), atypical mycobacteria, and certain systemic fungi—further emphasize the need for aerosol protection. Handle other infectious samples according to good laboratory procedures and methods to prevent spread of disease. Because spills may generate aerosols, observe proper safety precautions for aerosol containment.

Do not run toxic, pathogenic, or other hazardous materials in this rotor without taking all appropriate safety precautions. Biosafe containment should be used when Risk Group II materials (as identified in the World Health Organization *Laboratory Biosafety Manual*) are handled; materials of a higher group require more than one level of protection.

## Symmetrical and Balanced Loading

To ensure optimal performance and stability, the rotor must be loaded symmetrically. Two factors affect symmetric loading:

1. The buckets or carriers must be loaded symmetrically with respect to their pivotal axes.
2. The rotor should be loaded symmetrically with respect to its center of rotation.

For best results, load opposing buckets or carriers with the same type of labware containing the same amounts of fluid of equal density. Additionally, opposing buckets and their contents should weigh approximately the same (within 10 grams).

It is not necessary to completely fill all tubes or positions in buckets; however, partially filled buckets must be balanced with respect to the bucket pivotal axis. Each tube should be placed so that its weight is balanced by a tube in a diametrically opposite position across the pivotal axis in the same bucket. Opposing buckets should also be filled the same way.


## Using Buckets

You can load buckets before or after they are installed on the rotor yoke. In either case, we recommend filling the appropriate labware first and then loading the labware into the buckets to avoid tripping the imbalance detector.

- 1 Load the filled containers into the buckets (see *Tubes*, page 4 for tube information).

**NOTE** Tubes placed in opposing buckets must be filled to the same level with liquid of the same density. Fewer than maximum number of tubes can be centrifuged in the multi-tube adapters if they are arranged symmetrically in the buckets (see Figure 3).

**Figure 3** Arranging tubes in the Rotor

**NOTE** All buckets must be attached to the rotor for every run, whether loaded or empty.  
Never run the rotor with fewer than four positions filled


- 2 Attach each bucket to the yoke by aligning the grooves in the bucket sides with the pivot pins, then sliding the buckets down until the pivot pins are seated in the bucket pockets.

**CAUTION**

Attach all four buckets, loaded or empty, to the rotor yoke. If only two buckets are filled, place them in opposite positions on the yoke. All four positions must contain buckets during a run. Always have adapters installed in all four buckets, even if there are no tubes in the adapters.

- 3 Gently swing the buckets to ensure that they are properly seated on the pivot pins.
- 4 Make sure buckets are properly seated over the pivot pins on the rotor yoke as shown in Figure 4.

Figure 4 Seating Buckets on Yoke Pivot Pins

- |              |                                                |
|--------------|------------------------------------------------|
| 1. CORRECT   | 3. Pivot pins are seated in bucket pockets     |
| 2. INCORRECT | 4. Pivot pins are not seated in bucket pockets |

**NOTE** Buckets must be inserted at an angle, or the bucket bottom will rest on the rotor base and the pivot pins will not seat in the pin pockets.

## Operation

- 1 Refer to the centrifuge instruction manual for centrifuge operation.
- 2 See *Run Speeds*, page 9, for information about speed limitations.

## Removal and Sample Recovery

### CAUTION

If disassembly reveals evidence of leakage, you should assume that some fluid escaped the rotor. Apply appropriate decontamination procedures to the centrifuge and accessories.

- 1 Remove the buckets from the centrifuge.
- 2 Remove the tubes from the buckets.

**NOTE** If glass tubes break, remove the glass very carefully from the adapter and bucket. If all the glass particles are not contained in the bucket or adapter, be careful when examining or cleaning the centrifuge gasket and chamber as glass particles may be embedded in their surfaces.

- 3 If removing the rotor yoke, turn the T-handle wrench to the left (counterclockwise) to loosen the tie-down screw.

Lift the yoke straight up and off the drive shaft.

**NOTE** Refer to the centrifuge manual for proper installation and removal instructions. Failure to follow printed instructions may result in damage to the rotor and the centrifuge.

## Run Speeds

The centrifugal force at a given radius in a rotor is a function of speed. Comparisons of forces between different rotors are made by comparing the rotors' relative centrifugal fields (RCF). When rotational speed is adjusted so that identical samples are subjected to the same RCF in two different rotors, the samples are subjected to the same force. If RPM is entered, the RCF at each speed is automatically calculated by the centrifuge software; if the RCF is entered, the centrifuge calculates the equivalent revolutions per minute (RPM). *Do not select rotational speeds that exceed 4700 RPM.*

### Speed Derating

The maximum run speed listed in the rotor specifications is for operation when all conditions are within the standard specifications. Maximum allowable run speed for the rotor must be reduced in some circumstances.

1. The rotor can process solutions at the maximum rated speed if the solution density is 1.2 g/mL or less. When centrifuging solutions of densities greater than 1.2 g/mL, the maximum run speeds must be reduced according to the following equation to protect the rotor from excessive stresses due to the added tube load:

$$\text{reduced maximum speed} = (4700 \text{ rpm}) \sqrt{\frac{1.2 \text{ g/mL}}{\text{density of tube contents}}}$$

2. When some glass tubes are used, reduce rotor speed to prevent tube breakage. Because the strength of tubes can vary from lot to lot, and will depend on handling and usage, Beckman Coulter highly recommends that you pretest these tubes in the appropriate labware in the SX4700 rotor, using water samples, to determine optimal operating conditions.


**NOTE** Scratches (even microscopic ones) significantly weaken glass tubes.

## Care and Maintenance

---

### Maintenance

*Do not use sharp tools on the rotor, as they can scratch the anodized surface. Corrosion begins in scratches and may open fissures in the rotor with continued use.*

- 1 Periodically (at least monthly) inspect the rotor yoke and buckets, especially inside cavities, for rough spots or pitting, white powder deposits—frequently aluminum oxide—or heavy discoloration.

If any of these signs are evident, do not run the rotor.

Contact your Beckman Coulter representative for information about the Field Rotor Inspection Program and the rotor repair center.

- 2 Regularly check the condition of the tie-down screw (B31136); if it is worn or damaged, replace it.
- 3 Approximately once a week, and after cleaning and/or autoclaving, lubricate the rotor pins and pin sockets with a lubricant such as Spinkote (306812).

Store the rotor in a dry environment (not in the centrifuge). Refer to *Chemical Resistances* (IN-175) for the chemical compatibilities of rotor and accessory materials. Your Beckman Coulter representative provides contact with the Field Rotor Inspection Program and the rotor repair center.

### Cleaning

*Wash rotor components immediately if salts or other corrosive materials are used or if spillage has occurred. Do not allow corrosive materials to dry on the rotor.*

Under normal use, wash the rotor frequently (at least weekly) to prevent buildup of residues.

**NOTE** Do not wash the rotor components in a dishwasher. Do not soak the rotor in detergent solution for long periods, such as overnight.

## Rotor Components

- 1 Wash the rotor yoke and tie-down screw in a mild detergent, such as Solution 555 (339555), that won't damage the rotor.  
Dilute the detergent with water (10 parts water to 1 part detergent).  
The Rotor Cleaning Kit (339558) contains two plastic-coated brushes and two quarts of Solution 555 for use with rotors and accessories.
- 2 Thoroughly rinse the cleaned rotor components with distilled water.
- 3 Air-dry the rotor components upside down.  
*Do not use acetone to dry the rotor.*
- 4 Lubricate the rotor pins with Spinkote (306812).

Before reinstalling the rotor yoke, lightly lubricate the drive hole with Spinkote (306812) to prevent the rotor from sticking, as follows:

- 1 Apply the lubricant onto a swab.
- 2 Draw the coated swab through a paper towel to remove excess lubricant.
- 3 Lightly coat the inside of the drive hole with the lubricant remaining on the swab.

## Buckets

- 1 Wash the buckets with a mild detergent such as Solution 555 (diluted with 10 parts water to 1 part detergent) and a soft brush.
- 2 Thoroughly rinse with water and air-dry.

## Adapters

- 1 Use a mild detergent such as Solution 555, diluted 10 to 1 with water, and a soft brush to scrub the adapters.
- 2 Rinse thoroughly with water.
- 3 Air-dry upside down.

## Decontamination

If aluminum rotor components become contaminated with radioactive material, decontaminate them using a solution that will not damage the anodized surfaces. Beckman Coulter has tested a number of solutions and found two that do not harm anodized aluminum: RadCon Surface Spray or IsoClean Solution (for soaking),\* and Radiacwash.†

While Beckman Coulter has tested these materials and found that they do not damage components, no guarantee of decontamination is expressed or implied. Follow appropriate decontamination procedures as directed by your laboratory safety officer.

If the rotor or other components are contaminated with toxic or pathogenic materials, follow appropriate decontamination procedures as directed by your laboratory safety officer.

## Sterilization and Disinfection

- The rotor and all rotor components can be autoclaved at 121°C for up to an hour. Place the rotor yoke and buckets in the autoclave upside down.
- Ethanol (70%) or hydrogen peroxide (6%) may be used on all rotor components, including those made of plastic. Bleach (sodium hypochlorite) may be used, but may cause discoloration of anodized surfaces. Use the minimum immersion time for each solution, per laboratory standards.

\* In U.S., contact Nuclear Associates (New York); in Eastern Europe and Commonwealth States, contact Victoreen GmbH (Munich); in South Pacific, contact Gammasonics Pty. Ltd. (Australia); in Japan, contact Toyo Medic Co. Ltd. (Tokyo).

† In U.S., contact Biodex Medical Systems (Shirley, New York); internationally, contact the U.S. office to find the dealer closest to you.

 **CAUTION**

**Ethanol is a flammability hazard. Do not use it in or near operating centrifuges.**

While Beckman Coulter has tested these methods and found that they do not damage the rotor or components, no guarantee of sterility or disinfection is expressed or implied. When sterilization or disinfection is a concern, consult your laboratory safety officer regarding proper methods to use.

Refer to publication IN-192, included in each box of tubes or bottles, for tube sterilization and disinfection procedures.

## Tube Breakage

 **CAUTION**

**To reduce the potential for corrosion, clean buckets thoroughly immediately following a tube breakage.**

If a glass tube breaks, remove the glass very carefully from the adapter and bucket. Imbedded glass particles that remain in the bucket or adapters can cause tube failure during subsequent runs.

## Storage

When it is not in use, store the rotor in a dry environment (not in the centrifuge).

## Returning a Rotor

Before returning a rotor or accessory for any reason, prior permission must be obtained from Beckman Coulter, Inc. The authorization form may be obtained from your local Beckman Coulter sales office. The form, entitled *Returned Material Authorization (RMA)* for United States returns or *Returned Goods Authorization (RGA)* for international returns, should contain the following information:

- rotor type and serial number,
- history of use (approximate frequency of use),
- reason for the return,
- original purchase order number, billing number, and shipping number, if possible,
- name and email address of the person to be notified upon receipt of the rotor or accessory at the factory,
- name and email address of the person to be notified about repair costs, etc.

To protect our personnel, it is the customer's responsibility to ensure that all parts are free from pathogens and/or radioactivity. Sterilization and decontamination must be done before returning the parts. Smaller items (such as tubes, bottles, etc.) should be enclosed in a sealed plastic bag.

All parts must be accompanied by a note, plainly visible on the outside of the box or bag, stating that they are safe to handle and that they are not contaminated with pathogens or radioactivity. Failure to attach this notification will result in return or disposal of the items without review of the reported problem.

Use the address label printed on the authorization form when mailing the rotor and/or accessories.

Customers located outside the United States should contact their local Beckman Coulter office.

Supply List

Call Beckman Coulter Customer Service at 1-800-742-2345 (U.S.A. or Canada) or contact your local Beckman Coulter Representative for detailed information on ordering parts and supplies. For your convenience, a partial list is given below.

Replacement Rotor Parts

| Description                                | Part Number |
|--------------------------------------------|-------------|
| SX4700 rotor assembly (includes 4 buckets) | B30593      |
| Buckets (set of 2)                         | B30718      |
| Tie-down screw                             | B31136      |
| T-handle rotor wrench                      | 368246      |

Other

**NOTE** For MSDS information, go to the Beckman Coulter website at [www.beckmancoulter.com](http://www.beckmancoulter.com).

| Description                   | Part Number |
|-------------------------------|-------------|
| Adapters                      | see Table 1 |
| Spinkote lubricant (2 oz)     | 306812      |
| Silicone vacuum grease (1 oz) | 335148      |
| Rotor Cleaning Kit            | 339558      |
| Solution 555 (1 qt)           | 339555      |
| Rotor cleaning brush          | 339379      |

# Beckman Coulter, Inc.

## Benchtop Rotor Warranty

Subject to the conditions specified below and the warranty clause of the Beckman Coulter, Inc., terms and conditions in effect at the time of sale, Beckman Coulter agrees to correct either by repair or, at its election, by replacement, any defects of material or workmanship which develop within seven (7) years after delivery of a benchtop centrifuge rotor to the original buyer by Beckman Coulter or by an authorized representative, provided that investigation and factory inspection by Beckman Coulter discloses that such defect developed under normal and proper use. Should a Beckman Coulter centrifuge be damaged due to a failure of a rotor covered by this warranty, Beckman Coulter will supply free of charge all centrifuge parts required for repair.

### Conditions

Except as otherwise specifically provided herein, this warranty covers the rotor only and Beckman Coulter shall not be liable for damage to accessories or ancillary supplies including but not limited to (i) tubes, (ii) tube caps, (iii) tube adapters, or (iv) tube contents.

This warranty is void if the rotor has been subjected to customer misuse such as operation or maintenance contrary to the instructions in the Beckman Coulter rotor or centrifuge manual.

This warranty is void if the rotor is operated with a rotor drive unit or in a centrifuge unmatched to the rotor characteristics or operated in a Beckman Coulter centrifuge that has been improperly disassembled, repaired, or modified.


Thermoplastic rotors or components used in some benchtop centrifuges are warranted for one (1) year from date of purchase.

### Disclaimer

IT IS EXPRESSLY AGREED THAT THE ABOVE WARRANTY SHALL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND OF THE WARRANTY OF MERCHANTABILITY AND THAT NEITHER BECKMAN COULTER, INC. NOR ITS SUPPLIERS SHALL HAVE ANY LIABILITY FOR SPECIAL OR CONSEQUENTIAL DAMAGES OF ANY KIND WHATSOEVER ARISING OUT OF THE MANUFACTURE, USE, SALE, HANDLING, REPAIR, MAINTENANCE, OR REPLACEMENT OF THE PRODUCT.

Warranty-2

PN B30343AA

## Related Documents

### **Allegra X-5 Centrifuge (B29071)**

- Safety
- Introduction
- Description
- Operation
- Troubleshooting
- Care and Maintenance
- Installation

Available in hard copy or electronic pdf by request.

### **Additional References**

- Chemical Resistances for Beckman Coulter Centrifugation Products (IN-175)
- Use and Care of Tubes and Bottles (IN-192)

Available in hard copy or electronic pdf by request.

Available at [www.beckmancoulter.com](http://www.beckmancoulter.com)

[www.beckmancoulter.com](http://www.beckmancoulter.com)

© 2013 Beckman Coulter, Inc.  
All Rights Reserved

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