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1. Hawks RL, Chiang CN, editors. *Urine Testing for Drugs of Abuse*. NIDA Research Monograph 73. Rockville, MD: National Institute on Drug Abuse (NIDA), Department of Health and Human Services; 1986.
2. *Mandatory Guidelines for Federal Workplace Drug Testing Programs*. National Institute on Drug Abuse. Federal Register Vol. 69, No. 71, 2004:19650–61.
3. Baselt RC. *Disposition of Toxic Drugs and Chemicals in Man*, 5th ed. Foster City, CA: Chemical Toxicology Institute; 2000:211–15, 407–12, 589–93.
4. Rubenstein KE, Schneider RS, Ullman EF. Homogeneous enzyme immunoassay: a new immunochemical technique. *Biochem Biophys Res Commun* 1972;47:846–51.
5. US Department of Labor, Occupational Safety and Health Administration. 29 CFR Part 1910.1030. *Bloodborne Pathogens*.
6. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*, 5th ed. Washington, DC: US Government Printing Office, January 2007.
7. World Health Organization. *Laboratory Biosafety Manual*, 3rd ed. Geneva: World Health Organization, 2004.
8. Sewell DL, Bove KE, Callihan DR, et al. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Third Edition (M29-A3)*. Wayne, PA: Clinical and Laboratory Standards Institute, 2005.
9. Armbruster DA, Dasgupta A, et al. *Toxicology and Drug Testing in the Clinical Laboratory; Approved Guideline—Vol. 27, No. 15*. Wayne, PA: Clinical and Laboratory Standards Institute, 2007.
10. Hayes LW, Krasselt WG, Mueggler PA. Concentrations of morphine and codeine in serum and urine after ingestion of poppy seeds. *Clin Chem* 1987;33(6):806–8.

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| Configure assay parameters — General                 |                                        |                                                  |                               |
|------------------------------------------------------|----------------------------------------|--------------------------------------------------|-------------------------------|
| <input checked="" type="radio"/> General             | <input type="radio"/> Calibration      | <input type="radio"/> SmartWash                  | <input type="radio"/> Results |
| Assay: <b>Opi3Q</b>                                  | Type: <b>Photometric</b>               | Version: <b>†</b>                                |                               |
| Number: <b>2843</b>                                  | Assay availability: <b>Enabled</b>     | Run controls for onboard reagents by: <b>Lot</b> |                               |
| <input checked="" type="radio"/> Reaction definition | <input type="radio"/> Reagent / Sample | <input type="radio"/> Validity checks            |                               |
| Reaction mode: <b>Rate up</b>                        |                                        |                                                  |                               |
| Primary                                              |                                        | Secondary                                        | Read times                    |
| Wavelength: <b>340</b> / <b>††</b>                   |                                        | Main: <b>20 – 24</b>                             |                               |
| Last required read: <b>24</b>                        |                                        | Flex: <b>— —</b>                                 |                               |
| Absorbance range: <b>-0.1000 – 3.2000</b>            |                                        | Color correction: <b>— —</b>                     |                               |
| Sample blank type: <b>None</b>                       |                                        |                                                  |                               |

| O Reaction definition                |                              |          |          | ● Reagent / Sample |                  | O Validity checks |  |
|--------------------------------------|------------------------------|----------|----------|--------------------|------------------|-------------------|--|
| Reagent: <b>OPIM9</b>                | Reagent volume: <b>100</b>   |          | R1       | R2                 |                  |                   |  |
| Diluent: <b>None</b>                 | Water volume: <b>—</b>       |          |          |                    |                  |                   |  |
| Diluent dispense mode: <b>Type 0</b> | Dispense mode: <b>Type 0</b> |          | Type 0   | Type 0             |                  |                   |  |
| Dilution name: <b>Sample</b>         | Diluted sample               | Diluent  | Water    | Dilution factor    | Default dilution |                   |  |
| <b>STANDARD : 12.0</b>               | <b>—</b>                     | <b>—</b> | <b>—</b> | <b>1:1.00</b>      | <b>●</b>         |                   |  |

| O Reaction definition       |  | O Reagent / Sample |  | ● Validity checks |  |
|-----------------------------|--|--------------------|--|-------------------|--|
| Reaction check: <b>None</b> |  |                    |  |                   |  |
| Rate linearity %: <b>—</b>  |  |                    |  |                   |  |

| Configure assay parameters — Calibration     |                                              |                                 |                                       |
|----------------------------------------------|----------------------------------------------|---------------------------------|---------------------------------------|
| <input type="radio"/> General                | <input checked="" type="radio"/> Calibration | <input type="radio"/> SmartWash | <input type="radio"/> Results         |
| Assay: <b>Opi3Q</b>                          | Calibration method: <b>Linear</b>            |                                 |                                       |
| <input checked="" type="radio"/> Calibrators | <input type="radio"/> Volumes                | <input type="radio"/> Intervals | <input type="radio"/> Validity checks |
| Calibrator set: <b>OPQ</b>                   | Calibrator level: <b>(†††)</b>               | Concentration: <b>†</b>         |                                       |
|                                              | Blank: <b>OPQ1 (DOA MC Neg Cal)</b>          | <b>0</b>                        |                                       |
|                                              | Cal 1: <b>OPQ2 (Opi 300 Cal)</b>             |                                 |                                       |
| Replicates: <b>2</b> [Range 1 – 3]           |                                              |                                 |                                       |

| O Calibrators          |                    | ● Volumes   |                | O Intervals |          | O Validity checks |  |
|------------------------|--------------------|-------------|----------------|-------------|----------|-------------------|--|
| Calibrator: <b>OPQ</b> | Calibrator level   | Sample      | Diluted sample | Diluent     | Water    |                   |  |
|                        | Blank: <b>OPQ1</b> | <b>12.0</b> | <b>—</b>       | <b>—</b>    | <b>—</b> |                   |  |
|                        | Cal 1: <b>OPQ2</b> | <b>12.0</b> | <b>—</b>       | <b>—</b>    | <b>—</b> |                   |  |

| O Calibrators                     |  | O Volumes |  | ● Intervals |  | O Validity checks |  |
|-----------------------------------|--|-----------|--|-------------|--|-------------------|--|
| Calibration intervals:            |  |           |  |             |  |                   |  |
| Full interval: <b>312</b> (hours) |  |           |  |             |  |                   |  |
| Calibration type:                 |  |           |  |             |  |                   |  |
| Adjust type: <b>None</b>          |  |           |  |             |  |                   |  |

| O Calibrators                    |  | O Volumes    |  | O Intervals  |  | ● Validity checks |  |
|----------------------------------|--|--------------|--|--------------|--|-------------------|--|
| Blank absorbance range:          |  |              |  |              |  |                   |  |
| Span:                            |  | <b>Blank</b> |  | <b>Blank</b> |  |                   |  |
| Span absorbance range:           |  |              |  |              |  |                   |  |
| Expected cal factor:             |  | <b>0.00</b>  |  |              |  |                   |  |
| Expected cal factor tolerance %: |  | <b>0</b>     |  |              |  |                   |  |

| Configure assay parameters — SmartWash                        |                                   |                                            |                               |                                      |
|---------------------------------------------------------------|-----------------------------------|--------------------------------------------|-------------------------------|--------------------------------------|
| <input type="radio"/> General                                 | <input type="radio"/> Calibration | <input checked="" type="radio"/> SmartWash | <input type="radio"/> Results | <input type="radio"/> Interpretation |
| Assay: <b>Opi3Q</b>                                           |                                   |                                            |                               |                                      |
| COMPONENT                                                     | REAGENT / ASSAY                   | WASH                                       | Volume                        | Replicates                           |
| Cuvette                                                       | <b>Trig**</b>                     | <b>10% Detergent B</b>                     | <b>345</b>                    |                                      |
| **Not required for ARCHITECT Software version 7.00 and above. |                                   |                                            |                               |                                      |

| Configure assay parameters — Results |                                   |                                 |                                          |                                      |
|--------------------------------------|-----------------------------------|---------------------------------|------------------------------------------|--------------------------------------|
| <input type="radio"/> General        | <input type="radio"/> Calibration | <input type="radio"/> SmartWash | <input checked="" type="radio"/> Results | <input type="radio"/> Interpretation |
| Assay: <b>Opi3Q</b>                  |                                   | Assay number: <b>2843</b>       |                                          |                                      |
| Dilution default range:              |                                   | Result units: <b>—</b>          |                                          |                                      |
| Low-Linearity:                       |                                   |                                 |                                          |                                      |
| High-Linearity:                      |                                   |                                 |                                          |                                      |

| Configure assay parameters — Interpretation |                                   |                                 |                               |                                                 |
|---------------------------------------------|-----------------------------------|---------------------------------|-------------------------------|-------------------------------------------------|
| <input type="radio"/> General               | <input type="radio"/> Calibration | <input type="radio"/> SmartWash | <input type="radio"/> Results | <input checked="" type="radio"/> Interpretation |
| Assay: <b>Opi3Q</b>                         |                                   | Assay number: <b>2843</b>       |                               |                                                 |
| Name:                                       | Range:                            | Results review required:        |                               |                                                 |
| <b>Negative</b>                             | <b>&lt; 300</b>                   | <input type="checkbox"/>        |                               |                                                 |
| <b>Positive</b>                             | <b>&gt;= 300</b>                  | <input type="checkbox"/>        |                               |                                                 |

| Configure result units            |               |
|-----------------------------------|---------------|
| Assay: <b>Opi3Q</b>               |               |
| Version: <b>†</b>                 |               |
| Result units: <b>—</b>            |               |
| Decimal places: <b>0</b>          | [Range 0 – 4] |
| Correlation factor: <b>1.0000</b> |               |
| Intercept: <b>0.0000</b>          |               |

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| Configure assay parameters — General                 |                                    |                                                  |                               |
|------------------------------------------------------|------------------------------------|--------------------------------------------------|-------------------------------|
| <input checked="" type="radio"/> General             | <input type="radio"/> Calibration  | <input type="radio"/> SmartWash                  | <input type="radio"/> Results |
| Assay: <b>Opi1SQ</b>                                 | Type: <b>Photometric</b>           | Version: †                                       |                               |
| Number: <b>2841</b>                                  | Assay availability: <b>Enabled</b> | Run controls for onboard reagents by: <b>Lot</b> |                               |
| <input checked="" type="radio"/> Reaction definition |                                    |                                                  |                               |
| Reaction mode: <b>Rate up</b>                        |                                    | Read times                                       |                               |
| Primary                                              | Secondary                          |                                                  |                               |
| Wavelength: <b>340</b>                               | ††                                 | Main: <b>20 – 24</b>                             |                               |
| Last required read: <b>24</b>                        |                                    | Flex: — —                                        |                               |
| Absorbance range: <b>-0.1000 – 3.2000</b>            |                                    | Color correction: — —                            |                               |
| Sample blank type: <b>None</b>                       |                                    |                                                  |                               |

| Configure assay parameters — Results |                                   |                                 |                                          |
|--------------------------------------|-----------------------------------|---------------------------------|------------------------------------------|
| <input type="radio"/> General        | <input type="radio"/> Calibration | <input type="radio"/> SmartWash | <input checked="" type="radio"/> Results |
| Assay: <b>Opi1SQ</b>                 |                                   | Assay number: <b>2841</b>       |                                          |
| Dilution default range:              |                                   | Result units: <b>ng/mL</b>      |                                          |
| Low-Linearity:†† <b>10.3</b>         |                                   |                                 |                                          |
| High-Linearity: <b>1000.0</b>        |                                   |                                 |                                          |

| Configure assay parameters — Calibration               |                                              |                                 |                                       |
|--------------------------------------------------------|----------------------------------------------|---------------------------------|---------------------------------------|
| <input type="radio"/> General                          | <input checked="" type="radio"/> Calibration | <input type="radio"/> SmartWash | <input type="radio"/> Results         |
| Assay: <b>Opi1SQ</b> Calibration method: <b>Linear</b> |                                              |                                 |                                       |
| <input checked="" type="radio"/> Calibrators           | <input type="radio"/> Volumes                | <input type="radio"/> Intervals | <input type="radio"/> Validity checks |
| Calibrator set: <b>OPSQ</b>                            | Calibrator level: (†††)                      | Concentration:‡                 |                                       |
|                                                        | Blank: <b>OPSQ1 (DOA MC Neg Cal)</b>         | <b>0.0</b>                      |                                       |
|                                                        | Cal 1: <b>OPSQ2 (Opi Cal 150)</b>            |                                 |                                       |
|                                                        | Cal 2: <b>OPSQ3 (Opi Cal 300)</b>            |                                 |                                       |
|                                                        | Cal 3: <b>OPSQ4 (Opi Cal 500)</b>            |                                 |                                       |
|                                                        | Cal 4: <b>OPSQ5 (Opi Cal 1000)</b>           |                                 |                                       |
| Replicates: <b>2</b>                                   | [Range 1 – 3]                                |                                 |                                       |

| Configure assay parameters — Results |                                   |                                 |                                          |
|--------------------------------------|-----------------------------------|---------------------------------|------------------------------------------|
| <input type="radio"/> General        | <input type="radio"/> Calibration | <input type="radio"/> SmartWash | <input checked="" type="radio"/> Results |
| Assay: <b>Opi1SQ</b>                 |                                   | Assay number: <b>2841</b>       |                                          |
| Dilution default range:              |                                   | Result units: <b>ng/mL</b>      |                                          |
| Low-Linearity:†† <b>10.3</b>         |                                   |                                 |                                          |
| High-Linearity: <b>1000.0</b>        |                                   |                                 |                                          |

| Configure assay parameters — Calibration               |                                              |                                 |                                       |
|--------------------------------------------------------|----------------------------------------------|---------------------------------|---------------------------------------|
| <input type="radio"/> General                          | <input checked="" type="radio"/> Calibration | <input type="radio"/> SmartWash | <input type="radio"/> Results         |
| Assay: <b>Opi1SQ</b> Calibration method: <b>Linear</b> |                                              |                                 |                                       |
| <input checked="" type="radio"/> Calibrators           | <input type="radio"/> Volumes                | <input type="radio"/> Intervals | <input type="radio"/> Validity checks |
| Calibrator set: <b>OPSQ</b>                            | Calibrator level: (†††)                      | Concentration:‡                 |                                       |
|                                                        | Blank: <b>OPSQ1 (DOA MC Neg Cal)</b>         | <b>0.0</b>                      |                                       |
|                                                        | Cal 1: <b>OPSQ2 (Opi Cal 150)</b>            |                                 |                                       |
|                                                        | Cal 2: <b>OPSQ3 (Opi Cal 300)</b>            |                                 |                                       |
|                                                        | Cal 3: <b>OPSQ4 (Opi Cal 500)</b>            |                                 |                                       |
|                                                        | Cal 4: <b>OPSQ5 (Opi Cal 1000)</b>           |                                 |                                       |
| Replicates: <b>2</b>                                   | [Range 1 – 3]                                |                                 |                                       |

| Configure assay parameters — Calibration               |                                      |                                 |                                       |
|--------------------------------------------------------|--------------------------------------|---------------------------------|---------------------------------------|
| <input type="radio"/> General                          | <input type="radio"/> Calibration    | <input type="radio"/> SmartWash | <input type="radio"/> Results         |
| Assay: <b>Opi1SQ</b> Calibration method: <b>Linear</b> |                                      |                                 |                                       |
| <input checked="" type="radio"/> Calibrators           | <input type="radio"/> Volumes        | <input type="radio"/> Intervals | <input type="radio"/> Validity checks |
| Calibrator set: <b>OPSQ</b>                            | Calibrator level: (†††)              | Concentration:‡                 |                                       |
|                                                        | Blank: <b>OPSQ1 (DOA MC Neg Cal)</b> | <b>0.0</b>                      |                                       |
|                                                        | Cal 1: <b>OPSQ2 (Opi Cal 150)</b>    |                                 |                                       |
|                                                        | Cal 2: <b>OPSQ3 (Opi Cal 300)</b>    |                                 |                                       |
|                                                        | Cal 3: <b>OPSQ4 (Opi Cal 500)</b>    |                                 |                                       |
|                                                        | Cal 4: <b>OPSQ5 (Opi Cal 1000)</b>   |                                 |                                       |
| Replicates: <b>2</b>                                   | [Range 1 – 3]                        |                                 |                                       |

| Configure assay parameters — Calibration               |                                      |                                 |                                       |
|--------------------------------------------------------|--------------------------------------|---------------------------------|---------------------------------------|
| <input type="radio"/> General                          | <input type="radio"/> Calibration    | <input type="radio"/> SmartWash | <input type="radio"/> Results         |
| Assay: <b>Opi1SQ</b> Calibration method: <b>Linear</b> |                                      |                                 |                                       |
| <input checked="" type="radio"/> Calibrators           | <input type="radio"/> Volumes        | <input type="radio"/> Intervals | <input type="radio"/> Validity checks |
| Calibrator set: <b>OPSQ</b>                            | Calibrator level: (†††)              | Concentration:‡                 |                                       |
|                                                        | Blank: <b>OPSQ1 (DOA MC Neg Cal)</b> | <b>0.0</b>                      |                                       |
|                                                        | Cal 1: <b>OPSQ2 (Opi Cal 150)</b>    |                                 |                                       |
|                                                        | Cal 2: <b>OPSQ3 (Opi Cal 300)</b>    |                                 |                                       |
|                                                        | Cal 3: <b>OPSQ4 (Opi Cal 500)</b>    |                                 |                                       |
|                                                        | Cal 4: <b>OPSQ5 (Opi Cal 1000)</b>   |                                 |                                       |
| Replicates: <b>2</b>                                   | [Range 1 – 3]                        |                                 |                                       |

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| Configure assay parameters — General                 |                                        |                                       |                               |
|------------------------------------------------------|----------------------------------------|---------------------------------------|-------------------------------|
| <input checked="" type="radio"/> General             | <input type="radio"/> Calibration      | <input type="radio"/> SmartWash       | <input type="radio"/> Results |
| Assay: Opi2Q                                         | Type: Photometric                      | Version: †                            |                               |
| Number: 2844                                         | Assay availability: Enabled            |                                       |                               |
| Run controls for onboard reagents by: Lot            |                                        |                                       |                               |
| <input checked="" type="radio"/> Reaction definition | <input type="radio"/> Reagent / Sample | <input type="radio"/> Validity checks |                               |
| Reaction mode: Rate up                               | Primary Secondary                      | Read times                            |                               |
| Wavelength: 340 / ‡‡                                 | Main: 20 – 24                          | Flex: —                               |                               |
| Last required read: 24                               | Color correction: —                    |                                       |                               |
| Absorbance range: -0.1000 – 3.2000                   |                                        |                                       |                               |
| Sample blank type: None                              |                                        |                                       |                               |

| O Reaction definition         |                | ● Reagent / Sample    |                 | O Validity checks |  |
|-------------------------------|----------------|-----------------------|-----------------|-------------------|--|
| Reagent: OPIM9                |                | Reagent volume: 100   |                 | R1 R2             |  |
| Diluent: None                 |                | Water volume: —       |                 |                   |  |
| Diluent dispense mode: Type 0 |                | Dispense mode: Type 0 |                 | Type 0            |  |
| Dilution name Sample          | Diluted sample | Diluent Water         | Dilution factor | Default dilution  |  |
| STANDARD : 3.2                | —              | —                     | 1:1.00          | ●                 |  |

| O Reaction definition |  | O Reagent / Sample |  | ● Validity checks |  |
|-----------------------|--|--------------------|--|-------------------|--|
| Reaction check: None  |  |                    |  |                   |  |
| Rate linearity %: —   |  |                    |  |                   |  |

| Configure assay parameters — Calibration     |                                              |                                 |                                       |
|----------------------------------------------|----------------------------------------------|---------------------------------|---------------------------------------|
| <input type="radio"/> General                | <input checked="" type="radio"/> Calibration | <input type="radio"/> SmartWash | <input type="radio"/> Results         |
| Assay: Opi2Q                                 | Calibration method: Linear                   |                                 |                                       |
| <input checked="" type="radio"/> Calibrators | <input type="radio"/> Volumes                | <input type="radio"/> Intervals | <input type="radio"/> Validity checks |
| Calibrator set: DOAQ                         | Calibrator level: (†††)                      | Concentration: †                |                                       |
|                                              | Blank: DOAQ1 (DOA MC Neg Cal)                | 0                               |                                       |
|                                              | Cal 1: DOAQ2 (DOA MC I Cal 2)                |                                 |                                       |
| Replicates: 2 [Range 1 – 3]                  |                                              |                                 |                                       |

| O Calibrators    |  | ● Volumes        |        | O Intervals    |         | O Validity checks |  |
|------------------|--|------------------|--------|----------------|---------|-------------------|--|
| Calibrator: DOAQ |  | Calibrator level | Sample | Diluted sample | Diluent | Water             |  |
| Blank: DOAQ1     |  | 3.2              | —      | —              | —       | —                 |  |
| Cal 1: DOAQ2     |  | 3.2              | —      | —              | —       | —                 |  |

| O Calibrators              |  | O Volumes |  | ● Intervals |  | O Validity checks |  |
|----------------------------|--|-----------|--|-------------|--|-------------------|--|
| Calibration intervals:     |  |           |  |             |  |                   |  |
| Full interval: 312 (hours) |  |           |  |             |  |                   |  |
| Calibration type:          |  |           |  |             |  |                   |  |
| Adjust type: None          |  |           |  |             |  |                   |  |

| O Calibrators                      |  | O Volumes |  | O Intervals |  | ● Validity checks |  |
|------------------------------------|--|-----------|--|-------------|--|-------------------|--|
| Blank absorbance range: —          |  |           |  |             |  |                   |  |
| Span: Blank – Blank                |  |           |  |             |  |                   |  |
| Span absorbance range: —           |  |           |  |             |  |                   |  |
| Expected cal factor: 0.00          |  |           |  |             |  |                   |  |
| Expected cal factor tolerance %: 0 |  |           |  |             |  |                   |  |

| Configure assay parameters — SmartWash                        |                                   |                                            |                               |                                      |
|---------------------------------------------------------------|-----------------------------------|--------------------------------------------|-------------------------------|--------------------------------------|
| <input type="radio"/> General                                 | <input type="radio"/> Calibration | <input checked="" type="radio"/> SmartWash | <input type="radio"/> Results | <input type="radio"/> Interpretation |
| Assay: Opi2Q                                                  |                                   |                                            |                               |                                      |
| COMPONENT                                                     | REAGENT / ASSAY                   | WASH                                       | Volume                        | Replicates                           |
| Cuvette                                                       | Trig**                            | 10% Detergent B                            | 345                           |                                      |
| **Not required for ARCHITECT Software version 7.00 and above. |                                   |                                            |                               |                                      |

| Configure assay parameters — Results |                                   |                                 |                                          |                                      |
|--------------------------------------|-----------------------------------|---------------------------------|------------------------------------------|--------------------------------------|
| <input type="radio"/> General        | <input type="radio"/> Calibration | <input type="radio"/> SmartWash | <input checked="" type="radio"/> Results | <input type="radio"/> Interpretation |
| Assay: Opi2Q                         |                                   |                                 |                                          |                                      |
| Assay number: 2844                   |                                   |                                 |                                          |                                      |
| Dilution default range: —            |                                   |                                 |                                          |                                      |
| Low-Linearity: —                     |                                   |                                 |                                          |                                      |
| High-Linearity: —                    |                                   |                                 |                                          |                                      |

| Configure assay parameters — Interpretation |                                   |                                 |                               |                                                 |
|---------------------------------------------|-----------------------------------|---------------------------------|-------------------------------|-------------------------------------------------|
| <input type="radio"/> General               | <input type="radio"/> Calibration | <input type="radio"/> SmartWash | <input type="radio"/> Results | <input checked="" type="radio"/> Interpretation |
| Assay: Opi2Q                                |                                   |                                 |                               |                                                 |
| Assay number: 2844                          |                                   |                                 |                               |                                                 |
| Name:                                       | Range:                            | Results review required:        |                               |                                                 |
| Negative                                    | < 2000                            | <input type="checkbox"/>        |                               |                                                 |
| Positive                                    | >= 2000                           | <input type="checkbox"/>        |                               |                                                 |

| Configure result units     |               |
|----------------------------|---------------|
| Assay: Opi2Q               |               |
| Version: †                 |               |
| Result units: —            |               |
| Decimal places: 0          | [Range 0 – 4] |
| Correlation factor: 1.0000 |               |
| Intercept: 0.0000          |               |

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| Configure assay parameters — General      |                             |                   |           |
|-------------------------------------------|-----------------------------|-------------------|-----------|
| ● General                                 | ○ Calibration               | ○ SmartWash       | ○ Results |
| Assay: Opi6SQ                             | Type: Photometric           | Version: †        |           |
| Number: 2842                              | Assay availability: Enabled |                   |           |
| Run controls for onboard reagents by: Lot |                             |                   |           |
| ● Reaction definition                     | ○ Reagent / Sample          | ○ Validity checks |           |
| Reaction mode: Rate up                    |                             |                   |           |
| Primary                                   | Secondary                   | Read times        |           |
| Wavelength: 340 / ††                      | Main: 20 – 24               |                   |           |
| Last required read: 24                    | Flex: — —                   |                   |           |
| Absorbance range: -0.1000 – 3.2000        | Color correction: — —       |                   |           |
| Sample blank type: None                   |                             |                   |           |

| ○ Reaction definition         |                |                       |        | ● Reagent / Sample |                  | ○ Validity checks |  |
|-------------------------------|----------------|-----------------------|--------|--------------------|------------------|-------------------|--|
| Reagent: OPIM9                |                | Reagent volume: 100   | R1     | R2                 |                  |                   |  |
| Diluent: None                 |                | Water volume: —       |        |                    |                  |                   |  |
| Diluent dispense mode: Type 0 |                | Dispense mode: Type 0 | Type 0 | Type 0             |                  |                   |  |
| Dilution name Sample          | Diluted sample | Diluent               | Water  | Dilution factor    | Default dilution |                   |  |
| STANDARD : 3.2                | —              | —                     | —      | 1:1.00             | ●                |                   |  |

| ○ Reaction definition |  | ○ Reagent / Sample |  | ● Validity checks |  |
|-----------------------|--|--------------------|--|-------------------|--|
| Reaction check: None  |  |                    |  |                   |  |
| Rate linearity %: —   |  |                    |  |                   |  |

| Configure assay parameters — Calibration |                                |                  |                   |
|------------------------------------------|--------------------------------|------------------|-------------------|
| ○ General                                | ● Calibration                  | ○ SmartWash      | ○ Results         |
| Assay: Opi6SQ                            | Calibration method: Linear     |                  |                   |
| ● Calibrators                            | ○ Volumes                      | ○ Intervals      | ○ Validity checks |
| Calibrator set: DOASQ                    | Calibrator level: (†††)        | Concentration: † |                   |
|                                          | Blank: DOASQ1 (DOA MC Neg Cal) | 0.0              |                   |
|                                          | Cal 1: DOASQ2 (DOA MC I Cal 1) |                  |                   |
|                                          | Cal 2: DOASQ3 (DOA MC I Cal 2) |                  |                   |
|                                          | Cal 3: DOASQ4 (DOA MC I Cal 3) |                  |                   |
|                                          | Cal 4: DOASQ5 (DOA MC I Cal 4) |                  |                   |
| Replicates: 2                            | [Range 1 – 3]                  |                  |                   |

| ○ Calibrators     |               | ● Volumes        |        | ○ Intervals    |         | ○ Validity checks |  |
|-------------------|---------------|------------------|--------|----------------|---------|-------------------|--|
| Calibrator: DOASQ |               | Calibrator level | Sample | Diluted sample | Diluent | Water             |  |
|                   | Blank: DOASQ1 | 3.2              | —      | —              | —       | —                 |  |
|                   | Cal 1: DOASQ2 | 3.2              | —      | —              | —       | —                 |  |
|                   | Cal 2: DOASQ3 | 3.2              | —      | —              | —       | —                 |  |
|                   | Cal 3: DOASQ4 | 3.2              | —      | —              | —       | —                 |  |
|                   | Cal 4: DOASQ5 | 3.2              | —      | —              | —       | —                 |  |

| ○ Calibrators              |  | ○ Volumes |  | ● Intervals |  | ○ Validity checks |  |
|----------------------------|--|-----------|--|-------------|--|-------------------|--|
| Calibration intervals:     |  |           |  |             |  |                   |  |
| Full interval: 312 (hours) |  |           |  |             |  |                   |  |
| Calibration type:          |  |           |  |             |  |                   |  |
| Adjust type: None          |  |           |  |             |  |                   |  |

| ○ Calibrators                      |  | ○ Volumes |  | ○ Intervals |  | ● Validity checks |  |
|------------------------------------|--|-----------|--|-------------|--|-------------------|--|
| Blank absorbance range:            |  | —         |  | —           |  |                   |  |
| Span: Blank                        |  | Blank     |  | Blank       |  |                   |  |
| Span absorbance range:             |  | —         |  | —           |  |                   |  |
| Expected cal factor: 0.00          |  |           |  |             |  |                   |  |
| Expected cal factor tolerance %: 0 |  |           |  |             |  |                   |  |

| Configure assay parameters — SmartWash                        |                 |                 |           |                  |
|---------------------------------------------------------------|-----------------|-----------------|-----------|------------------|
| ○ General                                                     | ○ Calibration   | ● SmartWash     | ○ Results | ○ Interpretation |
| Assay: Opi6SQ                                                 |                 |                 |           |                  |
| COMPONENT                                                     | REAGENT / ASSAY | WASH            | Volume    | Replicates       |
| Cuvette                                                       | Trig**          | 10% Detergent B | 345       |                  |
| **Not required for ARCHITECT Software version 7.00 and above. |                 |                 |           |                  |

| Configure assay parameters — Results |               |             |                     |                  |
|--------------------------------------|---------------|-------------|---------------------|------------------|
| ○ General                            | ○ Calibration | ○ SmartWash | ● Results           | ○ Interpretation |
| Assay: Opi6SQ                        |               |             | Assay number: 2842  |                  |
| Dilution default range:              |               |             | Result units: ng/mL |                  |
| Low-Linearity:††                     | 44.1          |             |                     |                  |
| High-Linearity:                      | 6000.0        |             |                     |                  |

| Configure result units     |               |
|----------------------------|---------------|
| Assay: Opi6SQ              |               |
| Version: †                 |               |
| Result units: ng/mL        |               |
| Decimal places: 1          | [Range 0 – 4] |
| Correlation factor: 1.0000 |               |
| Intercept: 0.0000          |               |

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1. US Department of Labor, Occupational Safety and Health Administration. 29 CFR Part 1910.1030. *Bloodborne Pathogens*.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*, 5th ed. Washington, DC: US Government Printing Office, January 2007.
3. World Health Organization. *Laboratory Biosafety Manual*, 3rd ed. Geneva: World Health Organization, 2004.
4. Sewell DL, Bove KE, Callihan DR, et al. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Third Edition (M29-A3)*. Wayne, PA: Clinical and Laboratory Standards Institute, 2005.

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# ARCHITECT cSystems Assay Parameters

## Ethanol

| Configure assay parameters — General                 |                             |                                        |  |
|------------------------------------------------------|-----------------------------|----------------------------------------|--|
| <input checked="" type="radio"/> General             |                             | <input type="radio"/> Calibration      |  |
| <input type="radio"/> SmartWash                      |                             | <input type="radio"/> Results          |  |
| <input type="radio"/> Interpretation                 |                             |                                        |  |
| Assay: EtOHM                                         | Type: Photometric           | Version: †                             |  |
| Number: 2847                                         | Assay availability: Enabled |                                        |  |
| <input checked="" type="radio"/> Reaction definition |                             | <input type="radio"/> Reagent / Sample |  |
| <input type="radio"/> Validity checks                |                             |                                        |  |
| Reaction mode: Rate up                               |                             | Read times                             |  |
| Primary Secondary                                    |                             | Main: 20 – 24                          |  |
| Wavelength: 340 / ‡‡                                 |                             | Flex: — —                              |  |
| Last required read: 24                               |                             | Color correction: — —                  |  |
| Absorbance range: -0.1000 – 3.2000                   |                             |                                        |  |
| Sample blank type: None                              |                             |                                        |  |

| <input type="radio"/> Reaction definition |        | <input checked="" type="radio"/> Reagent / Sample |         | <input type="radio"/> Validity checks            |                 |
|-------------------------------------------|--------|---------------------------------------------------|---------|--------------------------------------------------|-----------------|
| Reagent: ETOHM                            |        | Reagent volume: 100                               |         | R1 R2                                            |                 |
| Diluent: None                             |        | Water volume: —                                   |         | 100                                              |                 |
| Diluent dispense mode: Type 0             |        | Dispense mode: Type 0                             |         | Type 0                                           |                 |
| Dilution name                             | Sample | Diluted sample                                    | Diluent | Water                                            | Dilution factor |
| STANDARD                                  | 12.0   | —                                                 | —       | —                                                | 1:1.00          |
| <input type="radio"/> Reaction definition |        | <input type="radio"/> Reagent / Sample            |         | <input checked="" type="radio"/> Validity checks |                 |
| Reaction check: None                      |        |                                                   |         |                                                  |                 |
| Rate linearity %: —                       |        |                                                   |         |                                                  |                 |

| Configure assay parameters — Calibration     |  |                                              |  |
|----------------------------------------------|--|----------------------------------------------|--|
| <input type="radio"/> General                |  | <input checked="" type="radio"/> Calibration |  |
| <input type="radio"/> SmartWash              |  | <input type="radio"/> Results                |  |
| <input type="radio"/> Interpretation         |  |                                              |  |
| Assay: EtOHM                                 |  | Calibration method: Linear                   |  |
| <input checked="" type="radio"/> Calibrators |  | <input type="radio"/> Volumes                |  |
| <input type="radio"/> Intervals              |  | <input type="radio"/> Validity checks        |  |
| Calibrator set: ETH                          |  | Calibrator level: (†††)                      |  |
| Blank: ETH1 (Ethanol Neg Cal)                |  | Concentration: †                             |  |
| Cal 1: ETH2 (Ethanol 100 Cal)                |  | 0.0                                          |  |
| Replicates: 2 [Range 1 – 3]                  |  |                                              |  |

| <input type="radio"/> Calibrators |  | <input checked="" type="radio"/> Volumes |  | <input type="radio"/> Intervals |  | <input type="radio"/> Validity checks |  |
|-----------------------------------|--|------------------------------------------|--|---------------------------------|--|---------------------------------------|--|
| Calibrator: ETH                   |  | Calibrator level                         |  | Sample                          |  | Diluted sample                        |  |
| Blank: ETH1                       |  | 12.0                                     |  | —                               |  | —                                     |  |
| Cal 1: ETH2                       |  | 12.0                                     |  | —                               |  | —                                     |  |

| <input type="radio"/> Calibrators |  | <input type="radio"/> Volumes |  | <input checked="" type="radio"/> Intervals |  | <input type="radio"/> Validity checks |  |
|-----------------------------------|--|-------------------------------|--|--------------------------------------------|--|---------------------------------------|--|
| Calibration intervals:            |  | Full interval: 312 (hours)    |  | Calibration type:                          |  | Adjust type: None                     |  |

| <input type="radio"/> Calibrators |  | <input type="radio"/> Volumes |  | <input type="radio"/> Intervals |  | <input checked="" type="radio"/> Validity checks |  |
|-----------------------------------|--|-------------------------------|--|---------------------------------|--|--------------------------------------------------|--|
| Blank absorbance range:           |  | Span: Blank – Blank           |  | Span absorbance range:          |  | Expected cal factor: 0.00                        |  |
| Expected cal factor tolerance %:  |  | 0                             |  |                                 |  |                                                  |  |

| Configure assay parameters — SmartWash |                 |                                      |        |                                            |
|----------------------------------------|-----------------|--------------------------------------|--------|--------------------------------------------|
| <input type="radio"/> General          |                 | <input type="radio"/> Calibration    |        | <input checked="" type="radio"/> SmartWash |
| <input type="radio"/> Results          |                 | <input type="radio"/> Interpretation |        |                                            |
| Assay: EtOHM                           |                 |                                      |        |                                            |
| COMPONENT                              | REAGENT / ASSAY | WASH                                 | Volume | Replicates                                 |
| Cuvette                                | Trig            | 10% Detergent B                      | 345    |                                            |

| Configure assay parameters — Results |  |                                          |  |
|--------------------------------------|--|------------------------------------------|--|
| <input type="radio"/> General        |  | <input type="radio"/> Calibration        |  |
| <input type="radio"/> SmartWash      |  | <input checked="" type="radio"/> Results |  |
| <input type="radio"/> Interpretation |  |                                          |  |
| Assay: EtOHM                         |  | Assay number: 2847                       |  |
| Dilution default range:              |  | Result units: mg/dL                      |  |
| Low-Linearity: ††                    |  | 10.0                                     |  |
| High-Linearity:                      |  | 600.0                                    |  |

| Configure result units     |                                 |
|----------------------------|---------------------------------|
| Assay: EtOHM               | Version: †                      |
| Result units: mg/dL        | Decimal places: 1 [Range 0 – 4] |
| Correlation factor: 1.0000 | Intercept: 0.0000               |

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- Baselt RC. *Disposition of Toxic Drugs and Chemicals in Man*, 5th ed. Foster City, CA: Chemical Toxicology Institute; 2000:323-6.
- Beutler HO: Ethanol. In: Bergmeyer HU, editor. *Methods of Enzymatic Analysis*, Vol. VI, 3rd ed. New York: Academic Press, 1984:598-606.
- Redetzki HM, Dees WL. Comparison of four kits for enzymatic determination of ethanol in blood. *Clin Chem* 1976;22, 83-6.
- US Department of Labor, Occupational Safety and Health Administration. 29 CFR Part 1910.1030. *Bloodborne Pathogens*.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*, 5th ed. Washington, DC: US Government Printing Office, January 2007.
- World Health Organization. *Laboratory Biosafety Manual*, 3rd ed. Geneva: World Health Organization. 2004
- Sewell DL, Bove KE, Callihan DR, et al. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Third Edition (M29-A3)*. Wayne, PA: Clinical and Laboratory Standards Institute, 2005.
- Guder WG, da Fonseca-Wollheim F, Heil W, et al. *The Quality of Diagnostic Samples*. Darmstadt, Germany: GIT Verlag; 2001:28, 50.
- Heise HA. Concentrations of alcohol in samples of blood and urine taken at the same time. *J For Sci* 1967;12, 454.
- Jones AW. Urine as a biological specimen for forensic analysis of alcohol and variability in the urine-to-blood relationship. *Toxicol Rev* 2006;25:15-35.

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1. Hawks RL, Chiang CN, editors. *Urine Testing for Drugs of Abuse*. NIDA Research Monograph 73. Rockville, MD: National Institute on Drug Abuse (NIDA), Department of Health and Human Services; 1986.
2. *Mandatory Guidelines for Federal Workplace Drug Testing Programs*. National Institute on Drug Abuse. Federal Register Vol. 69, No. 71, 2004:19650–61.
3. Pohland A, Boaz HE, Sullivan HR. Synthesis and identification of metabolites resulting from the biotransformation of *d,l*-methadone in man and in rat. *J Med Chem* 1971;14:194–7.
4. Baselt RC, Casarett LJ. Urinary excretion of methadone in man. *Clin Pharm Therap* 1972;13:64–70.
5. Baselt RC. *Disposition of Toxic Drugs and Chemicals in Man*, 5th ed. Foster City, CA: Chemical Toxicology Institute; 2000:523-7.
6. Rubenstein KE, Schneider RS, Ullman EF. Homogeneous enzyme immunoassay: a new immunochemical technique. *Biochem Biophys Res Commun* 1972;47:846-51.
7. US Department of Labor, Occupational Safety and Health Administration. 29 CFR Part 1910.1030. *Bloodborne Pathogens*.
8. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*, 5th ed. Washington, DC: US Government Printing Office, January 2007.
9. World Health Organization. *Laboratory Biosafety Manual*, 3rd ed. Geneva: World Health Organization, 2004.
10. Sewell DL, Bove KE, Callihan DR, et al. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Third Edition (M29-A3)*. Wayne, PA: Clinical and Laboratory Standards Institute, 2005.
11. Shoemaker MJ, Earley RJ, Gorsky JE, et al. *Urine Drug Testing in the Clinical Laboratory; Approved Guideline—Vol. 19, No. 6*. Wayne, PA: Clinical and Laboratory Standards Institute, 1999(17).

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| Configure assay parameters — General                                                                                                                                                                                                                      |       |                     |             |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------|---------------------|-------------|
| <input checked="" type="radio"/> General <input type="radio"/> Calibration <input type="radio"/> SmartWash <input type="radio"/> Results <input type="radio"/> Interpretation                                                                             |       |                     |             |
| Assay:                                                                                                                                                                                                                                                    | MethQ | Type:               | Photometric |
| Number:                                                                                                                                                                                                                                                   | 2846  | Assay availability: | Enabled     |
| <input checked="" type="radio"/> Reaction definition <input type="radio"/> Reagent / Sample <input type="radio"/> Validity checks                                                                                                                         |       |                     |             |
| Reaction mode: <b>Rate up</b><br>Primary      Secondary      Read times<br>Wavelength: 340 / <b>††</b> Main: 20 – 24<br>Last required read: 24      Flex: — —<br>Absorbance range: -0.1000 – 3.2000      Color correction: — —<br>Sample blank type: None |       |                     |             |

| Configure assay parameters — General                                                                                                                                                                                                                                                                                                                                       |  |  |  |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|--|--|
| <input type="radio"/> Reaction definition <input checked="" type="radio"/> Reagent / Sample <input type="radio"/> Validity checks                                                                                                                                                                                                                                          |  |  |  |
| Reagent: <b>METHM</b> Reagent volume: 100      R1      R2<br>Diluent: <b>None</b> Water volume: —<br>Diluent dispense mode: <b>Type 0</b> Dispense mode: <b>Type 0</b> Type 0<br>Dilution name      Sample      Diluted sample      Diluent      Water      Dilution factor      Default dilution<br><b>STANDARD</b> : 12.0      —      —      —      =      1:1.00      ● |  |  |  |
| <input type="radio"/> Reaction definition <input type="radio"/> Reagent / Sample <input checked="" type="radio"/> Validity checks                                                                                                                                                                                                                                          |  |  |  |
| Reaction check: <b>None</b><br>Rate linearity %: —                                                                                                                                                                                                                                                                                                                         |  |  |  |

| Configure assay parameters — Calibration                                                                                                                                                                |  |  |  |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|--|--|
| <input type="radio"/> General <input checked="" type="radio"/> Calibration <input type="radio"/> SmartWash <input type="radio"/> Results <input type="radio"/> Interpretation                           |  |  |  |
| Assay: <b>MethQ</b> Calibration method: <b>Linear</b>                                                                                                                                                   |  |  |  |
| <input checked="" type="radio"/> Calibrators <input type="radio"/> Volumes <input type="radio"/> Intervals <input type="radio"/> Validity checks                                                        |  |  |  |
| Calibrator set: <b>DOAQ</b> Calibrator level: ( <b>†††</b> )      Concentration:†<br>Blank: <b>DOAQ1</b> (DOA MC Neg Cal)      0<br>Cal 1: <b>DOAQ2</b> (DOA MC I Cal 2)<br>Replicates: 2 [Range 1 – 3] |  |  |  |
| <input type="radio"/> Calibrators <input checked="" type="radio"/> Volumes <input type="radio"/> Intervals <input type="radio"/> Validity checks                                                        |  |  |  |
| Calibrator: <b>DOAQ</b><br>Calibrator level      Sample      Diluted sample      Diluent      Water<br>Blank: <b>DOAQ1</b> 12.0      —      —      —<br>Cal 1: <b>DOAQ2</b> 12.0      —      —      —   |  |  |  |
| <input type="radio"/> Calibrators <input type="radio"/> Volumes <input checked="" type="radio"/> Intervals <input type="radio"/> Validity checks                                                        |  |  |  |
| Calibration intervals:<br>Full interval: <b>240</b> (hours)<br>Calibration type:<br>Adjust type: <b>None</b>                                                                                            |  |  |  |
| <input type="radio"/> Calibrators <input type="radio"/> Volumes <input type="radio"/> Intervals <input checked="" type="radio"/> Validity checks                                                        |  |  |  |
| Blank absorbance range: — —<br>Span: <b>Blank</b> — <b>Blank</b><br>Span absorbance range: — —<br>Expected cal factor: <b>0.00</b><br>Expected cal factor tolerance %: <b>0</b>                         |  |  |  |

| Configure assay parameters — SmartWash                                                                                                                                        |                 |                 |        |            |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------|-----------------|--------|------------|
| <input type="radio"/> General <input type="radio"/> Calibration <input checked="" type="radio"/> SmartWash <input type="radio"/> Results <input type="radio"/> Interpretation |                 |                 |        |            |
| Assay: <b>MethQ</b>                                                                                                                                                           |                 |                 |        |            |
| COMPONENT                                                                                                                                                                     | REAGENT / ASSAY | WASH            | Volume | Replicates |
| Sample Probe <sup>†††</sup>                                                                                                                                                   |                 | 0.5% Acid Wash  |        |            |
| Cuvette                                                                                                                                                                       | Trig            | 10% Detergent B | 345    |            |
| <sup>†††</sup> Sample probe Sample wash protocol is <b>Maximum wash</b> .                                                                                                     |                 |                 |        |            |

| Configure assay parameters — Results                                                                                                                                          |  |  |  |  |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|--|--|--|
| <input type="radio"/> General <input type="radio"/> Calibration <input type="radio"/> SmartWash <input checked="" type="radio"/> Results <input type="radio"/> Interpretation |  |  |  |  |
| Assay: <b>MethQ</b> Assay number: <b>2846</b>                                                                                                                                 |  |  |  |  |
| Dilution default range:<br>Low-Linearity:<br>High-Linearity:                                                                                                                  |  |  |  |  |

| Configure assay parameters — Interpretation                                                                                                                                   |  |        |     |                          |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|--------|-----|--------------------------|
| <input type="radio"/> General <input type="radio"/> Calibration <input type="radio"/> SmartWash <input type="radio"/> Results <input checked="" type="radio"/> Interpretation |  |        |     |                          |
| Assay: <b>MethQ</b> Assay number: <b>2846</b>                                                                                                                                 |  |        |     |                          |
| Name:                                                                                                                                                                         |  | Range: |     | Results review required: |
| <b>Negative</b>                                                                                                                                                               |  | <      | 300 | <input type="checkbox"/> |
| <b>Positive</b>                                                                                                                                                               |  | >=     | 300 | <input type="checkbox"/> |

| Configure result units |                        |
|------------------------|------------------------|
| Assay:                 | <b>MethQ</b>           |
| Version:               | <b>†</b>               |
| Result units:          | —                      |
| Decimal places:        | <b>0</b> [Range 0 – 4] |
| Correlation factor:    | <b>1.0000</b>          |
| Intercept:             | <b>0.0000</b>          |

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| Configure assay parameters — General                 |                                    |                                        |  |
|------------------------------------------------------|------------------------------------|----------------------------------------|--|
| <input checked="" type="radio"/> General             |                                    | <input type="radio"/> Calibration      |  |
| <input type="radio"/> SmartWash                      |                                    | <input type="radio"/> Results          |  |
| <input type="radio"/> Interpretation                 |                                    |                                        |  |
| Assay: <b>MethSQ</b>                                 | Type: <b>Photometric</b>           | Version: †                             |  |
| Number: <b>2845</b>                                  | Assay availability: <b>Enabled</b> |                                        |  |
| <input checked="" type="radio"/> Reaction definition |                                    | <input type="radio"/> Reagent / Sample |  |
| <input type="radio"/> Validity checks                |                                    |                                        |  |
| Reaction mode: <b>Rate up</b>                        |                                    |                                        |  |
| Primary                                              |                                    | Secondary                              |  |
| Wavelength: <b>340</b>                               |                                    | Main: <b>20 – 24</b>                   |  |
| Last required read: <b>24</b>                        |                                    | Flex: — —                              |  |
| Absorbance range: <b>-0.1000 – 3.2000</b>            |                                    | Color correction: — —                  |  |
| Sample blank type: <b>None</b>                       |                                    |                                        |  |

| Configure assay parameters — Results                 |  |                                          |  |
|------------------------------------------------------|--|------------------------------------------|--|
| <input type="radio"/> General                        |  | <input type="radio"/> Calibration        |  |
| <input type="radio"/> SmartWash                      |  | <input checked="" type="radio"/> Results |  |
| <input type="radio"/> Interpretation                 |  |                                          |  |
| Assay: <b>MethSQ</b>                                 |  | Assay number: <b>2845</b>                |  |
| Dilution default range:                              |  | Result units: <b>ng/mL</b>               |  |
| Low-Linearity:†† <b>10.0</b>                         |  |                                          |  |
| High-Linearity: <b>1000.0</b>                        |  |                                          |  |
| <input checked="" type="radio"/> Reaction definition |  |                                          |  |
| <input type="radio"/> Reagent / Sample               |  |                                          |  |
| <input type="radio"/> Validity checks                |  |                                          |  |
| Reagent: <b>METHM</b>                                |  |                                          |  |
| Diluent: <b>None</b>                                 |  |                                          |  |
| Diluent dispense mode: <b>Type 0</b>                 |  |                                          |  |
| Dilution name                                        |  |                                          |  |
| Sample                                               |  |                                          |  |
| sample                                               |  |                                          |  |
| Diluent                                              |  |                                          |  |
| Water                                                |  |                                          |  |
| Dilution factor                                      |  |                                          |  |
| 1:1.00                                               |  |                                          |  |
| Default dilution                                     |  |                                          |  |
| ●                                                    |  |                                          |  |
| Reaction check: <b>None</b>                          |  |                                          |  |
| Rate linearity %: —                                  |  |                                          |  |

| Configure assay parameters — Calibration     |  |                                              |  |
|----------------------------------------------|--|----------------------------------------------|--|
| <input type="radio"/> General                |  | <input checked="" type="radio"/> Calibration |  |
| <input type="radio"/> SmartWash              |  | <input type="radio"/> Results                |  |
| <input type="radio"/> Interpretation         |  |                                              |  |
| Assay: <b>MethSQ</b>                         |  | Calibration method: <b>Linear</b>            |  |
| <input checked="" type="radio"/> Calibrators |  | <input type="radio"/> Volumes                |  |
| <input type="radio"/> Intervals              |  | <input type="radio"/> Validity checks        |  |
| Calibrator set: <b>DOASQ</b>                 |  |                                              |  |
| Calibrator level: (†††)                      |  |                                              |  |
| Blank: <b>DOASQ1 (DOA MC Neg Cal)</b>        |  |                                              |  |
| Cal 1: <b>DOASQ2 (DOA MC I Cal 1)</b>        |  |                                              |  |
| Cal 2: <b>DOASQ3 (DOA MC I Cal 2)</b>        |  |                                              |  |
| Cal 3: <b>DOASQ4 (DOA MC I Cal 3)</b>        |  |                                              |  |
| Cal 4: <b>DOASQ5 (DOA MC I Cal 4)</b>        |  |                                              |  |
| Replicates: <b>2</b> [Range 1 – 3]           |  |                                              |  |

| Configure assay parameters — Calibration |  |                                          |  |
|------------------------------------------|--|------------------------------------------|--|
| <input type="radio"/> Calibrators        |  | <input checked="" type="radio"/> Volumes |  |
| <input type="radio"/> Intervals          |  | <input type="radio"/> Validity checks    |  |
| Calibrator: <b>DOASQ</b>                 |  |                                          |  |
| Calibrator level                         |  |                                          |  |
| Sample                                   |  |                                          |  |
| Diluted sample                           |  |                                          |  |
| Diluent                                  |  |                                          |  |
| Water                                    |  |                                          |  |
| Blank: <b>DOASQ1</b>                     |  |                                          |  |
| Cal 1: <b>DOASQ2</b>                     |  |                                          |  |
| Cal 2: <b>DOASQ3</b>                     |  |                                          |  |
| Cal 3: <b>DOASQ4</b>                     |  |                                          |  |
| Cal 4: <b>DOASQ5</b>                     |  |                                          |  |

| Configure assay parameters — Calibration   |  |                                       |  |
|--------------------------------------------|--|---------------------------------------|--|
| <input type="radio"/> Calibrators          |  | <input type="radio"/> Volumes         |  |
| <input checked="" type="radio"/> Intervals |  | <input type="radio"/> Validity checks |  |
| Calibration intervals:                     |  |                                       |  |
| Full interval: <b>240</b> (hours)          |  |                                       |  |
| Calibration type:                          |  |                                       |  |
| Adjust type: <b>None</b>                   |  |                                       |  |

| Configure assay parameters — Calibration  |  |                                                  |  |
|-------------------------------------------|--|--------------------------------------------------|--|
| <input type="radio"/> Calibrators         |  | <input type="radio"/> Volumes                    |  |
| <input type="radio"/> Intervals           |  | <input checked="" type="radio"/> Validity checks |  |
| Blank absorbance range:                   |  |                                                  |  |
| Span: <b>Blank</b> – <b>Blank</b>         |  |                                                  |  |
| Span absorbance range:                    |  |                                                  |  |
| Expected cal factor: <b>0.00</b>          |  |                                                  |  |
| Expected cal factor tolerance %: <b>0</b> |  |                                                  |  |

| Configure assay parameters — SmartWash     |                 |                                   |            |            |
|--------------------------------------------|-----------------|-----------------------------------|------------|------------|
| <input type="radio"/> General              |                 | <input type="radio"/> Calibration |            |            |
| <input checked="" type="radio"/> SmartWash |                 | <input type="radio"/> Results     |            |            |
| <input type="radio"/> Interpretation       |                 |                                   |            |            |
| Assay: <b>MethSQ</b>                       |                 |                                   |            |            |
| COMPONENT                                  | REAGENT / ASSAY | WASH                              | Volume     | Replicates |
| Sample Probe†††                            |                 | <b>0.5% Acid Wash</b>             |            |            |
| Cuvette                                    | <b>Trig</b>     | <b>10% Detergent B</b>            | <b>345</b> |            |

††† Sample probe Sample wash protocol is Maximum wash.

| Configure assay parameters — Results |  |                                          |  |  |
|--------------------------------------|--|------------------------------------------|--|--|
| <input type="radio"/> General        |  | <input type="radio"/> Calibration        |  |  |
| <input type="radio"/> SmartWash      |  | <input checked="" type="radio"/> Results |  |  |
| <input type="radio"/> Interpretation |  |                                          |  |  |
| Assay: <b>MethSQ</b>                 |  | Assay number: <b>2845</b>                |  |  |
| Dilution default range:              |  | Result units: <b>ng/mL</b>               |  |  |
| Low-Linearity:†† <b>10.0</b>         |  |                                          |  |  |
| High-Linearity: <b>1000.0</b>        |  |                                          |  |  |

| Configure result units                 |  |
|----------------------------------------|--|
| Assay: <b>MethSQ</b>                   |  |
| Version: †                             |  |
| Result units: <b>ng/mL</b>             |  |
| Decimal places: <b>1</b> [Range 0 – 4] |  |
| Correlation factor: <b>1.0000</b>      |  |
| Intercept: <b>0.0000</b>               |  |

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- Hawks RL, Chiang CN, editors. *Urine Testing for Drugs of Abuse*. NIDA Research Monograph 73. Rockville, MD: National Institute on Drug Abuse (NIDA), Department of Health and Human Services; 1986.
- Mandatory Guidelines for Federal Workplace Drug Testing Programs*. National Institute on Drug Abuse. Federal Register Vol. 69, No. 71, 2004:19650–61.
- Kalant OJ. *The Amphetamines: Toxicity and Addiction*. Thomas, Springfield. 1966:151.
- Morgan JP. Substance Abuse: *Clinical Problems and Perspectives*. Lowison JH and Ruiz P, editors. Williams & Wilkins, Baltimore, MD. 1981:167.
- Baselt RC. *Disposition of Toxic Drugs and Chemicals in Man*, 5th ed. Foster City, CA: Chemical Toxicology Institute; 2000:49–51, 527–31.
- Respiratory and cerebral stimulants. *AHFS Drug Information '87*. Bethesda, MD: American Society of Hospital Pharmacists, Inc; 1987:1105–7.
- Hawks RL, Chiang CN, editors. *Urine Testing for Drugs of Abuse*. NIDA Research Monograph 73. Rockville, MD: National Institute on Drug Abuse (NIDA), Department of Health and Human Services; 1986:35–6, 95–7.
- Rubenstein KE, Schneider RS, Ullman EF. Homogeneous enzyme immunoassay: a new immunochemical technique. *Biochem Biophys Res Commun* 1972;47:846–51.
- US Department of Labor, Occupational Safety and Health Administration. 29 CFR Part 1910.1030. *Bloodborne Pathogens*.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*, 5th ed. Washington, DC: US Government Printing Office, January 2007.
- World Health Organization. *Laboratory Biosafety Manual*, 3rd ed. Geneva: World Health Organization, 2004.
- Sewell DL, Bove KE, Callihan DR, et al. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Third Edition (M29-A3)*. Wayne, PA: Clinical and Laboratory Standards Institute, 2005.

- Armbruster DA, Dasgupta A, et al. *Toxicology and Drug Testing in the Clinical Laboratory; Approved Guideline—Vol. 27, No. 15*. Wayne, PA: Clinical and Laboratory Standards Institute, 2007.

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# ARCHITECT cSystems Assay Parameters

## Amphetamine/Methamphetamine (Amp5Q) 500 Качественный

| Configure assay parameters — General                                                                                                                                          |  |                                    |            |  |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|------------------------------------|------------|--|
| <input checked="" type="radio"/> General <input type="radio"/> Calibration <input type="radio"/> SmartWash <input type="radio"/> Results <input type="radio"/> Interpretation |  |                                    |            |  |
| Assay: <b>Amp5Q</b>                                                                                                                                                           |  | Type: <b>Photometric</b>           | Version: † |  |
| Number: <b>2879</b>                                                                                                                                                           |  | Assay availability: <b>Enabled</b> | Lot        |  |
| Run controls for onboard reagents by: *                                                                                                                                       |  |                                    |            |  |
| <input checked="" type="radio"/> Reaction definition <input type="radio"/> Reagent / Sample <input type="radio"/> Validity checks                                             |  |                                    |            |  |
| Reaction mode: <b>Rate up</b><br>Primary      Secondary      Read times                                                                                                       |  |                                    |            |  |
| Wavelength: <b>340</b> / ‡‡                                                                                                                                                   |  | Main: <b>20 – 24</b>               |            |  |
| Last required read: <b>24</b>                                                                                                                                                 |  | Flex: — —                          |            |  |
| Absorbance range: <b>-0.1000 – 3.2000</b>                                                                                                                                     |  | Color correction: — —              |            |  |
| Sample blank type: <b>None</b>                                                                                                                                                |  |                                    |            |  |

| Configure assay parameters — Reagent / Sample                                                                                     |            |                                            |         |                                                    |
|-----------------------------------------------------------------------------------------------------------------------------------|------------|--------------------------------------------|---------|----------------------------------------------------|
| <input type="radio"/> Reaction definition <input checked="" type="radio"/> Reagent / Sample <input type="radio"/> Validity checks |            |                                            |         |                                                    |
| Reagent: <b>AMPM9</b>                                                                                                             |            | Reagent volume: <b>100</b> <b>100</b>      |         |                                                    |
| Diluent: <b>None</b>                                                                                                              |            | Water volume: — —                          |         |                                                    |
| Diluent dispense mode: <b>Type 0</b>                                                                                              |            | Dispense mode: <b>Type 0</b> <b>Type 0</b> |         |                                                    |
| Dilution name                                                                                                                     | Sample     | Diluted sample                             | Diluent | Water                                              |
| <b>STANDARD</b>                                                                                                                   | <b>8.0</b> | —                                          | —       | —                                                  |
| Dilution factor: <b>1:1.00</b>                                                                                                    |            |                                            |         | Default dilution: <input checked="" type="radio"/> |

| Configure assay parameters — Validity checks                                                                                      |  |  |  |  |
|-----------------------------------------------------------------------------------------------------------------------------------|--|--|--|--|
| <input type="radio"/> Reaction definition <input type="radio"/> Reagent / Sample <input checked="" type="radio"/> Validity checks |  |  |  |  |
| Reaction check: <b>None</b>                                                                                                       |  |  |  |  |
| Rate linearity %: —                                                                                                               |  |  |  |  |

| Configure assay parameters — Calibration                                                                                                                                      |  |                                               |  |  |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|-----------------------------------------------|--|--|
| <input type="radio"/> General <input checked="" type="radio"/> Calibration <input type="radio"/> SmartWash <input type="radio"/> Results <input type="radio"/> Interpretation |  |                                               |  |  |
| Assay: <b>Amp5Q</b>                                                                                                                                                           |  | Calibration method: <b>Linear</b>             |  |  |
| <input checked="" type="radio"/> Calibrators <input type="radio"/> Volumes <input type="radio"/> Intervals <input type="radio"/> Validity checks                              |  |                                               |  |  |
| Calibrator set: <b>DOAQL</b>                                                                                                                                                  |  | Calibrator level: (†††)      Concentration: † |  |  |
|                                                                                                                                                                               |  | Blank: <b>DOAQL1 (DOA MC Neg Cal)</b> 0       |  |  |
|                                                                                                                                                                               |  | Cal 1: <b>DOAQL2 (DOA MC 1 Cal 1)</b>         |  |  |
| Replicates: <b>2</b> [Range 1 – 3]                                                                                                                                            |  |                                               |  |  |

| Configure assay parameters — Volumes                                                                                                             |            |                |         |       |
|--------------------------------------------------------------------------------------------------------------------------------------------------|------------|----------------|---------|-------|
| <input type="radio"/> Calibrators <input checked="" type="radio"/> Volumes <input type="radio"/> Intervals <input type="radio"/> Validity checks |            |                |         |       |
| Calibrator: <b>DOAQL</b>                                                                                                                         |            |                |         |       |
| Calibrator level                                                                                                                                 | Sample     | Diluted sample | Diluent | Water |
| Blank: <b>DOAQL1</b>                                                                                                                             | <b>8.0</b> | —              | —       | —     |
| Cal 1: <b>DOAQL2</b>                                                                                                                             | <b>8.0</b> | —              | —       | —     |

| Configure assay parameters — Intervals                                                                                                           |  |         |  |  |
|--------------------------------------------------------------------------------------------------------------------------------------------------|--|---------|--|--|
| <input type="radio"/> Calibrators <input type="radio"/> Volumes <input checked="" type="radio"/> Intervals <input type="radio"/> Validity checks |  |         |  |  |
| Calibration intervals:                                                                                                                           |  |         |  |  |
| Full interval: <b>312</b>                                                                                                                        |  | (hours) |  |  |
| Calibration type:                                                                                                                                |  |         |  |  |
| Adjust type: <b>None</b>                                                                                                                         |  |         |  |  |

| Configure assay parameters — Validity checks                                                                                                     |  |  |  |  |
|--------------------------------------------------------------------------------------------------------------------------------------------------|--|--|--|--|
| <input type="radio"/> Calibrators <input type="radio"/> Volumes <input type="radio"/> Intervals <input checked="" type="radio"/> Validity checks |  |  |  |  |
| Blank absorbance range: — —                                                                                                                      |  |  |  |  |
| Span: <b>Blank</b> — <b>Blank</b>                                                                                                                |  |  |  |  |
| Span absorbance range: — —                                                                                                                       |  |  |  |  |
| Expected cal factor: <b>0.00</b>                                                                                                                 |  |  |  |  |
| Expected cal factor tolerance %: <b>0</b>                                                                                                        |  |  |  |  |

| Configure assay parameters — SmartWash                                                                                                                                        |                 |                        |            |            |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------|------------------------|------------|------------|
| <input type="radio"/> General <input type="radio"/> Calibration <input checked="" type="radio"/> SmartWash <input type="radio"/> Results <input type="radio"/> Interpretation |                 |                        |            |            |
| Assay: <b>Amp5Q</b>                                                                                                                                                           |                 |                        |            |            |
| COMPONENT                                                                                                                                                                     | REAGENT / ASSAY | WASH                   | Volume     | Replicates |
| <b>R1</b>                                                                                                                                                                     | <b>ACETM</b>    | <b>Detergent A</b>     | <b>345</b> | <b>1</b>   |
| <b>R2</b>                                                                                                                                                                     | <b>ACETM</b>    | <b>Detergent A</b>     | <b>345</b> | <b>1</b>   |
| <b>Cuvette</b>                                                                                                                                                                | <b>Trig**</b>   | <b>10% Detergent B</b> | <b>345</b> |            |
| **Not required for ARCHITECT Software version 7.00 and above.                                                                                                                 |                 |                        |            |            |

| Configure assay parameters — Results                                                                                                                                          |  |                           |  |  |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|---------------------------|--|--|
| <input type="radio"/> General <input type="radio"/> Calibration <input type="radio"/> SmartWash <input checked="" type="radio"/> Results <input type="radio"/> Interpretation |  |                           |  |  |
| Assay: <b>Amp5Q</b>                                                                                                                                                           |  | Assay number: <b>2879</b> |  |  |
| Dilution default range:                                                                                                                                                       |  | Result units: —           |  |  |
| Low-Linearity:                                                                                                                                                                |  |                           |  |  |
| High-Linearity:                                                                                                                                                               |  |                           |  |  |

| Configure assay parameters — Interpretation                                                                                                                                   |               |                           |  |  |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------|---------------------------|--|--|
| <input type="radio"/> General <input type="radio"/> Calibration <input type="radio"/> SmartWash <input type="radio"/> Results <input checked="" type="radio"/> Interpretation |               |                           |  |  |
| Assay: <b>Amp5Q</b>                                                                                                                                                           |               | Assay number: <b>2879</b> |  |  |
| Name:                                                                                                                                                                         | Range:        | Results review required:  |  |  |
| <b>Negative</b>                                                                                                                                                               | < <b>500</b>  | <input type="checkbox"/>  |  |  |
| <b>Positive</b>                                                                                                                                                               | >= <b>500</b> | <input type="checkbox"/>  |  |  |

| Configure result units            |               |
|-----------------------------------|---------------|
| Assay: <b>Amp5Q</b>               |               |
| Version: †                        |               |
| Result units: —                   |               |
| Decimal places: <b>0</b>          | (Range 0 – 4) |
| Correlation factor: <b>1.0000</b> |               |
| Intercept: <b>0.0000</b>          |               |

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# ARCHITECT cSystems Assay Parameters

## Amphetamine/Methamphetamine 1000 Качественный

| Configure assay parameters — General                 |                                        |                                       |                               |
|------------------------------------------------------|----------------------------------------|---------------------------------------|-------------------------------|
| <input checked="" type="radio"/> General             | <input type="radio"/> Calibration      | <input type="radio"/> SmartWash       | <input type="radio"/> Results |
| Assay: <b>AmpQ</b>                                   | Type: <b>Photometric</b>               | Version: †                            |                               |
| Number: <b>2849</b>                                  | Assay availability: <b>Enabled</b>     |                                       |                               |
| Run controls for onboard reagents by: * <b>Lot</b>   |                                        |                                       |                               |
| <input checked="" type="radio"/> Reaction definition | <input type="radio"/> Reagent / Sample | <input type="radio"/> Validity checks |                               |
| Reaction mode: <b>Rate up</b>                        |                                        |                                       |                               |
| Primary                                              |                                        | Secondary                             | Read times                    |
| Wavelength: <b>340</b> / ††                          |                                        | Main: <b>20 – 24</b>                  |                               |
| Last required read: <b>24</b>                        |                                        | Flex: — —                             |                               |
| Absorbance range: <b>-0.1000 – 3.2000</b>            |                                        | Color correction: — —                 |                               |
| Sample blank type: <b>None</b>                       |                                        |                                       |                               |

| Configure assay parameters — Results |                                   |                                 |                                          |
|--------------------------------------|-----------------------------------|---------------------------------|------------------------------------------|
| <input type="radio"/> General        | <input type="radio"/> Calibration | <input type="radio"/> SmartWash | <input checked="" type="radio"/> Results |
| Assay: <b>AmpQ</b>                   |                                   | Assay number: <b>2849</b>       |                                          |
| Dilution default range:              |                                   | Result units: —                 |                                          |
| Low-Linearity:                       |                                   |                                 |                                          |
| High-Linearity:                      |                                   |                                 |                                          |

| Configure assay parameters — Interpretation |                                   |                                 |                               |
|---------------------------------------------|-----------------------------------|---------------------------------|-------------------------------|
| <input type="radio"/> General               | <input type="radio"/> Calibration | <input type="radio"/> SmartWash | <input type="radio"/> Results |
| Assay: <b>AmpQ</b>                          |                                   |                                 |                               |
| Assay number: <b>2849</b>                   |                                   |                                 |                               |
| Name:                                       | Range:                            | Results review required:        |                               |
| <b>Negative</b>                             | < 1000                            | <input type="checkbox"/>        |                               |
| <b>Positive</b>                             | >= 1000                           | <input type="checkbox"/>        |                               |

| Configure assay parameters — Calibration     |                                              |                                 |                                       |
|----------------------------------------------|----------------------------------------------|---------------------------------|---------------------------------------|
| <input type="radio"/> General                | <input checked="" type="radio"/> Calibration | <input type="radio"/> SmartWash | <input type="radio"/> Results         |
| Assay: <b>AmpQ</b>                           |                                              |                                 |                                       |
| Calibration method: <b>Linear</b>            |                                              |                                 |                                       |
| <input checked="" type="radio"/> Calibrators | <input type="radio"/> Volumes                | <input type="radio"/> Intervals | <input type="radio"/> Validity checks |
| Calibrator set: <b>DOAQ</b>                  | Calibrator level: (†††)                      | Concentration: ‡                |                                       |
|                                              | Blank: <b>DOAQ1 (DOA MC Neg Cal)</b>         | <b>0</b>                        |                                       |
|                                              | Cal 1: <b>DOAQ2 (DOA MC I Cal 2)</b>         |                                 |                                       |
| Replicates: <b>2</b> [Range 1 – 3]           |                                              |                                 |                                       |

| Configure assay parameters — Results |                                   |                                 |                                          |
|--------------------------------------|-----------------------------------|---------------------------------|------------------------------------------|
| <input type="radio"/> General        | <input type="radio"/> Calibration | <input type="radio"/> SmartWash | <input checked="" type="radio"/> Results |
| Assay: <b>AmpQ</b>                   |                                   |                                 |                                          |
| Assay number: <b>2849</b>            |                                   |                                 |                                          |
| Name:                                | Range:                            | Results review required:        |                                          |
| <b>Negative</b>                      | < 1000                            | <input type="checkbox"/>        |                                          |
| <b>Positive</b>                      | >= 1000                           | <input type="checkbox"/>        |                                          |

| Configure assay parameters — Interpretation |                                   |                                 |                               |
|---------------------------------------------|-----------------------------------|---------------------------------|-------------------------------|
| <input type="radio"/> General               | <input type="radio"/> Calibration | <input type="radio"/> SmartWash | <input type="radio"/> Results |
| Assay: <b>AmpQ</b>                          |                                   |                                 |                               |
| Assay number: <b>2849</b>                   |                                   |                                 |                               |
| Name:                                       | Range:                            | Results review required:        |                               |
| <b>Negative</b>                             | < 1000                            | <input type="checkbox"/>        |                               |
| <b>Positive</b>                             | >= 1000                           | <input type="checkbox"/>        |                               |

| Configure assay parameters — Calibration     |                                              |                                 |                                       |
|----------------------------------------------|----------------------------------------------|---------------------------------|---------------------------------------|
| <input type="radio"/> General                | <input checked="" type="radio"/> Calibration | <input type="radio"/> SmartWash | <input type="radio"/> Results         |
| Assay: <b>AmpQ</b>                           |                                              |                                 |                                       |
| Calibration method: <b>Linear</b>            |                                              |                                 |                                       |
| <input checked="" type="radio"/> Calibrators | <input type="radio"/> Volumes                | <input type="radio"/> Intervals | <input type="radio"/> Validity checks |
| Calibrator set: <b>DOAQ</b>                  | Calibrator level: (†††)                      | Concentration: ‡                |                                       |
|                                              | Blank: <b>DOAQ1 (DOA MC Neg Cal)</b>         | <b>0</b>                        |                                       |
|                                              | Cal 1: <b>DOAQ2 (DOA MC I Cal 2)</b>         |                                 |                                       |
| Replicates: <b>2</b> [Range 1 – 3]           |                                              |                                 |                                       |

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| Configure assay parameters — General                 |                                        |                                       |                               |
|------------------------------------------------------|----------------------------------------|---------------------------------------|-------------------------------|
| <input checked="" type="radio"/> General             | <input type="radio"/> Calibration      | <input type="radio"/> SmartWash       | <input type="radio"/> Results |
| Assay: <b>AmpSQ</b>                                  | Type: <b>Photometric</b>               | Version: <b>†</b>                     |                               |
| Number: <b>2848</b>                                  | Assay availability: <b>Enabled</b>     |                                       |                               |
| Run controls for onboard reagents by: <b>Lot</b>     |                                        |                                       |                               |
| <input checked="" type="radio"/> Reaction definition | <input type="radio"/> Reagent / Sample | <input type="radio"/> Validity checks |                               |
| Reaction mode: <b>Rate up</b>                        | Primary                                | Read times                            |                               |
| Wavelength: <b>340</b>                               | Secondary                              | Main: <b>20 – 24</b>                  |                               |
| Last required read: <b>24</b>                        |                                        | Flex: <b>— —</b>                      |                               |
| Absorbance range: <b>-0.1000 – 3.2000</b>            |                                        | Color correction: <b>— —</b>          |                               |
| Sample blank type: <b>None</b>                       |                                        |                                       |                               |

| <input type="radio"/> Reaction definition | <input checked="" type="radio"/> Reagent / Sample | <input type="radio"/> Validity checks |    |
|-------------------------------------------|---------------------------------------------------|---------------------------------------|----|
| Reagent: <b>AMPM9</b>                     | Reagent volume: <b>100</b>                        | R1                                    | R2 |
| Diluent: <b>None</b>                      | Water volume: <b>—</b>                            |                                       |    |
| Dilution name: <b>STANDARD</b>            | Dispense mode: <b>Type 0</b>                      | Type 0                                |    |
| Sample                                    | Diluted sample                                    | Default dilution                      |    |
| 8.0                                       | —                                                 | 1:1.00                                |    |

| <input type="radio"/> Reaction definition | <input type="radio"/> Reagent / Sample | <input checked="" type="radio"/> Validity checks |
|-------------------------------------------|----------------------------------------|--------------------------------------------------|
| Reaction check: <b>None</b>               |                                        |                                                  |
| Rate linearity %: <b>—</b>                |                                        |                                                  |

| Configure assay parameters — Calibration     |                                              |                                 |                                       |
|----------------------------------------------|----------------------------------------------|---------------------------------|---------------------------------------|
| <input type="radio"/> General                | <input checked="" type="radio"/> Calibration | <input type="radio"/> SmartWash | <input type="radio"/> Results         |
| Assay: <b>AmpSQ</b>                          | Calibration method: <b>Linear</b>            |                                 |                                       |
| <input checked="" type="radio"/> Calibrators | <input type="radio"/> Volumes                | <input type="radio"/> Intervals | <input type="radio"/> Validity checks |
| Calibrator set: <b>DOASQ</b>                 | Calibrator level: <b>(†††)</b>               | Concentration: <b>†</b>         |                                       |
| Replicates: <b>2</b>                         | Blank: <b>DOASQ1 (DOA MC Neg Cal)</b>        | <b>0.0</b>                      |                                       |
|                                              | Cal 1: <b>DOASQ2 (DOA MC I Cal 1)</b>        |                                 |                                       |
|                                              | Cal 2: <b>DOASQ3 (DOA MC I Cal 2)</b>        |                                 |                                       |
|                                              | Cal 3: <b>DOASQ4 (DOA MC I Cal 3)</b>        |                                 |                                       |
|                                              | Cal 4: <b>DOASQ5 (DOA MC I Cal 4)</b>        |                                 |                                       |

| <input type="radio"/> Calibrators | <input checked="" type="radio"/> Volumes | <input type="radio"/> Intervals | <input type="radio"/> Validity checks |
|-----------------------------------|------------------------------------------|---------------------------------|---------------------------------------|
| Calibrator: <b>DOASQ</b>          | Calibrator level                         | Sample                          | Diluted sample                        |
|                                   | Blank: <b>DOASQ1</b>                     | <b>8.0</b>                      | <b>—</b>                              |
|                                   | Cal 1: <b>DOASQ2</b>                     | <b>8.0</b>                      | <b>—</b>                              |
|                                   | Cal 2: <b>DOASQ3</b>                     | <b>8.0</b>                      | <b>—</b>                              |
|                                   | Cal 3: <b>DOASQ4</b>                     | <b>8.0</b>                      | <b>—</b>                              |
|                                   | Cal 4: <b>DOASQ5</b>                     | <b>8.0</b>                      | <b>—</b>                              |

| <input type="radio"/> Calibrators | <input type="radio"/> Volumes | <input checked="" type="radio"/> Intervals | <input type="radio"/> Validity checks |
|-----------------------------------|-------------------------------|--------------------------------------------|---------------------------------------|
| Calibration intervals:            |                               |                                            |                                       |
| Full interval: <b>312</b> (hours) |                               |                                            |                                       |
| Calibration type:                 |                               |                                            |                                       |
| Adjust type: <b>None</b>          |                               |                                            |                                       |

| <input type="radio"/> Calibrators         | <input type="radio"/> Volumes | <input type="radio"/> Intervals | <input checked="" type="radio"/> Validity checks |
|-------------------------------------------|-------------------------------|---------------------------------|--------------------------------------------------|
| Blank absorbance range: <b>—</b>          |                               |                                 |                                                  |
| Span: <b>Blank</b> — <b>Blank</b>         |                               |                                 |                                                  |
| Span absorbance range: <b>—</b>           |                               |                                 |                                                  |
| Expected cal factor: <b>0.00</b>          |                               |                                 |                                                  |
| Expected cal factor tolerance %: <b>0</b> |                               |                                 |                                                  |

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\*\*Not required for ARCHITECT Software version 7.00 and above.

| Configure assay parameters — Results |                                   |                                 |                                          |
|--------------------------------------|-----------------------------------|---------------------------------|------------------------------------------|
| <input type="radio"/> General        | <input type="radio"/> Calibration | <input type="radio"/> SmartWash | <input checked="" type="radio"/> Results |
| Assay: <b>AmpSQ</b>                  |                                   |                                 |                                          |
| Default dilution range:              |                                   |                                 |                                          |
| Low-Linearity: <b>††</b> 100.0       |                                   |                                 |                                          |
| High-Linearity: <b>2000.0</b>        |                                   |                                 |                                          |
| Assay number: <b>2848</b>            |                                   |                                 |                                          |
| Result units: <b>ng/mL</b>           |                                   |                                 |                                          |

| Configure result units            |               |
|-----------------------------------|---------------|
| Assay: <b>AmpSQ</b>               |               |
| Version: <b>†</b>                 |               |
| Result units: <b>ng/mL</b>        |               |
| Decimal places: <b>1</b>          | (Range 0 – 4) |
| Correlation factor: <b>1.0000</b> |               |
| Intercept: <b>0.0000</b>          |               |

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1. Hawks RL, Chiang CN, editors. *Urine Testing for Drugs of Abuse*. NIDA Research Monograph 73. Rockville, MD: National Institute on Drug Abuse (NIDA), Department of Health and Human Services; 1986.
2. *Mandatory Guidelines for Federal Workplace Drug Testing Programs*. National Institute on Drug Abuse. Federal Register Vol. 69, No. 71, 2004:19650–61.
3. Baselt RC. *Disposition of Toxic Drugs and Chemicals in Man*, 5th ed. Foster City, CA: Chemical Toxicology Institute; 2000:44–6, 60–1, 75–6, 108–9, 110–1, 507–9, 538–9, 541–2, 668–70, 689–91, 774–6.
4. Rubenstein KE, Schneider RS, Ullman EF. Homogeneous enzyme immunoassay: a new immunochemical technique. *Biochem Biophys Res Commun* 1972;47:846–51.
5. US Department of Labor, Occupational Safety and Health Administration. 29 CFR Part 1910.1030, *Occupational Exposure to Bloodborne Pathogens*.
6. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; January 2007.
7. World Health Organization. *Laboratory Biosafety Manual*. Geneva: World Health Organization. 2004
8. Sewell DL, Bove KE, Callihan DR, et al. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Third Edition (M29-A3)*. Wayne, PA: Clinical and Laboratory Standards Institute, 2005.
9. Shoemaker MJ, Earley RJ, Gorsky JE, et al. *Urine Drug Testing in the Clinical Laboratory; Approved Guideline—Vol. 19, No. 6*. Wayne, PA: Clinical and Laboratory Standards Institute, 1999(17).

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| Configure assay parameters — General                                                                                                                                                                                                                                                             |                                    |                   |  |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------|-------------------|--|
| <input checked="" type="radio"/> General <input type="radio"/> Calibration <input type="radio"/> SmartWash <input type="radio"/> Results <input type="radio"/> Interpretation                                                                                                                    |                                    |                   |  |
| Assay: <b>BarbQ</b>                                                                                                                                                                                                                                                                              | Type: <b>Photometric</b>           | Version: <b>†</b> |  |
| Number: <b>2851</b>                                                                                                                                                                                                                                                                              | Assay availability: <b>Enabled</b> |                   |  |
| <input checked="" type="radio"/> Reaction definition <input type="radio"/> Reagent / Sample <input type="radio"/> Validity checks                                                                                                                                                                |                                    |                   |  |
| Reaction mode: <b>Rate up</b><br>Primary      Secondary      Read times<br>Wavelength: <b>340</b> / <b>‡‡</b> Main: <b>20 – 24</b><br>Last required read: <b>24</b> Flex: <b>— —</b><br>Absorbance range: <b>-0.1000 – 3.2000</b> Color correction: <b>— —</b><br>Sample blank type: <b>None</b> |                                    |                   |  |

| Configure assay parameters — General                                                                                                                                                                                                                                                                                                                                                                              |  |  |  |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|--|--|
| <input type="radio"/> Reaction definition <input checked="" type="radio"/> Reagent / Sample <input type="radio"/> Validity checks                                                                                                                                                                                                                                                                                 |  |  |  |
| Reagent: <b>BARBM</b> Reagent volume: <b>100</b> R1      R2<br>Diluent: <b>None</b> Water volume: <b>—</b><br>Diluent dispense mode: <b>Type 0</b> Dispense mode: <b>Type 0</b> Type 0<br>Dilution name      Sample      Diluted sample      Diluent      Water      Dilution factor      Default dilution<br><b>STANDARD</b> : <b>12.0</b> —      —      —      — <b>1:1.00</b> <input checked="" type="radio"/> |  |  |  |
| <input type="radio"/> Reaction definition <input type="radio"/> Reagent / Sample <input checked="" type="radio"/> Validity checks                                                                                                                                                                                                                                                                                 |  |  |  |
| Reaction check: <b>None</b><br><br>Rate linearity %: <b>—</b>                                                                                                                                                                                                                                                                                                                                                     |  |  |  |

| Configure assay parameters — Calibration                                                                                                                                                                 |  |  |  |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|--|--|
| <input type="radio"/> General <input checked="" type="radio"/> Calibration <input type="radio"/> SmartWash <input type="radio"/> Results <input type="radio"/> Interpretation                            |  |  |  |
| Assay: <b>BarbQ</b> Calibration method: <b>Linear</b>                                                                                                                                                    |  |  |  |
| <input checked="" type="radio"/> Calibrators <input type="radio"/> Volumes <input type="radio"/> Intervals <input type="radio"/> Validity checks                                                         |  |  |  |
| Calibrator set: <b>DOAQ</b> Calibrator level: <b>(†††)</b> Concentration: <b>‡</b><br>Blank: <b>DOAQ1 (DOA MC Neg Cal)</b><br>Cal 1: <b>DOAQ2 (DOA MC I Cal 2)</b><br>Replicates: <b>2</b> [Range 1 – 3] |  |  |  |
| <input type="radio"/> Calibrators <input checked="" type="radio"/> Volumes <input type="radio"/> Intervals <input type="radio"/> Validity checks                                                         |  |  |  |
| Calibrator: <b>DOAQ</b><br>Blank: <b>DOAQ1</b> Sample      Diluted sample      Diluent      Water<br>Cal 1: <b>DOAQ2</b> <b>12.0</b> —      —      —<br>Cal 1: <b>DOAQ2</b> <b>12.0</b> —      —      —  |  |  |  |
| <input type="radio"/> Calibrators <input type="radio"/> Volumes <input checked="" type="radio"/> Intervals <input type="radio"/> Validity checks                                                         |  |  |  |
| Calibration intervals:<br>Full interval: <b>312</b> (hours)<br>Calibration type:<br>Adjust type: <b>None</b>                                                                                             |  |  |  |
| <input type="radio"/> Calibrators <input type="radio"/> Volumes <input type="radio"/> Intervals <input checked="" type="radio"/> Validity checks                                                         |  |  |  |
| Blank absorbance range: <b>—</b><br>Span: <b>Blank</b> —<br>Span absorbance range: <b>—</b><br>Expected cal factor: <b>0.00</b><br>Expected cal factor tolerance %: <b>0</b>                             |  |  |  |

| Configure assay parameters — SmartWash                                                                                                                                        |                 |                           |                        |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------|---------------------------|------------------------|
| <input type="radio"/> General <input type="radio"/> Calibration <input checked="" type="radio"/> SmartWash <input type="radio"/> Results <input type="radio"/> Interpretation |                 |                           |                        |
| Assay: <b>BarbQ</b>                                                                                                                                                           |                 |                           |                        |
| COMPONENT                                                                                                                                                                     | REAGENT / ASSAY | WASH                      | Volume      Replicates |
| Cuvette                                                                                                                                                                       | Trig            | <b>10% Detergent B***</b> | <b>345</b>             |
| ***Select "Detergent B" for software prior to Version 2.2.                                                                                                                    |                 |                           |                        |

| Configure assay parameters — Results                                                                                                                                          |  |  |  |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|--|--|
| <input type="radio"/> General <input type="radio"/> Calibration <input type="radio"/> SmartWash <input checked="" type="radio"/> Results <input type="radio"/> Interpretation |  |  |  |
| Assay: <b>BarbQ</b> Result units: <b>—</b>                                                                                                                                    |  |  |  |
| Assay defaults:<br>Low-Linearity:<br>High-Linearity:                                                                                                                          |  |  |  |

| Configure assay parameters — Interpretation                                                                                                                                   |                  |                          |  |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------|--------------------------|--|
| <input type="radio"/> General <input type="radio"/> Calibration <input type="radio"/> SmartWash <input type="radio"/> Results <input checked="" type="radio"/> Interpretation |                  |                          |  |
| Assay: <b>BarbQ</b>                                                                                                                                                           |                  |                          |  |
| Name:                                                                                                                                                                         | Range:           | Results review required: |  |
| <b>Negative</b>                                                                                                                                                               | <b>&lt; 200</b>  | <input type="checkbox"/> |  |
| <b>Positive</b>                                                                                                                                                               | <b>&gt;= 200</b> | <input type="checkbox"/> |  |

| Configure result units                 |  |
|----------------------------------------|--|
| Assay: <b>BarbQ</b>                    |  |
| Version: <b>†</b>                      |  |
| Result units: <b>—</b>                 |  |
| Decimal places: <b>0</b> [Range 0 – 4] |  |
| Correlation factor: <b>1.0000</b>      |  |
| Intercept: <b>0.0000</b>               |  |

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## Configure assay parameters — General

|                                                                                                                                                                                                                                                                                                  |                                    |                          |                   |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------|--------------------------|-------------------|
| <input checked="" type="radio"/> General <input type="radio"/> Calibration <input type="radio"/> SmartWash <input type="radio"/> Results <input type="radio"/> Interpretation                                                                                                                    |                                    |                          |                   |
| Assay: <b>BarbSQ</b>                                                                                                                                                                                                                                                                             |                                    | Type: <b>Photometric</b> | Version: <b>†</b> |
| Number: <b>2850</b>                                                                                                                                                                                                                                                                              | Assay availability: <b>Enabled</b> |                          |                   |
| <input checked="" type="radio"/> Reaction definition <input type="radio"/> Reagent / Sample <input type="radio"/> Validity checks                                                                                                                                                                |                                    |                          |                   |
| Reaction mode: <b>Rate up</b><br>Primary      Secondary      Read times<br>Wavelength: <b>340</b> / <b>‡‡</b> Main: <b>20 – 24</b><br>Last required read: <b>24</b> Flex: <b>— —</b><br>Absorbance range: <b>-0.1000 – 3.2000</b> Color correction: <b>— —</b><br>Sample blank type: <b>None</b> |                                    |                          |                   |

|                                           |               |                                                   |          |                                       |                 |
|-------------------------------------------|---------------|---------------------------------------------------|----------|---------------------------------------|-----------------|
| <input type="radio"/> Reaction definition |               | <input checked="" type="radio"/> Reagent / Sample |          | <input type="radio"/> Validity checks |                 |
| Reagent: <b>BARBM</b>                     |               | Reagent volume: <b>100</b>                        |          | R1 R2                                 |                 |
| Diluent: <b>None</b>                      |               | Water volume: <b>—</b>                            |          | —                                     |                 |
| Diluent dispense mode: <b>Type 0</b>      |               | Dispense mode: <b>Type 0</b>                      |          | <b>Type 0</b> <b>Type 0</b>           |                 |
| Dilution name                             | Sample        | Diluted sample                                    | Diluent  | Water                                 | Dilution factor |
| <b>STANDARD</b>                           | : <b>12.0</b> | <b>—</b>                                          | <b>—</b> | <b>—</b>                              | <b>1:1.00</b>   |
|                                           |               |                                                   |          | Default dilution                      |                 |
|                                           |               |                                                   |          | <input checked="" type="radio"/>      |                 |

|                                           |  |                                        |  |                                                  |  |
|-------------------------------------------|--|----------------------------------------|--|--------------------------------------------------|--|
| <input type="radio"/> Reaction definition |  | <input type="radio"/> Reagent / Sample |  | <input checked="" type="radio"/> Validity checks |  |
| Reaction check: <b>None</b>               |  |                                        |  |                                                  |  |
| Rate linearity %: <b>—</b>                |  |                                        |  |                                                  |  |

## Configure assay parameters — Calibration

|                                                                                                                                                                               |  |                                                                                                                                                                                                           |  |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|
| <input type="radio"/> General <input checked="" type="radio"/> Calibration <input type="radio"/> SmartWash <input type="radio"/> Results <input type="radio"/> Interpretation |  |                                                                                                                                                                                                           |  |
| Assay: <b>BarbSQ</b>                                                                                                                                                          |  | Calibration method: <b>Linear</b>                                                                                                                                                                         |  |
| <input checked="" type="radio"/> Calibrators <input type="radio"/> Volumes <input type="radio"/> Intervals <input type="radio"/> Validity checks                              |  |                                                                                                                                                                                                           |  |
| Calibrator set: <b>DOASQ</b>                                                                                                                                                  |  | Calibrator level: <b>(†††)</b> Concentration: <b>‡</b>                                                                                                                                                    |  |
| Replicates: <b>2</b> [Range 1 – 3]                                                                                                                                            |  | Blank: <b>DOASQ1 (DOA MC Neg Cal)</b><br>Cal 1: <b>DOASQ2 (DOA MC I Cal 1)</b><br>Cal 2: <b>DOASQ3 (DOA MC I Cal 2)</b><br>Cal 3: <b>DOASQ4 (DOA MC I Cal 3)</b><br>Cal 4: <b>DOASQ5 (DOA MC I Cal 4)</b> |  |

|                                   |                      |                                          |                |                                 |          |                                       |  |
|-----------------------------------|----------------------|------------------------------------------|----------------|---------------------------------|----------|---------------------------------------|--|
| <input type="radio"/> Calibrators |                      | <input checked="" type="radio"/> Volumes |                | <input type="radio"/> Intervals |          | <input type="radio"/> Validity checks |  |
| Calibrator: <b>DOASQ</b>          |                      |                                          |                |                                 |          |                                       |  |
|                                   | Calibrator level     | Sample                                   | Diluted sample | Diluent                         | Water    |                                       |  |
|                                   | Blank: <b>DOASQ1</b> | <b>12.0</b>                              | <b>—</b>       | <b>—</b>                        | <b>—</b> |                                       |  |
|                                   | Cal 1: <b>DOASQ2</b> | <b>12.0</b>                              | <b>—</b>       | <b>—</b>                        | <b>—</b> |                                       |  |
|                                   | Cal 2: <b>DOASQ3</b> | <b>12.0</b>                              | <b>—</b>       | <b>—</b>                        | <b>—</b> |                                       |  |
|                                   | Cal 3: <b>DOASQ4</b> | <b>12.0</b>                              | <b>—</b>       | <b>—</b>                        | <b>—</b> |                                       |  |
|                                   | Cal 4: <b>DOASQ5</b> | <b>12.0</b>                              | <b>—</b>       | <b>—</b>                        | <b>—</b> |                                       |  |

|                                   |  |                               |  |                                            |  |                                       |  |
|-----------------------------------|--|-------------------------------|--|--------------------------------------------|--|---------------------------------------|--|
| <input type="radio"/> Calibrators |  | <input type="radio"/> Volumes |  | <input checked="" type="radio"/> Intervals |  | <input type="radio"/> Validity checks |  |
| Calibration intervals:            |  |                               |  |                                            |  |                                       |  |
| Full interval: <b>312</b>         |  | (hours)                       |  |                                            |  |                                       |  |
| Calibration type:                 |  |                               |  |                                            |  |                                       |  |
| Adjust type: <b>None</b>          |  |                               |  |                                            |  |                                       |  |

|                                           |  |                               |  |                                 |  |                                                  |  |
|-------------------------------------------|--|-------------------------------|--|---------------------------------|--|--------------------------------------------------|--|
| <input type="radio"/> Calibrators         |  | <input type="radio"/> Volumes |  | <input type="radio"/> Intervals |  | <input checked="" type="radio"/> Validity checks |  |
| Blank absorbance range: <b>—</b> <b>—</b> |  |                               |  |                                 |  |                                                  |  |
| Span: <b>Blank</b> <b>—</b> <b>—</b>      |  |                               |  |                                 |  |                                                  |  |
| Span absorbance range: <b>—</b> <b>—</b>  |  |                               |  |                                 |  |                                                  |  |
| Expected cal factor: <b>0.00</b>          |  |                               |  |                                 |  |                                                  |  |
| Expected cal factor tolerance %: <b>0</b> |  |                               |  |                                 |  |                                                  |  |

## Configure assay parameters — SmartWash

|                                                                                                                                                                               |                 |                    |        |            |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------|--------------------|--------|------------|
| <input type="radio"/> General <input type="radio"/> Calibration <input checked="" type="radio"/> SmartWash <input type="radio"/> Results <input type="radio"/> Interpretation |                 |                    |        |            |
| Assay: <b>BarbSQ</b>                                                                                                                                                          |                 |                    |        |            |
| COMPONENT                                                                                                                                                                     | REAGENT / ASSAY | WASH               | Volume | Replicates |
| Cuvette                                                                                                                                                                       | Trig            | 10% Detergent B*** | 345    |            |
| ***Select "Detergent B" for software prior to Version 2.2.                                                                                                                    |                 |                    |        |            |

## Configure assay parameters — Results

|                               |  |                                   |  |                                 |  |                                          |  |                                      |  |
|-------------------------------|--|-----------------------------------|--|---------------------------------|--|------------------------------------------|--|--------------------------------------|--|
| <input type="radio"/> General |  | <input type="radio"/> Calibration |  | <input type="radio"/> SmartWash |  | <input checked="" type="radio"/> Results |  | <input type="radio"/> Interpretation |  |
| Assay: <b>BarbSQ</b>          |  | Result units: <b>ng/mL</b>        |  |                                 |  |                                          |  |                                      |  |
| Assay defaults:               |  |                                   |  |                                 |  |                                          |  |                                      |  |
| Low-Linearity:††              |  | <b>25.0</b>                       |  |                                 |  |                                          |  |                                      |  |
| High-Linearity:               |  | <b>1000.0</b>                     |  |                                 |  |                                          |  |                                      |  |

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1. Hawks RL, Chiang CN, editors. *Urine Testing for Drugs of Abuse*. NIDA Research Monograph 73. Rockville, MD: National Institute on Drug Abuse (NIDA), Department of Health and Human Services; 1986.
2. *Mandatory Guidelines for Federal Workplace Drug Testing Programs*. National Institute on Drug Abuse. Federal Register Vol. 69, No. 71, 2004:19650-61.
3. Baselt RC. *Disposition of Toxic Drugs and Chemicals in Man*, 5th ed. Foster City, CA: Chemical Toxicology Institute; 2000:25-7, 94-5, 153-6, 191-2, 194-6, 198-200, 245-8, 322-3, 360-2, 377-80, 398-9, 482-4, 499-500, 580-3, 617-9, 641-2, 719-20, 756-7, 802-3, 851-4.
4. Rubenstein KE, Schneider RS, Ullman EF. Homogeneous enzyme immunoassay: a new immunochemical technique. *Biochem Biophys Res Commun*, 1972;47:846-51.
5. US Department of Labor, Occupational Safety and Health Administration. 29 CFR Part 1910.1030. *Bloodborne Pathogens*.
6. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*, 5th ed. Washington, DC: US Government Printing Office, January 2007.
7. World Health Organization. *Laboratory Biosafety Manual*, 3rd ed. Geneva: World Health Organization, 2004.
8. Sewell DL, Bove KE, Callihan DR, et al. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Third Edition (M29-A3)*. Wayne, PA: Clinical and Laboratory Standards Institute, 2005.
9. Shoemaker MJ, Earley RJ, Gorsky JE, et al. *Urine Drug Testing in the Clinical Laboratory; Approved Guideline—Vol. 19, No. 6*. Wayne, PA: Clinical and Laboratory Standards Institute, 1999(17).
10. Borrey D, Meyer E, Duchateau L, et. al. Enzymatic hydrolysis improves the sensitivity of Emit screening for urinary benzodiazepines. *Clin Chem* 2002;48(11):2047-9.

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### Configure assay parameters — General

|                                                      |                                        |                                       |                               |                                      |
|------------------------------------------------------|----------------------------------------|---------------------------------------|-------------------------------|--------------------------------------|
| <input checked="" type="radio"/> General             | <input type="radio"/> Calibration      | <input type="radio"/> SmartWash       | <input type="radio"/> Results | <input type="radio"/> Interpretation |
| Assay: BenzQ                                         |                                        | Type: Photometric                     |                               | Version: †                           |
| Number: 2853                                         |                                        | Assay availability: Enabled           |                               |                                      |
| <input checked="" type="radio"/> Reaction definition | <input type="radio"/> Reagent / Sample | <input type="radio"/> Validity checks |                               |                                      |
| Reaction mode: Rate up                               |                                        |                                       |                               |                                      |
| Primary                                              |                                        | Secondary                             |                               | Read times                           |
| Wavelength: 340 / ‡‡                                 |                                        | Main: 20 – 24                         |                               |                                      |
| Last required read: 24                               |                                        | Flex: — —                             |                               |                                      |
| Absorbance range: -0.1000 – 3.2000                   |                                        | Color correction: — —                 |                               |                                      |
| Sample blank type: None                              |                                        |                                       |                               |                                      |

|                                                                                    |                                                   |                                       |  |  |
|------------------------------------------------------------------------------------|---------------------------------------------------|---------------------------------------|--|--|
| <input type="radio"/> Reaction definition                                          | <input checked="" type="radio"/> Reagent / Sample | <input type="radio"/> Validity checks |  |  |
| Reagent: BENZM                                                                     |                                                   |                                       |  |  |
| Diluent: None                                                                      |                                                   |                                       |  |  |
| Diluent dispense mode: Type 0                                                      |                                                   |                                       |  |  |
| Dispense mode: Type 0 Type 0                                                       |                                                   |                                       |  |  |
| Dilution name Sample Diluted sample Diluent Water Dilution factor Default dilution |                                                   |                                       |  |  |
| STANDARD : 12.0 — — — 1:1.00 ●                                                     |                                                   |                                       |  |  |

|                                           |                                        |                                                  |
|-------------------------------------------|----------------------------------------|--------------------------------------------------|
| <input type="radio"/> Reaction definition | <input type="radio"/> Reagent / Sample | <input checked="" type="radio"/> Validity checks |
| Reaction check: None                      |                                        |                                                  |
| Rate linearity %: —                       |                                        |                                                  |

### Configure assay parameters — Calibration

|                                              |                                              |                                 |                                       |                                      |
|----------------------------------------------|----------------------------------------------|---------------------------------|---------------------------------------|--------------------------------------|
| <input type="radio"/> General                | <input checked="" type="radio"/> Calibration | <input type="radio"/> SmartWash | <input type="radio"/> Results         | <input type="radio"/> Interpretation |
| Assay: BenzQ                                 |                                              | Calibration method: Linear      |                                       |                                      |
| <input checked="" type="radio"/> Calibrators | <input type="radio"/> Volumes                | <input type="radio"/> Intervals | <input type="radio"/> Validity checks |                                      |
| Calibrator set: DOAQ                         |                                              | Calibrator level: (†††)         |                                       | Concentration: †                     |
|                                              |                                              | Blank: DOAQ1 (DOA MC Neg Cal)   |                                       | 0                                    |
|                                              |                                              | Cal 1: DOAQ2 (DOA MC I Cal 2)   |                                       |                                      |
| Replicates: 2 [Range 1 – 3]                  |                                              |                                 |                                       |                                      |

|                                   |                                          |                                 |                                       |               |
|-----------------------------------|------------------------------------------|---------------------------------|---------------------------------------|---------------|
| <input type="radio"/> Calibrators | <input checked="" type="radio"/> Volumes | <input type="radio"/> Intervals | <input type="radio"/> Validity checks |               |
| Calibrator: DOAQ                  |                                          |                                 |                                       |               |
| Calibrator level                  |                                          | Sample                          | Diluted sample                        | Diluent Water |
| Blank: DOAQ1                      |                                          | 12.0                            | —                                     | —             |
| Cal 1: DOAQ2                      |                                          | 12.0                            | —                                     | —             |

|                                   |                               |                                            |                                       |
|-----------------------------------|-------------------------------|--------------------------------------------|---------------------------------------|
| <input type="radio"/> Calibrators | <input type="radio"/> Volumes | <input checked="" type="radio"/> Intervals | <input type="radio"/> Validity checks |
| Calibration intervals:            |                               |                                            |                                       |
| Full interval: 168                |                               | (hours)                                    |                                       |
| Calibration type:                 |                               |                                            |                                       |
| Adjust type: None                 |                               |                                            |                                       |

|                                    |                               |                                 |                                                  |
|------------------------------------|-------------------------------|---------------------------------|--------------------------------------------------|
| <input type="radio"/> Calibrators  | <input type="radio"/> Volumes | <input type="radio"/> Intervals | <input checked="" type="radio"/> Validity checks |
| Blank absorbance range: — —        |                               |                                 |                                                  |
| Span: Blank — Blank                |                               |                                 |                                                  |
| Span absorbance range: — —         |                               |                                 |                                                  |
| Expected cal factor: 0.00          |                               |                                 |                                                  |
| Expected cal factor tolerance %: 0 |                               |                                 |                                                  |

### Configure assay parameters — SmartWash

|                               |                                   |                                            |                               |                                      |
|-------------------------------|-----------------------------------|--------------------------------------------|-------------------------------|--------------------------------------|
| <input type="radio"/> General | <input type="radio"/> Calibration | <input checked="" type="radio"/> SmartWash | <input type="radio"/> Results | <input type="radio"/> Interpretation |
| Assay: BenzQ                  |                                   |                                            |                               |                                      |
| COMPONENT                     | REAGENT / ASSAY                   | WASH                                       | Volume                        | Replicates                           |
| R1                            | BILT0                             | Detergent A                                | 345                           | 1                                    |
| R1                            | ACETM                             | Detergent A                                | 345                           | 1                                    |
| R1                            | ALT00                             | Detergent A                                | 345                           | 1                                    |
| R1                            | AST00                             | Water                                      | 345                           | 1                                    |
| R1                            | LIP00                             | Water                                      | 345                           | 1                                    |
| R2                            | ACETM                             | Detergent A                                | 345                           | 1                                    |
| R2                            | ALT00                             | Detergent A                                | 345                           | 1                                    |
| R2                            | AST00                             | Water                                      | 345                           | 1                                    |
| R2                            | LIP00                             | Water                                      | 345                           | 1                                    |
| Cuvette                       | Trig                              | 10% Detergent B                            | 345                           |                                      |

### Configure assay parameters — Results

|                               |                                   |                                 |                                          |                                      |
|-------------------------------|-----------------------------------|---------------------------------|------------------------------------------|--------------------------------------|
| <input type="radio"/> General | <input type="radio"/> Calibration | <input type="radio"/> SmartWash | <input checked="" type="radio"/> Results | <input type="radio"/> Interpretation |
| Assay: BenzQ                  |                                   | Assay number: 2853              |                                          |                                      |
| Dilution default range:       |                                   | Result units: —                 |                                          |                                      |
| Low-Linearity:                |                                   |                                 |                                          |                                      |
| High-Linearity:               |                                   |                                 |                                          |                                      |

### Configure assay parameters — Interpretation

|                               |                                   |                                 |                               |                                                 |
|-------------------------------|-----------------------------------|---------------------------------|-------------------------------|-------------------------------------------------|
| <input type="radio"/> General | <input type="radio"/> Calibration | <input type="radio"/> SmartWash | <input type="radio"/> Results | <input checked="" type="radio"/> Interpretation |
| Assay: BenzQ                  |                                   | Assay number: 2853              |                               |                                                 |
| Name:                         |                                   | Range:                          | Results review required:      |                                                 |
| Negative                      |                                   | < 200                           | <input type="checkbox"/>      |                                                 |
| Positive                      |                                   | >= 200                          | <input type="checkbox"/>      |                                                 |

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| Configure assay parameters — General                 |                                        |                                       |                               |
|------------------------------------------------------|----------------------------------------|---------------------------------------|-------------------------------|
| <input checked="" type="radio"/> General             | <input type="radio"/> Calibration      | <input type="radio"/> SmartWash       | <input type="radio"/> Results |
| Assay: <b>BenzSQ</b>                                 |                                        | Type: <b>Photometric</b>              | Version: <b>†</b>             |
| Number: <b>2852</b>                                  | Assay availability: <b>Enabled</b>     |                                       |                               |
| <input checked="" type="radio"/> Reaction definition | <input type="radio"/> Reagent / Sample | <input type="radio"/> Validity checks |                               |
| Reaction mode: <b>Rate up</b>                        |                                        |                                       |                               |
| Primary                                              |                                        | Secondary                             | Read times                    |
| Wavelength: <b>340</b> / <b>††</b>                   |                                        | Main: <b>20 – 24</b>                  | Flex: <b>— —</b>              |
| Last required read: <b>24</b>                        |                                        | Flex: <b>— —</b>                      |                               |
| Absorbance range: <b>-0.1000 – 3.2000</b>            |                                        | Color correction: <b>— —</b>          |                               |
| Sample blank type: <b>None</b>                       |                                        |                                       |                               |

| Configure assay parameters — General      |                                                   |                                                |                 |
|-------------------------------------------|---------------------------------------------------|------------------------------------------------|-----------------|
| <input type="radio"/> Reaction definition | <input checked="" type="radio"/> Reagent / Sample | <input type="radio"/> Validity checks          |                 |
| Reagent: <b>BENZM</b>                     |                                                   | Reagent volume: <b>100</b> <b>R1</b> <b>R2</b> |                 |
| Diluent: <b>None</b>                      |                                                   | Water volume: <b>—</b>                         |                 |
| Diluent dispense mode: <b>Type 0</b>      |                                                   | Dispense mode: <b>Type 0</b> <b>Type 0</b>     |                 |
| Dilution name                             | Sample                                            | Diluted sample                                 | Dilution factor |
| <b>STANDARD</b>                           | <b>12.0</b>                                       | <b>—</b>                                       | <b>1:1.00</b>   |

| Configure assay parameters — General      |                                        |                                                  |  |
|-------------------------------------------|----------------------------------------|--------------------------------------------------|--|
| <input type="radio"/> Reaction definition | <input type="radio"/> Reagent / Sample | <input checked="" type="radio"/> Validity checks |  |
| Reaction check: <b>None</b>               |                                        |                                                  |  |
| Rate linearity %: <b>—</b>                |                                        |                                                  |  |

| Configure assay parameters — Calibration     |                                              |                                   |                                       |
|----------------------------------------------|----------------------------------------------|-----------------------------------|---------------------------------------|
| <input type="radio"/> General                | <input checked="" type="radio"/> Calibration | <input type="radio"/> SmartWash   | <input type="radio"/> Results         |
| Assay: <b>BenzSQ</b>                         |                                              | Calibration method: <b>Linear</b> |                                       |
| <input checked="" type="radio"/> Calibrators | <input type="radio"/> Volumes                | <input type="radio"/> Intervals   | <input type="radio"/> Validity checks |
| Calibrator set: <b>DOASQ</b>                 | Calibrator level: <b>(†††)</b>               |                                   | Concentration: <b>0.0</b>             |
| Replicates: <b>2</b> [Range 1 – 3]           | Blank: <b>DOASQ1</b> (DOA MC Neg Cal)        |                                   |                                       |
|                                              | Cal 1: <b>DOASQ2</b> (DOA MC I Cal 1)        |                                   |                                       |
|                                              | Cal 2: <b>DOASQ3</b> (DOA MC I Cal 2)        |                                   |                                       |
|                                              | Cal 3: <b>DOASQ4</b> (DOA MC I Cal 3)        |                                   |                                       |
|                                              | Cal 4: <b>DOASQ5</b> (DOA MC I Cal 4)        |                                   |                                       |

| Configure assay parameters — Calibration |                                          |                                 |                                       |
|------------------------------------------|------------------------------------------|---------------------------------|---------------------------------------|
| <input type="radio"/> Calibrators        | <input checked="" type="radio"/> Volumes | <input type="radio"/> Intervals | <input type="radio"/> Validity checks |
| Calibrator: <b>DOASQ</b>                 | Calibrator level                         | Sample                          | Diluted sample                        |
| Blank: <b>DOASQ1</b>                     | <b>12.0</b>                              | <b>—</b>                        | <b>—</b>                              |
| Cal 1: <b>DOASQ2</b>                     | <b>12.0</b>                              | <b>—</b>                        | <b>—</b>                              |
| Cal 2: <b>DOASQ3</b>                     | <b>12.0</b>                              | <b>—</b>                        | <b>—</b>                              |
| Cal 3: <b>DOASQ4</b>                     | <b>12.0</b>                              | <b>—</b>                        | <b>—</b>                              |
| Cal 4: <b>DOASQ5</b>                     | <b>12.0</b>                              | <b>—</b>                        | <b>—</b>                              |

| Configure assay parameters — Calibration |                               |                                            |                                       |
|------------------------------------------|-------------------------------|--------------------------------------------|---------------------------------------|
| <input type="radio"/> Calibrators        | <input type="radio"/> Volumes | <input checked="" type="radio"/> Intervals | <input type="radio"/> Validity checks |
| Calibration intervals:                   |                               |                                            |                                       |
| Full interval: <b>168</b>                |                               | (hours)                                    |                                       |
| Calibration type:                        |                               |                                            |                                       |
| Adjust type: <b>None</b>                 |                               |                                            |                                       |

| Configure assay parameters — Calibration  |                               |                                 |                                                  |
|-------------------------------------------|-------------------------------|---------------------------------|--------------------------------------------------|
| <input type="radio"/> Calibrators         | <input type="radio"/> Volumes | <input type="radio"/> Intervals | <input checked="" type="radio"/> Validity checks |
| Blank absorbance range: <b>—</b>          |                               |                                 |                                                  |
| Span: <b>Blank</b> — <b>Blank</b>         |                               |                                 |                                                  |
| Span absorbance range: <b>—</b>           |                               |                                 |                                                  |
| Expected cal factor: <b>0.00</b>          |                               |                                 |                                                  |
| Expected cal factor tolerance %: <b>0</b> |                               |                                 |                                                  |

| Configure assay parameters — SmartWash |                                   |                                            |                               |                                      |
|----------------------------------------|-----------------------------------|--------------------------------------------|-------------------------------|--------------------------------------|
| <input type="radio"/> General          | <input type="radio"/> Calibration | <input checked="" type="radio"/> SmartWash | <input type="radio"/> Results | <input type="radio"/> Interpretation |
| Assay: <b>BenzSQ</b>                   |                                   |                                            |                               |                                      |
| COMPONENT                              | REAGENT / ASSAY                   | WASH                                       | Volume                        | Replicates                           |
| R1                                     | BILT0                             | Detergent A                                | 345                           | 1                                    |
| R1                                     | ACETM                             | Detergent A                                | 345                           | 1                                    |
| R1                                     | ALT00                             | Detergent A                                | 345                           | 1                                    |
| R1                                     | AST00                             | Water                                      | 345                           | 1                                    |
| R1                                     | LIP00                             | Water                                      | 345                           | 1                                    |
| R2                                     | ACETM                             | Detergent A                                | 345                           | 1                                    |
| R2                                     | ALT00                             | Detergent A                                | 345                           | 1                                    |
| R2                                     | AST00                             | Water                                      | 345                           | 1                                    |
| R2                                     | LIP00                             | Water                                      | 345                           | 1                                    |
| Cuvette                                | Trig                              | 10% Detergent B                            | 345                           |                                      |

| Configure assay parameters — Results |                                   |                                 |                                          |
|--------------------------------------|-----------------------------------|---------------------------------|------------------------------------------|
| <input type="radio"/> General        | <input type="radio"/> Calibration | <input type="radio"/> SmartWash | <input checked="" type="radio"/> Results |
| Assay: <b>BenzSQ</b>                 |                                   | Assay number: <b>2852</b>       |                                          |
| Dilution default range:              |                                   | Result units: <b>ng/mL</b>      |                                          |
| Low-Linearity: <b>††</b> <b>18.6</b> |                                   |                                 |                                          |
| High-Linearity: <b>1000.0</b>        |                                   |                                 |                                          |

| Configure result units                 |  |
|----------------------------------------|--|
| Assay: <b>BenzSQ</b>                   |  |
| Version: <b>†</b>                      |  |
| Result units: <b>ng/mL</b>             |  |
| Decimal places: <b>1</b> [Range 0 – 4] |  |
| Correlation factor: <b>1.0000</b>      |  |
| Intercept: <b>0.0000</b>               |  |

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1. Hawks RL, Chiang CN, editors. *Urine Testing for Drugs of Abuse*. NIDA Research Monograph 73. Rockville, MD: National Institute on Drug Abuse (NIDA), Department of Health and Human Services; 1986.
2. *Mandatory Guidelines for Federal Workplace Drug Testing Programs*. National Institute on Drug Abuse. Federal Register Vol. 69, No. 71, 2004:19650–61.
3. Baselt RC. *Disposition of Toxic Drugs and Chemicals in Man*, 5th ed. Foster City, CA: Chemical Toxicology Institute; 2000:205–10.
4. Pearman K. Cocaine: a review. *J Laryngol Otol* 1979;93:1191–9.
5. Stewart DJ, Inaba T, Lucassen M, et al. Cocaine metabolism: cocaine and norcocaine hydrolysis by liver and serum esterases. *Clin Pharmacol Ther* 1979;25(4):464–8.
6. Hamilton HE, Wallace JE, Shimek EL, et al. Cocaine and benzoylecgonine excretion in humans. *J Forensic Sci* 1977;22(4):697–707.
7. Van Dyke C, Byck R, Barash PG, et al. Urinary excretion of immunologically reactive metabolite(s) after intranasal administration of cocaine, as followed by enzyme immunoassay. *Clin Chem* 1977;23(2):241–4.
8. Rubenstein KE, Schneider RS, Ullman EF. Homogeneous enzyme immunoassay: a new immunochemical technique. *Biochem Biophys Res Commun* 1972;47:846–51.
9. US Department of Labor, Occupational Safety and Health Administration. 29 CFR Part 1910.1030. *Bloodborne Pathogens*.
10. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*, 5th ed. Washington, DC: US Government Printing Office, January 2007.
11. World Health Organization. *Laboratory Biosafety Manual*, 3rd ed. Geneva: World Health Organization, 2004.
12. Sewell DL, Bove KE, Callihan DR, et al. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Third Edition (M29-A3)*. Wayne, PA: Clinical and Laboratory Standards Institute, 2005.
13. Armbruster DA, Dasgupta A, et al. *Toxicology and Drug Testing in the Clinical Laboratory; Approved Guideline—Vol. 27, No. 15*. Wayne, PA: Clinical and Laboratory Standards Institute, 2007.

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| Configure assay parameters — General                 |  |  |                                                   |  |                                 |                                                  |                               |  |                                      |  |
|------------------------------------------------------|--|--|---------------------------------------------------|--|---------------------------------|--------------------------------------------------|-------------------------------|--|--------------------------------------|--|
| <input checked="" type="radio"/> General             |  |  | <input type="radio"/> Calibration                 |  | <input type="radio"/> SmartWash |                                                  | <input type="radio"/> Results |  | <input type="radio"/> Interpretation |  |
| Assay: CocSQ1                                        |  |  | Type: Photometric                                 |  | Version: †                      |                                                  |                               |  |                                      |  |
| Number: 2880                                         |  |  | Assay availability: Enabled                       |  | Lot                             |                                                  |                               |  |                                      |  |
| Run controls for onboard reagents by: *              |  |  |                                                   |  |                                 |                                                  |                               |  |                                      |  |
| <input checked="" type="radio"/> Reaction definition |  |  | <input type="radio"/> Reagent / Sample            |  |                                 | <input type="radio"/> Validity checks            |                               |  |                                      |  |
| Reaction mode: Rate up                               |  |  | Primary                                           |  |                                 | Secondary                                        |                               |  | Read times                           |  |
| Wavelength: 340                                      |  |  | / ††                                              |  |                                 | Main: 18 – 21                                    |                               |  | Flux: — —                            |  |
| Last required read: 21                               |  |  |                                                   |  |                                 |                                                  |                               |  |                                      |  |
| Absorbance range: -0.1000 – 3.2000                   |  |  | Color correction: — —                             |  |                                 |                                                  |                               |  |                                      |  |
| Sample blank type: None                              |  |  |                                                   |  |                                 |                                                  |                               |  |                                      |  |
| <input type="radio"/> Reaction definition            |  |  | <input checked="" type="radio"/> Reagent / Sample |  |                                 | <input type="radio"/> Validity checks            |                               |  |                                      |  |
| Reagent: COCM9                                       |  |  | Reagent volume: 100                               |  |                                 | R1                                               |                               |  | R2                                   |  |
| Diluent: None                                        |  |  | Water volume: —                                   |  |                                 |                                                  |                               |  |                                      |  |
| Diluent dispense mode: Type 0                        |  |  | Dispense mode: Type 0                             |  |                                 | Type 0                                           |                               |  | Type 0                               |  |
| Dilution name: Sample                                |  |  | Dilution factor: 1:1.00                           |  |                                 |                                                  |                               |  | Default dilution                     |  |
| STANDARD : 9.2                                       |  |  |                                                   |  |                                 |                                                  |                               |  |                                      |  |
| <input type="radio"/> Reaction definition            |  |  | <input type="radio"/> Reagent / Sample            |  |                                 | <input checked="" type="radio"/> Validity checks |                               |  |                                      |  |
| Reaction check: None                                 |  |  |                                                   |  |                                 |                                                  |                               |  |                                      |  |
| Rate linearity %: —                                  |  |  |                                                   |  |                                 |                                                  |                               |  |                                      |  |

| Configure assay parameters — Calibration     |  |  |                                              |  |                                            |  |                                                  |  |                                      |  |
|----------------------------------------------|--|--|----------------------------------------------|--|--------------------------------------------|--|--------------------------------------------------|--|--------------------------------------|--|
| <input type="radio"/> General                |  |  | <input checked="" type="radio"/> Calibration |  | <input type="radio"/> SmartWash            |  | <input type="radio"/> Results                    |  | <input type="radio"/> Interpretation |  |
| Assay: CocSQ1                                |  |  | Calibration method: Linear                   |  |                                            |  |                                                  |  |                                      |  |
| <input checked="" type="radio"/> Calibrators |  |  | <input type="radio"/> Volumes                |  | <input type="radio"/> Intervals            |  | <input type="radio"/> Validity checks            |  |                                      |  |
| Calibrator set: DOASQ                        |  |  | Calibrator level: (†††)                      |  | Concentration: †                           |  |                                                  |  |                                      |  |
| Replicates: 2 [Range 1 – 3]                  |  |  | Blank: DOASQ1 (DOA MC Neg Cal)               |  | 0.0                                        |  |                                                  |  |                                      |  |
|                                              |  |  | Cal 1: DOASQ2 (DOA MC I Cal 1)               |  |                                            |  |                                                  |  |                                      |  |
|                                              |  |  | Cal 2: DOASQ3 (DOA MC I Cal 2)               |  |                                            |  |                                                  |  |                                      |  |
|                                              |  |  | Cal 3: DOASQ4 (DOA MC I Cal 3)               |  |                                            |  |                                                  |  |                                      |  |
|                                              |  |  | Cal 4: DOASQ5 (DOA MC I Cal 4)               |  |                                            |  |                                                  |  |                                      |  |
| <input type="radio"/> Calibrators            |  |  | <input checked="" type="radio"/> Volumes     |  | <input type="radio"/> Intervals            |  | <input type="radio"/> Validity checks            |  |                                      |  |
| Calibrator: DOASQ                            |  |  | Calibrator level                             |  | Sample                                     |  | Diluted sample                                   |  | Diluent                              |  |
|                                              |  |  | Blank: DOASQ1                                |  | 9.2                                        |  | —                                                |  | —                                    |  |
|                                              |  |  | Cal 1: DOASQ2                                |  | 9.2                                        |  | —                                                |  | —                                    |  |
|                                              |  |  | Cal 2: DOASQ3                                |  | 9.2                                        |  | —                                                |  | —                                    |  |
|                                              |  |  | Cal 3: DOASQ4                                |  | 9.2                                        |  | —                                                |  | —                                    |  |
|                                              |  |  | Cal 4: DOASQ5                                |  | 9.2                                        |  | —                                                |  | —                                    |  |
| <input type="radio"/> Calibrators            |  |  | <input type="radio"/> Volumes                |  | <input checked="" type="radio"/> Intervals |  | <input type="radio"/> Validity checks            |  |                                      |  |
| Calibration intervals:                       |  |  | Full interval: 312 (hours)                   |  |                                            |  |                                                  |  |                                      |  |
| Calibration type:                            |  |  | Adjust type: None                            |  |                                            |  |                                                  |  |                                      |  |
| <input type="radio"/> Calibrators            |  |  | <input type="radio"/> Volumes                |  | <input type="radio"/> Intervals            |  | <input checked="" type="radio"/> Validity checks |  |                                      |  |
| Blank absorbance range:                      |  |  | Span: Blank – Blank                          |  |                                            |  |                                                  |  |                                      |  |
| Span absorbance range:                       |  |  | Expected cal factor: 0.00                    |  |                                            |  |                                                  |  |                                      |  |
| Expected cal factor tolerance %:             |  |  | 0                                            |  |                                            |  |                                                  |  |                                      |  |

| Configure assay parameters — SmartWash                                     |                 |                                      |        |                                            |
|----------------------------------------------------------------------------|-----------------|--------------------------------------|--------|--------------------------------------------|
| <input type="radio"/> General                                              |                 | <input type="radio"/> Calibration    |        | <input checked="" type="radio"/> SmartWash |
| <input type="radio"/> Results                                              |                 | <input type="radio"/> Interpretation |        |                                            |
| Assay: CocSQ1                                                              |                 |                                      |        |                                            |
| COMPONENT                                                                  | REAGENT / ASSAY | WASH                                 | Volume | Replicates                                 |
| R1                                                                         | DIG00           | Detergent A                          | 345    | 1                                          |
| R1                                                                         | AMIK9           | Detergent A                          | 345    | 1                                          |
| R1                                                                         | VANCO           | Detergent A                          | 345    | 1                                          |
| R1                                                                         | GENT9           | Detergent A                          | 345    | 1                                          |
| R1                                                                         | TOBRA           | Detergent A                          | 345    | 1                                          |
| R1                                                                         | DGT08           | Detergent A                          | 345    | 1                                          |
| R2                                                                         | DIG00           | Detergent A                          | 345    | 1                                          |
| R2                                                                         | AMIK9           | Detergent A                          | 345    | 1                                          |
| R2                                                                         | VANCO           | Detergent A                          | 345    | 1                                          |
| R2                                                                         | GENT9           | Detergent A                          | 345    | 1                                          |
| R2                                                                         | TOBRA           | Detergent A                          | 345    | 1                                          |
| R2                                                                         | DGT08           | Detergent A                          | 345    | 1                                          |
| Cuvette                                                                    | Trig**          | 10% Detergent B                      | 345    |                                            |
| **Не требуется для программного обеспечения ARCHITECT, версия 7.00 и выше. |                 |                                      |        |                                            |

| Configure assay parameters — Results       |                                   |
|--------------------------------------------|-----------------------------------|
| <input type="radio"/> General              | <input type="radio"/> Calibration |
| <input checked="" type="radio"/> SmartWash | <input type="radio"/> Results     |
| <input type="radio"/> Interpretation       |                                   |
| Assay: CocSQ1                              | Assay number: 2880                |
| Dilution default range:                    | Result units: ng/mL               |
| Low-Linearity: ††                          | 40.0                              |
| High-Linearity:                            | 1000.0                            |

| Configure result units          |  |
|---------------------------------|--|
| Assay: CocSQ1                   |  |
| Version: †                      |  |
| Result units: ng/mL             |  |
| Decimal places: 1 [Range 0 – 4] |  |
| Correlation factor: 1.0000      |  |
| Intercept: 0.0000               |  |

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| Configure assay parameters — General      |                             |            |                       |
|-------------------------------------------|-----------------------------|------------|-----------------------|
| General                                   | Calibration                 | SmartWash  | Results               |
| Assay: CocQ                               | Type: Photometric           | Version: † |                       |
| Number: 2855                              | Assay availability: Enabled |            |                       |
| Run controls for onboard reagents by: Lot |                             |            |                       |
| Reaction definition                       |                             |            |                       |
| Reaction mode: Rate up                    | Primary                     | Secondary  | Read times            |
| Wavelength: 340                           | /                           | ‡‡         | Main: 18 – 21         |
| Last required read: 21                    |                             |            | Flex: — —             |
| Absorbance range: -0.1000 – 3.2000        |                             |            | Color correction: — — |
| Sample blank type: None                   |                             |            |                       |

| Reaction definition           | Reagent / Sample      | Validity checks |
|-------------------------------|-----------------------|-----------------|
| Reagent: COCM9                | Reagent volume: 100   | R1 R2           |
| Diluent: None                 | Water volume: —       |                 |
| Diluent dispense mode: Type 0 | Dispense mode: Type 0 | Type 0          |
| Dilution name: Sample         | Diluted sample        | Diluent         |
| STANDARD : 8.0                | Water                 | Dilution factor |
|                               |                       | 1:1.00          |

| Reaction definition  | Reagent / Sample | Validity checks |
|----------------------|------------------|-----------------|
| Reaction check: None |                  |                 |
| Rate linearity %: —  |                  |                 |

| Configure assay parameters — Calibration |                               |                  |         |
|------------------------------------------|-------------------------------|------------------|---------|
| General                                  | Calibration                   | SmartWash        | Results |
| Assay: CocQ                              | Calibration method: Linear    |                  |         |
| Calibrators                              |                               |                  |         |
| Calibrator set: DOAQ                     | Calibrator level: (†††)       | Concentration: † |         |
|                                          | Blank: DOAQ1 (DOA MC Neg Cal) | 0                |         |
|                                          | Cal 1: DOAQ2 (DOA MC I Cal 2) |                  |         |
| Replicates: 2                            | [Range 1 – 3]                 |                  |         |

| Calibrators      | Volumes          | Intervals | Validity checks |
|------------------|------------------|-----------|-----------------|
| Calibrator: DOAQ | Calibrator level | Sample    | Diluted sample  |
|                  | Blank: DOAQ1     | 8.0       |                 |
|                  | Cal 1: DOAQ2     | 8.0       |                 |

| Calibrators            | Volumes | Intervals                  | Validity checks |
|------------------------|---------|----------------------------|-----------------|
| Calibration intervals: |         |                            |                 |
|                        |         | Full interval: 312 (hours) |                 |
| Calibration type:      |         |                            |                 |
|                        |         | Adjust type: None          |                 |

| Calibrators                        | Volumes | Intervals                 | Validity checks |
|------------------------------------|---------|---------------------------|-----------------|
| Blank absorbance range:            |         |                           |                 |
|                                    |         | Span: Blank – Blank       |                 |
| Span absorbance range:             |         |                           |                 |
|                                    |         | Expected cal factor: 0.00 |                 |
| Expected cal factor tolerance %: 0 |         |                           |                 |

| Configure assay parameters — SmartWash |                 |                 |         |                |
|----------------------------------------|-----------------|-----------------|---------|----------------|
| General                                | Calibration     | SmartWash       | Results | Interpretation |
| Assay: CocQ                            |                 |                 |         |                |
| COMPONENT                              | REAGENT / ASSAY | WASH            | Volume  | Replicates     |
| R1                                     | DIG00           | Detergent A     | 345     | 1              |
| R1                                     | AMIK9           | Detergent A     | 345     | 1              |
| R1                                     | VANCO           | Detergent A     | 345     | 1              |
| R1                                     | GENT9           | Detergent A     | 345     | 1              |
| R1                                     | TOBRA           | Detergent A     | 345     | 1              |
| R1                                     | DGT0B           | Detergent A     | 345     | 1              |
| R2                                     | DIG00           | Detergent A     | 345     | 1              |
| R2                                     | AMIK9           | Detergent A     | 345     | 1              |
| R2                                     | VANCO           | Detergent A     | 345     | 1              |
| R2                                     | GENT9           | Detergent A     | 345     | 1              |
| R2                                     | TOBRA           | Detergent A     | 345     | 1              |
| R2                                     | DGT0B           | Detergent A     | 345     | 1              |
| Cuvette                                | Trig**          | 10% Detergent B | 345     |                |

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1. US Department of Labor, Occupational Safety and Health Administration. 29 CFR Part 1910.1030. *Bloodborne Pathogens*.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*, 5th ed. Washington, DC: US Government Printing Office, January 2007.
3. World Health Organization. *Laboratory Biosafety Manual*, 3rd ed. Geneva: World Health Organization, 2004.
4. Sewell DL, Bove KE, Callihan DR, et al. *Protection of Laboratory Workers from Occupationally Acquired Infections: Approved Guideline—Third Edition (M29-A3)*. Wayne, PA: Clinical and Laboratory Standards Institute, 2005.

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1. Hawks RL, Chiang CN, editors. *Urine Testing for Drugs of Abuse*. NIDA Research Monograph 73. Rockville, MD: National Institute on Drug Abuse (NIDA), Department of Health and Human Services; 1986.
2. *Mandatory Guidelines for Federal Workplace Drug Testing Programs*. National Institute on Drug Abuse. Federal Register Vol. 69, No. 71, 2004:19650–61.
3. Baselt RC. *Disposition of Toxic Drugs and Chemicals in Man*, 5th ed. Foster City, CA: Chemical Toxicology Institute; 2000:812–7.
4. Perez-Reyes M, Di Guiseppe S, Mason AP, Davis KH. Passive inhalation of marijuana smoke and urinary excretion of cannabinoids. *Clin Pharmacol Ther* 1983;34:36.
5. Ferslew KE, Manno JE, Manno BR. Determination of urinary cannabinoid metabolites following incidental exposure to marijuana smoke. *Res Commun Substance Abuse*, 1983:289.
6. Rubenstein KE, Schneider RS, Ullman EF. Homogeneous enzyme immunoassay: a new immunochemical technique. *Biochem Biophys Res Commun* 1972;47:846–51.
7. US Department of Labor, Occupational Safety and Health Administration. 29 CFR Part 1910.1030, *Occupational Exposure to Bloodborne Pathogens*.
8. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; January 2007.
9. World Health Organization. *Laboratory Biosafety Manual*. Geneva: World Health Organization. 2004.
10. Sewell DL, Bove KE, Callihan DR, et al. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Third Edition (M29-A3)*. Wayne, PA: Clinical and Laboratory Standards Institute, 2005.
11. Shoemaker MJ, Earley RJ, Gorsky JE, et al. *Urine Drug Testing in the Clinical Laboratory; Approved Guideline—Vol. 19, No. 6*. Wayne, PA: Clinical and Laboratory Standards Institute, 1999(17).
12. Blanc JA, Manneh VA, Ernst R, et al. Adsorption losses from urine-based cannabinoid calibrators during routine use. *Clin Chem* 1993;39(8):1705–12.
13. Roth KDW, Siegel NA, et al. Investigation of the effects of solution composition and container material type on the loss of 11-nor- $\Delta^9$ -THC-9-carboxylic acid. *J Anal Tox* 1996;20:291–300.

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## Configure assay parameters — General

|                                                                                                                                                                               |                     |                     |             |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|---------------------|-------------|
| <input checked="" type="radio"/> General <input type="radio"/> Calibration <input type="radio"/> SmartWash <input type="radio"/> Results <input type="radio"/> Interpretation |                     |                     |             |
| Assay:                                                                                                                                                                        | THC5Q               | Type:               | Photometric |
| Number:                                                                                                                                                                       | 2858                | Assay availability: | Enabled     |
| <input checked="" type="radio"/> Reaction definition <input type="radio"/> Reagent / Sample <input type="radio"/> Validity checks                                             |                     |                     |             |
| Reaction mode:                                                                                                                                                                | Rate up             |                     |             |
|                                                                                                                                                                               | Primary   Secondary | Read times          |             |
| Wavelength:                                                                                                                                                                   | 340 / ##            | Main:               | 20 – 24     |
| Last required read:                                                                                                                                                           | 24                  | Flex:               | — —         |
| Absorbance range:                                                                                                                                                             | -0.1000 – 3.2000    | Color correction:   | — —         |
| Sample blank type:                                                                                                                                                            | None                |                     |             |

|                                                                                                                                   |        |                 |                                                      |
|-----------------------------------------------------------------------------------------------------------------------------------|--------|-----------------|------------------------------------------------------|
| <input type="radio"/> Reaction definition <input checked="" type="radio"/> Reagent / Sample <input type="radio"/> Validity checks |        |                 |                                                      |
| Reagent:                                                                                                                          | THCM9  | Reagent volume: | 100   100                                            |
| Diluent:                                                                                                                          | None   | Water volume:   | — —                                                  |
| Diluent dispense mode:                                                                                                            | Type 0 | Dispense mode:  | Type 0   Type 0                                      |
| Dilution name                                                                                                                     | Sample | Diluted sample  | Diluent   Water   Dilution factor   Default dilution |
| STANDARD :                                                                                                                        | 12.0   | —               | — = 1:1.00   ●                                       |
| <input type="radio"/> Reaction definition <input type="radio"/> Reagent / Sample <input checked="" type="radio"/> Validity checks |        |                 |                                                      |
| Reaction check: None                                                                                                              |        |                 |                                                      |
| Rate linearity %: —                                                                                                               |        |                 |                                                      |

## Configure assay parameters — Calibration

|                                                                                                                                                                               |       |                     |                             |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------|---------------------|-----------------------------|
| <input type="radio"/> General <input checked="" type="radio"/> Calibration <input type="radio"/> SmartWash <input type="radio"/> Results <input type="radio"/> Interpretation |       |                     |                             |
| Assay:                                                                                                                                                                        | THC5Q | Calibration method: | Linear                      |
| <input checked="" type="radio"/> Calibrators <input type="radio"/> Volumes <input type="radio"/> Intervals <input type="radio"/> Validity checks                              |       |                     |                             |
| Calibrator set:                                                                                                                                                               | THCQ  | Calibrator level:   | (†††)   Concentration: ‡    |
|                                                                                                                                                                               |       | Blank:              | THCQ1 (DOA MC Neg Cal)      |
|                                                                                                                                                                               |       | Cal 1:              | THCQ2 (Cannabinoids 50 Cal) |
| Replicates:                                                                                                                                                                   | 2     | [Range 1 – 3]       |                             |

|                                                                                                                                                  |       |                |                 |
|--------------------------------------------------------------------------------------------------------------------------------------------------|-------|----------------|-----------------|
| <input type="radio"/> Calibrators <input checked="" type="radio"/> Volumes <input type="radio"/> Intervals <input type="radio"/> Validity checks |       |                |                 |
| Calibrator:                                                                                                                                      | THCQ  | Diluted sample | Diluent   Water |
| Calibrator level                                                                                                                                 |       | Sample         |                 |
| Blank:                                                                                                                                           | THCQ1 | 12.0           | — —             |
| Cal 1:                                                                                                                                           | THCQ2 | 12.0           | — —             |

|                                                                                                                                                  |  |  |  |
|--------------------------------------------------------------------------------------------------------------------------------------------------|--|--|--|
| <input type="radio"/> Calibrators <input type="radio"/> Volumes <input checked="" type="radio"/> Intervals <input type="radio"/> Validity checks |  |  |  |
| Calibration intervals:                                                                                                                           |  |  |  |
| Full interval: 312 (hours)                                                                                                                       |  |  |  |
| Calibration type:                                                                                                                                |  |  |  |
| Adjust type: None                                                                                                                                |  |  |  |

|                                                                                                                                                  |  |  |  |
|--------------------------------------------------------------------------------------------------------------------------------------------------|--|--|--|
| <input type="radio"/> Calibrators <input type="radio"/> Volumes <input type="radio"/> Intervals <input checked="" type="radio"/> Validity checks |  |  |  |
| Blank absorbance range: — —                                                                                                                      |  |  |  |
| Span: Blank — —                                                                                                                                  |  |  |  |
| Span absorbance range: — —                                                                                                                       |  |  |  |
| Expected cal factor: 0.00                                                                                                                        |  |  |  |
| Expected cal factor tolerance %: 0                                                                                                               |  |  |  |

## Configure assay parameters — SmartWash

|                                                                                                                                                                               |                 |                    |                     |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------|--------------------|---------------------|
| <input type="radio"/> General <input type="radio"/> Calibration <input checked="" type="radio"/> SmartWash <input type="radio"/> Results <input type="radio"/> Interpretation |                 |                    |                     |
| Assay:                                                                                                                                                                        | THC5Q           |                    |                     |
| COMPONENT                                                                                                                                                                     | REAGENT / ASSAY | WASH               | Volume   Replicates |
| Cuvette                                                                                                                                                                       | Trig            | 10% Detergent B*** | 345                 |
| *** Select "Detergent B" for software prior to version 2.2.                                                                                                                   |                 |                    |                     |

## Configure assay parameters — Results

|                                                                                                                                                                               |       |               |   |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------|---------------|---|
| <input type="radio"/> General <input type="radio"/> Calibration <input type="radio"/> SmartWash <input checked="" type="radio"/> Results <input type="radio"/> Interpretation |       |               |   |
| Assay:                                                                                                                                                                        | THC5Q | Result units: | — |
| Assay defaults:                                                                                                                                                               |       |               |   |
| Low-Linearity:                                                                                                                                                                |       |               |   |
| High-Linearity:                                                                                                                                                               |       |               |   |

## Configure assay parameters — Interpretation

|                                                                                                                                                                               |        |                          |  |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------|--------------------------|--|
| <input type="radio"/> General <input type="radio"/> Calibration <input type="radio"/> SmartWash <input type="radio"/> Results <input checked="" type="radio"/> Interpretation |        |                          |  |
| Assay: THC5Q                                                                                                                                                                  |        |                          |  |
| Name:                                                                                                                                                                         | Range: | Results review required: |  |
| Negative                                                                                                                                                                      | < 50   | <input type="checkbox"/> |  |
| Positive                                                                                                                                                                      | >= 50  | <input type="checkbox"/> |  |

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## Configure assay parameters — General

|                                                                        |                                        |                                       |                               |                                      |
|------------------------------------------------------------------------|----------------------------------------|---------------------------------------|-------------------------------|--------------------------------------|
| <input checked="" type="radio"/> General                               | <input type="radio"/> Calibration      | <input type="radio"/> SmartWash       | <input type="radio"/> Results | <input type="radio"/> Interpretation |
| Assay: <b>THC1SQ</b> Type: <b>Photometric</b> Version: <b>†</b>        |                                        |                                       |                               |                                      |
| Number: <b>2856</b> Assay availability: <b>Enabled</b>                 |                                        |                                       |                               |                                      |
| <input checked="" type="radio"/> Reaction definition                   | <input type="radio"/> Reagent / Sample | <input type="radio"/> Validity checks |                               |                                      |
| Reaction mode: <b>Rate up</b>                                          |                                        |                                       |                               |                                      |
| Primary Secondary Read times                                           |                                        |                                       |                               |                                      |
| Wavelength: <b>340</b> / <b>‡‡</b> Main: <b>20 – 24</b>                |                                        |                                       |                               |                                      |
| Last required read: <b>24</b> Flex: <b>— —</b>                         |                                        |                                       |                               |                                      |
| Absorbance range: <b>-0.1000 – 3.2000</b> Color correction: <b>— —</b> |                                        |                                       |                               |                                      |
| Sample blank type: <b>None</b>                                         |                                        |                                       |                               |                                      |

|                                                                                    |                                                   |                                       |  |  |
|------------------------------------------------------------------------------------|---------------------------------------------------|---------------------------------------|--|--|
| <input type="radio"/> Reaction definition                                          | <input checked="" type="radio"/> Reagent / Sample | <input type="radio"/> Validity checks |  |  |
| Reagent: <b>THCM9</b> Reagent volume: <b>100</b> <b>100</b>                        |                                                   |                                       |  |  |
| Diluent: <b>None</b> Water volume: <b>— —</b>                                      |                                                   |                                       |  |  |
| Diluent dispense mode: <b>Type 0</b> Dispense mode: <b>Type 0</b> <b>Type 0</b>    |                                                   |                                       |  |  |
| Dilution name Sample Diluted sample Diluent Water Dilution factor Default dilution |                                                   |                                       |  |  |
| <b>STANDARD : 12.0 — — — = 1:1.00 ●</b>                                            |                                                   |                                       |  |  |

|                                           |                                        |                                                  |
|-------------------------------------------|----------------------------------------|--------------------------------------------------|
| <input type="radio"/> Reaction definition | <input type="radio"/> Reagent / Sample | <input checked="" type="radio"/> Validity checks |
| Reaction check: <b>None</b>               |                                        |                                                  |
| Rate linearity %: <b>—</b>                |                                        |                                                  |

## Configure assay parameters — Calibration

|                                                                                        |                                              |                                 |                                       |                                      |
|----------------------------------------------------------------------------------------|----------------------------------------------|---------------------------------|---------------------------------------|--------------------------------------|
| <input type="radio"/> General                                                          | <input checked="" type="radio"/> Calibration | <input type="radio"/> SmartWash | <input type="radio"/> Results         | <input type="radio"/> Interpretation |
| Assay: <b>THC1SQ</b> Calibration method: <b>Linear</b>                                 |                                              |                                 |                                       |                                      |
| <input checked="" type="radio"/> Calibrators                                           | <input type="radio"/> Volumes                | <input type="radio"/> Intervals | <input type="radio"/> Validity checks |                                      |
| Calibrator set: <b>THC100SQ</b> Calibrator level: <b>(†††)</b> Concentration: <b>‡</b> |                                              |                                 |                                       |                                      |
| Blank: <b>THC100SQ1 (DOA MC Neg Cal)</b>                                               |                                              |                                 |                                       |                                      |
| Cal 1: <b>THC100SQ2 (Cannabinoids 20 Cal)</b>                                          |                                              |                                 |                                       |                                      |
| Cal 2: <b>THC100SQ3 (Cannabinoids 50 Cal)</b>                                          |                                              |                                 |                                       |                                      |
| Cal 3: <b>THC100SQ4 (Cannabinoids 100 Cal)</b>                                         |                                              |                                 |                                       |                                      |
| Replicates: <b>2</b> [Range 1 – 3]                                                     |                                              |                                 |                                       |                                      |

|                                                                |                                          |                                 |                                       |  |
|----------------------------------------------------------------|------------------------------------------|---------------------------------|---------------------------------------|--|
| <input type="radio"/> Calibrators                              | <input checked="" type="radio"/> Volumes | <input type="radio"/> Intervals | <input type="radio"/> Validity checks |  |
| Calibrator: <b>THC100SQ</b>                                    |                                          |                                 |                                       |  |
| Calibrator level Sample Diluted sample Diluent Water           |                                          |                                 |                                       |  |
| Blank: <b>THC100SQ1</b> <b>12.0</b> <b>—</b> <b>—</b> <b>—</b> |                                          |                                 |                                       |  |
| Cal 1: <b>THC100SQ2</b> <b>12.0</b> <b>—</b> <b>—</b> <b>—</b> |                                          |                                 |                                       |  |
| Cal 2: <b>THC100SQ3</b> <b>12.0</b> <b>—</b> <b>—</b> <b>—</b> |                                          |                                 |                                       |  |
| Cal 3: <b>THC100SQ4</b> <b>12.0</b> <b>—</b> <b>—</b> <b>—</b> |                                          |                                 |                                       |  |

|                                   |                               |                                            |                                       |  |
|-----------------------------------|-------------------------------|--------------------------------------------|---------------------------------------|--|
| <input type="radio"/> Calibrators | <input type="radio"/> Volumes | <input checked="" type="radio"/> Intervals | <input type="radio"/> Validity checks |  |
| Calibration intervals:            |                               |                                            |                                       |  |
| Full interval: <b>312</b> (hours) |                               |                                            |                                       |  |
| Calibration type:                 |                               |                                            |                                       |  |
| Adjust type: <b>None</b>          |                               |                                            |                                       |  |

|                                           |                               |                                 |                                                  |
|-------------------------------------------|-------------------------------|---------------------------------|--------------------------------------------------|
| <input type="radio"/> Calibrators         | <input type="radio"/> Volumes | <input type="radio"/> Intervals | <input checked="" type="radio"/> Validity checks |
| Blank absorbance range: <b>— —</b>        |                               |                                 |                                                  |
| Span: <b>Blank</b> <b>— —</b>             |                               |                                 |                                                  |
| Span absorbance range: <b>— —</b>         |                               |                                 |                                                  |
| Expected cal factor: <b>0.00</b>          |                               |                                 |                                                  |
| Expected cal factor tolerance %: <b>0</b> |                               |                                 |                                                  |

## Configure assay parameters — SmartWash

|                                                             |                                   |                                            |                               |                                      |
|-------------------------------------------------------------|-----------------------------------|--------------------------------------------|-------------------------------|--------------------------------------|
| <input type="radio"/> General                               | <input type="radio"/> Calibration | <input checked="" type="radio"/> SmartWash | <input type="radio"/> Results | <input type="radio"/> Interpretation |
| Assay: <b>THC1SQ</b>                                        |                                   |                                            |                               |                                      |
| COMPONENT                                                   | REAGENT / ASSAY                   | WASH                                       | Volume                        | Replicates                           |
| Cuvette                                                     | Trig                              | <b>10% Detergent B***</b>                  | <b>345</b>                    |                                      |
| *** Select "Detergent B" for software prior to version 2.2. |                                   |                                            |                               |                                      |

## Configure assay parameters — Results

|                                                 |                                   |                                 |                                          |                                      |
|-------------------------------------------------|-----------------------------------|---------------------------------|------------------------------------------|--------------------------------------|
| <input type="radio"/> General                   | <input type="radio"/> Calibration | <input type="radio"/> SmartWash | <input checked="" type="radio"/> Results | <input type="radio"/> Interpretation |
| Assay: <b>THC1SQ</b> Result units: <b>ng/mL</b> |                                   |                                 |                                          |                                      |
| Assay defaults:                                 |                                   |                                 |                                          |                                      |
| Low-Linearity:†† <b>10.0</b>                    |                                   |                                 |                                          |                                      |
| High-Linearity: <b>100.0</b>                    |                                   |                                 |                                          |                                      |

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## Configure assay parameters — General

|                                                      |  |                                        |  |                                       |  |                               |  |                                      |  |
|------------------------------------------------------|--|----------------------------------------|--|---------------------------------------|--|-------------------------------|--|--------------------------------------|--|
| <input checked="" type="radio"/> General             |  | <input type="radio"/> Calibration      |  | <input type="radio"/> SmartWash       |  | <input type="radio"/> Results |  | <input type="radio"/> Interpretation |  |
| Assay: <b>THC2SQ</b>                                 |  | Type: <b>Photometric</b>               |  | Version: <b>†</b>                     |  |                               |  |                                      |  |
| Number: <b>2857</b>                                  |  | Assay availability: <b>Enabled</b>     |  |                                       |  |                               |  |                                      |  |
| <input checked="" type="radio"/> Reaction definition |  | <input type="radio"/> Reagent / Sample |  | <input type="radio"/> Validity checks |  |                               |  |                                      |  |
| Reaction mode: <b>Rate up</b>                        |  | Primary                                |  | Secondary                             |  | Read times                    |  |                                      |  |
| Wavelength: <b>340</b> / <b>‡</b>                    |  |                                        |  |                                       |  | Main: <b>20 – 24</b>          |  |                                      |  |
| Last required read: <b>24</b>                        |  |                                        |  |                                       |  | Flex: <b>— —</b>              |  |                                      |  |
| Absorbance range: <b>-0.1000 – 3.2000</b>            |  |                                        |  |                                       |  | Color correction: <b>— —</b>  |  |                                      |  |
| Sample blank type: <b>None</b>                       |  |                                        |  |                                       |  |                               |  |                                      |  |

|                                           |  |                                                   |  |                                       |  |               |  |                 |  |
|-------------------------------------------|--|---------------------------------------------------|--|---------------------------------------|--|---------------|--|-----------------|--|
| <input type="radio"/> Reaction definition |  | <input checked="" type="radio"/> Reagent / Sample |  | <input type="radio"/> Validity checks |  |               |  |                 |  |
| Reagent: <b>THCM9</b>                     |  | Reagent volume: <b>100</b>                        |  | R1                                    |  | R2            |  |                 |  |
| Diluent: <b>None</b>                      |  | Water volume: <b>—</b>                            |  |                                       |  |               |  |                 |  |
| Diluent dispense mode: <b>Type 0</b>      |  | Dispense mode: <b>Type 0</b>                      |  | <b>Type 0</b>                         |  | <b>Type 0</b> |  |                 |  |
| Dilution name                             |  | Sample                                            |  | Diluent                               |  | Water         |  | Dilution factor |  |
| <b>STANDARD</b>                           |  | <b>: 6.0</b>                                      |  | <b>—</b>                              |  | <b>—</b>      |  | <b>1:1.00</b>   |  |
|                                           |  |                                                   |  |                                       |  |               |  | <b>•</b>        |  |

|                                           |  |                                        |  |                                                  |  |  |  |  |  |
|-------------------------------------------|--|----------------------------------------|--|--------------------------------------------------|--|--|--|--|--|
| <input type="radio"/> Reaction definition |  | <input type="radio"/> Reagent / Sample |  | <input checked="" type="radio"/> Validity checks |  |  |  |  |  |
| Reaction check: <b>None</b>               |  |                                        |  |                                                  |  |  |  |  |  |
| Rate linearity %: <b>—</b>                |  |                                        |  |                                                  |  |  |  |  |  |

## Configure assay parameters — Calibration

|                                              |  |                                                |  |                                                |  |                                       |  |                                      |  |
|----------------------------------------------|--|------------------------------------------------|--|------------------------------------------------|--|---------------------------------------|--|--------------------------------------|--|
| <input type="radio"/> General                |  | <input checked="" type="radio"/> Calibration   |  | <input type="radio"/> SmartWash                |  | <input type="radio"/> Results         |  | <input type="radio"/> Interpretation |  |
| Assay: <b>THC2SQ</b>                         |  | Calibration method: <b>Linear</b>              |  |                                                |  |                                       |  |                                      |  |
| <input checked="" type="radio"/> Calibrators |  | <input type="radio"/> Volumes                  |  | <input type="radio"/> Intervals                |  | <input type="radio"/> Validity checks |  |                                      |  |
| Calibrator set: <b>THC200SQ</b>              |  | Calibrator level: <b>(†††)</b>                 |  | Concentration: <b>†</b>                        |  |                                       |  |                                      |  |
| Replicates: <b>2</b> [Range 1 – 3]           |  | Blank: <b>THC200SQ1 (DOA MC Neg Cal)</b>       |  | Cal 1: <b>THC200SQ2 (Cannabinoids 50 Cal)</b>  |  |                                       |  |                                      |  |
|                                              |  | Cal 2: <b>THC200SQ3 (Cannabinoids 100 Cal)</b> |  | Cal 3: <b>THC200SQ4 (Cannabinoids 200 Cal)</b> |  |                                       |  |                                      |  |

|                                   |  |                                          |  |                                 |  |                                       |  |          |  |
|-----------------------------------|--|------------------------------------------|--|---------------------------------|--|---------------------------------------|--|----------|--|
| <input type="radio"/> Calibrators |  | <input checked="" type="radio"/> Volumes |  | <input type="radio"/> Intervals |  | <input type="radio"/> Validity checks |  |          |  |
| Calibrator: <b>THC200SQ</b>       |  | Calibrator level                         |  | Sample                          |  | Diluted sample                        |  | Diluent  |  |
| Blank: <b>THC200SQ1</b>           |  | <b>6.0</b>                               |  | <b>—</b>                        |  | <b>—</b>                              |  | <b>—</b> |  |
| Cal 1: <b>THC200SQ2</b>           |  | <b>6.0</b>                               |  | <b>—</b>                        |  | <b>—</b>                              |  | <b>—</b> |  |
| Cal 2: <b>THC200SQ3</b>           |  | <b>6.0</b>                               |  | <b>—</b>                        |  | <b>—</b>                              |  | <b>—</b> |  |
| Cal 3: <b>THC200SQ4</b>           |  | <b>6.0</b>                               |  | <b>—</b>                        |  | <b>—</b>                              |  | <b>—</b> |  |

|                                   |  |                               |  |                                            |  |                                       |  |  |  |
|-----------------------------------|--|-------------------------------|--|--------------------------------------------|--|---------------------------------------|--|--|--|
| <input type="radio"/> Calibrators |  | <input type="radio"/> Volumes |  | <input checked="" type="radio"/> Intervals |  | <input type="radio"/> Validity checks |  |  |  |
| Calibration intervals:            |  |                               |  |                                            |  |                                       |  |  |  |
| Full interval: <b>312</b> (hours) |  |                               |  |                                            |  |                                       |  |  |  |
| Calibration type:                 |  |                               |  |                                            |  |                                       |  |  |  |
| Adjust type: <b>None</b>          |  |                               |  |                                            |  |                                       |  |  |  |

|                                   |  |                               |  |                                 |  |                                                  |  |  |  |
|-----------------------------------|--|-------------------------------|--|---------------------------------|--|--------------------------------------------------|--|--|--|
| <input type="radio"/> Calibrators |  | <input type="radio"/> Volumes |  | <input type="radio"/> Intervals |  | <input checked="" type="radio"/> Validity checks |  |  |  |
| Blank absorbance range:           |  | <b>—</b>                      |  | <b>—</b>                        |  | <b>—</b>                                         |  |  |  |
| Span:                             |  | <b>Blank</b>                  |  | <b>—</b>                        |  | <b>—</b>                                         |  |  |  |
| Span absorbance range:            |  | <b>—</b>                      |  | <b>—</b>                        |  | <b>—</b>                                         |  |  |  |
| Expected cal factor:              |  | <b>0.00</b>                   |  | <b>—</b>                        |  | <b>—</b>                                         |  |  |  |
| Expected cal factor tolerance %:  |  | <b>0</b>                      |  | <b>—</b>                        |  | <b>—</b>                                         |  |  |  |

## Configure assay parameters — SmartWash

|                                                             |                 |                                   |        |                                            |  |                               |  |                                      |  |
|-------------------------------------------------------------|-----------------|-----------------------------------|--------|--------------------------------------------|--|-------------------------------|--|--------------------------------------|--|
| <input type="radio"/> General                               |                 | <input type="radio"/> Calibration |        | <input checked="" type="radio"/> SmartWash |  | <input type="radio"/> Results |  | <input type="radio"/> Interpretation |  |
| Assay: <b>THC2SQ</b>                                        |                 |                                   |        |                                            |  |                               |  |                                      |  |
| COMPONENT                                                   | REAGENT / ASSAY | WASH                              | Volume | Replicates                                 |  |                               |  |                                      |  |
| Cuvette                                                     | Trig            | 10% Detergent B***                | 345    |                                            |  |                               |  |                                      |  |
| *** Select "Detergent B" for software prior to version 2.2. |                 |                                   |        |                                            |  |                               |  |                                      |  |

## Configure assay parameters — Results

|                               |  |                                   |  |                                 |  |                                          |  |                                      |  |
|-------------------------------|--|-----------------------------------|--|---------------------------------|--|------------------------------------------|--|--------------------------------------|--|
| <input type="radio"/> General |  | <input type="radio"/> Calibration |  | <input type="radio"/> SmartWash |  | <input checked="" type="radio"/> Results |  | <input type="radio"/> Interpretation |  |
| Assay: <b>THC2SQ</b>          |  | Result units: <b>ng/mL</b>        |  |                                 |  |                                          |  |                                      |  |
| Assay defaults:               |  |                                   |  |                                 |  |                                          |  |                                      |  |
| Low-Linearity:††              |  | <b>10.0</b>                       |  |                                 |  |                                          |  |                                      |  |
| High-Linearity:               |  | <b>200.0</b>                      |  |                                 |  |                                          |  |                                      |  |

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| Assay Configuration: SmartWash Page |            |      |      |
|-------------------------------------|------------|------|------|
| Rgt Probe                           | Reagent    | Wash | Vol  |
|                                     | ____       | ____ | ____ |
| Cuvette                             | Assay Name | Wash | Vol  |
|                                     | ____       | ____ | ____ |
| Sample Probe                        | Wash       |      |      |
|                                     | ____       |      |      |

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| Assay Configuration: SmartWash Page |            |      |      |
|-------------------------------------|------------|------|------|
| Rgt Probe                           | Reagent    | Wash | Vol  |
| ____                                | ____       | ____ | ____ |
| Cuvette                             | Assay Name | Wash | Vol  |
| ____                                | ____       | ____ | ____ |
| Sample Probe                        | Wash       |      |      |
| ____                                | ____       |      |      |

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1. US Department of Labor, Occupational Safety and Health Administration. 29 CFR Part 1910.1030. *Bloodborne Pathogens*.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*, 5th ed. Washington, DC: US Government Printing Office, January 2007.
3. World Health Organization. *Laboratory Biosafety Manual*, 3rd ed. Geneva: World Health Organization, 2004.
4. Sewell DL, Bove KE, Callihan DR, et al. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Third Edition (M29-A3)*. Wayne, PA: Clinical and Laboratory Standards Institute, 2005.
5. Blanc JA, Manneh VA, et al. Adsorption losses from urine-based cannabinoid calibrators during routine use. *Clin Chem* 1993;39:1705–12.
6. Roth KDW, Siegel NA, et al. Investigation of the effects of solution composition and container material type on the loss of 11-nor- $\Delta^9$ -THC-9-carboxylic acid. *J Anal Tox* 1996;20:291–300.

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1. Hawks RL, Chiang CN, editors. *Urine Testing for Drugs of Abuse*. NIDA Research Monograph 73. Rockville, MD: National Institute on Drug Abuse (NIDA), Department of Health and Human Services;1986.
2. *Mandatory Guidelines for Federal Workplace Drug Testing Programs*. National Institute on Drug Abuse. Federal Register Vol. 69, No. 71, 2004:19650–61.
3. Baselt RC, Cravey RH. *Disposition of Toxic Drugs and Chemicals in Man*. Chemical Toxicology Institute, 4th ed. Foster City, CA;1995.
4. Rubenstein KE, Schneider RS, Ullman EF. Homogeneous enzyme immunoassay: a new immunochemical technique. *Biochem Biophys Res Commun* 1972;47:846–51.
5. US Department of Labor, Occupational Safety and Health Administration. 29 CFR Part 1910.1030. *Bloodborne Pathogens*.
6. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*, 5th ed. Washington, DC: US Government Printing Office, January 2007.
7. World Health Organization. Laboratory Biosafety Manual, 3rd ed. Geneva:World Health Organization, 2004.
8. Sewell DL, Bove KE, Callihan DR, et al. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Third Edition (M29-A3)*. Wayne, PA: Clinical and Laboratory Standards Institute, 2005.
9. Armbruster DA, Dasgupta A, et al. *Toxicology and Drug Testing in the Clinical Laboratory; Approved Guideline—Vol. 27, No. 15*. Wayne, PA: Clinical and Laboratory Standards Institute, 2007.

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| Configure assay parameters — General                 |                                        |                                       |                               |
|------------------------------------------------------|----------------------------------------|---------------------------------------|-------------------------------|
| <input checked="" type="radio"/> General             | <input type="radio"/> Calibration      | <input type="radio"/> SmartWash       | <input type="radio"/> Results |
| Assay: EcstSQ Type: Photometric Version: †           |                                        |                                       |                               |
| Number: 2861                                         |                                        | Assay availability: Enabled           |                               |
| Run controls for onboard reagents by: Lot            |                                        |                                       |                               |
| <input checked="" type="radio"/> Reaction definition | <input type="radio"/> Reagent / Sample | <input type="radio"/> Validity checks |                               |
| Reaction mode: Rate up                               |                                        |                                       |                               |
| Primary                                              |                                        | Secondary                             | Read times                    |
| Wavelength: 340 / ††                                 |                                        | Main: 20 – 24                         |                               |
| Last required read: 24                               |                                        | Flex: —                               |                               |
| Absorbance range: -0.1000 – 3.2000                   |                                        | Color correction: —                   |                               |
| Sample blank type: None                              |                                        |                                       |                               |

| O Reaction definition         |        |                       |         | ● Reagent / Sample |                 | O Validity checks |  |
|-------------------------------|--------|-----------------------|---------|--------------------|-----------------|-------------------|--|
| Reagent: ECSTM                |        | Reagent volume: 100   |         | R1                 | R2              |                   |  |
| Diluent: None                 |        | Water volume: —       |         |                    |                 |                   |  |
| Diluent dispense mode: Type 0 |        | Dispense mode: Type 0 |         | Type 0             | Type 0          |                   |  |
| Dilution name                 | Sample | Diluted sample        | Diluent | Water              | Dilution factor | Default dilution  |  |
| STANDARD                      | 8.0    | —                     | —       | —                  | 1:1.00          | ●                 |  |
| O Reaction definition         |        |                       |         | O Reagent / Sample |                 | ● Validity checks |  |
| Reaction check: None          |        |                       |         |                    |                 |                   |  |
| Rate linearity %: —           |        |                       |         |                    |                 |                   |  |

| Configure assay parameters — Calibration     |                                              |                                          |                                       |
|----------------------------------------------|----------------------------------------------|------------------------------------------|---------------------------------------|
| <input type="radio"/> General                | <input checked="" type="radio"/> Calibration | <input type="radio"/> SmartWash          | <input type="radio"/> Results         |
| Assay: EcstSQ Calibration method: Linear     |                                              |                                          |                                       |
| <input checked="" type="radio"/> Calibrators | <input type="radio"/> Volumes                | <input type="radio"/> Intervals          | <input type="radio"/> Validity checks |
| Calibrator set: ECSTSQ                       |                                              | Calibrator level: (†††) Concentration: ‡ |                                       |
| Replicates: 2 [Range 1 – 3]                  |                                              | Blank: ECSTSQ1 (DOA MC Neg Cal) 0.0      |                                       |
|                                              |                                              | Cal 1: ECSTSQ2 (Ecstasy 250 Cal)         |                                       |
|                                              |                                              | Cal 2: ECSTSQ3 (Ecstasy 500 Cal)         |                                       |
|                                              |                                              | Cal 3: ECSTSQ4 (Ecstasy 750 Cal)         |                                       |
|                                              |                                              | Cal 4: ECSTSQ5 (Ecstasy 1000 Cal)        |                                       |

| O Calibrators      |  |                  |  | ● Volumes |                | O Intervals |       | O Validity checks |  |
|--------------------|--|------------------|--|-----------|----------------|-------------|-------|-------------------|--|
| Calibrator: ECSTSQ |  | Calibrator level |  | Sample    | Diluted sample | Diluent     | Water |                   |  |
| Blank: ECSTSQ1     |  | 8.0              |  | —         | —              | —           | —     |                   |  |
| Cal 1: ECSTSQ2     |  | 8.0              |  | —         | —              | —           | —     |                   |  |
| Cal 2: ECSTSQ3     |  | 8.0              |  | —         | —              | —           | —     |                   |  |
| Cal 3: ECSTSQ4     |  | 8.0              |  | —         | —              | —           | —     |                   |  |
| Cal 4: ECSTSQ5     |  | 8.0              |  | —         | —              | —           | —     |                   |  |

| O Calibrators              |  |  |  | O Volumes |  | ● Intervals |  | O Validity checks |  |
|----------------------------|--|--|--|-----------|--|-------------|--|-------------------|--|
| Calibration intervals:     |  |  |  |           |  |             |  |                   |  |
| Full interval: 312 (hours) |  |  |  |           |  |             |  |                   |  |
| Calibration type:          |  |  |  |           |  |             |  |                   |  |
| Adjust type: None          |  |  |  |           |  |             |  |                   |  |

| O Calibrators                      |  |  |  | O Volumes |  | O Intervals |  | ● Validity checks |  |
|------------------------------------|--|--|--|-----------|--|-------------|--|-------------------|--|
| Blank absorbance range:            |  |  |  | —         |  |             |  |                   |  |
| Span: Blank                        |  |  |  | Blank     |  |             |  |                   |  |
| Span absorbance range:             |  |  |  | —         |  |             |  |                   |  |
| Expected cal factor: 0.00          |  |  |  |           |  |             |  |                   |  |
| Expected cal factor tolerance %: 0 |  |  |  |           |  |             |  |                   |  |

| Configure assay parameters — SmartWash                        |                                   |                                            |                               |                                      |  |
|---------------------------------------------------------------|-----------------------------------|--------------------------------------------|-------------------------------|--------------------------------------|--|
| <input type="radio"/> General                                 | <input type="radio"/> Calibration | <input checked="" type="radio"/> SmartWash | <input type="radio"/> Results | <input type="radio"/> Interpretation |  |
| Assay: EcstSQ                                                 |                                   |                                            |                               |                                      |  |
| COMPONENT                                                     | REAGENT / ASSAY                   | WASH                                       | Volume                        | Replicates                           |  |
| Cuvette                                                       | Trig**                            | 10% Detergent B                            | 345                           |                                      |  |
| **Not required for ARCHITECT Software version 7.00 and above. |                                   |                                            |                               |                                      |  |

| Configure assay parameters — Results |                                   |                                 |                                          |                                      |  |
|--------------------------------------|-----------------------------------|---------------------------------|------------------------------------------|--------------------------------------|--|
| <input type="radio"/> General        | <input type="radio"/> Calibration | <input type="radio"/> SmartWash | <input checked="" type="radio"/> Results | <input type="radio"/> Interpretation |  |
| Assay: EcstSQ                        |                                   |                                 | Assay number: 2861                       |                                      |  |
| Default dilution range:              |                                   |                                 | Result units: ng/mL                      |                                      |  |
| Low-Linearity: †† 22.0               |                                   |                                 |                                          |                                      |  |
| High-Linearity: 1000.0               |                                   |                                 |                                          |                                      |  |

| Configure result units          |  |
|---------------------------------|--|
| Assay: EcstSQ                   |  |
| Version: †                      |  |
| Result units: ng/mL             |  |
| Decimal places: 1 [Range 0 – 4] |  |
| Correlation factor: 1.0000      |  |
| Intercept: 0.0000               |  |

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**Microgenics Corporation**
**EC/REP** Microgenics GmbH
Abbott Laboratories Inc.