State of Minnesota

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SECRETARY OF STATE

APOSTILLE

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1. Country: United States of America

This Public Document: Filtek Bulk Fill Posterior Restorative

2. has been signed by Angela Maxcine Draper

3. acting in the capacity of Notary Public, State of Minnesota

4. bears the seal/stamp of Angela Maxcine Draper, Notary Public, State of Minnesota

CERTIFIED

5. at: St. Paul, Minnesota

6. Date: September 10, 2015

7. by: Secretary of State, State of Minnesota

8. File No: 8416323-1

9. Seal/Stamp:

10. Signature:

Steve Simon

Steve Simon
Secretary of State

Statement

To Whom It May Concern:

Subject: Filtek™ Bulk Fill Posterior Restorative

In my capacity as Senior Regulatory Affairs Specialist, 3M ESPE Dental Products, St. Paul, Minnesota, USA, I am providing the attached Instructions for Use, dated 2015-04, for subject product.

Sincerely,

Scott Erickson, RAC
Senior Regulatory Affairs Specialist
3M ESPE Dental Products


10 Sep 2015
Date

Attachment

STATE OF MINNESOTA)
COUNTY OF RAMSEY)^{ss}

On September 10, 2015, before me, a Notary Public within and for said county, personally appeared Scott Erickson known by me to be the person described in and who executed the foregoing instrument and acknowledged that he executed the same as his own free act and deed.

Angela Maxcine Draper, Notary Public
My Commission Expires, January 31, 2019

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ENGLISH

General Information

3M™ ESPE™ Filtek™ Bulk Fill Posterior Restorative material, is a visible-light activated, restorative composite optimized to create posterior restorations simpler and faster. This bulk fill material provides excellent strength and low wear for durability. The shades are semi-translucent and low stress curing, enabling up to a 5 mm depth-of-cure (refer to curing recommendations in table below). With excellent polish retention, Filtek Bulk Fill Posterior Restorative is also useful for anterior restorations that call for a semi-translucent shade. All shades are radiopaque. Filtek Bulk Fill Posterior Restorative is offered in A1, A2, A3, B1, and C2 shades. The fillers are a combination of a non-agglomerated/non-aggregated 20 nm silica filler, a non-agglomerated/non-aggregated 4 to 11 nm zirconia filler, an aggregated zirconia/silica cluster filler (comprised of 20 nm silica and 4 to 11 nm zirconia particles), and a ytterbium trifluoride filler consisting of agglomerate 100 nm particles. The inorganic filler loading is about 76.5% by weight (58.4% by volume). Filtek Bulk Fill Posterior Restorative contains ERGP-DMA, diurethane-DMA, and 1, 12-dodecane-DMA. Filtek Bulk Fill Posterior Restorative is applied to the tooth following use of a methacrylate-based dental adhesive, such as manufactured by 3M ESPE, which permanently bonds the restoration to the tooth structure. Filtek Bulk Fill Posterior Restorative is packaged in traditional syringes and single-dose capsules.

Indications:

Filtek Bulk Fill Posterior Restorative is indicated for use in:

- Direct anterior and posterior restorations (including occlusal surfaces)
- Base/liner under direct restorations
- Core build-ups
- Splinting
- Indirect restorations including inlays, onlays and veneers
- Restorations of deciduous teeth
- Extended fissure sealing in molars and premolars
- Repair of defects in porcelain restorations, enamel, and temporaries

Precautionary Information for Patients

This product contains substances that may cause an allergic reaction by skin contact in certain individuals. Avoid use of this product in patients with known acrylate allergies. If prolonged contact with oral soft tissue occurs, flush with large amounts of water. If allergic reaction occurs, seek medical attention as needed, remove the product if necessary and discontinue future use of the product.

Precautionary Information for Dental Personnel

This product contains substances that may cause an allergic reaction by skin contact in certain individuals. To reduce the risk of allergic response, minimize exposure to these materials. In particular, avoid exposure to uncured product. If skin contact occurs, wash skin with soap and water. Use of protective gloves and a no-touch technique is recommended. Acrylates may penetrate commonly used gloves. If product contacts glove, remove and discard glove, wash hands immediately with soap and water and then re-glove. If allergic reaction occurs, seek medical attention as needed. 3M ESPE MSDS information can be obtained from www.3MESPE.com or contact your local subsidiary.

Instructions for Use

Preparation

- Prophy:** Teeth should be cleaned with pumice and water to remove surface stains.
- Shade Selection:** Prior to isolation of tooth, select the appropriate shade(s) of Filtek Bulk Fill Posterior Restorative using a standard VITAPAN® classical shade guide.
Note: As Filtek Bulk Fill Posterior Restorative is semi-translucent, the location of the restoration, underlying tooth color or adjacent restorations may influence the final appearance of the restoration.
- Isolation:** A rubber dam is the preferred method of isolation. Cotton rolls and an evacuator can also be used.

Directions

Direct Restorations

4. Cavity Preparation:

- 4.1 Anterior restorations:** Use conventional cavity preparations for all Class III, IV and V restorations.
- 4.2 Posterior restorations:** Prepare the cavity. Line and point angles should be rounded. No residual amalgam or other base material should be left in the internal form of the

preparation that would interfere with light transmission and therefore, the hardening of the restorative material.

5. Placement of Matrix:

- 5.1 Anterior restorations:** Mylar strips and crown forms may be used to minimize the amount of material used.
- 5.2 Posterior restorations:** Place a thin dead-soft metal, or a pre-contoured-mylar or a pre-contoured-metal matrix band and insert wedges firmly. Burnish the matrix band to establish proximal contour and contact area. Adapt the band to seal the gingival area to avoid overhangs.

Note: The matrix may be placed following the enamel etching and adhesive application steps if preferred.

- 6. Pulp protection:** If a pulp exposure has occurred and the situation warrants a direct pulp capping procedure, use a minimum amount of calcium hydroxide on the exposure followed by an application of 3M™ ESPE™ Vitrebond™ or Vitrebond™ Plus Light Cure Glass Ionomer. Vitrebond or Vitrebond Plus liner/bases may also be used to line areas of deep cavity excavation.

- 7. Adhesive System:** To bond Filtek Bulk Fill Posterior Restorative to tooth structure, use of a 3M™ ESPE™ dental adhesive system (for example 3M™ ESPE™ Single Bond Universal) is recommended. Refer to adhesive system product instructions for full instructions and precautions for the products. After curing the adhesive, continue to maintain isolation from blood, saliva and other fluids and proceed immediately to placement of Filtek Bulk Fill Posterior Restorative.

Note: Follow the adhesive system instructions for use for recommended silane treatment during repair of ceramic restorations, followed by the adhesive application.

8. Delivery:

Follow the directions corresponding to the dispensing system chosen.

Dispensing the Composite:

- 8.1 Syringe:** Dispense the necessary amount of restorative material from the syringe onto the mix pad by turning the handle slowly in a clockwise manner. To prevent oozing of the restorative when dispensing is completed, turn the handle counterclockwise a half turn to stop paste flow. Immediately replace syringe cap. If not used immediately, the dispensed material should be protected from light.
- 8.2 Single-Dose Capsule:** Insert capsule into 3M™ ESPE™ Restorative Dispenser. Refer to separate restorative dispenser instructions for full instructions and precautions. Extrude restorative directly into cavity.

9. Placement:

- 9.1 Avoid intense light in the working field.** Exposure of paste to intense light may cause premature polymerization.
- 9.2 Capsule:** Start dispensing in the deepest portion of the preparation, holding the tip close to the preparation surface. Withdraw the capsule tip slowly as the cavity is filled, and avoid lifting the tip out of dispensed material while dispensing, to reduce voids. When dispensing has been completed, drag the capsule tip against the cavity wall while withdrawing from the operative field. For proximal areas, hold the tip against the matrix to aid material flow into the proximal box.
- 9.3 Slightly overfill the cavity** to permit extension of composite beyond cavity margins. Contour and shape with appropriate composite instruments.

- 10. Curing:** This product is intended to be cured by exposure to a halogen or LED light with a minimum intensity of 550 mW/cm² in the 400-500 nm range. Cure each increment by exposing its entire surface to a high intensity visible light source, such as a 3M ESPE curing light. Hold the light guide tip as close to the restorative as possible during light exposure. Use light cure chart to determine appropriate cure times and conditions for all shades.

Caries Classification	Increment Depth	All halogen lights (with output of 550-1000 mW/cm ²)	3M™ ESPE™ LED lights (with output 1000-2000 mW/cm ²)
Classes I, III, IV and V	4 mm	40 sec	20 sec
Class II	5 mm	20 sec occlusal, 20 sec buccal, 20 sec lingual	10 sec occlusal, 10 sec buccal, 10 sec lingual

Note: For class II restorations, remove the matrix band prior to the buccal and lingual curing steps.

- 11. Contouring:** Contour restoration surfaces with fine finishing diamonds, burs or stones. Contour proximal surfaces with Sof-Lex™ Finishing Strips, manufactured for 3M ESPE.
- 12. Adjust Occlusion:** Check occlusion with a thin articulating paper. Examine centric and lateral excursion contacts. Carefully adjust occlusion by removing material with a fine polishing diamond or stone.
- 13. Finish and Polishing:** Polish with the Sof-Lex™ Finishing and Polishing System.

Indirect Procedure for Inlays, Onlays or Veneers

1. Dental Operatory Procedure

- 1.1 Shade selection:** Choose the appropriate shade(s) of Filtek Bulk Fill Posterior Restorative prior to isolation.
- 1.2 Preparation:** Prepare the tooth.

- 1.3 Impressioning: After preparation is complete, make an impression of the prepared tooth by following the manufacturer's instructions of the impression material chosen. An impressioning material, such as manufactured by 3M ESPE, may be used.

2. Laboratory Procedure

- 2.1 Pour the impression of the preparation with die stone. Place pins at the preparation site at this time if a "triple tray" type of impression was used.
- 2.2 Separate the cast from the impression after 45 to 60 minutes. Place pins in die and base the cast as for a typical crown and bridge procedure. Mount or articulate the cast to its counter model on an adequate articulator.
- 2.3 If a second impression was not sent, pour a second cast using the same impression registration. This is to be used as a working cast.
- 2.4 Section out the preparation with a laboratory saw and trim away excess or, expose the margins so they can be easily worked. Mark the margins with a red pencil if needed. Add a spacer at this time if one is required.
- 2.5 Soak the die in water, then with a brush, apply a very thin coat of separating medium to the preparation, let it dry somewhat, then add another thin layer.
- 2.6 Add the first increment of composite to the floor of the preparation, stay short of the margins, and follow the cure recommendations described in the Direct Restoration section (Step 10).
- 2.7 Place and cure additional increments of composite. Allow for the last increment (incisal) to include the contact areas.
- 2.8 Place the die back into the articulated arch. Add the last increment of composite to the occlusal surface. Overfill very slightly mesially, distally, and occlusally. This will allow for the mesiodistal contacts and the proper occlusal contact when the opposing arch is brought into occlusion with the uncured increment. Light cure for only ten seconds, then remove the die to prevent adhering to adjacent surfaces. Finish the curing process following the cure times in the Direct Restoration section (Step 10).
- 2.9 With the occlusal contacts already established, begin removing the excess composite from around the points of contact. Develop the inclines and ridges as per remaining occlusal anatomy.
- 2.10 Care must be taken when removing the prosthesis from the die. Break off small amounts of the die from around the restoration, the die stone should break away cleanly from the cured restoration, until all of the restoration is recovered.
- 2.11 Using the master die, check the restoration for flash, undercuts, and fit. Adjust as necessary, and then polish as noted above in Direct Restorative steps 11-13.

3. Dental Operatory Procedure

- 3.1 Roughen the interior surfaces of the indirect restoration.
- 3.2 Clean the prosthesis in a soap solution in an ultrasonic bath and rinse thoroughly.
- 3.3 Cementation: Cement the prosthesis using a 3M™ ESPE™ resin cement system.

Storage and Use

1. This product is designed to be used at room temperature. If stored in cooler allow product to reach room temperature prior to use. Shelf life at room temperature is 36 months. Ambient temperatures routinely higher than 27°C/80°F may reduce shelf life. See outer package for expiration date.
 2. Do not expose restorative materials to elevated temperatures, or to intense light.
 3. Do not store materials in proximity to eugenol containing products.
- Disinfect this product using an intermediate level disinfection process (liquid contact) as recommended by the Centers for Disease Control and endorsed by the American Dental Association. Guidelines for Infection Control in Dental Health-Care Settings – *MMWR*, December 19, 2003:52(RR-17), Centers for Disease Control and Prevention.

Disposal

See the Material Safety Data Sheet (available at www.3MESPE.com or through your local subsidiary for disposal information).

Customer Information

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Caution: U.S. Federal Law restricts the device to sale or use on the order of a dental professional.

Warranty

3M ESPE warrants this product will be free from defects in material and manufacture. 3M MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the products for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M ESPE's sole obligation shall be repair or replacement of the 3M ESPE product.

Limitation of Liability

Except where prohibited by law, 3M ESPE will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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VITAPAN® is not a registered trademark of 3M or 3M Deutschland GmbH.
3M, ESPE, Filtek, Vitrebond and Sof-Lax are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada.
© 3M 2015. All rights reserved.