Niall Sheehy
Notary Public

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NIALL SHEEHY

NOTARY PUBLIC

6/7, Glenthworth St.,
Limerick
Ireland

Telephone +353 61 414355


Facsimile + 353 61 414738

Email nsheehy@sellors.ie

TO ALL TO WHOM THESE PRESENTS SHALL COME, I, Niall Sheehy, Notary Public, duly authorised, admitted, and sworn, practicing at 6/7 Glenthworth St., Limerick, Republic of Ireland, DO HEREBY CERTIFY that I am familiar with documents as issued from time to time by Cook Ireland Ltd and that the original document hereunto annexed namely, **Instructions for Use-Evolution Biliary Stent System Covered** upon which I have signed my name and affixed my Seal of Office is a true and faithful original and that the original appears to me to be genuine.

AND I HAVE CERTIFIED as herein at the request of Cook Ireland Ltd having its place of business at O'Halloran Road, National Technology Park, in the city of Limerick, Ireland, and for the benefit of the proper authorities in **Russian Federation** and for all to whom it shall come.

IN FAITH AND TESTIMONY whereof I have hereunto set my hand and affixed my Seal of Office.

Niall Sheehy
Notary Public

Date

NIALL SHEEHY
Notary Public
6/7 Glenthworth St., Limerick.
Commissioned For Life

APOSTILLE (Convention de La Haye du 5 octobre 1961)			
1. Country: Pays/País:		IRELAND	
This public document Le présent acte public / El presente documento público			
2. has been signed by a été signé par ha sido firmado por		Niall Sheehy	
3. acting in the capacity of agissant en qualité de quien actúa en calidad de		Notary Public	
4. bears the seal / stamp of est revêtu du sceau / timbre de y está revestido del sello / timbre de			
Certified Attesté / Certificado			
5. at à / en	Cork	6. the le / el día	16/08/2021
7. by par / por	Department of Foreign Affairs		
8. No sous no bajo el número	7742352021		
9. Seal / stamp: Sceau / timbre Sello / timbre	10. Signature: Signature: Firma: <i>Mary Kelleher</i>		
This Apostille only certifies the authenticity of the signature and the capacity of the person who has signed the public document. It does not certify the content of the document for which it was issued. To verify the issuance of this Apostille, see www.austlii.edu.au/au/other/dfat/ie			

327404

INSTRUCTIONS FOR USE

Evolution Biliary Stent System Covered

PURPOSE

These devices are intended for the palliation of malignant neoplasms in the biliary tree.

DEVICE DESCRIPTION

The Evolution Biliary Stent System is available in two versions:

- Evolution Biliary Stent System - Partially Covered (EVO-PC-8-9-6-B; EVO-PC-8-9-8-B; EVO-PC-10-11-4-B; EVO-PC-10-11-6-B; EVO-PC-10-11-8-B)
- Evolution Biliary Stent System - Full Covered (EVO-FC-8-9-6-B; EVO-FC-8-9-8-B; EVO-FC-10-11-4-B; EVO-FC-10-11-6-B; EVO-FC-10-11-8-B)

The Evolution Biliary Duct Stent System consists of the following two components:

- self-expanding metallic stents
- stent delivery system

Stent Delivery System (Introducer and Handle)

The stent delivery system comprises an introducer and a handle. The introducer consists of two co-axial catheters (i.e., an inner and outer catheter). A preloaded stent, placed between the inner and outer catheter, is released by retracting the outer catheter. Once exposed, the stent "self-expands" to conform to the surrounding anatomy. The movement of the inner catheter with respect to the outer catheter is controlled by the handle, which allows for desheathing (to deploy the stent) and resheathing (to recapture the stent during stent deployment). There are four radiopaque markers to aid in the deployment of the stent while using fluoroscopy. Three of these markers are on the inner catheter and one on the outer catheter. Inner catheter marker - Length: 3mm +/- 0.25mm. Outer catheters marker - Length: 1.02 mm (+/- 0.08 mm). The total length of the stent is indicated by radiopaque markers on the inner catheter, indicating the actual length of the stent in a collapsed state. The stent frame has a diameter of 0.180 mm. The introducer contains a zip port on the outer catheter to allow for the exit of a pre-positioned wire guide. The Cook Evolution Biliary Stent System features a single-handed deployment handle for retracting the outer catheter, which provides greater control to the user during stent deployment. The system also allows the user to recapture the stent if the user fails to place the stent in the desired location. This can be achieved by pressing a button on the side of the handle, which reverses the direction of the outer catheter movement. An internal gearing mechanism within the handle makes this possible. Once the stent is fully recaptured within the outer catheter, the stent may redeploy.



Figure. Handle. Length 298.29mm (+/- 0.5mm)

The differences in the design of the self-expanding metallic stents from the aforementioned versions are described below.

The Cook Evolution Biliary Stent System Partially Covered (EVO-PC-8-9-6-B; EVO-PC-8-9-8-B; EVO-PC-10-11-4-B; EVO-PC-10-11-6-B; EVO-PC-10-11-8-B)

Stent constructed of a single, woven, platinum-core nitinol radiopaque wire. The stent foreshortens due to its design. The stent has a flange at both ends to aid in mitigating migration after the stent has been placed in the bile duct. A silicone membrane is applied to the surface of the stent; this covers the complete body of the partially covered stent (not including the flanges at both ends). The covering is intended to reduce the risk of tissue ingrowth. To aid in visibility under fluoroscopy, the stent is manufactured with a platinum-cored nitinol wire. A lasso loop is included at the proximal end, which can be used to reposition the stent after placement or for acute removal of the stent, if desired. The proximal and distal ends of the stent are designed with rounded "crowns" to potentially reduce tissue trauma.


Figure 1: The Cook Evolution Biliary Stent System Partially Covered

The Cook Evolution Biliary Stent System - Full Covered (EVO-FC-8-9-6-B; EVO-FC-8-9-8-B; EVO-FC-10-11-4-B; EVO-FC-10-11-6-B; EVO-FC-10-11-8-B)

Stent constructed of a single, woven, platinum-core nitinol radiopaque wire. The stent foreshortens due to its design. The stent has a flange at both ends to aid in mitigating migration after the stent has been placed in the bile duct. A silicone membrane is applied to the surface of the stent; this covers the full stent in the fully covered version (except where the lasso passes through the silicone). The covering is intended to reduce the risk of tissue ingrowth. To aid in visibility under fluoroscopy, the stent is manufactured with a platinum-cored nitinol wire. A lasso loop is included at the proximal end, which can be used to reposition the stent after placement or for acute removal of the stent, if desired. The proximal and distal ends of the stent are designed with rounded "crowns" to potentially reduce tissue trauma.

Figure 2: The Cook Evolution Biliary Stent System Full Covered

DESIGN

1 Fr=0.33 mm

Catalogue number	Introduction System Diameter (Fr / mm)	Introduction System Working Length (m)	Stent length (mm)	Stent body diameter (mm)	Stent Flange Diameter (mm)
The Cook Evolution Biliary Stent System Partially Covered					
EVO-PC-8-9-6-B	8.5/ 2.8	2	60	8	9
EVO-PC-8-9-8-B	8.5/ 2.8	2	80	8	9
EVO-PC-10-11-4-B	8.5/ 2.8	2	40	10	11
EVO-PC-10-11-6-B	8.5/ 2.8	2	60	10	11
EVO-PC-10-11-8-B	8.5/ 2.8	2	80	10	11
The Cook Evolution Biliary Stent System Full Covered					
EVO-FC-8-9-6-B	8.5/ 2.8	2	60	8	9
EVO-FC-8-9-8-B	8.5/ 2.8	2	80	8	9
EVO-FC-10-11-4-B	8.5/ 2.8	2	40	10	11
EVO-FC-10-11-6-B	8.5/ 2.8	2	60	10	11
EVO-FC-10-11-8-B	8.5/ 2.8	2	80	10	11

Radial resistance at 37° C

Test Conditions: The test article was loaded into an introducer sheath and deployed. The test articles were tested using a low friction, stainless steel machine solutions radial expansion force gauge with a temperature within the tester maintained at 37 °C ± 2 °C.

Table: Radial force results at 3 mm diameter (Expansion).

Description Body diameter x Length (mm)	Expansion force per unit length (N/mm)	Pass/Fail
Evolution Biliary Stent System - Full Covered (10x40)	Range: 0.46 to 0.52 Mean: 0.49 SD: 0.02	Pass
Evolution Biliary Stent System - Full Covered (8x60)	Range: 0.41 to 0.53 Mean: 0.46 SD: 0.03	
Evolution Biliary Stent System - Partially Covered (8x60)	Range: 0.42 to 0.44 Mean: 0.43 SD: 0.01	

Table: Results of radial reinforcement 6 mm or 8 mm depending on the original stent diameter (Expansion)

Description Body diameter x Length (mm)	Expansion force per unit length (N/mm)	Pass/Fail
Evolution Biliary Stent System - Full Covered (8x60)	Range: 0.15 to 0.19 Mean: 0.17 SD: 0.01	Pass
Evolution Biliary Stent System - Full Covered (10x40)	Range: 0.16 to 0.18 Mean: 0.17 SD: 0.01	
Evolution Biliary Stent System - Partially Covered (10x40)	Range: 0.11 to 0.14 Mean: 0.12 SD: 0.01	

Axial/pivotal stent strength

Objective: The purpose of this test is to establish and compare the axial force of the devices using a 3-point bend test.

Table: Results of axial force testing

Test Article	Combined Loading/ Unloading Slope (N/mm)	
	Mean	SD
		Min Max

Evolution Biliary Stent (Fully Covered) 8 x 60 mm	0.25	SD = 0.01 Min = 0.23 Max = 0.27
Evolution Biliary Stent (Partially Covered) 8 x 60 mm	0.25	SD = 0.01 Min = 0.24 Max = 0.26
Evolution Biliary Stent (Fully Covered) 10 x 40 mm	0.23	SD = 0.01 Min = 0.22 Max = 0.24
Evolution Biliary Stent (Partially Covered) 10 x 40 mm	0.22	SD = 0.01 Min = 0.20 Max = 0.22

Stent deployment strength

Acceptance criteria: The stent deployment force must be below 63N when deployed from a full length introducer into a water bath heated to 37°C ±2°C in a simulated use model.

Result:

Stent Deployment Force:

The peak force range: 31.70N to 61.95N

Acceptance criteria: The stent deployment force must be below 18N when deployed from a 'shortie' introducer into a water bath heated to 37°C ±2°C.

Result:

Stent Deployment Force:

The peak force range: 11.25N to 17.35N

Connection strength

Strength Testing- Stent Grasping Feature.

Table : Acceptance criteria and results: strength testing – stent grasping feature

Test	Acceptance Criteria	Results	Pass/ Fail
Stent Grasping Feature - <i>Tensile test</i> .	The Stent Grasping Feature must withstand a tensile force of >20.4N when using a rat-tooth forceps at a tensile speed of 330mm/min.	All test articles met the acceptance criteria. Peak Force: Mean:63.216 N SD:3.968 N Min: 57.029 N Max: 71.130 N	Pass

Strength Testing – Inner Catheter Assembly

Table: Acceptance criteria and results: strength testing – inner catheter assembly.

Test	Acceptance Criteria	Results	Pass/ Fail
Yellow Marker (Stent Carrier) to Retrieval Loop Pusher – <i>Compression test.</i>	The Stent Carrier to Retrieval Loop Pusher joint must withstand a compressive force of >23.4N at 330mm/min crosshead speed (compressive force to be applied to Yellow Marker).	All test articles withstood a force of >23.4N. The minimum force tested to was 23.504 N.	Pass
RO Markers - <i>Tensile test.</i>	The RO marker bands must remain attached to Introduction System for a tensile force of >7N when a 330mm/min tensile force is applied.	All test articles met the acceptance criteria. Peak Force: Mean:26.591N SD:6.862 N Min: 18.697 N Max: 42.969 N	Pass
Retrieval Loop Suture to Retrieval Loop Shaft - <i>Tensile test.</i>	The tensile strength of the Retrieval Loop wire/suture to Retrieval Loop Pusher joint must be >23.4N at 330mm/min cross-head speed.	All test articles withstood a force of >23.4N. The minimum force tested to was 23.516 N.	Pass
Tip (Stent Carrier) to Retrieval Loop Shaft - <i>Tensile test.</i>	The Tip must remain attached to the Stent Carrier when a tensile force of >15N is applied at 330mm/min cross head speed.	All test articles met the acceptance criteria. Peak Force:	Pass
	The Stent Carrier must remain attached to the Retrieval Loop Pusher when a tensile force of >15N is applied at 330mm/min cross head speed.	Mean:32.576 N SD: 6.167 N Min: 22.878 N Max: 46.808 N	Pass
FLLA Socket to Retrieval Loop Shaft - <i>Tensile test.</i>	The strength of the Retrieval Loop shaft to FLLA Socket joint must be >81.9N at 330mm/min cross-head speed.	All test articles met the acceptance criteria. Peak Force: Mean:141.207 N SD:6.333 N Min: 122.317 N Max: 152.304 N	Pass

Strength Testing - Outer Sheath Assembly.**Table:** Acceptance criteria and results: strength testing – outer sheath assembly

Test	Acceptance Criteria	Results	Pass/ Fail
Zip Port (Proximal Flexor	The Distal Flexor to Proximal Flexor joint (Zip Port) must achieve a	All test articles withstood a force of >23.4N.	Pass

Test	Acceptance Criteria	Results	Pass/ Fail
to Distal Flexor joint) - <i>Compression test.</i>	compressive force strength of >23.4N at 330mm/min cross head speed.	The minimum force tested to was 23.503N.	
Zip Port (Proximal Flexor to Distal Flexor joint) - <i>Tensile test.</i>	The Distal Flexor to Proximal Flexor joint (Zip Port) must achieve a tensile strength of >23.4N at 330mm/min cross head speed.	All test articles met the acceptance criteria. Peak Force: Mean:45.827 N SD:2.055 N Min: 42.894 N Max: 48.580 N	Pass
Shuttle to Proximal Flexor - <i>Compression test.</i>	The Proximal Flexor to Handle Shuttle joint must achieve a compressive force strength of >81.9N at 330mm/min cross head speed.	All test articles withstood a force of >81.9N. The minimum force tested to was 82.002N.	Pass
Shuttle to Proximal Flexor - <i>Tensile test.</i>	The Proximal Flexor to Handle Shuttle joint must achieve a tensile strength of >81.9N at 330mm/min cross-head speed.	All test articles met the acceptance criteria. Peak Force Mean:108.663 N SD:0.939 N Min: 106.192 N Max: 110.007 N	Pass
Shuttle to Belt - <i>Tensile test.</i>	The Shuttle to Belt joint must achieve a tensile strength of >81.9N at 330mm/min cross-head speed.	All test articles met the acceptance criteria. Peak Force: Mean:446.834 N SD:100.183 N Min: 228.860 N Max: 587.805 N	Pass

Strength Testing – Lock Wire Assembly.

Table: Acceptance criteria and results: strength testing – lock wire assembly.

Test	Acceptance Criteria	Results	Pass/ Fail
Lock Wire Assembly (Stent Lock Wire to Lock Wire Anchor Assembly) - <i>Tensile test.</i>	The Lock Wire Anchor Assembly to Lock Wire joint force must be >35N at 330mm/min cross head speed.	All test articles met the acceptance criteria. Peak Force: Mean:67.400 N SD:2.499 N Min: 64.131 N Max: 73.755 N	Pass

All results indicated apply to all versions Evolution Biliary Stent System - Partially Covered and Evolution Biliary Stent System - Full Covered.

Silicone Coating Details of Evolution Biliary Stent System - Partially Covered:

Catalogue number	The minimum diameter of the silicone coating, mm	Maximum diameter of silicone coating, mm	Silicone coating thickness, mm
EVO-PC-8-9-6-B	45	55	0.06 (+0.03/-0.02 mm)
EVO-PC-8-9-8-B	65	75	0.06 (+0.03/-0.02 mm)
EVO-PC-10-11-4-B	25	35	0.06 (+0.03/-0.02 mm)
EVO-PC-10-11-6-B	45	55	0.06 (+0.03/-0.02 mm)
EVO-PC-10-11-8-B	65	75	0.06 (+0.03/-0.02 mm)

Silicone coating thickness of Evolution Biliary Stent System - Full Covered - 0.06 mm (+0.03/-0.02 mm)

For stents covered with silicone, the availability of information on the area of the stents is not provided for by the production.

Temperature conditions for transportation:

Products are stable at a temperature of transportation from -50 to +50 C.

INSTRUCTIONS FOR USE

Evolution Biliary Stent System - Partially Covered (EVO-PC-8-9-6-B; EVO-PC-8-9-8-B; EVO-PC-10-11-4-B; EVO-PC-10-11-6-B; EVO-PC-10-11-8-B)

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

IMPORTANT INFORMATION

Please review prior to use.

Warning: The safety and effectiveness of this device for use in the vascular system has not been established.

INTENDED USE

This device is used in palliation of malignant neoplasms in the biliary tree.

STENT DESCRIPTION

This flexible, self-expanding biliary metal stent is constructed of a single, woven, platinum core nitinol radiopaque wire. The partially covered stent has the potential to reduce tumor ingrowth through the stent. The stent has a flange at both ends to aid in preventing migration after the stent has been placed in the bile duct. There is a yellow grasping loop at the end of the stent, the purpose is to reposition distally/remove the stent during initial stent placement procedure in the event of incorrect placement. The stent

foreshortens due to its design. The total length of the stent is indicated by radiopaque markers on the inner catheter (**Fig. D, 2 & 4**), indicating the actual length of the stent in a collapsed state. The proximal and distal ends of the stent are designed with numerous rounded “crowns” to potentially reduce tissue trauma.

DELIVERY SYSTEM DESCRIPTION

The stent is mounted on an inner catheter, which accepts a .035-inch wire guide, and is constrained by an outer catheter. A pistol-grip delivery handle allows controlled stent deployment or recapture. The position of the red marker on top of the handle in relation to the “point-of-no-return” position on the handle label indicates the point after which it is not possible to recapture the deployed stent. There are four radiopaque markers to aid in the deployment of the stent while using fluoroscopy, and a yellow endoscopic marker for transpapillary placement.

CONTRAINDICATIONS

Those specific to ERCP and any procedure to be performed in conjunction with stent placement.

Additional contraindications include, but are not limited to: inability to pass the wire guide or stent through the obstructed area, biliary duct strictures of benign etiology, biliary obstruction preventing endoscopic cholangiography, concurrent perforated bile duct, those patients for whom endoscopic procedures are contraindicated, patients with coagulopathy, concurrent bile duct stones, very small intrahepatic ducts and any use other than those specifically outlined under Intended Use.

POTENTIAL ADVERSE EVENTS

Potential adverse events associated with ERCP include, but are not limited to: allergic reaction to contrast or medication, aspiration, cardiac arrhythmia or arrest, cholangitis, cholecystitis, cholestasis, hemorrhage, hypotension, infection, pancreatitis, perforation, respiratory depression or arrest, sepsis.

Additional adverse events that can occur in conjunction with biliary stent placement include, but are not limited to: allergic reaction to nickel, bile duct ulceration, death (other than due to normal disease progression), fever, inflammation, ingrowth due to tumor or excessive hyperplastic tissue, nausea, obstruction of the pancreatic duct, pain/discomfort, perforation, recurrent obstructive jaundice, stent migration, stent misplacement, stent occlusion, trauma to the biliary tract or duodenum, tumor overgrowth, vomiting.

PRECAUTIONS

Refer to the package label for the minimum channel size required for this device.

A complete diagnostic evaluation should be performed prior to placement to measure the stricture length and determine the proper stent length. The stent length chosen should allow for additional length on either side of the stricture.

If wire guide or stent cannot advance through obstructed area, do not attempt to place stent.

Stent should be placed using fluoroscopic and endoscopic monitoring.

The stent should only be placed with the Cook delivery system, which is provided with each stent.

This stent is intended for palliative treatment only. Alternate methods of therapy should be investigated prior to placement.

Long-term patency with this stent has not been established. Periodic evaluation of the stent is advised. Assessment must be made to determine the necessity of sphincterotomy or balloon dilation prior to stent placement. In the event sphincterotomy or balloon dilation is required, all appropriate cautions, warnings and contraindications must be observed.

The stent cannot be recaptured after the deployment threshold has been passed. The position of the red marker on the top of the handle in relation to point-of-no-return position on the handle label indicates when threshold has been passed.

Do not attempt to reload a deployed stent.

Do not use this device for any other purpose than the stated intended use.

Not compatible with THSF wire guide.

Use of this device is restricted to a trained healthcare professional.

Passing a second stent delivery system through a just deployed stent is not recommended and could cause the stent to dislodge. To reduce the risk of stent dislodgement ensure that the delivery system is completely recaptured following stent deployment.

WARNINGS

The safety and effectiveness of this device for use in the vascular system has not been established.

Placement of a partially covered stent across a branch duct or major bifurcation may result in obstruction.

This stent is not intended to be removed after initial stent placement and is considered a permanent implant. Attempts to remove the stent after placement may cause damage to the surrounding mucosa.

This device is not intended to be deployed through the wall of a previously placed or existing metal stent.

Doing so could result in difficulty or inability to remove introducer.

The stent contains nickel, which may cause an allergic reaction in individuals with nickel sensitivity.

This stent system is designed for single use only. Attempts to reprocess, resterilize and/or reuse may lead to device failure and/or transmission of disease.

After stent placement, additional methods of treatment such as chemotherapy and irradiation may increase the risk of stent migration due to tumor shrinkage, stent erosion, and/or mucosal bleeding.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Testing of overlapping stents has not been completed and is not recommended.

MRI INFORMATION

This symbol means the device is MR conditional

Non-clinical testing has demonstrated that the Evolution Biliary Stent is MR Conditional according to ASTM F2503. A patient with this stent can be scanned safely anytime after placement under the following conditions.

Static Magnetic Field

- Static magnetic field of 3 Tesla or less
- Highest spatial magnetic gradient field of 1600 Gauss/cm or less

Non-clinical evaluation was conducted in a MR system (Excite, General Electric Healthcare) with a maximum spatial magnetic gradient field of 1600 Gauss/ cm as measured with a gaussmeter in the position of the static magnetic field pertinent to the patient (i.e., outside of scanner covering, accessible to a patient or individual).

MRI-Related Heating

1.5 and 3.0 Tesla Systems: It is recommended to scan in normal operation mode (whole body averaged specific absorption rate (SAR) = 2.0W/ kg) ("Normal Operating Mode" is defined as the mode of operation of the MR system in which none of the outputs have a value that cause physiological stress to the patient) for 15 minutes of scanning (i.e., per scanning sequence).

1.5 Tesla Temperature Rise

In non-clinical testing, the Evolution Biliary Stent produced a maximum temperature rise of 2.1°C (for a single stent) during 15 minutes of MR imaging (i.e., for one scanning sequence) performed in a MR 1.5 Tesla System (Siemens Medical Solutions, Malvern, PA, Software Numaris/4) at an MR system reported

whole-body-averaged SAR of 2.9 W/kg (associated with a calorimetry measured whole body averaged value of 2.1 W/kg).

3.0 Tesla Temperature Rise

In non-clinical testing, the Evolution Biliary Stent produced a maximum temperature rise of 3.1°C (for a single stent) during 15 minutes of MR imaging (i.e., for one scanning sequence) performed in a MR 3 Tesla System (Excite, GE Electric Healthcare, Milwaukee, WI, Software 14X.M5) at an MR system reported whole-body-averaged SAR of 2.9 W/kg (associated with a calorimetry measured whole body averaged value of 2.7 W/kg).

Image Artifacts

MR image quality may be compromised if the area of interest is within the lumen or within approximately 7 mm of the position of the Evolution Biliary Stent, as found during non-clinical testing using the sequence: T1-weighted, spin echo and gradient echo pulse sequence in a 3.0 Tesla MR system (Excite, General Electric Healthcare, Milwaukee, WI). Therefore, it may be necessary to optimize MR imaging parameters for the presence of this metallic stent.

It is recommended that patients register the conditions under which the implant can safely be scanned with the MedAlert Foundation (www.medicalert.org) or an equivalent organization.

EQUIPMENT REQUIRED

Stent and Delivery System

.035 inch wire guide

Fusion Wire Lock

PROCEDURE PREPARATION

With duodenoscope in place, introduce a wire guide, floppy tip first and advance until it is fluoroscopically visualized in position through stricture.

INSTRUCTIONS FOR USE

1. Remove stent delivery system from package and backload device over a prepositioned wire guide, ensuring that the wire guide exits catheter at Zip port (**Fig. A1, A2**).
2. Unlock wire guide from Fusion Wire Guide Locking Device.
3. Introduce device in short increments into accessory channel of duodenoscope until Zip port is inside of the accessory channel, then relock wire guide. Continue advancing device in short increments.
4. With elevator open, advance device until endoscopically visualized exiting duodenoscope.
5. For transpapillary placement reference the yellow marker for accurate stent positioning. **Note:** Stents bridging papilla should extend beyond papilla and into duodenum approximately 0.5 cm after deployment.
6. Under fluoroscopic guidance, with elevator open, continue to advance the device in short increments until the stent is fluoroscopically visualized through the stricture. Fluoroscopically visualize radiopaque markers on inner catheter at either end of stent and position the inner radiopaque markers (**Fig. D,2 and D,4**) a minimum of 1 cm beyond the stricture so that it bridges stricture completely.
7. Confirm desired stent position fluoroscopically. To deploy stent, remove red safety guard from handle, (**Fig. B**) then squeeze the trigger. **Note:** Each full trigger squeeze will deploy stent by an equal amount.
8. If stent repositioning is required during deployment, it is possible to recapture stent. **Note:** It is not possible to recapture stent after passing point-of-no-return, indicated when red marker on top of handle has passed the point-of-no-return position on the handle label (**Fig. C**).

To fluoroscopically monitor the point-of-no-return (**Fig. D, 5**), the stent can be recaptured until the point where the outer catheter radiopaque marker (**Fig. D, 1**) is aligned with the second inner catheter radiopaque marker (**Fig. D, 3**).

Note: It is **not** possible to recapture stent after the outer radiopaque marker (**Fig. D, 1**) passes the second inner radiopaque marker (**Fig. D, 3**).

Note: The stent can be recaptured a maximum of three times.

9. To reposition the stent, the stent must first be recaptured and the elevator must be open. **Note:** Do not push forward on the delivery system with stent partially deployed. Push direction button on side of delivery handle to opposite side (**Fig. E1**). **Note:** Hold thumb on button when squeezing trigger for first time to recapture. Continue squeezing trigger as required to recapture stent by desired amount.

10. To resume deployment, push button to opposite side again and hold button for first stroke while squeezing trigger (**Fig. E2**).




11. When stent point-of-no-return has been passed, disconnect luer lock fitting to remove safety wire completely from delivery handle (**Fig. F**).






12. Continue deploying stent by squeezing trigger.

13. After deployment, fluoroscopically confirm full stent expansion. While maintaining wire guide position, push direction button on side of delivery handle to opposite side (**Fig. E1**). Squeeze the trigger to completely recapture the delivery system. **Unlock the wire guide from Fusion Wire Guide locking device.** The device can be safely removed with the elevator of the duodenoscope fully down.

To reposition distally/remove stent during initial stent placement procedure, grasp the yellow grasping loop at end of stent with forceps (**Fig. G**) and reposition distally/remove as required.

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

INSTRUCTIONS FOR USE

Evolution Biliary Stent System - Full Covered (EVO-FC-8-9-6-B; EVO-FC-8-9-8-B; EVO-FC-10-11-4-B; EVO-FC-10-11-6-B; EVO-FC-10-11-8-B)

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

IMPORTANT INFORMATION

Please review prior to use.

Warning: The safety and effectiveness of this device for use in the vascular system has not been established.

INTENDED USE

This device is used in palliation of malignant neoplasms in the biliary tree.

STENT DESCRIPTION

This flexible, self-expanding biliary metal stent is constructed of a single, woven, platinum core nitinol radiopaque wire. The fully covered stent has the potential to reduce tumor ingrowth through the stent. The stent has a flange at both ends to aid in preventing migration after the stent has been placed in the bile duct. There is a yellow grasping loop at the end of the stent, the purpose is to reposition distally/remove the stent during initial stent placement procedure in the event of incorrect placement. The stent foreshortens due to its design. The total length of the stent is indicated by radiopaque markers on the inner catheter (**Fig. D,2 & 4**), indicating the actual length of the stent in a collapsed state. The proximal and distal ends of the stent are designed with numerous rounded "crowns" to potentially reduce tissue trauma.

DELIVERY SYSTEM DESCRIPTION

The stent is mounted on an inner catheter, which accepts a .035-inch wire guide, and is constrained by an outer catheter. A pistol-grip delivery handle allows controlled stent deployment or recapture, and a yellow endoscopic marker for transpapillary placement. The position of the red marker on top of the handle in relation to the "point of no-return" position on the handle label indicates the point after which it is not possible to recapture the deployed stent. There are four radiopaque markers to aid in the deployment of the stent while using fluoroscopy.

CONTRAINDICATIONS

Those specific to ERCP and any procedure to be performed in conjunction with stent placement.

Additional contraindications include, but are not limited to: inability to pass the wire guide or stent through the obstructed area, biliary duct strictures of benign etiology, biliary obstruction preventing endoscopic cholangiography, concurrent perforated bile duct, those patients for whom endoscopic procedures are contraindicated, patients with coagulopathy, concurrent bile duct stones, very small intrahepatic ducts and any use other than those specifically outlined under Intended Use.

POTENTIAL ADVERSE EVENTS

Potential adverse events associated with ERCP include, but are not limited to: allergic reaction to contrast or medication, aspiration, cardiac arrhythmia or arrest, cholangitis, cholecystitis, cholestasis, hemorrhage, hypotension, infection, pancreatitis, perforation, respiratory depression or arrest, sepsis.

Additional adverse events that can occur in conjunction with biliary stent placement include, but are not limited to: allergic reaction to nickel, bile duct ulceration, death (other than due to normal disease progression), fever, inflammation, ingrowth due to tumor or excessive hyperplastic tissue, nausea, obstruction of the pancreatic duct, pain/discomfort, perforation, recurrent obstructive jaundice, stent migration, stent misplacement, stent occlusion, trauma to the biliary tract or duodenum, tumor overgrowth, vomiting.

PRECAUTIONS

Refer to the package label for the minimum channel size required for this device.

A complete diagnostic evaluation should be performed prior to placement to measure the stricture length and determine the proper stent length. The stent length chosen should allow for additional length on either side of the stricture.

If wire guide or stent cannot advance through obstructed area, do not attempt to place stent.

Stent should be placed using fluoroscopic and endoscopic monitoring.

The stent should only be placed with the Cook delivery system, which is provided with each stent.

This stent is intended for palliative treatment only. Alternate methods of therapy should be investigated prior to placement.

Long-term patency with this stent has not been established. Periodic evaluation of the stent is advised. Assessment must be made to determine the necessity of sphincterotomy or balloon dilation prior to stent

placement. In the event sphincterotomy or balloon dilation is required, all appropriate cautions, warnings and contraindications must be observed.

The stent cannot be recaptured after the deployment threshold has been passed. The position of the red marker on the top of the handle in relation to point-of-no-return position on the handle label indicates when threshold has been passed.

Do not attempt to reload a deployed stent.

Do not use this device for any other purpose than the stated intended use.

Not compatible with THSF wire guide.

Use of this device is restricted to a trained healthcare professional.

Passing a second stent delivery system through a just deployed stent is not recommended and could cause the stent to dislodge.

WARNINGS

Placement of a fully covered stent across a branch duct or major bifurcation may result in obstruction.

This stent is not intended to be removed after initial stent placement procedure and is considered a permanent implant. Attempts to remove the stent after placement may cause damage to the surrounding mucosa.

Testing of overlapping stents has not been completed and is not recommended.

This device is not intended to be deployed through the wall of a previously placed or existing metal stent. Doing so could result in difficulty or inability to remove introducer.

The stent contains nickel, which may cause an allergic reaction in individuals with nickel sensitivity.

This stent system is designed for single use only. Attempts to reprocess, resterilize and/or reuse may lead to device failure and/or transmission of disease.

After stent placement, additional methods of treatment such as chemotherapy and irradiation may increase the risk of stent migration due to tumor shrinkage, stent erosion, and/or mucosal bleeding.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

The safety and effectiveness of this device for use in the vascular system has not been established.

MRI INFORMATION

This symbol means the device is MR conditional.

Non-clinical testing has demonstrated that the Evolution Biliary Stent is MR Conditional according to ASTM F2503. A patient with this stent can be scanned safely anytime after placement under the following conditions.

Static Magnetic Field

- Static magnetic field of 3 Tesla or less
- Highest spatial magnetic gradient field of 1600 Gauss/cm or less

Non-clinical evaluation was conducted in a MR system (Excite, General Electric Healthcare) with a maximum spatial magnetic gradient field of 1600 Gauss/cm as measured with a gaussmeter in the position of the static magnetic field pertinent to the patient (i.e., outside of scanner covering, accessible to a patient or individual).

MRI-Related Heating

1.5 and 3.0 Tesla Systems: It is recommended to scan in normal operation mode (whole body averaged specific absorption rate (SAR) = 2.0W/kg) ("Normal Operating Mode" is defined as the mode of operation of the MR system in which none of the outputs have a value that cause physiological stress to the patient) for 15 minutes of scanning (i.e., per scanning sequence).

1.5 Tesla Temperature Rise

In non-clinical testing, the Evolution Biliary Stent produced a maximum temperature rise of 2.1°C (for a single stent) during 15 minutes of MR imaging (i.e., for one scanning sequence) performed in a MR 1.5 Tesla System (Siemens Medical Solutions, Malvern, PA, Software Numaris/4) at an MR system reported whole-body-averaged SAR of 2.9 W/kg (associated with a calorimetry measured whole body averaged value of 2.1 W/kg).

3.0 Tesla Temperature Rise

In non-clinical testing, the Evolution Biliary Stent produced a maximum temperature rise of 3.1°C (for a single stent) during 15 minutes of MR imaging (i.e., for one scanning sequence) performed in a MR 3 Tesla System (Excite, GE Electric Healthcare, Milwaukee, WI, Software 14X.M5) at an MR system reported whole-body-averaged SAR of 2.9 W/kg (associated with a calorimetry measured whole body averaged value of 2.7 W/kg).

Image Artifacts

MR image quality may be compromised if the area of interest is within the lumen or within approximately 7 mm of the position of the Evolution Biliary Stent, as found during non-clinical testing using the sequence: T1-weighted, spin echo and gradient echo pulse sequence in a 3.0 Tesla MR system (Excite, General Electric Healthcare, Milwaukee, WI). Therefore, it may be necessary to optimize MR imaging parameters for the presence of this metallic stent.

It is recommended that patients register the conditions under which the implant can safely be scanned with the MedicAlert Foundation (www.medicalert.org) or an equivalent organization.

EQUIPMENT REQUIRED

Stent and Delivery System.

035 inch wire guide

Fusion Wire Lock

PROCEDURE PREPARATION

With duodenoscope in place, introduce a wire guide, floppy tip first and advance until it is fluoroscopically visualized in position through stricture.

INSTRUCTIONS FOR USE

1. Remove stent delivery system from package and backload device over a prepositioned wire guide, ensuring that the wire guide exits catheter at Zip port (**Fig. A1, A2**).
2. Unlock wire guide from Fusion Wire Guide Locking Device.
3. Introduce device in short increments into accessory channel of duodenoscope until Zip port is inside of the accessory channel, then relock wire guide. Continue advancing device in short increments.
4. With elevator open, advance device until endoscopically visualized exiting duodenoscope.
5. For transpapillary placement reference the yellow marker for accurate stent positioning. **Note:** Stents bridging papilla should extend beyond papilla and into duodenum approximately 0.5 cm after deployment.
6. Under fluoroscopic guidance, with elevator open, continue to advance the device in short increments until the stent is fluoroscopically visualized through the stricture. Fluoroscopically visualize radiopaque markers on inner catheter at either end of the stent and position the inner radiopaque markers (**Fig. D,2 and D,4**) a minimum of 1 cm beyond the stricture so that it bridges stricture completely.
7. Confirm desired stent position fluoroscopically. To deploy stent, remove red safety guard from handle, (**Fig. B**) then squeeze the trigger. **Note:** Each full trigger squeeze will deploy stent by an equal amount.
8. If stent repositioning is required during deployment, it is possible to recapture stent. **Note:** It is **not** possible to recapture stent after passing point-of-no-return, indicated when red marker on top of handle has passed the point-of-no-return position on the handle label (**Fig. C**).

To fluoroscopically monitor the point-of-no-return (Fig. D, 5), the stent can be recaptured until the point where the outer catheter radiopaque marker (Fig. D, 1) is aligned with the second inner catheter radiopaque marker (Fig. D, 3).

Note: It is **not** possible to recapture stent after the outer radiopaque marker (Fig. D, 1) passes the second inner radiopaque marker (Fig. D, 3).

Note: The stent can be recaptured a maximum of three times.

9. To reposition the stent, the stent must first be recaptured and the elevator must be open. **Note:** Do not push forward on the delivery system with stent partially deployed. Push direction button on side of delivery handle to opposite side (Fig. E1). **Note:** Hold thumb on button when squeezing trigger for first time to recapture. Continue squeezing trigger as required to recapture stent by desired amount.

10. To resume deployment, push button to opposite side again and hold button for first stroke while squeezing trigger (Fig. E2).


11. When stent point-of-no-return has been passed, disconnect luer lock fitting to remove safety wire completely from delivery handle (Fig. F).






12. Continue deploying stent by squeezing trigger.

13. After deployment, fluoroscopically confirm full stent expansion. While maintaining wire guide position, push direction button on side of the delivery handle to opposite side (Fig. E1). Squeeze the trigger to completely recapture the delivery system. **Unlock the wire guide from Fusion Wire Guide locking device.** The device can be safely removed with the elevator of the duodenoscope fully down.

To reposition distally/remove stent during initial stent placement procedure, grasp the yellow grasping loop at end of stent with forceps (Fig. G) and reposition distally/remove as required.

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

MATERIALS

Evolution Biliary Stent System - Partially Covered (EVO-PC-8-9-6-B; EVO-PC-8-9-8-B; EVO-PC-10-11-4-B; EVO-PC-10-11-6-B; EVO-PC-10-11-8-B); Evolution Biliary Stent System - Fully Covered (EVO-FC-8-9-6-B; EVO-FC-8-9-8-B; EVO-FC-10-11-4-B; EVO-FC-10-11-6-B; EVO-FC-10-11-8-B)

Stent Components		Patient contact
Stent wire	Nickel-titanium (nitinol) wire with Platinum Core Nitinol alloy (nickel-titanium) Nickel 54.5 -57.0 Carbon 0.05 Max Oxygen 0.05 Max Hydrogen 0.005 Max Titan Decisive	Yes
Stent Cover	Silicone MED-4755	Yes
Stent Sub Assembly (includes stent coating)		
Stent Grasping Feature	Ultrahigh molecular weight polyethylene	Yes

	(UHMWPE) / Nitinol Nitinol alloy (nickel-titanium) Nickel 54.5 -57.0 Carbon 0.05 Max Oxygen 0.05 Max Hydrogen 0.005 Max Titan Decisive / PTFE Teflon	
Stent Cover (Second treatment)	Silicone MED-4755	Yes
Delivery system		
Second Assembly		
Location	Description	
Handle	EVO biliary nozzle	Lexan No
Introducer	Suture: white force fibre	Ultrahigh molecular weight polyethylene (UHMWPE) Yes
Introducer	Adhesive: Loctite 4061 (15544)	Adhesive (Cyanoacrylate) Yes
Introducer	Adhesive: Loctite 4161(26078)	Adhesive (Cyanoacrylate) No
Introducer	Retrieval Loop Assembly	Stainless steel AISI 304 (Alloy code: C Max 0.08 Cr 18 - 20 Fe 66.345 - 74 Mn Max 2 Ni 8 - 10.5 P Max 0.045 S Max 0.03 Cu Max 1) / Solder - Silver (alloy: 78% silver, 20% palladium and 2% platinum) Yes
Introducer	Tubing: Nylon Special Clear	Nylon 6 / Vestamid L2101F Yes
Introducer	Adhesive: Loctite 20268/26073 (14500)	Adhesive (Cyanoacrylate) No
Introducer	Cannula: Stainless Steel, Tapered	Stainless steel AISI 304 (Alloy code: C Max 0.08 Cr 18 - 20 Fe 66.345 - 74 Mn Max 2 Ni 8 - 10.5 P Max 0.045 S Max 0.03 Cu Max 1) No
Introducer	Tube: HDPE	High density polyethylene 276-73 Yes
Introducer	Tantalum Radiopaque Marker Band	Tantalum Grade R05200 Yes
Introducer	Adhesive	Loctite 4061 / Acrylic Yes

Introducer	Connector FLLA	Polycarbonate Iupilon	No
Introducer	Cannula: Stainless Steel, Flared	Stainless steel AISI 304 (Alloy code: C Max 0.08 Cr 18 - 20 Fe 66.345 - 74 Mn Max 2 Ni 8 - 10.5 P Max 0.045 S Max 0.03 Cu Max 1)	No
Handle	Adapter: Female Luerlock Adapter) (FLLA	Polycarbonate Iupilon	No
Handle	Belt Pulleys	Acetal	No
Handle	CAM Roller	Acetal	No
Handle	Handle Outer Shell Left	Iupilon polycarbonate / Acrylonitrile butadiene styrene blend (ABS POLIMAXX GA800)	No
Handle	Handle Outer Shell Right	Iupilon polycarbonate / Acrylonitrile butadiene styrene blend (ABS POLIMAXX GA800)	No
Handle	Idler Pin	Stainless steel AISI 304	No
Handle	RMN Pulley Stabiliser	Aluminum (Alloy A199.8 1080A: Silicon - 0.15%; Iron - 0.15%; Copper - 0.03%; Manganese - 0.02%; Magnesium - 0.02%; Zinc - 0.06%; Titanium -0.02%; Aluminum not less than 99.80%)	No
Handle	Optibelt Drive Belt	Polychloropen Neopren, Kevlar-49 Glasscord, Polyamide	No
Handle	SEMS Compression Spring	Stainless steel 304 with PTFE micron coating	No
Handle	PONR Label	Fasson PP Light Top	No
Introducer	EVO Biliary Introducer Tip	Pebax with with 40% Barium (Natural)	Yes
Introducer	Tubing	Nylon 12	Yes
Introducer	Tubing: PTFE Extruded	PTFE (Polytetrafluoroethylene), Natural	No
Introducer	Stainless Steel Crimping Band	Stainless steel 304	No
Introducer	Tubing: Nylon Special Clear	Nylon 12	Yes
Introducer	Tubing: Nylon Special Clear (Distal Outer Clear Tube)	Nylon 12, Natural Vestamid L2101F	Yes
Introducer	Adhesive: Loctite 4061 (15544)	Adhesive (Cyanoacrylate)	Yes
Introducer	Tubing: Braided Polyimide	Braided Polyimide Polyamide Elastomer - Blend of Nylon 12 and Polyester Block-Polyamide Copolymer	Yes

		(proprietary formula)	
Zip Port Sub Assembly			
Introducer	Tube	Nylon 12	No
Introducer	Adhesive: Loctite 4061 (15544)	Adhesive (Cyanoacrylate)	Yes
Introducer	EVO Biliary Shuttle Cap	Polycarbonate Iupilon	No
Introducer	Shuttle	Polycarbonate Lexan	No
Introducer	Tubing:	Nylon 12	No
Introducer	Flexor: Distal section	Pebax / PTFE / Stainless Steel Ring: Platinum with Iridium	Yes
Introducer	Flexor: Distal section	Westamide and Pebax / PTFE (Teflon) / 304 stainless steel Band: platinum with iridium	Yes
Introducer	Flexor: Hypotube Flexor, Braided (Proximal Flexor)	Nylon 6/304 Stainless Steel Wire Tube: PTFE (natural) Tecaflon PTFE	No
Introducer	Wire: Tibeta PTFE Coated (Lock Wire)	Ti-Beta PTFE coated (Teflon)	Yes
Introducer	Retaining Wire Housing	Polycarbonate Lexan	No
Introducer	Retaining Wire Anchor	Stainless steel 304	No
Introducer	Lock Wire Screw	Polycarbonate Lexan	No
Handle			
Various Components		PolyCarbonate Lexan / Acetal / Delrin / Stainless Steel 304	No

COMPLIANCE WITH RELEVANT REGULATIONS AND TECHNICAL STANDARDS

Cook Ireland Ltd manufacture all medical products in accordance with Medical Devices Directive 93/42/EEC.

All of Cook facilities are also certified to ISO 13485 (2016) *Quality Management System*.

DISPUTE RESOLUTION PROCEDURE

Please notify Cook for return authorization.

Cook Medical Europe Ltd

O'HALLORAN ROAD

NATIONAL TECHNOLOGY PARK

LIMERICK, IRELAND

TEL: +353 61 334440 FAX: +353 61 334441

WWW.COOKMEDICAL.COM

WARRANTY INFORMATION

Authorized representative on the territory of the Russian Federation - Company "SCHAG" (ZAO "SCHAG"). Russia, 119002, Moscow, Karmanitzki per. 9, floor 5, office 501

Tel. (495) 956-13-09;

Fax (495) 956-1310,

e-mail: info@schag.ru

The shelf life of products is two (2) year. Manufacturer's warranty extends to the quality of products.

SHELF LIFE

The "use by" date as denoted by symbols on the labeling. The shelf life is shown by indication the time of manufacture (factory symbol) to the stated expiration date (hourglass symbol). The difference between both dates is the shelf life.

Evolution Biliary Stent System - Partially Covered and Evolution Biliary Stent System - Full Covered have expiration date two (2) years.

STERILIZATION DETAILS

These products supplied sterile in peel open packages utilizing an EtO (Ethylene Oxide) validated method to a sterility assurance level (SAL) of 10^{-6} . These devices are intended for one-time use, sterile if unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile.

CONDITIONS OF USAGE

These devices are being used in the conditions of the patient care medical institutions.

STORAGE RULES

These products should be stored in a dry place away from extreme temperatures / data should be stored in a dark, dry, cool place, avoiding prolonged exposure to light.

TRANSPORTATION RULES


Transportation shall account for the type of transportation employed, the destination, the vendor used for transportation and an estimation of the frequency and harshness of any handling required during transportation.

LABELING

These products are supplied in the individual packaging, which should show no signs of damage or opening. Upon detecting the bends, breaks and other defects the device cannot be used.

The label of the device shows the following device parameters:

- the name of the medical device;
- trademark of the manufacturing company;
- catalog number;
- device diagram with dimensions;
- link to the glossary of symbols;
- method of sterilization;
- number of the series (lot);
- symbol "use before";
- symbol "date of manufacture";
- symbol "prohibition of reuse";
- the symbol "refer to the instructions for use";
- symbol "be careful!";
- symbol "to sale by or on the order of a physician";
- site of the manufacturing company;
- barcodes;
- the logo of the manufacturing company;
- name of the manufacturing company;
- address of the manufacturing company;
- symbol of the presence of the CE certificate.

REPAIRS

Devices are supplied sterile and for single-use only. Do not attempt to reuse or carry out maintenance.

DISPOSAL RULES

Upon completion of procedure, the devices are to be disposed of as per institutional guidelines for biohazardous medical waste. For EU countries - in accordance with Council Directive 93/42/EEC/ For Russian Federation – in accordance with national requirements, regulated by sanitary rules and norms 2.1.3684-21 of January 28, 2021 "Sanitary and epidemiological requirements for the maintenance of the territories of urban and rural settlements, for water bodies, drinking water and drinking water supply, atmospheric air, soils, living quarters, operation of industrial, public premises, organization and implementation of sanitary and anti-epidemic (preventive) measures". Epidemiologically hazardous waste - class B.

Maryna Zhukava

Regulatory Affairs Specialist

29-July-2021

Date

Niall Sheehy

NIALL SHEEHY
Notary Public
67 Glentworth St., Limerick.
Commissioned For Life

09-08-2021

Cook Ireland Ltd.

O'Halaran Road
National Technology Park
Limerick, Ireland

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