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**NIALL SHEEHY**

**NOTARY PUBLIC**

6/7, Glentworth St.,  
Limerick  
Ireland

Telephone +353 61 414355

Facsimile + 353 61 414738

Email [nsheehy@sellors.ie](mailto:nsheehy@sellors.ie)

TO ALL TO WHOM THESE PRESENTS SHALL COME, I, Niall Sheehy, Notary Public, duly authorised, admitted, and sworn, practicing at 6/7 Glentworth St., Limerick, Republic of Ireland, DO HEREBY CERTIFY that I am familiar with documents as issued from time to time by Cook Ireland Ltd and that the original document hereunto annexed namely, **Instructions for Use-EVOLUTION nitinol uncoated stents for maintaining patency of the gastrointestinal tract (GIT)** upon which I have signed my name and affixed my Seal of Office is a true and faithful original and that the original appears to me to be genuine.

AND I HAVE CERTIFIED as herein at the request of Cook Ireland Ltd having its place of business at O'Halloran Road, National Technology Park, in the city of Limerick, Ireland, and for the benefit of the proper authorities in **Russian Federation** and for all to whom it shall come.

IN FAITH AND TESTIMONY whereof I have hereunto set my hand and affixed my Seal of Office.


Niall Sheehy  
Notary Public

9<sup>th</sup> August 2021

Date

NIALL SHEEHY  
Notary Public  
6/7 Glentworth St., Limerick.  
Commissioned For Life

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This Apostille only certifies the authenticity of the signature and the capacity of the person who has signed the public document. It does not certify the content of the document for which it was issued. To verify the issuance of this Apostille, see [www.austlii.edu.au/au/other/dfat/ie](http://www.austlii.edu.au/au/other/dfat/ie)

327421

## INSTRUCTIONS FOR USE

### **EVOLUTION nitinol uncoated stents for maintaining patency of the gastrointestinal tract (GIT)**

#### **PURPOSE**

These devices are intended to maintain patency of the gastrointestinal tract.

#### **DEVICE DESCRIPTION**

Evolution Duodenal Stent System – Uncovered (EVO-22-27-12-D; EVO-22-27-6-D; EVO-22-27-9-D).

Evolution Colonic Stent System – Uncovered (EVO-25-30-10-C; EVO-25-30-6-C; EVO-25-30-8-C)

consists of the following two components:

- Self Expanding Metal Stent
- Stent Delivery System (this is made up of the handle and introducer)


#### **The Self Expanding Metal Stent:**

Self-expanding flexible stent is constructed of a single, woven, nitinol wire. The stent foreshortens due to its design. The total length of the stent is indicated by radiopaque markers (3 mm (+/- 0.25 mm)) on the inner catheter, indicating the actual length of the stent at nominal stent diameter.

Both ends of the stent are expanded: the increased diameter of the stent at the stent ends is intended to provide resistance to migration. To aid in visibility under fluoroscopy, there are four tantalum marker bands (0.991mm (+/- 0.127mm)) at either end of the stent. The proximal and distal ends of the stent are designed with numerous rounded “crown” to potentially reduce tissue trauma.

#### **The Stent Delivery System:**

The stent is mounted on an inner catheter, which accepts a 0.035 inch wire guide, and is constrained by an outer catheter. A pistol-grip delivery handle allows stent deployment or recapture. The introducer system consists of two catheters that operate co-axially with the preloaded stent placed between the inner catheter and outer catheter, and a handle to deliver the stent. The stent releases by retracting the outer catheter to expose the stent, then the stent self expands and conforms to the anatomy of the duodenum or colon.

**Figure:** Handle - Length 381.99mm (+/- 0.5mm)

**Figure:** Evolution Stent

Stent cell area:

Evolution Duodenal Stent System - Uncovered (EVO-22-27-12-D; EVO-22-27-6-D; EVO-22-27-9-D) - 3.82 mm<sup>2</sup> (± 10%)

Evolution Colonic Stent System - Uncovered (EVO-25-30-10-C; EVO-25-30-6-C; EVO-25-30-8-C) - 3.18 mm<sup>2</sup> (± 10%)

### **Construction**

1Fr = 0.33 mm

| Part Number    | Introduction System Diameter | Introduction System Length (cm) | Stent Length (mm) | Stent Body Diameter (mm) | Stent Flange Diameter (mm) |
|----------------|------------------------------|---------------------------------|-------------------|--------------------------|----------------------------|
| EVO-22-27-6-D  | 10 Fr/3.3mm                  | 230 cm ± 3cm                    | 60                | 22                       | 27                         |
| EVO-22-27-9-D  | 10 Fr/3.3mm                  | 230 cm ± 3cm                    | 90                | 22                       | 27                         |
| EVO-22-27-12-D | 10 Fr/3.3mm                  | 230 cm ± 3cm                    | 120               | 22                       | 27                         |
| EVO-25-30-6-C  | 10 Fr/3.3mm                  | 230 cm ± 3cm                    | 60                | 25                       | 30                         |
| EVO-25-30-8-C  | 10 Fr/3.3mm                  | 230 cm ± 3cm                    | 80                | 25                       | 30                         |
| EVO-25-30-10-C | 10 Fr/3.3mm                  | 230 cm ± 3cm                    | 100               | 25                       | 30                         |

Stent frame thickness - 0.200 mm (± 0.008 mm)

**Table:** Stent Dimensions

| Stent Length*                                                                       | Flange Length | Body Diameter | Flange Diameter | Location of marker band from end of stent |
|-------------------------------------------------------------------------------------|---------------|---------------|-----------------|-------------------------------------------|
| Evolution Duodenal Stent System - Uncovered (EVO-22-27-12-D; EVO-22-27-6-D; EVO-22- |               |               |                 |                                           |

| 27-9-D)                                                                                   |              |              |              |               |
|-------------------------------------------------------------------------------------------|--------------|--------------|--------------|---------------|
| 60 mm ± 3 mm                                                                              | 15 mm ± 2 mm | 22 mm ± 1 mm | 27 mm ± 1 mm | 5 mm ± 2.5 mm |
| 90 mm ± 3 mm                                                                              | 20 mm ± 2 mm |              |              |               |
| 120 mm ± 3 mm                                                                             | 20 mm ± 2 mm |              |              |               |
| Evolution Colonic Stent System - Uncovered (EVO-25-30-10-C; EVO-25-30-6-C; EVO-25-30-8-C) |              |              |              |               |
| 60 mm ± 3 mm                                                                              | 15 mm ± 2 mm | 25 mm ± 1 mm | 30 mm ± 1 mm | 5 mm ± 2.5 mm |
| 80 mm ± 3 mm                                                                              | 20 mm ± 2 mm |              |              |               |
| 100 mm ± 3 mm                                                                             | 20 mm ± 2 mm |              |              |               |

\* The dimension indicates the actual length of the stent at nominal stent diameter

**Table:** Dimension of the Delivery System

| Delivery System Length | Delivery System OD            | Width of Handle      |
|------------------------|-------------------------------|----------------------|
| 2300 mm $\pm$ 30 mm    | 3.40 mm (10 Fr) $\pm$ 0.08 mm | 382 mm $\pm$ 0.10 mm |

**Table:** Dimension of Outer catheter

| Stent length                                                                               | Outer Catheter Length | Distal Section (Clear) Length of the Outer Catheter |
|--------------------------------------------------------------------------------------------|-----------------------|-----------------------------------------------------|
| Evolution Duodenal Stent System - Uncovered (EVO-22-27-12-D; EVO-22-27-6-D; EVO-22-27-9-D) |                       |                                                     |
| 60 mm                                                                                      | 2,359 mm ± 2 mm       | 138 mm +/- 2 mm                                     |
| 90 mm                                                                                      |                       | 191 mm +/- 2 mm                                     |
| 120 mm                                                                                     |                       | 246 mm +/- 2 mm                                     |
| Evolution Colonic Stent System - Uncovered (EVO-25-30-10-C; EVO-25-30-6-C; EVO-25-30-8-C)  |                       |                                                     |
| 60 mm                                                                                      | 2,359 mm ± 2 mm       | 152 mm +/- 2 mm                                     |
| 80 mm                                                                                      |                       | 189 mm +/- 2 mm                                     |
| 100 mm                                                                                     |                       | 237 mm +/- 2 mm                                     |

**Table:** Dimension of Delivery system (Inner catheter)

|                                                                                            | Pushing catheter Length | Length of Tubing | Distance between yellow marker and tip | Distance between proximal ends of marker bands* |
|--------------------------------------------------------------------------------------------|-------------------------|------------------|----------------------------------------|-------------------------------------------------|
| Evolution Duodenal Stent System - Uncovered (EVO-22-27-12-D; EVO-22-27-6-D; EVO-22-27-9-D) |                         |                  |                                        |                                                 |
| 60 mm                                                                                      | 2,542 mm ± 2 mm         | 1,900 mm ± 10 mm | 123 mm-126 mm                          | 60 mm ± 3 mm                                    |
| 90 mm                                                                                      | 2,489 mm ± 2 mm         |                  | 176 mm-179 mm                          | 90 mm ± 3 mm                                    |
| 120 mm                                                                                     | 2,434 mm ± 2 mm         |                  | 231 mm-234 mm                          | 120 mm ± 3 mm                                   |
| Evolution Colonic Stent System - Uncovered (EVO-25-30-10-C; EVO-25-30-6-C; EVO-25-30-8-C)  |                         |                  |                                        |                                                 |

|        |                     |                      |                |                   |
|--------|---------------------|----------------------|----------------|-------------------|
| 60 mm  | 2,528 mm $\pm$ 2 mm | 1,900 mm $\pm$ 10 mm | 137 mm -140 mm | 60 mm $\pm$ 3 mm  |
| 80 mm  | 2,491 mm $\pm$ 2 mm |                      | 174 mm -177 mm | 80 mm $\pm$ 3 mm  |
| 100 mm | 2,443 mm $\pm$ 2 mm |                      | 222 mm -225 mm | 100 mm $\pm$ 3 mm |

\* This indicates stent nominal length.

The markers on the inner catheter are 3 mm (+/- 0.25 mm).

Tantalum markers are 0.991mm (+/- 0.127mm).

**Table: Dimension of Tip**

| Tip Length         | Tip Outer Diameter at Dimpled Portion | Tip Inner Diameter    |                      |
|--------------------|---------------------------------------|-----------------------|----------------------|
|                    |                                       | Proximal site         | Distal site          |
| 11 mm $\pm$ 0.3 mm | 3.12 mm $\pm$ 0.05 mm                 | 1.22 mm $\pm$ 0.05 mm | 1.0 mm $\pm$ 0.05 mm |

### Change of length

#### Stent Dimensions:

Objective: To evaluate the stent dimensions of the Evolution Colonic/ Duodenal Stent post 25 months real time aging and post deployment.

Number of Samples: Twenty-Two (22) test samples of 10 cm accelerated aged stents (25mm Body Diameter 30mm Flare Diameter, 10 cm long). If all samples meet the acceptance criteria, a sample size gives 90% confidence of 90% reliability of the population of all stents meeting the acceptance criteria. All samples were double ETO sterilised and have been subjected to 25 months real time aging.

#### Test conditions:

Following stent deployment/ simulated use testing, the stents were placed in a water bath which was maintained at 37°C  $\pm$  2°C. After 30 minutes the outer diameter of the stent body and two flares was measured.

#### Acceptance Criteria and Results:

**Table: Stent Dimensions (Post Aging) Acceptance Criteria & Results:**

| Testing          | Acceptance Criteria                                                                                                                                  | Result                                                                                                                                                                                                 |
|------------------|------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Stent Dimensions | Body Diameter: 22.8 – 27.3 mm<br>Flare Diameter: 27.55 – 32.55 mm<br><br>(After deployment allowing the stent to recover at 37°C +/-2°C for 30mins.) | <b>Pass</b><br><br><u>Body Diameter:</u><br>Min: 24.29 mm<br>Max: 25.11 mm<br><u>Flare 1 Diameter:</u><br>Min: 29.00 mm<br>Max: 29.98 mm<br><u>Flare 2 Diameter:</u><br>Min: 29.03 mm<br>Max: 30.25 mm |

**Analysis:** All the test samples met the predetermined acceptance criteria following 25 month real time aging testing. The results of this test support a 24 month shelf life for the Evolution Colonic/ Duodenal Stent System.

### **Stent deployment strength**

**Table:** Results Radial (compression) force.

| Description                                                                  | Radial (compression) force per unit length N/mm |                                       |
|------------------------------------------------------------------------------|-------------------------------------------------|---------------------------------------|
|                                                                              | At stent compressed to 15 mm diameter           | At stent compressed to 10 mm diameter |
| Cook Evolution Duodenal Stent (EVO-22-27-12-D; EVO-22-27-6-D; EVO-22-27-9-D) | 0.176 (Mean of 10 results)                      | 0.205 (Mean of 10 results)            |
| Cook Evolution Colonic Stent (EVO-25-30-10-C; EVO-25-30-6-C; EVO-25-30-8-C)  | 0.154 (Mean of 10 results)                      | 0.178 (Mean of 10 results)            |

**Table:** Results Expansion force.

| Description                                                                  | Expansion force per unit length N/mm  |                                       |
|------------------------------------------------------------------------------|---------------------------------------|---------------------------------------|
|                                                                              | At stent compressed to 15 mm diameter | At stent compressed to 10 mm diameter |
| Cook Evolution Duodenal Stent (EVO-22-27-12-D; EVO-22-27-6-D; EVO-22-27-9-D) | 0.092 (Mean of 10 results)            | 0.141 (Mean of 10 results)            |
| Cook Evolution Colonic Stent (EVO-25-30-10-C; EVO-25-30-6-C; EVO-25-30-8-C)  | 0.082 (Mean of 10 results)            | 0.125 (Mean of 10 results)            |

### **Connection strength**

| Testing                                | Acceptance Criteria | Pass $\checkmark$<br>Fail X  | Force (N)                   |
|----------------------------------------|---------------------|------------------------------|-----------------------------|
| Distal tip joint strength              | Peak force > 15 N.  | $\checkmark$                 | Min: 25.4<br>Max: 87.5      |
| Marker band joint strength             | Peak force > 5 N.   | $\checkmark$                 | Min: 20.2<br>Max: 44.8      |
| Proximal PEEK flare joint strength     | Peak force > 45 N.  | $\checkmark$                 | Min: 137<br>Max: 175        |
| PTFE / tantalum band crimp strength    | Peak force > 15 N.  | $\checkmark$                 | Min: 38.5<br>Max: 89.9      |
| Cannula to Nitinol wire joint strength | Peak force > 15 N.  | $\checkmark$                 | Min: 61.7<br>Max: 95.7      |
| Retrieval loop retention force         | Peak force > 30N.   | $\checkmark$                 | Min: 46.986<br>Max: 171.956 |
| Marker band joint                      | Peak force > 5 N.   | <b>Pass</b><br>Min = 7.239 N |                             |

|                                   |                    |                                                                      |                       |
|-----------------------------------|--------------------|----------------------------------------------------------------------|-----------------------|
| strength                          |                    | Max = 24.900 N<br>Mean = 15.994 N                                    |                       |
| Retrieval Loop Tensile Strength   | Peak force > 30N   | <b>Pass</b><br>Min = 99.845 N<br>Max = 181.821 N<br>Mean = 161.338 N |                       |
| Outer sheath flare joint strength | Peak force > 45 N. | √                                                                    | Min: 77.4<br>Max: 105 |
| PEEK / Polyimide joint strength   | Peak force > 15 N. | √                                                                    | Min: 27.3<br>Max: 198 |

|                           |                                              |                                                 |
|---------------------------|----------------------------------------------|-------------------------------------------------|
| Handle Mechanism Strength | Deployment / Desheathing: Peak force > 5Kg   | Min: 130 N / 13.26 kg<br>Max: 205 N / 20.9 kg   |
|                           | Recapturing / Resheathing : Peak force > 5Kg | Min: 99.5 N / 10.15 kg<br>Max: 172 N / 17.54 kg |

### **Fatigue testing**

#### **Sling Radial Fatigue Testing**

**Objective:** To evaluate the fatigue life of the duodenal and colonic stents by subjecting them to time-accelerated, physiologically modelled, diameter-controlled radial compression for a specified number of cycles. Test articles were submersed in distilled water maintained at 37°C+/-2°C.

Table: Sling Radial Fatigue Testing Acceptance Criteria & Results.

| Testing                      | Acceptance Criteria                                                | Pass √<br>Fail X |
|------------------------------|--------------------------------------------------------------------|------------------|
| Sling Radial Fatigue Testing | Stent must remain intact                                           | √                |
|                              | Stent wire must not break during the test                          | √                |
|                              | Marker band must remain on the stent and not detach from the stent | √                |
|                              |                                                                    |                  |

**Analysis:** The test samples met the predetermined acceptance criteria demonstrating an acceptable fatigue life of the stent. In each case, upon visual inspection there was no fracture of the stent wire, loss of stent marker bands or changes in the stent geometry.

Temperature conditions for transportation:

The products are stable at a temperature of transportation -50 to +50 C.

## **INSTRUCTIONS FOR USE**

### **Evolution Duodenal Stent System - Uncovered (EVO-22-27-12-D; EVO-22-27-6-D; EVO-22-27-9-D)**

**CAUTION:** U.S. Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

### **IMPORTANT INFORMATION**

Please review prior to use.

### **INTENDED USE**

This device is used for palliative treatment of duodenal or gastric outlet obstruction and duodenal strictures caused by malignant neoplasms.

### **STENT DESCRIPTION**

This flexible, self-expanding stent is constructed of a single, woven, nitinol wire. The stent has a flange at both ends to aid in preventing migration. The stent foreshortens due to its design. The total length of the stent is indicated by radiopaque markers on the inner catheter, indicating the actual length of the stent at nominal stent diameter. To aid in visibility under fluoroscopy, there are four tantalum marker bands at either end of the stent. The proximal and distal ends of the stent are designed with numerous rounded "crowns" to potentially reduce tissue trauma.

### **DELIVERY SYSTEM DESCRIPTION**

The stent is mounted on an inner catheter, which accepts a .035-inch wire guide, and is constrained by an outer catheter. An endoscopically and fluoroscopically visible yellow marker defines the proximal end of the stent when constrained in the catheter. A pistol-grip delivery handle allows controlled stent deployment or recapture. The position of the red marker on top of the handle in relation to the point-of-no-return position on the handle label indicates the point after which it is not possible to recapture the deployed stent.

### **CONTRAINDICATIONS**

Those specific to GI endoscopy and any procedure to be performed in conjunction with stent placement. Additional contraindications include, but are not limited to: enteral ischemia, suspected or impending perforation, intra-abdominal abscess/perforation, inability to pass wire guide or stent through obstructed area, patients for whom endoscopic procedures are contraindicated, coagulopathy/patients with elevated bleeding times, benign disease.

### **POTENTIAL ADVERSE EVENTS**

Potential adverse events associated with GI endoscopy include, but are not limited to: airway obstruction, allergic reaction to contrast or medication, aspiration, biliary obstruction, cardiac arrhythmia or arrest, cholangitis, fever, hemorrhage, hypotension, infection, perforation, reflux, respiratory depression or arrest.

Additional adverse events include, but are not limited to: allergic reaction to nickel, bowel impaction, death (other than due to normal disease progression), erosion of the luminal mucosa, foreign body sensation, inadequate expansion, intestinal perforation, nausea/vomiting, pain/discomfort, pancreatitis, pressure necrosis, septicemia, stent misplacement and/or migration, stent occlusion, tumor ingrowth or overgrown, ulcerations.

## **PRECAUTIONS**

Refer to product package label for the minimum channel size required for this device.

A complete diagnostic evaluation must be performed prior to use to determine proper stent size.

If wire guide or stent cannot advance through obstructed area, do not attempt to place stent.

Stent should be placed using fluoroscopic monitoring with endoscopy.

The stent should only be placed with the Cook delivery system, which is provided with each stent.

This device is intended for palliative treatment only. Alternate methods of therapy should be investigated prior to placement.

Long-term patency of this device has not been established. Periodic evaluation is advised.

Stent cannot be retrieved after the deployment threshold has been passed. The position of the red marker on the top of the handle in relation to the point-of-no-return position on the handle label indicates when threshold has been passed.

Do not use this device for any purpose other than stated intended use.

Use of this device is restricted to a trained healthcare professional.

## **WARNINGS**

The stent is not intended to be removed or repositioned after stent placement and is intended to remain in the body permanently. Attempts to remove or reposition stent after placement may cause damage to surrounding tissue or mucosa.

The stent contains nickel, which may cause an allergic reaction in individuals with nickel sensitivity.

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

After stent placement, alternative methods of treatment such as chemotherapy and radiation should not be administered as this may increase risk of stent migration due to tumor shrinkage, stent erosion, and/or mucosal bleeding.

## **MR INFORMATION**

This symbol means the device is MR Conditional.

Non-clinical testing has demonstrated that this stent is MR Conditional according to ASTM F2503. A patient with this stent can be scanned safely immediately after placement under the following conditions.

### **Static Magnetic Field**

- Static magnetic field of 3 Tesla or less
- Maximum spatial magnetic gradient of 720 Gauss/cm or less

Non-clinical evaluation was conducted in a 3 Tesla MR system (Excite, General Electric Healthcare) with a maximum spatial magnetic gradient field of 720 Gauss/cm as measured with a gaussmeter in the position of the static magnetic field pertinent to the patient (i.e., outside of scanner covering, accessible to a patient or individual).

### **MRI-Related Heating**

•1.5 and 3.0 Tesla Systems: It is recommended to scan in normal operation mode (whole body averaged specific absorption rate (SAR)  $\leq 2.0$  W/kg) ("Normal Operating Mode" is defined as the mode of operation of the MR system in which none of the outputs have a value that cause physiological stress to the patient) for 15 minutes of scanning (i.e., per scanning sequence).

•Non-clinical testing was conducted on the stent under the following conditions, and produced a maximum temperature rise of 3.7 °C:

- a maximum whole body average specific absorption rate (SAR) of 2.9 W/kg (corresponding to a calorimetry measured value of 2.1 W/kg) for 15 minutes of MR scanning in a 1.5 Tesla Magnetom (Siemens Medical Solutions, Malvern, PA, Software Numaris/4) MR scanner.

- a maximum whole body average specific absorption rate (SAR) of 2.9 W/kg (corresponding to a calorimetry measured value of 2.7 W/kg) for 15 minutes of MR scanning in a 3.0 Tesla Excite (GE Electric Healthcare, Milwaukee, WI, Software 14X.M5) MR scanner.

#### **Image Artifacts**

MR image quality may be compromised if the area of interest is within the lumen of the Evolution Stent or within approximately 6 mm of the position of the Evolution Stent as found during non-clinical testing using the sequences: T1-weighted, spin echo pulse sequence and gradient echo pulse sequence in a 3.0 Tesla Excite (GE Electric Healthcare, Milwaukee, WI, Software 14X.M5) MR system with body radiofrequency coil. Therefore, it may be necessary to optimize MR imaging parameters for the presence of this metallic stent.

Additional Information: The safety of performing an MRI procedure in a patient with overlapping duodenal stents or other MRI-conditional device(s) in direct contact with this device has not been determined. Performing MRI in such situations is not recommended.

#### **EQUIPMENT REQUIRED**

Stent and delivery system

.035 inch wire guide

#### **SELECTION OF STENT**

A stent length of at least 4 cm longer than the stricture is recommended (e.g. approximately 2 cm longer at either end of the stricture)

#### **PATIENT PREPARATION**


With duodenoscope in place, introduce a wire guide until it is fluoroscopically visualized in position through stricture.

#### **INSTRUCTIONS FOR USE**

1. Remove stent delivery system from the package.
2. Introduce delivery system in short increments over wire guide, into accessory channel, until it is endoscopically visualized exiting scope.
3. To ensure that stent will bridge stricture after deployment, fluoroscopically position radiopaque markers on inner catheter beyond extremities of stricture to be crossed. **Note:** Yellow marker on the delivery system located at the proximal end of stent can also serve as an endoscopic/fluoroscopic reference for positioning of the proximal (nearest user) end of the stent relative to the stricture.
4. Confirm desired stent position fluoroscopically and deploy stent by removing red safety guard from the handle. (**Fig. 1**)
5. Continue deploying stent by squeezing trigger. **Note:** Each trigger squeeze will deploy stent by an equal amount.
6. If stent repositioning is required during deployment, it is possible to recapture stent. **Note:** It is not possible to recapture stent after passing point-of-no-return, indicated when red marker on top of the handle has passed the point-of-no-return position on the handle label. The stent can be recaptured a maximum of five times. (**Fig. 2**)

7. To recapture stent, push direction button on side of delivery system to opposite side. (Fig. 3) **Note:** Hold thumb on button when squeezing trigger for first time to recapture.
8. Continue squeezing trigger as required to recapture stent by desired amount.
9. To resume deployment, push button to opposite side again and hold thumb on button for first stroke while squeezing trigger. (Fig. 4)
10. When stent point-of-no-return has been passed, pull safety wire out of delivery handle near wire guide port. (Fig. 5)
11. Continue deploying stent by squeezing trigger.
12. After deployment, fluoroscopically confirm full stent expansion. Once full expansion is confirmed, delivery system can be safely removed.

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

## **INSTRUCTIONS FOR USE**

**Evolution Colonic Stent System - Uncovered (EVO-25-30-10-C; EVO-25-30-6-C; EVO-25-30-8-C)**

**CAUTION:** U.S. Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

### **IMPORTANT INFORMATION**

Please review prior to use.

### **INTENDED USE**

This device is used for palliative treatment of colonic obstruction or colonic strictures caused by malignant neoplasms, and to relieve large bowel obstruction prior to colectomy in patients with malignant strictures.

### **STENT DESCRIPTION**

This flexible, self-expanding stent is constructed of a single, woven, nitinol wire. The stent has a flange at both ends to aid in preventing migration. The stent foreshortens due to its design. The total length of the stent is indicated by radiopaque markers on the inner catheter, indicating the actual length of the stent at nominal stent diameter. To aid in visibility under fluoroscopy, there are four main tantalum marker bands at either end of the stent. The proximal and distal ends of the stent are designed with numerous rounded "crowns" to potentially reduce tissue trauma.

### **DELIVERY SYSTEM DESCRIPTION**

The stent is mounted on an inner catheter, which accepts a .035-inch wire guide, and is constrained by an outer catheter. An endoscopically and fluoroscopically visible yellow marker defines the proximal end of the stent when constrained in the catheter. A pistol-grip delivery handle allows controlled stent deployment or recapture. The position of the red marker on top of the handle in relation to the point-of-no-return position on the handle label indicates the point after which it is not possible to recapture the deployed stent.

### **CONTRAINDICATIONS**

Those specific to GI endoscopy and any procedure to be performed in conjunction with stent placement. Additional contraindications include, but are not limited to: enteral ischemia, suspected or impending perforation, intra-abdominal abscess/ perforation, inability to pass wire guide or stent through obstructed area, patients for whom endoscopic procedures are contraindicated, significant coagulopathy, benign disease.

### **POTENTIAL ADVERSE EVENTS**

Potential adverse events associated with GI endoscopy include, but are not limited to: allergic reaction to contrast or medication, aspiration, cardiac arrhythmia or arrest, fever, hemorrhage, hypotension, infection, perforation, respiratory depression or arrest.

Additional adverse events include, but are not limited to: allergic reaction to nickel, bowel impaction, constipation, death (other than due to normal disease progression), diarrhea, erosion of the luminal mucosa, foreign body sensation, inadequate stent expansion, intestinal perforation, nausea/vomiting, pain/discomfort, peritonitis, pressure necrosis, septicemia, stent misplacement and/or migration, stent occlusion, symptoms of tenesmus or urgency/incontinence, tumor ingrowth or overgrowth, ulcerations.

### **PRECAUTIONS**

Page 11 of 18

Refer to product package label for the minimum channel size required for this device.  
A complete diagnostic evaluation must be performed prior to use to determine proper stent size.  
If wire guide or stent cannot advance through obstructed area, do not attempt to place stent.  
Stent should be placed endoscopically with fluoroscopic monitoring.  
The stent should only be placed with the Cook delivery system, which is provided with each stent.  
This device is intended for palliative treatment only. Alternate methods of therapy should be investigated prior to placement.  
Long-term patency of this device has not been established. Periodic evaluation is advised.  
Use of this device is restricted to a trained healthcare professional.  
Do not use this device for any purpose other than stated intended use.  
Stent cannot be retrieved after the deployment threshold has been passed. The position of the red marker on the top of the handle in relation to the point-of-no-return position on the handle label indicates when threshold has been passed.

### **WARNINGS**

The stent is not intended to be removed or repositioned after stent placement and is intended to remain in the body permanently. Attempts to remove or reposition stent after placement may cause damage to surrounding tissue or mucosa.

This stent contains nickel, which may cause an allergic reaction in individuals with nickel sensitivity.

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

To minimize pain and tenesmus, the stent end nearest to the anal canal/ anus should be placed 2 cm above the anal canal or 6 cm from the anus.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

After stent placement, alternative methods of treatment such as chemotherapy and radiation should not be administered as this may increase risk of stent migration due to tumor shrinkage, stent erosion, and/or mucosal bleeding.

The device should be used with caution and only after careful consideration in patients with:

- Patients with radiation colitis or proctitis.
- Patients with elevated bleeding times, coagulopathies.

### **MRI INFORMATION**

**This symbol means the device is MR conditional**

Non-clinical testing has demonstrated that the Evolution Colonic Stent is MR Conditional according to ASTM F2503. A patient with this stent can be scanned safely anytime after placement under the following conditions:

- Static magnetic field of 3.0 Tesla or less.
- Spatial magnetic gradient field of 720 Gauss/cm. Non-clinical evaluation was conducted in an MR system (Excite, General Electric Healthcare) with a maximum spatial magnetic gradient field of 720 Gauss/cm as measured with a gaussmeter in the position of the static magnetic field pertinent to the patient (i.e., outside of scanner covering, accessible to a patient or individual).

#### **MRI-Related Heating**

- 1.5 and 3.0 Tesla Systems: It is recommended to scan in normal operation mode (whole body averaged

specific absorption rate (SAR) = 2.0 W/kg) ("Normal Operating Mode" is defined as the mode of operation of the MR system in which none of the outputs have a value that cause physiological stress to the patient) for 15 minutes of scanning (i.e., per scanning sequence).

#### 1.5 Tesla Temperature Rise

In non-clinical testing, the Evolution Colonic Stent produced maximum temperature rises of 1.9 °C during 15 minutes of MR imaging (i.e., for one scanning sequence) performed in a 1.5 Tesla MR System (Magnetom, Siemens Medical Solutions, Malvern, PA, Software Numaris/4, Version Syngo MR 2002B DHHS Active-shielded) at an MR system reported whole-body-averaged SAR of 2.9 W/kg (associated with a calorimetry measured whole body averaged value of 2.1 W/kg).

#### 3.0 Tesla Temperature Rise

In non-clinical testing, the Evolution Colonic Stent produced maximum temperature rises of 2.7 °C during 15 minutes of MR imaging (i.e., for one scanning sequence) performed in a 3 Tesla MR System (Excite, HDx, General Electric Healthcare, Milwaukee, WI, Software 14X.M5) at an MR system reported whole-body-averaged SAR of 2.9 W/kg (associated with a calorimetry measured whole body averaged value of 2.7 W/kg).

#### **Image Artifacts**

MR image quality may be compromised if the area of interest is within the lumen or within approximately 3 mm of the position of the Evolution Colonic Stent, as found during non-clinical testing using the sequence: T1-weighted, spin echo and gradient echo pulse sequence in a 3.0 Tesla MR system (Excite, General Electric Healthcare, Milwaukee, WI, Software 14X.M5). Therefore, it may be necessary to optimize MR imaging parameters for the presence of this metallic stent.

Additional Information: The safety of performing an MRI procedure in a patient with overlapping colonic stents or other MRI-conditional device(s) in direct contact with this device has not been determined. Performing MRI in such situations is not recommended.

### EQUIPMENT REQUIRED

Stent and delivery system  
.035-inch wire guide

### SELECTION OF STENT

A stent length of at least 4 cm longer than the stricture is recommended (e.g. approximately 2 cm longer at either end of the stricture)

### PATIENT PREPARATION

With colonoscope in place, introduce a wire guide until it is fluoroscopically visualized in position through stricture.

### INSTRUCTIONS FOR USE

1. Remove stent delivery system from the package.
2. Introduce delivery system in short increments over wire guide, into accessory channel, until it is endoscopically visualized exiting the scope.
3. To ensure that stent will bridge stricture after deployment, fluoroscopically position radiopaque markers on inner catheter beyond extremities of stricture to be crossed. **Note:** Yellow marker on the delivery system located at proximal end of stent can also serve as an endoscopic/fluoroscopic reference for positioning of the proximal (nearest user) end of the stent relative to the stricture.
4. Confirm desired stent position fluoroscopically and deploy stent by removing red safety guard from the handle. (**Fig. 1**)
5. Continue deploying stent by squeezing trigger. **Note:** Each trigger squeeze will deploy stent by an

Page 13 of 18

equal amount.

6. If stent repositioning is required during deployment, it is possible to recapture stent. **Note:** It is not possible to recapture stent after passing point-of-no-return, indicated when red marker on top of the handle has passed the point-of-no-return position on handle label. (Fig. 2) The stent can be recaptured a maximum of five times.

7. To recapture stent, push direction button on side of delivery system to opposite side. (Fig. 3) **Note:** Hold thumb on button when squeezing trigger for first time to recapture.

8. Continue squeezing trigger as required to recapture stent by desired amount.


9. To resume deployment, push button to opposite side again and hold thumb on button for first stroke while squeezing trigger. (Fig. 4)

10. When stent point-of-no-return has been passed, pull safety wire out of delivery handle near wire guide port. (Fig. 5)

11. Continue deploying stent by squeezing trigger.

12. After deployment, fluoroscopically confirm full stent expansion. Once full expansion is confirmed, delivery system can be safely removed.

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

## MATERIALS

| Component Description                       | Material                                                                                                                                                                 | Patient contacting |
|---------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------|
| <b>Stent Components</b>                     |                                                                                                                                                                          |                    |
| Stent Nitinol Wire                          | Nitinol<br>Part            Mass%<br>Nickel        54.5 -57.0<br>Carbon        0.05 Max<br>Oxygen        0.05 Max<br>Hydrogen     0.05 Max<br>Titanium      Decisive      | Yes                |
| Stent Marker Bands                          | Tantalum Grade R05200                                                                                                                                                    | Yes                |
| Silicone                                    | Silicone        Grade        Vinyl        Special<br>Polydimethylsiloxane (MED-4755)                                                                                     | Yes                |
| <b>Introducer System</b>                    |                                                                                                                                                                          |                    |
| Introducer System Tip (Inner Catheter)      | Pebax with 40% Barium                                                                                                                                                    | Yes                |
| Polyimide Tubing (Inner Catheter)           | Braided Polyimide<br>Polyamide Elastomer - Blend of Nylon 12 and Polyester<br>Block-Polyamide Copolymer (proprietary formula)                                            | Yes                |
| Introducer Marker Bands (Inner Catheter)    | Tantalum Grade R05200                                                                                                                                                    | Yes                |
| Introducer Pusher (Inner Catheter)          | Inner: Nylon 12 with Bismuth<br>Outer: Nylon 12                                                                                                                          | Yes                |
| Introducer PEEK Tubing (Inner Catheter)     | Annealed Polyetheretherketone (PEEK)<br>381G                                                                                                                             | Yes                |
| Introducer PTFE Tubing (Inner Catheter)     | Polytetrafluoroethylene (PTFE) Teflon                                                                                                                                    | Yes                |
| Introducer Nylon Tubing (Inner Catheter)    | Nylon 6                                                                                                                                                                  | No                 |
| Retrieval Loop (Inner Catheter)             | Stainless Steel ASTM A313                                                                                                                                                | Yes                |
| Introducer support cannula (Inner Catheter) | Stainless Steel ASTM A313                                                                                                                                                | No                 |
| Nitinol Retention wire                      | Nitinol Wire<br>Part            Mass%<br>Nickel        54.5 -57.0<br>Carbon        0.05 Max<br>Oxygen        0.05 Max<br>Hydrogen     0.05 Max<br>Titanium      Decisive | Yes                |
| Outer Catheter                              | PTFE lining (MS # 122), Nylon 6 Outer surface,<br>Braid: Stainless Steel ASTM A313<br>Coil: Stainless Steel ASTM A313                                                    | Yes                |

| Component Description | Material                                                                                                                                                                                                                                                                                                                                    | Patient contacting |
|-----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------|
| Cannula               | Stainless Steel ASTM A313                                                                                                                                                                                                                                                                                                                   | No                 |
| Tubing                | Tantalum Grade R05200                                                                                                                                                                                                                                                                                                                       | Indirect           |
| <b>Handle</b>         |                                                                                                                                                                                                                                                                                                                                             |                    |
| Various components    | Acetal / PolyCarbonate GE Plastics Grade Lexan / 10% glass-filled polycarbonate / PolyCarbonate Grade Lexan / PolyCarbonate Grade Lexan / Acrylonitrile butadiene styrene blend Lanxess / Stainless Steel 304 with PTFE coating Teflon / Polychloroprene Neopren, Glass fibre Kevlar-49/steel tension members / Stainless Steel / Aluminium | No                 |
| <b>Consumables</b>    |                                                                                                                                                                                                                                                                                                                                             |                    |
| Adhesive              | Loctite 4061                                                                                                                                                                                                                                                                                                                                | Yes                |
| Adhesive              | Loctite 20268/ 26073/ 4013                                                                                                                                                                                                                                                                                                                  | No                 |

#### **COMPLIANCE WITH RELEVANT REGULATIONS AND TECHNICAL STANDARDS**

Cook Ireland Ltd manufacture all medical products in accordance with Medical Devices Directive 93/42/EEC.

All of Cook facilities are also certified to ISO 13485 (2016) *Quality Management System*.

#### **DISPUTE RESOLUTION PROCEDURE**

Please notify Cook for return authorization.

Cook Medical Europe Ltd

O'HALLORAN ROAD

NATIONAL TECHNOLOGY PARK

LIMERICK, IRELAND

TEL: +353 61 334440 FAX: +353 61 334441

WWW.COOKMEDICAL.COM

#### **WARRANTY INFORMATION**

Authorized representative on the territory of the Russian Federation - Company "SCHAG" (ZAO "SCHAG"). Russia, 119002, Moscow, Karmanitzki per. 9, floor 5, office 501

Tel. (495) 956-13-09;

Fax (495) 956-1310,

e-mail: [info@schag.ru](mailto:info@schag.ru)

The shelf life of products is two (2) year. Manufacturer's warranty extends to the quality of products.

#### **SHELF LIFE**

The "use by" date as denoted by symbols on the labeling. The shelf life is shown by indication the time of manufacture (factory symbol) to the stated expiration date (hourglass symbol). The difference between both dates is the shelf life.

Evolution Duodenal Stent System – Uncovered (EVO-22-27-12-D; EVO-22-27-6-D; EVO-22-27-9-D), Evolution Colonic Stent System – Uncovered (EVO-25-30-10-C; EVO-25-30-6-C; EVO-25-30-8-C) have a two (2) year shelf life.

### **STERILIZATION DETAILS**

These devices are supplied sterile in peel open packages utilizing an EtO (Ethylene Oxide) validated method to a sterility assurance level (SAL) of  $10^{-6}$ . These devices are intended for one-time use, sterile if unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile.

### **CONDITIONS OF USAGE**

These devices are being used in the conditions of the patient care medical institutions.

### **STORAGE RULES**

These devices should be stored in a dry location, away from temperature extremes / these devices should be stored in a dark, dry, cool place, with avoidance of extended exposure to light.

### **TRANSPORTATION RULES**



Transportation shall account for the type of transportation employed, the destination, the vendor used for transportation and an estimation of the frequency and harshness of any handling required during transportation.

### **PACKAGING AND LABELLING**

Supplied in the individual packaging, which should show no signs of damage or opening. Upon detecting the bends, breaks and other defects the device cannot be used.

The label of the device shows the following device parameters:

- the name of the medical device;
- trademark of the manufacturing company;
- catalog number;
- device diagram with dimensions;
- link to the glossary of symbols;
- method of sterilization;
- number of the series (lot);
- symbol "use before";
- symbol "date of manufacture";
- symbol "prohibition of reuse";
- the symbol "refer to the instructions for use";
- symbol "be careful!";
- site of the manufacturing company;
- barcodes;
- the logo of the manufacturing company;
- name of the manufacturing company;
- address of the manufacturing company;
- symbol of the presence of the CE certificate.

## REPAIRS

Devices are supplied sterile and for single-use only. Do not attempt to reuse or carry out maintenance.

## DISPOSAL RULES

Upon completion of procedure, the devices are to be disposed of as per institutional guidelines for biohazardous medical waste. For EU countries - in accordance with Council Directive 93/42/EEC/ For Russian Federation – in accordance with national requirements, regulated by sanitary rules and norms 2.1.3684-21 of January 28, 2021 "Sanitary and epidemiological requirements for the maintenance of the territories of urban and rural settlements, for water bodies, drinking water and drinking water supply, atmospheric air, soils, living quarters, operation of industrial, public premises, organization and implementation of sanitary and anti-epidemic (preventive) measures". Epidemiologically hazardous waste - class B.

Maryna Zhukava  
Regulatory Affairs Specialist

29-July-2021  
Date

## Cook Ireland Ltd.

O'Halaran Road  
National Technology Park  
Limerick, Ireland

Company Registration No. 230176  
VAT Registration No. IE-8203176.0

Niall Sheehy  
NIAL SHEEHY  
Notary Public  
67 Glentworth St., Limerick.  
Commissioned For Life

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