Directions for Use

For dental use only

Uniclip Plastic Dental Root-Canal Posts in individual packs

Manufactured by:


Maillefer Instruments Holding Sarl, Chemin du Verger, 3, CH-1338, Ballaigues, Switzerland, Tel: +41(0)21 843-97 66, Fax: +41(0) 21 843-9293

1. Description
2. Photos
3. Composition
4. Indication for use
5. Contraindications
6. Precautions
7. Adverse reactions
8. Step-by-step instructions
9. Removing and replacement
10. Disinfection, cleaning and sterilization
11. Transportation
12. Storage and shelf life
13. Environmental requirements
14. Conditions of use
15. Disposal
16. Manufacturer's warranties
17. Authorized Representative

1. Description

The Uniclip posts are plastic non-sterile dental endodontic posts for single use.


Post No./Size	108	208	308	110	210	310
Initial diameter (mm)	0.8	0.8	0.8	1.0	1.0	1.0
Head length (mm)	5	5.5	6	5	5.5	6
Body length (mm)	7.46	9.2	11	7.6	9.3	11
Color	white	yellow	red	blue	green	black

Pack size:

Uniclip plastic posts in individual packs x 100 pcs.; sizes 108, 208, 308, 110, 210, 310.

2. Photos

Uniclip plastic posts; sizes from left to right: 108, 208, 308, 110, 210, 310

3. Composition

Plastic – polystyrene PS (PG-33).

4. Indication for use

The Uniclip posts are used for coronal reconstruction of the root filled tooth.

5. Contraindications

None.


6. Precautions





The post must be cleaned with alcohol before insertion in the canal.

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8. Step-by-step instructions

	Choose the correct root post by placing the transparent drawing (with full size outlines of the different sized root posts) over the radiograph. Note the color of the indicated post and select a drill set of the same color. The length of the tapered section of the root canal should be $\frac{2}{3}$ of the total length. For curved canals, select a root post with a tapered section no longer than the straight section of the coronal portion of the canal.
	Phase 1 Preparatory drilling Using the special penetration drill (REF C 0212), work at low speed (800-1200 rpm) but high torque. The full length of the cutting section should penetrate the canal.
	Intermediate drilling If necessary, use the special precalibrating drill (REF C 0214). This drill is correct for all lengths of root post with a diameter of 1 mm at the tip.
	Phase 2 Precision drilling Prepare the canal to its final size using the special calibrating drill. Employ the same method as for the special penetration drill. The finished preparation will be 0.2 mm longer than the root post, which will prevent over-stressing during tightening.

	Phase 3 Inserting the plastic burn-out post Try each root post in its correctly sized canal, but without pushing in completely. The position of the tapered section of the root post should correspond to the final position of the cutting section of the precision calibrating drill. Push the root post in completely. You are now ready for the final test before taking the impression.
	Phase 4 Before taking the impression Check that the heads of all posts cannot interfere with each other during insertion. The impression tray must not touch the head of the post at any time.
	Phase 5 Taking the impression Inject the impression material into the retention boxes and around the posts, covering them completely.
	Phase 6 Trying in the post before luting Cut a flat all the way along the post, extending it by a groove running under the base and along the side of the retaining box. This allows excess cement to escape and ensure perfect sealing of the built up tooth.
	Reconstruction

	Build up the tooth in the usual way.
--	--------------------------------------

9. Removing and replacement

These procedures are not suitable for plastic posts.

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Step-by-step procedure for sterilization of Unclip posts:

Operation	Operating mode	Warning
DISINFECTION	Immersion in NaOCl solution (at least 2.5%) during 5 min at room temperature.	Do not use disinfecting solutions containing phenol or any products which are not compatible with the Unclip posts.

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MAILLEFER INSTRUMENTS HOLDING Sarl

A handwritten signature in blue ink, appearing to read "Grégory THOMAS".

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Directions for Use

For dental use only

Unimetric Titanium Posts in individual packs

Manufactured by:


**Maillefer Instruments Holding Sarl, Chemin du Verger, 3, CH-1338, Ballaigues,
Switzerland, Tel: +41(0)21 843-97 66, Fax: +41(0) 21 843-9293**

1. Description
2. Photos
3. Composition
4. Indication for use
5. Contraindications
6. Warnings
7. Precautions
8. Adverse reactions
9. Step-by-step instructions
10. Removing and replacement
11. Disinfection, cleaning and sterilization
12. Transportation
13. Storage and shelf life
14. Environmental requirements
15. Conditions of use
16. Disposal
17. Manufacturer's warranties
18. Authorized Representative

1. Description


The Unimetric posts are non-sterile titanium endodontic posts for single use. The Unimetric posts are available in two diameters (0.8 or 1.0 mm) and in various sizes with a long (sizes 108L, 208L, 308L, 110L, 210L, 310L) or a short head (sizes 008S, 108S, 208S, 308S, 110S, 210S, 310S).

Unimetric 0.8mm: initial diameter 0.8 mm, taper: 2.65°

Post No.	008S	108S	108L	208S	208L	308S	308L
Size							
Head length (mm)	2	2	2.5	2.25	2.85	2.5	3.2
Body length (mm)	6.25	7.5	7.5	9.25	9.25	11	11
Color	Purple	White	White	Yellow	Yellow	Red	Red

Unimetric 1.0mm: initial diameter 1.0 mm, taper: 2.65°

Post	110S	110L	210S	210L	310S	310L
Size						
Head length (mm)	2	2.5	2.25	2.85	2.5	3.2
Body length (mm)	7.5	7.5	9.25	9.25	11	11
Color	Blue	Blue	Green	Green	Black	Black

Unimetric pack size

- Unimetric titanium posts in individual packs x 25 pcs.; diameter 0.8 mm or 1.0 mm; sizes 008S, 108S, 108L, 208S, 208L, 308S, 308L, 110S, 110L, 210S, 210L, 310S, 310L.

Accessories:

- 1. Unimetric 0.8 mm Tray for sterilization.
- 2. Unimetric 1.0 mm Tray for sterilization.

2. Photos

	Unimetric titanium posts with a long (upper) or a short (lower) head. On the photo 110L (upper) and 110S (lower), diameter 1.0mm.
	Unimetric 0.8 mm Tray for sterilization
	Unimetric 1.0 mm Tray for sterilization

3. Composition

Posts: Titanium DIN 3.7165.

Trays: Polyphenylene Sulfide (PPS, Ryton R4), polystyrenecover is made of Polyetherimide (PEI, Ultem 1010).

4. Indication for use

Coronal restorations reinforced with a metallic post. Restorations on non-vital teeth requiring a post.






5. Contraindications
None known.

6. Warnings
None known.

7. Precautions
The post must be cleaned with alcohol before insertion in the canal. Avoid touching the post with your fingers after cleaning.

8. Adverse reactions
No adverse reactions have been reported.

9. Step-by-step instructions

	Take a root post of the estimated correct size from the Unimetric Set and compare it with the X-ray photo. The ideal size is the one which, when the head is level with the floor of the pulp cavity, has the tip of its screw thread 2/3 down the root canal. Choose a long or a short head, according to the final level of the reconstructed tooth. Note the color of the root post selected, and use a drill set of the same color.
	Using the penetration drill, work at low speed (800-1200 rpm). The full length of the cutting section should penetrate the canal. The final position of the head is indicated by either of the two rings on the drill: one representing a long head, the other a short head. Drill to a depth which will give a suitable final position for the selected root post.
	If necessary, use the special precalibrating drill to precalibrate the root canal. While drilling, continue until a resistance is felt.
	Give the canal its final size by using the precision calibrating drill and the same method as for the penetration drill. The finished hole will be 2/10 longer than the root post, which will prevent oversteering when tightening.
	Try each root post in its correctly sized canal, but without pushing in completely. Check that the heads of all posts cannot interfere with each other during insertion. If any post is too high, repeat operations no 2 & no 4.

 6  REF A 0022	Fill each finished root canal with zinc oxyphosphate cement, using a Lentulo.
 7  REF C 0201	Insert each root post using tweezers or the wrench. Do not screw tight! Also, do not simply pull the wrench clear, instead tilt sideways while maintaining the axial pressure. Wait for the cement to harden, then remove any excess.

10. Removing and replacement

Unimetric titanium posts could be removed or replaced if necessary.

The procedure for removing a metal post and its replacement is a standard procedure that is part of the common dentist training. The removing could be carried out using ultrasonic tips, for example Start-X # 4 (Maillefer, Switzerland). Using of operating microscope is strongly recommended for such procedure.

General steps of the removal:

- Dentist removes crowns, inlays, old fillings and other artificial constructs.
- Dentist disrupts the cement that holds the post in the tooth using an appropriate ultrasonic tip.
- Dentist pull the post with forceps.

11. Disinfection, cleaning and sterilization

Foreword: For hygiene and sanitary safety purposes, all instruments must be cleaned, disinfected and sterilized before each usage to prevent any contamination. This concerns the first use as well as the subsequent ones.

General recommendations:

- Use only disinfecting solutions which are approved for their efficacy and in accordance with the DFU of the disinfecting solution manufacturer. For all metal instruments, it is recommended to use anticorrosion disinfecting and cleaning agents.
- For your own safety, please wear personal protective equipment (gloves, glasses, mask).
- The user is responsible for the sterility of the product for the first cycle and each further usage as well as for the usage of damaged or dirty instruments where applicable after sterilization.
- The appearance of defects such as cracks, deformations (bent, twisted), corrosion, loss of color coding or marking, indicates that the devices are not able to fulfill the intended use with the required safety level.
- Use quality water in accordance with local regulations especially for the last rinsing step or with a washer-disinfector.

Step-by-step procedure for sterilization of implantable devices (posts) and accessories:

Operation	Operating mode	Warning
Automated cleaning with washer-disinfector	• Place the devices in a kit, support or container to avoid any contact between instruments or posts.	• Discard any instruments with large obvious defects (broken, bent or twisted). • Avoid any contact between instruments or posts when

	Put them in the washer-disinfector (at least 5 min at 90 °C).	placing kits, supports or container in the washer-disinfector. - Follow instructions and observe concentrations given by the manufacturer (see also general recommendations). - Use only approved washer-disinfector according to EN ISO 15883, maintain and calibrate it regularly.
OR		
Manual cleaning or assisted by an ultrasonic device	- Place the devices in a kit, support or container to avoid any contact between instruments or posts. - Immerse in the disinfecting solution with cleaning properties, assisted by an ultrasonic device if suitable.	- No visible impurities should be observed on the instruments. - Discard any instruments with large obvious defects (broken, bent or twisted). - Follow instructions and observe concentrations given by the manufacturer (see also general recommendations). - The disinfecting solution should be free of the following corrosion inhibitors: aldehydes and di- or triethanolamines.
Rinsing	Abundant rinsing (at least 1 min).	- Use quality water in accordance with local regulations. - If a disinfecting solution contains a corrosion inhibitor, it is recommended to rinse the instruments just before the autoclaving. - Dry on a single use non-woven cloth, or with a drying machine or filtered compressed air.
Inspection	- Inspect devices and sort out those with defects. - Assemble the devices (if applicable).	- Dirty instruments must be cleaned and disinfected again. - Discard instruments which show any deformations (bent, twisted), damages (broken, corroded) or defects (loss of color coding or marking) affecting the resistance, the safety or the performance of the instrument or posts.
Packaging	Place the devices in a kit, support or container to avoid any contact between instruments or posts and pack the devices in sterilization pouches.	- Avoid any contact between instruments or posts. Place the devices in kits, supports or containers. - Check the validity period of the sterilization pouch. - Use packaging which are

		resistant up to a temperature of 141°C and in accordance with EN ISO 11607.
Sterilization	• Steam sterilization at 134 °C during 18 min.	• The instruments, posts and the plastic supports must be sterilized according to the packaging labeling. • Use only autoclaves according to the requirements of EN 13060, EN 285. • Use a validated sterilization procedure according ISO 17665. • Follow the maintenance procedure for the autoclave device given by the manufacturer. • Use only this recommended sterilization procedure. • Control the efficiency (packaging integrity, no humidity, color change of sterilization indicators, physico-chemical integrators, digital records of cycles parameters). • Ensure traceability of procedure records.
Storage	• After sterilization, keep devices in sterilization packaging in a dry and clean environment.	• Sterility cannot be guaranteed if packaging is open, damaged or wet. • Check the packaging and the medical devices before using them (packaging integrity, no humidity and validity period).

12. Transportation

No special transportation conditions are required. The devices should be transported according to the rules applicable to each transport mode.

13. Storage and shelf life

The shelf life is unlimited. No special storage conditions are required. After sterilization, keep devices sealed in sterilization pouches in a dry and clean environment.

14. Environmental requirements applicable to use of product

The product does not have an adverse impact on the environment during its life cycle.

15. Conditions of use

Use at: +10 - +35. The devices are not subject to repair. No special maintenance is required.

16. Disposal

Dispose in accordance with the local regulations and the policies and procedures of your medical facility.

17. Manufacturer's warranties

The manufacturer is not responsible and shall not indemnify for any possible damage or accidents caused by:

- Use of instruments unrelated to the system and instruments which can affect normal operation.
- Failure to follow the instructions for use.

The user is responsible for testing of materials for fitness and use for any purpose not expressly specified in the instructions!

For product quality issues, please contact the Authorized Manufacturer's Representative.

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Directions for Use

For dental use only

Radix-Anker-Standard and Radix-Anker-Long Titanium Posts in individual packs

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1. Description

Radix-Anker-Standard and Radix-Anker-Long posts are non-sterile titanium endodontic posts for single use. The Radix-Anker-Long posts have a longer head as compared to the Radix-Anker-Standard.

Radix-Anker-Standard post sizes:

Post No.	1	2	3	3B
Size				
Head length, mm	2.75	3.85	4.95	4.95
Head diameter, mm	2.5	3.0	4.0	3.5
Body length, mm	6	7	9	9
Initial diameter, mm	1.15	1.35	1.6	1.6

Radix-Anker-Long post sizes:

Post No.	1	2	3
Size			
Head length, mm	4.5	6	7.5
Head diameter, mm	2.0	2.5	3.0
Body length	6	7	9
Initial diameter, mm	1.15	1.35	1.6

Pack size

- 1. Radix-Anker-Standard titanium posts in individual packs x 12 pcs.; sizes 1, 2, 3, 3B.2.
- 2. Radix-Anker-Long titanium posts in individual packs x 12 pcs.; sizes 1, 2, 3.

Accessories:

- 1. Radix-Anker-Standard Tray for sterilization.
- 2. Radix-Anker-Long Tray for sterilization.

2. Photos

	Radix-Anker-Standard titanium posts
---	-------------------------------------

	Radix-Anker-Long titanium posts
	Radix-Anker-Standard Tray for sterilization
	Radix-Anker-Long Tray for sterilization

3. Composition

Posts: titanium DIN 3.7165.
 Trays: Polyphenylene Sulfide (PPS, Ryton R4), polystyrenecover is made of Polyetherimide (PEI, Ultem 1010).

4. Indication for use

The Radix-Anker posts are used to support the coronal restoration in case of insufficient residual tooth structure (< 4 mm). This product should be used only in hospital environments, clinics or dental offices by qualified dental personnel.

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6. Warnings

To obtain improved distribution of forces on the dental reconstruction, natural dentine must be present over 2 or 3 mm all around the pivot.

7. Precautions

The post must be cleaned with alcohol before insertion in the canal. Avoid touching the posts with your fingers after cleaning. Dentsply Maillefer recommends the use of a rubber dam.



8. Adverse reactions

In the present technical state, no adverse reactions have been reported so far.

9. Step-by-step instructions

The following rules must be respected in preparing the anchorage cavity and in seating and inserting the Radix-Anker

1 	The crown part of the canal will remain open or will be opened by manual instruments down to penetration depth of the post.
2 	The canal cavity is initially drilled with the non-cutting tip penetration drill C0220.
3 	The seating for the Radix-Anker head is drilled with the root facer C0217. (The stability of Radix-Anker is thus ensured)

4 	Precision drilling is performed manually with the spiral bur C0216 and the mandrel C0219.
5 	Using the gauge C0191 the preparation for the Radix-Anker can be checked as well as its position. The cutting of threads in the root cavity will be realized by screwing the post in place alternately by a half-turn clockwise and a quarter-turn anti-clockwise.
6 	Unscrew the post to insert the sealing cement in the root canal and on the post itself.
7 	When the anchor post is reinserted, it must first be turned to the left until it engages with the pre-cut thread, and only then screwed to the right.
8 	Clean the head and root section, then apply a dentine adhesive. The composite material, resin or amalgam is applied for crown restorations, temporary crowns and fillings.

10. Removing and replacement

Radix Anker titanium posts could be removed or replaced if necessary. The procedure for removing a metal post and its replacement is a standard procedure that is a part of the common dentist training. The removing could be carried out using ultrasonic tips, for example Start-X # 4 (Maillefer, Switzerland). Using of operating microscope is strongly recommended for such procedure.

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Operation	Operating mode	Warning
Automated cleaning with washer-disinfector	• Place the devices in a kit, support or container to avoid any contact between instruments or posts. • Put them in the washer-disinfector (at least 5 min at 90 °C).	• Discard any instruments with large obvious defects (broken, bent or twisted). • Avoid any contact between instruments or posts when placing kits, supports or container in the washer-disinfector. • Follow instructions and observe concentrations given by the manufacturer (see also general recommendations).

		• Use only approved washer-disinfector according to EN ISO 15883, maintain and calibrate it regularly.
OR		
Manual cleaning or assisted by an ultrasonic device	• Place the devices in a kit, support or container to avoid any contact between instruments or posts. • Immerse in the disinfecting solution with cleaning properties, assisted by an ultrasonic device if suitable.	• No visible impurities should be observed on the instruments. • Discard any instruments with large obvious defects (broken, bent or twisted). • Follow instructions and observe concentrations given by the manufacturer (see also general recommendations). • The disinfecting solution should be free of the following corrosion inhibitors: aldehydes and di- or triethanolamines.
Rinsing	• Abundant rinsing (at least 1 min).	• Use quality water in accordance with local regulations. • If a disinfecting solution contains a corrosion inhibitor, it is recommended to rinse the instruments just before the autoclaving. • Dry on a single use non-woven cloth, or with a drying machine or filtered compressed air.
Inspection	• Inspect devices and sort out those with defects. • Assemble the devices (if applicable).	• Dirty instruments must be cleaned and disinfected again. • Discard instruments which show any deformations (bent, twisted), damages (broken, corroded) or defects (loss of color coding or marking) affecting the resistance, the safety or the performance of the instrument or posts.
Packaging	• Place the devices in a kit, support or container to avoid any contact between instruments or posts and pack the devices in sterilization pouches.	• Avoid any contact between instruments or posts. Place the devices in kits, supports or containers. • Check the validity period of the sterilization pouch. • Use packaging which are resistant up to a temperature of 141°C and in accordance with EN ISO 11607.
Sterilization	• Steam sterilization at 134 °C during 18 min.	• The instruments, posts and the plastic supports must be sterilized according to the packaging labeling. • Use only autoclaves according to the requirements of EN 13060, EN 285.

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13.Storage and shelf life

The shelf life is unlimited. No special storage conditions are required. After sterilization, keep devices in sterilization pouches in a dry and clean environment.

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MAILLEFER INSTRUMENTS HOLDING Sàrl

A handwritten signature in blue ink, appearing to be "Grégory THOMAS".

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Directions for Use

For dental use only

EasyPost Dental Glass Fiber Endodontic Posts in individual packs

Manufactured by:

Maillefer Instruments Holding Sarl, Chemin du Verger, 3, CH-1338, Ballaigues, Switzerland, Tel: +41(0)21 843-97 66, Fax: +41(0) 21 843-9293


1. Description
2. Photos
3. Composition
4. Indication for use
5. Contraindications
6. Precautions
7. Adverse reactions
8. Step-by-step instructions
9. Removal and replacement
10. Disinfection, cleaning and sterilization
11. Transportation
12. Storage and shelf life
13. Environmental requirements
14. Conditions of use
15. Disposal
16. Manufacturer's warranties
17. Authorized Representative

1. Description

EasyPost are non-sterile dental glass fiber endodontic posts for single use.

EasyPost sizes:

Post No. Size	1	2	3	4	5	6
OD0, mm	0.8	0.8	1.0	1.0	1.3	1.3
OD3, mm	1.35	1.47	1.67	1.83	2.04	2.22
L1, mm	20.0	20.0	20.0	20.0	20.0	20.0
L2, mm	6.00	7.25	7.2	9.00	8.00	10.00
L3, mm	13.77	12.52	12.51	10.71	11.63	9.63

Taper: 2.65°.

Pack size:

EasyPost glass fiber posts in individual packs x 10 pcs.; sizes 1, 2, 3, 4, 5, 6.

2. Photos

	EasyPost glass fiber post, size 1
	EasyPost glass fiber post, size 2
	EasyPost glass fiber post, size 3
	EasyPost glass fiber post, size 4
	EasyPost glass fiber post, size 5
	EasyPost glass fiber post, size 6

3. Composition

Material	Weight content (%)
Unidirectional glass fiber	81.3 ± 1
Bisphenol A epoxy resin	18.2 ± 1
Peroxide catalyst: 2,5-dimethyl-2,5-di (2-ethylhexanoylperoxy) hexane	0.5 ± 0.2
White polyurethane	0.7 ± 0.2

4. Indication for use

The EasyPost is used to support the coronal restoration in case of insufficient residual tooth structure (< 4 mm). This product should be used only in hospital environments, clinics or dental offices by qualified dental personnel.

5. Contraindications

Insufficient residual dentin (at least 2 mm of tooth structure is required around the preparation).

6. Precautions

- The post must be cleaned with alcohol before insertion in the canal.
- Avoid reprocessing the post more than once because numerous reprocessing cycles might reduce the useful life of the device and damage the marking.
- Avoid touching the posts with your fingers after cleaning.
- Shortening the post should be done outside of the mouth.
- Dentsply Maillefer recommends the use of a rubber dam.

7. Adverse reactions

No adverse reactions have been reported.

8. Step-by-step instructions

- 1) Root canal preparation. In case of insufficient residual tooth substance (< 4 mm), the EasyPost is used to support the coronal restoration.
- 2) Predrilling with a Largo Peeso Reamer, Ref. A 0009 no. 1 or no. 2, depending on the diameter (rotation speed 800-1200 rpm).
- 3) Precision drilling with the calibrating drill, Ref. C0601 no. 1, no. 2, no. 3, or no. 4, depending on the case (1000-2000 rpm).
- 4) Rinse and dry the canal.
- 5) Check positioning in mouth. Please keep a dentine collar of at least 2 mm around the preparation.
- 6) Draw the post outside the oral cavity using a diamond-tipped disc.
- 7) Clean and condition the dentine walls with EDTA during 1 minute.
- 8) Rinse with sodium hypochlorite then dry with air and large diameter paper points.
- 9) Clean the post with alcohol.
- 10) Apply adhesive to the post (silane does not need to be applied to the post): it is recommended to use a compomer or composite luting cement such as Dyract Cem in conjunction with Prime&Bond + Self Cure. The post as well as the space and access cavity must be treated with Prime&Bond + Self Cure. Please follow the manufacturer's instructions. The compomer or composite luting cement can be applied in the canal using a rotary paste filler (Lentulo).
- 11) Insert the post in the canal: please apply some pressure to ensure maximum adhesion.
- 12) Remove any excess.
- 13) Crown section: reconstruct a prosthetic false stump. We recommend using a composite core build-up material (e.g. Dentsply Core X).

9. Removal and replacement

Easypost glass fiber posts could be removed or replaced if necessary.

The procedure for removing fiber glass post and its replacement is a standard procedure that is a part of the common dentist training. The removing could be carried out using ultrasonic tips, for example Start-X # 3 or # 4 (Maillefer, Switzerland). Using of operating microscope is strongly recommended for such procedure.

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10. Disinfection, cleaning and sterilization

Foreword: For hygiene and sanitary safety purposes, all products and accessories must be cleaned, disinfected and sterilized before each usage to prevent any contamination. This concerns the first use as well as the subsequent ones.

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Operation	Operating mode	Warning
Automated cleaning with washer-disinfector	• Place the devices in a kit, support or container to avoid any contact between instruments or posts. • Put them in the washer-disinfector (at least 5 min at 90 °C).	• Discard any instruments with large obvious defects (broken, bent or twisted). • Avoid any contact between instruments or posts when placing kits, supports or container in the washer-disinfector. • Follow instructions and observe concentrations given by the manufacturer (see also general recommendations). • Use only approved washer-disinfector according to EN ISO 15883, maintain and calibrate it regularly.

OR		
Manual cleaning or assisted by an ultrasonic device	- Place the devices in a kit, support or container to avoid any contact between instruments or posts. - Immerse in the disinfecting solution with cleaning properties, assisted by an ultrasonic device if suitable.	- No visible impurities should be observed on the instruments. - Discard any instruments with large obvious defects (broken, bent or twisted). - Follow instructions and observe concentrations given by the manufacturer (see also general recommendations). - The disinfecting solution should be free of the following corrosion inhibitors: aldehydes and di- or triethanolamines.
Rinsing	Abundant rinsing (at least 1 min).	- Use quality water in accordance with local regulations. - If a disinfecting solution contains a corrosion inhibitor, it is recommended to rinse the instruments just before the autoclaving. - Dry on a single use non-woven cloth, or with a drying machine or filtered compressed air.
Inspection	- Inspect devices and sort out those with defects. - Assemble the devices (if applicable).	- Dirty instruments must be cleaned and disinfected again. - Discard instruments which show any deformations (bent, twisted), damages (broken, corroded) or defects (loss of color coding or marking) affecting the resistance, the safety or the performance of the instrument or posts.
Packaging	Place the devices in a kit, support or container to avoid any contact between instruments or posts and pack the devices in sterilization pouches.	- Avoid any contact between instruments or posts. Place the devices in kits, supports or containers. - Check the validity period of the sterilization pouch. - Use packaging which are resistant up to a temperature of 141°C and in accordance with EN ISO 11607.
Sterilization	Steam sterilization at 134 °C during 18 min.	- The instruments, posts and the plastic supports must be sterilized according to the packaging labeling. - Use only autoclaves according to the requirements of EN 13060, EN 285.

		• Use a validated sterilization procedure according ISO 17665. • Follow the maintenance procedure for the autoclave device given by the manufacturer. • Use only this recommended sterilization procedure. • Control the efficiency (packaging integrity, no humidity, color change of sterilization indicators, physico-chemical integrators, digital records of cycles parameters). • Ensure traceability of procedure records.
Storage	• After sterilization, keep devices in sterilization packaging in a dry and clean environment.	• Sterility cannot be guaranteed if packaging is open, damaged or wet. • Check the packaging and the medical devices before using them (packaging integrity, no humidity and validity period).

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12. Storage and shelf life

The shelf life is unlimited.

No special storage conditions are required.

Keep away from direct sunlight.

After sterilization, keep devices in sterilization pouches in a dry and clean environment.

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MAKE CHANGES TO THE FOLLOWING PARAGRAPHS OF THE DIRECTIONS FOR USE

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REMOVE FROM THIS SECTION ALL OTHER ITEMS AND LIST OF ACCESSORIES.

2. Photos

REMOVE FROM THIS PARAGRAPH ALL OTHER ITEMS.

3. Composition

Plastic – polystyrene PS (PG-33).

8. Step-by-step instructions

ADD INFORMATION:

Instruments for the posts placement are sold separately.

CHANGE THE TITLE OF PARAGRAPH 9 OF THE INSTRUCTIONS AND REPLACE ITS CONTENTS:

9. Disinfection.

Foreword: For hygiene and sanitary safety purposes, products must be disinfected before usage to prevent any contamination.

Operation	Operating mode	Warning
DISINFECTION	Immersion in NaOCl solution (at least 2.5%) during 5 min at room temperature.	Do not use disinfecting solutions containing phenol or any products which are not compatible with the Uniclip posts.

13. Conditions of use

The devices are regularly operated at temperature: 32 ° to 42 ° C. The devices are not subject to repair. No special maintenance is required.

Ballaigues, June 2, 2015

Gregory THOMAS

A handwritten signature in blue ink, appearing to be "G. Thomas", written over a horizontal line.

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AMEND SECTION:

Unimetric pack size:

1. Unimetric titanium posts in individual packs x 25 pcs.; diameter 0.8 mm or 1.0 mm; sizes 008S, 108S, 108L, 208S, 208L, 308S, 308L, 110S, 110L, 210S, 210L, 310S, 310L.

REMOVE ALL OTHER ITEMS FROM THIS SECTION.


AMEND SECTION:

Accessories:

1. Unimetric 0.8 mm Tray for sterilization.
2. Unimetric 1.0 mm Tray for sterilization.

REMOVE ALL OTHER ITEMS FROM THIS SECTION.

2. Photos

	Unimetric 0.8 mm Tray for sterilization
	Unimetric 1.0 mm Tray for sterilization

REMOVE ALL OTHER ITEMS FROM THIS PARAGRAPH.

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