
January 21, 2014

The attached copy of the Medtronic Arrive Instructions for Use is an extract of the original English and Russian language translations of the Instructions for Use.

Thank you very much.

Sincerely,


Laura Sobrin
Regulatory Affairs Manager, International
Oscor Inc.

US Notary Public

State of Florida
City of Palm Harbor
County of Pinellas

Sworn to and subscribed before me this 21st day of January 2014 by Laura Sobrin, I hereinto set my hand and seal.

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USA

Caution: Federal Law (USA)

restricts this device to sale by or
 on the order of a physician.

4875 Palm Harbor Blvd. Palm Harbor, FL 34683

P: 727.937.2511 | F: 727.934.9835

WWW.OSCOR.COM

The following list includes trademarks or registered trademarks of Medtronic in the United States and possibly in other countries. All other trademarks are the property of their respective owners. ARRIVE, Medtronic

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1 Description

The ARRIVE (10 Fr) Braided Transseptal Sheath offers percutaneous entry and provides a conduit to deliver diagnostic and therapeutic catheters to specific heart chambers and locations. It provides support for positioning and maintaining the position of catheters at specific locations within the heart.


Figure 1. Arrive Braided Transseptal Sheath

The sheath has a hemostasis valve at the hub to minimize blood loss during catheter introduction or exchange, and a sideport for fluid infusion, blood sampling, or pressure monitoring. The sideport has a three-way stopcock attached.

The sheath has a soft tip with a lubricious inner lumen and a radiopaque marker which aids in defining the location of the tip. The dilator is designed to conform to the inner diameter and curve of the sheath and has a tapered tip. Refer to the following table for available sizes:

Model	Inner Diameter	Curve	Overall Length	Usable Length
990061	10 Fr	55 °–120 °	63 cm	61 cm
990079	10 Fr	55 °–120 °	81 cm	79 cm

1.1 Contents of package

The ARRIVE Braided Transseptal Sheath is supplied sterile. The package contains the following items:

- 1 ARRIVE Braided Transseptal Sheath
- 1 dilator
- product documentation

2 Indications for use

The ARRIVE Braided Transseptal Sheath is indicated for the percutaneous introduction of various types of cardiovascular catheters to all heart chambers including the left atrium via transseptal puncture.

3 Contraindications

Contraindications for the ARRIVE Braided Transseptal Sheath include, but are not limited to, the following conditions:

- Known or suspected left atrial myxoma
- Previous intra-atrial septal patch or device placement
- Previous or systemic embolization from the left side of the heart

4 Warnings and precautions

For single use only – This sheath is intended only to be used once for a single patient. Do not reuse, reprocess, or resterilize these devices for purpose of reuse. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or create a risk of contamination of the device that could result in patient injury, illness, or death.

Sterile package inspection – Inspect the sterile packaging and sheath prior to use. If the sterile packaging or sheath is damaged, do not use the sheath. Contact your Medtronic representative.

Expiration date – Check to verify the device is within its expiration date. Do not use if the product date has expired.

Qualified users – This device should be used only by physicians thoroughly trained in percutaneous procedures.

Biohazard disposal – Discard all used devices and sterile components in accordance with hospital procedures.

Fluoroscopy required for sheath placement – Use of fluoroscopy during sheath insertion and placement is strongly advised. Inserting the sheath without fluoroscopy may result in damage to cardiac and vascular structures. Carefully consider potential risks before using the device in pregnant women.

Other catheters, devices, or wires – Avoid sheath entanglement with other catheters, devices, or wires. Such entanglement may necessitate surgical intervention.

Frequent aspiration and flushing – Aspirate and flush the sheath frequently to help minimize the potential for embolic events resulting from the introduction of air or the formation of clot within the sheath.

Sheath handling –

- Use extreme care when manipulating the sheath and/or dilator. Lack of careful attention can result in injury such as perforation or tamponade.
- Do not use excessive force to advance or withdraw the sheath, especially if resistance is encountered.
- Do not use the sheath if it is kinked, damaged, or cannot be straightened.
- Do not at any time preshape or bend the sheath. Bending or kinking can cause damage to the device.

Radiofrequency ablation – When using the sheath in the presence of RF ablation, care must be taken to assure all ablating elements are outside the sheath.

Back-bleeding – Assure the stopcock is in the closed position after flushing to prevent back-bleeding.

Remove catheters and dilators slowly – Remove devices from the sheath slowly. Rapid removal may damage the valve components resulting in blood flow through the valve and cause a vacuum allowing air to enter the sheath.

Sideport aspiration – Aspirate the sideport when withdrawing the catheter, probe, or dilator to remove fibrin deposits that may have accumulated in or on the tip of the sheath.

Sideport infusion – Infusion through the sideport and catheter should be done after all air is removed from the system.

5 Adverse events

Potential adverse events associated with percutaneous procedures include, but are not limited to, the following conditions

- Access site complications
- Air embolus
- Aortic puncture
- Arrhythmias
- Atrial septal defect
- AV fistula formation
- Coronary artery spasm or damage
- Death
- Device entrapment
- Dissection
- Hematoma
- Hemorrhage
- Myocardial infarction
- Nerve damage
- Pacemaker or defibrillator lead displacement
- Perforation or tamponade
- Pericardial effusion
- Pleural effusion
- Pseudoaneurysm
- Pulmonary edema
- Stroke
- Thromboembolic events
- Valve damage
- Vasovagal reaction
- Vessel spasm or trauma

6 Instructions for use

- Following standard percutaneous puncture, insert a guide wire through the introducer needle into the vessel.
- Advance the guide wire to the desired position.
Note: At no time should the guide wire be advanced or withdrawn when resistance is met.
- Hold the guide wire in place and remove the introducer needle. Do not withdraw the guide wire back into the cannula.
- Prior to use flush both dilator and sheath using heparinized saline.
- Wet the dilator shaft with sterile saline solution prior to insertion through the hemostatic valve.

6.1 Preparing and managing the transseptal sheath

- Assemble the dilator and sheath together until the dilator hub locks into the sheath hub.
Note: Any device/component inserted through the hemostatic valve of the sheath should be wetted and placed through the center of the valve to prevent tearing of the seal and leakage.
- Insert a 0.032" (0.035" maximum) guide wire through the vasculature and into the desired heart chamber using standard vascular access techniques.
- Thread the dilator sheath assembly over the guide wire using a slight twisting motion.
- Position the sheath in the right atrium to prepare for the transseptal puncture.
Note: If a transseptal needle is used, it is recommended to use together with a stylet in order to prevent skiving of the inner lumen of the dilator.
- Insert the properly prepared transseptal needle into the sheath dilator assembly.
Warning: Avoid unnecessary catheter exchanges to minimize sheath related embolic events.
- After a successful transseptal puncture advance the sheath into the left atrium. Under fluoroscopic guidance ensure the sheath is in correct position.
- Slowly remove the needle from the dilator.
- Slowly remove the dilator from the sheath.
Note: Remove the dilator slowly to reduce the possibility of creating a vacuum in the sheath. Blood should aspirate freely through the sideport. If not, withdraw the sheath as the sheath tip may be resting against the wall of the atrium or a pulmonary vein.
- Aspirate blood through the side arm and be sure the sheath is clear of air. Perform aspiration or spontaneous bleed back from the transseptal sheath prior to insertion of a catheter into the sheath.
Note: Perform this step following all catheter exchanges, preferably on the right side.
- Attach the sheath sideport to the monitoring line.

6.2 Following successful transseptal puncture

- Introduce the properly prepared catheter through the hemostatic valve into the left atrium.
- After removal of the sheath, use standard technique to achieve hemostasis.

7 Specifications

	Model 990061	Model 990079
Sheath		
Overall length	63 cm	81 cm
Usable length	61 cm	79 cm
Curve profile	55° - 120°	55° - 120°
Inner diameter	10 Fr	10 Fr
Outer diameter	13 Fr	13 Fr
Dilator		
Overall length	68 cm	86 cm
Useable length	66 cm	84 cm
Guide wire compatibility	0.032" (0.035" max)	0.032" (0.035" max)
Transseptal needle compatibility	71 cm	89 cm

8 Medtronic disclaimer of warranty

For complete warranty information, see the accompanying disclaimer of warranty document.

9 Service

Medtronic employs highly trained representatives and engineers located throughout the world to serve you and, upon request, to provide training to qualified hospital personnel in the use of Medtronic products. Medtronic also maintains a professional staff to provide technical consultation to product users. For more information, contact your local Medtronic representative, or call or write Medtronic at the appropriate telephone number or address listed on the back cover.

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© Medtronic, Inc. 2011
Medtronic Australasia Pty Ltd