NIDEK CO., LTD.

34-14 Maehama, Hiroishi-cho, Gamagori, Aichi 443-0038, Japan
Manufacturers, Exporters & Importers of Ophthalmic Instruments, and Opto-Electronics Instruments

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TEL +81-533-67-6611
FAX +81-533-67-6610
URL <http://www.nidek.co.jp>
<http://www.nidek.com>

To whom it may concern


July 14 2015

CERTIFICATE

It is hereby certified that the attached sheets are the copies of:

Operator's Manual for MP-3

NIDEK CO., LTD


Yoneji Mizuno
Senior Manager
Regulatory Affairs Dept.
NIDEK CO., LTD

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MICROPERIMETER MP-3

OPERATOR'S MANUAL

Eye & Health Care

NIDEK CO., LTD.

NIDEK CO., LTD.
(Manufacturer)

: 34-14, Maehama, Hiroishi-cho, Gamagori, Aichi 443-0038, Japan
Telephone: +81-533-67-6611
Facsimile: +81-533-67-6610

NIDEK CO., LTD.
(Tokyo Office)

: 3F Sumitomo Fudosan Hongo Bldg., 3-22-5, Hongo,
Bunkyo-Ku, Tokyo 113-0033, Japan
Telephone: +81-3-5844-2641
Facsimile: +81-3-5844-2642

NIDEK INCORPORATED
(United States Agent)

: 47651 Westinghouse Drive, Fremont, California 94539, U. S. A.
Telephone: +1-510-226-5700
Facsimile: +1-510-226-5750

NIDEK S.A.
(EU Authorized Representative)

: Europarc 13, rue Auguste Perret, 94042 Créteil, France
Telephone: +33-1-49 80 97 97
Facsimile: +33-1-49 80 32 08

MD VISION Ltd.
(Authorized Representative
in Russian Federation)

: 117312, Moscow, Gubkin str. 14, Russia
Telephone: +7-495-989-80-56
Facsimile: +7-495-989-80-56

March 2015
16202--P902-A0
Printed in Japan

© 2015 NIDEK CO., LTD

Before Use

This Operator's Manual contains information necessary for the operation of the NIDEK MICROPERIMETER MP-3.

This manual includes operating procedures, safety precautions, specifications, and information about accessories and maintenance. This manual is necessary for proper use. Especially, the safety precautions and operating procedures must be thoroughly understood prior to the operation of the instrument.

Keep this manual handy for reference.

If you encounter any problems or have questions about the instrument, please contact NIDEK or your authorized distributor.

Table of Contents

1 SAFETY PRECAUTIONS - - - 1

- 1.1 For Safe Use - - - 1
- 1.2 Usage Precautions - - - 2
- 1.3 Labels and Symbols - - - 8

2 BEFORE USE - - - 11

- 2.1 Outline of Device - - - 11
 - 2.1.1 Intended Use - - - 11
 - 2.1.2 Principle - - - 11
- 2.2 Configuration and Functions - - - 12
 - 2.2.1 Device description - - - 12
 - 2.2.2 Home screen - - - 17
 - 2.2.3 PATIENTS screen - - - 19
 - 2.2.4 Test screen - - - 21
 - 2.2.5 EXAMS screen - - - 28
 - 2.2.6 View screen - - - 30
- 2.3 Packed Contents - - - 36
- 2.4 Initial Use - - - 37

3 OPERATING PROCEDURE - - - 41

- 3.1 Operation Flow Chart - - - 41
- 3.2 Preparation for Test - - - 42
- 3.3 Basic Test Procedure - - - 44
 - 3.3.1 Procedures when operation errors occur - - - 47
 - 3.3.2 Shutdown - - - 50
- 3.4 Connecting Other Devices - - - 52
 - 3.4.1 Computer connection procedure - - - 52
- 3.5 Operation Procedures When the Device is Connected Over a LAN - - - 53
 - 3.5.1 LAN connection: Preparation for Test - - - 53
 - 3.5.2 LAN connection: Completion of test - - - 54
 - 3.5.3 Patient registration - - - 55
 - 3.5.4 LAN connection: Microperimetry test - - - 59
 - 3.5.5 Microperimetry test practice mode - - - 71
 - 3.5.6 LAN connection: Retinography - - - 73
 - 3.5.7 LAN connection: Fixation test - - - 75
 - 3.5.8 If any satisfactory color fundus image cannot be captured - - - 78
- 3.6 FollowUp Test - - - 80

4 DEVICE SETTINGS AND MAINTENANCE - - - 83

- 4.1 Troubleshooting - - - 83
 - 4.2 Error Messages and Remedies - - - 84
-

-
- 4.3 Attaching Chinrest Paper - - - 87
 - 4.4 Parameter Settings - - - 88
 - 4.4.1 Setting procedure - - - 88
 - 4.4.2 Exam setting - - - 90
 - 4.4.3 Application setting - - - 96
 - 4.4.4 System settings - - - 97
 - 4.5 Cleaning - - - 100
 - 4.5.1 Cleaning the objective lens - - - 100
 - 4.6 Consumables and Maintenance Parts List - - - 101

5 SPECIFICATIONS AND TECHNICAL INFORMATION - - - 103

- 5.1 Specifications - - - 103
- 5.2 Exam setting: Test configuration settings - - - 105
- 5.3 Glossary - - - 107
- 5.4 EMC (Electromagnetic Compatibility) - - - 108

6 WARRANTY - - - 111

SAFETY PRECAUTIONS

1.1 For Safe Use

1

BEFORE USE, READ THIS MANUAL.

Be sure to read the operator's manual prior to operation of the device to understand the safety precautions and operating procedures thoroughly.

In this manual, a signal word is used to designate the degree or level of safety alerting.

WARNING

Indicates a potentially hazardous situation which, if not avoided, might result in death or serious injury.

CAUTION

Indicates a potentially hazardous situation which, if not avoided, might result in minor or moderate injury or property damage accident.

Even situations indicated by  **WARNING**,  **CAUTION** may result in serious injury under certain conditions.

Safety precautions must be strictly followed at all times.

1.2 Usage Precautions

Safety of LED

CAUTION

- The device is a Class 1 LED product and the LED used for the device is safe in expected use conditions including situations such as looking into the LED using an optical system. However, it is recommended to observe the following precautions when using the device:
 - 1) Do not direct LED beams to human eyes when unnecessary.
 - 2) Do not look into the objective lens (LED beams aperture) for a long period of time for observation.
 - 3) If the device cannot be used properly or any device problem cannot be solved, immediately stop using the device, then contact NIDEK or your authorized distributor.Use of the device outside the scope of this manual may cause unanticipated adverse events or adverse device effects.

Before use

WARNING

- Be sure to use a grounded power outlet.
Electric shock or fire may result from device malfunction or electric leakage.

CAUTION

- Do not use this device for other than the intended purpose.
NIDEK is not responsible for accidents or malfunctions caused by misuse.
- Cautions for safety and operating procedures must be thoroughly understood before using this device. Do not use accessories that are not specified by NIDEK.
Use of the device outside the scope of the specified use may result in unexpected troubles and adverse events.
- Use of the device is limited to doctors or persons qualified by the law of each country.
Use of the device outside the scope of the specified use may result in unexpected troubles and adverse events. NIDEK is not responsible for accidents or malfunctions caused by misuse.
- Never disassemble nor touch the interior of the device or touch its interior.
There is no part within the device that requires servicing by the user.

⚠ CAUTION

- Be sure to install the device where the conditions below can be maintained.

In addition, use the device under the conditions below.

Ambient temperature: 10 to 35°C (50 to 95°F)

Humidity: 30 to 90% (non-condensing)

Atmospheric pressure: 800 to 1060 hPa

A location with low dust

A location not exposed to water

Protected from interference lights

Protected from exposure to strong electromagnetic waves

Level and stable surface free from vibration and bumping.

The room can be darkened to the degree that a newspaper can barely be read

If the device is not installed and used under the above conditions, the reliability of measurement is lowered, and malfunction may result. In addition, there is a possibility of injury if the device receives shock and falls down.

- Avoid installing the device where it is exposed to direct air-conditioning flow.

Changes in temperature may result in condensation inside the device or adversely affect the device performance.

- Be sure to use a (HOSPITAL GRADE) power outlet which meets the power specification requirements.

The device may not perform properly, or malfunction or fire may occur.

- Never use a power strip or extension cable to supply the device with power.

The electrical safety may be lowered.

- Do not use a power cord other than the one provided. In addition, do not connect the provided power cord to any other device.

Failure or fire may result.

- Do not place heavy objects on the power cord.

A damaged power cord may cause fire or electric shock.

- Before connecting the cable, turn OFF the power switch and disconnect the power cord from the power outlet.

Malfunction of the device may result.

- Install the device so that the outlet that the mains plug is inserted into is easily accessible during use. In addition, ensure that the power cord can be disconnected without the use of any tool.


Failure to do so may interfere with disconnecting the power from the input power source in case of an abnormality.

- Insert the plug properly into the connector with the correct orientation according to the indication. Do not apply undue force to make the connections.

- The device main body should be carried by two persons holding it at positions (A) and (B) (both right and left sides).

Hold it by the bottom of the base unit. Avoid lifting it by the forehead rest or the main unit.

If only one person carries the device, or areas other than the base are used for lifting and the device falls, there is a possibility of injury or malfunction.

During use

CAUTION

- Before use, perform visual and operation checks. If abnormal conditions are encountered, stop using the device.
If the device is used under abnormal conditions, intended results may not be achieved. Also unanticipated malfunctions or health hazards may occur due to improper diagnosis.
- Before and after use of the device, and before examining each patient, clean the chinrest and forehead rest with clean gauze or absorbent cotton. If necessary, dampen a cloth with rubbing alcohol and gently wipe them off.
If the chinrest paper is used, remove a single sheet from the stack of sheets.
- Take care not to place hands or fingers under (or between) the moving parts (image capturing unit, main unit, and chinrest). Be sure to also caution patients.
Hands or fingers may be pinched and may result in injury.
- Keep the objective lens free of fingerprints and dust. Also confirm that it is not dirty before use.
The quality of the captured image may be compromised.
- In the event of smoke or strange odors, immediately turn off the device and disconnect the power plug from the power outlet. After confirming that the device no longer generates smoke, contact NIDEK or your authorized distributor.
Use of the device under such abnormal conditions may cause fire or electric shock. In case of fire, use a dry chemical (ABC) extinguisher to extinguish the fire.
- Before performing the test, explain the purpose and method sufficiently to patients.
- When bringing the measuring unit near to the patient's face or moving it right and left, take care not to contact the patient's face.
Injury may occur.
- Instruct the patient to focus on the fixation target with their eyes wide open. Start the test after confirming that the instruction is properly followed by the patient.
If the instructions are not properly followed, proper test result may not be obtained.
- After the test is complete, and the patient leaving from the device, be sure not to stand up holding the chinrest.
The device main body may topple over resulting in injury.
- Never press the LCD touch-screen with hard point objects such as ball-point pen. Use the accessory touch pen. Never place magnetic objects near the LCD touch-screen.
The device malfunction may occur.
- According to the LCD touch-screen characteristics, some pixels may be constantly lit or not lit. This is not failure so the LCD touch-screen may continue to be used.
- Should the device fail, disconnect the power cord from the power outlet and contact NIDEK or your authorized distributor without touching the interior of the device.
- After a long period of disuse, check for any abnormality before use.

⚠ CAUTION

- If the device is connected to a computer that does not comply with IEC 60601-1 (except one that uses an AC adapter that meets the Class II requirements of IEC 60950-1), supply power to the device and computer through an isolation transformer.


Electric shock may result. Contact NIDEK or your authorized distributor for installing isolation transformers.

- When connecting to peripheral equipment such as a computer through LAN port via a medical facility network, insert or connect an isolation transformer between the medical electrical equipment and network devices (HUB etc.), or the network devices and other electrical equipment.

Electric shock may result. Contact NIDEK or your authorized distributor for installing isolation transformers.

- Use devices that comply with IEC 60601-1 in the patient environment. If any device that does not comply with IEC 60601-1 is to be used, use an isolation transformer or common protective grounding.

The figure below shows the volume of space (patient environment) in which contact can occur between the patient and any part of the device or between the patient and any other person(s) touching the device.

After use

CAUTION

- When the device is not in use, turn off the power switch and place the dust cover over the device main body.
Dust may affect the image capturing performance.
- Make sure that the power switch is turned off before connecting or disconnecting the power cord to or from the power outlet.
If the power cord is connected or disconnected with the power switch on, device malfunction may occur.
- If the device is not to be used for a long time, disconnect the power cord from the power outlet.
- Maintain the ambient temperature and humidity in the following ranges during transport and storage of the device in the packed condition.
Temperature: -30 to +60°C (-22 to 140°F) (during transport),
-10 to +55°C (14 to 131°F) (during storage)
Humidity: 10 to 95 % (non-condensing)
Atmospheric pressure: 700 to 1060 hPa
A location with low dust
A location not exposed to water
A location not exposed to direct sunlight
- To transport the device, set the device main body to packing mode and pack the device in the packing material in which the device was delivered. In addition, avoid vibration or shock to the device.
Excessive vibration or shock may reduce the device reliability.

Maintenance

CAUTION

- When replacing the xenon flash lamp, contact NIDEK or your authorized distributor.
Extremely high voltage part exists inside of the xenon flash lamp. When opening the lamphouse cover or replacing the lamp, electric shock may result.
For the replacement parts, "O Parts to be maintained by service personnel" (page 101).
- Maintenance of the device must be performed every 6 months.
For maintenance details, contact NIDEK or your distributor. If maintenance cannot be performed by the user, request maintenance from NIDEK or your distributor.
- Only service personnel trained by NIDEK can repair the device.
NIDEK will not be responsible for accidents caused by improper servicing.
- Before performing maintenance, wipe the surface of the device properly for disinfection.
- When the device is sent back to NIDEK for repair or maintenance, wipe the surfaces of the device (especially, the areas that come into contact with the patient) with a clean cloth dampened with rubbing alcohol.
- Do not use the device beyond its service life.
Even with proper maintenance and check, after time, the device reliability or safety may become degraded and fail to achieve the target values.

Disposal

CAUTION

- When disposing of the device and accessories, sort them by material and follow local ordinances and recycling regulations. Follow local governing ordinances and recycling plans regarding disposal or recycling of device components, particularly when disposing of the lithium ion battery, circuit board, plastic parts that contain brominated flame retardant, LCD, or power cord.

It is recommended to entrust the disposal to a designated industrial waste disposal contractor.

Inappropriate disposal may contaminate the environment.

- When disposing of packing materials, sort them by material and follow local ordinances and recycling regulations.

Inappropriate disposal may contaminate the environment.


1

1.3 Labels and Symbols

- Cautionary labels are provided on the device. If the labels are peeling off or the characters are wearing off, contact NIDEK or your authorized distributor.

	Indicates that the operator is advised to refer to the related instructions in the operator's manual.
	Indicates that the degree of protection against electric shock is of a Type B Applied Part. The applied parts are the forehead rest, chinrest, and response button. → "2.2.1 Device description" (page 12)
	Indicates that when the switch is pressed to this symbol side, power is not supplied to the device.
	Indicates that when the switch is pressed to this symbol side, power is supplied to the device.
	Indicates that the device must be supplied only with alternating current.
	CE Marking Indicates that the product conforms fully to the requirements of the Medical Device Directive (93/42/EEC). This device is classified as a Class IIa product according to the Medical Device Directive.
	Indicates that this product must be disposed of in a separate collection of electrical and electronic equipment in EU.
	Indicates the date of manufacture.
	Indicates the manufacturer.
	Indicates the port to which to connect the response button terminal.
	Indicates the knob adjusting the fundus observation light intensity.

SAFETY PRECAUTIONS: Labels and Symbols

SAFETY PRECAUTIONS: Labels and Symbols

BEFORE USE

2.1 Outline of Device

- THE NIDEK MICROPERIMETER, MP-3 measures the visual sensitivity of a specified area on the fundus and captures color fundus images. The fundus image overlaid with the retinal sensitivity mapping is displayed on the screen for fundus-image-correlated evaluation.

2

2.1.1 Intended Use

This device is indicated for measuring macular sensitivity, fixation stability and the locus of fixation, as well as providing color fundus imaging.

2.1.2 Principle

(1) Visual sensitivity mapping function

Visual sensitivity mapping can be displayed using static perimetry principle.

Stimuli and a fixation target are displayed by the built-in LCD projector. While focusing on the fixation target, the patient presses the response button to indicate that they saw a stimulus projected at the location and light intensity specified by the internal measurement program or the operator. The patient's response signals are computed by the device and measurement results are displayed on the LCD.

(2) Fundus photography function

Auto alignment is performed while observing on the LCD, the front of the patient's eye illuminated by the eye front illumination LED (infrared light). After the alignment is roughly complete, the mode automatically changes to the fundus observation mode.


The fundus of the patient's eye is illuminated by the fundus illumination LED (infrared light). After alignment and focusing are automatically performed on the fundus, white light from the xenon flash lamp is emitted on the fundus. The light reflected from the fundus is captured by the built-in color CCD camera for fundus image capture.

2.2 Configuration and Functions

- The MP-3 can be used over a LAN network with a computer running NAVIS-EX and provided NAVIS-EX Viewer. ➡ "During use" (page 4)

In this manual, the applicable computer is referred to as "PC", and network connection as "LAN connection".

2.2.1 Device description

1 LCD touch-screen

Various operation screens, image capturing screen, and test data are displayed.

Press displayed buttons to operate the device.

When the buttons are difficult to use by finger, use the touch pen.

The 10.4-inch color LCD monitor can be adjusted by pulling the lower part to tilt it at various angles.

If operator uses the device in a standing posture, tilt the screen at a suitable angle.

To return the screen to the original position, raise it as high as it goes (horizontal position), then slowly lower it.

When the screen is at the origin position, it is secured with a magnet.

2 Release button

Captures fundus images when the auto shot function is set to off (manual).

3 Joystick knob

Moves the measuring unit up and down. (Clockwise: Up, Counterclockwise: Down)

4 Joystick

Slightly moves the measuring unit left and right or forward and backward.

5 Coarse adjustment control

Roughly moves the main unit left and right or forward and backward.

6 Pilot lamp

Illuminates when the device is turned on (I).

Blinks after the device is turned off and when the device is in Sleep mode*1.

↩ "3.3.2 Shutdown" (page 50), ↩ "O System1/4" (page 97).


7 Measuring unit

8 Main unit

9 Base

2

*1The device automatically enters sleep mode when the device is left idle for a pre-determined period of time to save power. (The screen turns off.) However, the device will not enter sleep mode during testing. The preset times for sleep mode are 5, 10, and 15 minutes. (The factory setting is 15 minutes.)

**10 Forehead rest**

The patient's forehead is rested on the forehead rest to prevent their head from moving.
Clean it before each patient.

11 Objective lens**12 Chinrest****13 Eye front illumination LED**

Infrared LED used for auto alignment during test with illuminating patient's eye front

14 Eye level marker

Used as a guide for the patient's eye level.

15 Response button

The device records that the patient visually recognized a stimulus when the patient presses the response button during the visual sensitivity test.

Instruct patients to press the response button only when they clearly recognize stimuli.

When the button is pressed, a short beep sounds.


16 Response button holder

Attachable by the base built-in magnet. When moving the device, remove the holder.

17 Response button connector

Port to connect the Response button

When transporting the device, disconnect the response button from the port.

2

⚠ CAUTION

- Equipment connected to the analog or digital interfaces must be certified according to the representative appropriate national standards (such as IEC60601-1). Furthermore, all configurations must comply with the device standard IEC60601-1. Anyone who connects additional equipment to the signal input part or signal output part configures a medical system, and is therefore responsible that the system complies with the requirements of the system standard IEC 60601-1. If in doubt, consult NIDEK or your authorized distributor.
- For the external device to be connected over a LAN, use the EMC complaint device (CISPR22, 24)

18 Chinrest up/down buttons (▲ / ▼)

When the chinrest reaches the upper (or lower) limit of its movable range, the upper limit mark  (or ) is displayed on the screen.

The chinrest can be moved up and down after the eye side (right or left) is selected on each test screen.

19 Focus knob

Adjusts the focus of the fundus image. (Clockwise: to the +Diopter direction, Counterclockwise: to the - Diopter direction)

20 Eye front/fundus observation toggle button

Toggles between the eye front observation and fundus observation screens.

21 Fundus observation illumination intensity knob

Adjusts the intensity for fundus observation illumination (infrared light).

(Clockwise: Intensity is increased, Counterclockwise: Intensity is reduced.)

22 Power switch

Press the ON () to turn on power to the device.

Be sure to press the Shutdown button on the Home screen before turning off the device.

23 USB port

USB ports to connect other devices. Two ports are for the USB ports. (USB 2.0 supported)

Prior to use, consult NIDEK or your authorized distributor.

24 LAN port

Port to connect with a PC over a LAN.

The measurement data can be exchanged between the device and PC.

For use of the device, contact NIDEK or your authorized distributor.

25 Power inlet


The removable power cord is connected.

2.2.2 Home screen

The screen to select a test type.

Note

- Depending on the operation status, the test buttons that are not necessary to be used are displayed in gray. The buttons displayed in gray are inactive.


2

1 Test button

Press the button of the desired test.

When the device is connected to a PC over a LAN: When no patient is selected, the patient list screen is displayed.

	STATIC MICROPERIMETRY Microperimetry test (visual sensitivity mapping) Performs the test with the default test configuration. The name of the default test configuration is displayed on the button. Example to the left: Normal  "4.4.2 Exam setting" (page 90)
	STATIC MICROPERIMETRY PRACTICE Microperimetry test practice mode Response button practice function  "3.5.5 Microperimetry test practice mode" (page 71)
	RETINOGRAPHY Fundus photography Performs the test with the default test configuration. The name of the default test configuration is displayed on the button. Example to the left: Retinography  "4.4.2 Exam setting" (page 90)

	FIXATION Fixation test Performs the test with the default test configuration. The name of the default test configuration is displayed on the button. Example to the left: Fixation  "4.4.2 Exam setting" (page 90)
---	--

2 Test configuration area

Press the test configuration area to display the registered test configurations. Select the desired test configuration to be performed.

When the device is connected to a PC over a LAN: When no patient is selected, the patient list screen is displayed.

Test configurations can be added on the setting screen.  "4.4.2 Exam setting" (page 90)

3 Operation buttons (Home screen)

Buttons to display the patient list, set the test parameter, or shut down the device.

	Displays the confirmation message dialog box. Pressing the OK button on the confirmation message dialog box shuts down the device.
 *1	Displays the patient list screen.
 *1	Deselects the patient.
 *1	When a patient is selected: Displays the EXAMS screen. When no patient is selected: Displays the patient list screen.
	Displays the setting screen.

*1: The button becomes enabled only when the device is connected to a PC over a LAN.

2.2.3 PATIENTS screen

The screen to register and select patients.

This screen is displayed only when the device is connected to a PC over a LAN.

Note

- Depending on the operation status, the test buttons that are not necessary to be used are displayed in gray. The buttons displayed in gray are inactive.

1 Patient list

Displays the information of the registered or examined patients.

2 Patient list scroll buttons

Buttons to scroll the patient list.

	Displays the first page (the earliest patient).
	Moves up the list by page (displays the previous page).
	Moves down the list by page (displays the next page).
	Displays the last page (the most recent patient).

3 Operation buttons (Patient list screen)

Buttons to display Home screen, register a new patient, and select a registered patient.

	Displays the Home screen.
--	---------------------------

	Confirms the patient selected (touched) on the patient list. After the confirmation, the screen changes as follows: When the patient list screen is displayed by pressing the Patients button: The Home screen appears. When the patient list screen is displayed by pressing the test button/test configuration without selecting a patient: The Test screen appears.
	Deselects the patient.
	Displays the patient registration screen and registers a new patient. ➡ "3.5.3 Patient registration" (page 55)

4 Patient search button (PATIENTS screen)

Touch the field to display the keyboard button.

Input a patient name or ID on the keyboard screen to search.

The item to be searched is displayed in the patient search field and the applicable patients are displayed in the patient list.


To return to the patient list, press the button again and delete the item to be searched in the keyboard screen.

2.2.4 Test screen

The basic screen for performing tests.

Note

- Depending on the operation status, the test buttons that are not necessary to be used are displayed in gray. The buttons displayed in gray are inactive.

1 Test name

*A: **STATIC MICROPERIMETRY**: Microperimetry test

*B: **RETINOGRAPHY**: Fundus image capture

*C: **FIXATION**: Fixation test

*D: **STATIC MICROPERIMETRY PRACTICE**: Microperimetry test practice

2 Parameters

Displays various parameters for the test.

PARAMETERS	Parameter contents	*A	*B	*C	*D
Configuration	Test configuration name	○	○	○	○
Background	Background luminance	○	○ (4 asb only)	○	○
Fixation Target	Fixation target information (Shape, color, size)	○	○	○	○

BEFORE USE: Configuration and Functions

PARAMETERS	Parameter contents	*A	*B	*C	*D
Pattern	Stimuli arrangement pattern When Pattern type is set to Automatic Stimuli arrangement pattern name display When Pattern type is set to Semi-automatic or Manual The corresponding Semi-automatic / Manual display	○	—	—	○ (Manual only)
Stimulus	Stimulus size	○	—	—	○
Dynamic range	Dynamic range of stimuli luminance	○	—	—	○
Strategy	Threshold strategy	○	—	—	○
Options	Displayed when the following options are set for the selected test configuration  "O Microperimetry test settings" (page 91)  "O Fixation" (page 95)				
	Optional Color Fundus Image:Yes ----Color Fundus	○	—	○	—
	Refinement:Yes ----Refinement	○	—	—	○ (only Yes)
	Recheck:Yes ----Recheck	○	—	—	○ (only Yes)
	Pre-test:Yes ----Pre-test	○	—	—	—

3 Status

Displays the test status or result.

STATUS	Status / Result during the test	*A	*B	*C	*D
Tracking	Status of fundus tracking Success: Success (green) Failing: Failing (red)	○	—	○	○
Trigger	Response button status Released: Released (white) Pressed: Pressed (green)	○	—	—	○
Fixation Stability	Fixation status (fixation stability) Stable: Stable (green) Relatively Unstable: Relatively unstable (orange) Unstable: Unstable (red)  "3 Test information display area" (page 33)	○	—	○	○
Elapsed	Elapsed test time	○	—	○	○
Tracked	Total duration of the test (Excluding the time taken for alignment)	○	—	○	○
Remaining	Remaining time for the fixation test	—	—	○	—
Completion	Test progress Microperimetry test n/N N: Total number of stimuli to be projected n: Number of the completed stimuli Fixation test n%	○	—	○	○
False Positive	Result of the false positive test performed during the test: f / F F: Total number of times that the false positive test has performed. f: Number of times that the patient has pressed the response button during the false positive test. Ex.: 0/3 The test has been performed three times and the patient has not pressed the response button. (green) Ex.: 1/3 The test has been performed three times and the patient has pressed the response button once. (red)	○	—	—	○

4 Patient's eye display

Displays the image of the eye being tested.

5 Stimulus luminance scale

Indicates stimulus luminance levels and their corresponding display color for the microperimetry test.

Not displayed for fixation test and retinography.

6 Message field

Displays operating instructions for the test step by step from start to finish.

7 Operation buttons

Buttons to select the eye to be tested, move the chinrest up and down, set the test parameters, and perform operations such as cancellation of the test.

 Abort	Displays the confirmation message dialog box. Pressing the OK button in the confirmation message dialog box cancels the test and displays the Home screen.
 Right	Eye selection: Right eye
 Left	Eye selection: Left eye
 Up	Moves the chinrest up.
 Down	Moves the chinrest down.
 OK	Proceeds to the next procedure.
 Stop	Displays the confirmation message dialog box. Pressing the OK button in the confirmation message dialog box ends the test and displays the test result confirmation screen.
	On the image deletion confirmation screen: Displays the confirmation message dialog box. Pressing the OK button in the confirmation message dialog box displays the test result confirmation screen.
 Pause	Pauses the test. Pressing the Resume button resumes the test.
 Retake	Returns to the eye front alignment and conduct color fundus image capture.
 Take	After the microperimetry test or fixation test are canceled: Returns to the eye front alignment and captures the fundus image. When "Optional Color Fundus Image" parameter is set to Yes → "O Microperimetry test settings" (page 91) → "O Fixation" (page 95)

 Accept	Only when RETINOGRAPHY is performed: Captured image is sent to NAVIS-EX. When other tests are also performed: Displays the test result confirmation screen.
	On the test results confirmation screen: Sends the test result to NAVIS-EX and move to the View screen.
 Discard	Only when RETINOGRAPHY is performed: Displays the Home screen. When other tests are also performed: Displays the image deletion confirmation screen.
	On the test result confirmation screen: Deletes the test result and displays the Home screen.
 Settings	Displays the setting screen.

2

8 Eye indication

Indicates whether the eye being tested is right or left.

When the right or left is selected, icons are displayed and the icon of the tested eye is displayed in blue.

The indication can be toggled between R/L or OD/OS.

The default setting is OD/OS.

	Right eye: When the right eye is selected, the icon is displayed in blue.
	Left eye: When the left eye is selected, the icon is displayed in blue.

9 Patient information

Displays the ID, name, and age of the selected patient.

10 Auto function buttons

Auto functions for the test.

The functions are performed in the order starting from the top.

 Pupil Alignment	Eye front alignment Auto: Auto alignment Manual: Adjust with the joystick  "O Failure of auto alignment to the eye front" (page 47)
 Retina Alignment	Aligns to the fundus Auto: Auto alignment Manual: Adjust with the joystick  "O Failure of auto alignment to the fundus" (page 48)
 Retina Focusing	Adjusts focus on the retina Auto: Auto focus Manual: Adjusted with the focus knob.
 Position Check	Alignment automatic check Automatically checks the eye and the device position. Functions when the retina alignment is set to Auto.

 IR Control	Adjustment of the fundus observation illumination intensity Auto: Automatic illumination control Manual: The illumination intensity is manually adjusted using the fundus observation illumination intensity knob.
 Autoshot	Autoshot function Auto: Automatic color fundus image capture Manual: Images are captured by the release button. When "Optional Color Fundus Image" is set to No, the Auto shot button is not displayed.  "O Microperimetry test settings" (page 91)  "O Fixation" (page 95)

11 Process status

Indicates the status of each Auto function by color.

Gray *1,*2	Not yet performed
Green blinking icon*1	Being performed
Green*3	Performed The auto function button starts blinking when the button is pressed.
Red	Error is generated during processing.

*1: Each press of the auto function icon toggles between Auto and Manual.

*2: When the auto function name is displayed in gray, the function cannot be performed.

*3: The auto function button starts blinking when the button is pressed allowing each function to be manually performed. Pressing the button again switches to auto.

12 Auto indication

Indicates whether each function is set to Auto or not by color.

Blue	Auto The step is automatically performed.
Black	Manual The auto setting is canceled or cannot be performed.

13 Indicator

Indicates flash intensity, focus (Diopter), and fundus observation illumination intensity.

	Sets and displays the flash intensity for color fundus image capture. To reduce intensity: Press  . To increase intensity: Press  .
	Displays fundus focus (Diopter). - Diopter direction: Press  or turn the focus knob counterclockwise. + Diopter direction: Press  or turn the focus knob clockwise.

	Indicates the fundus observation illumination intensity (infrared light). Adjustable when IR Control is set to manual. To reduce intensity: Turn the fundus observation illumination intensity knob counterclockwise. To increase intensity: Turn the fundus observation illumination intensity knob clockwise.
---	--

14 Fixation target setting button

Pressing the button displays the fixation target settings at the bottom of the screen.

The shape, color, and size of the fixation target can be changed.

2

Shape	 Press and select from the drop-down menu.  "O Microperimetry test settings" (page 91)  "O Retinography" (page 93)  "O Fixation" (page 95)
Color	
Thickness	
Size Radius Size Distance	
 Done	Pressing this button completes the settings and closes the fixation target settings.

2.2.5 EXAMS screen

The screen to display the test results of the selected patients.

This screen is displayed only when the device is connected to a PC over a LAN.


FollowUp can be performed by selecting past test results.

➡ "3.6 FollowUp Test" (page 80)

How to display the EXAMS screen: ➡ "Home screen" (page 17)

Note

- Depending on the operation status, the test buttons that are not necessary to be used are displayed in gray. The buttons displayed in gray are inactive.

1 Patient information

2 Test results

Displays the test results (image and test conditions).

3 Test result selection button

Pressing the button displays the selected test results.

Pressing the button selects or deselects the icon.

Selected: The icon or the outline of the icon is displayed in blue.

Deselected: The icon or the outline of the icon is displayed in white.

	Displays the microperimetry test results.
	Displays the Fundus image capture results.
	Displays the fixation test results.
	Displays the test results of the right eye.
	Displays the test results of the left eye.

2

4 Test list scroll buttons

Buttons to scroll the test list.

	Displays the first page (the first test result).
	Moves up the test results by a page (displays the previous page).
	Moves down the test results by a page (displays the next page).
	Displays the last page (the most recent test result).


5 Operation buttons (EXAMS screen)

Buttons to view test results and perform FollowUp.

	Displays the Home screen.
	After selecting a test result, pressing this button displays the View screen.
	After selecting a test result, pressing this button performs FollowUp under the same conditions as the test result.

2.2.6 View screen

The screen to display the selected test result.
This screen is displayed only when the device is connected to a PC over a LAN.
Performs FollowUp under the same conditions as the displayed test result.






1 Test name



2 Test result image switch button

Pressing the button displays the desired test result image.
*A: **STATIC MICROPERIMETRY**: Microperimetry test
*B: **RETINOGRAPHY**: Fundus image capture
*C: **FIXATION**: Fixation test
Pressing the button selects or deselects the icon.
Selected: The icon or the outline of the icon is displayed in blue.
Deselected: The icon or the outline of the icon is displayed in white.

Button	Explanation	*A	*B	*C
 Color fundus image	Press either button to toggle between the color fundus image and IR fundus image. Selecting one icon deselects the other icon. The icon is not displayed for the data with no color fundus image.	○	—	○
 IR fundus image				
 Fixation target	Chooses whether or not to display the fixation target.	○	○	○

Button	Explanation	*A	*B	*C
 Grid display	Whether or not to display the grid 	○	○	○
 Fixation distribution color map	Whether or not to display the fixation distribution color map or fixation point Fixation distribution color map: High density Purple ↔ Low density Blue			
 Fixation point display	Fixation point display: Distribution is indicated by dots. Only one icon can be selected at a time. 	○	—	○

Button	Explanation	*A	*B	*C
 Fixation circle display	Whether or not to display the fixation circle and BCEA circle Fixation circle display: Fixation Stability result display Percentage of fixation points within the circles with 2° and 4° in diameter. <Display example> Circle at 2° (Φ) Percentage of fixation points 99.7% Circle at 4° (Φ) Percentage of fixation points 100.0% → "3 Test information display area" (page 33) → "5.3 Glossary" (page 107) BCEA display: Instability of the fixation point is indicated by size, area (angle of view), and major axis inclination of ellipses for three levels of standard deviation. <Display example> Ellipse 68.2% (1Std Div) Ellipse that represents instability of fixation point within 68.2% Area 0.2°² Ellipse area unit: (angle of view)² Axes Major 0.3° -minor 0.2° Major: long axis diameter (radius) Minor: short axis diameter (radius) Ellipse 95.4% (2Std Div) Ellipse that represents instability of fixation point within 95.4% Area 1.8°² Axes Major 0.9° -minor 0.7° Ellipse 99.6% (3Std Div) Ellipse that represents instability of fixation point within 99.6% Area 1.8°² Axes Major 0.9° -minor 0.7° Axes slope 0.0° Major axis inclination of ellipses: displayed in the range of -90.0° to 90.0° 0.0 is a horizontal status → "5.3 Glossary" (page 107) Only one icon can be selected at a time.    Circle at 4° (Φ) Percentage of fixation points 100.0%    Ellipse 99.6% (3 Std Dev) Area 1.8°² Axes Major 0.9° - minor 0.7° Axes slope 0.0°	○	—	○
 BCEA display				

Button	Explanation	*A	*B	*C
 Stimuli display	Whether or not to display the stimuli (microperimetry test results) 	○	—	—
 Fixation distribution graph	Whether or not to display the Fixation distribution graph or Fixation time profile graph Fixation time profile graph: [Distance from the fixation center] vs [Instability rate] Distribution in the horizontal and vertical directions are indicated in "X DISPLACEMENT" and "Y DISPLACEMENT". Min.: Lower left, Max.: Upper right, Std.: Standard deviation Fixation time profile graph: "Elapsed time" vs "Distance from the fixation center" Only one icon can be selected at a time.			
 Fixation time profile graph		○	—	○

2

3 Test information display area

PARAMETERS: Parameter settings for each test.

"2.2.4 Test screen" → "Parameters" (page 21)

RESULTS: Test results display (No display for retinography).

"2.2.4 Test screen" → "Status" (page 23)

Fixation Stability: → "5.3 Glossary" (page 107)

Stable (green): 75% or greater of all the fixation points are within the 2°(angle of view) circle.

Rel.Unstable (orange): Relatively Unstable; 75% or greater of all the fixation points exists within the 4°(angle of view) circle and less than 75% within the 2° circle.

Unstable (Red): Less than 75% of all the fixation points exists within the 4°(angle of view) circle.

COMMENT: Comments for patient information (comment can be edited only when the device is connected to the PC via LAN)

In Remarks under the COMMENT field, capture information when a color fundus image is captured is displayed.

Remarks: Tv=64. ISO=200. FL=09

Tv: shutter speed, ISO: ISO sensitivity, FL: flash level

When the color fundus image becomes out of focus during the microperimetry test and the Fixation test, the following message is displayed in the Remarks.

Warning: spherical correction modified for the fundus image acquisition.

4 Stimulus luminance scale

Indicates stimulus luminance levels and their corresponding display color for the microperimetry test. Not displayed for fixation test and retinography.

5 Operation buttons (View screen)

Buttons to view test results and perform FollowUp.

 Home	Home	Displays the Home screen.
 Back	Back	Displays the EXAMS screen.
 FollowUp	FollowUp	Performs FollowUp under the same conditions as the test result.
 Edit	Edit	Color fundus image editing Only available for the test data with color fundus images

6 Color fundus image editing

Pressing the Edit button displays the brightness and contrast adjustment buttons at the bottom of the View screen.

Brightness and contrast level: The level at the image capture is set as 0.

 Brightness	Brightness adjustment range: -20 to +20  : Reduces brightness.  : Increases brightness.	
 Contrast	Contrast adjustment range: -20 to +20  : Reduces contrast.  : Increases contrast.	
 Revert	Reverts the brightness and contrast to the previously saved settings.	
 Reset	Resets the brightness and contrast to 0.	
 Done	Finishes the image editing and saves the current settings. The brightness and contrast adjustment buttons disappear and operation buttons are displayed again.	

Note

- When the brightness and contrast are changed and saved, *** mark is attached to the right of the Color Fundus Image button . When the brightness and contrast is reset to the level at the image capture with the  button, *** mark disappears.

7 Patient information**8 LCD angle confirmation marks**

Indicates whether or not the LCD touch-screen is at an appropriate angle to observe the test result image.

A checkered appearance of the LCD angle confirmation marks indicates that the operator is not viewing the display at the correct angle for appropriate observation. Change the viewing angle so that the marks appear evenly before observing the test results.

 : Checkered pattern	Inappropriate viewing angle
 : Even color	Appropriate viewing angle

9 Magnification slider

Enlarges the displayed test result image.


Pressing and dragging the enlarged image assist the examiner to observe the test results in details.

	Pressing this button enlarges the image.
	Moving up the slider enlarges the image. Moving down the slider reduces the image to its original size.
	Pressing the button reduces the image to its original size.
	Pressing this button displays the image in the original size (fitting the frame).

BEFORE USE: Packed Contents

2.3 Packed Contents

- The following are included in the standard configuration. Check the contents before use.

Part name	Quantity	Appearance
Response button Response button holder	1 unit for each	
Touch pen Pen stand	1 unit for each	
Power cord	1 unit	
Chinrest paper Chinrest paper pins (2 units)	1 set for each	
Cap holder Objective lens cap (The objective lens cap is attached on the objective lens at the time of ship- ment.)	1 unit for each	
Blower Dust cover	1 unit for each	
Operator's manual Quick reference guide	1 volume for each	
"NAVIS-EX" installation CD Operator's manual	1 unit for each	

2.4 Initial Use


- Install the device where the conditions specified in "Before Use" are satisfied.

The following describes procedures to prepare for initial use of the device.

When connecting the device with a PC over a LAN, also see "3.4 Connecting Other Devices" (page 52).

- 1 Install the device on a stable place.
- 2 Attach the cap holder according to the procedure below.


- 1) Tie the straps(*A) on the cap holder(*B) and the objective lens cap(*A) as shown to the right.

2

BEFORE USE: Initial Use


- 2) Adhere the cap holder on the right or left side of the device.
- 3) The cap holder is adhesive-backed. Adhere it after removing the paper backing.
- 4) Place the objective lens cap on the objective lens when the device is not used, or the cap holder when the device is used.

3 Attach the holder for the response button to its receptacle on the side of the base.

Note


- The receptacle for the response button holder is magnetic, so the holder is attachable.


4 Connect the response button connector to the response button port.

Note

- To firmly connect the connector, perform the following.
 - (1) Align the guides (*D) of response button connector and port.
 - (2) Insert the connector straight in.
 - (3) Rotate the male screw of the terminal to completely insert the connector.

5 Connect the power cord to the inlet of the device.

- 6 Confirm that the power to the device is turned off (○). Plug the power cord to the power outlet.

 WARNING

- Be sure to use a grounded power outlet.
Electric shock or fire may result from device malfunction or electric leakage.


- 7 Turn on (I) power to the device.

Wait until the Home screen appears.


When the power is turned on, the initialization screen briefly appears after several tens of seconds then changes to the Home screen.

 Note

- The LCD remains black for about one minute until the Home screen is displayed. At this time the device is initializing. Do not turn off the power switch.

Initial screen

HOME screen

The device is now ready for use.

 Note


- Set the parameters according to the use of the device.  "4.4 Parameter Settings" (page 88)

BEFORE USE: Initial Use

OPERATING PROCEDURE

3.1 Operation Flow Chart

3

3.2 Preparation for Test

1 Turn on (|) power to the device.

Wait until the Home screen appears.

When the power is turned on, the initialization screen briefly appears after several tens of seconds. When the power is turned on, the measuring unit makes a slight movement in all directions to determine the initial position for the measuring unit.

Then the Home screen appears.

Note

- For one minute until the Home screen is displayed, the LCD becomes black. As the device is processing the startup, do not turn off the power off by mistake.

2 Remove the objective lens cap.

3 Perform a check before use.

Before using the device, check the following points:

- No error message appears.
- The objective lens is clean.
- The examiner operates the device with the buttons on the Home screen.

If any abnormalities occur, immediately stop using the device and perform remedies in accordance with "4.1 Troubleshooting" (page 83).

4 Conduct patient preparation.

1) Describe to the patient the purpose and procedures for the test.

Note


- Instruct the patient to focus on the fixation target projected in front of their eye during the test.
If the patient moves their face or eyes frequently, the test results may be affected. In addition, this cause the test to pause, resulting in longer test time.
- Instruct the patient to do the following when performing the microperimetry test.
Hold the Response button with either hand easy to press, and press down on the depressed portion with the thumb.
While focusing on the fixation target, the patient should press the response button when they see a stimulus around the target. (Do not press and hold the response button.)
- A short press produces a short beep and a long press produces a long beep (1.5 seconds). After a long beep sounds, instruct the patient not to press and hold the button.
After a long beep sounds, the button becomes inactive until released.

- 2) Wipe the forehead rest(*A) and chinrest(*B) with a clean cotton swab or gauze dampened with rubbing alcohol.
- 3) For the microperimetry test, instruct the patient to hold the response button.
- 4) Instruct the patient to place their chin fully onto the chinrest and lightly press their forehead against the forehead rest.

The chinrest height cannot be adjusted before the test selection.

If the chinrest height is not suitable and difficult to put on, conduct steps 1 to 4 in "3.3 Basic Test Procedure" (page 44) and adjust the chinrest height.

Note

- Do not allow extraneous illumination to come into the patient's field of vision.

Testing in an environment that is too bright may affect the test results due to the patient's pupil failing to dilate enough.

3.3 Basic Test Procedure

1 Select the test to be performed on the Home screen.

Details of each test

↪ "3.5.4 LAN connection: Microperimetry test" (page 59)

↪ "3.5.5 Microperimetry test practice mode" (page 71)

↪ "3.5.6 LAN connection: Retinography" (page 73)

↪ "3.5.7 LAN connection: Fixation test" (page 75)

2 Press the or button to select the eye to be tested.

The message below is displayed in the message field.

Select Eye
Press Right or Left to select the patient's eye

The main unit moves to the selected eye side.

CAUTION

- Before pressing the button, instruct the patient to keep their face, hand, or fingers away from the main unit and the measuring unit.

Hands or fingers may be pinched or the patient's face may contact the unit.

3 Press the button.

When the right or left is selected, Eye indication icons are displayed and the icon of the tested eye is displayed in blue. The message field displays the following messages.

Select Eye
Press Right or Left to select the patient's eye


Adjust Chin rest
Press Up or Down or use the chin rest controls to adjust the chin rest. When finished, press OK.


4 Use the chinrest up/down buttons (,) to align the height of the patient's eyes with the eye level marker^(*A).

Align the height while observing the patient.

The patient's eye front image is displayed on the screen.

- 5 Press the **▲Up** or **▼Down** button or chinrest up/down buttons (▲, ▼) to align the center of the alignment guide spots(*C) with the eye level marker(*B).



- 6 Press the **✓OK** button.

Auto alignment to the eye front starts.

When the alignment completes, the process proceeds to the next step.


For procedures for manual operation or when operation errors occur, refer to “3.3.1 Procedures when operation errors occur” (page 47).

Step	Message field	
 Pupil Align- ment	Automatic Pupil Alignment in progress...	Auto alignment to the eye front
		The screen changes to the fundus observa- tion screen.
 Retina Align- ment	Automatic Retina Alignment in progress...	Auto alignment to the fundus

OPERATING PROCEDURE: Basic Test Procedure

Step	Message field	
	↓	
 Retina Focusing	Automatic Focus in progress...	Auto focus to the fundus.
	↓	
	Reference fiber selection in progress...	Optical working dot positioning (automatic)
 Position Check	↓	
	Post Focus Retina Alignment in progress...	Sets positioning to the fundus Auto tracking to fundus alignment Functions when the fundus alignment is set to Auto.
 IR Control	—	Automatic illumination control of the fundus observation illumination (infrared light). Operates after auto (or manual) focusing. The illumination intensity can be manually adjusted using the fundus observation illumination intensity knob.

After the steps above are all complete, the fixation target is displayed on the fundus (*D) observation image.

The message below is displayed in the message field.

Move Fixation Target
 Drag the fixation target to change the field of view. When finished, press OK

- 7** Drag the fixation target on the fundus observation screen to display the desired fundus area to be tested.
- 8** Press the  button.
Start the selected test.
- 9** After the test is complete, check the test result.
- 10** Press the  button.
The test result is displayed on the View screen.
- 11** Press the  button.
The Home screen appears.
- 12** To conduct the next test, start from the (Step 4) in "3.2 Preparation for Test".
To finish the test, see  "3.3.2 Shutdown" (page 50).

3.3.1 Procedures when operation errors occur

This section describes the procedures for each auto function when an operation error occurs before the test.

When an error occurs, the auto setting is changed to the manual operation.

Note

- Failure of auto alignment to the eye front or fundus, may be caused by the following. To prevent or reduce the failure, properly instruct the patient.
 - Patient moved their face.
 - Patient moved their eyes frequently without gazing the fixation target.

○ Failure of auto alignment to the eye front

The message below is displayed in the message field, and auto alignment is canceled.

Manual Pupil Alignment

Use the joystick to align the patient's pupil.

When the solid circle matches the six dots on the live image, press OK.

Pupil alignment failed

Follow the message and perform the alignment in the procedures below.

- 1) Align the six alignment guide spots (*A) around the alignment circle (*B) using the coarse adjustment control, joystick, and joystick knob.

*B	Alignment circle (solid line) Displayed only for the manual alignment. Set the alignment position by aligning the guide spots.
*C	Required pupil diameter mark (dashed line) Displayed only for the manual alignment. The mark indicates the required pupil diameter of 4 mm.




When the main unit and measuring unit exceed the movable range, the limit mark (blue arrows) appears.

Operate the coarse adjustment control, joystick, and joystick knob in the direction indicated by the arrows.

▲▲ / ▼▼	Up / Down
◀◀ / ▶▶	Right / Left
▲▲▲ / ▼▼▼	Forward / Backward


- 2) Press the  button.

The screen changes to the fundus observation screen for auto alignment to the fundus.

○ Retina Alignment Failure of auto alignment to the fundus

The message below is displayed in the message field and auto alignment is canceled.

Follow the message and perform the alignment in the procedures below.

- 1) Move the measuring unit forward or backward using the coarse adjustment control and joystick to the position where the optical working dots can be seen most clearly.

The optical working dots are a pair of two blinking lights.

When the main unit and measuring unit exceed the movable range, the limit mark (blue arrows) appears.


"O Failure of auto alignment to the eye front" (page 47)

Operate the coarse adjustment control, joystick, and joystick knob in the direction indicated by the arrows.

- 2) Operate the coarse adjustment control, joystick, and joystick knob so that the right and left optical working dots (*D) are at symmetrical positions on the optical working dot charts (*E).

- 3) Press the button.


The auto focus to the fundus starts.

○ Retina Focusing Failure of auto focus to the fundus

When the auto focus cannot be performed, the message window appears on the screen.

The auto focus has failed or the refractive power of the eye is outside the auto focus range.

Note

- When the Split bar--Show parameter is set to "No", the focus bar is not displayed.
In this case, manually adjust the focus with the diopter correction lens or pause the test and set the parameter to ON again.


Parameter setting → "O System1/4" (page 97)

● To manually adjust the focus without diopter correction lens

- 1) Press the No button in the message dialog box.
Manually adjust the focus.
 - 2) Perform the following so that the upper and lower lines of the focus split (*F) form a single line.
The focus split is displayed on top of the focus bar (*G).
- Rotate the focus knob.
 - Press ◀ or ▶ of the focusing indicator (◀ ▶).


3) Press the [OK] button.

The fixation target positioning starts.

● To manually adjust the focus with diopter correction lens

- 1) Press the Yes button in the message dialog box.
The focus bar disappears from the fundus screen and the focus split is not displayed.
- 2) Observe the fundus image and adjust the focus as follows:
 - Rotate the focus knob.
 - Press ◀ or ▶ of the focusing indicator (◀ ▶).

3) Press the [OK] button.

The fixation target positioning starts.

3.3.2 Shutdown

○ Normal shut down

- 1 In the Home screen, press the  button.

The shutdown confirmation message appears.

- 2 Press the OK button.

- LCD touch-screen: Blacks out (back light is turned off.)
- The main and measuring units return to their standard positions.

- 3 Wait until the pilot LED starts to blink and the beep sounds.

Note

- Until the device can be turned off, the LCD becomes black. However, as the device is processing the startup, do not turn off the power off by mistake until the LED blinks and a beep sounds.

- 4 Turn off (○) power to the device.

Note

- Three beeps sound for three times when the pilot LED starts blinking. The beeps continue to sound until the power is turned off.

- 5 Place the objective lens cap on the objective lens.

- 6 Clean the forehead rest and chinrest and place the dust cover on the device.

Wipe with a clean cotton swab or gauze dampened with rubbing alcohol.

Keep the device clean for the next use.

Note

- Be sure to always place the dust cover on the device when the device is not used.

○ Finishing the test when the device is transported

When transporting the device, set the device to packing mode. In Packing mode, the measuring unit and chinrest are automatically set to the position specified for transport.

- 1 In the same manner as in "○ Normal shut down", turn off (○) power to the device.


- 2 While pressing and holding the chinrest down button (▼), turn on (I) power to the device.

Continue to hold the chinrest down button until a long beep sounds (about 10 seconds).

When the device is started and enters packing mode, a message in the right screen is displayed and the chinrest and the measurement part are positioned at the bottom.

Message:

"Please wait...
Do not turn off the power."

Note

- While the message "Please wait... Do not turn off the power." is displayed, do not turn off the power.

3 When the messages shown to the right appears, press the Shutdown button on the screen.

- LCD touch-screen: Blacks out (back light is turned off.)

Message:

"This device is now ready to be packed.
Please press the Shutdown button.
Shut down the system.
(After the pilot LED has blinked,
please turn off power.)"

4 Wait until the pilot LED starts to blink.

Note

- Until the device can be turned off, the LCD becomes black. However, as the device is processing the startup, do not turn off the power off by mistake until the LED blinks and a beep sounds.

5 Turn off (O) power to the device.

6 Place the objective lens cap on the objective lens.

7 Remove the response button, power cord, and various connection cables from the device.

8 Clean the forehead rest and chinrest.

Wipe with a clean cotton swab or gauze dampened with rubbing alcohol.

9 Pack the device in the packing material supplied at the time of purchase.

3.4 Connecting Other Devices


- The MP-3 is meant to be used over a LAN network with a PC running the NAVIS-EX and provided MP Viewer for NAVIS-EX (hereinafter, MP Viewer) installed.

For the installation method and use of NAVIS-EX and MP Viewer, refer to the operator's manual of each software.

CAUTION

- If the device is connected by LAN to other devices such as an external PC via a network of the medical facility, do not connect the device to a network that can connect to the Internet.
Be sure to configure the local network only with related devices such as the PC to which NAVIS-EX is installed. NIDEK will not assume responsibility or compensate for damages caused by any virus infection and development.
- When connecting the MP-3 over a LAN, use a network hub.
Otherwise proper communication may not be performed.
- Be sure to turn off power to a device when connecting a LAN cable.
Malfunction may result if a LAN cable is connected with the power turned on.

3.4.1 Computer connection procedure

Vertically insert the cables of each device in the correct orientation.

Remove the LAN cable (Fig. 1) while holding down the locking clip indicated by "*".

Note

- To connect the device to a network (LAN), set parameters of the device and PC after consultation with the network administrator of the facility. Only authorized service personnel are allowed to perform the connection procedures.
- When not set to be connected to LAN (→ "O System 2/4" (page 99)), do not connect the LAN cable to the PC. The operation may become slow.

CAUTION

- When NAVIS-EX is set to IP address mode and the MP-3 is connected over a LAN without using a network hub, turn on the MP-3 first and then the PC.
NAVIS-EX cannot be activated properly.


3.5 Operation Procedures When the Device is Connected Over a LAN

3.5.1 LAN connection: Preparation for Test

1 Turn on the computer.


2 Start up NAVIS-EX on the computer.

The login dialog box appears.

3 Enter a user name and password, then click the Login.

The NAVIS-EX Patient List screen appears on the computer screen.

NAVIS-EX Patient List screen

4 Turn on (I) power to the device.

Wait until the Home screen appears.

When the power is turned on, the initialization screen briefly appears after several tens of seconds then changes to the Home screen.

Note

- For approximately one minute until the Home screen is displayed, the LCD becomes black. However, as the device is processing the startup, do not turn off the power off by mistake.

5 Register / Select a patient ID to be tested from the Home screen.

1) Press the Patients () button.

The patient list screen appears.

2) Press the patient to be tested on the patient list screen.

The applicable patient is outlined in blue.

To register a new patient, press the New () button.  "3.5.3 Patient registration" (page 55)

3) Press the Select () button.

The Home screen appears.

Note

- Pressing the  button may display the message to the right.

The Image List screen of the same patient is opened or the same patient information is edited in NAVIS-EX. Display the Patient List screen in NAVIS-EX.

 This patient can't be selected since it has been locked by an other application.

6 Conduct patient preparation.

 "Preparation for Test" (page 42)  (Step4)

3.5.2 LAN connection: Completion of test

1 In the Home screen, press the button.

The shutdown confirmation message appears.

 Do you really want to shutdown?
OK Cancel

2 Press the OK button.

- LCD touch-screen: Blacks out (back light is turned off).
- The main and measuring units return to their reference positions.

3 Wait until the pilot LED starts to blink and the beep sounds.

Note

- Until the device can be turned off, the LCD becomes black. However, as the device is processing the startup, do not turn off the power off by mistake until the LED blinks and a beep sounds.

4 Turn off (○) power to the device.

Note

- Three beeps sound for three times when the pilot LED starts blinking. The beeps continue to sound until the power is turned off.

5 Shut down NAVIS-EX on the computer.

6 Shut down the computer and turn off power.

7 Place the objective lens cap on the objective lens.

8 Clean the forehead rest and chinrest and place the dust cover on the device.

Wipe with a clean cotton swab or gauze dampened with rubbing alcohol.

Keep the device clean for the next use.

Note

- Be sure to always place the dust cover on the device when the device is not used.

3.5.3 Patient registration

This section describes the patient registration procedures.

- 1 Press the New () button on the patient list screen.


The new patient registration screen appears.

- 2 Enter the patient information and press the Save () button.

The patient information is registered and the patient list screen appears.

To cancel the registration, press the Cancel

() button.


Note

- Once the patient is registered, they cannot be edited or deleted on the device. Edit or delete the data on NAVIS-EX.
- The patient can also be registered with NAVIS-EX. For the method of patient registration, refer to the NAVIS-EX Operator's Manual.
- If the patient information is not displayed on the upper right on the screen, the patient is not selected. Select the patient again.

○ Patient information entry items


Patient ID	Pressing the field displays the keyboard. Use the keyboard for entry.  • <i>Patient ID input: How to use the keyboard</i> (page 56)
First Name	Pressing each field displays the keyboard. Use the keyboard for entry.  • <i>Inputting last name, first name, comment, and history using keyboard</i> (page 57)
Middle Name	
Last Name	
Date of Birth	Pressing the Date of Birth field displays the date entry window.  • <i>Entering the date of birth</i> (page 57)
Race	Pressing the entry field displays a pull-down menu. Select the applicable race.
Sex	☐ Male ☐ Female ☐ N/S: Touch the corresponding field (check mark is placed).

History	There are two ways of entry. Entry using the keyboard Pressing the field displays the keyboard. Use the keyboard for entry.  • <i>Inputting last name, first name, comment, and history using keyboard" (page 57)</i> Selecting from the History list Pressing the Select ... () to the right of the entry field displays the History list. After selecting the disease name, press the OK button.
Comment	Pressing the field displays the keyboard. Use the keyboard for entry.  • <i>Inputting last name, first name, comment, and history using keyboard" (page 57)</i>

● Patient ID input: How to use the keyboard

1) Press the entry field.

The keyboard appears.

2) Use the keyboard for entry.

Shift () key: Switches the characters.

Lowercase characters / Numbers ↔ Uppercase characters / Symbols

( , ) keys: Moves the cursor. (The cursor can be moved by directly pressing the entered character.)

BACKSPACE  key: Deletes the input characters. Move the cursor and delete characters.

 button: Cancels the entry. Pressing this button closes the keyboard and cancels the entry.


3) Press the OK key.

Closes the keyboard and displays the entered information in each field.

- Inputting last name, first name, comment, and history using keyboard

- 1) Press the entry field.

The keyboard appears.

Alphabet lowercase characters / Numbers Alphabet uppercase characters / Symbols

Example: Last Name entry keyboard

3

- 2) Use the keyboard for entry.

When the "A" key is displayed in orange

Input alphabet, numbers, and symbols directly.

Shift () key: Pressing the key switches the characters.

Alphabet lowercase characters / Numbers ↔ Alphabet uppercase characters / Symbols

, keys: Moves the cursor. (The cursor can be moved by directly pressing the entered character.)

BACKSPACE key: Deletes the input characters. Move the cursor and delete characters.

button: Cancels the entry. Pressing this button closes the keyboard and cancels the entry.

- 3) Press the OK key.

Closes the keyboard and displays the entered information in each field.

- Entering the date of birth

- 1) Press the entry field.

The calendar to enter the date of birth is displayed.

- 2) Select the date of birth.

Select in the order of year(*A), month(*B), and day(*C).

To cancel the entry, press the button.

"March 4, 1954" is selected in the following procedures as an example.


- (1) Select the year of birth.

There are two ways to enter.

Select from forward/back

button: back one year

button: forward one year


OPERATING PROCEDURE: Operation Procedures When the Device is Connected Over a LAN

Search by year

Press the year button and select the desired year from the displayed list.

The example below shows how to select the year 1954.

(2) Select the month of birth.

There are two ways to enter.

Select from forward/back


< button: back one month

> button: forward one month

Selecting the month from the list

Press the month button and select the desired month from the displayed list.

The example below shows how to select "May".

3) Select the date of birth.

The calendar disappears and the date of birth is displayed in the Date of Birth field.

3.5.4 LAN connection: Microperimetry test

1 Select the microperimetry test "STATIC MICROPERIMETRY".

- When testing with the default test configuration

1) Press .

The STATIC MICROPERIMETRY screen appears.


- When testing with the registered Configuration

1) Touch the Configuration area.

The registered Configuration appear.

2) Press the desired Configuration.

Use the scroll bar ^(*A) to see all Configuration.
The STATIC MICROPERIMETRY screen appears.

2 Perform  (Step2) to  (Step6) in "Basic Test Procedure" (page 44).

3 Drag the fixation target on the fundus observation screen to display the desired fundus area to be tested.

4 Press the  button.

Perform the auto alignment to the fundus again.

For Manual, see "3.3.1 Procedures when operation errors occur" (page 47).


Step	Message field	
 Retina Alignment	Automatic Retina Alignment in progress...	Auto alignment to the fundus
	↓	
	Reference fiber selection in progress...	Optical working dot positioning (automatic)
	↓	
—	Automatic selection of Tracking Reference Image in progress...	Acquiring reference image for fundus tracking

OPERATING PROCEDURE: Operation Procedures When the Device is Connected Over a LAN

Step	Message field	
	↓	
—	Coordinate system calibration in progress...	Coordinate calibration A beep sounds. Instruct the patient to focus on the fixation target.

- When the patient's eye movement is stable

After the steps above (coordinate calibration) is complete, the Stimuli Pattern^(*E) is displayed on the acquired IR fundus image (reference image for fundus tracking).

The message below is displayed in the message field.

Move Stimuli Pattern
Drag the stimuli pattern to relocated the stimulation area.
When finished, press OK

↓
Proceed to Step 5.

- When the patient's eye movement is not stable

When the patient's eye movement is not stable, stop the steps above (coordinate calibration) and the Stimuli Pattern^(*E) is displayed on the acquired IR fundus image (reference image for fundus tracking).

◀ , ▶ buttons: 1/5 to 5/5 switching

The message below is displayed in the message field.

Select the tracking reference image and press OK.
Press Retake to acquire a new set of images.

Follow the message and perform either of the procedures below.

Select IR fundus image and press the Accept () button.

Conduct coordinate calibration for the selected image as the reference image for the fundus tracking.

Press the Discard () button.

Perform from the auto alignment to the coordinate calibration again.

When the Stimuli Pattern is displayed on the IR fundus image (reference image for fundus tracking) and the following message is displayed in the message field, proceed to Step 5.

Move Stimuli Pattern
Drag the stimuli pattern to relocated the stimulation area.
When finished, press OK

When the Pattern type is set to Automatic, the Stimuli Pattern is displayed as shown above. When the Pattern type is set to Semi-automatic or Manual, the operator sets the stimuli projection range and the stimuli number for each test.

🔗 "O Stimuli pattern (Pattern type): Semi-automatic/Manual" (page 69)

🔗 "O Microperimetry test settings" (page 91)

5 Drag the Stimuli Pattern on the fundus to the position to be tested.

Note

- When a part of Stimuli Pattern is outside the fundus area, the message below is displayed in the message field.

ex.) Eleven points of stimuli are outside the fundus area and cannot be projected.


There are 11 nominally unprojectable stimuli. Do you want to continue?

When the test is continued, the stimuli that cannot be projected is displayed as "×".

6 Inform the patient that the test will start. Instruct them to press the response button when they see a stimulus.

7 Press the button.

The stimuli is projected on the fundus, and the test starts.

*1: Indicates that Position Check is being performed. This procedure is not performed in manual mode.

During the test, the projection status of the stimuli is displayed in the center of the screen. The response button operation status is displayed in the status display on the left side of the screen.

	Pink: Points that have not been tested. The number indicates the luminance. Other colors: Points that have been tested or are being tested. The color corresponds with the color of the stimuli luminance scale. ● : Points that have been projected. ▲ : Points being tested. The luminance of the next projection is increased to the displayed luminance. ▼ : Points being tested. The luminance of the next projection is decreased to the displayed luminance. × : Projection cannot be performed. ➡ (Step5) Note The point being projected is blinking.
	STATUS Trigger Released: Response button is not pressed (white) Pressed: Response button is pressed (green)

- To pause or resume the test

- 1) Press the  button.

The test is paused.

The message below is displayed in the message field.

Exam Paused
Press Resume to restart.

- 2) Press the  button.

The test is resumed.

- To cancel the test

- 1) Press the  button.

The test is aborted.

The cancellation confirmation dialog box appears.

 **Do you really want to stop the exam?**

OK Cancel

- 2) Press the OK button.

The test results until cancellation are displayed.

Note

- When the color fundus image capture function is set to Yes, the  button is displayed under the screen. The color fundus image is captured by pressing the Take button.

- 3) To save the test result, press the  button.

To discard the test result, press the  button.

Accept: Saves the test result and displays the View screen.

The auto tracking is canceled and the main unit and measuring unit return to their reference positions.

Proceed to Step 11.

Discard: The discard confirmation dialog box is displayed.


OK: Discards the test result and displays the Home screen.

Cancel: Returns to the test result screen.

Are you sure you want to discard current examination?

● Procedure for when the position of the eye is largely shifted during the test

The test is canceled and the fundus alignment screen (manual mode) appears. The message below is displayed in the message field.

3

Manual Retina Alignment

Use the joystick to align the patient's retina.

When the two spots are focused and centered on the live image, press OK

Position check failed

- 1) Use the joystick to correct the alignment.

👉 "O Failure of auto alignment to the fundus" (page 48)

- 2) Press the  button.

The Test screen appears and the test continues.

- 3) Press the Retina Alignment button on the right of the screen to set the fundus alignment to Auto. When set to Auto, the right side of the button becomes blue(*A).

- 4) During the test, if the patient's eye deviated largely and the same message as above is displayed again, conduct above steps 1) to 3) to continue the test.

● Finishing the test

When the test finishes, the confirmation dialog box asking whether or not to perform the additional test appears.

Stimulation is completed. Do you want to add or recheck some stimuli?

- 1) When additional stimuli are not necessary, press "No". When additional stimuli are necessary, press "Yes".

No: No additional test is necessary.


When "Optional Color Fundus Image" is set to Yes, proceed to color fundus image capture.

When the parameter is set to No, finish the test and proceed to Step 11.

For "Optional Color Fundus Image" setting, see ➡ "O Failure of auto alignment to the fundus" (page 48).


When "Optional Color Fundus Image" is set to Yes, proceed to color fundus image capture.

When color fundus image capture is automatically performed (➡ "10 Auto function buttons" (page 25))

The retinography preview screen (simple display) is displayed.

Proceed to Step 8.

Retinography preview screen
(simple display)

When color fundus image capture is manually performed

Press the release button by following the message below.

Yes: Additional test is necessary.

The additional test screen appears.

➡ "O Additional test" (page 67)


After additional test, when "Optional Color Fundus Image" is set to Yes, proceed to color fundus image capture.

➡ "O Microperimetry test settings" (page 91)

➡ "O Fixation" (page 95)

The capture procedure is the same as "No": No additional test.

When "Optional Color Fundus Image" is set to No, finish the test and proceed to Step 11.

Additional test screen

- 8 When the captured image is satisfactory, wait until the preview screen appears. When the captured image is not satisfactory, press the Retake () button.

Note


- When the color fundus image cannot be captured satisfactorily, see "3.5.8 If any satisfactory color fundus image cannot be captured" (page 78).

- Wait: The captured image is satisfactory.

The retinography preview screen (detail display) is displayed.

Proceed to Step 9.

The IR fundus image and color fundus image from the test are overlaid and displayed with a checkerboard background.


Retinography preview screen
(detail display)

- Retake: The captured image is not satisfactory.

The fundus image capture is performed again and the retinography preview screen (simple display) is displayed.

- 9 Confirm that the images overlaid in the checked pattern (IR and color fundus images) are properly aligned by observing the blood vessel alignment.

- 10 After checking the alignment, proceed to the next step with the following buttons.

-  [Accept]: Alignment is proper.

The test result confirmation screen with the color fundus image overlaid is displayed.

- 1) To check the image, change the displayed position by dragging the image or use the magnification slider to enlarge the image for observation.

Test results confirmation screen

- 2) After the image is checked, proceed to the next operation with the following buttons.

-  [Accept]: Saves the test result and displays the View screen.

The auto tracking is canceled and the main unit and measuring unit return to their reference positions.

Proceed to Step 17.

-  [Discard]: The discard confirmation dialog box is displayed.

 "To cancel the test" (page 62) Step 3)


-  [Retake]: Capture an additional color fundus image.

[Retake] in  (Step 8): The captured image is not satisfactory.

-  [Manual Registration]: The alignment can be corrected.

The Manual Registration screen appears.

 "O Manual Registration" (page 68)


Manual Registration screen

-  [Stop]: Cancel the test.

The confirmation dialog box to cancel the test appears.

 "To cancel the test" (page 62)

- 11 The test result can be observed in details on the View screen.

 "View screen" (page 30)

- 12 Press the Home  button.

The Home screen appears.

○ Additional test

When the microperimetry test is complete, methods to retest a specified point or to test an additional point are as follows.


In Step 7 (page 61), the confirmation dialog box asking whether or not to add more stimuli appears.

Stimulation is completed. Do you want to add or recheck some stimuli?

1 Click the Yes button.

The additional test screen appears.

The area outlined with the dotted squares on the fundus image is enlarged and displayed in the lower left of the screen.

3

Note

- The zoom in the enlarged screen can be changed with the parameter setting.

High precision pointing view zoom factor: 1x / 2x / 4x—2x is the default.

➡ "4.4.3 Application setting" (page 96)

2 Drag the square(*A) or enlarged screen(*B) on the fundus image so that the center of the square(*C) is aligned with the position desired to be tested.

3 Slightly align the position by pressing ◀, ▶, ▲, or ▼(*D).

4 Select the initial luminance by pressing the stimuli scale.

When the desired luminance is not displayed, press ◀ or ▶(*E) to move the scale.

The selected luminance is outlined in blue(*F).

5 Press the Add vertex **+ Add/Recheck** button.

The added or recheck stimulus is displayed as a point to be projected (in pink) in the center of the square.

6 Repeat Steps 2 to 5 as many times as necessary to the desired number of points.

7 Press the **■** button.

Only the set points are tested.

When the test is complete, the confirmation dialog box displayed above appears again.

○ Manual Registration

When the color fundus is captured after the microperimetry test, this method corrects the alignment between the acquired IR and color fundus image. The manual registration can also be performed in the Fixation Test.

- 1** Press the Manual Registration  button in the retinography preview screen (detail display).

The Manual Registration screen appears.

- 2** For the IR fundus image and color fundus image, drag the square(*A) or enlarged screen(*B) on the fundus image so that the center of the square(*C) is aligned with a mark (Landmark) position of each image.

 "O Additional test" (page 67)

Guideline of the mark (Landmark): junction or crossing point of vessels.

- 3** Press the Add pair  button.

The first mark has been specified.

- 4** Repeat Steps 2 to 3 to specify the second mark.

To respecify the mark:




 [Undo] : Deletes the landmark.

 [Redo] : Displays the squares at the previous position.

- 5** Press the  button.

The retinography preview screen (detail display) appears.

The alignment of the images has been corrected.

○ Stimuli pattern (Pattern type): Semi-automatic/Manual

When microperipetry test is performed with the test configuration whose Pattern type is Semi-automatic or Manual, the operator sets the stimuli projection range and the stimuli number for each test.


● Semi-automatic

Specify the stimuli projection area and the number of stimuli to be projected for the test.

- 1) Perform (Step 1) to (Step 4) in "3.5.4 LAN connection: Microperimetry test" (page 59).

The stimuli projection setting screen is displayed.

The area in the dotted square on the fundus image is enlarged and displayed in the lower left of the screen.

3

Note

- The zoom in the enlarged screen can be changed with the parameter setting.

High precision pointing view zoom factor: 1x / 2x / 4x-----2x is the default.

➡ "4.4.3 Application setting" (page 96)

The message below is displayed in the message field.

Pattern Editing Wizard

Draw the pattern outline. When finished, press Accept pattern area.

- 2) Drag the square(*A) or enlarged screen(*B) on the fundus image so that the center of the square(*C) is aligned with the position to be the stimuli projection area.
- 3) Finely adjust the position by pressing , , , or (*D).
- 4) Press the Add vertex button.
- 5) Repeat Steps 2 to 4 to define the stimuli projection area.


To respecify the pair:

[Undo] : Deletes the added point.

[Redo] : Displays the square at the previous position.

- 6) Press the Accept Pattern Area button.

The stimuli projection area is defined. The stimuli projection points are displayed in the area.

- 7) Increase or decrease the total number of the stimuli to the desired number.

- button: Decrease the total number of stimuli to be projected.
+ button: Increase the total number of stimuli to be projected.

- 8) Select the initial luminance by pressing the stimuli scale.

When the desired luminance is not displayed, press ◀ or ▶ (*E) to move the scale.

The selected luminance is outlined in blue(*F).

- 9) Inform the patient that the test will start, and instruct them to press the response button when they see a stimulus.

- 10) Press the Accept density button.


The stimuli is projected on the fundus, and the test starts. ➡ (Step7) (page 61)

● Manual

Specify the stimuli for performing test.

- 1) Perform (Step 1) to (Step 4) in "3.5.4 LAN connection: Microperimetry test" (page 59).

The area in the dotted square on the fundus image is enlarged and displayed in the lower left of the screen.

Note

- The zoom in the enlarged screen can be changed with the parameter setting.

High precision pointing view zoom factor: 1x / 2x / 4x———2x is the default.

➡ "4.4.3 Application setting" (page 96)

The message below is displayed in the message field.

Pattern Editing

Edit pattern to change the stimulation area. When finished, press OK.


- 2) Drag the square(*A) or enlarged screen(*B) on the fundus image so that the center of the square(*C) is aligned with the position to be tested.

3) Finely adjust the position by pressing , , , or  (*D).

4) Select the stimuli luminance by pressing the stimuli scale.

When the desired luminance is not displayed, press  or  (*E) to move the scale.

The selected luminance is outlined in blue(*F).

5) Press the Add  button.

The added or recheck stimulus is displayed as a point to be projected (in pink) in the center of the square.

6) Repeat Steps 2 to 5 as many times as necessary for the desired number of points.

7) Inform the patient that the test will start, and instruct them to press the response button when they see a stimulus.

8) Press the  button.

The stimuli is projected on the fundus, and the test starts.  (Step7) (page 61)

3

3.5.5 Microperimetry test practice mode

The Microperimetry test is performed by the patient focusing at the fixation target and pressing the response button when they see a stimulus.

This is a function to practice the operation of the response button.

1 Press .

The screen moves to the practice screen (STATIC MICROPERIMETRY PRACTICE).

2 Perform (Step 2) to (Step 4) in "Basic Test Procedure" (page 44).

3 Drag the fixation target on the fundus observation screen to display the desired fundus area to be tested.


4 Press the  button.

Perform the auto alignment to the fundus again.

Step	Message field	
 Retina Alignment	Automatic Retina Alignment in progress...	Auto alignment to the fundus
	↓	
	Reference fiber selection in progress...	Optical working dot positioning (automatic)
	↓	
—	Automatic selection of Tracking Reference Image in progress...	Acquiring reference image for fundus tracking
	↓	

—	Coordinate system calibration in progress...	Coordinate calibration A beep sounds. Instruct the patient to focus on the fixation target.
---	--	--

After the steps above are complete, the stimuli projection setting screen is displayed.

The message below is displayed in the message field.

Pattern Editing

Edit pattern to change the stimulation area. When finished, press OK.

- 5 Drag the square(*A) or enlarged screen(*B) on the fundus image so that the center of the square(*C) is aligned with the position to be tested.

Note

- The zoom in the enlarged screen can be changed with the parameter setting.

High precision pointing view zoom factor: 1x / 2x / 4x——2x is the default.

🔗 "4.4.3 Application setting" (page 96)

- 6 Finely adjust the position by pressing ◀, ▶, ▲, or ▼ (*D).

- 7 Select the stimuli luminance by pressing the stimuli scale.

When the desired luminance is not displayed, press ◀ or ▶ (*E) to move the scale.

The selected luminance is outlined in blue (*F).

- 8 Press the Add + button.

The added or recheck stimulus is displayed as a point to be projected (in pink) in the center of the square.


- 9 Repeat Steps 5 to 8 as many times as necessary for the desired number of points.

- 10 Inform the patient that the test will start, and instruct them to press the response button when they see a stimulus.


- 11 Press the button.

Start the test.

The stimuli are projected with the specified brightness. The brightness is not changed.

	Pink: Points to be projected. The number indicates the luminance.
	Other than pink: Test finished
	● 14: Patient pressed the Response button. (The color corresponds with the color of the stimuli luminance scale.) ○ 0: Patient did not press the Response button. ×: Stimuli are outside the fundus area and cannot be projected.
	The point being projected is blinking.


When the test finishes, the confirmation dialog box asking whether or not to perform the additional test appears.

- 12** When additional stimuli are necessary, press the Yes button. To finish the test, press the No button.


Yes: The screen changes to the screen to the right.

Addition button:

Perform step 5 to step 10.

The finished test points can be specified again to perform the test.

No: Complete the test and return to the Home screen.

3.5.6 LAN connection: Retinography

- 1** Select "RETINOGRAPHY".


Select a test configuration in the same manner as "3.5.4 LAN connection: Microperimetry test" (page 59) (Step 1).


- 2** Perform (Step 2) to (Step 6) in "Basic Test Procedure" (page 44).

- 3** Drag the fixation target on the fundus observation screen to display the fundus area to be captured.


- 4** Press the  button.

Perform the auto alignment to the fundus again and capture a color fundus image with the auto shot function.

After the steps above are complete, the test result confirmation screen appears.

Test results confirmation screen

- 5** To check the captured image, change the displayed position by dragging the image or use the magnification slider to enlarge the image.

Note

- When the color fundus image cannot be captured satisfactorily, see "3.5.8 If any satisfactory color fundus image cannot be captured" (page 78).

- 6** Press the corresponding button.

[Accept]: Saves the test results and displays the View screen.

The auto tracking is canceled and the main unit and measuring unit return to their reference positions.

Proceed to Step 7.

[Discard]: Discards the test result. The discard confirmation dialog box is displayed.

OK: Discards the captured image and displays the Home screen.

Cancel: Displays the test result confirmation screen.

Are you sure you want to discard current examination?

- 7** The test result can be observed in detail on the View screen.

☞ "View screen" (page 30)

- 8** Press the Home  button.

The Home screen appears.

3.5.7 LAN connection: Fixation test

1 Select the fixation test "FIXATION".

Select a test configuration in the same manner as "3.5.4 LAN connection: Microperimetry test" (page 59). ➡ (Step 1)

2 Perform (Step 2) to (Step 6) in "Basic Test Procedure" (page 44).

3 Drag the fixation target on the fundus observation screen to display the desired fundus area to be tested.

4 Press the button.

Perform the following steps and start the fixation test.

Step	Message field	
 Retina Alignment	Automatic Retina Alignment in progress...	Auto alignment to the fundus
	↓	
	Reference fiber selection in progress...	Optical working dot positioning (automatic)
	↓	
—	Automatic selection of Tracking Reference Image in progress...	Acquiring reference image for fundus tracking
	↓	
—	Coordinate system calibration in progress...	Coordinate calibration A beep sounds. Instruct the patient to focus on the fixation target.
	↓	
—	Exam in progress...	Starting the fixation test
	↓	
—	Automatic Alignment Keeping in progress... ^{*1}	Executing the fixation test
	↓	
—	Exam in progress...	

*1: Indicates that Position Check is being performed. This procedure is not performed in manual mode.

● Finishing the test

When the specified test duration expires, the test is complete.


When "Optional Color Fundus Image" is set to Yes, proceed to color fundus image capture.

When the parameter is set to No, finish the test and proceed to Step 6.


For "Optional Color Fundus Image" setting, see ➡ "O Fixation" (page 95).

When "Optional Color Fundus Image" is set to Yes, the procedures to capture color fundus image are as follows.

When color fundus image capture is automatically performed (⇒ "10 Auto function buttons" (page 25))

The retinography preview screen (simple display) is displayed.

Retinography preview screen (simple display)

When color fundus image capture is manually performed

Press the release button by following the message below.

- To pause or resume the test

⇒ "3.5.4 LAN connection: Microperimetry test" (page 59) (Step 7)

- To cancel the test

⇒ "3.5.4 LAN connection: Microperimetry test" (page 59) (Step 7)

- Procedure for when the position of the eye is largely shifted during the test

⇒ "3.5.4 LAN connection: Microperimetry test" (page 59) (Step 7)

5 Perform Steps 8 to 10 in "3.5.4 LAN connection: Microperimetry test" (page 59).

Manual Registration "O Manual Registration" (page 68)

- 6** The test result can be observed in detail on the View screen.




 "2.2.6 View screen" (page 30)

- 7** Press the Home  button.



The Home screen appears.

3.5.8 If any satisfactory color fundus image cannot be captured

According to the alignment between the main body and the patient's eyes, a black shadow may be projected on the color fundus image. Use the table below to improve image capture quality.

Symptom	Cause and Remedy	
Satisfactory image		
Near the center of the image is dark	The patient's pupil may not be dilated enough. Improve the patient's pupil dilation. If the pupil is difficult to dilate, increase the flash intensity for image capture.	
White reflection on the image periphery	The main unit may be shifted from the patient's eye. Occurred on the upper part: move the main unit upward Occurred on the lower part: move the main unit downward Occurred on the right part: move the main unit rightward Occurred on the left part: move the main unit leftward	
White reflection on the lower part of the image.	- The patient's eyelid may be closing. Instruct the patient to open one's eye wide. If the patient cannot open their eye sufficiently, lift the eyelid taking care not to press against the eyeball. - The image may have been captured at the moment of blinking. Ask the patient to refrain from blinking and conduct capture as quickly as possible.	

OPERATING PROCEDURE: Operation Procedures When the Device is Connected Over a LAN

Symptom	Cause and Remedy
A circular or linear faded white shape appears in the image.	- The patient's eyelashes may be obscuring the required pupil diameter mark. Instruct the patient to open one's eye wide. If the patient cannot open their eye sufficiently, lift the eyelid taking care not to press against the eyeball. - The patient's tear may be shedding too much. Capture the image after wiping the patient's tear. 
The center of the image is white. The image becomes white overall.	The objective lens may be smudged. → "4.5.1 Cleaning the objective lens" (page 100) 

3.6 FollowUp Test


The FollowUp test can be performed under the same conditions as the past test results.

The FollowUp test may be performed on the test results for the microperimetry test and fixation test.

The FollowUp test cannot be performed on the fundus image capture.

- 1 Select the desired past test results on the Exam screen to perform FollowUp.
- 2 To check the past test results, press the View  button to display the results on the View screen. Check the test data.
If the checking is not necessary, perform Step 3 on the Exam screen.
- 3 Press the FollowUp  button.
The same Exam screen as the past test is displayed.
- 4 For the microperimetry test, instruct the patient to hold the response button.
- 5 Perform (Step 2) to (Step 5) in "Basic Test Procedure" (page 44).
- 6 Press the  button.

Start alignment and capture a reference image for tracking the fundus (IR fundus image).


After the image acquisition, the acquired image and the reference image for tracking of the past data are overlaid and displayed in an alternately checkered pattern with a checkerboard background.

 (Black and white): Follow-up tracking reference image

Acquired reference image for tracking

 (Purple): Baseline tracking reference image

Reference image for tracking of the past data

FollowUp reference image check screen
(Microperimetry test)

7 Confirm that the overlaid images (reference images for tracking) are properly aligned by observing the alignment of marks such as blood vessels.

8 Perform operation in the same manner as "3.5.4 LAN connection: Microperimetry test". ➡ (Step10)

- [Accept]: Alignment is proper.

Microperimetry test:

Microperimetry test is performed with the same fixation target positions and the same stimuli pattern as the past test data.

Fixation test

Microperimetry test is performed with the same fixation target positions as the past test data.

- [Manual Registration]: The alignment can be corrected.

The Manual Registration screen appears.

➡ "O Manual Registration" (page 68)

9 After the test is complete, the test result can be observed in details on the View screen.

➡ "2.2.6 View screen" (page 30)

For saved test results, a display indicates whether there is a FollowUp test.

History(*A): Follow up of 2014 /10/14 10:33 (1)

YYYY/MM/DDDD Time (Number of times FollowUp test was performed)

The test data on which the FollowUp test was conducted is identified by the "Baseline -" display.

If FollowUp is performed on FollowUp data, the previous FollowUp becomes the Baseline data.

ex) Select the data with perform the FollowUp test.

Baseline data of "Follow up of 2014 /10/14 10:33 (1)" and "Follow up of" 2015 /1/4 11:12 (2)" are the same.

Press the Home button.

The Home screen appears.

OPERATING PROCEDURE: FollowUp Test

DEVICE SETTINGS AND MAINTENANCE

4.1 Troubleshooting

- In the event that the device does not work correctly, attempt to correct the problem according to the following table before contacting NIDEK or your authorized distributor.

Symptom	Remedy
The pilot lamp does not illuminate even when power to the device is turned on.	The power cord may not be plugged into the power outlet. After checking the power inlet and power cable connection, firmly insert to the power outlet again.
The LCD touch-screen suddenly turns off.	If the pilot lamp is blinking, the device is in Sleep mode. Press any button to restore the device.
The main unit cannot be moved to the right or left. The main unit cannot be moved forward or backward.	The movement range of the measuring unit is designed to be limited in order to prevent it from contacting the patient's face. Operate the measuring unit again after moving it back to the operator's side.
The measuring unit cannot be operated by the joystick.	When the selection of the eye to be tested and the chinrest adjustment are not complete, the operation cannot be performed. Be sure to complete those procedures before operating the measuring unit.

If the symptom cannot be corrected by the above remedies, contact NIDEK or your authorized distributor.

4.2 Error Messages and Remedies

- If a response button error or test result data sending error appears on the screen, correct the error according to the instructions provided in the Cause and Remedy column.

When any error other than the response button error or test result data sending error occurs, turn off power to the device and contact NIDEK or your authorized distributor.

For details of turning off the power, see "3.3.2 Shutdown" (page 50).

If the error cannot be remedied, contact NIDEK or your authorized distributor. Prepare in advance the serial numbers of the main body and the PC, error information, and the symptom.

○ Response button error

Error related to the response button. When an error is generated, a message to the right is displayed.

 Please check the response button condition

Response button error generation status	Cause and Remedy
At startup	• The response button may not be properly connected or the response button may be broken. • Turn off the device power and firmly reconnect the response button connector to the response button port. • If the same message appears when the device power is turned on again, turn off the device and contact NIDEK or your authorized distributor.
During test	• Patient may be pressing the response button (approximately 45 seconds or more) • Instruct the patient not to press and hold the response button. Whether the response button is "pressed or released" can be checked in the status display on the left of the test screen.  "3 Status" (page 23) Trigger Released: white Pressed: green When "Released" is confirmed, the test can now be proceeded. (Touching the message deletes the display.) • When a status of pressing and holding the response button continues, the device malfunction is possible. Turn off power to the device and contact NIDEK or your authorized distributor.

Note

- Even if the response button error occurs, the fundus photography "RETINOGRAPHY" and the fixation test "FIXATION" can be performed. To perform the test, touch the message and delete the display beforehand.

○ Test result data sending error

If the test result data fails to be sent to NAVIS-EX (after pressing the Accept button on the test result screen), follow the instructions below.

If the same error message appears repeatedly after conducting the following procedures, turn off the device and contact NIDEK or your authorized distributor.

If the error occurs, a message to the right is displayed.

✕ Exam Saving failed! Data can not be sent to the server.

1 Press the OK button.

The message changes.

Depending on where the communication error occurs, it is possible that no message is displayed. In this case, return to the Home screen and proceed to Step 3.

✕ NAVIS-EX connection failed: please check connection settings.

2 Press the OK button.

The Home screen is displayed.

4

3 Press the Settings Settings button.

The System2/4 screen is displayed.


Pending Exams Management field

Some test data have been saved on this client but not sent to NAVIS-EX sever.

This could be due to a connection error between the device and the server.

To check and save the connection settings, press the Connect and send data button.

Pressing 'Delete data' pending test will be erased and all information will be lost.

4 In accordance with the message, perform the following operation.

- To send the test result data again

1) Check the LAN cable connection between the device and the PC.

2) Press the Connect and send data  (*A) button.

Success of the connection with NAVIS-EX:

Proceed to Step 3).

! The connection to Navis-EX database server has been successfully established!

Failure of the connection with NAVIS-EX:

Failed to connect to the Navis-EX database server: review the connection settings.

If the connection to NAVIS-EX fails repeatedly, press the OK button and then the Home  button. Then return to the Home screen and turn off power to the device and contact NIDEK or your authorized distributor.

3) Press the OK button.


The confirmation message of the test results to be sent is displayed.

<Test result data information>

Exam type

Patient

Eye: left/right eye

4) Press the OK button.

From the System2/4 screen, when the "Data storage field being pended" is not displayed, the test result is sent.]

 Note

- If data fails to be sent to NAVIS-EX, restart the PC running NAVIS-EX and return to step 3) above to perform the operation again.

5) Press the Home button to display the Home screen.

6) Press the Exams button to check the test result data in the Exam screen.

● To cancel the test result data

1) Press the Delete data button.

A message to the right is displayed.

Do you want to delete the pending exams from the client?
Pressing Yes all the locally stored data will be erased.

2) Press the Yes button and then the Home button to return to the Home screen.

4.3 Attaching Chinrest Paper


1 Pull out the two chinrest paper pins from the chinrest.

2 Remove a proper amount of chinrest paper from a new pad.

A new pad of chinrest paper cannot be attached to the chinrest at one time. Be sure to use a pad of 6 mm or less removed from the whole pad.

3 Insert the chinrest paper pins into the chinrest paper.

Insert the chinrest paper pins into the holes on the chinrest paper.

4 Attach the chinrest paper to the chinrest

1) While holding the chinrest paper pins and chinrest paper, insert a chinrest paper pin into one of the two holes on the chinrest.


2) Insert the other chinrest paper pin into the other hole on the chinrest.


4.4 Parameter Settings

4.4.1 Setting procedure

Change the settings of the device on the Settings screen.

- 1 To display the Settings screen, press Settings  in the lower right of the Home screen or the screen before the test starts.
- 2 Press Exams , Application , or System  in the upper left of the Settings screen to display the desired setting items.

Exams	Exam setting: Test configuration settings 
Application	Application setting: Display language, Right / Left eye indication, patient name display, zoom 

System	System settings (1/4 to 4/4) Setting items: 1/4: Observation camera, fundus image camera, beep sound, sleep time, LCD brightness, clock, focus bar, default flash intensity for fundus image capture, and chin-rest initialization 2/4: Network Confirmation items 3/4: Software version 4/4: Software license terms
	

4


3 Change the settings as desired.

When attempting to move another screen, the saving confirmation message appears.

Select save "Yes" or abandon "No".

 There are unsaved changes: would you like to save them?

4.4.2 Exam setting

1 Exam type

Press the "Exam type" field and select the desired test type from the displayed pull-down menu (Retinography, Fixation, Static micropertimetry).

2 Configuration

Press the "Configuration" field and select the desired test configuration from the displayed pull-down menu (registered test configuration list).

The selected test configuration is displayed in each setting field.

Note

- The test configuration named "FACTORY", is the pattern specified by the manufacturer. It cannot be deleted or changed. Only configurations created by the user can be deleted or changed. However, the "FACTORY" test configuration name can be copied to create a new test configuration.

• Procedures to create a new test configuration by copying the test configuration named, "FACTORY." (page 91)

3 Set as default [Set as default] button

The selected test configuration is set as the default. (The test configuration will be displayed on the Test button of the Home screen.)

DEFAULT is displayed as the name of test configuration set as the default.

Ex.: Normal [FACTORY] [DEFAULT]

4 Delete [Delete] button

The selected test configuration is deleted.

5 New [New] button

A new test configuration can be created.

- 1) Press the New button.

The Configuration field becomes blank.

- 2) Select the desired setting for each field.
- 3) Press the Save button.

6 [Save] button

The changed or newly created test configuration is saved and registered.

- Procedures when the settings are changed

- 1) After changing the settings, press the Save button.
The changed settings are saved.

- Procedures when a new test pattern is created

- 1) Press the Save button.

The keyboard appears.

 “● Patient ID input: How to use the keyboard” (page 56)

4

- Procedures to create a new test configuration by copying the test configuration named, "FACTORY".

- 1) Select the test configuration with the desired "FACTORY".
- 2) Change the desired settings with each settings.
- 3) Press the Save button.

The keyboard appears.

 “● Patient ID input: How to use the keyboard” (page 56)

○ Microperimetry test settings

Select from a pull-down menu.

General	
Exam type	Select Static microperimetry
Configuration*1	Select from the registered Configuration list.
Automatism	
Pupil Alignment	Choose whether or not to perform auto alignment to the eye front. Automatic / Manual
Retina Alignment	Choose whether or not to perform auto alignment to the fundus. Automatic / Manual
Retina Focusing	Choose whether or not to perform auto focusing to the fundus. Automatic / Manual

DEVICE SETTINGS AND MAINTENANCE: Parameter Settings

IR Control	Choose whether or not to perform automatic illumination control of the fundus observation illumination (infrared light). Automatic / Manual
Optional Color Fundus Image	Choose whether or not to perform Optional Color Fundus Image. Yes / No No: During the test, the Auto Shot button is not displayed on the test screen.
Auto Shot	Select auto shot function (color fundus image capture) Yes / No
Fixation Target*2	
Shape	Fixation target shape Single Cross ($\oplus$) / Circle ($\bigcirc$) / Four Crosses ($\oplus\oplus\oplus\oplus$)
Color	Fixation target color white / yellow / blue / red
Thickness (width)	Fixation target width unit: angle of view 0.10° / 0.20° / 0.30° / 0.40° / 0.50°
Size / Radius / Size Distance	Fixation target size unit: angle of view Options in the pull-down menu When Single Cross is selected Size: 1.00° to 20.00°(1°increments) When Circle is selected Radius: 0.5°and 1.00° to 20.00°(1°increments) When Four Crosses is selected Size (size of each cross): 1.00° to 20.00°(1°increments) Distance (distances between upper and lower, and right and left): 1.00° to 20.00°(1°increments)
Examination	
Pattern type	Stimuli arrangement Automatic / Semi-automatic / Manual
Pattern name*1	Stimuli arrangement pattern (Only when Pattern type is set to Automatic.) Select from the registered stimuli arrangement patterns.
Strategy*3	Threshold strategy 4-2 / 4-2-1
Duration	Stimulus duration 100ms / 200ms
Refinement*4	Stimuli can be added. (When Pattern type is set to Automatic or Semi-automatic) Yes / No
Recheck*5	The tested stimuli can be projected for recheck. (Only when Refinement is set to Yes.) Yes / No
Pre-test*6	Stimuli projection Pre-test (When Pattern type is set to Automatic or Semi-automatic) Yes / No
Perimetric scale	Microperimetry scale "□MP-1 scales" *7 not checked: 31.4asb-white on white-34dB (standard) Touch "□MP-1 scales" *7 to check and select: 4asb-white on white-20dB / 4asb-white on white-34dB

Stimulus	Stimuli size Goldmann I / Goldmann II / Goldmann III / Goldmann IV / Goldmann V
Attenuation	Initial setting for the stimuli luminance (Only when Pattern type is set to Automatic.) 0 dB to 34 dB (1 dB increments)

*1: When the registered stimuli arrangement pattern is edited or deleted with MP Viewer, the test configuration using the registered pattern is selected and displayed,  mark appears before the name and [OUTDATED] is displayed after the name as shown below.

For the test configuration in this case, the stimuli arrangement pattern before editing or deleting can be used. The stimuli arrangement pattern before editing and deleting cannot be used to create any new configuration.

ex): 

 new pattern [OUTDATED]

*2: The fixation target can be changed in the Test screen.  "14 Fixation target setting button" (page 27)

*3: Threshold strategy

4-2:

If the initial stimulus is not recognized by the patient, the subsequent stimuli at that location are increased by 4 dB in intensity until the patient recognizes the stimulus (first reversal). The subsequent stimuli are reduced by 2 dB in intensity until the patient no longer recognizes the stimulus. (second reversal)

If the initial stimulus is recognized by the patient, the subsequent stimuli at the location are reduced by 4 dB in intensity until the patient no longer recognizes the stimulus (first reversal). The subsequent stimuli are increased by 2 dB in intensity until the patient recognizes the stimulus (second reversal).

4-2-1:

After the second response reversal, stimuli are adjusted by 1 dB until the third response reversal occurs.

*4: The test can be continued by adding stimuli.  "O Additional test" (page 67)

*5: Specify the completed test point to restart the test.  "O Additional test" (page 67)

*6: The stimuli pattern is divided into quarters and one stimulus in each quadrant is projected to determine the initial luminance for the other stimuli in each quadrant.

*7: When the box is checked, the perimetric scale can be selected with the background luminance 4asb, equivalent to NIDEK MP-1.

○ Retinography

Select from a pull-down menu.

General	
Exam type	Retinography selection
Configuration	Select from the registered test configurations list.
Automatism	
Pupil Alignment	Choose whether or not to perform auto alignment to the eye front. Automatic / Manual
Retina Alignment	Choose whether or not to perform auto alignment to the fundus. Automatic / Manual
Retina Focusing	Choose whether or not to perform auto focusing. Automatic / Manual
IR Control	Choose whether or not to perform automatic illumination control of the fundus observation illumination (infrared light). Automatic / Manual
Auto Shot	Choose whether or not to automatically perform color fundus image capture. Yes / No

DEVICE SETTINGS AND MAINTENANCE: Parameter Settings

Fixation Target ^{*1}	
Shape	Fixation target shape Single Cross (+) / Circle (○) / Four Crosses ($\begin{smallmatrix} + & + \\ + & + \end{smallmatrix}$)
Color	Fixation target color white / yellow / blue / red
Thickness (width)	Fixation target width unit: angle of view 0.10° / 0.20° / 0.30° / 0.40° / 0.50°
Size / Radius / Size Distance	Fixation target size unit: angle of view Options in the pull-down menu When Single Cross is selected Size: 1.00° to 20.00°(1°increments) When Circle is selected Radius: 0.5° and 1.00° to 20.00°(1°increments) When Four Crosses is selected Size (size of each cross): 1.00° to 20.00°(1°increments) Distance (distances between upper and lower, and right and left crosses): 1.00° to 20.00°(1°increments)

^{*1}: The fixation target can be changed in the Test screen.  "14 Fixation target setting button" (page 27)

○ Fixation


Select from a pull-down menu.

General	
Exam type	Fixation selection
Configuration	Select from the registered Configuration list.
Automatism	
Pupil Alignment	Choose whether or not to perform auto alignment to the eye front. Automatic / Manual
Retina Alignment	Choose whether or not to perform auto alignment to the fundus. Automatic / Manual
Retina Focusing	Choose whether or not to perform auto focusing to the fundus. Automatic / Manual
IR Control	Choose whether or not to perform automatic illumination control of the fundus observation illumination (infrared light). Automatic / Manual
Optional Color Fundus Image	Choose whether or not to perform Optional Color Fundus Image. Yes / No No: During the test, the Auto Shot button is not displayed on the test screen.
Auto Shot	Select auto shot function (color fundus image capture) Yes / No
Background	
Luminance range	Background luminance "□MP-1 Luminances" *1 not checked: 31.4asb (standard) Touch "□MP-1 Luminances" *1 and check to select: 4asb
Fixation Target*2	
Shape	Fixation target shape Single Cross (+) / Circle (○) / Four Crosses (+ + +)
Color	Fixation target color white / yellow / blue / red
Thickness (width)	Fixation target width unit: angle of view 0.10° / 0.20° / 0.30° / 0.40° / 0.50°
Size / Radius / Size Distance	Fixation target size unit: angle of view Options in the pull-down menu When Single Cross is selected Size: 1.00° to 20.00°(1°increments) When Circle is selected Radius: 0.5° and 1.00 to 20.00°(1°increments) When Four Crosses is selected Size (size of each cross): 1.00° to 20.00°(1°increments) Distance (distances between upper and lower, and right and left crosses): 1.00° to 20.00°(1°increments)
Examination	
Duration	Test duration 5s to 3600s (5s increments)

*1: When the box is checked, the background luminance 4 asb, equivalent to NIDEK MP-1 can be selected.

*2: The fixation target can be changed in the Test screen.  "14 Fixation target setting button" (page 27)

4.4.3 Application setting

Select from a pull-down menu.

Note


- The value underlined in the table is the default value.

Application	
Language	Display language. Displayed in English The setting is default and cannot be changed by users.
Eye abbreviation format	Eye indication R-L / <u>OD-OS</u>
Patient name format	Specify the patient name format. Last-Middle-First / <u>First-Middle-Last</u>
High precision pointing view zoom factor	Zoom 1x / <u>2x</u> / 4x

1 [Save] button

The changed settings are saved.

4.4.4 System settings

System settings consist of four pages. The page can be changed with the (▲, ▼) (*A) button.

○ System1/4

When not instructed, change the item from a pull-down menu.

Note

- The value underlined in the table is the default value.

Date Time	
Date	yyyy/mm/dd Press the entry field and enter with the calendar. • "Entering the date of birth" (page 57)
Time	HH:MM:SS AM/PM Press the time field and enter with the input window. Change the hour, minute, or second, with the feed buttons (▼, ▲). Confirms the setting with the OK button. HH: 12:MM:SS AM up to 11:MM:SS PM MM: 00 to 59 SS: 00 to 59
Color Camera (Color camera: Fundus image capture settings)	
Resolution	Image size (Pixels) 1M / 3M / 5M / 8M / <u>12M</u>
Quality	Image quality <u>Fine</u> / Normal / Basic
ISO speed	ISO sensitivity 100 / <u>200</u> / 400
Contrast	Contrast -2 / -1 / <u>0</u> / +1 / +2
Sharpness	Sharpness 0 / +1 / <u>+2</u> / +3 / +4 / +5

DEVICE SETTINGS AND MAINTENANCE: Parameter Settings



Saturation	Saturation -2 / -1 / <u>0</u> / +1 / +2
White balance	White balance <u>Day white</u> / Color temperature
Red gain	White balance correction R value. (When White balance is set to Day white) -128 to <u>30</u> to 127 (1 increments)
Blue gain	White balance correction B value. (When White balance is set to Day white) -128 to <u>-15</u> to 127 (1 increments)
Color Temperature	Color Temperature (When White balance is set to Color Temperature) 4000 K to <u>5000 K</u> to 8000 K (100 K increments)
IR Camera (IR camera: observation image setting)	
Enhance IR	Enhancement of an observation image <u>Yes</u> / No
Split bar	
Show	Choose whether or not to use a focus bar. <u>Yes</u> / No
Display	
Backlight	LCD back light brightness Low / <u>Medium</u> / High
Sound	
Volume	Beep tone Off / Low / <u>High</u>  : Check the selected beep tone.
Flash	
Flash level	Default value of flash intensity for color fundus image capture When the ISO speed is set to 100: 10 to <u>13</u> to 16 (1 increments) When the ISO speed is set to 200: 6 to <u>9</u> to 12 (1 increments) When the ISO speed is set to 400: 2 to <u>5</u> to 8 (1 increments)
Chinrest	
Initialization	Initialization of the chinrest position (for each test) <u>Off</u> / On
Sleep	
Timeout	Sleep mode time 5 min / 10 min / <u>15 min</u>

1 [Save] button

The changed settings are saved.

○ System 2/4

Network setting for connecting with NAVIS-EX

Date Management	
NAVIS-EX	Choose whether or not to use NAVIS-EX. Select or deselect the box "□Use NAVIS-EX".
Server name/IP address	Server name/IP address entry Press the entry field and enter with the keyboard.  "• Patient ID input: How to use the keyboard" (page 56)
Port number	Port number entry Press the entry field and enter with the keyboard.  "• Patient ID input: How to use the keyboard" (page 56)
Test connection	Test for connection with NAVIS-EX Press the Test connection...  button to perform the test. Test connection succeeded:  Test connection failed:  If the test is failed, check the server name, port number, connection PC setting, and LAN cable connection.

4

1 [Save] button

The changed settings are saved.

1) Press the Save button.

A message is displayed.

When NAVIS-EX is connected:

Switching to NAVIS-EX mode requires restarting the system. The system will be shutdown. Press OK to proceed.

When NAVIS-EX is not connected (standalone):

Switching to stand-alone mode requires restarting the system. The system will be shutdown. Press OK to proceed.

2) Press the OK button.

The input contents are saved and the MP-3 shuts down.

○ System3/4

The device information such as serial number and software information.

○ System4/4

The license information of the software on the device.

4.5 Cleaning

If the cover or the LCD monitor becomes contaminated, wipe them with a soft cloth. For severe stains, soak the cloth in a neutral detergent, wring well, and then wipe. Finally dry with a soft, dry cloth.

CAUTION

- Never use an organic solvent such as paint thinner or alcohol.
- Wipe the LCD touch-screen gently with a soft cloth. Never press hard objects against or place magnetic objects near the LCD touch-screen.
The LCD touch-screen may become damaged or malfunction may occur. Malfunction of the device may result.
- Never use a sponge or cloth soaked in water.
If liquid enters the device, malfunction may result.

4.5.1 Cleaning the objective lens

Fingerprints or dust on the objective lens reduces the quality of the fundus images. Check the cleanliness of the objective lens before use. If the objective lens is contaminated, clean it.

- 1 Use a blower to blow off any dust from the objective lens.
- 2 Wrap the cleaning paper around a thin stick (or use a cotton swab) and dip the stick into the alcohol solution so that an appropriate amount soaks into the tip of the stick, then wipe the objective lens.

Note

- Do not use a hard stick such as one made of metal. Use something that will not scratch the glass.
- Starting from the center, gently wipe the lens by moving the stick outward in a spiral.

- 3 If the lens is still contaminated, use a new swab and repeat the above process until the lens is clean.

4.6 Consumables and Maintenance Parts List

Part name	Part number	Remarks
Chinrest paper	32903-M047	1 stack

○ Parts to be maintained by service personnel

Part name	Part number	Remarks
Xenon flash lamp	15411-E101	1 set

The service life of the xenon flash lamp is 15,000 times of emission with the maximum intensity use condition (level 17).

A lower intensity extends the service life. If images are captured only with the standard level 9, the service life may become 10 times or more longer than the above.

When the lamp is replaced with a spare, restock the spare.

DEVICE SETTINGS AND MAINTENANCE: Consumables and Maintenance Parts List

SPECIFICATIONS AND TECHNICAL INFORMATION

5.1 Specifications

• Microperimetry test	Visual field angle	40°
	Maximum stimulus luminance	10,000 asb
	Background luminance	31.4 asb / 4 asb
	Stimuli luminance threshold range	Standard: Background luminance: 31.4 asb 10031.4 asb (0 dB) to 35.4 asb (34dB) (Contrast: 319.5 to 0.13) Equivalent to NIDEK MP-1 "Background luminance: 4 asb" 404 asb (0dB) to 8 asb (20dB) (Contrast: 101 to 1.01) 1,004 asb (0 dB) to 4.4 asb (34dB) (Contrast: 251 to 0.1)
	Stimulus duration	100 msec / 200 msec
	Stimuli size	Goldman I / II / III / IV / V compatible
	Threshold strategy	4-2 / 4-2-1
	Static stimulus	
	White on White	
• Fundus photography	Picture angle	45°± 5 % (when the patient's eye refraction is 0D)
	Zoom	× 0.42
	Camera resolution	Center 60 lines/mm or greater Intermediary 40 lines/mm or greater Periphery 25 lines/mm or greater
	Pixel pitch on the fundus	4.4 μm or less
	Minimum pupil diameter required	Φ4 mm
• Fixation test	Fixation target (used in conjunction with the microperimetry test)	Shape: select from three types of shape
		Color: select from white / yellow / red / blue
• Working distance	45.7 mm (distance from the objective lens to the corneal surface)	
• Diopter correction range	-25D to +15D	

SPECIFICATIONS AND TECHNICAL INFORMATION: Specifications

• Fundus auto focus range	-12D to +15D	
• Other functions	Observation / display method	10.4-inch color LCD touch-screen
• Power input	Voltage, frequency Power consumption	AC 100 to 240 V $\pm$ 10%, 50/60 Hz 160VA
• Dimensions and mass	Dimensions Mass	334 (W) $\times$ 562 (D) $\times$ 560 (H) mm 36kg
• Environmental conditions (during use)	Temperature Humidity Atmospheric pressure Installation location Others	10 to 35°C (50 to 95°F) 30 to 90% (non-condensing) 800 to 1060 hPa Indoors A well ventilated place free from hazardous particles, smoke, or fumes
• Environmental conditions (during transport and storage)	Temperature Humidity Atmospheric pressure	-30 to +60°C (-22 to 140°F) (during transport), -10 to +55°C (14 to 131°F) (during storage) 10 to 95 % (non-condensing) 700 to 1060 hPa
• Others	Service life Unit per package	8 years from the date of initial operation * Proper maintenance, inspection, and consumable parts replacement are necessary.  "Maintenance" (page 6)  "4.6 Consumables and Maintenance Parts List" (page 101) 1 unit
• Classifications	Form of protection against electrical shock: Class I Form of protection against electrical shock (applied parts): Type B applied part Degree of protection against harmful ingress of water: IPX0 Degree of safety application in the presence of flammable anesthetics and/or flammable cleaning agents: This device is not suitable to be used in a potentially flammable environment. Degree of suitability for use in an oxygen rich environment: This device is not intended for use in an oxygen rich environment. Mode of operation: This device is classified as a continuous operating device. Conformity to electromagnetic compatibility standard: In compliance with IEC 60601-1-2:2007 - EMC	
Accessories		
• Standard accessories	Response button, Response button holder, touch pen, pen stand, and power cable Chinrest paper, chinrest paper pin, object lens cap, cap holder, blower brush, dust cover, operator's manual, MP Viewer for NAVIS-EX installation CD, MP viewer for NAVIS-EX operator's manual, USB license key, USB programme recovery key, Motorized Optical Table, PC cabinet	

5.2 Exam setting: Test configuration settings

The test configuration and stimuli arrangement pattern presets by NIDEK are as follows:

○ Test configuration

Test feature and configurations

AMD	Normally used for the patient with age-related macular degeneration (AMD). This configuration is used with the following settings. - Stimulus: Goldmann III, 200 msec - Threshold strategy: 4-2 - Fixation target: Single Cross, size 1°, thickness 0.2°, color red - Pattern: MP1Macula_20deg
Macular Hole	Normally used for the patient with macular hole. This configuration is used with the following settings. - Stimulus: Goldmann I, 200 msec - Threshold strategy: 4-2 - Fixation target: Circle, radius 1°, thickness 0.2°, color red - Pattern: MP1Macula_8deg
Maculopathy	Normally used for the patient with maculopathy. This configuration is used with the following settings. - Stimulus: Goldmann III, 200 msec - Threshold strategy: 4-2 - Fixation target: Circle, radius 1°, thickness 0.2°, color red - Pattern: MP1Macula_12deg
Medical	Normally used for the patient with age-related macular degeneration (AMD), diabetic retinopathy, retinal degeneration, or retinal vascular diseases or such. This configuration is used with the following settings. - Stimulus: Goldmann III, 200 msec - Threshold strategy: 4-2 - Fixation target: Circle, radius 1°, thickness 0.2°, color red - Pattern: MP1Medical_20deg
Normal	Normally used for the patient with no obvious eye diseases. This configuration is used with the following settings. - Stimulus: Goldmann III, 200 msec - Threshold strategy: 4-2 - Fixation target: Circle, radius 1°, thickness 0.2°, color red - Pattern: MP1Normal_20deg
Surgical	Normally used for the patient with macular hole, macular pucker, or retinal detachment. This configuration is used with the following settings. - Stimulus: Goldmann III, 200 msec - Threshold strategy: 4-2 - Fixation target: Circle, radius 1°, thickness 0.2°, color red - Pattern: MP1Surgical_20deg

○ Stimuli arrangement pattern (Pattern name)

Feature of stimuli arrangement patterns

Humphery_10-2 [FACTORY]	This pattern consists of 68 stimuli in the range of 20° angle of view.
MP1Macula_10deg [FACTORY]	This pattern consists of 40 stimuli in the range of 10° angle of view. This pattern is used to study the macular mid-section sensitivity in the specified range of angle of view.
MP1Macula_12deg [FACTORY]	This pattern consists of 45 stimuli in the range of 12° angle of view. This pattern is used for the configuration "Maculopathy".
MP1Macula_20deg [FACTORY]	This pattern consists of 76 stimuli in the range of 20° angle of view. This pattern is used for the configuration "AMD".
MP1Macula_8deg [FACTORY]	This pattern consists of 45 stimuli in the range of 8° angle of view. This pattern is used for the configuration "Macular Hole".
MP1Medical_20deg [FACTORY]	This pattern consists of 29 stimuli in the range of 20° angle of view. This pattern is used for the configuration "Medical".
MP1Normal_20deg [FACTORY]	This pattern consists of 33 stimuli in the range of 20° angle of view. This pattern is used for the configuration "Normal".
MP1Retina_40deg [FACTORY]	This pattern consists of 100 stimuli in the range of 40° angle of view. As this test range is large, this pattern is normally used to study the retinal overall sensitivity.
MP1Surgical_20deg [FACTORY]	This pattern consists of 33 stimuli in the range of 20° angle of view. This pattern is used for the configuration "Surgical".

5.3 Glossary

- The following terms are used in the device and operator's manual.

Parameter	Explanation
asb (apostilb)	A unit to indicate background luminance and stimuli luminance. A larger value indicates a greater brightness.
BCEA	BCEA (Bivariate Contour Ellipse Area) indicates instability of the fixation point by size, area (angle of view), and major axis inclination of ellipses for three levels of standard deviation. Ellipse 68.2% (1Std Div): Ellipse that represents instability of fixation point within 68.2% Ellipse 95.4% (2Std Div): Ellipse that represents instability of fixation point within 95.4% Ellipse 99.6% (3Std Div): Ellipse that represents instability of fixation point within 99.6% List of references Timberlake et al. <i>Retinal Location of the Preferred Retinal Locus Relative to the Fovea in Scanning Laser Ophthalmoscope Images, Optometry and Vision Science</i>, vol. 82, N. 3, pp. E177, 2005.
False positive	False positive test During the microperimetry test, do not project the stimuli and check the patient's response (press or do not press the response button).
Fixation Stability	Fixation stability is indicated by percentage of fixation points within the circle with 2° and 4° in diameter. List of references Fujii et al. <i>Patient selection for macular translocation using SLO, Ophthalmology</i>, vol 109, nr. 9, 1999
white on white	The background and the stimuli are both white (achromate).
Auto alignment	A function to automatically align (forward and backward, right and left, up and down) to the eye front and the fundus
Auto shot	A function to automatically capture a color image when alignment and focus are optimized on the fundus
Auto tracking	Fundus images during the test are automatically tracked, compared to, and positioned on (forward and backward, right and left, up and down) the beginning image at the regular time intervals.
Auto focusing	A function to adjust focus on the retina automatically

5.4 EMC (Electromagnetic Compatibility)

The device complies with the International Electrotechnical Commission standards (IEC 60601-1-2: 2007) for electromagnetic compatibility.

WARNING

- Use the specified accessories, optional accessories, and cables.
Otherwise increased emissions or decreased immunity of the device may result.
- Do not use the device near, on, or under other electronic equipment.

CAUTION

- The device needs special precautions regarding EMC.
The device needs to be installed and used in accordance with the EMC information provided in this manual.
- Portable and mobile RF communications equipment can affect the device.

○ Specified cable

Part name	Connector Shielded	Cable Shielded	Ferrite Core	Length (m)
Power cable	-	No	No	2.5
Carl code (response button)	No	No	Yes	0.3
LAN cable	Yes	No	No	2.0

○ Essential performance

Continuous operation (measurement status)

SPECIFICATIONS AND TECHNICAL INFORMATION: EMC (Electromagnetic Compatibility)

Guidance and manufacturer's declaration - electromagnetic emissions		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	*1	
Voltage fluctuations/ Flicker emissions IEC 61000-3-3	*2	

* 1 For the regions where the rated voltage is 220 V or greater, this device complies with class A. For the regions where the rated voltage is 127 V or less, this standard is not applicable.

* 2 For the regions where the rated voltage is 220 V or greater, this device complies with this standard. For the regions where the rated voltage is 127 V or less, this standard is not applicable.

Guidance and manufacturer's declaration - electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floor should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage, dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles < 5% U_T (> 95% dip in U_T) for 5 sec	<5% U_T (> 95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles < 5% U_T (> 95% dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE U_T is the a.c. mains voltage prior to application of the test level.			

SPECIFICATIONS AND TECHNICAL INFORMATION: EMC (Electromagnetic Compatibility)

Guidance and manufacturer's declaration - electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms ($V_1=3$)	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d=1.2 \sqrt{P}$ 150 kHz to 80 MHz $d=1.2 \sqrt{P}$ 80 MHz to 800 MHz $d=2.3 \sqrt{P}$ 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m ($E_1=3$)	
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.			
^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

Recommended separation distances between portable and mobile RF communications equipment and the device			
The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d=1.2 \sqrt{P}$	80 MHz to 800 MHz $d=1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d=2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

6.

WARRANTY

1. General. NIDEK warrants to the owner of a valid proof of purchase ("Consumer") that the products and all accessories provided by NIDEK in the sales package, but excluding immaterial goods, ("Product") are free from material defects in materials, construction and workmanship, pursuant to the following terms and conditions, when installed and used normally and in accordance with operation instructions. This warranty extends only to the Consumer of Products purchased and used within the sales area covered by, NIDEK other applicable agreements notwithstanding.

2. Warranty. During the warranty period, NIDEK authorized service representatives will repair or replace, at written notice by the Consumer, at NIDEK option, without charge, a Product which is not free from material defects in construction, material or workmanship ("Materially Defective Product") and therefore is unusable or defective. If NIDEK repairs this product, NIDEK may use new or refurbished replacement parts. If NIDEK chooses to replace this product, NIDEK may choose to replace it with a new or refurbished product of the same or similar design. NIDEK will return repaired or replacement Product to Consumer in working condition only, if not agreed otherwise in writing. Replaced parts, modules or equipment will be returned into NIDEK property. Repair or replacement of Product, at NIDEK option, is Consumer's exclusive remedy.

3. Duration of Warranty Period: The Warranty Period for the Product extends for one (1) s. The Limited Warranty period starts at the date of the delivery to the Consumer, at the latest with the date of issue of the delivery note (from Distributor to Consumer). The Warranty ends at the latest one (1) year from the date shown on the note of delivery to the Distributor. For consumable products carrying an expiration date, the Limited Warranty ends, at the latest, at the date of expiration printed on the product packaging.

4. Coverage of Warranty: The warranty covers exclusively and only defects in material, construction or workmanship of the Products and does among others not cover:

- a) Product that has been subjected to misuse, accident, incident, shipping or other physical damage, improper installation, abnormal operation or handling, neglect, inundation, fire, smoke, liquid intrusion in any form, chemical or electrolytic influences, normal indicated or excessive use; or
- b) Product that has been used by a non-authorized user, or
- c) Product that has been damaged due to repair, service, alteration, and use of parts other than NIDEK or modification by anyone other than an authorized service representative of, NIDEK or
- d) Consumable or Disposable Products or components,
- e) Product to the extent that the problem experienced is caused by radiation in any form, signal conditions, network reliability or cable or antenna systems; or
- f) Product to the extent that the problem is caused by use with electrical or electronic accessories not provided by NIDEK ; or
- g) Product whose warranty/quality sticker, product serial number label or electronic serial number information have been removed, altered or rendered illegible; or
- h) Charges for installation or set-up, adjustments of customer controls, and installation or repair of systems outside the unit, as well as
- i) Any other reasons, which NIDEK not responsible for


5. Place of Jurisdiction, Applicable Legislation: These Warranty Terms are subject to Japanese law. The courts at the domicile of NIDEK shall hold jurisdiction.

Registered No : 454 2015

NOTARIAL CERTIFICATE

On this 15th day of July, 2015, before me, Hidetoshi Somiya, a Notary in and for Tokyo Legal Affairs Bureau, personally appeared Miyuki Shinohara , with satisfactory evidence of her identification and of her being the agent of Yoneji Mizuno , Corporate Officer and Senior Manager, Regulatory Affairs Department of NIDEK Co.,Ltd, and who is authorized to sign the attached document on behalf of the said company, which is organized and existing according to the laws of Japan, stated in my presence that the principal, said Yoneji Mizuno had acknowledged himself to have signed the attached document.

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APOSTILLE

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1. Country: JAPAN
This public document
2. has been signed by **Hidetoshi Somiya**
3. acting in the capacity of Notary of the Tokyo Legal Affairs Bureau
4. bears the seal/stamp of **Hidetoshi Somiya, Notary**
Certified
5. at Tokyo
6. Jul.15.2015
7. by the Ministry of Foreign Affairs
8. 15-N^o030652
9. Seal/stamp:
10. Signature

A. Ogawa

Ayako OGAWA

For the Minister for Foreign Affairs

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