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COOK MEDICAL EUROPE LTD.

O'Halloran Road,  
National Technology Park  
Limerick, Ireland.

Company Registration No. 492272

VAT Registration No. 0776234K

## Plastic Endoprotheses Stents to Maintain Patency of the Gastrointestinal Tract

### INSTRUCTIONS FOR USE

#### Table of conversions

| Parameter | Conversion     |
|-----------|----------------|
| 1 French  | 1 mm = 1Fr / 3 |
| 1 Inch    | 25,4 mm        |

I certify that I am satisfied that this document  
is an original, duly signed by the person who  
had lawful authority to sign same.  
this 28 day of August 2015

Niall Sheehy, Notary Public,  
6/7 Glentworth St.  
Limerick.

2015

Fusion® Marathon™ Anti-Reflux Biliary Stent (FS-MAR-10-5; FS-MAR-10-6; FS-MAR -10-7;  
FS-MAR -10-8; FS-MAR -10-9; FS-MAR -10-10; FS-MAR -10-11; FS-MAR -10-12)

### DEVICE DESCRIPTION

| RPN          | Stent diameter,<br>Fr /mm | Stent length, cm | Pushing catheter<br>diameter, Fr /mm | Guiding catheter<br>diameter, Fr /mm |
|--------------|---------------------------|------------------|--------------------------------------|--------------------------------------|
| FS-MAR-10-5  | 10 / 3.3                  | 5                | 9 / 3                                | 5.5 / 1.8                            |
| FS-MAR-10-6  | 10 / 3.3                  | 6                | 9 / 3                                | 5.5 / 1.8                            |
| FS-MAR-10-7  | 10 / 3.3                  | 7                | 9 / 3                                | 5.5 / 1.8                            |
| FS-MAR-10-8  | 10 / 3.3                  | 8                | 9 / 3                                | 5.5 / 1.8                            |
| FS-MAR-10-9  | 10 / 3.3                  | 9                | 9 / 3                                | 5.5 / 1.8                            |
| FS-MAR-10-10 | 10 / 3.3                  | 10               | 9 / 3                                | 5.5 / 1.8                            |
| FS-MAR-10-11 | 10 / 3.3                  | 11               | 9 / 3                                | 5.5 / 1.8                            |
| FS-MAR-10-12 | 10 / 3.3                  | 12               | 9 / 3                                | 5.5 / 1.8                            |

Fusion® Marathon™ Anti-Reflux Biliary Stent is a stent with an ePTFE sock, 2cm long, attached to the distal or duodenal end. The purpose of the sock extension is to reduce duodenobiliary reflux through the stent. The product consists of a stent and introducer (comprised of a guiding catheter 5.5 Fr/ 1.8 mm and pushing catheter 9 Fr/ 3 mm). A Fusion wire guide is recommended for use with this product and is sold separately. The product code for the Fusion® Marathon™ Anti-Reflux Biliary Stent (FS-MAR) are FS-MAR-X-Y where X denotes the French (FR) size and Y denotes the length (cm).

### INTENDED USE

This device is intended for endoscopic placement of a preloaded biliary stent to drain obstructed biliary ducts and reduce duodenal content reflux.

### NOTES

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Do not use this device for any purpose other than the stated intended use.

Store in a dry location, away from temperature extremes.

### CONTRAINDICATIONS

Contraindications include those specific to ERCP.

Inability to pass the wire guide or stent through the obstructed area.

### POTENTIAL COMPLICATIONS

Potential complications associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Additional complications that may occur with biliary stent placement include, but are not limited to: trauma to the biliary tract or duodenum, obstruction of the pancreatic duct, stent migration.

### PRECAUTIONS

**U.S. Federal law restricts this device to sale by or on the order of a physician.**

Refer to package label for minimum channel size required for this device.

Do not use this device if stent cannot advance through strictured area.

Stent should be checked periodically for signs of biliary obstruction. This stent should not be left indwelling for more than three months or as directed by physician. The wire guide diameter and the inner lumen of the wire-guided device must be compatible.

### SYSTEM PREPARATION

Remove foam block from guiding catheter.

### INSTRUCTIONS FOR USE

**I. If using Intra Ductal Exchange (IDE) port and a pre-positioned short wire (see figure 1).**

1. Unlock short wire from the Wire Guide Locking Device, and advance introducer with preloaded stent over pre-positioned wire guide, ensuring that wire guide exits guiding catheter at IDE port.

2. Introduce device into endoscope accessory channel and relock wire guide.
  3. With elevator open advance device in short increments until it is endoscopically visualized exiting scope.
  4. Advance device into appropriate duct.
  5. Fluoroscopically monitor radiopaque band on tip of guiding catheter. When a sufficient length of guiding catheter is above majority of stricture, disconnect Luer lock fitting on introducer.  
**Important:** A sufficient portion of guiding catheter must be above stricture prior to stent positioning.
  6. Maintain position of guiding catheter and advance stent into desired position using pushing catheter.
  7. Fluoroscopically and endoscopically confirm desired stent position.
  8. Unlock wire guide from Wire Guide Locking Device.
  9. Retract wire guide until it exits guiding catheter at IDE port.
- Note: Fluoroscopically visualize radiopaque band at IDE port. When radiopaque distal tip of wire guide passes band, a disengagement from wire guide lumen will occur.**
10. Advance disengaged wire guide to maintain ductal access.
  11. Relock wire guide into Wire Guide Locking Device.
  12. Retract guiding catheter into endoscope while maintaining position of pushing catheter.
  13. Confirm position of stent once guiding catheter is removed.
  14. Remove pushing catheter.

**II. If using Proximal Wire Port (PWP) and a pre-positioned long wire guide (see figure 2).**

**Note: For best results, wire guide should be kept wet.**


1. Remove stylet.
2. Advance introducer with preloaded stent over pre-positioned wire guide and into accessory channel of endoscope until proximal end of wire guide exits PWP.


**Note: Ensure that wire guide does not exit at IDE port.**

**REFER TO STEPS 3-7 IN "SECTION I", THEN RESUME WITH STEP 3 BELOW:**

3. Retract and remove both guiding catheter and wire guide into endoscope while maintaining position of pushing catheter.
4. Confirm position of stent once guiding catheter and wire guide are removed.
5. Remove pushing catheter.

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

**This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease**

Oasis® One Action Stent Introduction System with Pre-Loaded Cotton Leung® Biliary Stent (OACL-10-12; OACL-10-15; OACL-10-5; OACL-10-7; OACL-10-9; OACL-11.5-12; OACL-11.5-15; OACL-11.5-5; OACL-11.5-7; OACL-11.5-9; OACL-7-12; OACL-7-15; OACL-7-5; OACL-7-7; OACL-7-9; OACL-8.5-12; OACL-8.5-15; OACL-8.5-5; OACL-8.5-7; OACL-8.5-9)

| RPN          | Stent diameter, Fr /mm | Stent length, cm | Pushing catheter diameter, Fr /mm | Guiding catheter diameter, Fr /mm |
|--------------|------------------------|------------------|-----------------------------------|-----------------------------------|
| OACL-10-12   | 10 / 3.3               | 12               | 10 / 3.3                          | 6 / 2                             |
| OACL-10-15   | 10 / 3.3               | 15               | 10 / 3.3                          | 6 / 2                             |
| OACL-10-5    | 10 / 3.3               | 5                | 10 / 3.3                          | 6 / 2                             |
| OACL-10-7    | 10 / 3.3               | 7                | 10 / 3.3                          | 6 / 2                             |
| OACL-10-9    | 10 / 3.3               | 9                | 10 / 3.3                          | 6 / 2                             |
| OACL-11.5-12 | 11.5 / 3.8             | 12               | 11 / 3.7                          | 6 / 2                             |
| OACL-11.5-15 | 11.5 / 3.8             | 15               | 11 / 3.7                          | 6 / 2                             |
| OACL-11.5-5  | 11.5 / 3.8             | 5                | 11 / 3.7                          | 6 / 2                             |
| OACL-11.5-7  | 11.5 / 3.8             | 7                | 11 / 3.7                          | 6 / 2                             |
| OACL-11.5-9  | 11.5 / 3.8             | 9                | 11 / 3.7                          | 6 / 2                             |
| OACL-7-12    | 7 / 2.3                | 12               | 8 / 2.7                           | 5 / 1.7                           |
| OACL-7-15    | 7 / 2.3                | 15               | 8 / 2.7                           | 5 / 1.7                           |

|              |           |    |           |         |
|--------------|-----------|----|-----------|---------|
| OACL-7-5     | 7 / 2.3   | 5  | 8 / 2.7   | 5 / 1.7 |
| OACL-7-7     | 7 / 2.3   | 7  | 8 / 2.7   | 5 / 1.7 |
| OACL-7-9     | 7 / 2.3   | 9  | 8 / 2.7   | 5 / 1.7 |
| OACL-8.5-12  | 8.5 / 2.8 | 12 | 8.5 / 2.8 | 5 / 1.7 |
| OACL -8.5-15 | 8.5 / 2.8 | 15 | 8.5 / 2.8 | 5 / 1.7 |
| OACL -8.5-5  | 8.5 / 2.8 | 5  | 8.5 / 2.8 | 5 / 1.7 |
| OACL -8.5-7  | 8.5 / 2.8 | 7  | 8.5 / 2.8 | 5 / 1.7 |
| OACL -8.5-9  | 8.5 / 2.8 | 7  | 8.5 / 2.8 | 5 / 1.7 |

Oasis® One Action Stent Introduction System with Pre-Loaded Cotton Leung® Biliary Stent includes an Introducer with four radiopaque bands 5 cm apart (comprised of 170 cm pushing catheter and 205 cm guiding catheter) and pre-loaded Cotton-Leung® stent with positioning sleeve attached on it. It is recommended that this product be used with a Tracer Metro wire guide (supplied separately). The product code for the Oasis® One Action Stent Introduction System with Pre-Loaded Cotton Leung® Biliary Stent (OACL) is OACL-X-Y, where X denotes the French (FR) size and Y denotes the length in centimeters (cm) of the stent.

#### **INTENDED USE**

This device (with preloaded stent, if applicable) is intended for endoscopic biliary stent placement to drain obstructed bile ducts.

#### **NOTES**

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

#### **CONTRAINDICATIONS**

Those specific to ERCP.

#### **POTENTIAL COMPLICATIONS**

Those associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Additional complications associated with biliary stent placement include, but are not limited to: trauma to the biliary tract or duodenum, obstruction of the pancreatic duct, stent migration.

## PRECAUTIONS

Refer to package label for minimum channel size required for this device.

Wire guide diameter and inner lumen of wire-guided device must be compatible.

Not compatible with pigtail stents.

Do not use this device if stent cannot advance through strictured area.

Preloaded stents should not be left indwelling for more than three months or as directed by a physician. Periodic evaluation is recommended.

## SYSTEM PREPARATION

1.If applicable, pre-load desired stent and positioning sleeve onto tip of introducer.

## INSTRUCTIONS FOR USE

### **I. If using Intra Ductal Exchange (IDE) port and a pre-positioned short wire. (See fig. 1)**

1.Unlock short wire from Wire Guide Locking Device and advance introducer with preloaded stent over pre-positioned wire guide ensuring that wire guide exits guiding catheter at IDE port.

2.Introduce device into endoscope accessory channel and relock wire guide.

3.With elevator open, advance device in short increments until endoscopically visualized exiting scope. **Note:** As front flap of stent enters accessory channel, keep positioning sleeve over back flap until completely in accessory channel. Slide positioning sleeve over pushing catheter, keeping it clear of accessory channel.

4.Advance device into appropriate duct.

5.Fluoroscopically monitor radiopaque band(s) on guiding catheter. When a sufficient length of guiding catheter is above majority of stricture, disconnect Luer lock fitting on introducer.

**Important:** A sufficient portion of guiding catheter must be above stricture prior to stent positioning.

6.Maintain position of guiding catheter and advance stent into desired position using pushing catheter.

7.Fluoroscopically and endoscopically confirm desired stent placement.

8.Unlock wire guide from Wire Guide Locking Device.

9.Retract wire guide until it exits guiding catheter at IDE port.

**Note:** Fluoroscopically visualize radiopaque band at IDE port. When radiopaque distal tip of wire guide passes band, a disengagement from wire guide lumen will occur.

10.Advance disengaged wire guide to maintain ductal access.

11.Relock wire guide into Wire Guide Locking Device.

12.Retract guiding catheter into endoscope while maintaining position of pushing catheter.

13.Confirm position of stent once guiding catheter is removed.

14.Remove pushing catheter.

### **II. If using Proximal Wire Port (PWP) and a pre-positioned long wire guide. (See fig. 2)**

**Note:** For best results, wire guide should be kept wet.

1.Remove stylet, if applicable.

2.Advance introducer with preloaded stent over pre-positioned wire guide and into accessory channel of endoscope until proximal end of wire guide exits PWP. **Note:** If applicable, ensure that wire guide does not exit at IDE port.


**REFER TO STEPS 3-7 IN "SECTION I", THEN RESUME WITH STEP 3 BELOW:**

3.Retract and remove both guiding catheter and wire guide into endoscope while maintaining position of pushing catheter.

4.Confirm position of stent once guiding catheter and wire guide are removed.

5.Remove pushing catheter.

**Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.**

Oasis® One Action Stent Introduction System with Pre-Loaded ST-2 Tannenbaum® Biliary Stent (OATS-8.5-15; OATS -11.5-15; OATS -8.5-12; OATS -8.5-5; OATS -8.5-7; OATS -8.5-9; OATS -10-5; OATS -10-7; OATS -10-9; OATS -10-12; OATS -10-15; OATS -11.5-7; OATS -11.5-9; OATS -11.5-5; OATS-11.5-12)

| RPN          | Stent diameter, Fr /mm | Stent length, cm | Pushing catheter diameter, Fr /mm | Guiding catheter diameter, Fr /mm |
|--------------|------------------------|------------------|-----------------------------------|-----------------------------------|
| OATS-8.5-15  | 8.5 / 2.8              | 15               | 8.5 / 2.8                         | 5 / 1.7                           |
| OATS -8.5-12 | 8.5 / 2.8              | 12               | 8.5 / 2.8                         | 5 / 1.7                           |
| OATS -8.5-5  | 8.5 / 2.8              | 5                | 8.5 / 2.8                         | 5 / 1.7                           |
| OATS -8.5-7  | 8.5 / 2.8              | 7                | 8.5 / 2.8                         | 5 / 1.7                           |
| OATS -8.5-9  | 8.5 / 2.8              | 9                | 8.5 / 2.8                         | 5 / 1.7                           |
| OATS -10-5   | 10 / 3.3               | 5                | 10 / 3.3                          | 6 / 2                             |
| OATS -10-7   | 10 / 3.3               | 7                | 10 / 3.3                          | 6 / 2                             |

|               |            |    |          |       |
|---------------|------------|----|----------|-------|
| OATS -10-9    | 10 / 3.3   | 9  | 10 / 3.3 | 6 / 2 |
| OATS -10-12   | 10 / 3.3   | 12 | 10 / 3.3 | 6 / 2 |
| OATS -10-15   | 10 / 3.3   | 15 | 10 / 3.3 | 6 / 2 |
| OATS -11.5-7  | 11.5 / 3.8 | 7  | 11 / 3.7 | 6 / 2 |
| OATS -11.5-9  | 11.5 / 3.8 | 9  | 11 / 3.7 | 6 / 2 |
| OATS -11.5-5  | 11.5 / 3.8 | 5  | 11 / 3.7 | 6 / 2 |
| OATS-11.5-12  | 11.5 / 3.8 | 12 | 11 / 3.7 | 6 / 2 |
| OATS -11.5-15 | 11.5 / 3.8 | 15 | 11 / 3.7 | 6 / 2 |

Oasis® One Action Stent Introduction System with Pre-Loaded ST-2 Tannenbaum® Biliary Stent includes Introducer with four radiopaque bands 5 cm apart (comprised of 170 cm pushing catheter and 205 cm guiding catheter) and preloaded ST-2 Soehendra Tannenbaum stent with positioning sleeve attached on it. It is recommended that this set be used with a Tracer Metro wire guide (supplied separately). The product code for the Oasis® One Action Stent Introduction System with Pre-Loaded ST-2 Tannenbaum® Biliary Stent (OATS) is OATS-X-Y, where X denotes the French (FR) size and Y denotes the length in centimeters (cm) of the stent.

#### **INTENDED USE**

This device (with preloaded stent, if applicable) is intended for endoscopic biliary stent placement to drain obstructed bile ducts.

#### **NOTES**

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

#### **CONTRAINDICATIONS**

Those specific to ERCP.

#### **POTENTIAL COMPLICATIONS**

Those associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Additional complications associated with biliary stent placement include, but are not limited to: trauma to the biliary tract or duodenum, obstruction of the pancreatic duct, stent migration.

## PRECAUTIONS

Refer to package label for minimum channel size required for this device.

Wire guide diameter and inner lumen of wire-guided device must be compatible.

Not compatible with pigtail stents.

Do not use this device if stent cannot advance through strictured area.

Preloaded stents should not be left indwelling for more than three months or as directed by a physician. Periodic evaluation is recommended.

## SYSTEM PREPARATION

1.If applicable, pre-load desired stent and positioning sleeve onto tip of introducer.

## INSTRUCTIONS FOR USE

### **I. If using Intra Ductal Exchange (IDE) port and a pre-positioned short wire. (See fig. 1)**

1.Unlock short wire from Wire Guide Locking Device and advance introducer with preloaded stent over pre-positioned wire guide ensuring that wire guide exits guiding catheter at IDE port.

2.Introduce device into endoscope accessory channel and relock wire guide.

3.With elevator open, advance device in short increments until endoscopically visualized exiting scope. **Note:** As front flap of stent enters accessory channel, keep positioning sleeve over back flap until completely in accessory channel. Slide positioning sleeve over pushing catheter, keeping it clear of accessory channel.

4.Advance device into appropriate duct.

5.Fluoroscopically monitor radiopaque band(s) on guiding catheter. When a sufficient length of guiding catheter is above majority of stricture, disconnect Luer lock fitting on introducer.

**Important:** A sufficient portion of guiding catheter must be above stricture prior to stent positioning.

6.Maintain position of guiding catheter and advance stent into desired position using pushing catheter.

7.Fluoroscopically and endoscopically confirm desired stent placement.

8.Unlock wire guide from Wire Guide Locking Device.

9.Retract wire guide until it exits guiding catheter at IDE port.

**Note:** Fluoroscopically visualize radiopaque band at IDE port. When radiopaque distal tip of wire guide passes band, a disengagement from wire guide lumen will occur.

10.Advance disengaged wire guide to maintain ductal access.

11.Relock wire guide into Wire Guide Locking Device.

12.Retract guiding catheter into endoscope while maintaining position of pushing catheter.

13.Confirm position of stent once guiding catheter is removed.

14.Remove pushing catheter.

### **II. If using Proximal Wire Port (PWP) and a pre-positioned long wire guide. (See fig. 2)**

**Note:** For best results, wire guide should be kept wet.

1.Remove stylet, if applicable.

2.Advance introducer with preloaded stent over pre-positioned wire guide and into accessory channel of endoscope until proximal end of wire guide exits PWP. **Note:** If applicable, ensure that wire guide does not exit at IDE port.


**REFER TO STEPS 3-7 IN "SECTION I", THEN RESUME WITH STEP 3 BELOW:**

3.Retract and remove both guiding catheter and wire guide into endoscope while maintaining position of pushing catheter.

4.Confirm position of stent once guiding catheter and wire guide are removed.

5.Remove pushing catheter.

**Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.**

Geenen® Pancreatic Stent Set (GEPD-10-12; GEPD-10-15; GEPD-10-3; GEPD -10-5; GEPD -10-7; GEPD -10-9; GEPD -11.5-10; GEPD -11.5-12; GEPD -11.5-6; GEPD -11.5-8; GEPD-3-10; GEPD -3-11; GEPD -3-12; GEPD -3-13; GEPD -3-14; GEPD -3-15; GEPD -3-3; GEPD -3-4; GEPD -3-5; GEPD -3-6; GEPD -3-7; GEPD -3-8; GEPD -3-9; GEPD -5-10; GEPD -5-11; GEPD-5-12; GEPD -5-13; GEPD -5-14; GEPD -5-15; GEPD -5-2; GEPD -5-3; GEPD -5-4;- GEPD 5-5; GEPD -5-6; GEPD -5-7; GEPD -5-8; GEPD -5-9; GEPD -7-10; GEPD -7-11; GEPD -7-12; GEPD -7-13; GEPD -7-14; GEPD -7-15; GEPD -7-3; GEPD -7-4; GEPD -7-5; GEPD -7-6; GEPD -7-7; GEPD -7-8; GEPD -7-9; GEPD -8.5-12; GEPD -8.5-3; GEPD -8.5-5; GEPD -8.5-7; GEPD -8.5-9)

| RPN        | Stent diameter, Fr /mm | Stent length, cm | Pushing catheter diameter, Fr /mm | Guiding catheter diameter, Fr /mm |
|------------|------------------------|------------------|-----------------------------------|-----------------------------------|
| GEPD-10-12 | 10 /3.3                | 12               | 10 /3.3                           | 6 /2                              |
| GEPD-10-15 | 10 /3.3                | 15               | 10 /3.3                           | 6 /2                              |
| GEPD-10-3  | 10 /3.3                | 3                | 10 /3.3                           | 6 /2                              |
| GEPD-10-5  | 10 /3.3                | 5                | 10 /3.3                           | 6 /2                              |
| GEPD-10-7  | 10 /3.3                | 7                | 10 /3.3                           | 6 /2                              |

|              |            |    |          |       |
|--------------|------------|----|----------|-------|
| GEPD-10-9    | 10 / 3.3   | 9  | 10 / 3.3 | 6 / 2 |
| GEPD-11.5-10 | 11.5 / 3.8 | 10 | 11 / 3.6 | 6 / 2 |
| GEPD-11.5-12 | 11.5 / 3.8 | 12 | 11 / 3.6 | 6 / 2 |
| GEPD-11.5-6  | 11.5 / 3.8 | 6  | 11 / 3.6 | 6 / 2 |
| GEPD-11.5-8  | 11.5 / 3.8 | 8  | 11 / 3.6 | 6 / 2 |
| GEPD-3-10    | 3 / 1      | 10 | 3 / 1    | N/A   |
| GEPD-3-11    | 3 / 1      | 11 | 3 / 1    | N/A   |
| GEPD-3-12    | 3 / 1      | 12 | 3 / 1    | N/A   |
| GEPD-3-13    | 3 / 1      | 13 | 3 / 1    | N/A   |
| GEPD-3-14    | 3 / 1      | 14 | 3 / 1    | N/A   |
| GEPD-3-15    | 3 / 1      | 15 | 3 / 1    | N/A   |
| GEPD-3-3     | 3 / 1      | 3  | 3 / 1    | N/A   |
| GEPD-3-4     | 3 / 1      | 4  | 3 / 1    | N/A   |
| GEPD-3-5     | 3 / 1      | 5  | 3 / 1    | N/A   |
| GEPD-3-6     | 3 / 1      | 6  | 3 / 1    | N/A   |
| GEPD-3-7     | 3 / 1      | 7  | 3 / 1    | N/A   |
| GEPD-3-8     | 3 / 1      | 8  | 3 / 1    | N/A   |
| GEPD-3-9     | 3 / 1      | 9  | 3 / 1    | N/A   |
| GEPD-5-10    | 5 / 1.7    | 10 | 5 / 1.7  | N/A   |
| GEPD-5-11    | 5 / 1.7    | 11 | 5 / 1.7  | N/A   |
| GEPD-5-12    | 5 / 1.7    | 12 | 5 / 1.7  | N/A   |
| GEPD-5-13    | 5 / 1.7    | 13 | 5 / 1.7  | N/A   |
| GEPD-5-14    | 5 / 1.7    | 14 | 5 / 1.7  | N/A   |
| GEPD-5-15    | 5 / 1.7    | 15 | 5 / 1.7  | N/A   |
| GEPD-5-2     | 5 / 1.7    | 2  | 5 / 1.7  | N/A   |
| GEPD-5-3     | 5 / 1.7    | 3  | 5 / 1.7  | N/A   |

|             |           |    |           |         |
|-------------|-----------|----|-----------|---------|
| GEPD-5-4    | 5 / 1.7   | 4  | 5 / 1.7   | N/A     |
| GEPD-5-5    | 5 / 1.7   | 5  | 5 / 1.7   | N/A     |
| GEPD-5-6    | 5 / 1.7   | 6  | 5 / 1.7   | N/A     |
| GEPD-5-7    | 5 / 1.7   | 7  | 5 / 1.7   | N/A     |
| GEPD-5-8    | 5 / 1.7   | 8  | 5 / 1.7   | N/A     |
| GEPD-5-9    | 5 / 1.7   | 9  | 5 / 1.7   | N/A     |
| GEPD-7-10   | 7 / 2.3   | 10 | 7 / 2.3   | N/A     |
| GEPD-7-11   | 7 / 2.3   | 11 | 7 / 2.3   | N/A     |
| GEPD-7-12   | 7 / 2.3   | 12 | 7 / 2.3   | N/A     |
| GEPD-7-13   | 7 / 2.3   | 13 | 7 / 2.3   | N/A     |
| GEPD-7-14   | 7 / 2.3   | 14 | 7 / 2.3   | N/A     |
| GEPD-7-15   | 7 / 2.3   | 15 | 7 / 2.3   | N/A     |
| GEPD-7-3    | 7 / 2.3   | 3  | 7 / 2.3   | N/A     |
| GEPD-7-4    | 7 / 2.3   | 4  | 7 / 2.3   | N/A     |
| GEPD-7-5    | 7 / 2.3   | 5  | 7 / 2.3   | N/A     |
| GEPD-7-6    | 7 / 2.3   | 6  | 7 / 2.3   | N/A     |
| GEPD-7-7    | 7 / 2.3   | 7  | 7 / 2.3   | N/A     |
| GEPD-7-8    | 7 / 2.3   | 8  | 7 / 2.3   | N/A     |
| GEPD-7-9    | 7 / 2.3   | 9  | 7 / 2.3   | N/A     |
| GEPD-8.5-12 | 8.5 / 2.8 | 12 | 8.5 / 2.8 | 5 / 1.7 |
| GEPD-8.5-3  | 8.5 / 2.8 | 3  | 8.5 / 2.8 | 5 / 1.7 |
| GEPD-8.5-5  | 8.5 / 2.8 | 5  | 8.5 / 2.8 | 5 / 1.7 |
| GEPD-8.5-7  | 8.5 / 2.8 | 7  | 8.5 / 2.8 | 5 / 1.7 |
| GEPD-8.5-9  | 8.5 / 2.8 | 9  | 8.5 / 2.8 | 5 / 1.7 |

Geenen® Pancreatic Stent Set (GEPD) comprises of a stent with 8 cm positioning sleeve attached on and a pushing catheter 170 cm in length. A guiding catheter 320 cm in length is also

Page 12 of 34

supplied with the 8.5Fr/ 2.8 mm, 10Fr/ 3.3 mm and 11.5Fr/ 3.8 mm stent sets. The Geenen® Pancreatic Stent Set is contoured for positioning in the pancreatic duct. The stent is non-pigtail stent, has offset flaps to help minimize migration, side holes to facilitate drainage and has a tapered tip to facilitate smooth cannulation. The product code for the Geenen® Pancreatic Stent/Stent Set is GEPD-X-Y where X denotes the French (FR) size and Y denotes the length in centimetres (cm).

### **INTENDED USE**

This device is used to drain obstructed pancreatic ducts.

### **NOTES**

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

For 8.5 and 10 FR Sof-Flex® stents, lubricate the guiding catheter with water-soluble lubricant.

Retain the positioning sleeve (if any)\* for use when introducing the stent flaps into the accessory channel.

### **CONTRAINDICATIONS**

Those specific to ERCP and any procedure to be performed in conjunction with stent placement. Inability to pass the wire guide or stent through the obstructed area.

### **POTENTIAL COMPLICATIONS**

Those associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Those associated with pancreatic stent placement include, but are not limited to: trauma to the pancreatic tract or duodenum, stent migration.

### **PRECAUTIONS**

Refer to package label for minimum channel size required for this device.

The positioning sleeve\* or the pigtail straightener\* is not intended for use in the accessory channel of the endoscope.

Care must be exercised when straightening the pigtail curls\* in order to avoid kinking or breaking the stent.

The tapered tip end of the stent\* must be positioned in the pancreatic duct while the other end remains in the duodenum.

This device should not be left indwelling for more than three months or as directed by a physician. Periodic evaluation is recommended.

Select the Cook stent introducer system (if not included) in the appropriate French size. Stent introducer system may be supplied separately.

## WARNINGS

MRI compatibility of pancreatic stents with radiopaque bands has not been established.

## INSTRUCTIONS FOR USE

### I. FOR NON-PIGTAIL STENTS:

1. Gently ensure full extension of all side flaps\*.
2. Load positioning sleeve onto duodenal flap end of stent.

#### For 7 FR and smaller stents:

- 3a. Introduce stent, tapered tip first, and positioning sleeve onto a pre-positioned wire guide.

#### For 8.5 FR and larger stents:

- 3b. Remove Tuohy-Borst adapter from end of guiding catheter\*, then introduce guiding catheter\* into accessory channel over a pre-positioned wire guide.
- 3c. Introduce stent, tapered tip first, and positioning sleeve onto guiding catheter\* and pre-positioned wire guide.

**THEN REFER TO STEPS 4-7 BELOW**

### II. FOR PIGTAIL STENTS:

#### For 7 FR and smaller stents:

1. Advance pigtail straightener\* along shaft of stent from the tapered tip end in order to straighten pigtail curl. (See fig. 1)
2. Introduce stent, tapered tip first, and pigtail straightener\* onto pre-positioned wire guide.
3. Advance pushing catheter over wire guide to advance pigtail stent into accessory channel.

#### For 8.5 FR and larger stents:


- 3a. Advance pigtail straightener\* along shaft of stent from the tapered tip end in order to straighten pigtail curl. (See fig. 1)
  - 3b. Introduce stent, pigtail first, and pigtail straightener\* onto guiding catheter\*.
- Advance guiding catheter over pre-positioned wire guide.

**THEN REFER TO STEPS 4-7 BELOW**

4. Advance guiding catheter\* and/or pushing catheter in 1-2 cm increments until stent is in desired position.
5. Fluoroscopically and endoscopically confirm desired stent position.
6. After confirming stent position, gently remove wire guide, then the guiding catheter\* from endoscope while maintaining position of stent with pushing catheter.
7. Gently remove pushing catheter from accessory channel.

**Upon completion of procedure, dispose of device(s) per institutional guidelines for biohazardous medical waste.**

\*if any

The Geenen® Pancreatic Stent Set (GPDS-5-15; GPDS-5-2; GPDS-5-12; GPDS-5-3; GPDS-5-5; GPDS-5-7; GPDS-5-9; GPDS-7-12; GPDS-7-5; GPDS-7-7; GPDS-7-9; GPDS-7-12).

| RPN       | Stent diameter, Fr /mm | Stent length, cm | Pushing catheter diameter, Fr /mm |
|-----------|------------------------|------------------|-----------------------------------|
| GPDS-5-15 | 5 / 1.7                | 15               | 5 / 1.7                           |
| GPDS-5-2  | 5 / 1.7                | 2                | 5 / 1.7                           |
| GPDS-5-12 | 5 / 1.7                | 12               | 5 / 1.7                           |
| GPDS-5-3  | 5 / 1.7                | 3                | 5 / 1.7                           |
| GPDS-5-5  | 5 / 1.7                | 5                | 5 / 1.7                           |
| GPDS-5-7  | 5 / 1.7                | 7                | 5 / 1.7                           |
| GPDS-5-9  | 5 / 1.7                | 9                | 5 / 1.7                           |
| GPDS-7-5  | 7 / 2.3                | 5                | 7 / 2.3                           |
| GPDS-7-7  | 7 / 2.3                | 7                | 7 / 2.3                           |
| GPDS-7-9  | 7 / 2.3                | 9                | 7 / 2.3                           |
| GPDS-7-12 | 7 / 2.3                | 12               | 7 / 2.3                           |

The Geenen® Pancreatic Stent Set (GPDS) comprises of a stent with 8 cm positioning sleeve attached on it and pushing catheter 170 cm in length. The Geenen® Pancreatic Stent Set is contoured for positioning in the pancreatic duct. The stent is non-pigtail stent, has offset flaps to help minimize migration, side holes to facilitate drainage and has a tapered tip in larger French sizes to facilitate smooth cannulation. The product code for the Geenen® Pancreatic Stent Set is GPDS-X-Y, where X denotes the French (FR) size and Y denotes the length in centimetres (cm).

#### **INTENDED USE**

This device is used to drain obstructed pancreatic ducts.

#### **NOTES**

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.  
Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.  
For 8.5 and 10 FR Sof-Flex® stents, lubricate the guiding catheter with water-soluble lubricant.  
Retain the positioning sleeve (if any)\* for use when introducing the stent flaps into the accessory channel.

### **CONTRAINDICATIONS**

Those specific to ERCP and any procedure to be performed in conjunction with stent placement.  
Inability to pass the wire guide or stent through the obstructed area.

### **POTENTIAL COMPLICATIONS**

Those associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Those associated with pancreatic stent placement include, but are not limited to: trauma to the pancreatic tract or duodenum, stent migration.

### **PRECAUTIONS**

Refer to package label for minimum channel size required for this device.

The positioning sleeve\* or the pigtail straightener\* is not intended for use in the accessory channel of the endoscope.

Care must be exercised when straightening the pigtail curls\* in order to avoid kinking or breaking the stent.

The tapered tip end of the stent\* must be positioned in the pancreatic duct while the other end remains in the duodenum.

This device should not be left indwelling for more than three months or as directed by a physician. Periodic evaluation is recommended.

Select the Cook stent introducer system (if not included) in the appropriate French size. Stent introducer system may be supplied separately.

### **WARNINGS**

MRI compatibility of pancreatic stents with radiopaque bands has not been established.

### **INSTRUCTIONS FOR USE**

#### **I. FOR NON-PIGTAIL STENTS:**

1. Gently ensure full extension of all side flaps\*.
2. Load positioning sleeve onto duodenal flap end of stent.

##### **For 7 FR and smaller stents:**

- 3a. Introduce stent, tapered tip first, and positioning sleeve onto a pre-positioned wire guide.

##### **For 8.5 FR and larger stents:**

- 3b. Remove Tuohy-Borst adapter from end of guiding catheter\*, then introduce guiding catheter\* into accessory channel over a pre-positioned wire guide.
- 3c. Introduce stent, tapered tip first, and positioning sleeve onto guiding catheter\* and pre-positioned wire guide.

**THEN REFER TO STEPS 4-7 BELOW**

#### **II. FOR PIGTAIL STENTS:**

**For 7 FR and smaller stents:**

1. Advance pigtail straightener\* along shaft of stent from the tapered tip end in order to straighten pigtail curl. (See fig. 1)
2. Introduce stent, tapered tip first, and pigtail straightener\* onto pre-positioned wire guide.
3. Advance pushing catheter over wire guide to advance pigtail stent into accessory channel.

**For 8.5 FR and larger stents:**


- 3a. Advance pigtail straightener\* along shaft of stent from the tapered tip end in order to straighten pigtail curl. (See fig. 1)
  - 3b. Introduce stent, pigtail first, and pigtail straightener\* onto guiding catheter\*.
- Advance guiding catheter over pre-positioned wire guide.

**THEN REFER TO STEPS 4-7 BELOW**

4. Advance guiding catheter\* and/or pushing catheter in 1-2 cm increments until stent is in desired position.
5. Fluoroscopically and endoscopically confirm desired stent position.
6. After confirming stent position, gently remove wire guide, then the guiding catheter\* from endoscope while maintaining position of stent with pushing catheter.
7. Gently remove pushing catheter from accessory channel.

**Upon completion of procedure, dispose of device(s) per institutional guidelines for biohazardous medical waste.**

\*if any

Cotton-Huibregtse® Biliary Stent Set (CHBS-10-10; CHBS-10-11; CHBS-10-12; CHBS-10-15; CHBS-10-3; CHBS-10-5; CHBS-10-7; CHBS-10-9; CHBS-11.5-10; CHBS-11.5-12; CHBS-11.5-15; CHBS-11.5-5; CHBS-11.5-7; CHBS-11.5-9; CHBS-7-10; CHBS-7-12; CHBS-7-15; CHBS-7-5; CHBS -7-3; CHBS -7-6; CHBS-7-7; CHBS-7-9; CHBS-8.5-12; CHBS-8.5-15; CHBS-8.5-5; CHBS-8.5-7; CHBS-8.5-9)

| RPN        | Stent diameter, Fr /mm | Stent length, cm | Pushing catheter diameter, Fr /mm | Guiding catheter diameter, Fr /mm |
|------------|------------------------|------------------|-----------------------------------|-----------------------------------|
| CHBS-10-10 | 10 / 3.3               | 10               | 10 / 3.3                          | 6 / 2                             |
| CHBS-10-11 | 10 / 3.3               | 11               | 10 / 3.3                          | 6 / 2                             |
| CHBS-10-12 | 10 / 3.3               | 12               | 10 / 3.3                          | 6 / 2                             |
| CHBS-10-15 | 10 / 3.3               | 15               | 10 / 3.3                          | 6 / 2                             |

|              |            |    |           |         |
|--------------|------------|----|-----------|---------|
| CHBS-10-3    | 10 / 3.3   | 3  | 10 / 3.3  | 6 / 2   |
| CHBS-10-5    | 10 / 3.3   | 5  | 10 / 3.3  | 6 / 2   |
| CHBS-10-7    | 10 / 3.3   | 7  | 10 / 3.3  | 6 / 2   |
| CHBS-10-9    | 10 / 3.3   | 9  | 10 / 3.3  | 6 / 2   |
| CHBS-11.5-10 | 11.5 / 3.8 | 10 | 11 / 3.7  | 6 / 2   |
| CHBS-11.5-12 | 11.5 / 3.8 | 12 | 11 / 3.7  | 6 / 2   |
| CHBS-11.5-15 | 11.5 / 3.8 | 15 | 11 / 3.7  | 6 / 2   |
| CHBS-11.5-5  | 11.5 / 3.8 | 5  | 11 / 3.7  | 6 / 2   |
| CHBS-11.5-7  | 11.5 / 3.8 | 7  | 11 / 3.7  | 6 / 2   |
| CHBS-11.5-9  | 11.5 / 3.8 | 9  | 11 / 3.7  | 6 / 2   |
| CHBS-7-10    | 7 / 2.3    | 10 | 7 / 2.3   | N/A     |
| CHBS-7-12    | 7 / 2.3    | 12 | 7 / 2.3   | N/A     |
| CHBS-7-15    | 7 / 2.3    | 15 | 7 / 2.3   | N/A     |
| CHBS-7-5     | 7 / 2.3    | 5  | 7 / 2.3   | N/A     |
| CHBS-7-3     | 7 / 2.3    | 3  | 7 / 2.3   | N/A     |
| CHBS-7-6     | 7 / 2.3    | 6  | 7 / 2.3   | N/A     |
| CHBS-7-7     | 7 / 2.3    | 7  | 7 / 2.3   | N/A     |
| CHBS-7-9     | 7 / 2.3    | 9  | 7 / 2.3   | N/A     |
| CHBS-8.5-12  | 8.5 / 2.8  | 12 | 8.5 / 2.8 | 5 / 1.7 |
| CHBS-8.5-15  | 8.5 / 2.8  | 15 | 8.5 / 2.8 | 5 / 1.7 |
| CHBS-8.5-5   | 8.5 / 2.8  | 5  | 8.5 / 2.8 | 5 / 1.7 |
| CHBS-8.5-7   | 8.5 / 2.8  | 7  | 8.5 / 2.8 | 5 / 1.7 |
| CHBS-8.5-9   | 8.5 / 2.8  | 9  | 8.5 / 2.8 | 5 / 1.7 |

This product typically contains a stent with positioning sleeve, guiding catheter length 320 cm and pushing catheter length 170 cm (except the CHBS-7-Y which does not have guiding catheter). The stent has a proximal and a distal flap. It also has a tapered tip at the proximal end

to facilitate smooth cannulation. The product code for the Cotton-Huibregtse® Biliary Stent Set is CHBS-X-Y, where X denotes the French (FR) size and Y denotes the length in centimeters (cm).

### **INTENDED USE**

This device is used to drain obstructed biliary ducts.

### **NOTES**

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Do not use this device for any purpose other than stated intended use.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Notify Cook for return authorization.

This device must be stored in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

For 10 FR Sof-Flex® stents, lubricate the guiding catheter with water-soluble lubricant.

Retain the positioning sleeve (if any)\* for use when introducing the stent flaps into the accessory channel.

### **CONTRAINDICATIONS**

Those specific to ERCP and any procedure to be performed in conjunction with stent placement. Inability to pass wire guide or stent through obstructed area.

### **POTENTIAL COMPLICATIONS**

Those associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Those associated with biliary stent placement include, but are not limited to: trauma to the biliary tract or duodenum, obstruction of the pancreatic duct, stent migration.

### **PRECAUTIONS**

Refer to package label for minimum channel size required for this device.

This stent must be placed under fluoroscopic monitoring.

The positioning sleeve\* or the pigtail straightener\* is not intended for use in the accessory channel of the endoscope.

Care must be exercised when straightening the pigtail curls\* in order to avoid kinking or breaking the stent.

The tapered tip end of the stent\* or side holes\* must be positioned in the common bile duct while the other end remains in the duodenum.

The longer flaps of the ST-2 Biliary Stent should be positioned in the duct while the shorter flaps remain in the duodenum. (See fig. 1)

This device should not be left indwelling for more than three months or as directed by a physician. Periodic evaluation is recommended.

Select the Cook stent introducer system (if not included) in the appropriate French size.

### **Warning**

For stents with radiopaque bands, MRI compatibility has not been established.

## INSTRUCTIONS FOR USE

### I. FOR NON-PIGTAIL STENTS:

1. Gently ensure full extension of all side flaps\*.
2. Load positioning sleeve onto duodenal flap end of stent.

For 5 to 7 FR stents:

- 3.a. Introduce stent, tapered tip first, and positioning sleeve onto a pre-positioned wire guide.

For 8.5 FR and larger stents:

- 3.b. Remove Tuohy-Borst adapter from end of guiding catheter, then introduce guiding catheter into accessory channel over a pre-positioned wire guide.
- 3.c. Introduce stent, tapered tip first, and positioning sleeve onto guiding catheter and pre-positioned wire guide.

THEN REFER TO STEPS 4-7 BELOW

### II. FOR PIGTAIL STENTS:

1. Advance pigtail straightener along shaft of stent to tapered tip end in order to straighten pigtail curl. (See fig. 2)
2. Introduce stent, tapered tip first, and pigtail straightener onto pre-positioned wire guide until straightener reaches second curl.
3. Advance pushing catheter over wire guide to advance pigtail stent into accessory channel.

THEN REFER TO STEPS 4-7 BELOW



4. Advance guiding catheter\* and/or pushing catheter in 1-2 cm increments until stent is in desired position.
5. Fluoroscopically and endoscopically confirm desired stent position.

**Note (for non-pigtail stents only):** Replace Tuohy-Borst adapter onto end of guiding catheter\*. Inject contrast, if desired, to fluoroscopically visualize stent position.

6. After confirming stent position, gently remove wire guide, then guiding catheter\* from endoscope while maintaining position of stent with pushing catheter.
7. Gently remove pushing catheter from accessory channel.

**Upon completion of procedure, dispose of introducer components per institutional guidelines for biohazardous medical waste.**

\* (if any)

Cotton-Leung® (Amsterdam) Biliary Stent Set (CLBS-10-10; CLBS-10-11; CLBS-10-12; CLBS-10-13; CLBS-10-14; CLBS-10-15; CLBS-10-16; CLBS-10-17; CLBS-10-18; CLBS-10-3; CLBS-10-4; CLBS-10-5; CLBS-10-6; CLBS-10-7; CLBS-10-8; CLBS-10-9; CLBS-11.5-10; CLBS-11.5-11; CLBS-11.5-12; CLBS-11.5-13; CLBS-11.5-14; CLBS-11.5-15; CLBS-11.5-16; CLBS-11.5-17; CLBS-11.5-18; CLBS-11.5-3; CLBS-11.5-5; CLBS-11.5-6; CLBS-11.5-7; CLBS-11.5-8; CLBS-11.5-9; CLBS-5-5; CLBS-5-7; CLBS-5-9; CLBS-5-12; CLBS-7-10; CLBS-7-11; CLBS-7-12; CLBS-7-13; CLBS-7-14; CLBS-7-15; CLBS-7-16; CLBS-7-17; CLBS-7-18; CLBS-7-3; CLBS-7-4; CLBS-7-5; CLBS-7-6; CLBS-7-7; CLBS-7-8; CLBS-7-9; CLBS-8.5-10; CLBS-8.5-11; CLBS-8.5-12; CLBS-8.5-13; CLBS-8.5-14; CLBS-8.5-15; CLBS-8.5-16; CLBS-8.5-17; CLBS-8.5-18; CLBS-8.5-3; CLBS-8.5-5; CLBS-8.5-6; CLBS-8.5-7; CLBS-8.5-8; CLBS-8.5-9)

| RPN          | Stent diameter,<br>Fr /mm | Stent<br>length, cm | Pushing catheter<br>diameter, Fr /mm | Guiding catheter<br>diameter, Fr /mm |
|--------------|---------------------------|---------------------|--------------------------------------|--------------------------------------|
| CLBS-10-10   | 10 / 3.3                  | 10                  | 10 / 3.3                             | 6 / 2                                |
| CLBS-10-11   | 10 / 3.3                  | 11                  | 10 / 3.3                             | 6 / 2                                |
| CLBS-10-12   | 10 / 3.3                  | 12                  | 10 / 3.3                             | 6 / 2                                |
| CLBS-10-13   | 10 / 3.3                  | 13                  | 10 / 3.3                             | 6 / 2                                |
| CLBS-10-14   | 10 / 3.3                  | 14                  | 10 / 3.3                             | 6 / 2                                |
| CLBS-10-15   | 10 / 3.3                  | 15                  | 10 / 3.3                             | 6 / 2                                |
| CLBS-10-16   | 10 / 3.3                  | 16                  | 10 / 3.3                             | 6 / 2                                |
| CLBS-10-17   | 10 / 3.3                  | 17                  | 10 / 3.3                             | 6 / 2                                |
| CLBS-10-18   | 10 / 3.3                  | 18                  | 10 / 3.3                             | 6 / 2                                |
| CLBS-10-3    | 10 / 3.3                  | 3                   | 10 / 3.3                             | 6 / 2                                |
| CLBS-10-4    | 10 / 3.3                  | 4                   | 10 / 3.3                             | 6 / 2                                |
| CLBS-10-5    | 10 / 3.3                  | 5                   | 10 / 3.3                             | 6 / 2                                |
| CLBS-10-6    | 10 / 3.3                  | 6                   | 10 / 3.3                             | 6 / 2                                |
| CLBS-10-7    | 10 / 3.3                  | 7                   | 10 / 3.3                             | 6 / 2                                |
| CLBS-10-8    | 10 / 3.3                  | 8                   | 10 / 3.3                             | 6 / 2                                |
| CLBS-10-9    | 10 / 3.3                  | 9                   | 10 / 3.3                             | 6 / 2                                |
| CLBS-11.5-10 | 11.5 / 3.8                | 10                  | 11 / 3.7                             | 6 / 2                                |
| CLBS-11.5-11 | 11.5 / 3.8                | 11                  | 11 / 3.7                             | 6 / 2                                |

|              |            |    |          |       |
|--------------|------------|----|----------|-------|
| CLBS-11.5-12 | 11.5 / 3.8 | 12 | 11 / 3.7 | 6 / 2 |
| CLBS-11.5-13 | 11.5 / 3.8 | 13 | 11 / 3.7 | 6 / 2 |
| CLBS-11.5-14 | 11.5 / 3.8 | 14 | 11 / 3.7 | 6 / 2 |
| CLBS-11.5-15 | 11.5 / 3.8 | 15 | 11 / 3.7 | 6 / 2 |
| CLBS-11.5-16 | 11.5 / 3.8 | 16 | 11 / 3.7 | 6 / 2 |
| CLBS-11.5-17 | 11.5 / 3.8 | 17 | 11 / 3.7 | 6 / 2 |
| CLBS-11.5-18 | 11.5 / 3.8 | 18 | 11 / 3.7 | 6 / 2 |
| CLBS-11.5-3  | 11.5 / 3.8 | 3  | 11 / 3.7 | 6 / 2 |
| CLBS-11.5-5  | 11.5 / 3.8 | 5  | 11 / 3.7 | 6 / 2 |
| CLBS-11.5-6  | 11.5 / 3.8 | 6  | 11 / 3.7 | 6 / 2 |
| CLBS-11.5-7  | 11.5 / 3.8 | 7  | 11 / 3.7 | 6 / 2 |
| CLBS-11.5-8  | 11.5 / 3.8 | 8  | 11 / 3.7 | 6 / 2 |
| CLBS-11.5-9  | 11.5 / 3.8 | 9  | 11 / 3.7 | 6 / 2 |
| CLBS-5-5     | 5 / 1.7    | 5  | 5 / 1.7  | N/A   |
| CLBS-5-7     | 5 / 1.7    | 7  | 5 / 1.7  | N/A   |
| CLBS-5-9     | 5 / 1.7    | 9  | 5 / 1.7  | N/A   |
| CLBS-5-12    | 5 / 1.7    | 12 | 5 / 1.7  | N/A   |
| CLBS-7-10    | 7 / 2.3    | 10 | 7 / 2.3  | N/A   |
| CLBS-7-11    | 7 / 2.3    | 11 | 7 / 2.3  | N/A   |
| CLBS-7-12    | 7 / 2.3    | 12 | 7 / 2.3  | N/A   |
| CLBS-7-13    | 7 / 2.3    | 13 | 7 / 2.3  | N/A   |
| CLBS-7-14    | 7 / 2.3    | 14 | 7 / 2.3  | N/A   |
| CLBS-7-15    | 7 / 2.3    | 15 | 7 / 2.3  | N/A   |
| CLBS-7-16    | 7 / 2.3    | 16 | 7 / 2.3  | N/A   |
| CLBS-7-17    | 7 / 2.3    | 17 | 7 / 2.3  | N/A   |
| CLBS-7-18    | 7 / 2.3    | 18 | 7 / 2.3  | N/A   |

|             |           |    |           |         |
|-------------|-----------|----|-----------|---------|
| CLBS-7-3    | 7 / 2.3   | 3  | 7 / 2.3   | N/A     |
| CLBS-7-4    | 7 / 2.3   | 4  | 7 / 2.3   | N/A     |
| CLBS-7-5    | 7 / 2.3   | 5  | 7 / 2.3   | N/A     |
| CLBS-7-6    | 7 / 2.3   | 6  | 7 / 2.3   | N/A     |
| CLBS-7-7    | 7 / 2.3   | 7  | 7 / 2.3   | N/A     |
| CLBS-7-8    | 7 / 2.3   | 8  | 7 / 2.3   | N/A     |
| CLBS-7-9    | 7 / 2.3   | 9  | 7 / 2.3   | N/A     |
| CLBS-8.5-10 | 8.5 / 2.8 | 10 | 8.5 / 2.8 | 5 / 1.7 |
| CLBS-8.5-11 | 8.5 / 2.8 | 11 | 8.5 / 2.8 | 5 / 1.7 |
| CLBS-8.5-12 | 8.5 / 2.8 | 12 | 8.5 / 2.8 | 5 / 1.7 |
| CLBS-8.5-13 | 8.5 / 2.8 | 13 | 8.5 / 2.8 | 5 / 1.7 |
| CLBS-8.5-14 | 8.5 / 2.8 | 14 | 8.5 / 2.8 | 5 / 1.7 |
| CLBS-8.5-15 | 8.5 / 2.8 | 15 | 8.5 / 2.8 | 5 / 1.7 |
| CLBS-8.5-16 | 8.5 / 2.8 | 16 | 8.5 / 2.8 | 5 / 1.7 |
| CLBS-8.5-17 | 8.5 / 2.8 | 17 | 8.5 / 2.8 | 5 / 1.7 |
| CLBS-8.5-18 | 8.5 / 2.8 | 18 | 8.5 / 2.8 | 5 / 1.7 |
| CLBS-8.5-3  | 8.5 / 2.8 | 3  | 8.5 / 2.8 | 5 / 1.7 |
| CLBS-8.5-5  | 8.5 / 2.8 | 5  | 8.5 / 2.8 | 5 / 1.7 |
| CLBS-8.5-6  | 8.5 / 2.8 | 6  | 8.5 / 2.8 | 5 / 1.7 |
| CLBS-8.5-7  | 8.5 / 2.8 | 7  | 8.5 / 2.8 | 5 / 1.7 |
| CLBS-8.5-8  | 8.5 / 2.8 | 8  | 8.5 / 2.8 | 5 / 1.7 |
| CLBS-8.5-9  | 8.5 / 2.8 | 9  | 8.5 / 2.8 | 5 / 1.7 |

This product typically contains a stent with positioning sleeve, guiding catheter length 320 cm and pushing catheter length 170 cm (except the CLBS-7-XX and CLBS-5-XX, which does not contain guiding catheter). The stent has a proximal and a distal flap. It also has a tapered tip in larger French sizes at the proximal end to facilitate smooth cannulation. The product code for the

Cotton-Leung® Biliary Stent Set is CLBS-X-Y, where X denotes the French (FR) size and Y denotes the length in centimeters (cm).

### **INTENDED USE**

This device is used to drain obstructed biliary ducts.

### **NOTES**

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Do not use this device for any purpose other than stated intended use.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Notify Cook for return authorization.

This device must be stored in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

For 10 FR Sof-Flex® stents, lubricate the guiding catheter with water-soluble lubricant.

Retain the positioning sleeve (if any)\* for use when introducing the stent flaps into the accessory channel.

### **CONTRAINDICATIONS**

Those specific to ERCP and any procedure to be performed in conjunction with stent placement. Inability to pass wire guide or stent through obstructed area.

### **POTENTIAL COMPLICATIONS**

Those associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Those associated with biliary stent placement include, but are not limited to: trauma to the biliary tract or duodenum, obstruction of the pancreatic duct, stent migration.

### **PRECAUTIONS**

Refer to package label for minimum channel size required for this device.

This stent must be placed under fluoroscopic monitoring.

The positioning sleeve\* or the pigtail straightener\* is not intended for use in the accessory channel of the endoscope.

Care must be exercised when straightening the pigtail curls\* in order to avoid kinking or breaking the stent.

The tapered tip end of the stent\* or side holes\* must be positioned in the common bile duct while the other end remains in the duodenum.

The longer flaps of the ST-2 Biliary Stent should be positioned in the duct while the shorter flaps remain in the duodenum. (See fig. 1)

This device should not be left indwelling for more than three months or as directed by a physician. Periodic evaluation is recommended.

Select the Cook stent introducer system (if not included) in the appropriate French size.

### **Warning**

For stents with radiopaque bands, MRI compatibility has not been established.

### **INSTRUCTIONS FOR USE**

#### **I. FOR NON-PIGTAIL STENTS:**

1. Gently ensure full extension of all side flaps\*.

2. Load positioning sleeve onto duodenal flap end of stent.

For 5 to 7 FR stents:

3.a. Introduce stent, tapered tip first, and positioning sleeve onto a pre-positioned wire guide.

For 8.5 FR and larger stents:

3.b. Remove Tuohy-Borst adapter from end of guiding catheter, then introduce guiding catheter into accessory channel over a pre-positioned wire guide.

3.c. Introduce stent, tapered tip first, and positioning sleeve onto guiding catheter and pre-positioned wire guide.

THEN REFER TO STEPS 4-7 BELOW

## II. FOR PIGTAIL STENTS:

1. Advance pigtail straightener along shaft of stent to tapered tip end in order to straighten pigtail curl. (See fig. 2)

2. Introduce stent, tapered tip first, and pigtail straightener onto pre-positioned wire guide until straightener reaches second curl.

3. Advance pushing catheter over wire guide to advance pigtail stent into accessory channel.

THEN REFER TO STEPS 4-7 BELOW

4. Advance guiding catheter\* and/or pushing catheter in 1-2 cm increments until stent is in desired position.

5. Fluoroscopically and endoscopically confirm desired stent position.



**Note (for non-pigtail stents only):** Replace Tuohy-Borst adapter onto end of guiding catheter\*. Inject contrast, if desired, to fluoroscopically visualize stent position.

6. After confirming stent position, gently remove wire guide, then guiding catheter\* from endoscope while maintaining position of stent with pushing catheter.

7. Gently remove pushing catheter from accessory channel.

**Upon completion of procedure, dispose of introducer components per institutional guidelines for biohazardous medical waste.**

\* (if any)

Solus® Double Pigtail Stent with Introducer (ZSS-10-10-RB; ZSS-10-12-RB; ZSS-10-15-RB; ZSS-10-1-RB; ZSS-10-3-RB; ZSS-10-4-RB; ZSS-10-5-RB; ZSS-10-7-RB; ZSS-10-9-RB)

| RPN          | Stent diameter, Fr /mm | Stent length, cm | Pushing catheter diameter, Fr /mm | Guiding catheter diameter, Fr /mm |
|--------------|------------------------|------------------|-----------------------------------|-----------------------------------|
| ZSS-10-10-RB | 10 / 3.3               | 10               | 10 / 3.3                          | 5 / 1.65                          |
| ZSS-10-12-RB | 10 / 3.3               | 12               | 10 / 3.3                          | 5 / 1.65                          |
| ZSS-10-15-RB | 10 / 3.3               | 15               | 10 / 3.3                          | 5 / 1.65                          |
| ZSS-10-1-RB  | 10 / 3.3               | 1                | 10 / 3.3                          | 5 / 1.65                          |
| ZSS-10-3-RB  | 10 / 3.3               | 3                | 10 / 3.3                          | 5 / 1.65                          |
| ZSS-10-4-RB  | 10 / 3.3               | 4                | 10 / 3.3                          | 5 / 1.65                          |
| ZSS-10-5-RB  | 10 / 3.3               | 5                | 10 / 3.3                          | 5 / 1.65                          |
| ZSS-10-7-RB  | 10 / 3.3               | 7                | 10 / 3.3                          | 5 / 1.65                          |
| ZSS-10-9-RB  | 10 / 3.3               | 9                | 10 / 3.3                          | 5 / 1.65                          |

#### **INTENDED USE**

This device is used to drain obstructed biliary ducts.

#### **NOTES**

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

#### **CONTRAINDICATIONS**

Those specific to ERCP and any procedure to be performed in conjunction with stent placement. Inability to pass wire guide or stent through obstructed area.

#### **POTENTIAL COMPLICATIONS**

Those associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Additional complications associated with biliary stent placement include, but are not limited to: trauma to biliary tract or duodenum, obstruction of pancreatic duct, stent migration.

#### **PRECAUTIONS**

Refer to package label for minimum channel size required for this device.

A complete diagnostic evaluation must be performed prior to use to determine proper stent size.

This stent must be placed under fluoroscopic monitoring.

This stent must only be placed using the Solus Stent Introduction System, which is included with this stent set.

This stent should not be left indwelling for more than three months or as directed by a physician. Periodic evaluation is recommended.

Pigtail straightener is used to straighten pigtail curls on stent and to aid introduction of stent into accessory channel. Pigtail straightener is not intended for use in accessory channel.

Assessment must be made to determine necessity of sphincterotomy prior to stent placement. If sphincterotomy is required, all appropriate cautions, warnings and contraindications must be observed.


### WARNING

MRI compatibility of this device has not been established.

### INSTRUCTIONS FOR USE

1. Advance pigtail straightener, if applicable, along shaft of stent to end. (See fig. 1) **Note:** Care must be exercised when straightening pigtail curls in order to avoid kinking or breaking the stent.
2. Lubricate guiding catheter with water-soluble lubricant. Back load stent and pigtail straightener onto tip of soft stent introducer. Advance preloaded stent and soft stent introducer over a pre-positioned wire guide.
3. Slide pigtail straightener, if applicable, back along stent, straightening second pigtail curl. Advance stent and remove pigtail straightener as stent enters accessory channel. **Caution:** Pigtail straightener is not intended for use in accessory channel.
4. Fluoroscopically monitor radiopaque bands on tip of guiding catheter. Radiopaque bands on catheter are arranged in 5 cm increments. When a sufficient length of guiding catheter (a minimum of one radiopaque band) is above majority of stricture, disconnect Luer lock fitting on introducer. **Caution:** A sufficient portion of guiding catheter must be above stricture prior to stent positioning.
5. Maintain position of guiding catheter and advance stent into desired position using pushing catheter.
6. Fluoroscopically and endoscopically confirm desired stent placement.
7. Retract and remove both guiding catheter and wire guide into endoscope while maintaining position of pushing catheter.
8. Confirm position of stent once guiding catheter and wire guide are removed.
9. Remove pushing catheter.

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

## **COMPLIANCE WITH RELEVANT REGULATIONS AND TECHNICAL STANDARDS**

Please notify Cook for return authorization.  
Cook Medical Europe Ltd  
O'HALLORAN ROAD  
NATIONAL TECHNOLOGY PARK  
LIMERICK, IRELAND  
TEL: +353 61 334440 FAX: +353 61 334441  
[WWW.COOKMEDICAL.COM](http://WWW.COOKMEDICAL.COM)

## **COMPLIANCE WITH RELEVANT REGULATIONS AND TECHNICAL STANDARDS**

Cook Ireland Ltd manufactures all medical products in accordance with Medical Devices Directive 93/42/EEC.

All of Cook facilities are also certified to ISO 13485 (2003/2012) *Quality Management System*.

## **USAGE CONDITIONS**

Plastic Endoprostheses Stents to Maintain Patency of the Gastrointestinal Tract are used in the patient care medical institutions.

## **SHELF LIFE**

The "use by" date as denoted by symbols on the labeling. The shelf life is shown by indication the time of manufacture (factory symbol) to the stated expiration date (hourglass symbol). The difference between both dates is the shelf life.

Plastic Endoprostheses Stents to Maintain Patency of the Gastrointestinal Tract have a three (3) year shelf life.

## **STERILIZATION DETAILS**

Plastic Endoprostheses Stents to Maintain Patency of the Gastrointestinal Tract are supplied sterile in peel open packages utilizing an EtO (Ethylene Oxide) validated method to a sterility assurance level (SAL) of  $10^{-6}$ . These devices are intended for one-time use, sterile if unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile.

## **STORAGE RULES**

These devices should be stored in a dark, dry, cool place; avoid extended exposure to light.

## **TRANSPORTATION RULES**

Transportation shall account for the type of transportation employed, the destination, the vendor used for transportation and an estimation of the frequency and harshness of any handling required during transportation.

### **PACKAGING AND LABELING**

Plastic Endoprotheses Stents to Maintain Patency of the Gastrointestinal Tract are supplied in the individual packaging, which should show no signs of damage or opening. Upon detecting the bends, breaks and other defects the device cannot be used.


The label of the device shows the following device parameters:

- company-manufacturer;
- manufacturing address;
- address of the EU representative;
- device name (trade mark);
- set components;
- size of the stent delivery system;
- size of endoscope channel;
- size of wire guide;
- sterilization method;
- serial number (LOT);
- marking about CE certification;
- marking on reusability;
- marking for "Attention. Consult accompanied documentation"
- shelf life

# Oasis® One Action Stent Introduction System with Pre-Loaded Cotton Leung® Biliary Stent

REF OACL-11.5-12

REF Test

Rx ONLY

|            |         |   |   |
|------------|---------|---|---|
| STERILE EO | 2018-08 |   |   |
| LOT Test   | 2015-08 | ⊗ | Ⓢ |

CE 0123


# Fusion® Marathon™ Anti-Reflux Biliary Stent

REF FS-MAR-10-5


REF Test

BIBI00001

Rx ONLY

|            |         |   |   |
|------------|---------|---|---|
| STERILE EO | 2018-05 |   |   |
| LOT Test   | 2015-05 | ⊗ | Ⓢ |

CE 0123


Geenen®

Pancreatic Stent Set

REF GPDP-11.5-10

REF Test

|            |         |   |   |         |
|------------|---------|---|---|---------|
| STERILE EO | 2018-07 |   |   | Rx ONLY |
| LOT Test   | 2015-07 | ② | ③ |         |

Manufactured by:  
Cook Ireland Ltd,  
O'Halloran Road, National Technology Park,  
Limerick, Ireland.

CE 0123


Geenen®

Pancreatic Stent Set

REF GPDS-5-5

REF Test

|            |         |   |   |         |
|------------|---------|---|---|---------|
| STERILE EO | 2018-05 |   |   | Rx ONLY |
| LOT Test   | 2015-05 | ② | ③ |         |

Manufactured by:  
Cook Ireland Ltd,  
O'Halloran Road, National Technology Park,  
Limerick, Ireland.

CE 0123


Geenen®

Pancreatic Stent Set

REF GPDS-7-12

REF Test

|            |         |         |   |    |
|------------|---------|---------|---|----|
| STERILE EO | 2018-07 | 2015-07 | 2 | 18 |
| LOT Test   |         |         |   |    |

COOK MEDICAL  
Manufactured by:  
Cook Ireland Ltd,  
O'Halloran Road, National Technology Park,  
Limerick, Ireland.

CE 0123


Cotton-Huibregtse®

Biliary Stent Set

REF CHBS-11.5-12

REF Test

|            |         |         |   |    |
|------------|---------|---------|---|----|
| STERILE EO | 2018-05 | 2015-05 | 2 | 18 |
| LOT Test   |         |         |   |    |

COOK MEDICAL  
Manufactured by:  
Cook Ireland Ltd,  
O'Halloran Road, National Technology Park,  
Limerick, Ireland.

CE 0123


## Cotton-Leung® (Amsterdam)

Biliary Stent Set

REF CLBS-11.5-7

REF Test


## Solus®

Double Pigtail Stent with Introducer

REF ZSS-10-4-RB

REF Test

## DISPOSAL RULES

Upon completion of procedure, the devices are to be disposed of as per institutional guidelines for biohazardous medical waste.

Regulatory Affairs Manager

COOK MEDICAL EUROPE LTD. 28 Aug 15

O'Halloran Road, Date  
National Technology Park  
Limerick, Ireland.

Company Registration No. 492272  
VAT Registration No. 9776234K

COOK MEDICAL EUROPE LTD.

O'Halloran Road,  
National Technology Park  
Limerick, Ireland.

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
Manufactured by  
Cook Ireland Ltd.  
O'Halloran Road, National Technology Park,  
Limerick, Ireland.

CE 0123


# Fusion® Marathon™

Anti-Reflux Biliary Stent

REF FS-MAR-10-5

REF Test

REF FS-MAR-10-5

Rx ONLY

|            |         |   |   |
|------------|---------|---|---|
| STERILE EO | 2018-05 |   |   |
| LOT Test   | 2015-05 | ⊗ | ⓘ |

Manufactured by  
Cook Ireland Ltd.  
O'Halloran Road, National Technology Park,  
Limerick, Ireland.


CE 0123

# Oasis® One Action Stent Introduction System with Pre-Loaded ST-2 Tannenbaum® Biliary Stent

REF OATS-10-5

REF Test

Rx ONLY

|            |         |   |   |
|------------|---------|---|---|
| STERILE EO | 2018-05 |   |   |
| LOT Test   | 2015-05 | ⊗ | Ⓜ |

www.cookmedical.com

**COOK MEDICAL**  
Manufactured by:  
Cook Ireland Ltd,  
O'Halloran Road, National Technology Park,  
Limerick, Ireland.


CE 0123

# Oasis® One Action Stent Introduction System with Pre-Loaded ST-2 Tannenbaum® Biliary Stent

REF OATS-11.5-9

REF Test

Rx ONLY

|            |         |   |   |
|------------|---------|---|---|
| STERILE EO | 2018-05 |   |   |
| LOT Test   | 2015-05 | ⊗ | Ⓜ |

www.cookmedical.com

**COOK MEDICAL**  
Manufactured by:  
Cook Ireland Ltd,  
O'Halloran Road, National Technology Park,  
Limerick, Ireland.

CE 0123


Geenen®

Pancreatic Stent Set

REF GEPD-11.5-10

REF Test

|            |         |   |   |          |
|------------|---------|---|---|----------|
| STERILE EO | 2018-07 |   |   | REF ONLY |
| LOT Test   | 2015-07 | ② | ③ |          |

COOK MEDICAL

Manufactured by:  
Cook Ireland Ltd.  
O'Halloran Road, National Technology Park,  
Limerick, Ireland.


CE 0123

Geenen®

Pancreatic Stent Set

REF GPDS-5-5

REF Test

|            |         |   |   |          |
|------------|---------|---|---|----------|
| STERILE EO | 2018-05 |   |   | REF ONLY |
| LOT Test   | 2015-05 | ② | ③ |          |

COOK MEDICAL

Manufactured by:  
Cook Ireland Ltd.  
O'Halloran Road, National Technology Park,  
Limerick, Ireland.


CE 0123

Geenen®

Pancreatic Stent Set

REF GPDS-7-12

REF Test

|            |         |   |   |         |
|------------|---------|---|---|---------|
| STERILE EO | 2018-07 |   |   | Rx ONLY |
| LOT Test   | 2015-07 | 2 | 1 |         |

Manufactured by:  
Cook Ireland Ltd.,  
O'Halloran Road, National Technology Park,  
Limerick, Ireland.

CE 0123


Cotton-Huibregtse®

Biliary Stent Set

REF CHBS-11.5-12

REF Test

|            |         |   |   |         |
|------------|---------|---|---|---------|
| STERILE EO | 2018-05 |   |   | Rx ONLY |
| LOT Test   | 2015-05 | 2 | 1 |         |

Manufactured by:  
Cook Ireland Ltd.,  
O'Halloran Road, National Technology Park,  
Limerick, Ireland.

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


COOK MEDICAL EUROPE

Date

O'Halloran Road, National Technology Park


Limerick, Ireland.

Company Registration No. 402272

VAT Registration No. 9776234K

I certify that I am satisfied that this document  
is an original, duly signed by the person who  
had lawful authority to sign same,  
this 25<sup>th</sup> day of August 2015

  
Niall Sheehy, Notary Public,  
677 Glentworth St.  
Limerick.

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