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Uno™ Neurovascular Embolization System

DEVICE DESCRIPTION

The Uno™ Neurovascular Embolization System is a self-expanding vessel occlusion device. The Uno Device is oval-shaped, comprised of nitinol and secured at both ends with platinum marker bands. The device incorporates a PTFE partial cover. The proximal marker band attaches to a wire that pushes the device through a commercially available microcatheter to the intended treatment site. The "delivery wire" detaches from the Uno Device by electrolytic means after deployment with the Reverse Medical Detachment System.

The Uno System is packaged as a single unit with the Uno Device, introducer sheath and detachable delivery wire. The system is provided sterile, non-pyrogenic, and is intended for single use only.

The Uno System is specifically designed for use with the Reverse Medical Detachment System. This unit is sold separately.

Figure 1: Reverse Medical Uno Neurovascular Embolization System


Table 1. MVP System Dimensions

Model	A. Diameter	B. Length	Target Vessel Diameter
	Unconstrained	Unconstrained	
Uno-3	5.3mm	12mm	1.5 – 3.0mm
Uno-5	6.5mm	12mm	3.0 – 5.0mm

Table 2. Recommended Microcatheter

Model	Rec. Microcatheter ID	Rec. Microcatheter Max. Length
Uno-3	0.021"	150cm
Uno-5	0.027"	150cm

INDICATIONS FOR USE

The Reverse Medical Uno Neurovascular Embolization System is indicated to obstruct blood flow in the neurovasculature in vessels between 1.5 – 5 mm.

The Reverse Medical Corporation Uno System is designed for elective closure of vessels in the neurovasculature such as closure of a parent artery due to aneurysms after test occlusion; closure of an artery due to uncontrollable bleeding in the internal carotid or intracranial vertebral artery; arteriovenous fistulas; or uncontrollable bleeding in the skull base.

CONTRAINDICATIONS

None known.

COMPATIBILITY

Refer to product label for device dimensions. Refer to labeling provided with other medical technologies to determine compatibility.

WARNINGS

- The Uno System is not for placement into intracranial aneurysms.
- The safety and effectiveness of the Uno System has not been established for cardiac uses (e.g., cardiac septal occlusion, patent ductus arteriosus, paravalvular leak closure).

- The Uno Neurovascular Embolization System should only be used by physicians who are familiar with angiographic and interventional procedures. It is important to read the instructions for use prior to using this product.
- Use prior to the expiration date printed on the product packaging label.
- The Uno System is provided sterile and non-pyrogenic unless the unit package is opened or damaged. Do not use if the packaging is breached or damaged.
- Inspect the Uno System prior to use for any irregularities or damage and discard if any inconsistencies are observed.
- The Uno System is intended for single use only. After use, dispose of the device in accordance with hospital and/or local government policy. The device should not be cleaned, re-processed, re-sterilized, or reused as these processes may damage the components and surface of the device and compromise performance. In addition, the risk of infection of a reprocessed device has not been qualified in order to reduce this risk to an acceptable level. Structural integrity and/or function may be impaired by reuse or cleaning.
- The Uno System should be advanced or manipulated under fluoroscopic guidance. Do not advance or withdraw the device when excessive resistance is met until the cause of resistance is determined.
- Do not use a power injection syringe to inject contrast solution with the device through a microcatheter.
- Do not advance the delivery wire after detaching it from the Uno Device.
- Only use the battery operated Reverse Medical Detachment System as the detachment system power source.
- Physicians must be prepared to deal with urgent situations which require removal of the device. This includes the availability of an on-site surgeon.
- Do not use with catheters that are contraindicated for embolic devices (coils and/or plugs).

PRECAUTIONS

- The Uno Device consists of a nickel-titanium alloy that is generally considered safe. However, in-vitro testing demonstrated that nickel is released from the device for 60 days minimum. Patients who are allergic to nickel may have an allergic reaction to the device, especially patients with a history of metal allergies. Certain allergic reactions can be serious; patients should be instructed to notify their physicians immediately if they suspect they are experiencing an allergic reaction such as difficulty breathing or inflammation of the face or throat. Some patients may develop an allergy to nickel if the device is implanted.
- Physicians should exercise clinical judgment in situations that involve use of anticoagulants or antiplatelet drugs before, during and/or after use of the device. Please note care in specific populations: pregnancy – care should be taken to minimize radiation exposure to the mother and fetus, and nursing mothers – there has been no quantitative assessment of the presence of leachables in breast milk.
- Verify Uno System compatibility when using other ancillary devices in intravascular procedures. Physician must be familiar with percutaneous, intravascular techniques and possible complications associated with the procedure.
- The Uno System should be hydrated prior to use.
- Exercise care in handling the Uno System to reduce the chance of accidental damage.
- There is a potential for higher push forces when used with slotted nitinol composite microcatheters.
- Verify that the diameter of any accessory device that is used is compatible with the Uno System prior to use.
- Potential adverse events that may occur during or after a procedure placing this device include, but are not limited to:
 - Air embolus
 - Allergic reaction/toxic effects
 - Bleeding
 - Death
 - Device migration
 - Fever
 - Foreign material embolic event
 - Hemolysis
 - Infection
 - Occlusion of unintended vessel
 - Peripheral embolism
 - Recanalization
 - Residual flow
 - Stroke/TIA
 - Surgical intervention
 - Vascular access site complication
 - Vessel trauma/perforation

MATERIALS REQUIRED FOR USE WITH THE UNO SYSTEM

- Delivery of the Uno System requires the use of a Microcatheter (see Table 2 for recommended size).
- Reverse Medical Detachment System (ED2-RM) is required to detach the Uno Device in the desired location. Accessories required:
 - 9-volt alkaline battery
 - 20 or 22 gauge sterile needle
 - Reverse Medical Sterile Cable Set (RMCS-2.75)
- Heparinized Saline Solution
- Rotating Hemostatic Valve (RHV)
- 3-Way Stopcock
- Saline Bag Cuff
- Access Sheath
- Guide Catheter compatible with a Microcatheter recommended in Table 2
- Guide Wire compatible with the Guide Catheter and microcatheter
- Infusion Stand

MRI INFORMATION

MR Conditional

The Uno Device was determined to be MR conditional.

Non-clinical testing demonstrated that the Uno Device is MR Conditional. A patient with this device can be scanned safely, immediately after placement under the following conditions:

- Static magnetic field of 3-Tesla or less
- Maximum spatial gradient magnetic field of 9,000-Gauss/cm or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 4-W/kg for 15 minutes of scanning (i.e., per pulse sequence)
- First Level Controlled Operating Mode for the MR system

MRI-RELATED HEATING

In non-clinical testing, the Uno Device produced the following temperature rise during MRI performed for 15-min of scanning (i.e., per pulse sequence) in the 3-Tesla (3-Tesla/128-MHz, Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI) MR system:

Highest temperature change +1.7° C

Therefore, the MRI-related heating experiments for the Uno Device at 3-Tesla using a transmit/receive RF body coil at an MR system reported whole body averaged SAR of 2.9 - W/kg (i.e., associated with a calorimetry measured whole body averaged value of 2.7-W/kg) indicated that the greatest amount of heating that occurred in association with these specific conditions was equal to or less than +1.7° C.

ARTIFACT INFORMATION

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the Uno™ Device. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary. The maximum artifact size (i.e., as seen on the gradient echo pulse sequence) extends approximately 5-mm relative to the size and shape of the Uno Device.

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	70-mm ²	51-mm ²	187-mm ²	115 mm ²
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

PREPARATION FOR USE

- In order to achieve optimal performance of the Uno System and to reduce the risk of thromboembolic complication, a continuous saline flush should be maintained between a) access sheath and guide catheter, b) Microcatheter and guide catheter, and c) Microcatheter and guidewire and Uno System.
- Place the appropriate guide catheter following recommended procedure. Connect a rotating hemostatic valve (RHV) to the hub of the guide catheter. Attach a 3-way stopcock to the side arm of the RHV. Connect a heparinized saline line for continuous flush.
- Check all fittings so air is not introduced into guide catheter or Microcatheter during continuous flush.

DIRECTIONS FOR USE

- Access the targeted vessel and perform an angiogram using standard technique to measure the vessel diameter at the desired occlusion site. Ensure the targeted vessel meets the recommended diameter (see Table 1).
- Ensure that the occlusion site is long enough to accommodate the implanted Uno Device without obstructing unintended vessels.
- Select a Microcatheter (see Table 2 for recommended size) and prepare according to the manufacturer's instructions for use.
Note: the delivery Microcatheter length must be less than or equal to 150cm.
- Insert the guidewire and advance the Microcatheter over the guidewire until the distal tip is near the leading edge of the occlusion site.
- Remove the guidewire from the Microcatheter.
- Remove the Uno System from the dispenser tube by gently removing the yellow introducer sheath from the two outer clips and pull the Uno System out of the Clear Housing Tube. Gently grasp the introducer sheath proximal to the Uno Device, slowly remove the Uno Device, delivery wire and introducer sheath simultaneously from the dispenser tube. See Figure 2.


Figure 2. Packaged Uno System

- Inspect the Uno System thoroughly for damage. If damage exists, replace with a new Uno System.
- Flush the Uno System by gently immersing a few centimeters of the introducer sheath and the Uno Device in heparinized saline.
- Draw the Uno Device into the introducer sheath by pulling the delivery wire until the distal tip of the Uno Device is just covered by the introducer sheath.
Note: If the Uno Device is pulled too far into the introducer sheath it may not flush correctly; two holes in the introducer sheath located 3cm and 6cm proximal to the tip facilitate air purging and flushing.
- Insert the Uno System, by inserting the distal end of the introducer sheath through the RHV and into the hub of the Microcatheter until the sheath is firmly seated. Tighten the RHV to prevent back flow of blood, but not so tight as to damage the Uno System during its introduction into the Microcatheter.

- Usually verify the flushing solution is infusing normally. Saline drip should be seen exiting the proximal end of the introducer sheath.
2. Loosen the RHV enough to advance the Uno Device and delivery wire as a system, but not enough to allow back flow of blood. Remove the introducer sheath over the delivery wire's proximal end.
 3. Advance the delivery wire and Uno Device until the distal platinum marker of the Uno Device is aligned with the Microcatheter distal marker band.
 4. Confirm that the marker bands are positioned optimally, and deploy the Uno Device by maintaining forward pressure on the delivery wire and unsheath the Uno Device by pulling the Microcatheter back slowly.
 15. If the device position is unsatisfactory:
 - Stabilize the delivery wire and advance the Microcatheter distally until the Uno Device is recaptured within the catheter.
 - Reposition and deploy the device, or remove the device from the patient.
 - The device should not be resheathed and deployed more than a maximum of 3 times.
 16. The Uno Device should be just beyond the microcatheter; the proximal marker of the Uno Device should be positioned no more than 2-3 mm beyond the tip of the microcatheter for proper detachment.
 17. If the device position is satisfactory, proceed with detachment:
 - a. Insert a sterile needle (20 G or 22G) at the patient's shoulder (M. Deltoides) or groin.
 - b. Insert the bayonet plug of the sterile Cable Set into the terminal of the Detachment System (for detailed instructions on connecting the Cable Set to the Detachment System, see instructions for Use of the Reverse Medical Detachment System)
 - c. Clip the connector end of the sterile black cable onto the sterile needle.
 - d. Clip the connector end of the sterile red cable onto the proximal end of the delivery wire. Ensure that the delivery wire is resting on a dry, clean surface.
 - e. Confirm under fluoroscopy that the Microcatheter tip is proximal to the detachment zone and the proximal marker band of the Uno Device.
 - f. Switch "ON" the Detachment System. (For detailed instructions on starting a detachment, see instructions for Use of the Reverse Medical Detachment System.)
 - g. An acoustic signal and "Detach" Indicator flashes on the Detachment System, which indicates the detachment of the Uno Device.
 - h. Verify successful detachment fluoroscopically by slowly pulling the delivery wire proximally to make sure the Uno Device does not move. In the unlikely event the Uno Device moves proximally, advance the delivery wire distally to re-establish the Uno Device and Microcatheter marker alignment. Start the Detachment System again to resume the detachment process. Verify Uno Device detachment as above.
 - i. When Uno Device detachment has been detected and fluoroscopically confirmed, disconnect the red cable from the delivery wire. Remove the needle from the patient arm or groin.
 - j. Switch "OFF" the Detachment System.
 - k. Hold the Microcatheter in place, retract the delivery wire into the Microcatheter, and remove the delivery wire from the patient. Complete the procedure following standard technique.
 18. After completing the procedure, withdraw and discard the microcatheter, Cable Set, and applicable device accessories. Store the Reverse Medical Detachment System in a clean, dry and secure place. The Detachment System may be cleaned with a soft cloth dampened with 70% alcohol.
 19. Dispose of the battery in accordance with hospital, administrative, and/or local government policy. Remove battery when the Detachment System is not in use.

STORAGE

Avoid exposure to water, sunlight, extreme temperatures and high humidity during storage. Store the Uno System under controlled room temperature. See the product label for the device shelf life. Do not use the device beyond the labeled shelf life.

MATERIALS

The Uno System does not contain any latex or PVC materials.

WARRANTY

Reverse Medical Corp. (Reverse) warrants that reasonable care has been used in the design and manufacture of this device. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular use. Handling, storage, cleaning, and sterilization of this device as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond Reverse's control directly affect the device and the results obtained from its use. Reverse's obligation under this warranty is limited to the repair or replacement of this device and Reverse shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this device. Reverse neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device. Reverse assumes no liability with respect to devices reused, reprocessed, or resterilized and makes no warranties, expressed or implied, including but not limited to merchantability or fitness for intended use, with respect to such device.

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