

ERMA Inc.	Technical file	Document number
	Medical device	Revision: 1 Date:2015-06-01
	Patho Cutter	Page 1/9

Instruction for use

Patho Cutter

Medical device, types R 22⁰, R 35⁰, I, II, and HP

	Date	Name and position	Signature
Prepared by	May 1, 2015	Hisao Nakashima	
Verified by	May 26 , 2015	Yutaka Namiki (General Manager)	
Approved by	June 1, 2015	Hiroshi Shimosaka (President)	

ERMA Inc.	Technical file	Document number
	Medical device	Revision: 1 Date: 2015-06-01
	Patho Cutter	Page 2/9

Table of contents		
Section	Description	Page No.
1	Administrative information	3
2	Device family information	3
3	Device component information	6
4	Packaging information	6
5	Manufacturing information	6
6	Sterilization	7
7	The procedure of recycling	7
8	Requirements for environmental protection	7
9	Device testing	7
10	Results of the risk analysis	8
11	Applicable technical standards	8
12	Revision history	8
13	Project marking	8
14	Information about the authorized representative in the Russian Federation	8
15	Marking	9

ERMA Inc.	Technical file	Document number
	Medical device	Revision: 1 Date: 2015-06-01
	Patho Cutter	Page 3/9

1. Administrative information

1.1. Company name:

Name: ERMA Inc.

Address: 2-31-6 Yushima, Bunkyo-ku, Tokyo 113-0034; Japan

Phone: 81-3-3818-62-81

Fax: 81-3-3813-7301

E-mail: trade@erma.co.jp

Website: www.erma.co.jp

Commercial reg. No.:

ISO cert. numbers: EN ISO 13485:2012 ISO

1.2. Device classification: "Other" than those covered by Annex II

1.3. Notified body: Kai Europe GMBH

Kottendorfer Strasse 5, 42697 Solingen, Germany

2. Device family information

2.1. General product description:

This product was developed as the flat wedged shape blade with the Very sharp blade tip. Normally, this product is used shaving the paraffin-fixed anatomy to the super thin section. This product is single use only.

Main components

1. Blade (Stainless steel: 30KH13 <equivalent to X30Cr13 of ISO)
2. Plastic case (polycarbonate)
3. Inner box for blade in plastic case

ERMA Inc.	Technical file	Document number
	Medical device	Revision: 1 Date:2015-65-01
	Patho Cutter	Page 4/9

2.2. Patho Cutter product range

Product	Design	Cat. No.	Product quantity
Patho Cutter R 35 ⁰	As per above picture	08-636-0	1 pcs/10 pcs/120 pcs
Patho Cutter R 22 ⁰	As per above picture	08-637-0	1 pcs/10 pcs/120 pcs
Patho Cutter I	As per above picture	08-635-0	1 pcs/10 pcs/120 pcs
Patho Cutter II	As per above picture	08-635-1	1 pcs/10 pcs/120 pcs
Patho Cutter HP	As per above picture	08-640-0	1 pcs/10 pcs/60 pcs

2.3.Intended use

Patho Cutter can be used for shaving the paraffin-fixed anatomy to the super thin section.

2.4. The method of use, instruction of use

Patho Cutter can be attached to the corresponding blade holder and then attached to microtome to shave the paraffin- fixed anatomy to the super thin section.

ERMA Inc.	Technical file	Document number
	Medical device	Revision: 1 Date:2015-06-01
	Patho Cutter	Page 5/9

2.5.Shelf life

Expiration date has not been set.

2.6.Preservation and Storage Method

The blade should be preserved in room temperature (1 to 30 degree C) without high humidity (more than 85%) in the plastic case. Direct sunlight and water soaking must be avoided. Preserve without load in clean and dry place. Don't preserve in the place where chemical is placed and gas occurs.

2.7. Usage Environment

The blade should be used in room temperature (1 to 30 degree C) without high humidity (more than 85%)

2.8. Durability

Under the normal usage, the average number of paraffin block to be produced per blade in good quality is as follows for each model.

	Hard Tissue	Soft Tissue
Patho Cutter R35	15.5	32.0
Patho Cutter R22	14.0	33.2
Patho Cutter I	10.2	22.4
Patho Cutter II	12.5	19.0
Patho Cutter HP	18.2	17.6

2.9.Hardness of material (HV)

R35	R22	I	II	HP
640	645	663	676	629

2.10. Application

Patho Cutter is applicable to any type of micrtome

2.11. Warranty Period

1 year after manufacturing date

ERMA Inc.	Technical file	Document number
	Medical device	Revision: 1 Date:2015-06-01
	Patho Cutter	Page 6/9

3. Device Component information

Component	Material
Blade	Special Stainless
Case	Plastic

4. Packaging – transport packaging

1 cartridge (11.2x2.9x3.2cm) (10.2x4.1x3.9cm) for HP only	Contains 50 blades
1 box (11.4 x 15 x 6.6cm) (20.9x10.3x8cm) for HP only	Contains 10 cartridge
1 shipping box (32 x 25.4 x 22.8cm) (33.5 x 22.3 x 18.4cm) for HP only	Contains 12 boxes 6 boxes for HP only

The export double carton and polyethelene packing material is used.

Products can be transported by any means of transport at a temperature of - 25 °C to + 50 °C and humidity from 10% to 95%.

ERMA Inc.	Technical file	Document number
	Medical device	Revision: 1 Date:2015-06-01
	Patho Cutter	Page 7/9

5. Manufacturing Process

The manufacture process of the Patho Cutter is as described as next.

- 1, Pressing; (Punching operation with machinery to make the shape of blade.)
- 2, Heat treatment
- 3, Edge grinding; (Making the ideal shape of the blade edge.)
- 4, Cleaning; (Remove remaining substances and oil.)
- 5, Plasma Ion process; (Plasma Ion Process to make the blade edge harden.
for Patho Cutter R35 and R22 only)
- 6, Coating (Chromium alloy coating and PTFE coating to reduce the friction during
slicing and to improve the durability of blade)
- 7, Inspection; (Quality control.)

6. Sterilization

Microtome blades are non-sterile product and single use only.

7. The procedure of recycling

All the materials of the prodcut are recyclable. But the product must be disposed of in accordance with national standards, taking into account a possible biohazard samples.

8. Requirements for environmental protection

Pollution-free materials are used.

The emission of green effect gas during production process is monitored and restricted by manufacturer's standard.

The amount of the waste disposal is reported by manufacturer's internal audit.

9. Device testing

The quality test for material is carried out to avoid the deformity of the product. During the manufacture ring process, the drop test, visual test and other quality tests are carried out to maintain the quality of blade. These tests are carried out based on the manufacturer's own standard.

ERMA Inc.	Technical file	Document number
	Medical device	Revision: 1 Date:2015-06-01
	Patho Cutter	Page 8/9

10. Results of the risk analysis

Patho Cutter is an IVD product and non-invasive as long as it is used for its own purpose. But the improvement for the structure of plastic case lessen the risk of injury caused by the misuse of blade.

11. Applicable technical standards

The safety and quality standard for Patho Cutter is based on the manufacture's criteria.

12. Revision history

Date	Notification No.	Item to be changed	Description
05/10/13	13-28&29&30&31	Patho Cutter I Patho Cutter II Patho Cutter R35 Patho Cutter R22 Patho Cutter HP	The packing method
05/08/14	14-59&60&61	Patho Cutter I Patho Cutter II Patho Cutter R35 Patho Cutter R22 Patho Cutter HP	The quantity of box per shipping box has been changed from 10 to 12

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14. Information about the authorized representative in the Russian Federation

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ERMA Inc.	Technical file	Document number
	Medical device	Revision: 1 Date:2015-06-01
	Patho Cutter	Page 9/9

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CERTIFICATE OF SIGNATURE

We hereby certify that the signature of the under-mentioned person, affixed to the attached document, appears to be genuine and authentic based on our record.

Title of attached document: INSTRUCTION FOR USE PATHO CUTTER

Date of attached document: JUNE 1, 2015

Name of person: HIROSHI SHIMOSAKA

Title: PRESIDENT

Name of company: ERMA INC.

THE TOKYO CHAMBER OF COMMERCE AND INDUSTRY

Satoru Ichinokawa
Authorized signatory

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APOSTILLE

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1. Country: JAPAN
This public document
2. has been signed by Yutaka Hokura
3. acting in the capacity of Notary of the Tokyo Legal Affairs Bureau
4. bears the seal/stamp of Yutaka Hokura, Notary
Certified
5. at Tokyo
6. June 24, 2015
7. by the Ministry of Foreign Affairs
8. 15-N^o028325
9. Seal/stamp:
10. Signature

A. Ogawa

Ayako OGAWA

For the Minister for Foreign Affairs

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