1. INDICATION FOR USE

Asahi Leocceed™-H series dialyzer (hereafter called Leocceed-H) is intended for use in hemodialysis (HD) for the treatment of patients suffering from acute or chronic renal failure.

2. CONTRAINDICATIONS

None known. In general, contraindications for HD are applicable.

3. ADVERSE REACTIONS

- Consult a responsible physician immediately whenever an adverse reaction occurs.
- Some patients may experience hypersensitivity reactions during treatment. Symptoms and signs include acute shortness of breath with wheezing; respiratory arrest; itching; flushing; mood discomfort; abnormal sweating; hives; generalized or localized skin redness; edema of the face, hands, or feet; hypertension; hypotension; elevated pulse rates; arrhythmia; ocular hyperemia; hypoesthesia; fever; leukopenia; and thrombocytopenia.
- Side effects such as hypertension, hypotension, headache, chest pain, stomachache, fatigue and nausea which may be associated with hypovolemia or hypervolemia can usually be avoided by careful management of the patient's fluid volume and electrolyte balance as well as the blood flow rate and transmembrane pressure (TMP).

4. WARNINGS

Follow these instructions and those of the HD machine supplier. Improper use could cause injury or death. The manufacturer will not be responsible for patient safety if these instructions and the HD machine supplier instructions are not followed. Make sure these instructions are available to the operator at the treatment site.

- Leocceed-H is intended for single use only. Reprocessing of Leocceed-H may lead to adverse reactions for the patient and/or device failure.
- The user must use Leocceed-H before the expiration date. Refer to the Leocceed-H product label for the expiration date.
- Leocceed-H must be kept in a clean, dry environment away from direct sunlight, heat, and humidity and at temperatures of 0-30°C (32-86°F).
- Leocceed-H is a plastic product, and must not be exposed to vibration or impact during transportation, storage or handling.
- Leocceed-H must not be exposed to organic solvents such as alcohol to prevent cracking or deformation.
- Inspect each Leocceed-H before use, and replace it if it shows signs of damage, damage or improper sealing of the sterilization bag, or misaligned or separated stoppers.
- Ensure that the dialysate and blood compartments of Leocceed-H are rinsed before use according to the priming procedure described below.
- Use aseptic technique throughout set-up, priming and treatment. Do not touch surfaces covered by protective stoppers.
- The saline must be a physiological solution that is sterile and non-pyrogenic.
- Due to its high ultrafiltration coefficient, Leocceed-H must only be used with HD machines equipped with a UF controller (Refer to the instructions of the HD machine).
- Use only physiologic fluids prescribed by a physician with the dialyzer. Fluids must conform to local regulations, standards, laws, or ISO water quality standards for dialysis. Refer to fluid labeling for complete instructions.
- Ensure that the sensors, detectors and alarm system of the HD machine are in proper operating condition including the venous air sensor and the blood leak detector in accordance with the instructions of the HD machine before starting treatment.
- Ensure that air is completely purged from Leocceed-H and blood lines during priming; failure in purging air may result in blood clotting during treatment. Further, ensure that no air enters the patient's blood vessel during blood return.
- During priming and treatment, verify that there is no leakage from the dialyzer, blood lines, blood ports, and their connecting parts. If leakage is observed, replace the dialyzer or the blood lines.
- If dialysate is to be used instead of saline, care must be taken to ensure its purity in accordance with the applicable water quality standards (for example, dialysate intended for infusion defined in ISO 11663:2009).
- In the event of a problem during treatment, such as blood leakage or coagulation, immediately discontinue the treatment under the direction of a physician and replace Leocceed-H with a new primed dialyzer.

- Leocceed-H has multiple connection points. Failure to make the proper connections can result in compromised treatment, injury, or death due to blood or therapy fluid leak.
- Leocceed-H may remove medications given in the arterial blood line. Give medications post dialyzer unless instructed otherwise.
- Do not permit the TMP to exceed 600 mmHg (80 kPa) with TMP being defined as per ISO 8637:2010.
- The recommended ranges of the treatment parameters specified in the APPENDIX should be considered when planning treatment for patients.
- Dispose the used Leocceed-H and the blood lines in strict accordance with official and institutional standards for medical waste disposal.

5. SPECIFICATIONS & PERFORMANCE

- Membrane: Polysulfone / PVP
- Housing: Polyethylene terephthalate glycol
- Headers: Acrylonitrile butadiene styrene
- Stoppers: Polyethylene
- Potting material: Polyurethane

Sterilized by gamma rays

Fluid pathways are sterile and non-pyrogenic.

Refer to the APPENDIX for design specifications and *in vitro* performance information regarding individual Leocceed-H models.

6. ANTICOAGULATION

Anticoagulant requirements will vary with the patient's condition; the type, quantity, and method of anticoagulant administration should be in accordance with the directions by a physician who is familiar with the patient's condition.

7. OPERATION

1. PRIMING

- 1) Place Leocceed-H vertically in the holder, with the red header (arterial end) down.
 - 2) Connect the arterial and venous blood lines to the blood inlet and outlet ports. Connect the arterial blood line to the saline bag for infusion. Rinse the hollow fibers and the blood lines with at least 500 mL saline at approx. 100 mL/min. Discard the first 500 mL primed saline. If any air is trapped in the dialyzer or blood lines, introduce additional saline for rinsing while gently tapping the dialyzer with one's hands or pinching and releasing the blood lines until the air is completely removed. Do not attempt to remove the air by tapping the dialyzer with hard instruments, such as forceps, because this may damage the dialyzer.
 - 3) Connect the dialysate lines to the dialysate ports. Allow the dialysate to flow at approx. 500 mL/min for 2 min or more. Ensure that the dialysate flows upwards to facilitate removal of air from the dialysate compartment. Optionally, the saline can be recirculated to remove air after the above priming procedure.
 - 4) Discard primed or recirculated saline before connecting to the patient. Fill the entire hemodialysis circuit with fresh saline prior to connecting to the patient. Do not infuse the primed or recirculated saline into the patient.
- Instead of saline, dialysate intended for infusion may be used for priming and ending the treatment.

2. STARTING TREATMENT

- 1) Connect the arterial blood line to the patient's arterial cannula. Start blood pump (approx. 100 mL/min) to displace saline.
- 2) Connect the venous blood line to the venous cannula.

3. ENDING TREATMENT

- 1) Adjust the UF flow rate and stop the blood pump.
- 2) Remove the arterial cannula from the patient and connect the arterial blood line to the saline bag. Feed appropriate volume of saline into the hollow fibers and the blood lines by starting the blood pump slowly to rinse back patient's blood.
- 3) Clamp the venous blood line and remove the venous cannula from the patient.

8. LIMITED LIABILITY

The manufacturer will not be responsible for any warranty given by independent distributors and will not be liable for any damage caused by the reuse, misuse, modification or improper handling of Leocceed-H, non-compliance with warnings and instructions supplied by the manufacturer, events occurring after the product is released, failure to ensure the product is in proper condition before use, or failure to use Leocceed-H with compatible medical devices.

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1. INDICATION FOR USE

Asahi Leocceed™-N series dialyzer (hereafter called Leocceed-N) is intended for use in hemodialysis (HD) for the treatment of patients suffering from acute or chronic renal failure.

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- Leocceed-N is intended for single use only. Reprocessing of Leocceed-N may lead to adverse reactions for the patient and/or device failure.

- The user must use Leocceed-N before the expiration date. Refer to the Leocceed-N product label for the expiration date.

- Leocceed-N must be kept in a clean, dry environment away from direct sunlight, heat, and humidity and at temperatures of 0-30°C (32-86°F).

- Leocceed-N is a plastic product, and must not be exposed to vibration or impact during transportation, storage or handling.

- Leocceed-N must not be exposed to organic solvents such as alcohol to prevent cracking or deformation.

- Inspect each Leocceed-N before use, and replace it if it shows signs of damage, damage or improper sealing of the sterilization bag, or misaligned or separated stoppers.

- Ensure that the dialysate and blood compartments of Leocceed-N are rinsed before use according to the priming procedure described below.

- Use aseptic technique throughout set-up, priming and treatment. Do not touch surfaces covered by protective stoppers.

- The saline must be a physiological solution that is sterile and non-pyrogenic.

- Due to its high ultrafiltration coefficient, Leocceed-N must only be used with HD machines equipped with a UF controller (Refer to the instructions of the HD machine).

- Use only physiologic fluids prescribed by a physician with the dialyzer. Fluids must conform to local regulations, standards, laws, or ISO water quality standards for dialysis. Refer to fluid labeling for complete instructions.

- Ensure that the sensors, detectors and alarm system of the HD machine are in proper operating condition including the venous air sensor and the blood leak detector in accordance with the instructions of the HD machine before starting treatment.

- Ensure that air is completely purged from Leocceed-N and blood lines during priming; failure in purging air may result in blood clotting during treatment. Further, ensure that no air enters the patient's blood vessel during blood return.

- During priming and treatment, verify that there is no leakage from the dialyzer, blood lines, blood ports, and their connecting parts. If leakage is observed, replace the dialyzer or the blood lines.

- If dialysate is to be used instead of saline, care must be taken to ensure its purity in accordance with the applicable water quality standards (for example, dialysate intended for infusion defined in ISO 11663:2009).

- In the event of a problem during treatment, such as blood leakage or coagulation, immediately discontinue the treatment under the direction of a physician and replace Leocceed-N with a new primed dialyzer.

- Leocceed-N has multiple connection points. Failure to make the proper connections can result in compromised treatment, injury, or death due to blood or therapy fluid leak.

- Leocceed-N may remove medications given in the arterial blood line. Give medications post dialyzer unless instructed otherwise.

- Do not permit the TMP to exceed 600 mmHg (80 kPa) with TMP being defined as per ISO 8637:2010.

- The recommended ranges of the treatment parameters specified in the APPENDIX should be considered when planning treatment for patients.

- Dispose the used Leocceed-N and the blood lines in strict accordance with official and institutional standards for medical waste disposal.

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Refer to the APPENDIX for design specifications and *in vitro* performance information regarding individual Leocceed-N models.

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7. OPERATION

1. PRIMING

- Place Leocceed-N vertically in the holder, with the red header (arterial end) down.

- Connect the arterial and venous blood lines to the blood inlet and outlet ports. Connect the arterial blood line to the saline bag for infusion. Rinse the hollow fibers and the blood lines with at least 500 mL saline at approx. 100 mL/min. Discard the first 500 mL primed saline. If any air is trapped in the dialyzer or blood lines, introduce additional saline for rinsing while gently tapping the dialyzer with one's hands or pinching and releasing the blood lines until the air is completely removed. Do not attempt to remove the air by tapping the dialyzer with hard instruments, such as forceps, because this may damage the dialyzer.

- Connect the dialysate lines to the dialysate ports. Allow the dialysate to flow at approx. 500 mL/min for 2 min or more. Ensure that the dialysate flows upwards to facilitate removal of air from the dialysate compartment. Optionally, the saline can be recirculated to remove air after the above priming procedure.

- Discard primed or recirculated saline before connecting to the patient. Fill the entire hemodialysis circuit with fresh saline prior to connecting to the patient. Do not infuse the primed or recirculated saline into the patient.

- Instead of saline, dialysate intended for infusion may be used for priming and ending the treatment.

2. STARTING TREATMENT

- Connect the arterial blood line to the patient's arterial cannula. Start blood pump (approx. 100 mL/min) to displace saline.

- Connect the venous blood line to the venous cannula.

3. ENDING TREATMENT

- Adjust the UF flow rate and stop the blood pump.

- Remove the arterial cannula from the patient and connect the arterial blood line to the saline bag. Feed appropriate volume of saline into the hollow fibers and the blood lines by starting the blood pump slowly to rinse back patient's blood.

- Clamp the venous blood line and remove the venous cannula from the patient.

8. LIMITED LIABILITY

The manufacturer will not be responsible for any warranty given by independent distributors and will not be liable for any damage caused by the reuse, misuse, modification or improper handling of Leocceed-N, non-compliance with warnings and instructions supplied by the manufacturer, events occurring after the product is released, failure to ensure the product is in proper condition before use, or failure to use Leocceed-N with compatible medical devices.

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