Biosealer CR 4

including

QUAD plus sealer
Automatic Sealing Handle
Bench Unit
Manual Handle
Ergonomic Manual Handle

User's Manual
and
Technical reference

Ljungberg & Kögel AB

This page is consciously left blank

1 Important

This user's manual is written for the person responsible for the operation of Biosealer CR 4. The operational methods and routines are developed and tested to ensure a reliable, safe and effective operation of CR4. It is important that the operator has studied and understood the contents of this manual before using CR4.

- This instrument is a sealing equipment using radio frequency for welding and the electrical emission at the operating frequency 40.68 MHz is high. Make sure that other instruments and equipment near the sealing unit can withstand this emission.
- Never touch the electrodes with your fingers during the sealing period! This may cause burn damages!
- Ensure that the PVC tube is dry on the surface.

As indicated by the -symbol, CR4 contains electronic components which are classified as dangerous for the environment and the equipment must be left to a recycling station when disposed. Local regulation decides how to recycle the CR4, but it is preferred to disassemble the CR4 in its main parts and recycle the metal parts, plastic part, electronics etc.

2 Warranty

- (1) The manufacturer, hereby guarantees the original buyer that Biosealer CR 4 is manufactured in a professional and quality manner, and will be free from all faults during a period of one year from the date of delivery from the manufacturer.
- (2) The warranty includes equipment or components that prove to have faults during the warranty period. The manufacturer will without cost for the customer, repair or replace the equipment that is faulty.
- (3) The warranty is not valid if the equipment has been repaired by anyone else than qualified personnel, that is approved by the manufacturer.
- (4) The warranty is not valid if the equipment has been changed in any way that according to the manufacturers opinion, affects the reliability or stability of the instrument.
- (5) The warranty is not valid when the serial number has been changed, crossed over or been removed, or if the fault has been caused by misuse or abnormal use.
- (6) In these cases the manufacturer or the manufacturer's representative will inform the customer about the decision, and if wished by the client will repair the equipment for normal rate. An estimated price can be given on request.

Abelko Innovation is committed to develop high-quality equipment and technical services to all our customers. We welcome any inputs on technical issues that are encountered so that they can be resolved quickly and in the most appropriate manner. Please submit your comments/feedbacks through your local distributors or alternatively email us directly at info@abelko.se

3 Table of contents

1	Important	3
2	Warranty	4
3	Table of contents	5
4	Properties.....	7
5	Description	7
5.1	Power Unit.....	7
5.2	Bench Unit.....	8
5.3	Handle Unit.	9
5.4	Automatic Sealing Handle (ASH).....	9
5.5	QUAD plus Unit.....	10
5.6	Ergonomic Sealing Handle (MSH-II).....	11
6	Installation	11
6.1	Check selected voltage.	11
7	Operation with the Bench Unit	12
7.1	Connection of Bench Unit	12
7.2	Connection of AC-line.....	12
7.3	Sealing	12
8	Operation with the Handle Unit.....	13
8.1	Connection of Handle Unit.....	13
8.2	Connection of AC-line.....	13
8.3	Sealing	13
9	Operation with the Ergonomic Handle Unit	14
9.1	Connection of Ergonomic Handle Unit	14
9.2	Connection of AC-line.....	14
9.3	Sealing	14
10	Operation with the QUAD plus Unit	15
10.1	Connection(s) of QUAD plus Unit	15
10.1.1	to primary Power Unit	15
10.1.2	to Bench Unit - option	15
10.1.3	to secondary Power Unit - option	15
10.2	Connection of AC-line.....	15
10.3	Sealing	16
11	Operation with the Automatic Sealing Handle (ASH).....	17
11.1	Connection of ASH.....	17
11.2	Connection of AC-line.....	17
11.3	Sealing	17
12	Adjustments	18
12.1	Adjustment of the sealing time.	18
13	Maintenance	18
13.1	Cleaning the electrodes.....	18
13.1.1	Bench Unit.....	18
13.1.2	Handle unit	19
13.1.3	QUAD plus unit.....	19
13.1.4	Automatic Sealing Handle	19
13.1.5	Ergonomic Handle	20
14	Technical description.....	21
14.1	Block diagram.....	21
14.2	Function	22
14.3	Power supply	22
14.4	Timer board	23
14.4.1	Layout.....	23
14.4.2	Circuit description	23
15	Trouble shooting.....	24
16	Technical data.....	25
16.1	Power Unit.....	25
16.2	Bench Unit.....	25
16.3	Handle Unit	25
16.4	Ergonomic Handle Unit.....	25

CR4 User's Manual

16.5	QUAD plus Unit	25
16.6	Automatic sealing Handle	25
17	Spare parts list Power Unit	26
18	Spare parts list Handle Unit	27
19	Spare parts list Bench Unit	28
20	Spare parts Automatic Sealing Handle.....	29
21	Spare parts Ergonomic Sealing Handle MSH-II	30
22	MDD Declaration	31

Rev history

2009-03-10 TI Ergonomic handle included
2009-08-21 TI MDD changed, spelling control whole doc.
2009-11-03 TI MSH II coax cable new art. No. Spelling corrected.
2010-04-08 TI Page 12, CR4 for tubes of up to 6,5 mm is on demand
2010-09-06 TI, MDD changed back to 93/42/EEC
2012-11-09 JA Added Abelko Innovation commitment concerning technical feedback
2014-05-05 DC Address in DoC updated
2015-04-13 DC Corrected cable length
2015-04-20 DC 1.74m cable for MSH changed to 1.75m


4 Properties

- **Model for continuous use.** Biosealer CR 4 is built for sealing PVC-tubes, especially for blood bags tubes or sets for plasmaferes. CR 4 has a powerful RF-unit (high frequency), which makes it suitable both for routine procedures at donation rooms, and repeated operations at preparation rooms without overheating.
- **Bench Unit and Handle Unit.** The Power Unit can either be connected to a Bench Unit or an automatic or manual Handle Unit. The Bench Unit automatically seals when a tube is in position. The Handle Unit is easy to move and therefore it is easy to reach the part of the tube, which shall be sealed. Both units has a 1,9 m cable. The electrodes are well protected by a plastic cover and are easy to clean.
- **QUAD plus Unit.** The Power Unit can also be connected to a QUAD plus Unit. The QUAD plus Unit automatically seals 1-4 seal when the tube is in position or on command. It is connected to the Power Unit by 1.5 m cables.
- **Wide and safe seal.** The RF gives a wide and safe seal of about 3mm with a marking in the middle for making it easier to pull the tubes apart.


5 Description

5.1 Power Unit

- (1) **AC-line indication.** When the power is turned on this indication is lit.
- (2) **Sealing indication.** When the electrodes are closed around the tube and sealing starts, this indication is lit and stays so until the sealing is completed.
- (3) **Indication of high temperature.** If the temperature inside the Power Unit exceeds the predetermined alarm level, this indication is lit.
- (4) **Adjustment of the sealing time.** When you seal other types of blood bags, the sealing time can be adjusted with this potentiometer. See part "Adjustments".
- (5) **Contact for the signal cable.** When using the Bench Unit or the QUAD plus Unit, connect the signal cable from the Bench Unit or the QUAD plus Unit to this contact.
- (6) **Contact for RF.** The coax cable from the Bench Unit, the Handle Unit or the QUAD plus Unit shall be connected to this contact.
- (7) **Coax cable.** This cable is used for connecting either the Bench Unit or the Handle Unit to the contact for RF (HF-OUT).


- (8) **AC-line inlet.** For the connection of the power cable.
- (9) **Power switch.** To turn off or on the power. When the power is on, the indication is lit.
- (10) **Power switch box.** There are two fuses in the box. Selected mains voltage is shown.
You can choose the mains voltage with the switch in the box. See part "Installation".


5.2 Bench Unit

- 11) **Electrodes.** Put the PVC-tube between the electrodes. Then the tube is automatically sealed.
- 12) **Indication of sealing.** Is lit when sealing starts and stays lit until the seal is completed.
- 13) **Signal cable.** Before turning the power on, the signal cable from the Bench Unit shall be connected to the Power Unit.


5.3 Handle Unit.

- (14) **Electrodes.** Put the PVC-tube between the electrodes and then press the Handle Unit to seal the tube.
- (15) **Sealing indication.** Is lit when sealing starts and is lit until the seal is completed.
- (16) **Protection cover.** Protects electrodes from damage and also prevent unintentional finger contact with electrodes.



5.4 Automatic Sealing Handle (ASH)

- (1) **Start button.** A press starts the seal process.
- (2) **Sealing indication.** A green lamp indicates when there is high frequency power between electrodes.
- (3) **Electrodes.** Transfers the high frequency energy to the PVC tube. The tube to be sealed must be placed between the electrodes.
- (4) **Protection cover.** Protects electrodes from damage and also prevent unintentional finger contact with electrodes.


5.5 QUAD plus Unit

- (17) **Electrodes.** Put the PVC-tube between the electrodes and the tube is automatically sealed if AUTO is selected.
- (18) **Signal cable.** Before turning the power on, the signal cable from the QUAD plus Unit shall be connected to the Power Unit.
- (19) **Coax cable.** Before turning the power on, the coax cable from the QUAD plus Unit shall be connected to the Power Unit.
- (20) **Knurled screw.** Remove when cleaning electrodes.
- (21) **Connection for primary Power Unit.**
- (22) **Connection for secondary Power Unit.**
- (23) **Connection for Bench Unit.**
- (24) **READY indication.** Lit when the unit is ready to perform a seal sequence. It is turned off during sealing or if a connected Bench Unit is sealing. When seal sequence is finished it is lit again. If the tube is not in position when the SEAL button is pressed (e.g. manual start of seal) it is flashing for two seconds.
- (25) **SEAL button.** Start seal sequence when pressed if the tube is in position.
- (26) **Sealing selection and indication.** Press to select number of seals that shall be performed and the corresponding indication is lit. During sealing it indicates which clamp is active.
- (27) **AUTO selection and indication.** When lit, seal starts when the tube is in position.

5.6 Ergonomic Sealing Handle (MSH-II)


- (1) **Sealing indication.** Is lit when sealing starts and is lit until the seal is completed.
- (2) **Electrodes.** Transfers the high frequency energy to the PVC tube. The tube to be sealed must be placed between the electrodes.
- (3) **Protection cover.** Protects electrodes from damage and also prevent unintentional finger contact with electrodes.
- (4) **Coaxial cable**

6 Installation

6.1 Check selected voltage.

- (1) The predetermined voltage setting is shown in the switch box. Check that shown voltage agrees with your local AC-line.
- (2) If you have to change the voltage switch, then first remove the power cable from the mains inlet and then the fuse box, by pressing the catch for the box.
- (3) Possible voltage selections is shown on the switch. Put the switch in the box with the selected voltage value on top. This value will show in the small window on the fuse box.


7 Operation with the Bench Unit

- Always turn the power off before you connect or disconnect a cable or both cables.
- Always use the coax cable that is delivered with the equipment. It is not allowed to extend the cable. If you need a longer cable, then order a longer cable, see spare parts list.


7.1 Connection of Bench Unit

- (1) Connect the Bench Unit to the RF-contact on the Power Unit with the coax cable.
- (2) Connect the signal cable from the Bench Unit to the contact on the Power Unit.

7.2 Connection of AC-line.

- (1) Connect the power cable to the AC-inlet on CR 4
- (2) Connect the power cable to the AC-line.
- (3) Turn the power on.


7.3 Sealing

- 1) Place the PVC-tube in the slot of the Bench Unit.
- 2) The sealing automatically starts, and the indication diodes are lit on the Power Unit and the Bench Unit.

Note !

Do not pull the PVC-tube when sealing.

- (14) The PVC-tube is sealed in appr. 1 sec. The electrodes automatically releases the tube and the indications is switched off. CR4 is adjusted for tubes of 4-5 mm thickness. CR4 for tubes of up to 6,5 mm is on demand. See "Adjustments" for sealing of other tubes.

Note !

When sealing the PVC-tube on two or more places the distance between the seals must be more than 50 mm. Shorter distance may cause leakage.


8 Operation with the Handle Unit

Always turn the power off before you connect or disconnect the coax cable.


8.1 Connection of Handle Unit

- (1) Connect the Handle Unit to the RF-contact on the Power Unit with the coax cable.


8.2 Connection of AC-line

- (1) Connect the power cable to the power inlet on CR 4.
- (2) Connect the power cable to the AC-line.
- (3) Turn the power on by the switch on the backside of the Power Unit

8.3 Sealing

- (1) Place the PVC-tube in the slot of the Handle Unit.
- (2) Press the Handle Unit. The sealing indication on both the Power Unit and Handle Unit is lit.
- (3) The PVC-tube is sealed in appr. 1 sec. The PVC-tube is sealed in appr. 1 sec. Wait a short moment (0.5 sec) after the indication is turned off to let the seal cool down before releasing the tube. CR 4 is adjusted for PVC-tubes of 4-5 mm thickness. See "Adjustments" for sealing other tubes.

Note !

When sealing the PVC-tube on two or more places the distance between the seals must be more than 50 mm. Shorter distance may cause leakage.

9 Operation with the Ergonomic Handle Unit


Always turn the power off before you connect or disconnect the coax cable.

9.1 Connection of Ergonomic Handle Unit

- (1) Connect the Handle Unit to the RF-contact on the Power Unit with the coax cable.


9.2 Connection of AC-line

- (1) Connect the power cable to the power inlet on CR 4.
- (2) Connect the power cable to the AC-line.
- (3) Turn the power on by the switch on the backside of the Power Unit

9.3 Sealing

- (1) Place the PVC-tube in the slot of the Handle Unit.
- (2) Press the Handle Unit. The sealing indication on both the Power Unit and Handle Unit is lit as long as the RF is on.
- (3) The PVC-tube is sealed in appr. 1 sec. Wait a short moment (0.5 sec) after the indication is turned off to let the seal cool down before releasing the tube.
- (4) CR 4 is adjusted for PVC-tubes of 4-5 mm thickness. See "Adjustments" for sealing other tubes.

Note !

When sealing the PVC-tube on two or more places the distance between the seals must be more than 50 mm. Shorter distance may cause leakage.

10 Operation with the QUAD plus Unit

- Always turn the power off before you connect or disconnect a cable or cables.
- Always use the coax cable that is delivered with the equipment. It is not allowed to extend the cable.

10.1 Connection(s) of QUAD plus Unit

10.1.1 to primary Power Unit


- (1) Connect the QUAD plus Unit primary RF and the primary Power Unit RF with the 1.5 m coax cable.
- (2) Connect the QUAD plus Unit primary Signal and the primary Power Unit SIGNAL with the 1.5 m signal cable.

10.1.2 to Bench Unit - option

- (3) The Bench Unit is connected to the QUAD plus Unit by 1.9 m coax and signal cables.


10.1.3 to secondary Power Unit - option

- (4) Connect the QUAD plus Unit secondary RF and the secondary Power Unit RF with the 1.5 m coax cable.
- (5) Connect the QUAD plus Unit secondary Signal and the secondary Power Unit SIGNAL with the 1.5 m signal cable.


10.2 Connection of AC-line.

- (1) Connect the power cable to the AC-inlet on CR 4
- (2) Connect the power cable to the AC-line.
- (3) Turn the power on. Always turn on the primary CR4 first and then the secondary CR4.
- (4) Always turn off the secondary CR4 first and then the primary CR4.

10.3 Sealing

- (1) Choose number of seals by pressing down button 1,2,3 or 4. The corresponding LED is lit.
- (2) To choose automatic start of seal, press down AUTO and the indication is lit.
- (3) Place the PVC-tube in the slot of the QUAD plus Unit.
- (4) The sealing automatically starts, SEALING diode on the Power Unit is lit once every seal and the READY-indication on the QUAD plus Unit is off until the sequence is completed and then lit again.
- (5) For manual start of seal (e.g. AUTO is not selected and lit) press SEAL (see 4).
- (6) For sealing with the Bench Unit see 7.3 . The READY-indication is off when the Bench Unit is sealing.
- (7) When the secondary Power Unit is connected to the QUAD plus Unit it will perform 2 seals simultaneously.

Note !

Do not pull the PVC-tube when sealing.

11 Operation with the Automatic Sealing Handle (ASH)


- Always turn the power off before you connect or disconnect a cable or cables.
- Always use the coax cable that is delivered with the equipment. It is not allowed to extend or shorten the cable.

11.1 Connection of ASH

- (1) Connect the signal cable from ASH to SIGNAL at Power Unit.
- (2) Connect the coax cable from ASH to RF OUT at Power Unit.

11.2 Connection of AC-line.

- (3) Connect the power cable to the AC-inlet on CR 4
- (4) Connect the power cable to the AC-line.
- (5) Turn the power on by the switch on the backside of the Power Unit

11.3 Sealing

- (1) Place the PVC-tube in the slot of the Handle Unit.
- (2) Press the start button. The electrodes are automatically pressing the tube and makes a seal. The sealing indication on both the Power Unit and ASH is lit.

Note !

Do not pull the PVC-tube when sealing.

- (3) The PVC-tube is sealed in appr. 1 sec. The electrodes are automatically released and the indicators are turned off. CR 4 is adjusted for PVC-tubes of 4-5 mm thickness. See "Adjustments" for sealing other tubes.

c

Note !


When sealing the PVC-tube on two or more places the distance between the seals must be more than 50 mm. Shorter distance may cause leakage.

12 Adjustments

12.1 Adjustment of the sealing time.

- (1) Remove the plastic cover for the potentiometer. It shows in the opening.
- (2) For thicker tubes than 4-5 mm, turn the potentiometer clockwise to increase the sealing time. The RF then will be turned on for a longer period.
- (3) If the sealing power is too strong and the sealed parts of PVC tube melts, turn the potentiometer counter clockwise and thereby reduce the sealing time. The RF is then turned on for a shorter period.
- (4) When the adjustment is completed replace the plastic cover.

Potentiometer for adjustment of sealing time

13 Maintenance



Always turn the power off before cleaning the electrodes !

13.1 Cleaning the electrodes.


13.1.1 Bench Unit

- (1) Remove the coax-cable from the Bench Unit.
- (2) Press gently on top of the cover and pull the lever for the locking device on the backside. This will release the cover locking.
- (3) Lift the cover.
- (4) Clean both electrodes with distilled water.
- (5) Dry up the electrodes carefully with a soft and lint-free cloth.
- (6) Put the cover back in position by pressing the cover while pushing the locking device.
- (7) Replace the coax cable.


13.1.2 Handle unit

- (1) Remove the coax-cable from the handle unit.
- (2) Loosen a bit the two screws holding the protection cover and tilt it down.
- (3) Clean both electrodes with distilled water.
- (4) Dry up the electrodes carefully with a soft and lint-free cloth.
- (5) Tilt up the protection cover and fasten it with the two screws.
- (6) Replace the coax cable.

13.1.3 QUAD plus unit

- (1) Remove all coax- and signal cables from the QUAD plus unit.
- (2) Unscrew the two knurled screws holding the cover and lift it.
- (3) Clean all electrodes with distilled water.
- (4) Dry up the electrodes carefully with a soft and lint-free cloth.
- (5) Put the cover back in position and fasten it with the two knurled screws.
- (6) Replace the cables.

13.1.4 Automatic Sealing Handle

- (1) Remove the signal and coax cable from the Power Unit.
- (2) Pull the protection cover forwards.
- (3) Clean both electrodes with distilled water.
- (4) Dry up electrodes with a soft and lint-free cloth.
- (5) Put the protection cover back again and reconnect the cables.


- Do not press the electrodes excessively in order to keep electrodes parallel.
- Check that the electrodes are completely dry.

13.1.5 Ergonomic Handle

(1) Remove the coaxial cable from the handle unit

(2) Remove the protective cover by first pressing the front down as shown in the figure.


(3) Then pull the cover forwards.

(4) Clean both electrodes with distilled water.

(5) Dry up the electrodes carefully with a soft and lint-free cloth.

(6) Replace the coax cable.

14 Technical description

CR 4 consists of a RF module operating at 40,68 MHz with a maximum output power of 90W. The RF-generator (HFG-01) is a complete module mounted on the backside of CR 4 to achieve the best possible cooling. HFG-01 is provided with current, 24 V, from a transformer with a rectifier, with a maximum possible current of 6 ADC.


Note!

Due to regulations of maximum emission (radiation) at other frequencies than 40,68 MHz the HFG-01 module is soldered up and shall be regarded as a component. It can not and it is not allowed to be adjusted or repaired by anyone else than the manufacturer. In the event of malfunction of HFG-01 it must be sent to the manufacturer for repair.

The voltage that occurs during sealing may cause burn damages if someone touches the fixed electrode during the sealing period. Therefore it is absolutely necessary to turn off the power before cleaning or service of the Handle Unit, Bench Unit or the QUAD plus Unit.

14.1 Block diagram

14.2 Function

The sealing function is controlled by the timer board mounted in the front on the Power Unit. The timer board consists of a microprocessor which governs all functions. When the Handle Unit is pressed, or the tube is placed between the electrodes on the Bench Unit, the connection from the central conductor in the coax cable to ground is released and a DC-voltage of 12 V is supplied from the timer board.

This voltage is recognized by the microprocessor (input PB0) on the timer card. The processor first checks if there is any leakage of current between the electrodes, if so this indicates moist on the tube. The limit is set to $5\mu\text{A}$, which corresponds to 2,5 V on input PB1 on the processor. If the value exceeds this limit, the seal is cancelled and an indication is given by the yellow diode "sealing" which is flashing. This means that the tube is damp and that there is a risk for sparking and leakage.


If no damp is found, the transistor T5 is turned on to activate the solenoid in the Bench Unit, if this unit is connected to the Power Unit. This signal is on during appr. 150 msec to achieve a compression of the tube, before the RF is turned on. Transistor T3 turns on the RF generator and keeps it on during the period determined by the timer potentiometer P1. Thereafter the RF is turned off, and appr. 250 msec later T5 releases the solenoid in the Bench or the QUAD plus unit.

The temperature on the RF-module is measured with a termistor mounted on the rear part above the RF-module and if the temperature exceeds predetermined value on potentiometer P2 the start of sealing is cancelled and the red diode "Temp alarm" is lit. The limit is set to 60°C at factory, but can be changed with potentiometer P2 according to the values in the circuit description.

If the electrodes are shortened during sealing the RF-generator immediately is stopped.


Note! Never adjust P2 to a lower value than 2,2V, which means 70°C , since this may cause over temperature in the RF-generator and damage the transistors.

14.3 Power supply


14.4 Timer board

14.4.1 Layout

14.4.2 Circuit description

15 Trouble shooting

Do not adjust the electrodes by yourself. They can be damaged by erroneous performed adjustment.

Sealing does not start and the diode "Sealing" is not lit when the Handle Unit is pressed.	• Test another Handle Unit. • Change the coax cable. • Change the timer board.
The diode "Sealing" is lit but not the lamp on the Handle Unit	• Test another Handle Unit. • Change the coax cable. • Control the cables to the RF-generator. • Measure the voltage on the 22000uF capacitor. It shall be 26-30VDC. • Change the RF-generator.
The diode "Temp Alarm" is lit and sealing does not start.	• The RF-generator is overheated. Let it cool down until the alarm is turned off. Control set alarm level on the timer board, see table in the circuit description above.
The diode "Sealing" flashes and sealing does not start.	• The tube is damp, dry it.
Sparks appear by the electrodes during sealing.	• The electrodes are not parallel. Contact service personnel.
The seal is done but with bad quality, leakage may appear.	• Check that the coax cable is the right type, the length must be 1.9m for sealing handle type 1 and Bench Unit. Ergonomic sealing handle has a special type coax cable.

16 Technical data

16.1 Power Unit

Voltage demands :	100/115/230/240 VAC 50/60 Hz
Power consumption:	220 VA
Net fuses:	2x T2.5 AL 250V
Internal fuses :	Protection for high current, type PTC
Sealing time:	Stepless adjustable 0.5 till 5 sec
Working frequency:	40.68 MHz
Output RF power	90W / 50Ω
Dimensions:	290(L) x 220 (W) x 88 (H) mm, (incl. RF-outlet)
Weight:	4.1 kg
Relative humidity incl. storage and transport:	10-95% not condensing
Working temperature:	10-40 °C
Storing-and transp. temp:	-40...+ 70 °C
Operation:	Continuously
Noice level	~60dB

Classification: Protection against electrical chock: class 1, type B

16.2 Bench Unit

Dimensions:	150 (L) x 60 (W) x 73(H) mm
Weight:	0.9 Kg
Length of cable:	1.9 m (4.2 or 9.1 m is available as option)

16.3 Handle Unit

Dimensions:	206 (L) x 27 (W) x 35 (H) mm
Weight:	0.3 kg
Length of cable:	1.9 m (4.2 or 9.1 m cable is available as option)
Manual pressing force	~11N

16.4 Ergonomic Handle Unit

Dimensions:	180 (L) x 35 (W) x 145 (H) mm
Weight:	0.220 kg
Length of cable:	Coiled cable 1.2m (Total length 1.75m)
Manual pressing force	~20N


16.5 QUAD plus Unit

Dimensions:	220 (L) x 319 (W) x 63 (H) mm
Weight:	3.9 kg
Length of cable:	1.5 m

16.6 Automatic sealing Handle


Mode of application:	Recommended max 1 seal/minute alternative 10 seals during 1 minute and after that standby in 5 minutes.
Dimensions:	188 (L) x 40 (W) x 40 (H) mm
Weight:	0.36 kg (0.5 kg incl. the cable)
Length of cable:	1,75 m

17 Spare parts list Power Unit





9-30401-00	Bottom plate (chassis)	9-30413-00	Ferrite core $\varnothing$ 29mm
9-30402-00	Rubber feet x 4	9-30414-00	HFG-01, Rf-generator
9-30403-00	Coverlid for timer adjust	9-30415-00	Ferrite core $\varnothing$ 13mm
9-30404-00	3-pol contact for Bench Unit	9-30416-00	PC-board timer/control
9-30405-00	Internal coaxial cable, BNC output	9-30417-00	Cover
9-30406-00	Diode bridge 20A	9-30418-00	Mains supply cable, Europe
9-30407-00	Mounting plate capacitor	9-30419-00	Connector socket incl. PTC fuse
9-30408-00	Capacitor 22000uF/40V	9-30420-00	Bleeding resistor
9-30409-00	Transformer RTR19340		
9-30410-00	Mains input module incl. filter		
9-30411-00	Cable HFG power supply		
9-30412-00	Cable HFG to timer board		

18 Spare parts list Handle Unit


9-32401-00	Cover
9-32402-00	BNC chassis female contact
9-32403-00	Switch mom close
9-32404-00	Internal coax cable
9-32405-00	Coil incl. delrin core and dim glowing lamp
9-32406-00	Dim glowing lamp
9-32408-00	Plate spring
9-32409-00	Handle
9-32410-00	Chassis
9-32411-00	Pin
9-32413-00	Rubber tile
9-32414-00	Electrode house complete
9-32415-00	Electrode cover
9-32416-19	Coax cable 1.9m
9-32416-42	Coax cable 4.2m
9-32416-91	Coax cable 9.1m

19 Spare parts list Bench Unit


9-33401-00	Cover	9-33411-00	Micro switch
9-33402-00	Signal cable	9-33412-00	Welding house chassis complete
9-33403-00	Locking device	9-33413-00	Pull back spring
9-33404-00	Chassis BNC contact	9-33414-00	Steering pin complete
9-33405-00	Hot moving electrode	9-33415-00	Pin compl. with cases and locking nut
9-33406-00	Fixed electrode	9-33416-00	Clamp spring
9-33407-00	Circuit board for connections	9-33417-00	Piston complete incl. coil and dim glowing lamp
9-33408-00	Solenoid	9-33418-00	Dim glowing lamp
9-33409-00	Chassis	9-33419-00	Welding house cover, micro switch side
9-33410-19	Coax cable 1.9m	9-33420-00	Welding house cover
9-33410-42	Coax cable 4.2m	9-33421-00	Coax cable internal
9-33410-91	Coax cable 9.1m	9-33422-00	RF-cable flexible

20 Spare parts Automatic Sealing Handle

9-34401-00	Cover with start button
9-34402-00	Chassis complete
9-34403-00	Coaxial/Signal Cable
9-34404-00	End plate
9-34405-00	Movable electrode
9-32406-00	Fixed electrode
9-32407-00	Electrode cover

21 Spare parts Ergonomic Sealing Handle MSH-II

9-35401-00	Right cover	9-35411-00	Pushing spring
9-35402-00	Left cover	9-35412-00	Piston complete incl. coil
9-35403-00	Handle	9-35413-00	Moving electrode
9-35404-00	Protection cover	9-35414-00	Fixed electrode
9-35405-00	Internal coax cable with SMA contact	9-34515-00	PC board complete incl. micro switch (Not in picture)
9-35406-00	Indicator lens	9-34516-00	3 Socket head cap screw for screen plates
9-35407-00	Screen plate bottom	9-35417-00	Screw Torx T25x8
9-35408-00	Screen plate top	9-35418-00	Screw Torx T30x8
9-35409-00	Pushing knob	9-35419-16	Coax cable coiled 1.2m, max 1.6m drawn out, total 1.75m (Not in picture)
9-35410-00	Pull back spring		

22 MDD Declaration

DECLARATION OF CONFORMITY

according to the

Medical Devices Directive 93/42/EEC

STANDARDS TO WHICH CONFORMITY IS DECLARED:

EN 60 601-1 (SEMKO certificate 971041/01)

EN 60 601-1-2 (SEMKO TR 9709076D.2)

Including: EN 55011 Kl B, EN 61000-4-2, EN 61000-4-3, EN 61000-4-4, EN 61000-4-5


Manufacturer: Abelko Innovation
Address: 972 54 Luleå
Telephone: +46 920-220360
Telefax: +46 920-220068

Marketing: Ljungberg&Kögel AB
Address: Box 1032, 251 10 Helsingborg
Telephone: +46 042-139860
Telefax: +46 042-132181

Type of Equipment: Tube-sealer Biosealer
Model: CR4
Product class: Class 1

I, the undersigned, hereby declare that the equipment specified above conforms to the above Directive and Standards.


Date of issue: 26 June 2014

/Thommy Lundström

Position/title: President

The document contains 34 bound numbered and sealed pages.

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A handwritten signature in blue ink, consisting of several loops and a long horizontal stroke.

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**User's Manual
and
Technical reference**

Ljungberg & Kögel AB

Rev 2015-04-14

1. Important

This user's manual is written for the person responsible for the operation of Biosealer CR4 AA. The operational methods and routines are developed and tested to ensure a reliable, safe and effective operation of CR4 AA. It is important that the operator has studied and understood the contents of this manual before using CR4 AA.

- This instrument is a sealing equipment using radio frequency for welding and the electrical emission at the operating frequency 40.68 MHz is high. Make sure that other instruments and equipment near the sealing unit can withstand this emission.
- Never touch the electrodes with your fingers during the sealing process! This may cause burn damages!
- Ensure that the PVC tube is dry on the surface.

As indicated by the -symbol, CR4 contains electronic components which are classified as dangerous for the environment and the equipment must be left to a recycling station when disposed. Local regulation decides how to recycle the biosealers, but it is preferred to disassemble the biosealer in its main parts and recycle the metal parts, plastic part, electronics etc.

2. Warranty

- 2.1. The manufacturer hereby guarantees the original buyer that Biosealer CR4 AA is manufactured in a professional and quality manner, and will be free from all faults during a period of one year from the date of delivery from the manufacturer.
- 2.2. The warranty includes equipment or components that prove to have faults during the warranty period. The manufacturer will without cost for the customer, repair or replace the equipment that is faulty.
- 2.3. The warranty is not valid if the equipment has been repaired by anyone else than qualified personnel, approved by the manufacturer.
- 2.4. The warranty is not valid if the equipment has been changed in any way that according to the manufacturers opinion, affects the reliability or stability of the instrument.
- 2.5. The warranty is not valid when the serial number has been changed, crossed over or been removed, or if the fault has been caused by misuse or abnormal use.
- 2.6. In these cases the manufacturer or a representative will inform the customer about the decision, and if wished by the client will repair the equipment for normal rate. An estimated price can be given on request.

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3. Table of contents

1. Important	2
2. Warranty	2
3. Table of contents	3
4. Properties.....	4
5. Description	4
5.1. Power Unit.....	4
5.2. Bench Unit.....	5
5.3. Handle Unit	5
5.4. Ergonomic Sealing Handle (MSH-III).....	6
6. Operation with the Bench Unit	7
6.1. Connection of Bench Unit	7
6.2. Sealing	7
7. Operation with the Handle Unit.....	8
7.1. Connection of Handle Unit	8
7.2. Sealing	8
8. Operation with the Ergonomic Handle Unit	9
8.1. Connection of Ergonomic Handle Unit MSH-III.....	9
8.2. Sealing	9
9. Adjustments	10
9.1. Change to manual timer	10
10. Maintenance	10
10.1. Cleaning the electrodes.....	10
10.1.1. Bench Unit.....	10
10.1.2. Handle unit	11
10.1.3. Ergonomic Handle MSH-III	11
11. Technical description.....	12
12. Trouble shooting.....	13
13. Technical data.....	14
13.1. Power Unit.....	14
13.2. Bench Unit.....	14
13.3. Handle Unit	14
13.4. Ergonomic Handle Unit.....	14
14. Spare parts list Power Unit	15
15. Spare parts list Handle Unit.....	16
16. Spare parts list Bench Unit	17
17. Spare parts list Ergonomic Sealing Handle MSH-III.....	18
MDD Declaration	19

Rev history

2011-11-08 TI Ergonomic handle MSH3 included
 2011-11-28 TI Coiled cable Ergo Handle changed part number
 2011-12-08 TI Spare parts added 2-pole contact
 2012-11-09 JA Added Abelko Innovation commitment concerning technical feedback
 2012-12-13 JA Minor changes in the chapter layout
 2012-12-13 DC Partlist updated with 5.1
 2013-01-05 DC Type in DoC updated to Biosealer
 2014-04-17 DC Address in DoC updated, Added explanation that CR4 is the same as CR4 AA
 2014-05-07 DC RoHS added in DoC
 2015-04-14 DC updated cable length
 2015-04-20 DC 1.74m cable for MSH changed to 1.75m

4. Properties

- 4.1. **Model for continuous use.** Biosealer CR4 AA is built for sealing PVC-tubes, especially blood bags tubes or sets for plasmaferes. CR4 AA has a powerful HF-unit (high frequency), which makes it suitable both for routine procedures at donation rooms and repeated operations at preparation rooms without overheating.
- 4.2. **Bench Unit and Handle Unit.** The Power Unit can either be connected to a handle Unit and a Bench Unit or two Handle Units. Bench Unit automatically seals when a tube is in position. The Handle Unit is easy to move and therefore it is easy to reach the part of the tube which shall be sealed. Both units has a 1,9 m cable. The electrodes of the Handle Units are well protected by a plastic cover and are easy to clean.
- 4.3. **Wide and safe seal.** The high frequency HF gives a wide and safe seal of about 3mm with a marking in the middle for making it easier to pull the tubes apart.

5. Description

5.1. Power Unit

5.1.1. Connectors and indicators

HANDLE: BNC contact for HF output to Ergonomic Handle or Manual Handle


HANDLE/BENCH UNIT: BNC contact for HF output to Ergonomic or Manual Handle or a Bench Unit.

SIGNAL CABLE: Signal output for controlling the Bench Unit electromagnet


STATUS: BI-colour LED indicating power on (fixed green), sealing (blinking orange) or some error (blinking red).

TEMP ALARM: Red LED if HF generator is overheated.

5.2. Bench Unit


5.3. Handle Unit.

5.4. Ergonomic Sealing Handle (MSH-III)

- (1) **Sealing indication.** Is blinking blue during sealing process and turns to green when ready.
- (2) **Electrodes.** Transfers the high frequency energy to the PVC tube. The tube to be sealed must be placed between the electrodes.
- (3) **Protection cover.** Protects electrodes from damage and also prevent unintentional finger contact with electrodes.
- (4) **Coaxial cable**

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
6.1. Connection of Bench Unit

- 6.1.1. Connect the Bench Unit to the HF-contact on the Power Unit with the coax cable.
- 6.1.2. Connect the signal cable from the Bench Unit to the signal contact on the Power Unit.


6.2. Sealing

- 6.2.1. Place the PVC-tube in the slot of the Bench Unit.
- 6.2.2. The sealing automatically starts, and the indication diodes are lit on the Power Unit and the Bench Unit.

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- 6.2.3. The PVC-tube is sealed in appr. 1-2 sec. The electrodes automatically release the tube and the indication is switched off. CR4 AA is adjusted for tubes of 4-5 mm thickness. CR4 AA for tubes of up to 6,5 mm is on demand. See chapter 9 for sealing of other tubes.

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

7.1. Connection of Handle Unit

- 7.1.1. Connect the Handle Unit to the HF-contact on the Power Unit with the coax cable.

7.2. Sealing

- 7.2.1. Place the PVC-tube in the slot of the Handle Unit.
- 7.2.2. Press the Handle Unit. The sealing indication on both the Power Unit and Handle Unit is lit.
- 7.2.3. The PVC-tube is sealed in appr. 1-2 sec. The PVC-tube is sealed in appr. 1 sec. Wait a short moment (0.5 sec) after the indication is turned off to let the seal cool down before releasing the tube. CR4 AA is adjusted for PVC-tubes of 4-5 mm thickness. See chapter 9 for sealing other tubes.

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8. Operation with the Ergonomic Handle Unit

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8.1. Connection of Ergonomic Handle Unit MSH-III


8.1.1. Connect the Handle Unit to the HF-contact on the Power Unit with the coax cable.

8.2. Sealing

8.2.1. Place the PVC-tube in the slot of the Handle Unit.

8.2.2. Press the Handle Unit. The sealing indication on the Handle Unit is blinking blue during seal process and turns to green when ready after about 1-2 sec.

8.2.3. CR4 AA is adjusted for PVC-tubes of 4-5 mm thickness. See chapter 9 for sealing other tubes.

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
9. Adjustments

9.1. Change to manual timer

Note! This is an operation which must be done by an authorized technician. We do not recommend manual timer settings since this can lead to overheating the HF generator or the electro magnet in the Bench Unit. It should only be used when sealing very special tubes. If in doubt contact Ljungberg&Kögel before thi settings.

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
10.1. Cleaning the electrodes.


10.1.1. Bench Unit

- 10.1.1.1. Remove the coax-cable from the Bench Unit.
- 10.1.1.2. Press gently on top of the cover and pull the lever for the locking device on the backside. This will release the cover locking.
- 10.1.1.3. Lift the cover. Clean both electrodes with distilled water or alcohol.
- 10.1.1.4. Dry up the electrodes carefully with a soft and lint-free cloth.
- 10.1.1.5. Put the cover back in position by pressing the cover while pushing the locking device.
- 10.1.1.6. Replace the coax cable.


10.1.2. **Handle unit**

- 10.1.2.1. Remove the coax-cable from the handle unit.
- 10.1.2.2. Loosen a bit the two screws holding the protection cover and tilt it down.
- 10.1.2.3. Clean both electrodes with distilled water.
- 10.1.2.4. Dry up the electrodes carefully with a soft and lint-free cloth.
- 10.1.2.5. Tilt up the protection cover and fasten it with the two screws.
- 10.1.2.6. Replace the coax cable.


10.1.3. **Ergonomic Handle MSH-III**

- 10.1.3.1. Remove the coaxial cable from the handle unit
- 10.1.3.2. Remove the protective cover by first pressing the front down as shown in the figure.

- 10.1.3.3. Then pull the cover forwards.

- 10.1.3.4. Clean both electrodes with distilled water.
- 10.1.3.5. Dry up the electrodes carefully with a soft and lint-free cloth.
- 10.1.3.6. Replace the coax cable.

11. Technical description

CR4 AA consists of a HF generator module operating at 40,68 MHz with a maximum output power of 100W. The HF generator (HFG-01) is a complete module mounted on the backside of CR4 AA to achieve the best possible cooling. HFG-01 is provided with current, 24 V, from a 200W medical approved power supply.

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The voltage that occurs during sealing may cause burn damages if someone touches the hot electrode during the sealing process. Therefore it is absolutely necessary to turn off the power before cleaning or service of the Handle Unit or Bench Unit.

11.1. Sealing Function

The sealing function is controlled by the timer board mounted on the HFG. The timer board consists of a microprocessor which governs all functions. When the Handle Unit is pressed, or when a tube is placed in the Bench Unit the connection from the central conductor in the coax cable to ground is released and a DC-voltage of 5 V is supplied from the timer board.

This voltage is recognized by the microprocessor on the timer board. The processor first checks if there is any leakage of current between the electrodes, if so this indicates damp on the tube. The limit is set to 6 μ A. If the value exceeds this limit, the seal is cancelled and an indication is given by the diode STATUS which is flashing red indicating "Moist on tube". There will also be an audible warning signal. This means that the tube is damp and that there is a risk for sparking and leakage.

If no damp is found the HF generator is turned on. The microprocessor has an automatic sensing algorithm which turns off the HFG when the seal is finished. After this there is a cooling period of 500 ms to let the seal cool down and then the seal is completed and the Status diode and the lamp on the handle are turned green and the indication on the Bench Unit is turned off and the electromagnet releases the tube.

The temperature of the HFG is measured with a termistor mounted on the rear part and if the temperature exceeds an alarm set level of 75 °C the start of sealing is cancelled and the red diode "Temp Alarm" is blinking.

If the electrodes are shortened during sealing the HF-generator immediately is stopped.

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Sealing does not start and the diode "Status" is not blinking when the handle is pressed.	• Test another handle. • Change the coax cable
The diode "Status" is blinking orange but not the lamp on the handle or the Bench Unit.	• Test another handle. • Test another Bench Unit • Change the coax cable • Change the HF-generator. (Only for service personnel).
The diode "Temp Alarm" is lit and sealing does not start.	• The HF-generator is overheated. Let it cool down until the alarm is turned off.
The diode "Status" flashes red and sealing does not start.	• The tube is damp, dry it.
Sparks appear by the electrodes during sealing.	• This can happen if the electrodes are overheated due to too many seals in a short time. Let the electrodes cool down.
Sparks appear by the electrodes during sealing.	• If there has been a leakage and flashes have created carbon particles on the electrodes, these particles must be removed completely before sealing can continue.
Sparks appear by the electrodes during sealing.	• It can indicate that the electrodes are not parallel. Contact service personnel.
The seal is done but with bad quality, leakage may appear.	• Check that the coax cable is the right type, the length must be 1.9 m for the Bench Unit and the Manual Handle. The Ergonomic Handle must have the coiled cable. <i>Note! The cables must not be shortened or lengthened.</i>
Short-circuit is shown on display or Status is blinking red during a seal, blue lamp is turned off.	• The electrodes are getting hot and the tubes melts easier which may cause the electrodes to shorten during sealing. Let the electrodes cool down by pausing for 10-15 minutes.

13. Technical data

13.1. Power Unit

Voltage demands:	90-250 VAC 50/60 Hz
Power consumption:	200 VA
Net fuses:	2x T 3,15A L 250V
Internal fuses:	Protection for high current, type PTC
Frequency:	40,68 MHz crystal controlled
Output HF power:	100W/50Ω maximum power
Seal time:	Automatic steppless or manually steppless adjustable 2-10 sec
Dimensions:	290x205x85 mm (LxWxH)
Weight:	2.1 Kg
Rel. Humidity incl. storage and transport:	10-95% not condensing
Working temperature:	10-40 °C
Storage- and transportation temperature:	-40- +70 °C
Operation:	Recommended max 1 seal each 3rd sec during continuous use or max 150 seals in sequence and after that 15 minutes rest.
Noice level	~60dB during normal use
Protective classification:	Protection against electrical shock: Class I type B. 

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

14. Spare parts list Power Unit





9-50401-00	Bottom plate (chassis)	9-50408-00	24V power contact incl cable
9-50402-00	Rubber feet x 4 (Not in picture)	9-50409-00	Bench unit internal signal cable
9-50403-00	Power supply 200W	9-50410-00	HF generator incl. timer board
9-50404-00	20-pole FFC cable	9-50411-00	Mains inlet
9-50405-00	Front board CR5	9-50412-00	Mains internal cable
9-50406-00	Signal contact bench unit	9-50413-00	Ferrite core $\varnothing$ 29mm
9-50407-00	Internal coax cable	9-50414-00	Cover (Not in picture)
		9-50415-00	2-pole input power contact

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
17. Spare parts list Ergonomic Sealing Handle MSH-III

9-35401-00	Right cover	9-35411-00	Pushing spring
9-35402-00	Left cover	9-35412-00	Piston complete incl. coil
9-35403-00	Handle	9-35413-00	Moving electrode
9-35404-00	Protection cover	9-35414-00	Fixed electrode
9-35405-00	Internal coax cable with SMA contact	9-34515-00	PC board complete incl. micro switch (Not in picture)
9-35405-01	Washer for SMA contact		
9-35406-00	Indicator lens	9-34516-00	3 Socket head cap screw for screen plates
9-35407-00	Screen plate bottom	9-35417-00	Screw Torx T25x8
9-35408-00	Screen plate top	9-35418-00	Screw Torx T30x8
9-35409-00	Pushing knob	9-35419-16	Coax cable coiled 1.2m, max 1.6m drawn out, total 1.75m (Not in picture)
9-35410-00	Pull back spring		

MDD Declaration

DECLARATION OF CONFORMITY

according to the
Medical Devices Directive 93/42/EEC
and
RoHS directive, 2011/65/EU

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Type of Equipment: Biosealer
Model: CR4 AA
Product class: Class 1

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Date of issue: 26 June 2014

A handwritten signature in black ink, appearing to read 'Tommy Lundström'.

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Ljungberg & Kögel AB

Important

This user's manual is written for the person responsible for the operation of Biosealer CR6. The operational methods and routines are developed and tested to ensure a reliable, safe and effective operation of CR6. It is important that the operator has studied and understood the contents of this manual before using CR6.

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Electrical safety

Only use battery chargers approved for CR6 or charger of another manufacturer recommended by the manufacturer of the CR6 and ensure that the chargers are marked with accurate AC-line voltage.

Batteries

The batteries used in CR6 contains lead and worn out batteries must be left to a recycling station.

Disposal of CR6

As indicated by the  and **Pb** -symbols Biosealer CR6 contains electronic components which are classified as dangerous for the environment and the equipment must be left to a recycling station when disposed. Local regulation decides how to recycle the CR6, but it is preferred to disassemble the CR6 in its main parts and recycle the metal parts, plastic part, electronics etc.

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Rev history:

2007-01-11: Warranty is nor valid for batteries.
2009-05-13: Ergonomic sealing handle included
2009-11-03 TI: P16 coax cable new art-no
2010-08-13 TI: P 5: Sealing can be done with or without the charger connected to mains power.
2012-11-09 JA: Added Abelko Innovation commitment concerning technical feedback
2014-05-05 DC: Address in DoC updated
2015-04-13 DC: Updated information about chargers
2015-04-20 DC: 1.74m cable for MSH changed to 1.75m

Table of contents:

Important	2
Electrical safety	2
Batteries	2
Disposal of CR6	2
Features	4
Function	5
Operation	5
Battery test	6
Temperature-alarm	6
Timer setting	6
Maintenance and service	6
Cleaning the electrodes	6
Charging batteries	7
Change of batteries	8
Technical description	9
Block diagram	9
Sealing function	10
Layout Timer board	11
Circuit description Timer board	11
Circuit diagram CR6	12
Trouble shooting	13
Spare parts list Power Unit	14
Spare parts list Manual handle	15
Spare parts Ergonomic Sealing Handle MSH-II	16
Spare parts list Automatic handle	17
Technical data	18
Warranty	18
MDD-declaration	19

Features

Biosealer CR6 is a battery operated sealing unit for PVC-tubes. With new and fully charged batteries it manages more than 750 seals (500 with automatic handle) on 4-5 mm PVC-tubes. The high capacity makes it suitable to use CR6 both in donation rooms and laboratories.

Biosealer CR6 is comprised of 3 units:

- Manual sealing handle type 1 or ergonomic sealing handle type MSH-II or automatic handle ASH.
- Power unit with batteries, RF generator and a control unit monitored by a microprocessor.
- Battery charger

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Function

Sealing is done in about 1 sec and the width of the seal is about 3 mm with a marking in the middle which makes it easier to pull the tubes apart. Before sealing, the microprocessor controlled timer board checks that the tube not is damp, to avoid sparks and thereby leakage. After each seal a multicoloured diode is lit. The diode indicates the remaining capacity of the batteries. This can also be checked before sealing with the push button at the front of the power unit.

Operation

1. Ensure that the power is off and thereafter connect the sealing handle.
2. Turn the power on and place the PVC-tube between the electrodes.
Manual handle: Press the handle. The sealing automatically starts and two indications are given. The yellow diode "sealing" on the power unit is lit and the orange lamp, indicating high frequency power, on the handle is lit. The intensity of the orange lamp is reduced during the sealing and is turned off when the seal is completed. Wait a short moment (0.5 sec) after the indication is turned off to let the seal cool down before releasing the tube. The yellow diode is also turned off when the seal is completed.
Automatic handle: Push the start button. The sealing automatically starts and two indications is given. The yellow LED "sealing" on the power unit is lit and the green LED, indicating high frequency power, on the handle is lit. The green and the yellow LED is turned off when the seal is completed.
3. Immediately after the sealing the multicoloured diode on the power unit is lit. It shows green when the remaining capacity of the battery is 70-100%, orange for 25-70% and red when the capacity is lower than 25%.
4. A new seal can be made immediately after the yellow LED "sealing" has been turned off.

Note:

- Never have a shorter distance than 5 cm between two seals. Otherwise the pressure in the tube may cause leakage.
- Never stretch the tube during sealing, since this may cause leakage.
- **Manual handle:** Continuous sealing with an interval of 1 seal/sec during several minutes can cause overheating of the power unit or the electrodes with a declined quality of the seals or cause "temp alarm" on the power unit. Always use the delivered standard coax cable 1.9 m . If a longer cable is needed , 4.2 m or 9.1m shall be used, see list of spare parts.

Note! Ergonomic sealing handle has a special coax cable not compatible with standard coax cable.

- **Automatic handle:** Continuous sealing with an interval of less then one minute or more then 10 seals a minute can cause overheating and weakening of the solenoid which may result that the tube is not squeezed together properly. Red LED on the handle indicates overheating. If overheating occurs stop sealing immediately. Wait 5 minutes after the red LED has switch off before start sealing.
- Sealing can be done with or without the charger connected to mains power.

Battery test

Battery test can be done at any time by pressing the button "batt test". The "batt status" indicator will show the battery capacity in the same way as after each seal with green light for 70-100%, orange light for 25-70% and red light when capacity is less than 25%.

Temperature-alarm

If the temperature of the high frequency generator exceeds predetermined alarm level, in factory set to 60 °C, the diode "Temp alarm" is lit. Sealing can not be done before the RF- generator has cooled and the diode is turned off.

Timer setting

If the RF is turned off before the seal is completed the timer must be increased. A thicker tube normally demands a longer seal-time. The time can be set between 0.5-5 seconds with the potentiometer under behind the lid to the right of LED "sealing". Avoid setting the timer too long since this will decrease battery capacity and heat up the electrodes.

Remove the lid

1. Adjust time with a small screw driver. Clockwise turn increase time.
2. Replace the lid when ready.

Maintenance and service

The maintenance that can be done by the operator is limited to cleaning. All other service must be done by an approved service personnel.

The power unit and the handle can be cleaned with a mild soap solution.


Always turn off the power unit and disconnect the battery charger from the AC-line before cleaning the equipment.

Cleaning the electrodes


Manual handle type1 : Remove the coaxial cable. Unscrew the 2 screws to the plastic cover and remove it. Clean the electrodes with distilled water and dry with a soft cloth. Replace the coax cable.

Ergonomic Handle :

- (1) Remove the coaxial cable from the handle unit
- (2) Remove the protective cover by first pressing the front down as shown in the figure.

- (3) Then pull the cover forwards.
- (4) Clean both electrodes with distilled water.
- (5) Dry up the electrodes carefully with a soft and lint-free cloth.
- (6) Replace the coax cable.

Automatic handle:

Remove cover and clean the electrodes with distilled water and dry with a soft cloth.

Charging batteries

Only battery chargers intended for charging of CR6 or charger of another manufacturer recommended and approved by the manufacturer of the CR6 are allowed to be used.

1. Check the type marking sign that the power of the charger corresponds to the local AC -line voltage.
2. Connect the charger to the AC-line and the low voltage connector to the input at the back side of the CR6 marked "Charger Inlet". (It is allowed to make the connections in the reverse order).
3. A battery which has a capacity of 25% is fully loaded within 8 hours.
4. Always have the charger connected to CR6 when it is not used mobile.
5. Sealing can be done when the charger is connected.
6. It lacks significance whether CR6 is turned on or off during charging.
7. Do not let the capacity of the batteries be reduced to more than 25%, since this will shorten the life length of the batteries.
8. The life length of the batteries depends mainly on the surrounding temperature and as a guideline the life length at 25 °C temperature is 3-5 years and at 35 °C temperature it is 1-2 years. Repeated discharging of the batteries below 25% reduces the life length.

Change of batteries


Change of batteries must be done by approved service personnel.

The battery charger and the handle shall be disconnected when the batteries are changed. The batteries used in CR6 contains lead and worn out batteries must be left to a recycling station according to local regulations.

When the battery capacity begins to be insufficient for sealing during one day, the batteries have to be changed.

1. Unscrew the 4 screws for the cover to the power unit and remove the cover.
2. Loosen the connections to the batteries and carefully note how they are placed.
3. Unscrew the 2 screws for the battery clamp and remove the batteries.
4. Replace the batteries with new batteries of same model, see technical data, and attach the battery clamp.
5. Connect the cables in the same way as they were before the changing of the batteries.
6. Attach the cover.

Batteries location and connection:

Wrong connection of the batteries may cause that the cables get over-heated and that the RF-generator or the control card is damaged.

Technical description


Biosealer CR6 consists of a radio frequency (RF) generator operating at 40,68 MHz with a maximum effect of 80W. RF-generator HFG-01 is a complete module mounted in the power unit. HFG-01 is provided with current (24 V) by batteries with serial coupling. The batteries do not need maintenance and are charged with 28 VDC from an external battery charger.

Note!

Due to regulations of maximum radiation at other frequencies than 40,68 MHz the HFG-01 module is soldered up and shall be regarded as a component. It can not and it is not allowed to be adjusted or repaired by anyone else than the manufacturer. In the event of malfunction of HFG-01 it must be sent to the manufacturer for repair.

The voltage that occurs during sealing may cause burn damages if someone touches the fixed electrode during the sealing period. Therefore it is absolutely necessary to turn off the power before cleaning or service of any of the units.

Block diagram

Sealing function

The sealing function is controlled by the timer board mounted in the front of the power unit. The timer board consists of a microprocessor which governs all functions. When the sealing handle is pressed, the connection from the central conductor in the coax cable to ground is released and a DC-voltage of 12 V is supplied from the timer board.

This voltage is recognized by the microprocessor (input PB0) on the timer card. The processor first checks if there is any leakage of current between the electrodes, if so this indicates damp on the tube. The limit is set to 5 μ A, which corresponds to 2,5 V at input PB1 on the processor. If the value exceeds this limit, the seal is cancelled and an indication is given by the yellow diode "sealing" which is flashing. This means that the tube is damp and that there is a risk for sparking and leakage.

If no damp is found, the transistor T4 turns on to activate the solenoid of automatic handle if connected. This signal is on for about 150 ms to get the tube squeezed together properly before the RF is turned on. The RF generator is turned on by transistor T3 and keeps it on during the period determined by the timer potentiometer P1. Thereafter the RF is turned off by T3 and 250 ms later the solenoid is turned off by T4.


About 100 msec before the RF is turned off, the battery capacity is measured with the input PB2 on the microprocessor. After completed sealing the diode for the battery is lit and during 3 sec the capacity is shown. Green light for 70-100% left, green/ red light for 25-70% and red light if the capacity is below 25%.

The capacity of the batteries also can be tested before sealing, by pushing the "batt test" button at the front of the power unit. The microprocessor then discharge the batteries with approximately 0,5 A during 0,5 sec with transistor T7 and resistor R19, and measures the voltage at input PB5 and indication is given in the same way as described above.


The temperature of the RF-module is measured with a termistor mounted on the rear part above the RF-module. If the temperature exceeds predetermined value set to 60 °C at factory the start of sealing is cancelled and the red diode "Temp alarm" is lit.

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
Circuit diagram CR6


Trouble shooting

Sealing does not start and the diode "Sealing" is not lit when the handle is pressed.	• Test another handle. • Change the coax cable • Change the timer board
The diode "Sealing" is lit but not the lamp on the handle.	• Test another handle. • Change the coax cable • Control the cables to the RF-generator. • Change the RF-generator.
The diode "Temp Alarm" is lit and sealing does not start.	• The RF-generator is overheated. Let it cool down until the alarm is turned off. • Control set alarm level on the timer board, see table in the circuit description above.
The diode "Sealing" flashes and sealing does not start.	• The tube is damp, dry it.
Sparks appear by the electrodes during sealing.	• The electrodes are not parallel. Contact service personnel.
The seal is done but with bad quality, leakage may appear.	• Check that the coax cable is the right type, the length must be 1.9m for bench unit and sealing handle typ1. Ergonomic sealing handle has a special type coax cable.


Note: Ljungberg & Kögel AB will supply more detailed service-instructions with component lists and trouble-shooting guides on request.

Spare parts list Power Unit


9-30601-00	BATTERY x 2
9-30602-00	BATTERY CLAMP
9-30603-00	BOTTOM PLATE POWER UNIT
9-30604-00	BATTERY CHARGER INLET
9-30605-00	CABLE COAX. STANDARD FOR HF UNIT
9-30606-00	RUBBER FEET x 4
9-30607-00	RUBBER FEET FOR BATTERY POSITIONING x 2
9-30608-00	HOLDER FOR PC-BOARD
9-30609-00	POWER PANEL SWITCH
9-30610-00	TIMER/LOGIC PC-BOARD
9-30611-00	COWER LID FOR TIME ADJUST
9-30612-00	RF-UNIT
9-30613-00	COVER POWER UNIT
9-30614-00	INTERNAL POWER SUPPLY CABLE (NOT IN PICTURE)
9-30615-0	INTERNAL SIGNAL CABLE WITH SOLENOID CONTACT (NOT IN PICTURE)

Spare parts list Manual handle


9-32401-00	Cover
9-32402-00	BNC chassis female contact
9-32403-00	Switch mom close
9-32404-00	Internal coax cable
9-32405-00	Coil incl. delrin core and dim glowing lamp
9-32406-00	Dim glowing lamp
9-32408-00	Plate spring
9-32409-00	Handle
9-32410-00	Chassis
9-32411-00	Pin
9-32413-00	Rubber tile
9-32414-00	Electrode house complete
9-32415-00	Electrode cover
9-32416-19	Coax cable 1.9m
9-32416-42	Coax cable 4.2m
9-32416-91	Coax cable 9.1m
9-32417-00	Rubber tile switch

Spare parts Ergonomic Sealing Handle MSH-II

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Spare parts list Automatic handle

9-34401-00	Cover with start button
9-34402-00	Chassis complete
9-34403-00	Coaxial/Signal cable
9-34404-00	Rear cover (not in picture)
9-34405-00	Moving electrode
9-32406-00	Fixed electrode
9-32407-00	Electrode cover

Technical data

Batteries:	2x12V 2Ah Lead acid
Internal fuse:	Overload protection type PTC, self recovery.
Effect consumption:	170W maximum effect
Sealing time:	Stepless adjustable 0,5-5 sec
Frequency:	40,68 MHz crystal controlled
Output effect:	90W/50Ω maximum effect
Dimensions:	266x119x88 mm (LxWxH) incl. rf-outlet
Weight:	3.1 Kg incl. batteries
Rel. Humidity incl. storage and transport:	10-95% not condensing
Working temperature:	10-40 °C
Storage- and transport- tation temperature:	-40- +50 °C. (The life length of the batteries is shortened at high temperatures).
Operation:	Manual handle: Recommended max 1 seal each 3:rd sec during continuous use. Automatic handle: Recommended max 1 seal each minute during continuous use or max 10 seals during 1 minute and then rest for 5 minutes.
Noice level	~60dB during normal operation
Manual pressing force	Manual sealing handle: ~11N Ergonomic sealing handle: ~20N
Symbol: 	Protection against electrical chock: class II type B.

Warranty

The manufacturer hereby guarantees the original buyer that Biosealer CR 6 is manufactured in a professional and quality manner, and will be free from all faults during a period of one year from the date of delivery from the manufacturer. The warranty includes equipment and components that proves to have faults during the warranty period. The manufacturer will without cost for the customer, repair or replace the equipment that is faulty. The warranty is not valid for the batteries since they are considered to be expendable supplies. The warranty is not valid if the equipment has been repaired by anyone else than qualified personnel, that is approved by the manufacturer. The warranty is not valid if the equipment has been changed in any way that according to the manufacturers opinion, affects the reliability or stability of the instrument. The warranty is not valid when the serial number has been changed, crossed over or been removed, or if the fault has been caused by misuse or abnormal use. In these cases the manufacturer or the manufacturer's representative will inform the customer about the decision, and if wished by the client will repair the equipment for normal rate. An estimated price can be given on request.

DECLARATION OF CONFORMITY

according to the
Medical Devices Directive 93/42/EEC

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
Type of Equipment: Tube-sealer
Model: CR6
Product class: Class 1

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/Thommy Lundström

Position/title: President

The document contains 20 bound numbered and sealed pages.

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Biosealer CR6-AA

Including

**Bench Unit
Manual Handle
Ergonomic Handle MSH-III**

User Manual

Ljungberg & Kögel AB

This page is consciously left empty

This u
opera
opera
manua

7
em
Ne

Only u
the mar

The bat
but sho

Bioseale

directive
disposed
CR6 AA

Manufac

Marketin

Abelko In
customers
quickly an
distributor

Important

This user's manual is written for the person responsible for the operation of Biosealer CR6-AA. The operational methods and routines are developed and tested to ensure a reliable, safe and effective operation of CR6-AA. It is important that the operator has studied and understood the contents of this manual before using CR6-AA.

This instrument is a sealing equipment using radio frequency for welding and the electrical emission, at the operating frequency 40,68 MHz, is high. Make sure that other instruments and equipment near the sealing unit can withstand this emission.

Never touch the electrodes with your fingers during the sealing process! This may cause burn damages!

Ensure that the tubing is dry on the surface.

Electrical safety

Only use the battery charger intended for CR6-AA or charger of another manufacturer recommended by the manufacturer of CR6-AA and make sure it is marked with the accurate mains voltage.

Battery

The battery used in CR6-AA contains NiMH and is not considered to be dangerous for the environment but should be left to a recycling centre for recycling according to local regulation

Environmental effects

Biosealer CR6-AA contains none or insignificant parts of the forbidden hazardous substances in RoHS

directive 2002/95/EG but as indicated by the  -symbol, it should be left to a recycling station when disposed. Local regulation decides how to recycle the CR6 AA, but it is preferred to disassemble the CR6 AA in its main parts and recycle the metal parts, plastic part, electronics etc

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Table of contents

Important	3
Electrical safety	3
Battery.....	3
Environmental effects.....	3
Manufacturer	3
Table of contents	4
Properties	5
Function	6
Description	6
Ergonomic handle (MSH-III)	6
Manual handle	6
Bench Unit (BU).....	7
Operation	8
Connecting the MSH-III	8
Sealing with manual handle.....	8
Connecting the manual handle	9
Sealing	9
Connecting the Bench Unit.....	9
Sealing with Bench Unit.....	9
Adjustments	10
Changing to manual timer.	10
High temp alarm	10
Charging the battery	10
Maintenance and service	12
Cleaning the electrodes.....	12
Change of batteries	13
Trouble shooting.....	14
Spare parts list Power Unit.....	15
Spare parts list Ergonomic Sealing Handle MSH-III	16
Spare parts list Handle Unit.....	17
Spare parts list Bench Unit	18
Technical data.....	19
Power Unit.....	19
Bench Unit.....	19
Handle Unit	19
MSH III.....	19
Limited Warranty.....	20
Declaration of Conformity.....	21

Revision history:

2012-11-09 JA	Added Abelko Innovation commitment concerning technical feedback
2012-12-13 DC	Partlist for MSH3 updated with 5.1
2013-01-13 DC	DoC updated from Tube-sealer to BioSealer
2013-11-20 DC	DoC updated
2014-04-17 DC	Address in DoC updated
2014-08-13 DC	Manual sealing handle added
2015-04-14 DC	Cable length and information about charger updated
2015-04-20 DC	1.74m cable for MSH changed to 1.75m

Properties

Biosealer CR6-AA is a battery operated sealing unit for PVC-tubes, especially suitable for tubing from blood bags and tubing from sets for plasmapheresis.

It can also seal PVC/EVA and some types of pure EVA tubes.

CR6-AA has automatic timer sensing when the seal is ready for optimal seal time and battery endurance independent of type of tubes and thickness.

Appliances are either manual handle or Bench Unit (BU)

With new and fully charged battery it can make 1500 seals (1200 with BU) 4-5 mm PVC tubes.


This high capacity makes it suitable to use CR6-AA either mobile or stationary in donation rooms and laboratories.

CR6-AA also works with older L&K handles and BU:s.

Biosealer CR6-AA is comprised of 3 units:

- Manual handle, ergonomic manual handle MSH-III or BU including associated cables.
- Power Unit including battery, high frequency generator and timer unit.
- Battery charger.

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
Function

Sealing is done in 1-2 sec depending on type of tube and the width of the seal is about 3 mm with a marking in the middle which makes it easier to pull the tubes apart. Before sealing, the microprocessor controlled timer board checks that the tube is not damp, to avoid sparks and thereby leakage. After each seal a multicolored diode is lit indicating remaining capacity of the battery.


Description

Ergonomic handle (MSH-III)


- (1) **Seal indicator.** Lamp blinking blue during seal process and turns green when ready.
- (2) **Electrodes.** Transfer the high frequency energy to the tube to be sealed
- (3) **Protective cover.** Protects the electrodes from damage and also prevent unintentional finger contact with electrodes.
- (4) **Coaxial cable**

Manual handle

Bench Unit (BU)

Operation

Always turn the power off before you connect or disconnect the coax cable!

Connecting the MSH-III

Connect the MSH-III to the power unit with the coax cable.

***Use only hand power for connection.
Pulling too hard with a tool may damage the connector.***


Turn on the power switch, the lamps for STATUS and BATTERY will turn green. If the BATTERY lamp is orange or red the battery must be recharged before operation, see charge of battery on page 8.

Sealing with manual handle

Place the PVC-tube in the slot of the handle unit and press the handle.

The sealing will start directly and a blue lamp on the handle will start blinking and STATUS will blink orange on the power unit.


When the seal is ready the lamp on the handle and STATUS will turn green.

If the tube is easy welded or a great number of seals has been done in a short time making the electrodes hot it can happen that the plastic melts that much that the electrodes make contact and if this happens the green lamp will not be lighted. This is not an error, the seal will be satisfactory.


Connecting the manual handle

Connect the Handle Unit to the HF-contact on the Power Unit with the coax cable.

Sealing

1. Place the PVC-tube in the slot of the Handle Unit.
2. Press the Handle Unit. The sealing indication on both the Power Unit and Handle Unit is lit.
3. The PVC-tube is sealed in appr. 1-2 sec. Wait a short moment (0.5 sec) after the indication is turned off to let the seal cool down before releasing the tube. The handle is adjusted for PVC-tubes of 4-5 mm thickness. When sealing other tubes, please contact Ljungberg & Kögel.

Connecting the Bench Unit

Connect the BU to the power unit with the coax cable and signal cable.

Turn on the power switch, the lamps for STATUS and BATTERY will turn green. If the BATTERY lamp is orange or red the battery must be recharged before operation, see charge of battery on page 8.

Sealing with Bench Unit

Hold the tube with both hands and place it in the opening on the BU. The sealing will start directly and STATUS will blink orange on the power unit and an orange lamp on the BU will be lighted during the time high frequency is on. The intensity will go down during the melting process and turns off when ready.

Check of battery capacity is done regularly and BATTERY is green when it is 50-100% left, orange for 25-50% and red when there is less than 25% left.

A new seal may be done immediate after STATUS is green on power unit.

Don't make seals with shorter distance than 5 cm since the pressure in the tube may cause leakage.

Don't stretch the tube during sealing since this may cause leakage.

Don't make more than 50 continuous seals without a pause. Try to keep a seal rate of max 33% duty cycle.


Adjustments

Changing to manual timer.

This is an operation that must be done by an authorized technician. We recommend not using manual timer since this may cause overheated HF generator or solenoid in the BU. It should only be used for very special tubes. Please contact Ljungberg&Kögel if in doubt.

If you have hard sealed tubes that needs longer timer i.e. more time for the high frequency, you can disconnect the automatic state by turning the potentiometer TP2 on the timer PC board clockwise with a small screwdriver. Turning it more than 20% turns on the manual timer and the timer will be set from 2 to 10 sec. (potentiometer 20-100%). To be able to locate TP2 you have to lift off the cover by unscrewing the 4 screws on the side.

High temp alarm

If the temperature of the high frequency generator exceeds alarm level, set to 75°C, the diode "TEMP ALARM" is lit. Sealing cannot be done before the RF- generator has cooled down and the diode is turned off.

Charging the battery

The "Battery charger for CR6-AA" is specially designed for charging the NiMH batteries powering CR6-AA. It is medical approved according to EN 60601-1 and UL recognized for US and Canada.

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Only battery charger type "Battery charger for CR6-AA" or another charger recommended by L&K may be used.

1. Check that the local mains supply corresponds to the charger type sign.
2. Connect the power supply cord to the charger and then to the AC-line and the low voltage connector to the input at the back side of the CR6-AA marked "BATTERY CHARGER INLET". It is ok to connect the charger to CR6-AA first and then to AC-line.

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3. A battery with a capacity of 25% is fully loaded within 6-8 hours.
4. Since this is a NiMH battery pack it should never be overcharged. The recommendation for longest working life is to use the sealer unconnected to mains until the BATTERY turns orange and then recharge it.
5. Sealing can be done when the charger is connected.
6. Sealing can be done with the needle in the arm with or without connected charger.
7. It lacks significance whether CR6-AA is turned on or off during charging.
8. Avoid to completely unloading the battery since this will shorten the life of the battery.
9. The working life of the battery depends on the amount of charging cycles and the surrounding temperature included the temperature rise that occurs when batteries are overcharged. Typical expected working life is 3 years with a charging frequency of 200 cycles/year and 25°C. Deep discharging and overcharging shortens working life.
10. If the sealer is not used for a while the battery needs maintenance recharging every 3 months.
11. If the battery has been completely discharged during some time it may need repeatedly charging and discharging to obtain full capacity.

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For the best seal performance it is important that the electrodes always are clean. After a leakage it is very important to clean the electrodes carefully.

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

Cleaning the electrodes


MSH-III

1. Remove the coaxial cable from the handle unit
2. Remove the protective cover by first pressing the front down as shown in the picture.
3. Then pull the cover forwards.
4. Clean both electrodes with distilled water or alcohol.
5. Dry up the electrodes carefully with a soft and lint-free cloth.
6. Replace the protective cover and the coax cable.


Handle unit

1. Remove the coax-cable from the handle unit.
2. Loosen a bit the two screws holding the protection cover and tilt it down.
3. Clean both electrodes with distilled water.
4. Dry up the electrodes carefully with a soft and lint-free cloth.
5. Tilt up the protection cover and fasten it with the two screws.
6. Replace the coax cable.

Bench Unit (BU)

1. Unplug the coaxial cable from the BU.
2. Turn BU upside down, see picture.
3. Clean both electrodes with distilled water or alcohol.
4. Dry up the electrodes carefully with a soft and lint-free cloth.
5. Replace the protective cover and the coax cable.

A more complete cleaning may be done by opening the cover. Unplug the coax cable and pull the locking clamp backwards.

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Change of batteries must be done by approved service personnel. The battery charger and the handle must be disconnected when the batteries are changed.

CR6-AA must be turned OFF. Always use battery packs from L&K

1. Unscrew the 4 screws holding the cover and remove it.
2. Loosen the battery cables (red and black).
3. Unscrew the 2 screws to the battery holder, now the battery pack can be removed.
4. Change battery pack and reassemble in reversed order.

The battery used in CR6-AA contains NiMh and is not considered to be dangerous for the environment but should be left to a recycling station for recycling.

Trouble shooting

The blue lamp on the ergonomic handle stops blinking but the green lamp is not lightened.	This may occur if the welding makes the plastic melts until the electrodes makes contact. This is no error, the seal will be satisfactory.
Sealing does not start and the lamp "Status" is not blinking when the handle is pressed	Test another handle. Test another coax cable.
The lamp "Status" is blinking orange but not the lamp on the handle or the Bench Unit.	Test another handle. Test another Bench Unit. Change the coax cable Change the HF-generator. (Only for service personnel).
The lamp "Temp Alarm" is lit and sealing does not start.	The HF-generator is overheated. Let it cool down until the alarm is turned off.
The lamp "Status" flashes red and sealing does not start.	The tube is damp, dry it.
Sparks appear by the electrodes during sealing.	This can happen if the electrodes are overheated due to too many seals in a short time. Let the electrodes cool down.
	If there has been a leakage and flashes have created carbon particles on the electrodes, these particles must be removed completely before sealing can continue.
	It can indicate that the electrodes are not parallel. Contact service personnel.


The electrodes must be adjusted by an authorized service technician. They can be damaged by erroneous performed adjustment.


More detailed troubleshooting, training of technician and service instructions can be obtained by Ljungberg & Kögel on demand.

Spare parts list Power Unit


9-31601-00	Battery pack
9-31602-00	Battery holder
9-31603-00	HFG incl. timer board
9-31604-00	Front board
9-31605-00	Internal coaxial cable
9-31606-00	Ferrite ring
9-31607-00	C91 contact, (signal to BU)
9-31608-00	Inlet for charger including cable
9-31609-00	On/Off switch
9-31610-00	20-pole FFC cable
9-31611-00	Overlay front
9-31612-00	Overlay back
9-31613-00	Rubber feet x 4, not in picture
9-31614-00	Chassis
9-31615-00	Cover, not in picture

Spare parts list Ergonomic Sealing Handle MSH-III





9-35401-00	Right cover	9-35411-00	Pushing spring
9-35402-00	Left cover	9-35412-00	Piston complete incl. coil
9-35403-00	Handle	9-35413-00	Moving electrode
9-35404-00	Protection cover	9-35414-00	Fixed electrode
9-35405-00	Internal coax cable with SMA contact	9-34515-00	PC board complete incl. micro switch (Not in picture)
9-35405-01	Washer for SMA contact	9-34516-00	3 Socket head cap screw for screen plates
9-35406-00	Indicator lens	9-35417-00	Screw Torx T25x8
9-35407-00	Screen plate bottom	9-35418-00	Screw Torx T30x8
9-35408-00	Screen plate top	9-35419-16	Coax cable coiled 1.2m, max 1.6m drawn out, total 1.75m (Not in picture)
9-35409-00	Pushing knob		
9-35410-00	Pull back spring		

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Spare parts list Bench Unit

9-33401-00	Cover	9-33411-00	Micro switch
9-33402-00	Signal cable	9-33412-00	Welding house chassis complete
9-33403-00	Locking device	9-33413-00	Pull back spring
9-33404-00	Chassis BNC contact	9-33414-00	Steering pin complete
9-33405-00	Hot moving electrode	9-33415-00	Pin compl. with cases and locking nut
9-33406-00	Fixed electrode	9-33416-00	Clamp spring
9-33407-00	Circuit board for connections	9-33417-00	Piston complete incl. coil and dim glowing lamp
9-33408-00	Solenoid	9-33418-00	Dim glowing lamp
9-33409-00	Chassis	9-33419-00	Welding house cover, micro switch side
9-33410-19	Coax cable 1.9m	9-33420-00	Welding house cover
9-33410-42	Coax cable 4.2m	9-33421-00	Coax cable internal
9-33410-91	Coax cable 9.1m	9-33422-00	HF-cable flexible

Technical data

Power Unit

Battery:	NiMh battery pack 24V/ 3,8Ah, 170x34x68 (LxWxH) mm.
Internal fuse:	Overload protection type PTC, self-recovery, and thermo fuse.
Power output:	100W /50 Ohm.
Seal time:	Steppless automatic 1-2 sec, or manual steppless setting 2-10 sec.
Frequency:	40,68 MHz.
Dimensions:	266x119x88 mm (LxWxH) incl. HF contact.
Weight:	2,15 Kg incl. battery.
Rel. humidity incl. storage and transport:	10-95% not condensing
Working temperature:	10-40 °C
Storage- and transportation temperature:	-40- +50 °C. (The working life of the batteries is shortened at high temperatures).
Operation:	Recommended max 1 seal each 3rd sec during continuous use or max 50 seals in sequence and after that 15 minute rest.
Noice level:	~60dB during normal operation
Manual pressing force:	Manual sealing handle: ~11N Ergonomic sealing handle: ~20N

Protection against electrical chock: Class II type B 

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Bench Unit

Dimensions:	150 (L) x 60 (W) x 73(H) mm
Weight:	0.9 Kg
Length of cable:	1.9 m (4.2 m or 9.1 m is available as option)

Handle Unit

Dimensions:	206 (L) x 27 (W) x 35 (H) mm
Weight:	0.3 kg
Length of cable:	1.9 m (4.2 or 9.1 m cable is available as option)

MSH III

Dimensions:	180 (L) x 35 (W) x 145 (H) mm
Weight:	0.220 kg
Length of cable:	Coiled cable 1.2m (Total length 1.75m)

Limited Warranty

The manufacturer hereby guarantees the original buyer that Biosealer CR6-AA is manufactured in a professional and quality manner, and will be free from all faults during a period of one year from the date of delivery from the manufacturer. The warranty includes equipment and components that proves to have faults during the warranty period.

The manufacturer will without cost for the customer, repair or replace the equipment that is faulty. The warranty is not valid for the battery pack since this is considered to be expendable supply. The warranty is not valid if the equipment has been repaired by anyone else than qualified personnel, that is approved by the manufacturer. The warranty is not valid if the equipment has been changed in any way that according to the manufacturers opinion, affects the reliability or stability of the instrument. The warranty is not valid when the serial number has been changed, crossed over or been removed, or if the fault has been caused by misuse or abnormal use. In these cases the manufacturer or the manufacturers representative will inform the customer about the decision, and if wished by the client will repair the equipment for normal rate. An estimated price can be given on request.

Includ

Manu
Addre
Telep
Telef

Mark
Add
Tele
Tele

Typ
Mo
Pro

I, t
an

Declaration of Conformity

according to the
Medical Devices Directive, 93/42/EEC
and
RoHS directive, 2011/65/EU

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Type of Equipment: Biosealer
Model: CR6-AA
Product class: Class 1

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
/ Thommy Lundström

Position/title: CEO

The document contains 22 bound numbered and sealed pages.

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Ljungberg & Kögel AB

Important

This user's manual is written for the person responsible for the operation of Biosealer CR6-PS. The operational methods and routines are developed and tested to ensure a reliable, safe and effective operation of CR6-PS. It is important that the operator has studied and understood the contents of this manual before using CR6-PS.

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Electrical safety

Only use the battery charger intended for CR6-PS or charger of another manufacturer recommended by the manufacturer of the CR6-PS and make sure it is marked with the accurate mains voltage.

Disposal of CR6-PS

As indicated by the  -symbol, Biosealer CR6-PS contains electronic components which are classified as dangerous for the environment and the equipment must be left to a recycling station when disposed. Local regulation decides how to recycle the CR6-PS, but it is preferred to disassemble the CR6-PS in its main parts and recycle the metal parts, plastic part, electronics etc. The batteries must be left to a recycling station and recycled according to local regulation.

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Revision history:

2006-04-21 TI	Trouble shooting p10 completed with check length of cable.
2007-01-16 TI	Warrant is not valid for battery pack
2008/01/25 TI	Change charger connected to mains power
2008-08-08 TI	New change charger connected
2009-11-03 TI	Change back to charger can be connected, Ergonomic Sealing Handle included, P6 timer adjust
2010-09-06 TI	MDD changed
2012-11-09 JA	Added Abelko Innovation commitment concerning technical feedback
2014-05-05 DC	Address in DoC updated
2015-04-13 DC	Updated information about battery charger
2015-04-20 DC	Cable length for MSH updated

Table of contents:

Important	2
Electrical safety	2
Disposal of CR6-PS	2
Features	4
Function.....	5
Operation.....	5
Battery test.....	5
Temperature-alarm	6
Timer adjust	6
Charging batteries	6
Maintenance and service	7
Cleaning the electrodes	7
Change of batteries	8
Technical description	8
Block diagram	9
Sealing Function.....	9
Trouble shooting	10
Spare parts Power Unit.....	11
Power cord sets for battery charger	11
Spare parts Handle Unit	12
Spare parts Ergonomic Sealing Handle MSH-II	13
Technical data	14
Warranty.....	14
MDD Declaration.....	16


Features

Biosealer CR6-PS is a battery operated sealing unit for PVC-tubes using environmentally approved cadmium-free NiMH batteries. With new and fully charged batteries it manages more than 500 seals on 4-5 mm PVC-tubes. The high capacity makes it suitable to use CR6-PS both in donation rooms and laboratories.

Biosealer CR6-PS is comprised of 3 units:

- Manual Handle Unit type 1 or Ergonomic Sealing Handle.
- Power Unit with batteries, RF generator and a control unit monitored by a microprocessor.
- Battery charger for NiMH batteries

Optional there is a strap with fastenings for mobile use.

Ergonomic handle:

Func

Sealing
which
checks
diode i
sealing

1. Ens
2. Tur
sea
Un
the
yel
3. Im
the
lov
4. A

Note

- No
ca
- N
- C
P
th
- A
E
- S
d
T
c

Batt
indi
oran

Function

Sealing is done in about 1 sec and the width of the seal is about 3 mm with a marking in the middle which makes it easier to pull the tubes apart. Before sealing, the microprocessor controlled timer board checks that the tube not is damp, to avoid sparks and thereby leakage. After each seal a multicoloured diode is lit. The diode indicates the remaining capacity of the batteries. This can also be checked before sealing with the push button at the front of the Power Unit.

Operation

1. Ensure that the power switch on the front is in off position and connect the Handle Unit.
2. Turn the power switch on, place the PVC-tube between the electrodes and press the handle. The sealing automatically starts and two indications is given. The yellow diode "SEALING" on the Power Unit is lit and the orange lamp, indicating high frequency power, on the handle is lit. The intensity of the orange lamp is reduced during the sealing and is turned off when the seal is completed. The yellow diode is also turned off when the seal is completed
3. Immediately after the sealing the multicoloured diode on the Power Unit is lit. It shows green when the remaining capacity of the battery is 70-100%, orange for 25-70% and red when the capacity is lower than 25%.
4. A new seal can be made immediately after the diode "sealing" has been turned off.

Note:

- Never have a shorter distance than 5 cm between two seals. Otherwise the pressure in the tube may cause leakage.
- Never stretch the tube during sealing, since this may cause leakage.
- Continuous sealing with an interval of 1 seal/sec during several minutes can cause overheating of the Power Unit or the fixed electrode with a declined quality of the seals or cause "TEMP ALARM" on the Power Unit.
- Always use the delivered standard coax cable 2.3 m for the metal handle and the coiled 2.15m for the Ergonomic Handle. The cable may not be shortened.
- Sealing can be done with or without the charger connected to mains power. It is recommended to disconnect the charger from mains power if there is long periods, 10-30 minutes, between sealing. The reason is that the battery may be overcharged otherwise. In continuous usage it is ok to have the charger connected.

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Temperature-alarm

If the temperature of the high frequency generator exceeds alarm level, set to 60 °C, the diode "TEMP ALARM" is lit. Sealing can not be done before the RF- generator has cooled down and the diode is turned off.

Timer adjust

If the RF-power stops before the seal is complete the sealing time needs to be adjusted. A ticker tube usually needs a longer sealing time. The sealing time can be set to 0.5-5 sec with the potentiometer behind the cover "TIMER ADJ". Avoid setting the timer too long since this will decrease battery capacity and heat up the electrodes.

1. Remove the cover.
2. Adjust the timer with a little screwdriver, turning clockwise will increase the time.
3. Put the cover back when adjustments are completed.

Charging batteries

The battery charger for CR6-PS is specially designed for charging the NiMH batteries powering CR6-PS.

It has a controlled charging procedure which is vital to avoid overcharge of batteries since this will shorten the life. The controller is using both ΔV sensing and timeout limits.

(ΔV sensing is a method to recognize when battery is fully charged since at this point the voltage starts to decrease instead of increase depending on temperature rise in battery).

More information about the charger can be obtained from the manufacturer of the CR6-PS.

Note! Only battery charger intended for the CR6-PS or a charger of another manufacturer recommended by the manufacturer of the CR6-PS may be used

1. Check that the local mains supply corresponds to the charger type sign.
2. Connect the power supply cord to the charger and then to the AC-line and the low voltage connector to the input at the back side of the CR6-PS marked "BATTERY CHARGER INLET".
3. A battery with a capacity of 25% is fully loaded within 1-2 hours.
4. Sealing can be done when the charger is connected.
5. It lacks significance whether CR6-PS is turned on or off during charging.
6. Avoid to completely unload the batteries since this will shorten the life of the batteries.
7. The life of the batteries depends on the amount of charging cycles and the surrounding temperature included the temperature rise that occurs when batteries are overcharged. Typical expected life time is 3 years with a charging frequency of 200 cycles/year.

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

(6) Replace the coax cable.

Change of batteries

Change of batteries must be done by approved service personnel.
The battery charger and the handle must be disconnected when the batteries are changed.
Always use battery packs from Ljungberg & Kögel

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1. Unscrew the 4 screws holding the back plate and remove it.
2. Lift out the battery pack and loosen the connections to the battery.
3. Replace with a new one of the same model and connect the cables, red to red and black to black, and put the cables back as they were before.
4. Put the back plate back and fasten it with the 4 screws.


Technical description

Biosealer CR6-PS consists of a radio frequency (RF) generator operating at 40,68 MHz with a maximum effect of 90W. RF-generator HFG-02 is a complete module mounted in the Power Unit. HFG-02 is provided with current (24 V) by a NiMH battery pack. The batteries do not need maintenance and are charged with a special NiMH charger.

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Block diagram**Sealing Function**

The sealing function is controlled by the timer board mounted on the HFG-02. The timer board consists of a microprocessor which governs all functions. When the Handle Unit is pressed, the connection from the central conductor in the coax cable to ground is released and a DC-voltage of 12 V is supplied from the timer board.

This voltage is recognized by the microprocessor (input PB0) on the timer card. The processor first checks if there is any leakage of current between the electrodes, if so this indicates damp on the tube. The limit is set to 5 μ A, which corresponds to 2,5 V at input PB1 on the processor. If the value exceeds this limit, the seal is cancelled and an indication is given by the yellow diode "SEALING" which is flashing. This means that the tube is damp and that there is a risk for sparking and leakage.

If no damp is found, the transistor T3 turns on the RF generator and keeps it on during the period determined by the timer potentiometer P1. Thereafter the RF is turned off.

About 100 msec before the RF is turned off, the battery capacity is measured with the input PB2 on the microprocessor. After completed sealing the diode for the battery is lit and during 3 sec the capacity is shown. Green light for 70-100% left, orange light for 25-70% and red light if the capacity is below 25%.

The capacity of the batteries also can be tested before sealing, by pushing the "BATT TEST" button at the front of the Power Unit. The microprocessor then discharge the batteries with approximately 0,5 A during 0,5 sec with transistor T7 and resistor R19, and measures the voltage at input PB5 and indication is given in the same way as described above.

The temperature of the RF-module is measured with a termistor mounted on the rear part above the RF-module and if the temperature exceeds programmed set level 60 °C the start of sealing is cancelled and the red diode "TEMP ALARM" is lit.

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
Trouble shooting

Sealing does not start and the diode "Sealing" is not lit when the handle is pressed.	• Test another handle. • Change the coax cable
The diode "SEALING" is lit but not the lamp on the handle.	• Test another handle. • Change the coax cable • Change the RF-generator. (Only for service personnel).
The diode "TEMP ALARM" is lit and sealing does not start.	• The RF-generator is overheated. Let it cool down until the alarm is turned off.
The diode "SEALING" flashes and sealing does not start.	• The tube is damp, dry it.
Sparks appear by the electrodes during sealing.	• The electrodes are not parallel. Contact service personnel.
The seal is done but with bad quality, leakage may appear.	• Check that the coax cable is the right type, the length must be 2.3m for CR6PS.

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Spare parts Power Unit


9-81601-00	HFG-02 RF-generator
9-81602-00	Battery pack NiMH 24V/1100mA incl. wires and connectors
9-81603-00	Timer PC- board
9-81604-00	Front plate
9-81605-00	Back plate
9-81606-00	Cover upper half
9-81607-00	Cover bottom half
9-81608-00	Strap
9-81609-00	Strap holder incl. 2 screws
9-81610-00	ITW Fixlock ETG20
9-81611-00	Front and back plate countersunk screws x 4
9-81612-00	Rubber feet x 4 (Not in picture)
9-81613-00	On/off -switch (Not in picture)
9-81614-00	Battery. test button (Not in picture)
9-81615-00	Charger inlet (Not in picture)
9-81616-00	IC CPU ST6 CR-6PS (Not in picture)
9-81617-00	IC LM358 op-amp. (Not in picture)

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
9-35401-00	Right cover	9-35411-00	Pushing spring
9-35402-00	Left cover	9-35412-00	Piston complete incl. coil
9-35403-00	Handle	9-35413-00	Moving electrode
9-35404-00	Protection cover	9-35414-00	Fixed electrode
9-35405-00	Internal coax cable with SMA contact	9-34515-00	PC board complete incl. micro switch (Not in picture)
9-35406-00	Indicator lens	9-34516-00	3 Socket head cap screw for screen plates
9-35407-00	Screen plate bottom	9-35417-00	Screw Torx T25x8
9-35408-00	Screen plate top	9-35418-00	Screw Torx T30x8
9-35409-00	Pushing knob	9-35419-17	Coax cable coiled 1.2m, max 1.7m drawn out, total 2.15m. (Not in picture)
9-35410-00	Pull back spring		

Technical data

Batteries:	Battery pack 24V NiMH (Ljungberg&Kögel)
Battery safety:	Overload protection type PTC, self recovery, and thermo fuse.
Effect consumption:	120W maximum effect
Sealing time:	Stepless adjustable 0,5-5 sec
Frequency:	40,68 MHz crystal controlled
Output effect:	90W/50Ω maximum effect
Dimensions:	167x141x35 mm (LxWxH) incl. RF-outlet
Weight:	1.3 Kg incl. batteries
Rel. Humidity incl. storage and transport:	10-95% not condensing
Working temperature:	10-40 °C
Storage- and transport- temperature:	-40- +50 °C. (The life length of the batteries is shortened at high temperatures).
Operation:	Recommended max 1 seal each 3:rd sec during continuous use.
Noice level:	~60dB during normal operation
Manual pressing force:	Manual sealing handle: ~11N Ergonomic sealing handle: ~20N
Symbol: 	Protection against electrical chock: class II type B.

This device complies with requirements in MDD 93/42 EEC

Manufactured by: Abelko Innovation
Industrivägen 17
972 54 Luleå
Marketing: Ljungberg&Kögel AB
Box 1032
251 10 Helsingborg

Warranty

The manufacturer hereby guarantees the original buyer that Biosealer CR6-PS is manufactured in a professional and quality manner, and will be free from all faults during a period of one year from the date of delivery from the manufacturer. The warranty includes equipment and components that proves to have faults during the warranty period. The manufacturer will without cost for the customer, repair or replace the equipment that is faulty. The warranty is not valid for the battery pack since this is considered to be expendable supply. The warranty is not valid if the equipment has been repaired by anyone else than qualified personnel, that is approved by the manufacturer. The warranty is not valid if the equipment has been changed in any way that according to the manufacturer's opinion, affects the reliability or stability of the instrument. The warranty is not valid when the serial number has been changed, crossed over or been removed, or if the fault has been caused by misuse or abnormal use. In these cases the manufacturer or the manufacturer's representative will inform the customer about the decision, and if wished by the client will repair the equipment for normal rate. An estimated price can be given on request.

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
Type of Equipment: Biosealer
Model: CR6-PS
Product class: Class 1

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
Position/title: President

1100

The document contains 18 bound numbered and sealed pages.

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Biosealer CR6-PS AA

Including automatic timer

User's manual
and
Technical reference

Ljungberg & Kögel AB

Important

This user's manual is written for the person responsible for the operation of Biosealer CR6-PS AA. The operational methods and routines are developed and tested to ensure a reliable, safe and effective operation of CR6-PS AA. It is important that the operator has studied and understood the contents of this manual before using CR6-PS AA.

- This instrument is a sealing equipment using radio frequency for welding and the electrical emission, at the operating frequency 40,68 MHz, is high. Make sure that other instruments and equipment near the sealing unit can withstand this emission.
- Never touch the electrodes with your fingers during the sealing period! This may cause burn damages!

Electrical safety

Only use the battery charger intended for CR6-PS AA or charger of another manufacturer recommended by the manufacturer of CR6-PS AA and make sure it is marked with the accurate mains voltage.

Disposal of CR6-PS AA

As indicated by the -symbol, Biosealer CR6-PS AA contains electronic components and must be left to a recycling station when disposed. Local regulation decides how to recycle the biosealers, but it is preferred to disassemble the biosealer in its main parts and recycle the metal parts, plastic part, electronics etc.

Environmentally effects

Biosealer CR6-PS AA contains none or insignificant parts of the forbidden hazardous substances in RoHS directive 2002/95/EG. The batteries must be left to a recycling station and recycled according to local regulation

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Table

Important
Electrical
Disposal
Environment
Manufact
Features..
Function..
Operation
Manual se
Temperatu
Battery tes
Charging t
Maintenance
Cleaning
Change of
Technical c
Block di
Sealing F
Trouble sho
Spare parts
Power cord
Spare parts
Spare par
Technical d
Warranty ...
MDD Decla

Revision history:
2006-04-21 TI
2007-01-16 TI
2008/01/25 TI
2008-08-08 TI
2009-11-03 TI
2010-09-06 TI
2011-05-23 TI
2011-08-08 TI
2011-11-28 TI
2012-11-09 JA
2012-12-04 TI
2012-12-13 DC
2013-01-15 DC
2014-04-17 DC
2014-05-07 DC
2014-06-02 DC
2015-04-14 DC
2015-04-20 DC

Table of contents:

Important	3
Electrical safety	3
Disposal of CR6-PS AA	3
Environmentally effects	3
Manufacturer	3
Features	3
Function	5
Operation	6
Manual setting of timer	6
Temperature alarm	6
Battery test	7
Charging the batteries	7
Maintenance and service	7
Cleaning the electrodes	8
Change of batteries	8
Technical description	9
Block diagram	9
Sealing Function	10
Trouble shooting	10
Spare parts Power Unit	11
Power cord sets for battery charger	12
Spare parts Metal Handle	12
Spare parts Ergonomic Sealing Handle MSH-III	13
Technical data	14
Warranty	15
MDD Declaration	15
	16

Revision history:

2006-04-21 TI	Trouble shooting p10 completed with check length of cable.
2007-01-16 TI	Warrant is not valid for battery pack
2008/01/25 TI	Change charger connected to mains power.
2008-08-08 TI	New change charger connected
2009-11-03 TI	Change back to charger can be connected, Ergonomic Sealing Handle included, P6 timer adjust,
2010-09-06 TI	MDD changed
2011-05-23 TI	AA introduced + some other changes
2011-08-08 TI	MDD changed
2011-11-28 TI	MSH-III spare parts coax cable better specified
2012-11-09 JA	Added Abelko Innovation commitment concerning technical feedback
2012-12-04 TI	Check of charging led
2012-12-13 DC	Partlist for MSH3 updated with 5.1
2013-01-15 DC	Title update to contain AA
2014-04-17 DC	Address in DoC updated
2014-05-07 DC	RoHS added in DoC
2014-06-02 DC	REACH added to DoC
2015-04-14 DC	Note about charger added
2015-04-20 DC	Cable length for MSH updated


Features

Biosealer CR6-PS AA is a battery operated sealing unit for PVC, PVC/EVA and EVA-tubes using environmentally approved cadmium-free NiMH batteries. With new and fully charged batteries it manages more than 500 seals on 4-5 mm PVC-tubes. The high capacity makes it suitable to use CR6-PS AA both in donation rooms and laboratories.

Biosealer CR6-PS AA is comprised of 3 units:

- Manual metal handle unit with straight coax cable or an ergonomic handle with coiled coax cable.
- Power Unit with batteries and a carry strap, RF generator and a control unit monitored by a microcontroller.
- Battery charger for NiMH batteries

Optional there is a carry bag for mobile use.

Sealing
makes i
that the
The dio
the push

1. E
2. T
se
3. T
se
rel
4. O
ora
ha
5. Im
sho
who
6. A n

Note:

- Never ha
cause lea
- Never str
- **Do not re**
- Continuo
the Power
the Power
- Always us
- Sealing ca
disconnect
The reason
charger co

Note! Manual
except in spec

Function

Sealing is done in 1-2 sec and the width of the seal is about 3 mm with a marking in the middle which makes it easier to pull the tubes apart. Before sealing, the microcontroller controlled timer board checks that the tube is not damp to avoid sparks and thereby leakage. After each seal a multicolored diode is lit. The diode indicates the remaining capacity of the batteries. This can also be checked before sealing with the push button at the front of the Power Unit.

Operation

1. Ensure that the power switch on the front is in OFF position and connect the Handle Unit.
2. Turn the power switch ON, place the PVC-tube between the electrodes and press the handle. The sealing automatically starts and the yellow diode "SEALING" on the Power Unit is lit.
3. The high frequency is indicated with a blinking blue lamp on the ergonomic handle. When the seal is ready the blinking blue lamp become green indicating seal ready and the handle can be released.
4. On the metal handle an orange lamp is indicating high frequency power. The intensity of the orange lamp is reduced during the sealing and is turned off when the seal is completed. *Wait half a second to let the plastic stiffen before releasing the handle.*
5. Immediately after the sealing the multicolored diode BATT COND on the Power Unit is lit. It shows green when the remaining capacity of the battery is 70-100%, orange for 25-70% and red when the capacity is lower than 25%.
6. A new seal can be made immediately after the diode "SEALING" has been turned off.

Note:

- Never have a shorter distance than 5 cm between two seals. Otherwise the pressure in the tube may cause leakage.
- Never stretch the tube during sealing, since this may cause leakage.
- *Do not release the handle while welding is in progress since this may cause leakage.*
- Continuous sealing with an interval of 1 seal/sec during several minutes may cause overheating of the Power Unit or the electrodes with a declined quality of the seals or cause "TEMP ALARM" on the Power Unit.
- Always use the delivered coax cables. *The cables may not be shortened or lengthened.*
- Sealing can be done with or without the charger connected to mains power. It is recommended to disconnect the charger from mains power if there is long periods, 10-30 minutes, between sealing. The reason is that the battery may be overcharged otherwise. In continuous usage it is ok to have the charger connected.

Manual setting of timer

Note! Manual setting of timer will shut off the automatic timer function. It is not recommended except in special cases. Consult L&K's representative if in doubt.

When sealing thick or stiff PVC or EVA tubes it may be necessary to increase the seal time more than what you get with the autotimer. The timer can be manually set with the potentiometer TIMER ADJ with a small screwdriver. The potentiometer set to >20% clockwise will disconnect the auto timer function and give a fixed timer. The potentiometer in the middle will give 5 sec and the max clockwise will give 10 sec.

1. Remove the plastic cover TIMER ADJ on the frontplate.
2. Adjust the timer with a little screwdriver, turning clockwise will increase the time.
3. Put the cover back when adjustments are completed.

Temperature alarm

If the temperature of the high frequency generator exceeds alarm level, set to 75°C, the diode "TEMP ALARM" is lit. Sealing can not be done before the RF- generator has cooled down and the diode is turned off.

Battery test

Battery test can be done at any time by pressing the button "BATT TEST". The "BATT STATUS" indicator will show the battery capacity in the same way as after each seal with green light for 70-100%, orange light for 25-70% and red light when capacity is less than 25%.

Charging the batteries

The "Battery charger for CR6-PS AA" is specially designed for charging the NiMH batteries powering CR6-PS AA. It is medical approved according to EN 60601-1 and UL recognized for US and Canada.

It has a controlled charging procedure which is vital to avoid overcharge of batteries since this will shorten the working life. The controller is using both DeltaV sensing and timeout limits.

(DeltaV sensing is a method to recognize when battery is fully charged since at this point the voltage starts to decrease instead of increase depending on temperature rise in battery).

Note! Only battery charger type "Battery charger for CR6-PS AA" or another charger recommended by L&K may be used.

1. Check that the local mains supply corresponds to the charger type sign.
2. Connect the power supply cord to the charger and then to the AC-line and the low voltage connector to the input at the back side of the CR6-PS AA marked "BATTERY CHARGER INLET". It is ok to connect the charger to CR6-PS AA first and then to AC-line.
3. A battery with a capacity of 25% is fully loaded within 1-2 hours.
4. Since this is a NiMH battery pack it should never be overcharged. The recommendation for longest working life is to use the sealer unconnected to mains until the BATT COND led turns orange and then recharge it.
5. Sealing can be done when the charger is connected.
6. It lacks significance whether CR6-PS AA is turned on or off during charging.

7. Avoid
8. The w
tempe
expec
discha
9. If the
- 10.If the
discha

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Maintena
approved

Alv
cle

The Powe

Cleaning

Metal ha
Clean the
cable.

Ergonom

- (1) Remo
- (2) Remo
the fro
- (3) Then
- (4) Clean
- (5) Dry u
and li
- (6) Repla

7. Avoid to completely unload the battery since this will shorten the life of the battery.
8. The working life of the battery depends on the amount of charging cycles and the surrounding temperature included the temperature rise that occurs when batteries are overcharged. Typical expected working life is 3 years with a charging frequency of 200 cycles/year and 25°C. Deep discharging and overcharging shortens working life.
9. If the sealer not is used for a while the battery needs maintenance recharging every 3 months.
10. If the battery has been completely discharged during some time it may need repeatedly charging and discharging to obtain full capacity.

Check that the charging led on the charger will be orange or red when the charging is started. If it is green directly after start the mains cable must be unplugged and reconnected once again.

Maintenance and service

Maintenance to be done by the operator is limited to cleaning. All other service must be done by an approved service personnel.

Always turn off the Power Unit and disconnect the battery charger from the mains supply before cleaning the equipment.



The Power Unit and the handle can be cleaned with a mild soap solution.

Cleaning the electrodes

Metal handle: Remove the coaxial cable. Unscrew the 2 screws holding the plastic cover and remove it. Clean the electrodes with distilled water and dry with a soft cloth. Replace the plastic cover and the coax cable.

Ergonomic Handle :

- (1) Remove the coaxial cable from the handle unit
- (2) Remove the protective cover by first pressing the front down as shown in the picture.
- (3) Then pull the cover forwards.
- (4) Clean both electrodes with distilled water or alcohol.
- (5) Dry up the electrodes carefully with a soft and lint-free cloth.
- (6) Replace the protective cover and the coax cable.

Change of batteries

When the battery capacity begins to be insufficient for one days sealing, the batteries has to be changed.

Change of batteries must be done by approved service personnel.
The battery charger and the handle must be disconnected when the batteries are changed.
CR6-PS AA must be turned OFF.
Always use battery packs from L&K

1. Unscrew the 4 screws holding the back plate and remove it.
2. Lift carefully out the battery pack and loosen the connections to the battery.
3. Replace with a new one of the same model and connect the cables, red to red and black to black, and put the cables back as they were before.
4. Put the back plate back and fasten it with the 4 screws.

Technical description

Biosealer CR6-PS AA consists of a high frequency (HF) generator operating at 40,68 MHz with a maximum effect of 100W. HF-generator HFG-02 is a complete module mounted in the Power Unit. HFG-02 is provided with current (24 V) by a NiMH battery pack. The batteries do not need maintenance and are charged with a special NiMH charger.

Note!

Due to regulations of maximum radiation at other frequencies than 40,68 MHz the HFG-02 module is soldered up and shall be regarded as a component. It can not and it is not allowed to be adjusted or repaired by anyone else than the manufacturer. In the event of malfunction of HFG-02 it must be sent to the manufacturer for repair.

The voltage that occurs during sealing may cause burn damages if someone touches the fixed electrode during the sealing process. Therefore it is absolutely necessary to turn off the power before cleaning or service of any of the units.

Block d

E
b

Sealing

The sealing
of a micro
the central
the timer b

This voltage
there is any
set to 5µA,
yellow diode
sparking ar


If no damp
auto adjust
500msec ar

The battery
the battery i
25-70% and

The capacit
the front of

The tempera
module and

Block diagram

Sealing Function

The sealing function is controlled by the timer board mounted on the HFG-02. The timer board consists of a microcontroller which governs all functions. When the Handle Unit is pressed, the connection from the central conductor in the coax cable to ground is released and a DC-voltage of 5V is supplied from the timer board.

This voltage is recognized by the microcontroller on the timer board and the processor first checks if there is any leakage of current between the electrodes, if so this indicates damp on the tube. The limit is set to 5 μ A, and if the value exceeds this limit, the seal is cancelled and an indication is given by the yellow diode "SEALING" which is flashing. This means that the tube is damp and that there is a risk for sparking and leakage.

If no damp is found, the HF generator is started and keeps on during a time that is determined by the auto adjust function in the microcontroller. Thereafter the RF is turned off and there is cooling time of 500msec and after that the green light turns on the ergonomic handle.

The battery capacity is measured once a minute with a load of 100 Ohm and after each seal the diode for the battery is lit and during 3 sec the capacity is shown. Green light for 70-100% left, orange light for 25-70% and red light if the capacity is below 25%.

The capacity of the batteries can also be tested before sealing, by pushing the "BATT TEST" button at the front of the Power Unit.

The temperature of the RF-module is measured with a termistor mounted on the rear part above the RF-module and if the temperature exceeds programmed set level 75 °C the start of sealing is cancelled and

the red diode "TEMP ALARM" is lit. If the electrodes are shortened during sealing the RF-generator is immediately stopped.

Trouble shooting


Sealing does not start and the diode "SEALING" is not lit when the handle is pressed.	• Test another handle. • Change the coax cable
The diode "SEALING" is lit but not the lamp on the handle.	• Test another handle. • Change the coax cable • Change the RF-generator. (Only for service personnel).
The diode "TEMP ALARM" is lit and sealing does not start.	• The RF-generator is overheated. Let it cool down until the alarm is turned off.
The diode "SEALING" flashes and sealing does not start.	• The tube is damp, dry it.
Sparks appear by the electrodes during sealing.	• This can happen if the electrodes are overheated due to too many seals in a short time. Let the electrodes cool down.
Sparks appear by the electrodes during sealing.	• If there has been a leakage and flashes have created carbon particles on the electrodes, this particles must be removed completely before sealing can continue.
Sparks appear by the electrodes during sealing.	• It can also indicate that the electrodes are not parallel. Contact service personnel.
The seal is done but with bad quality, leakage may appear.	• Check that the coax cable is the right type, the cable must be marked CR6-PS AA.

Note: L&K will supply more detailed service-instructions with component lists and trouble-shooting guides on request.

9-81601-
9-81602-
9-81603-
9-81604-
9-81605-
9-81606-
9-81607-
9-81608-
9-81609-
9-81610-
9-81611-
9-81612-
9-81613-
9-81614-
9-81615-

9-81601-5
9-81601-5
9-81601-5
9-81601-5

Spare parts Power Unit


9-81601-00	HFG-02 RF-generator
9-81602-00	Battery pack NiMH 24V/1100mA incl. wires and connectors
9-81603-01	Timer PC-board
9-81604-01	Front plate
9-81605-00	Back plate
9-81606-00	Cover upper half
9-81607-00	Cover bottom half
9-81608-00	Strap
9-81609-00	Strap holder incl. 2 screws
9-81610-00	ITW Fixlock ETG20
9-81611-00	Front and back plate countersunk screws x 4
9-81612-00	Rubber feet x 4 (Not in picture)
9-81613-00	On/off-switch (Not in picture)
9-81614-00	Battery. test button (Not in picture)
9-81615-00	Charger inlet (Not in picture)

Power cord sets for battery charger

9-81601-51	Power cord set (EUR)
9-81601-52	Power cord set (UK)
9-81601-53	Power cord set (US/JPN)
9-81601-54	Power cord set (AUS)

Spare parts Metal Handle

9-32401-00	Cover
9-32402-00	BNC chassis female contact
9-32403-00	Switch normally closed
9-32404-00	Internal coax cable
9-32405-00	Coil incl. delrin core and dim glowing lamp
9-32406-00	Dim glowing lamp
9-32408-00	Plate spring
9-32409-00	Handle
9-32410-00	Chassis
9-32411-00	Pin
9-32413-00	Rubber tile stop
9-32414-00	Electrode house complete
9-32415-00	Electrode cover
9-32416-23	Coax cable 2.3m
9-32417-00	Rubber tile switch

Spare parts Ergonomic Sealing Handle MSH-III

9-35401-00	Right cover	9-35411-00	Pushing spring
9-35402-00	Left cover	9-35412-00	Piston complete incl. coil
9-35403-00	Handle	9-35413-00	Moving electrode
9-35404-00	Protection cover	9-35414-00	Fixed electrode
9-35405-00	Internal coax cable with SMA contact	9-34515-00	PC board complete incl. micro switch (Not in picture)
9-35405-01	Washer for SMA contact	9-34516-00	3 Socket head cap screw for screen plates
9-35406-00	Indicator lens	9-35417-00	Screw Torx T25x8
9-35407-00	Screen plate bottom	9-35418-00	Screw Torx T30x8
9-35408-00	Screen plate top	9-35419-17	Coax cable coiled 1.2m, max 1.7m drawn out, total 2.15m. (Not in picture)
9-35409-00	Pushing knob		
9-35410-00	Pull back spring		

Technical data

Batteries:	Battery pack 24V/1100mAh NiMH (Ljungberg&Kögel)
Battery safety:	Overload protection type PTC, self recovery, and thermo fuse.
Effect consumption:	120W maximum effect
Sealing time:	Steppless Automatic or manual steppless adjustable 2-10 seconds
Frequency:	40.68 MHz crystal controlled
Output effect:	100W/50Ω maximum effect
Dimensions:	167x141x35 mm (LxWxH) incl. RF-outlet
Weight:	1.3 Kg incl. batteries
Rel. Humidity incl. storage and transport:	10-95% not condensing
Working temperature:	10-40 °C
Storage- and transportation temperature:	-40- +50 °C. (The working life of the batteries is shortened at high temperatures).
Operation:	Recommended max 1 seal each 3:rd sec during continuous use.
Noice level:	~60dB during normal operation
Manual pressing force:	Manual sealing handle: ~11N Ergonomic sealing handle: ~20N
Protection against electrical chock:	Class II type B 

This equipment complies with the demands in medical directive MDD 93/42 EEC

Manufactured by: Abelko Innovation
972 54 Luleå
Sweden

Warranty

The manufacturer hereby guarantees the original buyer that Biosealer CR6-PS AA is manufactured in a professional and quality manner, and will be free from all faults during a period of one year from the date of delivery from the manufacturer. The warranty includes equipment and components that proves to have faults during the warranty period. The manufacturer will without cost for the customer, repair or replace the equipment that is faulty. The warranty is not valid for the battery pack since this is considered to be expendable supply. The warranty is not valid if the equipment has been repaired by anyone else than qualified personnel, that is approved by the manufacturer. The warranty is not valid if the equipment has been changed in any way that according to the manufacturers opinion, affects the reliability or stability of the instrument. The warranty is not valid when the serial number has been changed, crossed over or been removed, or if the fault has been caused by misuse or abnormal use. In these cases the manufacturer or the manufacturer's representative will inform the customer about the decision, and if wished by the client will repair the equipment for normal rate. An estimated price can be given on request.

MDD Decl

Manufa
Address:
Telepho
Telefax

Market
Address
Teleph
Telefax

Type o
Model
Produ

I, the u
Regul

Date

Signa

Posi

DECLARATION OF CONFORMITY

according to the
Medical Devices Directive 93/42/EEC
and
RoHS directive, 2011/65/EU
and
REACH regulation, 1907/2006

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Manufacturer: Abelko Innovation
Address: Industrivägen 17, 972 54 Luleå
Telephone: +46 920-220360
Telefax: +46 920-220068

Marketing: Ljungberg&Kögel AB
Address: Box 1032 251 10 Helsingborg
Telephone: +46 042-139860
Telefax: +46 042-132181

Type of Equipment: Biosealer
Model: CR6-PS AA
Product class: Class 1

I, the undersigned, hereby declare that the equipment specified above conforms to the above Directives, Regulations and Standards.

Date of issue: 26 June 2014


Signature:

/ Thommy Lundström

Position/title: President

The document contains 16 bound numbered and sealed pages.

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