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Supra Phob

Hydrophobic Aspheric Optics Foldable Intraocular Lens

Patient Name

Date of Implant

Implanted Eye

L:

R:

Hospital

Lens Identification Label

Supra Phob

Preceded System
Instruction for Use

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Step 2:

Slightly move the lens protector Cap forward from the injector.

Step 3:

Check whether the lens placed properly in the cartridge, and to ensure that posterior side of the lens facing downwards.

Step 4:

Apply a thin layer of viscoelastic from the barrel to the tip of the cartridge to serve as a lubricant for the free movement of lens during injections.

Step 5:

Gently remove the lens protector cap and close the cartridge properly.

Step 6:

Make sure that the lens inside the cartridge should be properly loaded or neither the haptic nor optic is caught between the shutters.

Step 7:


Ensure that there is no gap between the shutters, and then gently push the plunger towards the lens.

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Mfgd. by
ALFA APMARQUE OCULAR DEVICES (P) LTD.
B-11, 8/1, 233, N. of 4th A, VILLAGE OF MAHARAJA,
MUMBAI-400 002, INDIA.
www.alfamarque.com

EMERGO EUROPE
Meenestraat 15,
2513 BH The Hague,
The Netherlands.

REF. : SPNT200 POWER: 23.0D
Ø B : 6.00mm Ø T : 13.00mm
B. No. : 3716 ⚙ : 12-2016
S. No. : 1112319 033
Supra Phob

CE 1023

REF. : SPNT200 POWER: 23.0D
Ø B : 6.00mm Ø T : 13.00mm
B. No. : 3716 ⚙ : 12-2016
S. No. : 1112319 033
Supra Phob

CE 1023

REF. : SPNT200 POWER: 23.0D
Ø B : 6.00mm Ø T : 13.00mm
B. No. : 3716 ⚙ : 12-2016
S. No. : 1112319 033
Supra Phob

CE 1023

PRODUCT INFORMATION

Supra Phob

Hydrophobic Foldable Intraocular Lens with Aspheric Optics

DESCRIPTION OF THE DEVICE

Intra Ocular Lenses are optical implants for the human crystalline lens in the visual correction of aphakia. The lens optic is cut out from HYDROPHOBIC ACRYLIC BLANKS with UV blocking properties.

MODE OF ACTION

Intra Ocular Lens (IOL) functions as a refracting medium to replace the natural lens in the correction of aphakia. The IOL can be placed in Posterior chamber. The following indications and contra indications are based on research of medical literature and are to be used only as guides. The list is indicative and not to be viewed as complete or comprehensive.

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POTENTIAL HAZARDS

The relevant hazards and side effects of Intra Ocular lens essentially the same as (but not limited to) those of routine cataract extraction. These include the following -

- Intra Ocular Foreign Body Debris
- Lens Implant Lux Anterior Chamber
- Secondary Cataract Formation
- Malpositioned Lens
- Aphakic Glaucoma
- Temporary Flat Anterior Chamber
- Retinal Detachment
- Lens Implant Lux Posterior Chamber
- Endothelial Corneal Dystrophy
- Vitreous Hemorrhage into the Anterior Chamber
- Infection
- Pupillary Block
- Corneal Dystrophy

A' CONSTANT

A' Constant value is to be located on the Box Cover Label. This number is given only as a guideline and not based on any clinical data. It is advised that the surgeon calculates his own value for each style and specifications of Intra Ocular Lens.

PACKAGE

The Lens is supplied in a sterile Package. Validity of sterilization applies as long as the sealed inner peel-pouch is not disturbed or damaged. Any damage to the peel-pouch or any accidental opening of the peel-pouch is to be declared as "NOT STERILE".

DIRECTIONS FOR USE

PLEASE examine the peel-pouch prior to opening to assure sterility. Any Damage to the peel-pouch should be viewed very seriously and the IOL inside might not be sterile. After proper examination, peel the pouch and remove the Lens case. Hold and lift the lens using smooth edged forceps. Use BSS (or) Saline to clear off the static charges that might have developed on the lens surface. Do not touch the central optic portion of the lens with forceps. Do not use rubber gloves dusted with talc powder. It may cause irritation.

REPORTING

Adverse reactions and potentially sight threatening complications that may reasonably be regarded as lens related

and that were not previously expected in nature, severity or degree of incidence should be reported to Appasamy Ocular Devices (P) Ltd. This information is requested from the all implanting surgeons in order to document potentially long term effects of IOL implants.

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RETURN GOODS POLICY

Appasamy Ocular Devices (P) Ltd., accepts returned lenses for exchange only. No cash refunds will be issued. To return lenses, you must first obtain a return authorization number from customer services department. No returned goods will be accepted without proper authorization number. Returned lenses should be shipped by traceable method. No credit will be given to lost or damaged lenses in shipment. Lenses will be replaced as long as they are returned within six months of their original invoice date.

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1. The lens is re-sterilized.
2. Lens is altered in any manner.
3. Lens is repackaged by any one other than Appasamy Ocular Devices (P) Ltd.

DISPOSING OF NON-STERILE OR CONTAMINATED MEDICAL DEVICES
There are no specific guidelines for disposing of non-sterile or contaminated medical devices. Follow standard procedures for disposing outdated or contaminated products.

PRECAUTIONS

- Do not sterilize the IOLs
- Do not Re-use the product
- Do not store above 45°C (Fragile, avoid shocks)
- Use only sterile BSS irrigating solution for rinsing / soaking
- Do not use if the pouch is opened
- Handle gently to avoid lens surface and haptic damage
- Use as soon as the IOL is out of lens case.

Mfgd by:
A&D APPASAMY OCULAR DEVICES (P) LTD.
R.S. No. 9/12 & 3, N.H.45-A, Villupuram Main Road,
Vadungulam, Puducherry - 605 102, INDIA.

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
Supra Phob-bby

Blue Blocking Yellow
Preloaded System

- ⚠ Attention, see instructions for use ⓘ
- ⚠ Achtung, Gebrauchsanweisung einsehen ⓘ
- ⚠ Atención, ver instrucciones de uso ⓘ
- ⚠ Attention, voir mode d'emploi ⓘ
- ⚠ Attenzione, vedere le istruzioni per l'uso ⓘ

STERILE EO

Do not resterilize

Do not reuse

Store away from sunlight

Store in dry place


T° < 40°C / 104°F

Do not use if package damaged

Contains : One (1) Hydrophobic Single Piece Biconvex Foldable Lens One (1) STERILE EO Delivery System

REF. : SupraPhob BBY POWER: 20.00D
Ø B : 6.00mm Ø T : 13.00mm
S. No. : P1310121 091
B. No. : 4587
A. CONSTANT : 118.0
Supra Phob - BBY

S/n P1310121 091

Optic 6.00 mm Power 20.00 D
Length 13.00 mm
Exp 9-2018 Ref SUPRAPHOB BBY

Patient Lens Implant
Identification Card

EU Representative :
EMERGO EUROPE
Molenstraat 15,
2513 BH The Hague,
The Netherlands

Supra Phob-bby

Blue Blocking Yellow

Preloaded System

Hydrophobic Foldable Intraocular Lens
with
Aspheric Optics

CE 1023

Patient Name

Date of Implant

Implanted Eye

L:

R:

Hospital

Lens Identification Label

REF. : SupraPhob BBY POWER: 20.00D
Ø_B : 6.00mm Ø_T : 13.00mm
B. No.: 4587 : 09-2018
S. No.: P1310121 091
Supra Phob - BBY CE 1023

REF. : SupraPhob BBY POWER: 20.00D
Ø_B : 6.00mm Ø_T : 13.00mm
B. No.: 4587 : 09-2018
S. No.: P1310121 091
Supra Phob - BBY CE 1023

REF. : SupraPhob BBY POWER: 20.00D
Ø_B : 6.00mm Ø_T : 13.00mm
B. No.: 4587 : 09-2018
S. No.: P1310121 091
Supra Phob - BBY CE 1023

REF. : SupraPhob BBY POWER: 20.00D
Ø_B : 6.00mm Ø_T : 13.00mm
B. No.: 4587 : 09-2018
S. No.: P1310121 091
Supra Phob - BBY CE 1023


REF. : SupraPhob BBY POWER: 20.00D
Ø_B : 6.00mm Ø_T : 13.00mm
B. No.: 4587 : 09-2018
S. No.: P1310121 091
Supra Phob - BBY CE 1023

Supra Phob

Preloaded System
Instruction for Use

Step 1:


Peel the blister cover and take the injector carefully.

Step 2:

Close the cartridge by pressing the wings of the cartridge. The Lens protecting cover will automatically come out. Click sound indicates the proper closure of cartridge.

Step 3:

Apply the disco elastic from the front barrel of the cartridge upto the pusher in the injector before implantation of the IOL.

Step 4:

Then gently push the plunger towards the lens.

Step 5:

The lens moves towards the cartridge. The lens will completely come out of the cartridge tip end.

Mfgd. by:

APPASAMY OCULAR DEVICES (P) LTD.
R.S.No. 9/1,2&3, N.H.45-A, Villupuram Main Road,
Vadamangalam, Puducherry - 605 102. INDIA

EC REP

EMERGO EUROPE
Molenstraat 15,
2513 BH The Hague,
The Netherlands.

CE 1023

PRODUCT INFORMATION

Supra Phob-bby

Blue Blocking Yellow
Preloaded System

Hydrophobic Foldable Intraocular lens with Aspheric Optics

DESCRIPTION OF THE DEVICE

Intra Ocular Lenses are optical implants for the human crystalline lens in the visual correction of aphakia. The lens optic is let out from HYDROPHOBIC ACRYLIC BLANKS with UV blocking properties.

MODE OF ACTION

Intra Ocular Lens (IOL) functions as a refracting medium to replace the natural lens in the correction of aphakia. The IOL can be placed in Posterior chamber. The following indications and contra indications are based on research of medical literature and are to be used only as guides. The list is indicative and not be viewed as complete or comprehensive.

YELLOW LENSES

Throughout life, the crystalline lens protects the retina from hazardous UV light. Cataract extraction removes this protective barrier and may accelerate retinal photo toxicity to control the amount of radiation transmitted to the retina after cataract surgery. IOLs were created with the capability of blocking UV radiation. All IOLs today block short wavelength UV radiation up to 400nm.

Blue-absorbing chromophores added to IOL material to expand protection and decrease the amount of UV-blue energy (400- to 500-nm range) reaching the retina. It is thought that adding a blue-blocking chromophore increased protection from photic retinopathy, because it simulated the properties of the natural crystalline lens.

Yellow IOLs approximate the UV and blue-light attenuating properties of the human crystalline lens. This IOL contains a covalently bound chromophore, that providing increased protection against harmful UV light while allowing optimal blue light transmission, improving contrast sensitivity in low light conditions especially in eyes at risk of age-related macular degeneration.

INDICATIONS

- Monocular Cataract
- Mature Cataract
- Congenital Cataract
- Occupational Needs
- Traumatic Cataract

CONTRAINDICATIONS

- Chronic Severe Uveitis
- Epithelial Dystrophy
- Proliferative Diabetic Retinopathy
- Choroidal Hemorrhage
- Microphthalmos
- Aniridia
- Concomitant Severe Eye Disease
- Glaucoma Problem
- Rubella Cataract
- Massive Vitreous Loss
- In Cataracts Present in Children

POTENTIAL HAZARDS

The relevant hazards and side effects of Intra Ocular lens essentially the same as (but not limited to) those of routine cataract extraction. These include the following:-

- Intra Ocular Foreign Body Debris
- Temporary Corneal Edema
- Secondary Cataract Formation
- Malpositioned Lens
- Aphakic Glaucoma
- Temporary Flat Anterior Chamber
- Retinal Detachment
- Lens Implants loop Amputation
- Endothelial Corneal Dystrophy
- Vitreous Herniation into the Anterior Chamber
- Infection
- Pupilary Block
- Corneal Dystrophy

'A' CONSTANT

'A' Constant value is to be located on the Box Cover Label. This number is given only as a guideline and not based on any clinical data. It is advised that the surgeon calculates his own value for each style and specifications of Intra Ocular Lens.

In Collaboration with

GANTEC

CORPORATION U.S.A.
2354 Hassell Road, Suite D,
Hoffman Estates, IL 60195

PACKAGE

The Lens is supplied in a Sterile Package. Validity of sterilization applies as long as the sealed inner peel-pouch is not disturbed or damaged. Any damage to the peel-pouch or any accidental opening of the peel-pouch is to be declared as "NOT STERILE".

DIRECTIONS FOR USE

PLEASE examine the peel-pouch prior to opening to assure sterility. Any Damage to the peel-pouch should be viewed very seriously and the IOL inside might not be sterile. After proper examination, peel the pouch and remove the Lens case. Hold and lift the lens using smooth edged forceps. Use BSS (or) Saline to clear off the static charges that might have developed on the lens surface. Do not touch the central optic portion of the lens, with forceps, it will leave permanent forceps marks on the visual axis. Do not use rubber gloves dusted with talc powder. It may cause irritation.

REPORTING

Adverse reactions and potentially sight threatening complications that may reasonably be regarded as lens related and that were not previously expected in nature severity or degree of incidence should be reported to Appasamy Ocular Devices (P) Ltd. This information is requested from the all implanting surgeons in order to document potentially long term effects of IOL implants.

PATIENT IDENTITY CARD

The packaging contains product identification stickers for maintaining a record of the IOL implantation. Surgeons are requested to give the "Patient ID Card" to the patient after implantation and advise them to carry the card at all times.

RETURN GOODS POLICY

Appasamy Ocular Devices (P) Ltd., accepts returned lenses for exchanges only. No cash refunds will be issued. To return lenses, you must first obtain a return authorization number from customer services department. No returned goods will be accepted without proper authorization number. Returned lenses should be shipped by traceable method. No credit will be given to lost or damaged lenses in shipment. Lenses will be replaced as long as they are returned within six months of their original invoice date.

WARNING

Intra Ocular Lens supplied is void of all warranties expressed or implied, if


1. The lens is re sterilized.
2. Lens is altered in any manner.
3. Lens is repackaged by any one other than Appasamy Ocular Devices (P) Ltd.

DISPOSING OF NON STERILE OR CONTAMINATED MEDICAL DEVICES

There are not specific guidelines for disposing of Non sterile or contaminated medical devices. Follow standard procedures for discarding outdated or contaminated products.

PRECAUTIONS

- ⊗ Do not re sterilize the IOLs
- ⊗ Do not Re-use the product
- ⚠ $T < 40^{\circ}\text{C} / 104^{\circ}\text{F}$ Do not store above 45°C (Fragile, avoid shocks)
- ⚠ Use only sterile BSS irrigating solution for rinsing / soaking
- ⊗ Do not use if the pouch is opened
- ⚠ Handle gently to avoid lens surface and haptic damage
- ⚠ Use as soon as the IOL is out of lens case.

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

A. CONSTANT : 118.00

Supra PhobToric

Foldable Hydrophobic Acrylic Intraocular Lens

Supra PhobToric

- Double Muscle Dual Haptic - For better Centration & stability
- Anterior Aspheric Toric - Better Functional vision
- Astigmatism corrected - Sharp Vision

Preloaded System

Here

Open

SN AT1405847 001

Ref SP-TORIC T3 Power20.00 Cyl 1.50 D

Optic 6.00 mm Length 13.00 mm

Exp 04-2019

SUPRAPHOB TORIC

Foldable Hydrophobic Acrylic Intraocular Lens

Supra Phob **Toric**

Preloaded System

Patient Lens Implant
Identification Card

CAUTION Patient Lens Implant
Identification Card
The lens is for identification

PATIENT'S NAME

:

SURGEON'S NAME

:

SURGERY DATE

:

RIGHT

☐

LEFT

☐

Supra Phob Toric

Foldable Hydrophobic Acrylic Intraocular Lens Preloaded System

DESCRIPTION:

Intraocular lenses are optical implants for the human crystalline lens in the visual correction of aphakia. The lens optic is lathe cut from Hydrophobic Acrylic blanks with UV blocking property with more than 90% light transmittance. The lens has a good memory and preloaded in an injecting system for easy delivery. After surgical insertion into the eye, the lens perfectly unfolds to restore the original shape. The dual haptics provide for proper positioning and fixation of the IOL optic within the Capsular bag.

MODE OF ACTION

Supraphob Toric IOLs are intended to be positioned in the posterior chamber of the eye. Replacing the natural crystalline lens and which allows the lens to function as a refractive medium in the correction of aphakia. These IOL's have a biconvex-Anterior Toric optic with cylinder axis marks to denote the flat meridian (plus cylinder axis). Alignment of the toric IOL axis marks with the post-operative steep corneal meridian allows the lens to correct astigmatism.

INDICATIONS

The Supraphob Toric Posterior chamber intraocular lenses are intended for primary implantation in the capsular bag of the eye for visual correction of aphakia and pre-existing corneal astigmatism secondary to removal of a cataractous lens in adult patients with or without presbyopia, who desire improved uncorrected distance vision, reduction of residual refractive cylinder and increased spectacle independence for distance vision.

'A' CONSTANT

'A' constant value is printed on the secondary packing. This number is given only as a guideline and not based on any clinical data. It is advised that the surgeon calculates his own value for each style and specifications of Intra Ocular Lens.

PACKAGE

The lens is supplied in a Pre loaded Blister Sterile Package. Validation of Sterilization applies as long as the sealed inner Blister pack is not disturbed or damaged. Any damage to the Blister pack or any accidental opening of the Blister pack is to be declared as "NOT STERILE".

DIRECTIONS FOR USE

1. Examine the label on the unopened package for model, power (spherical equivalent and cylinder) and expiration date.
2. After opening the secondary package, verify the lens information (model, power, Spherical power, Cylinder power and serial number) is consistent with information on outer package labeling.
3. This device is sterile until the primary blister is obtained. Inspect the blister carefully for tears, cuts, punctures or other signs that the blister has been opened or damaged.
4. Apply the Viscoclastic tip of the cartridge for a smooth delivery of the lens as well back side of the lens. Refer the illustration.
5. Apply a gentle pressure to the injector to deliver the lens comfortably.

REPORTING

Adverse reactions and potentially sight threatening complications that may reasonably be regarded as lens related and that were not previously expected in nature severity or degree of incidence should be reported to Appasamy Ocular Devices (P) Ltd, PONDY. This information is requested from all implanting surgeons in order to document potentially long-term effects of IOL implants.

PATIENT IDENTITY CARD

The packing contains product identification stickers for maintaining a record of the IOL implantation and identification. Surgeons are requested to give the "Patient ID Card" to the patients after implantation and advise them to carry the card at all times.

RETURN GOODS POLICY

Appasamy Ocular Devices (P) Ltd accepts returned lenses for exchanges only. No cash refunds will be issued. To return lens, you must first obtain a return authorization number from customer services department. No returned goods will be accepted without proper authorization number. Returned lenses should be shipped by traceable method. No credit will be given to lost or damaged lenses in shipment. Lenses will be replaced as long as they are returned within six months of their original invoice date.

SELECTION PLACEMENT OF THE Supraphob TORIC IOL

The astigmatism to be corrected should be determined from keratometry and biometry. The size and location of the surgical incision may affect the amount and axis of corneal astigmatism. In order to optimize IOL selection and axis placement. We provide a web-based tool (www.appasamy.com) for the surgeon. Pre-operative keratometry and biometry data, incision location, and the surgeon's estimated surgically induced corneal astigmatism are used to determine the appropriate Supraphob Toric IOL model spherical equivalent lens power, and axis of placement in the eye. For optimum results, the surgeon must ensure the correct placement and orientation of the lens within the capsular bag. The Anterior surface of the IOL is marked with a thin line at the haptic/optic junction that identify the flat meridian of the Supraphob Toric optic. This identification from an imaginary line representing the plus cylinder axis. The Supraphob Toric IOL cylinder axis marks should be aligned with the post-incision steep control meridian (indented axis of placement).

WARNINGS

1. This lens should not be implanted if the posterior capsule is ruptured OR the zonules are damaged OR a primary posterior capsulotomy is planned.

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
A. CONSTANT : 118.00

Supra Phob Rezen

- Diffractive - Refractive Optics
- Advanced 360° Square Edge
- Aberration Controlled
- Enhanced Contrast
- Pre-loaded System

Here

Open


SN RN1305003 033

Ref SPNT200MF Power21.00 Adl 3.50 D
Optic 6.00 mm Length 13.00 mm
Exp 04-2018

SUPRAPHOB REGEN

Patient Lens Implant Identification Card

Supra Phob Regen

- Diffractive - Refractive Optics
- Advanced 360° Square Edge
 - Aberration Controlled
- Enhanced Contrast • Preloaded System

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EC REP

EMERGO EUROPE
Molenstraat 15,
2513 BH The Hague,
The Netherlands.

Marketed by:

 **APPASAMY ASSOCIATES**
Empowering vision®
Chennai, India

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GANTEC CORPORATION U.S.A