DIRECTIONS

for use of frames for sight correction glasses

1. PURPOSE: Metal, plastic, combined, semi-rimless and rimless frames for sight correction glasses produced by the company "ARTS OPTICAL COMPANY LIMITED" are used to fix lenses in a correct position in relation to the eyes.

2. MAIN TECHNICAL CHARACTERISTICS AND STRUCTURE

Frames for sight correction glasses consist of a frontal part, temples, nose pads, hinges. Hinges can be single and spring.

Frames for sight correction glasses are characterized by 3 main geometric parameters: width (a), width of a nose bridge (b) and length of temples up to the curve (L).

The parameters that are used for frames produced by the company "ARTS OPTICAL COMPANY LIMITED" are indicated in the table below:

Frames for sight correction glasses		Distance between vertical tangents and a bezel groove in the rim (a), in mm.	Minimum distance between vertical tangents and the basis of a bezel groove in the rims (b), in mm.	Length of a temple from the axis of the hinge to the beginning of the curve (L), in mm	Weight of a frame without lenses, in grammes
Metal frames for sight correction glasses	men	56±0.2	18±0.2	110±1.0	19.0±1.0
	women	54±0.2	16±0.2	90±1.0	19.0±1.0
	unisex	56±0.2	18±0.2	110±1.0	19.0±1.0
Plastic frames for sight correction glasses	men	56±0.2	18±0.2	110±1.0	20.5±1.0
	women	54±0.2	16±0.2	90±1.0	20.2±1.0
	unisex	56±0.2	18±0.2	110±1.0	20.4±1.0
Combined frames for sight correction glasses	men	56±0.2	18±0.2	110±1.0	19.9±1.0
	women	54±0.2	16±0.2	90±1.0	19.4±1.0
	unisex	56±0.2	18±0.2	110±1.0	19.8±1.0
Semi-rimless frames for sight correction glasses	men	56±0.2	18±0.2	110±1.0	19.0±1.0
	women	54±0.2	16±0.2	90±1.0	19.0±1.0
	unisex	56±0.2	18±0.2	110±1.0	19.0±1.0
Rimless frames for sight correction glasses	men	n/a	n/a	110±1.0	18.6±1.0
	women	n/a	n/a	90±1.0	18.4±1.0
	unisex	n/a	n/a	110±1.0	18.6±1.0

surface of frames should not contain any sharp edges and burrs, megascopic blisters, scratches, bubbles, solder beadings. Temples should turn around the axis smoothly, without jamming.

Decorative furnishing (if applicable) is made of metal (monel KA500).

Decorative covering are not applicable.

Frames for sight correction glasses are in compliance with the norms: ISO EN 12870.

Frames for sight correction glasses are chosen taking into account anthropometric data and information of a patient.

3. CONTRAINDICATIONS FOR USE can be:

- individual intolerance to sight correction when wearing glasses;
- impossibility to correct sight when wearing glasses.

4. MATERIALS USED:

Plastic frames are produced of cellulose acetate Parade 209389. Metal, combined and semi-rimless frames are made of copper-nickel alloy (monel KA500), their temple tips are made of cellulose acetate Parade 209389, their nose pads are made of silicon LSR50 re.

Semi-rimless frames are made of copper-nickel alloy (monel KA500), their temple tips are made of cellulose acetate Parade 209389, their nose pads are made of silicon LSR50 re, fixing nylon thread is made of nylon polymer Nilor.

Non-lenses are made of acrylic Polycryl Clear.

Thread for semi-rimless frames is made of nylon polymer Nilor.

All materials correspond to biological security requirements.

5. MARKING:

Frames for sight correction glasses produced by the company "ARTS OPTICAL COMPANY LIMITED" have marking in temples. There is the name of the producer in the right temple. There are technical parameters in the left temple: horizontal width of lenses to be fixed in (a), width of the nose bridge (b) and full length of temples; code or name of model; colour code or covering code.

6. PREPARATION TO WORK OF THE ITEM:

When being inserted into it, an optical frame must not lose its initial shape. The rims have to be symmetrical. Splits, scratches, gaps between optical lenses and rims of the frontal part of the frame, as well as other defects of assembled glasses are not acceptable.

Optical lenses inserted into a frame must not move or fall down from the rims of the frame during normal exploitation of glasses.

7. STORAGE AND USAGE RULES:

- it is not allowed to make significant mechanical efforts to a frame;
- the frame is to be kept in its case or in another package that excludes possibility of mechanical damage;
- whole average service life period of a frame is not less than 3 years;
- storage conditions: in heated accommodations, at ambient temperature 5°C to $+50^{\circ}\text{C}$, at relative humidity of 80%.
- when storing, it is not allowed to use aggressive substances (acids, alkalis) or organic solvents or their vapours.

8. USAGE CONDITIONS

Frames for sight correction glasses are used in all climatic zones, at any time of day and night.

9. USAGE PRECAUTIONS

Use optical frames with care when going in for sport and when doing quick movements with the head, leans, jumps, and strikes.

10. EVENTUAL SIDE EFFECTS WHEN USING

No side effects have been revealed when using the medical item.

11. INFORMATION ABOUT EVENTUAL INFLUENCE OF USAGE OF THE MEDICAL ITEM ON ABILITY TO DRIVE A VEHICLE AND ON CONTROL MECHANISMS

Frames for sight correction glasses ensure possibility to drive a vehicle and control mechanisms for drivers having problems with sight (if frames are used with production of sight correction glasses).

12. STERILITY

Not sterile.

13. MAINTENANCE AND REPAIR

Frames for sight correction glasses are to be repaired in special laboratories having licences for production, maintenance and repair.

14. WARRANTY

Shelf life of frames: for plastic and combined frames – 5 years; for metal frames – 10 years.


Whole service life period of frames is not less than 3 years in case of no more than 30000 movements of temples.

15. TRANSPORTATION

Transportation of frames for sight correction glasses can be made by any means of transport, in covered vehicles, according to freight rules applicable for such a vehicle.

16. UTILIZATION

Frames for sight correction glasses are to be utilized in the same way as domestic waste.

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