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Head of Global Regulatory Affairs
Mark Smyth *ysmyth*

Hear now. And always

Cochlear®

Contents

Get started	5
How to get Custom Sound® 4.2 Service Pack 3	5
Symbols used in this guide	6
Introduction	6
System requirements	8
Log in to Custom Sound® software	9
Create a clinician	10
Edit a clinician	11
Set clinician preferences	11
Navigate Custom Sound® software.....	14
Menu bar	14
Toolbars	16
Status bar	18
Keyboard shortcuts	21
Work with recipients	24
Create a recipient	24
Edit a recipient	26
Export a recipient	26
Import a recipient	27
Import recipient data	27
Archive or restore a recipient	28
Connect a sound processor	30
Connect a remote assistant	31
Start a session	31
Program recipients.....	32
Obtain impedance measurements.....	33
Measure impedances	33
View impedance measurements	34
Flag or unflag electrodes	35
Obtain NRT® measurements	36
Perform AutoNRT®	36
Work with MAPs.....	38
Create or open a MAP	38
Upgrade or convert a MAP	39
Compare MAPs	40
Archive or restore a MAP	40
Set levels.....	42
Work with the Set Levels screen	42
Select channels	43
Configure live program settings	43
Talk to a recipient	44
Work with Implant ID	44
Use HearingMentor	45
Set T- and C-levels	45
Select NRT®/objective levels	46
Use the behavioural method	47

Use the NRT®/objective offset method.....	47
Use the NRT®/objective preset method.....	48
Measure individual channels.....	49
Sweep levels.....	49
Shift levels.....	50
Tilt levels.....	50
Adjust MAP parameters.....	51
Adjust basic MAP parameters.....	51
Adjust advanced MAP parameters.....	52
Modify channel gains.....	53
Reorder channels.....	54
Double channels.....	55
Add notes and generate reports.....	55
Create additional MAPs.....	56
Predict T- and C-levels.....	56
Create progressive MAPs.....	56
Program an acoustic component.....	57
Set acoustic parameters.....	57
Determine power requirements.....	59
Estimate compliance levels.....	59
Check battery suitability.....	60
Optimise the power level.....	60
Apply a manual power level.....	63
Determine compatibility of MAPs.....	63
Perform bilateral balance.....	64
Balance T- and C-levels.....	64
Configure and write programs.....	65
Configure a program.....	66
Select a previous environment.....	75
Write to a sound processor.....	75
Finalise a session.....	76
View session history.....	77
View usage data.....	78
View CP900 usage data.....	78
Work with databases.....	79
Create or edit a connection.....	79
Connect to a server.....	80
Create a database.....	80
Upgrade a database.....	81
Back up or restore a database.....	81
Delete a database.....	82
Appendix 1: Cochlear™ Nucleus® CR120/CR220 Intraoperative Remote Assistant.....	83
Appendix 2: Remote Assistant Fitting.....	84
Appendix 3: Working with Nucleus® Fitting Software.....	90
Other information.....	93
Certification, applied standards and copyright information.....	93
Table of symbols.....	93
Legal statement.....	93

Get started

How to get Custom Sound® 4.2 Service Pack 3

Cochlear invites you to register to download **Custom Sound Suite 4.2 Service Pack**

3. To download, you need to register by going to the following location:

<https://software.cochlear.com/MZ0gMW>

Custom Sound 4.2 Service Pack 3 is also available in CD format if preferred.

Symbols used in this guide

 **Note:** Important information or advice. Can save inconvenience.

Introduction

Custom Sound® 4.2 software is intended to be used by cochlear implant and hearing care professionals to create and modify hearing profiles in order to provide comfortable and usable stimulation to cochlear implant recipients.

The cochlear implant system components compatible with Custom Sound 4.2 Software are listed below:

Cochlear Nucleus® Implants	Cochlear Nucleus® Sound Processors									
	CP920 sound processor	CP910 sound processor	CP810 sound processor	CP802 sound processor	Freedom® Hybrid™	Freedom®	ESPrnt™ 3G	ESPrnt™	SPrint™	Spectra
Nucleus® CI551 cochlear implant			☑	☑		☑				
Nucleus® CI532 cochlear implant	☑	☑	☑	☑	☑	☑				
Nucleus® CI522 cochlear implant	☑	☑	☑	☑	☑	☑				
Nucleus® CI513 cochlear implant	✓	✓	✓	☑	☑	✓				
Nucleus® CI512 cochlear implant	✓	✓	✓	☑	☑	✓				
Nucleus® ABI541 auditory brainstem implant			☑			☑				
Nucleus® CI422 cochlear implant	☑	☑	☑	☑	☑	☑				
Nucleus® Hybrid™ L24 cochlear implant	☑	☑	☑		☑	☑				
Nucleus® Hybrid™ S12 cochlear implant			☑		☑	☑				
Nucleus® Hybrid™ S8 cochlear implant			☑		☑	☑				
Nucleus® Freedom® CI24RE Contour Advance™ cochlear implant	✓	✓	✓	☑	☑	✓	✓			

	CP920 sound processor	CP910 sound processor	CP810 sound processor	CP802 sound processor	Freedom® Hybrid™	Freedom®	ESPr™ 3G	ESPr™	SPrint™	Spectra
Nucleus® Freedom® CI24RE Straight cochlear implant	✓	✓	✓	☑	☑	✓	✓			
Nucleus® Freedom® CI24RE Contour® cochlear implant	✓	✓	✓	☑	☑	✓	✓			
Nucleus® 24 ABI	✓	✓				✓			✓	
Nucleus® 24 Double Array	✓	✓	✓	☑		✓	✓	✓	✓	
Nucleus® CI24R (CA) cochlear implant	✓	✓	✓	☑	☑	✓	✓	✓	✓	
Nucleus® CI24R (CS) cochlear implant	✓	✓	✓	☑	☑					
Nucleus® CI24R (ST) cochlear implant	✓	✓	✓	☑	☑	✓				
Nucleus® CI24M cochlear implant	✓	✓	✓	☑	☑	✓				
Nucleus® 22 cochlear implant						✓	✓*	✓*		✓
Cochlear Nucleus® Remote Assistants	Cochlear Nucleus® Sound Processors									
Nucleus® CR110 remote assistant			✓							
Nucleus® CR230 remote assistant	✓	✓	✓							

✓ Available

☑ Available as a regional option only

* ESPr™ 3G for Nucleus® 22 or ESPr™ 22

Note:

- CI551, CI522, ABI541, CI422 and Nucleus® Hybrid™ implant types are available as regional options only.
- CP802 and Freedom® Hybrid™ sound processors, and the CP810 sound processor for Hybrid™ and ABI541, are available as regional options only.
- Not all products are approved in all markets.
- The CR210 remote control is not compatible with Custom Sound software.

System requirements

The minimum computer system requirements are listed below:

Operating system	Windows® 8 (or 8.1), Windows® 7, Windows Vista® (Service Pack 1 or higher), Windows® XP (Service Pack 3)
Processor	Intel® Pentium® 4 2.4 GHz
RAM	1 GB or greater
Hard drive	7 GB or greater with at least 5 GB of available space
CD-ROM	Required for software installation only
Screen resolution	1024 x 768
Serial port*	1 x port supporting 115 k baud (PPS)
USB port*	1 x port supporting USB 1.0

* Bilateral programming requires the simultaneous connection of two separate programming systems compatible with the sound processor types being used. Depending upon your programming equipment configuration, you may require 2 x PCI bus slots, 2 x serial ports or 2 x USB 1.0 ports. For more information please contact Cochlear.

Any IBM compatible computer system connected to the Cochlear™ system, when used in programming mode, should meet FCC 47CFR Part 15 Class B and CISPR22 Class B emissions requirements. Consult the documentation provided by the computer system manufacturer to ensure compliance, prior to connecting any Cochlear device.

Log in to Custom Sound® software

The Login window displays when Custom Sound® software is launched, and allows you to log in to either the application or training mode. Training mode simulates connection to a sound processor and allows you to learn the features of Custom Sound without requiring physical hardware.

To log in to Custom Sound software:

1. Select your clinician name from the **Clinician** drop-down list.
2. Type your password in the **Password** field.
A password is only required if one has been set up in the Preferences window.
3. Click **Login**.
The Recipient screen displays.

To log in to training mode:

1. Select your clinician name from the **Clinician** drop-down list.
2. Type your password in the **Password** field (if required).
3. Select the **Training Mode** check box.
4. Select the types of sound processors you wish to train on from the **Processor** drop-down lists.
Up to two CP900 series, CP800 series, Freedom® Hybrid™ and/or Freedom® sound processors can be selected.
5. Click **Login**.
The Recipient screen displays, and the title bar indicates the software is operating in training mode. The status bar indicates the software is using simulated hardware, and if there is a programming interface attached Custom Sound software does not connect to it.

Note:

- 'Training Mode' is automatically added to the title of MAPs created in training mode.
- When not in training mode, MAPs created in training mode can be opened and written to the recipient's sound processor.

Create a clinician

The username entered when the Custom Sound® software database is installed is automatically created as a clinician with clinic administrator rights. Users with clinic administrator rights can create additional clinicians and edit preferences for all clinicians as required. Creating individual users allows each clinician to set their own preferences. The recipient database is shared across all users.

To create a clinician:

1. Click the **Tools** menu and select **My Clinic**.
The Clinic window displays a list of registered clinicians.
2. Click **Create New Clinician**.
The Preferences window displays with the **Personal** tab selected.
3. Type the clinician's name in the **Username** field.
4. Type the clinician's email address in the **Email** field (optional).
The email address automatically populates the Feedback form when providing feedback about Custom Sound software to Cochlear™.
5. To give the clinician administrator rights, select the **Clinic Administrator** check box.
6. To set a password for the clinician, click **Set Password**, type the password details and click **OK**.
7. To set additional preferences for the clinician, click each of the tabs in the Preferences window and select the desired settings.
Clinicians can set their own preferences at any time.
8. Click **OK** to save the changes and close the Preferences window.
The clinician displays in the Registered Clinicians list.
9. Create additional clinicians as desired, or click the **Close** button  to close the Clinic window.

 **Note:** If Quick/Standard Programming has been installed as a regional option, **Quick Programming** or **Standard Programming** can be selected in the **Personal** tab of the Preferences window. Standard programming provides access to all clinical features in Custom Sound software, while quick programming provides access to the most used clinical features only.

Edit a clinician

Existing clinicians can be edited as desired by clinic administrators only.

To edit a clinician:

1. Click the **Tools** menu and select **My Clinic**.
The Clinic window displays a list of registered clinicians.
2. Select the desired clinician and click **Edit Selected Clinician**.
The Preferences window displays with the **Personal** tab selected.
3. Edit the existing clinician details as desired.
Clinicians can be marked as inactive by selecting the **This Account is Disabled** check box, and a new password can be entered should a clinician forget their password.
4. To set additional preferences for the clinician, click each of the tabs in the Preferences window and select the desired settings.
5. Click **OK** to save the changes and close the Preferences window.
6. Click the **Close** button  to close the Clinic window.

To edit the clinic name, click **Edit Clinic Details** in the Clinic window and overwrite the existing details.

 **Note:** If Quick/Standard Programming has been installed as a regional option, the default user interface for the clinic can be changed. Click **Edit Clinic Details** in the Clinic window and select **Quick Programming** or **Standard Programming**. The selected option applies by default when creating a clinician, and can be edited in the Preferences window for each clinician as required.

Set clinician preferences

Each clinician can set their own default preferences that apply each time Custom Sound® software is launched.

To set clinician preferences:

1. Click the **Tools** menu and select **My Preferences**.
The Preferences window displays with the **Personal** tab selected.
2. Click each tab in turn and select the desired settings.
3. Click **OK** to save the changes and close the Preferences window.

 **Note:** The **OK** button will be available only if changes have been made to the clinician preferences.

The following settings can be selected from the tabs in the Preferences window:

The following settings can be modified:

Personal		
	Username and Password	The username and password (optional) are used to log in to Custom Sound software and allow individual preferences to be set for each clinician.
	User Interface	If Quick/Standard Programming has been installed as a regional option, quick or standard programming can be selected. The user interface can be set by clinic administrators only and cannot be changed during a session. If Quick Programming is selected, limited parameters can be set in the Preferences window.
General		
	Levels and Stimulations	The specified parameters apply by default when setting T- and C-levels using any programming method, or when performing a sweep of the channels.
	Global Modifier	The specified parameters apply by default when shifting or tilting T- and C-level profiles.
	Miscellaneous	Use spacebar to stop: all current stimulations can be stopped by pressing the spacebar on the keyboard.
		Show dialogue when importing existing data: a prompt displays when importing existing recipients, enabling personal data to be excluded from the import.
		Delay auto compliance until going live: voltage compliance is automatically estimated when a MAP goes live, and is not tested when a MAP is opened or changed.
		Auto flag faulty electrodes: short- and open-circuit electrodes are automatically flagged when performing an implant test (postoperative mode only).
Progressive MAPs		
	Behavioural	When creating progressive MAPs, the C-level profile on each additional MAP is adjusted by the specified amount.
	NRT®/Objective Preset	When using the NRT/objective preset method, the T- and C-level profiles on each additional MAP are offset from the objective measurement profile by the specified amounts.
Display		
	Data Grid Advanced Properties	The selected options display in the Data Grid on the Set Levels screen.

Recipient Information Display	The selected options display in the Recipient list on the Recipient screen.
MAP Display	Show parent MAP levels: when creating a new MAP based on a saved MAP, markers display in the Channel Grid indicating the original T- and C-levels. Show levels as a tooltip hint: when the cursor is placed over a channel in the Channel Grid, the T- and C-level values display as a tooltip. Show levels on channel marker labels: the T- and C-level values of each selected channel display in the Channel Grid.
Miscellaneous	Show implant communication messages: when Custom Sound software detects an error in communicating with the implant, a warning message displays on the status bar.
Streamlined	
Default Channels	Behavioural channels: the selected channels display as measurable when using the Behavioural method. T-offset channel: the selected channel displays as measurable when using the NRT/objective offset method.
NRT	
General	Number of electrodes: the specified number of electrodes is selected by default when performing AutoNRT®. Include NRT analysis data in CDX files: NRT measurements are included in recipient files when exporting recipients.
Display	Show traces during AutoNRT: the NRT Traces window displays automatically when performing AutoNRT. Show NRT/objective markers: NRT/objective level markers display in the Channel Grid.
Intraoperative Options	Perform electrode conditioning: electrodes are conditioned at the specified current level prior to measuring the NRT level. Perform enhanced interpolation: automatic interpolation of NRT levels is enabled.
Postoperative Options	Starting current level: each channel starts measuring the NRT level from the specified current level. Auto increase waiting time (secs): the AutoNRT algorithm waits the specified amount of time before increasing the stimulus required to perform the next measurement.
Reports	
Reports	Reports can be viewed in the Report viewer or the web browser. When the Report viewer is used, notes added to the report have text wrapping applied based on the number of characters specified.

The Preferences window can also be accessed from the Programming screen by clicking the **My Preferences** button  on the Programming toolbar.

Navigate Custom Sound® software

The following tools can be used to navigate Custom Sound® software:

- [The menu bar at the top of the screen](#)
- [The Programming and Set Levels toolbars](#)
- [The status bar at the bottom of the screen](#)
- [Keyboard shortcuts available throughout the application](#)

Menu bar

The following actions can be selected from the menu bar:

Menu	Action
File	
Import AutoNRT® and impedance measurements	Import AutoNRT and impedance measurements that have been exported from a remote assistant.
Exit	Exit the Custom Sound® software.
Session	
Intraoperative Testing	Start an intraoperative testing session for the selected recipient.
Program a Recipient	Start a programming session for the selected recipient.
Session History	Access the Session History window to view details of previous sessions.
End Session	End the current session.
Recipient	
Import Recipients	Import recipients from a file.
Create	Access the Recipient Details window to create a new recipient.
Edit	Access the Recipient Details window to edit an existing recipient.
Archive	Archive the selected recipient(s).
Talk to Recipient	Talk to a recipient during a programming session using a previously saved MAP.
Export	Export the selected recipient(s) to a file.
Export Anonymously	Export the selected recipient(s) to a file, excluding personal details.
Manage	Access the Manage Recipients window to archive or restore a recipient.
Implant	
Auto Power Setting	Access the Auto Power Setting window to perform the measurements required to optimise the power level.
Impedance Measurement History	Access the Impedance Measurements window to view impedance measurements from previous implant tests.

MAP		
	Save	Save the currently selected MAP.
	Save All	Save all open MAPs.
	Close	Close the currently selected MAP.
	Close All	Close all open MAPs.
	Convert MAP	Convert an existing MAP for use with a new sound processor.
	Upgrade MAP	Upgrade an existing MAP for use with a Freedom® sound processor. Note: MAP upgrades are not available between Freedom® Hybrid™ sound processors and CP900 sound processors.
	Compare MAPs/Compare with Another MAP	Access the MAP Comparison window to view differences between MAPs.
	Manage	Access the Manage MAPs window to archive or restore existing MAPs.
	Link Doubled Channels	Link the T- and C-levels on adjacent channels that have the same active and indifferent electrodes.
Processor		
	Connect	Connect Custom Sound software to the sound processor attached to the Programming Pod, Portable Programming System (PPS) or Clinical Programming System (CPS).
	Reset	Reset the connected sound processor and clear all programming slots.
	SPrint™ Master Reset	Upgrade the SPrint supervisor software and clear all programming slots.
Tools		
	My Preferences	Access the Preferences window to set clinician preferences.
	My Clinic	Access the Preferences window to create or edit clinicians and set clinician preferences (clinic administrators only).
	Disable/Enable Pod Communication	Disable/enable communication between Custom Sound software and the attached Programming Pod.
Help		
	Help	Access the Help file to view information on how to use Custom Sound software.
	Provide Feedback	Access the Feedback form to provide feedback about Custom Sound software to Cochlear™.
	About	Access the About Custom Sound window to view information about the current version of Custom Sound software.

Toolbars

The Programming toolbar is available on the Programming screen, regardless of the task being performed. The Set Levels toolbar is available on the Set Levels screen only.

	Tooltip	Action
Programming toolbar		
	Edit the Selected Recipient	Access the Recipient Details window to edit an existing recipient.
	Previous Task	Return to the previous task in the Programming Tasks or IntraOp Tasks navigation bar.
	Next Task	Progress to the next task in the Programming Tasks or IntraOp Tasks navigation bar.
	Connect to an Attached Processor	Connect Custom Sound® software to the sound processor attached to the Programming Pod, Portable Programming System (PPS) or Clinical Programming System (CPS).
	My Preferences	Access the Preferences window to set clinician preferences.
	Talk to Recipient	Talk to a recipient during a programming session using a previously saved MAP.
	Stop All Stimulations	Stop all current stimulations.
	End Programming Session	End the current session.
Set Levels toolbar		
	Save MAP	Save the currently selected MAP.
	MAP Parameters	Access the MAP Parameters window to set basic and advanced MAP parameters.
	Acoustics	Access the Acoustics window to prescribe acoustic levels (Freedom® Hybrid™ and CP900 series sound processors only).
	Undo	Undo the previous action, up to a maximum of 20 actions. Previous actions can be selected from the drop-down list.
	Redo	Redo the previous action, up to a maximum of 20 actions. Previous actions can be selected from the drop-down list.
	Go Live	Measure T- and C-levels using live voice testing.
	Hug T-Profile	Set C-levels to follow the T-level profile with a dynamic range of 2.

	Predict Levels	Predict T- and C-levels when creating a new MAP from an existing MAP using a different stimulation rate or pulse width.
	Make All Channels Measurable	Enable T- and C-levels to be measured on all channels.
	Make Only Selected Channels Measurable	Enable T- and C-levels to be measured on the selected channels only. The levels on the non-selected channels are automatically interpolated.
	Reverse Electrode Order	Reverse the tonotopic order of channels by changing the channel-to-electrode configuration for all channels (ABI541 and Nucleus® 24 ABI implants only).
	Title/Notes	Access the Title and Notes window to enter a title and/or notes.
	Modify Gains	Access the Gain Shapers window to adjust high and/or low frequency gains.
	Configure Live Program Settings	Access the Live Program Settings window to adjust the volume, sensitivity and environment settings that apply when going live.
	Create Progressive MAPs	Create additional MAPs from an initial MAP, with the C-levels adjusted on each progressive MAP.
	MAP Report	Access the Report viewer to print the MAP details (saved MAPs only).
	Estimate Compliance Levels	Estimate the maximum C-level achievable for each channel in a MAP.
	Measure Skin Flap and Optimise Power Level	Estimate the skin flap thickness and optimise the power level for the MAP (Freedom® sound processors with Nucleus® 22 implants only).
	Close MAP	Close the currently selected MAP.

Status bar

The status bar provides information on the current status of the recipient, the connected sound processor and programming interface, and the implant. The status bar is colour coded to differentiate between right ear (red) and left ear (blue) implants and sound processors.

The following icons may display on the status bar:

	Tooltip or message	Status
Auto stimulate status		
	Auto stimulate	Automatic stimulation is currently enabled when setting T- and C-levels.
Programming interface status		
	(Programming interface, com #)	The specified programming interface is connected.
	(Remote assistant type)	The specified remote assistant is connected. The right-click menu allows you to: • upgrade the software for the remote assistant • export the recipient's AutoNRT® and impedance data.
	Updating device software...	The remote assistant software is being updated.
Sound processor status		
	(Sound processor type, serial number)	The specified sound processor is connected. The right-click menu allows you to connect or reset a sound processor. Note: the sound processor icon will vary according to the sound processor that is connected.
	(Sound processor type, serial number)	The specified Freedom® sound processor is using a Nucleus® 22 coil, and is configured for Nucleus 22 use only.
	No processor connected	No sound processor is connected.
	Writing...	Programs are being written to the connected sound processor.

		Compliance levels are being estimated for the connected sound processor.
	Getting battery suitability...	Battery suitability is being checked for the connected sound processor.
	Implant ID enabled for MAPs on this processor	Implant ID is enabled for the connected sound processor.
Program status		
	(MAP title, environment, volume and sensitivity settings)	The programming slot on the connected sound processor contains the specified program. The right-click menu allows you to open the MAP, go live with the program, or erase the program from the programming slot. A red dot displayed on the icon indicates the presence of MAPs that have been modified using Remote Assistant Fitting.
		The programming slot on the connected sound processor is empty.
	Corrupt MAP in slot	The MAP in the programming slot is corrupt.
		The MAP in the programming slot is not in the database.
		The specified program is currently in live mode.
Implant status (excluding Nucleus® 22 implants)		
	Implant detected	Communication with the implant has been established.
	No communication with implant	Communication with the implant cannot be established. This may indicate the coil has come off the recipient's implant or cables are disconnected or faulty.
	Incorrect implant type detected	The selected sound processor belongs to a different implant type.
	Telemetry being affected by other processors	Telemetry is being affected by another sound processor in too close proximity. For bilateral recipients, the sound processor for the contralateral ear may need to be switched off.
	Could not read implant ID. Please remove coil from implant, wait 5 seconds then reconnect coil.	The implant ID could not be read during the implant ID check.

Recipient status		
	Recipient (<i>name surname</i>) is not currently open	The recipient of the connected sound processor is in the database. Click to start a programming session.
	Recipient of this processor is unknown. Click to register	The recipient of the connected sound processor is not in the database. The Unknown Recipient prompt displays automatically and allows you to register the recipient or reset the sound processor. If a session is currently open, you can also associate the programs on the sound processor with the open recipient.
	Recipient (<i>name surname</i>)	The recipient of the connected sound processor is in the database, and a programming session is open.
	Recipient (<i>name surname</i>) is not currently open	The recipient of the connected sound processor differs to the currently open programming session. Click to start a programming session for the connected recipient.

Keyboard shortcuts

The following keyboard shortcuts are available in Custom Sound® software:

Shortcut	Action
Stop stimulation	
Esc	Stop stimulation.
Spacebar	Stop stimulation (where clinician preference is selected).
General	
K	Shift the keyboard focus to the Data Grid and the Channel Grid.
F1	Access the Help file.
F3	Connect/disconnect a sound processor.
F4	Enable talk over mode.
F5	Go live.
F6	Enable/disable Programming Pod communication.
F9	Access the Preferences window.
Ctrl+S	Save the currently selected MAP in the Set Levels screen.
Shift+Ctrl+S	Save all open MAPs in the Set Levels screen.
Ctrl+F4	Close the currently selected MAP in the Set Levels screen.
Shift+Ctrl+F4	Close all open MAPs in the Set Levels screen.
Ctrl+P	Print a report in the Report viewer.
Ctrl+Z	Undo.
Ctrl+Y	Redo.
Select channels in the Channel Grid	
←	Move the channel focus to the next left measurable channel.
→	Move the channel focus to the next right measurable channel.
Ctrl+←	Move the channel focus to the next left interpolated channel.
Ctrl+→	Move the channel focus to the next right interpolated channel.
V	Toggle between select and deselect a channel.
X	Toggle between activate and deactivate a channel.
I	Toggle between measure and interpolate a channel.
Ctrl+D	Toggle between double and unlink a channel.

Select mode for measuring T-, C- and objective measurement levels

T	Switch to T mode (all keyboard actions will be performed on T-levels).
C	Switch to C mode (all keyboard actions will be performed on C-levels).
O	Switch to O mode (all keyboard actions will be performed on objective measurement levels).
B	Switch to T & C mode (all keyboard actions will be performed on both T- and C-levels).
M	Toggle between T, C and T & C modes.

Measure T-, C- and objective measurement levels (depending on selected mode)

<i>n</i> Enter	Set the T- and/or C-level to <i>n</i> .
U <i>n</i> Enter	Set the up step size to <i>n</i> .
D <i>n</i> Enter	Set the down step size to <i>n</i> .
↑	Increase the T- and/or C-level or objective measurement level by the up step size.
↓	Decrease the T- and/or C-level or objective measurement level by the down step size.
Page Up <i>n</i>	Increase the T- and/or C-level by <i>n</i> .
Page Down <i>n</i>	Decrease the T- and/or C-level by <i>n</i> .
Z	Toggle between automatic stimulation on and off.
S <i>n</i> Enter	Set the number of stimuli to <i>n</i> .
F2	Produce one stimulus at T- or C-level on the focused channel.
G <i>n</i>	Produce <i>n</i> stimuli at T- or C-level on the focused channel.
J	Produce the set number of stimuli at T- or C-level on the focused channel.

Sweep levels (depending on selected mode)

R	Sweep right (ascending frequency) for the set number of channels.
Ctrl+R <i>n</i> Enter	Sweep right (ascending frequency) for <i>n</i> channels.
L	Sweep left (descending frequency) for the set number of channels.
Ctrl+L <i>n</i> Enter	Sweep left (descending frequency) for <i>n</i> channels.
N <i>n</i> Enter	Set the number of channels to be stimulated to <i>n</i> .
N S Enter	Set the number of channels to be stimulated to selected.
N A Enter	Set the number of channels to be stimulated to all.
P <i>n</i> Enter	Set the percentage of the dynamic range to <i>n</i> .

Shift levels (depending on selected mode)

Ctrl+↑	Increase global T- and/or C-levels.
Ctrl+↓	Decrease global T- and/or C-levels.
↑	Increase global C-levels in live.
↓	Decrease global C-levels in live.
H P Enter	Set the shift mode to percentage of the dynamic range.
H C Enter	Set the shift mode to current level.
A A Enter	Set the channel selection to all.
A S Enter	Set the channel selection to selected.
Q n Enter	Set the shift step size to <i>n</i> .

Tilt levels (depending on selected mode)

[	Left tilt global T- and/or C-levels.
]	Right tilt global T- and/or C-levels.
F n Enter	Set the tilt value to <i>n</i> .

Use the Data Grid

Shift+↑ / Shift+↓	Move up/down one row in the Data Grid.
Shift+← / Shift+→	Move left/right one column in the Data Grid.
↑ / ↓	Select a value from the drop-down list in the selected cell (Active Electrode, Stimulation Mode/Indifferent Electrode, Gain and Pulse Width fields only).

Work with recipients

The Recipient screen displays a list of all current recipients contained in the Custom Sound® software database.

The Recipient screen displays automatically when you log in to Custom Sound software.

Recipients can be sorted by clicking the header of the column you wish to sort by, or filtered using the **Search** field. Valid search criterion include first name, last name, date of birth, national ID number and cross reference ID number.

The Recipient screen allows you to:

- [Create a new recipient](#)
- [Edit an existing recipient](#)
- [Export a recipient](#)
- [Import a recipient](#)
- [Archive or restore a recipient](#)

Note: The default columns displayed in the Recipient list can be changed in the **Display** tab of the Preferences window (F9).

Create a recipient

In order to perform intraoperative testing or program a recipient's sound processor, the recipient must first be created in Custom Sound® software.

To create a recipient:

1. Click **Create** in the Recipient screen.

Alternatively, click the **Recipient** menu and select **Create**.

The Recipient Details window displays with the **General** tab selected. Mandatory fields are indicated by a **Required** symbol .

2. Type the recipient's name in the **First Name** and **Last Name** fields, and enter additional personal details as desired.

To enter a date of birth, select the check box in the **Date of Birth** field and either type the date or select it from the drop-down calendar.

3. Click **Add**.

The Implant Details window displays.

4. Type the implant name in the **Implant Type** field, or select the implant from the drop-down list.

Alternatively, click the **Search for an Implant** button . The Implant Search window displays a list of available implants, and can be filtered by clearing the check boxes corresponding to the implant type, implant series or sound processor you wish to exclude. To view an image of each implant type, select the **Show Implant Images** check box. Select the implant in the Filtered Results list, and click **OK**.

Note: Ensure the correct implant is selected. Incorrect implant selection may result in excessive or unsafe charge being delivered to the recipient.

5. Select the implant side from the **Ear** drop-down list.

If more than one implant is entered per side, the first implant will be marked as explanted.

6. Enter additional implant details as desired, and click **OK**.

7. Click the **Address** tab and enter the recipient's address details as desired.

8. Click the **Sound Processors** tab and click **Add**.

The Sound Processor Details window displays.

9. Select the type of sound processor from the **Model** drop-down list.

10. Enter additional sound processor details as desired, and click **OK**.

11. Click the **Acoustics** tab, and click the points on the left and/or right audiogram that correspond to the recipient's pure tone audiometric data.

The right-click menu allows you to delete a point, clear all points, or mark a point as no response.

12. Click the **Notes** tab and type additional information about the recipient as desired.

13. To print a report of the recipient details, click **Report** and use the Print toolbar on the Report viewer.

14. Click **OK** to save the changes and close the Recipient Details window.

The recipient is added to the list on the Recipient screen, and the Programming screen displays in postoperative mode.

 **Note:** For Freedom® Hybrid™ sound processors and CP900 series sound processors in Hybrid mode, the audiogram must be entered in the **Acoustics** tab to enable acoustic levels to be prescribed.

Edit a recipient

Existing recipients can be edited as required.

To edit a recipient:

1. Select the desired recipient in the Recipient list.
2. Click **Edit**.
Alternatively, click the **Recipient** menu and select **Edit**.
The Recipient Details window displays.
3. Edit the existing recipient details as desired.
4. Click **OK** to save the changes and close the Recipient Details window.

The Recipient Details window can also be accessed from the Programming screen by clicking the **Edit the Selected Recipient** button  on the Programming toolbar.

Export a recipient

A recipient's details can be exported to a file, which can then be imported into another Custom Sound® software database. The Anonymous Export function removes the recipient's personal details from the exported file.

To export a recipient:

1. Select the desired recipient in the Recipient list.
2. Right-click on the recipient and select **Export** or **Export Anonymously**.
Alternatively, click the **Recipient** menu and select **Export** or **Export Anonymously**.
The Save As window displays.
3. Navigate to the desired file location and enter an appropriate filename.
4. Click **Save**.

The recipient details are saved to the specified location as a Cochlear Data Exchange (.cdx4) file.

Note: By default, NRT® data is not exported with the recipient's details. To include NRT data, access **Tools > My Preferences**, click the **NRT** tab and select the **Include NRT analysis data in CDX files** check box.

Import a recipient

Recipient files that have been exported from another Custom Sound® software database can be imported into the Recipient list.

To import a recipient:

1. Click the **Recipients** menu and select **Import Recipients**.

The Open window displays.

2. Navigate to the .cdx, .cdx2, .cdx3, .cdx3a or .cdx4 file you wish to import.

3. Click **Open**.

The file is imported into Custom Sound software. A new recipient is created, or, if the recipient already exists, the current details are overwritten.

Note: To view a prompt that allows you to include or exclude recipient details when importing a file, access **Tools > My Preferences**, click the **General** tab and select the **Show dialogue when importing existing data** check box.

Import recipient data

Recipient AutoNRT® and impedance data can be imported into Custom Sound® software. Prior to importing the data, it must have been:

- exported from a remote assistant
- saved to a location it can be imported from.

Note: this option is not available from a CR110 or CR230 remote assistant.

To export AutoNRT and impedance data from a remote assistant:

1. [Connect the recipient's remote assistant](#).
2. Right-click on the remote assistant icon in the [status bar](#).
3. Select **Export AutoNRT and Impedance measurements**.

The Export Remote Assistant AutoNRT and impedance measurements to... window will display to allow you to select the file location.

4. Navigate to the location you wish to save the .crf file.
5. Click **Save**.

The .crf file will be saved to the selected location.

To import recipient AutoNRT and impedance data:

1. Click the **File** menu and select **Import AutoNRT and impedance measurements**.

The Open Cochlear Remote File window will display to allow the selection of the data file exported from the remote assistant.

2. Navigate to the .crf file you wish to import.
3. Click **Open**.

The file is imported into Custom Sound software.

The data will be imported and stored in the Custom Sound software database. Data that corresponds to an implant known to Custom Sound software will be associated with the relevant implant. If the implant does not exist, the data will be saved and matched to an implant either when a new implant with the correct Implant ID is created, or an existing implant has its Implant ID changed to match the Implant ID in the imported data.

After the data import process has completed, Custom Sound software will display the results depending on the success of the import:

- If no new valid data has been saved to the database, Custom Sound software will display a message in the status bar to indicate the reason why no changes have been made, or
- When data has been successfully associated with implants, a window will display the recipients, implants, and any additional information relating to the imported data.

Archive or restore a recipient

Recipients who are no longer active can be archived from the Recipient list. An archived recipient is not deleted from the database, and can be restored at a later date if required.

To archive a recipient:

1. Select the desired recipient in the Recipient list.
 2. Right-click on the recipient and select **Archive**.
- Alternatively, click the **Recipient** menu and select **Archive**.
- The recipient is removed from the Recipient list.

To restore an archived recipient:

1. Click **Manage** in the Recipient screen.

Alternatively, click the **Recipient** menu and select **Manage**.

The Manage Recipients window displays.

2. Select the recipient you wish to restore in the Archived list.

Recipients can be sorted by clicking the header of the column you wish to sort by, or filtered using the **Search** field.

3. Click **Restore**.

The recipient is removed from the Archived list to the Active list, and is restored to the Recipient screen.

4. Click the **Close** button  to close the Manage Recipients window.

Recipients can also be archived in the Manage Recipients window by selecting the recipient in the Active list and clicking **Archive**.

Connect a sound processor

A sound processor is connected to a computer via the appropriate programming interface:

- **Programming Pod:** for use with CP900 series, CP800 series, Freedom® Hybrid™ and Freedom® sound processors.
- **Portable Programming System (PPS):** for use with ESprit™ series, SPrint™ and Spectra sound processors.
- **Clinical Programming System (CPS):** for use with ESprit series, SPrint and Spectra sound processors.

To connect a sound processor:

1. Connect the Programming Pod, PPS or CPS to the computer on which Custom Sound® software is running.
Programming Pods must be connected to a powered USB port and not a USB port on the monitor or keyboard.
2. Connect the sound processor to the programming interface.
Note: Ensure you have the correct sound processor for the recipient.
3. If Custom Sound software does not automatically detect the sound processor, click the **Connect to an Attached Processor** button  on the Programming toolbar.
Alternatively, click the **Processor** menu and select **Connect**.
The status bar indicates the sound processor is connected, and displays the programs that have been written to each programming slot (if any).
4. If the firmware of the connected sound processor is not supported by the current version of Custom Sound software, the firmware is automatically updated. If any existing programs in the sound processor are incompatible with the updated firmware, the Update Processor prompt displays:
 - Click **Erase Programs and Update** to erase the unsupported programs and update the firmware.
 - Click **Disconnect** to disconnect the sound processor from Custom Sound without updating the firmware.

To reset a sound processor and clear all current programs from the programming slots, use one of the following options:

- Click the **Processor** menu, select **Reset** > (*Sound Processor*) and click **Yes** to confirm.
- Right-click on the **Sound Processor** button  on the status bar, select **Reset** and click **Yes** to confirm.

To disable communication between Custom Sound software and an attached sound processor, click the **Tools** menu and select **Disable Pod Communication**.

Note: To configure the port used to connect a PPS, access the Login window, click **Hardware Options** and select the desired parameters.

Connect a remote assistant

Custom Sound® software allows you to apply updated software to a connected remote assistant when required.

To connect a remote assistant:

1. Connect the remote assistant via a USB cable to the computer on which Custom Sound software is running.

The status bar indicates the remote assistant is connected.

2. If the software of the connected remote assistant is not supported by the current version of Custom Sound software, the Update Remote Assistant prompt displays:
 - Click **Update** to update the remote assistant software.
 - Click **Do Not Update** to close the prompt without updating the software.

If you choose not to update the remote assistant software when prompted, the software can be updated at any time by right-clicking on the **Remote Assistant** button  on the status bar and selecting **Upgrade Software**.

Start a session

The Recipient screen allows you to select a recipient in the Recipient list and start an intraoperative or postoperative session.

To start an intraoperative testing session for the selected recipient, use one of the following options:

- Click **Intraoperative Testing**.
- Click the **Session** menu and select **Intraoperative Testing**.

The Programming screen displays, with the available intraoperative tasks listed in the IntraOp Tasks navigation bar.

To start a postoperative programming session for the selected recipient, use one of the following options:

- Double-click the recipient in the Recipient list.
- Click **Program a Recipient**.
- Click the **Session** menu and select **Program a Recipient**.
- Click the **Recipient** button  in the status bar.

The Programming screen displays, with the available programming tasks listed in the Programming Tasks navigation bar.

To close a session at any time, click the **Session** menu and select **End Session**, or click the **End Programming Session** button  on the Programming toolbar.

Note: If the Meningitis Vaccination Message has been installed as a regional option, the Meningitis Vaccination Message window displays the first time an intraoperative or postoperative session is opened for the recipient. Review the information provided, and click **Continue Session** to acknowledge the message and start the session.

Program recipients

The Programming screen allows you to perform intraoperative testing or program a recipient's sound processor during an initial activation or follow-up session.

The Programming screen displays when you create or select a recipient and start a session.

The Programming screen is colour coded to differentiate between right (red) and left (blue) ear implants, and contains the following areas:

- **IntraOp Tasks/Programming Tasks navigation bar:** provides links to the screens used to test an implant and/or program a sound processor. The **Red**  and **Blue**  arrows indicate the task that is currently selected for the corresponding implant and sound processor.
- **Recipient Summary navigation bar:** provides links to details of the recipient's implant(s) and of previous programming sessions.
- **MAP Summary navigation bar:** displays details of the currently selected MAP.
- **Implant tabs and Sound Processor drop-down lists:** enable selection of the left or right implant and the corresponding sound processor.
- **Programming toolbar:** allows various actions to be applied to the currently open session.
- **Status bar:** displays the programming slots for each sound processor and provides information on the status of tasks as they are performed.

The various programming screens open in the main display area as each task is selected.

The Programming screen allows you to:

- [Obtain impedance measurements](#)
- [Obtain NRT measurements](#)
- [Create a new MAP, and convert, upgrade or open an existing MAP](#)
- [Measure threshold \(T\) and comfort \(C\) levels, and set additional MAP parameters](#)
- [Perform bilateral balancing](#)
- [Set environment options and write programs to the sound processor](#)
- [Generate reports and finalise a programming session](#)

Obtain impedance measurements

The Impedance screen allows an implant test to be performed. An implant test measures the impedance for each electrode and detects any short- or open-circuit electrodes.

The Impedance screen displays automatically when you start a session. To return to the Impedance screen at anytime, click **Measure Impedances** in the Programming Tasks or IntraOp Tasks navigation bar.

Custom Sound® software measures impedances using the following stimulation modes:

- **Common Ground (CG):** one electrode is designated as active while the remaining electrodes together form the indifferent 'electrode'.
- **Monopolar (MP1, MP2 and MP1+2):** one or both of the extracochlear indifferent electrodes is used in combination with intracochlear active electrodes.

An implant test detects the following short- or open-circuit electrodes:

- **Intracochlear electrodes:** detects short- and open-circuit electrodes in Common Ground mode, and open-circuit electrodes in Monopolar modes.
- **Extracochlear electrodes:** detects open-circuit electrodes in Monopolar modes, but does not detect short-circuit electrodes between extracochlear electrodes.

Note:

- The impedance measurement function is not available for Nucleus® 22 implants.
- Impedances cannot be measured using an ESPrit™ series sound processor. Where the implant type supports impedance measurement, an alternate sound processor can be used to perform the implant test.

Measure impedances

An implant test can be performed during an intraoperative or postoperative session. By default, Custom Sound® software measures the impedance of each unflagged electrode in CG, MP1, MP2 and MP1+2 stimulation modes.

To measure impedances:

1. (Intraoperative testing only) Place the coil in a sterile bag on the skin flap covering the implant.
2. To change the default measurement settings, click **Customise** in the Impedance screen.
The Impedance Measurement Options window displays.
3. Select or clear the desired check boxes to specify which stimulation modes and electrodes to test, and click **OK**.
4. Click **Measure**.

The electrodes for which successful impedance measurements were obtained in all stimulation modes display in green . Any electrodes for which short- or open-circuits were detected display in red , and are automatically flagged in postoperative mode only.

Flagged electrodes will display in brown in subsequent programming sessions.

Note: To disable the automatic flagging of electrodes in postoperative mode, access **Tools > My Preferences**, click the **General** tab and clear the **Auto Flag Faulty Electrodes** check box.

View impedance measurements

Impedance measurements are stored in Custom Sound® software and can be viewed at any time.

To view impedance measurements from the last implant test:

1. Click **Details** in the Impedance screen.

Alternatively, click the **Last Implant Test** link  in the Recipient Summary navigation bar.

The Impedance Measurement Details window displays the impedance measurements for each electrode in each stimulation mode. Short- or open-circuit electrodes are indicated by a **Cross** symbol  and flagged electrodes are indicated by a **Flag** symbol .

2. To flag or unflag an electrode, select or clear the corresponding check box in the **Flagged** column.

The Add Note to Electrode window displays. Type any notes next to the date stamp and click **Save and Close**. Electrode notes can be viewed by placing the cursor over the corresponding **Flag** symbol .

3. To print a report of the measurement details, click **Report** and use the Print toolbar on the Report viewer.
4. Click **Close** or **OK** to close the Impedance Measurement Details window.

Note: Automatically flagged electrodes will display with a yellow flag . Manually flagged electrodes will display with a red flag .

To view impedance measurements from historical implant tests:

1. Click the **Implant** menu and select **Impedance Measurement History**.

The Impedance Measurements window displays a list of all impedance measurements obtained for the implant. Impedance measurements for each stimulation mode (by default, CG, MP1, MP2 and MP1+2) are grouped by the date and time they were obtained.

2. Select a group of measurements.

Alternatively, clear the **Group by Date** check box and select individual measurements performed on different dates. Hold down the **Ctrl** key to select multiple measurements (up to a maximum of four).

3. Click **View Measurement Details** to open the Impedance Measurement Details window and view the impedance measurements for each electrode.

Alternatively, click **Print a Report of the Measurement** to open the Report viewer and print the report.

4. Click the **Close** button  to close the Impedance Measurements window.

Flag or unflag electrodes

Manually flagged electrodes are not used when measuring impedances in postoperative mode or when creating new MAPs. Automatically flagged electrodes are retested when an impedance test is performed.

Electrodes can be manually flagged or unflagged prior to repeating an implant test.

To flag or unflag electrodes:

1. Click the **Implant** link  in the Recipient Summary navigation bar.
The Implant Details window displays.
2. Click the **Electrodes** tab.
3. Select or clear the check box in the **Flagged** column for the desired electrode.
The electrode is marked as flagged or not flagged respectively.
4. Type any additional notes for the electrode in the **Electrode Notes** field.
5. Click **OK** to close the Implant Details window.

Obtain NRT® measurements

Neural Response Telemetry (NRT®) records neural activity within the cochlea in response to electric stimulation from the implant. The AutoNRT® screen allows T-NRT levels to be measured.

To access the AutoNRT screen, click **Perform AutoNRT** in the Programming Tasks or IntraOp Tasks navigation bar.

Note:

- The AutoNRT function is not available for ABI541, Nucleus® 24 or Nucleus® 22 implant types.
- AutoNRT cannot be performed using an ESPr™ series sound processor. Where the implant type supports AutoNRT, an alternate sound processor can be used to obtain NRT measurements.

Perform AutoNRT®

AutoNRT® can be performed during an intraoperative or postoperative session. An implant test must be performed prior to obtaining NRT® measurements and any flagged electrodes will not be measured.

To perform AutoNRT:

1. Select the number of electrodes you wish to measure in the Number of Electrodes box on the AutoNRT screen.

Alternatively, click the Channel Numbers in the Channel Grid to select the channels you wish to measure.

By default, Custom Sound® software sets the number of electrodes to 5 (where available) and selects an even distribution along the array.

2. (Postoperative mode only) Select the starting current level (5 to 245 CL) and step size (3 to 6).
3. To view NRT traces while testing is performed, click **Show NRT Traces**.

The NRT Traces window displays in front of the AutoNRT screen.

4. Click **Start**.

Starting at the specified current level, the NRT level gradually increases by the step size indicated until the threshold is detected for each selected channel in turn. If you do not wish the NRT levels to increase automatically, clear the **Auto Increase** check box and click **Increase** until the thresholds are detected (postoperative mode only).

An information message displays when the measurements are complete. The NRT values are automatically saved and are available for use when programming MAPs.

5. To measure additional electrodes, select a greater value in the Number of Electrodes box and repeat the measurement process.

To stop AutoNRT at any time, click **Stop**. Click **Continue** when you are ready to continue, or click **Skip Electrode** to skip the currently selected electrode and measure the next one.

Note:

- (Intraoperative mode only) By default, each electrode is conditioned at 230 current levels prior to measuring the NRT level.

- (Intraoperative mode only) Enhanced interpolation can be selected as a clinician preference, enabling automatic interpolation of NRT levels when the distance between two measured electrodes is less than 8 channels. Measurements are obtained 6 CL above and below the interpolated thresholds. If a neural response is detected above the threshold and none is detected below, the interpolated value is selected. If no response is detected above the threshold or one is detected below, standard AutoNRT is performed.
- The default parameters for performing AutoNRT can be changed in the **NRT** tab of the Preferences window (**F9**).

Work with MAPs

The MAP screen displays a list of all MAPs for the currently selected implant, and can be used to create a new MAP or to open an existing MAP.

To access the MAP screen, click **Open or Create MAP** in the Programming Tasks navigation bar.

MAPs can be sorted by number, title or creation date by clicking on the header of the column you wish to sort by.

Note: The **Last Modified By** column displays the Cochlear™ software that was last used to update the MAP. MAPs that have been created or modified using Nucleus® Fitting Software or Remote Assistant Fitting can be modified using Custom Sound® software.

MAPs created or modified using Custom Sound software can be modified using Nucleus Fitting Software or Remote Assistant Fitting, although certain parameters apply.

To view a report of an existing MAP, right-click on the MAP in the MAP list and select **Report**. The report includes MAP details, channel details, acoustic parameters (where appropriate) and software used and can be printed using the Print toolbar.

The MAP screen allows you to:

- [Create a new MAP or open an existing MAP](#)
- [Upgrade or convert a MAP for use with a new sound processor](#)
- [Compare MAPs](#)
- [Archive or restore a MAP](#)

Create or open a MAP

The MAP screen allows you to create a new MAP or open an existing MAP. Where changes are made to an existing MAP, a new MAP will be automatically created allowing the original MAP to be preserved.

To create a new MAP:

1. Select the sound processor and the desired MAP parameters from the drop-down lists in the Create a New MAP area.

It is recommended the default parameters are used during an initial activation.

2. Click **Create**.

The MAP is created using the specified parameters and opens in the Set Levels screen. Each MAP is assigned a number, incremented by one for each successive MAP.

Note: For CP900 series sound processors, the **Hybrid Mode** check box will be available on the **Open or Create MAP** screen. The check box will be selected if an [audiogram](#) has been created for the recipient. Refer to [Program an acoustic component](#) for further information on acoustic programming.

To open an existing MAP, use one of the following options:

- Double-click the MAP in the MAP list.
- Right-click on the MAP and select **Open**.
- Select the MAP (or hold down the **Ctrl** key and select multiple MAPs) and click **Open**.

- Right-click on the **Program**  in the status bar and select **Open MAP**.

The selected MAP or MAPs open in the Set Levels screen.

To save a MAP at any time, click the **Save MAP** button  on the Set Levels toolbar, or click the **MAP** menu and select **Save** or **Save All**.

To close a MAP, click the **Close MAP** button  on the Set Levels toolbar, or click the **MAP** menu and select **Close** or **Close All**.

Upgrade or convert a MAP

An existing MAP can be upgraded or converted for use with a different sound processor type. Custom Sound® software provides the following options:

- **Upgrade MAP:** the frequency table and as many parameters as possible are set to the defaults of the new sound processor.
- **Convert MAP:** where possible, the frequency table of the original MAP is copied to the new MAP. If the frequency table cannot be copied, Custom Sound software will adjust the frequency boundaries or use the default table for the new sound processor.

To upgrade or convert a MAP:

1. Select the desired MAP in the MAP list.
2. Right-click on the MAP, point to **Upgrade MAP** or **Convert MAP** and select the sound processors you wish to change from and to.

Alternatively, click the **MAP** menu, point to **Upgrade MAP** or **Convert MAP** and select the sound processors you wish to change from and to.

The values of the original MAP are copied to a new MAP. Parameter values that are invalid for the new MAP are set to the default values for the new sound processor type, and details of the changes display in the MAP Summary window.

3. To print a report of the changes, click **Print Summary** and use the Print toolbar on the Report viewer.
4. To add details of the changes to the MAP notes, select **Add Summary to the MAP Notes of the New MAP**, and click **OK**.

MAP notes are saved with the MAP and added to MAP reports.

5. Click the **Close** button  to close the MAP Summary window.

The new MAP displays in the Set Levels screen.

Note:

- MAP conversions and upgrades are not available between Freedom® Hybrid™ sound processors and ESPr™ series or SPrint™ sound processors.
- MAP conversions are available between Freedom® Hybrid™ sound processors and CP900 sound processors.
- MAP upgrades are not available between Freedom® Hybrid™ sound processors and CP900 sound processors.
- The MAP Summary window is only available for Nucleus 24 and Nucleus 22 implant types.

Compare MAPs

The MAP Compare function allows you to compare up to four different MAPs, including those created for different sound processors and different implants.

To compare MAPs created for the same implant:

1. Hold down the **Ctrl** key and select the MAPs you wish to compare in the MAP list.
2. Right-click on one of the highlighted MAPs and select **Compare MAPs**.

Alternatively, click the **MAP** menu and select **Compare MAPs**.

The MAP Comparison window displays the differences between the MAPs, including the T- and C-levels, pulse width, gains and frequency boundaries for each channel.

3. To test a MAP using live voice testing, select the sound processor from the **Processor** drop-down list and click the **Go Live** button  corresponding to the desired MAP. Click the **Stop** button  to stop live voice testing.
4. To print the MAP comparison details, click **Print**. Type any notes to be added to the report in the MAP Comparison Report window, and click **Report**.

The MAP comparison details display in the Report viewer, and can be printed using the Print toolbar.

To compare MAPs created for bilateral implants:

1. Select one of the MAPs you wish to compare in the MAP list.
 2. Right-click on the MAP and select **Compare with Another MAP**.
- Alternatively, click the **MAP** menu and select **Compare with Another MAP**.

The Select MAPs window displays a list of available MAPs for each implant.

3. Select the check boxes of the required MAPs and click **Compare**.

The MAP Comparison window displays the differences between the MAPs.

Note: Custom Sound® software also allows you to compare programs from previous sessions. See [View session history](#).

Archive or restore a MAP

MAPs that are no longer required can be archived from the MAP list. An archived MAP is not deleted from the database, and can be restored at a later date if required.

To archive a MAP:

1. Select the desired MAP in the MAP list.
2. Right-click on the MAP and select **Archive**.

The MAP is removed from the MAP list.

To restore an archived MAP:

1. Click the **MAP** menu and select **Manage**.
The Manage MAPs window displays.
2. Select the MAP you wish to restore in the Archived list.
3. Click **Restore**.

The MAP is removed from the Archived list to the Active list, and is restored to the MAP screen.

4. Click the **Close** button  to close the Manage MAPs window.

MAPs can also be archived in the Manage MAPs window by selecting the MAP in the Active list and clicking **Archive**.

Set levels

The Set Levels screen is used to set threshold (T) and comfort (C) levels for the currently selected MAP, and to adjust MAP parameters as required.

The Set Levels screen displays automatically when you create or open a MAP. To return to the Set Levels screen at any time, click **Set Levels** in the Programming Tasks navigation bar.

The Set Levels screen contains the following areas:

- **Set Levels toolbar:** allows various actions and parameters to be applied to the currently selected MAP.
- **Programming tab:** contains the controls used to set T- and C-levels and check battery suitability. The **Expand**  and **Contract**  arrows allow the various panels to be expanded or collapsed.
- **HearingMentor™ tab:** contains an inbuilt help system for trouble-shooting hearing performance.
- **Channel Grid:** displays channels for each electrode along the array. Each active channel displays the following markers:
 -  C-level
 -  T-level
 -  NRT/objective level

The markers can be adjusted by clicking and dragging them to the desired current level.

Channels corresponding to flagged electrodes are automatically disabled and indicated by a **Flag** symbol . Electrode notes can be viewed by placing the cursor over the symbol. The channels on which high and low frequency tones will be presented are indicated by a **Bell** symbol .

- **Data Grid:** displays the parameters and values that apply to individual channels. The C-level, T-level and dynamic range values display by default, and additional parameters can be viewed by right-clicking on the Data Grid and selecting from the list of options. Values can be typed directly in the Data Grid cells or selected from drop-down lists.
- **Message box:** displays information and warning messages relating to the current MAP.

The Set Levels screen allows you to:

- [Set T- and C-levels using a variety of programming methods](#)
- [Adjust MAP parameters](#)
- [Create additional MAPs from an initial MAP](#)
- [Program an acoustic component \(CP900 series and Freedom™ Hybrid™ sound processors only\)](#)
- [Determine voltage compliance and battery suitability](#)

Work with the Set Levels screen

The following actions can be applied at anytime while working in the Set Levels screen:

- [Select channels to be measured or interpolated](#)
- [Configure the volume, sensitivity and environment settings that apply to live voice testing](#)
- [Talk to a recipient](#)
- [Work with Implant ID](#)
- [Use HearingMentor™ to trouble-shoot hearing performance](#)

Select channels

The currently selected channel displays in gold in the Channel Grid on the Set Levels screen. Multiple channels can be selected by clicking the Channel Numbers, which remain highlighted until the channels are de-selected. The right-click menu displays a list of actions that can be applied to the selected channel or channels.

A narrow bar indicates T- and C-levels are automatically interpolated for that channel. Interpolation is only available in Monopolar stimulation modes and between channels that have the same pulse width. To make an interpolated channel measurable use one of the following options:

- Click the **Make All Channels Measurable** button  on the Set Levels toolbar to access the entire array.
- Highlight specific Channel Numbers in the Channel Grid and click the **Make Only Selected Channels Measurable** button  on the Set Levels toolbar. The non-selected channels are interpolated.
- Double-click an interpolated channel in the Channel Grid to make it measurable (and vice-versa).

Configure live program settings

Custom Sound® software allows you to adjust the volume, sensitivity and environment settings that apply to live voice testing. By default, Custom Sound software sets the ADRO setting and volume and sensitivity levels based on the Everyday environment that was last written to the sound processor. If more than one Everyday environment exists, the settings from the first program are selected. If no Everyday environment exists, the settings from the Noise, Focus or Music environment are selected in that order. Where no previous environments exist, the default Everyday settings are selected.

 **Note:** For CP900 series sound processors, Custom Sound software sets the configuration based on the first program that was last written to the sound processor. If more than one Automatic program exists, the settings from the first program are selected. If no Automatic program exists, the settings from the program in the first slot are selected. Where no previous programs exist, the default Automatic settings are selected.

To configure live program settings:

1. Click the **Configure Live Program Settings** button  on the Set Levels toolbar.
The Live Program Settings window displays.
2. Drag the **Volume** slider to the desired value (0 to 9, or 1 to 10 for CP900 series and CP800 series sound processors).
The specified volume level enables stimulation at C-level.
3. Drag the **Sensitivity** slider to the desired value (0 to 20, or 0 to 31.5 for ESPrit™ series sound processors).
The sensitivity level determines the minimum input signal level required for stimulation. A higher sensitivity setting requires a lower sound processor level to cause stimulation (and vice-versa).
4. Click the **ADRO®** option button to enable the ADRO option.
Alternatively, click the **None** option button to go live with ADRO disabled.

 **Note:** For CP900 series sound processors several parameters can be enabled or disabled, including Audibility (Boost soft sounds, Soften loud sounds) and Noise Reduction (Background, Wind).

5. Click the **Close** button  to close the Live Program Settings window.

The selected settings apply when you go live. The volume and sensitivity levels and the program configuration can be adjusted for individual programs in the Write to Processor screen (see [Configure a program](#)).

Note:

- For ESprit series sound processors, only the sensitivity level can be modified in the Live Program Settings window.
- For Spectra sound processors, the live program settings cannot be configured.

Talk to a recipient

Custom Sound® software allows you to talk to a recipient during a programming session using a previously saved MAP.

To talk to a recipient:

1. Click the **Talk to Recipient** button  on the Programming toolbar.
2. Select the desired MAP in the MAP list and click **Select** (if required).
By default, the MAP in programming slot 1 on the sound processor is used (where available). To use a different MAP, click **Select a Different MAP for Talk Over**, select the desired MAP and click **Select**.
Talk over mode is activated using the selected MAP and the live program settings.
3. To adjust the volume or sensitivity, drag the **Volume** or **Sensitivity** slider to the desired value.
The volume and sensitivity settings reflect the values in the Live Program Settings window, and vice versa.
4. Click **Done** to stop using talk over mode.

Note: For ESprit™ series, SPrint™ and Spectra sound processors, if you select a different MAP to that in programming slot 1, talk over runs the selected MAP in programming slot 1 and the original MAP is overwritten.

Work with Implant ID

The Implant ID feature in the CP900 series, CP800 series, Freedom® Hybrid™ and Freedom® sound processors allows the sound processor to be associated with a specific implant or implant type. This reduces the risk of unintended stimulation due to mixing up sound processors between recipients or between ears. Implant ID is not available for other sound processors.

A padlock icon in the Implant tab on the Programming screen indicates the Implant ID status:

- **CI500 series, CI422, Nucleus® Hybrid™ implants, Freedom® (CI24RE) implants and ABI541 implant types:** the padlock icon displays in the Implant tab after you go live for the first time in a programming session, and the Implant ID updated message displays on the status bar.
- **Nucleus® 24 and Nucleus® 22 implant types:** the padlock icon displays in the Implant tab when you first open a programming session.

Once you write a program to the sound processor, the status bar will also display the padlock icon, indicating that Custom Sound® software has programmed the sound processor with Implant ID. Implant ID will not be activated in the sound processor until a program has been written.

The padlock icons are colour-coded to indicate the level of protection:

-  **Gold:** the sound processor that is programmed for use with this specific implant will not stimulate another implant. Gold padlock protection is not available for Nucleus 24, Nucleus 22 or Freedom® (CI24RE) implant types prior to serial number 1020050437005.

-  **Green:** a sound processor that is programmed for use with a specific type of implant will not stimulate another implant which is a different type.
-  **Grey:** Implant ID has been disabled by a clinician.
- **No icon:** Implant ID is not available.

Where Implant ID is disabled or not available, a sound processor that is placed on the wrong implant could result in unintended stimulation.

If you experience problems with Implant ID, the auto correct option allows you to clear any previously recorded ID value and re-read the ID of that implant (excluding Nucleus 24 and Nucleus 22 implant types). Should the problem persist, Implant ID can be disabled for all implant types excluding Nucleus 22 implants.

To auto correct or disable Implant ID:

1. Ensure the coil is on the recipient's implant.
2. Click the **Implant** link  in the Recipient Summary navigation bar.
The Implant Details window displays.
3. Click the **Implant ID** tab.
4. Click **Auto Correct** or **Disable Implant ID** as required.
5. Click **OK** to apply the selected option and close the Implant Details window.

To re-enable Implant ID, re-access the Implant Details window and click **Enable Implant ID** on the **Implant ID** tab.

Use HearingMentor

HearingMentor provides a knowledge base of commonly reported sound quality symptoms, together with recommended actions.

To use HearingMentor:

1. Click the **HearingMentor** tab.
2. Type a symptom in the **Search** field, or select an option in the list of symptoms.
The Select an Action box displays a list of suggested actions.
3. Select the desired action.
The Details box indicates how the selected action will change the sound quality of the MAP and advises whether the change can be automatically applied by Custom Sound® software.
4. Click **Apply** (where enabled).
Changes are made to the MAP as detailed, and a counter appears alongside the selected action indicating the number of times the action has been applied.
Where Apply is not enabled, the suggested action can be performed manually.

To undo or redo an action applied through HearingMentor, click the **Undo**  or **Redo**  button on the Set Levels toolbar.

Set T- and C-levels

Custom Sound® software includes three streamlined programming methods designed to simplify programming without compromising outcomes. Cochlear™ recommends using one of these methods at the initial activation.

The streamlined programming methods are:

- 
Behavioural method: T-levels are measured on selected channels using psychophysics and C-levels are measured simultaneously using live voice testing. This method is appropriate for recipients who can make reliable behavioural responses to sound.
- 
NRT/objective offset method: the T-level profile is offset from the objective measurement profile and a single offset channel is measured using psychophysics. C-levels are measured using live voice testing. This method is appropriate for recipients who may give limited behavioural responses.
- 
NRT/objective preset method: the T- and C-level profiles are offset from the objective measurement profile and additional MAPs are created with the C-levels set progressively higher on each MAP. This method is appropriate for recipients who may not give any reliable behavioural responses.

In addition to the streamlined programming methods, individual channels can be measured as required.

To pause or stop at any time while setting T- and C-levels, click the **Stop All Stimulations** button  on the Programming toolbar.

Note:

- Streamlined programming is only available for MAPs created using Monopolar stimulation modes.
- Streamlined programming is not available for ABI541, Nucleus® 24 ABI and Nucleus® 22 implant types.

Select NRT/objective levels

NRT® or objective measurements obtained via Custom Sound® software or Custom Sound® EP software can be imported into the Set Levels screen and used as a guide for setting T- and C-levels. NRT/objective measurements must be selected in order to use an NRT/objective streamlined programming method.

To select NRT/objective levels:

- Click **Select NRT/Objective Levels** or **Import NRT/Objective Levels** in the Method panel.
The Select NRT/Objective Levels window displays.
- Select the type of measurement (for example, NRT) from the **Measurement Type** drop-down list.
A list of available measurements displays, detailing the date of measurement and the electrodes measured.
- Select the measurements you wish to use.
Multiple measurements can be selected by holding down the **Ctrl** key and clicking the required measurements.
- Click **OK**.
The NRT/objective levels for the measured channels display in the Channel Grid.

Note:

- For Nucleus® 24 implant types, a prompt displays when **Import NRT/Objective Levels** is selected. Click **From Custom Sound EP** to access the Select NRT/Objective Levels window; or, click **From AGF File** to open an externally saved .agf file from NRT 3.x.
- For Nucleus® 22 implants, NRT levels cannot be measured. The Method panel provides the option to **Select Objective Levels** only.

Use the behavioural method

The behavioural method allows T- and C-levels to be set based on behavioural responses from the recipient. T-levels on the streamlined programming channels (by default, channels 22, 16, 11, 6, and 1 where available) are measured using psychophysics and C-levels are measured simultaneously using live voice testing. Where streamlined programming is supported, the behavioural method is selected by default when a new MAP is created.

To use the behavioural method:

1. Select the step sizes (1 to 10 CL) and the number of stimulations (1 to 999) in the Streamlined Programming panel.
2. Click the **Increase** button  in the Set Ts area.
The T-level on the first streamlined channel increases by the step size indicated, and automatic stimulation occurs based on the number of stimulations specified.
3. Continue to increase the T-level until the threshold is reached.
A smaller step size can be selected and the **Decrease** button  used as the threshold is neared. To manually stimulate a channel at anytime, click the **Stimulate** button .
4. Click each of the remaining streamlined channels in the Channel Grid in turn, and set the T-level for each.
The T-levels for the non-measured channels are automatically interpolated.
5. Click the **Go Live** button .
6. Talk to the recipient and click the **Increase** button  in the Live Cs area.
The C-levels increase globally by the step size indicated.
7. Continue to adjust the C-levels until ongoing conversational speech is comfortably loud.
8. Click the **Stop** button  to stop live voice testing.

Note:

- The default channels for streamlined programming can be changed in the **Streamlined** tab of the Preferences window (F9).
- The default parameters for setting T- and C-levels can be changed in the **General** tab of the Preferences window.

Use the NRT®/objective offset method

The NRT®/objective offset method offsets the T-level profile from the objective measurement profile, and enables a single offset channel to be measured using psychophysics. C-levels are measured using live voice testing.

To use the NRT/objective offset method:

1. Click **Change Programming Method** in the Method panel.
The Select Programming Method window displays.
2. Click the **NRT/Objective Offset** option button, and click **OK**.
If NRT/objective levels have not been selected, the Select NRT/Objective Levels window displays. Select the measurement type and the measurements you wish to use, and click **OK**.

The T-level on the offset channel (by default, channel 11 where available) is set to 80 or 100 CL depending on the type of implant. The remaining T-levels are interpolated to keep the difference between the starting T-level and the objective measurement level equal to that of the offset channel.

3. Select the offset channel in the Channel Grid.
4. Select the step sizes (1 to 10 CL) and the number of stimulations (1 to 999) in the Streamlined Programming panel.
5. Click the **Increase** button  in the Set Ts area.
The T-levels increase globally by the step size indicated, and automatic stimulation of the offset channel occurs based on the number of stimulations specified.
6. Continue to increase the T-levels until the threshold is reached.
A smaller step size can be selected and the **Decrease** button  used as the threshold is neared.
7. Click the **Go Live** button .
8. Talk to the recipient and click the **Increase** button  in the Live Cs area.
The C-levels increase globally by the step size indicated.
9. Continue to adjust the C-levels, watching for any indications that sound is too loud.
10. Click the **Stop** button  to stop live voice testing.

 **Note:** The default T-offset channel can be changed in the **Streamlined** tab of the Preferences window (F9).

Use the NRT®/objective preset method

The NRT®/objective preset method offsets the T- and C-level profiles from the objective measurement profile. Additional MAPs are created with the C-levels set progressively higher on each MAP.

To use the NRT/objective preset method:

1. Click **Change Programming Method** in the Method panel.
The Select Programming Method window displays.
2. Click the **NRT/Objective Preset** option button, and click **OK**.
If NRT/objective levels have not been selected, the Select NRT/Objective Levels window displays. Select the measurement type and the measurements you wish to use, and click **OK**.
Custom Sound® software sets T-levels 40 CL below the corresponding objective measurement levels and automatically interpolates any non-measured channels. The C-levels on the initial MAP are set 10 CL above the T-levels. Additional MAPs are automatically created with the C-levels set 10 CL higher in each progressive MAP.
3. Click the **Go Live** button  in the Streamlined Programming panel of each MAP in turn, watching for any indications that sound is too loud.
4. Click the **Stop** button  for each MAP to stop live voice testing.

Note:

- For ESPrit™ series sound processors, only two MAPs are created.

- The default amounts by which the T- and C-levels are offset can be changed in the **Progressive MAPs** tab of the Preferences window (F9).

Measure individual channels

When you do not wish to use a streamlined programming method or streamlined programming is unavailable (for example, for ABI541, Nucleus® 24 ABI or Nucleus® 22 implants), T- and C-levels can be measured on each individual channel.

To measure individual channels:

1. Ensure the behavioural method is selected in the Streamlined Programming panel.
2. Click the **Make All Channels Measurable** button  on the Set Levels toolbar to access the entire array (if required).
3. Select the channel you wish to measure in the Channel Grid.
4. Click the **C** or **T** option button in the Set Ts and Cs panel based on which level you wish to set.
5. Select the step sizes (1 to 10 CL) and the number of stimulations (1 to 999).
6. Click the **Increase**  button.
The T- or C-level on the selected channel increases by the step size indicated, and automatic stimulation occurs based on the number of stimulations specified.
7. Continue to click the **Increase**  and **Decrease**  buttons until the desired level is reached, adjusting the step sizes as required.
8. If you do not wish to measure every channel and the MAP supports interpolation, highlight the Channel Numbers that have been measured in the Channel Grid and click the **Interpolate** button .
The T- and/or C-levels between the measured channels are automatically interpolated.
9. To test the T- and C-levels using live voice testing, click the **Go Live** button  on the Set Levels toolbar.

Click the **Stop All Stimulations** button  on the Programming toolbar to stop live voice testing.

The **T & C** option button allows you to adjust the T- and C-levels for the selected channel simultaneously, keeping the dynamic range constant.

Note:

- Interpolation is only available for MAPs created using Monopolar stimulation modes and between channels that have the same pulse width.
- Interpolation is not available for ABI541 or Nucleus® 24 ABI implants due to the variation in T- and C-levels from one electrode to the next.

Sweep levels

The Sweep function allows the sequential stimulation of each active channel along the electrode array in order to check for channels that are too loud or too soft. Levels can be balanced by sweeping channels in small groups and confirming equal perception of loudness with the recipient.

To sweep levels:

1. Click the **C** or **T** option button in the Sweep panel based on which levels you wish to sweep.

Alternatively, select a percentage from the drop-down list to sweep at a percentage of the dynamic range.

2. Click the **Selected** option button and highlight the desired Channel Numbers in the Channel Grid.
Alternatively, click the **Channels** option button and select a number of channels from the drop-down list.
3. Select the number of sweep cycles you wish to perform (1 to 20).
4. Select the channel you wish to start from in the Channel Grid, and click the **Sweep Left**  or **Sweep Right**  button.

The specified channels are sequentially stimulated in the selected direction.

To pause or stop the sweep at any time, click the **Stop All Stimulations** button  on the Programming toolbar.

Shift levels

The Shift function allows selected T- and/or C-levels to be increased or decreased simultaneously.

To shift levels:

1. Click the **C**, **T** or **T & C** option button in the Shift panel based on which levels you wish to shift.
2. Click the **CL** option button to shift the levels by an absolute amount.
Alternatively, click the **% of DR** option button to shift the levels by a percentage of the dynamic range for each channel.
3. Click the **All** option button to shift the whole profile.
Alternatively, click the **Selected** option button and highlight the desired Channel Numbers in the Channel Grid.
4. Select the amount (1 to 10 CL) or the percentage (1 to 30 %) the levels will be adjusted by.
5. Click the **Increase**  or **Decrease**  button.

The specified levels shift in the selected direction by the amount indicated.

Tilt levels

The Tilt function allows the T- and/or C-level profile to be rotated left or right around a central channel.

To tilt levels:

1. Click the **C**, **T** or **T & C** option button in the Tilt panel based on which levels you wish to tilt.
2. Select the amount the levels will be adjusted by (1 to 10).
3. Click the **Tilt Left**  or **Tilt Right**  button.

The specified levels rotate in the selected direction by the amount indicated.

 **Note:** By default, a non-linear tilt is enabled. To enable a linear tilt, access **Tools > My Preferences**, click the **General** tab and select **Linear** from the **Tilt Profile Type** drop-down list.

Adjust MAP parameters

MAP parameters and channel configuration can be adjusted at the discretion of the clinician. Custom Sound® software allows you to:

- [Adjust basic MAP parameters](#)
- [Adjust advanced MAP parameters](#)
- [Modify channel gains](#)
- [Reorder channels when the tonotopic order of one or more channels is atypical](#)
- [Double-channel MAP](#)
- [Add notes and generate a MAP report](#)

Adjust basic MAP parameters

Basic MAP parameters are applied when a MAP is created.

To adjust basic MAP parameters:

1. Click the **MAP Parameters** button  on the Set Levels toolbar.
The MAP Parameters window displays with the **Basic** tab selected.
2. Adjust the desired fields:
 - **Strategy:** the sound coding strategy defines how the sound processor will analyse acoustic signals and code them for delivery to the implant.
 - **Stimulation Mode:** the stimulation mode determines the location of the indifferent electrode relative to the active electrode, and is applied here at MAP level. The stimulation mode/indifferent electrode and the active electrode can be modified on a per channel basis using the Data Grid on the Set Levels screen.
 - **Channel Rate:** the channel rate (or stimulation rate) determines the frequency in Hertz of the biphasic current pulse delivered to a channel. The channel rate applies to all channels in a MAP.
 - **Maxima:** the maxima is a parameter of the SPEAK, ACE™ and MP3^{co}™ sound coding strategies, and refers to the frequency ranges in an audio signal that contain the greatest amounts of energy. The value specifies the number of maxima that will be selected for any given signal.
 - **Pulse Width:** the value specifies the duration of the pulse in microseconds, and is applied here at MAP level. The pulse width can be modified on a per channel basis using the Data Grid on the Set Levels screen.
 - **Frequency Table:** the frequency table defines the frequency range (or bandwidth) that is assigned to each active channel in the MAP. A given channel receives stimulation when its bandwidth has sufficient energy in it to be selected as a maxima. The frequency table depends upon the sound coding strategy and the number of channels in the MAP.
 - **Amplitude Scale:** the amplitude scale defines the amount of charge delivered for each stimulation, and the type (current level, stimulus level or ESprit™ stimulus level) can be selected for Nucleus® 22 implants only.
3. Click the **Close** button  to save the changes and close the MAP Parameters window.

Adjust advanced MAP parameters

Advanced MAP parameters are applied by default when a MAP is created.

To adjust advanced MAP parameters:

1. Click the **MAP Parameters** button  on the Set Levels toolbar.
The MAP Parameters window displays with the **Basic** tab selected.
2. Click the **Advanced** tab, and adjust the desired fields.
3. Click the **Close** button  to close the MAP Parameters window.

The following parameters are available in the Advanced tab based on the type of sound processor:

CP900 series, CP800 series, Freedom® Hybrid™, Freedom® and SPrint™ sound processors

CP900 series, CP800 series, Freedom® Hybrid, Freedom® and SPrint™ sound processors	
Signal Preconditioning	
	Volume Adjustment (% DR): specifies the maximum percentage of the dynamic range by which C-levels can be adjusted when modifying the volume level.
Tones and Indicators	
	Tone Level (% DR): specifies the loudness of tones as a percentage of the dynamic range.
	Tone Channel: specifies the channel used for tones to indicate successful key press actions (for example, when changing volume levels).
	Low Tone Channel: specifies the channel used for tones to indicate unsuccessful key press actions or warnings (for example, low battery warnings).
Analysis	
	T-SPL (dB): specifies the SPL for a speech signal at which T-level stimulation starts (excluding CP802 and SPrint sound processors).
	C-SPL (dB): specifies the SPL for a speech signal at which C-level stimulation starts (excluding CP802 and SPrint sound processors).
	Loudness Growth: (previously known as Q-value) controls the steepness of the loudness growth function.
	Base Level: specifies the minimum level at which electric stimulation begins (SPrint sound processors only).
Stimulation	
	Jitter: specifies the percentage by which the stimulation rate will be varied in a positive or negative direction for SPEAK MAPs (250 Hz) only. For SPrint sound processors, Jitter can also be applied for ACE™ and CIS MAPs in a negative direction only.
Advanced RF	
	RF Link Integrity Setting: a trouble-shooting parameter to be adjusted under the supervision of a Cochlear representative (CP900 series and CP800 series sound processors only). Please contact your Cochlear representative if you require further information.

ESPr™ series sound processors

ESPr™ series sound processors

Processor Control

Low Tone Battery Alert (mV): specifies the level at which the sound processor will warn that the battery is low.

Low Tone Battery Cut Out Alert (mV): specifies the level at which the sound processor will shut down.

Tones and Indicators

Tone Level (% DR): specifies the loudness of tones as a percentage of the dynamic range.

Tone Duration (msec): specifies the duration of tones.

Tone Channel: specifies the channel used for tones to indicate warnings (for example, low battery warnings).

Analysis

Base Level: specifies the minimum level at which electric stimulation begins.

Extended Base Level (dB): specifies the maximum level to which the base level can rise.

Loudness Growth: (previously known as Q-value) controls the steepness of the loudness growth function.

Stimulation

Jitter: specifies the percentage by which the stimulation rate will be varied in a positive or negative direction for SPEAK MAPs (250 Hz), and in a negative direction only for ACE and CIS MAPs.

Spectra sound processors

Spectra sound processors

Analysis

Base Level: specifies the minimum level at which electric stimulation begins.

Loudness Growth: (previously known as Q-value) controls the steepness of the loudness growth function.

Modify channel gains

In addition to the user-controlled microphone sensitivity that determines the overall gain applied to the input signal, each channel in a MAP includes an adjustable gain control. Channel gains can be modified during live voice testing to improve the perceived sound quality of a MAP.

To modify individual channel gains using the Data Grid:

1. Right-click on the Data Grid and select **Show Gain**.
The gain value displays for each channel.
2. Click the desired cell and press the **Down (↓)** arrow key to select a value from the drop-down list.

To modify global channel gains using the Gain Shapers window:

1. Click the **Modify Gains** button  on the Set Levels toolbar.
The Gain Shapers window displays.

2. Click the desired button in the Actions box:

-  Increase low frequency gains only
-  Increase low and high frequency gains
-  Increase high frequency gains only
-  Increase all frequency gains
-  Centre all frequency gains
-  Decrease all frequency gains
-  Decrease low frequency gains only
-  Decrease low and high frequency gains
-  Decrease high frequency gains only

The gain is applied to each channel output and affects the stimulation levels sent to the implant.

To view the gain value of each channel in the Gain Shapers graph, select the **Show Values** check box.

3. Click the **Close** button  to close the Gain Shapers window.

Note:

- For SPEAK, ACE™ and MP3™ strategies, channel gains are applied before the maxima are selected and can influence the channels selected for stimulation. Reducing the gain on a given channel de-emphasises that channel's output, making it less likely to be selected. Conversely, increasing the gain emphasises the output.
- For CIS strategies, all of the activated channels are stimulated during each frame, regardless of channel amplitude. Channel gains do not influence the channels selected for stimulation, but do increase or decrease stimulation levels.

Reorder channels

Channels may need to be reordered when the tonotopic order of one or more channels is atypical.

The tonotopic order of electrodes usually follows an orderly progression of pitch from high to low from base to apex. Channels that are not in the correct pitch order will not match up correctly with the frequency table used in the MAP, and can result in poor performance. Where channels are clearly out of pitch order and there are no electrode problems, the order of the channels can be changed.

Channel reordering is particularly important for ABI541 and Nucleus® 24 ABI implants as the pitch order is more complex with these devices.

To reorder a channel:

1. Right-click on the Data Grid and select **Show Active Electrode**.
2. Click the **Make All Channels Measurable** button  on the Set Levels toolbar to make all channels measurable (if required).
3. Select the channel you wish to move, hold down the **Ctrl** key and drag the channel to the desired location.

The Data Grid confirms the active electrode of the reordered channel is now stimulating the new channel location. The channel configuration of the adjacent channels change by one, and the affected Channel Numbers display in red. The psychophysical data of the reordered channel remains unchanged.

 **Note:** Pitch order can be confirmed with the recipient by sweeping the channels at C-level after they have been balanced for loudness. See [Sweep levels](#).

Double channels

Doubling a channel allocates two channels to the same electrode pair, and links the channels to ensure they keep the same T- and C-levels.

To double channels:

1. Right-click on the Data Grid and select **Show Active Electrode**.
2. Click the **Make All Channels Measurable** button  on the Set Levels toolbar to access the entire array (if required).
3. Highlight the Channel Numbers you wish to exclude from the MAP in the Channel Grid.
4. Right-click on the Channel Grid and select **Disable Selected Channels**.
5. Right-click on each channel you wish to double and select **Double Channel (#)**.
The selected channel doubles to the next basal channel, and the levels of the two channels are linked. Changes in the channel-to-electrode allocation are reflected in the Data Grid, and the affected Channel Numbers display in red.

Add notes and generate reports

A title and additional notes can be added to a MAP, and a report of the MAP details can be generated and printed.

To add a title and notes to a MAP:

1. Click the **Title/Notes** button  on the Set Levels toolbar.
The Title and Notes window displays.
2. Type a title for the MAP and any additional notes in the **Title** and **Notes** fields.
Notes will be added to MAP reports.
3. Click the **Close** button  to save the changes and close the Title and Notes window.

To generate a MAP report:

1. Click the **Save MAP** button  on the Set Levels toolbar to save the MAP (if required).
2. Click the **MAP Report** button  on the Set Levels toolbar.
The Report viewer displays the details of the MAP, including sound processor and implant details, channel details and MAP parameters.
3. Use the Print toolbar to print the MAP report.

Create additional MAPs

Custom Sound® software allows you to create MAPs from an initial MAP using the following methods:

- [Adjust the stimulation rate and pulse width, and predict the T- and C-levels](#)
- [Create progressive MAPs with the C-levels automatically adjusted](#)

Predict T- and C-levels

A saved MAP can be used to create an additional MAP with a different stimulation rate and pulse width, and the T- and C-levels can be predicted for the new MAP. The predict T- and C-levels function can only be used with Monopolar stimulation modes, and is not available for Nucleus® 22 implants.

To predict T- and C-levels:

1. Open the initial saved MAP.
2. Click the **MAP Parameters** button  on the Set Levels toolbar.
The MAP Parameters window displays.
3. Select a new value from the **Channel Rate** and/or **Pulse Width** drop-down lists.
4. Click the **Close** button .
The new MAP is created and the T- and C-levels are reset.
5. Select a single channel in the Channel Grid (typically Channel 11) and measure the T- and C-levels.
6. Click the **Predict Levels** button  on the Set Levels toolbar.
Using the measured channel as a base, the shape of the T- and C-level profiles on the original MAP is applied to the new MAP on channels with the same stimulation mode and pulse width as the measured channel.

For more information on measuring T- and C-levels, see [Measure individual channels](#).

Create progressive MAPs

After creating an initial MAP, additional MAPs can be created with the C-levels automatically adjusted in each progressive MAP.

To create progressive MAPs:

1. Open the initial MAP.
2. Click the **Create Progressive MAPs** button  on the Set Levels toolbar.
The initial MAP defaults to programming slot 1, and three additional MAPs are created for the remaining programming slots:
 - Programming slot 2: MAP 1 C-levels + 5 CL
 - Programming slot 3: MAP 1 C-levels + 10 CL
 - Programming slot 4: MAP 1 C-levels - 5 CL

Note:

- For ESPr™ series sound processors, a progressive MAP is created for programming slot 2 only.
- For Spectra sound processors, progressive MAPs cannot be created because there are no additional programming slots in the sound processor.

- The default values by which the C-levels are shifted can be changed in the **Progressive MAPs** tab of the Preferences window (F9).

Program an acoustic component

Hybrid™ recipients use acoustic stimulation (sound processor and acoustic component) and electric stimulation (sound processor and cochlear implant) to access the entire sound spectrum.

When creating a MAP for a Freedom® Hybrid™ sound processor or a CP900 series sound processor (when the acoustic component is required), the acoustic component needs to be programmed in addition to the electric component. Custom Sound® software automatically creates a non-overlapping MAP, where low frequency information is presented via acoustic stimulation only and high frequency information is presented via modified frequency-to-electrode allocation tables.

Custom Sound software provides three prescriptive methods which specify the target gain and maximum output level for each frequency band of amplification, and can be used as a baseline for programming the acoustic component:

- CHP (Cochlear Hybrid Prescription)
- NAL-RP (National Acoustic Laboratories-Revised, Profound)
- DSL (Desired Sensation Levels)

A choice of two fitting procedures can be used with the prescriptive methods:

- **WDRC (Wide Dynamic Range Compression):** designed to amplify low-level input signals more than high-level input signals to provide audibility, while ensuring comfortable volume levels at all times. WDRC may be suitable for recipients with loudness tolerance or loudness recruitment issues. WDRC is the default setting for CHP.
- **Linear:** provides the same amount of amplification (gain) to all input sound levels until the maximum power output is achieved. Linear gain may be suitable for recipients who do not like to hear too much background noise in quiet, or when the gain for soft sounds increases the potential for feedback (particularly combined with large vents). Linear is the default setting for NAL-RP and DSL prescriptions.

Set acoustic parameters

Custom Sound® software allows you to set the acoustic parameters for a Freedom® Hybrid™ or a CP900 series sound processor based on the recipient's pure tone audiometric data. Gain and maximum power output (MPO) can be programmed for the low frequency range of hearing (125 to 1800 Hz).

To set the acoustic parameters:

1. Click the **Acoustics** button  on the Set Levels toolbar.
The Acoustics window displays.
2. If the recipient's audiogram has not already been entered, click **View Recipient's Audiogram**.
The Recipient Details window displays with the **Acoustics** tab selected. Select the points on the left and/or right audiogram that correspond to the recipient's pure tone audiometric data, and click **OK**. Click **Recalculate** in the Set Prescription panel.
 **Note:** If an existing audiogram has not been reviewed in the past twelve months, a new audiogram should be entered.
3. Select **Show ZCC Values** or **Show Real Ear Values** from the drop-down list.

The Frequency Graph displays four curves: Target Gain and Aided Gain (60 dB SPL input) and Target MPO and Aided MPO (OSPL 90). The lower frequency boundary for electric stimulation is automatically set based on the recipient's audiogram.

4. Select a prescription method (DSL, NAPL-RP or CHP) from the **Method** drop-down list in the Set Prescription panel, and click the **WDRC** or **Linear** option button.

By default, CHP WDRC is selected.

5. Select the vent size (0 to 3 mm) in the **Vent Size** field.

Custom Sound software automatically recalculates the gain and MPO levels.

6. Click the **MPO** or **Gain** option button in the Set Gains/MPO panel based on the amplification characteristics you wish to set.

7. Select the step size (1 to 3), and click the **Go Live** button .

8. Select a channel in the Frequency Graph, and use the left **Increase**  and/or **Decrease**  buttons to adjust the level on the selected channel only.

Alternatively, use the right buttons to adjust the overall gain or MPO.

9. Click the **Stop** button  to stop live voice testing.

To balance loudness:

1. Go live using acoustic stimulation only by removing the coil from the recipient's implant. Obtain a loudness rating of 'comfortable but slightly soft' by adjusting the overall gain only.
2. Go live using electric stimulation only by clearing the **Hybrid Mode** check box. Obtain a loudness rating of 'comfortable but slightly soft'.
3. Go live using both electric and acoustic stimulation by selecting the **Hybrid Mode** check box. Confirm an overall loudness rating of comfortable.

Note:

- The boundary of where acoustic stimulation stops and electric stimulation begins is determined by the limit of recipient's residual hearing (defined as those frequencies where thresholds ≥ 90 dBHL).
- Additional acoustic channels can be added if required by right-clicking on the lowest deactivated channel (i.e. in the shaded grey area) and selecting **Enable Channel**.
- The timing of the onset of the acoustic stimulation relative to the electric stimulation can be adjusted by setting the **Delay** time in the Miscellaneous panel. Adjustment to the delay is not recommended at this stage.

Determine power requirements

Depending on the type of sound processor and the type of implant, Custom Sound® software allows you to:

- [Estimate the maximum C-level for a channel to remain in voltage compliance](#)
- [Determine the battery types suitable for use with a MAP](#)
- [Optimise the power level of the sound processor in order to maximise battery life](#)

Wherever possible, Custom Sound software performs these actions automatically.

Estimate compliance levels

Once impedance measurements have been performed, Custom Sound® software allows you to view the voltage compliance level for each channel in a MAP. Where a channel is out of compliance, the maximum voltage available from the implant is not sufficient to generate the desired current level on that channel.

Compliance levels are automatically estimated when you do one of the following for the first time (or for the first time after changing the pulse width, stimulation mode or stimulation rate):

- Go live with a MAP
- Click a **Battery** icon  in the Battery Suitability box
- Write a program to the sound processor

To estimate compliance levels at any time:

1. Ensure the coil is on the recipient's implant.
2. Click the **Estimate Compliance Levels** button  on the Set Levels toolbar.

Red markers display in the Channel Grid, indicating the compliance levels for the battery type selected in the Battery Suitability box.

Note: If impedance measurements have not been performed, the red markers will not be displayed.

Where the C-level on a channel exceeds the compliance level, the channel is deemed out of compliance and displays in red.

3. To display the compliance levels for another battery type, click the corresponding check box in the Battery Suitability box.

Note:

- For ABI541 and Nucleus® 24 ABI implants, compliance levels can only be estimated by clicking the **Estimate Compliance Levels** button .
- Compliance levels cannot be estimated for Nucleus® 22 implants.
- Compliance levels cannot be estimated for ESprit™ series, SPrint™ and Spectra sound processors.
- If you wish compliance levels to be estimated by default when a MAP is first opened or created, access **Tools > My Preferences**, click the **General** tab and clear the **Delay Auto Compliance Until Going Live** check box.

Check battery suitability

The Battery Suitability box displays the battery types applicable to the connected sound processor, and allows you to check which battery types are suitable for use with a MAP. A question mark next to a battery icon indicates battery suitability needs to be checked, and this information updates whenever the basic MAP parameters are changed. The battery types may include (from left to right): 2 Zinc Air, 3 Zinc Air, Rechargeable and 2 AAA.

To check battery suitability:

1. Click a **Battery** icon  in the Battery Suitability box.

Custom Sound® software performs the required measurements, and the Estimated Minimum Battery Life window displays the estimated hours for each battery type. Battery types that may not be suitable are highlighted in red or orange, and additional information displays in the Message box. Estimated hours are not available when the red marker is displayed.

2. Click the **Close** button  to close the Estimated Minimum Battery Life window.

The Battery Suitability box displays a green, orange or red marker for each battery type:

-  **Green:** the battery type is suitable for use with the MAP.
-  **Orange:** the battery type may result in sound degradation or sound processor cut-outs towards the end of battery life.
-  **Red:** the battery type cannot deliver enough power for the MAP, which may result in unreliable operation or intermittency.

Note:

- Battery suitability cannot be checked for ESPr™ series, SPrint™ and Spectra sound processors.
- For Nucleus® 22 implants, you are prompted to confirm the recipient's coil is on the implant when you click a battery icon. Alternatively, if a skin flap measurement has not already been taken, you are prompted to estimate the skin flap thickness and optimise the power level.

Optimise the power level

Custom Sound® software allows an automatic or manual power setting to be applied to a MAP. It is recommended Auto Power is used wherever possible, enabling the power level for the sound processor to be automatically optimised.

The type of power setting currently selected displays in the Battery Suitability box (where available):

-  Auto Power (default)
-  Manual Power

Once determined, the calculated power level displays in the **Power** tab of the MAP Parameters window, or can be viewed by placing the cursor over the power icon in the Battery Suitability box.

The power level can be optimised for the following sound processors only:

CP900 series, CP800 series, Freedom® Hybrid™ and Freedom® sound processors (excluding Nucleus® 22 implants)

Custom Sound software automatically optimises the power level when you:

- Click a **Battery** icon  in the Battery Suitability box
- Write a program to the sound processor

If the coil is not on the implant when power optimisation occurs, the power level is set to the previously saved level (where available) or the default level.

Note: When the power level is optimised, a stimulation occurs which may be heard by the recipient.

Freedom sound processors (Nucleus 22 implants)

Custom Sound software prompts you to measure the skin flap and optimise the power level when you first do one of the following:

- Go live with a MAP
- Click a **Battery** icon  in the Battery Suitability box
- Write a program to the sound processor

The calculated power level is used for all MAPs created or opened in the same programming session.

To optimise the power level:

1. Ensure the coil is on the recipient's implant.
2. Click the **Measure Skin Flap and Optimise Power Level** button  on the Set Levels toolbar.

Alternatively, click **Yes** when the Auto Skin Flap Measurement prompt displays.

Custom Sound software measures the thickness of the skin flap and optimises the power level. The skin flap measurement displays in the Message box, and will be included in the Session History report.

If you choose not to measure the skin flap, the power level is set to the previously saved level (where available) or the default level.

ESPril™ 3G sound processors (Nucleus 24 implant types)

Custom Sound software prompts you to perform an implant measurement when you first do one of the following for the initial MAP only:

- Go live with a MAP
- Write a program to the sound processor

The implant measurement is saved for the recipient until such time as a new measurement is taken. Custom Sound software uses the measurement when automatically calculating the optimised power level for each new MAP. You cannot go live or write programs to the sound processor until the implant measurement is performed, or a manual power setting is selected.

To obtain the implant measurement and optimise the power level:

1. Ensure the coil is on the recipient's implant.
2. Click the **Implant** menu and select **Auto Power Setting**.
Alternatively, click **OK** when the Optimise Power prompt displays.
The Auto Power Setting window displays.
3. Confirm the correct sound processor and coil type are selected, and click **Measure**.
An information message displays when measurement is complete.
4. Click **Yes, Accept and Close** to accept the measurement.
5. Click the **MAP Parameters** button  on the Set Levels toolbar, and click the **Power** tab.

Alternatively, click the **Power** icon  in the Battery Suitability box.

The auto power level is shown. The power level is calculated to ensure all channels remain in compliance, and will be recalculated if the T- and C-levels are adjusted.

6. Click **Voltage**.
The Voltage Compliance window displays the channel voltage and the auto power level, and can be used to determine if or how you can reduce the power level for the MAP.

Note: For ESPrIt 3G sound processors with Nucleus 24 implants only:

- Power levels can only be optimised when impedance measurements have been obtained. Impedances can be measured using an alternate sound processor type.
- Power levels cannot be optimised for MAPs created using Pseudomonopolar or Variable stimulation modes.

ESPrIt 3G and ESPrIt™ 22 sound processors (Nucleus 22 implants)

Custom Sound software prompts you to measure the skin flap and optimise the power level when you first do one of the following:

- Go live with a MAP
- Write a program to the sound processor

You cannot go live or write programs to the sound processor until the skin flap measurement is performed, or a manual power setting is selected.

To optimise the power level:

1. Ensure the coil is on the recipient's implant.
2. Click the **Implant** menu and select **Auto Power Setting**.
Alternatively, click **OK** when the Optimise Power prompt displays.
The Auto Power Setting window displays.
3. Confirm the correct sound processor and coil type are selected, and click **Measure**.
An information message displays when measurement is complete.
4. Click **Yes, Accept and Close** to accept the measurement.
The power level is automatically optimised.

Apply a manual power level

When it is not possible to use an automatic power setting, a manual power level can be applied to a MAP.

To apply a manual power level:

1. Click the **MAP Parameters** button  on the Set Levels toolbar.

Alternatively, click the **Power** icon  in the Battery Suitability box.

The MAP Parameters window displays.

2. Click the **Power** tab (if required).
3. Click the **Manual** option button and type a power level in the corresponding field.
4. Click the **Close** button  to close the MAP Parameters window.

Determine compatibility of MAPs

The message box on the Set Levels screen may display messages regarding the compatibility of the MAP to other Cochlear™ software (Nucleus® Fitting Software or Remote Assistant Fitting).

For the MAP to be modified by the recipient using Remote Assistant Fitting or by a Clinician using Nucleus Fitting Software, it must be within compliance levels. If the MAP is not compatible with other Cochlear software, updates can only be made using Custom Sound® software.

[Estimating compliance levels](#) will indicate where the MAP is not compatible with other Cochlear software. If modifications will be made to the MAP using Remote Assistant Fitting or Nucleus Fitting Software, [adjust the MAP parameters](#) to allow compatibility.

The following criteria must be met to allow modifications to a MAP by Remote Assistant Fitting or Nucleus Fitting Software:

Criteria to allow MAP adjustments by Remote Assistant Fitting or Nucleus® Fitting Software	Remote Assistant Fitting	Nucleus Fitting Software
The calculated power level for the MAP must not be above the compliance threshold for standard battery operation. This criteria applies whether the MAP uses Auto or Manual power.	✓	
Impedances and link constants are available to be written with the MAP.	✓	
MAP stimulation and channel modes are all equal to each other and must be either Monopolar 1, Monopolar 2 or Monopolar 1+2.	✓	✓
MAP pulse width and channel pulse widths are all equal to each and the pulse width is less than or equal to 100us.	✓	✓
All enabled channels have an active electrode number (AE) equal to the channel number.	✓	✓
The frequency table is equal to the number of enabled channels in the MAP.	✓	✓
The MAP strategy must be either ACE, SPEAK or MP3 sm ™.	✓	✓
The MAP must have 12 or more enabled channels.	✓	✓

Perform bilateral balance

The Bilateral Balance screen allows the loudness of bilateral MAPs to be balanced simultaneously.

To access the Bilateral Balance screen, click **Bilateral Balance** in the Programming Tasks navigation bar.

The Bilateral Balance screen displays the currently open MAP for each implant. To display another MAP, click the **Select MAP** button  for the left or right side, select the desired MAP from the MAP list and click **Select**.

Balance T- and C-levels

The loudness of bilateral MAPs can be balanced using simultaneous live voice testing.

To balance T- and C-levels:

1. Click the **Go Live** button  for the right or left side in the Bilateral Balance screen.
Alternatively, click the central **Go Bilateral Live** button to go live on both MAPs simultaneously.
2. Click the **C, T** or **T & C** option button in the desired Shift box based on which levels you wish to shift.
You can shift levels on the right or left side, or use the Bilateral Shift box to adjust both MAPs simultaneously.
3. Select the amount (1 to 10 CL) the levels will be adjusted by.
4. Click the **Increase**  or **Decrease**  button.
The specified levels shift in the selected direction by the amount indicated.
5. Continue to adjust the levels for each side until the desired levels are reached.
If required, the frequency table, sensitivity and volume for each side can be adjusted in the corresponding Parameters box.
6. Click the **Stop** button  for each side to stop live voice testing.

 **Note:** Click the **Stop** button  at any time if the recipient experiences discomfort.

Configure and write programs

The Write to Processor screen is used to apply listening environments to MAPs and to write programs to the sound processor.

To access the Write to Processor screen, click **Write to Processor** in the Programming Tasks navigation bar.

The Write to Processor screen contains the following areas:

- **Programs:** displays up to four Program boxes depending on the type of sound processor, each comprising a MAP and an environment.
- **Processor:** indicates the MAP number and the corresponding environment that will be written to each programming slot on the sound processor.

For CP900 series sound processors SmartSound® iQ provides two options for configuring programs:

-  **SCAN:** provides adjustments in response to an automatic classification of the listening environment.
-  **Program 2:** allows full customised configuration for microphone directionality, loudness control and noise reduction. The custom option also allows for the selection of remote assistant icons.

For all other sound processors, the SmartSound® environments contain the recommended settings for different listening situations, and up to four environments can be applied to a single MAP:

 **Everyday:** for typical listening situations.

 **Noise:** for listening situations with significant background noise.

 **Focus:** for focusing on a person or small group and removing background noise.

 **Music:** for listening to live music or music from an audio source.

When the Write to Processor screen is first accessed during a session, the Program boxes automatically populate the following:

- The MAP that is currently selected in the Set Levels screen.
- The environments that were last written to the sound processor, irrespective of the MAP or MAPs they were applied to. Where no previous environments exist, the default environments are used.
- The volume and sensitivity settings only from the Live Program Settings window.

When creating progressive MAPs (including MAPs created using the NRT®/objective preset method), the Program boxes automatically populate the progressive MAPs, and assign the environments that were last written to the sound processor in conjunction with previous progressive MAPs.

Where no previous environments exist, the Everyday environment or the Automatic environment (CP900 series sound processors) is assigned to each MAP.

To populate the Program boxes with the programs currently in the sound processor, click **Read** and click **Yes** to confirm.

The Write to Processor screen also displays the power setting for each MAP and the battery suitability for the sound processor (where available). Click a **Power** icon  or a **Battery** icon  to view the estimated hours for each battery type in the Estimated Minimum Battery Life window. If compliance levels and battery suitability have not previously been checked, Custom Sound® software will perform these measurements when an icon is clicked. For more information, see [Determine power requirements](#).

Configure a program

Each program comprises a MAP and a corresponding environment, and can be configured as desired. In addition, the sound processor can be configured for CP900 series and CP800 series sound processors.

To configure a program:

1. Ensure the desired MAP is displayed in the Program box.

To display another MAP, click the **Select MAP** button , select the desired MAP in the MAP list and click **Select**.

2. To apply a different environment, select the desired option from the **SmartSound Environment** drop-down list.

 **Note:** The **SmartSound Environment** drop-down list option is not available for CP900 series sound processors.

3. To change the settings for the environment, click the **Configure Environment** button .

4. Select the desired settings, and click **OK**.

The settings available vary based on the type of sound processor.

5. Select values in the **Volume**  and **Sensitivity**  fields (where available).
6. (CP900 series and CP800 series sound processors only) To change the settings for the sound processor, click the **Configure Processor** button .

 **Note:** The Configure Processor button will vary based on the type of sound processor.

The Processor Configuration window displays.

7. Select the desired settings, and click **OK**.

The selected settings apply to all programs. For bilateral recipients with two CP900 series or CP800 series sound processors, the selected settings apply to both sound processors.

To reset the environment or sound processor settings to the Custom Sound® software defaults, click **Reset** in the Environment Configuration or Processor Configuration window. For Freedom® Hybrid™ sound processors, Freedom® sound processors, CP900 series and CP800 series sound processors in the Processor Configuration window only, click **Defaults > Adult Defaults** or **Defaults > Paediatric Defaults** to select the default settings recommended for adults or paediatrics respectively.

To enable an empty programming slot, select **Empty** from the **MAP** drop-down list in the Program box.

The following settings are available in the Environment Configuration and Processor Configuration windows based on the type of sound processor:

CP920 and CP910 sound processors

CP920 and CP910 sound processors		
Environment configuration	Default setting	
	SCAN 	Program 2 
Microphone Directionality		
Standard: Zoom and Beam™ microphone directionality are disabled.		✓
Zoom: enables the microphones to be focused in the forward direction.		
Beam™: enables an adaptive beamformer designed to improve the signal-to-noise ratio of speech coming from one direction in a noisy environment.		
Audibility		
Boost soft sounds (ADRO®): if selected, ADRO® enables a digital pre-processing algorithm designed to make soft sounds more audible in quiet and maintain comfort for loud sounds. If Custom programming is being used, the Whisper™ check box is also available. Whisper™ enables a fast-acting compression scheme for improved access to soft and moderately loud sounds.	✓	✓
Soften loud sounds (ASC): if selected, Auto Sensitivity Control is enabled. If Custom programming is being used, the Less (60dB) check box is available with the Auto Sensitivity Control breakpoint set to 60 decibels.		✓
Noise Reduction		
Background: enables automatic background noise reduction.	✓	✓
Wind: enables automatic wind noise reduction. Comfort is improved when wind is automatically detected, by switching to a low noise microphone setting and suppressing loud channels.	✓	✓
Remote Assistant Icon		
If Custom programming is being used, an icon can be selected for the remote assistant that will be associated with the program. Either a numbered icon can be selected, or one of 10 icons depicting a scene are available.  Note: the choice of icon does not dictate the program type or place any restrictions on configuration.		✓

Processor configuration		Default setting	
		Adults	Paediatrics
Recipient-adjustable settings			
Processor Button Lock: locks the buttons on the sound processor.			
Telecoil Mixing Ratio: specifies the ratio of the telecoil input in relation to the microphone input.		3:1	1:1
Accessory Mixing Ratio: specifies the ratio of the accessory input in relation to the microphone input.		2:1	1:1
Processor Indicator Lights: displays a light on the sound processor to indicate functionality or fault conditions. The options available in the list are: - Child - indicator lights are displayed for microphone input, auxiliary input, alarms and processor auto-off when the sound processor is on or off the recipient's head. - Monitor - indicator lights are displayed for auxiliary input when the sound processor is off the recipient's head, alarms and processor auto-off. - Adult - indicator lights are displayed for alarms when the sound processor is off the recipient's head and processor auto-off.		Monitor	Child
Private Tones: delivers a private tone when a warning is signalled or a key is pressed on the sound processor or remote assistant.  Note: the choice of icon does not dictate the program type or place any restrictions on configuration.		✓	✓
Clinician-adjustable settings			
Allow Telecoil: enables the telecoil to be activated on the sound processor. When selected, Telecoil can be turned on or off by the recipient using the sound processor or the remote assistant.		✓	✓
Allow Auto Telecoil: enables the sound processor to automatically activate the telecoil in response to an incoming signal. When selected, Auto Telecoil can be turned on or off by the recipient using the remote assistant. Allow Telecoil must be selected to allow Auto Telecoil to be enabled.		✓	
Auto Processor Off: enables the sound processor to turn off automatically when the coil is off the implant for more than two minutes.		✓	
Processor Interface: specifies the user interface option for the sound processor. When Simple is selected, volume and sensitivity cannot be adjusted using the controls on the sound processor. When Advanced (Volume) or Advanced (Sensitivity) is selected, volume or sensitivity respectively can be adjusted via the sound processor in addition to the remote assistant.		Simple	Simple

Remote assistant settings		
Allow Volume Control: enables the volume to be adjusted using the controls on the remote assistant. When cleared, the volume is applied at the level set in the Program box and cannot be adjusted by the remote assistant.	✓	✓
Allow Sensitivity Control: enables the sensitivity to be adjusted using the controls on the remote assistant. When cleared, the sensitivity is applied at the level set in the Program box and cannot be adjusted using the remote assistant.		
Loudness Control Preference: specifies the loudness control option for the remote assistant. Available options are Volume or Sensitivity.	Volume	Volume
Allow Hearing Profile Adjustments: enables the recipient to use Remote Assistant Fitting to adjust a MAP created by their clinician using Custom Sound®. The MAP must meet criteria to be adjusted via Remote Assistant Fitting.		

CP810 and CP802 sound processors

CP810 and CP802 sound processors

Environment configuration

Default setting

CP810 sound processors only

SmartSound options

ADRO®: enables a digital pre-processing algorithm designed to make soft sounds more audible in quiet and maintain comfort for loud sounds.

✓

✓

✓

✓

Autosensitivity™: enables the automatic adjustment of microphone sensitivity based on the noise floor of the surrounding environment.

✓

✓

✓

Zoom: enables the microphones to be focused in the forward direction.

✓

Beam™: enables an adaptive beamformer designed to improve the signal-to-noise ratio of speech coming from one direction in a noisy environment.

✓

Whisper™: enables a fast-acting compression scheme for improved access to soft and moderately loud sounds.

✓

None: enables SmartSound options to be disabled.

Parameters

Auto-sens. Brkpt-SPL (dB): specifies the sound pressure level at which the Autosensitivity breakpoint is set. When the noise floor is above or below the breakpoint, the Autosensitivity control gradually reduces or increases the sensitivity respectively.

57

Processor configuration

Default setting

Adults

Paediatrics

Recipient-adjustable settings

Processor Button Lock: locks the buttons on the sound processor.

Telecoil Mixing Ratio: specifies the ratio of the telecoil input in relation to the microphone input (CP810 sound processors only).

3:1

1:1

Accessory Mixing Ratio: specifies the ratio of the accessory input in relation to the microphone input (CP810 sound processors only).

2:1

1:1

Indicator Light - Sound: displays a light on the sound processor to indicate an incoming audio signal is being received.

✓

Indicator Light - Key Presses and Alerts: displays a light on the sound processor when a warning is signalled or a key is pressed on the sound processor or remote assistant. The light is only applicable when the coil is on the implant.

✓

Private Tones: delivers a private tone when a warning is signalled or a key is pressed on the sound processor or remote assistant.

✓

✓

 **Note:** Recipients can change the recipient-adjustable parameters at any time using the remote assistant.

Clinician-adjustable settings

Allow Telecoil: enables the telecoil to be activated on the sound processor. When selected, Telecoil can be turned on or off by the recipient using the sound processor or the remote assistant (CP810 sound processors only).	✓	✓
Allow Auto Telecoil: enables the sound processor to automatically activate the telecoil in response to an incoming signal. When selected, Auto Telecoil can be turned on or off by the recipient using the remote assistant. Allow Telecoil must be selected to allow Auto Telecoil to be enabled (CP810 sound processors only).	✓	
Auto Processor Off: enables the sound processor to turn off automatically when the coil is off the implant for more than two minutes.	✓	
Processor Interface: specifies the user interface option for the sound processor. When Simple is selected, volume and sensitivity cannot be adjusted using the controls on the sound processor. When Advanced (Volume) or Advanced (Sensitivity) is selected, volume or sensitivity respectively can be adjusted via the sound processor in addition to the remote assistant.	Simple	Simple

Remote assistant settings

Allow Volume Control: enables the volume to be adjusted using the controls on the remote assistant. When cleared, the volume is applied at the level set in the Program box and cannot be adjusted by the remote assistant.	✓	✓
Allow Sensitivity Control: enables the sensitivity to be adjusted using the controls on the remote assistant. When cleared, the sensitivity is applied at the level set in the Program box and cannot be adjusted using the remote assistant.	✓	✓
Allow Hearing Profile Adjustments: enables the recipient to use Remote Assistant Fitting to adjust a MAP created by their clinician using Custom Sound®. The MAP must meet criteria to be adjusted via Remote Assistant Fitting.		
Allow Creation of New Hearing Profile: enables the recipient to use Remote Assistant Fitting to create a brand new MAP using their remote assistant.		

Freedom Hybrid and Freedom sound processors

Environment configuration

Default setting

SmartSound options

ADRO: enables a digital pre-processing algorithm designed to make soft sounds more audible in quiet and maintain comfort for loud sounds.

✓

✓

✓

✓

Autosensitivity: enables the automatic adjustment of microphone sensitivity based on the noise floor of the surrounding environment.

✓

✓

Beam: enables an adaptive beamformer designed to improve the signal-to-noise ratio of speech coming from one direction in a noisy environment.

✓

Whisper: enables a fast-acting compression scheme for improved access to soft and moderately loud sounds.

✓

None: enables SmartSound options to be disabled.

Note: For Freedom Hybrid sound processors, ADRO, Autosensitivity and Whisper apply to electric processing only, while Beam applies to both electric and acoustic processing.

Parameters

Adults

Paediatrics

Auto-sens. Brkpt-SPL (dB): specifies the sound pressure level at which the Autosensitivity breakpoint is set. When the noise floor is above or below the breakpoint, the Autosensitivity control gradually reduces or increases the sensitivity respectively.

57

57

Telecoil Mixing: enables the sound processor to receive a mix of telecoil and microphone input.

Telecoil Always On: enables the telecoil on the sound processor to be on at all times.

Sensitivity Mode: when Manual is selected, enables the sensitivity to be adjusted using the controls on the sound processor. When Fixed is selected, the sensitivity is applied at the level set in the Program box and cannot be adjusted by the recipient.

Manual

Manual

Volume Control: enables the volume to be adjusted using the controls on the sound processor. When cleared, the volume is applied at the level set in the Program box and cannot be adjusted by the recipient.

✓

✓

Accessory Mixing Ratio: specifies the ratio (1:1 to 10:1) of the accessory input in relation to the microphone input.

2:1

1:1

Telecoil Mixing Ratio: specifies the ratio (1:1 to 10:1) of the telecoil input in relation to the microphone input.

3:1

1:1

Tones and indicator lights		Adults	Paediatrics
Private Tones: delivers a private tone when a warning is signalled.		✓	✓
Private Key Press Notification: delivers a private tone when a key is pressed on the sound processor.		✓	✓
Public Tones: delivers a public tone when a warning is signalled.			✓
Public Key Press Notification: delivers a public tone when a key is pressed on the sound processor.			✓
Indicator Light - Sound: displays a light on the sound processor to indicate an incoming audio signal is being received (Freedom sound processors only).			✓
Indicator Light - Help: displays a light on the sound processor when a warning is signalled (Freedom sound processors only).			✓

ESPr™ series sound processors

ESPr™ series sound processors

Environment configuration

Default setting

SmartSound options

Autosensitivity: enables the automatic adjustment of microphone sensitivity based on the noise floor of the surrounding environment.

✓

Parameters

Auto-sens. Brkpt-SPL (dB): specifies the sound pressure level at which the Autosensitivity breakpoint is set. When the noise floor is above or below the breakpoint, the Autosensitivity control gradually reduces or increases the sensitivity respectively.

21*, -15**

Base Level Breakpoint (dB): specifies the noise level in decibels.

10*, -18**

Telecoil/Lapel Mic Mixing: enables the sound processor to receive a mix of telecoil/lapel microphone and microphone input.

Sensitivity Mode: when Manual is selected, enables the sensitivity to be adjusted using the controls on the sound processor. When Fixed is selected, the sensitivity is applied at the level set in the Program box and cannot be adjusted by the recipient.

Manual

Volume Control: enables the volume to be adjusted using the controls on the sound processor. When cleared, the volume is applied at the level set in the Program box and cannot be adjusted by the recipient. Volume Control can only be selected when Sensitivity Mode is set to Fixed or Autosensitivity is enabled.

*ESPr, **ESPr™ 3G and ESPr™ 22

SPrint™ sound processors

SPrint sound processors		
Environment configuration		Default setting
		 
SmartSound options		
ADRO: enables a digital pre-processing algorithm designed to make soft sounds more audible in quiet and maintain comfort for loud sounds.		☒ ☒
Autosensitivity: enables the automatic adjustment of microphone sensitivity based on the noise floor of the surrounding environment.		☐ ☒
Parameters		
Volume Control: enables the volume to be adjusted using the controls on the sound processor. When cleared, the volume is applied at the level set in the Program box and cannot be adjusted by the recipient.		☒

Spectra sound processors

Spectra sound processors	
Environment configuration	Default setting
	
Parameters	
Auto-sens. Brkpt-SPL (dB): specifies the sound pressure level at which the Autosensitivity breakpoint is set. When the noise floor is above or below the breakpoint, the Autosensitivity control gradually reduces or increases the sensitivity respectively.	

Select a previous environment

 **Note:** This option is not available for CP900 sound processors.

Environments that have been written to the sound processor in a previous session can be selected and applied to a new program or programs. The environments retain all previous configuration settings and volume and sensitivity settings.

To select a previous environment for a single program:

1. Click the **Select a Previously Used Environment** button  in the desired Program box.
The Select Environment window displays a list of all environments written to the sound processor, grouped by the date and time of the session.
2. Click the desired environment.
The corresponding environment, volume and sensitivity settings display.
3. Click **Use Environment**.
The environment displays in the selected Program box.

To select all environments from a previous session:

1. Click the **Select a Previously Used Environment** button  in any Program box.
The Select Environment window displays a list of all environments written to the sound processor, grouped by the date and time of the session.
2. Click the sound processor for the desired session.
The environments written to the programming slots display. Place the cursor over each programming slot to view the corresponding environment, volume and sensitivity settings.
3. Click **Use all Environments**.
The environments display in the Program boxes in the order they were previously written.

 **Note:** The environment icons will vary depending on the type of sound processor and environments being used.

Write to a sound processor

Once programs have been configured they can be written to the recipient's sound processor.

To write to a sound processor:

1. Click the **Go Live**  and **Stop**  buttons to test each program (optional).
2. Click **Write**.
The programs are written to the programming slots on the sound processor. The status bar indicates progress and displays the MAP number(s) and corresponding environments when writing is complete.

To write the same programs to a second sound processor, disconnect the current sound processor and connect the other. The Write to Processor screen retains the program configurations, which can be written to the second sound processor.

Finalise a session

Once programs have been written to the sound processor, the Finalise Programming screen allows you to add session notes and generate reports prior to closing the programming session.

To access the Finalise Programming screen, click **Finalise Programming** in the Programming Tasks navigation bar.

The following reports can be generated from the Finalise Programming screen:

- **Recipient Take Home Report:** provides details of each program written to the sound processor, including the MAP number and title, the environment, suitable battery types and any MAP notes.
- **Session Report:** provides details of the recipient's programming session, including the program details and the MAP report for each MAP written to the sound processor.

The Current Session Details area displays the environment and MAP number for each program written to the sound processor during the session:

- To view the volume and sensitivity settings for an environment, place the cursor over the desired **Environment** icon.
- To view the settings for a CP900 series or a CP800 series sound processor in the Processor Configuration window, click **Sound Processor Configuration**.

To finalise the programming session:

1. Type any notes about the session in the Session Notes area.
Notes are included on the session report and saved against the session history.
2. To generate reports, select the **Recipient Take Home Report** and/or **Session Report** check boxes.
3. Click **End Session**.
The session closes, and the selected reports display in the Report viewer.
4. Use the Print toolbar to print the reports.
5. Click the **Close** button  to close the Report viewer and return to the Recipient screen.

Note: Additional reports are available throughout Custom Sound® software:

- Recipient details report (see [Create a recipient](#))
- MAP report (see [Add notes and generate reports](#))
- MAP upgrade/conversion report (see [Upgrade or convert a MAP](#))
- MAP comparison report (see [Compare MAPs](#))
- Program comparison report (see [View session history](#))
- Historical session and recipient take home reports (see [View session history](#))

View session history

The details of each programming session are saved in Custom Sound® software, and previous sessions can be viewed in the Session History window.

To view session history:

1. Click the **Session** menu and select **Session History**.

Alternatively, click the **Last Session History** link  in the Recipient Summary navigation bar.

The Session History window displays the previous sessions by date and time, the clinician name, the software and software version and the session notes. For bilateral recipients, session details for both the right and left implants display.

2. Select a session in the Session History list.

The Sound Processor panel displays the environment and MAP number for each program written to the sound processor during the selected session.

3. Select from the following actions in the Sound Processor panel:

- To view the volume and sensitivity settings for an environment, place the cursor over the desired **Environment** icon.
- To view the settings for an environment in the Environment Configuration window, click the desired **Environment** icon.
- To open a MAP in the Set Levels screen, click the desired **MAP** link.
- To view the settings for a CP900 series or CP800 series sound processor in the Processor Configuration window, ensure the correct sound processor is selected and click **Sound Processor Configuration**.
- To restore the programs from the selected session to the connected sound processor, click **Write to Processor**.

4. To add notes or edit existing notes for the selected session, type the desired text in the Notes panel.
5. To print a report for the selected session, click **Print Report** or **Print Recipient Take Home Report** in the Session panel and use the Print toolbar on the Report viewer.

To compare programs:

1. Hold down the **Ctrl** key and select the sessions you wish to compare in the Session History list.
2. Click **Compare Programs** in the Session panel.

The Select Programs window displays a list of available programs for each programming session.

3. Select the check boxes of the required programs and click **Compare**.

The Program Comparison window displays the differences between the MAPs and the environments.

4. To test a program using live voice testing, select the sound processor from the **Processor** drop-down list and click the **Go Live** button  corresponding to the desired program. Click the **Stop** button  to stop live voice testing.

5. To print the program comparison details, click **Print**. Type any notes to be added to the report in the Program Comparison Report window, and click **Report**.

The program comparison details display in the Report viewer, and can be printed using the Print toolbar.

Click the **Close** button  to save any changes and close the Session History window.

View usage data

A recipient's usage data can be displayed via the data viewer.

View CP900 usage data

1. Select the recipient from the list of recipients on the Recipient screen.
2. Click the View Usage Data button (at the bottom left of the screen), or click Recipient in the dropdown menu (at the top of the screen) and select View Usage Data.

The data viewer can be accessed at any point when the recipient's record is open.

The Data Viewer panel displays pie charts showing daily averages for Time on Air, Scenes and Program usage.

 **Note:** For a bilateral recipient two sets of data will be displayed.

3. To view Loudness data, click on the Show button to the left of the Scenes graph and select Loudness. The Loudness data will display.
4. To view Accessory usage, click on the Show button to the left of the Program usage graph and select Accessory usage. The Accessory usage data will display.

Trend graphs: For Time on Air, Scenes and Loudness, a trend graph is displayed under the pie chart display data for the last 5 periods. If the mouse is rolled over each of these graphs, they will display the length of each period in days.

 **Note:** Each period corresponds to the dates the sound processor is connected to Custom Sound.

Volume or Sensitivity can be selected by moving the rollover switch under the Program usage pie chart. This controls whether volume or sensitivity information is shown when the mouse is rolled over the chart.

Accessory Usage: If the mouse is rolled over the Accessory usage pie chart, the number of uses of each accessory during the last period is displayed

 **Note:** The sound processor does not need to be connected to see previous periods. Up to five periods are shown in the Data viewer.

Work with databases

A database is automatically created when Custom Sound® Suite is first installed. The database is shared across Custom Sound® software and Custom Sound® EP software, and contains the records for all recipients created in either application. The Cochlear™ Database Manager allows you to manage databases and to select the database you wish to connect to.

To access the Cochlear Database Manager, navigate to the Windows® **Start** menu > **Programs** > **Cochlear** and select **Cochlear Database Manager**.

The Cochlear Database Manager contains the following tabs:

- **Connections:** displays the list of connections to existing databases. When you log on to Custom Sound software or Custom Sound EP software, the software connects to the database specified by the current connection. The connection that is currently in use is indicated by a **Check** symbol .
- **Databases:** displays the list of existing databases. The Databases tab allows you to connect to the database server, create new databases and manage existing databases.

Create or edit a connection

The Cochlear™ Database Manager allows you to create a new connection to an existing database, and to edit or delete existing connections.

To create a connection:

1. Click **Add** in the **Connections** tab.
The Add Connection window displays. Mandatory fields are indicated by a **Required** symbol .
2. Type a name for the connection in the **Connection Name** field.
3. Type the name of the server on which the database resides in the **Server\Instance Name** field.
4. If a password is required to log on to the database server, select the **Use a Specific Username and Password** option button and type the username and password in the appropriate fields.
Otherwise, retain the **Use Windows Integrated Security** option button.
5. Select the database you wish to connect to from the **Database** drop-down list.
6. To adjust the amount of time that Custom Sound® software or Custom Sound® EP software will wait for a response from the database server, type the number of seconds in the **Connection Timeout** field.
7. Click **OK**.
The connection is added to the Connections list.
8. To select the new connection as the default, right-click on the connection and select **Set as Active Connection**.
The connection is marked with a **Check** symbol .

To edit a connection, right-click on the connection and select **Edit Connection**. Edit the existing details as required, and click **OK**.

To delete a connection, right-click on the connection, select **Delete Connection** and click **Yes** to confirm the deletion. Deleting a connection does not delete the database it connects to.

Connect to a server

In order to create or manage databases, you must first connect to the server on which the databases reside.

To connect to a server:

1. Click the **Databases** tab in the Cochlear™ Database Manager.
2. Click **Connect**.

The Connect window displays.

3. Type the name of the server on which the databases reside in the **Server\Instance Name** field, or retain the default setting.
4. Retain the **Use Windows Integrated Security** option button.

Alternatively, select the **Use a Specific Username and Password** option button and type the username and password for the database server in the appropriate fields.

5. Click **OK** to save the changes and close the Connect window.

The databases that exist on the server display in the Database list.

Once the connection to the server has been established, the Cochlear Database Manager allows you to:

- [Create a new database](#)
- [Upgrade an existing database](#)
- [Backup or restore an existing database](#)
- [Delete an existing database](#)

Create a database

The Cochlear™ Database Manager allows you to create a new database as desired. By default, a connection to the new database is automatically created, but is not selected as the default connection.

To create a database:

1. Click **Create** in the **Databases** tab.

The Create Database window displays.

2. Type a name for the database in the **Database** field.
3. Type the clinic and clinician names in the appropriate fields.

The clinician name is used to automatically create a clinician with clinic administrator rights in Custom Sound® software.

4. Type a name for the connection in the **Connection Name** field.
5. Retain the **Use Windows Integrated Security** option button.

Alternatively, select the **Use a Specific Username and Password** option button and type the username and password for the database server in the appropriate fields.

6. Click **Create**.

The database displays in the Database list, and a connection is created in the **Connections** tab.

If you do not wish to automatically create a new connection when the database is created, clear the **Create a New Connection for This Database** check box.

Upgrade a database

An existing database may need to be upgraded when a new version of the Custom Sound® Suite software is released. By default, the Cochlear™ Database Manager backs up the existing database prior to performing the upgrade.

To upgrade a database:

1. Select the database in the Database list.
2. Click **Upgrade**.
3. Type a file location and filename for the backup file, or retain the default file path.
4. Click **Upgrade**.

The database is upgraded.

If you do not wish to back up the database prior to upgrading, clear the **Backup the Database Before Upgrading** check box.

 **Note:** When you upgrade the Custom Sound Suite Software from 2.0 or later, the installer provides the option to upgrade the database. When the upgrade option in the installer is selected, the Custom Sound Suite database that is currently in use is upgraded. Any additional databases or databases created using an earlier version of the software can be upgraded using the Cochlear Database Manager.

Back up or restore a database

The Cochlear™ Database Manager allows you to back up a database to an external file. A database can be restored from a previous backup if required.

To back up a database:

1. Select the database in the Database list.
 2. Click **Backup**.
- Alternatively, right-click on the database and click **Backup**.
3. Type a file location and filename for the backup file, or retain the default file path.
 4. Click **Backup**.

A backup of the database is saved to the specified location.

To restore a database:

1. Click **Restore** in the **Databases** tab.
 2. Type a new name for the database in the **Database** field.
 3. Type the file path of the file you wish to restore in the **Backup File** field.
- Alternatively, click **Browse**, navigate to the desired file, and click **OK**.
4. Click **Restore**.

The restored database displays in the Database list.

Delete a database

Databases that are no longer required can be deleted from the Database list. A deleted database cannot be restored. A database can only be restored from a backup file, and it is recommended you back up the database prior to deletion.

To delete a database:

1. Select the database you wish to delete in the Database list.
2. Click **Delete**.
Alternatively, right-click on the database and select **Delete Database**.
3. Click **Yes** to confirm the deletion.

Appendix 1: Cochlear™ Nucleus® CR120/CR220 Intraoperative Remote Assistant

This section provides clinicians who use Custom Sound® software with an introduction to using the intraoperative remote assistant.

Import recipient data

Recipient AutoNRT® and impedance data can be imported into Custom Sound software. Prior to importing the data, it must have been:

- exported from a remote assistant
- saved to a location it can be imported from.

To download AutoNRT and impedance data from a remote assistant:

1. Connect the recipient's remote assistant.
2. Right-click on the remote assistant icon in the status bar.
3. Select **Export AutoNRT and Impedance measurements**.

The Export Remote Assistant AutoNRT and impedance measurements to... window will display to allow you to select the file location.

4. Navigate to the location you wish to save the .crf file.
5. Click **Save**.

The .crf file will be saved to the selected location.

To import recipient AutoNRT and impedance data:

1. Click the **File** menu and select **Import AutoNRT and impedance measurements**.

The Open Cochlear Remote File window will display to allow the selection of the data file exported from the remote assistant.

2. Navigate to the .crf file you wish to import.
3. Click **Open**.

The file is imported into Custom Sound Software.

The data will be imported and stored in the Custom Sound Software database. Data that corresponds to an implant known to Custom Sound Software will be associated with the relevant implant. If the implant does not exist, the data will be saved and matched to an implant either when a new implant with the correct Implant ID is created, or an existing implant has its Implant ID changed to match the Implant ID in the imported data.

After the data import process has completed, Custom Sound Software will display the results depending on the success of the import:

- If no new valid data has been saved to the database, Custom Sound Software will display a message in the status bar to indicate the reason why no changes have been made, or
- When data has been successfully associated with implants, a window will display the recipients, implants, and any additional information relating to the imported data.

Appendix 2: Remote Assistant Fitting

This section provides clinicians who use Custom Sound® software with an introduction to remote assistant fitting.

Work with MAPs

The MAP screen displays a list of all MAPs for the currently selected implant, and can be used to create a new MAP or to open an existing MAP.

To access the MAP screen, click **Open or Create MAP** in the Programming Tasks navigation bar.

MAPs can be sorted by number, title or creation date by clicking on the header of the column you wish to sort by.

Note: The **Last Modified By** column displays the Cochlear™ software that was last used to update the MAP. MAPs that have been created or modified using Remote Assistant Fitting can be modified using Custom Sound Software.

MAPs created or modified using Custom Sound Software can be modified using Remote Assistant Fitting, although certain parameters apply.

To view a report of an existing MAP, right-click on the MAP in the MAP list and select **Report**. The report includes MAP details, channel details, acoustic parameters (where appropriate) and software used and can be printed using the Print toolbar.

Compatibility of MAPs

The message box on the Set Levels screen may display messages regarding the compatibility of the MAP for use with Remote Assistant Fitting.

For the MAP to be modified using Remote Assistant Fitting it must be within compliance limits. If the MAP is not within compliance limits, updates can only be made using Custom Sound Software.

Estimating compliance levels in Custom Sound Software will indicate where the MAP is not compatible with Remote Assistant Fitting. If modifications will be made to the MAP using Remote Assistant Fitting, adjust the MAP parameters to allow compatibility.

Criteria to allow MAP modifications by Remote Assistant Fitting	
	The calculated power level for the MAP must not be above the compliance threshold for standard battery operation. This criteria applies whether the MAP uses Auto or Manual power.
	Impedances and link constants are available to be written with the MAP.
	MAP stimulation and channel modes are all equal to each other and must be either Monopolar 1, Monopolar 2 or Monopolar 1+2.
	MAP pulse width and channel pulse widths are all equal to each other and the pulse width is less than or equal to 100µs.
	Channel interphase gaps are all equal to each other.
	All enabled channels have an active electrode number (AE) equal to the channel number.
	The frequency table is equal to the number of enabled channels in the MAP.
	The MAP strategy must be either ACE, SPEAK or MP3 ⁰⁰⁰ ™.
	The MAP must have 12 or more enabled channels.

Configure a program

Each program comprises a MAP and a corresponding environment, and can be configured as desired. In addition, the sound processor can be configured for CP900 series and CP800 series sound processors.

To configure a program:

1. Ensure the desired MAP is displayed in the Program box.

To display another MAP, click the **Select MAP** button , select the desired MAP in the MAP list and click **Select**.

2. To apply a different environment, select the desired option from the **SmartSound Environment** drop-down list.

 **Note:** The **SmartSound Environment** drop-down list option is not available for CP900 series sound processors.

3. To change the settings for the environment, click the **Configure Environment** button .
4. Select the desired settings, and click **OK**.

The settings available vary based on the type of sound processor.

5. Select values in the **Volume**  and **Sensitivity**  fields (where available).
6. (CP900 series and CP800 series sound processors only) To change the settings for the sound processor, click the **Configure Processor** button .

 **Note:** The Configure Processor button will vary based on the type of sound processor.

The Processor Configuration window displays.

7. Select the desired settings, and click **OK**.

The selected settings apply to all programs. For bilateral recipients with two CP900 series or CP800 series sound processors, the selected settings apply to both sound processors.

To reset the environment or sound processor settings to the Custom Sound® software defaults, click **Reset** in the Environment Configuration or Processor Configuration window. For Freedom® Hybrid™ sound processors, Freedom® sound processors, CP900 series and CP800 series sound processors in the Processor Configuration window only, click **Defaults > Adult Defaults** or **Defaults > Paediatric Defaults** to select the default settings recommended for adults or paediatrics respectively.

To enable an empty programming slot, select **Empty** from the **MAP** drop-down list in the Program box.

The following settings are available in the Environment Configuration and Processor Configuration windows based on the type of sound processor:

CP920 and CP910 sound processors

CP920 and CP910 sound processors		
Environment configuration	Default setting	
	SCAN	Program 2
		
Microphone Directionality		
Standard: Zoom and Beam™ microphone directionality are disabled.		✓
Zoom: enables the microphones to be focused in the forward direction.		
Beam™: enables an adaptive beamformer designed to improve the signal-to-noise ratio of speech coming from one direction in a noisy environment.		
Audibility		
Boost soft sounds (ADRO®): if selected, ADRO® enables a digital pre-processing algorithm designed to make soft sounds more audible in quiet and maintain comfort for loud sounds. If Custom programming is being used, the Whisper™ check box is also available. Whisper™ enables a fast-acting compression scheme for improved access to soft and moderately loud sounds.	✓	✓
Soften loud sounds (ASC): if selected, Auto Sensitivity Control is enabled. If Custom programming is being used, the Less (60dB) check box is available with the Auto Sensitivity Control breakpoint set to 60 decibels.		✓
Noise Reduction		
Background: enables automatic background noise reduction.	✓	✓
Wind: enables automatic wind noise reduction. Comfort is improved when wind is automatically detected, by switching to a low noise microphone setting and suppressing loud channels.	✓	✓
Remote Assistant Icon		
If Custom programming is being used, an icon can be selected for the remote assistant that will be associated with the program. Either a numbered icon can be selected, or one of 10 icons depicting a scene are available.  Note: the choice of icon does not dictate the program type or place any restrictions on configuration.		✓
Processor configuration	Default setting	
	Adults	Paediatrics
Recipient-adjustable settings		
Processor Button Lock: locks the buttons on the sound processor.		

Telecoil Mixing Ratio: specifies the ratio of the telecoil input in relation to the microphone input.	3:1	1:1
Accessory Mixing Ratio: specifies the ratio of the accessory input in relation to the microphone input.	0:1	1:1
Processor Indicator Lights: displays a light on the sound processor to indicate functionality or fault conditions. The options available in the list are: - Child - indicator lights are displayed for microphone input, auxiliary input, alarms and processor auto-off when the sound processor is on or off the recipient's head. - Monitor - indicator lights are displayed for auxiliary input when the sound processor is off the recipient's head, alarms and processor auto-off. - Adult - indicator lights are displayed for alarms when the sound processor is off the recipient's head and processor auto-off.	Monitor	Child
Private Tones: delivers a private tone when a warning is signalled or a key is pressed on the sound processor or remote assistant.  Note: the choice of icon does not dictate the program type or place any restrictions on configuration.	✓	✓
Clinician-adjustable settings		
Allow Telecoil: enables the telecoil to be activated on the sound processor. When selected, Telecoil can be turned on or off by the recipient using the sound processor or the remote assistant.	✓	✓
Allow Auto Telecoil: enables the sound processor to automatically activate the telecoil in response to an incoming signal. When selected, Auto Telecoil can be turned on or off by the recipient using the remote assistant. Allow Telecoil must be selected to allow Auto Telecoil to be enabled.	✓	
Auto Processor Off: enables the sound processor to turn off automatically when the coil is off the implant for more than two minutes.	✓	
Processor Interface: specifies the user interface option for the sound processor. When Simple is selected, volume and sensitivity cannot be adjusted using the controls on the sound processor. When Advanced (Volume) or Advanced (Sensitivity) is selected, volume or sensitivity respectively can be adjusted via the sound processor in addition to the remote assistant.	Simple	Simple
Remote assistant settings		
Allow Volume Control: enables the volume to be adjusted using the controls on the remote assistant. When cleared, the volume is applied at the level set in the Program box and cannot be adjusted by the remote assistant.	✓	✓
Allow Sensitivity Control: enables the sensitivity to be adjusted using the controls on the remote assistant. When cleared, the sensitivity is applied at the level set in the Program box and cannot be adjusted using the remote assistant.		
Loudness Control Preference: specifies the loudness control option for the remote assistant. Available options are Volume or Sensitivity.	Volume	Volume
Allow Hearing Profile Adjustments: enables the recipient to use Remote		

Assistant Fitting to adjust a MAP created by their clinician using Custom Sound®. The MAP must meet criteria to be adjusted via Remote Assistant Fitting.		
---	--	--

CP810 and CP802 sound processors

CP810 and CP802 sound processors					
Environment configuration		Default setting			
					
		CP810 sound processors only			
SmartSound options					
ADRO®: enables a digital pre-processing algorithm designed to make soft sounds more audible in quiet and maintain comfort for loud sounds.		✓	✓	✓	✓
Autosensitivity™: enables the automatic adjustment of microphone sensitivity based on the noise floor of the surrounding environment.		✓	✓	✓	
Zoom: enables the microphones to be focused in the forward direction.			✓		
Beam™: enables an adaptive beamformer designed to improve the signal-to-noise ratio of speech coming from one direction in a noisy environment.				✓	
Whisper™: enables a fast-acting compression scheme for improved access to soft and moderately loud sounds.					✓
None: enables SmartSound options to be disabled.					
Parameters					
Auto-sens. Brkpt-SPL (dB): specifies the sound pressure level at which the Autosensitivity breakpoint is set. When the noise floor is above or below the breakpoint, the Autosensitivity control gradually reduces or increases the sensitivity respectively.		57			
Processor configuration		Default setting			
		Adults	Paediatrics		
Recipient-adjustable settings					
Processor Button Lock: locks the buttons on the sound processor.					
Telecoil Mixing Ratio: specifies the ratio of the telecoil input in relation to the microphone input (CP810 sound processors only).		3:1		1:1	
Accessory Mixing Ratio: specifies the ratio of the accessory input in relation to the microphone input (CP810 sound processors only).		2:1		1:1	
Indicator Light - Sound: displays a light on the sound processor to indicate an incoming audio signal is being received.				✓	
Indicator Light - Key Presses and Alerts: displays a light on the sound processor when a warning is signalled or a key is pressed on the sound processor or remote assistant. The light is only applicable when the coil				✓	

	is on the implant.		
	Private Tones: delivers a private tone when a warning is signalled or a key is pressed on the sound processor or remote assistant.	✓	✓
 Note: Recipients can change the recipient-adjustable parameters at any time using the remote assistant.			
Clinician-adjustable settings			
	Allow Telecoil: enables the telecoil to be activated on the sound processor. When selected, Telecoil can be turned on or off by the recipient using the sound processor or the remote assistant (CP810 sound processors only).	✓	✓
	Allow Auto Telecoil: enables the sound processor to automatically activate the telecoil in response to an incoming signal. When selected, Auto Telecoil can be turned on or off by the recipient using the remote assistant. Allow Telecoil must be selected to allow Auto Telecoil to be enabled (CP810 sound processors only).	✓	
	Auto Processor Off: enables the sound processor to turn off automatically when the coil is off the implant for more than two minutes.	✓	
	Processor Interface: specifies the user interface option for the sound processor. When Simple is selected, volume and sensitivity cannot be adjusted using the controls on the sound processor. When Advanced (Volume) or Advanced (Sensitivity) is selected, volume or sensitivity respectively can be adjusted via the sound processor in addition to the remote assistant.	Simple	Simple
Remote assistant settings			
	Allow Volume Control: enables the volume to be adjusted using the controls on the remote assistant. When cleared, the volume is applied at the level set in the Program box and cannot be adjusted by the remote assistant.	✓	✓
	Allow Sensitivity Control: enables the sensitivity to be adjusted using the controls on the remote assistant. When cleared, the sensitivity is applied at the level set in the Program box and cannot be adjusted using the remote assistant.	✓	✓
	Allow Hearing Profile Adjustments: enables the recipient to use Remote Assistant Fitting to adjust a MAP created by their clinician using Custom Sound®. The MAP must meet criteria to be adjusted via Remote Assistant Fitting.		
	Allow Creation of New Hearing Profile: enables the recipient to use Remote Assistant Fitting to create a brand new MAP using their remote assistant.		

Appendix 3: Working with Nucleus® Fitting Software

This section provides clinicians who use Custom Sound® Software with an introduction to Nucleus® Fitting Software, and to describe the interoperability between the two applications.

Nucleus Fitting Software is intended to be used by cochlear implant and hearing care professionals to create and modify hearing profiles in order to provide comfortable and usable stimulation to cochlear implant recipients.

The cochlear implant system components compatible with Nucleus® Fitting Software are listed below.

For additional information on using Nucleus Fitting Software, refer to the Nucleus Fitting Software eHelp.

	Cochlear Nucleus® Sound Processors				
	CP800 Series		Freedom®	CP900 Series	
	CP802 Sound Processor	CP810 Sound Processor	Freedom® Sound Processor	CP910 Sound Processor	CP920 Sound Processor
Cochlear Nucleus® Implants					
Nucleus® CI522 Cochlear Implant	✓	✓	✓	✓	✓
Nucleus® CI513 Cochlear Implant	✓	✓	✓	✓	✓
Nucleus® CI512 Cochlear Implant	✓	✓	✓	✓	✓
Nucleus® CI422 Cochlear Implant	✓	✓	✓	✓	✓
Nucleus® CI24RE (H) Cochlear Implant	✓	✓	✓	✓	✓
Nucleus® CI24RE (CA) Cochlear Implant	✓	✓	✓	✓	✓
Nucleus® CI24RE (ST) Cochlear Implant	✓	✓	✓	✓	✓
Nucleus® CI24RE (CS) Cochlear Implant	✓	✓	✓	✓	✓
Nucleus® CI24M Cochlear Implant	✓	✓	✓	✓	✓
Nucleus® CI24R (ST) Cochlear Implant	✓	✓	✓	✓	✓
Nucleus® CI24R (CS) Cochlear Implant	✓	✓	✓	✓	✓
Nucleus® CI24R (CA) Cochlear Implant	✓	✓	✓	✓	✓

Cochlear Nucleus® Remote Assistants					
Nucleus® CR110 Remote Assistant		✓			
Nucleus® CR230 Remote Assistant		✓		✓	✓

 **Note:** Nucleus Fitting Software does not support auditory brainstem or double array implant types.

Not all products are approved in all markets.

Connecting to Nucleus® Fitting Software

A recipient's sound processor and implant are connected to Nucleus® Fitting Software via a programming pod, in the same way they are connected to Custom Sound® software.

1. Connect the Programming Pod.

The programming pod connection status changes to 'connected' when successfully connected.

2. Connect the sound processor you wish to fit during the session to the Programming Pod system.

The sound processor connection status changes to 'connecting' while the device information is read.

The status changes to 'connected' when successfully connected.

3. Place the coil on the client's head.

The coil connection status changes to 'testing' while automatic tests are run (e.g. an impedance check, power checks). An impedance check is an implant test that measures the impedance for each electrode and detects any short- or open-circuit electrodes.

This testing may take a short time to complete.

The status changes to 'connected' when successfully connected.

MAP parameters in Nucleus® Fitting Software

If a MAP is created using Nucleus Fitting Software it uses the global default MAP parameters:

- Strategy = ACE™
- Rate = 900 Hz
- Pulse width = 25 µs
- Maxima = 8
- Mode = MP1+2.

If a MAP is created using Custom Sound Software with non-default parameters, Nucleus Fitting Software will preserve these parameters, although adjustments cannot be made.

Compatibility of MAPs

All MAPs that have been created using Nucleus Fitting Software can be opened and modified using Custom Sound software.

A MAP created using Custom Sound Software can be modified using Nucleus Fitting Software, although it must be within compliance limits. If the MAP is not compatible, Nucleus Fitting Software will automatically turn down the Master Volume to ensure compliance.

If you know that modifications will be made to a MAP using Nucleus Fitting Software, consider the criteria below when adjusting the MAP parameters.

Criteria to allow MAP adjustments by Nucleus® Fitting Software

MAP stimulation and channel modes are all equal to each other and must be either Monopolar 1, Monopolar 2 or Monopolar 1+2.
MAP pulse width and channel pulse widths are all equal to each and the pulse width is less than or equal to 100µs.
All enabled channels have an active electrode number (AE) equal to the channel number. MAPs with advanced programming (e.g. double channel mapping, reordered channel mapping cannot be programmed using Nucleus Fitting Software).
Non-default frequency allocation tables are not used. The frequency table is equal to the number of enabled channels in the MAP.
The MAP strategy must be either ACE, SPEAK or MP3 ⁰⁰⁰ ™.
The MAP must have 12 or more enabled channels.

Electrode flagging

Electrodes will display in Nucleus Fitting Software with a state of either **On**, **Off** or **Automatically** turned off. The table below displays the action in Custom Sound software and the corresponding state in Nucleus Fitting Software.

Custom Sound Software	Nucleus® Fitting Software
Flagged during an impedance measurement	Automatically turned off
Manually flagged	Electrode is turned off
Channel is disabled	Electrode is turned off

When you turn an electrode **Off** in Nucleus Fitting Software, it displays as manually flagged in Custom Sound.

Other information

Certification, applied standards and copyright information

The Custom Sound® software fulfils the essential requirements listed in Annex 1 of the EC Directive 90/385/EEC on Active Implantable Medical Devices.

The year in which authorization to affix the CE mark was granted was:

2009: System 4 and System 5

2013: System 6

2014: Custom Sound Suite 4.2

Table of symbols

Symbol	Description
	Manufacturer.
	Date of manufacture.
	Authorised representative in the European Community.
	Refer to instruction manual.
	Specific warnings or precautions associated with the device, which are not otherwise found on the label.
	CE mark approved by notified body 0197.

Legal statement

The statements made in this version of the Help are believed to be true and correct as of the date of publication. However, specifications are subject to change without notice.

Copyright 2014 Cochlear Ltd. All rights reserved.

Cochlear, the elliptical logo, Nucleus, ACE, AutoNRT, Autosensitivity, Beam, Whisper, Contour, Contour Advance, Freedom, Hybrid, ESPrit, SPrint, MP3^{ooo}, Hear now And Always and Custom Sound are either trademarks or registered trademarks of Cochlear Ltd.

ADRO is a registered trademark of Dynamic Hearing Pty Ltd. IBM is a registered trademark of IBM Corporation. Intel and Pentium are registered trademarks of Intel Corporation. Microsoft and Windows are registered trademarks of Microsoft Corporation.

Cochlear Ltd (ABN 96 002 618 073) 1 University Avenue, Macquarie University, NSW 2109, Australia

Tel: +61 2 9428 6555 Fax: +61 2 9428 6352

Cochlear Ltd (ABN 96 002 618 073) 14 Mars Road, Lane Cove, NSW 2066, Australia

Tel: +61 2 9428 6555 Fax: +61 2 9428 6352

Cochlear Americas 13059 E Peakview Avenue, Centennial, CO 80111, USA

Tel: +1 303 790 9010 Fax: +1 303 792 9025

Cochlear Canada Inc 2500-120 Adelaide Street West, Toronto, ON M5H 1T1, Canada

Tel: +1 416 972 5082 Fax: +1 416 972 5083

Cochlear AG EMEA Headquarters, Peter Merian-Weg 4, 4052 Basel, Switzerland

Tel: +41 61 205 0404 Fax: +41 61 205 0405

ECREP Cochlear Deutschland GmbH & Co. KG Karl-Wiechert-Allee 76A, 30625 Hannover, Germany

Tel: +49 511 542 770 Fax: +49 511 542 7770

Cochlear Europe Ltd 6 Dashwood Lang Road, Bourne Business Park, Addlestone, Surrey KT15 2HJ, United Kingdom

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Tel: +90 216 538 5900 Fax: +90 216 538 5919

Cochlear (HK) Limited Room 1204, 12/F, CRE Building, No 303 Hennessy Road, Wanchai, Hong Kong SAR

Tel: +852 2530 5773 Fax: +852 2530 5183

Cochlear Korea Ltd 1st floor, Cheongwon building, 828-5, Yuksam dong, Kangnam gu, Seoul, Korea

Tel: +82 2 533 4663 Fax: +82 2 533 8408

Cochlear Limited (Singapore Branch) 6 Sin Ming Road, #01-16 Sin Ming Plaza Tower 2, Singapore 575585

Tel: +65 6553 3814 Fax: +65 6451 4105

Cochlear Medical Device (Beijing) Co Ltd

Unit 2208 Germdale Tower B, 91 Jianguo Road, Chaoyang District, Beijing 100022, P.R. China

Tel: +86 10 5909 7800 Fax: +86 10 5909 7900

Cochlear Medical Device Company India Pvt. Ltd.

Ground Floor, Platina Building, Plot No C-59, G-Block, Bandra Kurla Complex, Bandra (E), Mumbai – 400 051, India

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Cochlear Limited, 1 University Avenue, Macquarie University, NSW 2109,