
Operational documentation

for light therapy device Bioptron 2 with floor stand (PAG-880-SET),
Bioptron Pro 1 with table stand (PAG-990), Bioptron Pro 1 with floor
stand (PAG-991)

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


Document History

Revision	CO/NCR/CAPA	Description of Change	Implementation Date
003	CAPA 54	DEKRA audit: minor non-conformity 3	June 12 th , 2013
004	A0233	Remove label "Keep dry" and harmonization Type Label	November 8 th , 2013
005	A0250	Re-labeling of packagings (all Set devices)	January 09 th , 2014
006	A0208	New BPro1 (Project Beam)	May 8 th , 2014

Statement of Confidentiality

The information contained in this document is considered proprietary to Bioptron AG, and may not be disclosed by others without prior written consent of Bioptron AG. The distribution of this document shall be limited to only those persons with the direct need to evaluate its content.

Pro1 (PAG-990)


BIOPTRON[®]
LIGHT THERAPY SYSTEM By Depter Group

BIOPTRON AG, 8832 Wollerau, Switzerland

1 Piece

BIOPTRON Pro 1

100-240 V~, 50/60 Hz

Part no.: PAG-990
Weight: gross 6.05 kg net 4.60 kg
Volume: 0.059 m³

Production code: 006-2014-20


CE 0124

MADE IN SWITZERLAND

2014

Packaging Label PAG-990

Label on Packaging

"Storage Transport Conditions"

Type: BIOPTRON Pro1 Power: 100-240 V~, 50/60 Hz, 75 W Code: BIOPTRON AG, Made in Switzerland          	Type Label (303176) Type: Model: Power Supply Information: Serial Number: TechnicalCode-Year-Week-Day(1-7)ConsecutiveNumber(1-999) Protection Class: Label "Keep dry" Label "Warning Light Source" Label "WEEE Disposal" Label "Isolation class" Label "CE-Conformity Assessment Body" Label "ETL" (for Canada & US) Label "C-Tick" (for Australia & New Zealand)
	Label Quality Control Passed (302810)
	Label "Read IFU" (303817)

Pro1 Floor Stand (PAG-991-FS)

BIOPTRON  BIOPTRON AG, 8832 Wollerau, Switzerland BPro1 Floor Stand CE  2014 Piece: 1 Part no.: PAG-991-FS Weight: gross 6.86 kg net 6.02 kg Volume: 0.024 m³ Rev: 06  7 612675 001379	Packaging Label Floor Stand
BIOPTRON Pro 1 Floor Stand BIOPTRON AG 8832 Wollerau Switzerland CE  2014 333443	Type Label Floor Stand

Pro1 Table Stand (PAG-990-TS)

BIOPTRON 

BIOPTRON AG, 8832 Wollerau, Switzerland

BIOPTRON Pro1 Table Stand


CE 
2014

Piece: 1
Part no.: PAG-990-TS
Weight: gross 1.90 kg net 1.30 kg
Volume: 0.019 m³
Revision: 00

Packaging Label Table Stand
(EAN code)

Guarantee Label Table Stand
(EAN code)

Certificate


BIOPTRON Floor Stand, B2

General

Customer item Number PAG-883-Y (PAG-880-SET-2)


Technical Information

Item Number	000556
Device description	Accessory to B2 Device (PAG-880/ PAG-880-SET)
Material	Steel with coating CEWEPOL type 15 83936, thermoplastic Cycloy
Weight Net kg.	8.00 (without device)
Dimensions while folded, cm	69 x 20 x 10
Dimensions while unfolded, cm (w x h x l)	72 x 87-147 x 80
Range of height adjustment of the stand, cm	from 87 to 147
The angles of rotation of the rod with the device	от 0° до 360°
Color	RAL9016 (traffic white) / RAL9017 (traffic black)
EC directive	93/42/EEC
Device type	Medical device class I/ supplement
Labelling	CE
GMDN Code	35239

Logistics Information


Manufacturer	BIOPTRON AG, Sihleggstrasse 23, 8832 Wollerau; Switzerland
Customer	Home Art & Sales Services; Sihleggstrasse 23, 8832 Wollerau, Switzerland
Net Weight kg	8.00
Gross Weight kg	8.80
Volume m ³	0.02
Cardbox dimension cm (l x w x h)	70 x 21 x 11
Customs Tariff Code	7326.9098

Date:

17.05.2016

Anouar Hicham
 Head of RA & QA Departments
 BIOPTRON AG

17.05.2016

Tim Scheffele
 General Manager
 BIOPTRON AG

BIOPTRON AG

Sihleggstrasse 23, CH-8832 Wollerau
 Phone: +41 43 900 98 00 Fax: +41 43 888 28 99
 INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU.com

Bank: UBS AG, CH-6900 Lugano
 SWIFT: UBSWCHZH80A, BC 247

IBAN CHF: CH94 0024 7247 5714 9301 B
 IBAN EUR: CH24 0024 7247 5714 9362 L
 IBAN USD: CH57 0024 7247 5714 9361 C

Certificate


BIOPTRON Floor Stand, BPro1

General

Customer item Number PAG-991-FS


Technical Information

Item Number	001379
Device description	Accessory to BPro1 Device (PAG-991)
Material	Steel with coating CEWEPOL type 15 83936, thermoplastic Cycloy
Weight Net kg.	5.98 (without device)
Dimensions while folded, cm	98 x 18 x 10
Dimensions while unfolded, cm (w x h x l)	57.5 x 61.8 -81 x 69.5
Range of height adjustment of the stand, cm	from 61.8 to 81
The angles of rotation of the rod with the device	120° from the vertical axis of the device
Color	RAL9016 (traffic white) / RAL7035 (light grey)
EC directive	93/42/EEC
Device type	Medical device class I/ supplement
Labelling	CE
GMDN Code	35239

Logistics Information

Manufacturer	BIOPTRON AG, Sihleggstrasse 23, 8832 Wollerau; Switzerland
Customer	Home Art & Sales Services; Sihleggstrasse 23, 8832 Wollerau, Switzerland
Net Weight kg	5.98
Gross Weight kg	6.9
Volume m ³	0.024
Cardbox dimension cm (l x w x h)	100 x 20 x 12
Customs Tariff Code	7326.9098

Date:

17.05.2016

 Anouar Hicham
 Head of RA & QA Departments
 BIOPTRON AG

 Tim Scherf
 General Manager
 BIOPTRON AG

BIOPTRON AG

Sihleggstrasse 23, CH-8832 Wollerau
 Phone: +41 43 888 28 00 Fax: +41 43 888 28 99
 INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

Bank: UBS AG, CH-6900 Lugano
 SWIFT: UBSWCHZH80A, BC 247

IBAN CHF: CH94 0024 7247 5714 9301 B
 IBAN EUR: CH24 0024 7247 5714 9362 L
 IBAN USD: CH57 0024 7247 5714 9361 C

Certificate


BIOPTRON Table Stand Unit, BPro1

General

Customer item Number PAG-990-TS

Technical Information

Item Number	001461
Device description	Accessory to BPro1 Device (PAG-990)
Material	Thermoplastic Cyclooy
Weight Net kg.	1.20 (without device)
Dimensions, cm (l x w x h)	38,9 x 25,7 x 27,7
Range of adjustment of the stand	Not adjustable
Color	White
EC directive	93/42/EEC
Device type	Medical device class I/ supplement
Labelling	CE
GMDN Code	35239

Logistics Information

Manufacturer	BIOPTRON AG, Sihleggstrasse 23, 8832 Wollerau; Switzerland
Customer	Home Art & Sales Services; Sihleggstrasse 23, 8832 Wollerau, Switzerland
Net Weight kg	1.2
Gross Weight kg	2.1
Volume m ³	0.112
Cardbox dimension cm (l x w x h)	60 x 54 x 36

Date:

17.05.2016

Anouar Hicham
 Head of RA & QA Departments
 BIOPTRON AG

Tim Scheifele
 General Manager
 BIOPTRON AG

BIOPTRON AG

Sihleggstrasse 23, CH-8832 Wollerau
 Phone: +41 43 999 28 00 Fax: +41 43 888 28 99
 INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

Bank: UBS AG, CH-6900 Lugano
 SWIFT: UBSWCHZH80A, BC 247

IBAN CHF: CH94 0024 7247 5714 9301 B
 IBAN EUR: CH24 0024 7247 5714 9362 L
 IBAN USD: CH57 0024 7247 5714 9361 C


Item 102939


Content

102939
PAG-990-TBAG Transport Bag BPro1

Technical Information

Applied Quality Standard	IEC 60529:1989/A1:1999 EN 60529:1991/AC:1993
Handle	25mm PP Ribbon
Handle color	White
Technical Description	IP Bag
Weight Net	-
Outer dimensions (w x h x d)	1000mm x 720mm
Material	0.2mm Clear PVC
Color	Clear/transparent


Labeling / Positioning	
------------------------	---

Size: 10 cm x 13 cm

102939

BIOPTRON
LIGHT THERAPY SYSTEM
www.bioptron.com



Folded dimensions	160mm x 260mm (lxb)
-------------------	---------------------

Logistic Information

Country of Origin CN
Customs Tariff Code -

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
Logistic Information

Country of Origin	CN
Customs Tariff Code	

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 TSM SUCCESS MANUAL® designed and copyright by ergo use SWISS

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Anouar Hicham
Head of RA & QA Departments
BIOPTRON AG

Tim Scheifele
General Manager
BIOPTRON AG

MedAll

Table 1: Emission

Guidance and manufacturer's declaration – electromagnetic emissions		
The light therapy device MedAll is intended for use in the electromagnetic environment specified below. The customer or the user of the light therapy device MedAll should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The light therapy device MedAll uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The light therapy device MedAll is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class C (and Class A)	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

Table 2: Immunity (all devices)

Guidance and manufacturer's declaration – electromagnetic immunity			
The light therapy device MedAll is intended for use in the electromagnetic environment specified below. The customer or the user of the light therapy device MedAll should assure that it is used in such an environment.			
Immunity test standard	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±8 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient / burst IEC 61000-4-4	±2 kV for power supply lines	±4 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line to line ±2 kV line to earth	±1 kV line to line	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply lines IEC 61000-4-11	<5% U_n (0,5 cycle) 40 % U_n (5 cycles) 70 % U_n (25 cycles) <5% U_n for 5s	<5% U_n (0,5 cycle) 40 % U_n (5 cycles) 70 % U_n (25 cycles) <5% U_n for 5s	Mains power quality should be that of a typical commercial or hospital environment.
Note: U_n is the a.c. mains voltage prior to application of the test level.			
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	200 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Table 3: Immunity (not life-supporting devices)

Guidance and manufacturer's declaration – electromagnetic immunity			
The light therapy device MedAll is intended for use in the electromagnetic environment specified below. The customer or the user of the light therapy device MedAll should assure that it is used in such an environment.			
Electromagnetic environment – guidance			
Portable and mobile RF communications equipment should be used no closer to any part of the light therapy device MedAll, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.			
Immunity test standard	IEC 60601 test level	Compliance level	Recommended separation distance ^a
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	V1 = 10 Vrms (estimated value) 150 kHz to 80 MHz	$d = 0.35 \sqrt{P}$ 150 kHz to 80 MHz
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 800 MHz	E1 = 10 V/m 80 MHz to 800 MHz	$d = 0.35 \sqrt{P}$ 80 MHz to 800 MHz
Radiated RF IEC 61000-4-3	3 V/m 0.8 to 2.5 GHz	E2 = 10 V/m 800 MHz to 2.7 GHz	$d = 0.7 \sqrt{P}$ 0.8 to 2.7 GHz
Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 			
Note 1: At 80MHz and 800MHz, the higher frequency range applies. Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people			
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the MedAll is used exceeds the applicable RF compliance level above, the light therapy device MedAll should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the light therapy device MedAll. ^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m. ^c Possible shorter distances outside ISM bands are not considered to have a better applicability of this table.			

Table 4: Recommended separation distances (not life-supporting devices)

Recommended separation distances between portable and mobile RF communications equipment and the light therapy device MedAll			
The light therapy device MedAll is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the light therapy device MedAll can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the light therapy device MedAll as recommended below, according to the maximum output power of the communication equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150kHz to 80MHz $d = 0.35 \sqrt{P}$	80MHz to 800MHz $d = 0.35 \sqrt{P}$	800MHz to 2500MHz $d = 0.7 \sqrt{P}$
0.01	0.04 m	0.04 m	0.07 m
0.1	0.12 m	0.12 m	0.22 m
1	0.35 m	0.35 m	0.7 m
10	1.11 m	1.11 m	2.2 m
100	3.5 m	3.5 m	7 m
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
Note 1: At 80MHz and 800MHz, the separation distance for the higher frequency range applies.			
Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
Note 3: An additional factor of 10/3 is used in calculating the recommended separation distance to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas.			

BIOPTRON Pro 1 (BPro1)

Table 1: Emission

Guidance and manufacturer's declaration – electromagnetic emissions		
The therapy lamp BPro1 is intended for use in the electromagnetic environment specified below. The customer or the user of the therapy lamp BPro1 should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The therapy lamp BPro1 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The therapy lamp BPro1 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class C (and Class A)	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

Table 2: Immunity (all devices)

Guidance and manufacturer's declaration – electromagnetic immunity			
The therapy lamp BPro1 is intended for use in the electromagnetic environment specified below. The customer or the user of the therapy lamp BPro1 should assure that it is used in such an environment.			
Immunity test standard	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient / burst IEC 61000-4-4	±2 kV for power supply lines	±4 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line to line ±2 kV line to earth	±1.5 kV line to line	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply lines IEC 61000-4-11	<5% U_r (0,5 cycle) 40 % U_r (5 cycles) 70 % U_r (25 cycles) <5% U_r for 5s	<5% U_r (0,5 cycle) 40 % U_r (5 cycles) 70 % U_r (25 cycles) <5% U_r for 5s	Mains power quality should be that of a typical commercial or hospital environment.
Note: U_r is the a.c. mains voltage prior to application of the test level.			
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	200 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Table 3: Immunity (not life-supporting devices)

Guidance and manufacturer's declaration – electromagnetic immunity			
The therapy lamp BPro1 is intended for use in the electromagnetic environment specified below. The customer or the user of the therapy lamp BPro1 should assure that it is used in such an environment.			
Electromagnetic environment – guidance			
Portable and mobile RF communications equipment should be used no closer to any part of the therapy lamp BPro1, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.			
Immunity test standard	IEC 60601 test level	Compliance level	Recommended separation distance ^c
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	V1 = 10 Vrms (estimated value) 150 kHz to 80 MHz	$d = 0.35 \sqrt{P}$ 150 kHz to 80 MHz
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 800 MHz	E1 = 10 V/m 80 MHz to 800 MHz	$d = 0.35 \sqrt{P}$ 80 MHz to 800 MHz
Radiated RF IEC 61000-4-3	3 V/m 0.8 to 2.5 GHz	E2 = 10 V/m 800 MHz to 2.7 GHz	$d = 0.7 \sqrt{P}$ 0.8 GHz to 2.7 GHz
where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).			
Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b			
Interference may occur in the vicinity of equipment marked with the following symbol: 			
Note 1: At 80MHz and 800MHz, the higher frequency range applies.			
Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people			
a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the therapy lamp BPro1 is used exceeds the applicable RF compliance level above, the therapy lamp BPro1 or the device, which contains it should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device containing the therapy lamp BPro1.			
b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.			
c Possible shorter distances outside ISM bands are not considered to have a better applicability of this table.			

Table 4: Recommended separation distances (not life-supporting devices)

Recommended separation distances between portable and mobile RF communications equipment and the therapy lamp BPro1			
The therapy lamp BPro1 is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the therapy lamp BPro1 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the therapy lamp BPro1 as recommended below, according to the maximum output power of the communication equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 0.35 \sqrt{P}$	80 MHz to 800 MHz $d = 0.35 \sqrt{P}$	800 MHz to 2500 MHz $d = 0.7 \sqrt{P}$
0,01	0.04 m	0.04 m	0.07 m
0,1	0.12 m	0.12 m	0.22 m
1	0.35 m	0.35 m	0.7 m
10	1.11 m	1.11 m	2.2 m
100	3.5 m	3.5 m	7.0 m
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
Note 3: An additional factor of 10/3 is used in calculating the recommended separation distance to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas.			

BIOPTRON 2 (B2)

Table 1: Emission

Guidance and manufacturer's declaration – electromagnetic emissions		
The light therapy device B2 is intended for use in the electromagnetic environment specified below. The customer or the user of the light therapy device B2 should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The light therapy device B2 uses RF energy only for its internal function. Therefore, its RF emission is very low and is not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The light therapy device B2 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

Table 2: Immunity (all devices)

Guidance and manufacturer's declaration – electromagnetic immunity			
The light therapy device B2 is intended for use in the electromagnetic environment specified below. The customer or the user of the light therapy device B2 should assure that it is used in such an environment.			
Immunity test standard	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient / burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for signal lines	±2 kV for power supply lines ±2 kV for signal lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line to line ±2 kV line to earth	±1 kV line to line ±2 kV line to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply lines IEC 61000-4-11	<5% U_r (0,5 cycle) 40 % U_r (5 cycles) 70 % U_r (25 cycles) <5% U_r for 5s	<5% U_r (0,5 cycle) 40 % U_r (5 cycles) 70 % U_r (25 cycles) <5% U_r for 5s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the light therapy device B2 requires continued operation during power mains interruptions, it is recommended that the light therapy device B2 is powered from an UPS.
Note: U_r is the a.c. mains voltage prior to application of the test level.			
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Table 3: Immunity (not life-supporting devices)

Guidance and manufacturer's declaration – electromagnetic immunity			
The light therapy device B2 is intended for use in the electromagnetic environment specified below. The customer or the user of the light therapy device B2 should assure that it is used in such an environment.			
Electromagnetic environment – guidance			
Portable and mobile RF communications equipment should be used no closer to any part of the light therapy device B2, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.			
Immunity test standard	IEC 60601 test level	Compliance level	Recommended separation distance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	V1 = 10 Vrms 150 kHz to 80 MHz	$d = 0.35 \sqrt{P}$ 150 kHz to 80 MHz
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 800 MHz	E1 = 10 V/m 80 MHz to 800 MHz	$d = 0.35 \sqrt{P}$ 80 MHz to 800 MHz
Radiated RF IEC 61000-4-3	3 V/m 800 MHz to 2.5 GHz	E2 = 7 V/m 800 MHz to 2.7 GHz	$d = \sqrt{P}$ 800 MHz to 2.7 GHz
where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).			
Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b			
Interference may occur in the vicinity of equipment marked with the following symbol: 			
Note 1: At 80MHz and 800MHz, the higher frequency range applies.			
Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the light therapy device B2 is used exceeds the applicable RF compliance level above, the light therapy device B2 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the light therapy device B2.			
^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.			

Table 4: Recommended separation distances

Recommended separation distances between portable and mobile RF communications equipment and the light therapy device B2			
The light therapy device B2 is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the light therapy device B2 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the light therapy device B2 as recommended below, according to the maximum output power of the communication equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150kHz to 80MHz $d = 0.35 \sqrt{P}$	80MHz to 800MHz $d = 0.35 \sqrt{P}$	800MHz to 2700MHz $d = \sqrt{P}$
0.01	3.5 cm	3.5 cm	10 cm
0.1	7.8 cm	7.8 cm	32 cm
1	35 cm	35 cm	1 m
10	78 cm	78 cm	3.2 m
100	3.5 m	3.5 m	10 m
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
Note 1: At 80MHz and 800MHz, the separation distance for the higher frequency range applies.			
Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

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GUARANTEE CERTIFICATE

EN

This BIOPTRON product is manufactured with the latest production methods and quality components to ensure trouble-free performance.

The Guarantee

With this Consumer Guarantee, BIOPTRON AG guarantees this product to be free of defects in materials and workmanship at the time of its original purchase or at the time it was taken on hire purchase terms by the consumer from the retailer for the period of ONE year. If during this period of guarantee the product proves defective due to improper materials or workmanship, the Retailer, Authorized Service Agents or Authorized Servicing Dealers will, without charge for labour or parts, repair or replace this product or its defective parts on the conditions explained hereafter.

Conditions

- I. The guarantee is valid when presented together with the original invoice/receipt issued to the consumer by the retailer, and if the guarantee card states a) the purchaser's name, b) the retailer's name and address, c) the model name and serial number of the purchased product and d) the date of purchase of the product.
BIOPTRON AG reserves the right to refuse guarantee service if this information has been removed or changed after the original purchase of the product by the consumer from the retailer.
- II. This product will not be considered defective in materials or workmanship if it needs to be adapted, changed or adjusted in order to conform to the national or local technical or safety standards in any country other than the one which the product was originally designed and manufactured for.
This guarantee does not cover reimbursement for a) such adaptations, changes or adjustments or attempts to do so, whether properly performed or not, or b) any damage resulting thereof.
- III. This guarantee covers none of the following:
 - A. Periodic check-ups, maintenance and repair or replacements of parts due to normal wear and tear;
 - B. Home service transport costs and other costs and risks of transport relating directly or indirectly to the guarantee of this product.
 - C. Damage to this product resulting from:
 1. Abuse and misuse, including but not limited to:
 - a) the use of this product for other than its normal purposes or not in accordance with BIOPTRON's instructions for the proper use and maintenance of this product; and
 - b) the installation of this product in a manner inconsistent with the technical or safety standards in the country where the product is used;
 2. Repairs done by unauthorized service stations;
 3. Accidents, acts of God or any cause beyond the control of BIOPTRON AG, including but not limited to lightning, water, fire, public disturbances and improper ventilation.
- IV. This guarantee only covers BIOPTRON units which were imported and sold by an authorized BIOPTRON retailer. The guarantee claim must be made via your retailer.
- V. This guarantee does not affect the consumer's statutory rights under applicable national legislation in force, nor the consumer's rights against the retailer arising from the sales/purchase contract.
In the absence of applicable national legislation, this guarantee will be the consumer's sole and exclusive legal remedy, and neither BIOPTRON AG nor its Subsidiaries or Distributors shall be liable for any incidental or consequential damages for breach of any express or implied guarantee of this product.

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BIOPTRON 
LIGHT THERAPY SYSTEM By Zepher Group

© COPYRIGHT BIOPTRON AG - ZEPHER GROUP - Released on 2014-02 - Code 302087 Rev03

Legalization of signature

The authorized person of legalization of the canton of Schwyz, Switzerland, certifies, that Mr. Anouar Hicham, born on 21.04.1969, from Germany, in Herbrechtingen Germany, identified by Identity Card no. LFGHJHY7V of Germany, issued on 10.12.2010, and Mr. Timothy David Scheifele, born on 12. Mai 1963, from Zürich ZH, in Baar ZG, identified by Identity Card no. E1664844 of Switzerland, issued on 16.10.2013, have acknowledged the foregoing signatures as their own signatures. Mr. Anouar Hicham and Mr. Timothy David Scheifele are known personally to the authorized person of legalization as well as their signatures.

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© COPYRIGHT BIOPTRON AG – ZEPTER GROUP – Released on 2014-03 – Code 303704-RU
© COPYRIGHT BIOPTRON AG - ZEPTER GROUP - Published on 2013-03 - Code: 303202 Rev.04