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CALIFORNIA COPY CERTIFICATION BY DOCUMENT CUSTODIAN

I, PRIYA NAVKAR, hereby swear (or affirm) that the attached reproduction of Alligator IFU is a true, correct and complete photocopy of a document in my possession.

Description of Original Document

Priya Navkar

Signature of Custodian of Original Document

9775 Toledo Way Irvine CA 92618

Kansas

State of California

County of Orange

Subscribed and sworn to (or affirmed) before me on this 18 day of May, 2015, by

Date Month Year

PRIYA NAVKAR

Name of Custodian of Original Document

proved to me on the basis of satisfactory evidence to be the person(s) who appeared before me.

Signature

Joanne Charles

Signature of Notary Public

Please Notary Seal Above

OPTIONAL

Though the information in this section is not required by law, it may prove valuable to persons relying on the document and could prevent fraudulent removal and reattachment of this form to another document.

Description of Attached Document Copy

Title or Type of Document: Notarized for Russia

Document Date: _____ Identifying No.: _____ No. of Pages: _____

Signer(s) or Issuing Agency: _____

Capacity Claimed by Custodian

☐ Individual ☐ Attorney ☐ Trustee ☐ Business Proprietor or Manager

☐ Corporate Officer — Title: _____

☐ University or School Officer — Title: _____

☐ Governmental Officer or Agent — Title: _____

☐ Other: _____

Custodian Is Representing: _____

Navan

INSTRUCTIONS FOR USE

AlligatorTM Retrieval Device

70730-001 Rev. 02/13

This page is intentionally left blank.

TABLE OF CONTENTS

Alligator™ Retrieval Device

1. Introduction	1
2. Safety	2
3. Components	3
4. Assembly	4
5. Operation	5
6. Maintenance	6
7. Troubleshooting	7
8. Specifications	8
9. Warranty	9
10. Contact Information	10
11. Appendix A: Diagrams	11
12. Appendix B: Parts List	12
13. Appendix C: Safety Instructions	13
14. Appendix D: User Manual	14
15. Appendix E: Technical Specifications	15
16. Appendix F: Troubleshooting Guide	16
17. Appendix G: Maintenance Schedule	17
18. Appendix H: Warranty Claim Form	18
19. Appendix I: Contact Information	19
20. Appendix J: Glossary	20
21. Appendix K: Index	21
22. Appendix L: Revision History	22
23. Appendix M: Change Log	23
24. Appendix N: Revision History	24
25. Appendix O: Change Log	25
26. Appendix P: Revision History	26
27. Appendix Q: Change Log	27
28. Appendix R: Revision History	28
29. Appendix S: Change Log	29
30. Appendix T: Revision History	30
31. Appendix U: Change Log	31
32. Appendix V: Revision History	32
33. Appendix W: Change Log	33
34. Appendix X: Revision History	34
35. Appendix Y: Change Log	35
36. Appendix Z: Revision History	36

TABLE OF CONTENTS

Alligator™ Retrieval Device

English	4
Français.	5
Deutsch.	6
Italiano	7
Español	8
Svenska.	9
Nederlands.	10
Português	11
Suomi.	12
Dansk	13
Ελληνικά	14
Česky	15
Magyar	16
Русский.	17
Polski	18
Türkçe.	19
Norsk	20
Slovenčina	21
Română.	22
Български	23
한국어.	24
عربي.	25
Symbol Glossary	26

English Instructions for Use

EN

Alligator™ Retrieval Device

CAUTION


Federal (USA) law restricts this device to sale, distribution, and use by or on the order of a physician. This device should be used only by physicians with a thorough understanding of angiography and/or percutaneous neurointerventional procedures.

This device should be used only by physicians trained in percutaneous, intravascular techniques and procedures at medical facilities with the appropriate fluoroscopic equipment.

CONTENTS

One Alligator™ Retrieval Device (ARD), one introducer and one device handle, as an assembly.

Figure 1

DESCRIPTION

The Alligator™ retrieval device (ARD) is a retriever with grasping jaws attached to the tip of a flexible wire. This device is designed to be used in conjunction with commercially available 3F microcatheter with minimum ID (inside diameter) of 0.021 in. (0.53mm). The grasping jaws and the distal end of the ARD are made of radio opaque alloy facilitating fluoroscope visualization.

INDICATIONS FOR USE

This device is intended for use in the peripheral and neurovasculature for foreign body retrieval.

CONTRAINDICATIONS

This device is not intended for:

- Removal of foreign objects entrapped by tissue growth
- Removal of implanted pacing leads

POTENTIAL COMPLICATIONS

Potential complications include, but are not limited to, the following:

- Vessel perforation or dissection
- Vasoospasm
- Thromboemboli
- Neurological deficits including stroke and death
- Hemorrhage

DIRECTIONS FOR USE

1. Using standard interventional techniques, position the tip of the microcatheter proximal to the object to be retrieved.
2. Select the appropriate size Alligator™ retrieval device based on vessel diameter. In general, the labeled diameter of the Alligator™ retrieval device should approximate that of the vessel in which it will be used.
3. Secure the introducer and the device handle and gently withdraw the ARD from packaging tube until the grasping jaws are inside the introducer. The device handle may now be removed and discarded.
4. Place the tip of introducer into the luer hub of the microcatheter and gently advance the Alligator™ retrieval device into the microcatheter. Make sure that the tip of the introducer is held firmly and deeply in the hub of the microcatheter before advancing the Alligator™ retrieval device.
5. Remove introducer after the ARD is inserted into the microcatheter by pulling it back over the wire.
6. Continuous heparinized saline flush during use of the ARD is recommended.
7. Under fluoroscopic guidance, advance the ARD grasping jaws to the tip of the microcatheter. Do not advance the grasping jaws beyond the tip of the microcatheter.
8. Open the jaws of the Alligator™ retrieval device by holding the ARD wire in place and gently withdrawing the microcatheter. This will unsheath the grasping jaws.
9. Advance the Alligator™ retrieval device gently in the vessel under fluoroscopic guidance until the tips of the jaws are past the proximal edge of the foreign object.
10. Holding the ARD wire still, gently advance the microcatheter over the wire until the radiopaque marker at the tip of the microcatheter passes the base of the grasping jaws.
11. Maintain tension on the grasping jaws (and retain the hold on the foreign body) by gently pulling the proximal end of the ARD. Do not pull so hard that the jaws are completely drawn into the microcatheter or the ARD will release the foreign object. The ARD cannot be pulled back entirely into the microcatheter when object is grasped.
12. Gently retract the microcatheter and ARD at the same time. If the foreign object is not well visualized on fluoroscopy, then a gentle contrast injection through the guiding catheter may help opacify the vessel and confirm that the ARD has grasped the foreign object.

- Avoid rotating the ARD during use. Rotating the ARD, particularly when it is sheathed within the microcatheter, may twist and entangle the ARD grasping jaws.
- When advancing an open (unsheathed) ARD across a vessel bifurcation, confirm that all jaws remain in the same vessel before closing the jaws with the microcatheter.
- Do not use the ARD as an access wire. Instead, gain access with a regular guidewire and microcatheter, and then unsheath the ARD close to the foreign body.
- This device is supplied STERILE for single use only. Do not reprocess or re-sterilize. Reprocessing and re-sterilization increase the risks of patient infection and compromised device performance.

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CALIFORNIA ALL-PURPOSE ACKNOWLEDGMENT

CIVIL CODE § 1189

A notary public or other officer completing this certificate verifies only the identity of the individual who signed the document to which this certificate is attached, and not the truthfulness, accuracy, or validity of that document.

State of California

County of Orange

On 19th May 2015 before me,

Joanne Charles, Notary Public

Here Insert Name and Title of the Officer

personally appeared


Priya Nankar

Name(s) of Signer(s)

who proved to me on the basis of satisfactory evidence to be the person(s) whose name(s) is/are subscribed to the within instrument and acknowledged to me that he/she/they executed the same in his/her/their authorized capacity(ies), and that by his/her/their signature(s) on the instrument the person(s), or the entity upon behalf of which the person(s) acted, executed the instrument.

I certify under PENALTY OF PERJURY under the laws of the State of California that the foregoing paragraph is true and correct.

WITNESS my hand and official seal.

Signature

Joanne Charles
Signature of Notary Public

Place Notary Seal Above

OPTIONAL

Though this section is optional, completing this information can deter alteration of the document or fraudulent reattachment of this form to an unintended document.

Description of Attached Document

Title or Type of Document: Alligator Notarized - RV Document Date: _____

Number of Pages: _____ Signer(s) Other Than Named Above: _____

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Micro Therapeutics, Inc.
Dbc Ev3 Neurovascular, 9775 Toledo Way,
Irvine, CA 92618, USA
Regulatory Affairs Product Associate

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