
Confirm

Christoph Bosshard

CEO

Oertli Instrumente AG

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
VC821100

CAUTION:

U.S. Federal Law restricts this device to sale by or on the order of a physician.

Document no.: VV016039E

Version: 01

Publication date: 14.02.2014

Original service manual

© Oertli Instrumente AG. All rights reserved.

Oertli, the Oertli logo, CataRhex, the CataRhex logo and ParaProg are registered trademarks of Oertli Instrumente AG.

This service manual shall not be multiplied in part or in full by electronic, (photo) mechanical means or in any other way, nor used without authorization for competitive purposes without the consent of Oertli Instrumente AG.

Subject to modifications.

Contact and service:

Oertli Instrumente AG

Hafnerwisenstrasse 4

CH – 9442 Berneck, Switzerland

Tel: +41 (0)71 747 42 00

Fax: +41 (0)71 747 42 90

E-mail: info@oertli-instruments.com

Website: www.oertli-instruments.com

Contents

1	About this document.....	4
1.1	Warnings and symbols in this manual.....	4
2	Safety instructions.....	4
3	Block diagram CataRhex 3.....	5
4	Installation.....	5
4.1	General remarks.....	5
4.2	Power supply.....	5
5	Device composition.....	6
5.1	Block diagram control unit.....	6
5.2	Logic PCB.....	7
5.3	Power PCB.....	8
5.4	Pedal PCB.....	9
6	Service mode operation.....	9
6.1	General information regarding service mode.....	9
6.2	Display software (SW) version and serial number.....	10
6.3	Setting default values.....	10
6.4	Save and import surgeon data.....	10
6.5	Checking and calibrating the pedal.....	11
6.6	Calibration of peristaltic sensor.....	12
6.7	Button and sound test.....	12
6.8	Operating voltages.....	12
6.9	Temperature.....	12
6.10	Installation of HFDS (additional functions).....	13
6.11	VIT-PN corrective (only after replacement of Vit-PN module).....	13
7	Installation and uninstallation of components.....	14
7.1	Replacement of fuses.....	14
7.2	Replacement of control panel.....	14
7.3	Open housing.....	15
7.4	Replacement of peristaltic pump.....	18
7.5	Replacement of Vit-PN module.....	19
7.6	Replacement of coupling plate.....	21
7.7	Replacement of system board.....	22
7.8	Replacement of power supply unit.....	22
7.9	Replacement of clamp for I/V pole.....	23
8	Software update.....	24
8.1	Device.....	24
8.2	Pedal.....	25
9	Installing replacement parts of pedal.....	25
9.1	Replacement of cover plate.....	25
9.2	Replacement of pedal electronics board.....	25
9.3	Replacement of pedal cable.....	27
9.4	Replacement of pedal sensor.....	28
9.5	Replacement of rocker.....	29
9.6	Replacement of pedal handle.....	30
9.7	Replacement of rubber mat.....	30
9.8	Replacement of rubber buffer.....	31
9.9	Set of replacement screws.....	32
9.10	Set of small replacement parts.....	33
9.11	Thread locking.....	33
10	Order number for replacement parts.....	33
10.1	Control unit.....	33
10.2	Pedal.....	33
11	Messages, warnings and error messages / fault correction.....	34
12	Authorized service points.....	36
12.1	Sending defective device, instrument or part to Oertli Instrumente AG.....	36
13	Functional test.....	37

1 About this document

The service manual is supplied with the device. The service manual is a supplement to VV016038 CataRhex 3 instructions for use and describes safe and appropriate service and maintenance of the device.
The technical data of the unit can be found in the technical data chapter in the Instructions for Use.

1.1 Warnings and symbols in this manual

Warning sign/word	Danger level	Consequences of non-compliance
 DANGER!	Imminent danger	Death, serious injury
 WARNING!	Possible danger	Death, serious injury
 CAUTION!	Possible dangerous situation	Slight injury
NOTE!	Possible dangerous situation	Material damage

Symbol	Meaning
	Safety sign
	Electrostatic sensitive devices
	Information that makes working with the device easier

2 Safety instructions

DANGER!
Unplug mains (the power cord) before opening the unit!


CAUTION!
Repairs to the device may only be performed by service technicians authorized by Oertli and based on the last valid version of the service manual. Improper repairs compromise the safety of staff and patients and will void any warranty.

NOTICE
This unit contains **Electrostatic Sensitive Devices (ESD)**. When opening the unit, it is very important that you and the unit are grounded at all times.

3 Block diagram CataRhex 3

4 Installation

4.1 General remarks


When placing the device, ensure clearance of at least 20cm around the ventilation openings. Do not cover the device when in operation.

4.2 Power supply


Type	Value
Supply system	Voltage: 100 V-240 VAC Frequency: 50 Hz-60 Hz Output: 150 VA

5 Device composition

5.1 Block diagram control unit


5.2 Logic PCB


Number	Abbreviation	Function
JP1	-	Connection control unit
JP2	DIAGNOSE	RS232 interface (only for manufacturer)
JP3	CAN1	CAN connection pedal
JP4	OCDS1	Empty
JP5	CAN2	CAN reserve
JP6	USB PROGRAMM	Remains open (only for manufacturer)
JP7	FORCE SENSOR	
JP8	-	Connection control unit
JP9	RESERVE	Empty
JP10	USB	USB connection
JP11	ENCODER	Connection encoder of peristaltic pump
JP12	VACTEST	Remains open (only for manufacturer)
JP13	SWITCH	Connection of tubing set detection switch

5.3 Power PCB

Number	Abbreviation	Function
JP1	IRR VALVE	Connection of irrigation valve
JP2	PHACO	Connection of phaco handpiece
JP3	POWER	Power supply for power supply unit
JP4	PUMP	Connection of peristaltic pump
JP6	DIA	Connection of diathermy handpiece
JP7	VIT PN	Connection of Vit-PN module
JP9	VIT	Remains open

5.4 Pedal PCB

Number	Abbreviation	Function
JP1	SUPPLY	CAN, +24V of the CataRhex 3
JP4	Diagnose	RS232 interface: remains open (only for manufacturer)
JP5	Sensor	Vertical angle sensor
JP6	Switch	Left and right directional switch

6 Service mode operation

6.1 General information regarding service mode

After starting (RESET or network ON) you access the service mode by simultaneously pressing the AUX and VIT keys for at least three seconds.

To navigate, use the buttons below the display.

Lower line with symbols

These can be used to move the cursor.

- Direction arrows (▲) (▼) Use these direction arrows, to jump to the next or previous menu and to move the cursor up or down.
- Direction arrows (◀) (▶) Use these direction arrows to move the cursor left or right.
- Enter button (↵) Use the enter button, to confirm a selected function.

The function selected with the cursor will flash.

To exit the service mode, press the RESET button at any time.


The service mode is divided into the following different menus:

SURG-MEM INTERNAL	- Resetting of surgeon data for one or all surgeon memories. - Internal copying of surgeon data.
SURG-MEM USB	- Saving and retrieving data with USB stick. - Copying of surgeon data from USB stick.
PEDAL	Calibration of the pedal.
SENSOR SET	Calibration of the force sensor.
BUTTONS	Testing the buttons on the control panel and sound system.
VOLTAGES	Display of current device voltages.
TEMPERATURE	Temperature value of the sensor in the device.
ADDITIONAL FUNCTIONS	Generating request code, installing and uninstalling additional functions.

6.2 Display software (SW) version and serial number

The SW version and the device serial number will be displayed after entering the service mode.

Display:

SNXXYYNNNN Serial number
SWCXX.YY SW version application programme
BOOTYYY SW version boot programme

6.3 Setting default values

All the limit values (stored values), pedal control and ParaProg settings can be reset to the values recommended by the manufacturer in the SURG-MEM INTERNAL menu.

- In the SURG-MEM INTERNAL menu, select the surgeon storage location first or choose ALL for all surgeons by keeping the AUX key depressed while pressing the keys below "SELECT SURGEON".
- Confirm "SET ALL PARAMETER DEFAULT?" by pressing the Enter button (↵).
- Confirm "CONFIRM DEFAULT?" by pressing the Enter button (↵).

6.4 Save and import surgeon data

It is possible to store all surgeon and system parameters on a USB stick in the SURG-MEM USB menu. The USB stick must be formatted as described in the section "Software update / device".

Possibilities:

"WRITE ALL TO USB"	All data (settings, ParaProg, pedal settings, system parameters) from all surgeons will be stored on the USB stick.
"READ ALL FROM USB"	All surgeon data (settings, ParaProg, pedal settings) of all surgeons will be restored from the USB stick (with the surgeon's name on the USB stick).
"READ SURGEON FROM USB"	All surgeon data (settings, ParaProg, pedal settings) of only one surgeon will be copied from the USB stick to a specified surgeon memory in the CataRhex 3 device. The surgeon name in CataRhex 3 will not be changed. It is very important that the software version in the destination unit is equal to or higher than the software version in the unit containing the source data!

Possible messages on the display	Error / action
Insert memory stick	Insert USB stick. The USB stick must be formatted as described in the section "Software update / device".
File not found	Surgeon data file with the name CR3S_DOC.BIN is not on the USB stick.
Copy data to USB (.)	Wait until data are written.
Copy data from USB (.)	Wait until data are read.
Memory stick removed; data saving incomplete	Do not remove the USB stick while writing.
Memory stick removed; data load incomplete	Do not remove the USB stick while reading.  CAUTION! Surgeon data on the device may be inconsistent! Repeat surgeon data read.
File on memory stick is corrupted.	Surgeon data file is corrupted. Write the data from the source device to the USB stick again.
Data saved; remove memory stick.	Data transfer procedure was successful. Remove the USB stick!
Data read; remove memory stick.	Data access procedure was successful. Remove the USB stick!
USB ERROR	An error occurred. Please repeat the surgeon data access or transfer procedure.

In the menu SURG-MEM INTERNAL with the function "COPY SURGEON MEMORY" it is possible to copy all surgeon settings (settings, ParaProg, pedal settings) from one surgeon memory to another selectable memory position. The surgeon name from the destination will not be changed.

6.5 Checking and calibrating the pedal

The pedal settings can be checked and re-calibrated in the PEDAL menu. As part of this process, the mechanical positions are synchronised with the sensor values. These values are then stored in the pedal. This means that a calibrated pedal can also be operated with a different CataRhex 3 unit.

The CataRhex 3 unit will automatically detect whether another type of pedal is connected. SwissTech pedals can be used with CataRhex 3 as well. OS3 and faros pedals cannot be used with CataRhex 3.

Display:

SERVICE MODE	PEDAL
VERTICAL POS 0	0% NoSwitch HW:001 SW:007
SET PEDAL?	YES NO
▼ ▲	◀ ▶

Pedal position (e.g. POS 0)

- 0 = Neutral
- 1 = IRR
- 2 = Aspiration
- 3 = Instrument active
- 4 = Reflux

0%, vertical deflection

0 ... 100% (per pedal position)

NoSwitch, horizontal deflection

NoSwitch = horizontal central
LEFT = horizontal left
RIGHT = horizontal right

Pedal calibration:

- Move the cursor to "SET PEDAL?"
- Move the cursor to "REALLY CHANGE SETTINGS?"
- Follow the instructions displayed.

Select "YES" and press the Enter button (↵).
Select "YES" and press the Enter button (↵).

When entering settings, ensure that the 100% end value can be read off for the horizontal and vertical deflection in each area (otherwise repeat the calibration).

6.6 Calibration of peristaltic sensor

The vacuum for the peristaltic pump is measured with the force sensor.

The offset adjustment for the force sensor can be made in the SENSOR SET menu.

Display:

SERVICE MODE	SENSOR SET
ACTUAL OFFSET: 0333	MEASURED OFFSET: 0338
CHANGE OFFSET:	YES NO
▼ ▲	◀ ▶

ACTUAL OFFSET: **0333**
MEASURED OFFSET: **0338**

Saved offset value (latest adjustment)
Measured force sensor offset value

Sensor calibration:

- Move the cursor to "CHANGE OFFSET?"
- Remove the tubing system.
- Press the Enter button (↵) for "AFTERWARDSSTORE OFFSET:"

Select "YES" and press the Enter button (↵).
Select "STORE"

CAUTION!

A poorly adjusted sensor can lead to poor chamber stability!

6.7 Button and sound test

All buttons on the control panel, the tube switches and the loudspeaker can be tested in the "BUTTONS / SOUND" menu.

- When actuating the buttons, the following applies:
 - If a button is pressed, the name of the button will appear on the display.
 - If no button or more than one button is pressed, "NO OR MANY KEYS" is shown on the display.
- When the button below "SOUND" is actuated, a clear, continuous, loud acoustic signal can be heard as long as the button is pressed.
- To exit the "BUTTONS / SOUND" menu, press the Enter button (↵)
- When inserting the tubing system, the following applies:
 - If the tube system is correctly inserted, "BOTH TUBE SWITCHES" appears on the display.
 - If the tube system is only partially inserted, "TUBE SWITCH 1" or "TUBE SWITCH 2" appears on the display.
 - If the tube system is not inserted, "NO SWITCH" appears on the display.

6.8 Operating voltages

The operating voltages can be checked in the "VOLTAGES" menu.

The following applies:

+ 3.3V ±0.2V
+ 5.0V ±0.3V
+12.0V ±0.5V
- 12.0V ±0.5V
+24.0V ±1.0V

6.9 Temperature

The temperature in the device can be controlled in the "TEMPERATURE" menu.

The temperature sensor is located on the logic print.

The temperature in the device must not exceed 55°C.

6.10 Installation of HFDS (additional functions)

In the "ADDITIONAL FUNCTIONS" menu, request codes for the requirement of installation keys for additional functions can be generated. Moreover, the desired additional functions can be released using the installation key in this menu. A function that has been released can also be uninstalled again.

Display:

SERVICE MODE	ADDITIONAL FUNCTIONS
GENERATE REQUEST CODE	INSTALL FUNCTIONS
UNINSTALL FUNCTIONS	
▼ ▲	◀ ▶

Generate a request code:

- Move the cursor to "GENERATE REQUEST CODE" and press the Enter button (↵).
- Select the desired function and press the button below "GENERATE".
- The request code will appear on the display.
- Insert a USB stick. The USB stick must be formatted as described in the section "Software update / device".
- Press the button underneath "TO USB". The request code will be stored on the USB stick.
- Wait until "Data saved; remove memory stick" appears on the display.
- The request code is now saved on the USB stick in the file "CR3SCODE.BIN".
- Send the request code to Oertli.

Installation of additional functions:

- Move the cursor to "INSTALL FUNCTIONS" and press the Enter button (↵).
- If the installation key received from Oertli is in the file named CR3SK_XXYYNNNN, where XXYYNNNN is the serial number of the device, use a USB stick with this file. The USB stick must be formatted as described in the section "Software update / device".
- Move the cursor to "INSTALL FUNCTIONS" and press the Enter button (↵).
- The installation key will be read from the USB stick and appear on the display.
- Press the button below "INSTALL". The installation key will be checked and, provided it is valid, installed. Otherwise "INVALID" will appear on the display.

Uninstallation of additional functions:

- Move the cursor to "UNINSTALL FUNCTIONS" and press the Enter button (↵).
- Select the desired function and press the button underneath "REMOVE".
- Confirm uninstallation by pressing the button below "CONFIRM".

6.11 VIT-PN corrective (only after replacement of Vit-PN module)

The correction factor for the built-in VIT-PN unit can be set in the "VIT-DC CALIBRATION" menu. To access this menu, keep the VIT key depressed after start-up.

Display:

SERVICE MODE	VIT-DC CALIBRATION
VIT-DC CORRECTION 120%	
CHANGE CORRECTION?	YES NO
▼ ▲	◀ ▶ ↵ EXIT

Change correction:

- Move the cursor to "CHANGE CORRECTION?"
- The Vit-PN pump will start.
- Change the correction value with the left / right buttons.
- Move the cursor to "Finish?"
- The Vit-PN pump will stop.
- Press the button below "EXIT" to leave the menu.

Select "YES" and press the Enter button (↵).

Select "YES" and press the Enter button (↵).

7 Installation and uninstallation of components

NOTE!

The installation of components is done in the reverse order of their uninstallation. Each section specifically indicates any points which are important to observe.

NOTE!


A function test must be carried out on the device after every intervention (see Chapter 12).

7.1 Replacement of fuses

~~VX520010~~ Replacement fuses

Quantity	Description
10	Fuses T 3,15 AH, 250V

Work steps:

DANGER!

Take the electric plug out of the socket!

1. The fuse holder [1] can be pulled out by pressing the tabs [2] at the top and bottom of the fuse holder.
2. Replace defective fuses [3] with new ones of the same type. You will find the values indicated below the mains plug on the device.

DANGER!

Do not use any makeshift fuses!

DANGER!

Do not short-circuit the fuse holder!


3. Insert the fuse holder.

7.2 Replacement of control panel

~~VX210106~~ Control panel



Quantity	Description
1	Control panel
2	M4x8 countersunk screw, TORX

Work steps:


DANGER!
Take the electric plug out of the socket!

1. Unscrew the two screws on the upper cover plate and remove the upper cover plate.

2. Disconnect both flat ribbon cables from the logic print.
3. Pull out both flat ribbon cables through the front frame of the housing.

4. Fold up the control unit.
5. Undo the two countersunk head screws on the torque hinge and take off the control unit.
6. Carry out assembly of the new control unit in reverse order of dismantling.

No operating unit can be installed which has already been used with another control device.

7.3 Open housing


DANGER!
When assembling the supporting structure and front frame, the grounding connection between the supporting structure and front frame must be plugged in!

Quantity	Description
1	Upper cover plate with bottle / bag holder
2	M3x8 oval head screw, TORX, ecosyn-fix

Follow Step 1.

Quantity	Description
1	Under cover plate
3	M3x8 oval head screw, TORX, ecosyn-fix


Follow step 2.

DANGER!

Take the electric plug out of the socket!


1. Unscrew the two screws on the upper cover plate and remove the upper cover plate.


2. Unscrew the three screws on the lower cover plate and remove the lower cover plate.


3. Disconnect the "CAN" and "DIAGNOSE" plugs from the logic print.

4. Undo six screws on the supporting structure and carefully remove the supporting structure from the front frame.

5. Disconnect the cable from the power supply on the power print.


6. Unplug the grounding cable from the supporting structure of the front frame.

7.4 Replacement of peristaltic pump


Quantity	Description
1	Pump motor
4	M3x12 cylinder bolt

Quantity	Description
1	Pump wheel including cylinder and screw

Work steps:

1. Open housing (refer to Chapter 7.3).
2. Unplug the "PUMP" cable from the power print, then disconnect the "ENCODER" cable at the logic print.


3. Turn the pump wheel until the stud points to the side.
4. Undo stud with Allen wrench.


The stud is stuck and can only be undone by applying great force or heat!


5. Remove pump wheel.


6. Undo the four screws on the front frame and remove the pump motor.

7. Build in new pump motor paying attention to cable orientation.
8. Fix pump motor with four screws on the front frame.

9. Align motor axis so that the surface area points to the side.

10. Attach pump wheel with pump rollers to motor axis.

NOTE!

The distance between the pump rollers and the front frame must be set to 0.2 mm!

11. Attach pump wheel with stud and secure with Loctite 243.


7.5 Replacement of Vit-PN module

Quantity	Description
1	Vit-PN module
4	M3x8 oval head screw, TORX


Work steps:


1. Open housing (refer to Chapter 7.3).
2. Disconnect "Vit-PN" cable from the power print.

3. Unplug both the black and green tubes on the Vit-PN unit by pressing the colour ring towards the plug.

4. Undo the four screws on the Vit-PN unit.

5. Remove the four screws from the Vit-PN unit.
6. Carry out assembly in reverse order of dismantling. Connect the black tube to the black connection (above), connect the green tube to the green connection (below).
7. Service mode: set supplied Vit-PN correction factor.

7.5.1 Test Vit-PN

The function test is done without an attached instrument.


1. Select the VIT function.
2. Press the pedal briefly to switch the valves.
3. Air must flow out of the lower, green connector.

 If air flows out of the black connector, the tubes were swapped between the Vit-PN module and the connectors. In that case, the stripper will not work properly. Change tubing accordingly.


7.6 Replacement of coupling plate

Quantity	Description
1	Coupling plate with force sensor
3	M3x8 countersunk screw, TORX
1	M3x6 oval head screw, TORX
1	M3 lamellar disc

Work steps:

1. Open housing (refer to Chapter 7.3).
2. Disassemble Vit-PN unit (refer to Chapter 7.5).
3. Pull off the connectors [1] for the force sensor, tube detection switches and IRR valve on the system board.
4. Remove grounding cable [2].
5. The entire coupling plate can be removed by undoing three screws [3].

DANGER!


The grounding cable must be fixed with the lamellar disc and overhead screw.

6. When putting in the coupling plate in place, ensure that the switch lever [4] can move freely; it must not be impeded by any part of the housing.


7.7 Replacement of system board

Quantity	Description
1	System board (contains Logic PCB and Power PCB)
4	M3x8 oval head screw, TORX

Work steps:

1. Open housing (refer to Chapter 7.3).
2. Unplug all cables on the logic PCB. Refer to Section 5.2 for an overview of the interfaces.
3. Unplug all cables on the power PCB. Refer to Section 5.3 for an overview of the interfaces.
4. Undo the four screws of the holding sheet.

5. Take off the holding sheet.
6. Take out the system board.
7. Carry out assembly in reverse order of dismantling. Make sure the printed circuit boards are correctly slotted in the circuit board guide.


No system unit can be installed which has already been used with another control unit.


7.8 Replacement of power supply unit

Quantity	Description
1	Power supply
3	M3x6 oval head screw, TORX
1	M3x8 cheese head screw
1	M3 spring ring


Work steps:

1. Open housing (refer to Chapter 7.3).
2. Disconnect mains switch cable from power supply.

3. Undo the four screws on the power supply unit and take out the power supply unit.

4. The assembly of the components is done in the reverse order of disassembly.

DANGER!

Fix grounding wire with cheese-head screw and spring ring during assembly.

7.9 Replacement of clamp for I/V pole


VX210111 Clamp for I/V pole

Quantity	Description
1	Clamp for I/V pole

8 Software update

8.1 Device

The service interface is located below the display.

- 1 Service interface
- 2 Connections for vitrectomy handpiece
- 3 Connection for diathermy handpiece
- 4 Connection for phaco handpiece

New software can be loaded onto the device using a conventional USB memory stick. The stick must be properly formatted.

The connector with the name "Service" is to be used as the interface. This is a standard USB interface.

DANGER!

This connector is not designed for connecting conventional components like printers, etc. Connecting such components poses a risk to the safety of the patient!

The USB controller in the device only supports USB memory devices formatted in FAT12, FAT16 or FAT32 (FAT being the name of the file system used) with a sector size of 512 bytes. In addition, only mass storage devices which support SCSI or AT instructions are acceptable.

1. In order to format the USB mass storage device correctly, a command window must be opened in Windows. Start → Execute → type in "cmd".
2. In the command window, type in the following line: `FORMAT "Name" /FS:FAT32 /A:512`.
"Name" is the designation of the drive for the USB mass storage device (e.g. H:).

On USB mass storage devices larger than 2GB this formatting procedure is not possible for the whole device. In that case, the device has to be divided into different partitions where the primary partition has to be smaller than 2GB.

Please note that executing the format command is usually only possible with administrator privileges.
The USB memory will not be recognized if a sector size other than 512 byte is used.

Once the USB memory is formatted, it can be used for any number of updates, as long as the file system is not altered.

The software (CR3AS_xxxx.BIN) can now be copied to the USB memory stick using a standard Explorer window. Care should be taken that all of the older software on the memory stick is deleted in advance and that the software is copied to the main directory.


Procedure:

1. Insert the USB stick into the "Service" connector. Ensure that it is the right way round.
2. When starting up, hold the "Reset" key depressed until the programme opens the boot loader (approx. 4 sec).

The update procedure will then start automatically. It is important not to interrupt this process (by pulling out the memory stick)!

A new programme will be installed during the update. It is subsequently tested for proper execution. If, for some reason, the procedure cannot be completed successfully, the partially programmed firmware on the device will be deleted. In this situation, the device cannot be operated until new software has been installed.

Function display during the update:

Following a successful software update:

Boot-Loader Software update OK Remove Memory stick
--

Following an unsuccessful software update:

Boot-Loader WARNING: DEVICE IS NOW WITHOUT FIRMWARE Remove Memory stick

or

Boot-Loader WARNING: FIRMWARE UPDATE NOT COMPLETED. PLEASE REPEAT FIRMWARE UPDATE!
--

8.2 Pedal


New software is installed by the manufacturer.

9 Installing replacement parts of pedal

9.1 Replacement of cover plate

Quantity	Description
1	Cover plate
9	M3x8 oval head screw, TORX

Work steps:


1. Unscrew the 9 oval head screws and remove the cover plate.
2. Mount the new cover plate.


9.2 Replacement of pedal electronics board

Quantity	Description
1	Pedal PCB
1	Seal


Work steps:


1. Unscrew and remove the cover plate (see Chapter 9.1).
2. Remove seal.

3. Unplug all cables.
4. Remove the old print.

5. Install the new print.
6. Plug all cables back in.


7. Apply the new seal.
8. Mount the cover plate.


9.3 Replacement of pedal cable

Quantity	Description
1	Cable for pedal
1	Seal
1	TEMPO SIL2 white


Work steps:

1. Unscrew and remove the cover plate (see Chapter 9.1).
2. Remove the seal.
3. Unplug the CAN cable, and remove the green rubber seal.

4. Insert the new CAN cable.
5. Pre-fill the pouch with TEMPO SIL2.


6. Press the CAN cable into the channel and fill the pouch to above the rim with TEMPO SIL2.
7. Place adhesive tape on the sealing surface and wait until TEMPO SIL2 has cured (approx. 2 minutes).
8. Remove adhesive tape and use a sharp razor blade to remove excess TEMPO SIL2.
9. Reconnect the CAN cable.
10. Apply the new seal and mount the cover plate.

9.4 Replacement of pedal sensor

Quantity	Description
1	Hall effect sensor, wired
1	Seal
1	TEMPO SIL2 white
2	M3x12 oval head screw, TORX

Work steps:

1. Unscrew and remove the cover plate (see Chapter 9.1)
2. Remove the seal.
3. Unplug the defective hall effect sensor and horizontal switch, and remove the green rubber seal.
4. Remove the defective hall effect sensor.
5. Mount the new hall effect sensor.
6. Pre-fill the pouch with TEMPO SIL 2.
7. Press the hall effect sensor cable and horizontal switch cable into the sealing compound.
8. Fill the pouch to above the rim with TEMPO SIL 2.
9. Place the adhesive tape on the sealing surface and wait until TEMPO SIL 2 has cured (approx. 2 minutes).
10. Remove adhesive and use a sharp razor blade to remove excess TEMPO SIL 2.
11. Plug in the hall effect sensor and horizontal switch to pedal print.
12. Apply the new seal and mount the cover plate.
13. Perform pedal check (see Chapter 6.5).


CAUTION!

A new sealing mat must be used after repair on the open pedal!


9.5 Replacement of rocker

Quantity	Description
1	Rocker
1	Wave vertical deflection
1	M4x16 grub screw


Work steps:


1. Unscrew and remove the cover plate (see Chapter 9.1).
2. Remove the seal.
3. Unplug the hall effect sensor and horizontal switch, and remove the green rubber seal.
4. Undo the grub screw.

5. Push out the wave vertical deflection component from the base body.

6. Raise the slotted lever and carefully position the rocker.
7. Take out the horizontal switch cable through the opening in the base body and remove the old rocker.

8. Assemble the new rocker in reverse order.


Make sure the cables are not kinked!

9.6 Replacement of pedal handle

Quantity	Description
1	Pedal handle
4	M5x16 countersunk screw, TORX, black

Work steps:


1. Undo the 4 countersunk head screws and remove the old pedal grip.
2. Insert the new pedal grip and fasten it with the 4 new countersunk head screws.

9.7 Replacement of rubber mat

Quantity	Description
1	Rubber mat

Work steps:


1. Remove the old or defective rubber mat.
2. Thoroughly clean the rocker attachment with medical benzene.
3. Apply the new rubber mat to the rocker attachment.

9.8 Replacement of rubber buffer

Quantity	Description
4	Rubber buffer with thread



Work steps:

1. Remove the old or defective rubber buffers.
2. Adhere the new rubber buffers on the base body.

9.9 Set of replacement screws

VX950045 Set of replacement screws

Quantity	Description	Picture
10	M4x8 countersunk bolt Torx	
20	M5x16 countersunk bolt Torx	
50	M3x8 eco-fix Torx	
10	M3x12 cylinder bolt Torx	
10	M3x6 grub screw I-6kt	
10	M3x8 countersunk bolt Torx	
10	M3 lamellar disc	
10	M3x6 eco-fix Torx	
10	M3x8 cylinder bolt Torx	
10	M3 spring ring	
10	M3x12 eco-fix Torx	

9.10 Set of small replacement parts

Quantity	Description
5m	Rubber sponge sealing ring d=3mm
5m	Sealing strip
20	Cable strap
10	Cable strap holder
10	Rubber buffer 12.7 x 3.5mm black

9.11 Thread locking

With the exception of the screws used to secure the electronic print, Loctite must be applied to all screws in order to guarantee long-term stability. The following Loctite thread lockers are recommended.

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10 Order number for replacement parts

10.1 Control unit

Replacement fuses
Upper cover plate
Lower cover plate
Pump motor
Pump wheel
VIT-PN module
Coupling plate with force sensor
System board
Control Panel
Power supply unit
Clamp for I/V pole
Set of small replacement parts
Set of replacement screws

10.2 Pedal

Pedal PCB
Cable for Pedal
Pedal sensor
Pedal handle
Rubber mat for pedal
Rubber buffer for pedal
Sealing cover for pedal
Pedal rocker

11 Messages, warnings and error messages / fault correction

Message	Explanation
SYSTEM READY	Successful completion of the system check when starting up the system.
SELECT SURGEON	Asks you to select the appropriate surgeon memory.
CONNECT PEDAL	The pedal is not connected correctly or not at all.
INSERT TUBE	The tubing system is not inserted correctly or not at all (2 recognition switches).
PREOP FILL RINSE 100%	PREOP is running, progress displayed as a percentage.
PREOP DISCONTINUED!	PREOP has been interrupted.
PHACO TEST DISCONTINUED!	The phaco test has been interrupted.
PHACO TEST ACTIVE	The phaco test is running, nothing else can be activated.
PHACO TEST O.K.	The phaco test has been successfully completed.
PLEASE REPEAT PHACO TEST	The phaco test could not be carried out due to activated pedal. The phaco test must be restarted.
TEST HANDPIECE!	The handpiece shows low performance.
HANDPIECE DEFECTIVE!	The handpiece can't be operated.
CHECK TIP!	The tip has not been screwed on properly.
STORING VALUES	Values are being stored (operator has to keep the key pressed).
VALUES STORED	Values have been stored (operator can let go of the key).
NOT STORED	Values have not been stored (operator did not hold the key down long enough).
TEMPERATURE TOO HIGH!	The temperature in the device is too high. Wait until the device has cooled down, the message will disappear.
NOT AVAILABLE	HFDS function not available.
SERIAL NUMBER IS MISSING	Serial number not programmed in the device.
SERIAL NUMBER DOES NOT MATCH	Invalid serial number programmed in the device.

Error Message	Cause	Corrective actions
ERROR 1 CALL SERVICE	Internal operating values (+3.3V, +5V, +12V, -12V, +24V) not within tolerance range	• Check the operating voltages in service mode. • Remove the cable between the power supply and the system board: check the voltage on the power supply. • Replace the system board or the power supply.
ERROR 2 ADJUST DIA	DIA adjustment has not been adjusted correctly or not at all	Replace the system board and notify the manufacturer.
ERROR 3 ADJUST PEDAL	Pedal is not adjusted correctly or not at all (impossible reference values for pedal positions)	The pedal has to be readjusted as described in the Section "Service mode operation / pedal calibration".
ERROR 4 GRAPHIC DISPLAY	Graphic display is not ready or does not respond	• Check the cable from the system board to the control panel. • Replace the control panel.
ERROR 8 PROGRAM FAILURE	Undefined state	• Check for the presence of extreme electromagnetic radiation from a nearby non-Oertli appliance. • Please notify the manufacturer if the message is displayed again.
ERROR 9 NV-RAM FAILURE	NV-RAM defective or missing	Replace the system board and notify the manufacturer.
ERROR 10 PERISTALTIC PUMP FAILURE	Pump does not react, pump does not rotate (interruption in the pump control circuit)	Check the encoder cable and the motor cable to the pump motor.
ERROR 11 OVERLOAD	U_{pow} too low (PTC2 has responded)	• Remove the continued short circuit in the PHACO, DIA or CAPS instruments. • Wait until the affected circuit has cooled down, the message will disappear.
ERROR 13 CHECK PUMP	The pump does not react nor rotate (interruption in the pump control circuit or resistance at the edge of the pump too strong)	Check tubing system. The message will disappear after approx. 10 seconds.
ERROR 15 FORCE SENSOR	Force sensor has not been plugged in or is defective	• Check the force-sensor cable from the system board to the coupling plate. • Replace the coupling plate with the force sensor.
ERROR 17 ADJUST FORCE SENSOR	Sensor has been adjusted incorrectly or not at all	The force sensor has to be readjusted as described in the section "Service mode operation / sensor calibration".
Display remains dark	No supply voltage	• Check the fuses at the mains power socket. • Check the voltage on the power supply.
	Display defective	• Check the cable from system board to the control panel. • Replace the control panel.
CONNECT PEDAL message remains, even when the pedal is connected	Pedal defective or another type of pedal is connected	• Check / replace the pedal cable. • Replace the pedal. • Use only CataRhex 3 or SwissTech pedals.
The PED sign will be displayed even if the pedal is not operated	Pedal is misadjusted	The pedal has to be readjusted as described in the Section "Service mode operation / pedal calibration".
Limit value cannot be reached by moving the pedal outwards	Pedal is misadjusted	The pedal has to be readjusted as described in the Section "Service mode operation / pedal calibration".

12 Authorized service points

Oertli Instrumente AG

Hafnerwisenstrasse 4

CH – 9442 Berneck, Switzerland

Tel.: +41 (0)71 747 42 00

E-mail: info@oertli-instruments.com

Website: www.oertli-instruments.com

Contact your local distribution partner or Oertli Instrumente AG for information on authorized service points in your country.

12.1 Sending defective device, instrument or part to Oertli Instrumente AG

1. Clean and sterilize instruments before returning them to the service point.
2. Send device, instrument or part to the service point.

13 Functional test

Customer:

Serial numbers:

Article N°	Serial N°	SW version

All tests must be done with fluid!

Peristaltic pump-system:

- Select I/A mode, press PREOP button
- Build occlusion with doctors I/A values:
- Build occlusion with 50ml/600mmHg
- Keep occluded and release pedal

Connect I/A-tubing ports together!

- 0-100% approx. in 1 min ☐
- Build-up time OK ☐
- in < 2.5 sec ☐
- Vents to atmospheric pressure ☐

- Activate reflux (heel position)
- Activate reflux & build occlusion
- Press the IRR button (continuous irrigation)
- Stops when reflux is activated

Disconnect I/A-tubing!

- Aspiration flows back ☐
- Pump stops at 100 mmHg ☐
- Irrigation continuous when pedal released ☐
- ☐

Phaco:

- Select PHACO, press "Test"
- Go for maximum phaco power
- Activate PULSE, BURST, CMP
- Build occlusion with surgeon values
- Check occlusion mode

(Mount tip, test-chamber & tubing)

- Test OK ☐
- Reaches 100% ☐
- Function OK ☐
- Release occlusion; chamber stays stable ☐
- Function OK ☐

Dia/Caps:

- Activate DIA; multimeter or instrument
- Activate CAPS; multimeter or instrument

- Function OK ☐
- Function OK ☐

Vit:

- Activate VIT with Instrument

- Function OK ☐

Enter service mode

- Pedal:
- Sensors (remove cassette / tubing):
- Buttons:

- All positions reach 100% ☐
- Sensor 0 ± 10, only adjust if necessary! ☐
- Push all buttons ☐

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
Caution!

The tubing set must be checked for completeness, possible damage and leaks prior to assembly.

Caution!

Products with visible damage or chipping must be sorted out and replaced. Do not continue to use products damaged in this way!

Warning!

In accordance with processing instruction TN999042, the tubing set must be cleaned, disinfected and sterilized prior to each application.

Warning!

Tubing sets intended for single-use must not be re-used.

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
Intended Use

I/A tubes with a coupling part for use with the faros device or CataRhex series. Suction and rinsing instrument for all surgical applications.

This hose has to be used exclusively with the CataRhex SwissTech, CataRhex 3 or faros units.

Delivery scope

The tubing set comprises:

- | | |
|-----------------|---|
| 1 Coupling part | 4 Irrigation tube |
| 1a Handle | 5 Irrigation port |
| 2 Pump tube | 5a Protection cap (single-use only) |
| 2a Tube stop | 7 Aspiration port |
| 2b Tube stop | 8 Infusion port |
| 3 Drainage tube | 9 Rinsing tube (designed for multiple use only) |

Application

To ensure correct assembly and dismantling of the tubing set, please follow the device instruction manual.

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
- 

Caution!
U.S. Federal Law restricts this device to sale by or on the order of a physician.
- Rx Only

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Warning!

In accordance with processing instruction TN999042, the diathermy, glaucoma and capsulotomy tips must be cleaned, disinfected and sterilized prior to each application.

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Caution!

Products with visible damage, deformation, corrosion, chipping or discoloration at their tips (usually black) must be sorted out and replaced. Do not continue to use products damaged in this way!

Make sure the power cable and the plug-in connector are undamaged.

Intended Use

Diathermy tips for use with a diathermy handpiece for diathermy and coagulation.

Diathermy, glaucoma and capsulotomy tips have to be used exclusively with compatible Oertli instruments.

Product description

Correct geometry of the chosen products:

VE201726 (Standard capsulotomy tip)

VE201751 (abee® glaucoma tip)

Application

Capsulotomy tips and the tip for glaucoma surgery (VE201751) may only be used if the anterior chamber is filled with a fluid for visco-surgery (e.g. Ophtalin, Healon, Viscoat, adatoceel, Ocucoat).

Sufficient physiological saline solution needs to be present at the points of use for the satisfactory function of the eraser tip (VE201721 – VE201722).

To ensure correct use of diathermy tips, please follow the device instruction manual.

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Reusable / autoclavable

Warning!

In accordance with processing instruction TN999042, the phaco handpiece must be cleaned, disinfected and sterilized prior to each application.

The maximum number of sterilization cycles is defined in processing instruction TN999042. This maximum number of sterilization cycles must not be exceeded.

Warning!

Tips intended for single-use must not be reused.

Warning!

The phaco tip vibrates at ultrasound frequency (28 kHz, max. 120 µm deflection) and can result in burns through friction if it is pressed against the incision.

Warning!

To reduce the risk of tissue burns, always fit the infusion sleeve over the tip, refrain from holding the instrument at an angle in the incision, keep the ultrasound energy as low as possible and preferably apply the ultrasound intermittently instead of continuously.

EXCEPTION: Bimanual phaco with CMP.

Caution!

Products with visible damage, deformation, corrosion or chipping must be sorted out and replaced. Do not continue to use products damaged in this way!

Make sure the power cable and the plug-in connector are undamaged and the lumen are not blocked.

Intended Use

The phaco handpiece has to be used exclusively for cataract surgery using the ultrasound emulsification method in conjunction with plug-in connector compatible Oertli equipment.

Application

Warning!

The infusion line must always remain connected to the handpiece in order to avoid excessive warming.

EXCEPTION: Bimanual phaco with CMP.

To ensure correct use of the phaco handpiece, please follow the device instruction manual.

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
Definitions

For single use only, do not reuse!

Content is sterile if package is unopened and undamaged (sterilized with ethylene oxide).

Do not use if package is damaged.

Application

WARNING!

Single use only – do not reuse!

WARNING!

Adequate cleaning and reconditioning of products with small lumen (inner diameter and length) cannot be guaranteed and there is an infection hazard for the patient.

- For use with Oertli® OS4, OS3, faros and CataRhex 3 console and in combination with trocar system of same size only.
- Observe color coding, connect green Luer with green console port (PN2: OS3, faros: right; OS4, CataRhex 3: bottom), black with black (PN1; OS3, faros: left; OS4, CataRhex 3: top)
- Do not cover the venting holes on the cutter body.
- Do not operate dry, always aspirate with liquid.
- Remove elongation sleeve A as desired.

WARNING!

The safety instructions in the device instruction manual must be observed.

WARNING!

Hazard: Use of instruments and consumables not manufactured by Oertli Instrumente AG.

Impact: Interference of system performance and potential hazard for user and patient

- Only use consumables and instruments in conjunction with Oertli® devices that are manufactured by Oertli Instrumente AG.

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