
NOTARIAL CERTIFICATE

BE IT KNOWN that I, Esmond Geoffrey John Jenkins of Old Manor Barn,
Shorncliffe, Cirencester in the County of Gloucestershire, England, a duly authorised
Notary Public, CERTIFY AS FOLLOWS:

1. That the paper writings annexed hereto is an original document produced to
me and signed by Clara Gladys Abbott the Principal Regulatory Specialist of
Abbott Diabetes Care Limited registration documents who is authorised to
sign Registration documents for the company and who is resident at 22 Old
School Close, Caversfield, Bicester in the County Oxfordshire OX27 8TW
and who is identified to me by Irish Passport number PC3375506 and UK
Driving licence number ABBOT852250CG9AE 36

SIGNED and sealed Bampton Law LLP Solicitors Market Square, Bampton in the
Country of Oxfordshire, England on 6th October 2017

Esmond Geoffrey John Jenkins

Notary Public
England and Wales

Protocol 17- 482

Esmond Jenkins Notary Public England and Wales
Clanfield House, Market Square
Bampton, Oxfordshire OX18 2JJ
Tel. 01993 852222

English version of ART36195-001-Instructions for Use

CLARA ABBOTT
06 OCT 17

INSTRUCTIONS FOR USE
Blood Glucose and Ketone Monitoring System
FreeStyle Optium Neo with Accessories

Manufacturing company:

(ABBOTT DIABETES CARE LIMITED, UK)

Corporate manufacturers:

1. Abbott Diabetes Care Ltd., RANGE ROAD WITNEY OXON OX29 0YL, UK
2. Abbott Diabetes Care Inc., 1360 South Loop Road, Alameda, CA 94502, USA
3. Asahi Polyslider Co., Ltd., Nakanoshima Mitsui Building 3-3-3, Nakanoshima Kita- ku, Osaka 533-0005, Japan
4. RR Donnelley Global Turnkey Solution Sp.z o.ul. Zakladowa 90/92 92-402 Lodz Poland
5. Owen Mumford Ltd., Brook Hill, Woodstock, OX20 1TU, UK
6. Flextronics Technology (Shenzhen) Co., Ltd. Building 1, Yusheng Industrial Park, 467 Xixiang Section, National Highway 107, Xixiang, Baoan District, 518126 Shenzhen, Guangdong, PEOPLE'S REPUBLIC OF CHINA
7. Facet Technologies, LLC, 101 Liberty Industrial Parkway, McDonough, GA, 30253, USA.

1. Name of the medical device:

Blood Glucose and Ketone Monitoring System - FreeStyle Optium Neo with Accessories:

I. Blood Glucose and Ketone Monitoring System - FreeStyle Optium Neo, including:

1. FreeStyle Optium Neo Blood Glucose and Ketone Monitoring System.
2. FreeStyle Optium Blood Glucose Test Strips (10, 25, 50, 100 pcs) manufactured by Abbott Diabetes Care Ltd., UK, Registration certificate No. 2012/13309 (if required).
3. FreeStyle Optium Blood β -Ketone Test Strips (10, 25, 50, 100 pcs) manufactured by Abbott Diabetes Care Ltd., UK, Registration certificate No. 2012/13309 (if required).
4. FreeStyle lancing device II -1pcs (if required).
5. Thin lancets (10, 25,50,100 pcs) (if required).
6. User's Manual - 3 pcs.
7. USB Cable (if required).

II. Accessories:

Carry case (1 pcs).

2. Purpose

FreeStyle Optium Neo Blood Glucose and Ketone Monitoring System is designed for stand-alone use at in-vitro diagnostics by the patient or medical staff as an aid in the management of diabetes. Used for measuring glucose in fresh whole blood samples taken from fingers, forearm, upper arm

or the base of the thumb. Use for measuring ketone (β -hydroxybutyrate) in fresh whole blood samples from fingers only. The system may also aid in the management of diabetes by providing the user with suggested recommendations to insulin dose(s) based on healthcare professional entered data.

1 Indications

The FreeStyle Optium Neo Blood Glucose and Ketone Monitoring System is for use outside the body only (in vitro diagnostic use) for self testing or professional use as an aid in the management of diabetes.

Used for measuring glucose in fresh whole blood samples taken from fingers, forearm, upper arm or the base of the thumb. Used for measuring ketone (β -hydroxybutyrate) in fresh whole blood samples from fingers only.

The system may also aid in the management of diabetes by providing the user with suggested recommendations to insulin dose(s) based on healthcare professional entered data.

2 Contraindications:

The FreeStyle Optium Neo Blood Glucose and Ketone Monitoring System is not for use in diagnosis of diabetes mellitus, but is to be used as an aid in monitoring the effectiveness of diabetes control programs.

3 Side Effects:

Bruising may occur at the test site. If bruising occurs, consider selecting another site.

3. Main Principle of Medical Device, Functioning and Specifics

The FreeStyle Optium Neo Blood Glucose and Ketone Monitoring System together with compatible glucose and ketone level test strips is based on amperometric analysis, with glucose level measured by reaction of blood glucose with glucose dehydrogenase (GDH), and β -hydroxybutyrate (β -OHB) is measured by its reaction with hydroxybutyrate dehydrogenase (HBDH) in blood samples or reference solutions by means of electric current. The device measures electric current resulted from reaction of glucose, ketones and special ferment applied to the test strip and than converts the measurements into a clinical indication. The value of the generated current is proportional to the analyzed substance (glucose or ketone) level in blood sample or reference solution.

4. Specifications

The FreeStyle Optium Neo Blood Glucose and Ketone Monitoring System must comply with following requirements:


Parameters	FreeStyle Optium Neo
Size	5.97 cm (w) x 8.68 cm (l) x 0.87 cm (d)
Weight	33 to 37 g (1.2 to 1.3 oz.) including batteries
Power Source	Two CR 2032 lithium (coin cell) batteries
Battery Life	Up to 3000 tests
Memory	Up to 1000 events, including blood glucose, blood ketone, and control solution results, insulin doses, and other meter information

Compatible Blood Glucose Test strips	FreeStyle Optium Blood Glucose Test Strips
Compatible Blood Ketones Test Strips	FreeStyle Optium Blood β -Ketone Test Strips
Assay Method	Amperometric electrochemistry
Assay Range	For blood glucose testing 1.1 - 27.8 mmol/L For blood ketone testing 0.0 - 8.0 mmol/L
Test Time	5 seconds - glucose 10 seconds - ketone
Hematocrit Range	30%-60% - glucose 30%-60% - ketone
Calibration Position	YSI x 1.12 (plasma equivalent) - glucose
Sample Type	Glucose – Capillary Ketone – Capillary and Venous
Meter Storage Temperature Range	-20 °C to 60 °C (-4 °F to 140 °F)
Humidity Range	10%-90%
Altitude**	7,200ft (2,195m) - glucose 7,200ft (2,195m) - ketone
Sample Type	Glucose – Capillary blood Ketone – Capillary and venous blood
Sample Sites	Glucose - Fingertip, base of thumb, forearm or upper arm Glucose - Fingertip (capillary sample) or venous blood (by health care professionals in hospital environment)
Sample Volume	$\geq 0.6 \mu\text{l}$ – glucose $1.5 \mu\text{l}$ - ketone
Re-Application	Glucose – 5 seconds Ketone – 10 seconds

* Testing for glucose was completed at 15-65% for glucose testing, however the claim will remain 30-60%.

** Testing completed at 10,000ft (3,148m) however the product claim will remain 7,200ft (2,195m).

5. Information on appearance and design (design, main properties and materials used for medical device components)

Material: Housing - PC-ABC Carbotex

Black colorant - Carbon Black

Display - Lexan polycarbonate matte film

Test strip window PC-ABC Carbotex

Size – 6.0 / 8.6 ± 5% cm

Weight with battery: max. 37 g ±5%

Dialog window: 4.1/4.3 ±1% cm

Mediated contract through patient's arm skin


All materials are non-toxic and comply with Regulation (EC) 1907/2006 by European Parliament and Council of Europe dated December 18, 2008 on Registration, Evaluation and Authorization of Chemicals (REACH).

6. Scope of supply:

Blood Glucose and Ketone Monitoring System - FreeStyle Optium Neo with Accessories:

- I. Blood Glucose and Ketone Monitoring System - FreeStyle Optium Neo, including:
 1. FreeStyle Optium Neo Blood Glucose and Ketone Monitoring System.

2. FreeStyle Optium Blood Glucose Test Strips (10, 25, 50, 100 pcs) manufactured by Abbott Diabetes Care Ltd., UK, Registration certificate No. 2012/13309 (if required).
 3. FreeStyle Optium Blood β -Ketone Test Strips (10, 25, 50, 100 pcs) manufactured by Abbott Diabetes Care Ltd., UK, Registration certificate No. 2012/13309 (if required).
 4. FreeStyle lancing device II -1pcs (if required).
 5. Thin lancets (10, 25,50,100 pcs) (if required).
 6. User's Manual - 3 pcs.
 7. USB Cable (if required).
- II. Accessories: Carry case, 1 pcs.

1

FreeStyle Optium Neo Overview

Display window

- Device main screen (shown here)
- Shows readings and other important info


USB port

- Insert cable here to transfer data when loading readings into PC (requires data management system)

Down button

Operating space

- Apply blood or reference solution on the operating space located at the end of test strip

2

Symbols shown on the screen

Symbol	Meaning	Symbol	Meaning
	Log		Insulin dose administered during breakfast
	Device is ready for sampling		Insulin dose administered during lunch
	Result of the reference solution analysis		Insulin dose administered during dinner
	Insulin level		Test in the fasted state
	Insulin		Ketones
	Morning dose of long-acting insulin		Battery low
	Evening dose of long-acting insulin		Connected to PC
	Insulin dose setting		Setup mode
	Low glucose level or level pattern		High glucose level or level pattern

3

Quick guide – blood glucose measuring

- 1. Wash and dry your hands before and after the measuring**

- 2. Insert the test strip.**

- 3. Puncture your skin at analysis site and apply blood on the strip.**

- 4. See the result.**
Additional information on test procedure is available in Section 9.

4

Readings

The FreeStyle Optium Neo Blood Glucose and Ketone Monitoring System is for use outside the body only (in vitro diagnostic use) for self-testing or professional use as an aid in the management of diabetes.

Use for measuring glucose in fresh whole blood samples taken from fingers, forearm, upper arm or the base of the thumb. Use for measuring ketone (β hydroxybutyrate) in fresh whole blood samples from fingers only.

The system may also aid in the management of diabetes by providing the user with suggested recommendations to insulin dose(s) based on healthcare professional entered data.

IMPORTANT:

- Use only FreeStyle Optium blood glucose test strips and FreeStyle Optium blood β -Ketone test strips. Other test strips may produce inaccurate results.
- See test strip instructions for use for more information about sample types.
- Read the instructions in this Owner's Setup Guide. Failure to follow instructions may cause incorrect results. Practice the testing procedures before using the meter.
- Follow your healthcare professional's advice when testing blood glucose levels and blood ketone levels.
- Observe caution when using around children. Small parts may constitute a choking hazard.

5

Contraindications

The FreeStyle Optium Neo Blood Glucose and Ketone Monitoring System is not for use in diagnosis of diabetes mellitus, but is to be used as an aid in monitoring the effectiveness of diabetes control programs.

6

Side Effects

Bruising may occur at the test site. If bruising occurs, consider selecting another site.

7

Getting to Know Your Meter

To turn your meter on:

- Press , or
- Insert a strip

To turn your meter off:

- Press and hold  for 3 seconds, or
- Do nothing for 2 minutes

Checking Your Meter Screen Every Time You Turn It On

Meter Off Screen

Your meter screen should be fully black when powered off. Each time you turn on your meter, a white start-up test screen will appear for 1 second.

If you see any white segments in the black off screen, or any black segments in the white test screen, there may be a problem with the meter. Contact Customer Care.

Start-up Test Screen

Note: If the meter battery is low, the  will appear in both the meter off screen and start-up test screen.


1. Start with meter off (no test strip inserted)

- Press  to turn on the meter

2. Press and hold the time (11:50) on the screen for 3 seconds until the screen changes.

12 h
(e.g. 1:23 PM)

24 h
(e.g. 13:23)

3. Set Time Format (12 hour or 24 hour clock)

- Press  or  to change
- Press  to continue

Note: Your meter can display either a 12 h (1:23 PM) or 24 h (13:23) time format. If you prefer the 12 h format, there is no "AM". If setting a PM time, continue to press until you see "PM".

4. Set Hours


- The hours blink. Press  or  to set the hours
- Press  to continue

5. Set Minutes


- The minutes blink. Press  or  to set the minutes
- Press  to continue

Setting Up the Date


6. Set Month

- The month blinks. Press  or  to set the month
- Press  to continue

7. Set Day

- The day blinks. Press  or  to set the day
- Press  to continue

8. Set Year

- Press  or  to set the year
- Press  to save. Setup is complete

Testing Your Blood Glucose or Blood Ketone

Performing a Blood Glucose or Blood Ketone Test

- 1. Open the foil test strip package at the notch and tear down to remove the test strip**

- 2. Insert the test strip into the meter until it stops.**

This will turn on the meter.

Notes:

- Be sure to check that your meter screen is working properly each time you turn your meter on. If you see any white segments in the black off screen, or any black segments in the white test screen, there may be a problem with the meter. (See Section 8, Getting to Know Your Meter, for more details)
- The meter turns off after 3 minutes of inactivity. Remove and reinsert the unused test strip to restart the meter.

The  blinks, indicating the meter is ready for you to apply a sample to the test strip.

Note: KET will appear on the screen if you have inserted a purple blood ketone strip.

- 3. Obtain a blood sample.**

Use the lancing device to obtain a blood sample. (See lancing device instructions for use for more information.)

- 4. Apply blood to the test strip.**

Bring the blood drop to the white area at the end of the test strip. The blood is drawn into the test strip.

Hold blood to test strip until you see 3 short lines on the meter screen. This means you have applied enough blood.

Notes:

- If you are testing blood glucose, you will see a 5-second countdown. If you are testing blood ketone, you will see a 10-second countdown.
- Do not remove the test strip from the meter during the countdown.
- If the countdown does not start, you may not have applied enough blood to the test strip. See test strip instructions for use for re-application instructions. If the countdown still does not start, remove the used strip and discard it properly. Start a new test with a new test strip.

Blood Glucose
Result example

Blood Ketone
Result example

5. View the result.

The test is complete when the result appears on the meter screen (examples shown). The result is stored in memory.

Note: If you see , it means the meter can suggest a new meal insulin dose. Ignore this if you're not about to eat a meal and take your meal insulin dose.

6. Press and hold  to turn off the meter. Discard the used test strip properly.

The meter displays blood glucose results in mmol/L. The unit of measurement is preset. You cannot change this setting.

IMPORTANT: The meter displays results from 1.1 - 27.8 mmol/L. Low or high blood glucose results can indicate a potentially serious medical condition.

Low Blood Glucose Results

What You See	What It Means	What To Do
 Solid red arrow  A blinking red arrow LO Appears when result is lower than 1.1 mmol/L.	Appears when result is lower than 3.9 mmol/L or the target set on the meter by your healthcare professional. A pattern of low glucose has developed. If 2 low results occur within the past 5 days AND both are within the same 3-hour time period, the meter will display a blinking ↓ Severe low blood glucose <i>or</i> There may be a problem with the test strip.	Follow your healthcare professional's advice to treat low blood glucose. Follow your healthcare professional's advice to treat low blood glucose. Repeat the test with a new test strip. If the result is LO, contact your healthcare professional immediately.

Note: If you see the error messages E-3 or E-4, consult the Error Messages section in this Owner's Setup Guide.

IMPORTANT: Contact your healthcare professional if you have symptoms that do not match your test result, and you have followed the instructions in this Owner's Setup Guide.

High Blood Glucose Results

What You See	What It Means	What To Do
 Solid yellow arrow  A blinking yellow arrow  KET HI Appears when result is higher than 27.8 mmol/L.	Appears when result is higher than 13.3 mmol/L or the target set on the meter by your healthcare professional. A pattern of high glucose has developed. If 3 high results occur within the past 5 days AND all are within the same 3-hour time period, the meter will display a blinking ↑. Blood glucose level is higher than or equal to 13.3 mmol/L. Severe high blood glucose <i>or</i> There may be a problem with the test strip.	Follow your healthcare professional's advice to treat high blood glucose. Follow your healthcare professional's advice to treat high blood glucose. Check blood ketone if checking ketones is part of your diabetes management program. Repeat the test with a new test strip. If the result is HI, contact your healthcare professional immediately.

11

Understanding Blood Glucose Test Results and Patterns

The meter displays ketone results in mmol/L, from 0.0 - 8.0 mmol/L. The unit of measurement is preset. You cannot change this setting.

12


IMPORTANT: Follow your healthcare professional's advice before you make any changes to your diabetes management program.

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Repeat a blood ketone test using a new blood ketone test strip when:

- HI appears on the display
- Your result is unusually high
- You question your result
- Your blood ketone result is 0.0 mmol/L, but your blood glucose is higher than 16.7 mmol/L

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View Logbook

Your meter logbook can store up to 1,000 events – including blood glucose, ketone, and control solution results, insulin doses, and other meter information.

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What You See	What It Means	What To Do
Result is between 0.6 – 1.5 mmol/L.	Result is between 0.6 – 1.5 mmol/L.	Follow your healthcare professional's advice.
Result is higher than 1.5 mmol/L.	You may be at risk of developing diabetic ketoacidosis (DKA) ²⁻⁶	Contact your healthcare professional immediately.
Appears when result is higher than 8.0 mmol/L.	Very high blood ketone or There may be a problem with the test strip.	Repeat the test with a new test strip. If the result is HI, contact your healthcare professional immediately.

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


Insulin Dose

Units of
Insulin

2. Press  to scroll through to view your logbook events.

Fasting result

Control Solution Result

Viewing Your Blood Glucose Averages

1. While on the Home screen,

press  to view your logbook events.

Press $\vee$ or $\wedge$ to scroll through 7-, 14- and 30-day averages.

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A control solution test should be performed when you are not sure of your results and want to confirm that your meter and test strips are working properly.

IMPORTANT:

- Use only MediSense glucose and ketone control solutions with the meter.
- Control solution results should fall within the control solution range printed on the test strip instructions for use.
- Check that the LOT number printed on the test strip foil packet and instructions for use match.
- Do not use control solution past its expiration date. Discard control solution 3 months after opening or on the expiration date printed on the bottle, whichever comes first. (Example: open April 15, discard July 15; write the discard date on the side of the bottle.)
- The control solution range is a target range for control solution only, not for blood glucose levels.
- Replace the cap securely on the bottle immediately after use.

IMPORTANT: (continued)

- Do not add water or other liquid to control solution.
- Control solution results do not reflect your blood glucose.
- Contact Customer Care for information on how to obtain control solution.

Performing a Control Solution Test

1. Open the foil test strip packet at the notch and tear down to remove the test strip.

2. Insert the test strip until the meter turns on.

Notes:

- Be sure to check that your meter screen is working properly each time you turn your meter on. If you see any white segments in the black off screen, or any black segments in the white test screen, there may be a problem with the meter. (See Section 8, Getting to Know Your Meter, for more details.)
- The meter turns off after 3 minutes of inactivity. Remove and reinsert the unused test strip to restart the meter.

The  blinks, indicating the meter is ready for you to apply a sample to the test strip.

Note: KET will appear on the screen if you have inserted a purple blood ketone strip.

3. Press and hold the down arrow for 3 seconds to **mark the test as a control solution test**. The  appears. The meter is now ready for you to apply control solution to the test strip.

IMPORTANT: The test result will be saved to memory as a blood result if not marked as a control solution test. This may affect your blood glucose averages.

4. Apply control solution to the test strip. Shake the control solution bottle to mix the solution.

Apply a drop of control solution to the white area at the end of the test strip in the area shown. The control solution is drawn into the test strip.

5. Hold the control solution to the test strip until:

- You see 3 short lines on the meter screen. This means you have applied enough control solution and the meter is reading the control solution.

Notes:

- If you are testing with a blood glucose strip, you will see a 5-second countdown. If you are testing with a blood ketone strip, you will see a 10-second countdown.
- Do not remove the test strip from the meter during the countdown.
- If the countdown does not start, remove and discard the used test strip, turn off the meter and try again with a new strip.

Examples:

**6. View the result.**

The test is complete (examples shown) when the result appears on the meter screen. The result is stored in memory as a control solution result.

Compare the control solution result to the range printed on the blood glucose or blood ketone test strip instructions for use. The result should fall within the range.

***Note:** KET appears with the result if performing a ketone control solution test.*

Out of Range Control Solution Results:

- Repeat the test if control solution results are outside the range printed on the test strip instructions for use.
- Stop using the meter if control solution results are consistently outside the range printed on the test strip instructions for use. Contact Customer Care.

Transferring meter data to a computer requires a compatible data management system. You will also need a micro USB cable to connect your meter's USB port to your computer. For more information, please contact Customer Care.

WARNING: To avoid the possibility of electric shock, never perform a blood glucose test while the meter is connected to the computer.

15

Error Messages

Message	What It Means	What To Do	Message	What It Means	What To Do
E-1	The temperature is too hot or too cold for the meter to work properly	1. Move the meter and test strips to a location where the temperature is within the test strip operating range. (See test strip instructions for use for the appropriate range.) 2. Wait for the meter and test strips to adjust to the new temperature. 3. Repeat the test using a new test strip. 4. If the error reappears, contact Customer Care.	E-3	Blood drop is too small or Incorrect test procedure or There may be a problem with the test strip	1. Turn off the meter. 2. Repeat the test using a new test strip. 3. If the error reappears, contact Customer Care.
E-2	Meter Error	1. Turn off the meter. 2. Repeat the test using a new test strip. 3. If the error reappears, contact Customer Care.	E-4	The blood glucose level may be too high to be read by the system or There may be a problem with the test strip	1. Repeat the test using a new test strip. 2. If the error reappears, contact Customer Care.
			E-5	Blood was applied to the test strip too soon	1. Review the testing instructions. 2. Repeat the test using a new test strip. 3. If the error reappears, contact Customer Care.

Message	What It Means	What To Do
E-6	Meter error	1. Check that you are using the correct strip for this meter. (See test strip instructions for use to verify your strip is compatible with this meter.) 2. Repeat the test using a test strip for use with your meter. 3. If the error reappears, contact Customer Care.
E-7	No coding required or Test strip may be damaged, used, or the meter does not recognize it	1. Check that you are using the correct test strip for this meter. (See test strip instructions for use to verify your strip is compatible with this meter.) 2. Repeat the test using a test strip for use with your meter. 3. If the error reappears, contact Customer Care.
E-9	Meter error	1. Turn off the meter. 2. Repeat the test using a new test strip. 3. If the error reappears, contact

16 Taking Care of Your Meter

Changing the battery

You will see this screen when your batteries are low.




Note: Your meter settings and logbook information will be saved when you change the batteries.

IMPORTANT: After you first see this warning, you can perform approximately 28 tests before you need to change the batteries.

WARNING: Batteries should be kept away from small children. If swallowed, contact a healthcare professional immediately.

Step	Action
	1. Turn meter over and slide open the battery door on the side as shown.

Step	Action
	2. Remove the old batteries.
	3. Install new batteries with (+) facing up. Note: The meter uses 2 replaceable CR 2032 coin cell batteries.
	4. Slide the door back into place until it clicks. Note: The next time you turn your meter on, it may prompt you to reset the time and date. (Refer to Section 9, Setting Up the Meter).

Note: When you no longer need the meter, remove the batteries and dispose of batteries and meter in compliance with your local government regulations.

The European Battery Directive requires separate collection of spent batteries, aiming to facilitate recycling and to protect the environment. The batteries in this product should be removed and disposed in accordance with local regulations for separate collection of spent batteries.

Cleaning Your Meter

Step	Action
1	Clean the outside of the meter using a damp cloth and: • Mild detergent/soap and water, or • 70% isopropyl alcohol, or • A mixture of 1 part household bleach, 9 parts water
2	Allow the meter to dry.

IMPORTANT: Do not place the meter in water or other liquids. Avoid getting dust, dirt, blood, control solution, water, or any other substance in the meter's test strip port and USB port.

1. Test strip is inserted in the strip port and nothing happens.

What It Means

Test strip is not inserted properly or fully into the meter

No batteries are installed; Batteries installed incorrectly

Dead batteries

Meter may be plugged into a computer (PC appears on meter screen)

Problem with the test strip

Problem with the meter

What To Do

1. With the contact bars (3 black lines) facing up, insert the test strip into the meter until it stops. This turns on the meter.
2. If the meter still does not turn on, contact Customer Care.

Refer to Section 17, Taking Care of Your Meter, on how to properly install batteries.

Dead batteries

Meter may be plugged into a computer (PC appears on meter screen)

Try a new test strip.

Contact Customer Care.

2. The test does not start after applying the blood sample.

Blood sample is too small

Sample applied after meter turns off

Problem with meter or test strip

1. See test strip instructions for use for re-application instructions.
2. Repeat the test using a new test strip.
3. If the test still does not start, contact Customer Care.

1. Review the testing instructions.
2. Repeat the test using a new test strip.
3. If the test still does not start, contact Customer Care.

1. Repeat the test using a new test strip.
2. If the test still does not start, contact Customer Care.

Introduction

This feature allows you to log insulin doses so they are recorded in the logbook. You can enable this feature at any time.

Set Up for Insulin Dose Logging

1. While on the Home screen, press  and hold for 3 seconds until the  appears. Insulin Dose Logging is now enabled.


Note: To turn off this feature, repeat this step.

How to Use Insulin Dose Logging

For long-acting insulin and/or insulin for breakfast, lunch, or dinner.

1. While on the Home screen, press 

morning	Long-acting insulin
evening	Long-acting insulin
breakfast	Meal insulin
lunch	Meal insulin
dinner	Meal insulin

3. Press  or  to enter the actual dose amount.

4. Press ☐ to log dose.
☒ means you took dose

2. Press or to enter the actual dose amount.

How to Log Additional Rapid-Acting Insulin Doses (e.g. snacks, bedtime correction, etc.)

1. While on the Home screen, press
 and hold  for 3 seconds until the screen changes.

3. Press ☐ to log dose.
☒ means you took dose

Assay method

Amperometry

Automatic shutoff

At least two minutes of inactivity

Battery life

Up to 3000 tests

Measurement range

For blood glucose testing 1.1 - 27.8 mmol/L
 For blood ketone testing 0.0 - 8.0 mmol/L

Memory

Up to 1000 events, including blood glucose, blood ketone, and control solution results, insulin doses, and other meter information

Minimum computer requirements

System must only be used with EN60950-1 rated computers. Use a USB certified cable

Operating relative humidity

10% to 90% (non-condensing)

Operating temperature

Meter: 10 °C to 50 °C (50 °F to 122 °F)
 System: See test strip instructions for use

Power source

Two CR 2032 lithium (coin cell) batteries

Size

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Storage temperature

-20 °C to 60 °C (-4 °F to 140 °F)

System altitude


See test strip instructions for use

Weight

33 to 37 g (1.2 to 1.3 oz.) including batteries

Note: For test strip specifications, see test strip instructions for use.

Electromagnetic Compatibility (EMC): FreeStyle Optium Neo meter has been tested for both electrostatic discharge and radio frequency interference. Emissions are low and unlikely to interfere with other nearby electronic equipment. To limit radio frequency interference do not use the FreeStyle Optium Neo meter near cellular or cordless telephones, radio transmitters or other electrical or electronic equipment that are sources of electromagnetic radiation, as these may interfere with the proper operation of the meter. Avoid use of the device in very dry environments, as electrostatic discharges from synthetic materials (e.g., carpets) could cause damage.

FreeStyle Optium Neo meter

Housing - PC-ABC Carbotex
Black colorant - Carbon Black
Button - colorless PC-ABC Carbotex
Main spring – AISI 304
size 9.5 cm $\pm$ 2 / 1.7 cm $\pm$ 2 weight 10,0 g $\pm$ 2%
Mediated contract through patient's arm skin

Principle of use of FreeStyle lancing device II together with Thin lancets.

FreeStyle lancing device II is intended for use with Thin lancets

Overview of the lancing device (refer to figures)

- (A) Depth adjusting disk
- (E) Lancet holder
- (I) Transparent cap (for sampling from alternative sites)
- (B) Depth indicator window
- (F) Puncture device cap
- (C) Puncture button
- (G) Lancet cap
- (D) Arming/triggering slider
- (H) Lancet

Important!

This FreeStyle lancing device II is intended for use by a single patient only and it is not suitable for use by health care professionals on multiple patients. Use of one device on multiple patients may result in infection and thus strictly prohibited.

Use of the FreeStyle lancing device II:

1. Remove cap from the lancing device (F) Break the cap off by bending it.
2. Insert the lancet (H) into the lancet holder (E). Insert the new lancet into the blue holder cylinder tightly (this may result in arming the device).
3. With a rotation movement, remove the lancet cap (G). Holding the puncture device firmly with one hand, unscrew the lancet cap with another hand.

4. Put on the puncture device cap (F)

Put the cap on so it fastened with a click. Take care not to touch the exposed needle of the lancet.

Note. Use dark cap for fingertip blood sampling and transparent cap for sampling from alternative sites.

5. Set the depth (A)

Puncture device can be adjusted for 9 different options of depth (including intermediate ones). Level 1 is for minimal depth, and 5 is for maximal depth. Minimal value is recommended for finger puncturing. For blood sampling from alternative sites, use transparent cap (I) with fixed depth value.

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Pull the grey slider back to click (the device may have been armed at stage 2 already).

Now, glucose level can be measured.

7. Choose an appropriate sampling site

a. Fingertip

- Press the puncture device (with dark cap on) to the side surface of your fingertip.
- Press the puncture button (C).
- If necessary, squeeze your finger to get a blood drop.
- More detailed info on analysis with use of your glucose level control system is given in the User's Manual.

b. Puncture in alternative sampling sites

More detailed info on alternative sampling sites (aside from fingertip) is given in the test strip instruction manual.

Warning!

Perform analysis with fingertip blood sample. Don't use alternative sampling sites in following cases:

- you suspect you have low glucose level (hypoglycemia) or low hypoglycemia insensitivity;
- you suspect that your glucose level changes quickly;
- results of analysis from samples taken from alternative sites don't match your well-being;
- less than two hours passed after meal, insulin administration or physical exercise.

If the set lacks transparent cap you need, contact the Customer Service.

- In order to make "fresh" blood get closer to the skin surface at the sampling site, rub this spot intensively for few second until you feel the warmth.

- Press the puncture device (with transparent cap on) tightly to the sampling site.
- Press the puncture button (C). Don't lift the device.

- Keep pressing the device against your skin, gradually increasing the pressure strength, for few seconds before the blood drop appears.
- Carefully lift the puncture device, try not to blur the blood sample. Note. Avoid puncturing moles, spots above veins, bones and sinews. Bruising may occur at the test site. In this case, consider choosing another sampling site.

8. Removing the lancet

- Remove the puncture device cap (F).
- Put the triggering slider down (D) to disconnect the lancet from the device.
- The lancet must be removed and disposed in according to local regulations on wastes disposal.

Taking Care of Your Meter

Clean the outside of the puncture device with an isopropyl (or soap water) alcohol cloth. If needed, the cap can be removed, washed with warm water and thoroughly rinsed. Keep the puncture device in the carry case or clean dry place.

Important!

Do not place the housing with puncture device in water or other liquids.

Thin lancets: AISI 304 material, PC-ABC Carbotex coating, no colorant

Thin lancets are used together with FreeStyle lancing device II in accordance with the procedure given above.

The devices makes short contact with capillary blood and immediate contact with skin - at the moment of puncture (automatic activation)

All finished products are non-toxic and comply with Regulation (EC) 1907/2006 by European Parliament and Council of Europe dated December 18, 2008 on Registration, Evaluation and Authorization of Chemicals (REACH). Devices are composed of stainless steel sharp and Polyethylene mold. Thin Lancets are irradiated for sterilization and the needle remains sterile until the plasticmolded seal of the twist-off cap is broken.

Carry case

Polyethylene terephthalate canvas Density of the polyethylene terephthalate canvas:

125 +/- 10 g/m²

Mediated contract through patient's arm skin

All materials are non-toxic and comply with Regulation (EC) 1907/2006 by European Parliament and Council of Europe dated December 18, 2008 on Registration, Evaluation and Authorization of Chemicals (REACH).

Abbott

Connecting cable

The USB cable is a yellow PVC and metal wire. All materials are non-toxic and comply with Regulation (EC) 1907/2006 by European Parliament and Council of Europe dated December 18, 2008 on Registration, Evaluation and Authorization of Chemicals (REACH).

Information on product materials contacting with human body (material, brand)

Parts description	Material
Top housing	PC-ABS Plastic acrylonitrile butadiene styrene, impact resistant industrial thermoplastic commercial resin based on acrylonitril copolymer with butadiene and styrene (name of the plastic is formed by first letters of the monomer names) is non-toxic and complies with Regulation (EC) 1907/2006 by European Parliament and Council of Europe dated December 18, 2008 on Registration, Evaluation and Authorization of Chemicals (REACH). Device makes short contact with intact skin.
Bottom housing	Housing - PC-A BC Carbotex black colorant – Carbon Black, impact resistant industrial thermoplastic commercial resin based on acrylonitril copolymer with butadiene and styrene (name of the plastic is formed by first letters of the monomer names) is non-toxic and complies with Regulation (EC) 1907/2006 by European Parliament and Council of Europe dated December 18, 2008 on Registration, Evaluation and Authorization of Chemicals (REACH). Device makes short contact with intact skin.
Battery door	Housing - PC-A BC Carbotex black colorant - Carbon Black, impact resistant industrial thermoplastic commercial resin based on acrylonitril copolymer with butadiene and styrene (name of the plastic is formed by first letters of the monomer names) is non-toxic and complies with Regulation (EC) 1907/2006 by European Parliament and Council of Europe dated December 18, 2008 on Registration, Evaluation and Authorization of Chemicals (REACH). Device makes short contact with intact skin.
Strip port	PC-ABC Carbotex strip window. Used in medical devices, non-toxic and complies with Regulation (EC) 1907/2006 by European Parliament and Council of Europe dated December 18, 2008 on Registration, Evaluation and Authorization of Chemicals (REACH). Device makes short contact with intact skin.
Back labels	Thermo and Lexan (polycarbonate high-duty resin) is non-toxic and complies with Regulation (EC) 1907/2006 by European Parliament and Council of Europe dated December 18, 2008 on Registration, Evaluation and Authorization of Chemicals (REACH). Device makes short contact with intact skin.
Buttons	PC-ABC Carbotex colorless

 Abbott

Parts description	Material
	Main spring -A1S1304 are non-toxic and comply with Regulation (EC) 1907/2006 by European Parliament and Council of Europe dated December 18, 2008 on Registration, Evaluation and Authorization of Chemicals (REACH). Device makes short contact with intact skin.
Display	Display - Lexan polycarbonate matte film E-paper display and e-ink are non-toxic and comply with Regulation (EC) 1907/2006 by European Parliament and Council of Europe dated December 18, 2008 on Registration, Evaluation and Authorization of Chemicals (REACH). Device makes short contact with intact skin.
FreeStyle Lancing Device II	Housing - PC-ABC Carbotex Black colorant - Carbon Black Button - colorless PC-ABC Carbotex Main spring – AISI 304 size 9.5 cm ± 2 / 1.7 cm ± 2 weight 10,0 g $\pm 2\%$ Device makes short contact with capillary blood and mediated skin contact during puncturing (automatic activation)
Thin lancets	Thin lancets are made of AISI 304, covered with PC-ABC Carbotex, no colorant. All materials are non-toxic and comply with Regulation (EC) 1907/2006 by European Parliament and Council of Europe dated December 18, 2008 on Registration, Evaluation and Authorization of Chemicals (REACH). Thin Lancets are irradiated for sterilization and the needle remains sterile until the plasticmolded seal of the twist-off cap is broken. Device makes short contact with capillary blood and mediated skin contact during puncturing (automatic activation).
Carry case	Polyethylene terephthalate canvas. Density of the polyethylene terephthalate canvas: 125 +/- 10 g/m2 Mediated contract through patient's arm skin

Environmental Protection

The FreeStyle Optium Neo Blood Glucose and Ketone Monitoring System complies with the EU Directive 2002/95/CE on Restriction of Hazardous Substances in electrical and electronic equipment: lead (Pb), cadmium (Cd), mercury (Hg), hexavalent chromium (Cr VI), polybrominated biphenyls (PBB) and polybromodiphenylether (PBDE)

FreeStyle Optium Blood Glucose Test Strips (10, 25, 50, 100 pcs), UK, Registration certificate No. 2012/13309 (if required).

- Basis: thermoplastic film of Coveme polyester.
- Conductive electrodes layer: conductive carbonaceous paste C2130514 D10 for chemical system support.
- Insulation layer: operation ink Erson E6110-116
- Ferments layer (concentration per strip):
- Glucose dehydrogenase (pseudomonadaceae species GDH-NAD) > 0.03 units
 - NAD+ (as potassium salt) > 0.1 µg
 - Phenanthroline quinone ≥ 0.02 µg
 - Nonreactive ingredients ≥ 16.3 µg
- Upper ribbon:
 - Transparent Mylan A (FLEXPET) polyester ribbon
 - White Melinex (Dupont) polyester separator
 - Adhesive layer based on Aroset PS 5308 and PS 5307 acrylic polymer
 - Colorant: black QLT4726Quartz Ultra C

Package: Thermally fixed PET 23-16 film, polyethylenterephthalate, continuous adhesive layer structure, with no spaces and breaches, no adhesion defects. Glue layer (width) 3.0-5.0 mm Gluing strength (opening force) 40-80 N.

When opening the package, the separation is continuous, uniform, homogenous, no peeling effect.

No material ruptures.

There are no colorant spots under the sample.

Contact with damaged skin

Specifications

Item	Value
Calibration	By plasma
Blood sample size	0.6 µl
Measurement time	5 seconds
Additional time for blood sample application	5 seconds
Sufficient volume confirmation signal	Yes
Measurement range	1.1 to 27.8 mmol/l
Measurement method	Amperometry
Hematocrit	30-60%
Storage temperature	4-30 °C
Shelf life	18 months

Performance characteristics

Performance characteristics of the FreeStyle Optium Blood Glucose Test Strips were estimated during laboratory and clinical tests.

Analysis limits: 1.1–2.8 mmol/l (20–500 mg/dl)

Analysis time: 5 seconds

Calibration reference

FreeStyle Optium Blood Glucose Test Strips were calibrated by glucose analyzer YSI (by YSI Incorporated)

Readings of the glucose levels in whole blood obtained from YSI analyzer are multiplied by 1.12 (plasma equivalent) to get glucose levels equal to plasma for calibration of FreeStyle Optium Blood Glucose Test Strips.

Repeatability

Repeatability tests demonstrated that results variability is within 3.0-4.4 %.

See Table 1 and Table 2.

Accuracy

Measurements of the capillary blood glucose level, obtained by specially trained staff, were compared to the results obtained with YSI glucose analyzer.

Accuracy of the system results obtained at glucose concentration < 5.55 mmol/l (100 mg/dl) with the FreeStyle Optium Neo Blood Glucose and Ketone Monitoring System.

Within ± 0.28 mmol/l (5 mg/dl)

105 / 154 (68.2 %)

Within ± 0.56 mmol/l (10 mg/dl)

149 / 154 (96.8 %)

Within ± 0.83 mmol/l (15 mg/dl)

154 / 154 (100.0 %)

Accuracy of the system results obtained at glucose concentration ≥ 5.55 mmol/l (100 mg/dl) with the FreeStyle Optium Neo Blood Glucose and Ketone Monitoring System.

Within ± 5 %

410 / 632 (64.9 %)

Within ± 10 %

581 / 632 (91.9 %)

Within ± 15 %

625 / 632 (98.9 %)

Accuracy of the system results obtained at glucose concentration from 1.51 mmol/l (29 mg/dl) to 24,3 mmol/l (438 mg/dl) with the FreeStyle Optium Neo Blood Glucose and Ketone Monitoring System.

Within ± 0.83 mmol/l (15 mg/dl) or 15%

779 / 786 (99.1 %)

Study of the fingertip capillary blood glucose sample results from 165 patients demonstrated following results: 98% of results were within ± 0.83 mmol/l (15 mg/dl) from reference values obtained from the YSI glucose analyzer, with glucose concentrations below 5.55 mmol/l (100 mg/dl), and 98% of results were within ± 15 % from reference values obtained from the YSI glucose analyzer, with glucose concentrations below 5.55 mmol/l (100 mg/dl).

Accuracy of the system results obtained at glucose concentration < 5.55 mmol/l (100 mg/dl) with the FreeStyle Optium Neo Blood Glucose and Ketone Monitoring System, with such meters as: Optium, Optium Xceed and Optium Easy.

Within ± 0.28 mmol/l (5 mg/dl)

125 / 222 (56.3 %)

Within ± 0.56 mmol/l (10 mg/dl)

194 / 222 (87.4 %)

Within ± 0.83 mmol/l (15 mg/dl)

219 / 222 (98.6 %)

Accuracy of the system results obtained at glucose concentration ≥ 5.55 mmol/l (100 mg/dl) with the FreeStyle Optium Neo Blood Glucose and Ketone Monitoring System, with such meters as: Optium, Optium Xceed and Optium Easy.

Within $\pm 5\%$	Within $\pm 10\%$	Within $\pm 15\%$
455 / 751 (60.6 %)	655 / 751 (87.2 %)	730 / 751 (97.2 %)

Accuracy of the system results obtained at glucose concentration from 1.8 mmol/l (33 mg/dl) to 27.4 mmol/l (494 mg/dl) with the FreeStyle Optium Neo Blood Glucose and Ketone Monitoring System, with such meters as: Optium, Optium Xceed and Optium Easy.

Within ± 0.83 mmol/l (15 mg/dl) or 15%	949 / 973 (97.5 %)
--	--------------------

Study of the fingertip capillary blood glucose sample results from 120 patients demonstrated following results: 95% of results were within ± 15 from reference values obtained from the YSI glucose analyzer, with glucose concentrations below 5.55 mmol/l (100 mg/dl), and within ± 0.83 mmol/g (15 mg/dl) with glucose concentrations below 5.55 mmol/l (100 mg/dl). 100% of results were within ± 0.83 mmol/l (15 mg/dl) from reference values obtained from the YSI glucose analyzer, with glucose concentrations below 5.55 mmol/l (100 mg/dl), and 94% of results were within ± 15 from reference values obtained from the YSI glucose analyzer, with glucose concentrations below 5.55 mmol/l (100 mg/dl).

Table 1 - Result repeatability estimation

	Low value	Moderate low value	Moderate high value	High value
Mean value, mmol/l	2.63 (47.3)	4.82 (86.9)	8.09 (145.8)	18.48 (332.9)
Standard deviation, mmol/l (mg/dl)	0.10 (1.9)	0.15 (2.7)	0.25 (4.5)	0.55 (10.0)
Variation coefficient, %	4.0	3.1	3.1	3.0

Table 2 - Intralaboratory repeatability

	Low value	Average value	High value
Mean value, mmol/l	2.40 (43.2)	5.09 (91.7)	16.23 (292.5)
Standard deviation, mmol/l (mg/dl)	0.11 (1.9)	0.18 (3.3)	0.50 (9.0)

 **Abbott**

Variation coefficient, %

4.4

3.6

3.1

FreeStyle Optium Blood β -Ketone Test Strips (10, 25, 50, 100 pcs), Registration certificate No. 2012/13309

Determination principle

When applying blood on the test strip, blood β -hydroxybutyrate reacts with test strip reagents, generating a weak electric current. Current rate is measured, and the result is shown on the FreeStyle Optium Neo Blood Glucose and Ketone Monitoring System display. Current rate depends on amount of β -hydroxybutyrate in the blood sample.

Measurement method Amperometry

Basis thermoplastic film of Coveme polyester.

- Conductive electrodes layer: conductive carbonaceous paste C2130514 D10 for chemical system support
- Insulation layer: operation ink Erson E6110-116
- Ferments layer (concentration per strip):
- β -oxybutyrate dehydrogenase (pyocyanic MMW peptide) ≥ 0.03 units
- NAD (free acid form) $\geq 1.67 \mu\text{g}$
- Phenanthroline quinone $\geq 0.29 \mu\text{g}$
- Non-reacting components $\geq 19.51 \mu\text{g}$

Upper ribbon:

- Transparent Mylan A (FLEXPET) polyester ribbon
- White Melinex (Dupont) polyester separator
- Adhesive layer based on Aroset PS 5308 and PS 5307 acrylic polimer
- Colorant: black QLT4726Quartz Ultra C
- Package: Thermally fixed PET 23-16 film

Polyethyleneterephthalate, continuous adhesive layer structure, with no spaces and breaches, no adhesion defects.

Glue layer (width) 3.0-5.0 mm

Gluing strength (opening force) 40-80 N.

When opening the package, the separation is continuous, uniform, homogenous, no peeling effect. No material ruptures. There are no colorant spots under the sample.

Contact with damaged skin

Determined values range: 0.0 - 8.0 mmol/l

Analysis time: 10 seconds

Additional time for blood sample application: 30 seconds

Results repeatability: With ketone concentration below 1.5 mM, results repeatability does not exceed 0.15 mM.
With ketone concentration above 1.5 mM, results repeatability does not exceed 10% of the variation coefficient

Blood sample volume: 1.5 μ l

Hematocrit value range 30-60%

Re-Application for the ketone level measurement: after 30 seconds

Calibration reference:

FreeStyle Optium Blood β -Ketone Test Strips are calibrated to show β -hydroxybutyrate concentration in plasma with Randox analyzer (RB1007).

Precision

Studies demonstrated that the measuring inaccuracy is within 3.1–3.8%. The results were obtained with laboratory study of the venous whole blood sample (n = 20 for each level). See Table 1.

Accuracy

Accuracy determination revealed comparability of the results obtained by trained operators and users with no special skills. The accuracy was studied in three clinics and four hospitals by comparing the whole blood β -hydroxybutyrate level measurements with measurements obtained from study of blood plasma in reference laboratory devices. See Table 2.

These studies demonstrated that data obtained with FreeStyle Optium Neo, FreeStyle Optium, Optium and Optium Xceed systems are comparable enough with the results obtained by reference laboratory method.

At least 95% of the results are within A and B areas of the Natrass inaccuracy grid.

Table 1 – Repeatability			
	Low	Average	High
Mean value, mmol/l	0.34	2.36	6.32
Standard deviation, mmol/l	0.03	0.09	0.2
Variation coefficient, %	--	3.8	3.1

Table 2 – Accuracy	
Samples quantity	288
Ketone concentration range, mmol/l	0.07-5.2
Slope	1.06
Ordinate axis crossing point, mmol/l	0.07
Correlation coefficient r	0.98

Conditions of Practical Application

Diabetes mellitus is a chronic disorder of glucose metabolism characterized by high bloodglucose levels. Clinically, it is associated with serious long-term complications. Currently, diabetes is incurable and requires lifelong management including significant ongoing lifestyle changes, self monitoring and the use of oral and injectable medication. There are two maintypes of diabetes mellitus: type 1 diabetes, characterized by a complete inability of the body to produce insulin, and type 2 diabetes, caused by a combination of insulin resistance and relative insulin deficiency.

Based on evidence from landmark randomized clinical trials, the aim of modern diabetes management is to achieve tight glycemic control in order to reduce the risk of serious long-term complications.

For patients with type 1 diabetes mellitus (T1DM), treatment involves self-monitoring of blood glucose (SMBG) and daily insulin administration. Although the majority of patients with type 2 diabetes mellitus (T2DM) are managed with lifestyle changes and oral agents (OAs) many people will require insulin at some stage. Diabetes has a significant impact on the patient, requiring lifetime management, regular healthcare appointments and substantial lifestyle changes.

The monitoring of blood glucose concentrations has always been a valuable tool for use in the diagnosis, treatment, and control of diabetes mellitus. Monitoring of blood glucose at home by patients with Type 1 and Type 2 diabetes (insulin dependent and non-insulin dependent diabetes mellitus) has become a routine part of diabetes therapy. The Diabetes Control and Complications Trial (DCCT) 5 recommended that Type 1 diabetics self test at least 3-4 times a day, while Type 2 diabetics should self-test at least once a day.

Blood Glucose Monitoring for Insulin dependent patients

For patients treated with insulin there is general agreement that self monitoring of blood glucose levels is a cornerstone of modern diabetes care. Regular Self-monitoring of blood glucose (SMBG) provides the patient with a real-time measurement of their blood glucose level. This information obtained at regular intervals can be used to:

- To assess the patient's current glycemic status
- To determine the impact of changes in lifestyle and medical therapies
- As an early warning system for the detection of low and high blood glucose levels
- To confirm the diagnosis of hypoglycemia
- To provide information to HCPs to optimize therapy (including basal insulin)
- To calculate an appropriate rapid acting insulin dose

There is consensus that SMBG is an integral part of modern Intensive Insulin Therapy. Self-monitoring of blood glucose (SMBG) is widely recognized as a core component of effective diabetes self-management particularly for patients using insulin. Ideally, the FreeStyle Optium Neo Blood Glucose and Ketone Monitoring System and test strips are chosen to suit individual needs.

According to International Guidelines, SMBG should be performed at least three times per day for patients using multiple daily injections. However, the optimal frequency of self monitoring should depend on the characteristics of an individual's blood glucose control, the insulin treatment regimen and personal preferences, for example the desired lifestyle.

Ketone (β -oxybutyrate) testing provides an additional tool for good diabetes management. When the body cannot process glucose for energy (this happens with diabetes) or if there is little or no glucose, it must get energy from another source. When this happens, the body processes the fat from the fat stores (in the body) for energy. During this process, ketones form in the blood. Blood ketones are another indicator that helps patients and their healthcare professionals determine if the patient's diabetes is out of control. Ketone testing is recommended.

- During illness (i.e., cold, flu or an infection)
- When blood glucose results are higher than the range set by the healthcare professional
- When patients are under a lot of stress
- During pregnancy

When a blood glucose result is 13.3 mmol/L (240 mg/dL) or higher, the meter will display "KET". It is suggested that the user (or healthcare professional) should then check the level of ketones in the blood, and seek medical advice.

The Diabetes Monitoring System is also capable of measuring β -hydroxybutyrate (β -OHB), which is the predominant ketone body in diabetic ketoacidosis (DKA). The American Diabetes Association recommends that diabetics should test the level of β -OHB in the blood when blood glucose levels are elevated (>240 mg/dL or >13.3 mmol/L). DKA is a serious complication of diabetes that requires immediate attention; any delay in the diagnosis and treatment of DKA may lead to an increase in morbidity and mortality.

Insulin Therapy and challenges with insulin dose estimations. Insulin can only be administered by subcutaneous injection or infusion. There are multiple different insulin preparations and regimens used in clinical practice. For patients who require insulin, accurate dosing is essential to maintain glucose levels within their target range in order to prevent disease related complications. At the same time, inaccurate insulin dosing can lead to problems with low blood glucose levels (hypoglycemia). For all patients with T1DM and an increasing proportion of patients with T2DM, modern management involves the use Multiple Daily Injection (MDI) therapy. This regimen requires patient education and intensive self monitoring of blood glucose.

Blood Glucose Monitoring for Insulin dependent patients

For patients treated with insulin there is general agreement that self monitoring of blood glucose levels is a cornerstone of modern diabetes care.¹ Regular Self-monitoring of blood glucose (SMBG) provides the patient with a real-time measurement of their blood glucose level. This information obtained at regular intervals can be used to:

- To assess the patient's current glycemic status

¹ ADA 2014 "Standards of medical care in diabetes 2014" Diabetes care 37 Suppl 1:S11-S24

- To determine the impact of changes in lifestyle and medical therapies
- As an early warning system for the detection of low and high blood glucose levels
- To confirm the diagnosis of hypoglycemia
- To provide information to HCPs to optimize therapy (including basal insulin)
- To calculate an appropriate rapid acting insulin dose

There is consensus that SMBG is an integral part of modern Intensive Insulin Therapy. Self-monitoring of blood glucose (SMBG) is widely recognized as a core component of effective diabetes self-management particularly for patients using insulin. Ideally, the FreeStyle Optium Neo Blood Glucose and Ketone Monitoring System and test strips are chosen to suit individual needs.²

According to International Guidelines, SMBG should be performed at least three times per day for patients using multiple daily injections. However, the optimal frequency of self monitoring should depend on the characteristics of an individual's blood glucose control, the insulin treatment regimen and personal preferences, for example the desired lifestyle⁷.

Ketone (β -oxybutyrate) testing provides an additional tool for good diabetes management. When the body cannot process glucose for energy (this happens with diabetes) or if there is little or no glucose, it must get energy from another source. When this happens, the body processes the fat from the fat stores (in the body) for energy. During this process, ketones form in the blood. Blood ketones are another indicator that helps patients and their healthcare professionals determine if the patient's diabetes is out of control. Ketone testing is recommended³.

- During illness (i.e., cold, flu or an infection)
- When blood glucose results are higher than the range set by the healthcare professional
- When patients are under a lot of stress
- During pregnancy

When a blood glucose result is 13.3 mmol/L (240 mg/dL) or higher, the meter will display "KET". It is suggested that the user (or healthcare professional) should then check the level of ketones in the blood, and seek medical advice.

The Diabetes Monitoring System is also capable of measuring β -hydroxybutyrate (β -OHB), which is the predominant ketone body in diabetic ketoacidosis (DKA). The American Diabetes Association recommends that diabetics should test the level of β -OHB in the blood when blood

glucose levels are elevated (>240 mg/dL or >13.3 mmol/L). DKA is a serious complication of diabetes that requires immediate attention; any delay in the diagnosis and treatment of DKA may lead to an increase in morbidity and mortality.⁴

Insulin Therapy and challenges with insulin dose estimations. Insulin can only be administered by subcutaneous injection or infusion. There are multiple different insulin preparations and regimens used in clinical practice. For patients who require insulin, accurate dosing is essential to maintain glucose levels within their target range in order to prevent disease related complications. At the same time, inaccurate insulin dosing can lead to problems with low blood glucose levels (hypoglycemia). For all patients with T1DM and an increasing proportion of patients with T2DM, modern management involves the use Multiple Daily Injection (MDI) therapy. This regimen requires patient education and intensive self-monitoring of blood glucose.⁵

Methods and means of sterilization

Clean and sterilize measurement window and its cover after used and at least once a day.

Cotton buds, swabs or fabric cloth slightly dampen in disinfecting solution can be used for sterilization. Use mixture of propanol-1, propanol-2 and glutaric dialdehyde, also known in some countries as Bacillol plus. Wait until wiped surfaces dry completely. Be careful, so that no liquid got into the device.

CAUTION!!! Never use organic solvents, don't spray the device, don't put it in disinfecting solutions.

⁴

Ketoacidosis (DKA), Living Well with Diabetes. American Diabetes Association, 1701 North Beauregard Street, Alexandria, VA 22311 1-800-DIABETES Copyright 1995-2013, American Diabetes Association.

⁵

NICE U.K. (2004). NICE Type 1 Diabetes. London: National Institute for Clinical Excellence.

Symbol	Meaning	Symbol	Meaning
	Refer to the Instructions for Use		Warning
	Temperature Limits		Expiry date
	Manufacturer		In-vitro diagnostics medical device
	Lot code		Catalog No.
	For single use only		For recycling
	Manufacture date		Serial No.
	Don't drink		Radiation-sterilized (lancets only)
	CE Symbol		
	Don't use if the package is damaged.		

Registration certificate No. dated by

 Don't expose to direct sunlight

An evaluation was completed to ensure the effect of short-term extreme temperature on the test strips which are contained within the FreeStyle Optium Neo Blood Glucose and Ketone Monitoring System was within specification. The results from the transit storage testing met the acceptance criteria as shown in tables below, and therefore demonstrate acceptable performance of the G3cP test strips when exposed to short term extremes of temperature. Therefore the acceptance criteria for transit storage testing have been met.

Mean YSI Reference		Mean Bias Difference from Control (22°C)	Stress Temperature (°C)			
mmol/l	mg/dl		-20	40	60	-20 and 60
2.6	46	mmol/l	0.00	-0.08	-0.25	-0.24
		mg/dl	-0.1	-1.5	-4.5	-4.4
6.4	114	%	-1.8	-2.5	-7.9	-8.5
22.4	405	%	0.0	-2.9	-8.6	-8.8

KEEP AWAY FROM CHILDREN!

Parameter	Value
Storage temperature	The FreeStyle Optium Neo Blood Glucose and Ketone Monitoring System: - 20°C to 60°C FreeStyle Optium Blood β -Ketone Test Strips: 18 °C to 30 °C FreeStyle Optium Blood Glucose Test Strips: 15 °C to 40 °C
Operating temperature	The FreeStyle Optium Neo Blood Glucose and Ketone Monitoring System - 10°C to 50°C FreeStyle Optium Blood β -Ketone Test Strips: 4 °C to 30 °C FreeStyle Optium Blood Glucose Test Strips: 4 °C to 30 °C
Operating humidity	10 to 90 % (without condensing)
Pressure	50 hPa to 106 hPa

The storage temperature range of the **FreeStyle Optium Neo Blood Glucose and Ketone Monitoring System** is 4 °C to 30 °C, for the set includes FreeStyle Optium Blood Glucose Test Strips with limited shelf life.

Manufacturer's Warranties

Abbott Diabetes Care Inc. ("Abbott") warrants that the the FreeStyle Optium Neo Blood Glucose and Ketone Monitoring System has no shelf life provided that the recommendations by manufacturer on application and maintenance service (if required) are observed.

Limited Warranty

Abbott Diabetes Care Inc. ("Abbott") warrants that the FreeStyle Optium Neo Blood Glucose and Ketone Monitoring System shall be free from defects in material and workmanship for a period of two (2) years from the original date of purchase or delivery date, whichever is later,

provided it is not modified, altered or misused. Under this Limited Warranty, if the FreeStyle Optium Neo Blood Glucose and Ketone Monitoring System is defective in material or workmanship, Abbott's sole obligation is to replace the FreeStyle Optium Neo Blood Glucose and Ketone Monitoring System, free of charge, with a same or similar device as determined by Abbott in its sole discretion. Abbott warrants the performance of the device in accordance with its specifications if it is used as directed and provided that the failure to perform or misperformance of the device has not been caused in whole or in part by the use test strips that are not the test strips manufactured by Abbott. This Limited Warranty covers only the FreeStyle Optium Neo Blood Glucose and Ketone Monitoring System and shall not apply to auxiliary equipment or disposable accessories, extends only to the original purchaser, and is not assignable or transferable. Please see the User's Manual/Owner's Setup Guide for Instructions for Use, distributor and manufacturer information. For warranty service, contact Customer Care for assistance and/or instructions for obtaining a replacement device. See User's Manual/Owner's Setup Guide for your Customer Care contact information. Abbott may require as a condition to obtaining warranty service that you return your FreeStyle Optium Neo Blood Glucose and Ketone Monitoring System, postage prepaid, to an address specified by Customer Care. To the extent possible under law, these are Abbott's complete express warranties for FreeStyle Optium Neo Blood Glucose and Ketone Monitoring System and state your exclusive remedies. Abbott make no other warranties, express or implied, and Abbott excludes and disclaims any other warranties including, but not limited to the implied warranties of merchantability and fitness for a particular purpose and any other remedies including, but not limited to compensations under contract or tort laws. Abbott does not warrant that operation of the FreeStyle Optium Neo Blood Glucose and Ketone Monitoring System will be uninterrupted or error free.

Limitation of Liability

To the extent possible under law, Abbott will not be liable for any incidental or consequential damages resulting directly or indirectly from use of the FreeStyle Optium Neo Blood Glucose and Ketone Monitoring System or its failure to perform in accordance with specifications. Some jurisdictions do not allow the exclusion or limitation of other express or implied warranties or incidental or consequential damages, so the above limitations or exclusions may not apply to you.

Service Center

Abbot Laboratories OOO, 125171, Moscow, Leningradskoe shosse, 16A, bldg. 1

Tel/Fax: Tel/Fax: +7 (495) 258 42 80 www.abbottdiabetes.ru

Claim Department

Abbot Laboratories OOO, 125171, Moscow, Leningradskoe shosse, 16A, bldg. 1

Tel/Fax: Tel/Fax: +7 (495) 258 42 80 www.abbottdiabetes.ru

Disposal Instructions for the Russian Federation

Disposal of the CR 2032 lithium batteries -

Spent batteries may not be disposed along with domestic waste.

Disposal of the FreeStyle Optium Neo Blood Glucose and Ketone Monitoring System, FreeStyle lancing device II within the Russian Federation is performed in accordance with Sanitary Rules and Regulations 2.1.7.2790-10, medical devices fall within Hazard Class A. They pertain to epidemiologically safe waste close by composition to solid domestic waste and should be disposed by an accredited organization in accordance with current Russian legislation.

Disposal of FreeStyle Optium Blood Glucose Test Strips, FreeStyle Optium Blood β -Ketone Test Strips, Thin lancets in the Russian Federation is performed in accordance with Sanitary Rules and Regulations 2.1.7.2790-10, these medical devices fall within Hazard Class B. They pertain to epidemiologically hazardous waste (infected and potentially infected waste).

Authorized representative of the manufacturer in the Russian Federation: Abbot Laboratories OOO, 125171, Moscow, Leningradskoe shosse, 16A, bldg. 1 Tel/Fax: +7 (495) 258 42 80 www.abbottdiabetes.ru

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
infab pyx to

NOTARIAL CERTIFICATE

LET IT KNOWN that I, Esmond Geoffrey John Jenkins of Old Manor Barn,
Shorncliffe, Cirencester in the County of Gloucestershire, England, a duly authorised
Notary Public, CERTIFY AS FOLLOWS:

That the paper writings annexed hereto is an original document produced to
me and signed by **Archana Sabananthan** the UK Regulatory Project
Manager and duly authorised signatory of Abbott Diabetes Care Limited and
who is identified to me by Australian Passport number N7480193.

SIGNED and sealed Bampton Law LLP Solicitors Market Square, Bampton in the
Country of Oxfordshire, England on 7th August 2017

Esmond Geoffrey John Jenkins

Notary Public
England and Wales


Protocol 17-389

Esmond G J Jenkins - Notary Public

Old Manor Barn
Shorncliffe
Cirencester, County of Gloucestershire
GL7 6DE, England

In this Owner's Guide only principles of use of FreeStyle Lancing Device II are specified.

This FreeStyle Lancing Device II is intended for use together with FreeStyle lancets or Thin lancets. To get additional information on your glucose monitoring system and instruction for measuring blood glucose please refer to the User's Manual.

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Important!

This FreeStyle Lancing Device II is intended for use by a single patient only and it is not suitable for use by health care professionals on multiple patients. Use of one device on multiple patients may result in infection and thus strictly prohibited.

Use of the FreeStyle Lancing Device II:

1. Remove cap from the lancing device (F). Break the cap off by bending it.

Archana Sabananthan
Regulatory Project Manager
Abbott Diabetes Care Ltd

ABBOTT

07/AUG/2017

2. Insert the lancet (H) into the lancet holder (E). Insert the new lancet into the blue holder cylinder tightly (this may result in arming the device).
3. With a rotation movement, remove the lancet cap (G). Holding the puncture device firmly with one hand, unscrew the lancet cap with another hand.
4. Put on the puncture device cap (F)
Put the cap on so it fastened with a click. Take care not to touch the exposed needle of the lancet.

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Puncture in alternative sampling sites

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- In order to make "fresh" blood get closer to the skin surface at the sampling site, rub this spot intensively for few second until you feel the warmth.
- Press the puncture device (with transparent cap on) tightly to the sampling site.
- Press the puncture button (C). Don't lift the device.

- Keep pressing the device against your skin, gradually increasing the pressure strength, for few seconds before the blood drop appears.

- Carefully lift the puncture device, try not to blur the blood sample. Note. Avoid puncturing moles, spots above veins, bones and sinews. Bruising may occur at the test site. In this case, consider choosing another sampling site.

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