Medtronic

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I, Sowmya Shetty, hereby swear (or affirm) that the attached Instructions for Use – Endurant II/Endurant IIs Stent Graft System is a true copy of the original document.

S Shetty
Sowmya Shetty
3850 Brickway Blvd. BWB2
Santa Rosa, CA 95403

17 Mar 2015
Date


COUNTRY: USA

State of California
ss.
County of Sonoma

Subscribed and sworn to (or affirmed) before me on this March 17, 2015, by Sowmya Shetty proved to me on the basis of satisfactory evidence to be the person who appeared before me.

A notary public or other officer completing this certificate verifies only the identity of the individual who signed the document to which this certificate is attached, and not the truthfulness, accuracy, or validity of that document.

Sandra Frieze, Notary Public
Sandra Frieze, Notary Public

Endurant™ II

Stent Graft System
Стент-графт система
支架系统
Sustav stent-grafta
Systém stentgraftu
Stentgraftsystem
Stentgraftsystem
Stenttisiirrajestelmä
Système d'endoprothèse
Gefäßprothesensystem
Σύστημα μασχέματος stent
Sztentgraftrendszer
Sistema di endoprotesi
Stentgraftsystem
System stentgraftu
Sistema de endoprótese
Sistem de grefa de tip stent
Система стент-графта
Stent graft system
Système stentgraftu
Sistema de endoprótesis
Stentgraftsystem
Stent Graft Sistemi

Instructions for Use • Инструкции за употреба • 使用説明
Upute za upotrebu • Návod k použití • Brugsanvisning
Gebruiksaanwijzing • Käyttöohjeet • Mode d'emploi
Gebrauchsanweisung • Οδηγίες χρήσης • Használati útmutató
Instruzioni per l'uso • Bruksanvisning • Instrukcja użytkowania
Instruções de utilização • Instrucţiuni de utilizare
Инструкция по эксплуатации • Uputstva za upotrebu
Pokyny na používanie • Instrucciones de uso • Bruksanvisning
Kullanım Talimatları

2014/JUL/21 at 4:29 p.m. Doc number: M052196T001 [multi23]

multilang_23:20140310, alldocMod:20140310, alldocStyle:20131127, master_pageset:20131127,
languageLookup:20130312, table:20140623, characterJoin:20080301, 3of9:20071016, 3of9_svg:
20071610
printspec:a6

Trademarks may be registered and are the property of their respective owners.

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Stent Graft System

TABLE OF CONTENTS

1. Device Description	8
2. Indications for Use	10
3. Contraindications	11
4. Warnings and Precautions	11
5. Adverse Events	12
6. Patient Selection and Treatment	13
7. Patient Counseling Information	13
8. How Supplied	13
9. Clinical Use Information	13
10. Implant Instructions	15
11. Bail-Out Techniques	23
12. Follow-up Imaging Recommendations	23
13. Additional Surveillance and Treatment	24

1. Device Description

The Endurant™ II/Endurant™ IIs stent graft system is designed for the endovascular repair of infrarenal abdominal aortic or aortoiliac aneurysms. When placed within the target lesion, the stent graft provides a permanent, alternative conduit for blood flow within the patient's vasculature by excluding the lesion from blood flow and pressure.

The stent graft system is comprised of 2 main components: the implantable stent graft and the disposable delivery system. The stent graft is preloaded into the delivery system and advanced to the aneurysm using fluoroscopic guidance. Upon deployment, the stent graft self-expands to conform to the shape and size of the seal zones above and below the aneurysm.

1.1. Stent Graft

The Endurant™ II/Endurant™ IIs stent graft (Figure 1) has 2 basic configurations: a bifurcated configuration and a limb configuration. Additional configurations include iliac extension, aortic extension, abdominal tube, and aorto-uni-iliac (AUI). After placement of the bifurcated or AUI device, limbs and additional stent grafts are introduced separately into the vessel and mated with the implanted configuration.

All configurations are composed of nitinol stents sewn to a fabric graft with nonresorbable sutures. Radiopaque markers are sewn onto the stent graft to aid in visualization and to facilitate accurate placement. The nitinol stents may also be visible under fluoroscopy.

Stent graft components should be oversized to be larger than the measured vessel inner diameter (Section 9.2). Table 1 contains a summary of the stent graft materials.


Figure 1. Stent Graft Configurations and Locations of RO Markers

1. Radiopaque Marker
2. "e" Marker
3. Radiopaque Gate Marker
4. Endurant II Aortic Extension/Abdominal Tube Configuration
5. Endurant II Bifurcated Configuration
6. Endurant IIS Bifurcated Configuration
7. Endurant II Iliac Extension Configuration
8. Endurant II Limb Configuration
9. Endurant II Aorto-Uni-Iliac Configuration
10. Overlap Marker

Note: This and all other product graphics appearing in this manual are not drawn to scale.

Table 1. Stent Graft Materials

Component	Material
Stents	Nickel-Titanium (Nitinol) Alloy
Button Radiopaque Markers	Platinum-Iridium Alloy
"e" Radiopaque Marker	Platinum
Contralateral Gate Marker	Platinum-Iridium Alloy
Graft Material	Polyester
Suture	Polyester and Polyethylene

The Endurant II/Endurant IIS stent graft system does not contain natural rubber latex. However, during the manufacturing/assembly process, it may have incidental contact with latex-containing products.

1.1.1. Bifurcated Configuration

The bifurcated stent graft is available in 2 configurations, the Endurant II bifurcated configuration and the Endurant IIS bifurcated configuration. The Endurant II bifurcated configuration is an aortoiliac stent graft that is available in 3 lengths. The Endurant IIS bifurcated configuration is an aortic configuration available in a single, shorter length (Figure 1). The proximal end of both bifurcated configurations deploy into the proximal neck and upper section of the aneurysm. The proximal end of the bifurcated configuration is composed of nitinol stents sewn to a fabric graft. The suprarenal portion of the proximal end is not covered with graft fabric (Figure 1). The suprarenal stent also contains anchor pins to fix the stent graft in place inside the aorta.

The aortic section bifurcates distally into 2 smaller tubes: an ipsilateral leg and a shorter contralateral leg. In the Endurant II bifurcated configuration, all stents on the ipsilateral leg are sewn to the outside of the fabric creating a smooth inner lumen. In the Endurant IIS bifurcated configuration, the 4 distal stents are sewn to the inside of the ipsilateral leg graft fabric. For all sizes, the stents on the contralateral leg are sewn to the inside of the graft fabric (Figure 1).

1.1.2. Limb Configuration

The proximal end of the limb configuration deploys within the legs of the bifurcated configuration, while the distal end deploys into the iliac artery. The proximal end of the limb configuration has an open web configuration (Figure 1), which contains no graft material in its stent valleys.

Note: A limb device is implanted on both the ipsilateral and contralateral legs of an Endurant IIS bifurcated configuration. See Deploy Limb Stent Graft into Ipsilateral Leg (Endurant IIS bifurcated configuration only) (Section 10.2.12).

1.1.3. Iliac Extension Configuration

An iliac extension configuration is available if additional distal stent graft length is needed. It has an open web configuration on its proximal end (Figure 1).

Note: An appropriately sized limb configuration can be used as an iliac extension configuration.

1.1.4. Aortic Extension and Abdominal Tube Configurations

Aortic extension and abdominal tube configurations are available if additional proximal stent graft length is needed. The aortic extension and abdominal tube stent grafts have a bare proximal suprarenal stent with anchor pins (Figure 1).

1.1.5. Aorto-Uni-Iliac (AUI) Configuration

The proximal end of the aorto-uni-iliac (AUI) configuration is deployed into the proximal neck and upper section of the aneurysm. All stents on the proximal aortic end of the AUI configuration are sewn to the outside of the graft fabric. The proximal stent (suprarenal) of the aortic section is not covered with fabric. As such, this bare stent design allows the AUI stent graft to be fixed above the renal arteries without obstructing them with graft fabric. Refer to Figure 1 for a pictorial representation of the proximal configuration. The suprarenal stent includes anchor pins to help fix the AUI device in place. The suprarenal stent is joined to the proximal edge of the graft by ultra high molecular weight polyethylene suture.

Distally, the aortic section tapers down to a smaller diameter tube. In the distal end of the tapered AUI device, the stents are sewn to the inside of the graft fabric.

Note: The Talent™ occluder system (available separately) can be used to assist in stent graft implantation.

1.2. Delivery System

The Endurant™ II delivery system that delivers all stent graft configurations consists of a single-use, disposable catheter with an integrated handle to provide controlled deployment. It is available in 14, 16, 18, and 20 French graft cover diameters and a working length of 57 cm ± 2 cm. The catheter assembly is flexible and compatible with a 0.035 in (0.88 mm) guidewire. There are 2 types of delivery systems: the aortic (Figure 2) and the iliac (Figure 3) delivery systems. The aortic delivery system delivers the bifurcated components, aortic extension, AUI, and abdominal tube stent graft configurations. The iliac delivery system delivers the limb and iliac extension stent graft configurations. The aortic delivery system features a tip capture mechanism, which is not present in the iliac delivery system.


Figure 2. Aortic Delivery System

- | | |
|--------------------|-----------------|
| 1. Rear Handle | 7. Graft Cover |
| 2. Back-End Wheel | 8. Markerband |
| 3. Screw Gear | 9. Spindle |
| 4. External Slider | 10. Sleeve |
| 5. Trigger | 11. Tapered Tip |
| 6. Front Grip | |


Figure 3. Iliac Delivery System

- | | |
|--------------------|----------------|
| 1. Rear Handle | 5. Front Grip |
| 2. Screw Gear | 6. Graft Cover |
| 3. External Slider | 7. Markerband |
| 4. Trigger | 8. Tapered Tip |

2. Indications for Use

The Endurant II/Endurant IIS stent graft system is indicated for the endovascular treatment of infrarenal abdominal aortic or aortoiliac aneurysms in patients with the following characteristics:

- adequate iliac or femoral access vessel morphology that is compatible with vascular access techniques, devices or accessories
- proximal neck length of ≥10 mm with insignificant calcification, or insignificant thrombus with ≤80° infrarenal and ≤45° suprarenal angulation and a vessel diameter approximately 10% to 20% smaller than the labeled Endurant II/Endurant IIS stent graft diameter.
- proximal neck length of ≥15 mm with insignificant calcification, or insignificant thrombus with ≤75° infrarenal and ≤60° suprarenal angulation and a vessel diameter approximately 10% to 20% smaller than the labeled Endurant II/Endurant IIS stent graft diameter.
- distal fixation length of ≥15 mm
- aortic neck diameters with a range of 19 to 32 mm
- iliac diameters with a range of 8 to 25 mm
- morphology suitable for aneurysm repair
- one of the following:
 - aneurysm diameter >5 cm

- aneurysm diameter of 4 to 5 cm, which has also increased in size by 0.5 cm in the last 6 months
- aneurysm that is at least 1.5 times the diameter of the normal infrarenal aorta

3. Contraindications

The Endurant II/Endurant IIS stent graft system is contraindicated in:

- patients who have a condition that threatens to infect the graft
- patients who are sensitive to or have allergies to the device materials listed in Table 1

Also consider patient selection information (Section 4.2).

4. Warnings and Precautions

Caution: Read all instructions carefully. Failure to properly follow the instructions, warnings, and precautions may lead to serious consequences or injury to the patient.

4.1. General

- The Endurant II/Endurant IIS stent graft system should only be used by physicians and teams trained in vascular interventional techniques, including training in the use of the device. Specific training expectations are described in Physician Training Requirements (Section 8.1).
- Always have a vascular surgery team available during implantation or reintervention procedures in the event that conversion to open surgical repair is necessary.

4.2. Patient Selection

- Inappropriate patient selection may result in poor device performance or device performance not otherwise in accordance with the specifications.
- Do not use the Endurant II/Endurant IIS stent graft system in patients unable to undergo, or who will not be compliant with, the necessary preoperative and postoperative imaging and implantation procedures (Section 9 to Section 12).
- The Endurant II/Endurant IIS stent graft system is not recommended in patients who cannot tolerate contrast agents necessary for intraoperative and postoperative follow-up imaging.
- The Endurant II/Endurant IIS stent graft system is not recommended in patients exceeding weight or size limits necessary to meet imaging requirements.
- Key anatomic elements that may affect successful exclusion of the aneurysm include severe proximal neck angulation with short proximal aortic neck ($>75^\circ$ with neck length <15 mm or $>60^\circ$ with neck length <10 mm); and thrombus or calcium formation at the arterial implantation sites, specifically the proximal aortic neck and distal iliac artery interface; narrowing of the aorta at the bifurcation point. Irregular calcification or plaque may compromise the fixation and sealing of the implantation sites. Necks exhibiting these key anatomic elements may be more conducive to graft migration. The use of a bifurcated device in a patient with a narrowing of the distal aorta may result in reduced flow through the limbs. The AUI stent graft is recommended for patients with significant distal aorta narrowing.
- Iliac conduits may be used to ensure the safe insertion of the delivery system if the patient's access vessels, as determined by treating physician, preclude safe insertion of the delivery system.
- The long-term safety and effectiveness of the Endurant II/Endurant IIS stent graft system has not been established.
- The safety and effectiveness of the Endurant II/Endurant IIS stent graft system has not been evaluated in patients who:
 - are less than 18 years old
 - are pregnant or lactating
 - have an aneurysm that is:
 - suprarenal
 - juxtarenal or pararenal
 - isolated iliofemoral
 - mycotic
 - inflammatory
 - pseudoaneurysm
 - have a dominant patent inferior mesenteric artery and an occluded or stenotic celiac or superior mesenteric artery
 - have an untreated thoracic aneurysm >4.5 cm in diameter
 - requires emergent aneurysm treatment, eg, trauma or rupture
 - have a history of bleeding diathesis or coagulopathy
 - have had a myocardial infarction (MI) or cerebrovascular accident (CVA) within 3 months prior to implantation
 - have a reversed conical neck, which is defined as a >4 mm distal increase over a 10 mm length
 - have a known hypersensitivity or contraindication to anticoagulants, antiplatelets, or contrast media, which is not amenable to pre-treatment
 - have significant (typically $>25\%$ of vessel circumference of aortic neck and iliac artery, or $>50\%$ of the length of the iliac artery) aortic mural thrombus at either the proximal or distal attachment location that would compromise bilateral fixation and seal of the device
 - have ectatic iliac arteries requiring bilateral exclusion of hypogastric blood flow
 - have arterial access site that is not expected to accommodate the diameter of the device (14 to 20 Fr) due to size or tortuosity
 - have active infection at the time of the index procedure documented by pain, fever, drainage, positive culture, or leukocytosis (WBC $>11,000/\text{mm}^3$) that is treated with antimicrobial agents (nonprophylactic)
 - have congenital degenerative collagen disease
 - have a creatinine >2.0 mg/dl (or >182 $\mu\text{mol/L}$)
 - are on dialysis
 - have connective tissue disorder
- All patients should be advised that endovascular treatment requires lifelong, regular follow-up to assess their health and the performance of their endovascular graft. Patients with specific clinical findings (eg, endoleaks, enlarging aneurysms, or changes in the structure or position of the endovascular graft) should receive enhanced follow-up. Specific follow-up guidelines are described in Section 12.
- Patients experiencing reduced blood flow through the graft limb or leaks may be required to undergo secondary interventions or surgical procedures.
- Intervention or conversion to standard open surgical repair following initial endovascular repair should be considered for patients experiencing enlarging aneurysms or endoleak. An increase in aneurysm size or persistent endoleak may lead to aneurysm rupture.

4.3. Before the Implant Procedure

- Preoperative planning for access and placement should be performed before opening the device packaging.
- Carefully inspect the Endurant II/Endurant IIS stent graft system packaging and device for damage or defects prior to use. Do not use product if any sign of damage or breach of the sterile barrier is observed. Do not attempt to resterilize the delivery system or stent graft.
- Do not bend, kink, or otherwise alter the Endurant II/Endurant IIS delivery system prior to implantation because it may cause deployment difficulties.

- To reduce the risk of thrombotic complications, an additional bolus of IV heparin should be administered before inserting the device.

4.4. During the Implant Procedure

- Exercise care in handling and delivery technique to help prevent vessel rupture.
- Studies indicate that the danger of micro-embolization increases with increased procedure duration.
- Renal complications may occur:
 - from an excess use of contrast agents
 - as a result of embolic or misplaced stent graft
- Do not deploy the stent grafts in a location that will cause an endoleak or occlude arteries necessary to supply blood flow to organs or extremities. This could necessitate surgical removal of the device.
- Use fluoroscopic guidance to advance the delivery system, detect kinking, or assess alignment problems with the stent graft devices. Do not use excessive force to advance or withdraw the delivery system when resistance is encountered. If the delivery system kinks during insertion, do not attempt to deploy the stent graft. Remove the device and insert a new system.
- Do not continue to torque the delivery system if the tip is not rotating with the delivery system.
- Exercise particular care in difficult areas, such as areas of stenosis, intravascular thrombosis, or in calcified or tortuous vessels. Balloon angioplasty at the site of a narrowed or stenotic vessel may be considered prior to attempting to gently reintroduce the catheter delivery system.
- An inadequate seal zone may result in increased risk of leakage into the aneurysm or migration of the stent graft.
- Systemic anticoagulation should be used during the implantation procedure based on hospital/physician protocol. If heparin is contraindicated, an alternative anticoagulant should be considered.
- Stent grafts cannot be replaced or drawn back into the delivery system, even if it is only partially deployed.
- If the graft cover is accidentally withdrawn, the device will prematurely deploy and may be incorrectly positioned.
- For the limb stent graft overlap criteria with Endurant II ipsilateral leg only, please refer to Table 8. As noted in Table 8, for the limb stent graft configurations that have an overlap criteria of 3 stents only, do not overlap more than 3 stents.
- When deploying the stent graft, be sure to hold the front grip stationary.
- If a compliant balloon catheter is used, do not over-inflate or inflate outside the graft material. Follow all manufacturer instructions regarding catheter operation.
- High pressure injections of contrast media made at the edges of the stent graft immediately after implantation can cause endoleaks.

4.5. Treatment and Follow-up

- Any endoleak left untreated during the implantation procedure must be carefully monitored after implantation.
- All patients with endovascular aneurysm repair should undergo periodic imaging to evaluate the stent graft, aneurysm size, and occlusion of vessels in the treatment area. Significant aneurysm enlargement (>5 mm), the appearance of a new endoleak, evidence of perigraft flow, change in aneurysm pulsatility, or migration resulting in an inadequate seal zone should prompt further investigation and may indicate the need for additional intervention or surgical conversion.
- Additional treatment including endovascular treatment or surgical conversion should be strongly considered in the following cases:
 - aneurysm growth >5 mm, with or without endoleak, since last follow-up
 - change in aneurysm pulsatility, with or without growth or endoleak
 - persistent endoleak, with or without aneurysm growth
 - stent graft migration resulting in an inadequate seal zone
 - decrease in renal function due to renal artery occlusion (migration or poor placement)
- Following endovascular aneurysm repair (EVAR), spinal cord ischemia (SCI) may result in a rare complication of paraplegia or paraparesis. Cerebrospinal fluid (CSF) drain is advised if spinal cord ischemia is suspected.

4.6. MRI Safety Information

Nonclinical testing has demonstrated that the Endurant II/Endurant IIS stent graft system is MR Conditional. It can be scanned safely in both 1.5 T and 3.0 T MR systems only, with using only the specific testing parameters (Section 9.6). Additional MRI safety information is found in Section 9.6.

5. Adverse Events

5.1. Potential Adverse Events

Adverse events that may occur or require intervention include, but are not limited to:

- amputation
- anesthetic complications and subsequent attendant problems (eg, aspiration)
- aneurysm enlargement
- aneurysm rupture and death
- aortic damage, including perforation, dissection, bleeding, rupture and death
- arterial or venous thrombosis or pseudoaneurysm
- arteriovenous fistula
- bleeding, hematoma, or coagulopathy
- bowel complications (eg, ileus, transient ischemia, infarction, necrosis)
- cardiac complications and subsequent attendant problems (eg, arrhythmia, myocardial infarction, congestive heart failure, hypotension, hypertension)
- claudication (eg, buttock, lower limb)
- death
- edema
- embolization (micro and macro) with transient or permanent ischemia or infarction
- endoleak
- fever and localized inflammation
- genitourinary complications and subsequent attendant problems (eg, ischemia, erosion, fistula, incontinence, hematuria, infection)
- hepatic failure
- impotence

- infection of the aneurysm, device access site, including abscess formation, transient fever, and pain
- lymphatic complications and subsequent attendant problems (eg, lymph fistula)
- neurologic local or systemic complications and subsequent attendant problems (eg, confusion, stroke, transient ischemic attack, paraplegia, paraparesis, paralysis)
- occlusion of device or native vessel
- pulmonary complications and subsequent attendant problems
- renal complications and subsequent attendant problems (eg, artery occlusion, contrast toxicity, insufficiency, failure)
- stent graft: improper placement; incomplete deployment; migration; suture break; occlusion; infection; stent fracture; graft twisting or kinking; insertion and removal difficulties; graft material wear; dilatation; erosion; puncture and perigraft flow
- surgical conversion to open repair
- vascular access site complications, including infection, pain, hematoma, pseudoaneurysm, arteriovenous fistula, dissection
- vascular spasm or vascular trauma (eg, iliofemoral vessel dissection, bleeding, rupture, death)
- vessel damage
- wound complications and subsequent attendant problems (eg, dehiscence, infection, hematoma, seroma, cellulitis)

6. Patient Selection and Treatment

6.1. Individualization of Treatment

Each Endurant II/Endurant IIS stent graft system must be ordered in the appropriate size to fit the patient's anatomy. Proper sizing of the device is the responsibility of the physician. The stent graft should be oversized to be larger than the vessel inner diameter (aortic configurations are recommended to be oversized 10 to 20%; iliac configurations are recommended to be oversized 10 to 25%). Refer to Section 9.2 for further details. The stent graft configurations cover aortic diameters ranging from 18 to 32 mm and iliac diameters from 8 to 25 mm. The recommended overall length of the stent graft including multiple deployed devices should extend from the lowest renal artery to just above the internal iliac or hypogastric artery. All lengths and diameters of the stent graft devices necessary to complete the procedure should be available to the physician, especially when pre-operative case planning measurements (treatment diameters/lengths) are not certain. Use of this approach allows for greater intraoperative flexibility to achieve optimal procedural outcomes.

Medtronic may consult with physicians to determine proper stent graft dimensions based on the physician's assessment of the patient's anatomical measurements. The benefits and risks previously described should be carefully considered for each patient before use of the stent graft system.

Note: Due to the nature of the design and the flexibility of the Endurant II/Endurant IIS stent graft system, the overall length of each stent graft may be shorter when deployed.

Caution: Vessel over-distension and damage, or partial stent graft infolding, may be caused by excessive oversizing of the stent graft in relation to the diameter of the blood vessel.

7. Patient Counseling Information

The physician should review the following risks and benefits when counseling the patient about this endovascular device and procedure:

- patient age and life expectancy
- risks and benefits related to open surgical repair
- risks and benefits related to endovascular repair
- risks related to noninterventional treatment or medical management
- risks of aneurysm rupture compared to endovascular repair
- possibility that subsequent endovascular or open surgical repair of the aneurysm may be required
- the long-term safety and effectiveness of the Endurant II/Endurant IIS stent graft system has not been established
- long-term, regular follow-up is needed to assess patient health status and stent graft performance
- patients with specific clinical findings (eg, endoleaks, enlarging aneurysms) should be monitored closely
- symptoms of aneurysm rupture

Medtronic recommends that the physician disclose to the patient, in written form, all risks associated with treatment using the Endurant II/Endurant IIS stent graft system. Details regarding risks occurring during and after implantation of the device are provided in Adverse Events (Section 5).

8. How Supplied

8.1. Sterility

Each stent graft configuration (bifurcated, AUI, limb, aortic and iliac extension, and abdominal tube) is individually contained within a delivery system. It is sterilized using electron beam and is supplied sterile for single use only.

- Do not reuse or attempt to resterilize.
- If the device is damaged or the integrity of the sterile barrier has been compromised, do not use the product and contact your Medtronic representative for return information.

8.2. Contents

- One Endurant II/Endurant IIS Stent Graft System
- One Instructions for Use

8.3. Storage

Store the system at room temperature in a dark, dry place.

9. Clinical Use Information

9.1. Physician Training Requirements

All physicians should be trained in the use of the Endurant II/Endurant IIS stent graft system before using it.

Caution: The Endurant II/Endurant IIS stent graft system should only be used by physicians and teams trained in vascular interventional techniques and in the use of this device.

The following are the knowledge and skill requirements for physicians using the Endurant II/Endurant IIS stent graft system:

- natural history of abdominal aortic aneurysms (AAA), aortoiliac aneurysms, and comorbidities associated with AAA repair
- radiographic, fluoroscopic, and angiographic image interpretation
- appropriate use of radiographic contrast material
- arterial cutdown, arteriotomy, and repair or percutaneous access and closure techniques

- nonselective and selective guidewire and catheter techniques
- embolization
- angioplasty
- endovascular stent placement
- snare techniques
- techniques to minimize radiation exposure
- device selection and sizing

9.2. Recommended Device Sizing

The Endurant II/Endurant IIs stent graft systems are available in the sizes described in Table 2 through Table 7. For questions about this device sizing, refer to contact information in the back of this Instructions for Use.

Table 2. Sizing Chart — Endurant II Bifurcated Configuration

OD (Fr)	Proximal x Distal Diameter (mm x mm)	Covered Length (mm)	Vessel Inner Diameter (mm)
20	36x20	145, 186	29-32
	36x16		
	32x20		
	32x16		
18	28x20	124, 145, 186	26-28
	28x16		
	28x13		
	25x16		
	25x13		23-25
	23x16		
	23x13		
	23x13		
16	25x16	103	21-22
	25x13		
	23x16		
	23x13		

Table 3. Sizing Chart — Endurant IIs Bifurcated Configuration

OD (Fr)	Proximal x Distal Diameter (mm x mm)	Covered Length (mm)	Vessel Inner Diameter (mm)
20	36x14	103	29-32
	32x14		
	28x14		
	25x14		
18	25x14	103	23-25
	25x13		
	23x14		
	23x13		

Table 4. Sizing Chart — Limb Configuration

OD (Fr)	Proximal x Distal Diameter (mm x mm)	Covered Length (mm)	Vessel Inner Diameter (mm)
16	16x28	82, 93, 124, 156, 199	23-25
	16x24		
	16x20		
	16x16		
	16x13	156, 199	19-22
	16x10		
	16x10		
	16x10		
14	16x16	82, 93, 124	15-18
	16x13		
	16x10		
	16x10		

Table 5. Sizing Chart — Iliac Extension Configuration

OD (Fr)	Proximal x Distal Diameter (mm x mm)	Covered Length (mm)	Vessel Inner Diameter (mm)
18	28x28	82	23-25
16	24x24		19-22
14	20x20		15-18
	13x13		10-11
	10x10		8-9

Table 6. Sizing Chart — Aortic Extension and Abdominal Tube Configuration

OD (Fr)	Proximal x Distal Diameter (mm x mm)	Covered Length (mm)	Vessel Inner Diameter (mm)
20	36x36	49, 70	29-32
	32x32		26-28
	28x28		23-25
	25x25		21-22
18	23x23		19-20
	23x23		19-20

Table 7. Sizing Chart — Aorto-Uni-Iliac (AUI) Configuration

OD (Fr)	Proximal x Distal Diameter (mm x mm)	Covered Length (mm)	Vessel Inner Diameter (mm)
20	36x14	102	29-32
	32x14		26-28
	28x14		23-25
	25x14		21-22
18	23x14		19-20
	23x14		19-20

9.3. Device Inspection

Inspect the device and packaging to verify that no damage or defects exist. If the 'Use by' date has elapsed, the device is damaged, or the sterile barrier has been compromised, do not use the device and contact a Medtronic representative for return or replacement information.

9.4. Additional Required Equipment

- Additional Endurant II or Endurant IIs stent graft systems of various lengths and diameters
- Fluoroscope with digital angiographic capabilities (C-arm or fixed unit). Fluoroscopic imaging and the ability to record and recall all imaging.
- Assorted guidewires of adequate length
- Heparinized saline solution

9.5. Additional Recommended Equipment

- Introducer sheaths

14 Instructions for Use English

- Power injector
- Radiopaque ruler with centimeter increments
- Assorted balloon catheters
- Compliant balloon catheters
- Radiopaque contrast media
- Sterile silicone lubricant or sterile mineral oil
- Interventional snare devices
- Endovascular coils and vascular plugs

9.6. MRI Information

Nonclinical testing has demonstrated that the Endurant II/Endurant IIS stent graft is MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 or 3.0 tesla only
- Maximum spatial gradient magnetic field of 2500 gauss/cm or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 4 W/kg (First Level Controlled Mode)

Under the scan conditions defined above, the Endurant II/Endurant IIS stent graft is expected to produce:

- a maximum temperature rise of 1.00°C after 15 minutes of continuous scanning in a 1.5 tesla scanner
- and a maximum temperature rise of 3.27°C after 15 minutes of a continuous scanning in a 3.0 tesla scanner

The image artifact extends approximately 5 mm and 8 mm from the device, both inside and outside the device lumen when scanned in nonclinical testing using the sequence: spin echo and gradient echo, respectively, in a 3.0 tesla Siemens TrioTim (VB 13 Software) MR system with a whole body coil.

10. Implant Instructions

10.1. Vascular Access and Device Preparation

Correct sizing of the aorta and iliac vessels must be determined before implantation of the aortic and iliac stent graft configurations using contrast-enhanced computer-aided tomography (CT), as well as angiograms of both the iliac arteries and aorta. 3D imaging may also be beneficial. Refer to Recommended Device Sizing (Section 9.2). These images should be available for review during the procedure. Vascular instruments and other surgical supplies necessary to gain access to the artery should also be available.

To reduce the risk of thromboembolism, it is recommended that the patient be heparinized for the duration of the procedure.

Caution: Do not retract the graft cover of the delivery system until it is accurately placed within the vasculature and ready for deployment.

Caution: Never advance or retract equipment from the vasculature without the use of fluoroscopy.

10.1.1. Vascular Access

1. Following aseptic procedure, perform a vascular access at the femoral arteries.
2. Place a guidewire in the ipsilateral femoral artery and advance it above the renal arteries.
3. From the contralateral side femoral artery, place a second guidewire directed to the abdominal aorta.
4. Over this guidewire, place an angiography catheter above the renal arteries.
5. Take an angiogram.

Note: An additional incision might be necessary to access the common iliac artery.

10.1.2. Device Preparation

1. View the delivery system under fluoroscopy to visualize the radiopaque markers on the stent graft. The radiopaque markers indicate the position of the proximal and distal edges of the graft material.
2. Turn the graft cover to align the radiopaque gate marker on the short leg of the bifurcated configuration with the patient's contralateral iliac artery.
3. Flush the guidewire lumen with heparinized saline solution.
4. Prior to insertion into the vessel, activate the hydrophilic coating by wiping the outer surface of the graft cover with a sterile gauze, saturated in saline, until the graft cover is slippery to touch.

10.2. Delivery Procedure

Medtronic recommends using an appropriate caliber introducer sheath to perform diagnostic tests. No sheath is necessary for the introduction of the delivery system or deployment of the stent graft.

Caution: Do not remove the guidewire while the delivery system is in the patient.

Warning: To prevent thrombotic problems, a second bolus of IV heparin is recommended before inserting the device.

10.2.1. Introduction of Bifurcated Configuration

Warning: Do not advance the delivery system without placing a guidewire.

1. Slowly insert the delivery system.
2. Advance over the guidewire so that the proximal most stents and the radiopaque markers are visualized in the target proximal aortic neck (Figure 4).
3. Inject contrast media through an angiographic (pigtail) catheter into the abdominal aorta and mark the position of the target location, either on the imaging screen or on the patient's body.
4. Adjust the position of the bifurcated stent graft configuration so that the top edge of the graft fabric is just below the lowest renal artery. (Note: The edge of the graft fabric is 0.5 to 1.0 mm above the top edge proximal radiopaque markers.)

Note: If the top edge of the graft fabric is to be placed very close to the renal arteries, contrast media may be injected to identify the location of the lower renal artery and verify the position before full deployment.

Caution: Once proximal position has been identified, do not move the patient or imaging equipment, as it may compromise accuracy of stent graft placement.

Caution: The angiographic catheter can be removed prior to deployment. However, if the angiographic catheter is not removed until after deployment, ensure that the tip is straightened (such as with a pigtail catheter) with a guidewire before removal so that the stent graft is not pulled down.

Caution: When aligning the position of stent graft, be sure the fluoroscope is angled perpendicular to the center line of the infrarenal aorta to avoid parallax or other sources of visualization error. Some cranial caudal angulation of the image intensifier (I-I) tube may be necessary, especially if there is anterior angulation of the aneurysm neck.


Figure 4. Introduction of Aortic Delivery System

10.2.2. Confirm Position

1. Ensure that the distal portion of the contralateral stub leg is above the aortic bifurcation and within the aneurysmal sac, and not within the iliac vessel.
2. Rotate the handle until the radiopaque marker on the distal stent of the contralateral stub leg is aligned with the contralateral iliac artery.

Note: When attempting to rotate the system, if the tip does not rotate with the handle, pull back the system and reposition until the intended position is achieved.

10.2.3. Deploy Proximal End of Bifurcated Configuration

1. With 1 hand on the front grip, hold the delivery system stationary.
2. With the other hand, slowly withdraw the graft cover by rotating the external slider counterclockwise (in the direction of the slider arrow), until the constrained suprarenal stent is exposed and 2 to 3 of the covered stents have been fully deployed (Figure 5).
3. Use angiography to verify position of the bifurcated configuration in relation to the renal arteries.
4. If needed, gently push the entire delivery system proximally or pull distally until the proximal end of the graft material is even with the distal edge of the lowest renal artery.

Note: In the unlikely event of delivery system failure that results in partial stent graft deployment due to graft cover severance, the "handle disassembly" technique may permit successful deployment of the stent graft. Refer to Section 11, Bail-Out Techniques.

Caution: Do not rotate the graft cover during deployment as this may torque the device and cause it to rotate during deployment.

Caution: If the graft cover is accidentally withdrawn, the stent graft will prematurely deploy and may be incorrectly positioned.

Warning: Failure to properly align the radiopaque markers may result in improper deployment of the stent graft.


Figure 5. Deploying Proximal End of Bifurcated Configuration

10.2.4. Deploy Contralateral Leg of Bifurcated Configuration

Continue holding the front grip of the delivery system stationary and then rotate the slider handle counterclockwise, stopping immediately after the contralateral leg is released from the graft cover or delivery sheath (Figure 6).


Figure 6. Deploy the Contralateral Leg

10.2.5. Release Proximal End of Suprarenal Stent

1. Use angiography to verify the position of the bifurcated configuration in relation to the renal arteries.
2. Continue to hold the delivery system stationary with 1 hand on the front grip.
3. With the other hand, rotate the back-end wheel clockwise, moving the tapered tip forward to release the proximal end of the suprarenal stent (Figure 7).
4. Observe the release of the suprarenal stent under fluoroscopy and continue turning the back-end wheel until it is completely clear of the delivery system spindle.

Note: In the unlikely event that the proximal end of the suprarenal stent cannot be released, refer to Section 11, Bail-Out Techniques.

Caution: In the unlikely event that the back-end wheel separates during wheel rotation, remove the wheel. Manually advance the exposed tabs forward on the screw gear until all of the suprarenal stents release from the spindle. Refer to Section 11, Bail-Out Techniques.


Figure 7. Release Proximal End of Suprarenal Stent

10.2.6. Deploy Distal End of Bifurcated Configuration

Either continue to rotate the external slider counterclockwise or while holding the front grip of the delivery system stationary, use thumb to pull the trigger on the external slider and pull it back until the bifurcated stent graft is completely deployed.

Note: Retract the graft cover past the flexible stent stop tip (approximately 10 mm) to ensure that the graft cover edge does not disturb the graft position during forward advancement of the catheter for tip recapture.

Caution: When using the trigger to rapidly deploy the stent graft, the delivery system must remain stationary. Do not rotate the delivery system during stent graft deployment.


Figure 8. Deploy Distal End of the Bifurcated Configuration

1. Endurant II bifurcated configuration
2. Endurant IIs bifurcated configuration

10.2.7. Recapture Spindle in Tapered Tip (Endurant II Bifurcated Configuration Only)

Note: For the Endurant IIs bifurcated configuration, leave the delivery system in situ while deploying the limb stent graft into the contralateral leg. Skip to Deploy Limb Stent Graft into Contralateral Leg (Section 10.2.9).

1. Continue to hold the delivery system stationary with 1 hand on the front grip.
2. Confirm the spindle has fully separated from the suprarenal stent; gently torque the delivery system if it has not fully separated.
3. Gently rotate the delivery system while pushing the entire delivery system approximately 3 cm proximally so that the tapered tip and spindle are completely clear of the suprarenal stent.
4. With the other hand, rotate the back-end wheel counterclockwise to recapture the spindle in the tapered tip (Figure 9).
5. Observe the recapture of the spindle within the sleeve of the tapered tip under fluoroscopy.
6. Continue turning the back-end wheel counterclockwise until the spindle has been completely recaptured and the back-end wheel is at the bottom (Figure 9).

Note: When pushing the delivery system forward, be careful not to displace the distal end of the ipsilateral limb.

Note: Ensure that the suprarenal stent is fully disengaged from the spindle before pushing the delivery system forward.

Note: If the spindle catches on the suprarenal stent during advancement, completely advance the back-end wheel clockwise. Using a gentle in-and-out motion with the delivery system, rotate the delivery system until the spindle slips past the suprarenal stent. Then continue with the withdrawal process.

Caution: Stop rotating the back-end wheel when the bottom of the back-end screw gear is reached.

Warning: Failure to adequately advance the delivery system to recapture the spindle can result in the trapping of a suprarenal apex within the tapered tip sleeve. This will alter the proximal landing zone during delivery system withdrawal.


Figure 9. Recapture Spindle in Tapered Tip

10.2.8. Remove Delivery System (Endurant II Bifurcated Configuration Only)

1. Continue to hold the delivery system stationary with 1 hand on the front grip and the other hand on the external slider.
2. Gently torque and withdraw the delivery system until the spindle is retracted into the fabric portion of the stent graft.
3. Pull back the external slider trigger and hold it stationary while bringing the front grip to the slider (Figure 10).
4. Use continual fluoroscopy and watch the top of the bifurcated stent graft while slowly pulling back the tapered tip into the graft cover of the delivery system.
5. Gently remove the delivery system. Use fluoroscopy to ensure that the bifurcated configuration does not move during withdrawal.

Note: Maintain vessel access and wire placement until all stent graft configurations are in place.


Figure 10. Remove Delivery System

10.2.9. Deploy Limb Stent Graft Into Contralateral Leg

1. Prepare the iliac stent graft system as described in Device Preparation (Section 10.1.2).
2. On the patient's contralateral side, insert a guidewire through the contralateral leg and aortic neck of the previously placed bifurcated configuration.
3. Place the delivery system over the guidewire and into the contralateral leg of the bifurcated stent graft.
4. Insert the limb stent graft into the contralateral leg of the bifurcated stent graft. The proximal radiopaque marker of the limb should be aligned to the radiopaque marker at the bifurcation of the bifurcated stent graft configuration (Figure 11).


Figure 11. Introduction of the Iliac Delivery System

1. Endurant II bifurcated configuration
2. Endurant IIs bifurcated configuration
5. Ensure there is a 3-stent overlap (Figure 12).
6. With 1 hand on the front grip, hold the delivery system stationary.
7. With the other hand, slowly withdraw the graft cover by rotating the external slider counterclockwise.
8. At any point, pull the slider trigger and pull the external slider all the way back to finish deploying the limb configuration.
9. Remove the delivery system.

Note: In the unlikely event of delivery system failure that results in partial stent graft deployment, the "handle disassembly" technique may permit the successful deployment of the stent graft configuration. Refer to Section 11, Bail-Out Techniques.

Caution: Do not torque the iliac delivery system while in the patient.


Figure 12. Deploy Limb Configuration

1. Endurant II bifurcated configuration
2. Endurant IIs bifurcated configuration

10.2.10. Recapture Spindle in Tapered Tip for Delivery System in the Ipsilateral Leg (Endurant IIs Bifurcated Configuration Only)

1. Hold the delivery system stationary with 1 hand on the front grip.
2. Confirm the spindle has fully separated from the suprarenal stent; gently torque the delivery system if it has not fully separated.
3. Gently rotate the delivery system while pushing the entire delivery system approximately 3 cm proximally so that the tapered tip and spindle are completely clear of the suprarenal stent.
4. With the other hand, rotate the back-end wheel counterclockwise to recapture the spindle in the tapered tip (Figure 13).
5. Observe the recapture of the spindle within the sleeve of the tapered tip under fluoroscopy.
6. Continue turning the back-end wheel counterclockwise until the spindle has been completely recaptured and the back-end wheel is at the bottom (Figure 13).

Note: When pushing the delivery system forward, be careful not to displace the distal end of the ipsilateral limb.

Note: Ensure that the suprarenal stent is fully disengaged from the spindle before pushing the delivery system forward.

Note: If the spindle catches on the suprarenal stent during advancement, completely advance the back-end wheel clockwise. Using a gentle in-and-out motion with the delivery system, rotate the delivery system until the spindle slips past the suprarenal stent. Then continue with the withdrawal process.

Caution: Stop rotating the back-end wheel when the bottom of the back-end screw gear is reached.

Warning: Failure to adequately advance the delivery system to recapture the spindle can result in the trapping of a suprarenal apex within the tapered tip sleeve. This will alter the proximal landing zone during delivery system withdrawal.


Figure 13. Recapture Spindle in Tapered Tip

10.2.11. Remove Delivery System (Endurant IIs Bifurcated Configuration Only)

1. Continue to hold the delivery system stationary with 1 hand on the front grip and the other hand on the external slider.
2. Gently torque and withdraw the delivery system until the spindle is retracted into the fabric portion of the stent graft.
3. Pull back the external slider trigger and hold it stationary while bringing the front grip to the slider (Figure 14).
4. Use continual fluoroscopy and watch the top of the bifurcated stent graft while slowly pulling back the tapered tip into the graft cover of the delivery system.
5. Gently remove the delivery system. Use fluoroscopy to ensure that the bifurcated configuration does not move during withdrawal.

Note: Maintain vessel access and wire placement until all stent graft configurations are in place.


Figure 14. Remove Delivery System

10.2.12. Deploy Limb Stent Graft Into Ipsilateral Leg (Endurant IIs Bifurcated Configuration Only)

1. Prepare the iliac stent graft system as described in Device Preparation (Section 10.1.2).
2. On the patient's ipsilateral side, track the delivery system over the guidewire and into the ipsilateral leg of the previously placed bifurcated configuration. The limb stent graft has 2 markers on its proximal edge, 2 markers on its distal edge, and 1 overlap marker approximately 25 mm distal of the proximal markers.
3. Position the device. Overlap criteria between the limb stent graft and the ipsilateral leg of the bifurcated stent graft is dependent on limb selection. See Table 8 for recommended device overlap.

Table 8. Recommended Device Overlap — Limb Stent Graft and Ipsilateral Leg of Endurant IIs Bifurcated Stent Graft

Proximal diameter	Distal diameter	Length	Overlap	Deployment reference
16	10	82	3 stents only	Refer to step a
		93		
	13	82		
		93		
	20	82		
		93		
	24	82		
		93		
16	10	82	3 to 5 stents	Refer to step b
		93		
		124		
	13	124		
		156		
		199		
	16	82		
		93		
		124		
		156		
	20	124		
		156		
		199		
	24	124		
		156		
		199		
	28	124		
		156		
		199		

- a. 3 stent overlap. Align the overlap marker on the limb stent graft with the 2 markers on the distal end of the ipsilateral leg of the Endurant IIs bifurcated stent graft (Figure 15).

Warning: For the limb stent graft overlap criteria with Endurant IIs ipsilateral leg only, please refer to Table 8. As noted in Table 8, for the limb stent graft configurations that have an overlap criteria of 3 stents only, do not overlap more than 3 stents.


Figure 15. Introduction of the Iliac Delivery System — 3 Stent Overlap

- b. 3 to 5 stent overlap. For the limb stent grafts that can have an overlap between 3 to 5 stents, minimum 3 stent overlap can be achieved by following instruction in step a (Figure 15) or the maximum 5 stent overlap can be achieved by aligning the proximal edge markers on the limb stent graft to the bifurcation marker on the Endurant IIS bifurcated stent graft (Figure 16).


Figure 16. Introduction of the Iliac Delivery System — 5 Stent Overlap

4. With 1 hand on the front grip, hold the delivery system stationary.
 5. With the other hand, slowly withdraw the graft cover by rotating the external slide counterclockwise.
 6. At any point, pull the slider trigger and pull the external slider all the way back to finish deploying the limb configuration.
 7. Remove the delivery system (Section 10.2.11).
- Note:** In the unlikely event of delivery system failure that results in partial stent graft deployment, the "handle disassembly" technique may permit the successful deployment of the stent graft configuration. Refer to Section 11, Bail-Out Techniques.
- Caution:** Do not torque the iliac delivery system while in the patient.

10.2.13. Iliac or Aortic Extension Stent Graft Configurations

1. If an aortic extension stent graft configuration is needed, ensure that there is a minimum 3-stent overlap between the aortic extension stent graft and the bifurcated stent graft configuration.
2. Follow the bifurcated stent graft deployment process, except rotate the handle to open the extension component entirely before releasing the proximal end of the suprarenal stent of the aortic configuration.
3. If an iliac extension stent graft configuration is needed, ensure that there is a minimum 3-stent overlap between the iliac extension and the configuration it is being mated with. This is achieved by aligning the limb overlap marker with the most distal marker on the configuration with which the extension is being mated.
4. Follow the limb stent graft configuration deployment process (Section 10.2.8).

10.2.14. Abdominal Tube Stent Graft Configuration

If an abdominal tube stent graft configuration is needed, follow the bifurcated stent graft deployment process, with the following exception: Rotate the handle to open the abdominal tube stent graft entirely before releasing the proximal end of the suprarenal stent.

10.2.15. AUI Stent Graft Configuration

Follow the bifurcated stent graft configuration deployment steps, except rotate the handle to open the AUI stent graft covered portion completely before releasing the proximal end of the suprarenal stent of the AUI stent graft. Follow the steps described in Section 10.2.7 to Section 10.2.8 to remove the delivery system.

In the event that an additional distal stent graft configuration is required, use the limb stent graft as the AUI distal extension. Place the delivery system over the existing guidewire and follow the limb stent graft deployment steps in Section 10.2.8. To ensure proper docking between the 2 stent grafts, align the overlap radiopaque marker on the limb stent graft with the distal radiopaque markers on the AUI stent graft to ensure a 3 stent overlap.

An occluder may be used to block flow through the contralateral iliac artery. Refer to Section 10.2.16 and the Talent endoluminal occluder system Instructions for Use.

10.2.16. Occluder Deployment

The Talent occluder system (packaged separately) is available to be used with Endurant II/Endurant IIS stent graft system and is typically used in conjunction with the AUI stent graft component. The Talent occluder system is closed at both ends to stop retrograde blood flow into the aneurysm sac.

For details on use and implantation of the Talent occluder system, refer to the appropriate sections of the Talent endoluminal occluder system Instructions for Use.

10.2.17. Smoothing Stent Graft Fabric and Modeling Stent Graft

The Reliant™ stent graft balloon catheter (packaged separately) can be used to assist in stent graft implantation by modeling the covered portion of the stent graft and removing wrinkles and folds from the graft material, as needed. Use the balloon catheter to model the proximal and distal seal zones as well as any overlapping connection (or junction) areas between the stent graft components. Sub-optimal expansion of the self-expanding stent graft components may also be improved by use of the balloon catheter. Refer to the Reliant stent graft balloon catheter IFU for specific instructions.

Note: The Reliant stent graft balloon catheter is recommended for use with the Endurant II/Endurant IIS stent graft system. Data is not available for the use of other balloon catheters in remodeling stent grafts.

Caution: Over inflation of the balloon can cause graft tears or vessel dissection or rupture.

Warning: When expanding a vascular prosthesis, there is an increased risk of vessel injury or rupture, and possible patient death, if the balloon's proximal and distal radiopaque markers are not completely within the covered (graft fabric) portion of the prosthesis.

Warning: Do not use the Reliant stent graft balloon catheter in the treatment of dissections.

10.2.18. Seal Entry Sites

1. Remove ancillary devices prior to repairing the entry site.
2. Repair the entry site with standard closure technique.

10.2.19. Verify Placement and Seal

1. At the completion of the procedure, perform angiography to assess the stent graft for proximal and distal endoleaks and to verify position of the implanted stent graft in relation to the aneurysm and renal arteries.
2. Leaks at the attachment or connection sites should be treated using the balloon to remodel the stent graft against the vessel wall.
3. Major leaks that cannot be corrected by re-ballooning may be treated by adding aortic or iliac extension components to the previously placed stent graft components.

Caution: Any leak left untreated during the implantation procedure must be carefully monitored after implantation.

11. Bail-Out Techniques

In the unlikely event of delivery system failure, the following bail-out techniques may be used.

11.1. Screw Gear Handle Disassembly

If partial stent graft deployment due to graft cover severance occurs, the screw gear handle disassembly technique may permit successful deployment of the stent graft.

1. Pull back the trigger and fully retract the slider.
2. Stabilize the delivery system.
3. Insert the tips of a pair of hemostats into each of the screw gear handle disassembly ports on the front grip.
4. Disengage the front grip from the screw gear by pressing the tips of the hemostats into the handle disassembly ports simultaneously advancing the front grip away from the screw gear.
5. Advance the front grip until it fully clears the screw gear.
6. Separate the screw gear halves in order to identify the location of graft cover severance.
7. Manually retract the graft cover with your fingers or with hemostats until the stent graft is fully deployed.
8. Follow the instructions for tip capture deployment and delivery system removal.

11.2. Ballooning

If the captured proximal tip of the suprarenal stent cannot be deployed and the back-end wheel section still works, the ballooning technique may permit successful deployment of the suprarenal stent.

1. Use a compliant or semi-compliant balloon (Reliant balloon recommended).
2. Insert the balloon and move it to the bifurcated configuration's aortic section.
3. Inflate the balloon inside the stent graft to vessel size to stabilize stent graft.
4. Follow the instructions for tip capture deployment and delivery system removal.

11.3. Back-End Handle Disassembly

If no or partial deployment of the proximal end of the suprarenal stent occurs due to back-end wheel failure, the back-end handle disassembly technique may permit the successful deployment of the suprarenal stent.

1. Use hemostats to depress the exposed tabs to disassemble the back-end wheel.
2. Insert the tips of hemostats into each of the rear handle disassembly ports.
3. Disengage the rear handle by pressing the tips of the hemostats into the handle disassembly ports simultaneously retracting the rear handle from the delivery system.
4. Stabilize the delivery system.
5. Manually push the exposed tabs of the back-end T-tube to release the suprarenal stent from the spindle.
6. Manually pull back the exposed tabs of the back-end T-tube to recapture the tapered tip after deployment.
7. Follow the instructions for delivery system removal.
8. Hold the exposed tabs of the back-end T-tube so that it remains retracted and the tapered tip recaptured during delivery system removal.

11.4. Snare the Tapered Tip

If the back-end handle disassembly technique is unsuccessful due to an excessively high deployment force, the snare the tapered tip technique may permit successful deployment of the suprarenal stent.

1. Use a snare device.
2. Advance the snare device to the delivery system tapered tip section through upper torso access (i.e. brachial).
3. Utilize fluoroscopy to snare the edge of the delivery system tapered tip.
4. Stabilize the delivery system, especially the back-end section.
5. Pull the snare device to separate the suprarenal stent from the tip capture.
6. Manually pull back the back-end T-tube to recapture the tapered tip after deployment.
7. Follow the instructions for delivery system removal.
8. Ensure that the back-end T-tube remains retracted and the tapered tip recaptured during delivery system removal.

12. Follow-up Imaging Recommendations**12.1. General**

Current imaging of stent graft patients includes abdominal X-ray and CT, with and without contrast medium. Alternative imaging modalities such as magnetic resonance imaging should be used in patients with impaired renal function or intolerance to contrast media.

Imaging should be decided based upon the physician's clinical assessment of the patient pre- and post-implantation of the stent graft. After endovascular graft placement, patients should be regularly monitored for perigraft flow, aneurysm growth or changes in the structure or position of the endovascular graft. At a minimum, annual imaging is required, including 1) abdominal radiographs to examine device integrity (stent fracture, separation between bifurcated device and proximal cuffs or limb extensions, if applicable), and 2) contrast and non-contrast CT to examine aneurysm changes, perigraft flow, patency, tortuosity and progressive disease. If renal complications or other factors

preclude the use of image contrast media, abdominal radiographs and duplex ultrasound may provide similar information.

12.2. X-ray

Abdominal X-rays should be used to assess the presence of stent graft fracture. Four-view kidney, ureter, bladder (KUB) X-rays should be taken. Posterior/anterior (PA) and lateral images are recommended for visualization of the stent graft. Ensure the entire device is captured on images for device assessment.

12.3. CT with Contrast

Contrast-enhanced CT should be used to assess stent graft fixation, deformation, apposition to the vessel wall at proximal and distal fixation sites, stent graft migration, stent graft patency, AAA size, occlusion of branch vessels, and endoleak (including source and type if present).

A pre-contrast scan of 5 mm thick slices is suggested to determine if there are calcifications or areas where metal artifacts may be misinterpreted as endoleak. An arterial phase with <3 mm slice thickness and overlapping images with coverage from the celiac artery to the external iliac is recommended. In aneurysms that are not shrinking and have no apparent endoleak or fixation problems, a late venous phase scan may be performed. The venous phase scan may also be performed with thicker collimation (5 mm). It is recommended that the source data set be archived in case specialized evaluation is needed later (volume measurements, 3-dimensional reconstruction, or computer-aided measurement software). If the aneurysm is not shrinking by more than 5 mm within the first year, volume measurements may be obtained as a more sensitive indicator of AAA size using 3-dimensional software. Patients who are allergic to contrast should be pre-medicated 12-24 hours prior to receiving the drug.

12.4. Noncontrast CT

For patients with impaired renal function or those who are allergic to contrast medium, a spiral CT without contrast may be considered to assess stent graft fixation, deformation, apposition to the vessel wall at proximal and distal fixation sites, stent graft migration, occlusion of vessels, and size of the AAA diameter and volume measurements.

12.5. Duplex Ultrasound

For patients with impaired renal function or those who are allergic to contrast medium, a color-duplex ultrasound may be considered to assess size of AAA diameter, endoleaks, and stent graft occlusion and stenosis.

12.6. MRI or MRA

Patients with impaired renal function, ie, renal insufficiency, may also be considered for magnetic resonance imaging or angiography (MRI, MRA) in facilities that have expertise in this area. Artifact may occur related to the stent, and care should be used to insure adequate imaging of the outer aneurysm wall to assess AAA size. Volume measurement may be helpful if the aneurysm is not clearly shrinking. If there are concerns regarding imaging of calcified areas, fixation sites, or the outer wall of the aneurysm sac, adjunctive CT without contrast may be needed.

12.7. Imaging Tests

Refer to the table below for the recommended follow-up imaging schedule after stent graft implant.

Table 9. Imaging Recommendations

Imaging Test	Initial and Annual Follow-up for Up to Five Years
CT ^a or MRA ^b with contrast	X
Abdominal X-rays (4-View, KUB)	X
^a CT evaluation may include "3-phase technique," volume studies, 3-D reconstruction, or computer-aided measurements.	
^b MRA may be used in patients with impaired renal failure or contrast intolerance.	

12.8. Supplemental Imaging

Note: Additional radiological imaging may be necessary to further evaluate the stent graft in situ based on findings revealed by one of the surveillance programs. The following recommendations may be considered.

- If there is evidence of poor or irregular position of the stent graft, severe angulation, kinking or migration of the stent graft on abdominal X-rays, a spiral CT should be performed to assess aneurysm size and the presence or absence of an endoleak.
- If a new endoleak or increase in AAA size is observed by spiral CT, adjunctive studies such as 3-D reconstruction or angiographic assessment of the stent graft and native vasculature may be helpful in further evaluating any changes of the stent graft or aneurysm.
- Spiral CT without contrast, MRI or MRA may be considered in select patients who cannot tolerate contrast media or who have renal function impairment. For centers with appropriate expertise, gadolinium or CO₂ angiography may be considered in patients with renal function impairment requiring angiographic assessment.

13. Additional Surveillance and Treatment

Additional endovascular repair or open surgical aneurysm repair should be considered for patients with evidence of suboptimal stent graft fixation, proximal endoleak, distal endoleak, junction endoleak, unknown origin of persistent perigraft flow, or increase in AAA size > 5mm.

Disclaimer of Warranty

ALTHOUGH THE MEDTRONIC VASCULAR ENDURANT II/ENDURANT II's STENT GRAFT AND DELIVERY SYSTEM, HEREAFTER REFERRED TO AS THE "PRODUCT", HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, MEDTRONIC, INC., MEDTRONIC VASCULAR, INC. AND THEIR RESPECTIVE AFFILIATES, (COLLECTIVELY "MEDTRONIC") HAVE NO CONTROL OVER THE CONDITIONS UNDER WHICH THIS PRODUCT IS USED. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESSED AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

The exclusions and limitations set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

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© 2014 Medtronic, Inc.