13 "10", 2016
R. Broekhuizen

INSTRUCTIONS FOR USE
OF MEDICAL DEVICE

"Polymeric medical devices for collection, transportation, analysis and storage of biomaterials
with accessories"

***The undersigned authority legalizes the
signature of the signing proxy(ies).***

THE NETHERLANDS CHAMBER OF COMMERCE AND INDUSTRY

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R.D.M. van den Dool

Electronic Certifying Stamp

CONTINUOUS INNOVATION FOR PATHOLOGY

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APOSTILLE

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- 1 Country: THE NETHERLANDS
This public document
2. has been signed by **R.D.M. van den Dool**
- 3 acting in the capacity of official of the
Chamber of Commerce and Industry and
- 4 bears the seal/stamp of aforesaid Chamber of
Commerce and Industry

Certified

- 5 in 's-GRAVENHAGE
- 6 on **17-10-2016**
- 7 by the court registrar
- 8 no **2016-9660**
- 9 Seal/stamp:

10. Signature:

E. van Houte

E. van Houte

1. Prescription

Clinical pathology. Histological processing, embedding, archiving, analysis and transportation of biological materials.

2. Description and principle of operation

2.1 Method principle

Tissue-processing Tissue-Tec® III Uni-Cassette®

Design styles: White, Aqua, Blue, Yellow, Green, Gold, Red, Lilac, Orange, Pink, Tan, Orange Burn

Tissue-processing Tissue-Tec® III Uni-Cassette® is a polymeric device with holes; it consists of a container and an attached cover. The cover features a double system of plastic push locks. Biomaterials are placed inside. Locking of a cover creates an enclosed space with holes, which do not hinder circulation of solutions.

Pic. 1. Aqua Tissue-processing Tissue-Tec® III Uni-Cassette®

Tissue-processing Stacks Tissue-Tec® III Uni-Cassette®

Design styles: Orange, White, Aqua, Blue, Green, Red, Pink, Gold, Lilac, Orange Burn, Tan, Yellow

Tissue-processing Stacks Tissue-Tec® III Uni-Cassette® is a polymeric device with holes; it consists of a container and an attached cover. The cover features a double system of plastic push locks. Biomaterials are placed inside. Locking of a cover creates an enclosed space with holes, which do not hinder circulation of solutions. Specific stacking of tissue-processing cassettes in tubes allows automatic identification labeling of cassette surface with cassette-labeling printers.

Pic. 2. Red Tissue-processing Stacks Tissue-Tec® III Uni-Cassette®

Stacks Tissue-Tec® III Biopsy Uni-Cassette®

Design styles: Grey, Red, Blue, Yellow, Pink, Green, Orange, White

Stacks Tissue-Tec® III Biopsy Uni-Cassette® is a polymeric device with holes; it consists of a container and an attached cover. The cover features a double system of plastic push locks. Biomaterials are placed inside. Locking of a cover creates an enclosed space with holes, which do not hinder circulation of solutions. Specific stacking of biopsy cassettes in tubes allows automatic identification labeling of cassette surface with cassette-labeling printers.

Pic. 3. Orange Stacks Tissue-Tec® III Biopsy Uni-Cassette®

Tissue-Tec® III Biopsy Uni-Cassette®

Design styles: White, Blue, Yellow, Green, Red, Orange, Pink, Grey

Tissue-Tec® III Biopsy Uni-Cassette® is a polymeric device with holes; it consists of a container and an attached cover. The cover features a double system of plastic push locks. Biomaterials are placed inside. Locking of a cover creates an enclosed space with holes, which do not hinder circulation of solutions.

Pic. 4. Orange Tissue-Tec® III Biopsy Uni-Cassette®

Super Mega-Cassette

Design styles: White, Blue, Dark Blue, Pink, Yellow, Green, Dark Green, Red, Peach, Purple, Turquoise, Brown, Tan, Orange, Grey

Super Mega-Cassette is a polymeric device with holes; it consists of a container and an attached cover. The cover features a double system of plastic push locks. Biomaterials are placed inside. Locking of a cover creates an enclosed space with holes, which do not hinder circulation of solutions.

Pic. 5. Turquoise Super Mega-Cassette

Embedding Cassette Tissue-Tek® II Processing

Design styles: White, Aqua, Red, Blue, Yellow, Green, Gold, Tan, Lilac, Orange Burn, Orange, Clear, Pink, Grey

Embedding Cassette Tissue-Tek® II Processing is a container for preparing a paraffin block with biomaterials placed inside. When used in combination with a Processing/Embedding Cassette Metal Cover, it forms an enclosed space with holes, which do not hinder circulation of solutions.

Pic. 6. White Embedding Cassette Tissue-Tek® II Processing

Stacks Embedding Cassette Tissue-Tek® II Processing

Design styles: White, Aqua, Red, Blue, Yellow, Green, Gold, Tan, Lilac, Orange Burn, Orange, Clear, Pink, Grey

Stacks Embedding Cassette Tissue-Tek® II Processing is a container for preparing a paraffin block with biomaterials placed inside. When used in combination with a Processing/Embedding Cassette Metal Cover, it forms an enclosed space with holes, which do not hinder circulation of solutions. Specific stacking of processing/embedding cassettes in tubes allows automatic identification labeling of cassette surface with cassette-labeling printers.

Pic. 7. White Stacks Embedding Cassette Tissue-Tek® II Processing

Embedding Cassette Metal Cover

Processing/Embedding Cassette Metal Cover is used in combination with Processing/Embedding Cassettes and Stacks Processing/Embedding Cassettes. It forms an enclosed space with holes, which do not hinder circulation of solutions.

Biomaterials to be processed are placed inside.

Pic. 8. Processing/Embedding Cassette Metal Cover

2.2 Materials:

Name of a medical device

Tissue-processing Tissue-Tec® III Uni-Cassette®

Tissue-processing Stacks Tissue-Tec® III Uni-Cassette®

Stacks Tissue-Tec® III Biopsy Uni-Cassette®

Tissue-Tec® III Biopsy Uni-Cassette®

Super Mega-Cassette

Embedding Cassette Tissue-Tek® II Processing

Stacks Embedding Cassette Tissue-Tek® II Processing

Processing/Embedding Cassette Metal Cover

Material

POM (polyoxymethylene)

POM (polyoxymethylene)

POM (polyoxymethylene)

POM (polyoxymethylene)

POM (polyoxymethylene)

POM (polyoxymethylene)

POM (polyoxymethylene)

Stainless steel

According to general information on material properties, POM (polyoxymethylene) is resistant to weak acids, weak and strong alkaline solutions, organic solvents, such as gasoline, kerosene, oils and alcohols. Stainless steel is rather resistant to weak acids, coloring agents, and organic solvents.

3. Analytical characteristics

Tissue-processing Tissue-Tec® III Uni-Cassette®

Design styles: White, Aqua, Blue, Yellow, Green, Gold, Red, Lilac, Orange, Pink, Tan, Orange Burn

Tissue-processing Tissue-Tec® III Uni-Cassette® is used to arrange a tissue sample of a human, animal or plant and subsequent processing of this material with chemically active substances for dewatering, defatting and soaking of the sample with paraffin environment. Cassette containers are used on the next stage of material preparation for preparing a paraffin block with biomaterials placed inside as block containers, as well as on the following stages for sample microtomy and storage of a finished paraffin block with biomaterials placed inside.

Holes in the cassette cover are arranged all over the effective area, in four uniform-sized rectangular areas. Holes within each area are stacked on top of each other in relation to the long side of holes, i.e. arranged in columns. Each area has two columns with seven holes per column.

Basic holes in the cassette container are arranged all over the effective area of its bottom, in four uniform-sized areas. Holes within each area are stacked on top of each other in relation to the long side of holes, i.e. arranged in columns. Each area has two columns with seven holes per column. Each area also features 5 small-sized accessory holes, located in the point of conjunction with a short vertical side of the cassette

The cover is attached to the container with the help of two thin connecting strips made of the material identical to the cassette. Connecting strips remain intact in the course of cassette transportation and storage. However, their integrity may be breached during first-time use (when the user tries to lock the cassette).

Materials used with tissue-processing cassettes of all the above-mentioned types: biological and autopsy material of human and animal origin.

External dimensions	40 ± 1 x 25 ± 1 x 6.0 ± 0.5 mm
Hole size	5.2 ± 0.3 x 0.9 ± 0.1 mm
Cassettes feature angled surface for labeling	45 degrees angle
Capacity (volume)	3.56 ml ± 0.5 ml
Number of holes	Functional holes in the cover: 56 In the container: 56 basic holes and 20 accessory holes Total: 132

Tissue-processing Stacks Tissue-Tec® III Uni-Cassette®

Design styles: Orange, White, Aqua, Blue, Green, Red, Pink, Gold, Lilac, Orange Burn, Tan, Yellow

Tissue-processing Stacks Tissue-Tec® III Uni-Cassette® is used to arrange a tissue sample of a human, animal or plant and subsequent processing of this material with chemically active substances for dewatering, defatting and soaking of the sample with paraffin environment. Cassette containers are used on the next stage of material preparation for preparing a paraffin block with biomaterials placed inside as block containers, as well as on the following stages for sample microtomy and storage of a finished paraffin block with biomaterials placed inside.

Holes in the cassette cover are arranged all over the effective area, in four uniform-sized rectangular areas. Holes within each area are stacked on top of each other in relation to the long side of holes, i.e. arranged in columns. Each area has two columns with seven holes per column.

Basic holes in the cassette container are arranged all over the effective area of its bottom, in four uniform-sized areas. Holes within each area are stacked on top of each other in relation to

the long side of holes, i.e. arranged in columns. Each area has two columns with seven holes per column. Each area also features 5 small-sized accessory holes, located in the point of conjunction with a short vertical side of the cassette

The cover is attached to the container with the help of two thin connecting strips made of the material identical to the cassette. Connecting strips remain intact in the course of cassette transportation and storage. However, their integrity may be breached during first-time use (when the user tries to lock the cassette).

Materials used with tissue-processing cassettes of all the above-mentioned types: biological and autopsy material of human and animal origin.

External dimensions	$40 \pm 1 \times 25 \pm 1 \times 6.0 \pm 0.5$ mm
Hole size	$5.2 \pm 0.3 \times 0.9 \pm 0.1$ mm
Cassettes feature angled surface for labeling	45 degrees angle
Capacity (volume)	3.56 ml $\pm$ 0.5 ml
Number of holes	Functional holes in the cover: 56 In the container: 56 basic holes and 20 accessory holes Total: 132

Stacks Tissue-Tec® III Biopsy Uni-Cassette®

Design styles: Grey, Red, Blue, Yellow, Pink, Green, Orange, White

Stacks Uni-Cassette® Biopsy is used to arrange a tissue sample of a human, animal or plant and subsequent processing of this material with chemically active substances for dewatering, defatting and soaking of the sample with paraffin environment. Cassette containers are used on the next stage of material preparation for preparing a paraffin block with biomaterials placed inside as block containers, as well as on the following stages for sample microtomy and storage of a finished paraffin block with biomaterials placed inside.

Holes in the cassette cover are arranged all over the effective area, in four uniform-sized rectangular areas. Holes within each area are stacked on top of each other in relation to the long side of holes, i.e. arranged in columns. Each area has two columns with seven holes per column.

Basic holes in the cassette container are arranged all over the effective area of its bottom, in four uniform-sized areas. Holes within each area are stacked on top of each other in relation to the long side of holes, i.e. arranged in columns. Each area has two columns with seven holes per column. Each area also features 5 small-sized accessory holes, located in the point of conjunction with a short vertical side of the container.

The cover is attached to the container with the help of two thin connecting strips made of the material identical to the cassette. Connecting strips remain intact in the course of cassette transportation and storage. However, their integrity may be breached during first-time use (when the user tries to lock the cassette).

Materials used with tissue-processing cassettes of all the above-mentioned types: biological and autopsy material of human and animal origin.

External dimensions	$40 \pm 1 \times 25 \pm 1 \times 6.0 \pm 0.5$ mm
Cassettes feature angled surface for labeling	45 degrees angle
Hole size	$0.9 \pm 0.1 \times 0.9 \pm 0.1$ mm
Capacity (volume)	3.56 ml $\pm$ 0.5 ml
Number of holes	Functional holes in the cover: 168 In the container: 168 basic holes and 20 accessory holes

Tissue-Tec® III Biopsy Uni-Cassette®

Design styles: White, Blue, Yellow, Green, Red, Orange, Pink, Grey

Tissue-Tec® III Biopsy Uni-Cassette® is used to arrange a tissue sample of a human, animal or plant and subsequent processing of this material with chemically active substances for dewatering, defatting and soaking of the sample with paraffin environment. Cassette containers are used on the next stage of material preparation for preparing a paraffin block with biomaterials placed inside as block containers, as well as on the following stages for sample microtomy and storage of a finished paraffin block with biomaterials placed inside.

Holes in the cassette cover are arranged all over the effective area, in four uniform-sized rectangular areas. Holes within each area are stacked on top of each other in relation to the long side of holes, i.e. arranged in columns. Each area has two columns with seven holes per column.

Basic holes in the cassette container are arranged all over the effective area of its bottom, in four uniform-sized areas. Holes within each area are stacked on top of each other in relation to the long side of holes, i.e. arranged in columns. Each area has two columns with seven holes per column. Each area also features 5 small-sized accessory holes, located in the point of conjunction with a short vertical side of the container.

The cover is attached to the container with the help of two thin connecting strips made of the material identical to the cassette. Connecting strips remain intact in the course of cassette transportation and storage. However, their integrity may be breached during first-time use (when the user tries to lock the cassette).

Materials used with tissue-processing cassettes of all the above-mentioned types: biological and autopsy material of human and animal origin.

External dimensions	$40 \pm 1 \times 25 \pm 1 \times 6.0 \pm 0.5$ mm
Cassettes feature angled surface for labeling	45 degrees angle
Hole size	$0.9 \pm 0.1 \times 0.9 \pm 0.1$ mm
Capacity (volume)	$3.56 \text{ ml} \pm 0.5 \text{ ml}$
Number of holes	Functional holes in the cover: 168 In the container: 168 basic holes and 20 accessory holes Total: 356

Super Mega-Cassette

Design styles: White, Blue, Dark Blue, Pink, Yellow, Green, Dark Green, Red, Peach, Purple, Turquoise, Brown, Tan, Orange, Grey

Super Mega-Cassettes are used to arrange a tissue sample of a human, animal or plant and subsequent processing of this material with chemically active substances for dewatering, defatting and soaking of the sample with paraffin environment. Cassette containers are used on the next stage of material preparation for preparing a paraffin block with biomaterials placed inside as block containers, as well as on the following stages for sample microtomy and storage of a finished paraffin block with biomaterials placed inside.

Super Mega-Cassettes are supplied with covers that are not initially attached to the container. The cassette is assembled in the course of operation.

Holes in the cassette cover are arranged all over the effective area, in four uniform-sized rectangular areas. Holes within each area are stacked on top of each other in relation to the long side of holes, i.e. arranged in columns. Each area has three columns with 11 holes per column.

Materials used with tissue-processing cassettes of all the above-mentioned types: biological and autopsy material of human and animal origin.

External dimensions	75 ± 1 x 52 ± 1 x 18.0 ± 0.5 mm
Hole size	5.2 ± 0.3 x 0.9 ± 0.1 mm
Cassettes feature angled surface for labeling	45 degrees angle
Capacity (volume)	46.5 ml ± 1 ml
Number of holes	Functional holes in the cover: 132 In the container: 132 Total: 264

Embedding Cassette Tissue-Tek® II Processing

Design styles: White, Aqua, Red, Blue, Yellow, Green, Gold, Tan, Lilac, Orange Burn, Orange, Clear, Pink, Grey

Embedding Cassette Tissue-Tek® II Processing are used for material preparation for preparing a paraffin block with biomaterials placed inside, as well as on the following stages for sample microtomy and storage of a finished paraffin block with biomaterials placed inside. In combination with Processing/Embedding Cassette Metal Covers, Stacks Processing/Embedding Cassettes are used to arrange a tissue sample of a human, animal or plant and subsequent processing of this material with chemically active substances for dewatering, defatting and soaking of the sample with paraffin environment.

Holes are arranged all over the effective area of the cassette container bottom. Holes are stacked on top of each other, i.e. arranged in columns. In total, the cassette has five columns with 7 holes per column and a missing central hole in the central row.

External dimensions	40 ± 1 x 25 ± 1 x 6.0 ± 0.5 mm
Cassettes feature angled surface for labeling	45 degrees angle
Circular holes	D = 2.4 ± 0.1 mm
Capacity (volume)	3.7 ml ± 0.5 ml
Number of holes	34

Stacks Embedding Cassette Tissue-Tek® II Processing

Design styles: White, Aqua, Red, Blue, Yellow, Green, Gold, Tan, Lilac, Orange Burn, Orange, Clear, Pink, Grey

Stacks Embedding Cassette Tissue-Tek® II Processing are used for material preparation for preparing a paraffin block with biomaterials placed inside, as well as on the following stages for sample microtomy and storage of a finished paraffin block with biomaterials placed inside. In combination with Processing/Embedding Cassette Metal Covers, Stacks Processing/Embedding Cassettes are used to arrange a tissue sample of a human, animal or plant and subsequent processing of this material with chemically active substances for dewatering, defatting and soaking of the sample with paraffin environment.

Holes are arranged all over the effective area of the cassette container bottom. Holes are stacked on top of each other, i.e. arranged in columns. In total, the cassette has five columns with 7 holes per column and a missing central hole in the central row.

External dimensions	40 ± 1 x 25 ± 1 x 6.0 ± 0.5 mm
Cassettes feature angled surface for labeling	45 degrees angle
Circular holes	D = 2.4 ± 0.1 mm
Capacity (volume)	3.7 ml ± 0.5 ml

Number of holes	34
-----------------	----

Processing/Embedding Cassette Metal Cover

In combination with Processing/Embedding Cassettes and Stacks Processing/Embedding Cassettes, Processing/Embedding Cassette Metal Covers are used to arrange a tissue sample of a human, animal or plant and subsequent processing of this material with chemically active substances for dewatering, defatting and soaking of the sample with paraffin environment.

Holes are arranged all over the effective area of the cassette cover. Holes are stacked on top of each other, i.e. arranged in columns. In total, the cassette has five columns: two leftmost and rightmost columns have 7 holes per column, and a central row has 6 holes.

Dimensions	$34 \pm 1 \times 17 \pm 1 \times 1.5 \pm 0.3$ mm
Circular holes	$D = 2.4 \pm 0.1$ mm
Features fixturing lobes on short edges	
Number of holes	34

4. Precautions

Safety regulations and personal hygiene should be observed when handling polymeric medical devices for collection, transportation, analysis and storage of biomaterials with accessories.

Ensure draft and well-functioning ventilation when using the device. Use disposable rubber gloves. Keep away from open flame. Storage and intake of food at the workstation is forbidden.

Avoid contact with skin, spillage or accidental swallowing of reference solutions. Induce vomiting in case of swallowing. Immediately consult a doctor if feeling unwell.

Observe security requirements and legislative instruments of the country of use when handling polymeric medical devices for collection, transportation, analysis and storage of biomaterials with accessories.

5. Preparation

All the aforementioned polymeric medical devices for collection, transportation, analysis and storage of biomaterials with accessories are non-sterile disposable medical devices.

Tissue-processing Tissue-Tec® III Uni-Cassette®, Tissue-processing Stacks Tissue-Tec® III Uni-Cassette®, Stacks Tissue-Tec® III Biopsy Uni-Cassette®, Tissue-Tec® III Biopsy Uni-Cassette®, Super Mega-Cassettes, Embedding Cassette Tissue-Tek® II Processing, Stacks Embedding Cassette Tissue-Tek® II Processing and Processing/Embedding Cassette Metal Cover are medical devices that do not require special preparation for operation.

Disinfection is possible as needed.

Polymeric medical devices for collection, transportation, analysis and storage of biomaterials with accessories are disinfection-resistant. Disinfection is conducted by using chemical substances, autoclaving or air.

Chemical disinfection is conducted by means of chlorine substances or other disinfectant detergents. Please note that chemical substances may be toxic upon contact with skin or if inhaled; use personal protective equipment (gloves, apron and mask). Chemical disinfectants may be supplied ready to use or require precise dilution to specific concentration. Rinse off all visible dirt from the devices with water and a detergent before disinfection; use additional accessories (brushes, sponges) if possible. After that, put devices in a disinfection solution for a specified time. Duration of exposure is determined on the basis of a detergent's manual. After disinfection, rinse devices with flowing water and dry.

Autoclave disinfections is conducted at $+110\text{ }^{\circ}\text{C} \pm 2\text{ }^{\circ}\text{C}$, 0.05 mPa (0.5 kgf/cm) for 20 minutes (maximum deviation – 5 minutes).

Dry-air sterilizer disinfection is conducted at $+120\text{ }^{\circ}\text{C} \pm 4\text{ }^{\circ}\text{C}$ for 45 minutes (maximum deviation – 5 minutes).

Medical devices are resistant to chemical sterilization (with chemical solutions).

Hydrogen peroxide solution 6% may be used to sterilize medical devices; duration of sterilization exposure varies depending on the temperature: $+18\text{ }^{\circ}\text{C} \pm 2\text{ }^{\circ}\text{C}$ – 360 minutes ± 5 minutes; $+50\text{ }^{\circ}\text{C} \pm 2\text{ }^{\circ}\text{C}$ – 180 minutes ± 5 minutes.

Devices should be fully submerged in hydrogen peroxide solution 6% during the sterilization exposure; all ports and cavities should be filled with the solution.

6. Analysis

Tissue-processing Tissue-Tec® III Uni-Cassette®, Tissue-Tec® III Biopsy Uni-Cassette®, Super Mega-Cassette

Histological examination:

- Put the dissected material in the container of a tissue-processing cassette, a biopsy cassette or a super mega-cassette (use biopsy plates or bags as needed).
- Secure cover of the tissue-processing cassette or the biopsy cassette on the container (lock the cassette).
- Label one of three labeling surfaces with a laboratory permanent marker.
- Process the material in a conventional manner (manually or automatically).

Histological embedding:

- After material processing has been completed, open cover of the tissue-processing cassette, the biopsy cassette or the super mega-cassette: disconnect it from the container and remove.
- Embed histological material in a conventional manner using container of the tissue-processing cassette or the biopsy cassette for paraffin block.

Tissue-processing Stacks Tissue-Tec® III Uni-Cassette®, Stacks Tissue-Tec® III Biopsy Uni-Cassette®

Histological examination:

- Use cassette printer to label containers of tissue-processing stacks cassettes or biopsy stacks cassettes.
- Put the dissected material in the container of a marked tissue-processing cassette, a biopsy cassette or a super mega-cassette (use biopsy plates or bags as needed).
- Secure cover of the tissue-processing cassette or the biopsy cassette on the container.
- Lock cover of the tissue-processing or the biopsy cassette until it clicks.
- Process the material in a conventional manner (manually or automatically).

Histological embedding:

- After material processing has been completed, open cover of the tissue-processing cassette, the biopsy cassette or the super mega-cassette: disconnect it from the container and remove.
- Embed histological material in a conventional manner using container of the tissue-processing cassette or the biopsy cassette for paraffin block.

Embedding Cassette Tissue-Tek® II Processing

Histological examination:

- Put the dissected material in the container of a tissue-processing cassette or a biopsy cassette (use biopsy plates or bags as needed).

- Secure Processing/Embedding Cassette Metal Cover on the Processing/Embedding Cassette.
- Label the device with a laboratory permanent marker.
- Process the material in a conventional manner (manually or automatically).

Histological embedding:

- After material processing has been completed, open cover of the Processing/Embedding Cassette: disconnect it from the container and remove.
- Embed histological material in a conventional manner using container of the Processing/Embedding Cassette for paraffin block.

Stacks Embedding Cassette Tissue-Tek® II Processing

Histological examination:

- Use cassette printer to label containers of stacks processing/embedding cassettes.
- Put the dissected material in the container of a tissue-processing cassette or a biopsy cassette (use biopsy plates or bags as needed).
- Secure Processing/Embedding Cassette Metal Cover on the Processing/Embedding Cassette.
- Label the device with a laboratory permanent marker.
- Process the histological material in a conventional manner (manually or automatically).

Histological embedding:

- After material processing has been completed, open cover of the Processing/Embedding Cassette: disconnect it from the container and remove.
- Embed histological material in a conventional manner using container of the Stacks Processing/Embedding Cassette for paraffin block.

7. Safety requirements

All cassette covers are tightly connected with the container in the course of operation, which excludes the possibility of spontaneous unlocking with further loss of diagnostic material. Polymeric medical devices for collection, transportation, analysis and storage of biomaterials with accessories belong to "B" Class according to Sanitary Regulations and Standards SanPiN 2.1.7.2790-10 based on their epidemiological, toxicological and radiation hazard level.

8. Storage and handling conditions

Use polymeric medical devices for collection, transportation, analysis and storage of biomaterials with accessories at +10-35 °C.

Keep devices away from direct sunlight, sources of heat and open flame. Do not expose devices to rough mechanical impact.

Average lifetime – 5 years.

9. Requirements for the protection of the environment. Disposal conditions

Polymeric medical devices for collection, transportation, analysis and storage of biomaterials with accessories are clean and environmentally friendly products.

Detoxify waste before disposal. The detoxification method is selected on the basis of capabilities of the facility, where the medical device is used. Dispose of the medical advice in accordance with norms, regulations and legislative instruments of the country of use.

10. Packaging:

Medical devices are transported in carton boxes and polyethylene bags.

Name	Design styles	Shipping package	Mass of one MD (g)	Package mass (kg)
Tissue-processing Tissue-Tec® III Uni-Cassette®	White, Aqua, Blue, Yellow, Green, Gold, Red, Lilac, Orange, Pink, Tan, Orange Burn	Carton box; dimensions: 350 x 150 x 250 mm; 500 pcs	3.4 ± 0.5	1.8
		Carton box; dimensions: 350 x 290 x 250 mm; 2,000 pcs (4 boxes; each box holds 500 pcs)		7.2
Tissue-processing Stacks Tissue-Tec® III Uni-Cassette®	Orange, White, Aqua, Blue, Green, Red, Pink, Gold, Lilac, Orange Burn, Tan, Yellow	Carton box; dimensions: 350 x 150 x 250 mm; 400 pcs (10 transparent plastic tubes; each tube holds 40 cassettes. Every 200 cassettes are put in a polyethylene bag)	3.4 ± 0.5	1.65
Stacks Tissue-Tec® III Biopsy Uni-Cassette®	Grey, Red, Blue, Yellow, Pink, Green, Orange, White	Carton box; dimensions: 350 x 150 x 250 mm; 400 pcs (10 transparent plastic tubes; each tube holds 40 cassettes. Every 200 cassettes are put in a polyethylene bag)	3.4 ± 0.5	1.65
Tissue-Tec® III Biopsy Uni-Cassette®	White, Blue, Yellow, Green, Red, Orange, Pink, Grey	Carton box; dimensions: 350 x 150 x 250 mm; 500 pcs	3.4 ± 0.5	1.9
		Carton box; dimensions: 350 x 290 x 250 mm; 2,000 pcs (4 carton boxes; each box holds 500 pcs)		7.2
Super Mega-Cassette	White, Blue, Dark Blue, Pink, Yellow, Green, Dark Green, Red,	Carton box; dimensions: 295 x 170 x 325 mm; 100 pcs	26 ± 1.2	2.8
		Carton box;		22.4

	Peach, Purple, Turquoise, Brown, Tan, Orange, Grey	dimensions: 700 x 630 x 350 mm; 800 pcs (8 boxes; each box holds 100 pcs)		
Processing/Embedding Cassette Tissue-Tek® II Processing	White, Aqua, Red, Blue, Yellow, Green, Gold, Tan, Lilac, Orange Burn, Orange, Clear, Pink, Grey	Carton box; dimensions: 350 x 290 x 250 mm; 2,000 pcs (4 polyethylene bags; each bag holds 500 pcs)	2.5 ± 0.5	6.0
Stacks Processing/Embedding Cassette Tissue-Tek® II Processing	White, Aqua, Red, Blue, Yellow, Green, Gold, Tan, Lilac, Orange Burn, Orange, Clear, Pink, Grey	Carton box; dimensions: 320 x 180 x 260 mm; 800 pcs (20 transparent plastic tubes inside the box; each tube holds 40 cassettes)	2.5 ± 0.5	2.4
Processing/Embedding Cassette Metal Cover		Polyethylene bag; 25 pcs	2.5 ± 0.5	70

11. Labeling

Device labeling contains the following information: reference number, name, packaging, lot number (LOT), informational signs*, product information (color / volume / size / number of glasses / shape / presence of cover / presence of calibrated scale), production date, storage conditions, manufacturer and supplier details, contact information for reclamations.

Reuse is FORBIDDEN

WARNING! READ AND ACKNOWLEDGE ACCOMPANYING DOCUMENTS

In vitro diagnostic medical device

The product meets the essential requirements of EC directives and harmonized standards of the European Union; the product conformity to directives was assessed

12. Manufacturer warranty

Concerning polymeric devices for tissue processing please refer to 68A Prospekt Bolshoy V.O., Saint Petersburg, 199106, Russian Federation; tel.: +7 (812) 305 0606; e-mail: info@biovitrum.ru

13. International certificates and declarations

Polymeric medical devices for collection, transportation, analysis and storage of biomaterials with accessories have the following certificates:

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