I, **Michelle Metzger** hereby swear (or affirm) that the attached **Instructions for Use** for **Kitewire Deep Peripheral Guidewire** is an unaltered copy of the original document.

Country: Russia

Michelle Metzger
Michelle Metzger

9.22.15
Date

A notary public or other officer completing this certificate verifies only the identity of the individual who signed the document, to which this certificate is attached, and not the truthfulness, accuracy, or validity of that document.

State of California }

County of Sonoma }

Subscribed and sworn to (or affirmed) before me this 22nd day of September, 2015, by **Michelle Metzger** proved to me on the basis of satisfactory evidence to be the person who appeared before me.

Candi L. Nance

Candi L. Nance, Notary Public
Commission Number: 2125324
Commission Expiration Date: August 29, 2019

Description of Attached Document(s):

Document Description: Instructions for Use for Kitewire Deep Peripheral Guidewire

Number of pages of document: (32) Pages

Peripheral Guidewire 0.014" / 0.36mm

EN	Peripheral Guidewire 0.014" / 0.36mm Instructions for use	Page 5	FI	Perifeerinen ohjainlanka 0,014" / 0,36 mm Käyttöohjeet	Sivu 18
DE	Peripherer Führungsdraht, 0,014" / 0,36 mm Gebrauchsanleitung	Seite 6	CS	Periferní vodičí drát 0,014" / 0,36mm Návod k použití	Strana 19
FR	Guide périphérique de 0,36 mm / 0,014" Mode d'emploi	Page 7	SK	Periférny vodiaci drôt 0,014" / 0,36 mm Návod na použitie	Strana 20
IT	Filoguida periferico da 0,36 mm (0.014") Istruzioni per l'uso	Pagina 8	HU	Perifériás vezetődrót 0,014" / 0,36 mm Használati útmutató	21. oldal
ES	Guía periférica de 0,014" / 0,36 mm Instrucciones de uso	Página 9	KO	말초 가이드와이어 0.014" / 0.36mm 사용 설명서	22 페이지
NL	Perifere voerdraad 0,36 mm / 0,014 inch Gebruiksaanwijzing	Pagina 10	TR	Periferal Kılavuz Tel 0,014" / 0,36 mm Kullanım talimatları	Sayfa 23
DA	Perifer ledetråd 0,014" / 0,36 mm Brugsanvisning	Side 11	ET	Perifeerne juhtetraat 0.014" / 0,36 mm Kasutusjuhend	Lk 24
PL	Prowadnik obwodowy 0,014 cala / 0,36 mm Instrukcja użytkowania	Strona 12	SR	Periferna vodič-žica od 0,014 inča / 0,36 mm Upustvo za upotrebu	Stranica 25
PT	Fio-guia periférico de 0,014" / 0,36 mm Instruções de utilização	Página 13	LT	Periferinė kreipiamoji viela (0,014 colio / 0,36 mm) Naudojimo instrukcijos	26 puslapis
ZH	外周血管内导丝 0.014" / 0.36mm 使用说明	第 14 页	HR	Žica vodičica za periferno postavljanje 0,36 mm / 0.014" Upute za uporabu	Stranica 27
EL	Περιφερικό οδηγό σύρμα 0,014" / 0,36mm Οδηγίες χρήσης	Σελίδα 15	RO	Fir de ghidaj periferic 0,014" / 0,36 mm Instrucțiuni de utilizare	Pagina 28
NO	Perifer ledevaier 0,014" / 0,36 mm Bruksanvisning	Side 16	BG	Периферен водач 0,014" / 0,36 mm Инструкции за употреба	Страница 29
SV	Perifer ledare 0,36 mm / 0,014" Bruksanvisning	Sidan 17	RU	Периферический проводник 0,014 дюйма / 0,36 мм Инструкция по применению	Страница 30

BRIVANT LIMITED
Parkmore West Business Park,
Galway,
County Galway,
Ireland.

Tel: +353 91 385037
Fax: +353 91 766598

0050

Invatec S.p.A.
Via Martiri della Libertà, 7
25030 Roncadelle (BS), Italy
Tel: +39 030 2589311
Fax: +39 030 2589312

INVATEC
Innovative Technologies
We make ideas come alive

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DEVICE NAME

Kitewire™ Deep is a peripheral guidewire.

CAUTION

This device should be used only by physicians trained in angiography and percutaneous transluminal coronary angioplasty (PTCA), and / or percutaneous transluminal angioplasty (PTA).

Do not use if package is damaged.

Do not re-sterilize or re-use.

Do not wipe with alcohol.

Single Use Only.

- Reuse of single-use devices creates a potential risk of patient or user infections. Contamination of the device may lead to injury, illness or death of the patient.
- Cleaning, disinfection and sterilization may compromise essential material and design characteristics leading to device failure.

See Instructions for Use

Sterilized with ethylene oxide gas

Use the guidewire prior to the "Use before" date on the package label, preceded by the symbol .

Refer to the instructions supplied with any interventional devices to be used in conjunction with the Kitewire™ Deep for their intended uses, contraindications, and potential complications.

DESCRIPTION

The Kitewire™ Deep is a steerable guidewire available in several lengths and diameters. The distal tip is shapeable. Refer to the product label for product specifications (e.g. wire length, diameter and length of tip radiopacity).

STORAGE. Store in a cool, dark, dry place.

INDICATIONS FOR USE

The Kitewire™ Deep guidewires are intended for use in the coronary and peripheral vasculature.

CONTRAINDICATIONS

The Kitewire™ Deep is not intended for use in the cerebral vasculature. Patients judged not acceptable for percutaneous intervention (PCI).

PRECAUTIONS

- Failure to follow the instructions may compromise guidewire performance and result in complications.
- Prior to use, confirm compatibility of guidewire outer diameter with the balloon catheter.
- The tip section of the guidewire has a proper orientation for shaping. Identify the flexing plane before shaping. Shape in the same plane as that for flexure.
- Guidewire advancement, withdrawal, and torquing should be monitored by fluoroscopy.

DIRECTIONS FOR USE**A. Over-the-wire PTCA/PTA Balloon Catheters**

1. Prepare the balloon catheter according to manufacturer's directions.
2. If desired, shape the guidewire tip using standard technique. When preparing to shape the tip section, hold the tip at the middle joint and between the two fingers. Gently "brush" the spring coil to identify the plane of flexure and complete the shaping procedure.
3. Moisten the guidewire with sterile saline to increase the surface lubricity.
4. Place the guidewire into the dilatation catheter and advance it to the catheter lumen.
5. Remove the guidewire introducer by sliding it over the proximal end of the guidewire. Secure the guidewire within the haemostatic valve being careful not to over-tighten the compression fitting.
6. A torque application device may be applied to the proximal end of the guidewire.
7. Advance the guidewire into the target vessel and across the lesion using fluoroscopy to facilitate proper guidewire placement.
8. Holding the guidewire in position, advance the balloon catheter over the guidewire and into the target lesion.

9. Complete the procedure and remove the guidewire and balloon catheter according to procedural protocol.

B. Rapid Exchange Systems

1. If desired, shape the guidewire tip using standard technique. When preparing to shape the tip section, hold the tip at the middle joint and between the two fingers. Gently "brush" the spring coil to identify the plane of flexure and complete the shaping procedure.
2. Moisten the guidewire with sterile saline to increase surface lubricity.
3. Introduce the guidewire through the guiding catheter Y-adapter using a guidewire introducer.
4. Remove the guidewire introducer by sliding it over the proximal end of the guidewire. Secure the guidewire within the haemostatic valve being careful not to over-tighten the compression fitting.
5. Advance the guidewire into the target vessel and across the lesion using fluoroscopy to facilitate proper guidewire placement. A torque application device may be applied to the proximal end of the guidewire. Manifold flushing of the catheter lumen during guidewire advancement is recommended to remove any residual air.
6. Place the rapid-exchange balloon catheter on the guidewire and advance it into position according to manufacturer's instructions.
7. Holding the guidewire in position, advance the balloon catheter over the guidewire and into the target lesion.
8. Complete the procedure and remove the guidewire and balloon catheter according to procedural protocol.

WARNING

A guidewire is a delicate instrument and must not be advanced, withdrawn, or torqued if resistance is met. Guidewire manipulations must always be observed under fluoroscopy.

If the guidewire is removed and is to be re-inserted, it must be inspected for signs of damage (weakened or kinked segments) prior to re-introduction. Do not re-introduce if guidewire is weakened or kinked.

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I hereby certify that the multi-page document attached hereto is an "Appendix to the Instruction for Use for Kitewire Deep product."

Michael Dunning
QA/RA Manager

MARTIN J. KEARNS
NOTARY PUBLIC
1 DEVON PLACE
THE CRESCENT
GALWAY H91 N5H2
IRELAND
Telephone 091-589219
COMMISSIONED FOR LIFE

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
Martin J. Kearns

MARTIN J. KEARNS
NOTARY PUBLIC
1 DEVON PLACE
THE CRESCENT
GALWAY, H91 N5H2
IRELAND
Telephone 091-589219
COMMISSIONED FOR LIFE

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Martin J. Kearns

MARTIN J. KEARNS
NO LARY PUBLIC
1 DEVON PLACE
THE CRESCENT
GALWAY H91 N5H2
IRELAND
Telephone 091-589219
COMMISSIONED FOR LIFE

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Martin J. Kearns

MARTIN J. KEARNS
NOTARY PUBLIC
1 DEVON PLACE
THE CRESCENT
GALWAY H91 N5H2
IRELAND
Telephone 091-589219
COMMISSIONED FOR LIFE

- name of the manufacturer;
- manufacturing site;
- code UPN (unique product number);
- code CFN (customer facing number);
- date and number of Registration Certificate;
- sign of conformity in Russia ;
- warning and informing symbols and inscriptions.

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UTILIZATION

After use, this products may be a potential biohazard (hazard class: B). Handle and dispose of all such devices in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

The product dispose of in accordance with applicable law and the instructions of the medical institution

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STERILIZATION

The packaged guide wire shall be sterilised on a contract basis by a validated Ethylene Oxide (ETO) gas sterilisation cycle. ETO residuals will be less than 4mg per device after 2X sterilisation in accordance with ETO limits for limited exposure devices stated in ISO 10993-7:2008. The sterilization process will achieve a Sterility Assurance Level (SAL) of 10^{-6} .

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OPERATING CONDITIONS:

Do not use if package is opened or damaged.
Do not re-sterilize or re-use
Do not wipe with alcohol. Single use only. The product is intended for use by trained medical personnel in the operating room: at room

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Martin J. Kearns

MARTIN J. KEARNS
NOTARY PUBLIC
1 DEVON PLACE
THE CRESCENT,
GALWAY H91 N5H2
IRELAND
Telephone 091-589219
COMMISSIONED FOR LIFE

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The translation of this text was performed by me, the translator Mamedov Timur Djavanshirovich. The requirements regarding the text of the translation (maximum correctness, literate statement) have been explained to me.

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
I, Akimov Gleb Borisovich, Notary of the city of Moscow certify the authenticity of signature, made by the translator Mamedov Timur Djavanshirovich in my presence. His identity established.

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It is registered in the register under №
Fee collected: 100 rubles

Notary: /signature/

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Seal: Notary of Moscow Akimov G.B.
ITN 770400047857*1D-1*

Seal: Notary of Moscow Akimov G.B.
ITN 770400047857*1D-1*

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Stitched up, numbered and sealed () pages

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THE MATERIALS FROM WHICH IT IS MADE OR ITS MEDICAL DEVICE BASIC PARTS

denomination	materials type (trade mark)
Guidewire body	302 Stainless steel
Coil	Platinum Tungsten Alloy (Pt 92%/W8%)
Solder	Tin - 95% Silver - 5%
Polymer Jacket	Thermoplastic Polyurethane Estane ®88A
Coating	Distal polymer jacket coated with Harland Hydrophilic coating (Lubricent® Coating) Guidewire shaft coated with PTFE coating (Whitford Xylan 1770/G5841)
Tubes	Petrothene LR734001 & Marlex® HHM 5202BN Polyethylene
Clips for tubes	Petrothene LR734001& Marlex® HHM 5202BN Polyethylene & Westlake EN1817 Low Density Polyethylene
Retention clips for tubes	Thermoplastic elastomer Santoprene 8000 & REMAFIN-WHITE 06411
Materials primary sterile packaging /	A layer of polyethylene - Polyethylene 35709-G (12 PET/ 50 PE-EVA-E)
	A layer of high-density polyethylene - Tyvek® 1073B HCW CR27
	The ends of the package are sealed by heating

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
Weight of the Guidewire	1.18g (for 195cm) 1.54g (for 250 cm) 1.85g (for 300cm)
Surface	Guidewire will not have any foreign bodies or irregularities on the outer surface
Max Diameter	0.014" (0.36 mm)
Usable Length	195, 250, 300 cm $\pm$ 3cm
Coil length	3 cm ($\pm$ 3mm)
Joint length (Solder)	Distal Joint 0.5mm-0.8mm Proximal Joint 1mm-2mm
Location of tip radiopacity	The distal part of the reduced section is covered with a 3cm platinum tungsten spring coil to provide greater visibility on x-ray equipment.
Polymer jacket length	30 cm $\pm$ 1.5cm
Hydrophilic Coating length	30 cm ($\pm$ 1.5 cm) (polymer jacketed portion of guide wire)
Hydrophilic coating thickness	The diameter of the product part with applied hydrophilic coating should be no more than of the prescribed limits for shaft diameter
Tensile strength	Min. 1.1lbs (7,58 kPa)
Torque strength	Min 10 turns to failure of Guidewire
Torque response	Results should be comparable with the predicate devices
Tip flexibility	Results should be comparable within the range with predicate devices:

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MARTIN J. KEARNS
NOTARY PUBLIC
1 DEVON PLACE
THE CRESCENT
GALWAY H91 N5H2
IRELAND
Telephone 091-589219

			Force/ (g)	
	½ Span (12.5mm)	Min	0.065	
		Max	0.848	
	¾ Span (18.75mm)	Min	0.234	
		Max	4.087	
PTCA/PTA Balloon Catheters compatible with 0.014” (0,36 mm) guidewires				
PTCA/PTA Balloon Catheters compatible with 0.014” (0,36 mm) guidewires				

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INFORMATION ABOUT STICKER

The final packaging for the Russian market is accompanied with additional sticker.

Devices supplied in packaging marked with the following possible information:

- name of medical device in Russian language;
- name of the manufacturer;
- manufacturing site;
- code UPN (unique product number);
- code CFN (customer facing number);
- date and number of Registration Certificate;
- sign of conformity in Russia ;
- warning and informing symbols and inscriptions.

Sample of the sticker

	Product name	
	Name of the manufacturer	Code UPN
	Manufacturing site	Code CFN
	See IFU, Date and Number of RC	

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SHELF LIFE:

The shelf life of the product is 4 years.

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WARRANTY

The product and each component of its system have been designed, manufactured, tested and packaged with all reasonable care. The warnings contained in the instructions for use are expressly to be considered as an integral part of this provision. Company warranties the product until the expiry date indicated on the same. The warranty is valid provided that the use of the product was consistent with the instructions for use. Brivant is not liable for any direct, indirect, incidental or consequential damages caused by the product. Except in the case of fraud or grave fault on Brivant's part, compensation of any damage to the buyer will not, in any event, be greater than the invoice price of the disputed products. The guarantee contained in this provision incorporates and substitutes the legal guarantees for defects and compliance, and excludes any other possible liability of Brivant, however originating, from this product supplied. These limitations of liability and warranty are not intended to contravene any mandatory provisions of law applicable. If any clause of the disclaimer is considered by a competent court to be invalid or to be in conflict with the applicable law, the remaining part of it shall not be affected and shall remain in full force and effect. The invalid clause shall be substituted by a valid clause which best reflects Brivant's legitimate interest in limiting its liability or warranty. No person has any authority to bind Brivant to any warranty or liability regarding the product.

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Contact for all questions, including complaints, repair and service of medical devices the Authorized Representative in Russia:

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