Robert C. Meagher
Senior Director
World Wide Regulatory Affairs
Carestream Health, Inc.
150 Verona Street
Rochester, New York 14608
Telephone 585-627-6528

"26" August 2016 г.

Carestream

USER MANUAL

I. DRYVIEW 6950 Laser Imaging System with accessories consisting of:

1. DRYVIEW 6950 Laser Imaging System – 1pc.
2. Operational documentation - 1pc.

II. Accessories:

1. SORTER/DV6950, TX65 – 1pc.
2. Siesmic Kit /DV6950, TX65 – 1pc.
3. Filter is common with DV6950/TX65– 1pc.

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PHOTO, DIAGRAM, DRAWING MEDICAL DEVICE

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NOTICE

All rights reserved. No part of this manual may be reproduced or copied in any form by any mean graphic, electronic, or mechanical, including photocopying, typing, or information retrieval systems without written permission.

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SAFETY MANUAL SAFETY AND RELATED INFORMATION

The information contained herein is based on the experience and knowledge relating to the subject matter gained by the manufacturer prior to publication.

No patent license is granted by this information.

The manufacturer reserves the right to change this information without notice and makes no warranty, express or implied, with respect to this information. The manufacturer shall not be liable for any loss or damage, including consequential or special damages, resulting from the use of this information, even if loss or damage is caused by the manufacturer's negligence or other fault.

Safety, Warnings, and Cautions

Please read and understand all instructions before using this product.

DANGER

- This equipment is operated with hazardous voltage which can shock, burn, or cause death.
- Remove wall plug before servicing equipment. Never pull on cord to remove from outlet. Grasp plug and pull to disconnect. Do not attempt to service or repair the laser imager yourself to avoid exposure to dangerous voltage, laser beam, or other danger. Always call an authorized service provider for any service or repair.
 - Do not operate equipment with a damaged power cord.
 - Do not use an extension cord to power this equipment.
 - Do not operate equipment with any of the safety interlocks overridden.
 - Position the power cord so it will not be tripped over or pulled.
 - Connect this equipment to a grounded wall outlet.
 - A power cord is provided with this equipment. All countries must use an agency-approved power cord with plug type suitable for the country of use. Contact a qualified dealer for help.
 - Do not operate equipment with the covers open.

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**DANGER**

This equipment contains moving parts that may be accessible to the user. Loose clothing, jewelry or long hair may cause personal injury or damage to the equipment.

**DANGER**

This equipment is not contained in a sealed cabinet. Do not use this equipment in locations where it can come in contact with liquids, including body fluids.

**DANGER**

This equipment is not a mobile device. Casters shall only be used by qualified service personnel during product installation. Operators are not allowed to relocate/move the equipment.

**Caution**

The operator must not touch or have contact simultaneously with the patient and the laser imaging system.

**Caution**

Do not use a cell phone within 2.0 m (6.6 ft) of a laser imager. This proximity includes any imager behind a wall adjacent to your location.

**Caution**

Do not use a microwave oven within 4.0 m (13.1 ft) of a laser imager. Electromagnetic radiation from a microwave oven is only an issue if after the oven door is closed and latched, the seal does not maintain an electromagnetic tight fit between the oven door and oven main housing. Determining if the seal has an electromagnetic tight fit requires special detection equipment.

**Caution**

Do not use in the presence of flammable anesthetics, oxygen, or nitrous oxide. This equipment does not have a gas-sealed electronics enclosure and could ignite any flammable or explosive gases present in its environment.

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**Caution**

This equipment uses a DICOM network port, and is intended to connect to other medical devices over the network. It is not intended to be connected directly to other medical devices. Only qualified personnel may provide installation and service.

**Caution**

This device should not be used in close contact with MRI devices, due to possible very high magnetic fields near an MRI unit. The magnetic field in the area where this equipment is installed must be less than 50 G.

**Caution**

Do not substitute or modify any part of this equipment.

**Caution**

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the User Guide and other User Documentation, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

**Caution**

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

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**Caution**

Do not use isopropyl alcohol to clean the exterior surfaces of the laser imager.

**Caution**

In the U.S., exhausted filters are considered to be non-hazardous waste according to the US Environmental Protection Agency Resource Conservation Recovery Act (RCRA). Municipality owned and licensed solid waste management facilities are an appropriate disposal option. Contact your local or state solid waste authorities to determine if additional disposal requirements apply. In other regions, contact local or regional solid waste authorities for proper disposal guidance.

**Caution**

Lithium batteries should only be replaced by an authorized service provider. The laser imager uses a lithium battery series CR 2032 to power the clock and calendar circuitry. THERE IS A DANGER OF EXPLOSION IF THE BATTERY IS REPLACED INCORRECTLY. The battery must be replaced only with the same or equivalent type. The battery should be replaced every 6 - 8 years. The U.S. EPA's RCRA does not regulate disposal of this lithium battery. Users should discard spent batteries in municipal trash unless their community offers a battery collection program. In other regions, contact local or regional solid waste authorities for proper disposal guidance.

**Laser Warning**

The equipment uses an invisible laser beam with a maximum power of 120 milliwatts. Laser radiation may be present when the machine operates without the rear cover installed. Covers with this label may only be removed by an authorized service provider. USE OF CONTROLS OR ADJUSTMENTS, OR PERFORMANCE OF PROCEDURES OTHER THAN THOSE SPECIFIED HEREIN, MAY RESULT IN EYE DAMAGE.

Safety Labels

Safety labels are attached to the laser imager in compliance with international standards.

English Text on Labels

Some names on the labels are shortened and left in English. Below is a key to understand the meanings of the shortened words on the labels:

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1. **Class 3B invisible laser radiation.** This label states: "When open and interlocks defeated, avoid exposure to the beam."
2. **Hazard symbol**

Table 1: Laser specifications

Type	Scanning (moving) laser beam emitting from a diode
Wavelength	810 ±10 nanometers
Maximum power	120 mW
Beam divergence from Laser Diode	Minimum: 5 °, maximum: 32 °

High Voltage Warning

This warning label indicates that high voltage is present under panels or enclosures where labels are attached. These panels may only be removed by an authorized service provider.

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

Item	Label	Description
1	FCC compliance	Describes compliance, if applicable for the country of installation.

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2	Product	States that the imager is a Laser Imaging Printer.
3	Agency labels and Class 1 Laser Safety	• High voltage. Indicates that high voltage is present under panels where the label is attached. Only an authorized service provider should attempt access. • Static Sensitive Equipment. Identifies static-sensitive components. Connect a personal grounding strap to the appropriate ground before servicing this laser imager. These panels may only be removed by an authorized service provider. • Radio Frequency Energy. Indicates that the laser imager can radiate radio frequency energy. If not installed and used in accordance with the instructions, the laser imager may cause harmful interference to radio communications. • Class 1 Laser. Indicates that the laser imager complies with IEC requirements for Class 1 Laser systems.
4	Grounding reliability	States that grounding reliability can only be achieved when the equipment is connected to an equivalent receptacle marked "Hospital Only" or "Hospital Grade."
5	Japanese import license	Allows importation into Japan.

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6	Serial plate	Shows the serial number and model number of the imager along with other important data items.
7	Power cord inlet	Covers the power cord inlet when shipped from manufacturing. Shows the voltage at which the laser imager must be operated. The label is removed or moved during installation.
8	Safety sign ISO 7010-M002	Indicates that the operator must read the user documentation.

Hot and Sensitive Surface Labels

Figure 2: Drum labels: Hot Surface and Sensitive Surface (No Sharp Objects in this Vicinity)

This label indicates that you must use care where the label is installed:

- Hot surface! Take care to avoid possible burns.
- Sensitive coating! Do not use abrasive or sharp objects around the drums.

Safety and Health Compliance

This equipment has been tested for and complies with the following Safety and Emissions Standards. Certificates of Compliance and Declarations of Conformity have been issued as shown below.

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United States

• 21 CFR 1040.10 Class I

Code of Federal Regulations Title 21 Food and Drugs

Chapter I Food and Drug Administration, Department of Health and Human Services

Volume 8 - Parts 800 to 1299

Subchapter J - Radiology Health

Part 1040 - Performance Standards for Light Emitting Products

Section 10 - Laser Products

• UL 60601-1:2003 - Medical Electrical Equipment. General Requirement for Safety.

• ANSI/AAMI ES60601-1:2005 - Medical Electrical Equipment. Part 1: General requirements for basic safety and essential performance.

Canada

CSA C22.2 NO 60601-1 CAN/CSA:2008 - Medical electrical equipment - Part 1: General requirements for basic safety and essential performance.

Europe

EN 60601-1:2006 - Medical electrical equipment - Part 1: General requirements for basic safety and essential performance.

International

• IEC 60601-1:1988 - Medical Electrical Equipment. General Requirement for Safety.

• IEC 60601-1:2005 - Medical electrical equipment - Part 1: General requirements for basic safety and essential performance.

• IEC 60825-1:2007 - Safety of laser products - Part 1: Equipment classification and requirements.

EMC Standards

United States

• FCC Rules and Regulations, Title 47, Part 15, Subpart B, Class A: Radio Frequency Devices: Unintentional Radiators.

• This equipment has been tested and been found to comply with the limits for a Class A digital device pursuant to part 15 of the FCC rules. Those limits are designed to provide reasonable protection against harmful interference in a commercial or light industrial installation.

• FCC Rules and Regulations, Title 47, Part 15, Subpart C, Radio Frequency Devices:

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Code of Federal Regulations Title 21 Food and Drugs (раздел 21 Свода

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Part 1040 - Performance Standards for Light Emitting Products (часть 1040 —

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Intentional Radiators. "FCC ID: U726950"

- RF Exposure Guidance: This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter, except in accordance with FCC multi-transmitter product procedures.

Canada

- CAN/CSA-C22.2 NO. 60601-1-2-08 Medical Electric Equipment - Part 1-2: General requirements for safety - Collateral standard: Electromagnetic compatibility - Requirements and Tests.

Intentional Radiation "IC: 7027A-6950"

- This device complies with Industry Canada licence-exempt RSS standard(s). Operation is subject to the following two conditions: (1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.

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• Under Industry Canada regulations, this radio transmitter may only operate using an antenna of a type and maximum (or lesser) gain approved for the transmitter by Industry Canada. To reduce potential radio interference to other users, the antenna type and its gain should be so chosen that the equivalent isotropically radiated power (e.i.r.p.) is not more than that necessary for successful communication.

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- This Class A digital apparatus complies with Canadian ICES-003 (A).

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- This Class A digital apparatus meets all requirements of the Canadian Interference-Causing Equipment Regulations.

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Europe

EN60601-1-2: Medical Electrical Equipment - Part 1-2: General requirements for safety - Collateral standard: Electromagnetic compatibility - Requirements and Tests.

Rest of World

IEC 60601-1-2: Medical Electrical Equipment - Part 1-2: General requirements for safety - Collateral standard: Electromagnetic compatibility - Requirements and Tests.

Guidance and Manufacturer's Declaration for Electromagnetic Emissions

The system is intended for use in the electromagnetic environment specified below. The customer or user of the system should ensure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment— Guidance
RF emissions: - EN55011 - CISPR 11	Group 1	The system uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions: - EN55011 - CISPR 11	Class A	The system is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

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Harmonics emissions: • EN61000-3-2 • IEC 61000-3-2	Class A	
The system is suitable for use everywhere, including those establishments directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.		
Voltage fluctuations and flicker emissions: • EN61000-3-3 • IEC 61000-3-3	Complies	

Guidance and Manufacturer's Declaration for Electromagnetic Immunity

The system is intended for use in the electromagnetic environment specified below. The customer or user of the laser imager should ensure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment—Guidance
Electrostatic discharge (ESD): • EN61000-4-2 • IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete , or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.

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Electrical fast transient/burst:	±2 kV for power supply lines • EN61000-4-4 ±1 kV for input/output lines • IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment
Surge:	±1 kV differential mode ±2 kV common mode • EN61000-4-5 • IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment—Guidance
Voltage dips, short interruptions and voltage variations on power supply lines: • EN61000-4-11 • IEC 61000-4-11	<5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 sec.	<5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 sec.	Mains power quality should be that of a typical commercial or hospital environment. If the user of the laser imager requires continued operation during power mains interruptions, it is recommended that the laser imager be powered from an uninterruptible power supply or a battery.

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Power frequency (50/60 Hz) magnetic field: • EN61000-4-8 • IEC 61000-4-8	3 A/m	3 A/m	Mains power quality should be that of a typical commercial or hospital environment.
--	-------	-------	---

* UT is the a.c. mains voltage prior to application of the test level

Guidance and Manufacturer's Declaration for Electromagnetic Immunity

The system is intended for use in the electromagnetic environment specified below. The customer or user of the laser imager should ensure that it is used in such an environment. Portable and mobile RF communications equipment should be used no closer to any part of the system, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.

Immunity Test	IEC 60601 Test Level	Comp liance Level	Electromagnetic Environment— Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz-80 MHz	3 Vrms	$d = 1.17 \sqrt{P}$
Radiated RF C 61000-4-3	3 v/m 80 MHz-2.5 GHz	3 v/m	$d = 1.17 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.33 \sqrt{P}$ 800 MHz to 2.5 GHz where d is the recommended separation distance in meters (m) P is the maximum output rating of the transmitter in watts (W) according to the transmitter manufacturer. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey*, should be less than the compliance level in each frequency range*. Interference may occur in the vicinity of equipment marked with the following symbol:

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
Note

At 80 MHz and 800 MHz, the higher frequency range applies. These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

* Field strengths from fixed transmitters, such as base station for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the laser imager is used exceeds the applicable RF compliance level above, the laser imager should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the laser imager.

† Over the frequency range 150 kHz–80 MHz, field strengths should be less than 3 V/m.

Additional Guidance and Manufacturer's Declaration—Electromagnetic Emissions/Immunity

Electromagnetic Compatibility Precautions

Medical electrical equipment requires special precautions regarding electromagnetic compatibility (EMC). Medical equipment must be installed and put into service according to the EMC information provided in this document.

Communications Equipment

Portable and mobile radio frequency (RF) communications equipment can affect medical electrical equipment EMC performance.

Replacement of Cables, Accessories, or Transducers

The use of cables, accessories or transducers other than those specified in this document with the exception of transducers or cables sold by the manufacturer of the equipment as replacement parts for internal components, may result in increased emissions or decreased immunity of the medical equipment.

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Other Equipment

The System should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the System should be observed to verify normal operation in the configuration in which it will be used.

Shielded Locations

The System is fully compliant with the requirements of IEC 60601-1-2 without being located in a shielded room.

External Cables Listing

Description	Length	Shielded?
Cable, Category 5E, Shielded RJ45 to RJ45, Patch Cord, 4 Pair	2.1 m (7.0 ft)	Yes
Cord Assembly, Power, China	2.5 m (8.2 ft)	No
Cord Assembly, Power, Continental Europe Schuko	2.5 m (8.2 ft)	No
Cord Assembly, Power, India/South Africa	2.5 m (8.2 ft)	No
Cord Assembly, Power, North American	3.0 m (10.0 ft)	No
Cord Assembly, Power, UK, Fused	2.5 m (8.2 ft)	No

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Recommended Separation Distance Between Portable and Mobile RF Communications Equipment and the System

The system is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the system can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communication equipment (transmitters) and the system as recommended below, according to the maximum output of the communications equipment.

Rated Maximum Output Power of Transmitter (P)	Separation Distance (d) According to Frequency of Transmitter		
	Meters		
Watts	150 kHz-80 MHz d = 1.17 √P	80 MHz-800 MHz d = 1.17 √P	800 MHz- 2.5GHz d = 2.33 √P
0.01	0.12	0.12	0.24
0.10	0.37	0.37	0.74
1.00	1.17	1.17	2.33
10.00	3.70	3.70	7.37
100.00	11.70	11.70	23.30

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note
At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

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EU Directives

- 93/42/EEC Title: Council Directive Concerning Medical Devices.
- 1999/5/CE Title: Council Directive Concerning Radio Equipment and Telecommunications Terminal Equipment.

Authorized European Representative

Carestream Health France
1, rue Galilee
93192 NOISY-LE-GRAND CEDEX
FRANCE

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APPLICATION: Overview


Figure 1: Imager features

- 1 **Touch screen.** Your interface to the imager, the touch screen provides basic functionality, configuration/setup features, and status information.
- 2 **Power switch.** Press the power switch to turn the imager on or off. Check that **Ready** displays on the touch screen before switching the power off

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- 3 Right cover.** Protects sensitive electronic equipment. The right cover is only accessed by service personnel.
- 4 Front door.** Covers the three film supplies and two filters. Check that **Ready** displays on the touch screen before opening the door.
- 5 Filter area.** The two filters are at the bottom of the imager, inside the front door.
- 6 Left cover.** Covers the left side of the imager. You might remove the left cover to clear an occasional film jam. The left cover is interlocked.
- 7 Top cover or film sorter (paper tray).** If the imager does not include the optional sorter, the processed film exits to the top cover. The film sorter receives the processed film and transports the film to the selected film bin on top of the imager. Each of the five sorter bins can hold 50 sheets of film.

Important

To raise the sorter to access the top cover of the imager, remove any films that are resting on a sorter tray, and then grip either sorter handle with both hands on one side. Gently rotate up until the sorter is at an approximate right angle to the imager top. With the sorter fully raised, you can remove the top cover if needed.

Caution

To return the sorter to the operating position, grip one handle with both hands, and gently rotate the sorter down to rest on the imager top cover. Do not let the sorter fall down freely

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
- 8 **Siesmic Kit** (steel anchor bolts D 8 mm, L 80 mm - 4 pcs., Mounting brackets), via which the camera is attached to the floor surface in areas of seismic activity

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
- 1 **Processor area and drums.** Provides the heat that processes the image on the film.
 - 2 **Exposure transport.** Moves the film past the scanning laser beam.
 - 3 **Pickup assembly (upper, middle, and lower).** Lifts a single sheet of film from the supply cartridge and feeds it into the rollers.
 - 4 **Rollback assembly (upper, middle and lower).** Rolls the film cartridge cover back so the pickup assembly can lift the film. When the imager is not printing, the cover is closed over the film cartridge to protect the film from light.
 - 5 **Vertical transport.** Orients the film while moving it (the film) from the supply to the imaging portion of the imager.
-
- 6 **Accumulator.** Receives the film before it is imaged. The film is sent from the accumulator up to the optics module, and then to the processor assembly where the heat is applied to process the image.
 - 7 **Charcoal filters.** Absorb potential odors that may be caused by heat processing.
 - 8 **Feed rollers.** Move the film through the imager.
 - 9 **Film supply (upper, middle, and lower).** Holds the film cartridges. You load the film into the film supplies from the front of the imager (the front door).
 - 10 **Optics module.** Writes the image onto the film while the film is moved through the exposure transport.
 - 11 **Densitometer.** The densitometer receives the film after it passes through the processor and performs a density check to confirm image quality when a calibration is run.

How the Imager Works

The imager receives, processes, manages, and prints the images on film. As images are received for printing, they are placed in a sequential print queue and printed in order. The imager can continue to accept incoming print jobs even if temporarily unable to print (if the film supply is empty, etc.).

Print Sequence

Each time the imager receives a print request, it determines the requested film size and type, and then selects the appropriate film cartridge. The following print sequence occurs. The arrows show the film path.

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

Figure 3: Print sequence and film path

1. Suction cups in the pickup areas lift a single sheet of film out of a cartridge and feed the film into the transport rollers.
2. The transport rollers move the film down into the accumulator area.
3. The film reverses direction in the accumulator and moves up during imaging (as the optics module writes the image onto film). Then the film moves up to the processor.
4. As the film passes over the two processor drums, the heat generated by the drums develops the film.
5. The exit rollers move the developed film through the cooling section and the turnaround to the sorter (if applicable) or to the top cover.

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Supported Film Sizes

The imager is available in different configurations; the following film sizes can be supported:

- 35 x 43 cm (14 x 17 in.)
- 35 x 35 cm (14 x 14 in.)
- 28 x 35 cm (11 x 14 in.)
- 25 x 30 cm (10 x 12 in.)
- 20 x 25 cm (8 x 10 in.)

Clear, blue, and mammography film types may be available.

Mammography Imaging

To manage and support mammography imaging, the imager provides these features.

- You can load mammography film in any of the film supplies.
- The imager can be configured to print mammography images exclusively on mammography film. If a mammography image request is received, and:
 - a cartridge with mammography film is not loaded, and/or
 - the requested mammography film size is not loaded, then the imager will request the appropriate film size and type. You must load the requested size and type to print the mammography job.
- Several mammography-specific test prints are provided, including a test print for the routine calibration of mammography film.

Automatic Image Quality and Processing

An internal densitometer enables the imager to automatically adjust image processing parameters using Automatic Image Quality Control (AIQC) to produce an optimal image. The imager adjusts these parameters each time it prints a calibration film.

A calibration film (a series of gray bands of increasing density) is printed when:

- A film cartridge is inserted in the imager with film of a new manufacturing lot number.
- You request a calibration film at the touch screen or the Web Portal.
- A film cartridge is inserted into the imager for which a current calibration is not stored.

When the imager is calibrating the film, the calibration symbol appears on the film supply icon.

If the results of a film calibration are slightly outside of the target range, the cartridge enters AIQC manual mode, and the symbol (below) appears on the film supply icon. You can still print with this film cartridge, but you should verify the image quality on the printed film.

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
Figure 4: Film cartridge in manual mode

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Monitor the Imaging System from a Remote Location

You can view and manage the imager features and connections directly at the imager OR from across the room, across the medical center, or across town. The Web Portal enables you to view the imager features, connections, status, etc., from a networked personal computer. You can also check film count and the loaded film size and type.

Related topics:

Access the Web Portal

Enhanced Serviceability with Remote Monitoring

Smart Link is the remote technology that acts as a bridge between your imager and the service Customer Success Network. Smart Link technology, with Remote Management Services, is designed to enhance efficient system serviceability and simplify analytical and service processes through the monitoring of your equipment. Additionally, this simplifies the service process by providing qualified service personnel with faster, easier, and more complete access to the operational history of each system.

Remote Management Services provides the following for the imager:

- Ability to monitor and diagnose error conditions without introducing “downtime”
- Firewall-safe, Internet transmission of machine data, while complying with patient confidentiality regulations
- Minimal effort to set up communication
- Flexibility to facilitate specific configurations at each site

Installation, Setup

Imager installation and setup must be performed by a qualified service provider. Contact a qualified service provider with any questions.

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Basic Operating Tasks

During normal operation, the imager receives and automatically prints images sent by modalities over a network. Very little interaction is required. You can do the following:

- Turn the power on (I) and off (O)
- Load the film cartridges
- Monitor the touch screen for status and operating condition, and respond when appropriate

Sometimes it will be necessary to turn off an alarm, perform preventive maintenance, replace the filters, and perform other corrective actions such as a restart (as indicated on the touch screen).

You also might adjust image quality, use the mammography test prints, or do other troubleshooting tasks.

You can do all of these tasks at the imager touch screen, and some of the tasks can be done from a remote personal computer using the Web Portal.

Related topics:

Access the Web Portal

Imager Power

Power Switch Location

The power switch is located on the right side of the imager, near the front and top.

Start the Imager

Press the power switch to On (I). Wait as the imager warms up. The touch screen shows the progress as the imager becomes ready to print.

The warm-up and self-test period might last up to 35 minutes, depending on how long the imager has been powered off, the imager voltage, the ambient temperature, and whether the imager is in Power Save mode.

When the imager reaches operating temperature, the touch screen changes to show that the system is ready to print, and the imager prints any images that were received during the warm-up period.

During warm-up, you can use the touch screen to check the contents of the film supplies and open the front door to load new film cartridges if necessary.

Emergency Shutdown or Power Loss

In the event of a power loss, or if an emergency requires an immediate shutdown of the imager, films in process will not be completed. However, when power is restored, the imager will restart. After warming up, the imager automatically:

- ejects any jobs that were in process when the power was interrupted (the imager might eject blank films that had already left the film cartridge)

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- reprints those jobs that were not completed before the power was interrupted

Shut down or Restart the Imager

If you can, wait until any films in progress have finishing printing (the imager status, in the upper left corner of the touch screen, states *Printing* when prints are in progress).

Then press the power switch to Off (O). The imager will shut down.

To restart the imager, press the power switch to On (I). Wait as the imager warms up.

After warming up, the imager automatically ejects any partially-printed films (blank films might be ejected if they had already left the film cartridge), and reprints any films that were in process when power was interrupted.

Power Save

The imager has a Power Save option, which lets you reduce the power consumption of the imager by turning off power, fans, and motors during periods of inactivity. You must schedule the Power Save mode. To learn more, touch **Utilities > Power Save**, and then touch the ? button to review the online Help.

Related topics:

Touch Screen Help

Touch Screen Overview

Use the touch screen to monitor imager operations. For example, check that the imager is ready to print or check how much film is remaining in each cartridge.

The touch screen also displays information on how to resolve any errors that might occur.

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

Figure 5: Example touch screen

Item Feature	Item Feature Description	Item Feature Description
1	Imager status	Reports the status. Example messages: Ready, Self Test, Warming, Printing Disabled, Not Ready If the imager is not ready to print, check that the film cartridges are properly inserted, all doors/covers are closed, and check for any errors on the touch screen. If the imager has been off, it must warm up before it will be ready to print. The approximate number of minutes until the imager is fully warmed up displays in the imager status location on the touch screen.
2	Print jobs	Reports the number of print requests in the various queues.

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**Important**

To avoid exposing the film to light, do not open the imager front door until the **Pause** (||) symbol is off. Wait until the symbol is off to replace a cartridge.

Figure 6: Pause symbol illuminated (on)

- During most normal operations, including when the imager is idle, the symbol is off.
- The symbol is on while the imager is processing images, test prints, and during calibration.

If the **Pause** symbol is on, you can press the symbol to temporarily pause printing. When you press it, the symbol changes from blue to orange as jobs in progress finish printing. When the **Pause** symbol turns off, you can open the front door.

**Note**

After you press the **Pause** symbol, printing is suspended for approximately 5 minutes, or until you open and close the front door.

Displays the film count, type, and size. When film count drops below 8, the icon changes from blue to yellow. When film is depleted, the icon reports **Empty**.

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5	Help	Displays general online Help.
6	View Jobs	Displays the print jobs: queued to print, waiting for film, queued for retry, and unprintable. View the lists of jobs in each category, jobs by modality, etc.
7	Utilities	Displays the options for user operations, system information, configuration, and service features.
8	Messages	When the imager has messages to report (such as error messages), they appear in this area. When a message appears, touch it to learn more.

Clean the Touch Screen

Use a soft cloth to wipe the touch screen free of fingerprints or smudges. If needed, you can use a standard screen or glass cleaner.

Touch Screen Functions

Use the following table to navigate the touch screen. All functions begin at the **Main Screen**.

How To	Touch
Add or configure a modality	Utilities > Configuration > Configure Modality
Adjust volume	Utilities > User Operations > Speaker Volume
Assign sorter bins	Utilities > User Operations > Sorter Bin
Back up the system	Utilities > Configuration > Backup
Calibrate film	Film supply icon > Calibrate Film Supply

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Calibrate touch screen	Utilities > User Operations > Calibrate Touch
Change between metric and US (imperial)	Utilities > Configuration > Film Display Units
Change date/time	Utilities > Configuration > Date & Time
Change filters	Utilities > User Operations > Change Filter
Change the language	Utilities > Configuration > Change Language
Change the service phone number	Utilities > Service > Service Number
Change the brightness level of the touch screen	Utilities > User Operations > Display Brightness
Check values and print a QC mammography wedge test	Utilities > User Operations > Mammo QC
Configure the imager	Utilities > Configuration > Installation Assistant
Configure the imager on the network	Utilities > Configuration > Network
Disable printing	Utilities > Service > Printing Enabled
Enable printing	Utilities > Service > Printing Disabled
Log a filter change	Utilities > User Operations > Change Filter > Log Change Or Utilities > Service > Preventive Maintenance

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Log that the preventive maintenance has been completed	Utilities > Service > Preventive Maintenance
How To	Touch
Ping	Utilities > Service > Ping
Reset or modify modality settings	Utilities > Configuration > Reset Modality
Restore from backup file	Utilities > Configuration > Restore
View on the reverse side of the film	Utilities > Configuration > Reverse View
Run mammo test patterns	Film supply icon > Mammo Test Patterns
Run test patterns	Film supply icon > Test Patterns
Schedule power save	Utilities > User Operations > Power Save Schedules
Start power save	Utilities > User Operations > Power Save
Stop or start Remote Management Services	Utilities > Configuration > Remote Management
Switch the imager offline	Utilities > User Operations > Imager Online
Update software	Utilities > Service > Update Software
View DICOM associations	Utilities > View System Information > Diagnostic Data
View error history	Utilities > View System Information > Error History

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
View imager system information	Utilities > View System Information > System Info
View job lists	View Jobs
View network information	Utilities > View System Info > System Info > Network
View print counts	Utilities > View System Info > Print Counts
View print counts by modality	Utilities > View System Info > Modality Counts
View prints to filter change due	Utilities > View System Info > Print Counts
View the number of prints until the preventive maintenance is due	Utilities > View System Info > Print Counts
View software versions	Utilities > View System Info > System info > Versions

Touch Screen Help

Information about imager functions is provided through the touch screen Help. Help screens are available in two forms: general Help and screen-specific Help.

The **Main Screen** and most other screens contain a ? (**Help**) button located in the upper-right corner of the screen.

General Help

General Help provides basic information about the imager and its functions.


To access the general Help, touch the ? button on the **Main Screen**. On the screen that displays, touch a topic heading, and then touch the specific topic.

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Screen-specific Help


Screen-specific Help provides information and/or instructions concerning the screen on which the ? button was touched.

To access screen-specific help, touch the ? button on the current screen.

Work with Film Cartridges

The film supply icons on the touch screen correlate with the position in the imager (upper, middle, and lower).

The icons display the size of the film loaded, the number of film sheets remaining, and in some cases, status.

Example Film Supply Icon	Description
	The film count is 125, which indicates a full cartridge of 8 x 10 in. film. Ready to print from this film supply.
	Six sheets of film 8 x 10 in. film remaining. The yellow color indicates that film count is low.
	This cartridge with 68 sheets of 8 x 10 in. film has been calibrated, but the film is out of the calibration range. The icon indicates that you are printing in Automatic Image Quality Control (AIQC) manual mode. You can still print with this film cartridge, but you should verify the image quality on the printed film.

Important

To avoid exposing the film to light, do not open the imager front door until the **Pause** (||) symbol is off (**not** illuminated). Wait until the symbol is off to replace a cartridge.

Figure 7: Pause symbol illuminated (on)

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
If the **Pause** symbol is on, you can press the symbol to temporarily pause printing. Any jobs in progress finish printing, then the film cartridge cover closes. When the **Pause** symbol turns off, you can open the front door.

Load and Unload Film Cartridges

The film supply icons on the touch screen display the number of film sheets remaining. As a film cartridge runs low on film (fewer than 8 sheets remaining), the icon changes from blue to yellow. When a film cartridge is empty, the film supply icon displays **Empty**.


Note

Make sure the **Pause** (||) symbol is **off** before you open the imager front door. The following shows an illuminated **Pause** (||) symbol:

1. Open the imager front door.

H245_0005

2. To remove a film cartridge, hold the edges of the cartridge and lift the cartridge out of the film supply.
3. Discard the empty film cartridge.

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

4. To load a new or different cartridge, insert the cartridge into the film supply with the cartridge label facing up and the perforations facing in toward the imager.

5. Close the front door.


6. Check that the touch screen indicates the new film count. A new film cartridge contains 125 sheets.

Calibrate the Imager for the Loaded Film


The imager automatically calibrates film with a new media lot number. Occasionally, it may be necessary to calibrate the imager manually. In that case, the film supply icon reports the need to calibrate:

1. Touch the film supply icon.
2. Touch **Calibrate Film Supply**. Touch **Yes** to continue.

When the calibration is complete, the film supply icon displays the print count again. If the results of a film calibration are slightly outside of the target range, the cartridge enters Automatic Image Quality Control (AIQC) manual mode, and this symbol appears on the film supply icon:

You can still print with this film cartridge, but you should verify the image quality on the printed film. If the images are not acceptable, insert a different film cartridge.

Make a Test Print (SMPTE or Mammography)

You can initiate a test print to verify imager operation or mammography image quality.

1. Touch the film supply icon for the supply where you want to make a test print (upper, middle, or lower film supply):

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

2. Touch your choice:

» **Test Patterns**—Density Test (SMPTE), TG18-PQC

» **Mammo Test Patterns**—Mid-Density FlatField, TG18-PQC, TG18-QC

Note

Mammo Test Patterns are only available when mammography film is loaded in the film supply.

3. Confirm that the film supply and film size are correct, and then:

» Increase or decrease density

» Increase or decrease number of copies

4. Touch **Print**.

Note

Requirements for mammography quality control vary by region.

Manage Print Jobs

Most of the time, it is not necessary to manage the print jobs. Intervention may be needed in less common circumstances.

If a message appears on the main screen, for example, *2 jobs waiting for media* or *Queued print jobs need the following media*, touch the message to learn how to proceed and/or to silence the alarm. Press X to close the message.

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
You might view the print jobs at the touch screen when:


- you want to delete a print job
- a job cannot print, for example, because the requested film size is not available
- you want to view the queued print jobs, to see where your job is ordered in the queue

Print jobs are arranged in queues:

- **Queued to print**—Provides a list of the jobs waiting to be printed.
- **Waiting for media**—A different size or type of media is needed for these jobs.
- **Queued for retry**—These jobs did not print successfully on the previous attempt. The imager will automatically try again to print them.
- **Unprintable jobs**—These jobs cannot be printed on this imager. (Normally, unprintable jobs are automatically deleted from the imager.)

To view the print queues:

1. From the **Main Screen**, touch **View Jobs**.
2. Touch a list to view the jobs. The number of jobs for each list displays. Select a modality to view, or view the resulting list for all modalities.

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
3. Delete the job or all jobs as needed.

Filter Maintenance

The imager has two charcoal filters that must be maintained as part of preventive maintenance. The filters reduce potential odors from the imager during film processing.

Caution

In the U.S., exhausted filters are considered to be non-hazardous waste according to the U.S. Environmental Protection Agency Resource Conservation Recovery Act (RCRA). Municipally owned and licensed solid waste management facilities are an appropriate disposal option. Contact your local or state solid waste authorities to determine if additional disposal requirements apply. In other regions, contact local or regional solid waste authorities for proper disposal guidance.


The filters must be replaced after 15,000 prints. Keep two new filters in inventory to replace the used filters when needed.

Replace the Filters

The filters must be replaced after 15,000 prints.

Prerequisites:

The touch screen displays **Filter Change Due**.

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
Note

Make sure the Pause (||) symbol is **off** before you open the imager front door. The following shows an illuminated **Pause (||)** symbol:

1. Open the front door.

H245_0006

2. Remove the filter by lifting it up and pulling it forward.
3. The second filter is behind the first. Reach in and pull the second filter forward and out of the imager.
4. Install the new filters.
5. Close the front door.

Postrequisites:

Go to **Utilities > User Operations > Change Filter** and touch **Log Filter Change**. You must log the change to clear the **Filter Change Due** message and to reset the imager counter (so it can properly alert you next time the filters must be changed).

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
Access the Web Portal

Prerequisites:

A personal computer (desktop or laptop) connected to the network

You can use the Web Portal to monitor the imager from a personal computer. The PC can be in any location. It is not necessary to use the Web Portal, but you may find it convenient.

1. On a desktop or laptop computer, start WINDOWS INTERNET EXPLORER (IE).

Note

Only IE 6, 7, and 8 have been qualified with the imager. See the note below if using IE 8 or 9.

2. In the address field, type **http://<imager's IP address>**

Note

You can obtain the IP address from the imager touch screen. Touch **Utilities > Configuration > Network**.

3. Click **Go**.

The Web Portal opens to the Device Status (Home) screen.

- The main window shows general status, the number of print jobs queued, the number of jobs waiting for film, and the film count for each cartridge.
- The center panel displays the screens where you view and perform tasks. Online Help is available by selecting **Documentation** from the left panel, and the Help icon provides context-sensitive Help for fields and pages.
- The left panel also provides the links to all other screens.

Note

If you are using WINDOWS INTERNET EXPLORER 8 or 9, place the INTERNET EXPLORER (IE) window into compatibility view. After you have opened the IE window, click the Compatibility View toolbar button. This will correct some potential viewing issues. If the icon is not on the toolbar, select **Compatibility View** from the **Tools** menu.

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Troubleshooting Errors with the Imager

In the event of any failure or error in the printer, the display will show the headline "ALERT" or "FAULT" and a message explaining the nature of the error

Note

Some errors result in an audible alarm, alerting you that the imager needs attention. This is helpful if you do not work next to the imager. However, you can change the volume of the alarm or mute it if the alarm is unnecessary for you. Go to **Utilities > Speaker Volume** to change the setting.


Error Correction

When an error occurs, you are typically alerted with a message and alarm.

When an error occurs:

- An alarm sounds
- The correction procedure appears on the touch screen

Figure 8: Error description and resolution appear on the touch screen

1. If an alarm is sounding, stop the alarm at the touch screen. Touch **Silence Alarm**.
2. Correct the error:
 - » View the error correction procedure on the touch screen. Follow the on-screen

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directions.

» If a **Current Errors** button displays, touch it to view the error correction procedure.

If there are multiple errors, use the arrow buttons to scroll to view the error you want. Touch to select it.

3. When finished, touch **Close** to return to the **Main Screen**.

Look Up an Error to Resolve It

If you have the 5-digit error code, you can view the error correction steps from the online Help.

1. To access the online Help from the **Main Screen**, touch the ? button.

2. From the main menu, touch **Error Resolution > Look Up an Error to Resolve It > Error Codes Listing**.

A list of all error codes and resolutions displays.

3. Scroll to locate and view the error correction procedure on the touch screen. Follow the on-screen directions.

Labeled Areas to Check for Film Jams

Under the left panel, the imager has numbered labels (1–4) to indicate the location to check for film jams. These areas are referenced in the jam resolution instructions.

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

H245_0007

Label	Imager Area
1	Film supplies
2	Lower film path
3	Upper film path/processor drums
4	Cooling section/densitometer/exit area

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Manually Remove Film or Close a Cartridge

For some film jams, you can remove the film by turning a knob to move the film out of the imager. It also might be necessary to raise the cooling cover or the processor drum cover to remove jammed film underneath.

Caution


Always:

- make sure the **Pause** (||) symbol is **off** before opening any imager covers.
- use the blue-colored knobs and handles to safely open or lift covers on the imager.

1. Raise the sorter (if installed). Remove any films that are resting on a sorter tray, and grip a handle on either side of the sorter with both hands. Rotate the sorter up until it is at an approximate right angle to the imager top.

2. Remove the imager left cover:

- a. If the **Pause** (||) symbol is illuminated, press it and wait until it goes off.
- b. Open the front door and turn the knurled knob on the inside left.
- c. Remove the top cover and turn the knurled knob on the top left.
- d. Lift the left cover from the imager.

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

3. Following the error resolution instructions on the touch screen, locate the jammed film and manually remove it as needed by turning the appropriate knobs:


Knobs

Description

1	Processor drums	Turn to remove film from the drum area.
---	-----------------	---

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
2	Exposure transport	Turn to remove film from the exposure transport area.
3	Rollback knob	If the film becomes misplaced around the film cartridge, you can turn the knob to manually open or cover the film cartridge.
4	Vertical transport	Turn to move film up through the vertical transport.
5	Cooling area	Turn to move the film through the cooling area.
6	Exit roller	Turn to move the film through the turnaround area and exit area.

4. To check for film under the cooling cover (at the top of the imager), unscrew the knurled knobs and rotate the cooling cover up.



Hot Surface: The processor drums and rollers are hot. Take care when opening the processor drum cover.

5. To check for film under the processor drum cover, disengage the cover: push back on the blue lever on the left side of the cover, and grip the blue handle to rotate the cover up.

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
6. When the film has been located and removed, close the imager covers.

Touch Screen is Not Functional


If the touch screen is not responding, use the power switch on the imager to turn power off, and then on. If the touch screen is still non-responsive, contact a qualified service provider.

Information for Support

If you cannot correct a condition and need help, call a qualified service provider. To locate the phone number, touch **Utilities > Service > Service Number**.

Have the following information ready when you call:


- Model number and/or serial number
- Imager ID/K#
- Any error code and error message that is displayed on the touch screen

Note

You can locate the imager identifying information on the serial plate on the back of the imager, or go to **Utilities > System Information** on the touch screen.

Have access to the imager when you call; the Support Center personnel may ask you to perform some troubleshooting techniques to resolve the issue.

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About Configuring Modalities or Optimizing Printed Images

The imager provides a set of image processing parameters that is likely to produce acceptable images from each image source (modality).

From the touch screen, you can accept the system recommendations, refine the parameters by providing more information about the modality, or manually set up different parameters to meet individual or site preferences.

You can optimize the printed images using multiple approaches:

1. First Print—From a new modality, send an image to print at the imager. The imager recognizes that the modality is a new image source, attempts to identify the modality, and offers a set of image processing parameters based on what it can tell about the modality. Try the recommended parameters first; print an image and check if the image look is optimal. Follow the instructions on the touch screen as you go through this process.

2. Reset Modality—To help the imager apply the best parameters for the modality, you may select or confirm the modality type, manufacturer, and model. Touch **Utilities > Configuration > Reset Modality**, and follow the instructions on the screen.

Note

This approach is beneficial if the images are not sufficiently optimized after the **First Print** process. The **Reset Modality** process lets you refine the imager's understanding of the modality. For example, the imager might partially recognize a modality by type and manufacturer but not the model. In that case, you can specify the model, which will allow the imager to refine the parameters it applies to the images from that modality.

As instructed on the touch screen, send another image to print after resetting the parameters. If the image is still not optimized, you can change individual parameters (see step 3), or contact a qualified service provider for help.

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Important

Manually altering individual parameters is not recommended unless you have in-depth knowledge of image processing. Contact a qualified service provider for help as needed.

3. Configure Modality—Manually adjust individual processing parameters such as density, contrast, and smoothing. Specify if the settings at the imager should override processing parameters set up at the modality.

Touch Utilities > Configuration > Configure Modality, and alter parameters as needed.

Send another image to print after altering any parameters. If the image is still not optimized, contact a qualified service provider for help.

Film Technical Information

This section describes the characteristics of Laser Imaging Film, not the operation of the imager. The Laser Imaging Film is a high-resolution, infrared-sensitive, photothermographic film designed specifically for the imager.

Spectral Sensitivity of the Film

The Laser Imaging Film is infrared sensitive and has been sensitized to the infrared laser diode of the imagers. When handled according to instructions on the daylight-load film package, safelights are not needed. If you remove undeveloped film from the daylight-load package, you will need a darkroom setting and a green safelight.

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
Film Image Quality

The Laser Imaging Film delivers diagnostic-quality, continuous-tone images along with sharp alphanumeric and optimum contrast. This high-quality, silver-based film provides health care providers with the same diagnostic information they are accustomed to viewing—including spatial resolution, contrast, and gray levels. Because it is a totally dry imaging process, there is no image quality variability due to wet chemistry.

Undeveloped Film Handling and Storage


To achieve consistent results up to the expiration date indicated on the film package, the Laser Imaging Film must be stored in a cool, dry place (4–24 °C, or 39–75 °F) and properly shielded from x-rays, gamma rays, or other penetrating radiation.

The film can withstand short-term temperature spikes (up to 35 °C, or 95 °F) for several hours without any significant effect on film quality or performance. Temperatures above 35 °C (95 °F) will gradually diminish shelf life.

Hands must be clean, dry, and free of lotions. Handle film carefully by the edges to avoid physical strains such as pressure, creasing, or buckling.

Developed Film Handling and Archival

Handling the Laser Imaging Film requires reasonable care. Spills, humidity, and other moisture typically have no significant effect on developed films. However, prolonged

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exposure to intense light or excessive heat (equal to or greater than 54.4 °C or 130.0 °F) for more than three hours may cause some gradual darkening of images. Leaving films in vehicles in hot climates for extended periods of time is not recommended.

For best results, store the film in sleeves when not being reviewed. The Laser Imaging Film can be left on a light box for more than 24 hours. In extreme cases in which the light boxes are exceptionally hot (equal to or greater than 49 °C or 120 °F), the manufacturer recommends removing them prior to eight hours of continuous exposure.

Take care when using spotlight viewing for more than 30 seconds because temperatures near the light source may exceed 82.2 °C (180.0 °F). Use in slide projectors is not recommended due to the high temperatures generally found in these devices.

With dry technology, a small amount of final development occurs when the film exits the imager and is initially exposed to ambient or view-box lighting. This is virtually undetectable and has no effect on image quality (typically 0.02 change in density). This small density increase is uniform and permanent upon full exposure of the film under normal handling conditions (room light or view box).

The Laser Imaging Film has been tested and can be archived for more than 100 years when stored at American National Standards Institute (ANSI) recommended storage conditions (equal to or less than 25 °C or 77 °F). Processing film should be stored within the temperature range of 16–27 °C, or 60–80 °F, and at 30–50 % RH. Developed films may be stored at higher temperatures; however, that will reduce the number of years the film can be stored. For example, storing films at a constant elevated temperature of 32.2 °C (90.0 °F) may reduce archive capability to 30 years.

Exposing Film to Moisture

The Laser Imaging Films typically withstand humidity, spills, and other forms of water without any significant effect on image quality or film integrity. If needed, film can be cleaned with a clean, damp cloth.

Odor Dissipation

Dry technology eliminates nearly all processing odors. While some low-level odors are produced during the development process, they pose no known adverse health risks. Processing odor levels are further reduced by non-hazardous, disposable filters in the imager. The filters trap most low-level odors and prevent them from dissipating into the work environment. To help maintain optimum performance, the filters require periodic replacement. The imager requires no special venting.

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Heat Dissipation

The imager uses controlled heat to develop the Laser Imaging Film. The heat has virtually no effect on the air temperature of the work area. The amount of heat dissipated into an area during a day is typically less than the heat generated by two to four 100-watt light bulbs.

Specifications

Equipment Specifications

Unpacked	
Height	115.0 cm (45.3 in.) without sorter
	159.0 cm (62.6 in.) with sorter
Width	62.6 cm (24.7 in.)
Depth	68.6 cm (27.0 in.) without sorter
	71.8 cm (28.3 in.) with sorter
Weight	147 kg (325 lb) without sorter
	177 kg (390 lb) with sorter
 Note The approximate weight of a 35 x 43 cm (14 x 17 in.) film cartridge is 3.5 kg (7.7 lb).	

Note

The unpacked height is for the imager, out of the crate, with the film exit trays installed (the maximum height of the imager). The packed height is for the imager in the crate without trays installed.

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Power supply - 200-240V
Power (max) - 6A
Mains frequency - 50-60 Hz
Safety class - I

Imager Placement Requirements

Follow the placement requirements for the imager location.

- Allow 61 cm (24 in.) clearance around the top, sides, and back of the imager. This space is required for service access.
- The floor must be level (within 1°) and capable of supporting the imager's weight.
- Place the imager in an area with good ventilation. A small, confined room is not recommended.
- Avoid placing the imager in direct or excessive sunlight (for example, near a large window).

Environmental Requirements

Follow the environmental requirements for the imager location.

Temperature

- Operating: 15 to 33 °C (59 to 91 °F)
- Storage: -35 to 60 °C (-31 to 140 °F)

Relative Humidity

- Operating: 20–80 % RH, non-condensing
- Storage: 10–90 % RH, non-condensing

Altitude

-30 m (-100 ft) to 3,000 m (9,800 ft) above sea level

Environmental Effects

Acoustical noise:

- Less than or equal to 55 dB-A at 1 m from the front of the imager during idle or standby
- Less than or equal to 75 dB-A momentary at 1 m during normal operation

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- Less than or equal to 65 dB-A at 1 m during full printing cycle

Power Requirements

- 90–132 V (ac); 50/60 Hz; maximum 12 amps
- 180–264 V (ac); 50/60 Hz; maximum 6 amps

Connect the equipment to a power source that is suitable for the voltage and current ratings shown on the rating label. The single-phase power source, with grounding, must be provided within 2.5 m (8.0 ft) of the imager.

The wire must be insulation-rated for 600 V (ac). A dedicated line is recommended.

Power Cords

A power cord is provided with this equipment. All countries must use an agency-approved power cord with a plug type suitable for the country of use. Contact a qualified dealer for help.

Network Requirements

The imager receives digital images from medical imaging devices (modalities) over a 10/100/1000 Base-T Ethernet network.
Network Connection - Built-in DICOM server.

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USER MANUAL CONVENTIONS

Note

Notes provide additional information, such as expanded explanations, hints, or reminders.

Important

Important notices highlight critical policy information that affects how you use this manual and this product.

Caution

Cautions point out procedures that you must follow precisely to avoid injury to yourself or others, damage to the system or any of its components, or loss of data or corruption of files in software applications.

DANGER

Danger identifies procedures that you must follow precisely to avoid injury to yourself or others.

Laser Warning

Laser warnings warn personnel that access to laser radiation is possible and all personnel must avoid direct exposure to the beam.

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THE MATERIALS FROM WHICH IT IS MADE OR ITS MEDICAL DEVICE BASIC PARTS

Product name	Material
case	Polymer polycarbonate 104-701, polycarbonate dye dark gray 1-5%
display	Polyethylene terephthalate CAS No. 25038-59-9

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Film Recycling

According to the U.S. Environmental Protection Agency (EPA) standards, the Laser Imaging Film is not considered hazardous and requires no special disposal procedures. However, the film does contain silver and polyester that may be recovered by using one of several recycling processes.

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CUSTOMER SERVICE SUPPORT

LLC " Carestream Health " Moscow, street 4th Backbone, b. 11
Tel. : 8 (495) 660-56-90
Fax: 8 (495) 660-56-91

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CHARACTERISTICS

Parametr / element	Value / characteristic
Capture unit, pcs.	3 (upper, middle, lower)
Number of bins	5
Quality laser printing	650 pixels per inch The pixel size of 39 microns
Grayscale	14- bit
Productivity per hour films:	
Film size 35 x 43 cm	before 160
Film size 20 x 25 cm	before 250
prepress time, s	60
Film for laser exposure	
• Lifetime (100+) years film archivability	
Pictures with the standard optical density, Dmax	3.3 - blue substrate 3.2 - a transparent substrate

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Robert C. Meagher
Senior Director
World Wide Regulatory Affairs
Carestream Health, Inc.
150 Verona Street
Rochester, New York 14608
Telephone 585-627-6528
STATE OF NEW YORK
COUNTY OF MONROE
Signed and attested before me on (August 26, 2016) by

CHRISTINE A. BENOIT-AVERILL
NOTARY PUBLIC-STATE OF NEW YORK
No. 01BE6320594
Qualified in Monroe County
My Commission Expires March 09, 2019

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