Additional Operational Documentation

Hemodialyzer Leocceed, modifications Leocceed-18N, Leocceed-18H

(Supplement for the Russian Federation only)

The authorized representative of the manufacturer in the Russian Federation:

Delrus Inc.

Russia, 620086, Sverdlovsk region, Ekaterinburg, Posadskaya str., 23

Content of supply

Each carton contains the following contents:


- Thirty (30) each dialyzers
- Five (5) each fiber packaging trays (to cushion and protect the product)
- One (1) Instructions for Use document (operational documentation)

Weight of the carton is 4.1 kg

Construction of dialyzer

The Leocceed-18H and Leocceed-18N Dialyzers utilize the components and materials identified in item 5 Specification and Performance (IFU). The housing is an extruded tube that is cut to length for the dialyzer housing. The blood compartment caps, dialysate compartment caps, and port covers are all injection molded components. The fiber bundles are constructed of spun fiber made from Polysulfone that are bound in bundles of 13,000 – 13,850 fibers per bundle.

The Dialyzers are assembled using automated equipment into their finished form. The Dialyzers are sterilized using gamma irradiation to a Sterility Assurance Level (SAL) of 10^{-6} and are validated to have a sterile, non-pyrogenic fluid pathway. The sterile barriers for the product are the four port caps: two blood port caps and two dialysate port caps.

Shelf-life

The Leocceed Dialyzers are validated for a labeled shelf life of 36 months from the date of manufacture.

Disposal conditions

Dispose the used Leoced and the blood lines in strict accordance with official and institutional standards for medical waste disposal.

The substance hazard category of dialyzers according to Russia relevant regulation is recommended as the category B - Epidemiological hazardous waste.

Transport and storage conditions

The information relating to suitable environmental condition and proper handling during transport and storage of the Leoced hemodialyzers is provided in the IFU and product labeling in form of international symbols and accompanying text (labeling is a part of operational documentation).


The information is summarized below:

Instruction for use:	Leoced must be kept in a clean, dry environment away from direct sunlight, heat and humidity and at temperatures of 0-30°C. Leoced is a plastic product, and must not be exposed to vibration or impact during transportation, storage and handling.
Labeling (on boxes):	Fragile handle with care This way up Keep dry Temperature limits 0-30°C
Labeling (on product):	Temperature limits 0-30°C
Relative Humidity:	There are no limits identified for relative humidity
Atmospheric Pressure:	There are no limits identified for atmospheric pressure

Approved:

Laura F. Plath
Regulatory Affairs Manager
NxStage Medical, Inc.

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