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


 MATRI™
BONE

Instructions for Use of Medical Device

**Mineralized Matrix Collagen Matri™ BONE,
design variants**

«APPROVED»
BIOM'UP SA, France
/ COO
(position)
Jean-Luc / IMBERT
(name)
/ J. Imbert
(signature)
Stamp
biom'up
BIOM'UP SA
8, Allée Irène Joliot Curie - Bât. M7
69800 SAINT PRIEST
SIRET : 481 014 041 00034
N°TVA : FR 04 481 014 041

1. Name of medical device

Mineralized Matrix Collagen Matri™ BONE, design variants:

1. **Mineralized Matrix Collagen Matri™ BONE** with particle size range 80 - 200 µm, dimensions 12 x 12 x 12 mm (MAB121212S);
2. **Mineralized Matrix Collagen Matri™ BONE** with particle size range 80 - 200 µm, dimensions 24 x 12 x 12 mm (MAB241212S);
3. **Mineralized Matrix Collagen Matri™ BONE** with particle size range 80 - 200 µm, dimensions 30 x 30 x 6 mm (MAB303006S);
4. **Mineralized Matrix Collagen Matri™ BONE** with particle size range 80 - 200 µm, dimensions 35 x 60 x 6 mm (MAB356006S);
5. **Mineralized Matrix Collagen Matri™ BONE** with particle size range 0.5 - 1 mm, dimensions 12 x 12 x 12 mm (MAB121212B);
6. **Mineralized Matrix Collagen Matri™ BONE** with particle size range 0.5 - 1 mm, dimensions 24 x 12 x 12 mm (MAB241212B);
7. **Mineralized Matrix Collagen Matri™ BONE** with particle size range 0.5 - 1 mm, dimensions 30 x 30 x 6 mm (MAB303006B);
8. **Mineralized Matrix Collagen Matri™ BONE** with particle size range 0.5 - 1 mm, dimensions 35 x 60 x 6 mm (MAB356006B);

2. Details regarding manufacturer of medical device

Manufacturer/Developer:

BIOM'UP SA

8, Allée Irène Joliot-Curie, 69800 Saint-Priest, France

Tel.: +33 (0)4 86 57 36 10, Fax: +33 (0)4 37 69 00 84

3. Intended use of the medical device assigned by the manufacturer

Mineralized Matrix Collagen Matri™ BONE, design variants (hereafter Matri™ BONE) is used for filling different types of bone lesions resulting from surgery, trauma, or pathology.

Indications:

In orthopaedics it is indicated in case of:

- bone sampling (samples of cortico-cancellous bone, ligament graft, etc.);
- bone tumour resection;
- osteosynthesis;
- arthrodesis, spondylosyndesis;
- pseudarthrosis;
- osteotomy (sternotomy, etc.)
- any other mild bone lesion.

In the maxillo-facial and dental fields, the product is indicated for bone augmentation and bone reconstruction:

- filling of mild bone defects in maxillo-facial surgery and stomatology;
- reconstruction of the alveolar ridge;
- used in combination with guided tissue regeneration products or guided bone regeneration products (such as Cova™ MAX) as an implant for filling bone and sinus defects;
- filling of peri-implantation zones in combination with guided bone regeneration products.

Contraindications:

The product should not be implanted on an infected, inflamed, necrotic site or one that is likely to become so. The product should not be in contact with circulating cerebrospinal fluid. This is a pork product: it should not

be used in patients with a known or suspected sensitivity to these by-products. No controlled study has been conducted on the use of the product in children, pregnant or lactating women. Therefore, any use in children or women during pregnancy or lactation must be preceded by a very careful risks/benefits evaluation.

This product should not be used in the following cases:

- metabolic diseases (diabetes, osteomalacia, hyperparathyroidism);
- osteoporosis;
- severe renal or hepatic dysfunction;
- high-dose corticotherapy.

Interactions with other agents:

To avoid any interaction, contact of the product with iodine solutions and in general, any other non-physiological product, must be avoided. Collagen decreases the adhesiveness of certain cements. The product should not be used on sites where prosthetic material is fixed by cement.

Possible side effects:

No allergic reaction associated with this product has yet been observed. However, in light of its composition, allergic or inflammatory reactions cannot be completely excluded. The risks related to the intervention remain unchanged.



Application Conditions:



This product must be manipulated and implanted by qualified personnel according to the general principles of sterility and pre-medication.

4. Description of medical device

The product is a haemostatic, resorbable, osteoconductive bone substitution matrix that is comprised of hydroxyapatite, calcium phosphate and porcine a-cellular collagen. Current knowledge supports the theory that the porcine species is not affected by prion diseases. Nevertheless, our extraction process includes steps recognized for inactivating viruses. The mineralized part of the biomaterial has optimized porosity which provides the osteoconductive properties required for bone filling. Collagen provides its haemostatic properties and offers a growth matrix for the osteoblasts. Bone remodeling is optimum after 3-4 months.

5. Technical description of medical device

№	Version	Model		Weight		Young's Modulus	Maximum compressive strength
				product in the package	product without packaging		
1	Mineralized Matrix Collagen Matri TM BONE	Particle size range - 80 - 200 µm, Dimensions 12 x 12 x 12 mm (MAB121212S)		38,9 g	0,2 g	1930 MPa	1510 N/cm ²
2	Mineralized Matrix Collagen Matri TM BONE	Particle size range - 0.5 - 1 mm, Dimensions 12 x 12 x 12 mm (MAB121212B)		39,0 g	0,2 g	1950 MPa	1530 N/cm ²
3	Mineralized Matrix Collagen Matri TM BONE	Particle size range - 80 - 200 µm, Dimensions 24 x 12 x 12 mm (MAB241212S)		39,1 g	0,4 g	1970 MPa	1550 N/cm ²
4	Mineralized Matrix Collagen Matri TM BONE	Particle size range - 0.5 - 1 mm, Dimensions 24 x 12 x 12 mm (MAB241212B)		39,1 g	0,5 g	1970 MPa	1550 N/cm ²

	Mineralized Matrix Collagen Matri™ BONE	Particle size range - 80 - 200 µm, Dimensions 30 x 30 x 6 mm (MAB303006S)		39,2 g	0,6 g	2000 MPa	1570 N/cm ²
6	Mineralized Matrix Collagen Matri™ BONE	Particle size range - 0.5 - 1 mm, Dimensions 30 x 30 x 6 mm (MAB303006B)		39,2 g	0,5 g	2000 MPa	1570 N/cm ²
7	Mineralized Matrix Collagen Matri™ BONE	Particle size range - 80 - 200 µm, Dimensions 35 x 60 x 6 mm (MAB356006S)		39,9 g	1,2 g	2100 MPa	1610 N/cm ²
8	Mineralized Matrix Collagen Matri™ BONE	Particle size range - 0.5 - 1 mm, Dimensions 35 x 60 x 6 mm (MAB356006B)		39,9 g	1,2 g	2100 MPa	1610 N/cm ²


Matri™ BONE before hydration	Matri™ BONE after hydration
	

Table of general technical characteristics

Parameter	Specification
Performance	The mineralized part of the biomaterial has optimized porosity which provides the osteoconductive properties required for bone filling. Collagen provides its haemostatic properties and offers a growth matrix for the osteoblasts. Bone remodelling is optimum after 3-4 months.
Description	Collagen foam with bone substitute material
Defaults	No surface heterogeneity
Color	White to light yellow
Homogeneity	Homogen foam with an even spread of bone substitute
Composition of the medical device	• 10% collagen • 90% biphasic calcium phosphate (40% tricalcium phosphate ceramic + 60% synthetic hydroxyapatite)
Integrity	Dense foam
Mechanical properties	Foam can be pressed after hydration without significant bone substitute loss
Hydration time	< 240 seconds
Swelling	> 3% of weight
Smell	No specific chemical odour
pH	from 6.5 to 7.5
Loss on drying	≤10%
Bacterial Endotoxin	≤ 20 UI/device
Complete resorption time	occurs within 13 weeks

6. Instructions for use

For better handling, the product must be held with gloves or clean dry utensils.

According to basic surgical methods:

- the product can be cut with sterile surgical scissors to suit the desired size;
- the implantation site should not be rinsed;
- the product can be placed dry for better haemostasis;
- according to the indications, the product can be mixed with blood, bone marrow, or even stem cells for large defects;
- for orthopaedic treatment of deep, extensive bone defects and segmented defects of over 1 to 2 cm, mixing the matrix with a platelet concentrate or autologous bone shavings is recommended;
- the matrix is placed in contact with the healthy bone for optimum osteoconduction; in orthopaedics contact must be made with bone that has been well cleaned;
- the defect should be entirely filled; the product should be compacted without being squeezed down in the cavity that is to be filled;
- ideally, the product should be soaked in blood from healthy adjacent bone. Excess blood should be removed with a compress, so that if necessary, additional product can be introduced for a more complete filling;
- the product is intended to be left on the site. It should be removed in case of post-operative infection;
- the surgical intervention is terminated in the usual manner.

In orthopaedics it is strongly recommended to place a drain at the site, but the product and the drains must not be in contact.

In dental and maxillo-facial surgery, the muco-periosteal flap must cover the graft site and be sutured. If the site cannot be completely covered, a resorbable membrane (such as Cova™ MAX) must be used as a complement. In case of sinus filling, since the sinus cavity is not vascularized, soaking the product in blood is highly recommended. Compaction of the product is essential in this indication. The use of guided tissue regeneration membranes is also highly recommended to prevent any invagination of soft tissue in the sinus and to facilitate handling of the sinus membrane.

Considering that the product is osteo-conductive, it cannot regenerate bone by itself. The surgical technique should enable the product to be invaded by osteoblasts and endogenous growth factors.

7. Safety requirements

Precautions before use:

- Do not use after the expiry date indicated on the package;
- Ensure that the packaging of each product is intact before use;
- Do not use if the package or the product is damaged;
- Do not use if the product has been exposed to extreme heat;
- Each product is sterile and intended for a single use. The product should not be used or sterilized a second time. The consequences of re-sterilizing and/or re-using the product were not evaluated during the designing. Therefore, the mechanical strength and resorption capacity of the matrix are not guaranteed in case of re-sterilization and/or re-use.

Precautions for use/warning:

- The product must be manipulated and implanted by qualified personnel according to the general principles of sterility and pre-medication.
- After opening the sachet, the product must be used immediately on a site that has been prepared for the implantation.
- If the product is cut to the desired size, the strips from the cutting must absolutely not be re-used.
- Regardless of the reason for removal, a product that has been removed from an implant site should never be put back.

Inadequate regeneration is possible in certain unfavourable cases, notably in case of significant bone lesions or in regions of low bone regeneration. To avoid inadequate regeneration, please follow the indications and the instructions for use of the product closely.

- The product is intended for use as treatment of mild bone defects and presents no resistance to stress. In orthopaedics therefore, the matrix must, if necessary, be used in combination with an appropriate osteosynthetic material.
- In exceptional cases when the product has been preserved at extreme temperatures over long periods, its appearance and mechanical properties might be modified. These modifications do not affect patient security or haemostatic and bone reconstruction performances.

8. Labelling, operating conditions, transportation and storage

Labelling

Sterilization using
radiation

Do not re-sterilize

Do not reuse

Do not use if packaging
is damaged

Caution! Refer to the accompanying
documentation

Temperature Range

It conforms to European
Directive 93/42/EEC

The following Symbols are marked on transport package:

Avoid exposure of sunlight

Protect from moisture

Fragile, handle with care

Operation

After opening the primary packaging, the product is intended to be used immediately. There is no specific preparation of the device.

Transportation

This medical device is transported by all kinds of transport in covered vehicles in accordance with the requirements and regulations of the transportation of goods, which are applicable in each type of transport. During transportation, this medical device is stable at temperature from + 2° C to + 25° C, with a relative humidity of < 60% +/- 5%. Keep the device in sheltered from direct sunlight, moisture, heat and freezing place.

Storage

During storage, this medical product is stable at temperature from + 2° C to + 25° C, with a relative humidity of < 60% +/- 5%. Keep the device in sheltered from direct sunlight, moisture, heat and freezing place.

9. Warranty

BIOM'UP undertakes the sale and delivery of products complying with the requirements of current regulations applicable in France for this category of products. Products should be corresponded to the technical documentation of the manufacture.

BIOM'UP guarantees the replacement without charge, in the same quantities, of the product recognized defective by BIOM'UP, exclusive of any compensation or damages and on condition that the customer notify BIOM'UP of the defect immediately on its discovery, by any written means within the legal period of guarantee (expiry date).

BIOM'UP does not guarantee under any circumstances any indirect and/or intangible damage that could be caused by the products, by the customer or any third party for any other use not in accordance with the manufacturer's requirements.

10. Shelf life

Shelf life of the medical device is 36 months. This medical device must not be used after the date printed on the package.




11. Procedure for disposal or destruction of the medical device

Disposal of the medical device must be carried out in accordance with accepted standards of medical practice and local, state and federal laws and regulations.

12. Reclamation

Reclamations are to be sent to manufacturer's authorized official in the territory of the Russian Federation:
Limited Liability Company Technical-commercial Firm «MEDTRONIK» (LLC TKF «MEDTRONIK»),
121359, Moscow, Marshal Timoshenko str, h. 19, bldg. 1A,
Tel.: +7-499-1496924, Fax: +7-499-1413139, medtronik@medtronik.ru

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