
TRI HAWK
FOR THE BUR CONNOISSEUR

Xenopus Inc. Technical Information about the "Tri Hawk Dental Bur" For Registration in the Russian Federation

Written by: Scott Macdonald
Position: Senior Quality and Engineering Manager
Date: July 2016

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Table of contents

1. Letter of qualification for Scott Macdonald
2. Letter explaining diamond bur numbering system
3. Letter explaining presence of Strauss Company
4. Instructions for use FG carbide dental burs
5. Instructions for use FG diamond dental rotary instruments
6. Technical information FG carbide dental burs and diamond dental rotary instruments

XENOPUS


XENOPUS INC.
PO Box 619, 2 Prospect Road
Morrisburg, Ontario, Canada K0C 1X0
Tel: 613-543-4500 Fax: 613-543-4500
www.trihawk.com

To whom it may concern,

July 29, 2016


This letter is written to inform you that Scott Macdonald is the Senior Quality and Engineering Manager of Xenopus Inc. His responsibilities include the management of the Quality Control Department, which include the verification of product specifications based on ISO standards, the application of those standards, and the testing of product to ensure standards are met. Scott Macdonald is our expert in technical matters and as such is authorized to sign technical documents on behalf of Xenopus Inc.

Yours sincerely,

Joanne Kydd

C.O.O. Xenopus Inc.

Scott Macdonald

Senior Quality & Engineering Mgr.

XENOPUS

XENOPUS INC.
PO Box 619, 2 Prospect Road
Morrisburg, Ontario, Canada K0C 1X0
Tel: 613-543-4500 Fax: 613-543-4501
www.trihawk.com

10 June 2016

To whom it may concern,

Based on review of Conclusion No. 13/EB-16-299 the following comments are given:

Under section 8.2.1: The number codes for diamond instruments: The current numbering is correct. The first 2 digits of the number code, using the generic 80, correctly identify the material of the working part as diamond. The generic 80 designation is used throughout the industry and is correct based on Table 1 of ISO 6360-1/2004 where it states that "*Number 80 may be used for general information on diamond grit*". The generic number is used because diamond grit falls also under the designation of abrasive instruments under ISO 6360-6/2004 which states specifically in section 4.1 "*For diamond rotary instruments, the designation of grit size shall be used as given in ISO 7711-3*" and as such ISO 6360-6/2004 is the standard used to determine the numbers in position 10, 11, and 12 of the 15-digit number with reference to 7711-3/2004.

In addition as stated in section 5.3 of ISO 6360-1/2004 "*the third digit of the 15-digit overall number indicates... the binding of the grit of diamond and abrasive instruments.*" As such, for diamond instruments, the 3rd digit should be "6" read from table 3 in Section 5.3 of ISO 6360-1/2004 as plated metallic binding and not from Table 2 of the same standard.

Sincerely,

Scott Macdonald
Senior Quality & Engineering Mgr,
Xenopus Inc.

XENOPUS

XENOPUS INC.
PO Box 619, 2 Prospect Road
Morrisburg, Ontario, Canada K0C 1X0
Tel: 613-543-4500 Fax: 613-543-4501
www.trihawk.com

10 June 2016

To whom it may concern,

Based on review of Conclusion No. 13/EB-16-299 the following comments are given:

Under Section 8.3.2:

The company Strauss & Co Industrial Diamonds Ltd has been included in the Xenopus Inc. Technical File for the "Tri Hawk: Dental Diamond Bur", according to Directive 93/42/EEC concerning medical devices.

Strauss & Co Industrial Diamonds Ltd is an emergency supplier of FG Diamond Rotary Instruments to Xenopus Inc., in the event Xenopus Inc. is unable to meet production for an order. Although Strauss & Co Industrial Diamonds Ltd. supplies diamond rotary instruments infrequently to Xenopus Inc. We have been required by EEC regulations and by our notified body to include the Strauss Technical File in our diamond technical file.

Sincerely,

Scott Macdonald
Senior Quality & Engineering Mgr.,
Xenopus Inc.

Instructions for Use for FG Carbide Bur

In order for the **TRI HAWK bur** to perform as designed:

1. Prescribed use:

- a. Tri Hawk carbide bur is to be used by certified dental doctors for their own professional use in dental clinics, dental hospitals and dental offices.
- b. Tri Hawk carbide bur are to be used to cut hard materials in the mouth including but not limited to tooth structure, bone, and implant materials.
- c. Tri Hawk dental burs are single use instruments. Use the bur ONCE and throw it away just as you would a needle or scalpel. All brands of burs are compromised after the first use.
- d. There are no known side effects or contraindications to the proper and correct use of these instruments.

Note: Multiple usages and sterilization drastically reduce cutting performance and greatly increase the likelihood of breakage/fracture for ALL burs

2. Proper and safe use: The device is to be used on the instruction of and by Dental Doctor or licensed practitioner.

- a. **Clean and sterilize** burs in accordance with the directions below before use
- b. Dispose of used instruments after first use as per disposal instructions below.
- c. Always wear gloves when handling contaminated instruments.
- d. Eye protection must be worn to protect against ejected particles.
- e. Surgical masks must be worn to avoid inhalation of aerosol or dust generated.
- f. Burs are to be used in air powered turbine or electric contra-angle red ring 1:5 hand pieces at room temperature.
- g. Maintain handpieces in good working order and correctly lubricated.

Note: Ensure sure that your hand piece has been properly calibrated, balanced and maintained

- h. Ensure that the bur is securely gripped in the handpiece collet.
- i. Insert bur all the way to the end in hand piece chuck. Burs over 23mm (0.906 inches) in length should only be used per hand piece manufacturer's instructions. Burs must be securely tightened in hand piece chuck
- j. Do not use worn-out burs.
- k. Do not apply excessive pressure on the bur. Forces should not exceed 2 Newtons (approx. 0.45 pounds force). Higher forces do not result in a quicker material removal, but increase the thermal load on the tooth.
- l. Avoid removing the bur at too sharp an angle to avoid leverage and breakage.
- m. Spray cooling with a minimum flow of 50 ml/min must be provided. Additional cooling is required if burs longer than 21 mm or heads bigger than 2 mm are being used.
- n. The burs should be set to the working speed according to the graphical diagram below. It should be set before applying the bur to the object.

Carbide bur - guidance for maximum RPM:

Instrument head diameter 1/10 mm	Maximum permissible speed (RPM)	Recommended operational speed (RPM)
007 - 010	450,000	100,000 - 350,000
011 - 014	450,000	70,000 - 300,000
015 - 018	450,000	55,000 - 200,000
019 - 023	300,000	40,000 - 150,000
024 - 027	160,000	35,000 - 120,000
028 - 031	140,000	30,000 - 95,000

3. Product information:

a. Instrument Lengths

Head length (mm)	Overall length +/- 0.5 (mm)		
	Short	Regular	Extra Long
≤ 5.5	16.5	19	25
>5.5-8	18.5	21	27
>8-9	19.5	22	28
>9-10	21.5	24	30

Warning: Do not extend bur from chuck.

- b. **Nominal diameter:** Nominal diameter is designated in the last 3 digits of the ISO number on the label in tenths of a millimeter. Ex: a Tri Hawk 330 catalogue part has an ISO part number of 500-314-237-006-008 thus the part has a 0.8 mm nominal head diameter or 8 tenths of a millimeter.
- c. **Material and material hardness:** All Tri Hawk carbide dental burs are made from High carbon steel with a minimum hardness of 250 HV5 and have a tungsten carbide cutting tip.
- d. **Run-out:** Maximum radial run-out for all carbide burs is 0.05mm.
- e. **Number of blades and blade configurations:** number of blades are based on the type of toothing and the nominal diameter. See the table below for details of what Tri Hawk supplies.

Number of Blades

Type of toothing	ISO Number Designation	Nominal Diameter (1/10 mm)		
		004-020	021-023	024+
Right Helicoidal	006	6	8	10
Right Helicoidal with Left Cross Cut	007	6	8	10
Finishing bur, Straight, Medium	071	12	N/A	N/A

4. **Cleaning & Sterilization:** Instruments must be sterilized before first use and sterilized in an autoclave sterilizer.
 - a. **Prepare burs for sterilization:** Package the burs in containers dedicated for the selected autoclave machine. The autoclave manufacturers' instructions must be followed. Do not exceed the maximum recommended load.
 - b. **Inspection after preparation:** Visually inspect the burs after cleaning to ensure contamination has been removed.
 - c. **Sterilization by autoclave before first use:** Follow the autoclave manufacturer's instructions to sterilize the burs. Autoclave the instruments for a holding time not less than six minutes at a temperature of between 130°C and 135°C. Other validated time/temperature regimes may be used. The holding time is the minimum time for which the minimum temperature is sustained. Run the instruments through the drying cycle incorporated in the autoclave. Refer to the autoclave manual for available cycles and proper operation.
5. **Disposal:** Used burs should be considered to be contaminated and appropriate precautions should be taken when handling and disposing of a used bur. Used burs must be discarded in a "sharps" container.
6. **Storage:** The burs should be stored in the dry sterilization container (bur stand or pouch) until required up to the stated expiration date on the original package. Storage should be in dry, clean conditions between 0°C and 50°C. Staining or corrosion may occur if the burs are kept in a wet environment, therefore, to be avoided.
7. **Scope of Delivery:** Each bur is individually packaged and shipped inside a sealable plastic bag in quantities specified by each customer. **Each plastic bag is shipped with THIS "instructions for use" form.**

8. Transport: Burs are to be transported in their original (factory) packaging, placed in a sealable plastic bag. Burs must be transported between -50°C and 50°C. **All shipments must contain a copy of these instructions.**
9. If a manufacturing defect is found the burs will be reimbursed. A defect in the product can only be determined through testing of the returned product (a minimum of 20 burs in original packaging is required for official testing) by Xenopus Inc.. All potentially defective products in the Russian Federation must be returned via 000 «Белгио-Русс Технолоджи» (125466, Москва, ул. Юровская д. 92, пом. VI, ком. 1) Xenopus Inc.'s authorised representative. A manufacturing defect is limited to defects that are caused by a production process or faulty raw material. Xenopus Inc. takes no responsibility for misuse of its products nor problems caused by other equipment (e.g. turbine/contra-angle hand pieces).
10. Definition of Packaging Symbols:

Lot number: 

Expiry Date: 

Single use: 

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Note: It remains the responsibility of the facility and the trained personnel to ensure that these processes are performed using the appropriate equipment, materials and personnel. This process may require continuous validation and monitoring to ensure required results are achieved.

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Manufactured by:
Xenopus Inc.
2 prospect Rd. Morrisburg,
Ontario, Canada, K0C 1X0

Instructions for Use for FG Diamond Rotary instruments

In order for the *TRI HAWK diamond rotary instrument* to perform as designed:

1. **Prescribed use:**

- a. Tri Hawk diamond rotary instruments are to be used by certified dental doctors for their own professional use in dental clinics, dental hospitals and dental offices.
- b. Tri Hawk diamond rotary instruments are to be used to cut hard materials in the mouth including but not limited to tooth structure, bone, and implant materials.
- c. Tri Hawk dental instruments are single use instruments. Use the instrument ONCE and throw it away just as you would a needle or scalpel. All brands of instruments are compromised after the first use.
- d. There are no known side effects to the use of these instruments.

Note: Multiple usages and sterilization drastically reduce cutting performance and greatly increase the likelihood of breakage/fracture for ALL instruments

2. **Proper and safe use:** The device is to be used on the instruction of and by Dental Doctor or licensed practitioner.

- a. **Clean and sterilize** instruments in accordance with the directions below before use
- b. Dispose of used instruments after first use as per disposal instructions below.
- c. Always wear gloves when handling contaminated instruments.
- d. Eye protection must be worn to protect against ejected particles.
- e. Surgical masks must be worn to avoid inhalation of aerosol or dust generated.
- f. Instruments are to be used in air powered turbine or electric contra-angle red ring 1:5 hand pieces at room temperature.
- g. Maintain handpieces in good working order and correctly lubricated.

Note: Ensure sure that your hand piece has been properly calibrated, balanced and maintained

- h. Ensure that the instrument is securely gripped in the handpiece collet.
- i. Insert instrument all the way to the end in hand piece chuck. Instruments over 23mm (0.906 inches) in length should only be used per hand piece manufacturer's instructions. Instruments must be securely tightened in hand piece chuck
- j. Do not use worn-out instruments.
- k. Do not apply excessive pressure on the instrument. Forces should not exceed 2 Newtons (approx. 0.45 pounds force). Higher forces do not result in a quicker material removal, but increase the thermal load on the tooth.
- l. Avoid removing the instrument at too sharp an angle to avoid leverage and breakage.
- m. Spray cooling with a minimum flow of 50 ml/min must be provided. Additional cooling is required if instruments longer than 21 mm or heads bigger than 2 mm are being used.
- n. The instruments should be set to the working speed according to the graphical diagram below. It should be set before applying the instrument to the object.

Diamond rotary instrument - guidance for maximum RPM:

Instrument head diameter 1/10 mm	Maximum permissible speed (RPM)	Recommended operational speed (RPM)
007 - 010	450,000	100,000 - 220,000
011 - 014	450,000	70,000 - 220,000
015 - 018	450,000	55,000 - 160,000
019 - 023	300,000	40,000 - 120,000
024 - 027	160,000	35,000 - 110,000
028 - 031	140,000	30,000 - 95,000

3. Product information:

a. Instrument Lengths

Head length (mm)	Overall length +/- 0.5 (mm)		
	Short	Regular	Extra Long
≤ 5.5	16.5	19	25
>5.5-8	18.5	21	27
>8-9	19.5	22	28
>9-10	21.5	24	30

Warning: Do not extend instrument from chuck.

- b. Nominal diameter: Nominal diameter is designated in the last 3 digits of the ISO number on the label in tenths of a millimeter. Ex: a Tri Hawk 165-014C catalogue part has an ISO part number of 806-314-165-534-**014** thus the part has a 1.4 mm nominal head diameter or 14 tenths of a millimeter.
- c. Material and material hardness: All Tri Hawk diamond rotary instruments are made from stainless steel with a minimum hardness of 250 HV5 and use a metal plated metallic binding.
- d. Run-out: Maximum radial run-out is determined based on the grit size used. See the table below for determining maximum radial run-out:
- e. Designation, colour code, grit sizes: The table below shows the grit size and designation and colour designation as well as Tri Hawk catalogue abbreviation.

Designation, colour code, grit sizes, run-out

4.

Designation	Catalogue Abbreviation	Grit size range µm	Median grain size µm	Grit designation	Maximum Radial Run-out mm
ultra-fine	UF	4 to 14	8	M 4 to M 14	0.07
extra-fine	EF	10 to 36	25	M 10 to M 36	0.07
fine	F	27 to 76	46	M 27 to D 76	0.08
medium	M	64 to 126	107	D 64 to D 126	0.10
coarse	C	107 to 181	151	D 107 to D 181	0.12
very coarse	XC	151 to 213	181	D 151 to D 213	0.14

4. Cleaning & Sterilization: Instruments must be sterilized before first use and sterilized in an autoclave sterilizer.
 - a. Prepare instruments for sterilization: Package the instruments in containers dedicated for the selected autoclave machine. The autoclave manufacturers' instructions must be followed. Do not exceed the maximum recommended load.
 - b. Inspection after preparation: Visually inspect the instruments after cleaning to ensure contamination has been removed.
 - c. Sterilization by autoclave before first use: Follow the autoclave manufacturer's instructions to sterilize the instruments. Autoclave the instruments for a holding time not less than six minutes at a temperature of between 130°C and 135°C. Other validated time/temperature regimes may be used. The holding time is the minimum time for which the minimum temperature is sustained. Run the instruments through the drying cycle incorporated in the autoclave. Refer to the autoclave manual for available cycles and proper operation.
5. Disposal: Used instruments should be considered to be contaminated and appropriate precautions should be taken when handling and disposing of a used instrument. Used instruments must be discarded in a "sharps" container.
6. Storage: The instruments should be stored in the dry sterilization container (instrument stand or pouch) until required up to the stated expiration date on the original package. Storage should be in dry, clean conditions between 0°C and 50°C. Staining or corrosion may occur if the instruments are kept in a wet environment, therefore, to be avoided.
7. Scope of Delivery: Each instrument is individually packaged and shipped inside a sealable plastic bag in quantities specified by each customer. **Each plastic bag is shipped with THIS "instructions for use" form.**
8. Transport: Instruments are to be transported in their original (factory) packaging, placed in a sealable plastic bag. Instruments must be transported between -50°C and 50°C. **All shipments must contain a copy of these instructions.**
9. If a manufacturing defect is found the burs will be reimbursed. A defect in the product can only be determined through testing of the returned product (a minimum of 20 burs in original packaging is required for official testing) by Xenopus Inc.. All potentially defective products in the Russian Federation must be returned via 000 «Белгио-Русс Технолоджи» (125466, Москва, ул. Юровская д. 92, пом. VI, ком. 1) Xenopus Inc.'s authorised representative. A manufacturing defect is limited to defects that are caused by a production process or faulty raw material. Xenopus Inc. takes no responsibility for misuse of its products nor problems caused by other equipment (e.g. turbine/contrangle hand pieces).
10. Definition of Packaging Symbols:

Lot number: LOT Expiry Date:  Single use:  Conformité Européenne: 

Note: It remains the responsibility of the facility and the trained personnel to ensure that these processes are performed using the appropriate equipment, materials and personnel. This process may require continuous validation and monitoring to ensure required results are achieved.

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Manufactured by:
Xenopus Inc.
2 prospect Rd. Morrisinstrumentg,
Ontario, Canada, K0C 1X0

Item number	Item name**	material	sterility	re-usable / single use	Package*	Technical Info	Intended use
001-010C-SS	round coarse short diamond FG bur 806-313-001-534-010, 001-010C-SS	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 1.0mm, Head diameter: 1.0mm Total length: 16.5mm +/-0.5mm Coarse grit: diamond particle size 150-180µmm	Preparing small single surface cavities in children
001-014C-SS	round coarse short diamond FG bur, 806-313-001-534-014, 001-014C-SS	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 1.4mm, Head diameter: 1.4mm Total length: 16.5mm +/-0.5mm Coarse grit: diamond particle size 150-180µmm	Preparing small single surface cavities in children
001-018C-SS	round coarse short diamond FG bur, 806-313-001-534-018, 001-018C-SS	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 1.8mm, Head diameter: 1.8mm Total length: 16.5mm +/-1.5mm Coarse grit: diamond particle size 150-180µm	Preparing small single surface cavities in children
107-011C-SS	flat-end cylinder coarse short diamond FG bur, 806-313-107-534-011, 107-011C-SS	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 5.0mm, Head diameter: 1.1mm Total length: 16.5mm +/-0.5mm Coarse grit: diamond particle size 150-180µmm	Specifically used when flat surfaces are demanded. Sharpening angles, defining and flattening lateral or bottom walls in children
140-014C-SS	chamfer coarse short diamond FG bur, 806-313-140-534-014, 140-014C-SS	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 7.0mm, Head Diameter: 1.4mm Total length: 18.5mm +/-0.5mm Coarse grit: diamond particle size 150-180µmm	Provides cavities with roundish angles and heavy chamfers with a rounded general profile in children
141-014C-SS	chamfer coarse short diamond FG bur, 806-313-141-534-014, 141-014C-SS	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 8.0mm, Head Diameter: 1.4mm Total length: 19.5mm +/-0.5mm Coarse grit: diamond particle size 150-180µmm	Provides cavities with roundish angles and heavy chamfers with a rounded general profile in children
166-016C-SS	needle coarse short diamond FG bur, 806-313-166-534-016, 166-016C-SS	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 10.0mm, Head Diameter: 1.6mm Total length: 21.5mm +/-0.5mm Coarse grit: diamond particle size 150-180µmm	Provides a delicate definition for or either proximal or buccal sides. Chamfering or slicing the dental tissue as well as when grinding the crown surface for children
173-016C-SS	flat-end taper coarse short diamond FG bur, 806-313-173-534-016, 173-016C-SS	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 10.0mm, Head Diameter: 1.6mm Total length: 21.5mm +/-0.5mm Coarse grit: diamond particle size 150-180µmm	Designed to provide expulsive walls in operative dentistry. Useful in placing shoulders and defining the terminal line for crowns.
199-016C-SS	round-end taper coarse short diamond FG bur, 806-313-199-534-016, 199-016C-SS	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 10.0mm, Head Diameter: 1.6mm Total length: 21.5mm +/-0.5mm Coarse grit: diamond particle size 150-180µmm	Designed for preparing crown abutments with chamfered shoulders, large reductions of tooth structures and large cavity preparations for children
249-014C-SS	flame coarse short diamond FG bur, 806-313-249-534-014, 249-014C-SS	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 8.0mm, Head Diameter: 1.4mm Total length: 19.5mm +/-0.5mm Coarse grit: diamond particle size 150-180µmm	Diamond point provides sharp chamfers in prosthetic dentistry for children
289-014C-SS	torpedo coarse short diamond FG bur, 806-313-289-534-014, 289-014C-SS	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 8.0mm, Head Diameter: 1.4mm Total length: 19.5mm +/-0.5mm Coarse grit: diamond particle size 150-180µmm	useful when producing chamfered terminal margins. Provides 130-degree shoulders during the contouring of crowns for children.

**SS = 313 = Short shank = 16.5mm total length

314 = regular shank burs (no suffix) = 19.5mm

SU and Surg = 316 = Surgical shank = 25mm

* consult technical file for full details of packaging

Verstion 2 by S.M. May 2016

257-020C-SS	football coarse short diamond FG bur, 806-313-257-534-020, 257-020C-SS	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 4.5mm, Head Diameter: 2.0mm Total length: 16.5mm +/-0.5mm Coarse grit: diamond particle size 150-180µm	used for lingual grinding on incisors and cuspids and reducing the side for crown coverage in prosthetic dentistry for children
001-010C	round coarse diamond FG bur, 806-314-001-534-010, 001-010C	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 1.0mm, Head Diameter: 1.0mm Total length: 19.0mm +/-0.5mm Coarse grit: diamond particle size 150-180µm	Preparing small single surface cavities
001-014C	round coarse diamond FG bur, 806-314-001-534-014, 001-014C	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 1.4mm, Head Diameter: 1.4mm Total length: 19.0mm +/-0.5mm Coarse grit: diamond particle size 150-180µm	Preparing small single surface cavities
001-018C	round coarse diamond FG bur, 806-314-001-534-018, 001-018C	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 1.8mm, Head Diameter: 1.8mm Total length: 19.0mm +/-0.5mm Coarse grit: diamond particle size 150-180µm	Preparing small single surface cavities
001-022F	round fine diamond FG bur, 806-314-001-514-022, 001-022F	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 2.2mm, Head Diameter: 2.2mm Total length: 19.0mm +/-0.5mm Fine grit: diamond particle size 50µm	Preparing small single surface cavities
010-016C	inverted cone coarse diamond FG bur, 806-314-010-534-016, 001-016C	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 1.6mm, Head Diameter: 1.6mm Total length: 19.0mm +/-0.5mm Coarse grit: diamond particle size 150-180µm	Inverted cone burs are used primarily for producing undercuts in preparations where an increased strength of bond is required. They are also useful in the flattening of pulpal and gingival walls.
032-014C	double inverted cone coarse diamond FG bur, 806-314-032-534-014, 032-014C	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 2.6mm, Head Diameter: 1.4mm Total length: 19.0mm +/-0.5mm Coarse grit: diamond particle size 150-180µm	Inverted cone burs are used primarily for producing undercuts in preparations where an increased strength of bond is required. They are also useful in the flattening of pulpal and gingival walls.
255-016C	inverted cone coarse diamond FG bur, 806-314-255-534-016, 255-016C	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 6.5mm, Head Diameter: 1.6mm Total length: 21.0mm +/-0.5mm Coarse grit: diamond particle size 150-180µm	Inverted cone burs are used primarily for producing undercuts in preparations where an increased strength of bond is required. They are also useful in the flattening of pulpal and gingival walls.
289-012C	torpedo coarse diamond FG bur, 806-314-289-534-012, 289-012C	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 8.0mm, Head Diameter: 1.2mm Total length: 21.0mm +/-0.5mm Coarse grit: diamond particle size 150-180µm	useful when producing chamfered terminal margins. Provides 130-degree shoulders during the contouring of crowns
289-012F	torpedo fine diamond FG bur, 806-314-289-514-012, 289-012F	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 8.0mm, Head Diameter: 1.2mm Total length: 21.0mm +/-0.5mm Fine grit: diamond particle size 50µm	useful when producing chamfered terminal margins. Provides 130-degree shoulders during the contouring of crowns
289-014C	torpedo coarse diamond FG bur, 806-314-289-534-014, 289-014C	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 8.0mm, Head Diameter: 1.4mm Total length: 22.0mm +/-0.5mm Coarse grit: diamond particle size 150-180µm	useful when producing chamfered terminal margins. Provides 130-degree shoulders during the contouring of crowns
298-014F	gingival curitage fine diamond FG bur, 806-314-298-524-014, 298-014F	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 8.0mm, Head Diameter: 1.4mm Total length: 21.0mm +/-0.5mm Fine grit: diamond particle size 50µm	useful when producing chamfered terminal margins. Provides 130-degree shoulders during the contouring of crowns

298-014C	gingival curettage coarse diamond FG bur, 806-314-298-534-014, 298-014C	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 8.0mm, Head Diameter: 1.4mm Total length: 21.0mm +/-0.5mm Coarse grit: diamond particle size 150-180µm	Used for Gingival curettage preparations and gingival margin preparations.
298-016C	gingival curettage coarse diamond FG bur, 806-314-298-534-016, 298-016C	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 8.0mm, Head Diameter: 1.6mm Total length: 21.0mm +/-0.5mm Coarse grit: diamond particle size 150-180µm	Used for Gingival curettage preparations and gingival margin preparations.
298-018C	gingival curettage coarse diamond FG bur, 806-314-298-534-018, 298-018C	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 8.0mm, Head Diameter: 1.8mm Total length: 20.0mm +/-0.5mm Coarse grit: diamond particle size 150-180µm	Used for Gingival curettage preparations and gingival margin preparations.
540-010F	interproximal/separation fine diamond FG bur, 806-314-540-514-010, 540-010F	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 4.0mm, Head Diameter: 1.0mm Total length: 19.0mm +/-0.5mm Fine grit: diamond particle size 50µm	Used for small interproximal cavity preparations and for interproximal reduction
109-011C	flat end cylinder coarse diamond FG bur, 806-314-109-534-011, 109-011C	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 5.0mm, Head Diameter: 1.1mm Total length: 19.0mm +/-0.5mm Coarse grit: diamond particle size 150-180µm	Specifically used when flat surfaces are demanded. Sharpening angles, defining and flattening lateral or bottom walls
111-014C	flat end cylinder coarse diamond FG bur, 806-314-111-534-014, 111-014C	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 8.5mm, Head Diameter: 1.4mm Total length: 22.0mm +/-0.5mm Coarse grit: diamond particle size 150-180µm	Specifically used when flat surfaces are demanded. Sharpening angles, defining and flattening lateral or bottom walls
130-012C	pointed cylinder coarse diamond FG bur, 806-314-130-534-012, 130-012C	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 8.0mm, Head Diameter: 1.2mm Total length: 21.0mm +/-0.5mm Coarse grit: diamond particle size 150-180µm	Diamond point provides sharp chamfers in prosthetic dentistry.
140-014C	chamfer(round end cylinder) coarse diamond FG bur, 806-314-140-534-014, 140-014C	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 7.0mm, Head Diameter: 1.4mm Total length: 21.0mm +/-0.5mm Coarse grit: diamond particle size 150-180µm	Provides cavities with roundish angles and heavy chamfers with a rounded general profile
141-014C	chamfer(round end cylinder) coarse diamond FG bur, 806-314-141-534-014, 141-014C	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 8.0mm, Head Diameter: 1.4mm Total length: 22.0mm +/-0.5mm Coarse grit: diamond particle size 150-180µm	Provides cavities with roundish angles and heavy chamfers with a rounded general profile
233-012M	pear medium diamond FG bur, 806-314-233-524-012, 233-012M	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 3.2mm, Head Diameter: 1.2mm Total length: 19.0mm +/-0.5mm Medium grit: diamond particle size 107-120µm	may be used for much of the preparation of conservative to moderately large occlusal cavities in anterior teeth
249-012C	flame coarse diamond FG bur, 806-314-249-534-012, 249-012C	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 8.0mm, Head Diameter: 1.2mm Total length: 21.0mm +/-0.5mm Coarse grit: diamond particle size 150-180µm	Diamond point provides sharp chamfers in prosthetic dentistry.
249-014F	flame fine diamond FG bur, 806-314-249-514-014, 249-014F	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 8.0mm, Head Diameter: 1.4mm Total length: 21.0mm +/-0.5mm Fine grit: diamond particle size 50µm	Diamond point provides sharp chamfers in prosthetic dentistry.

**SS = 313 = Short shank = 16.5mm total length

314 = regular shank burs (no suffix) = 19.5mm

315 = long shank burs (no suffix) = 25mm

* consult technical file for full details of packaging

Version 2 by S.M. May 2016

249-014C	flame coarse diamond FG bur, 806-314-249-534-014, 249-014C	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 8.0mm, Head Diameter: 1.4mm Total length: 21.0mm +/-0.5mm Coarse grit: diamond particle size 150-180µm	Diamond point provides sharp chamfers in prosthetic dentistry.
249-016C	flame coarse diamond FG bur, 806-314-249-534-016, 249-016C	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 8.0mm, Head Diameter: 1.6mm Total length: 21.0mm +/-0.5mm Coarse grit: diamond particle size 150-180µm	Diamond point provides sharp chamfers in prosthetic dentistry.
257-020F	football fine diamond FG bur, 806-314-257-514-020, 257-020F	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 4.5mm, Head Diameter: 2.0mm Total length: 19.0mm +/-0.5mm Fine grit: diamond particle size 50µm	used for lingual grinding on incisors and cuspids and reducing the side for crown coverage in prosthetic dentistry.
257-020C	football coarse diamond FG bur, 806-314-257-534-020, 257-020C	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 4.5mm, Head Diameter: 2.0mm Total length: 19.0mm +/-0.5mm Coarse grit: diamond particle size 150-180µm	used for lingual grinding on incisors and cuspids and reducing the side for crown coverage in prosthetic dentistry.
277-018C	football coarse diamond FG bur, 806-314-277-534-018, 277-018C	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 4.0mm, Head Diameter: 1.8mm Total length: 19.0mm +/-0.5mm Coarse grit: diamond particle size 150-180µm	used for lingual grinding on incisors and cuspids and reducing the side for crown coverage in prosthetic dentistry.
277-023C	football coarse diamond FG bur, 806-314-277-534-023, 277-023C	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 5.0mm, Head Diameter: 2.3mm Total length: 19.0mm +/-0.5mm Coarse grit: diamond particle size 150-180µm	used for lingual grinding on incisors and cuspids and reducing the side for crown coverage in prosthetic dentistry.
035-033C	barrel coarse diamond FG bur, 806-314-035-534-033, 035-033C	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 4.4mm, Head Diameter: 3.3mm Total length: 19.0mm +/-0.5mm Coarse grit: diamond particle size 150-180µm	Large cavity and crown preparations
037-037C	barrel coarse diamond FG bur, 806-314-037-534-037, 037-037C	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 7.0mm, Head Diameter: 3.7mm Total length: 21.0mm +/-0.5mm Coarse grit: diamond particle size 150-180µm	Large cavity and crown preparations
165-014C	needle coarse diamond FG bur, 806-314-165-534-014, 165-014C	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 8.0mm, Head Diameter: 1.4mm Total length: 22.0mm +/-0.5mm Coarse grit: diamond particle size 150-180µm	Provides a delicate definition for or either proximal or buccal sides. Chamfering or slicing the dental tissue as well as when grinding the crown surface
166-014F	needle fine diamond FG bur, 806-314-166-514-014, 166-014F	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 10.0mm, Head Diameter: 1.4mm Total length: 21.0mm +/-0.5mm Fine grit: diamond particle size 50µm	Provides a delicate definition for or either proximal or buccal sides. Chamfering or slicing the dental tissue as well as when grinding the crown surface
166-014C	needle coarse diamond FG bur, 806-314-166-534-014, 166-014C	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 10.0mm, Head Diameter: 1.4mm Total length: 21.0mm +/-0.5mm Coarse grit: diamond particle size 150-180µm	Provides a delicate definition for or either proximal or buccal sides. Chamfering or slicing the dental tissue as well as when grinding the crown surface
166-016C	needle coarse diamond FG bur, 806-314-166-534-016, 166-016C	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 10.0mm, Head Diameter: 1.6mm Total length: 21.0mm +/-0.5mm Coarse grit: diamond particle size 150-180µm	Provides a delicate definition for or either proximal or buccal sides. Chamfering or slicing the dental tissue as well as when grinding the crown surface

**SS = 313 = Short shank = 16.5mm total length
314 = regular shank burs (no suffix) = 19.5mm
SU and Surg = 316 = Surgical shank = 25mm

* consult technical file for full details of packaging

Version 2.0.0.0

172-018C	flat-end taper coarse diamond FG bur, 806-314-172-534-018, 172-018C, 172-018C	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 8.0mm, Head Diameter: 1.8mm Total length: 22.0mm +/-0.5mm Coarse grit: diamond particle size 150-180µm	Designed to provide expulsive walls in operative dentistry. Useful in placing shoulders and defining the terminal line for crowns.
173-016F	flat-end taper fine diamond FG bur, 806-314-173-514-016, 173-016F	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 10.0mm, Head Diameter: 1.6mm Total length: 24.0mm +/-0.5mm Fine grit: diamond particle size 50µm	Designed to provide expulsive walls in operative dentistry. Useful in placing shoulders and defining the terminal line for crowns.
173-016C	flat-end taper coarse diamond FG bur, 806-314-173-534-016, 173-016C	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 10.0mm, Head Diameter: 1.6mm Total length: 24.0mm +/-0.5mm Coarse grit: diamond particle size 150-180µm	Designed to provide expulsive walls in operative dentistry. Useful in placing shoulders and defining the terminal line for crowns.
197-025C	round-end taper coarse diamond FG bur, 806-314-197-534-025, 197-025C	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 7.0mm, Head Diameter: 2.5mm Total length: 21.0mm +/-0.5mm Coarse grit: diamond particle size 150-180µm	Designed for preparing crown abutments with chamfered shoulders, large reductions of tooth structures and large cavity preparations
198-014C	round-end taper coarse diamond FG bur, 806-314-198-534-014, 198-014C	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 8.0mm, Head Diameter: 1.4mm Total length: 21.0mm +/-0.5mm Coarse grit: diamond particle size 150-180µm	Designed for preparing crown abutments with chamfered shoulders, large reductions of tooth structures and large cavity preparations
198-018C	round-end taper coarse diamond FG bur, 806-314-198-534-018, 198-018C	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 8.0mm, Head Diameter: 1.8mm Total length: 21.0mm +/-0.5mm Coarse grit: diamond particle size 150-180µm	Designed for preparing crown abutments with chamfered shoulders, large reductions of tooth structures and large cavity preparations
199-014C	round end taper coarse diamond FG bur, 806-314-199-534-014, 199-014C	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 10.0mm, Head Diameter: 1.4mm Total length: 24.0mm +/-0.5mm Coarse grit: diamond particle size 150-180µm	Designed for preparing crown abutments with chamfered shoulders, large reductions of tooth structures and large cavity preparations
199-016F	round end taper fine diamond FG bur, 806-314-199-514-016, 199-016F	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 10.0mm, Head Diameter: 1.6mm Total length: 24.0mm +/-0.5mm Coarse grit: diamond particle size 50µm	Designed for preparing crown abutments with chamfered shoulders, large reductions of tooth structures and large cavity preparations
199-016C	round-end taper coarse diamond FG bur, 806-314-199-534-016, 199-016C	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 10.0mm, Head Diameter: 1.6mm Total length: 24.0mm +/-0.5mm Coarse grit: diamond particle size 150-180µm	Designed for preparing crown abutments with chamfered shoulders, large reductions of tooth structures and large cavity preparations
198-018XC	round end taper super coarse diamond FG bur, 806-314-544-018, 198-018XC	Shank: AISI 420F Stainless Steel; Diamond: 100% Natural Diamonds; Plating binder: Nickel Sulphamate	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 8.0mm, Head Diameter: 1.8mm Total length: 21.0mm +/-0.5mm Coarse grit: diamond particle size 180 - 250µm	Used for grinding zirconium, E-max and other hard ceramics and for large reduction preparations.
FG 2SS	round short carbide FG bur, 500-313-001-006-010, FG 2SS	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 1.0mm, Head Diameter: 1.0mm Total length: 16.5mm +/-0.5mm Number of blades: 6	Round carbide burs are popular for the preparation of conservative single surface cavities in children
FG 4SS	round short carbide FG bur, 500-313-001-006-014, FG 2SS	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 1.4mm, Head Diameter: 1.4mm Total length: 16.5mm +/-0.5mm Number of blades: 6	Round carbide burs are popular for the preparation of conservative single surface cavities in children

**SS = 313 = Short shank = 16.5mm total length
314 = regular shank burs (no suffix) = 19.5mm
SU and Surg = 316 = Surgical shank = 25mm

* consult technical file for full details of packaging

Version 2 by S.M. May 2016

FG 33 1/2SS	inverted cone short carbide FG bur, 500-313-011-006-006, FG 33 1/2SS	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 0.6mm, Head Diameter: 0.6mm Total length: 16.5mm +/-0.5mm Number of blades: 6	inverted cone carbide burs are used primarily for producing undercuts at the junction of the pulpal floor and lateral walls in occlusal (class I) cervical (class IV) cavities and the occlusal locks of the class II cavities in children
FG 35SS	inverted cone short carbide FG bur, 500-313-011-006-010, FG 35SS	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 1.0mm, Head Diameter: 1.0mm Total length: 16.5mm +/-0.5mm Number of blades: 6	inverted cone carbide burs are used primarily for producing undercuts at the junction of the pulpal floor and lateral walls in occlusal (class I) cervical (class IV) cavities and the occlusal locks of the class II cavities in children
FG 330SS	Pear short carbide FG bur, 500-313-237-006-008, 330SS	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 2.0mm, Head Diameter: 0.8mm Total length: 16.5mm +/-0.5mm Number of blades: 6	may be used for much of the preparation of conservative to moderately large occlusal cavities in anterior teeth in children
FG 331SS	Pear short carbide FG bur, 500-313-237-006-010, FG 331SS	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 2.0mm, Head Diameter: 1.0mm Total length: 16.5mm +/-0.5mm Number of blades: 6	may be used for much of the preparation of conservative to moderately large occlusal cavities in anterior teeth in children
FG 245SS	Amalgam reduction short carbide FG bur, 500-313-236-006-008, FG 245SS	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 3.0mm, Head Diameter: 0.8mm Total length: 16.5mm +/-0.5mm Number of blades: 6	may be used for much of the preparation of conservative to moderately large occlusal cavities in anterior teeth in children
FG 56SS	straight fissure short carbide FG bur, 500-313-109-006-008, FG 56SS	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 4.2mm, Head Diameter: 0.8mm Total length: 16.5mm +/-0.5mm Number of blades: 6	used wherever straight parallel sides and flat floors are desired. Some applications include establishing initial outline form, gaining access to carious dentine and preparation of lateral walls in children
FG 57SS	straight fissure short carbide FG bur, 500-313-109-006-010, FG 57SS	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 4.2mm, Head Diameter: 1.0mm Total length: 16.5mm +/-0.5mm Number of blades: 6	used wherever straight parallel sides and flat floors are desired. Some applications include establishing initial outline form, gaining access to carious dentine and preparation of lateral walls in children
FG 556SS	straight fissure crosscut short carbide FG bur, 500-313-109-007-008, FG 556SS	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 4.2mm, Head Diameter: 0.8mm Total length: 16.5mm +/-0.5mm Number of blades: 6	used wherever straight parallel sides and flat floors are desired. Some applications include establishing initial outline form, gaining access to carious dentine and preparation of lateral walls in children
FG 557SS	straight fissure crosscut short carbide FG bur, 500-313-109-007-010, FG 557SS	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 4.2mm, Head Diameter: 1.0mm Total length: 16.5mm +/-0.5mm Number of blades: 6	used wherever straight parallel sides and flat floors are desired. Some applications include establishing initial outline form, gaining access to carious dentine and preparation of lateral walls in children
FG 170SS	tapered fissure short carbide FG bur, 500-313-170-006-010, FG 170SS	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 4.2mm, Head Diameter: 1.0mm Total length: 16.5mm +/-0.5mm Number of blades: 6	the removal of enamel, the access to the decayed dentine, preparation of composite, gold or ceramic resin, the precise forming of the cervical cavities in children
FG 700SS	tapered fissure short cross cut carbide FG bur, 500-313-170-007-010, FG 700SS	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 4.2mm, Head Diameter: 1.0mm Total length: 16.5mm +/-0.5mm Number of blades: 6	the removal of enamel, the access to the decayed dentine, preparation of composite, gold or ceramic resin, the precise forming of the cervical cavities in children
FG 1/2	round carbide FG bur, 500-314-001-006-006, FG 1/2	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 0.6mm, Head Diameter: 0.6mm Total length: 19.0mm +/-0.5mm Number of blades: 6	Round carbide burs are popular for the preparation of conservative single surface cavities.

**SS = 313 = Short shank = 16.5mm total length
314 = regular shank burs (no suffix) = 19.5mm
SU and Surg = 316 = Surgical shank = 25mm

* consult technical file for full details of packaging

Version 2 by S.M. May 2016

FG 2	round carbide FG bur, 500-314-001-006-010, FG 2	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 1.0mm, Head Diameter: 1.0mm Total length: 19.0mm +/-0.5mm Number of blades: 6	Round carbide burs are popular for the preparation of conservative single surface cavities.
FG 4	round carbide FG bur, 500-314-001-006-014,, FG 4	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 1.4mm, Head Diameter: 1.4mm Total length: 19.0mm +/-0.5mm Number of blades: 6	Round carbide burs are popular for the preparation of conservative single surface cavities.
FG 6	round carbide FG bur, 500-314-001-006-018, FG 6	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 1.8mm, Head Diameter: 1.8mm Total length: 19.0mm +/-0.5mm	Round carbide burs are popular for the preparation of conservative single surface cavities.
FG 8	round carbide FG bur, 500-314-001-006-023, FG 8	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 2.3mm, Head Diameter: 2.3mm Total length: 19.0mm +/-0.5mm Number of blades: 8	Round carbide burs are popular for the preparation of conservative single surface cavities.
FG 33 1/2	inverted cone carbide FG bur, 500-314-011-006-006, FG 33.5	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 0.6mm, Head Diameter: 0.6mm Total length: 19.0mm +/-0.5mm Number of blades: 6	inverted cone carbide burs are used primarily for producing undercuts at the junction of the pulpal floor and lateral walls in occlusal (class I) cervical (class IV) cavities and the occlusal locks of the class II cavities.
FG 35	inverted cone carbide FG bur, 500-314-011-006-006, FG 35	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 1.0mm, Head Diameter: 1.0mm Total length: 19.0mm +/-0.5mm Number of blades: 6	inverted cone carbide burs are used primarily for producing undercuts at the junction of the pulpal floor and lateral walls in occlusal (class I) cervical (class IV) cavities and the occlusal locks of the class II cavities.
FG 36	inverted cone carbide FG bur, 500-314-011-006-012, FG 36	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 1.2mm, Head Diameter: 1.2mm Total length: 19.0mm +/-0.5mm Number of blades: 6	inverted cone carbide burs are used primarily for producing undercuts at the junction of the pulpal floor and lateral walls in occlusal (class I) cervical (class IV) cavities and the occlusal locks of the class II cavities.
FG 37	inverted cone carbide FG bur, 500-314-011-006-014, FG 37	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 1.4mm, Head Diameter: 1.4mm Total length: 19.0mm +/-0.5mm Number of blades: 6	inverted cone carbide burs are used primarily for producing undercuts at the junction of the pulpal floor and lateral walls in occlusal (class I) cervical (class IV) cavities and the occlusal locks of the class II cavities.
FG 330	pear carbide FG bur, 500-314-237-006-008, FG 330	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 2.0mm, Head Diameter: 0.8mm Total length: 19.0mm +/-0.5mm Number of blades: 6	may be used for much of the preparation of conservative to moderately large occlusal cavities in anterior teeth
FG 331	pear carbide FG bur, 500-314-237-006-010, FG 331	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 2.0mm, Head Diameter: 1.0mm Total length: 19.0mm +/-0.5mm Number of blades: 6	may be used for much of the preparation of conservative to moderately large occlusal cavities in anterior teeth
FG 332	pear carbide FG bur, 500-314-237-006-012, FG 332	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 2.0mm, Head Diameter: 1.2mm Total length: 19.0mm +/-0.5mm Number of blades: 6	may be used for much of the preparation of conservative to moderately large occlusal cavities in anterior teeth
FG 245	amalgam reduction carbide FG bur, 500-314-236-006-008, FG 245	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 3.0mm, Head Diameter: 0.8mm Total length: 19.0mm +/-0.5mm Number of blades: 6	may be used for much of the preparation of conservative to moderately large occlusal cavities in anterior teeth.

**SS = 313 = Short shank = 16.5mm total length
 314 = regular shank burs (no suffix) = 19.5mm
 SU and Surg = 316 = Surgical shank = 25mm

* consult technical file for full details of packaging

Version 2 by S.M. May 2016

FG 56	straight fissure carbide FG bu, 500-314-109-006-008,	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 4.2mm, Head Diameter: 0.8mm Total length: 19.0mm +/-0.5mm Number of blades: 6	used wherever straight parallel sides and flat floors are desired. Some applications include establishing initial outline form, gaining access to carious dentine and preparation of laterall walls.
FG 57	straight fissure carbide FG bur, 500-314-109-006-010, FG 57	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 4.2mm, Head Diameter: 1.0mm Total length: 19.0mm +/-0.5mm Number of blades: 6	used wherever straight parallel sides and flat floors are desired. Some applications include establishing initial outline form, gaining access to carious dentine and preparation of laterall walls.
FG 58	straight fissure carbide FG bur, 500-314-109-006-012, FG 58	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 4.2mm, Head Diameter: 1.2mm Total length: 19.0mm +/-0.5mm Number of blades: 6	used wherever straight parallel sides and flat floors are desired. Some applications include establishing initial outline form, gaining access to carious dentine and preparation of laterall walls.
FG 556	straight fissure crosscut carbide FG bur, 500-314-109-007-008, FG 556	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 4.2mm, Head Diameter: 0.8mm Total length: 19.0mm +/-0.5mm Number of blades: 6	used wherever straight parallel sides and flat floors are desired. Some applications include establishing initial outline form, gaining access to carious dentine and preparation of laterall walls.
FG 557	straight fissure crosscut carbide FG bur, 500-314-109-007-010, FG 557	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 4.2mm, Head Diameter: 1.0mm Total length: 19.0mm +/-0.5mm Number of blades: 6	used wherever straight parallel sides and flat floors are desired. Some applications include establishing initial outline form, gaining access to carious dentine and preparation of laterall walls.
FG 558	straight fissure crosscut carbide FG bur, 500-314-109-007-012, FG 558	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 4.2mm, Head Diameter: 1.2mm Total length: 19.0mm +/-0.5mm Number of blades: 6	used wherever straight parallel sides and flat floors are desired. Some applications include establishing initial outline form, gaining access to carious dentine and preparation of laterall walls.
FG 170	tapered fissure carbide FG bur, 500-314-170-006-010, FG 170	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 4.2mm, Head Diameter: 1.0mm Total length: 19.0mm +/-0.5mm Number of blades: 6	the removal of enamel, the access to the decayed dentine, preparation of composite, gold or ceramic resin, the precise forming of the cervical cavities
FG 171	tapered fissure carbide FG bur, 500-314-170-006-012, FG 171	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 4.2mm, Head Diameter: 1.2mm Total length: 19.0mm +/-0.5mm Number of blades: 6	the removal of enamel, the access to the decayed dentine, preparation of composite, gold or ceramic resin, the precise forming of the cervical cavities
FG 700	tapered fissure cross cut carbide FG bur, 500-314-170-007-010, FG 700	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 4.2mm, Head Diameter: 1.0mm Total length: 19.0mm +/-0.5mm Number of blades: 6	the removal of enamel, the access to the decayed dentine, preparation of composite, gold or ceramic resin, the precise forming of the cervical cavities
FG 701	tapered fissure cross cut carbide FG bur, 500-314-170-007-012, FG 701	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 4.2mm, Head Diameter: 1.2mm Total length: 19.0mm +/-0.5mm Number of blades: 6	the removal of enamel, the access to the decayed dentine, preparation of composite, gold or ceramic resin, the precise forming of the cervical cavities
Talon 08	Talon 08 carbide FG bur, 500-314-139-006-008, Talon 08	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 4.2mm, Head Diameter: 0.8mm Total length: 19.0mm +/-0.5mm Number of blades: 6	Used for cutting all types of metal including: non-precious alloys, precious metal, and amalgam.
Talon 10	Talon 10 Carbide FG bur, 500-314-139-006-010, Talon 10	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 4.2mm, Head Diameter: 1.0mm Total length: 19.0mm +/-0.5mm Number of blades: 6	Used for cutting all types of metal including: non-precious alloys, precious metal, and amalgam.

**SS = 313 = Short shank = 16.5mm total length
314 = regular shank burs (no suffix) = 19.5mm
SU and Surg = 316 = Surgical shank = 25mm

* consult technical file for full details of packaging

Version 2 by S.M. May 2016

Technical information about the Tri Hawk Dental Rotary Instruments

Talon X10	Talon 10 cross cut Carbide FG bur, 500-314-139-007-010, Talon X10	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 4.2mm, Head Diameter: 1.0mm Total length: 19.0mm +/-0.5mm Number of blades: 6	Used for cutting all types of metal including: non-precious alloys, precious metal, and amalgam.
Talon 10SS	Talon 10 short Carbide FG bur, 500-313-139-006-010, Talon 10SS	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 4.2mm, Head Diameter: 1.0mm Total length: 16.5mm +/-0.5mm Number of blades: 6	Used for cutting all types of metal including: non-precious alloys, precious metal, and amalgam.
Talon 10SU	Talon 10 surgical Carbide FG bur, 500-316-139-006-010, Talon 10SU	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 4.2mm, Head Diameter: 1.0mm Total length: 25.0mm +/-0.5mm Number of blades: 6	Used for cutting all types of metal including: non-precious alloys, precious metal, and amalgam. Is excellent for removing wisdom teeth, and separating roots.
Talon 12	Talon 12 Carbide FG bur, 500-314-139-006-012, Talon 12	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 4.2mm, Head Diameter: 1.2mm Total length: 19.0mm +/-0.5mm Number of blades: 6	Used for cutting all types of metal including: non-precious alloys, precious metal, and amalgam.
Talon X12	Talon 12 cross cut Carbide FG bur, 500-314-139-007-012, Talon X12	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 4.2mm, Head Diameter: 1.2mm Total length: 19.0mm +/-0.5mm Number of blades: 6	Used for cutting all types of metal including: non-precious alloys, precious metal, and amalgam.
Talon 12SS	Talon 12 short Carbide FG bur, 500-313-139-006-012, Talon 12SS	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 4.2mm, Head Diameter: 1.2mm Total length: 16.5mm +/-0.5mm Number of blades: 6	Used for cutting all types of metal including: non-precious alloys, precious metal, and amalgam.
Talon 12SU	Talon 12 surgical Carbide FG bur, 500-316-139-006-012, Talon 12SU	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 4.2mm, Head Diameter: 1.2mm Total length: 25.0mm +/-0.5mm Number of blades: 6	Used for cutting all types of metal including: non-precious alloys, precious metal, and amalgam. Is excellent for removing wisdom teeth, and separating roots.
Talon 14SU	Talon 14 surgical Carbide FG bur, 500-316-139-006-012, Talon 14SU	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 4.2mm, Head Diameter: 1.4mm Total length: 25.0mm +/-0.5mm Number of blades: 6	Used for cutting all types of metal including: non-precious alloys, precious metal, and amalgam. Is excellent for removing wisdom teeth, and separating roots.
FG 7901	needle finishing carbide Fg bur, 500-314-207-071-009, FG 7901	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 3.5mm, Head Diameter: 0.9mm Total length: 19.0mm +/-0.5mm Number of blades: 12	Ideal burs for contouring and finishing of interproximal margins, occlusal margins, and facial surfaces. It can also be used for subgingival work.
FG 7902	needle finishing carbide Fg bur, 500-314-207-071-010, FG 7902	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 3.5mm, Head Diameter: 1.0mm Total length: 19.0mm +/-0.5mm Number of blades: 12	Ideal burs for contouring and finishing of interproximal margins, occlusal margins, and facial surfaces. It can also be used for subgingival work.
FG 7903	needle finishing carbide Fg bur, 500-314-207-071-012, FG 7903	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 3.5mm, Head Diameter: 1.2mm Total length: 19.0mm +/-0.5mm Number of blades: 12	Ideal burs for contouring and finishing of interproximal margins, occlusal margins, and facial surfaces. It can also be used for subgingival work.
FG 7404	Egg finishing carbide FG bur, 500-314-277-071-014, FG 7404	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 3.5mm, Head Diameter: 1.4mm Total length: 19.0mm +/-0.5mm Number of blades: 12	Egg shaped finishing burs are ideal for the lingual and occlusal surfaces of anterior and posterior teeth.

**SS = 313 = Short shank = 16.5mm total length

314 = regular shank burs (no suffix) = 19.5mm

Technical information about the Tri Hawk Dental Rotary Instruments

FG 7406	Egg finishing carbide FG bur, 500-314-277-071-018, , FG 7406	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 3.5mm, Head Diameter: 1.8mm Total length: 19.0mm +/-0.5mm Number of blades: 12	Egg shaped finishing burs are ideal for the lingual and occlusal surfaces of anterior and posterior teeth.
FG 2SU	round surgical carbide FG bur, 500-316-001-006-010, FG 2SU	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 1.0mm, Head Diameter: 1.0mm Total length: 25.0mm +/-0.5mm Number of blades: 6	Used to creat deep access into the pulp chamber.
FG 4SU	round surgical carbide FG bur, 500-316-001-006-014, FG 4SU	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 1.4mm, Head Diameter: 1.4mm Total length: 25.0mm +/-0.5mm Number of blades: 6	Used to creat deep access into the pulp chamber.
FG 6SU	round surgical carbide FG bur, 500-316-001-006-018, FG 6SU	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 1.8mm, Head Diameter: 1.8mm Total length: 25.0mm +/-0.5mm Number of blades: 6	Used to creat deep access into the pulp chamber.
FG 8SU	round surgical carbide FG bur, 500-316-001-006-023, FG 8SU	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 2.3mm, Head Diameter: 2.3mm Total length: 25.0mm +/-0.5mm Number of blades: 8	Used to creat deep access into the pulp chamber.
FG 57Surg	straight fissure surgical carbide FG bur, 500-316-109-006-010, Fg 57Surg	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 4.2mm, Head Diameter: 1.0mm Total length: 25.0mm +/-0.5mm Number of blades: 6	Longer shank provides greater visibility, also primarily used in root canal prepartions, tooth extractions and root separations.
FG 58Surg	straight fissure surgical carbide FG bur, 500-316-109-006-012, FG 58Surg	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 4.2mm, Head Diameter: 1.2mm Total length: 25.0mm +/-0.5mm Number of blades: 6	Longer shank provides greater visibility, also primarily used in root canal prepartions, tooth extractions and root separations.
FG 557Surg	straight fissure surgical crosscut carbide Fg bur, 500-316-109-007-010, FG 557Surg	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 4.2mm, Head Diameter: 1.0mm Total length: 25.0mm +/-0.5mm Number of blades: 6	Longer shank provides greater visibility, also primarily used in root canal prepartions, tooth extractions and root separations.
FG 558Surg	straight fissure surgical crosscut carbide FG bur, 500-316-109-007-012, FG 558Surg	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 4.2mm, Head Diameter: 1.2mm Total length: 25.0mm +/-0.5mm Number of blades: 6	Longer shank provides greater visibility, also primarily used in root canal prepartions, tooth extractions and root separations.
FG 171Surg	tapered fissure surgical carbide FG bur, 500-316-170-006-012, FG 171Surg	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 4.2mm, Head Diameter: 1.2mm Total length: 25.0mm +/-0.5mm Number of blades: 6	Longer shank provides greater visibility, also primarily used in root canal prepartions, tooth extractions and root separations.
FG 172Surg	tapered fissure surgical carbide FG bur, 500-316-170-006-014, FG 172Surg	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 4.2mm, Head Diameter: 1.4mm Total length: 25.0mm +/-0.5mm Number of blades: 6	Longer shank provides greater visibility, also primarily used in root canal prepartions, tooth extractions and root separations.
FG 701Surg	tapered fissure surgical cross cut carbide FG bur, 500-316-170-007-012, FG 701Surg	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 4.2mm, Head Diameter: 1.2mm Total length: 25.0mm +/-0.5mm Number of blades: 6	Longer shank provides greater visibility, also primarily used in root canal prepartions, tooth extractions and root separations.
FG 702SU	tapered fissure surgical cross cut carbide FG bur, 500-316-170-007-014, FG 702Surg	Shank: 1095 carbon steel Carbide: C2 ANSI 6-10% cobalt, remainder Tungsten Carbide	non-sterile	Single Use	foil on one side of the bur and clear Cellophane on the other	Head Length: 4.2mm, Head Diameter: 1.4mm Total length: 25.0mm +/-0.5mm Number of blades: 6	Longer shank provides greater visibility, also primarily used in root canal prepartions, tooth extractions and root separations.

**SS = 313 = Short shank = 16.5mm total length

314 = regular shank burs (no suffix) = 19.5mm

SU and Surg = 316 = Surgical shank = 25mm

* consult technical file for full details of packaging

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