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Release Date: October 2014

Warning

- Before using the system, check local regulations. If any local legislation is violated, use cannot be authorized.
 - In the United States, federal law restricts prescription medical devices to sale by or on the order of a physician, or properly licensed practitioner. Viora makes no representations regarding federal, state, or local laws or regulations that might apply to the use and operation of this system.
-

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Table of Contents

1.	Before You Start	9
1.1.	Conventions Used in the Manual.....	9
1.2.	Explanation of the Symbols Used in the System.....	11
2.	Introduction	13
2.1.	Regulating data for the construction of a medical device	14
2.1.1.	System Console	17
2.1.2.	Handpieces.....	17
2.1.3.	Disposable Sterilized Tip and Cover – V-FR Handpiece	17
2.1.4.	Filters and Spot for V-IPL.....	17
2.1.6.	Display Unit.....	18
2.2.	System Specifications	19
2.3.	Handpiece Specifications.....	19
2.3.1.	RF – V-ST Specifications	20
2.3.2.	RF – V-FR Specifications	20
2.3.3.	RF Specifications	20
2.3.4.	RF - V-Form - BC Applicator Specifications	20
2.3.5.	RF - V-Form - FC Applicator Specifications.....	20
2.3.6.	V-IPL Specifications	21
3.	Safety.....	22
3.1.	Technical requirements for medical device	23
3.2.	System Classifications.....	37
3.3.	Electrical Requirements	37
3.4.	Environmental Requirements.....	37
3.5.	Space Requirements	38
3.6.	Safety When Operating the System.....	38
3.7.	System Safety Features	38
3.8.	Potential Fire Hazards	39
3.9.	Patient and Personnel Safety	39
3.10.	Remote Interlock Connector	40

3.11.	Foot Pedal	40
3.12.	Safety in Maintaining the System.....	40
3.13.	Transport and Storage Conditions.....	41
3.14.	RF Based Technology Safety.....	41
3.14.1.	RF – V-ST Safety.....	41
3.14.2.	V-Form Handpiece Safety.....	41
3.14.3.	RF – V-FR Safety.....	41
3.15.	Light based Technology Safety	42
3.15.1.	Protective Eyewear	42
3.16.	Data for development of medical device.....	43
3.17.	Data for the manufacturing of medical device.....	44
4.	Labels.....	45
4.1.	System Labels.....	45
4.2.	Handpiece Labels.....	46
4.2.1.	RF – V-ST Handpiece and box	46
4.2.2.	RF – V-FR Handpiece and box	46
4.2.3.	RF – V-Form Handpiece and box.....	47
4.2.4.	V-IPL Handpiece and box	47
5.	System Installation	48
5.1.	The using data for medical device.....	48
5.2.	Unpacking.....	49
5.3.	Assembling the Cradle.....	51
5.4.	Connecting the Interlock	51
5.5.	Connecting the Foot Pedal.....	52
5.6.	Filling the Water Container	53
5.7.	Connecting the Power Cable	54
5.8.	Connecting the Handpieces	54
5.8.1.	Connecting the V-FR, V-ST, V-Form Handpieces	54
5.8.2.	Connecting the V-IPL Handpieces.....	56
5.8.3.	Circulating Water in V-IPL Handpieces.....	58
5.9.	Wheel Brakes.....	59
5.10.	Moving the System within the Facility	59
5.11.	Moving the System to another Facility.....	59

6	The data for operating medical device	60
6.1.	System On/Off	60
6.2.	Login/Logout	61
6.3.	Set time and Date	62
6.4.	Add Configuration.....	63
6.5.	Adding a New License	64
6.6.	Main Menu.....	65
6.7.	Set Up Menu.....	66
6.8.	Treatment Tips Screen	68
6.9.	RF – ST Treatment Screens	69
6.9.1.	ST Treatment Menu.....	69
6.9.2.	ST Basic Treatment Screen	69
6.9.3.	ST Advanced Treatment Screen	71
6.9.4.	V-ST Output Power	72
6.10.	FR Treatment Screens	73
6.10.1.	FR Treatment Menu	73
6.10.2.	FR Basic Treatment Screen	73
6.10.3.	FR Advanced Treatment Screen	74
6.10.4.	V-FR Output Power	75
6.11.	BC/FC Treatment Screens	77
6.11.1	BC or FC Treatment Screen	77
6.12.	V-IPL Treatment Screens.....	80
6.12.1.	Filter Confirmation Screen.....	80
6.12.2	V-IPL Setting Screen	82
6.12.3.	V-IPL Treatment Screen.....	83
6.14.	System Messages.....	85
7.	Clinical Guidelines for RF Based Technology.....	88
7.1.	Contraindications	88
7.2.	Possible Side Effects	89
7.3.	Realistic Expectations	89
7.4.	Consult Meeting.....	90
7.5.	Treatment Follow-Up.....	90
7.6.	V-ST Clinical Guide	91
7.6.1.	Indications of Use	91

7.6.2.	Test Procedure.....	91
7.6.3	Pre-Treatment Preparation.....	91
7.6.4.	Treatment Guidelines	92
7.6.5.	Treatment Schedule.....	92
7.6.6.	Post-Treatment Care.....	92
7.7.	V-FR Clinical Guide.....	92
7.7.1	Indications of use.....	92
7.7.2.	Pre-Treatment Preparation.....	93
7.7.3.	Test Procedure	93
7.7.4.	Treatment Guidelines	95
7.7.5.	Treatment Schedule.....	96
7.7.6.	Post-Treatment Care.....	96
7.8.	FC and BC Clinical Guide	97
7.8.1.	Indications of use.....	97
7.8.2.	Pre-Treatment Preparation.....	97
7.8.3.	Test Procedure	97
7.8.4.	Treatment Guidelines	100
7.8.5.	Treatment Schedule.....	101
7.8.6.	Post-Treatment Care.....	101
8.	Clinical Guidelines for Light Based Technology	101
8.1.	Contraindications	101
8.2.	Possible Side Effects	103
8.3.	Realistic Expectations	104
8.4.	Consult Meeting.....	104
8.5.	Test Procedure	105
8.6.	Treatment Follow-Up.....	106
8.7.	V-IPL Clinical Guide	106
8.7.1.	Indications of use.....	106
8.7.2.	Pre-Treatment Preparation.....	107
8.7.3.	Treatment Parameters.....	107
8.7.4.	Treatment Schedule.....	110
8.7.5.	Post-Treatment Care.....	110
8.8.1	Indications of use.....	111
8.8.2	Pre-Treatment Preparation.....	111

8.8.3.	Treatment Parameters.....	112
8.8.4.	Treatment Schedule.....	112
8.8.5.	Post-Treatment Care.....	113
9.	Data for the technical maintenance.....	114
9.1.	Cleaning the System.....	114
9.2.	Cleaning the V-ST, V-Form Electrodes.....	114
9.3.	Cleaning the V-IPL Handpiece.....	114
9.3.1.	Cleaning the Filters.....	114
9.3.2.	Cleaning the Spots.....	115
9.5.	Replacing the filters in the BC/FC applicators.....	116
9.6.	Replacing the Water Filter and Deionized Water.....	117
10.	Data for the service.....	118
10.1.	Fault Messages.....	118
10.2.	Failure Messages.....	122
10.3.	Operation troubleshooting.....	122
11.	About Viora.....	123
12.	Data for recycling.....	129
13.	System composition.....	129

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2. Introduction

V20 is a light based and Radiofrequency (RF) Non-Surgical medical aesthetic platform system.

The V20 platform, developed by Viora, utilizes radiofrequency (RF) for aesthetic and dermatological procedures, with applications for the face, neck, décolleté and hands. The system is based on the well-established technology of bipolar RF, which is extensively used for the coagulation and ablation of soft tissues. Over the past decade, RF technology has become widely used in the aesthetics and dermatology field for skin tightening and resurfacing.

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The V20 is in Class IIb according to Directive 93/42/EEC of 14 June 1993, amended by Council Directive 2007/47/EC of 5 September 2007, concerning medical devices. Annex IX, Rule 9 states "all active therapeutic devices intended to administer or exchange energy are in Class IIb".

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- Front connector for V-IPL handpieces
- Rear connector for V-FR, V-ST, V-Form handpieces
- RF Generator
- Switch Module
- Water Cooling system
- Display Unit
- An Emergency Button that completely stops V20 operation
- LED indicator used for V-IPL treatments
- Service panel

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1. RF V-ST
2. RF V-FR
3. V-Form
4. V-IPL

The V-FR, V-ST, V-Form handpieces are attached with a cable to the connector located on the rear panel of the main console.

The V-IPL handpieces are attached with a cable to the connector on the front panel of the main console.

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2.1.4 Filters and Spot for V-IPL

The V-IPL handpiece allows the operator to change the optical filters. The system will recognize the specific filter automatically.

Five (5) optical filters in **different** wavelength are used.

V-IPL handpiece has 2 options of spot size: the large spot is 6.4cm^2 and the small spot is 2.4cm^2 . The material of the **light** guide is Sapphire cooled by two TEC (Thermo Electro Cooling) units. The **operator** can change the spot size as necessary, according to the treatment area. Spot size will be automatically recognized by the system.

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1. Treatment type
2. Treatment settings
3. Pulse / Time counter
4. Treatment time

2.2 System Specifications

Supported Handpieces	RF – V-ST
	RF – V-FR
	V-Form
	V-IPL
Maximum Power	1000 watts
Electrical requirements	100-240 V~ 10A 50/60 Hz 100-120 V~ 10A 50/60 Hz Single Phase
Environmental Requirements	Operating temperature 18°-24°C (64°-75°F) and relative humidity from 10% to 80%, non- condensing
User Interface	10.4 inch touch screen
Weight	70 kg
Dimensions	112 (H) x 47 (W) x 67 (D) cm

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3 Safety

To ensure dependable, consistent and safe treatments, the following safety measures must be taken into consideration when operating and maintaining the system. The person operating the system must be familiar with the safety regulations provided in this chapter. The user **must** read the user manual completely before operating the V20 Platform. It is important that you read the **safety** chapter carefully read before you attempt to install or operate the system.

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User interface
On/Off switch on power inlet
RF release Trigger on RF H/P (V-FR, V-ST, V- Form), and on Foot Pedal when RF H/P (V-FR, V- ST, V-Form/FC) is connected
Light release Trigger on V-IPL, and on Foot Pedal when V-IPL H/P is connected
GUI - Graphical User Interface:
LCD
Size: 10.4"
Resolution: 800X600
Communication: LVDS 3.4.3.4.1.4 Back light: LED 3.4.3 Touch method: Resistive
Controller:
CPU: Freescale IMX53
Memory: Micro SD card 8 GB

Main CPU
Working Voltage: 3.3V
Consumption: 71 mA @ 3V
Clock: external 25 MHz

External software plug
Type: USB Dongle
Use: Technician, Production, SW update (for GUI and MAIN)
Memory: 1 GB

RF Board
Inductor Flat
Inout 48V
Output: up to

Pulse	10-
Pre-Pulse	80msec
RF	0.8. 1.7.
Calibrati	via SW

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Fans

	Back Fan:	
	Dimensions:	80X80X25mm.
	CFM:	37
	Operating Voltage:	12V
	Consumption:	1.4W
RF Fan:	Noise:	30dBA
	Type:	Heat Sink mounted
	Fan Dimensions:	80X80X25mm
	CFM:	27.16
	Operating Voltage:	12V
	Consumption:	1.32W
	Noise:	<27.4dbA

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3.1.1. Water System Specifications

Water tank

Materials:	PMMA
Dimension:	160x120x140mm
Capacity:	2000cc

Pump

Operating Voltage:	24V
Pressure:	30 PSI
Controlled:	By flow sensor

Flow Sensor

Materials:	Acrylic, Brass, 316SS
Close On-Off:	Differential
Connections:	1/4" female NPT standard
Close On-Off:	Differential
Minimum Flow rate:	1 liter \ minutes

Radiator

Weld Per:	AA8381
Pressure test:	150 PSI N2 Under water
Apply Label:	Mark date code and serial number
Materials:	Corrugated copper fin, Stainless steel 316l tube 0.71mm
Cover:	1.5 mm alum. 5052-H32
End Plates:	2.0 mm stainless steel 304

Fan of Radiator

Dimensions:	120X120X38mm
CFM:	151.90
Rated Voltage:	24V
Operating Voltage:	12V ~ 27.6V
Consumption:	12W
Noise:	47.50dBA
Rate:	2,800 RPM

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Dimension:	Diameter 50mm, Length 175mm
Connections:	1/4" female NPT standard
Max Operating Temp:	55°C
Max Operating Pressure:	6 Bar
Water Filter:	50 Micron

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Type: EEPROM
Size: 1kbit

Flash lamp

Type: Xenon
Cooled By: Dionzed water
Fill pressure: 700-1000 TORR
ID / OD: 6 mm / 7 mm

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Light Guide

Material: Sapphire
Surface Area: 15x16 mm², 40x16 mm²
Surfaces quality: 80-50 per MIL-O-13830B.
Surface AR coated with: <2% reflectivity at 415 nm
• 1200 nm spectral range.

Reflector

Material: Alum 6061
Surface roughness before coating: 80-50
Surface Finish: Metal or dielectric reflecting Coating per MIL-C-5541DClass 2, with reflectivity higher than 92% at 450-1200 nm wavelength range

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For V-IPL H/P: 415-1200 nm

3.1.11. Switches

The V20 console has an ON/OFF switch (on the inlet)

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The V20 has an optional footswitch, with the same functionality as the trigger for V-FR, V-ST and V-Form H/Ps.

V-IPL H/P has a trigger. Combining with Footswitch is must for releasing Light pulses.

3.1.12. Indicators

The V20 is provided with an LCD Graphic display, showing the different parameters and operation modes of the system.

The V20 system has a speaker for audio indications.

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3.3. Electrical Requirements

The V20 Platform is preset to accommodate the local line voltage according to the customer's order. Accordingly, the system will require a separate line supply of one of the following (appropriate to the customer's country):

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3.5. Space Requirements

The V20 system is designed for installation in a clinical environment. The base dimensions of the system are 112 cm (H) x 47cm (W) x 67cm (D). A clearance of 30 cm must be available in all directions to permit appropriate ventilation.

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Warning

- It is forbidden to connect any third party equipment to the system without

3.7. System Safety Features

The V20 Platform incorporates the following safety features:

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When working with the V-IPL handpiece, you must press both the foot pedal and trigger to emit energy.

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3.9. Patient and Personnel Safety

When handling the V20 Platform, the primary concern should be the safety of both operator and patient. Carefully read the following safety guidelines:

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All personnel in the treatment room during treatment must know how to shut down V20 in an emergency.

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3.10. Remote Interlock Connector

V20 has an Interlock connector that enables you to connect to an external Interlock at the entrance to the treatment room. When this option is used, the Interlock prevents V20 from operating when the entrance door is open.

3.11. Foot Pedal

V20 is supplied with a foot pedal to emit pulses on V-ST, V-FR and V-Form handpieces

In order to emit pulses with the V-IPL handpieces, you need to press both the trigger and the foot pedal simultaneously.

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The packaging protects the V20 Platform from breakage or damage according to the company's packaging instructions.

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3.14.3 RF – V-FR Safety

The V-FR handpiece will not operate until a V-FR tip is connected and recognized by V20.

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- Only personnel essential to the session and who are fully acquainted with the required safety procedures should be allowed entry into the room during treatment.
- All personnel in the treatment room during treatment must know all of the required safety procedures.
- Never look directly into the V-IPL handpiece even when wearing protective eyewear.
- Never aim the V-IPL handpiece at anything other than the treatment area.
- Never hold the V-IPL handpiece close to your eyes or those of the patient.
- Do not expose any skin to the light pulse except the test patch and the treatment area.
- When using the V-IPL handpiece apply ultrasound gel to the treatment area.

Warnings

Anyone working with light pulses must be aware of the possible dangers. During treatment, do not allow anyone who is not essential to the treatment into the treatment room.

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Warnings

- Always wear protective eyewear when applying a light-based technology.
- Different protective eyewear is indicated for use with V-IPL handpieces. Make sure to use the appropriate type.
- The maximum ocular hazard distance (OHD) for the V-IPL handpiece is 10 meter.

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Handpiece	Wavelength (nm)	Optical Density for Patient	Optical Density for operator and
V-IPL	415-1200	4	3

Table 3.1: Protective eyewear

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The V20 System meets the requirements of the following standards:

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g. EN 60601-1-6 (2010) - Medical Electrical Equipment - Part 1-6: General requirements for basic safety and essential performance - collateral standard: usability. Uncontrolled release of RF or Light pulse, due to SW corruption, shall be prevented by a designated safety mechanism in the main board. SW-activated safety mechanisms are described in the V20 SWRS.

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The V20 plastic will be manufactured in the following technologies:

- d. High pressure Molds.
- e. CNC

The V20 metal sheet will be manufactured in the following technologies:

- f. Punch
- g. Laser Cut.

Assembly process

The V20 system will be assembled according to the relevant procedure developed by Viora.

Testing Process

RF board, Switch Module, Main board, UI Board, Applicators Board shall be tested in a designated tester/Jig for each board. Assembly using customary techniques. Water system liquid shall be tested in a designated tester/Jig. Assembly using customary techniques. Each system will be tested under the Acceptance Test Procedure.

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	Protective eyewear	1	VTR9301-1
	Patient goggles	1	VTR9301-2
	Filter Cover	1	MF42049-0

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The V20 Platform was carefully packed to ensure the safe delivery and arrival of the system and its components.

Make sure that you have all of the contents. Refer to -Equipment List.

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Figure 5.5: Connecting the Interlock

5.5 Connecting the Foot Pedal

Connect the foot pedal cable to the connector on the rear panel


Figure 5.6: Connecting the Foot Pedal

Fasten the screws to secure it into place.

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- Take the water filling bag out of the water filling kit.
- Take the “S” shape hook, and hang the bag on the system’s handle.

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- On the rear panel, open the flap of the water tank.
- Connect the pipe to the system’s water tank connection.
- Make sure the pipe is properly connected by pulling it and verify it is stays connected.

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5.8 Connecting the Handpieces

Warning

Make sure that power of the V20 Platform is turned off whenever you connect or disconnect any of the handpieces.

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Figure 5.15: Connecting the V-IPL Handpieces

Figure 5.16: Connecting the V- IPL Handpieces

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5.8.3 Circulating Water in V-IPL Handpieces

Note

- User shall operate the water circulation with V-IPL handpieces before start working with these handpieces.
 - User shall check water level indicator and make sure water level is above 70% of its capacity before starting the water circulation operation.
-

- Connect the handpiece (V-IPL) to the system.
- Turn the system on, and enter system setup → Maintenance → Water Circulation Screen.
- Press Start to begin water circulation.
- Follow the screen and wait until the system finishes circulation.

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screen. To turn the system off:

Step 1. Press the Login/Logout button on the Main menu screen. Step 2. Press the ON/OFF button located at the back of the system.

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6.2 Login/Logout

The login screen protects against use of V20 by unauthorized personnel. When treatment is completed, press the logout button to prevent unauthorized users from activating the system.

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6.6 Main Menu

The main menu identifies which handpieces are connected to V20 and fades out the disconnected handpieces.

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6.9.2 ST Basic Treatment Screen


Figure 6.8: ST Basic Treatment

- A. Mode** - The system has four different RF treatment modes that can be controlled by the operator. The modes are based on the following integrated technologies:

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- B. **Energy** - To control the energy level, press the energy (+) or (-) keys. The smaller (+) / (-) keys will change energy by 2J. The large (+) / (-) keys will change energy by 10J.
- C. **Hour counter** - Counts total treatments time usage on the handpiece, per hours.
- D. **Minute counter** - Counts time usage for current treatment. Press the  key to zero the counter. If a pulse is not released within 3 minutes, the counter will automatically pause and the system will move to 'Standby' mode.
- E. **Cooling** - The RF electrodes are cooled during treatment. The cooling icon  indicates that cooling is activated. If a pulse is not released within 60 seconds, cooling will stop and the system will move to Standby mode. In order to resume treatment, press the  key.
- F. **Home** - Pressing the  key, the operator opens the main menu.
- G. **Menu** - Pressing the  key, the operator opens the treatment menu.
- H. **Basic / Advanced Mode** - The operator's current mode of operation is seen in the the colored key with the arrow facing down (basic or advanced  Basic > Advanced ).
- I. **Standby / Ready** -  indicates that the system is ready for work and pulses can be emitted by pressing the trigger on the handpiece or using the foot pedal. The system will move to  mode if:
 - Operator changes mode of treatment
 - Operator changes energy by 10J
 - The pulse interval or duration in 'advanced' mode is changed
 - A pulse is not emitted for at least 3 minutes

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- J. **Pulse Duration** -  Determines the length of the pulse during which energy is released ("On" time).
- K. **Pulse Interval** -  Determines the pause between pulses ("Off" time).

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6.10.2 FR Basic Treatment Screen



Figure 6.12: FR Basic Treatment

- A. **Program** – Indicates the selected program. Program may be changed through the treatment screen as well as through the treatment menu screen.
- B. **Energy** - Controls the desired RF energy using the energy (+) / (-) keys.
- C. **Tip pulse counter** - Counts total pulses count with inserted tip.
- D. **Pulse counter** - Counts emitted pulses for current treatment. Press the key to zero the counter. If a pulse is not released within 3 minutes, the counter will automatically pause and the system will move to 'Standby' mode.
- E. **Cooling** - Indicates whether the cooling feature is enabled. The operator may switch the cooling feature on and off as desired by pressing the cooling key.

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- F. **Home** - Pressing the  key opens the main menu.
- G. **Menu** - Pressing the  key opens the treatment menu.
- H. **Basic / Advanced Mode** - The operator's current mode of operation is seen in the colored key with the arrow facing down (basic or advanced ).
- I. **Standby / Ready** -  indicates that the system is ready for work and pulses can be emitted by pressing the trigger on the handpiece or foot pedal. System will move to  mode if:
- Operator changes energy of treatment.
 - A pulse is not emitted for at least 3 minutes.

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6.11.1 BC or FC Treatment


Figure 6.16: FC Treatment

- A. **Mode** - The system has four different RF treatment modes that can be controlled by the operator. The modes are based on the following integrated technologies:

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
Figure 6.20: V-IPL Setting

- A. **Skin Type** – Choose the skin type according to the Fitzpatrick
- Q. **Hair Color** – Choose the hair color according to the descriptions on screen.
- R. **Hair Thickness** – Choose the hair thickness.
- S. **Home** - Pressing the  key opens the main menu.
- T. **Change Filter** - Pressing the  key opens the filter confirmation screen.
- U. **Setting** - The colored key with the arrow facing down shows the operator's current filter in use.
- V. **Treat** - Pressing the  key opens the treatment screen.
- W. **Totals Pulse Counter** - Counts the pulse usage for the total treatments on the handpiece.

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
- H. **Home** - Pressing the  key opens the main menu.
- I. **Setting** - Pressing the  **Settings** key opens the setting menu.
- J. **Filter** - The colored key with arrow facing down shows the current filter in use.
- K. **Total duration bar** – indicates the total pulse duration (on and off time) only in multiple pulse.

L **Standby / Ready** -  indicates that the system is ready to be used and pulses can be emitted by pressing the trigger on the handpiece **and** stepping

on the foot pedal. The system will move to  mode if:

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- **Consumable time:** The system is set to remind the operator of the handpiece running time. A reminder message will appear on the screen alerting the operator of the remaining time or usage on the system and will refer the user to the system information screen for more details.
- **Cooling:** The system will give an alert when the electrodes should be cooled down during the ST treatment. A progress bar will indicate when the cooling process is complete.
- **Contact Cooling:** During IPL treatments, the system will give an alert when the handpieces should be cooled down during the treatment. The message will disappear when the cooling process is completed.
- **Water system cooling:** During IPL treatments, the system will give an alert when the water system should be cooled down during the treatment. The message will disappear once the cooling process is completed.
- **V-FR tip warnings**
 - Treatment tips should only be used for a **single treatment** application on a single patient to avoid cross contamination.
 - System will alert the operator if an **empty tip** is inserted.

- Every 20 pulses the operator should stop treatment and **clean** the V-FR tip of any debris that may affect the treatment. System will go into 'Standby' mode and the trigger/ foot pedal will be disabled until the "OK" key is pressed. To clean the tip, use a stiff-bristled brush dipped in alcohol. Make sure the tip is completely dry before resuming treatment.

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7.7.3 Test Procedure

Before receiving a patient for a V20 treatment session, a test procedure must be performed during the initial consult meeting. The test procedure should be performed on a non-exposed site on or near the treatment area. If more than one area is to be treated, a test should be performed on all treatment sites.

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- For facial treatments, perform the test on concealed areas such as the back of the neck. For body treatments perform the test on the target area.
- Choose the relevant treatment program (Shallow, Medium or Deep).
- Set the energy to 0.25 J.
- Apply one pulse. If no heat sensation was reported by the patient, increase the energy level by 0.25 J and apply an additional pulse in an adjacent spot.
- Continue in the same manner until the patient reports a noticeable heat sensation. Note the energy level you have reached (XJ)
- Release additional 2-3 pulses at the same energy (x J)
- Increase the energy level by 0.25J (X+0.25J) and release 2-3 pulses
- Increase the energy level again by 0.25J (X+0.50J) and release additional 2-3 pulses
- Never release pulses in the same spot
- Wait 10 minutes to allow a pattern to appear. The desired end point is local erythema and edema.
- If the skin does not show signs of erythema and edema after 10 minutes, increase the RF energy by an additional 0.25J (X+0.75J) and release an additional pulse adjacent to the location of the previous pulse.
- Wait an additional 10 minutes to allow the end point to appear.
- Treatment area must be examined 1-2 days (Fitzpatrick I–III) and 5-7 days (Fitzpatrick IV–VI) after test procedure. A grid-like pattern matching the electrode pins should be visible on the treated area.

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- Make sure that the tips are cleaned every 20 pulses while the system is in *Standby* mode. A note will appear on the screen every 20 pulses.
 - Undamaged – ensure that no pins are missing or damaged.
 - Examine the treated area for the desired end point; erythema and edema on the treatment area.
 - Edema will normally appear within several minutes after the treatment.
 - At the end of the treatment, examine the area to make sure there are no untreated patches of skin. If such patches exist follow the below instructions:
 - Set the system to *Standby* Mode.
 - Clean the untreated patches and the tip with a dry lint-free gauze pad.
 - Pass over the untreated patches and avoid treated areas.
 - At the end of every treatment session remove the treatment tip and clean the handpiece according to the cleaning guidelines. Insert the V-FR tip cover when the system is not in use in order to protect the handpiece.
-

Warning

- **Do not** treat over the eyelids or the lips.
 - **Do not** treat different patients with the same tip, in order to prevent cross- contamination.
 - **Do not** treat with an expired tip.
 - **Do not** use the V-FR Tip if the sterile package has been damaged
 - Note that in darker skin types erythema is not easily noticeable; therefore, the end point of edema will be more prominent in these patients.
 - Before cleaning the tip make sure to set the system to *Standby* mode
-

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7.8.4 Treatment Guidelines

1. Make sure there is **full contact** between the electrodes and the skin without putting unnecessary pressure on the applicator.
 2. Pulses are released when the trigger on the applicator is pressed down.
 3. If the trigger is held down, the system will continue to release consecutive pulses until the press button is released.
 4. Once the pulse and vacuum stop slide the applicator to the next spot on the treatment area.
 5. Overlap no more than 10% of the pervious pulse.
-

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7.8.6 Post-Treatment Care

If patient notes excessive heat sensation or strong erythema, the area may be cooled with the following: chilled Aloe Vera gel, cold (not frozen) packs or cold compresses.

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- History of Keloid scarring.
- Diabetes.
- Koebnerization of skin, such as Vitilgo (Hypopigmentation) and psoriasis.
- Treatment area exclusion criteria:
 - Inflammation at the treatment area
 - Tattoos
 - Moles, nevi
 - Permanent make up
 - Silicone implants in the treatment area.
 - Any inflammatory skin condition, active infection or herpes simplex in the treatment area.
- Effective long term hair removal with light based devices requires the presence of a hair follicle for photothermolysis to commence. Waxing, plucking, electrolysis, or any other method of long term hair removal that affects the follicle should be avoided prior to and in between photoepilation treatments.
- Patients should be warned of potential harm that may arise when treating postinflammatory hyperpigmentation (PIH).

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8.7.3 Treatment Parameters

8.7.3.1 Step-by-Step Operation

- Insert the required filter for the specific treatment into the designated slot in the V-IPL handpiece
- Connect the V-IPL handpiece to the system.
- Select the treatment from the main menu.
- Confirm the selected filter as it appears in the "Filter Confirmation" screen
- Choose the parameters per patient's characteristics. The selected parameters will be highlighted.

Note

The operator must choose all parameters before proceeding to the next screen. 'Treat' key will be activated only after choosing all parameters.

- Press "Treat".
- Review the parameters in the "Treatment Screen" and change if required.
You may change the pulse structure to 'Rapid' or 'Multiple' as necessary, per type of treatment.
- Change to 'Ready' mode and start the treatment.
- The LED located at the back of the handpiece, as well as on the system, will light up as an indication of the released pulse.
- In Standby mode, the LED light will be OFF.
- Upon entrance to the Ready Mode, a beeping sound is heard to indicate that the system is ready to emit a pulse.
- In Ready mode, the LED light will blink.
- Once a pulse has been emitted a beeping sound is heard to indicate the pulse is complete. The handpiece can then be safely moved to an adjacent spot on the treatment area.
- For every pulse released, the number of pulses indicated on the display screen will increase by one.
- When the treatment is completed, place the V-IPL handpiece in its cradle.

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8.8.2 Treatment Parameters

8.8.3.2 Step-by-Step Operation

- Connect the required spot to the Nd:YAG handpiece
- Connect the handpiece to the system
- Confirm the selected spot size appears in the “Spot Confirmation” screen.
- Choose the treatment screen, review the duration and fluence. Change

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Skin Rejuvenation

- Number of Treatments: 4-6

- Treatment Intervals: 2-4 weeks
- Maintenance: 1 treatment every 3-6 months, or as required

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The V20 system requires standard maintenance and cleaning procedures to ensure optimal system performance. However, note that the system should be technically serviced only by authorized service technicians.

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9.2 Cleaning the V-ST, V-Form Electrodes

V- ST, V-Form electrodes must be thoroughly cleaned after every treatment.

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9.3.1 Cleaning the Filters

Warning

- Filter should be constantly cleaned before and after every treatment.
 - Damaged or broken filters should be replaced.
 - System should be turned off before cleaning the filters
-

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Warning

- Never wash the 6.4cm² spot in water to avoid water penetrating the handpiece
 - System should be turned off before cleaning the 6.4cm² spot
-

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The V20 system is equipped with a self-testing software routine that continuously monitors system operation. The following troubleshooting guide does not attempt to list all possible system failures.

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Message on screen	Probable Cause	Corrective Action	Problem solved	Continuedaction
Bad coupling	Bad coupling between the electrodes and the skin (V-FR, V-ST, V-Form & V-	Make sure that all electrodes are in full contact with the skin	Yes	Resume work
			No	Contact the local distributor
Overheating	System is overheated	Wait a few minutes until system has cooled down and the message disappears from the screen.	Yes	Resume work
			No	Contact the local distributor
No Handpiece Connected	Improper connection of the handpiece	Disconnect and reconnect the handpiece	Yes	Resume work
			No	Contact the local distributor
No V-Form applicator	Only V-Form handpiece is connected to	Connect the	Yes	Resume work

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Filter is Disconnected	The filter is not connected to the V- IPL handpiece	Connect a Filter	Yes	Go To Filter confirmation
			No	1. Connect another Filter 2. Contact the local distributor
Spot size has been changed	When Nd:YAG laser Treatment is selected and spot	Confirm the change	Yes	Go to Setting screen
			No	NA
Spot is Disconnected	When Nd:YAG laser Treatment is selected from the main menu and there is no spot connected	Connect a Spot	Yes	Go to Spot Confirmation
			No	- Connect another spot - Contact the local distributor
Foot Pedal is Disconnected	Foot Pedal is disconnected in V- Laser and V- IPL handpieces	Connect the foot pedal	Yes	Resume work
			No	Contact the local distributor
Foot Pedal Pressed	If the pedal is pressed and user want to proceed to Ready Mode (all handpieces) If the pedal is pressed more than 15 second in V-FR, V-ST and V- Form handpieces	Release the Foot Pedal	Yes	Resume work
			No	1. Reconnect the Pedal (Confirm that the system is off) 2. Contact the local distributor
Handpiece trigger is Pressed	If the trigger is pressed and user want to go Ready Mode (all handpieces except	Release the Trigger	Yes	Resume work
			No	Contact the local distributor
Interlock is disconnected	Interlock is disconnected in V- Laser and V- IPL handpieces	Connect the Interlock	Yes	Resume work
			No	Contact the local distributor
Configuration or license code is incorrect	User tried to insert Incorrect code	Verify the code is entered correctly	Yes	Resume work
			No	Contact the local distributor
Configuration or	User tried to insert	Enter a new code	Yes	Resume work

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Message on screen	Corrective Action	Problem solved	Continuedaction
System is not turned on	A. Ensure proper connection to mains outlet B. Make sure the emergency switch is not pressed	Yes	Resume work
		No	Contact the local distributor
No sound	Makesurethevolume setting on the setup screen is not set to minimum.	Yes	Resume work
		No	Contact the local distributor
Noisy system	Reboot the system.	Yes	Resume work

		No	Contact the local distributor
Water leak from the system	Contact the local distributor	-	-
		-	•
Weak Vacuum in V- Form applicators (BC or FC) - Vacuum is not working properly	Replace the filter in the BC or FC applicator	Yes	Resume work
		No	Contact the local distributor

Table 10.3: Operation troubleshooting

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These technologies have revolutionized non-invasive treatments to improve control and increase opportunity for practitioners, while enhancing patient experience. In addition, this in- house expertise has been at the forefront of implementing new developments in response to market demand, improving upon treatments using existing technologies such as laser and IPL. By formulating exclusive combination protocols and combining different proprietary technologies, Viora achieves enhanced clinical results and greater customer and patient satisfaction.

APPENDIX A – EMC Tables

Table 1 – Guidance and MANUFACTURER’S declaration – ELECTROMAGNETIC EMISSIONS –
for all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacturer’s declaration – electromagnetic emissions		
The V20 is intended for use in the electromagnetic environment specified below. The customer or the user of the V20 should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The V20 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	The V20 is suitable for use in all establishments other than domestic, and may be used in domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, provided the following warning is heeded: Warning: This equipment/system is intended for use by healthcare professionals only. This equipment/system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the [ME EQUIPMENT or ME SYSTEM] or
Harmonic emissions IEC 61000-3-2	Compliance	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

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Guidance and manufacturer’s declaration –electromagnetic immunity			
The V20 is intended for use in the electromagnetic environment specified below. The customer or the user of the V20 should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_T (>95 % dip in U_T) for 0,5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 sec	>95 % dip for 10 ms 60 % dip for 100 ms 30 % dip for 500 ms 95 % dip for 5000 ms	Mains power quality should be that of a typical commercial or hospital environment. If the user of the V20 requires continued operation during power mains interruptions, it is recommended that the V20 be powered from an uninterruptible power

			$\left[\frac{1}{d}\right] \sqrt{\frac{P}{f}}$$\left[\frac{1}{d}\right] \sqrt{\frac{P}{f}}$ 800 MHz to 2,5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the V20 is used exceeds the applicable RF compliance level above, the [ME EQUIPMENT or ME SYSTEM] should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the V20. ^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than $[V/1]$ V/m.			

Table 4 – Recommended separation distances between portable and mobile RF communications equipment and the ME EQUIPMENT or ME SYSTEM – for ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING			
Recommended separation distances between portable and mobile RF communications equipment and the V20			
The V20 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the V20 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the V20 as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter [W]	Separation distance according to frequency of transmitter [m]		
	150 kHz to 80 MHz $\left[\frac{1}{\sqrt{P}} \right]$	80 MHz to 800 MHz $\left[\frac{1}{\sqrt{P}} \right]$	800 MHz to 2,5 GHz $\left[\frac{1}{\sqrt{P}} \right]$
0,01	0.12	0.12	0.23
0,1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.69	3.69	7.38
100	11.67	11.67	23.33
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

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The V20 Fractional treatment involves use of disposable tips composed of:

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Operation documentation of the manufacturer for medical device:

V₃₀

CE₀₃₄₄

Copyright © 2014 Viora. All rights reserved.

Viora reserves the right to make changes to its products or specifications to improve performance, reliability or manufacturability. Information furnished by Viora is believed to be accurate and reliable. However, Viora assumes no responsibility for its use. No license is granted by its implication or otherwise under any patent or patent rights of Viora.

No part of this document may be produced or transmitted in any form or by any means, electronic or mechanical, for any purpose, without the express written permission of Viora.

Viora has patents and pending patent applications, trademarks, copyrights, or other intellectual property rights covering subject matter in this document. The furnishing of this document does not give you any license to these patents, trademarks, copyrights, or other intellectual property rights except as expressly provided in any written agreement from Viora.

Specifications are subject to change without notice.

P/N: MK-124 EN F

Release Date: October 2014

Warning

- Before using the system, check local regulations. If any local legislation is violated, use cannot be authorized.
- In the United States, federal law restricts prescription medical devices to sale by or on the order of a physician, or properly licensed practitioner. Viora makes no representations regarding federal, state, or local laws or regulations that might apply to the use and operation of this system.

Manufacturer's Address

Viora Ltd.
3 Maskit Street
Herzliya, Israel 46733
Tel: +972 9955 1344
Fax: +972 9955 1345
Website:
www.vioramed.com E-mail:
info@vioramed.com

Manufacturer's Authorized Representative

Viora UK Ltd.
450 Bath Road,
Longford, Heathrow,
UB7 0EB, UK
Tel: +44 (0)208 757 8750
Fax: +44 (0)208 7578767

This device complies with the European Medical Device Directive 93/42/EEC

Table of Contents

1.	Before You Start	9
1.1.	Conventions Used in the Manual.....	9
1.2.	Explanation of the Symbols Used in the System.....	11
2.	Introduction	13
2.1.	Regulating data for the construction of a medical device	14
2.1.1.	System Console	17
2.1.2.	Handpieces.....	17
2.1.3.	Disposable Sterilized Tip and Cover – V-FR Handpiece	17
2.1.4.	Filters and Spot for V-IPL	17
2.1.5.	Spots for V-Laser Nd:YAG	17
2.1.6.	Display Unit.....	18
2.2.	System Specifications	19
2.3.	Handpiece Specifications	19
2.3.1.	RF – V-ST Specifications	20
2.3.2.	RF – V-FR Specifications	20
2.3.3.	RF Specifications	20
2.3.4.	RF - V-Form - BC Applicator Specifications	20
2.3.5.	RF - V-Form - FC Applicator Specifications.....	20
2.3.6.	V-IPL Specifications	21
2.3.7	V-Laser Nd:YAG Specifications	21
3.	Safety.....	23
3.1.	Technical requirements for medical device	24
3.2.	System Classifications.....	40
3.3.	Electrical Requirements	40
3.4.	Environmental Requirements.....	41
3.5.	Space Requirements	41
3.6.	Safety When Operating the System.....	41
3.7.	System Safety Features	42

3.8.	Potential Fire Hazards	42
3.9.	Patient and Personnel Safety	43
3.10.	Remote Interlock Connector	43
3.11.	Foot Pedal	43
3.12.	Safety in Maintaining the System.....	43
3.13.	Transport and Storage Conditions.....	44
3.14.	RF Based Technology Safety.....	44
3.14.1.	RF – V-ST Safety.....	44
3.14.2.	V-Form Handpiece Safety.....	44
3.14.3.	RF – V-FR Safety	44
3.15.	Light based Technology Safety	45
3.15.1	Protective Eyewear	45
3.16.	Data for development of medical device.....	46
3.17.	Data for the manufacturing of medical device.....	47
4.	Labels.....	48
4.1.	System Labels.....	48
4.2.	Handpiece Labels.....	49
4.2.1.	RF – V-ST Handpiece and box	49
4.2.2.	RF – V-FR Handpiece and box	49
4.2.3.	RF – V-Form Handpiece and box.....	50
4.2.4.	V-IPL Handpiece and box	50
4.2.5.	V-Laser Nd:YAG Handpiece	51
5.	System Installation	52
5.1.	The using data for medical device.....	52
5.2.	Unpacking.....	53
5.3.	Assembling the Cradle.....	55
5.4.	Connecting the Interlock	55
5.5.	Connecting the Foot Pedal.....	56
5.6.	Filling the Water Container	57
5.7.	Connecting the Power Cable.....	58
5.8.	Connecting the Handpieces	58
5.8.1.	Connecting the V-FR, V-ST, V-Form Handpieces	58
5.8.2.	Connecting the V-IPL & V-Laser Nd:YAG Handpieces.....	60

5.8.3.	Circulating Water in V-IPL & V-Laser Nd:YAG Handpieces.....	62
5.9.	Wheel Brakes.....	63
5.10	Moving the System within the Facility.....	63
5.11.	Moving the System to another Facility.....	63
6	The data for operating medical device	64
6.1.	System On/Off	64
6.2.	Login/Logout	65
6.3.	Set time and Date	66
6.4.	Add Configuration.....	67
6.5.	Adding a New License	68
6.6.	Main Menu.....	69
6.7.	Set Up Menu.....	70
6.8.	Treatment Tips Screen	72
6.9.	RF – ST Treatment Screens	73
6.9.1.	ST Treatment Menu.....	73
6.9.2.	ST Basic Treatment Screen	73
6.9.3.	ST Advanced Treatment Screen	75
6.9.4.	V-ST Output Power	76
6.10.	FR Treatment Screens	77
6.10.1.	FR Treatment Menu	77
6.10.2.	FR Basic Treatment Screen	77
6.10.3.	FR Advanced Treatment Screen	78
6.10.4.	V-FR Output Power.....	79
6.11.	BC/FC Treatment Screens	81
6.11.1	BC or FC Treatment Screen	81
6.12.	V-IPL Treatment Screens.....	84
6.12.1.	Filter Confirmation Screen.....	84
6.12.2	V-IPL Setting Screen	86
6.12.3.	V-IPL Treatment Screen.....	87
6.13.	V-Laser Nd:YAG Treatment Screens.....	90
6.13.1.	Spot Confirmation Screen	90
6.14.	V-Laser Nd:YAG Treatment Screen.....	91
6.15.	System Messages.....	93
7.	Clinical Guidelines for RF Based Technology.....	95

7.1.	Contraindications	95
7.2.	Possible Side Effects	96
7.3.	Realistic Expectations	96
7.4.	Consult Meeting.....	97
7.5.	Treatment Follow-Up.....	97
7.6.	V-ST Clinical Guide	98
7.6.1.	Indications of Use	98
7.6.2.	Test Procedure	98
7.6.3.	Pre-Treatment Preparation.....	98
7.6.4.	Treatment Guidelines	99
7.6.5.	Treatment Schedule.....	99
7.6.6.	Post-Treatment Care.....	99
7.7.	V-FR Clinical Guide.....	99
7.7.1.	Indications of use.....	99
7.7.2.	Pre-Treatment Preparation.....	100
7.7.3.	Test Procedure	100
7.7.4.	Treatment Guidelines	102
7.7.5.	Treatment Schedule	103
7.7.6.	Post-Treatment Care.....	103
7.8.	FC and BC Clinical Guide	104
7.8.1.	Indications of use.....	104
7.8.2.	Pre-Treatment Preparation.....	104
7.8.3.	Test Procedure	104
7.8.4.	Treatment Guidelines	107
7.8.5.	Treatment Schedule	107
7.8.6.	Post-Treatment Care.....	108
8.	Clinical Guidelines for Light Based Technology	108
8.1.	Contraindications	108
8.2.	Possible Side Effects	110
8.3.	Realistic Expectations	111
8.4.	Consult Meeting.....	111
8.5.	Test Procedure	112
8.6.	Treatment Follow-Up.....	112
8.7.	V-IPL Clinical Guide	113

8.7.1.	Indications of use.....	113
8.7.2.	Pre-Treatment Preparation.....	113
8.7.3.	Treatment Parameters.....	114
8.7.4.	Treatment Schedule.....	117
8.7.5.	Post-Treatment Care.....	117
8.8.	V-Laser Nd:YAG Clinical Guide.....	118
8.8.1	Indications of use.....	118
8.8.2	Pre-Treatment Preparation.....	119
8.8.3.	Treatment Parameters.....	119
8.8.4.	Treatment Schedule.....	120
8.8.5.	Post-Treatment Care.....	121
9.	Data for the technical maintenance.....	122
9.1.	Cleaning the System.....	122
9.2.	Cleaning the V-ST, V-Form Electrodes.....	122
9.3.	Cleaning the V-IPL Handpiece	122
9.3.1.	Cleaning the Filters.....	123
9.3.2.	Cleaning the Spots	123
9.4.	Cleaning the V-Laser Nd:YAG Handpiece	124
9.4.1.	Cleaning the Spots	124
9.5.	Replacing the filters in the BC/FC applicators.....	124
9.6.	Replacing the Water Filter and Deionized Water	126
10.	Data for the service.....	126
10.1.	Fault Messages	126
10.2.	Failure Messages.....	126
10.3	Operation troubleshooting	129
11.	About Viora.....	130
12.	Data for recycling.....	137
13.	System composition.....	137

1. Before You Start

Warning

- Read this manual to become familiar with all safety requirements and operating procedures before attempting to operate the system.
 - Any intense radio frequency (RF) device can cause injury if used improperly.
 - Any light based device can cause injury if used improperly.
 - High voltage is present inside the system. Always be aware of the possible dangers and take proper safeguards as described in this manual.
 - The system and its components must be serviced only by qualified personnel.
 - Only a trained person can upgrade the software version.
 - To avoid risk of electric shock, this equipment must only be connected to a supply main with protective earth.
 - No modification of this equipment is allowed.
-

1.1. Conventions Used in the Manual

The following conventions in the form of notes, cautions and warning are used in this manual:

Warning

A **Warning** alerts the operator to the possibility of injury, death, or serious adverse reactions associated with the use or misuse of the system.

Caution

A **Caution** alerts the operator to the possibility of a problem with the device associated with its use or misuse. Such problems include system malfunction, failure, and damage to the system or other property. The caution statement includes a precautionary action that should be taken to avoid the hazard.

Note

The content of this **Note** offers general information that is important to keep in mind.

Note

In many places throughout this manual you will see words or phrases in the text which are bolded in a sans-serif typeface, **such as this**. These are words and phrases that you will always see on the LCD display of the control panel.

1.2.Explanation of the Symbols Used in the System

Symbol	Description
	Attention
	CE Compliance
	CSA Compliance Symbol (244308 CSA master contract number, project 2137900)
	Dispose of according to WEEE regulations
	Fus
	Type BF Equipment
	Type B Equipment
	High frequency isolated patient circuit
	Manufacturer
	Non-ionizing electromagnetic radiation
	CAUTION: Federal law restricts this device to sale by or on the order of a physician
	Refer to instruction manual

	Danger: Invisible Laser Radiation
	LS EQUIPMENT - HAZARD symbol
	Emergency Stop
	Authorized Representative in the European Union
 10.08.2012	Discard Product After Expiration Date

2. Introduction

V30 is a light based and Radiofrequency (RF) Non-Surgical medical aesthetic platform system.

The V30 platform, developed by Viora, utilizes radiofrequency (RF) for aesthetic and dermatological procedures, with applications for the face, neck, décolleté and hands. The system is based on the well-established technology of bipolar RF, which is extensively used for the coagulation and ablation of soft tissues. Over the past decade, RF technology has become widely used in the aesthetics and dermatology field for skin tightening and resurfacing.

The platform utilizes V-IPL for a variety of aesthetic skin treatments such as long-term hair reduction, skin photorejuvenation, treatment of vascular lesions and solar lentigos and treatment of mild to moderate inflammatory acne on the face, neck, décolleté and body.

The advanced V-Laser Nd:YAG is used for vascular lesions, long-term hair removal and nail fungi clearance (Onychomycosis).

The V30 is in Class IIb according to Directive 93/42/EEC of 14 June 1993, amended by Council Directive 2007/47/EC of 5 September 2007, concerning medical devices. Annex IX, Rule 9 states "all active therapeutic devices intended to administer or exchange energy are in Class IIb".

Note

Only personnel who have undergone and completed a thorough training course are permitted to operate the system.

The intended users are aesthetic physicians or dermatologists. Eligibility of users shall be in accordance with local market regulations.

2.1. Regulating data for the construction of a medical device

Figure 2.1: System Console (Front)

Figure 2.2: System Console (Side)


Figure 2.3: Connector Cover

Figure 2.4: V-ST Handpiece

Figure 2.5: V-FR Hand Piece Tip

Figure 2.6: FR

Figure 2.7: V-FR Tip Cover

Figure 2.8: V-IPL Handpiece

Figure 2.9: V-IPL Filter

Figure 2.10: V-IPL Small Spot


Figure 2.11: V-IPL Filter Cover


Figure 2.12: V-Laser Nd:YAG Handpiece

Figure 2.13: V-Laser Nd:YAG Spots

**Figure 2.15: V-Form
Handpiece**

**Figure 2.16:
BC Applicator**

**Figure 2.17:
FC Applicator**

2.1.1 System Console

The system console houses the following components:

- Front connector for V-IPL and V-Laser Nd:YAG handpieces
- Rear connector for V-FR, V-ST, V-Form and V-Glide handpieces
- RF Generator
- Switch Module
- Water Cooling system
- Display Unit
- An Emergency Button that completely stops V30 operation
- LED indicator used for V-IPL and V-Laser Nd:YAG treatments
- Service panel

2.1.2 Handpieces

The system can accommodate four handpieces as follows:

1. RF V-ST
2. RF V-FR
3. V-Form
4. V-Glide
5. V-Laser Nd:YAG Long Pulse 1064nm
6. V-IPL

The V-FR, V-ST, V-Form and V-Glide handpieces are attached with a cable to the connector located on the rear panel of the main console.

The V-IPL and V-Laser Nd:YAG 1064nm handpieces are attached with a cable to the connector on the front panel of the main console.

2.1.3 Disposable Sterilized Tip and Cover – V-FR Handpiece

A disposable tip, containing a matrix of bi-polar electrode pins attaches to the end of the handpiece and is placed on the skin during treatment. The tip form is 13 rows and 8 columns of pins.


Figure 2.18: V-FR Tip Matrix

A thermoelectrically cooled bracket surrounding the grid area is designed to cool surrounding skin and reduce pain during treatment.

A V-FR tip cover is used when the V-FR handpiece is not in use to protect it.

2.1.4 Filters and Spot for V-IPL

The V-IPL handpiece allows the operator to change the optical filters. The system will

recognize the specific filter automatically.

Five (5) optical filters in different wavelength are used.

V-IPL handpiece has 2 options of spot size: the large spot is 6.4cm^2 and the small spot is 2.4cm^2 . The material of the light guide is Sapphire cooled by two TEC (Thermo Electro Cooling) units. The operator can change the spot size as necessary, according to the treatment area. Spot size will be automatically recognized by the system.

Figur 2.19: V-IPL

Filter

Figure 2.20: V-IPL

2.4 cm² Spot

2.1.5 Spots for V-Laser Nd:YAG

The V-Laser Nd:YAG handpiece allows the operator to change the spot (Light guide) which the system can automatically recognize.

Four (4) spot size options are available.

The material of the Light Guide is Sapphire cooled by two TEC (Thermo Electro Cooling) units.

Figure 2.21: V-Laser Nd:YAG Spot

2.1.6 Display Unit

The display unit - touch screen primarily presents the following information:

- 1. Treatment type
- 2. Treatment settings
- 3. Pulse / Time counter
- 4. Treatment time

2.2 System Specifications

Supported Handpieces	RF – V-ST
	RF – V-FR
	V-Form
	V-Glide
	V-IPL
	V-Laser Nd: YAG
Cooling System	DI Water cooling system
Maximum Power	1000 watts
Electrical requirements	100-240 V~ 10A 50/60 Hz 100-120 V~ 10A 50/60 Hz Single Phase
Environmental Requirements	Operating temperature 18°-24°C (64°-75°F) and relative humidity from 10% to 80%, non- condensing
User Interface	10.4 inch touch screen
Weight	70 kg
Dimensions	112 (H) x 47 (W) x 67 (D) cm

2.3 Handpiece Specifications

2.3.1 RF – V-ST Specifications

Treatment area	Frequency Mode	Frequency (Hz)	Energy (J/cm ³)	Pulse Duration
8mm x 8mm	Mode 1	0.8MHz	Up to 130 J/cm ³	50-200 ms
	Mode 2	1.7MHz		
	Mode 3	2.45MHz		
	Mode 4	0.8MHz, 1.7MHz, 2.45MHz		

2.3.2 RF – V-FR Specifications

Treatment area	Frequency (nm)	Energy/Power (Joule)	Pulse Duration (ms)
7.3mm x 21.8mm	1 MHz	0-10 J	10-100 ms

2.3.3 RF Specifications

2.3.4 RF - V-Form - BC Applicator Specifications

Treatment area	Frequency mode	Frequency (Hz)	Power (W)	Switching Pulse Duration (ms)
35.2 cm ³	Mode 1	0.8MHz	Up to 50 W	200 msec
	Mode 2	1.7MHz		
	Mode 3	2.45MHz		
	Mode 4	0.8MHz, 1.7MHz, 2.45MHz		

2.3.5 RF - V-Form - FC Applicator Specifications

Treatment area	Frequency mode	Frequency (Hz)	Power (W)	Switching Pulse Duration (ms)
0.637 cm ²	Mode 1	0.8MHz	Up to 25 W	200 msec
	Mode 2	1.7MHz		
	Mode 3	2.45MHz		
	Mode 4	0.8MHz, 1.7MHz, 2.45MHz		

2.3.6 V-IPL Specifications

Application	Filter	Fitzpatrick Type	Spot Size (cm ²)	Wavelength (nm)	Pulse Structure	Fluence (J/cm ²)	Pulse Duration (ms)	Pulse Interval (ms)	Repetition Rate
Hair removal	HR - 630	I - V	2.4, 6.4	630 - 1200	Single	5-25	20-50	-	0.4
					Multiple	5-25	5-15	15-75	0.4
					Rapid	5-8	7	-	2
	HR - 570	I - IV	2.4, 6.4	570 - 1200	Single	5-26	20-50	-	0.4
					Multiple	5-26	7-15	15-80	0.4
					Rapid	5-8	7	-	2
Skin rejuvenation	SR - 580	I - V	2.4, 6.4	580 - 1200	Single	5-28	10-30	-	0.4
					Multiple	5-28	4-8	5-40	0.4
					Rapid	5-8	7	-	2
Vascular and pigmented lesions	VP - 530	I - V	2.4, 6.4	530 - 1200	Single	5-30	10-30	-	0.4
					Multiple	5-30	4-8	5-40	0.4
					Rapid	5-8	7	-	2
Acne clearance	Acne - 415	I - V	2.4, 6.4	415 - 1200	Single	5-10	15-50	-	0.4
					Multiple	5-10	4-12	5-20	0.4
					Rapid	5-8	7	-	2

Note

Lifecycle of the V-IPL handpiece is 100,000 pulses

2.3.7 V-Laser Nd:YAG Specifications

Application	Spot Size (mm)	Wavelength (nm)	Fluence (J/cm ²)	Pulse Structure	Pulse Duration (ms)	Pulse Interval (ms)	Repetition Rate
Hair removal	9.5	1064	20-50	Single	10-25	-	Up to 1 Hz
				Multiple	5-10	5-15	
Vascular	2x 4 3.5	1064	20-500	Single	10-25	-	Up to 1 Hz
Onychomycosis	3.5 5.5			Multiple	5-10	5-15	

Note

Lifecycle of the Nd:YAG handpiece is 200.000 pulses

3 Safety

To ensure dependable, consistent and safe treatments, the following safety measures must be taken into consideration when operating and maintaining the system. The person operating the system must be familiar with the safety regulations provided in this chapter. The user **must** read the user manual completely before operating the V30 Platform. It is important that you read the **safety** chapter carefully read before you attempt to install or operate the system.

If any sections of the safety chapter are unclear, or if any questions arise prior to, or during installation and/or operation of the handpiece, contact your local distributor.

Viora assumes no liability whatsoever for any damage or injury caused as a result of actions taken which are not in strict accordance with the instructions provided in this manual.

Improper use may result in injury and the primary consideration when operating the system must be to maximize safety for both patient and operator.

Warning

- It is important that the operators of this system familiarize themselves with all safety requirements and operating procedures prior to operating the system.
- Any intense RF device may cause injury if used improperly.
- Any light based device can cause injury if used improperly.
- Failure of HF equipment could result in an unintended increase of output power.
- High voltage is present inside the system.
- Always be aware of the possible dangers and take proper safeguards as described in Section 3.5 of this user manual.
- For complete contact information please refer to page 2 of this user manual.

3.1. Technical requirements for medical device

Type: IEC, Medical Grade, illuminated
Rating: 10A
Fuse: 2 X (250 VAC, 10A, 5X20mm)
Line Filter
Grade: Medical
Rating: 10A
Power Cable
Compliance with standard UL817
Power supply
Grade: Medical
Input Voltage: 90-254 VAC, 50-60Hz
Output Voltage: 12/12/24/24 V
Output Power: up to 400W

User interface
On/Off switch on power inlet
RF release Trigger on RF H/P (V-FR, V-ST, V- Form), and on Foot Pedal when RF H/P (V-FR, V- ST, V-Form/FC) is connected
Light release Trigger on V-IPL and V-Laser Nd:YAG H/P, and on Foot Pedal when V-IPL or V-Laser Nd:YAG H/P is connected
GUI - Graphical User Interface:
LCD
Size: 10.4'
Resolution: 800X600
Communication: LVDS 3.4.3.4.1.4 Back light: LED 3.4.3 Touch
method:Resistive
Controller:
CPU: Freescale IMX53
Memory: Micro SD card 8 GB

Main CPU
Working Voltage: 3.3V
Consumption: 71 mA @ 3V
Clock: external 25 MHz
External software plug
Type: USB Dongle
Use: Technician, Production, SW update (for GUI and MAIN)
Memory: 1 GB

RF Board
Inductor Flat
Inout 48V
Outout: up to

Pulse	10-
Pre-Pulse	80msec
RF	0.8. 1.7.
Calibrati	via SW

Capacitors Charger

Grade: Medical

Input Voltage: 90-254 VAC, 50-60Hz

Output Voltage: Up to 500 vDC

Output Power: up to 1000J

Capacitors Bank

Capacitance: 26mf

Voltage: 500 vDC

Switch Module

Input Voltage: 24 vDC

Switching by: IGBT 400 A

Simmer Trigger: up to 1.5 kV

Unload Resistor: 250 ohm / 200 W

Fan:	Dimensions:	80X80X15mm
	CFM:	32.8
	Operating Voltage:	24V
	Consumption:	2.76W
	Noise:	34dBA

Fans

Back Fan:	Dimensions:	80X80X25mm.
	CFM:	37
	Operating Voltage:	12V
	Consumption:	1.4W
	Noise:	30dBA

RF Fan:	Type:	Heat Sink mounted
	Fan Dimensions:	80X80X25mm
	CFM:	27.16
	Operating Voltage:	12V
	Consumption:	1.32W
	Noise:	<27.4dbA

V-FR V-ST and V-Form H/P Connectors

No. of:	1
Type:	Female (socket)
Type:	Circular

V-IPL and V-Laser Nd:YAG H/P Connectors

No. of:	1
Type:	Female (socket)
Type:	Circular

Materials

Chassis:	Galvanized Sheet
Covers and panels:	ABS
H/P cradle:	ABS PC and TPR80

3.1.1. Water System Specifications

Water tank

Materials: PMMA
Dimension: 160x120x140mm
Capacity: 2000cc

Pump

Operating Voltage: 24V
Pressure: 30 PSI
Controlled: By flow sensor

Flow Sensor

Materials: Acrylic, Brass, 316SS
Close On-Off: Differential
Connections: 1/4" female NPT standard
Close On-Off: Differential
Minimum Flow rate: 1 liter \ minutes

Radiator

Weld Per: AA8381
Pressure test: 150 PSI N2 Under water
Apply Label: Mark date code and serial number
Materials: Corrugated copper fin, Stainless steel 316l tube 0.71mm
Cover: 1.5 mm alum. 5052-H32
End Plates: 2.0 mm stainless steel 304

Fan of Radiator

Dimensions: 120X120X38mm
CFM: 151.90
Rated Voltage: 24V
Operating Voltage: 12V ~ 27.6V
Consumption: 12W
Noise: 47.50dBA
Rate: 2,800 RPM

System Temperature Sensor:

Temperature upper limit: 40°C

Deionizer Cartridge with Water Filter

Dimension:	Diameter 50mm, Length
175mm	
Connections:	1/4" female NPT standard
Max Operating Temp:	55°C
Max Operating Pressure:	6 Bar
Water Filter:	50 Micron

3.1.2. Cart Specifications

Wheels

No. of:	4
Type:	with brakes
Diameter:	100mm

3.1.3. V-ST H/P Specification

Operation parameters

See section 3.13

Materials

Shell material: ABS PC
Electrodes material: Electrolytic tough pitch Copper,
UNS C11000, Hard Crom plated

Cable

Type: Servo, Shielded
RF Wires: 2X20 AWG, twisted pair, shielded
Signal Wires: 4X18 AWG
Rated Voltage: 1000V
Tested: 4000V
Rating: UL

TEC

Function: Electrodes cooling
Size: 15x30X3.5mm
Power: 8W (1units)
Voltage: 4V
Current: 2A
Delta Temp: 67 Celsius max

Temperature Sensor

Range: -55° to 150° Celcius
Accuracy: ±0.25° Celcius

CPU unit

Working Voltage: 3V
Consumption: 250 mA @ 3V
Clock: internal 8 MHz

Cooling fan

Dimensions: 25X25X10mm.
CFM: 3
Operating Voltage: 12V
Consumption: 0.45W
Noise: 16dBA

Memory chip

Type: EEPROM
Communication: SPI
Size: 1K

Treatment Area: 8mmX8mm

3.1.4. V-FR H/P Specification

Operation parameters:

Materials

H/P Shell material: PC

Cable

Type:	Servo, Shielded
RF Wires:	2X20 AWG, twisted pair, shielded
Signal Wires:	4X18 AWG
Rated Voltage:	1000V
Tested:	4000V
Rating:	UL

Disposable Tip - Regular

Electrode Matrix	
No. of electrodes:	56 (7 X 8)
Electrode's pitch:	1.3mm
Single electrode diameter:	0.3mm
Area of treatment:	14X15.5 mm
PCB material:	FR-4 2.4 mm thick

Connection to H/P	
Connector style:	contact
Identification:	EEPROM
Shell Material:	Santoprene TPV 281-87 med
Sterilization:	EtO Gas.

Vacuum Pump

Type:	Diaphragm
Motor:	Brushless
Voltage:	12V

Flow Rate: 1.5 L/Min @ 0 mbar, 0.6 @ -300 mbar
Max Vacuum: -300 mbar
Consumption: ~220mA @ 12V

Solenoid

Function: build and release vacuum at Tip
Way: 3
Operation Voltage: 3V

TEC

Function: Cooling the skin via the electrodes' array
Size: 15X15 mm
Voltage: 3.8V
Current: 4A
DT max: 66 Celsius

CPU unit

Working Voltage: 3V
Consumption: 250 mA @ 3V
Clock: internal 8MHz

Fan

Dimensions: 30X30X10mm.
CFM: 4.7
Operating Voltage: 12V
Consumption: 1.0 W
Noise: <23dBA

Temperature Sensor

Range: -55° to 150° Celcius
Accuracy: ±0.25° Celcius

Memory chip

Type: EEPROM
Communication: SPI
Size: 1K

3.1.5. V-Form H/P Specification

Combine vacuum and RF technology. Handpiece with 2 applicators: BC and FC.

Materials

Shell material: PC-LEXAN 500 + 10% FG

Cable

Type: 20/4 24/4 Cable 80C 300V
RF wires: 2X20 AWG , twisted pair , shielded
Signal wires: 4X24 AWG
Rated Voltage: 600V
Tested: 4000Vrms
Rating: UL
12V wires: 2X20 AWG , singles

Connector

No. of: 1
Type: Male
Type: Circular

Electrodes

No. of: 4 electrodes for FC applicator & 6 electrodes for BC
Materials: ABS polilac, plating : chrome

CPU unit

Type : PIC24FJ128GB206
Memory chip: Flash
Working Voltage: 3.3V
Consumption: 250 m
Clock: External 8 MHz

HP screen

Type : OEL Display Module

Size: 1.27"
Resolution: 128X96

Pump:

Type: Univision Technology Inc Vacuum pump
Flow: 11 l/min
Motor type: DC brush
Nominal speed: up to 600 rpm
Wight: 190g

Solenoid

Configuration: 3W2P
Power: 2W

IR Thermometer

Voltage : 3V
Diameter : 4.7 mm
Accuracy of +/-0.1°C

Filter

Operation range: 0-100 Kpa
Filtering accuracy: 10 µm
Temperature Range: 0-60°C
Filter area: 2.8 cm²

3.1.6. V-IPL H/P Specification

Operation parameters

See section 3.13

Materials

Shell material: Lexan 500 10% GF

Umbilical Cable

Type: Servo, Shielded
Lamp Wires: 2X12 AWG shielded
TEC Wires: 2X18 AWG
Applicator Board Cable: 4X28 AWG shielded
Trigger Cable: 2X22 AWG shielded
Rated Voltage: 1000V
Tested: 4000V
Rating: UL
Pipe Hot Water:
Color: Red
Material: Polyurethane
Diameter: 6/4mm
Pipe Cool Water:
Color: Blue
Material: Polyurethane
Diameter: 6/4mm

TEC

Function: Light Guide cooling
Size: 30x30X3.5mm
Power: 48W (1 units)
Voltage: 12V
Max Current: 6A
Delta Temp: 67 Celsius max

Temperature Sensor

Range: -55° to 150° Celsius
Accuracy: ±0.25° Celcius

CPU unit

Working Voltage: 3.3V
Consumption: 20 mA @ 3V
Clock: internal 8 MHz

Memory chip

Type: EEPROM
Size: 1kbit

Flash lamp

Type: Xenon
Cooled By: Dionized water
Fill pressure: 700-1000 TORR
ID / OD: 6 mm / 7 mm

Optical Tube

Dimension: Ø 10X Ø 9X L = 64 mm
Material: Fused Silica
Clear Optical: High transparency at 415-1200 nm
Operation temperature: 20-150 C

Light Guide

Material: Sapphire
Surface Area: 15x16 mm², 40x16 mm²
Surfaces quality: 80-50 per MIL-O-13830B.
Surface AR coated with: <2% reflectivity at 415 nm - 1200 nm spectral range.

Reflector

Material: Alum 6061

Surface roughness before coating: 80-50

Surface Finish: Metal or dielectric reflecting Coating per MIL-C-5541D Class 2, with reflectivity higher than 92% at 450-1200 nm wavelength range

Durability: Per MIL-M-13508C, including humidity
Operating temperature: 20-150° C

Optical Filters

No. of applications: 4
No. of Filters: **5**
Dimension: 48x18 mm
Types:

Acne 415-1200nm
VP 530-1200nm
SR 580-1200nm
HR 570-1200nm
HR 630-1200nm

3.1.7. V-Laser Nd:YAG H/P Specification

Operation parameters

Materials

Shell material: Lexan 500 10% GF

Umbilical Cable

Type: Servo, Shielded
Lamp Wires: 2X12 AWG shielded
TEC Wires: 2X18 AWG
Applicator Board Cable: 4X28 AWG shielded
Trigger Cable: 2X22 AWG shielded
Rated Voltage: 1000V
Tested: 4000V
Rating: UL
Pipe Hot Water:

Color: Red
Material: Polyurethane
Diameter: 6/4mm

Pipe Cool Water: Color: Blue
Material: Polyurethane
Diameter: 6/4mm

TEC

Function: Light Guide cooling
Size: 15x20X3.5 mm

Power: 9W (1units)
Voltage: 3.6V
Max Current: 3.7A

Delta Temp: 65 Celsius max

Temperature Sensor

Range: -55° to 150° Celcius
Accuracy: ±0.25° Celcius

CPU unit

Working Voltage: 3.3V
Consumption: 20 mA @ 3V
Clock: internal 8 MHz

Memory chip

Type: EEPROM
Size: 1kbit

Flash lamp

Type: Xenon
Cooled By: Dionzed water
Fill pressure: 1000-1150 TORR
ID / OD: 5 mm / 7 mm

Light Guide

Material: Sapphire
No. of Spot: 4
Surfaces quality: 80-50 per MIL-O-13830B.
Coating: AR@1064 nm T=99.5% normal incidence angle.
Damage threshold: >50J/cm2 for pulse 10 msec.
Spot size:
2x4 mm²
3.5 mm (Diameter)
5.5 mm (Diameter)
9.5 mm (Diameter)

Laser Crystal

Material:	Nd doped
YTTRIUM ALUMINUM GARNET	
Concentration:	1.2% +0/-0.1%
Dimension:	110 mm, Ø 8 mm
Surfaces quality:	10/5
Transferred wave front distortion:	□/16 per inch @
633 mm	
Wavelength:	1064 nm
Transmission:	60% +/-3%

3.1.8. Foot Pedal

Type:	Momentary close
Cable Termination:	DB9 Male
Sealing Rating:	IPX2

3.1.9. Shelf Life and Life Time

Console life time:	5 years
V-FR Tip shelf life:	2 years

3.1.10. Operation Profile

Operation limitations

V-IPL H/P work pulses limit:	100,000
V-Laser Nd:YAG H/P work pulses limit:	200,000
Disposable (regular) Tip Pulse limit:	600 pulses

RF Frequencies of operation

V-FR H/P:	1 MHz
V-ST & V-Form H/Ps:	0.8 MHz, 1.7 MHz, 2.45 MHz

RF Power output

For V-ST H/P:	0-130 J/cm ³
For V-FR H/P:	0-10 Joule
For BC applicator:	up to 50 W
For FC applicator:	up to 25 W
Pre-Pulse ST:	8W
Pre-Pulse FR:	4W
V-FR pre heating (Smart Heating) pulse:	9W or 5 J/cm ³

RF pulse width:

For V-ST H/P:	50-200msec
For V-FR H/P:	10-100msec
For V-Form H/P:	200msec
Pre Pulse width for all RF H/Ps:	80msec
Pre heating pulse for FR H/P:	100msec

V-IPL and V-Laser Nd:YAG Repetition rate

V-IPL H/P: up to 0.4Hz

V-Laser Nd:YAG H/P: up to 1 Hz

V-IPL and V-Laser Nd:YAG Power output

For V-IPL H/P:	5-35 J/cm ²
For V-Laser Nd:YAG H/P:	5-500 J/cm ²

V-IPL and V-Laser Nd:YAG pulse width:

For V-IPL H/P:	4-50 msec
For V-Laser Nd:YAG H/P:	10-95 msec

V-IPL and V-Laser Nd:YAG wave length:

For V-IPL H/P:	415-1200 nm
For V-Laser Nd:YAG H/P:	1064 nm

3.1.11. Switches

The V30 console has an ON/OFF switch (on the inlet)

V-ST, V-FR and V-Form H/P have a trigger for releasing RF pulses.

The V30 has an optional footswitch, with the same functionality as the trigger for V-FR, V-ST and V-Form H/Ps.

V-IPL and V-Laser Nd:YAG H/P has a trigger. Combining with Footswitch is must for releasing Light pulses.

3.1.12. Indicators

The V30 is provided with an LCD Graphic display, showing the different parameters and operation modes of the system.

The V30 system has a speaker for audio indications.

The V-FR handpiece is provided with a LED indicator to indicate RF Pulse release.

The V-Form H/P is provided with a LCD H/P screen

3.1.13. Physical Dimensions

System dimensions: 112(H)x47(W)x67(D) Cm

3.2. System Classifications

- Electric shock protection: Class I, Type BF
- Not suitable for use in presence of flammable substances
- Protection against ingress of liquids: Ordinary equipment
- Mode of operation – continuous

3.3. Electrical Requirements

The V30 Platform is preset to accommodate the local line voltage according to the customer's order. Accordingly, the system will require a separate line supply of one of the following (appropriate to the customer's country):

100-240 V~ 50/60 Hz

100-120 V~ 50/60 Hz

Warning

- To protect against risk of fire, replace the fuse with the same type and rating only.
 - Make sure to use a power cable with proper grounding.
 - If there is no grounding, user should contact a certified electrician in order to add grounding.
-

3.4. Environmental Requirements

- Electronic parts can be damaged by corrosive materials. The system must be operated in a non-corrosive environment.
- For optimal operation of the system, maintain a room temperature between 18°-24°C (64°-75°F) and relative humidity from 10% to 80%, with no condensation.
- Atmospheric Pressure Range: 70kPa – 106kPa.
- If the system has been stored in an environment cooler than the above temperatures, allow the system to warm up to room temperature before first activating the system.
- Metallic dust is potentially harmful to electrical equipment.

3.5. Space Requirements

The V30 system is designed for installation in a clinical environment. The base dimensions of the system are 112 cm (H) x 47cm (W) x 67cm (D). A clearance of 30 cm must be available in all directions to permit appropriate ventilation.

3.6. Safety When Operating the System

Keep all system covers closed at all times and do not touch any of the system's inner components. Only authorized service technicians are permitted to remove the panels. Always make sure that the V-FR tip is replaced before a new treatment session begins. Under no circumstances, should a V-FR tip be used for more than one patient. Avoid operating the system in the presence of potentially explosive or flammable materials.

Always make sure that all operating personnel have been properly trained and are familiar with the system's controls and are aware of the proper steps on how to instantly shut down the system when required.

Only aesthetic physicians and dermatologists are allowed to operate the system.

No part of the system should be immersed in liquid.

Medical electrical equipment can be affected by, and experience interference from, portable and mobile RF communications equipment.

Only original and specified accessories and cables should be used. Using any other equipment may result in increased emissions or decreased durability of the system.

Warning

- It is forbidden to connect any third party equipment to the system without
-

3.7. System Safety Features

The V30 Platform incorporates the following safety features:

Default treatment parameters are automatically suggested when the operator selects a treatment.

The handpiece is equipped with a trigger. If the trigger is not pressed, no energy will be emitted.

When working with the V-Laser Nd:YAG and V-IPL handpiece, you must press both the foot pedal and trigger to emit energy.

When working with RF-based handpieces, the system automatically tests for proper coupling prior to any pulse by releasing a short sub pulse at low voltage. If coupling between electrodes and skin is insufficient the RF energy will be disabled.

The handpiece and treatment tip must be properly connected to the system in order for the treatment mode to be activated. When handpiece or tips are not well connected, the system will alert the operator.

A V-IPL pulse cannot be emitted unless the handpiece, filter and spot are all properly connected.

The system can be switched to a 'Standby' mode, disabling the release of energy.

When changing treatment parameters, the system will automatically move to 'Standby' mode.

An Emergency Button to completely stop the system's operation.

Warning

Objects such as superficial metal or conductive implants, rods, plates or pins in the treatment areas are contraindicated to the treatment.

In case of improper coupling using an RF handpiece, an alert is indicated on screen and by audio signal the trigger should be released immediately and the handpiece should be re-oriented on the skin. In rare cases, superficial burns may occur from the heat generated by the RF electrodes.

3.8. Potential Fire Hazards

Do not use the system in the presence of explosive or flammable materials.

Do not use flammable substances when preparing the skin for treatment.

If alcohol is used for disinfecting the hand piece, it must be completely dry before the system is used.

3.9. Patient and Personnel Safety

When handling the V30 Platform, the primary concern should be the safety of both operator and patient. Carefully read the following safety guidelines:

Before beginning any treatment course, the operator should complete a detailed patient history to verify that the patient is eligible for the treatment and is in healthy condition.

The operator should provide the patient with information about the treatment, the treatment protocol, expected results and any risks or side effects associated with the treatment.

All personnel in the treatment room during treatment must know how to shut down V30 in an emergency.

Only trained operators are allowed to perform treatments.

A test procedure must be performed before any treatment course and on any new anatomical area, in order to evaluate the treatment parameters.

Warning

- Improper treatment may cause local burns.
 - Untrained operators should not operate the system
-

3.10. Remote Interlock Connector

V30 has an Interlock connector that enables you to connect to an external Interlock at the entrance to the treatment room. When this option is used, the Interlock prevents V30 from operating when the entrance door is open.

3.11. Foot Pedal

V30 is supplied with a foot pedal to emit pulses on V-ST, V-FR and V-Form handpieces

In order to emit pulses with the V-Laser Nd:YAG and V-IPL handpieces, you need to press both the trigger and the foot pedal simultaneously.

Warning

- 3.11.1. When the foot pedal or trigger is stuck, operator should manually release the pedal or trigger.
 - 3.11.2. Move to 'standby' mode before connecting or disconnecting the foot pedal.
-

3.12. Safety in Maintaining the System

Only authorized service technicians are qualified to open the system's panels to service the system.

Shut down the system and disconnect the power cable before you perform any of the maintenance procedures in chapter System

3.13. Transport and Storage Conditions

The optimal transport and storage conditions of the system are:

Ambient temperature: -5°C to 65°C.

Relative humidity: from 10% to 80%, non-condensing.

The packaging protects the V30 Platform from breakage or damage according to the company's packaging instructions.

Atmospheric Pressure Range: 50kPa – 106kPa.

3.14. RF Based Technology Safety

Any intense RF device may cause injury if used improperly.

RF energy is emitted into the skin via bi-polar electrodes. After applying ultrasound gel (when using the V-ST handpiece) or glycerin oil (when using V-Form or V-Glide handpieces), full coupling of both electrodes to the skin is necessary before pressing the trigger button, to ensure the RF energy is safely conducted. If coupling with the skin is broken due to incomplete contact or removal of the handpiece from the treatment area while the trigger is pressed, a beeping sound will indicate improper coupling.

Warnings

If there is improper coupling, release the trigger immediately and re-orient the electrodes on the skin. In rare cases, superficial burns can occur from the heat generated by the RF electrodes.

3.14.1. RF – V-ST Safety

Apply ultrasound gel to the treatment area to ensure proper contact.

3.14.2. V-Form Handpiece Safety

When operating the handpieces apply glycerin oil to the treatment area to ensure proper contact.

3.14.3 RF – V-FR Safety

The V-FR handpiece will not operate until a V-FR tip is connected and recognized by V30.

Warnings

V-FR tip replacements should only be purchased through Viora's authorized local distributor.

V-FR Tips are sterile. Individual tips should be discarded if they come in contact with an unsterile surface after removed from package (floor, hands)

Treatment tips should only be used for a **single treatment** application on a single patient to avoid cross contamination.

3.15. Light based Technology Safety

Light-based handpieces emit lasers or light pulses that can cause injury if used incorrectly. Guidelines to ensure optical safety are as follows:

Only personnel essential to the session and who are fully acquainted with the required safety procedures should be allowed entry into the room during treatment.

All personnel in the treatment room during treatment must know all of the required safety procedures.

Never look directly into the V-IPL or V-Laser Nd:YAG handpiece even when wearing protective eyewear.

Never aim the V-IPL or V-Laser Nd:YAG handpiece at anything other than the treatment area.

Never hold the V-IPL or V-Laser Nd:YAG handpiece close to your eyes or those of the patient.

Do not expose any skin to the light pulse except the test patch and the treatment area.

When using the V-IPL handpiece apply ultrasound gel to the treatment area.

Warnings

Anyone working with light pulses must be aware of the possible dangers.

During treatment, do not allow anyone who is not essential to the treatment into the treatment room.

Reflective objects such as jewelry, watches and surgical instruments or mirrors should be kept away from the flash lamp to avoid any light reflections.

IR emitted from the device may cause eye injury. Avoid eye exposure.

After 100K pulses or 3 years, whichever comes first, a Short Wavelength Boundary Test for the V-IPL handpiece needs to be performed. The handpiece should not deviate from the Set Value by more than 5%. The user is responsible for performing the test

Caution

UV is emitted from this device. Eye or skin irritation may result. Use of controls or adjustments or performance of procedures other than those specified herein may result in hazardous radiation exposure

3.15.1. Protective Eyewear

Protective eyewear should be appropriate to wavelength when there is exposure to Laser or light pulses.

Warnings

- Always wear protective eyewear when applying a light-based technology.
- Different protective eyewear is indicated for use with V-IPL and V-Laser Nd:YAG handpieces. Make sure to use the appropriate type.
- Use protective eyewear for V-Laser Nd:YAG according to local regulation.
- The maximum nominal ocular hazard distance (NOHD) for the V-Laser Nd:YAG handpiece is 836 meter.
- The maximum ocular hazard distance (OHD) for the V-IPL handpiece is 10 meter.

Viora protective eyewear should be used by patient physicians and staff as follows:

Handpiece	Wavelength (nm)	Optical Density for Patient	Optical Density for operator and
V-IPL	415-1200	4	3
V-Laser Nd:YAG	1064	>7	>7

Table 3.1: Protective eyewear

If for some reason the patient cannot wear the protective eyewear, provide patients with eye protection that completely blocks light from their eyes.

3.16. Data for development of medical device

The V30 System meets the requirements of the following standards:

- IEC 60601-1:2005 + CORR.1 (2006) + CORR.2 (2007) - Medical electrical equipment - Part 1: General requirements for basic safety and essential performance
- IEC 60601-2-2 (2009) - Medical electrical equipment, Part 2-2: Particular requirements for the safety of high frequency surgical equipment.
- IEC60601-1-2 (2007) - Medical electrical equipment - Part 1-2: Collateral Standard: Electromagnetic compatibility – Requirements and tests.
- IEC60825-1 – (2008) - Safety of laser products - part 1: equipment classification and requirements.
- IEC 60601-2-22 (2007) - Medical electrical equipment - part 2-22: particular requirements for basic safety and essential performance of surgical, cosmetic, therapeutic and diagnostic laser equipment.
- IEC 60601-2-57 (2011) - Particular Requirements for the Basic Safety and Essential Performance of non-laser light source equipment.

g. EN 60601-1-6 (2010) - Medical Electrical Equipment - Part 1-6: General requirements for basic safety and essential performance - collateral standard: usability. Uncontrolled release of RF or Light pulse, due to SW corruption, shall be prevented by a designated safety mechanism in the main board. SW-activated safety mechanisms are described in the V30 SWRS.

3. 17. Data for the manufacturing of medical device

The system's PCBs, as assembled by the Sub-contractor, include use of the following technologies:

- a. Surface Mount Technology.
- b. Reflow soldering.
- c. Through hole.

The V30 plastic will be manufactured in the following technologies:

- d. High pressure Molds.
- e. CNC

The V30 metal sheet will be manufactured in the following technologies:

- f. Punch
- g. Laser Cut.

Assembly process

The V30 system will be assembled according to the relevant procedure developed by Viora.

Testing Process

RF board, Switch Module, Main board, UI Board, Applicators Board shall be tested in a designated tester/Jig for each board. Assembly using customary techniques.

Water system liquid shall be tested in a designated tester/Jig. Assembly using customary techniques. Each system will be tested under the Acceptance Test Procedure.

4 Labels

4.1 System Labels

Figure 4.1: System Labels


Figure 4.2: CAUTION Neutral Fusing

Figure 4.3: CAUTION V-Laser Nd:YAG

Figure 4.4: CAUTION V-IPL

4.2 Handpiece Labels

4.2.1 RF – V-ST Handpiece and box

Figure 4.5: V-ST Handpiece Label


Figure 4.6: V-ST Handpiece box Label

4.2.2 RF – V-FR Handpiece and box

Figure 4.7: V-FR Handpiece Label


Figure 4.8: V-FR Handpiece Box Label

4.2.3 RF – V-Form Handpiece and box

Figure 4.11: V-Form Handpiece Label

Figure 4.12: V-Form Handpiece Set Box Label

4.2.4 V-IPL Handpiece and box

Figure 4.13: V-IPL Handpiece Label


Figure 4.14: V-IPL Handpiece Set Box Label

Figure 4.15: Optical Radiation Aperture Label

4.2.5 V-Laser Nd:YAG Handpiece


Figure 4.16: V-Laser Nd:YAG Handpiece Label
Box Label

Figure 4.17: V-Laser Nd:YAG Handpiece

Figure 4.18: V-Laser Nd:YAG Aperture Label

5 System Installation

5.1 The using data for the medical device

Component	Description	Quantity	Part No.
General	System console	1	VTL40000
	RF and V-IPL cradle	1	AS42103
	V-Laser Nd:YAG and V- IPL cradle	1	AS42102
	Screws for cradles	8	SC40046
	Allen key	1	MS40037
	Interlock Plug	1	HR40053
	Foot pedal	1	AS40030
	Power cable	1	VRT3001-UK
	User manual	1	MK-124
	Water filling kit	1	AS40011
	Connector cover	1	BM40017
V-ST	Handpiece	1	V-ST2001
V-FR	Handpiece	1	V-FR2002
	V-FR tips	1 box (10 tips)	FRTIP1000-2
	V-FR Tip cover	1	VTO-1001-8
	Cleaning Brush	1	VMD9101-2
V-Form	Handpiece Set	1	AS45930
	Handpiece	1	AS45900
	FC applicator	1	AS45910
	BC applicator	1	AS45920
	Filters	24	BM45008
	Hex key	1	MS45001
V-IPL	Handpiece set	1	VTL43000
	Handpiece	1	VTL43010
	V-IPL Optic Set	1 case with 5 filters and 1 spot	VTL43020

V-Laser Nd:YAG	Protective eyewear	1	VTR9301-1
	Patient goggles	1	VTR9301-2
	Filter Cover	1	MF42049-0
	Handpiece set	1	VTL44000
	Handpiece	1	VTL44010
	V-Laser Nd:YAG Optic set	1 case with 4 spots	VTL44020
	Protective eyewear	1	VTLC4000
	Patient goggles	1	VTR9301-2
	Danger Hazard sticker	1	DC40000

5.2 Unpacking

Warning

Only Viora-certified personal can unpack the device.

Note

Please do not place heavy objects on top of the system's package.

The V30 Platform was carefully packed to ensure the safe delivery and arrival of the system and its components.

Make sure that you have all of the contents. Refer to 5.15.1 Equipment List.

Make sure that there are no visible signs of damage to the system or to the package contents. If you are missing some of the contents or there are visible signs of damage, contact your local distributor.

- Open the packaging box by releasing the harnesses
- Remove the foam on the side and bottom.
- Attach the triangle foam to the bottom of the packing (ramp) using the Velcro.


Figure 5.1: Attach triangle sponge

- d. Release the front wheel brakes by pushing the clips up and turn them to be facing the ramp.


Figure 5.2: Release front wheels

- e. Carefully roll out the system using the handle.


Figure 5.3: roll out the system

5.3 Assembling the Cradle

- Connect the cradle to the cradle holder on the side of the system.
- Use the Allen key to fasten the screws.


Figure 5.4: Assembling the Cradle

5.4 Connecting the Interlock

Insert the interlock in to the interlock port.


Figure 5.5: Connecting the Interlock

5.5 Connecting the Foot Pedal

- Connect the foot pedal cable to the connector on the rear panel


Figure 5.6: Connecting the Foot Pedal

- Fasten the screws to secure it into place.

Note

- In order to avoid a potential accident, please store the foot pedal in a safe place when not in use
- The foot pedal must not get wet. Do not keep any liquids in the vicinity of the foot pedal.

5.6 Filling the Water Container

Warning

Use only deionized water when filling the water tank.

- Take the water filling bag out of the water filling kit.
- Take the "S" shape hook, and hang the bag on the system's handle.


Figure 5.7: Hanging the water filling bag

- On the rear panel, open the flap of the water tank.
- Connect the pipe to the system's water tank connection.
- Make sure the pipe is properly connected by pulling it and verify it is stays connected.


Figure 5.8: Connecting the pipe to the water tank

- Pour deionized water into the water bug to the maximum and wait until water meter on the rear panel is between the minimum and maximum, above 70% of its capacity.
- Remove the pipe and close the flap.

5.7 Connecting the Power Cable

- Connect the power cable to the power connection port on the system's rear panel.


Figure 5.9: Connecting the Power Cable

Plug the electrical cable into an appropriate electrical outlet. Refer to chapter 3.2 Electrical Requirements on Safety.

5.8 Connecting the Handpieces

Warning

Make sure that power of the V30 Platform is turned off whenever you connect or disconnect any of the handpieces.

5.8.1 Connecting the V-FR, V-ST, V-Form and V-Glide Handpieces

- Remove the handpiece from the package.
- Place the handpiece on the cradle.
- Connect the handpiece connector at the end of the cable to the connector on the rear panel by rotating it until it is securely locked.

Figure 5.10: Connecting the RF handpieces

5.8.1.1 Connecting the V-FR Tip

- To connect the V-FR tip to the handpiece, open the supplied box and remove one sterile V-FR tip.
- Peel off the seal and place the tip on the handpiece, without removing the tip from the blister package. Tips can be connected directly to the handpiece.

Figure 5.11: Connecting the V-FR Tip

5.8.1.2 Connecting the V-FR Tip Cover

The system is supplied with a black V-FR tip cover. When the V-FR handpiece is not in use, use the cover to protect it.


Figure 5.12: Connecting the V-FR Tip cover

5.8.1.3 Connecting the BC & FC Applicators On V-Form

- To connect the BC/FC applicator to the handpiece, align the applicator cam to the handpiece groove and push the applicator into the handpiece until locking sound is heard.
- To release the applicator from the handpiece, simply press the release button on the applicator and remove the applicator away from the handpiece.

Figure 5.13: Connecting the FC Applicator

Figure 5.14: Releasing the FC Applicator

5.8.2 Connecting the V-IPL & V-Laser Nd:YAG Handpieces

- Remove the handpiece from the package.
- Place the handpiece in the cradle.
- Connect the handpiece connector at the end of the cable to the connector on the front panel by rotating it until it is securely locked.

Note

Connector cover is supplied in order to protect the connector. In the first installation, please remove before connecting the V-IPL or V-Laser Nd:YAG handpieces

Figure 5.15: Connecting the V-Laser Nd:YAG /V- IPL Handpieces

Figure 5.16: Connecting the V-Laser Nd:YAG /V- IPL Handpieces

5.8.2.1 Connecting V-IPL Filter

Warning

Before changing the filter, make sure the system is on 'Standby' Mode.

- Remove the filter from the package.
- Clean the filter as described in Section 9.3.1 of this user manual.
- Insert the filter into the designated slot on the front of the handpiece.
- Make sure the filter is inserted securely.

Note

When not in use, filter cover should be inserted into the handpiece in order to protect the opening.

5.8.2.2 Connecting the 2.4cm² V-IPL Spot

Warning

When turning the system on make sure it has recognized the small spot as displayed on the top menu bar.

- Remove the 2.4cm² spot from the package.
- Clean the spot as described in Section 9.3.2 of this user manual.

- Connect the small spot to the permanent spot.
- Make sure the spot is securely connected.

5.8.2.3 Connecting V-Laser Nd:YAG Spots

- Remove a spot from the package.
- Clean the spot as described in Section 9.4.1 of this user manual.
- Connect the spot to the front of the handpiece.
- Make sure the spot is securely connected.


Figure 5.17: Connecting a V-Laser Nd:YAG Spot

5.8.3 Circulating Water in V-IPL & V-Laser Nd:YAG Handpieces

Note

- User shall operate the water circulation with V-IPL and V-Laser handpieces before start working with these handpieces.
- User shall check water level indicator and make sure water level is above 70% of its capacity before starting the water circulation operation.

- Connect the handpiece (V-IPL or V-Laser Nd:YAG) to the system.
- Turn the system on, and enter system setup → Maintenance → Water Circulation Screen.
- Press Start to begin water circulation.
- Follow the screen and wait until the system finishes circulation.

5.9 Wheel Brakes

The system is equipped with wheel brakes. To complete installation – please make sure that the wheel brakes are secured and locked.

5.10 Moving the System within the Facility

- Make sure the system is turned off and the power cable is disconnected.
- The system is equipped with wheel brakes. Make sure they are released to allow movement.
- Use the handle located at the back of the system to carefully move it around.

5.11 Moving the System to another Facility

Warning

Moving the system to another facility shall be done only by Viora certified personal.

- Connect the drainage pipe to the water bag.
- Connect the drainage pipe to the drainage plug on the rear panel of the system.
- Turn the system on, and enter system setup → Maintenance → Water Drainage Screen.
- Press Start to begin water drainage.
- Follow the screen and wait until the system finishes drainage.
- Turn off the system and remove the drainage pipe from the system.
- The system is equipped with wheel brakes. Make sure they are released to allow movement.
- Use the handle located at the back of the system to carefully move it around.
- Slowly load the system onto the vehicle. Make sure the system is secured in place and proper padding is used to protect the system while traveling.
- To unload the system in the new facility, carefully remove it from the vehicle and use the handle on the back of the system to place it in its new location.
- Reinstall the system.

Note

It is recommended that the system only be moved when it is packaged in its original case and coverings.

6 The data for the operating medical device

Warning

The operator should not use the handpiece if dropped on the floor

6.1 System On/Off

To turn the system on:

- Press the ON/OFF button located at the rear panel of the system.
- Following start up, the system will display the login

screen. To turn the system off:

- Step 1. Press the Login/Logout button on the Main menu screen. Step 2. Press the ON/OFF button located at the back of the system.

Note

In an emergency, press the Emergency Stop button.

6.2 Login/Logout

The login screen protects against use of V30 by unauthorized personnel. When treatment is completed, press the logout button to prevent unauthorized users from activating the system.

Note

It is forbidden to reveal the V30 password to any unauthorized person.


Figure 6.1: Login Screen

- Enter the login code supplied by your local distributor. The default code is 1234. Press the 'backspace' button  if you enter a wrong digit.
- Press **Login**. The system will verify the code and display the main menu screen.

6.3 Set time and Date


Figure 6.2 Set Time and Date screen

- A. Press the day button
- B. Press the minus/plus button until you reach the desired number.
- C. Repeat the same action with month/Year/Hour/Minutes buttons. Once date and hour are set, press OK button to confirm.

6.4 Add Configuration

The configuration codes are generated by Viora and provided to the local distributor upon the purchase of each system.

A configuration code determines the configuration of the actual device which then reflects the licensing codes to be provided to the user.


Figure 6.3: Add Configuration

- A. Type the 30 digits configuration code as provided to you by the local distributor.
- B. Press the OK button to confirm.

6.5 Adding a New License

The Licensing code determines the allowance of usage per application. Each license allows use for a specific application, for a pre-determined period of time and/or application usage.


Figure 6.4: Add License

- A. Type the 30 digits license code provided to you by the local distributor.
- B. Press the OK button to confirm.

6.6 Main Menu

The main menu identifies which handpieces are connected to V30 and fades out the disconnected handpieces.

Press the treatment selection key for the desired type of treatment.


Figure 6.5: Main Menu Screen

Note

Please make sure that the handpiece in use corresponds with the displayed treatment screen.

6.7 Set Up Menu

- **Treatment preferences** - Allows the operator to set a default mode for ST or FR treatment screen. The operator may always switch between basic and advanced mode through the treatment screen or by resetting the default preferences.
- **System information** - Shows indications for:
 - Serial numbers of the system and connected handpieces
 - Software versions
 - Usage data for connected handpieces and tips
- **Event History** - Shows history of system events such as self-test failures.
- **Treatment History** - Shows details about past treatments:
 - Treatment date and time
 - Duration
 - Handpiece used
 - Total pulses used
 - Last energy used
 - Last duration used
- **License** - - Shows details about:
 - Active licenses - Allows the operator to view all active licenses on the system.
 - Licenses History - Shows details (indication) about all the licenses that were on the systems (including expired codes).
 - Add License allows operator to add a new license.
 - Configuration History - Shows details of previous configuration.
 - Add Configuration- add a new configuration code

> Active Licenses

> Licenses History

>> Add License

> Configuration History

> Add Configuration

Figure 6.6: Software licensing configuration, through Set-up Menu

- **Maintenance** - Shows details about:
 - Water Drainage - Allows operator to drainage the water from the system.
 - Water Circulation - Allows operator to circulate water in system and in V-IPL and V-Laser Nd:YAG handpieces
 - Water Filter Status – Shows next Water Filter replacement date.
 - Reset Water Filter Status – Reset the next Water Filter replacement date.

6.8 Treatment Tips Screen

- From the treatment menu, press the  key to access treatment tips.
- Press the key for the desired category.
- Press **Close** to return to the treatment menu.


Figure 6.7: Treatment Tips Screen

6.9 RF – ST Treatment Screens

6.9.1 ST Treatment Menu

The ST Treatment menu offers four different treatment programs. One program is for body treatment and three are for facial treatments (forehead, cheeks, neck & décolleté). Choose the desired treatment area by pressing the touch screen.

6.9.2 ST Basic Treatment Screen


Figure 6.8: ST Basic Treatment

A. **Mode** - The system has four different RF treatment modes that can be controlled by the operator. The modes are based on the following integrated technologies:

Mode 1– Low RF frequency of 0.8 MHz (suitable for deep treatments).


Mode 2 – Medium RF frequency of 1.7 MHz (suitable for mid-layer treatments). Mode 3 – High RF frequency of 2.45 MHz (suitable for superficial treatments). Mode 4 – Combines all three frequencies for multi-layer treatments.

To choose the desired treatment mode, press the mode number

Note

Changing modes will move the system to **Standby**. Press the **Ready / Standby** key to resume work.

- B. **Energy** - To control the energy level, press the energy (+) or (-) keys. The smaller (+) / (-) keys will change energy by 2J. The large (+) / (-) keys will change energy by 10J.
- C. **Hour counter** - Counts total treatments time usage on the handpiece, per hours.
- D. **Minute counter** - Counts time usage for current treatment. Press the  key to zero the counter. If a pulse is not released within 3 minutes, the counter will automatically pause and the system will move to 'Standby' mode.
- E. **Cooling** - The RF electrodes are cooled during treatment. The cooling icon  indicates that cooling is activated. If a pulse is not released within 60 seconds, cooling will stop and the system will move to Standby mode. In order to resume treatment, press the  key.
- F. **Home** - Pressing the  key, the operator opens the main menu.
- G. **Menu** - Pressing the  key, the operator opens the treatment menu.
- H. **Basic / Advanced Mode** - The operator's current mode of operation is seen in the the colored key with the arrow facing down (basic or advanced  Basic  Advanced).
- I. **Standby / Ready** -  indicates that the system is ready for work and pulses can be emitted by pressing the trigger on the handpiece or using the foot pedal. The system will move to  mode if:
- Operator changes mode of treatment
 - Operator changes energy by 10J
 - The pulse interval or duration in 'advanced' mode is changed
 - A pulse is not emitted for at least 3 minutes

To resume work, press the **Standby** key. Wait for the cooling to restart (indicated by the progress bar).

To release an RF pulse, press the trigger located on the handpiece OR the foot pedal.

Note

- Pressing the trigger or foot pedal continuously, will emit 5 continuous RF pulses followed by a warning screen. Operator must press **OK** in order to resume work.
- Pressing the trigger or foot pedal accidentally with no contact on the skin, will cause a bad coupling alert on the screen.

Note

Follow the treatment movements as per the arrow indications on the image.

6.9.3 ST Advanced Treatment Screen


Figure 6.9: ST Advanced Treatment Screen

The advanced treatment screen is for experienced users and allows control of additional treatment parameters.

- J. **Pulse Duration** - Determines the length of the pulse during which energy is released ("On" time).
- K. **Pulse Interval** - Determines the pause between pulses ("Off" time).

Note

- In the advanced mode, pressing the trigger or foot pedal continuously, will emit 5 consequent RF pulses followed by a warning screen. Operator must press 'OK' in order to resume work.
- Pressing the trigger or foot pedal accidentally with no contact on the skin, will cause bad coupling alert on the screen.

Warning

Changing pulse duration is directly related to amount of energy released onto the treatment area.

6.9.4 V-ST Output Power


Figure 6.10: V-ST Output Power


Figure 6.11: V-ST Power Stability

6.10 FR Treatment Screens

6.10.1 FR Treatment Menu

The FR Treatment menu offers three different treatment programs (shallow, medium and deep). Upon choosing the treatment program, the relevant treatment screen opens.

6.10.2 FR Basic Treatment Screen


Figure 6.12: FR Basic Treatment

- A. **Program** – Indicates the selected program. Program may be changed through the treatment screen as well as through the treatment menu screen.
- B. **Energy** - Controls the desired RF energy using the energy (+) / (-) keys.
- C. **Tip pulse counter** - Counts total pulses count with inserted tip.
- D. **Pulse counter** - Counts emitted pulses for current treatment. Press the  key to zero the counter. If a pulse is not released within 3 minutes, the counter will automatically pause and the system will move to 'Standby' mode.
- E. **Cooling** - Indicates whether the cooling feature is enabled. The operator may switch the cooling feature on and off as desired by pressing the cooling key.

Note

Cooling reduces pain and makes the treatment more comfortable for the patient; therefore it is recommended to be used at all times, and as per treatment protocol.

- F. **Home** - Pressing the  key opens the main menu.
- G. **Menu** - Pressing the  **FR Menu** key opens the treatment menu.
- H. **Basic / Advanced Mode** - The operator's current mode of operation is seen in the the colored key with the arrow facing down (basic or advanced  **Basic**  **Advanced**).
- I. **Standby / Ready** -  indicates that the system is ready for work and pulses can be emitted by pressing the trigger on the handpiece or foot pedal. System will move to  mode if:
- Operator changes energy of treatment.
 - A pulse is not emitted for at least 3 minutes.

To resume work, press the **Standby** key.

To release an RF pulse, press the trigger located on the handpiece **or** the foot pedal.

6.10.3 FR Advanced Treatment Screen


Figure 6.13: FR Advance Treatment Screen

The advanced treatment screen is for experienced users and allows control of additional treatment parameters.

- J. **Pulse Duration** -  Determines the length of the pulse during which energy is released ("On" time).
- K. **Vacuum** -  Indicates whether vacuum feature is enabled. Operator may switch the vacuum feature on and off as desired by pressing the vacuum key.

Note

Vacuum ensures optimal coupling between skin and electrodes and is recommended for use at all times, per treatment protocols.

- L. **Smart Heat** -  Provides volumetric heating by releasing a low energy pulse for a more efficient fractional RF delivery due to lower skin impedance. Operator may switch the smart heat feature on and off as desired by pressing the Smart Heat key.

6.10.4 V-FR Output Power


Figure 6.14: V-FR Output Power


Figure 6.15: V-FR Power Stability

6.11 BC/FC Treatment Screens

6.11.1 BC or FC Treatment


Figure 6.16: FC Treatment

A. **Mode** - The system has four different RF treatment modes that can be controlled by the operator. The modes are based on the following integrated technologies:

Mode 1 – Low RF frequency of 0.8MHz (suitable for deep treatments).

Mode 2 – Medium RF frequency of 1.7 MHz (suitable for mid-layer

treatments). Mode 3 – High RF frequency of 2.45 MHz (suitable for superficial treatments). Mode 4 – Combines all three frequencies for multi layer treatments.

To choose the desired treatment mode, press the mode number

Note

Changing modes will move the system to **Standby**. Press the **Ready / Standby** key to resume work.

- B. **RF Energy level** -Press the desire number 1-4 while 1 is the lowest energy level and 4 is the highest level
- C. **Vacuum level** - Press the desire number 1-4 while 1 is the lowest vacuum level , 4 is the highest level and off is no vacuum .
- D. **Countdown clock** –press the clock button to insert the desire treatment time after setting the time the system will show the countdown left until the end of treatment
- E. **Temperature indicator** - displays the real time measured temperature
- F. **Applicator indicator**
- G. **Home** - Pressing the  key, operator opens the main menu.
- H. **Standby / Ready** -  indicates that the system is ready for work and pulses can be emitted by pressing the trigger on the handpiece or the foot pedal. System will move to  mode if:
 - Operator changes mode of treatment
 - Operator changes energy level
 - Operator changes vacuum level
 - Pulse is not emitted for 3 minutes
 or longer To resume work, press the **Standby** key.

To release an RF pulse, press the trigger located on the handpiece OR the foot pedal.


Figure 6.17: ST set Countdown clock screen

- I. **Delete button**
- J. **Set clock minutes** - 0-9 button to set the treatment duration
- K. **Confirmation button** – pressing on the OK button will set the time and back to treatment screen
- L. **Cancel button** - back to treatment screen

Note

- Pressing the trigger or foot pedal continuously, will emit 5 consequent RF pulses followed by a warning screen. Operator must press **OK** in order to resume work.
- Pressing the trigger or foot pedal accidentally with no contact on the skin, will cause bad coupling alert on the screen.


M. **Mode** - The system has four different RF treatment modes that can be controlled by the operator. The modes are based on the following integrated technologies:

Mode 1 – Low RF frequency of 0.8MHz (suitable for deep treatments).

Mode 2 – Medium RF frequency of 1.7 MHz (suitable for mid-layer

treatments). Mode 3 – High RF frequency of 2.45 MHz (suitable for

superficial treatments). Mode 4 – Combines all three frequencies for multi layer treatments.

To choose the desired treatment mode, press the mode number

Note

Changing modes will move the system to **Standby**. Press the **Ready / Standby** key to resume work.

N. **Energy** - To control the energy level, press the energy (+) or (-) keys. The smaller (+) / (-) keys will change energy by 2J. The large (+) / (-) keys will change energy by 10J.

O. **Minute counter** - Counts time usage for current treatment. Press the  key to zero the counter. When a pulse is not released 2 minutes or longer, counter will automatically pause and system will move to 'Standby' mode.

P. **Home** - Pressing the  key, operator opens the main menu.

6.12 V-IPL Treatment Screens

6.12.1 Filter Confirmation Screen

Upon choosing the V-IPL treatment on the main menu, the Filter Confirmation screen will open. It reflects the filter that has been connected and is currently in use while all other filters are grayed out.

Press the blinking **Confirm** key to enter the V-IPL setting screen.

Note

Each filter has a designated color both on the filter itself and on the top bar on the screen.

Filter Confirmation

> HR 570nm

> HR 630nm

> Acne 415nm

> VP 530nm

> SR 580nm

(Confirm >)

Figure 6.19: V-IPL Filter Confirmation Screen

6.12.2 V-IPL Setting Screen

Every filter has a unique setting screen according to clinical application and safety considerations.


Figure 6.20: V-IPL Setting

- A. **Skin Type** – Choose the skin type according to the Fitzpatrick
- Q. **Hair Color** – Choose the hair color according to the descriptions on screen.
- R. **Hair Thickness** – Choose the hair thickness.
- S. **Home** - Pressing the  key opens the main menu.
- T. **Change Filter** - Pressing the  Change Filter key opens the filter confirmation screen.
- U. **Setting** - The colored key with the arrow facing down shows the operator's current filter in use.
- V. **Treat** - Pressing the  key opens the treatment screen.
- W. **Totals Pulse Counter** - Counts the pulse usage for the total treatments on the handpiece.

Note

The operator must choose all parameters before advancing to the next screen. 'Treat' key will be activated only after choosing all parameters.

6.12.3.V-IPL Treatment Screen


Figure 6.21: V-IPL Treatment

A. **Pulse Structure** – When the operator presses the trigger and foot pedal once, there are 3 options for the pulse structure:

Single – One pulse

Multiple – 3 consecutive pulses

Rapid – One short pulse

Note

Single pulse is the default pulse structure.

- B. **Fluence** – Control the desired fluence using the (+) / (-) keys.
- C. **Duration** - Determines the length of the pulse during which energy is released (on time).
- D. **Interval** - Determines the delay between pulses (off time).

Warning

Changing pulse duration is directly related to the amount of energy released into the treatment area.

Note

For some pulse structures, **Duration** and **Interval** may be disabled.

- E. **Totals Pulse Counter** - Counts the pulse usage for the total treatments on the handpiece.
- F. **Pulse counter** - Counts emitted pulses for the current treatment. Press the  key to reset the counter to zero.
- G. **Cooling** - Indicates whether the cooling feature is enabled. The operator may switch the cooling feature on and off as desired by pressing the cooling key.

Note

- Cooling reduces pain and makes the treatment more comfortable to the patient. Therefore, it is recommended that you always use cooling.
- Cooling feature on is the default.
- Operator must always use ultrasound gel on the treatment area to prevent burns.

- H. **Home** - Pressing the  key opens the main menu.

- I. **Setting** - Pressing the  Settings key opens the setting menu.
- J. **Filter** - The colored key with arrow facing down shows the current filter in use.
- K. **Total duration bar** – indicates the total pulse duration (on and off time) only in multiple pulse.

L Standby / Ready -  indicates that the system is ready to be used and pulses can be emitted by pressing the trigger on the handpiece **and** stepping on the foot pedal. The system will move to  mode if:

- The operator changes:
 - Fluence
 - Pulse Structure
 - Duration
 - Interval
- The operator switches the cooling feature on or off
- A Pulse is not emitted for at least 3 minutes.

To resume work, press the **Standby** key.

To release a V-IPL pulse, press the trigger located on the handpiece **and** the foot pedal simultaneously.

6.13.V-Laser Nd:YAG Treatment Screens

6.13.1.Spot Confirmation Screen

Upon choosing the V-Laser Nd:YAG treatment on the main menu, the Spot Confirmation screen will open. It reflects the spot that has been connected and currently in use, while all other spots are grayed out.

Press the blinking **Confirm** key to enter the V-Laser Nd:YAG treatment screen.


Figure 6.22: V-Laser Nd:YAG Spot Confirmation Screen

6.14.V-Laser Nd:YAG Treatment Screen


Figure 6.23: V-Laser Nd:YAG Treatment Screen

Z. **Pulse Structure** – When the operator presses the trigger and foot pedal together, there are two options for the pulse structure:

Single – One pulse

Multiple – 3 consecutive pulses

Note

Single pulse is the default pulse structure.

- AA. **Duration** – Determines the length of the pulse during which energy is released (On time).
- BB. **Repetition** - Determines the number of pulses emitted per second.
- CC. **Fluence** – Controls the desired fluence using the (+) / (-) keys.
- DD. **Totals Pulse Counter** - Counts total pulse usage on the handpiece, per treatment.
- EE. **Pulse counter** – Counts the emitted pulses for the current treatment. Press the  key to reset the counter to zero.
- FF. **Cooling** - Indicates whether the cooling feature is enabled. The operator may switch the cooling feature on and off as desired, by pressing the cooling key.

Note

Cooling reduces pain and makes the treatment more comfortable for the patient. Therefore it is recommended that you always use cooling.

- GG. **Home** - Pressing the  key opens the main menu.
- HH. **Change Spot** - Pressing the  key opens the spot confirmation menu.
- II. **Treatment** - The colored key with the arrow facing down shows operator's current treatment in use.
- JJ. **Total duration bar** – Indicates the total pulse duration (on and off time) only while in multiple pulse.
- KK. **Standby / Ready** -  Indicates that the system is ready to work and pulses can be emitted by pressing the trigger on the handpiece **and** stepping on the foot pedal. System will move to  mode if:
- Operator changes:
 - Fluence
 - Pulse Structure
 - Duration
 - Interval
 - The operator switches the cooling feature to on or off
 - A pulse is not emitted for 3 minutes or longer.

To resume work, press the **Standby** key.

To release a Laser pulse, press the trigger located on the handpiece **and** the foot pedal simultaneously.

6.15. System Messages

The system continuously performs self-test procedures relating to the handpieces attached and the internal operation. Status or error messages will appear on the screen and display the relevant information when required. To continue, the operator must press the 'OK' key or wait for the progress bar to disappear.

The following screen categories may appear:

- **Consumable time:** The system is set to remind the operator of the handpiece running time. A reminder message will appear on the screen alerting the operator of the remaining time or usage on the system and will refer the user to the system information screen for more details.
- **Cooling:** The system will give an alert when the electrodes should be cooled down during the ST treatment. A progress bar will indicate when the cooling process is complete.
- **Contact Cooling:** During IPL and Nd:YAG treatments, the system will give an alert when the handpieces should be cooled down during the treatment. The message will disappear when the cooling process is completed.
- **Water system cooling:** During IPL and Nd:YAG treatments, the system will give an alert when the water system should be cooled down during the treatment. The message will disappear once the cooling process is completed.
- **V-FR tip warnings**
 - Treatment tips should only be used for a **single treatment** application on a single patient to avoid cross contamination.
 - System will alert the operator if an **empty tip** is inserted.
 - Every 20 pulses the operator should stop treatment and **clean** the V-FR tip of any debris that may affect the treatment. System will go into 'Standby' mode and the trigger/ foot pedal will be disabled until the "OK" key is pressed. To clean the tip, use a stiff-bristled brush dipped in alcohol. Make sure the tip is completely dry before resuming treatment.

Warning

- V-FR tip replacements should only be purchased through Viora's authorized local distributor.
-
- **Low impedance warning:** During FR and ST treatments, a warning will appear when the lower threshold of the impedance is detected. The operator should check the integrity of the tip and its electrodes and change tip if required.
 - **High impedance warning (bad coupling):** During the RF Treatments, the operator should make sure there is good contact between the skin and the electrodes.
 - **Water Filter expiry date:** When user checks the Water Filter status in Maintenance (Setup) screen, a warning will appear when the Water Filter expiry date has passed.

Figure 6.24: Warning Screen

Figure 6.25: Note Screen

7 Clinical Guidelines for RF Based Technology

7.1 Contraindications

The following conditions are contraindicated for treatment:

- **Pregnancy** and breast feeding.
- Treatment over tattoo or permanent makeup.
- Any **Permanent implant** in the designated treatment area (including but not limited to: silicon or injected chemical material, metal plates or other metal material).
- Any implanted **electronic device** anywhere in the body (such as cardiac pacemaker, defibrillator and cochlear implant).
- Severe concurrent conditions (such as **cardiac disorders**).
- **Blood coagulopathy** or excessive bleeding or bruising or use of blood thinning medications (anticoagulants), whether prescribed or over-the-counter.
- Active **skin disease** in the treatment area (such as psoriasis, sores, eczema and any type of rash).
- Tendency to **skin disorders** (such as keloids and impaired wound healing process) and extremely dry and fragile skin.
- **Skin cancer** (active or in the past), pre-malignant moles, malignancy (active or recent) or history of any kind of **cancer**.
- **Immunosuppressive diseases** (such as HIV Positive) or the use of immunosuppressive medications.
- **Endocrine disorder** (such as diabetes and thyroid disease).
- History of disease which may be stimulated by heat, such as recurrent Herpes in the treatment area.
- **Plastic surgery** such as face lifting or eyelid surgery within the past 12 months.
- Use of **Accutane** within the past 6 months.
- **Facial controlled abrasion** (such as dermabrasion, facial resurfacing, or medium / deep chemical peeling) within the past 3 months.
- **Botox injections** in the past 5-7 days
- **Chemical peel** or natural fillers (such as hyaluronic acid) in the past 2-4

weeks Additional Contraindications when treating with **V-FR** handpiece

- **Augmentation techniques** with injected bio-materials (such as injections of Botulinum toxin, collagen protein or fat) within the past 6 months.
- **Natural filler** (such as hyaluronic acid) in the treatment area within the past 6 months.
- **Extremely tanned skin** (including tanning beds or tanning creams) within the past 2 weeks.
- Any use of **NonSteroidal Anti-Inflammatory Drugs (NSAIDs)**, such as Nurofen, Advil, and Motrin) 7 days before and 7 days after every treatment.

Caution

Though the treatment is not advised for any of the above conditions, patients may be treated in some cases after consulting with their primary physician and providing a written consent for the treatment at their physician's discretion

7.2 Possible Side Effects

Improper use of the system may result in side effects. Although these effects are highly rare and expected to pass, any adverse reaction should be immediately reported to the patient's physician and to Viora's clinical team.

Appearance of side effects may occur during the treatment or shortly after. These side effects may include the following:

- Discomfort
- Appearance of any pigmentation (hyper-pigmentation or hypo-pigmentation)
- Changes in skin texture, such as burns, blisters and crust
- Excessive skin edema and /or erythema.
- Urticaria (hives)
- Scarring

The following adverse effects relate only to Fractional RF treatment:

- Post treatment mild or strong pain
- In-grown hair
- Skin infection

Caution

In the event side effects persist or are more severe than the above, the patient should consult a physician and notify Viora immediately.

Note

- Some dark-skinned patients (including Asian skin) may have a delayed response of one-to-two days after treatment and should be evaluated post-test procedure accordingly.
 - Proper use of the system by following the treatment protocol and operation instructions will prevent or minimize adverse effects.
-

7.3 Realistic Expectations

It is extremely important for the treatment's success and patient's satisfaction to set realistic expectations and advise the client of the following:

- Shortly after the treatment, erythema and edema are likely to appear in the treated zone. The edema is likely to clear up within a few hours while erythema may persist up to 12-24 hours following treatment.
- The degree of response to the treatment, and the number of treatment sessions required will vary among patients and will depend on the clinical and physiological condition at the start of the treatment regimen.
- The outcomes of the treatment will appear progressively, while temporary improvement in the skin will be observed immediate after each treatment; the maximum long-term result can be observed 2-4 months after the last treatment.
- Maintenance sessions may be required.

Additionally, the following should be explained to the patient prior to FR treatment:

- Light ablation may be detected with gradual improvement in skin texture/laxity.
- After the treatment and throughout the healing stage, small pin-point crusts may appear in the area of each ablated spot, at the contact area of the electrode pins with the skin. These crusts will undergo natural exfoliation.

7.4 Consult Meeting

The first visit of any patient should be dedicated to reviewing treatment possibilities, contraindications and performing a test procedure in order to establish that the patient is an appropriate candidate for the treatment, and obtain the optimal treatment parameters.

Upon the patient's initial visit, practitioners should follow the below steps during a consult meeting:

- Establish patient's needs and review symptoms.
- Set patient's expectations and desired end point for the treatment.
- Complete patient's medical history and review contraindications.
- Review with the patient the application method and operation, expected treatment results and possible side/adverse effects during or post treatment.
- Prior to commencing treatment sessions advise the patient to:
 - Avoid skin irritation or exposure to sun and using sunscreen during daylight hours.
 - Stop using any irritant topical agents for at least 2-3 days prior to treatment session.
- Sign the patient on an informed consent form including the above information.
- Establish a schedule of treatments.
- Take photos, measurements, etc.
- Perform a test procedure.

7.5 Treatment Follow-Up

Post-treatment follow-ups are very important to determine the most suitable regimen for each patient.

- The patient must be evaluated 2-3 days post treatment to make sure no side/adverse effects appear.
- Always inquire and update the patient history portfolio (including medical history and contraindications).
- Following completion of treatment course, patient should return after 4-6 weeks for evaluation.
- Prior to any maintenance treatment, perform a new test procedure in order to choose and verify the treatment parameters.

7.6 V-ST Clinical Guide

7.6.1 Indications of Use

The V-ST Handpiece is indicated for skin tightening on facial and body areas.

The V-ST handpiece utilizes proprietary CORE technology, bi-polar RF energy in three RF frequencies: 0.8, 1.7 and 2.45 MHz, as well as a combination of all three frequencies within a single pulse.

7.6.2 Test Procedure

The first clinical visit of any patient should be dedicated to reviewing treatment possibilities, contraindications and performing a test procedure in order to establish whether the patient is an appropriate candidate for the treatment, as well as testing for the optimal treatment parameters. In addition, performing a test procedure minimizes the possibilities for side effects and adverse response to the treatment. At the start of the treatment course and at every subsequent session, it is important to inquire and update the patient history portfolio, being mindful of new prescriptions, aesthetic or medical procedure, etc.

7.6.3 Pre-Treatment Preparation

The following steps should be taken before any treatment or test is performed.

- Ensure that the V-ST handpiece is properly inserted into the handpiece.
Make sure the electrodes are completely clean and dry. Left over residue on the handpiece may damage the system and cause superficial burns on the skin
- Use a mild cleansing solution (not alcohol based) to clean the treatment area.
- Make sure the treatment area is dry and free of any moisture, liquids or lint.
- If there is excessive hair on the treatment area, it is recommended to wax or shave it as it may interfere with the passage of the energy.

Note

Areas with metal dental implants (such as dental braces, caps and crowns) may be more sensitive to RF energy. In this case, the practitioner may isolate the area using gauze, dental rolls or a tongue depressor.

7.6.4 Treatment Guidelines

- Apply a layer of ultrasound gel (the layer should be half the height of the electrodes).
- Set the RF mode according to the treatment area and as indicated in the treatment protocol.
- Start with the minimal RF intensity indicated on the energy bar. Selecting the appropriate RF energy level is crucial for safe and effective operation.
- Place the handpiece on the treatment area in a perpendicular orientation and apply slight pressure, making sure both electrodes are in contact with the skin surface.
- Examine the skin's response to the chosen parameters before increasing the RF energy level.
- According to the practitioner's discretion and specific patient's reaction, the treatment parameters can be modified during subsequent sessions.

7.6.5 Treatment Schedule

Note that treatment sessions may be considerably longer and depend on the treatment area, physiological attributes and patient's response.

- Number of Treatments: 3-6
- Treatment Intervals: 2-4 weeks
- Maintenance: 1 treatment every 3-6 months, or as required.

7.6.6 Post-Treatment Care

- If patient notes excessive heat sensation or experiences strong erythema, the area may be cooled with the following; chilled Aloe Vera Gel, cold (not frozen) packs or cold compresses.
- In any case of rare blistering, refer the patient to a doctor.
- If adverse effects appear, promptly stop the treatment and apply cooling gel. Immediately refer to a doctor.
- Treatment should be resumed only after the patient has been thoroughly evaluated by a doctor and the side effects have subsided. Reduce the treatment parameters in subsequent treatments.
- Following the skin tightening treatment, patients should take special care of the treated area, making sure to avoid the use of hot water and discontinue use of abrasive or harsh products for two days.

7.7 V-FR Clinical Guide

7.7.1 Indications of use

The V-FR Handpiece is indicated for dermatological procedures requiring skin rejuvenation by means of pin-point skin ablation and coagulation.

The V-FR handpiece uses Viora's innovative SVC™ technology to give physicians greater control by delivering unprecedented penetration depth control, to determine the desired biological response. By delivering bi-polar RF energy in a pixilated fashion to the skin, isolated spots of skin are heated, leaving healthy surrounding tissue. This allows a faster healing process and extra cellular matrix remodeling mainly driven by fibroblast activity.

Warning

Dermatologists and aesthetic physicians may treat patients only after passing appropriate training and they must adhere to all local and state regulations regarding the use of this system.

7.7.2 Pre-Treatment Preparation

The following steps should be taken before any treatment or test is performed.

- Hair in the treatment area should be shaved. Hair may disrupt electrodes' full coupling with the skin and as result may cause burns.
- Use a mild cleansing solution to clean the treatment area. Make sure the treatment area is dry and free of any moisture, liquids or lint.
- Apply a topical anesthetic material (water-based only, containing at least 15% active ingredients) on the entire treatment area. Follow manufacturer guidelines when applying material.
- At the end of anesthetic activity (per manufacturer guidelines) wipe off the excess material from the skin with a tissue.
- Wash the skin with soap and lukewarm water and make sure the skin is completely clean of any residue.
- Once the skin is clean, sterilize the treatment area using an alcohol pad (at least 70% ethanol).
- Make sure the skin is completely dry before initiating the treatment.
- There is no need for the use of any medium solutions during the treatment (such as ultrasound gel or coupling gel).
- Before beginning the treatment make sure the handpiece is clean, dry and intact.
- Open the sealed pack of the treatment tip and check the integrity of the pins making sure they are all intact and that the tip's frame is undamaged.
- Connect the treatment tip to the handpiece and ensure the tip is secured in place.

Note

Areas with metal dental implants (such as dental braces, caps and crowns) may be more sensitive to RF energy. In this case, the practitioner may isolate the area using gauze, dental rolls or a tongue depressor.

7.7.3 Test Procedure

Before receiving a patient for a V30 treatment session, a test procedure must be performed during the initial consult meeting. The test procedure should be performed on a non-exposed site on or near the treatment area. If more than one area is to be treated, a test should be performed on all treatment sites.

The test procedure serves several purposes:

- Establishing the client is an appropriate candidate for the treatment:
 - For skin types I–III wait 1-2 days before assessing the skin response.
 - For skin types IV–VI wait 5-7 days before assessing the skin response.
- Providing the client with an opportunity to experience the FR treatment.
- Finding the correct parameters according to the skin's reaction and the patient's tolerance to the test procedure
- Minimizing possibilities for adverse effects.

Warning

- Patients with skin types V–VI may be more sensitive to treatment. Lower energy levels should be applied.
-

The test procedure should be conducted as follows:

- For facial treatments, perform the test on concealed areas such as the back of the neck. For body treatments perform the test on the target area.
- Choose the relevant treatment program (Shallow, Medium or Deep).
- Set the energy to 0.25 J.
- Apply one pulse. If no heat sensation was reported by the patient, increase the energy level by 0.25 J and apply an additional pulse in an adjacent spot.
- Continue in the same manner until the patient reports a noticeable heat sensation. Note the energy level you have reached (X)
- Release additional 2-3 pulses at the same energy (x J)
- Increase the energy level by 0.25J (X+0.25J) and release 2-3 pulses
- Increase the energy level again by 0.25J (X+0.50J) and release additional 2-3 pulses
- Never release pulses in the same spot
- Wait 10 minutes to allow a pattern to appear. The desired end point is local erythema and edema.
- If the skin does not show signs of erythema and edema after 10 minutes, increase the RF energy by an additional 0.25J (X+0.75J) and release an additional pulse adjacent to the location of the previous pulse.
- Wait an additional 10 minutes to allow the end point to appear.
- Treatment area must be examined 1-2 days (Fitzpatrick I–III) and 5-7 days (Fitzpatrick IV–VI) after test procedure. A grid-like pattern matching the electrode pins should be visible on the treated area.

Note

- The test procedure is to be executed and evaluated prior to beginning the complete treatment schedule.
- a) For the complete treatment protocol please refer to the treatment protocol attached to this user manual.

At the start of the treatment course, and at every subsequent session, it is important to inquire and update the patient history portfolio, being mindful of new prescriptions, aesthetic or medical procedures, etc.

7.7.4. Treatment Guidelines

- After completing the pre-treatment guidelines, set the treatment parameters according to the test results.
- Place the tip in full contact with the skin at a 90 degree angle to the surface of the target region.
- Apply one pulse by pressing down the trigger on the handpiece or the foot pedal.

Note

Pressing the trigger/ foot pedal continuously will change the system into 'Standby' mode. Only one pulse at a time is possible when treating with the V-FR hand piece.

- The LED located at the back of the handpiece will light up as an indication of the released pulse. The light has no bearing on the actual treatment and is only a visual indication that a pulse has been released.
- On Standby mode, the LED light will be OFF
- On Ready mode, the LED light will blink.
- Once the pulse is completed a beeping sound is heard to indicate the pulse is completed. The handpiece can then be safely moved to an adjacent spot in the treatment area.
- For every pulse released, the number of pulses indicated on the display screen will increase by one.
- Do not leave untreated spaces between the pulses and do not overlap between the pulses.
- Pass over the treatment area once.
- In case the system recognizes insufficient coupling, the RF energy will be disabled and the system will not emit a pulse. This will be accompanied by a message on screen and an audio indication. The system will not count this as a pulse.
- Throughout the treatment, make sure that the electrode pins are:
 - Clean and have no cellular remains (appearing as black material).

- Make sure that the tips are cleaned every 20 pulses while the system is in *Standby* mode. A note will appear on the screen every 20 pulses.
- Undamaged – ensure that no pins are missing or damaged.
- Examine the treated area for the desired end point; erythema and edema on the treatment area.
- Edema will normally appear within several minutes after the treatment.
- At the end of the treatment, examine the area to make sure there are no untreated patches of skin. If such patches exist follow the below instructions:
 - Set the system to *Standby* Mode.
 - Clean the untreated patches and the tip with a dry lint-free gauze pad.
 - Pass over the untreated patches and avoid treated areas.
- At the end of every treatment session remove the treatment tip and clean the handpiece according to the cleaning guidelines. Insert the V-FR tip cover when the system is not in use in order to protect the handpiece.

Warning

- **Do not** treat over the eyelids or the lips.
- **Do not** treat different patients with the same tip, in order to prevent cross- contamination.
- **Do not** treat with an expired tip.
- **Do not** use the V-FR Tip if the sterile package has been damaged
- Note that in darker skin types erythema is not easily noticeable; therefore, the end point of edema will be more prominent in these patients.
- Before cleaning the tip make sure to set the system to *Standby* mode

7.7.5. Treatment Schedule

Note that treatment sessions may be considerably longer and depend on the treatment area, physiological attributes and patient's response.

- Number of Treatments: 3-4
- Treatment Intervals: 4-6 weeks
- Maintenance: 1 treatment every 6-12 months, or as required.

7.7.6. Post-Treatment Care

- Emollient cream should be applied to the treatment area.
- If a patient notes excessive heat sensation or strong erythema or any discomfort, the area may be cooled with the following; chilled Aloe vera gel, cold (not frozen) packs or cold compresses.
- If blistering or ulcerated skin occurs (rarely), mild hydrocortisone ointment may be applied or a prescribed antibiotic ointment.

- Usually 1-3 days after the treatment small scabs (less than 1 mm in diameter) will appear and will possibly remain for several days. Do not scratch or touch the scabs and allow them to shed off naturally from the skin.
- Advise the patient of the following precautions during the first 2 days after FR treatment:
 - Avoid using hot water massages.
 - Avoid the use of abrasive or harsh products.
 - Avoid any mechanical or thermal damage to the treated area.
 - Avoid any type of skin tanning (by sun exposure, tanning beds or artificial sunless tanning lotions) during the entire course of treatment.
 - Keep the skin clean in order to avoid contamination or infection.
 - Use a high-factor sunscreen (at least 30 SPF) to protect the treated area for at least 1 month after every treatment.
- Moisturizing cream may be applied several hours after treatment (~½ a day).
- Make-up may be applied 24 hours after treatment.
- Avoid using any kind of soap on the treatment area up to 24 hours after the treatment.

7.8. FC and BC Clinical Guide

7.8.1 Indications of use

The V-Form Handpiece, with the BC and FC applicators, is indicated for temporarily reduction of cellulite and Body contouring via temporary circumferences reduction.

7.8.2. Pre-Treatment Preparation

The following steps should be taken before any treatment or test is performed.

- Use a mild cleansing solution (not alcohol based solution) to clean the treatment area.
- Make sure the treatment area is dry and free of any moisture, liquids or lint.
- Apply a very thin layer (<1 mm) of glycerin oil onto the treatment area to assure good coupling and conductivity.
- Place the applicator perpendicular to the treatment area.

7.8.3 Test Procedure

The first visit of any patient to a salon, spa or clinic should be dedicated to reviewing treatment possibilities, contraindications and performing a test spot in order to establish that the patient is an appropriate candidate for the treatment and test for the optimal treatment parameters. In addition, performing a test spot procedure minimizes the possibilities for side effects and adverse response to the treatment.

At the start of the treatment course and at every subsequent session, it is important to inquire and update the patient history portfolio, being mindful of new prescriptions, aesthetic or medical procedure, etc. Selecting the correct treatment parameters is essential for safe, accurate and efficient treatment.

The test procedure is divided into 2 steps;

A. Step 1 – Finding the patient comfort level

1. Insert the BC/FC applicator (depending on the treatment area).
2. Set the treatment parameters as follows:
 - RF Mode - **MODE IV**
 - RF Intensity – **Level 2**
 - Vacuum Intensity – **Level 1**
3. Identify a vector zone of approximately 4 cm x 24 cm in size (enough for four foot prints of the applicator), in the less visible site of the desired treatment area.
4. Divide the area into 4 squares as below:

I	II	III	IV
Vacuum 1	Vacuum 2	Vacuum 3	Vacuum 4

5. Perform 3 stacked pulses on the first square in the test area.
6. If no adverse effects appear and patient comfort allows, increase the vacuum intensity by 1 level and repeat the same procedure in the next square.
7. Continue in the same manner in square 3 (vacuum level 3) and square 4 (vacuum level 4).
8. Repeat the procedure 5 times passing through all squares and adjusting vacuum levels accordingly.
9. Continuously monitor the skin temperature (using the IR thermometer screen on the handpiece) not to exceed 42°C.
10. Ask the patient to assess the discomfort during the test procedure using the scale below and record it on the patient sheet:

Discomfort Scale:

0 – No Discomfort 1 – Mild 2 – Moderate 3 – Significant 4 - Unbearable

11. The highest vacuum level applied that generated moderate discomfort will be the selected vacuum level for the next phase.


Figure 7.1: IR Thermometer screen

B. Step 2 – Validating the parameters

1. Select another zone in the same treatment area (wide enough for 2 lines of passes) and apply 2-3 pulses in each spot until the skin temperature reaches an endpoint of 39-42⁰C.
2. Once endpoint is reached, continue in the same manner and maintain the temperature for 7 minutes performing a 'mini' treatment.

Note

In the parameter validation step, it is recommended to work on more sensitive areas within the treatment area; if such exist (this includes areas such as the flanks in the abdomen area). Ask for patient feedback during treatment.

3. Treatment area should be assessed 48-hours post test procedure for any local side effects.
4. If the patient developed side effects, the chosen vacuum level should be 1 below the initial patient's indication.
5. Make a record of the patient's subjective response as well as the skin's response on the patient's treatment sheet.

Note

- The operator may increase the parameters in subsequent treatments according to patient's comfort level and skin's response.
 - Make sure to follow clean up procedures for consumables and systems' at the end of the treatment session.
-

7.8.4 Treatment Guidelines

1. Make sure there is **full contact** between the electrodes and the skin without putting unnecessary pressure on the applicator.
 2. Pulses are released when the trigger on the applicator is pressed down.
 3. If the trigger is held down, the system will continue to release consecutive pulses until the press button is released.
 4. Once the pulse and vacuum stop slide the applicator to the next spot on the treatment area.
 5. Overlap no more than 10% of the previous pulse.
-

Note

Sliding the applicator before the vacuum pulse ends may cause bruising.

6. Clean the applicator and tip thoroughly at the end of every treatment session using an alcohol based solvent.
7. Let the handpiece and applicators dry completely before using them again.

7.8.5 Treatment Schedule

Note that treatment sessions may be considerably longer and depend on the treatment area, physiological attributes and patient's response.

- Number of Treatments: 6-8
- Treatment Intervals: 1 week
- Maintenance: 1 treatment every 3-4 months, or as required.

7.8.6 Post-Treatment Care

If patient notes excessive heat sensation or strong erythema, the area may be cooled with the following; chilled Aloe Vera gel, cold (not frozen) packs or cold compresses.

- If rare blistering occurs, mild hydrocortisone ointment may be applied or a prescribed antibiotic ointment.
- If adverse effects appear, promptly stop the treatment and apply cooling gel.
- Treatment should be resumed only after the patient has been thoroughly evaluated and the side effects have subsided. Reduce the treatment parameters in subsequent treatments.

8 Clinical Guidelines for Light Based Technology

8.1 Contraindications

The following conditions present exclusion criteria and should be strictly observed, evaluated and confirmed at each visit before the start of the treatment course and prior to every treatment session:

- Prolonged exposure to sun or artificial tanning during the past 4 weeks prior to treatment.
- Use of photosensitive medication or herbs (such as St. John's Wort) within 3 month prior to treatment.
- Recent or history of cancer, particularly skin cancer.
- Pregnancy, expectation of pregnancy, postpartum or nursing.
- Use of oral isotretinoin (Accutane) within the last 6 months.
- History or current use of Isotretinoin
- Use of Tretinoin (Retin A) within the last 2 weeks.
- Epilepsy.
- Fragile and dry skin.
- Hormonal disorders (such as diabetes), unless under control.
- Use of anticoagulants.
- Undiagnosed lesions.
- Diseases which may be stimulated by light.

- History of Keloid scarring.
- Diabetes.
- Koebnerization of skin, such as Vitiligo (Hypopigmentation) and psoriasis.
- Treatment area exclusion criteria:
 - Inflammation at the treatment area
 - Tattoos
 - Moles, nevi
 - Permanent make up
 - Silicone implants in the treatment area.
 - Any inflammatory skin condition, active infection or herpes simplex in the treatment area.
- Effective long term hair removal with light based devices requires the presence of a hair follicle for photothermolysis to commence. Waxing, plucking, electrolysis, or any other method of long term hair removal that affects the follicle should be avoided prior to and in between photoepilation treatments.
- Patients should be warned of potential harm that may arise when treating postinflammatory hyperpigmentation (PIH).

Note

- Contraindications should be thoroughly evaluated and confirmed at each patient visit!
 - Patients with Hirsutism may have an underlying endocrine disorder. Causes of hirsutism should be sought and treated by a physician; in this case the patient should supply written approval from their physician.
-

8.2 Possible Side Effects

Improper use of the system may result in side effects. Although these effects are highly rare and expected to eventually heal, any adverse reaction should be immediately reported to the patient's physician and to Viora's clinical team.

Appearance of side effects may occur during the treatment or shortly after. These side effects may include the following:

- **Pain**
Pain sensation is individual and may vary between patients. When a light pulse is triggered, the absorption of light and its transformation to heat may cause some degree of discomfort to the patient.
- **Blistering**
Inappropriate use of energy level may result in the formation of a crust or blister. In such cases, refer the patient to a physician.
- **Pigmentation Changes**
There may be a change of pigmentation in the treated area. Most cases of hypopigmentation or hyperpigmentation occur in people with darker skin types, or when the treated area has been exposed to sunlight following phototherapy treatment. In some patients, hyperpigmentation occurs despite protection from the sun. This usually fades within three to six months, but the change of pigment may be permanent in rare cases. In any case, it is necessary to refer the patient to a physician.
- **Scarring**
There is a small risk of scarring such as enlarged hypertrophic scars or in very rare cases, abnormal, large, raised keloid scars. To reduce the probability of scarring, refer the patient to a doctor.
- **Swelling**
Immediately after treatment, skin may temporarily swell and may be noticed at the nose and cheeks. Swelling usually subsides within three to seven days. In other cases of rare swelling, refer the patient to a doctor.
- **Paradoxical unwanted hair growth**

Warning

In any case of rare side effects refer the patient to a physician immediately.

8.3 Realistic Expectations

Setting realistic expectations is highly important to the treatment success and patient satisfaction and these should be set with the patient before the start of the treatment course.

The practitioner should thoroughly review treatment procedure, protocol, possible side effects, and the treatment regimen, number of sessions, intervals and post-treatment care.

Note

It is recommended to take pictures at the beginning of the treatment course and to log each treatment and note progress in a treatment log.

In addition, the Patients should be informed of the following:

- There may be some discomfort/pain associated with the treatments.
- Transient erythema/edema may appear immediately following treatment.
- There is a small risk of adverse reaction such as changes in the texture and pigmentation of the skin, which are usually transient.

Treatment results vary from person to person and there is a small probability that improvement will not be seen.

It is highly recommended that the patient sign an informed consent form after reviewing all relevant information. Since this form has juridical value and might be different from country to country, it is recommended to consult with a legal advisor. In general, all consent forms should include the statements below:

1. The patient received an explanation about the treatment, results and potential side effects.
2. The patient agrees to perform the treatments under these circumstances.

8.4 Consult Meeting

The first visit of any patient should be dedicated to reviewing treatment possibilities; contraindications and performing a test procedure in order to establish that the patient is an appropriate candidate for the treatment, and determining the optimal treatment parameters.

Upon the patient's initial visit, practitioners should follow the below steps during a consult meeting:

- Establish patient's needs and review symptoms.
- Set patient's expectations and desired end point for the treatment.
- Complete patient's medical history and review contraindications.
- Review with the patient the application method and operation, expected treatment results and possible side/adverse effects during or post treatment.
- Prior to commencing treatment sessions advise the patient to:
 - Avoid skin irritation or exposure to sun (4 weeks prior to treatment) and using sunscreen during daylight hours.

- Stop using any irritant topical agents for at least 2-3 days before treatment session.
- Sign the patient on an informed consent form including the information above.
- Establish a schedule of treatments.
- Take photos, measurements, etc.
- Perform a test procedure.

8.5 Test Procedure

Before receiving a patient for a treatment session, a test procedure must be performed during the initial consult meeting. The test procedure should be performed at a non-exposed site of the treatment area. If more than one area is to be treated, a test should be performed on each treatment site.

The test procedure serves several purposes:

- Establishing the client is an appropriate candidate for the treatment:
 - For skin types I–III wait 1-2 days before assessing the skin response.
 - For skin types IV–VI wait 5-7 days before assessing the skin response.
- Providing the client with an opportunity to experience the treatment.
- Finding the correct parameters according to patient's tolerance to the test procedure
- Minimizing possibilities for adverse effects.
- Follow the test procedure protocol provided with the system.

8.6 Treatment Follow-Up

Post-treatment follow-ups are very important to determine the most suitable regimen for each patient.

- The patient must be evaluated 2-3 days post treatment to make sure no side/adverse effects appear.
- Always inquire and update the patient history portfolio (including medical history and contraindications).
- Following completion of treatment course, patient should return after 4-6 weeks for evaluation.
- Prior to any maintenance treatment, perform a new test procedure in order to choose and verify the treatment parameters
- Patients should observe the following instructions and preferably be given an information sheet that details all post-care instructions:
 - Exposure to sunlight – Advise patients to use sunscreen with a **high SPF** and protect the treated area by **avoiding direct exposure** to sunlight and the use of tanning beds.
 - Maintain **body temperature** – Refrain from hot baths, sauna, steam etc., as well as excessive exercise for at least 48 hours.
 - Avoid picking or scratching treated area.
 - Avoid **exfoliation** and peels for 7 days.

- Refrain from wearing tight clothing.
- Home care depilation methods should be suspended.

8.7 V-IPL Clinical Guide

8.7.1 Indications of use

The V-IPL handpiece utilizes Intense Pulse Light technology which is based on the principle of Selective Photothermolysis. According to this principle, pulsed light parameters (wavelength, pulse duration, energy fluence and pulse configuration) are chosen to selectively target light sensitive molecules while minimizing damage to surrounding tissues. The V-IPL handpiece has five interchangeable filters and two different spot sizes, ensuring complete treatment coverage for all skin types and conditions. Applications for V-IPL indicated for::

- Moderate inflammatory acne (acne vulgaris)
- Benign pigmented epidermal and cutaneous lesions including warts, scars, striae, dyschromia, hyperpigmentation, melasma, ephelides (freckles), lentigines, nevi, melasma, and café-au-lait macules
- Benign cutaneous vascular lesions including port wine stains, hemangiomas, facial, truncal and leg telangiectasias, rosacea, erythema of rosacea, angiomas and spider angiomas, poikiloderma of Civatte, leg veins and venous malformations
- Removal of unwanted hair to effect stable long-term or permanent hair reduction

8.7.2 Pre-Treatment Preparation

- Before starting treatment make sure the handpiece, light guide surface and filter are completely clean and dry. Left over residue on the light guide may damage the system and cause superficial burns on the skin.
- Use a mild cleansing solution (not alcohol based) to clean the treatment area.
- Cover moles and beauty spots with a non-flammable piece of white pad.
- Apply a layer of ultrasound gel (~2 mm) to the treatment area (use clear and chilled ultrasound gel).
- Prior to turning the system ON ensure the emergency stop button is not pressed.
- Make sure operator and patient wear safety eyewear (supplied).
- Verify on the display screen that the actual filter and spot size are connected as required.
- Place the handpiece in full contact with the skin and press the button located on the handpiece.

Warning

- Cleaning the handpiece should always be done while the system is turned OFF by pressing the On/Off switch located on the rear panel of the system.

-
- Always apply a layer of ultrasound gel to the treatment area to prevent skin burns.
-

8.7.3 Treatment Parameters

8.7.3.1 Step-by-Step Operation

- Insert the required filter for the specific treatment into the designated slot in the V-IPL handpiece
- Connect the V-IPL handpiece to the system.
- Select the treatment from the main menu.
- Confirm the selected filter as it appears in the “Filter Confirmation” screen
- Choose the parameters per patient’s characteristics. The selected parameters will be highlighted.

Note

The operator must choose all parameters before proceeding to the next screen. ‘Treat’ key will be activated only after choosing all parameters.

- Press “Treat”.
- Review the parameters in the “Treatment Screen” and change if required.
You may change the pulse structure to ‘Rapid’ or ‘Multiple’ as necessary, per type of treatment.
- Change to ‘Ready’ mode and start the treatment.
- The LED located at the back of the handpiece, as well as on the system, will light up as an indication of the released pulse.
- In Standby mode, the LED light will be OFF.
- Upon entrance to the Ready Mode, a beeping sound is heard to indicate that the system is ready to emit a pulse.
- In Ready mode, the LED light will blink.
- Once a pulse has been emitted a beeping sound is heard to indicate the pulse is complete. The handpiece can then be safely moved to an adjacent spot on the treatment area.
- For every pulse released, the number of pulses indicated on the display screen will increase by one.
- When the treatment is completed, place the V-IPL handpiece in its cradle.

8.7.3.2 Example 1 – Hair Removal – Single Pulse

- Insert the 570nm filter into the V-IPL handpiece.
- Connect the V-IPL handpiece to the system.

- Confirm the selected filter as it appears in the “Filter Confirmation” screen.
- Select the patient’s characteristics. The Fitzpatrick scale is described in Table 1.


Figure 8.1: Example 1 – Setting screen

Reaction to Sun Exposure	Typical Features	Skin Type
Always burns, rarely tans	Pale skin, light eyes, light hair	Type I
Normally burns, sometimes tans	Fair skin, light eyes	Type II
Sometimes burns, usually tans	Darker white skin	Type III
Burns minimally, tans easily	Light brown skin	Type IV
Rarely burns, easily tans darkly	Brown skin	Type V
Never burns, always tans darkly	Dark brown or black skin	Type VI

Table 8.1: Fitzpatrick scale

Note

- Patients with freckled skin are classified 1 level higher than their base skin type on the Fitzpatrick scale.

- Proper maintenance of the handpiece and system will guarantee a long lifetime for the system and safe treatments for patients.

- In the treatment screen, review the duration and fluence. Change parameters as necessary.


Figure 8.2: Example 1 – Treatment screen

- To choose either **multiple** or **rapid** pulse structure, press the relevant selection on the left side, then review and edit the parameters, if necessary.
 - Adjust the duration by using the (+) / (-) buttons. Adjust the fluence by using the (+) / (-) buttons
- Change to 'Ready' mode and start the treatment.
- When the treatment is completed, place the V-IPL handpiece in its cradle.

8.7.4 Treatment Schedule

Note that treatment sessions may be considerably longer and depend on the treatment area, physiological attributes and patient's response.

Long-Term Hair removal

- Number of Treatments: 6-10
- Treatment Intervals: 4-10 weeks (as per treatment area)
- Maintenance: 1 treatment every 6 months, or as required

Skin Rejuvenation

- Number of Treatments: 4-6
- Treatment Intervals: 2-4 weeks
- Maintenance: 1 treatment every 3-6 months, or as required

Pigmented Lesions

- Number of Treatments: 2-6
- Treatment Intervals: 3-6 weeks
- Maintenance: As required

Vascular Lesions

- Number of Treatments: 4-8
- Treatment Intervals: 3-6 weeks
- Maintenance: As required

Acne clearance

- Number of Treatments: 8-10
- Treatment Intervals: 3-7 days
- Maintenance: As required

8.7.5 Post-Treatment Care

Should the patient experience any heat or burning sensation it is recommended to apply ice packs or other cooling means to the treatment area immediately following the treatment. This will also assist with sensitive skin that tends to swell. Chemical cold packs are not recommended if their temperature is below 4°C. In case of blisters or lacerated skin, refer the patient to a doctor.

Patients should observe the following instructions and ideally be given a document to take home, covering all post care instructions:

- Hydration - Patients should keep their body hydrated and drink plenty of water.
- Exposure to sunlight - patients are strongly advised to wear high SPF sunscreen (at least SPF 30) and protect the treated area by avoiding direct exposure to sunlight and avoid the use of tanning beds for the entire treatment course. The use of tanning beds or direct exposure to sunlight after treatment sessions may enhance melanin regeneration, which may result in hyperpigmentation.

- Maintain body temperature – Refrain from raising body temperatures by hot baths, sauna, steam etc., as well as avoid excessive exercise (as this heats up the body) for at least 48 hours to avoid adverse reactions.
- Treatment area – avoid the following:
 - Picking or scratching the treated area.
 - Exfoliating or peeling for 7 days.
 - Wearing tight clothing.
 - Swimming – avoid swimming in highly chlorinated water.
 - Hygiene – treatment area should remain clean and dry.
- Home care depilation methods (such as waxing, plucking, threading or creams) should be suspended for the entire duration of the photo-epilation treatment schedule as these interrupt with the hair's natural growth cycle and may affect treatment results.
- Makeup – makeup can be applied immediately after treatment. However, patients are advised to stop makeup application if the treatment site scales or exfoliates.

8.8 V-Laser Nd:YAG Clinical Guide

8.8.1 Indications of use

The V-Laser Nd:YAG handpiece emits long pulse light at 1064 nm and is indicated for the coagulation and hemostasis of soft tissue. The laser technology is based on the principle of Selective Photothermolysis. According to this principle, pulsed light parameters (wavelength, pulse duration, energy fluence and pulse configuration) are chosen to selectively target light sensitive molecules while minimizing damage to surrounding tissues. The V-Laser Nd:YAG includes four spot sizes for full customization, according to customers' needs.

The V-Laser Nd:YAG is indicated for:

- Removal of unwanted hair, for stable long-term, or permanent, hair reduction through selective targeting of melanin in hair follicles
- Removal or lightening of unwanted hair (with and without adjuvant preparation)
- Treatment of pseudofolliculitis barbae (PFB)
- Benign vascular lesions such as, but not limited to treatment of: port wine stains, hemangiomas. Warts, superficial and deep telangiectasias (venulectasias), reticular veins (0.1-4.0 mm dia.) of the leg, rosacea, venus lake, leg veins, spider veins, poikiloderma of civatte and angiomas
- Benign cutaneous lesions, such as, but not limited to: warts, scars, striae.
- Benign pigmented lesions such as, but not limited to: lentigos (age spots), solar lentigos (sun spots), café-au-lait macules, seborrheic, nevi and nevus of Ota, chloasma, verrucae, skin tags
- Pigmented lesions to reduce lesion size, for patients with lesions that would potentially benefit from aggressive treatment, and for patients with lesions that have not responded to other laser treatments

- The non-ablative treatment of facial wrinkles, such as, but not limited to: periocular wrinkles and perioral wrinkles
- Laser skin resurfacing procedures for the treatment of: acne scars and wrinkles
- Reduction of red pigmentation in hypertrophic and keloid scars where vascularity is an integral part of the scar
- Temporary increase of clear nail in patients with onychomycosis (e.g., dermatophytes *Trichophyton rubrum* and *T. mentagrophytes*, and/or yeasts *Candida albicans*, etc.).

8.8.2 Pre-Treatment Preparation

- Before starting treatment make sure the handpiece, filter and lightguide surface are completely clean and dry. Left over residue on the light guide and the filter may damage the system and cause superficial burns on the skin.
- Use a mild cleansing solution (not alcohol based) to clean the treatment area.
- Apply a layer of ultrasound gel (~2 mm) to the treatment area (use clear and chilled ultrasound gel).
- For Skin Rejuvenation application, cover moles and beauty spots with a non-flammable piece of white pad.
- Prior to turning the system ON ensure the emergency stop button is not pressed.
- Make sure operator and patient wear safety eyewear (supplied).
- Verify that on the display screen is that the actual filter and spot size are connected as required.
- Place the handpiece in full contact with the skin and press the button located on the handpiece.

8.8.3 Treatment Parameters

8.8.3.2 Step-by-Step Operation

- Connect the required spot to the Nd:YAG handpiece
- Connect the handpiece to the system
- Confirm the selected spot size appears in the "Spot Confirmation" screen.
- Choose the treatment screen, review the duration and fluence. Change parameters as necessary.


Figure 8.3: Treatment screen

- On the treatment screen, review the duration and fluence. Change parameters as necessary.
- To choose **multiple** pulse structure, press the relevant selection on the left side, then review and edit the parameter if necessary.
- Adjust the duration by using the (+) / (-) buttons. Adjust the fluence by using the (+) / (-) buttons. Change to 'Ready' mode and start the treatment,
- Change to 'Ready' mode and start the treatment.
- The LED located on the system, will light up as an indication of the released pulse.
- In Standby mode, the LED light will be OFF.
- In Ready mode, the LED light will blink.
- Upon entrance to the Ready Mode, a beeping sound is heard to indicate that the system is ready to emit a pulse.
- Once a pulse has been emitted a beeping sound is heard to indicate the pulse is complete. The handpiece can then be safely moved to an adjacent spot on the treatment area.
- For every pulse released, the number of pulses indicated on the display screen will increase by one.
- When the treatment is completed, place the handpiece in its cradle.

8.8.4 Treatment Schedule

Note that treatment sessions may be considerably longer and depend on the treatment area, physiological attributes and patient's response.

Long Term Hair removal

- Number of Treatments: 8-12
- Treatment Intervals: 4-10 weeks (depending on treatment area)
- Maintenance: 1 treatment every 6 months, or as required

Skin Rejuvenation

- Number of Treatments: 4-6
- Treatment Intervals: 2-4 weeks

- Maintenance: 1 treatment every 3-6 months, or as required

Vascular Lesions

- Number of Treatments: 4-8
- Treatment Intervals: 6-12 weeks
- Maintenance: As required

Onychomycosis

- Number of Treatments: 4-6
- Treatment Intervals: 7-10 days
- Maintenance: As required

8.8.5 Post-Treatment Care

Should the patient experience any heat or burning sensation it is recommended to apply ice packs or other cooling means to the treatment area immediately following the treatment. This will also assist with sensitive skin that tends to swell. Chemical cold packs are not recommended if their temperature is below 4°C. In case of blisters or lacerated skin, refer the patient to a doctor.

Patients should observe the following instructions and ideally be given a document to take home, covering all post care instructions:

- Hydration - Patients should keep their body hydrated and drink plenty of water.
- Exposure to sunlight - patients are strongly advised to use sunscreen and protect the treated area by avoiding direct exposure to sunlight and the use of tanning beds for the entire treatment course. The use of tanning beds or direct exposure to sunlight after treatment sessions may enhance melanin regeneration, which may result in hyper-pigmentation.
- Maintain body temperature – Refrain from raising body temperature by hot baths, sauna, steam etc., as well as avoid excessive exercise (as this heats up the body) for at least 48 hours to avoid adverse reactions.
- Treatment area – avoid the following:
 - Picking or scratching the treated area.
 - Exfoliating or peeling for 7 days.
 - Wearing tight clothing.
 - Swimming – avoid swimming in highly chlorinated water.
 - Hygiene – treatment area should remain clean and dry.
- Home care depilation methods (such as waxing, plucking, threading or creams) should be suspended for the entire duration of the photo-epilation treatment schedule as these interrupt with the hair's natural growth cycle and may affect treatment results.

Makeup – makeup can be applied immediately after treatment. However, patients are advised to stop makeup application if the treatment site scales or exfoliates

9 Data for the technical maintenance of medical device.

The V30 system requires standard maintenance and cleaning procedures to ensure optimal system performance. However, note that the system should be technically serviced only by authorized service technicians.

This section details the general routine maintenance procedures that may be performed by trained and approved operators.

Warning

DISCONNECT THE POWER SUPPLY CABLE BEFORE PERFORMING ANY MAINTENANCE OR SERVICE PROCEDURES TO AVOID ELECTRIC SHOCK!

- For continued protection against fire, replace the fuse only with one of the same type and rating
-

9.1 Cleaning the System

- It is recommended to clean the system at least once a week.
- Turn the system OFF.
- Wipe all external surfaces with a soft, damp, non-abrasive cloth.
- To clean the LCD screen, use a suitable solution for computer screens.

9.2 Cleaning the V-ST, V-Form Electrodes

V-ST, V-Form electrodes must be thoroughly cleaned after every treatment.

- Turn the system OFF.
- Wipe the V-ST electrodes with a dry cloth to remove excess ultrasound gel.
- Wipe the BC/FC applicators with either wet wipes (water or alcohol based) to remove residual glycerin oil.
- Using a cotton pad and an Isopropyl-based alcohol solvent (70%), carefully wipe the electrodes, making sure all corners are thoroughly cleaned.

Note

Gel accumulation on V-ST electrodes may cause a spark and potentially a superficial burn in the next treatment performed.

9.3 Cleaning the V-IPL Handpiece

9.3.1 Cleaning the Filters

Warning

- Filter should be constantly cleaned before and after every treatment.
- Damaged or broken filters should be replaced.
- System should be turned off before cleaning the filters

To clean the filter:

- Turn the system OFF.
- Remove the filter by pulling it from the top of the handpiece.
- Using the supplied cleaning cloth, carefully clean the surface of the filter (both sides).
- Using the supplied cleaning cloth and an Isopropyl-based alcohol solvent (70%), carefully wipe the surface of the filter (both sides).
- Carefully wipe the surface of the filter (both sides) until completely dried.
- Safely store the filter in the supplied padded case.

9.3.2 Cleaning the Spots

Spots should be cleaned of dirt and gel residue after each treatment.

Remaining residues on the light guide may be heated and harm the patient during treatment.

Warning

- Never wash the 6.4cm² spot in water to avoid water penetrating the handpiece
- System should be turned off before cleaning the 6.4cm² spot

To clean the 6.4cm² spot:

- Turn the system OFF.
- Make sure the 2.4cm² spot is disconnected.
- Using a soft tissue, wipe the surface of the spot to remove excess ultrasound gel, if used.
- Using the supplied cleaning cloth and an Isopropyl-based alcohol solvent (70%), carefully wipe the surface of the spot.
- Wipe the surface of the spot until a completely

dried. To clean the 2.4cm² spot:

- Disconnect the 2.4cm² spot from the 6.4cm² spot.
- Using a soft tissue, wipe the surface of the spot (both sides) to remove excess ultrasound gel.
- Using the supplied cleaning cloth and an Isopropyl-based alcohol solvent (70%), carefully wipe the surface of the spot (both sides).
- Wipe the surface of the spot (both sides) until completely dried.
- Safely store the spot in the padded case supplied.

9.4 Cleaning the V-Laser Nd:YAG Handpiece

9.4.1 Cleaning the Spots

Spots should be cleaned of dirt and gel after each treatment.

Remaining residue on the spot may be heated and harm the patient during treatment.

Warning

- Spot should not be cleaned right after treatment.
 - Operator must wait for the surface to cool down.
 - System should be turned off before cleaning the spots and the protective window
-

To clean the spots:

- Turn the system OFF.
- Remove the spot from the slot.
- Using a cotton pad and an Isopropyl-based alcohol solvent (70%), carefully wipe the surface of the spot (both sides).
- Using a cotton pad wipe the surface of the spot (both sides) until completely dried.
- Safely store the spot in the padded case

supplied. To clean the protective window

- Turn the system OFF.
- Remove the spot from the front of the handpiece.
- Using a cotton pad and an Isopropyl-based alcohol solvent (70%), carefully wipe the protective window on the front of the handpiece.
- Using a cotton pad wipe the protective window until completely dried.

9.5 Replacing the filters in the BC/FC applicators

1. Disconnect the applicator from the handpiece.


Figure 9.1: BC/FC applicators

2. By using the Hex key, Unscrew the filter plug


Figure 9.2 unscrewing the filter

3. Disconnect the filter from the filter plug.

Figure 9.3 disconnecting the filter

4. Insert the new filter into the filter house.

Figure 9.4 Filter house

5. Place the filter plug back in the applicator


Figure 9.5 replacing the filter plug

6. By using the Hex key, screw in the filter plug.

9.6 Replacing the Water Filter and Deionized Water

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Note

- Water Filter and deionized water should be replaced once a year.
- In clinics that work in high volume, you may need to replace Water Filter and deionized water sooner.

To view next Water Filter replacement date, see Maintenance in section 6.7

10 Data for the service of a medical device

The V30 system is equipped with a self-testing software routine that continuously monitors system operation. The following troubleshooting guide does not attempt to list all possible system failures.

Contacts in Russian Federation:

ООО «AlfaSpa Med», Russia, 109240, Moscow, Moscvoretskay nab. 2a, build 1.

Tel: 8(495)775-26-50, fax: 8(495)937-88-95.

e-mail: alfaspa@alfaspa.ru, www.alfaspa.ru

10.1 Fault Messages

The following table offers a list of the most probable system errors and instructions on how to solve them. If the corrective actions listed in the table do not solve the problem, please contact your local distributor for further assistance.

Message on screen	Probable Cause	Corrective Action	Problem solved	Continued action
Bad coupling	Bad coupling between the electrodes and the skin (V-FR, V-ST, V-Form & V-	Make sure that all electrodes are in full contact with the skin	Yes	Resume work
			No	Contact the local distributor
Overheating	System is overheated	Wait a few minutes until system has cooled down and the message disappears from the screen.	Yes	Resume work
			No	Contact the local distributor
No Handpiece Connected	Improper connection of the handpiece	Disconnect and reconnect the handpiece	Yes	Resume work
			No	Contact the local distributor
No V-Form applicator	Only V-Form handpiece is connected to	Connect the	Yes	Resume work

connected. Please connect V-Form applicator to resume	the system and no applicator is connected to the handpiece	desired applicator to the handpiece	No	Contact the local distributor
No Tip Connected	When FR Treatment is selected and there is no treatment tip connected	1. Connect the V-FR Tip 2. If there is already a tip connected, please replace it with a new one	Yes	Resume after a new tip has been connected
			No	Contact the local distributor
FR Tip Failure	The V-FR handpiece doesn't recognize the Tip	Replace the Tip	Yes	Resume work after new tip is
			No	Contact the local distributor

Filter is Disconnecte d	The filter is not connected to the V- IPL handpiece	Connect a Filter	Yes	Go To Filter confirmation
			No	1. Connect another Filter 2. Contact the local distributor
Spot size has been changed	When Nd:YAG laser Treatment is selected and spot	Confirm the change	Yes	Go to Setting screen
			No	NA
Spot is Disconnecte d	When Nd:YAG laser Treatment is selected from the main menu and there is no spot connected	Connect a Spot	Yes	Go to Spot Confirmation
			No	- Connect another spot - Contact the local distributor
Foot Pedal is Disconnecte d	Foot Pedal is disconnected in V- Laser and V- IPL handpieces	Connect the foot pedal	Yes	Resume work
			No	Contact the local distributor
Foot Pedal Pressed	If the pedal is pressed and user want to proceed to Ready Mode (all handpieces) If the pedal is pressed more than 15 second in V-FR, V-ST and V- Form handpieces	Release the Foot Pedal	Yes	Resume work
			No	1. Reconnect the Pedal (Confirm that the system is off) 2. Contact the local distributor
Handpiece trigger is Pressed	If the trigger is pressed and user want to go Ready Mode (all handpieces except	Release the Trigger	Yes	Resume work
			No	Contact the local distributor
Interlock is disconnecte d	Interlock is disconnected in V- Laser and V- IPL handpieces	Connect the Interlock	Yes	Resume work
			No	Contact the local distributor
Configuration or license code is incorrect	User tried to insert Incorrect code	Verify the code is entered correctly	Yes	Resume work
			No	Contact the local distributor
Configuration or	User tried to insert	Enter a new code	Yes	Resume work

License code has already been entered into this system	Used code		No	Contact the local distributor
Water System Failure – Failure code 401	Handpiece is not connected properly to svstem	Reconnect the V-IPL or V-Laser Handpiece	Yes	Resume work
			No	Contact the local distributor
General Main Board Failure – Failure code 106	General failure	Reboot the system.	Yes	Resume work
			No	Contact the local distributor

Table 10.1: Fault Messages

10.2 Failure Messages

The following table offers a list of Technical failure message that will appear under various, faulty conditions. The message will be accompanied by Error Number. Upon receiving one of the following warning messages, please contact your local distributor for more information and repair:

Message on screen	Error Number
General Main Board Failure	101 – 108
General RF Board Failure	201 - 204
Switch Module Failure	301 - 306
Water System Failure	401 - 404
General Handpiece Failure	501 - 513

Table 10.2: Failure Messages

10.3 Operation troubleshooting

The following table displays a list of operation failures that may occur without warning under various faulty conditions and instructions on how to solve them. If the corrective actions listed in the table do not solve the problem, please contact your local distributor for further assistance.

Message on screen	Corrective Action	Problem solved	Continued action
System is not turned on	A. Ensure proper connection to mains outlet B. Make sure the emergency switch is not pressed	Yes	Resume work
		No	Contact the local distributor
No sound	Make sure the volume setting on the setup screen is not set to minimum.	Yes	Resume work
		No	Contact the local distributor
Noisy system	Reboot the system.	Yes	Resume work

		No	Contact the local distributor
Water leak from the system	Contact the local distributor	-	-
		-	-
Weak Vacuum in V- Form applicators (BC or FC) - Vacuum is not working properly	Replace the filter in the BC or FC applicator	Yes	Resume work
		No	Contact the local distributor

Table 10.3: Operation troubleshooting

11 About Viora

Viora, a leading provider of medical aesthetic solutions, is passionate about helping customers deliver the most effective results to their patients. Since 2005, Viora has sought to improve upon treatment experiences by enhancing existing technology to give customers the best non-invasive solutions available.

Viora’s non-invasive and minimally invasive products treat a wide range of symptoms and are the perfect solutions to address today’s most commonly requested treatments, such as skin tightening, body & facial contouring, cellulite reduction, skin rejuvenation, acne clearance, skin resurfacing, stretch mark reduction, hair removal, skin lightening and more.

Viora’s skilled and knowledgeable clinical and technical teams, working with clinicians around the world have produced innovative, proprietary technologies such as CORE™ and SVC™.

These technologies have revolutionized non-invasive treatments to improve control and increase opportunity for practitioners, while enhancing patient experience. In addition, this in-house expertise has been at the forefront of implementing new developments in response to market demand, improving upon treatments using existing technologies such as laser and IPL. By formulating exclusive combination protocols and combining different proprietary technologies, Viora achieves enhanced clinical results and greater customer and patient satisfaction.

APPENDIX A – EMC Tables

Table 1 – Guidance and MANUFACTURER’S declaration – ELECTROMAGNETIC EMISSIONS →
for all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacturer’s declaration – electromagnetic emissions		
The V30 is intended for use in the electromagnetic environment specified below. The customer or the user of the V30 should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The V30 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	The V30 is suitable for use in all establishments other than domestic, and may be used in domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, provided the following warning is heeded:  Warning: This equipment/system is intended for use by healthcare professionals only. This equipment/system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the [ME EQUIPMENT or ME SYSTEM] or
Harmonic emissions IEC 61000-3-2	Compliance	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Table 2 – Guidance and MANUFACTURER’S declaration – electromagnetic IMMUNITY – for all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacturer’s declaration – electromagnetic immunity			
The V30 is intended for use in the electromagnetic environment specified below. The customer or the user of the V30 should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_T (>95 % dip in U_T) for 0,5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 sec	>95 % dip for 10 ms 60 % dip for 100 ms 30 % dip for 500 ms 95 % dip for 5000 ms	Mains power quality should be that of a typical commercial or hospital environment. If the user of the V30 requires continued operation during power mains interruptions, it is recommended that the V30 be powered from an uninterruptible power

			supply.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE - U_T is the AC mains voltage prior to application of the test level.

**Table 3 – Guidance and MANUFACTURER’s declaration – electromagnetic IMMUNITY –
For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING**

Guidance and manufacturer's declaration – electromagnetic immunity			
The V30 is intended for use in the electromagnetic environment specified below. The customer or the user of the V30 should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 TEST LEVEL	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 Vrms 150 kHz to 80 MHz 3 V/m 80 MHz to 2,5 GHz	3 Vrms 150 kHz to 80 MHz 3 V/m 80 MHz to 2,5 GHz	Portable and mobile RF communications equipment should be used no closer to any part of the V30, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = \sqrt{\frac{P}{k}}$ $d = \sqrt{\frac{P}{k}}$

			$[-] \sqrt{\quad} \quad [-] \sqrt{\quad}$ 800 MHz to 2,5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the V30 is used exceeds the applicable RF compliance level above, the [MEEQUIPMENT or MESYSTEM] should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the V30. ^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V1] V/m.			

Table 4 – Recommended separation distances between portable and mobile RF communications equipment and the ME EQUIPMENT or ME SYSTEM – for ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING			
Recommended separation distances between portable and mobile RF communications equipment and the V30			
The V30 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the V30 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the V30 as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter [W]	Separation distance according to frequency of transmitter [m]		
	150 kHz to 80 MHz $[\text{---}] \sqrt{\text{---}}$	80 MHz to 800 MHz $[\text{---}] \sqrt{\text{---}}$	800 MHz to 2,5 GHz $[\text{---}] \sqrt{\text{---}}$
0,01	0.12	0.12	0.23
0,1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.69	3.69	7.38
100	11.67	11.67	23.33
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

12. Data for recycling of medical device

The V30 Fractional treatment involves use of disposable tips composed of:

- Rigid PCB matrix.
- Electronic equipment (small IC, capacitors)
- Biocompatible frame made from plastic.

	Dispose of according to WEEE regulations
---	--

13. System composition

I. Main composition

1. Device Viora V30 -1 pcs.
2. Power cable – 1 pcs.
3. Manual– 1 pcs.

II. Accessories

1. Nozzle V-ST for lifting -1-10PCs
2. Nozzle V-FR for fractional procedure -1-10 pcs.
3. Disposable sterile pipette nozzle V-FR - 1-10 pcs.
4. Box for disposable sterile pipette nozzle V-FR - 1-10 pcs.
5. Protective cap for Nozzle V-FR - 1-5 pcs.
6. Nozzle V-IPL- 1-10 pcs.
7. V-IPL Filter HR-630- 1-10 pcs.
8. V-IPL Filter HR-570- 1-10 pcs.
9. V-IPL Filter SR-580- 1-10 pcs.
10. V-IPL Filter VP-530- 1-10 pcs.
11. V-IPL Filter AC-415- 1-10 pcs.
12. Applicator- Spot 2,4cm² for Nozzle V-IPL- 1-10 pcs.
13. Box for filters V-IPL- 1-10 pcs.
14. Protective cap for Nozzle V-FR - 1-10 pcs.
15. Nozzle V-Laser Nd:YAG - 1-10 pcs.
16. Applicator- spot V-Laser Nd:YAG 2x4mm²- 1-10 pcs.

- 17.** Applicator- spot V-Laser Nd:YAG Ø 3.5- 1-10 pcs.
- 18.** Applicator- spot V-Laser Nd:YAG Ø 5.5- 1-10 pcs.
- 19.** Applicator- spot V-Laser Nd:YAG Ø 9.5- 1-10 pcs.
- 20.** Box for applicator- spot V-Laser Nd:YAG- 1-10 pcs.
- 21.** Pedal- 1-5 pcs.
- 22.** Set for pouring water (bag -1 pcs., tube-1 pcs., drain tube -1 pcs., metal ring-1 pcs.)-1-5set.
- 23.** The security key- 1-5 pcs.
- 24.** Mount for Nozzles- 1-10 pcs.
- 25.** Protective cap for front connector- 1-5 pcs.
- 26.** Brush for Nozzles- 1-5 pcs.
- 27.** Box for Nozzles- 1-10 pcs.
- 28.** Goggles- 1-10 pcs.

Numbered, bound and sealed
Viora Ltd

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WARNING

INVISIBLE LASER RADIATION
AVOID EYE OR SKIN
EXPOSURE TO DIRECT OR
SCATTERED RADIATION
CLASS 4 LASER PRODUCT

IEC 60825-1,2007

Wavelength: 1064nm
Maximum Output: 500 J/cm²
Pulse Duration: 5-25 msec

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RISK GROUP 3

WARNING

The light emitted may result
in eye injury. Do not look
at the light source.

CAUTION

IR emitted from this device
may cause eye irritation.
Do not stare at the light source.

IEC 60601-2-57

Wavelength: 415-1200 nm
Maximum Output: 35 J/cm²
Pulse Duration: 4-50 msec.

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