Cory Marsh
Manager, Regulatory Operations
Merit Medical Systems, Inc.
30 June 2020

Merit Laureate Hydrophilic Guide Wire

Operational Documentation

MERIT LAUREATE

HYDROPHILIC GUIDE WIRE

INSTRUCTIONS FOR USE

Carefully read all warnings, precautions and directions prior to use. Failure to do so may result in the improper use of this device which could cause the following complications:

- Shearing of the hydrophilic guide wire
- Release of plastic pieces or fragments from the hydrophilic guide wire which may need to be retrieved from the vasculature.
- Vessel trauma

Description:

Merit Medical hydrophilic guide wires are constructed from a high quality, steerable, metallic (nitinol) core wire with a polymer coating. A hydrophilic coating is applied over the radiopaque polymer jacket. Guide wires are supplied sterile and non-pyrogenic.

Indications for use:

The Merit Hydrophilic Guide Wire is intended to be used in the peripheral vascular system to facilitate the placement of devices during diagnostic and interventional procedures.

Contraindications:

These guide wires are not intended for Percutaneous Transluminal Coronary Angioplasty use.

Warnings/Adverse Reactions:

- Inspect wire for damage prior to use, do not use a wire that has been bent, kinked, or damaged. Use of a damaged wire may result in vessel damage or wire fragment release into the vessel.
- Do not reshape the hydrophilic wire by any means. Attempting to reshape the wire may cause damage to the wire.
- Do not manipulate or withdraw the wire through a metal entry needle or a metal dilator, or use this wire with devices which contain metal parts such as atherectomy catheters or laser catheters or metal torque devices. This may result in destruction and/or separation of the outer polyurethane coating requiring retrieval. A plastic entry needle is recommended when using this wire for initial placement, or a catheter, introducer sheath or vessel dilator should replace the needle as soon as the guide wire has been inserted into the vessel.
- Never advance the guide wire against resistance without first determining the reason for the resistance under fluoroscopy. If resistance occurs and the cause of resistance cannot be determined, remove the guide wire and device as a unit. Excessive force against resistance may result in damage to the wire and/or to the vessel.
- When manipulating, advancing, exchanging, or withdrawing a catheter over the wire, secure and maintain the guide wire in place under fluoroscopy to avoid unexpected guide wire advancement; otherwise damage to the vessel wall by the wire's tip may occur.
- The hydrophilic guide wire should be used only by a physician, who is well trained in manipulation and observation of guide wires under fluoroscopy.

Other potential adverse reactions which may result from the improper use of a guide wire include, but

are not limited to:

- | | |
|---|-----------------------|
| • Thrombus | • Emboli |
| • Arterial or venous vessel wall damage | • Plaque dislodgment |
| • Hematoma at the puncture site | • Infection |
| • Vessel perforation | • Vessel spasm |
| • Hemorrhage | • Vascular thrombosis |

Precautions:

- When using a drug or a device concurrently with the wire, the operator should have a full understanding of the properties/characteristics of the drug or device so as to avoid damage to the hydrophilic guide wire.
- Use care when manipulating this guide wire through a tightened Hemostasis valve.

Cautions:

- At least 5 cm of the wire should protrude from the device hub at all times to prevent the wire from sliding entirely into the device due to the low sliding friction of this wire.
- Merit Medical hydrophilic guide wires are packaged in a plastic hoop fitted with a luer hub. This packaging is provided to facilitate compliance with the manufacturer recommended guidelines that the wire must be flushed with saline or heparinized saline prior to use (See instructions for use).
- Contents of unopened, undamaged package are sterile and non-pyrogenic.

Preparation for use:

1. The surface of the hydrophilic guide wire is not lubricious unless it is wet. Before attempting to remove the guide wire from its dispenser, inject sterile heparinized saline solution into the luer lock hub end of the dispenser to fill the dispenser coil. This will completely cover the guide wire surface, activate the hydrophilic coating, and make the guide wire very lubricious.

Warning: Failure to hydrate the dispenser hoop prior to guide wire removal may result in guide wire damage and/or difficult removal from the dispenser.

2. After hydrating the guide wire, gently grasp the J-straightener device and pull from the dispenser, once the straightener is separated from the dispenser, continue to remove the wire from the hoop.
3. If guide wire is not properly hydrated, it will be difficult to remove from the dispenser. Inject additional heparinized saline solution into dispenser and repeat step # 2.

Instructions for use:

1. Fill concurrent device with heparinized saline solution before and during use to ensure smooth movement of the hydrophilic guide wire within the device.
2. Use of sterilized gauze moistened with heparinized saline solution and/or a non metal torque device will facilitate handling of the wire.
3. Insert the guide wire into the device and advance to the desired location.

Warning: If movement of the wire within the device becomes diminished, remove guide wire and reactivate the hydrophilic coating by wetting its entire surface with a heparinized saline solution.

4. Wipe the guide wire with a 4x4 gauze moistened with heparinized saline solution to remove excess blood from the guide wire surface.

Warning: Do not use dry gauze as this may damage the guide wire surface resulting in increased resistance when the wire is reinserted into the device.

5. Re-hydrate the guide wire prior to reinsertion into any device or placement into a patient. If additional resistance is felt after re-hydration, replace guide wire.

6. Use of alcohol, antiseptic solutions, or other solvents must be avoided.

Warning: These solutions may adversely affect the surface of the hydrophilic guide wire.

7. After cleaning the wire, place into the saline filled hoop, proximal end first. The wire may also be placed in a guide wire basin and completely covered with heparinized saline solution.

Warning: Hydrophilic guide wires must be kept hydrated through out the entire procedure. Re-hydrate as necessary when the surface starts to dry out.

Px Only -Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

REUSE PRECAUTION STATEMENT


For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

ANNEX TO INSTRUCTIONS FOR USE

DEVICE SPECIFICATIONS

The wires are offered in the following configurations:

- 0.018" (0.46 mm) outer diameter: standard shaft, straight and angled tip configurations with lengths ranging from 80cm – 260cm.
- 0.025" (0.64 mm) outer diameter: standard shaft, straight and angled tip configurations with lengths ranging from 150cm – 260cm.
- 0.035" (0.89 mm) outer diameter: standard and stiff shafts, straight and angled tip configurations with lengths ranging from 80cm – 260cm.
- 0.038" (0.97 mm) outer diameter: standard and stiff shafts, straight and angled tip configurations with lengths ranging from 80cm – 260cm.

Figured Jacketed & Coated Product

MATERIAL OF CONSTRUCTION

In the figure below, the Laureate guidewire is presented to clarify the various components and materials.

COMPONENT	Material Specification
Core Wire	Nitinol
Jacket	Tungsten filled polyurethane
Coating	Hydrophilic Co-polymer

PACKAGING

The Merit Laureate is sold as a stand-alone device. The wire is placed inside a multiple loop flush dispenser, also referred to as a hoop. The dispenser hoop will have a female Luer connection which facilitates solution flushing through the hoop to hydrate the guide wire. The wire is placed so that the distal end of the wire comes out of the outer portion of the dispenser hoop. A wire holder will be placed on the other end of the hoop to introduce the wire into the catheter. The hooped guide wire assembly is placed in a Tyvek peel pouch that is sealed using a validated heat seal process. Five pouches are placed in a white box and completed with IFU. These white boxes can be placed inside a shipper carton/box

PRODUCT SHELF LIFE

The Merit Hydrophilic Guide Wire is labeled for a three (3) year shelf life.

Merit guarantees the compliance of the product with the requirements of regulatory documents during this period keeping the storage and transportation conditions and if the package is not damaged.

TRANSPORTATION, STORAGE AND USAGE

Transportation should be conducted in accordance with the agreed contract between the two sides
Storage conditions

Parameter	Value
Temperature	From +5°C to +27°C
Relative humidity	From 10% to 90%
Pressure	84.0 – 106.7 kPa (630 – 800 mm Hg)

Transport environment conditions

Parameter	Value
Temperature	From -35°C to +55°C (environment)
Relative humidity	From 10% to 90%
Pressure	84.0 – 106.7 kPa (630 – 800 mm Hg)

Usage:

Temperature Range: 25 ± 10 ° C

Relative air humidity: <80%

Atmospheric pressure: 84.0 -106.7 kPa (630-800 mm Hg).

Contact with the internal environment of the body at a temperature of 32°C to +42°C

Use Conditions – in a medical and preventive treatment facility

DISPOSAL

Dispose the devices according to local environmental requirements.

Class of hazard for medical wastes by reference to specification of morphological composition (in accordance with Russian Sanitary Regulations and Standards 2.1.7.2790-10) Class B

WARRANTY

WARRANTY

Limited Warranty and Disclaimer

Merit represents and warrants that the products shall materially conform to the specifications at the time of shipment to Customer and will remain materially free of defects in materials and workmanship for the period specified on the written limited warranty separately provided with each product. This warranty shall not apply to products that have been (i) modified, changed, repaired, refurbished, reprocessed, or altered by anyone other than Merit; (ii) subjected to misuse, mishandling, re-use, re-sterilization, accident, abuse, neglect or tampering; (iii) damaged by excessive physical, environmental or electrical stress; or that have had a serial number altered, defaced or removed; (iv) combined with any other product; (v) used outside the approved "Indications for Use" as cleared by the relevant competent authority, used contrary to the use outlined in the product specifications or in an application or environment for which such product was not designed or contemplated. Any technical product advice furnished by Merit is provided without duty or compensation and Merit assumes no obligation or liability, all such advice being given and accepted at Customer's risk. EXCEPT AS EXPRESSLY SET FORTH HEREIN MERIT MAKES NO REPRESENTATION OR WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, TO CUSTOMER, ITS END-USERS OR ANY THIRD PARTIES WITH RESPECT TO THE PRODUCTS. MERIT DISCLAIMS AND EXCLUDES ANY AND ALL OTHER EXPRESS, IMPLIED, AND STATUTORY WARRANTIES, INCLUDING WITHOUT LIMITATION, WARRANTIES AGAINST INFRINGEMENT, AND THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

Limitation of Liability

UNDER NO CIRCUMSTANCES SHALL MERIT BE LIABLE TO CUSTOMER OR ANY OTHER PERSON OR ENTITY FOR ANY PUNITIVE, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, WHETHER BASED UPON LOST GOODWILL, LOST RESALE PROFITS, WORK STOPPAGE, PRODUCT FAILURE, IMPAIRMENT OF OTHER PRODUCTS OR OTHERWISE AND WHETHER IN A CONTRACT, TORT OR OTHER ACTION FOR OR ARISING OUT OF BREACH OF WARRANTY, BREACH OF CONTRACT, DELAY, NEGLIGENCE, STRICT LIABILITY, OR OTHERWISE.

List of Product Label Part Numbers

	This symbol means that the device fully complies with European Directive 93/42/EEC.
	Manufacturer
	Authorized Representative

	Do not withdraw through metal devices.
	Single use
	Reference Number
	Lot Number
	Sterilization: Ethylene-Oxide Gas
	Caution: U. S. federal law restricts this device to sale by or on the order of a physician.
	Use by
	Caution: Consult Accompanying Documents
	Do not use, if package is damaged
	Non-pyrogenic

Customer Service

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095 U.S.A.
+ 1 801 253 1600

Merit Technologies LLC
119017 Russia, Moscow
1st Kazachii per. h.7 floor 1 office 6 ,
+ 7 495 221 89 02

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