August 1, 2016

Warranty Card (draft)

Medical Device
«Flexible fiber-optic nasopharyngolaryngoscopes»

Registration certificate No. _____

Model _____

Serial No. _____

Manufacturer: HUGER Medical Instrument Co., Ltd.

Address: No. 1328, Yaobei Road, Jiuting, Songjiang District, 201615, Shanghai China

An organization authorized by the manufacturer to accept and discharge consumers' claims regarding inadequate quality of goods: Limited Liability Company "Healthy World" 190000, Saint Petersburg, Pochtamtskaya ul., 10, lit. A, Pom. 1-N, tel: +7 (812) 385-74-55.

Manufacture date _____

Sale date _____

Buyer _____

Seller's representative _____

Signature

Buyer's

Address/Tel. _____

Warranty liability:

Seller undertakes to repair or replace the product within the warranty period charge-free, provided the Buyer follows storage, transportation and operation terms.

The warranty period is 12 months.

The warranty period is calculated from the sale date specified herein.

Signature:

Name: Fulong Jin

No. 1328, Yaobei Road, Jiuting, Songjiang District, 201615 Shanghai China

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
User Manual

Flexible fiber-optic nasopharyngolaryngoscopes

No. 1328, Yaobei Road, Jiuting, Songjiang District, 201615, Shanghai China

Phone: +86-21-67691764 Fax: +86-21-67691721

E-mail: inform@huger.cn Website: <http://www.huger.cn>

README

This manual extends to the model range of flexible fiber-optic nasopharyngolaryngoscopes (hereinafter referred to as the nasopharyngolaryngoscopes), manufactured by (HUGER Medical Instrument Co., Ltd.)

Area of use – referral medical and preventive treatment facilities.

Nasopharyngolaryngoscopes are manufactured in the following embodiments: FN-32A, FN-38A, FL-39A, FN-50A.

Nasopharyngolaryngoscopes are designed for visual examination and treatment of structures inside of nasal passages, including sinus openings, pharyngeal end of ear trumpet, larynx and vocal cords. Don't use nasopharyngolaryngoscopes for tasks, which are not included into the manual.

Please, read carefully this manual before commencement of work. The manual contains information on proper use and precautions while working with nasopharyngolaryngoscopes. Nasopharyngolaryngoscopes are very fragile, so proper treatment and maintenance according with the present manual, will considerably decrease the necessity of costly repair and extend the life of your nasopharyngolaryngoscope.

This manual describes recommended procedure of preparation and check-out of the nasopharyngolaryngoscope before commencement of work. The manual doesn't describe the procedure of endoscopic examination.

If you have any questions related to the content of this manual and safety regulations of nasopharyngolaryngoscopes, please contact the representative of Huger company.

PRECAUTIONS

Please, check the availability of nasopharyngolaryngoscope and all its accessories according to the completing units according to Section 3. If nasopharyngolaryngoscope its accessories are not available, or in case of any defect, get in touch with the representative of huger company. For examination of purpose and name of the nasopharyngolaryngoscope accessories, read section "Technical specifications", "Main specifications" and "Name of parts". Read carefully about the procedure of preparation, treatment and check-out of the nasopharyngolaryngoscope.

Nasopharyngolaryngoscope should be cleaned and disinfected before first usage.

Nasopharyngolaryngoscope and its accessories should be taken out from the case for transportation and it should be kept in accordance with requirements of Section 7. Transportation case is not designed for the storage of the nasopharyngolaryngoscope. Keep the transportation case for transportation of the nasopharyngolaryngoscope.

BEFORE COMMENCEMENT OF WORK

WARNING

Nasopharyngolaryngoscope is exacting product. Its construction includes many developments for patient's security. In particular, the system of distal end bend is designed in such way to provide mild reaction and maximum bend angle of the distal end while applying the optimal effort to the control handle of the distal end bend angle. Application of excessive force to the control handles can cause nasopharyngolaryngoscope breakdown and patient's injury. Before introduction of the nasopharyngolaryngoscope into patient's body, make sure that the control handle of the distal end is blocked in "F" position (without restraint), and distal end moves without resistance. In case of any resistance while introducing or checking out of the control mechanism of the distal end, stop using the nasopharyngolaryngoscope. Contact the representative of HUGER.

WARNING

Leakage test is required after every usage of the nasopharyngolaryngoscope. If water leakage is detected, don't use the nasopharyngolaryngoscope and contact the representative of HUGER. Otherwise, you may have serious problems and breakdown of the nasopharyngolaryngoscope is possible.

SYMBOLS

	Marking CE on this nasopharyngolaryngoscope indicates that this nasopharyngolaryngoscope was tested and complies with the requirements of Directive 93/42/ECC MDD
	Type BF
	Read User manuals
	Manufacturing date
	Manufacturer
	Authorized representative in the European Union
	Series number
IPX7	Degree of protection from water penetration

CONTENTS

1. PURPOSE AND MAIN TECHNICAL SPECIFICATIONS	5
2. NAMES OF PARTS	9
3. DELIVERY SET	12
4. OPERATION CONDITIONS	13
5. CHECK-OUT AND PRESTARTING PROCEDURES	13
6. APPLICATION OF THE NASOPHARYNGOLARYNGOSCOPE	16
7. CARE AND MAINTENANCE	17
8. TROUBLESHOOTING	23
9. MANUFACTURER WARRANTY	24
10. DISPOSAL CONSIDERATIONS	24
 ANNEX 1	 25

1. PURPOSE AND MAIN TECHNICAL SPECIFICATIONS

1.1 Purpose

Nasopharyngolaryngoscopes are designed for visual examination and treatment of structures inside of nasal passages, including sinus openings, pharyngeal end of ear trumpet, larynx and vocal cords.

Indications for use: Necessity of visual examination and treatment of structures inside of nasal passages, including sinus openings, pharyngeal end of ear trumpet, larynx and vocal cords.

Contraindications: n/a

Side effects: possible nasal bleeding when insertion of Nasopharyngolaryngoscope through the nose.

Nasopharyngolaryngoscope can be totally dipped into the water for cleaning and disinfection.

1.2 Main technical specifications

Item		Parameters				Units
		FN-32A	FN-38A	FL-39A	FN-50A	
Direction and observation angle		0±3	0±3	0±3	0±3	grade
Angle of vision		120±5	120±5	120±5	120±5	grade
Range of working distances		2-50	2-50	2-50	2-50	mm
Working section length		300±5	300±5	600±5	450±5	mm
Total length		No more than 600	No more than 600	No more than 900	No more than 750	mm
Working section width		3,2 _{-0,1}	3,8 _{-0,1}	4,1 _{-0,1}	5,0 _{-0,1}	mm
Distal end bend angle		Upwards 130±5, Downwards 130±5	Upwards 130±5, Downwards 130±5	Upwards 130±5, Downwards 130±5	Upwards 180±5, Downwards 130±5	grade
Width of forceps aperture		Non	Non	1,2 ^{+0,05}	2,0 ^{+0,05}	mm
Resolving power	At working distance 2 mm	4,2				mm ⁻¹
	At working distance 9 mm	3,1				mm ⁻¹
	At working distance 50 mm	0,3				mm ⁻¹
Net weight		No more than 800				g
Illumination of the object while using portable lantern at the calculated working distance (9,0±0,1) mm		No more than 3000				lux

Time of operation with fully charged battery		No more than 90	minute
Power voltage of the portable lantern		3,7±0,1	V
Battery	Voltage	3,7±0,1	B
	Capacity	3000±50	mA•h
	Time of full charge	No more than 4	hour
Electric power supply of the battery charger		220/50	V/Hz
Output voltage of the battery charger		4,2±0,1	V
Power consumption while battery charging		No more than 4,5	W
Maximum pressure of tester for leakage test		34±0,5	kPa
Cleaning brush	Length	1200±10	mm
	Diameter	1,8±0,1	mm
Mouthpiece	Ширина	50±0,5	mm
	Width	27±0,5	mm
	Length	40±0,5	mm
	Inner dimension of the mouthpiece window	(30 x 17)±0,5	mm
Length of irrigator of channels		1200±10	mm
Diameter of adapter branch for processing		5±0,1	mm
Length of light cable		1800±30	mm
Connector of the light cable connected to the nasopharyngolaryngoscope		STORZ	-
Connector of the light cable connected to the external lantern		STORZ	-
Case for transportation	length	600±3	mm
	width	470±3	mm
	height	160±3	mm

1.3 Security specifications

Category of product CE	IIa
Classification of protection against electric shock hazard	Internal power supply
Degree of protection against electric shock hazard	Working part of BF type
Degree of protection against explosion hazard	Out of service with deflagrating mixtures, anesthetics with air, oxygen or laughing gas.
Degree of protection against liquid penetration	Ordinary (without protection against water penetration), working part IPX7
Category according to disinfection methods	Methods of disinfection of equipment are recommended by the manufacturer.
Working hours	Continuous running
Operating life	5 years

1.4 Used materials.

o	Name	Used material	Material that contacts with the medical staff and (or) patient
1	Nasopharyngolaryngoscope*	Aluminum alloy 5251, Stainless steel AISI 316, polyvinylchloride SG-3	Stainless steel AISI 316, polyvinylchloride SG-3
1.1	Cover	Powder polyether paint PPGBH-635	Powder polyether paint PPGBH-635
2	Cleaning brush	polyvinylchloride SG-3, Stainless steel AISI 316, nylon K180	polyvinylchloride SG-3, Stainless steel AISI 316, nylon K180
3	Mouthpiece	Silicone LSG1001-08A/B	Silicone LSG1001-08A/B
4	Tester for leakage test	Polyvinylchloride SG-3, Stainless steel AISI 316 AISI 316	Polyvinylchloride SG-3, Stainless steel AISI 316 AISI 316
5	Irrigator of channels	Stainless steel AISI 316, silicone LSG1001-08A/B	Stainless steel AISI 316, silicone LSG1001-08A/B
6	Adapter for processing	ABS AF 312, silicone LSG1001-08A/B	ABS AF 312, silicone LSG1001-08A/B
7	Exhaustion valve	Polyvinylchloride SG-3	Polyvinylchloride SG-3
8	Forceps aperture valve	Polyvinylchloride SG-3	Polyvinylchloride SG-3
9	ETO cap	Aluminum alloy 5251	No contact
9.1	Cover	Powder polyether paint PPGBH-635	Powder polyether paint PPGBH-635
10	Eye shade	ABS AF 312	ABS AF 312
11	Eyeglass lead	ABS AF 312	ABS AF 312
12	Light cable	Stainless steel AISI 316, polyvinylchloride SG-3	Stainless steel AISI 316, polyvinylchloride SG-3
13	Safety cap of the light cable	polyvinylchloride SG-3	polyvinylchloride SG-3
14	Adapter connector for connection of the light cable to the external lantern	Stainless steel AISI 316	Stainless steel AISI 316
15	Adapter for the light cable with STORZ connector	Stainless steel AISI 316	Stainless steel AISI 316
16	Adapter for the light cable with WOLF connector	Stainless steel AISI 316	Stainless steel AISI 316
17	Portable lantern	Stainless steel AISI 316	Stainless steel AISI 316
18	Charging device	ABS AF 312	ABS AF 312
19	Battery	Aluminum alloy 5251, steel A572	No contact
20	Case for transportation	ABS AF 312	ABS AF 312
21	Keys from case for transportation	ABS AF 312, steel A572	ABS AF 312, steel A572

*- embodiment is chosen by the Customer from represented models: FN-32A, FN-38A, FL-39A, FN-50A.

1.5 Precautions.

- The nasopharyngolaryngoscope doesn't contain parts that can be fixed by the operator. Maintenance should be performed by technician authorized by the manufacturer.
- The nasopharyngolaryngoscope shouldn't be used at home.
- The nasopharyngolaryngoscope can't be used together with high-frequency surgical devices.
- The nasopharyngolaryngoscope should be checked out every time before using, it doesn't require regular maintenance.
- The nasopharyngolaryngoscope should be cleaned and disinfected before purpose use.
- Because of high intensity of light emission in the process of work, the temperature at the distal end of the nasopharyngolaryngoscope can exceed 41°C. To prevent from too high temperature of local tissues, distal end shouldn't stay at the same place for a long time.
- Do not bend nasopharyngolaryngoscope with force, make sure that you pull, rotate or turn gently the working part, as it may cause the risk of patient's injury.
- Before getting clear picture don't control pulls, don't switch on exhaust or move forward/backward the nasopharyngolaryngoscope, otherwise it may result in patient's injury.
- Don't press the distal end of the working part so that the surface of the objective lens is pressed into tissues, otherwise it will be impossible to get oriented.
- Don't bend and fold the section in a different way, as the nasopharyngolaryngoscope will be damaged.
- Don't press and squeeze bend section with force; otherwise the cover of section will be quickly worn out, deformed or torn off.
- Don't twist the working part, otherwise it will be damaged.
- When life expectancy of the nasopharyngolaryngoscope and its accessories comes to an end, all the equipment must be disposed according to effective local regulations or disposal requirements of the hospital.
- While using the equipment avoid excessive vibration and falls.
- In the process of nasopharyngolaryngoscope functioning doesn't look at the distal end, it may injure your eyes.
- It is strongly recommended to prepare spare parts for clinical use in case of nasopharyngolaryngoscopes breakdown.

2. NAME AND PURPOSE OF NASOPHARYNGOLARYNGOSCOPE'S PARTS

2.1 Name of nasopharyngolaryngoscope's parts


Fig. 1 Name of parts


Fig. 2 Name of parts

2.2 Purpose of the nasopharyngolaryngoscope's parts

- 1) Eyepiece is the element of optical system of the nasopharyngolaryngoscope, facing the user's eye, is designed for image viewing.
- 2) Control handle of distal end angle is used for regulation of working part angle. If the handle is turned "Upwards", flexible section will turn upwards too; if control handle is turned "Backwards", the working part will turn downwards.
- 3) Exhaust valve is designed to regulate aspiration. When exhaust valve is pressed the aspiration is performed (with connected aspirator, which is not supplied with the nasopharyngolaryngoscope).
- 4) Valve of the instrumental channel is designed to open/close instrumental channel of the nasopharyngolaryngoscope. When open it allows fluent introducing of the instrument in the process of functioning. When closed it provides hermiticity of the nasopharyngolaryngoscope. For connection inset the valve into the inlet of the instrumental channel. (for FL-39A , FN-50A)
- 5) Instrumental channel is used for endoscopic examination: collection of material for analysis and work with endoscopic instruments. (for FL-39A , FN-50A)
- 6) Working part serves to control bend of introduced part; distal end of the nasopharyngolaryngoscope can be moved by means of angular change.
- 7) Distal end is working part section situated on the side of the nasopharyngolaryngoscope opposite from the observer.
- 8) Connector for leakage test serves for connection of the tester for leakage test.
- 9) Light cable is designed for transmittance of light stream from external lantern (not delivered with nasopharyngolaryngoscope) to the nasopharyngolaryngoscope. Light cable has STORZ connector from both sides.
- 10) Cleaning brush is designed for cleaning of the instrumental channel from impurities after purpose use of the nasopharyngolaryngoscope. (for FL-39A , FN-50A)
- 11) Mouthpiece is used to direct distal end of the nasopharyngolaryngoscope to the glottal aperture in the process of functioning. Introduction of the mouthpiece in the patient's mouth before the procedure, prevent from bite and/or damage of the insertion tube of the nasopharyngolaryngoscope. Don't use mouthpiece with damages, deformation or other defects.
- 12) Tester for leakage test is designed to perform leakage test of the nasopharyngolaryngoscope. Tester presents a simple device with bulb for pressure pumping, manometer for control and strong flexible pipe with connector at the end. Tester allows detecting defects in flexible structures of the nasopharyngolaryngoscope at early stages, which enables to take timely correctional measures or replacement of defective nasopharyngolaryngoscope.
- 13) Irrigator of channels is designed for injection of solution with the help of syringe. (for FL-39A , FN-50A)
- 14) Adapter for procession is designed for automatic washing of the nasopharyngolaryngoscope with the help of washing machine. (for FL-39A , FN-50A)
- 15) ETO cap is designed for connection to the outlet for leakage test. Placement of ETO cap is necessary for connection of the tester for leakage test performance and for its sterilization. ETO cap should be removed in the process of patient's examination

and submersion. To place ETO cap into the outlet for leakage test, twist it clockwise to the click, to disconnect it twist it off pulling simultaneously.

- 16) Eye shade is designed for isolation for isolation of the observer from contact with metallic parts of the nasopharyngolaryngoscope.
- 17) Eyepiece lead is designed for protection of the nasopharyngolaryngoscope eye glass from dirt and dust.
- 18) Protective cap of the light cable is designed for protection of the light cable from dirt and dust.
- 19) Adaptor for connection of the light cable to external lantern Huger (not included into delivery set)
- 20) Adapter for the light cable with STORZ connection is designed for connection of the light cable with STORZ connector to the nasopharyngolaryngoscope.
- 21) Adapter for the light cable with WOLF connection is designed for connection of the light cable with WOLF connector to the nasopharyngolaryngoscope.
- 22) Portable lantern serves as source of light while using nasopharyngolaryngoscope; it is connected to the outlet of the lantern.
- 23) Charging device is designed for charging of battery.
- 24) Battery is designed for electric power supply of the lantern.
- 25) Case for transportation is designed for transportation of the nasopharyngolaryngoscope.
- 26) Keys from the case for transportation are designed for opening/closing of the case for transportation.

3. DELIVERY SET

No.	Name	Quantity, pieces
1	Flexible fiber-optic nasopharyngolaryngoscope*	1
2	Accessories	
2.1	Cleaning brush	No more than 5
2.2	Mouthpiece	No more than 5
2.3	Tester for leakage test	No more than 3
2.4	Irrigator of channels	No more than 5
2.5	Adapter for processing	No more than 5
2.6	Exhaust valve	No more than 20
2.7	Instrumental channel valve	No more than 20
2.8	ETO cap	1
2.9	Eye shade	1
2.10	Eyeglass lead	1
2.11	Light cable	No more than 3
2.12	Protective cap of the light cable	No more than 3
2.13	Adaptor for connection of the light cable to the external lantern	1
2.14	Adapter for the light cable with STORZ connection	1
2.15	Adapter for the light cable with WOLF connection	1
2.16	Portable lantern	1
2.17	Charging device	1
2.18	Battery	1
2.19	Case for transportation	1
2.20	Keys from the case for transportation	No more than 2
3	User manual	1

*- embodiment is chosen by the Customer from represented models: FN-32A, FN-

38A,
FL-39A, FN-50A.

4. OPERATION CONDITIONS

4.1 Operation conditions of the nasopharyngolaryngoscope:

Temperature	10 – 35°C
Relative humidity	30% – 85%
Atmosphere pressure	700 hPa – 1060 hPa

5. CHECK-OUT AND PRESTARTING PROCEDURES

5.1 Preparation of the nasopharyngolaryngoscope

Take the nasopharyngolaryngoscope out of the case for transportation, take off the eye glass lead and put on the eye shed.

5.2 Preparation and check-out of the lantern

5.2.1 Use of the nasopharyngolaryngoscope with portable lantern, included into delivery set.

5.2.1.2 Put battery into portable lantern observing polarity according to marking.

5.2.1.3 Connect portable lantern to the nasopharyngolaryngoscope using connector for the lantern. 5.2.1.4 Check out the lantern by pressing start button.

5.2.2 Use of the nasopharyngolaryngoscope with external lantern (not included into delivery set)

5.2.2.1 Connect light cable to the nasopharyngolaryngoscope with the help of adapter with STORZ connector, and the other end to the lantern with STORZ connector, if external lantern of Huger company is used (not included into delivery set), use adaptor for connection of the light cable to the external lantern).

5.2.2.2 If light cable with WOLF connector is used, connection to the nasopharyngolaryngoscope is done with the help of adaptor with WOLF connector.

 **Warning:** Working part of the external lantern must be BF or CF type.

5.3 Check-out of exhaust

- (1) To check out the exhaust connect exhaust valve to the aspirator, (not included into delivery set) and set suction of the aspiration no less than 50 kPa.
- (2) Connect exhaust valve to the aspirator and to the aspiration channel on the handle of the nasopharyngolaryngoscope.
- (3) Switch on aspirator. Plunge distal end into the water, regulation of exhaust is done by means of exhaust valve. (see fig.3)


Fig. 3

5.4 Check-out of endoscopic system

Before each application nasopharyngolaryngoscope should be checked out in accordance with the following procedure. If any malfunction or failure is suspected, stop using nasopharyngolaryngoscope and contact the representative of Huger.

1. Check-out of the working part of the nasopharyngolaryngoscope

- (1) Perform visual examination of the surface of the working part for any dents, prominent parts or other disorders.
- (2) Sweep the whole working part with your fingertips to exclude any prominent parts, inner deformations or other changes.

2. Check-out of the bend section of the working part

Try to manipulate control handle of distal end angle up to possible deviations on every direction. Make sure that the bend section of the working part bends easily and correctly, reaching maximum bends. At the same time check external surface of the bend section of the working part for holes, disruptions or other disorders (see fig.4)

- ▲ Doesn't press or twist bend section of the working part in the process of functioning.

Fig. 4

5.5 Check-out of the optical system

Turn focusing wheel on the eye glass to focus the image. Check out visualization of the object through the eye glass at the distance of 7 mm from objective lens (see fig. 5)


Fig. 5

5.6 Check-out of other parts

- (1) Check out your light cable for any damages (for example, distortions, cracks, and bends).
- (2) Check out light guide connector, control section for loose connections.

▲ Warning: Don't use nasopharyngolaryngoscope and contact the representative of Huger, if slightest defects are determined.

▲ Because of high intensity of light emission in the process of work, the temperature at the distal end of the nasopharyngolaryngoscope can exceed 41°C. To prevent from too high temperature of local tissues, distal end shouldn't stay at the same place for a long time

5.7 Preparation before commencement of work

1. Apply coat of silicone polish (lens cleaner, anti weeping agent).
Apply silicone polish (lens cleaner) on the edge of clean cloth and wipe slightly objective lens and lens of light guide. Remove the excess.

2. Application of lubricant

Lubricate working part with water soluble medical lubricant, accurately applying it to area of the distal end, avoiding application of the lubricant to the distal end.

▲ Don't use olive oil, petroleum or Vaseline based lubricants. Otherwise you can damage nasopharyngolaryngoscope material.

6. APPLICATION OF THE NASOPHARYNGOLARYNGOSCOPE

A person, working with the nasopharyngolaryngoscope must be familiar with clinical endoscopy. This manual describes basic manipulations with nasopharyngolaryngoscope.

6.1 Handling nasopharyngolaryngoscope

Hold nasopharyngolaryngoscope by the handle with left arm and the working part of the nasopharyngolaryngoscope by the right hand. Exhaust valve is regulated with forefinger of the left hand, and control handle of the distal end angle with the thumb

6.2 Introduction and examination

1. Regulate focus with focusing wheel.

Turn regulation wheel to focus the image.

2. Introduce the nasopharyngolaryngoscope.

2.1 The procedure is done in the sitting posture. A patient is sat down into an armchair with a slightly tilted-back head. Then the nasal cavity is anesthetized. It can be done with either a lidocaine-based gel, which is used for the endoscope lubrication, or a spray – anesthetic agent, which is applied directly in the nasopharynx. Then a nasopharyngolaryngoscope is carefully putted into the nose.

Attention! Do not be nervous and do not move during the survey. It is required to relax and breathe calmly and quietly.

2.2 Putting-in via mouth cavity

Put the operating part into the mouth cavity of the patient. Ask the patient in advance to bite the mouthpiece, then put the operating part into the patient's throat under the control of nasopharyngolaryngoscope. When the tip passes through the throat, the patient shall slightly bite the mouthpiece in order to hold it in position during the whole procedure. Exclude biting of the operating part of the nasopharyngolaryngoscope by the patient. Always use a mouthpiece.

3. Regulate lighting while using external lantern.

Regulate lightning on the lantern to achieve comfortable level of lightning. While lighting with xenon lamp, heating of the distal end of the nasopharyngolaryngoscope is possible. There is the possibility of risk of tissues burn as the result of long-term exposure to extensive irradiation. Always use lightning of minimal intensity, but it should be enough for adequate lightning.

4. Control of distal end angle.

Control handle of the distal end angle provides bending of the bend section of the working part of the nasopharyngolaryngoscope. Bending upwards/downwards and fixation of bend position is controlled with the thumb of the hand (see figure on page 9).

If control of the distant end angle is impossible, or there is any breakdown of the control process, please stop examining immediately. Put the handle in

middle position, and then under direct vision slowly take out the nasopharyngolaryngoscope.

5. Use of aspiration function.

By means of exhaust you can extract liquid and excretions. Connect aspirator (not included into the delivery set) to the exhaust valve. Turn on aspirator. Control of the exhaust is realized with the help of exhaust valve.

Avoid suction of mucous into aspiration channel, in this case stop aspiration immediately. Pull the nasopharyngolaryngoscope a bit; prevent the mucous from the contact of the aspiration channel.

6. Use of irrigation channel.

Connect syringe to the end of irrigation channel. Put irrigator of the channel through instrumental channel and inject the solution in the right place (see fig.6)


Fig 6.

6.3 Extraction of the nasopharyngolaryngoscope

Make sure that bend section of the working part is straightened (control handle of the distal end must be in the middle position).

Under direct vision slowly extract the nasopharyngolaryngoscope.

7. CARE AND MAINTENANCE

7.1 Summary

Before sterilization or disinfection the nasopharyngolaryngoscope must be thoroughly washed. Thorough washing can eliminate microscopic and organic materials, which affect disinfection/sterilization.

Nasopharyngolaryngoscopes are totally enclosed, they can be plunged into solution for cleaning and disinfection.

7.2 Cleaning of nasopharyngolaryngoscopes

After each examination, perform the following procedure:

1. Wipe the working part with the cloth and remove mucus.

2. Dip the working part into cleaning solution, wipe slowly with the cloth or sponge.

3. Cleansing, disinfection and sterilization of the exhaust valve.

According to annex 1 of this manual one exhaust valve is used for 6 usages.

4. Cleansing of aspiration channel and instrumental channel

Put the brush through the body of the exhaust valve so that the brush came out from the distal end of the nasopharyngolaryngoscope. Put the brush through the channel several times.

- Put the brush into the outlet of the instrumental channel so that the brush came out of the distal end of the nasopharyngolaryngoscope. Put the brush through the channel several times.

⚠ Please, insert the brush in direction shown on fig. 7.

Don't try to put the brush in opposite direction until head of the brush comes out of the distal end of the nasopharyngolaryngoscope.


Fig. 7

Close the hole of the aspiration channel with the finger, dip the distal end into the water for aspiration for 10 seconds, and aspirate water and air through the channel in turn for its washing.

5. Connect irrigator to the exhaust valve (as it is shown on fig. 8), wash all the channels with clean water or disinfection solution.


Fig. 8

6. Plunge nasopharyngolaryngoscope into cleaning solution, wipe all external surfaces and wash with water.
7. Wipe all external surfaces of the nasopharyngolaryngoscope.
8. Install valve of the instrumental channel.
9. Wipe the working part of the nasopharyngolaryngoscope with cloth moisture with alcohol.
10. Leakage test.

Insert ETO cap into the outlet for leakage test. Connect tester for leakage test to the nasopharyngolaryngoscope. Pump it till indicator arrow is in the test zone, pressure 18 kPa (don't exceed this level), plunge the nasopharyngolaryngoscope into the water (don't plunge manometer). Turn control handle three times passing the distal end up and down. Watch the nasopharyngolaryngoscope within 3 minutes. If air bubbles don't appear, nasopharyngolaryngoscope is leak proof. (Notice: at the beginning of immersion of the nasopharyngolaryngoscope into the water bubbles can be observed in the air, but not in the water, it's quite natural). Watch for 30 seconds. If there are no bubbles, nasopharyngolaryngoscope is leak proof.

⚠ If indicator arrow moves, it means that hermetic state is disturbed, take the nasopharyngolaryngoscope out of the water and contact HUGER representative to conduct repair works.

⚠ Leakage test must be performed after each use of the nasopharyngolaryngoscope. Don't use the nasopharyngolaryngoscope if hermetic state is disturbed. Contact the representative of HUGER for technical support. Otherwise serious breakdown of the nasopharyngolaryngoscope is possible.

11. In the process of cleaning of the nasopharyngolaryngoscope it can be entirely plunged into the water. Watch for 30 seconds. If there are no bubbles,

the nasopharyngolaryngoscope is leak proof; it can be totally plunged into the water for washing and disinfection.

7.3 Disinfection of the nasopharyngolaryngoscope

For high level disinfection we recommend to use "Sidex OPA":

Disinfection of the nasopharyngolaryngoscope must be done under following conditions:

Pressure: up to 1.5 hPa Temperature: up to 40°C

⚠ Disinfection under below mentioned conditions can cause serious breakdown of the nasopharyngolaryngoscope, which must be forbidden.

- Pressure exceeding 1.5 hPa and temperature more than 40°C.
- Washing and disinfection by ultrasound.
- Disinfection by boiling in water.
- Drying in dry-heat sterilizer for disinfection.
- Disinfection by steam.
- Disinfection with liquid cresol.
- Use of benclor, which was not diluted for disinfection.

⚠ in case of application of other disinfection methods that are not recommended, please, get in touch with HUGER for references.

1. Plunge nasopharyngolaryngoscope into disinfectant.
2. Time for soaking must be no less than 20 minutes.
3. After disinfection plunge nasopharyngolaryngoscope into the water. Wash it thoroughly and wipe all the surfaces with aseptic cloth without nap.
4. Wash all the channels with water and dry with the help of syringe and irrigator of channels.
5. Dry all the channels.

7.4 Washing with adapter for processing

1. Insert ends of the adapter for processing into instrumental channel, into the outlet for leakage test and into the tube of exhaust valve.
2. Connect adapter for processing to the washing machine (not delivered with the nasopharyngolaryngoscope).
3. Wash the nasopharyngolaryngoscope in the washing machine.

7.5 Sterilization of the nasopharyngolaryngoscope

Before sterilization the nasopharyngolaryngoscope must be thoroughly cleaned and dried, its ETO cap should be placed too.

Sterilization must be done with the help of ethylene oxide with the following parameters: temperature: $(55 \pm 5)^\circ \text{C}$, pressure: 1–1.7 kg/cm², humidity: $(50 \pm 10)\%$, gas concentration: (600 ± 20) mg/l, duration: 4 hours.

Results of sterilization depend on different factors, such as packing of the nasopharyngolaryngoscope or method of placement into sterilizer. Please, control the efficiency of sterilization with the help of biological and chemical indicators.

7.6 Storage and transportation of the nasopharyngolaryngoscope

(1) Before storage of the nasopharyngolaryngoscope it must be dried thoroughly. Give consideration to drying of the distal end, lens and light guide. Apply silicone resin to the edge of clean cloth and wipe delicately the surface of lens. It will prevent from film formation on the lens from remaining drops of water.

 Don't wipe off silicon resin.

(2) Storage area must be dry, clean, well ventilated with constant normal temperature. Avoid direct sunlight, high temperature and humidity.

(3) Nasopharyngolaryngoscope should be kept with straightened working part. If it is necessary to fold up the working part, you can do it to the same extent as when it is transported in the case.

(4) Don't use case for transportation for storage of the nasopharyngolaryngoscope. The case for transportation is designed for transportation of the nasopharyngolaryngoscope only.

(5) Transportation of the nasopharyngolaryngoscope must be done only in original case for transportation.

7.7 Precautions for storage and transportation

Give consideration to permissible limits of storage and transportation conditions:

1. Storage (Temperature / Humidity / Pressure)

10°C ~ 35°C / 30% ~ 85% / 700 hPa ~ 1060 hPa

2. Transportation (Temperature / Humidity / Pressure)

-50°C ~ +50°C / 10% ~ 95% / 500 hPa ~ 1060 hPa

7.8 Repair of the nasopharyngolaryngoscope

(1) In case of necessity of repair of the nasopharyngolaryngoscope, it should be delivered to the representative of HUGER in original package together with description of the breakdown or damage of the nasopharyngolaryngoscope.

(2) Minor problems in functioning of the nasopharyngolaryngoscope can be eliminated by the user. See section 8 "Troubleshooting". All the rest of repair works must be done by the representative of HUGER. Manufacturer bares no responsibility for injury of the patient/user, damage of the nasopharyngolaryngoscope, which was caused by the repair performed by unauthorized persons.

WARNING

For the purposes of infection control and security of persons, who will contact the nasopharyngolaryngoscope, before delivery to HUGER, it must be properly cleaned and disinfected at the high level. If nasopharyngolaryngoscope was used in patients with contagious diseases (hepatitis, HIV positive patients), please inform about Service Department of HUGER.

7.9 Battery

Power supply of the portable lantern is done from lithium-ion battery.

Power supply form the battery can be done in certain period of time. If you want to decrease the level of lightning from the portable lantern, take out the battery from the lantern and put in into charging device, observing polarity according to marking (time of full battery charging is 4 hours). The process of charging is inducted by yellow light-emitting diode, upon completion of charging light-emitting diode goes out. If you need to continue working, it is necessary to insert new battery.

 It is allowed to use lithium-ion battery supplied by the manufacturer.

8. Troubleshooting

Read carefully enumerated situations. If the problem can't be solved stop using nasopharyngolaryngoscope and deliver it to HUGER.

Problem		Possible reason	Solution
Image	Blurred image.	Objective lens doesn't correspond to the sight of the researcher.	Turn the wheel of eye glass focusing to focus the image.
		Objective lens is dirty.	Plunge the tip into the water to remove mucus.
		Color palette is disturbed.	Stop working; deliver nasopharyngolaryngoscope to HUGER for repair.
Exhaust	No aspiration or weak suction.	Exhaust valve is blocked.	Remove the valve. Wipe port with the cloth.
		Suction channel is blocked.	Disconnect suction pipe. Set exhaust valve loose and make higher pressure to remove remnants from the tube.
		Valve of the instrumental valve is damaged	Replace the valve of the instrumental channel.
	Exhaust valve stuck.	Exhaust valve is dirty.	Remove the valve and wash it. Wipe port with the cloth moisture with alcohol. Put the valve back.
	Leakage of water or air from the valve of instrumental channel.	Valve is worn out or absent.	Install new valve.
Distal end bend	Irregularity of the distal end.	Volume of movements is less than indicated in specifications.	Deliver nasopharyngolaryngoscope to the representative of for repair. (if resistance occurs in the process of movement of control handles, don't apply force, otherwise it may result in premature wear and damage of the nasopharyngolaryngoscope).
		Failures occur while bending.	
		Strong resistance.	
		Cover of the bend section is damaged.	Stop using the equipment. Deliver

9 MANUFACTURER WARRANTY

The manufacturer guarantees the compliance of the nasopharyngolaryngoscope with all the requirements of the present regulatory document while observing working, transportation and storage conditions by the customer.

Guarantee period of the nasopharyngolaryngoscope is 12 months from the date of sale.

In the matter of current or warranty repair contact the representative of Huger: Limited Liability Company "HEALTHY WORLD" 190000, Saint Petersburg, Pochtamtskaya ulitsa, 10, lit. A, pom. 1-H, phone: 7 (812) 385-74-55

10. DISPOSAL CONSIDERATIONS

Nasopharyngolaryngoscopes should be prepared and delivered to special licensed institutions dealing with disposal of electronic devices. Disposal should be performed in accordance with regulations effective on the territory of the Russian Federation.

Annex 1

1) Symbols

Following symbols are used in annex II

 Indicates dangerous condition disregard of which can cause death or serious injury.

 Indicates dangerous condition disregard of which can cause injury of medium or slight severity level, as well as the condition, disregard of which can cause breakdown of the equipment.

Note: additional information.

2) Purpose use

Exhaust valve is inserted into aspiration port of the nasopharyngolaryngoscope. Don't use the valve with any other purpose.

3) Operation procedure

-  The valve was not disinfected and sterilized before shipment. ETO gas sterilization or with high pressure/temperature is required after washing. The process of washing, disinfection and sterilization has been described above. Each valve is designed for 6 usages, if you determined the following troubles, immediately replace the valve:
 - The valve is faintly fixed in the port.
 - The solution injected with the syringe, is aspirated in the suction cup without using the valve.
 - The valve disconnects in the process of functioning.
 - Air aspirates from the end of nasopharyngolaryngoscope without using the valve.
- Don't use the valve with cracked or torn body, as well as stop using the valve if it can be easily withdrawn from the port, otherwise aspiration will be constant.
- If in the process of working it is impossible to stop aspiration, disconnect the tube of the valve from aspirator.

3.1) Connection of the aspiration valve to the nasopharyngolaryngoscope

⚠ Insert aspiration valve tight into the aspiration port of the nasopharyngolaryngoscope. If the valve was inserted in the wrong way, so that you can see a gap between the base of the valve and upper part of the aspiration port, the valve can be disconnected from the nasopharyngolaryngoscope and/or secretion of the patient can leak out or splash up from the port, seriously infringing infection control.

- Insert the exhaust valve into the nasopharyngolaryngoscope port, give consideration to the direction of the valve body. Press the valve with your finger to fix it with the click. (see Figure 1).


Figure 1. Connect the exhaust valve to the nasopharyngolaryngoscope

Check out if there is a good contact of the valve base with the port of the nasopharyngolaryngoscope. Incorrect installation implies gap between the base of the valve and upper part of the port (see Figure 2).

Figure 2. Correct installation

Incorrect installation

Note: Sometimes valve closes with a snap before complete installation. Press the valve to insert in completely into the cylinder.

3.2) Check-out before using

- Connect the tube of the exhaust valve to the aspirator. Turn on the aspirator.
- Plunge the distal end of the working part into aseptic water and press the button of the exhaust valve. Make sure that the water is constantly supplied into the aspirator bottle. Set the valve loose and make sure there is no aspiration. After check-out of aspiration take out the distal end of the water and aspirate air within 2 or 3 seconds to get the water out of the channel.

3.3) Operation procedure

Press and hold the button of the exhaust valve for activation of aspiration.

3.4) Removal of the exhaust valve from the nasopharyngolaryngoscope (Figure 3)


Figure 3. Removal of the exhaust valve from the nasopharyngolaryngoscope

- Revolve the valve tube clockwise.
- Remove the valve from the nasopharyngolaryngoscope.

4) Washing, disinfection and sterilization

- Always put on personal protective equipment from dangerous chemical and potentially infected materials. Personal protective equipment include glasses, mask, moisture protective suit, chemically stable gloves of suitable size, long enough to close all the skin.
- The results of sterilization depend on different factors, such as package of the instruments, method of placement and loading of the instrument into sterilizer. Please, control efficiency of sterilization with the help of biological and chemical indicators.

Note: Use cleaning solution in accordance with recommendations of the manufacturer of this solution.

4.1) Dismantling of the exhaust valve (Figure 4)

- Disconnect piston with the spring from the frame.
- Take spring out if the piston.
- Remove the valve from the nasopharyngolaryngoscope frame.

4.2) Cleaning (Figure 5)

- Insert the brush for cleaning of the channel outlet and rub inside several times.
- Insert the brush for cleaning of the channel into the plug of the tube and rub inside several times.
- Rub external surfaces of the frame and piston with the toothbrush in washing solution.

- Clean inner parts of the spring and the valve with the help of brush for cleaning of the channel outlet. Rub external surfaces with the toothbrush.
- Put all parts into ultrasound washing machine at 40 kHz and clean them with ultrasound within 5 minutes.
- Wash all the parts under running water.

Note: if some blood or mucus is left inside, it can become hard while soaking in disinfectant solution or while autoclaving.


Figure 4. Dismantling of the valve


Figure 5. Cleaning

4.3 Disinfection

- Prepare disinfectant solution according to recommendations of the solution manufacturer.
- If the valve is dismantled, plunge four parts into the disinfectant solution. If there are any air bubbles on the external and/or internal surfaces, wipe them with the cloth without nap or with 30 ml syringe.
- Clean the parts at the temperature and within the period recommended by the manufacturer of disinfectant solution.
- Take the parts out of the disinfectant solution and cautiously wash in aseptic water.
- Take the parts out of aseptic water and dry them thoroughly with aseptic cloth without nap.

4.4) Sterilization

When cleaning is done, apply ETO gas sterilization according to the parameters mentioned above in this manual.

5) Assembling and check-out of the aspiration valve

5.1) Assembling of aspiration valve

- Insert the valve into the frame, as it's shown on figure 6.
- Install the spring on the top of the frame.
- Insert the piston into the elements according to figure 6 and press several times.

5.2) Check-out after assembling (Figure 6)


Figure 6. Assembling of the exhaust valve

- Make sure that the spring is correctly installed into the frame.
- Check out of connection of the piston to the spring; make sure that the gap between the spring and the piston is small.
- Make sure that the valve is fixed in the frame in the right way.
- Examine visually for cracks on the spring and the valve.

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