**FRESENIUS
MEDICAL CARE**

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RUSSIAN FEDERATION

Instruction for use

proVencare
Double Lumen Catheter Sets

For registration purposes in the Russian Federation:

Bad Homburg, 20.12.2016

Place and date

Bettina Kriener / Product Center Manager
ESM & GxP Systems
Fresenius Medical Care AG & Co. KGaA
61346 Bad Homburg - Germany

**FRESENIUS
MEDICAL CARE**

Fresenius Medical Care
AG & Co. KGaA
61346 Bad Homburg v.d.H
Germany

Esther Nèu / Registration Manager
ESM & GxP Systems
Fresenius Medical Care AG & Co. KGaA
61346 Bad Homburg - Germany

Instructions for use

proVencare

Double Lumen Catheter Sets

GENERAL NOTES

Refer to product or carton label for:

	Manufacturing date		Catalogue number
	Use before		Batch code
	Ethylene oxide sterilization		Caution! Please refer to application instructions
	Catheter length		Units
	Outer Diameter		Latex free
	Single Use Only		Do not use if package is damaged

1. Definition of expressions used in this manual

WARNING! is used to alert the operator **NOT** to take a specific action that may cause a hazard which, if ignored, could result in an adverse reaction, injury or death.

CAUTION! is used to alert the operator to **TAKE** a specific action to protect against a possible hazard which, if ignored, could have an adverse reaction on the patient or equipment.

NOTE! is used to remind the operator of normal treatment functions and what is the suitable action in a particular situation.

2. General warnings and cautions

The medical techniques and procedures described in these instructions do not represent all medically acceptable protocols, nor are they intended as a substitute for the physician's experience and judgment in treating any specific patient. The following procedures are offered to clinicians only as a guide. Apply a catheter maintenance protocol and catheter insertion site maintenance protocol consistent with accepted standards of practice.

Read all instructions about the proVencare catheters carefully before using them.

To ensure safe and effective treatment, the following information must be observed.

Please ensure that you are completely familiar with the complications listed below prior to inserting a proVencare catheter:

air embolism, arterial puncture (e.g. A. femoralis, A. carotis, A. subclavia), atrial or ventricular perforation, bacteraemia, brachial plexus injury, cardiac arrhythmia, cardiac tamponade secondary to vessel wall injury, catheter embolism, central venous thrombosis, dysrhythmias, endocarditis, exit-site infection, exsanguination, femoral nerve injury, haematoma formation, haemorrhage, haemothorax, laceration of the vessel, luminal thrombosis, nerve damage, perforation of the vessel, pleural or mediastinal injury, pneumothorax, retroperitoneal bleed, right atrial puncture, septicaemia, subcutaneous haematoma, subcutaneous tunnel infection, superior vena cava puncture, thoracic duct laceration, thrombosis, trauma to major vessel or right atrium, vascular thrombosis.

In order to insert the catheter in the subclavian, jugular (internal or external) or femoral vein, standard operating-room protocol(s) should be used.

CAUTION!

1. For single use only. **DO NOT** reuse, reprocess or re-sterilise. Reuse, reprocessing or re-sterilisation may create a risk of contamination of the device and cause patient infection or cross infection, including, but not limited to the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

2. Contents sterile and non-pyrogenic in unopened, undamaged package.

3. **DO NOT** store with chemicals. Failure of the catheter can result.

4. **DO NOT** use catheter or accessories when the expiry date indicated on the labelling has been exceeded. The manufacturer will not be liable for any damage caused by use of expired products.

5. The devices should only be used by authorized personnel who are totally familiar with the procedure.

6. Aseptic technique, including the use of a face mask, hand washing and the wearing of sterile gloves must be used throughout all procedures in order to reduce the risk of disease transmission.

7. Before every use the correct position of the catheter must be checked.

8. After use, this product may be a potential hazard. Handle and dispose of it in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

9. The operating room is the preferred location for catheter placement. Percutaneous procedures require confirmation of correct placement by fluoroscopy or chest X-ray.

10. Avoid air embolism by keeping the catheter clamped at all times when not connected to syringes or blood lines.

11. Use only protection caps and other accessories provided by the manufacturer.

12. Luer connectors of proVencare catheters comply with international standard ISO 594. Safe and tight connection, however, can only be assured when the matching bloodline connector is also in compliance with this standard. Use of bloodline connectors deviating from the above standard may cause problems during connection such as stress cracks and leaks which are associated with risks of air embolism and blood loss.

13. Protection caps and tubing should be connected without applying excessive force.

14. An over-tightened catheter connection can crack luer connectors.

15. Place clamps only in the centre of the extension lines. Repeated clamping particularly near or on the hard plastic portion by the connectors may cause tubing fatigue and possible disconnections.

16. With each change in tubing connection, purge air from the catheter and all connecting tubing.

17. Before initiating treatment, aspirate indwelling citrate or heparin based locking solution and discard.

18. The catheter must be re-filled with indwelling citrate or heparin based locking solution after each use in order to maintain patency.

19. **DO NOT** infuse against a closed clamp or forcibly infuse a blocked catheter: back pressure will cause catheter failure.

20. Tape all non-locking connections to prevent accidental separation.

21. The catheter must always be secured against dislocation by additional use of medical tape or plasters.

22. Tape the caps between treatments to safeguard them against accidental removal.

23. Wound dressing must be kept dry. If adhesion of dressing is compromised by profuse perspiration or accidental wetting, the dressing must be changed by the medical or nursing staff under aseptic conditions.

24. A clotted catheter should be replaced if the clot cannot be aspirated or dissolved (see also Troubleshooting: Thrombus formation).

25. Patients requiring ventilator support are at increased risk of pneumothorax during subclavian vein cannulation.

26. Remove the catheter as soon as it is no longer required.

Instructions for use proVencare Double Lumen Catheter Sets

WARNING!

DO NOT expose extension tubes to sharp edges or objects such as metal clamps and **DO NOT** use scissors to remove dressings as extension tubes can be damaged and fatal blood loss or air embolism can occur.

DO NOT use alcohol, drugs containing alcohol, acetone based products or ointments. As an exception to the general exclusion of alcohol-containing disinfectants, 2% chlorhexidine gluconate solution in 70% isopropyl alcohol may be applied on the catheter or at the catheter entry site by wiping with pre-saturated disinfection wipes but limited to a total of 30 repetitions during the use of the catheter. Other forms of or more often application must not be used as failure of the catheter will result. Iodine based solutions without alcohol are the recommended antiseptic solutions that can be used for care of this catheter.

DO NOT re-sterilise the catheter by any method. The manufacturer will not be liable for any damage caused by reuse of this product.

DO NOT damage the catheter when suturing.

3. Specifications

Product description and materials:

The catheter is provided within a set. The Set includes additional equipment for insertion of the catheter using the Seldinger technique. Basic technical parameters of the set as well as used materials are represented in the following tables (see last page).

Priming volumes:

See preprint at the guiding insert of the clamps on every catheter.

Flow path:

In single-lumen catheters the blood flow path is indicated by a white clamp.

In multi-lumen catheters the arterial side is indicated by the red clamp and the venous side is indicated by the blue clamp.

In triple-lumen catheters the third, central lumen with a slide clamp can be used for guidewire placement, blood sampling and central venous pressure monitoring.

4. Indications

proVencare catheters are designed to be used for short-term vascular access for haemodialysis, haemofiltration and other extracorporeal blood cleaning processes. They can be used in the subclavian, jugular or femoral veins.

5. Contraindications

These catheters are only to be used as indicated. The catheter is not to be inserted in veins demonstrating thrombosis.

6. Adverse reactions

Special attention must be paid to patients with known hypersensitivity reactions to products sterilised with ethylene oxide, to components itemised in the material list (see last page) or other allergic reactions with other substances. The consultation of a physician is indispensable to judge the risk and take appropriate action if there is a suspicion of a hypersensitivity reaction

7. Directions for use

7.1 Guidelines for insertion of catheters

NOTE! The physician must use his or her own judgment, experience and discretion in determining which insertion technique should be used for each individual patient. The selection of the appropriate catheter length is at the sole discretion of the physician. To achieve proper tip placement, proper catheter length selection is important.

CAUTION! The location of the vessels should be detected by appropriate technique each time before insertion. Routine X-ray should always follow the initial insertion of this catheter to confirm proper placement prior to use.

7.2 Insertion site

WARNING! DO NOT place the catheter into or allow it to remain in the right atrium or right ventricle.

Tip Placement

(internal jugular / subclavian position)

7.3 Preparation

1. Provide a sterile operative field.
2. Using standard aseptic technique, prepare the insertion site and drape the prepared area.
3. Flush all lumens of the catheter with sterile saline and clamp the extension lines, leave the syringes at the tubing extensions.

WARNING! In High-Flow catheters: The venous clamp must not be closed in order to prevent damage to the inner dilator.

WARNING! To prevent air embolism, keep the catheter clamped at all times when not attached to a syringe or bloodlines.

4. For subclavian and jugular insertion: Place the patient in the Trendelenburg position, turning the head to the side opposite the insertion site.

For femoral insertion: Place the patient flat on their back.

Flex the patient's knee on the same side as the insertion site, while abducting and rotating the leg outwards.

5. A local anaesthetic is administered to the insertion site.

Instructions for use

proVencare

Double Lumen Catheter Sets

6. Connect a syringe to the Seldinger needle and insert it into the target vein. Aspirate a small amount of blood to ensure the needle is correctly positioned in the vein.

CAUTION! If arterial blood is aspirated, remove the needle and apply immediate pressure to the site. Ensure that the bleeding has stopped and that no haematomas have developed before attempting to cannulate the vein again.

7. Remove the syringe and hold the Seldinger needle closed with the thumb to prevent blood loss and air embolism.

8. Insert slowly the guidewire through the needle and into the vessel.

9. Advance the guidewire to the appropriate location in the vessel.

CAUTION!

- The guidewire should thread easily into the vein without resistance.
- If during insertion of the guidewire a resistance is noticed or further movement is not possible, the guidewire must be removed together with the combined product. If withdrawal cannot be easily accomplished, a chest X-ray should be obtained and further consultation requested.
- To avoid injury, do not insert the guidewire or dilators forcibly.
- Also, cardiac arrhythmia may result if the guidewire is allowed to pass into the right atrium. It is therefore recommended to monitor the procedure either by ECG or sonography.

WARNING! DO NOT pull back the guidewire through the Seldinger needle as this may sever the tip of the wire. If the guidewire must be removed, withdraw needle and guidewire.

10. While holding the guidewire in place, remove the Seldinger needle without accidentally altering the guidewire position.

11. Secure the guidewire in place.

12. Create a small incision at the insertion site.

13. Place the dilator on the distal end of the guidewire. Following this, dilate the subcutaneous tissue and vein wall, so that the catheter can easily be inserted into the target vein.

14. Withdraw the dilator.

15. Remove the syringe from the extension which is foreseen for the guidewire and set the catheter on the distal end of the guidewire.

For double-lumen catheters the venous lumen (blue) is used.

For triple-lumen catheters the middle lumen (white) is used.

16. By gentle rotation advance the catheter carefully through the subcutaneous tissue into the target vein. Insert the catheter to the appropriate depth. At the same time ensure that the guidewire is free by withdrawing it in small steps. For best safety and patient comfort the catheter should be inserted all the way down to the guide mark.

17. Withdraw the guidewire and close the clamp. For catheters with inner dilator / stylette, the inner dilator / stylette should be removed together with the guidewire.

WARNING! To prevent air embolism, expel all air from syringes before injecting any solutions.

18. Test the flow: Attach a 20 ml syringe, unclamp the catheter extension line and aspirate blood and any air present. The syringe should fill up to 20 ml within 6 seconds, which corresponds to a blood flow of 200 ml/min. This test should be performed on all lumens.

19. When satisfactory flow is achieved, flush each lumen with at least 20 ml of saline and close the clamp.

CAUTION! Clamp each lumen under positive pressure. Failure to clamp lumens can lead to air embolism or extensive blood loss.

20. Put caps provided by the manufacturer on the connectors of the extension lines.

21. Sew the catheter to the skin using the suture wings. Do not sew the catheter tubes.

22. Cover the entry wound with a tight dressing.

CAUTION! Correct catheter placement should be verified through the use of fluoroscopy or X-ray.

23. Enter the catheter length and catheter lot number in the patient's record.

24. If the catheter is not to be used immediately the catheter lumen must be locked (see chapter 7.5.)

7.4 Guidelines for catheter and exit-site care

WARNING! DO NOT use alcohol-based solutions to disinfect, wash, swab or soak the catheter. As exception to the general exclusion of alcohol-containing disinfectants, 2% chlorhexidine gluconate solution in 70% isopropyl alcohol may be applied on the catheter or at the catheter entry site by wiping with pre-saturated disinfection wipes but limited to a total of 30 repetitions during the use of the catheter. Other forms of or more often application must not be used as this can damage it. It is recommended to use a povidone iodine-based or peroxide solution.

- Keep the catheter and the exit-site dry.
- Never leave the catheter open. Always seal with syringes or caps.
- Ensure that the suture is correctly placed. The catheter should not be able to move in and out of the exit-site.
- Exchange the caps after each treatment.
- Cover the extension lines between treatments.

7.5 Guidelines for application of catheter lock solutions

To reduce the risk of a catheter occlusion between treatments keep the lumens of the catheter filled with a sufficient concentration and volume of a citrate or heparin based locking solution. In most cases, a concentration of 5,000 units/ml heparin is appropriate. Refer to the priming volume imprinted on the clamp-inserts of the catheter. The valid local regulations regarding locking solutions are to be followed. If applicable refer to the instructions for use of the lock solution.

1. After extracorporeal treatment, flush each lumen with 20 ml sterile saline and clamp the extensions.

2. Fill syringes according to the priming volumes of each lumen with the locking solution.

3. Remove the caps (or syringes) from the extensions and fit the prepared syringe to the corresponding Luer-Lock connector.

Instructions for use proVencare Double Lumen Catheter Sets

4. Open the clamps on the extensions. Instill the locking solution by flushing quickly to ensure that the locking solution reaches the end of the lumens. Instilling too slowly may cause the locking solution to exit the catheter from the proximal inlet holes, leaving the distal hole unprotected against thrombus formation.
5. Clamp the extension lines under positive pressure.
6. Remove the syringes and fit new sterile caps onto the Luer-Lock connectors of the extensions.

WARNING! To further prevent thrombus formation at the tip of the catheter, the clamp on the extension tube should not be opened once the lumen has been primed. If the clamp was opened, the priming volume of that lumen would increase slightly as a result of the extension tube returning to its "normal", unclamped state. This would create a vacuum at the tip, causing blood to be sucked up into the distal portion of the catheter, ultimately resulting in thrombus formation.

7.6 Guidelines for catheter removal

CAUTION! Only a qualified, licensed physician should remove the catheter.

1. **For subclavian and jugular positions:** Place the patient in the Trendelenburg position, turning the head to the side opposite the insertion site.

For femoral position: Place the patient flat on the back. Flex the patient's knee on the same side as the insertion site, while abducting and rotating the leg outwards.

2. Remove dressing and suture of the catheter.
3. Compress the venous insertion site and pull the catheter gently out of the vessel.
4. Apply appropriate wound dressing.

NOTE! Because of the risk of bleeding, the patient should, if possible, sit up rather than lie down after the catheter removal (not applicable with femoral catheters).

7.7 Troubleshooting

One-way obstruction

Suspect one-way obstruction when a lumen can be flushed but cannot be aspirated. This is usually caused by tip malposition with the catheter sucked against the vessel wall.

One of the following adjustments may eliminate the obstruction (sutures may have to be removed):

- Reposition the patient.
 - Have the patient hold his arm above his head and cough.
- Not applicable with femoral catheters.
- Flush with saline to push the catheter away from the wall.
 - Rotate the catheter. Not applicable if the catheter is sutured to skin or in a subcutaneous tunnel.
 - Withdraw the catheter slightly to reposition the tip. Not applicable if the catheter is sutured to skin or in a subcutaneous tunnel.

WARNING! DO NOT insert the catheter further into the vein.

Kinking

- Check catheter position.
- Check sutures at the exit-site.
- Does the catheter move in and out of the exit-site?
- Does it bend when the patient moves or lies on the same side as the catheter site?

Symptoms of kinking at the tip include decreased blood flow and increased venous resistance. An X-ray can confirm internal kinking.

Vessel collapse or spasm

Start and adjust pump speed slowly.

Thrombus formation

1. If catheter is blocked, first check to ensure it is not kinked.
2. If a thrombus is present, attempt to aspirate the clot gently with a syringe. If blood flow is still low, or blood cannot be withdrawn at all, the physician may choose to dissolve the clot with a chemical intervention using an approved thrombolytic agent.

WARNING! DO NOT use excessive force when flushing an obstructed lumen. Over pressurising the catheter will damage it and can cause trauma to the patient.

NOTE!

- A catheter could be blocked by a substance other than a blood clot.
- The cleared catheter may be used immediately.

Infection

Any signs of infection have to be considered as an indication for catheter removal.

8. Warranty Conditions

The manufacturer warrants that the catheter has been manufactured in accordance with its specifications and in compliance with the Medical Device Directive 93/42/EEC. Products with a manufacturing defect will be replaced if the defect is reported with the lot number.

The manufacturer shall not be liable for any misuse, improper handling, non-compliance with instructions for use or cautionary notes and for any damage incurred (except damage of the physical integrity) subsequent to the manufacturer's delivery of the medical device.

9. Maintenance

For single use only.

Do not use if the product or package is damaged.

After use, this product may be a potential hazard.

Handle and dispose of it in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

Maintenance is not provided!

10. Transportation and storage

Transportation is conducted by all types of enclosed transport (in rail carriages, containers, heated pressurized compartments, by plane, in closed cars, underdeck, etc.) according to the rules of carriage effective for the transport of the relevant type.

Transportation and storage temperature: +5 to +30°C.

Shelf-life is 5 years.

11. Authorized representative of the manufacturer in the Russian Federation:

Closed Joint-Stock Company (ZAO) "Fresenius SP", room 09, building 6, 2, Ugreshskaya str., 115088 Moscow.
Tel/Fax: 8(495)789-64-54

Basic technical parameters

FMC Article №	FMC Product №	Catheter mm	Flow venous lumen (ml/min)	Flow arterial lumen (ml/min)	Extension Line	Guide Wire (mm)	Dilator (mm)	Seldinger needle (mm)	LL-Protection cap (mm)	Wound dressing
FDC-607,5P	F00002309	2,2 x 75	136	153	Straigh	0,71 x 500	2,3 x 150	1,25 / 70	2 pieces	1 piece
FDC-610P	F00002310	2,2 x 100	131	141	Straigh	0,71 x 500	2,3 x 150	1,25 / 70	2 pieces	1 piece
FDC-612,5P	F00002311	2,2 x 125	126	128	Straigh	0,71 x 500	2,3 x 150	1,25 / 70	2 pieces	1 piece
FDC-810P	F00002312	2,7 x 100	182	199	Straigh	0,89 x 700	0,9 x 150	1,25 / 70	2 pieces	1 piece
FDC-812,5P	F00002313	2,7 x 120	178	191	Straigh	0,89 x 700	0,9 x 150	1,25 / 70	2 pieces	1 piece
FDC-815P	F00002314	2,7 x 150	173	183	Straigh	0,89 x 700	0,9 x 150	1,25 / 70	2 pieces	1 piece
FDC-1115	F00002315	3,7 x 150	323	352	Straigh	0,97 x 700	3,7 x 150	1,5 / 70	2 pieces	1 piece
FDC-1120	F00002317	3,7 x 200	309	322	Straigh	0,97 x 700	3,7 x 150	1,5 / 70	2 pieces	1 piece
FDC-1125	F00002318	3,7 x 250	296	292	Straigh	0,97 x 700	3,7 x 150	1,5 / 70	2 pieces	1 piece
FDC-1115J	F00002322	3,7 x 150	323	352	Curved	0,97 x 700	3,7 x 150	1,5 / 70	2 pieces	1 piece
FDC-1120J	F00002324	3,7 x 200	309	312	Curved	0,97 x 700	3,7 x 150	1,5 / 70	2 pieces	1 piece
FDC-1215	F00002319	4 x 150	394	406	Straigh	0,97 x 700	4 x 150	1,5 / 70	2 pieces	1 piece
FDC-1220	F00002320	4 x 200	384	388	Straigh	0,97 x 700	4 x 150	1,5 / 70	2 pieces	1 piece
FHC-1115	F00002330	3,8 x 150	324	329	Straigh	0,81 x 700	4 x 150	1,5 / 70	2 pieces	1 piece
FHC-1120	F00002331	3,8 x 200	315	315	Straigh	0,81 x 700	4 x 150	1,5 / 70	2 pieces	1 piece
FHC-1125	F00002332	3,8 x 250	297	302	Straigh	0,81 x 800	4 x 150	1,5 / 70	2 pieces	1 piece
FHC-1115J	F00002333	3,8 x 150	324	329	Curved	0,81 x 700	4 x 150	1,5 / 70	2 pieces	1 piece
FHC-1120J	F00002334	3,8 x 200	315	315	Curved	0,81 x 700	4 x 150	1,5 / 70	2 pieces	1 piece
FHC-1315	F00002335	4,3 x 150	366	382	Straigh	0,97 x 700	4,3 x 150	1,5 / 70	2 pieces	1 piece
FHC-1317,5	F00002336	4,3 x 175	357	374	Straigh	0,97 x 700	4,3 x 150	1,5 / 70	2 pieces	1 piece
FHC-1320	F00002337	4,3 x 200	349	363	Straigh	0,97 x 700	4,3 x 150	1,5 / 70	2 pieces	1 piece
FHC-1325	F00002338	4,3 x 250	316	352	Straigh	0,97 x 800	4,3 x 150	1,5 / 70	2 pieces	1 piece
FHC-1315J	F00002339	4,3 x 150	366	382	Curved	0,97 x 700	4,3 x 150	1,5 / 70	2 pieces	1 piece
FHC-1320J	F00002340	4,3 x 200	349	363	Curved	0,97 x 700	4,3 x 150	1,5 / 70	2 pieces	1 piece

Material list

Product	Description	Material	Abbreviation
Catheter	Catheter tube	Polyurethane	PUR
	Extension line	Polyurethane	PUR
	Connector	Polyvinyl chloride	PVC
Protection Cap	Protection Cap	Polyethylene	PE
Dilator	Dilator	Polyethylene	PE
Inner Dilator	Inner dilator tube	Tetrafluoroethylene / Hexafluoroisopropylene-copolymer	FEP
	Inner dilator hub	Polypropylene	PP
Guidewire	Guidewire	Stainless steel	
	Dispenser	Polypropylene	PP
Seldinger needle	Needle tube	Stainless steel	
	Needle hub	Polycarbonate	PC
Wound dressing	Film	Polyurethane	PUR

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