Jose Oswaldo Caputto Flosi
Jose Oswaldo Caputto Flosi
International Sales Manager
FANEM LTDA

Eng. José O. C. Flosi
INTERNATIONAL SALES MANAGE
FANEM LTDA

USER'S MANUAL

BILITRON® MODEL 3006

THE SUPER LED PHOTOTHERAPY


Technical Standard - NBR IEC 60601.1
NBR IEC 60601 2.50
NBR IEC 60601 1.2

Revision: 07/07
Edition: 05/07
Ref.: 006.410.550

Jose Osvaldo Caputto Flosi
International Sales Manager
FANEM LTDA

Eng. José O. C. Flosi
INTERNATIONAL SALES MANAGER
FANEM LTDA

Technical information

BILITRON® MODEL 3006

THE SUPER LED PHOTOTHERAPY

Technical Standard - NBR IEC 60601.1
NBR IEC 60601.2.50
NBR IEC 60601.1.2

Revision: 07/07
Edition: 08/07
Ref.: 006.410.650

TABLE OF CONTENTS

1 – Introduction.....	01
2 – Technical Characteristics.....	02
3 – Parts, Pieces and Accessories.....	05
4 – Precautions, Restrictions and Warnings.....	07
5 – Assembling the Equipment	10
6 – Operating the Equipment	10
7 – Preventive Maintenance, Corrective and Conservation.....	14
8 – Warranty and Technical Assistance.....	17

USER MANUAL
BILITRON® MODEL 3006
THE SUPERLED PHOTOTHERAPY

1 - Introduction

Phototherapy is commonly used in the treatment of neo-natal hyperbilirubinemia, by exposing the patient to radiation, focused on the blue visible light spectrum, for a period of time to be determined by the attending physician.

Ever since its very beginning in the end of the 50's, and later by the marked evolution of the types of lamps being used, the introduction of phototherapy in the treatment of neo-natal hyperbilirubinemia has minimized risks without, however, eliminating the adverse collateral effects caused by the ultra-violet and infra-red light on the patient's body. The large amount of different type lamps, and more recently, of conventional LEDS which, in order to obtain greater treatment efficiency, require many light sources emanating from the same device, which renders it cumbersome and heavy.

The Model 3006 BILITRON® may be considered a 7th. Generation phototherapy equipment available worldwide in the market. The treatment of hyperbilirubinemia through the use of light started with the use of low-intensity fluorescent lamps. Later on, fluorescent lamps spreading a blue spectrum were used. Still later, halogen lamps were introduced, as well as a variation of optical fibers. The LEDS (light-emitting diodes) were also introduced for phototherapy purposes, but, due to a low light-emission volume, hundreds of LEDS are needed to get satisfactory results.

The coming to the scene of the Super LED, with a high degree of blue spectrum radiation, and bearing a physical-chemical composition different from that of conventional LEDS, allowed the making of equipment made up of few LEDS, which made this equipment smaller and lighter.

The Super LEDS issue a safe type of radiation (according to IEC 825, Class 2 norms), completely eliminating ultra violet and infrared filters (see emissions graph), with the advantage of drastically reducing risks of burning, erythema and the unfelt loss of body fluids. It is then an innovating, micro-processed type of electronic photo-therapy, capable of eliminating undesirable radiations from the ultraviolet and from the infrared spectrums, incorporating an ellipse bearing a larger area with smaller losses at the edges, allowing a more anatomic application on the bodies of the patients.

The microprocessor electronics that has been introduced in this type of phototherapy allows the use of resources never used before, such as radiometers containing optical probes (the probe is optional), time of treatment, hour counter, luminous intensity adjustment (the first phototherapy with this kind of irradiation adjustment), a calendar clock, a distance vs. radiation converter, a data memory bank to issue a report and to communicate with the computer/printer or to communicate with the series 3620 radiometer.

With the intent to preserve its efficiency, and taking into account the nature of the task which the Model 3006 BILITRON® purports to accomplish, some special care is needed for the handling and the operation of the equipment. This User Manual provides general instructions for the installation, the use and the maintenance of the equipment by the operator, together with flaw diagnosing. Fanem® will not accept any liability if the user fails to follow instructions, and if s/he fails to abide by the warning and the maintenance recommendations contained in this manual, or if the user makes any repairs using non authorized components. The calibration and repairs must be executed only by qualified personnel. Any additional information are provided by your local distributor.

This manual should be read and perfectly understood, and it must be made readily available for everybody needing to work with the Model 3006 BILITRON® phototherapy unit. When not in use, the User Manual should be kept together with the BILITRON® unit. If there is any information that you are unable to understand, please get in touch with Fanem®'s representatives to get further details.

USER MANUAL

BILITRON® MODEL 3006

THE SUPERLED PHOTOTHERAPY

2 – Technical Characteristics

2.1 Input Voltage and Power

Model	Voltage (Volts) V~	Frequency Hertz (Hz)	Power VoltAmpere (VA)	Current I (A)	Fanem Ref. Code
3006 BTB	127/220/240	50/60	50	1	006.000.900
3006 BTI	127/220/240	50/60	50	1	006.001.900
3006 BTP	127/220/240	50/60	50	1	006.002.900

Voltage for 240V~ to be consulted.

2.2 Spectral Irradiation Graph



Figure 2. Spectral Irradiation from the Bilitron 3006 power source in the 400nm to 500nm interval.

2.3 Radiometer calibration curve (620 and 2620)

USER MANUAL

BILITRON® MODEL 3006

THE SUPERLED PHOTOTHERAPY

2.4 Dimensions and weight

♦ Dimensions

Model	Width (cm)	Depth (cm)	Height (cm)
3006	11,6	23,0	5,0
Arm for the 3006 BTB	3,5	10,0	30,0
Pedestal for the 3006 BTP	50,0	50,0	100,0 ~140,0
Alternates for the 3006 BTP	3" with a brake (2 units)		
	3" without a brake (1 unit)		

♦ Weight

Model	Net Weight (kg)	Weight When Packaged (kg)
3006 BTB	9,0	14,0
3006 BTI	1,0	3,0
3006 BTP	2,0	4,0

2.5 Symbols Used

	Type B Equipment Class II Isolation		Regular Equipment - Not Water Proof
	Attention: Hot Surface		Attention: Risk of Electric Shock
	Attention: Check Manual		Mandatory Grounding
	OFF (No Power Being Supplied)		ON (Electric Power Being Fed)
	Patient Eye Protector		

2.5.1 - Symbols - Packaging

	Fragile		This Side Up
	Shield Against Sunlight		Shield Against Rain
	Packing Limit		Attention: Check Accompanying Documents

USER MANUAL
BILITRON® MODEL 3006
THE SUPERLED PHOTOTHERAPY

2.6 Specifications

Protection fuse:

1A Type "T" - 3cm to 127V~ (for 3006-BTP models – 2 units)

0,5 Type "T" - 3cm to 220V~ (for 3006-BTP models – 2 units)

Average irradiance at 30 cm:

Luminous focus fulcrum..... 40 $\mu\text{W}/\text{cm}^2 \cdot \text{nm}$ ($\pm 10\%$)

Focus dimensions 250 x 230 mm

Irradiation adjustment band 20% to 100% of maximum irradiation ($\pm 10\%$)

Irradiation adjustment variation at every 10%

Radiometer:

Measurement band 0 a 200 $\mu\text{W}/\text{cm}^2 \cdot \text{nm}$

Resolution 0,1 $\mu\text{W}/\text{cm}^2 \cdot \text{nm}$

Reading precision $\pm 1 \mu\text{W}/\text{cm}^2 \cdot \text{nm}$

Noise < 52 dBA (45dBa environment)

Environmental Conditions:

Operating temperature tier..... 18°C to 28°C (room temperature)

Operating humidity tier..... 10% to 95% (non condensing)

Environmental Conditions for Storage and Transportation:

Room temperature 5°C to 55°C

Relative humidity 30% to 75% (non condensing)

Note.: These conditions apply to equipment that have been properly packaged.

USER MANUAL

BILITRON® MODEL 3006

THE SUPERLED PHOTOTHERAPY

3 - Parts, Pieces and Accessories



The Model 3006 BILITRON® phototherapy unit is made up of strong-impact proof plastic box, containing a module with 5 Superleds, one alpha-numeric display and a key board to program its functions, a micro ventilator for the Superleds module, and a micro processed controller.

The equipment is available in three types of assemblies, and it may be made ready to be installed in incubators, heated cradles or Mobile Spindles.


The Model 3006 BTI BILITRON® is assembled on rubber supports to be used inside of incubators.

The Model 3006 BTI BILITRON® is mounted using an articulate spindle system which supports the headstock, made of steel and painted with epoxy paint, which allows it to move to take advantage better of the light falling on the patient, and is allow the operator to adjust the focus over the patient. It is set directly on its adaptor block inside the Heated Cradles.

The Model 3006 BTP BILITRON® is mounted over a base supported by three 3" turning wheels, to make it easier to move around. It has a column with an adjustable spindle height, to allow it to be better positioned in the apparatus. There is an adjustable and articulate spindle on the BILITRON® body, allowing rotation adjustment to better focus the light over the patient.

AC Adaptor - Ref.: 006.054.900

Optical Kit - Ref.: 620.058.600
(optional)

USER MANUAL

BILITRON® MODEL 3006


THE SUPERLED PHOTOTHERAPY

Presentations:

1 - 3006 BTI
Ref.: 006.001.900

2 - 3006 BTP
Ref.: 006.000.900

3 - 3006 BTB
Ref.: 006.002.900

Item	Code	Description	Obs.
1	006.000.900	Bilitron® phototherapy unit <i>with an arm</i> and with a spindle to be assembled over the Cradle or over the incubator	BTB
2	006.001.900	Bilitron® <i>armless</i> phototherapy unit to be assembled inside the incubator.	BTI
3	006.002.900	Bilitron® phototherapy unit mounted on a Movable Spindle with turning supports	BTP
4	006.054.900	AC Adaptor for the Bilitron®	Included with Bilitron® BTB and BTI models
5	620.058.600	Optical Kit	optional
6	006.132.600	Articulate Arm to be used in Cradles or Incubators	Included with 006.000.900
7	006.133.600	Movable Spindle Set	Included with 006.002.900

3.4 - Power Source Module - Ref.: 006.053.900

A Set containing 5 Super LEDS, mounted over an aluminum dissipater, and bearing 5 directed light diffusers to better take advantage of the light.

The average life cycle of the Power Source Module is of around 20,000 hours. It is recommended, however, that the irradiation be monitored with the FANEM Radiation Monitor, according to normal maintenance procedures.

Attention: Always use the original FANEM Power Source Model, because these modules are controlled and selected individually according to the original project characteristics.

3.5 - Micro Ventilator Kit - Ref.: 006.251.025

Ventilator kit, located on the upper portion of the Box. It gets turned on whenever the equipment is powered up, causing the dissipation of the heat generated by the power source module, and consequently extending its useful life cycle.

3.6 - Optical Kit - Ref.: 620.058.600

Extension probe bearing a head containing an optical sensor, used to measure the irradiance issued by the Model 3006 BILITRON® phototherapy unit.

Attention: The optical kit is optional.

USER MANUAL

BILITRON® MODEL 3006

THE SUPERLED PHOTOTHERAPY

3.7 - Eyes Protector for Neonatal Phototherapy

Single part with design that it involves the head of the new born, confectioned in antiallergic material, with adjustment through Velcro® thus keeping the head of the new born at the original position. It internally contains a cushion projected ergonomically to avoid the eyes exposition of the new born from the entrances of light of the phototherapy.

Order Information:

Code: 006.103.023 – Regular: Occipital-frontal circumference 33-38cm (packing with 20 units)

Code: 006.102.023 – Premie: Occipital-frontal circumference 25-33cm (packing with 20 units)

Code: 006.101.023 – Micro: Occipital-frontal circumference 20-25cm (packing with 20 units.)

*** Accessories of exclusive use with the equipment.**

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Attention: This chapter of the User Manual contains extremely important information to insure the safety and the physical integrity of the patient, of the user and of the equipment. Read it CAREFULLY!

- ♦ It is recommended that the BILITRON® unit is used at a distance of not less than 30 cm between the light source and the patient.

Attention: Minimum distance of 30 cm between the light source and the patient.

Attention: Do not use the Bilitron® unit when the Power Source Module is in direct contact with the body of the patient.

- ♦ Check as to whether the network to which the equipment will be connected has the necessary characteristics to support the electric voltage and power outputs of the equipment, which are shown on a label attached to the device.
- ♦ The power cord must be plugged into a grounded outlet, set permanently in the wall, according to the norms and legislation in effect for low-voltage electric installations, and electric standards set by law for Health Establishments.

Attention: Do not use extensions or multiple plugs. If the ground is not appropriate, do not use the equipment.

Attention: In order to prevent damage to eyesight, the RN will need to use eye shades.

- ♦ The BILITRON® unit must be used according to the ideal comfort norms established for the nursery (between 23°C to 27°C).
- ♦ A phototherapy unit that is improperly used may expose the newly born to serious risks. This unit must be operated exclusively by trained and qualified personnel who are familiar with the risks and the benefits deriving from its use.
- ♦ This unit may not be used if any of its functions is not working properly. Technical servicing will be required in this case.

USER MANUAL

BILITRON® MODEL 3006

THE SUPERLED PHOTOTHERAPY

- ♦ It is unadvisable to leave the patient unattended when under the equipment.
- ♦ The temperature of the patient must be constantly monitored by means of a clinical thermometer.
- ♦ Thermotherapy devices (Heated Cradles, Incubators, Heated Mattresses, Radiation Heaters, etc), when coupled to the BILITRON® unit, may raise the patient's body temperature to dangerous levels. It is recommended that these equipment should be operated in a skin-to-skin controlled mode (self-controlled), or else in a manual mode. The heat rate of the heater must be reduced according to the temperature of the patient's body, and it should be monitored through a clinical thermometer.
- ♦ Always use original lamps. The use of different types of lamps may change the irradiation and temperature conditions regarding the levels that have been specified for the equipment.
- ♦ We recommend the usage of a FANEM Radiation Monitor to make sure that the lamp is emitting the perfect amount of radiation, in order to insure the miniam radiation of $30 \mu\text{W}/\text{cm}^2 \text{ nm}$ at a distance of 300 mm from a light source at the highest point of radiation (center of the luminous focus).
Note: $4 \mu\text{W}/\text{cm}^2 \text{ nm}$ is the value that is accepted internationally regarding adequate treatment (this is the least value according to Mims (1)).

Attention: Do not place the Bilitron® unit next to heat sources when it is used in Heated Cradles. Keep at least 40.0 cm away from the heat source.

Attention: Whenever the Bilitron® phototherapy equipment is used conjoined with a Heated Cradle, make sure that the phototherapy equipment is kept away from the area of the heat being irradiated by the cradle.

- ♦ Do not place any cloths or towels over the body of the Bilitron® unit. This may increase the temperature, and it will cause the irradiating source to be turned off.

Attention: Do not use the phototherapy equipment close to flammable anesthetics, combustible gases or cleaning agents which may provoke combustion.

Attention: The equipment is not adequate to be used with flammable anesthetic mixtures.

- ♦ The normal life cycle of the original Fanem lamps is of around 20,000 hours. However, it is recommended that the irradiation be constantly monitored by using the FANEM Radiation Monitor, pursuant the usual maintenance procedures, in order to get a better measure of the real efficacy of the lamps. Therefore, the irradiation sources (Super Leds) must be replaced whenever they undergo a loss of 25% of its total irradiance in regards to total bilirubin-Ebi.
- ♦ Do not activate the command keys with your nails, and do not use sharp pointed instruments to activate them.
- ♦ Use only original FANEM replacement parts and accessories in order to insure better performance and the safety of the equipment.

4.1 – Physiologic Effects

Phototherapy is commonly used in the treatment of neo-natal hyperbilirubinemia by exposing the patient to a blue spectrum concentrated radiation of visible light, for a period of time to be determined in each case by the attending physician.

The BILITRON® phototherapy unit allows a quick reduction in the serical bilirubin levels , thus reducing the need to use blood transfusions, and reducing treatment time.

- ♦ Even if transitory, collateral effects may be observed, as result of phototherapy. Among them, we could mention the increase in the flow of peripheral blood with the resulting dilation of blood vessels, erythema and an increase in the unfelt loss of body fluids, as well as changes in the intestinal tract.
- ♦ The hydric balance of the patient may be changed through the use of phototherapy.

USER MANUAL
BILITRON® MODEL 3006
THE SUPERLED PHOTOTHERAPY

- ♦ When newborn babies are exposed to the BILITRON® unit, as well as to any other type of phototherapy, they need adequate hydric support during treatment, as well as eyesight protection, besides routine assistance by nurses and doctors.
- ♦ Even if infrequent, some newborn babies exposed to the BILITRON® unit for over 4 days may present a tanned hue on their skins. These patients are normally premature babies, or babies showing serical bilirubin levels in excess of 2 mg%. When the BILITRON® unit is removed, the tanned hue goes away gradually without leaving any sequels.
- ♦ In newborn babies whose skin shows jaundice, it is usual that after of few hours of exposure to the BILITRON® unit, the area being illuminated appears to be lighter as well as whitish in color. This is due to the intense effect of the light over the bilirubinic pigment deposited on the skin.
- ♦ When newborn babies present a serious risk of developing high seric bilirubin levels, as well as in babies in whom the bilirubin concentration keeps on rising in spite of phototherapy, we recommend the use of two BILITRON® unit units.

Double exposure to phototherapy allows a quick reduction of the levels of seric bilirubin, reducing the need of blood transfusions and cutting down on the length of the treatment.

- ♦ The operator's eyesight may be harmed if s/he stays for a long time in the area where the patient is.
- ♦ Patients lying next to the Phototherapy equipment may need to be protected with armored folding screens, protective eyeglasses, etc.
- ♦ During the treatment, Bilirubin photo isomers may have a toxic effect.
- ♦ The seric Bilirubin levels in the patient should be monitored regularly.
- ♦ Liquid infusions and drugs in general should not be stored inside the area affected by radiation.
- ♦ Do not use plates or reflection sheets on the phototherapy equipment, since they may provoke dangerous body temperatures on the patient, and, together with the irradiating heat unit, they may cause serious injuries to the patient.
- ♦ The light source irradiated by the Super Led has extremely **low levels** of infrared rays, which are responsible for the heating of the bodies of the newborns. It is, therefore, recommended that their temperatures are checked as a matter of routine through the use of a calibrated clinical thermometer.

♦ **Environmental conditions:**

When used in open ambiances, and in unheated cradles, it is recommended that the temperature of the newborn babies be monitored to prevent hypothermia.

- ♦ Environmental conditions, such as the movement and the flow of air for instance, may affect the thermal balance of the patient.

4.1.1 – Suggestion to Improve the Efficiency of Phototherapy

- ♦ In order to improve the efficacy of the treatment, the newborn **must be unclothed**, being protected by eyeglasses only. The greater the are exposed to the light, the better will be the results.
- ♦ If, however, it is decided to use diapers, these diapers should be **as small as possible**, so as to cover only the smallest portion possible of the newborn's midriff.
- ♦ Because they prevent the access of the light to the skin of the newborn, large sized diapers reduce considerably the effectiveness of phototherapy.
- ♦ A larger amount of irradiance, allied to a larger bodily surface, will allow to attain more efficient clinical results. This will allow to significantly reduce the time for the phototherapy treatment, and consequently the early release of the patient and a resulting downsize in hospital costs.

4.2 – Electromagnetic compatibility and immunity

We are dealing with the ability of a piece of equipment and/or of a system to function in an electromagnetic environment, without inducing intolerable electromagnetic disturbances for any object or being inside of this ambience, and, on the other hand, allowing it to function without being degraded by electromagnetic disturbances.

This equipment has been protected and tested, and it complies with the following electromagnetic compatibility norms:

- ♦ EN 60601-1-2 (1993)

USER MANUAL
BILITRON® MODEL 3006
THE SUPERLED PHOTOTHERAPY

- ♦ NBR IEC 60601-1-2 (OCT/1997)
- ♦ CISPR11 (AUG/1999)
- ♦ IEC 61000-4-3 (1998)
- ♦ IEC 61000-4-4 (1995)
- ♦ IEC 61000-4-5 (1995)
- ♦ IEC 61000-4-6 (1996)
- ♦ IEC 61000-4-11 (1996)

It is within the parameters that have been specified regarding RF emissions, Immunity, electrostatic discharges, Irradiated and Transient Electromagnetic Radiofrequency Fields (Voltage Sprays and Bursts).

USER MANUAL

BILITRON® MODEL 3006


THE SUPERLED PHOTOTHERAPY

5 – Assembling the equipment


As mentioned previously, the model 3006 BILITRON® – FANEM phototherapy unit offers three assembly configurations.

The equipment comes ready to be installed in incubators (3006 BTI), Heated cradles (3006 BTB) or Movable Spindles (3006 BTP).


In case of configurations using movable spindles, the column shall be assembled over the vase over rotating wheels, using Allen M6 x 10 type screws. Place the pivot of the articulate spindle over the hole on the column, and tighten the holding screw.


3006 BTP

3006 BTI

3006 BTB

Attention: Before using the BILITRON® phototherapy unit, it should be cleaned and disinfected, as per instructions in chapter 7 of this manual, and/or pursuant the protocols that have been adopted by the Hospital Control and Infection Commission.

6 – Operating the Equipment


Attention: This chapter of the User Manual conveys extremely important information to insure the safety and the integrity of the patient, of the user and of the equipment. READ IT CAREFULLY!

- ♦ The plug of the power cable shall be connected to a grounded power outlet, set permanently in the wall, according to the norms and the laws in force dealing with low-voltage electric installations, and the laws pertaining to electric installations within Health Establishments.

Attention: Do not use multiple type extensions or plugs. If the ground is not adequate, do not use the equipment.

- ♦ Connect the power cable plug into the power grid, and turn on the general switch located on the side of the equipment. Check to ascertain whether the lamp got turned on.
- ♦ Position the phototherapy unit in the incubator or in the cradle where the patient lies. Adjust the height and the angle of incidence of the light over the patient.

Attention: In order to prevent damages to eyesight, the Newborn will need to have his/her eyes protected.

Attention: The equipment will reach its normal temperature and irradiation level conditions in less than 15 minutes.

USER MANUAL

BILITRON® MODEL 3006

THE SUPERLED PHOTOTHERAPY

- ♦ The angle of incidence of the light may be adjusted by moving the body of the reflector in the desired direction. The head may be rotated in order to better adjust the focus in the Articulate Spindle 3006 BTP or 3006 BTB models.


Attention: The Bilitron® 3006 BTP model may be only transported with the column in the lowered position.

- ♦ When coupled to the Heated Cradle, the unit will perform the same movements mentioned above, except for height adjustments. This movement will be compensated by the articulate spindle, which will bring the focus closer or it will move it away from the patient.

6.1 – Light Focus – how to adjust it

In order to take better advantage of the light falling on the patient, it is recommended that the focus of light be adjusted until it takes on an elliptical shape. That is, **the light shall range from the thorax to the start of the thighs**. See chart I.

Table I

Distance between patient and reflector D (mm)	Radiation $\mu\text{W}/\text{cm}^2 \cdot \text{nm}$ (on the center of the luminous focus)		
		A (mm)	B (mm)
300	35 - 45	250	230
400	20 - 35	320	270
500	15 - 25	400	370

Obs.: The format of the ellipse will depend on the angle of incidence of the light over the patient.


Figure 1

The Phototherapy unit when used on patients weighting in excess of 2500 grams should be placed some 40 cm away from the baby's body so as to allow the luminous focus to cover the whole body and the start of the thighs of the newborn. When the baby has been born prematurely and underweight, this height may be eventually reduced, provided the luminous focus covers the whole thorax down to the start of the thighs. The intensity of the radiation shall increase as the distance of the light source decreases from the patient, and vice-versa. See Table I and Figure 1.

- ♦ Always keep the focus of light pointed over the body and the start of the thighs of the newborn baby.


Attention: A minimum distance of 300 mm is recommended between the patient and the source of light.

Attention: Do not use the Bilitron® unit when the Power Source is directly in contact with the body of the patient.


USER MANUAL **BILITRON® MODEL 3006** **THE SUPERLED PHOTOTHERAPY**

6.2- Micro-processed Controller Unit

Alphanumeric Display: This display informs the user about the following items:

- 1- A calendar clock displaying the month and the year, the hours, the minutes and the seconds.
- 2- Overall usage time of the lamps (hour counter).
- 3- Radiation level measured in $\mu\text{W}/\text{cm}^2 \cdot \text{nm}$ – when used together with an optical sensor probe (optional).
- 4- Time of exposure.
- 5- Level of radiation computed for a given distance of the light focus (conversion).
- 6- Power level of the lamps when adjusted by the operator.
- 7- Sending of data to the printer or to the computer.

USER MANUAL


BILITRON® MODEL 3006

THE SUPERLED PHOTOTHERAPY

6.3 – Programming Mode

On the upper, front part of the headstock there are four keys used for function programming:

- 1- The MENU key
- 2- The UP key
- 3- The DOWN key
- 4- The CONFIRM key

Adjusting the calendar clock:

Press the MENU key until you see the function "System Config". Press the UP key until you see "Adjust Clock" displayed. Press the CONFIRM key until the digit starts blinking.

Adjust the desired value by activating the UP key. Press the CONFIRM key for the next digit, and make the necessary adjustments successively until all the values have been set. Press the MENU key to get back to the start screen.

Adjusting the Intensity of the Radiation:

The intensity of the light may be adjusted from 20% to 100% of nominal strength.

Press the MENU key until the function "Adjust Radiation" comes up. Press the DOWN or UP key. The display will show the level of the power being used expressed in percentages. Adjust the amount desired by using the UP and DOWN keys. Await the automatic return to the "Adjust Radiation" screen.

Time of Exposure:

Press the MENU key until the word "Therapy" is displayed, the total number of hours used will be informed on the second display line. It will be shown in hours.

To zap the counter, press the MENU key until you see "Zap Therapy" displayed. Press the CONFIRM key to zap the counter.

Lamp Counter:

Press the MENU key until you see the Word "Lamps" displayed.

You will see the lamp usage time displayed in terms of hours on the second line.

The counter showing the number of lamp usage hours will only be turned off when the general switch is activated while the CONFIRM key stays pressed down.

Note: This operation must be accomplished with the Source Module is replaced by a technician authorized to do so by Fanem Ltda.

Level of irradiation:

The Model 3006 BILITRON® phototherapy unit may also be used as a radiometer (*optional item*).

Plug the optical head into the outlet corresponding to the radiometer located on the right side of the BILITRON® box.

Press the MENU key until the value for radiation, expressed in $\mu\text{W}/\text{cm}^2 \cdot \text{nm}$ is displayed.

Note: When the optical sensor has not been coupled, the display will show "Optical Sensor".

Press the CONFIRM key to store the value of the irradiation shown in the display (* after R).

Press the UP key to store the value of the stored irradiation value, together with the date and the time. The file number will be shown very briefly.

The memory has the capacity to store up to 100 irradiation readings.

The file may be viewed on the "Report Radiation" screen. Press the UP or DOWN keys. The file number will be displayed on the first line, and the irradiation amount will be displayed on the second line. If you press the CONFIRM key, you will see the time and the date when the value was stored.

Data Transmission:

Connect a male DB9 cable into the BILITRON® DB9 connector.

In order to print the data regarding the hour counter, the time of exposure, the level of radiation and the calendar clock, press the MENU key until you see the information "Print File" displayed. Press the CONFIRM key to cause the transmission of the data to a printer or to a computer.

USER MANUAL

BILITRON® MODEL 3006

THE SUPERLED PHOTOTHERAPY

7 - Preventive Maintenance, Corrective and Conservation

This section contains cleaning and maintenance instructions. When necessary, disassembly instructions have been provided. The type of maintenance that is now shown in this section must be accomplished by qualified technicians only.

You should routinely inspect the BILITRON® phototherapy unit, and you should replace the accessories before putting the unit into service.

Attention: In order to start maintenance or cleaning procedures, make sure that the BILITRON® unit has been disconnected from the electrical power grid.

7.1 - Replacing the Power Source Module

The average life cycle for the Super LED is of around 20,000 hours, and it is recommended that the irradiation levels be checked periodically, as shown in the usual maintenance procedures.


Attention: Lamps must be replaced whenever they loose 25% of their total irradiance for bilirubin-Ebi.

Use only the original FANEM® Power Source Model, Ref. : 006.053.900

Use a star type screwdriver to replace the Power Source Module. Release the cover screw, open the access cover by pushing down on the latch, remove the Source Module, disconnect the cable from the socket and replace the module by a new one.

Note: the cover is not part of the Source Module.

Attention: Never replace the Source Module when the device is turned on.

Important Notice

Even if similar, and with the same technical specifications, the Super LEDS that have been specified as replacement parts must comply with a series of technical requirements specified in the project of the equipment, such as minimum irradiations, temperature tiers, etc., so as to achieve the maximum degree of safety and to attain the desired clinical effects. Therefore, to not use non original components. The Super LEDS sold by Fanem are 100% tested according to project requirements.

USER MANUAL
BILITRON® MODEL 3006
THE SUPERLED PHOTOTHERAPY

7.2 – Cleaning and Conservation

When the equipment is first delivered, and whenever the BILITRON® unit is not being used, or when a newborn baby is removed from the unit, when s/he gets a medical discharge, or yet, pursuant the procedures laid down by the Hospital Infection Control Commission (CCIH), it is recommended to carry out the cleaning procedure as determined by the disinfection protocol.

Attention: Do not use disinfectants containing abrasive elements or sodium hydrochloride, because these materials may cause the deterioration of the materials used in phototherapy.

In order to always get the best performance out of the BILITRON® unit over time, we recommend that at every 500 hours of use, the external diffusers of the equipment be cleaned by using a dry piece of cloth, or by softly blowing air through them in order to remove the accumulated dust. Never use inadequate chemical products, detergents or soaps to clean the equipment due to the risk of damaging its components.

We recommend the use of Quaternary Ammonia for perfect cleanliness. Since this product is a germicide to be used on plain surfaces, it will not affect the materials used to build the BILITRON® unit, thus allowing the simultaneous cleaning of its surface, fighting bacteria (the instructions found on the label of the product should be followed to the letter in order to achieve perfect asepsis).

Note: Use a cloth moistened with the enough of the product so as to prevent the liquid from flowing inside the equipment.

7.3 - Batteries

This equipment contains two 1.2 Vcc batteries to keep the clock running. These batteries have an average life cycle of 5 years, and they must be replaced by an authorized technician.

Attention: Do not use regular or alkaline batteries. Use only the FANEM® batteries (Ref.: 000.103.026).

Attention: the disposal of the **batteries** should comply with the laws in effect in the country.

7.4 – Protection Fuses

To replace the fuses protecting the equipment (3006-BTP model), the following procedure should be followed:

- ♦ Unplug the power cable.
- ♦ Make sure that the equipment and its accessories have been properly de-energized.
- ♦ With the help of a screw driver, press and rotate the fuse cover to the left.
- ♦ The fuse box cover will be disconnected, and the fuse will found plugged into this cover.

Go ahead with the replacement of the defective fuse by a new fuse, according to its electric specifications, and assemble the set again.

Use	Fuse
127V~ Grid	1A – Type T (3 cm)
220V~ Grid	0,5A – Type T (3 cm)

USER MANUAL

BILITRON® MODEL 3006

THE SUPERLED PHOTOTHERAPY

7.5 – Replacement Parts

Consult Section 3 of this manual for eventual replacement parts: Parts, Pieces and Accessories, together with their respective references.

In order to get the schematics, replacement parts, components or other additional references, please get in touch directly with Fanem Ltda or its business representatives.

The role and the safety of this equipment will only remain under warranty if the inspection, maintenance and repair services are performed by Fanem's Technical Assistance or by individuals who have been properly trained and qualified by Fanem Ltda.

Fanem Ltda. Does not accept any responsibility for damages that may occur to the equipment, and consequently to the patient as function of improper maintenance, which have not been performed by our Technical Assistance Department, or when used as non original replacements, spare parts/accessories.

The materials employed in the design of replacement parts and accessories, and consumption or wear and tear items aim to insure the perfect operation of the equipment according to its original characteristics as well as safety regarding toxicity, flammability and bio-compatibility of the materials employed.

Attention: Use only original FANEM replacement parts.

7.6 – Maintenance Chart

INTERVENTION	TERM	EXECUTING PARTY
Lamp replacement	Constant checking using the radiometer. Replacement when the irradiation is less than 25% Ebi (30 $\mu\text{W}/\text{cm}^2 \text{ nm}$ @ 300 mm of the source – Perpendicular Incidence when the adjustment is at its maximum).	User/Technician
Cleaning	At least at every 90 days	User/Technician
Cleaning the Micro-ventilator	At least at every 90 days	Technician
Routine Calibration /Inspection	6 months	Technician

7.7 – Disposal of

If the equipment or its component parts need to be discarded, and if there is no specific destination that may have been defined by the client, the item being discussed shall be forwarded to the maker or to its legal representative to be disposed of, pursuant the domestic laws that may be in force.

Attention: The disposal of the batteries must be accomplished pursuant the laws in force in the country.

Note: The equipment and/or its component parts must be kept clean and aseptic.

Non compliance with these terms exempts the vendor from any liability regarding possible impacts to the environment and/or people.

USER MANUAL
BILITRON® MODEL 3006
 THE SUPERLED PHOTOTHERAPY

7.8 - Trouble Shooting

Flaw	Possible Causes	Corrective Measures
1) The lamp does not light up	a) Burnt fuse b) Bad cable connection c) Bad lamp connection d) One or more lamps are burned	a) Replace fuse (see item 2.7) b) Connect the grid correctly. c) Connect the lamp correctly. d) Exchange the FANEM® Power Source by a new one (Ref.: 006.053.900).
2) Irradiance is too low	a) the BILITRON® is too far away from the patient b) The collimator is dirty and/or dusty c) The grid voltage is too low d) The lamp is at the limit of its useful life cycle e) Radiometer reading error f) Intensity adjustment less than 100%	a) Check for correct height pursuant the manual. b) Clean up the collimator, pursuant item 7 of the manual. c) Check to see whether the voltage is less than 10% of nominal. (Ref.: 006.053.900) Replace the FANEM® Power Source Model by a new one (Ref.: 006.053.900) e) Check the position of the sensor in the Radiometer, and check its battery. f) Adjust a new value for the effulgence.
3) The BILITRON® unit does not stay at the proper height to perform its function	a) The joints of the articulate spindle are worn down	a) Check to see whether the joint screws are tight. b) Exchange the spindle
4) The focus of the incident light is not perfect (according to the manual)	a) The cover of the Super Led Module is badly positioned.	a) Remove the cover and put it back, checking the way it fits.

If any of the above suggested solutions fail to resolve the problem, call FANEM® ' authorized technician or its closest local representative.

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- ♦ Like all FANEM® 's equipment, this one also carries a full warranty for one (01) year against possible manufacture defective parts (see annexed warranty).
- ♦ For all and any type of maintenance, whether covered or not by the warranty, always look for Technical Assistance that has been authorized by Fanem. Do not allow third, unqualified parties, to damage or to modify the original characteristics of your equipment.
- ♦ The term of warranty for the lamps starts as of their delivery.
- ♦ Always use Fanem®'s original parts.
- ♦ Registered in ANVISA under no. 10.224.620.049
- ♦ Technician Responsible
Engineer Orlando Rossi Filho
CREA/SP 98.435/D
- ♦ Legal Representative

Cinterqual LTDA
 Av. Ressano, 41
 1070-234
 Lisboa – Portugal
 Phone –351-21-386 73 21

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
Eng. José O. C. Flosi
INTERNATIONAL SALES MANAGER
FANEM LTDA

Technical information

REVERSE AND REFLECTED MICROPROCESSED PHOTOTHERAPY OF SUPERLEDS BILITRON BED™. MODEL 4006

THE SUPERLED PHOTOTHERAPY

TÜVRheinland

OCP 0004

Technical Standard – NBR IEC 60601-1
NBR IEC 60601-2-50
NBR IEC 60601-1-2

Review: 00/09
Edition: 06/09

REF: 406.100.550

USER MANUAL
REVERSE AND REFLECTED MICROPROCESSED PHOTOTHERAPY OF SUPERLEDS
BILITRON BED™ - MODEL 4006
THE SUPERLED PHOTOTHERAPY

CONTENTS

1 - Introduction	3
2 - Technical Features	3
2.1 Feeding Voltage and Power	
2.2 - Chart of Spectral Irradiance	4
2.3 - Radiometer Responsiveness Curve	4
2.4 - Dimensions and Weigh	5
2.5 - Symbols	5
2.5.1 - Symbols - Packaging	5
2.6 - Specifications	6
3 - Pieces, Parts and Fixtures	7
3.1 - Reflector Arch - REF: 406.104.321	9
3.2 - Clear-gel Mattress - REF: 406.106.321	9
3.3 - Radiation Module - REF: 406.101.600	9
3.4 - Micro-fans - REF: 000.641.026	9
3.5 - Optical Probe - REF: 620.058.600	9
3.6 - Patient's Skin Temperature Sensor - REF: 092.059.600	10
3.7 - Patch Kit	10
3.8 - Eye Protection for Neonatal Phototherapy	10
4 - Precautions, Restrictions and Warnings	10
4.1 - Physiological Effects	12
4.1.1 - Suggestion for Improving the Efficacy of Phototherapy	12
4.2 - Electromagnetic Compatibility and Immunity	12
4.3 - Special Care with the Clear-gel Mattress	15
5 - Installation of Equipment	16
6 - Operation of Equipment	17
6.1 - Lateral Panel I	18
6.2 - Microprocessed Control Unity	18
6.3 - Functions of the control unit	18
6.4 - Protection against overheating	22
7 - Preventive and Corrective Maintenance and Conservation	23
7.1 - Replacement of the Radiation Module	23
7.2 - Cleaning and Conservation	23
7.3 - Batteries	24
7.4 - Protection Circuit Breakers	24
7.5 - Replacement Parts	24
7.6 - Maintenance Table	25
7.7 - Disposal	25
7.8 - Diagnosis of Failures	26
8 - Warranty and Technical Assistance	26

USER MANUAL
REVERSE AND REFLECTED MICROPROCESSED PHOTOTHERAPY OF SUPERLEDS
BILITRON BED™ - MODEL 4006
THE SUPERLED PHOTOTHERAPY

1 - Introduction

Phototherapy is commonly used in the treatment of neonatal hyperbilirubinemia, subjecting the patient to a radiation exposure with a wavelength at the blue band of the spectrum of visible light, for a timeframe determined by the attending physician, in accordance to the circumstances.

The introduction of phototherapy in the treatment of neonatal bilirubinemia, since its inception in the late 1950s and throughout the major evolution in the types of lamps applied, minimized the risks but did not eliminate the adverse effects of ultraviolet and infrared rays on the patient's body. The treatment of bilirubinemia with light began with low-intensity fluorescent lights; and further employed fluorescent lights with emission at the blue band of the spectrum. Later, the halogeneous lights and a variation of optical fibers were introduced. The LEDs (light-emitting diodes) have also been introduced in phototherapy; however, due to the low emission of light of conventional LEDs, hundreds of them are required to obtain a satisfactory result.

The emergence of the Super LED with high-intensity radiation at the blue band of the spectrum, with a physical-chemical composition different from the one of the conventional LEDs, enables the building of equipments with fewer LEDs, rendering them more compact and lighter. The Super LEDs emit safe radiation (according to rule IEC 825 Class 2), completely eliminating the ultraviolet and the infrared filters (see chart of emissions), with the advantage of dramatically decreasing the risk of sunburn, rashes, and imperceptible water loss.

The BILITRON BED™ Model 4006, which employs Super LEDs as radiation source, might be considered the 7th generation of phototherapy in the world's market. This is an electronic phototherapy, microprocessed, capable of eliminating unwanted radiation in the spectrum of the ultraviolet, as well as from the infrared, incorporating a larger area ellipse with less loss at the edges, making the application more anatomic to the patients' bodies.

The BILITRON BED™ 4006 employs 17 Super LEDs and uses the reflector arch, which reflects to the thoracic region the rays that did not reach the back region, increasing the efficiency of the phototherapy.

The microprocessed electronic present in this phototherapy allows resources that were not previously feasible, such as: Radiometer with built-in optical probe (optional probe), built-in patient's skin temperature sensor (optional sensor), treatment time, hour meter, switch for light intensity, and clock / calendar.

Only qualified personnel should operate this product.

In order to preserve its efficiency and taking into account the nature of the works for which the BILITRON BED™ 4006 is proposed, certain precautions are necessary for the handling and the operation of the product. This User Manual provides general instructions for the installation, the use and the maintenance by the operator and diagnosis for failures. FANEM™ will not be liable if the User does not operate the unit in accordance with the instructions, does not follow the maintenance warnings and recommendations in this manual or if the User performs any repair with non-authorized components. Only qualified personnel should perform Calibration and repair. Any additional information might be secured at your local distributor.

This manual should be read, fully understood and be readily accessible to anyone who works with the Phototherapy BILITRON BED™ Model 4006 unit. When not in use, this User Manual should be kept together with the BILITRON BED™. If there is any information that you do not understand, please contact the representatives of FANEM™ for further clarifications


2 - Technical Features

2.1 - Feeding Voltage and Power

Input (Volts rms) V ~	Frequency Hertz (Hz)	Volt Ampere Power (VA)	Current I (A)	REF: Fanem Code
100-240	50/60	125	1	406.106.900 406.107.900


USER MANUAL
REVERSE AND REFLECTED MICROPROCESSED PHOTOTHERAPY OF SUPERLEDS
BILITRON BED™ - MODEL 4006
THE SUPERLED PHOTOTHERAPY

2.2 - Chart of Spectral Irradiance

[Legend]
Standardized Spectral Irradiance
Wavelength (nm)

2.3 - Radiometer Responsiveness Curve

[Legend]
Calibration Curve of the Radiometer
Wavelength (nm)

USER MANUAL **REVERSE AND REFLECTED MICROPROCESSED PHOTOTHERAPY OF SUPERLEDS** **BILITRON BED™ - MODEL 4006** **THE SUPERLED PHOTOTHERAPY**

2.4 - Dimensions and Weight

♦ **Dimensions**

Bilitron Bed™ 4006	Width (cm)	Length (cm)	Height (cm)
Unit only	40.0	59.5	41.5
Equipment with cart	46.8	83.8	120

Cart wheels for the Bilitron Bed™ 4006, with cart	4 "with brake
---	---------------

♦ **Weight**

Bilitron Bed™ 4006	Net Weight (kg)	Packaged Weight (kg)
Unit only	13.4	14.9
Equipment with cart	31.9	33.9

2.5 - Symbols

	BF Type Equipment Insulation Class I	IPX4
	Caution: Hot Surface	 Caution: Risk of Electric Shock
	Caution: See manual	Mandatory Grounding
	OFF (No Feeding Electrical Power)	 ON (With Feeding Electrical Power)
	Eye protection of patient	

2.5.1 - Symbols - Packaging

	Fragile		This Way Up
	Protect from Sunlight		Protect from Rain
	Stacking Limit		Caution: Consult accompanying documents

USER MANUAL **REVERSE AND REFLECTED MICROPROCESSED PHOTOTHERAPY OF SUPERLEDS** **BILITRON BED™ - MODEL 4006** **THE SUPERLED PHOTOTHERAPY**


2.6 - Specifications

Protection Circuit Breaker:

1.2 A Type "F" - 2 x 0.5 cm (2 units)

Characteristics of Irradiance:

Maximum Total Irradiance for Bilirubin (max)..... Ebi 47.5 $\pm$ 5 μ W/cm².nm (7.1 $\pm$ 0.7 mW/cm²)
 Average Total Irradiance for Bilirubin..... 32.6 μ W/cm².nm (4.9 mW/cm²)
 Dimensions of the elliptical spot 320 x 200 mm
 Irradiation adjustment range 0% to 100% of maximum radiation
 Variation of the irradiation adjustment..... each 10%
 Distance between the radiation source and the actual surface 12cm

Radiometer:

Range of measures 0 to 200 μ W/cm².nm
 Resolution 0.1 μ W/cm².nm
 Accuracy of reading $\pm$ 2 μ W/cm².nm

Noise < 52 dBA (site 45dBA)

Environmental Conditions:

Operation's Temperature range..... 18 °C / 64 °F to 28 °C / 82.4 °F (ambient)
 Operation's Humidity Range 10% to 95% (non-condensing)

Environmental Conditions for Transportation and Storage

Environmental Temperature 5 °C / 41 °F to 55 °C / 131 °F
 Relative Humidity 30% to 75% (non-condensing)

NOTE. Conditions specified with the property packed equipment.

USER MANUAL
REVERSE AND REFLECTED MICROPROCESSED PHOTOTHERAPY OF SUPERLEDS
BILITRON BED™ - MODEL 4006
THE SUPERLED PHOTOTHERAPY

3 - Pieces, Parts and Fixtures

BILITRON BED™ Model 4006 Phototherapy contains a module with 17 Super LEDs, and is equipped with two micro-fans for cooling and a microprocessed controller with alphanumeric display with keys for the programming of the functions.


The equipment has with the following items:

- Acrylic-reflector arch (REF: 406.104.321)
- Highly transparent gel mattress (REF: 406.106.321)
- Feeding Cable (REF: 006.161.600)


The BILITRON BED™ 4006 might be presented in two types of assembly: as equipment for use on tables or benches, and as equipment with cart.

The BILITRON BED™ 4006 with cart has 4-4" wheels, for easy transportation, and a column with a wireframe basket to place objects.


Two optional items might be used in the BILITRON BED™ 4006: Optical Probe (REF: 620.058.600) for measuring the irradiance (Radiometer), and patient's Skin Temperature Sensor (REF: 092.059.600).

Feeding Cable - REF: 006.161.600

Optical Probe - REF: 620.058.600 (optional)

Skin Temperature Sensor - REF: 092.059.600
(optional)

USER MANUAL


REVERSE AND REFLECTED MICROPROCESSED PHOTOTHERAPY OF SUPERLEDS

BILITRON BED™ - MODEL 4006


THE SUPERLED PHOTOTHERAPY

Presentations:

**1 - Bilitron Bed™ 4006 Phototherapy,
for use on bench
REF: 406.107.900**

**2 - Bilitron Bed™ 4006
Phototherapy, mounted on optional
cart and drawers**

Item	REF:	Description	Note
1	406.107.900	Phototherapy Bilitron Bed™ 4006 Phototherapy for use on bench	-
2	406.106.900	Bilitron Bed™ 4006 Phototherapy cart with fixed height cart	optional
3	006.161.600	Feeding cable for Bilitron Bed™ 4006	Supplied
4	406.106.321	Clear-Gel Mattress	Supplied
5	620.058.600	Optical Sensor Probe	optional
6	092.059.600	Patient's Skin Temperature Sensor	optional
7	406.104.321	Acrylic Reflector Arch	supplied
8	406.108.900	Bilitron Bed™ 4006 Phototherapy with cart with ergonomic column	optional
9	406.102.300	Wireframe basket for the Bilitron Bed™ 4006 cart (Not applicable, with drawer (s))	supplied
10	090.079.600	Chest with two drawers for the Bilitron Bed™ 4006 cart.	optional
11	090.078.600	Chest with one large drawer for the Bilitron Bed™ 4006 cart	optional

USER MANUAL
REVERSE AND REFLECTED MICROPROCESSED PHOTOTHERAPY OF SUPERLEDS
BILITRON BED™ - MODEL 4006
THE SUPERLED PHOTOTHERAPY

3.1 – Reflector Arch - REF: 406.104.321


Built in white acrylic and internally coated with reflective film to use the peripheral light that would normally be lost, therefore providing a better treatment efficacy. The function of the reflector arch is to reflect to the frontal surface of the patient the radiation that did not reach the back region.

Fixture for exclusive use with the equipment

3.2 - Clear-gel Mattress - REF: 406.106.321

Completely transparent, it provides the desired passage of light rays that are beneficial for the treatment, distributed in a large body surface area. It has an oval concavity that provides more comfort to the patients and facilitates their positioning.

Fixture for exclusive use with the equipment

Caution: The Clear-Gel Mattress requires special care to ensure the quality of the care and to avoid accidents. Refer to the item 4 of this manual.

3.3 Radiation Module- REF: 406.101.600

A set containing 17 Super LEDs with collimators, mounted on an aluminum dissipater with two micro-fans to cool of the dissipater.

The average life cycle of the radiation module is of, approximately, 10,000 hours; however, we recommend the continuous monitoring of the irradiance with the FANEM Radiation Monitor, as per the regular maintenance procedures.

Caution: Always use the original FANEM, Radiation Module because the light sources are controlled and individually selected as per their original features of the project.

3.4 – Micro-fans - REF: 000.641.026

Two micro-fans, located on the bottom of the box, operate constantly while the machine is on, yielding the dissipation of the heat generated by the radiation module and thus prolonging its life cycle.

Caution: Do not obstruct the airway of the micro-fans, as the module might burn.

3.5 - Optical Probe - REF: 620.058.600

Extension Probe, bearing a head with an optical sensor, used to measure the irradiance emitted by the BILITRON BED™ Model 4006 Phototherapy. The value of the irradiance measured is displayed on the alphanumeric device.

Caution: The optical probe is an optional item.

Fixture for exclusive use with the equipment

USER MANUAL
REVERSE AND REFLECTED MICROPROCESSED PHOTOTHERAPY OF SUPERLEDS
BILITRON BED™ - MODEL 4006
THE SUPERLED PHOTOTHERAPY

3.6 - Patient's Skin Temperature Sensor - REF: 092.059.600

Sensor used to measure the temperature of the patient's skin. The temperature value is displayed on the alphanumeric device.

Caution: The patient's skin temperature sensor is an optional item.

Fixture for exclusive use with the equipment

3.7 - Patch Kit

Antiallergic patch especially designed for the perfect fixation of the infant's skin temperature sensor.

Patch Kit	REF:
50 units	086.076.600
100 units	086.068.600

Caution: The patch kit is an optional item.

Fixture for exclusive use with the equipment

3.8 - Eye Protection for Neonatal Phototherapy

A single piece, with design that wraps around the entire head of the newborn, made of antiallergic material, adjusted by Velcro® thus preserving its original positioning. Contains, internally, a photo-opaque pad, ergonomically designed to prevent the exposure of the eyes of newborn to the light of the phototherapy.

Ordering Information:

REF: 150.079.600 - Large: Head circumference 30-38 cm (packages with 10 und.)

REF: 150.078.600 - Average: Head circumference 24-33 cm (packages with 10 und.)

REF: 150.077.600 - Small: Head circumference 20-28 cm (packages with 10 und.)

** Fixture for exclusive use with the equipment.*

Caution: - The manufacturer recommends single use.
- The Product supplied is not sterilized

4 - Precautions, Restrictions, and Warnings

Caution: This chapter of the User Manual contains extremely important information to ensure the safety and the integrity of the patient, of the User and of the Equipment. **READ CAREFULLY !**

- ♦ Check if the network to which the equipment will be connected has features to withstand the electrical conditions of voltage and Power of the equipment, which are indicated on the label affixed to the appliance.

USER MANUAL
REVERSE AND REFLECTED MICROPROCESSED PHOTOTHERAPY OF SUPERLEDS
BILITRON BED™ - MODEL 4006
THE SUPERLED PHOTOTHERAPY

- ♦ The plug of the feeding cable must be connected to a grounded outlet, permanently set on the wall, according to the current rules and legislations of low-power electrical installations for Health Care Facilities.

Caution: Do not use extensions or multiple sockets. If the grounding is not perfect, do not use the equipment.

Caution: To avoid damages to the eye, the newborn will need eye protection.

- ♦ The BILITRON BED™ 4006 should be used according to the ideal standards for environmental comfort in nurseries, i.e., from 23 to 27 ° C / 73.4 to 80.6 ° F.
- ♦ A misused unit of phototherapy might cause serious risks to the newborn. Only trained and qualified personnel, who know the risks and the benefits of its use, should operate this unit.
- ♦ Do not use the phototherapy unit if any of the functions of the unit does not operate properly . Please, call a service technician.
- ♦ It is not advisable to leave the patient unattended under the equipment.
- ♦ Constantly monitor the patient's temperature using a clinical thermometer, or with the BILITRON BED™ 4006 patient's skin temperature sensor (optional).
- ♦ Thermotherapy devices (thermal blankets, Heated Mattresses, Radiation Heaters, etc.) used with the BILITRON BED™, might increase the temperature of the patient's body to dangerous levels. We recommend that these equipments operate in the skin-controlled modality (servo controlled), or in the manual modality, where the heat power of the heater must be reduced according to the patient's body temperature, as measured by a clinical thermometer.
- ♦ Always use the original light sources of FANEM™. The use of different types of sources might alter the conditions of irradiance and temperature specified for the equipment.
- ♦ We recommend the use of a FANEM Radiation Monitor to confirm the perfect condition of the radiation sources and to allow adjustment of the desired level of irradiance.
Note: The internationally accepted minimum value for a treatment is of 4 $\mu\text{W}/\text{cm}^2 \text{ nm}$ (minimum value in accordance with Mims (1)), and this equipment offers far greater values and therefore more efficacy .
- ♦ Do not place cloths or towels on top of the BILITRON BED™ box, as this might cause an increase in temperature and the consequent shutdown of the radiating source.

Caution: Do not block the nozzles of the air circulation outlet on the back of BILITRON BED™, as this might cause an increase in temperature and the consequent shutdown of the radiating source.

Caution: Do not use the phototherapy equipment in the presence of flammable anesthetics, oxidizing gases or cleaning agents that might cause combustion.

Caution: This equipment is not suitable for use with flammable anesthetic mixes.

- ♦ Typically, the life cycle of the original sources is of 10,000 hrs; however, we recommend constant monitoring of the irradiance with the FANEM Radiation Monitor, according to the regular maintenance procedures for a better assessment of the actual efficacy of the light sources. Thus, the sources of irradiance (Super LEDs) should be replaced whenever they reach a loss of 25% of the total irradiance for Bilirubin-Ebi.
- ♦ Do not activate the control keys with your nails, and do not use sharp objects to trigger them.
- ♦ Use only FANEM original parts and fixtures to ensure a better performance and the safety of the equipment.

USER MANUAL
REVERSE AND REFLECTED MICROPROCESSED PHOTOTHERAPY OF SUPERLEDS
BILITRON BED™ - MODEL 4006
THE SUPERLED PHOTOTHERAPY

4.1 - Physiological Effects

Phototherapy is commonly used in the treatment of neonatal hyperbilirubinemia, subjecting the patient to a radiation exposure with a wavelength at the blue band of the spectrum of visible light, for a timeframe determined by the attending physician, in accordance to the circumstances.

The BILITRON BED™ 4006 Phototherapy allows the rapid reduction in serum bilirubin, thereby reducing the need for whole blood exchange transfusion and reducing the duration of treatment.

- ♦ However, we might observe the occurrence of transient effects might as a result of use of phototherapy. Among them are, mainly, the increase in peripheral blood flow with consequent vasodilatation, erythema and increased imperceptible water loss and alteration of bowels movement.
- ♦ The hydrous balance of the patient might be altered with the use of phototherapy.
- ♦ Therefore, newborns undergoing BILITRON BED™ Phototherapy, as well as using any other type of phototherapy, need, during the treatment, the proper hydrous aid, and adequate eye protection, in addition to the routine care of the nurses and of the physician.
- ♦ Although infrequent, some newborns undergoing BILITRON BED™ Phototherapy for more than 4 days might display a tanned skin. These patients are usually premature newborns and / or those whose direct bilirubin serum levels are above 2 mg%. With the withdrawal of BILITRON BED™ Phototherapy the color of the skin will gradually fade, without sequelae.
- ♦ In newborns, whose skin is very icterical, is common that after a few hours of Biltron Bed™ Phototherapy the illuminated area becomes more pallid and whitish. This is due to the intense action of light on the bilirubin pigment deposited in the skin.
- ♦ The eyes of the "operator" might suffer damages arising out of the long permanence in the area of the patient.
- ♦ Patients next to Phototherapy equipment might need protection, such as armored screens, goggles, etc..
- ♦ During treatment, the photoisomers of the bilirubin might cause toxic effects.
- ♦ Serum bilirubin of the patient should be measured regularly.
- ♦ Liquid infusions and drugs in general, should not be stored within the area of radiation.
- ♦ Do not use in the phototherapy equipment, reflection plates or sheets in addition to the original reflector arch of the equipment, as it might lead to dangerous body temperatures to the patient, and together with the unit of irradiant heat might cause serious damage to the patient.
- ♦ The source of light radiated from the SuperLED has extremely low levels of infrared rays, which are responsible for the warming-up of the body of the newborn, so we recommend the routine assessment of the temperature of the newborn with a verified clinical thermometer or with the original BILITRON BED™ patient's skin temperature sensor (optional).

Environmental Conditions:

- ♦ When used in open environments, we recommend the assessment of the newborn's temperature to prevent hypothermia.
- ♦ Environmental conditions such as, for instance, air movements and flows might affect the thermal balance of the patient.

4.1.1 - Suggestion for Improving the Efficacy of Phototherapy

- ♦ To improve the efficacy of the treatment, the newborn should be naked, wearing, only, the eye protection. The larger the area of skin exposed to the light, the better the result.
- ♦ However, if you choose to use diapers, it should be as small as possible, to cover only a small part of the perineum of the Newborn.
- ♦ Large diapers large, by preventing the action of light on the skin of newborns, greatly reduce the efficacy of the phototherapy.
- ♦ Higher irradiance, combined with greater body surface exposition, will achieve better clinical efficacy. This will allow a significant reduction in the length of the phototherapy treatment, and, therefore, allow early hospital discharge and reduction in hospital costs.

4.2 - Electromagnetic Compatibility and Immunity

It is the capability of the equipment and/or the system to operate in an electromagnetic environment without introducing intolerable electromagnetic disturbances to anything in this environment and, on the other hand, to operate without degradation, in the presence of an electromagnetic disturbance.

USER MANUAL

REVERSE AND REFLECTED MICROPROCESSED PHOTOTHERAPY OF SUPERLEDS

BILITRON BED™ - MODEL 4006

THE SUPERLED PHOTOTHERAPY

This equipment was designed, tested and meets the following EMC standards for electromagnetic compatibility:

- ♦ EN 60601-1-2
- ♦ NBR IEC 60601-1-2
- ♦ CISPR11
- ♦ IEC 61000-3-2
- ♦ IEC 61000-3-3
- ♦ IEC 61000-4-2

If within the parameters recommended for RF Emissions, Immunity, Electrostatic Discharge, Radio Frequency Electromagnetic Fields Irradiated, and Transients (Power Blasts and Outbreaks).

Caution: Mobile and portable RF communication equipment might affect Electrical-Medical Equipment.

Caution: The use of non-original fixtures, transducers, sensors and network cables might result in additional emissions or in a decrease of the immunity of the equipment.

Guidelines and Manufacturer's Statement - Electromagnetic Emissions

BILITRON BED™ 4006 Phototherapy is destined for the use in the following electromagnetic environment. We recommend that the client or the User of the Phototherapy ensure that it is used in such an environment.

Emission essays	Compliance	Electromagnetic environment – Standards
CISPR 11 RF Emissions	Group 1	Phototherapy uses RF energy only for its internal functions. However, its RF emissions are very low and are unlikely to cause any interference in nearby electronic equipment.
CISPR 11 RF Emissions	Class B	Phototherapy is suitable for the use in all facilities other than the households, and might be used in residential facilities and in those directly connected to the public net of low power electrical energy that feeds the facilities for household use, provided that the following warning is understood: Warning: This equipment is intended for use only by health care professionals. This equipment might cause radio interference or disrupt the operation of nearby equipment. It might be necessary to adopt procedures for mitigation, such as the realignment or the relocation of the Phototherapy or the local shielding.
IEC 61000-3-2 Harmonic Emissions	Class A	
Emissions due to Power fluctuation / flicker IEC 61000-3-3	according	

USER MANUAL
REVERSE AND REFLECTED MICROPROCESSED PHOTOTHERAPY OF SUPERLEDS
BILITRON BED™ - MODEL 4006
THE SUPERLED PHOTOTHERAPY

Guidance and Manufacturer's Statement - Electromagnetic Immunity

BILITRON BED™ 4006 Phototherapy destined for the use in the following electromagnetic environment.
 We recommend that the client or the User of the Phototherapy ensure that it is used in such an environment.

Immunity essay	ABNT NBR IEC 60601 Essay Level	Level of Compliance	Electromagnetic Environment – Standards
IEC 61000-4-6 RF- Conducted	3 V ms 150 KHz to 80 MHz	10V	Portable and mobile RF communication equipment should not be used in the proximities of any Bilitron BED™ 4006 Phototherapy including cables, bearing a separation distance less than the one that is recommended, calculated from the equation applicable to the transmitter frequency. Recommended Separation Distance $d = 0.35 \cdot P^{1/2}$ $d = 0.35 \cdot P^{1/2} \text{ 80MHz to 800 MHz}$ $d = 0.7 \cdot P^{1/2} \text{ 800 MHz to 2.5 MHz}$ where P is the nominal maximum outlet power of the transmitter in watts (W), according to the manufacturer of the transmitter, d is the recommended separation distance in meters (m) We recommend that the field intensity established by the RF transmitter, as determined by an on-site electromagnetic inspection, be smaller than the compliance level in each frequency band. ^a Interference might occur around the equipment marked with the following symbol: 
IEC 61000-4-3 RF- Conducted	3 V ms 80 MHz to 2.5 GHz	10 V / m	

NOTE 1: In 80 MHz and 800 MHz, the higher range frequency applies.

NOTE 2: These standards might not apply to all the situations. Electromagnetic propagation is affected by the absorption and the reflection of structures, objects and people.

^a The intensities of the field established by fixed transmitters, such as radio-base stations, telephone (cell/cordless) and mobile terrestrial radios, amateur radio, AM and FM – TV broadcast theoretically, might not be predicted accurately. To assess the electromagnetic environment due to fixed RF transmitters, we recommend that an electromagnetic inspection of the site. If the measure of the field intensity when the BILITRON BED™ 4006 Phototherapy is used exceeds the above-mentioned level of compliance, the BILITRON BED™ 4006 Phototherapy should be observed to see if its operation is normal. If abnormal performance is observed, additional procedures might be needed, such as the realignment or the relocation of the Phototherapy BILITRON BED™ 4006.

^b Above the frequency range of from 150 kHz to 80 MHz, the field intensity should be lower than 10 V / m.

USER MANUAL

REVERSE AND REFLECTED MICROPROCESSED PHOTOTHERAPY OF SUPERLEDS

BILITRON BED™ - MODEL 4006

THE SUPERLED PHOTOTHERAPY

Recommended separation distances between the portable and the mobile RF communications and the BILITRON BED™ 4006 Phototherapy

BILITRON BED™ 4006 Phototherapy is destined for the use in electromagnetic environments in which radiated RF disturbances are controlled. The User might help in the prevention of electromagnetic interference by maintaining a minimum distance between the portable and the mobile RF communication equipment (transmitters) and the BILITRON BED™ 4006 Phototherapy as recommended below, according to the maximum outlet power of the communication equipment.

Nominal outlet power of the transmitter W	Separation distance according to the frequency of the transmitter m		
	150 kHz to 80 MHz $d = 0,35 P^{1/2}$	80 MHz to 800 MHz $d = 0,35 P^{1/2}$	800 MHz to 2,5 GHz $d = 0,7 P^{1/2}$
0,01	0,04	0,04	0,07
0,1	0,11	0,11	0,22
1	0,35	0,35	0,7
10	1,11	1,11	2,21
100	3,50	3,50	7,00


For transmitters with nominal maximum outlet power not listed above, the separation distance d in meters (m) might be determined using the equation applicable to the frequency of the transmitter, where P is the nominal maximum outlet power of the transmitter.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range, applies.

NOTE 2: These standards might not apply to all the situations. Electromagnetic propagation is affected by the absorption and the reflection of structures, objects and people.

4.3 - Special Care with the Clear-gel Mattress

- ♦ The **Clear-gel Mattress** is a very important and sensible part of the BILITRON BED®. Treat it as special, and with great care.
- ♦ To keep it clean, use a soft cloth with **warm water and mild soap**. Do not use chemicals or alcohol to clean it.
- ♦ The use of **iodine solutions** employed in routine procedures for patients, might yield a "yellow / amber" color to the gel mattress, or even the prolonged use of the equipment might cause a **photosensitivity reaction** in the mattress, without a significant loss of radiance.
- ♦ Do not fold or roll up the gel mattress,, or it might lose its original shape, irreversibly, or rupture its protective film, thereby exposing the gel. It is important that the mattress be handled, always, in a horizontal position, as shown in the following figure.

Caution: Always handle the gel mattress in a horizontal position and never bend it or leave it upright, or the shape of the mattress might change, irreversibly.


- ♦ Due to the characteristics of the gel, if the mattress remains in an upright position, the gel will accumulate, slowly, at the bottom, as a result of the gravitational action, deforming the mattress. Therefore, you should keep the mattress, always, in a horizontal position.

USER MANUAL
REVERSE AND REFLECTED MICROPROCESSED PHOTOTHERAPY OF SUPERLEDS
BILITRON BED™ - MODEL 4006
THE SUPERLED PHOTOTHERAPY


- ♦ Do not use pointed or sharp objects in contact with the mattress. Beware of nails, rings, and bracelets while handling of the mattress, as you might puncture it.
- ♦ If the mattress is accidentally punctured, do not try to reuse it, because of the contamination of the gel, or as patient might swallow the gel. Please contact our Technical Support, so we can replace the defective mattress with a new product.

5 - Installation of Equipment

FANEM BILITRON BED™ Model 4006 Phototherapy, comes in two assembly configurations, as mentioned earlier.

BILITRON BED™ for use on bench

BILITRON BED™ with optional cart and drawers

Caution: Before using BILITRON BED™ Phototherapy, you must clean the equipment, and disinfect it as instructed under the item 7 of this manual and / or according to the protocols used by the Control and Hospital Infection Commission of the service.

USER MANUAL
REVERSE AND REFLECTED MICROPROCESSED PHOTOTHERAPY OF SUPERLEDS
BILITRON BED™ - MODEL 4006
THE SUPERLED PHOTOTHERAPY

6 - Operation of the Equipment


Caution: This chapter of the User Manual contains extremely important information to ensure the safety and the integrity of the patient, of the User, and of the Equipment. **READ CAREFULLY !**

- ♦ The plug of the feeding cable must be connected to a grounded outlet, permanently set on the wall, according to the current rules and legislations of low-power electrical installations for Health Care Facilities.

Caution: Do not use extensions or multiple sockets. If the grounding is not perfect, do not use the equipment.

- ♦ Connect the plug of the power cable to the electrical net, and turn on the main switch on the lateral side of the equipment. See if the lamps are on.
- ♦ Remove the reflector arch.
- ♦ Set the desired level of radiation (from 0 to 100%). We recommend the use of the Fanem (REF: 620.058.600) optical probe to assess and adjust the irradiance.
- ♦ Position the patient on the Clear-gel Mattress in the direction indicated in the marks of the equipment. As indicated in the following figure, the patient's head should be on the side where there are 2 Super LEDs at the extremities, which favor the treatment at the area of the head. Fix the skin temperature sensor (REF: 092.059.600) on the patient (optional). Reposition the reflector arch.

Caution: To avoid damage to the eye, the newborn will need eye protection.

Caution: The equipment has, in the Clear-gel Mattress, an oval concave area for optimal patient positioning on the actual area of treatment. We recommend the routine check of the position of the patient to ensure optimal positioning.

Caution: Equipment reaches its stable conditions of temperature and irradiance levels in less than 15 minutes.

USER MANUAL

REVERSE AND REFLECTED MICROPROCESSED PHOTOTHERAPY OF SUPERLEDS


BILITRON BED™ - MODEL 4006

THE SUPERLED PHOTOTHERAPY

6.1 – Lateral Panel

On the right side of the BILITRON BED™ is an access panel to the following functions:


- Input for the feeding cable
- On/off switch
- Circuit Breakers
- Input for the optical probe
- Input of the patient's skin temperature sensor
- RS-232 Interface

6.2 – Microprocessed Control Unity

Alphanumeric Display: informs the following items to the User:


- 1- Calendar clock with day, month, year, hours, minutes.
- 2- Patient's Skin temperature - when using the optional patient's skin sensor).
- 3- Level of radiation in $\mu\text{W}/\text{cm}^2 \cdot \text{nm}$ - when using the optical probe (optional).
- 4- Total time of use of the lamps (hour's meter).
- 5- Time of treatment.
- 6- Power level the lamp when adjusted by the operator.

6.3 - Functions of the control unit

On the front of Bilitron Bed™, there are four keys for programming the functions, namely:

- 1- ENTER key
- 2- DOWN key
- 3- UP key
- 4- MENU key



USER MANUAL
REVERSE AND REFLECTED MICROPROCESSED PHOTOTHERAPY OF SUPERLEDS
BILITRON BED™ - MODEL 4006
THE SUPERLED PHOTOTHERAPY

The MENU allows the switching in between the various system functions. The ENTER key is used to select a function or to confirm a change you made. The UP and DOWN keys allow you to change the value of a parameter.

The main screen Bilitron Bed™ might be one of the two shown in the figure. If the optional of the patient's skin temperature sensor is not enabled, the main screen displays the date, the time, and the level of radiation (from 0 to 100%). If the optional temperature sensor of the patient is on, it will display the date, the time and the temperature provided by the sensor.

Typically, when the display is on a screen that is not the main screen and no key is pressed, the system automatically returns to the main screen after a few seconds.

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The functions in the alphanumerical display are organized as follows:

- **Main screen:** indicates the current date and time and the level of radiation (or patient's skin temperature sensor, if this option is enabled)
- **Level of radiation:** indicates the level of radiation (0 to 100%)
- **Time of treatment.** indicates the length of the phototherapy, in hours and minutes
- **Total time of use:** indicates the total time of use of the Super LEDs in hours and minutes
- **Measurement of irradiance:** displays the value of irradiance (in $\mu\text{W}/\text{cm}^2 \cdot \text{nm}$) measured by the optical probe (optional)
- **Configurations:** configurations and system options
 - **Choose language:** allows you to select between Portuguese, English or Spanish
 - **Set Date and Time:** adjusts the current date and time
 - **Restart treatment time**

Monitoring the temperature of the patient:


If you purchased the Bilitron BED™ with the option of the patient's skin temperature sensor, the value read by the sensor will be displayed on the main screen when the sensor is plugged into its respective socket on the lateral side of the equipment.

Caution: We recommend the routine assessment of the temperature of the newborn with a verified clinical thermometer, or with the BILITRON BED™ patient's skin temperature sensor

USER MANUAL
REVERSE AND REFLECTED MICROPROCESSED PHOTOTHERAPY OF SUPERLEDS
BILITRON BED™ - MODEL 4006
THE SUPERLED PHOTOTHERAPY

If the temperature sensor is disconnected from the lateral side, then the screen will display horizontal lines instead of the value of temperature, as follows .

If you want to use the BILITRON BED™ patient's skin temperature sensor in equipment purchased without this feature, ask to an authorized FANEM technician to carry on the procedure to enable this option.

Adjustment of the level of radiation:

The light intensity might be adjusted from 0 to 100% of the total power.

Press MENU up to the "Level Radiation" function. On the screen, the indication of the power level will appear. Adjust to your desired level with the UP or DOWN keys.

The level of radiation might also be adjusted on the irradiance-measuring screen. See the section "Measure of irradiance."

Time of treatment.

Press the MENU key up to "treatment." The length of the treatment will appear on the second line of the display, in hours and minutes.

To restart the treatment clock , press MENU up to "Settings - Press ENTER. Press the ENTER key and then press the MENU key up to "Restart time of treatment?". Press the ENTER key to start this function. Press ENTER to restart the time of treatment or MENU to cancel it.

USER MANUAL
REVERSE AND REFLECTED MICROPROCESSED PHOTOTHERAPY OF SUPERLEDS
BILITRON BED™ - MODEL 4006
THE SUPERLED PHOTOTHERAPY

Total time of use of Super LEDs:

Press MENU up to "Total time of use." The total time of use of the lamps will be informed on the second line of the display, in hours and minutes. The counting of the total usage time will take place, only, if the radiation level is set to a value above 0%.

When the total time reaches 10,000 hours, the screen of "Total time of use" will flash the message "Check Super LEDs," indicating that you must assess the irradiance of the Super LEDs to determine the need to change the Radiation Module.

Although the average life cycle of the Radiation Module is of, approximately, 10,000 hours, we recommend the constant monitoring of the irradiance with FANEM Radiation Monitor, in accordance with the regular maintenance procedures to better assess the actual efficacy of the lamps. Thus, the sources of irradiance (Super LEDs) should be replaced whenever they reach a loss of 25% of the total irradiance for Bilirubin-Ebi.

The radiation module should be replaced by a FANEM Authorized Service. When this module is replaced, the technician must restart the count of the total time of use.

Measurement of irradiance:

Connect the optical probe (optional item, REF: 620,058,600) in the corresponding socket on the right side of the Bilitron BED™ box.

Press MENU up to "Measure the irradiance?". Press the ENTER key to enter this function. The first line will show the value of irradiance measured by the optical probe in $\mu\text{W}/\text{cm}^2\text{nm}$. The second line displays the radiation level in percentage, in which the radiation module is set. While the value of irradiance is monitored by the optical probe, the level of radiation might be changed through the UP or DOWN keys.

If you purchased the BILITRON BED™ with the optical probe, the radiometer will bear the calibration, already, carried out by the factory

Caution: For a reliable measure of irradiance with the BILITRON BED™, radiometer, employ the same optical probe used to calibrate the radiometer, always. Otherwise, the value of irradiance measured might not be correct.

If the BILITRON BED™ optical prober was purchased later, or if it was replaced, request an authorized FANEM Technician to carry out the radiometer calibration procedure, so the reading of the irradiance is reliable.

USER MANUAL
REVERSE AND REFLECTED MICROPROCESSED PHOTOTHERAPY OF SUPERLEDS
BILITRON BED™ - MODEL 4006
THE SUPERLED PHOTOTHERAPY

Language:


To select the language in which the information is displayed on the display, press MENU up to "Settings - Press ENTER. Press the ENTER key. The "Choose Language" will appear with the current language displayed on the second line of the display. Press the UP or DOWN keys to switch between the possible languages and the ENTER key to select the desired language. The other screens are displayed in the selected language.

Setting the clock / calendar:

Press MENU up to "Settings - Press ENTER. Press ENTER and then press the MENU key up to "Set Clock." Press the ENTER key to enter the adjustment mode.

The digit with flashing cursor might be changed with the up or down keys. Press the ENTER key to move the flashing cursor to the next digit and make adjustments, and so on to set all values. Press ENTER when the flashing cursor is on the last digit, to confirm the settings, or MENU to exit the adjustment mode without saving the changes.

6.4 - Protection against overheating

The BILITRON BED™ 4006 has the following protection devices against overheating:

- Two internal safety thermostats interrupt the supply of electrical power to the module when the radiation temperature inside of the module exceeds the maximum safety value. In this case, the lights go off and come back on, only, when the temperature of the radiation module is below the maximum safety value.
- There is one interruption detector for each micro-fan. In the event of obstruction of rotation, or failure of any one of micro-fans, the lamps, automatically, turn off to avoid overheating and the following message will appear on the display "Fan 1 interrupted " or "Fan 2 interrupted." The lamps light are turned on, again, when the faulty micro-fan is no longer obstructed.

USER MANUAL
REVERSE AND REFLECTED MICROPROCESSED PHOTOTHERAPY OF SUPERLEDS
BILITRON BED™ - MODEL 4006
THE SUPERLED PHOTOTHERAPY

7 - Preventive and Corrective Maintenance and Conservation

This section provides instructions for cleaning and maintenance. Where necessary, the instructions for the disassembly are provided. Maintenance that is not provided in this section should be performed only by the qualified technical service.

Routinely inspect the BILITRON BED™ Phototherapy and replace the assembly before putting the unit into use.

Caution: To initiate maintenance or cleaning procedures, make sure that the BILITRON BED™ is disconnected from the electrical net.

7.1 - Replacement of the Radiation Module

The average life cycle of the "SuperLED" of BILITRON BED™ is of, approximately, 10,000 hours, and we recommend the regular assessment of the irradiance levels, as per the regular maintenance procedures.

Caution: The lamps should be replaced when they reach a loss of 25% of its total irradiance for Bilirubin-Ebi.

Use, only, the original FANEM™ Radiation Module - REF: 406.101.600.
To replace the Radiation Module, contact the Fanem Service Center.

Important Note

Although similar, and bearing the same technical specifications, the Super LEDs specified for replacement must meet a series of technical requirements specified by the equipment design, as well as minimum irradiances, temperature ranges, etc., to obtain the complete safety and the clinical outcomes desired. So, do not use non-original components. The Super LEDs marketed by FANEM™ are 100% tested according to design requirements.

7.2 - Cleaning and Conservation

We recommend that upon the receipt of the initial equipment and whenever the BILITRON BED™ is not in use, or when a newborn is removed from the unit, or discharged, or, still in accordance with the Control and Hospital Infection Commission's procedures - CHIC. Perform the cleaning procedure followed the protocol of disinfection.

Caution: Do not use disinfectants that contain abrasives or sodium hypochlorite, as these might deteriorate the materials of phototherapy.

For better performances of the BILITRON BED™ over time, we recommend the carrying out of the cleaning of the radiation-emitting surface (external side) every 100 hours of machine use, using a dry cloth, or a light air jet, only to remove the accumulated dust. Never use inappropriate chemicals, detergents or soap products to clean, otherwise you may damage the components. Keep the Clear-gel Mattress, the acrylic crib and the reflector arch always clean. To clean, use a soft cloth with warm water and soap in an amount that suffices, but that does not allow the dripping of liquid into the equipment.

For a perfect cleaning, we recommend products based on quaternary ammonium that, as a germicidal product for fixed surface, it does not affect BILITRON BED™, while allowing the cleaning and the disinfection of the surface (strictly follow the instructions on the product label to ensure a perfect disinfection).

Note. Use a cloth soaked with sufficient, that does not allow the dripping of the product into the equipment.

USER MANUAL
REVERSE AND REFLECTED MICROPROCESSED PHOTOTHERAPY OF SUPERLEDS
BILITRON BED™ - MODEL 4006
THE SUPERLED PHOTOTHERAPY

7.3 - Batteries

This equipment has 2 -1.2 VDC batteries to operate the clock / calendar. These batteries have an average life cycle of 5 years, and should be replaced by an authorized technical service.

Caution: Do not use common or alkaline battery. Use only FANEM™ battery - REF: 000.103.026.

Caution: Disposal of **batteries** should be carried out according to the current legislation of the country.

7.4 - Protection Circuit Breakers

To replace the protection circuit breakers of the equipment, follow these steps:

- ◆ Remove the plug of the power cable from the power socket.
- ◆ Make sure that the equipment and its fixtures are not energized.
- ◆ With the aid of a screwdriver, press and rotate the front of the circuit breaker support to the left.
- ◆ The cover of the circuit breaker support will disconnect and the circuit breaker is attached to this cover.

Replace the defective circuit breaker with a new circuit breaker, according to its electrical specification, and reassemble the set.

Application	Fuse
100 - 240 V ~	1.2 A - Type "F" (2 cm)

7.5 - Replacement Parts

When you need replacement parts, refer to section 3 of this manual : Pieces, , Parts and Fixtures, with their respective references.

To obtain blueprints, parts, components or other references, please contact, directly, FANEM Ltda. or its sales representatives.

The function and the safety of this equipment is guaranteed only if the verification services, maintenance and repairs are carried out by the Fanem Technical Assistance or by individuals duly trained and qualified by FANEM Ltda.

Fanem Ltda. is not liable for any damage that might occur in equipment, as well as for consequences to the patient due to improper maintenance, not made by our Technical Assistance, or as a result of the use of non-original spare parts / fixtures.

Materials used in the design of the parts and the fixtures, and consumer and wear and tear items, are aimed to ensure the perfect operation of the equipment in accordance with its original features, as well as safety in respect to toxicity, flammability and biocompatibility of the materials used.

Caution: Use only FANEM™ original parts.

USER MANUAL **REVERSE AND REFLECTED MICROPROCESSED PHOTOTHERAPY OF SUPERLEDS** **BILITRON BED™ - MODEL 4006** *THE SUPERLED PHOTOTHERAPY*

7.6 - Maintenance Table

INTERVENTION	TERM	PERFORMER
Replace Lamps	Constant monitoring with radiometer. Replace with irradiance lower than 25% Ebi	Technical
Cleaning	Every 90 days (minimum).	User / Technical
Cleaning of Micro-fans	Every 90 days (minimum).	Technical
Evaluation / Routine Monitoring	6 months	Technical

7.7 - Disposal

If there is a need to dispose of the equipment or parts, and these do not have a specific destination defined by the customer, the item at hand should be sent to the manufacturer or to its representative for the due disposal arrangements according to the national legislation.

Caution: Disposal of **batteries** should be carried out according to the current legislation of the country.

Note. The equipment and / or parts thereof should be sent in clean and aseptic conditions.
The fulfillment of these conditions does not exempt the supplier of responsibility for possible impacts to the environment and / or persons.

USER MANUAL
REVERSE AND REFLECTED MICROPROCESSED PHOTOTHERAPY OF SUPERLEDS
BILITRON BED™ - MODEL 4006
THE SUPERLED PHOTOTHERAPY


7.8 - Diagnosis of Failures

Failure	Possible Causes	Correction
All LEDs are off and there is no specific message on the display.	a) the level of radiation at 0% b) burning of circuit breaker c) Poor connection of power cable d) Acting up of the safety thermostat due to overheating.	a) Adjust the level of radiation b) Replace the circuit breaker (see item 7.4) c) Connect to the network correctly. d) Ensure that the holes of the airflow, on the back of Bilitron Bed™ are not obstructed. If the problem persists, contact the Authorized Technical Service
All LEDs are off and display message: "Fan 1 interrupted" or "Fan 2 interrupted"	Failure in one of the micro-fans	Check if at the bottom of Bilitron Bed™, there is an object obstructing the rotation of the micro-fans. If the problem persists, contact the Authorized Technical Service to repair the defective micro-fans
One or more LEDs are off.	Burning of LEDs	Contact the Authorized Technical Service to replace the LEDs
The irradiance is very low	a) gel mattress or acrylic crib is dirty and / or dusty b) radiating surface is dirty and / or dusty c) Mains Power below 100 V ~. d) The lamps have already exceeded their life cycle e) Reading Error of the radiometer f) Adjust the low level of radiation	a) Clean the mattress and / or acrylic crib, according to item 7 of the manual. b) Clean the surface of radiation, according to item 7 of the manual. c) Check that the Power has fallen to below 100 V ~. d) Replace the FANEM™ Radiation Module with a new one - REF: 406.101.600. e) Check the position of the Radiometer Optical Probe. Compare the reading with standard radiometer. f) Set a new value for light intensity
Instead of displaying the value of the patient's skin temperature sensor, horizontal bars are displayed	a) patient's skin temperature sensor is disconnected b) Patient's skin temperature sensor with broken threads or broken bolts.	a) Connect the patient's skin temperature sensor b) Contact the Authorized Technical Assistance
Value of the radiometer reading does not change, even when placing it in different parts of the Mattress.	a) optical probe is disconnected b) non-calibrated optical probe or broken threads	a) Connect the optical probe b) Contact the Authorized Technical Assistance

If you cannot solve the issue with these procedures, call an authorized FANEM™ technician or your local dealer.

8 - Warranty and Technical Assistance

- Like all the equipments under the FANEM™ brand, this also receives a (01) one-year Warranty against possible manufacturing defects of the parts manufacturing (see the attached warranty).
- For all types of maintenance, within or after the expiration of the warranty, always look for an authorized Fanem Technical Assistance, do not allow any third party without the adequate technical skills, to damage or change the original characteristics of the equipment.
- Always use FANEM™ original parts.
- Registration ANVISA N°. 10.224.620.068
- Responsible Technician I
Eng Orlando Rossi Filho
CREA/SP 98.435/D

Jose Osvaldo Caputto Flosi
 International Sales Manager
 FANEM LTDA

9 - References

- MAISELS MJ: Neonatal Jaundice in Neonatology Pathophysiology and Management of the Newborn Avery GS (ed). JB Lippincott Co 2nd ed, page 473,1981

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

CE
0434

TABLE OF CONTENTS

1 – Introduction.....	01
2 – Technical Characteristics.....	02
3 – Parts and Accessories.....	05
4 – Precautions, Restrictions and Warnings.....	07
5 – Equipment Assembly.....	14
6 – Equipment Operation.....	17
7 – Preventive and Corrective Maintenance and Conservation.....	25
8 – Warranty and Technical Assistance.....	29

USER MANUAL
BILITRON™ SKY MODEL 5006
THE SUPERLED PHOTOTHERAPY

1 - Introduction

Phototherapy is commonly used in the treatment of neonatal hyperbilirubinemia, subjecting the patient to radiation exposure centered on the blue range of the visible light spectrum for a certain amount of time, defined by the medical protocol in service.

The introduction of phototherapy in the treatment of neonatal hyperbilirubinemia took place in the late 1950s and, later, presented a major evolution in the types of lamps applied, increasingly minimizing the risks without, however, eliminating the adverse side effects of ultraviolet and infrared rays over the patient's body. This type of light treatment began with low intensity fluorescent lamps, and, later, fluorescent lamps with emission in the blue band of the spectrum were introduced. Subsequently, halogen lamps were introduced and a variation with optical fibers. The LEDs (light emitting diodes) have also been introduced in phototherapy; however, due to low emission of light from conventional LEDs, hundreds of them are required to obtain a satisfactory result.

The emergence of SuperLED with high intensity radiation in the blue band of the spectrum, with differentiated physical and chemical composition from conventional LEDs, enabled the construction of equipments with a few LEDs, making them more compact and lighter. SuperLEDs emit safe radiation, completely eliminating the need for ultraviolet and infrared filters completely (see chart of emissions), with the advantage of dramatically decreasing the risk of burns, rashes and insensible loss of water dramatically.

BILITRON™ SKY Model 5006 is part of the 7th phototherapy generation from FANEM™. It is, therefore, an innovative electronic phototherapy, micro-processed, capable of eliminating unwanted radiation in the ultraviolet spectrum as well as in the infrared, incorporating a larger elliptical area with less loss on the edges, providing a more anatomic application to patients' bodies and increasing the effectiveness of the treatment.

The micro-processed electronics introduced in this phototherapy enables features that have never been used before, such as: Radiometer with a built-in optical probe (*optional probe*); treatment time, hour counter, light intensity adjustment, calendar clock, data memory for reporting emission and communication with printers and computers.

User's Responsibility: This product shall only be used by properly qualified personnel in compliance with local laws.

In order to preserve its efficiency and considering the nature of the work that BILITRON™ SKY Model 5006 is intended to perform, it is necessary to adopt special care during handling and operation. This User Guide provides instructions on the installation, use, operator maintenance and failure diagnosis. FANEM™ will not be responsible if the user does not operate the unit in accordance with its instructions, follow the maintenance warnings and recommendations provided in this manual or make any repair with unauthorized components. Worn or damaged parts must be replaced before use. Calibration, verification and possible repairs of the equipment must be performed periodically and must be carried out only by trained personnel. The equipment should not be changed without prior formal approval of FANEM LTDA. Any additional information is accessible through your local distributor.

This manual must be read and understood perfectly, remaining readily accessible to every professional working with the Phototherapy BILITRON™ SKY Model 5006 unit. When not in use, the User Manual should be stored with the equipment. If there is any information you do not understand, please contact FANEM™ representatives for additional information.

USER MANUAL

BILITRON™ SKY MODEL 5006


THE SUPERLED PHOTOTHERAPY

2 - Technical Characteristics

2.1 Power and Supply Voltage

Bilitron Sky 5006	Voltage (Volts) V~	Frequency Hertz (Hz)	Power Watts (W)	Current I (A)	Ref. Fanem Code
With rubber feet	100 - 240	50/60	60	1	506.000.900
With Pedestal	100 - 240	50/60	60	1	506.001.900

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
[Comprimento de Onda] – Wavelength

Spectral irradiance of the Bilitron Sky 5006 central radiance source between 400nm and 500nm, measured at 30cm distance

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[Comprimento de Onda] - Wavelength

USER MANUAL **BILITRON™ SKY MODEL 5006** THE SUPERLED PHOTOTHERAPY

2.4 Dimensions and Weight

♦ **Dimensions**

Bilitron Sky 5006	Width (cm)	Depth (cm)	Height (cm)
With rubber feet	36.0	27.0	6.5
In pedestal	50.0	50.0	100.0 - 161.0
Optional Caster for Bilitron Sky 5006 In pedestal		3" with brake (3 units)	

♦ **Weight**

Bilitron Sky 5006	Net Weight (Kg)	Packed Weight (Kg)
With rubber feet	1.8	4.0
In pedestal	10.0	16.0

2.5 Symbolology

	Type B Equipment Insulation Class I (Power Supply) Insulation Class II (Equipment)		Common Equipment - With No Protection Against Water Penetration
	Caution: Hot Surface		Caution: Risk of Electric Shock
	Caution: Consult Manual		Grounding Required
	OFF (With no Electric Supply Voltage)		ON (With Electric Supply Voltage)
	Patient Eye Protection	Continuous Operation Equipment	

2.5.1 - Symbolology - Packaging

	Fragile		This Way Up
	Protect from Sun Exposure		Protect from Rain
	Limit Stack		Caution: Consult Following Documents

USER MANUAL
BILITRON™ SKY MODEL 5006
THE SUPERLED PHOTOTHERAPY

2.6 Specifications

Average Irradiance at 40 cm:

Nominal Ebi.....	30 to 45 $\mu\text{W}/\text{cm}^2.\text{nm}$ ($\pm 10\%$)
Effective surface.....	320 x 220 mm
Adjustment range of irradiation.....	10% to 100% of max irradiation
Variation of irradiation adjustment.....	10%

Radiometer:

Measure range	0 to 200 $\mu\text{W}/\text{cm}^2.\text{nm}$
Resolution	0.1 $\mu\text{W}/\text{cm}^2.\text{nm}$
Reading accuracy	$\pm 1 \mu\text{W}/\text{cm}^2.\text{nm}$

Noise..... < 60 dBA (ambient 45dBA)

Environmental Conditions:

Operation Temperature range	18°C to 28°C (ambient)
Operation humidity range	10% to 95% (non condensing)

Environmental Conditions for Shipping and Storage:

Ambient Temperature	5°C to 55°C
Relative Humidity	30% to 75% (non-condensing)

Note.: Specified conditions with properly packed equipment.

USER MANUAL

BILITRON™ SKY MODEL 5006

THE SUPERLED PHOTOTHERAPY

3 - Parts and Accessories

The Phototherapy BILITRON™ SKY Model 5006 is assembled in a high impact plastic box containing three modules with 5 SuperLEDs, totalizing 15 SuperLEDs acting in the blue range, and three additional white LEDs for auxiliary examination lighting. It also contains an alphanumeric display and keyboard for programming functions, three micro-fans for SuperLEDs modules, RS-232 output and a micro-processed controller.


There are two assembly options:

- ♦ BILITRON™ SKY 5006 with rubber feet: assembled with rubber feet for use over incubators.
- ♦ BILITRON™ SKY 5006 in pedestal: assembled on a base with 3 casters of 2" for easier transport, column with adjustable height rod to allow a better positioning in devices such as Heated Cribs, Hospital Cribs, Incubators, etc. Arm with adjustable articulation and BILITRON™ body with rotation adjustment to best fit the spotlight on the patient.

Presentations:


1 - Phototherapy BILITRON™ SKY 5006

With rubber feet
REF: 506.000.900

2 - Phototherapy BILITRON™ SKY 5006

In Pedestal
REF: 506.001.900

USER MANUAL
BILITRON™ SKY MODEL 5006
THE SUPERLED PHOTOTHERAPY

3.1 - Radiation Source Module - REF: 006.061.900

The set contains 5 "SuperLEDs" assembled on an aluminum heatsink and 5 lenses for light collimation. The BILITRON™ SKY Model 5006 consists of Radiation Source modules.

The average life expectancy of each Radiation Source Module is about 20,000 hours; however, the manufacturer recommends constant monitoring of the Irradiance with FANEM Radiation Monitor, according to the usual maintenance procedures, since the emitted irradiance decreases according to the use of the radiation source.

Caution: Always use original FANEM™ radiation source modules. They are controlled and selected individually according to their original project characteristics.

3.2 - Micro-fan Kit - REF: 006.251.025

The Fan Assembly, located at the top of the box. It acts consistently with the powered equipment, providing dissipation of heat generated by its respective radiation source module, thus, prolonging its durability.

Caution: Do not obstruct the air passage of the micro-fans, under the risk of burning the radiation module.

3.3 - Optical Sensor Probe - REF: 620.058.600

Extension probe, equipped with an optical sensor head, used to measure the irradiance emitted by the phototherapy BILITRON™ SKY Model 5006.

Caution: The Optical Sensor Probe is an optional item.

3.7 - Eye Shield for Neonatal Phototherapy

Developed in hypoallergenic material with adjustable Velcro® maintaining its original position. It contains an internal photo-opaque pad ergonomically designed to prevent the exposure of the newborn's eyes to light leaks from the phototherapy process.

Information to order:

REF: 006.103.023 - Large: Circumference of the head 33-38cm (Kit with 20 units)

REF: 006.102.023 - Medium: Circumference of the head 24-33cm (Kit with 20 units)

REF: 006.101.023 - Small: Circumference of the head 20-28cm (Kit with 20 units)

** Accessory for exclusive use with the equipment.*

Caution: - Medical product for single use
- Non sterilized

USER MANUAL
BILITRON™ SKY MODELO 5006
THE SUPERLED PHOTOTHERAPY

4 - Precautions, Restrictions and Warnings

Caution: This chapter of the User Manual contains extremely important information to ensure the safety and integrity of the patient, user and equipment. Read it CAREFULLY!

- ♦ The use of BILITRON™ at a distance less than 30 cm between the radiation source and the patient is not recommended.

Caution: Minimum distance of 30 cm between the radiation source and the patient.

Caution: Do not use BILITRON™ with the Raditaion Module in direct contact with the patient's body.

- ♦ Verify if the network to which the equipment will be connected has features to withstand the electrical voltage and power conditions of the equipment indicated on the label attached to the appliance.
- ♦ The power supply cord plug must be connected into a grounded outlet, permanently fixed on the wall, in accordance with the standards and laws in effect for low voltage electrical installations and electrical laws for Health Care Establishments.

Caution: Do not use extension cords or multiple outlets. If appropriate grounding is not available, do not use the equipment.

Caution: In order to prevent eye injuries, the NB will need an eye protection.

- ♦ BILITRON™ must be used according to the ideal standards of environmental comfort to the nursery, i.e., (23 to 27°C).
- ♦ Incorrect use of the phototherapy unit may submit the newborn to serious risks. This unit must be operated only by trained and qualified personnel under the supervision of qualified medical professionals that are familiar with the risks and benefits of its use.
- ♦ This unit may not be used if any of its functions is not working properly. Qualified technical service should be required in this case.
- ♦ It is not advisable to leave a patient without constant medical care under the equipment.
- ♦ The positioning of the eye shield must be checked periodically.
- ♦ Continuously monitor the patient's temperature using a clinical thermometer.
- ♦ It should be checked regularly if there is any irritation in the patient.
- ♦ Thermal treatment devices, (Heated Crib, Incubators, Heated Mattresses, Radiation Heaters, Thermal Blankets, etc.) coupled to BILITRON™ can elevate the patient's body temperature to dangerous levels. It is recommended that these devices operate in the mode controlled by the skin (servo controlled) or, alternatively, in manual mode, with reduced heating power of the heater accompanied by constant measurements of the patient's body temperature, monitored using a clinical thermometer.
- ♦ Always use original FANEM™ radiation sources. The use of different types of sources may affect the irradiance and temperature conditions specified for the equipment.
- ♦ The patient's Bilirubin levels must be monitored according to the medical protocol in service.
- ♦ We recommend using a FANEM Radiation Monitor to confirm the perfect radiation conditions of the sources and to allow adjustment of the desired irradiance level, ensuring a minimum radiation of 10 $\mu\text{W}/\text{cm}^2 \text{ nm}$ at 300 mm from the light source in the most intense point of radiation (center of the light spot).

USER MANUAL

BILITRON™ SKY MODEL 5006

THE SUPERLED PHOTOTHERAPY

Note: The minimum internationally acceptable value for an appropriate treatment is $4 \text{ uW/cm}^2 \text{ nm}$ (minimum value according to Mims (1)), given the fact that this device offers much higher values, and, therefore, more effectiveness.

Caution: Do not place BILITRON™ near heat sources when applied in the Heated Crib. Maintain a minimum distance of 40.0 cm from the irradiant source of heat.

Caution: When using BILITRON™ Phototherapy equipment along with the Heated Crib, ensure that the phototherapy equipment does not invade the crib's irradiated heat area.

- ◆ Do not place cloths or towels over the BILITRON™ body. This may cause a temperature rise and consequent deactivation of the irradiating source by its overheating protection.

Caution: Do not block the exit holes of air circulation on the upper side of BILITRON™, since it may cause a temperature rise and consequent deactivation of the irradiating source by its overheating protection.

Caution: Do not use the phototherapy equipment in the presence of flammable anesthetics, combustible gases or cleaning agents which may cause combustion.

Caution: Equipment not suitable for the use in presence of flammable anesthetic mixtures.

- ◆ Usually, the durability expectancy of the original FANEM sources is approximately 20,000 hours, subject to a number of constructive and operational parameters. However, the constant monitoring of irradiance with FANEM Radiation Monitor is recommended, in accordance with usual maintenance procedures, to better assess the actual efficiency of light sources. Thus, even though the requirements of the technical standard IEC recommend that sources of irradiance (lamps) must be replaced when upon reaching a loss of 25% of total irradiance for bilirubin-Ebi, due to their superior irradiance characteristics, the BILITRON™ SuperLEDs continue to emit irradiances at highly effective levels.
- ◆ Do not press the command keys with finger nails or sharp objects.
- ◆ The blue light emitted by phototherapy equipment can disrupt clinical observations because it can mask changes in skin color, such as cyanosis. To carry out clinical observations on the patient, the phototherapy unit must be disconnected.
- ◆ Use only original FANEM parts and accessories to ensure better performance and safety of the equipment.

4.1 – Physiological Effects

Phototherapy is commonly used in the treatment of neonatal hyperbilirubinemia, subjecting the patient to a radiation exposure concentrated in the blue range of the spectrum of visible light, for a time to be determined, as appropriate, by the attending doctor.

BILITRON™ Phototherapy allows rapid reduction in serum bilirubin levels, reducing the need for blood transfusion and reducing the duration of treatment.

- ◆ Although transient, side effects can be observed due to the use of phototherapy. Among them, increase in peripheral blood flow with consequent vasodilation, erythema and increase of insensible water loss and alteration of intestinal transit stand out.
- ◆ The patient's water balance can be changed with the use of phototherapy.
- ◆ Therefore, newborns undergoing BILITRON™, as well as using any other type of phototherapy, need proper treatment of water support and eye protection, in addition to routine care of nurses and doctors.

USER MANUAL
BILITRON™ SKY MODELO 5006
THE SUPERLED PHOTOTHERAPY

- ◆ Although infrequent, some newborns undergoing BILITRON™ for more than four days can show tanned skin color. These patients are usually premature newborns and/or those with serum levels of direct bilirubin superior than 2 mg/dl. With the withdrawal of BILITRON™, the tanned color of the skin gradually fades without sequel.
- ◆ In newborns, whose skin is very jaundiced, it is common, after a few hours of Bilirubin phototherapy, that the illuminated area presents a paler and chalky color. This is due to intense action of light on the bilirubinic pigment deposited in the skin.
- ◆ For newborns at risk of developing elevated serum levels of bilirubin and for those whose bilirubin concentration continues to rise despite the use of phototherapy, we recommend using two BILITRON™.
Double phototherapy allows rapid decrease in serum levels of bilirubin, reducing the need for blood transfusion and reducing the duration of treatment.
- ◆ The eyes of the "operator" may be harmed in cases of long stay in the area of the patient. In this case, unplug the phototherapy device during the procedures, if necessary.
- ◆ Patients near Phototherapy equipment may need protection such as screens or shields, eyeglasses for protection, etc., because the exposure to light can cause transient effects, e.g. nausea or moderate dizziness.
- ◆ During treatment, Bilirubin photoisomers may cause toxic effects.
- ◆ Bilirubin serum levels of the patient must be measured regularly.
- ◆ Liquid infusions and drugs in general should not be stored within the area of radiation.
- ◆ Do not use plates or sheets of reflection in the phototherapy equipment, because they can cause dangerous body temperatures to the patient and, together with the unit of irradiant heat, they can cause serious injuries to the patient.
- ◆ The irradiated light source from SuperLED has extremely **low** levels of infrared rays, which are responsible for heating the body of the newborn, therefore, it is recommended to periodically check the newborn temperature with an inspected clinical thermometer.

Environmental Conditions:

- ◆ When used outdoors and in unheated cribs, it is recommended to check the temperature of the newborn to prevent hypothermia.
- ◆ Environmental conditions such as, for instance, moving and air flow, may affect the patient's thermal balance.

4.1.1 – Suggestion to Improve Phototherapy Effectiveness

- ◆ In order to improve treatment effectiveness, the newborn **should be naked**, only wearing the eye protection. The larger the area of the skin exposed to light is, the better the result is.
- ◆ However, in case of choosing the use of diapers, it should be the **smallest possible**, only covering a small portion of the Newborn's perineum.
- ◆ Given the fact that large diapers obstruct the action of the light on the newborn's skin, they also decrease the phototherapy effectiveness to a great extent.
- ◆ Increased irradiance combined with the largest body surface will allow a better clinical efficacy. This will enable significantly the reduction of time for phototherapy treatment and, therefore, allow early hospital discharge and reduction on hospital costs.

USER MANUAL

BILITRON™ SKY MODEL 5006

THE SUPERLED PHOTOTHERAPY

4.2 – Electromagnetic Compatibility and Immunity

It is the ability of the equipment and/or system to function in an electromagnetic environment without introducing intolerable electromagnetic disturbances to anything in this environment and, on the other hand, to function without degradation in the presence of an electromagnetic disturbance.

This equipment was designed, tested and meets the following electromagnetic compatibility standards.

- | | |
|---------------------|------------------|
| ♦ EN 60601-1-2 | ♦ IEC 61000-4-3 |
| ♦ NBR IEC 60601-1-2 | ♦ IEC 61000-4-4 |
| ♦ CISPR11 | ♦ IEC 61000-4-5 |
| ♦ IEC 61000-3-2 | ♦ IEC 61000-4-6 |
| ♦ IEC 61000-3-3 | ♦ IEC 61000-4-8 |
| ♦ IEC 61000-4-2 | ♦ IEC 61000-4-11 |

The equipment is within the parameters established for RF Emissions; Immunity; Electrostatic Discharge; Irradiated Electromagnetic Fields of Radiofrequency; and Transients (Bursts and Voltage Outbreaks).

Caution: Mobile and portable communication equipment of RF may affect Medical Electrical Equipment.

Caution: The use of accessories, transducers, sensors and network cables that are not original may result in increased Emissions or in a decrease of the equipment Immunity.

Guidelines and manufacturer statement – electromagnetic emissions

BILITRON™ SKY 5006 Phototherapy is intended to be used in electromagnetic environment specified below.
It is recommended that the client or user of BILITRON™ SKY 5006 phototherapy ensures that it is used in such environment.

Emission testing	Compliance	Electromagnetic environment - guidelines
RF Emissions CISPR 11	Group 1	BILITRON™ SKY 5006 Phototherapy uses RF energy only for its internal functions. However, its RF emissions are very low and it is unlikely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class A	BILITRON™ SKY 5006 Phototherapy is suitable for use in all establishments other than domestic and can be used in residential establishments and those directly connected to the public distribution of low voltage electricity which supply buildings for domestic use, provided that the following warning is understood: Warning: This equipment is intended only for the use of health professionals. This equipment may cause radio interference or disrupt equipment operations nearby. You may need to adopt mitigation procedures, such as reorientation or relocation of BILITRON™ SKY Phototherapy or local shielding
Harmonics Emissions IEC 61000-3-2	Class A	
Emissions due to fluctuation of voltage/scintillation IEC 61000-3-3	In accordance	

USER MANUAL

BILITRON™ SKY MODELO 5006

THE SUPERLED PHOTOTHERAPY

Guidelines and manufacturer statement – Electromagnetic Immunity

BILITRON™ SKY 5006 Phototherapy is intended to be used in electromagnetic environment specified below.
It is recommended that the client or user of BILITRON™ SKY 5006 phototherapy ensures that it is used in such environment.

Immunity Testing	Level of Testing ABNT NBR IEC 60601	Compliance Level	Electromagnetic Environment – Guidelines
Conducted RF IEC 61000-4-6	3 V ms 150 KHz to 80 MHz	10 V	Portable and mobile RF communication equipment should not be used near any part of BILITRON™ SKY 5006 Phototherapy, including cables, with separation distance less than the recommended, calculated from the equation applicable to the transmitter frequency.
Conducted RF IEC 61000-4-3	3 V ms 80 MHz to 2.5 GHz	10 V/m	Recommended Separation Distance $d = 0.35 \cdot P^{1/2}$ $d = 0.35 \cdot P^{1/2}$ 80MHz to 800 MHz $d = 0.7 \cdot P^{1/2}$ 800MHz to 2.5 MHz where P is the maximum output power of the transmitter in watts (W), according to the transmitter manufacturer and d is the recommended separation distance in meters (m) It is recommended that the field intensity established by the RF transmitter, as determined by an electromagnetic inspection on site, ^a is less than the compliance level in each rate of frequency ^b Interference may occur around the equipment marked with the following symbol: 

NOTE 1: In 80 MHz and 800 MHz, it is applied a higher rate of frequency.

NOTE 2: These guidelines may not be applicable in all situations. Electromagnetic propagation is affected by the absorption and reflection of structures, objects and people.

^a Field strengths set by fixed transmitters, such as base radio stations, telephone (cellular/cordless) and land mobile radios, amateur radio, AM/FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. In order to assess the electromagnetic environment due to fixed RF transmitters, it is recommended an electromagnetic inspection of the place. If the measure of the field intensity at the place where BILITRON™ SKY 5006 Phototherapy is used exceeds the level of compliance used above, BILITRON™ SKY 5006 Phototherapy should be observed to determine if the operation is Normal. If an abnormal performance is observed, additional procedures may be needed, such as reorienting or relocating BILITRON™ SKY 5006 Phototherapy.

^b Above the rate of frequency of 150 kHz to 80 MHz, the field intensity should be less than 10 V/m.

USER MANUAL
BILITRON™ SKY MODEL 5006
 THE SUPERLED PHOTOTHERAPY

Distances of separation recommended between portable and mobile RF communication equipment and BILITRON™ SKY 5006 Phototherapy			
BILITRON™ SKY 5006 Phototherapy is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The user can help to prevent electromagnetic interference by maintaining a minimum distance between the RF communication equipment, portable and mobile (transmitters), and BILITRON™ SKY 5006 Phototherapy as recommended below, according to the maximum output power of communication equipment.			
Maximum output power of the transmitter W	Distance of separation according to the transmitter frequency m		
	150 kHz to 80 MHz $d = 0.35 P^{1/2}$	80 MHz to 800 MHz $d = 0.35 P^{1/2}$	800 MHz to 2,5 GHz $d = 0.7 P^{1/2}$
0.01	0.04	0.04	0.07
0.1	0.11	0.11	0.22
1	0.35	0.35	0.7
10	1.11	1.11	2.21
100	3.50	3.50	7.00
For transmitters with maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the transmitter frequency, where P is the maximum output power of the transmitter.			
NOTE 1: In 80 MHz and 800 MHz, it is applied a higher rate of frequency.			
NOTE 2: These guidelines may not be applicable in all situations. Electromagnetic propagation is affected by the absorption and reflection of structures, objects and people.			

USER MANUAL
BILITRON™ SKY MODELO 5006
 THE SUPERLED PHOTOTHERAPY

Electromagnetic Immunity			
BILITRON™ SKY 5006 Phototherapy is intended to be used in electromagnetic environment specified below. It is recommended that the client or user of BILITRON™ SKY 5006 phototherapy ensures that it is used in such environment.			
Immunity Testing	Level of Testing NBR IEC 60601	Compliance Level	Electromagnetic Environment - Guidelines
Electrostatic Discharge IEC 61000-4-2 (ESD)	± 6 kV by contact	± 6 kV	The floor should be made of wood, concrete or ceramic. If the floor is made of synthetic material, the relative humidity should be at least 30%.
	± 8 kV by air	± 8 kV	
Electrical Transients Rapid/Pulse (Burst) IEC 61000-4-4	± 2 kV in the supply lines	± 2 kV in the supply lines	The quality of power supply should be the standard of a typical commercial or hospital environment.
	± 1 kV in the input and output lines	± 1 kV in the input and output lines	
Outbreaks IEC 61000-4-5	± 1 kV line to line	± 1 kV line to line	The quality of power supply should be the standard of a typical commercial or hospital environment.
	± 2 kV line to ground	± 2 kV line to ground	
Voltage drops, short interruptions and voltage variation in the input supply lines IEC 61000-4-11	< 5% U_T (>95% voltage drop in U_T) in 0.5 cycle	< 5% U_T (>95% voltage drop in U_T) in 0.5 cycle	The quality of power supply should be the standard of a typical commercial or hospital environment. If the BILITRON™ SKY 5006 Phototherapy user requires continuing operation during Power interruption, it is recommended that BILITRON™ SKY 5006 Phototherapy is supplied by a uninterrupted power supply or battery.
	40% U_T (60% voltage drop in U_T) in 5 cycles	40% U_T (60% voltage drop in U_T) in 5 cycles	
	70% U_T (30% voltage drop in U_T) in 25 cycles	70% U_T (30% voltage drop in U_T) in 25 cycles	
	< 5% U_T (>95% voltage drop in U_T) in 5 seconds	< 5% U_T (>95% voltage drop in U_T) in 5 seconds	
Magnetic field in the frequency of supply (50/60Hz) IEC 61000-4-8	3 A/m	3 A/m	Magnetic fields at the frequency of the supply should be at characteristic levels of a typical commercial or hospital environment.
NOTE: U_T is the AC supply voltage before the testing level application			

Caution: Phototherapy equipment have the characteristic of being used on incubators, hospital cribs, heated cribs, etc., i.e., they are used next to other health products, therefore, it is recommended that this equipment are observed in order to verify its normal operation according to the configuration used.

USER MANUAL


BILITRON™ SKY MODEL 5006

THE SUPERLED PHOTOTHERAPY


5 – Equipment Assembly

FANEM BILITRON™ SKY Model 5006 Phototherapy comes in two assembly configurations, as mentioned earlier.

BILITRON™ SKY Model 5006 with rubber feet The equipment comes ready for installation in Incubators. For BILITRON™ SKY Model 5006 configuration in pedestal, for use in Heated Cribs, Hospital Cribs, etc., see assembly instructions below.

**BILITRON™ SKY 5006
in pedestal**

**BILITRON™ SKY 5006
with rubber feet**

Caution: Before assembling or using BILITRON™ phototherapy, read carefully the entire User Manual.


Before using BILITRON™ phototherapy, the cleaning procedure should be performed, followed by an initial disinfection of the equipment as instructed in chapter 7 of this manual, and/or according to the protocols adopted by the Commission of Control and Hospital Infection of the service.

This appliance should only be installed, handled and maintained by qualified professionals.

USER MANUAL
BILITRON™ SKY MODELO 5006
THE SUPERLED PHOTOTHERAPY

5.1 – Pedestal Assembly - Bilitron® Sky 5006 Configuration – Pedestal

Step 1: Assembly of the Support Base


Attach the side brackets "A" and "B" in the center support rivets

Note 1: For attachment of the parts, 6 flat washers REF: 000.922.001, 6 spring washers REF: 000.939.001 and 6 Allen M6 screws REF: 000.142.015 must be used.

Note 2: Use an Allen M6 wrench for screws attachment (Allen wrench not included).

Step 2: Column Assembly


Attach the column in the center arm of the support base.

Note 1: For attachment of the parts, 4 flat washers REF: 000.922.001, 4 spring washers REF: 000.939.001 and 4 Allen M6 screws REF: 000.142.015 must be used.

Note 2: Use an Allen M6 wrench for screws attachment (Allen wrench not included).


USER MANUAL
BILITRON™ SKY MODEL 5006
THE SUPERLED PHOTOTHERAPY

Step 3: Articulated Arm Assembly

1 – Attach the horizontal rod in the attachment elbow of the vertical rod.

2 – Attach the horizontal arm at the attachment elbow using 1 Phillips M4 flat head screw REF.: 000.423.015 provided with the equipment.

3 – Attach the articulator in the support on the back of Bilitron™ Sky and connect the parts through the attachment screw-nut of the rod.

USER MANUAL
BILITRON™ SKY MODELO 5006
THE SUPERLED PHOTOTHERAPY

6 – Equipment Operation

Caution: This chapter of the User Manual contains extremely important information to ensure the safety and integrity of the patient, user and equipment. Read it CAREFULLY!

- ♦ Check the power cord, the radiation source and display for apparent signs of damage. If damaged, DO NOT place the unit in use and contact the nearest Authorized Technical Service from FANEM™.
- ♦ Make sure the fins of air circulation of the equipment are not covered or obstructed.
- ♦ The plug of the power supply cord must be connected into a grounded outlet, permanently fixed on the wall, in accordance with the standards and laws in effect for low voltage electrical installations and electrical laws for Health Care Establishments.

Caution: Do not use extension cords or multiple outlets. If there is not a perfect grounding, do not use the equipment.

- ♦ Connect the power cord plug into the outlet, and turn the switch on the back of the equipment, check if the radiation source is lit.
- ♦ When using BILITRON™ SKY 5006 pedestal version, check that the light-generating unit is firmly locked and attached to the pedestal set.
- ♦ Place phototherapy in the Incubator or Crib, where the patient is, adjust the height and angle of incidence of light on the patient.

Caution: In order to prevent eye injuries, the NB will need an eye shield.

Caution: The device reaches its stable conditions of temperature and irradiance levels in less time than 15 minutes.

- ♦ The angle of light incidence can be adjusted by moving the body of the reflector in the desired direction, the cap can be rotated to best adjust the light spot in the pedestal version.

Caution: In the BILITRON™ SKY 5006 in pedestal configuration, only carry the equipment with the column in the lowest height position.

USER MANUAL

BILITRON™ SKY MODEL 5006

THE SUPERLED PHOTOTHERAPY

6.1 –Light spot - “How to adjust”

For the best use of radiation on the patient, it is recommended that the light spot is adjusted so that the radiation falls from the thorax to the root of the thighs, see Table I.

Table I



Distance from the reflector to the patient d (mm)	Radiation $\mu\text{W}/\text{cm}^2 \cdot \text{nm}$ (in the center of the light spot)		
Distance	Radiation ($\pm 10\%$)	A (mm)	B (mm)
300	45 – 65	300	160
400	30 – 45	320	220
500	20 – 30	400	280

Image 01

Phototherapy on patients weighing over 2500 g should be positioned about 40 cm of them, allowing a light spot that covers the entire torso and root of the thighs of the newborn. In premature newborns and very low birth weight, this height can eventually be smaller, as long as the incidence of the light spot is always observed, covering the entire thorax and root of the thighs. The shorter the distance from the light source to the patient, the greater the radiation intensity will be, and vice versa, see Table I and Image 01.

♦ Always keep the light spot on the torso and root of the newborn's thighs. When the baby is positioned in Trendelenburg or Proclive, move the phototherapy equipment in order to keep the spotlight on the torso and the root of the thighs.


Caution: A minimum distance of 300 mm from the radiation source to the patient is recommended.

Caution: Do not use BILITRON® with the Radiation Source Module in straight contact with the patient's body.

USER MANUAL **BILITRON™ SKY MODELO 5006** THE SUPERLED PHOTOTHERAPY


SPECTRAL IRRADIANCE X DISTANCE

6.2- Micro-processed Controller Unit

Alphanumeric Display: This display informs the user the following items:

- 1- Calendar clock with Day, month and year, hours, minutes and seconds.
- 2- Total time of radiation source use (hour counter).
- 3- Level of radiation in $\mu\text{W}/\text{cm}^2 \cdot \text{nm}$ – when the probe with optical sensor is used (optional).
- 4- Treatment time.
- 5- Level of proportional Power of the light source when adjusted by the operator.
- 6- Data sent to printer or computer.

USER MANUAL


BILITRON™ SKY MODEL 5006

THE SUPERLED PHOTOTHERAPY

6.3 – Program Mode

On the front panel of the head, there are keys to program the functions, among them:

- 1- MENU key
- 2- INCREASE key
- 3- DECREASE key
- 4- ENTER key

The MENU key allows alternating among the various functions of the system. The ENTER key is to select a function or confirm a change made. The INCREASE and DECREASE keys allow you to change the value of a parameter.

05 / 21 / 11
11:03:02

BILITRON SKY™ main screen displays the current date and time.


Usually, when the display is on a screen that is not the main one and no key is pressed, the system automatically returns to the main screen after a few seconds.

Map of the functions of the control unit:

The functions on the alphanumeric display are organized as follows:

- **Main screen:** indicates current date and time
- **Radiation setting:** allows setting the level of radiation (from 10 to 100%)
- **Sensor – Calibration:** allows the calibration of the optical probe (optional) and consequent irradiance measure
- **Reports >:** allows saving files from 00 to 08
- **Radiation <:** allows saving irradiance readings
- **Print File:** allows sending files to the printer or computer
- **System Configuration:**
 - **Delete File:** allows deleting data/files already saved
 - **Set Clock:** allows setting current date and time
 - **Language:** allows selecting among Portuguese, English or Spanish
 - **Record:** allows activating or deactivating a record
- **Therapy:** indicates the phototherapy treatment time in hours
- **Reset Therapy:** allows resetting the phototherapy treatment time
- **SuperLED:** indicates the total time of SuperLEDs use, in hours and minutes

Clock/Calendar Setting:

Press MENU until the function “**System Config.**” is shown.

Press INCREASE until the function “**Set Clock**” is shown.

USER MANUAL

BILITRON™ SKY MODELO 5006

THE SUPERLED PHOTOTHERAPY

Press ENTER to enter in the setting mode.


The digit with flashing cursor can be changed by the INCREASE key. Press ENTER to move the flashing cursor to the next digit and make settings, and so on until all values are set.

Press ENTER when the flashing cursor is on the last digit to confirm settings.

Radiation Level Setting:


The light intensity can be set from 0 to 100% of total power.

Press MENU until the function "Set Radiation" is shown.

Press INCREASE to increase the value of light intensity. The level of power is shown on the display "Set P = 0 to 100%".

Press DECREASE to decrease the value of light intensity. The level of power is shown on the display "Set P = 0 to 100%".

Caution: BILITRON™ SKY has an additional light source which consists of 3 white LEDs Light of Day. This additional source remains on continuously when the equipment is energized, and it is more noticeable when the radiation level is set to very low values.

This lighting with white LEDs should be used for routine tests and procedures for newborns.

USER MANUAL
BILITRON™ SKY MODEL 5006
THE SUPERLED PHOTOTHERAPY


Treatment Time:

Press MENU until the function **"Therapy"** is shown.

Therapy
00013 Hs

The Treatment Time, total of hours used, will be informed on the second line of the, in hours.


To reset the counter, press MENU until the function **"Reset Therapy"** is shown.

Press ENTER. The hour counter will reset the Treatment Time.

Total Time of SuperLEDs use:

Press MENU until the function **"SuperLed"** is shown.

SuperLed
00013 Hs

The Total Time of SuperLEDs use, overall total of hours used by the SuperLEDs modules, will be informed on the second line of the display, in hours.

Although the average durability expectancy of the Radiation Module is around 20,000 hours, it is recommended constant monitoring of the irradiance with FANEM Radiation Monitor as the usual maintenance procedures to better assess the real effectiveness of the light source.

Thus, despite the requirements of IEC technical standard recommend that the sources of irradiance (SuperLEDs) should be replaced when they reach a loss of 25% of total irradiance for bilirubin-Ebi, due to its superior characteristics of irradiance, BILITRON™ SuperLEDs continue sending irradiances at highly effective levels.

The Radiation Module must be replaced by Authorized Technical Service from FANEM™. When this module is replaced, the responsible technician should restart the count of the total time of use.

Caution: The "Reset" of the Total Time of Use of the SuperLEDs should be done every time the LED Modules were replaced.

The replacement of SuperLED Modules should be done by FANEM™ Technical Assistance.

USER MANUAL


BILITRON™ SKY MODELO 5006

THE SUPERLED PHOTOTHERAPY

Irradiance Measure:

BILITRON SKY™ 5006 Phototherapy can also operate as a radiometer. For such thing, a FANEM™ Optical Probe must be used. (**Optical Probe – Optional Item - REF: 620.058.600**).

SENSOR

Connect the Cord of the Optical Probe (**Optional Item - REF: 620.058.600**) in the plug on the back part of BILITRON SKY™ body.

Press MENU until the screen with the irradiance value is shown in $\mu\text{W}/\text{cm}^2\text{nm}$.

27.0
 $\mu\text{W}/\text{cm}^2\text{nm}$

In case BILITRON SKY™ has been acquired with an Optical Probe, the radiometer Will come calibrated from the factory, and when the sensor is disconnected from the phototherapy, the value 0.0 $\mu\text{W}/\text{cm}^2\text{nm}$ will be shown.


To store in the memory the measurement of irradiance shown on the display, first press ENTER. An alert will appear (* after the R).

Press INCREASE so that the value is stored. The date and time that the measurement was made will also be memorized. Each value stored with its date and time establishes a file. The memory can store up to 100 irradiation files.

To view the files, press MENU until the screen "**Relat. > Radiat. <**" is shown.

Press INCREASE or DECREASE to view them. The numbers of the file will be on the display in the first row and the value of irradiance in the second row.

Press ENTER to show date and time of the corresponding file.

If BILITRON SKY™ has been acquired without the optical sensor, instead of presenting the value of irradiance in $\mu\text{W}/\text{cm}^2\text{nm}$, the message "**Sensor Not Calibrated**" will be shown. When the optical sensor is acquired, or if it replaces the existing one, ask a FANEM™ Authorized Technical Service to perform the calibration procedure of the radiometer, so that the reading of irradiance is reliable and so, it can be used.


Caution: For a reliable measurement of the irradiance with a radiometer from BILITRON SKY™, always use the same optical probe which the radiometer was calibrated with. Otherwise, the measured value of irradiance is probably not correct.

USER MANUAL

BILITRON™ SKY MODEL 5006

THE SUPERLED PHOTOTHERAPY

Language Selection:

Press MENU until the function "**System Config.**" is shown.

Press INCREASE until the screen indicating "**Language**" is shown on the first line and the current language on the second line.

Press ENTER to alternate among the possible languages.

Once the desired language is shown, press MENU to exit.
The other screens will be shown in the selected language.


Data Transmission:

BILITRON SKY™ is able to generate a report of information containing the total time of hours, the time of treatment, the current date and time and the stored irradiance files. This report is transmitted through the serial port (**RS 232**) of the device, and can be sent to a printer or computer containing a serial interface.

Connect BILITRON SKY™ to the printer or computer via a FANEM™ serial cable (supplied on request).

In the computer, a TELNET standard software may be used (such as Windows HyperTerminal, for example). The serial interface configuration of the equipment should be:

Bits per second: 9600 bps
Bits of data: 8
Parity: None
Bits of stop: 1
Flow control: Hardware

Press MENU until the function "**Print File**" is shown.

Press ENTER so that the report is transmitted.

BILITRON SKY™ 5006 has internal thermostats of safety which interrupt the supply of electrical power to the Radiation Module whether the inner temperature of the module exceed a maximum value of safety. In this case, the light source goes out and comes back on only when the temperature of the Radiation Module is below the maximum value of safety.

This section provides instructions for cleaning and maintenance. Where necessary, the disassembly instructions are provided. Maintenance that is not provided in this section should be performed only by qualified technical service.

Routinely, inspect BILITRON™ Phototherapy and replace accessories before putting the unit into service.

Caution: In order to initiate maintenance or cleaning procedures, make sure that BILITRON™ is unplugged.

The expected durability of "SuperLED" is around 20,000 hours, and it is recommended periodical assessment of the irradiance levels emitted according to the usual procedures of maintenance.

Caution: Despite the requirements of the technical standard IEC recommend that sources of irradiance (SuperLEDs) should be replaced when they reach a loss of 25% of total irradiance for bilirubin-Ebi, due to its superior characteristics of irradiance, the BILITRON™ SuperLEDs continue sending irradiances at highly effective levels.


USER MANUAL

BILITRON™ SKY MODEL 5006

THE SUPERLED PHOTOTHERAPY

Only use FANEM™ original Radiation Source Module. **Radiation Source Module - REF: 006.053.900**

How to replace the Radiation Source Module: With the assistance of a Philips type wrench, release the screws of the cover on the bottom. Open the access cover of the respective module and remove the Power Supply Module. Disconnect the cable from the socket and replace it with a new module.

Caution: Never change the Radiation Source Module with the device plugged on.

Important Observation

Although they are similar and have the same technical specifications, the SuperLEDs specified for replacement must meet a series of technical requirements specified by the equipment project, as well as minimum irradiations, temperature ranges, etc., for the achievement of total safety and clinical outcomes which are desired. Therefore, do not use non-original parts. 100% of SuperLEDs marketed by FANEM™ are tested according to the project requirements.

7.2 – Cleaning and Conservation

It is recommended by the initial receiving of the equipment and always that BILITRON™ is not being used for a long time or even according to procedures of the Commission of Control of Hospital Infection - HICC. Perform the cleaning procedure followed by the disinfection protocol.

Caution: Do not use cleaners or disinfectants containing abrasives or sodium hypochlorite, since these products can deteriorate the materials of phototherapy.

For the best performance of BILITRON™ over time, it is recommended that every 500 hours of machine use, a cleaning in the diffusers should be done (outside). Use a dry cloth or a gentle stream of air, only to remove accumulated dust. Never use inappropriate chemicals, detergents or scouring to do the cleaning, given the risk of damaging components.

For a perfect cleaning, we recommend Quaternary Ammonium which, as a germicidal product for fixed surface, does not affect BILITRON™ materials, allowing, simultaneously, superficial cleaning, and the fight against bacteria (strictly follow the instructions on the product label to ensure a perfect asepsis).

Note: Use a cloth soaked with the product enough so it does not drip liquid into the equipment.

USER MANUAL

BILITRON™ SKY MODELO 5006

THE SUPERLED PHOTOTHERAPY

7.3 - Batteries

This equipment has 2 batteries of 1.2 VDC for the clock maintenance. These batteries have average durability expectation of 5 years. They should be replaced by authorized technical service.

Caution: Do not use common or alkaline batteries. Only use FANEM™ battery - REF: 000.103.026.

Caution: The disposal of batteries must be carried out according to current legislation in the country.

7.4 - Switching Power Supply

This equipment works with 12 V $\pm$ and, for such, uses a switching Power supply with the following specifications:

Input: 100 – 240V~ (50/60Hz)

Output: 13.0 V $\pm$ / 3.7A max.

Caution: Only use FANEM® switching power supply – REF: 000.161.900.

Do not use other switching Power supplies which may have different electric characteristics and may cause the burning of phototherapy.

7.5 - Spare Parts

For possible spare parts, refer to section 3 of this manual: Parts and Accessories, with its respective references.

To obtain devices, parts, components or other additional references, contact FANEM LTDA. directly or its sales representatives.

The function and safety of this equipment are only guaranteed if the verification, maintenance and repair services are performed by the Technical Assistance from FANEM™ or by properly trained and qualified personnel by FANEM LTDA.

FANEM LTDA assumes no liability for any damage that may occur in equipment, as well as consequences to the patient due to improper maintenance, not made by our Technical Assistance, or when spare parts/accessories used are not original from the factory.

The materials used in the design of parts and accessories and consumer and abrasion items aim to ensure the perfect operation of equipment according to their original characteristics, as well as safety in respect to toxicity, flammability and biocompatibility of materials used.

Caution: Use only original FANEM™ parts.

USER MANUAL


BILITRON™ SKY MODEL 5006

THE SUPERLED PHOTOTHERAPY

7.5 – Chart of Maintenance

INTERVENTION	DEADLINE	EXECUTER
Source Change (SuperLED Module)	Constant verification with radiometer. Change with irradiance lower than 10 $\mu\text{W}/\text{cm}^2 \text{ nm}$ @ 300 mm of the source. (Perpendicular Incidence – at max adjustment of the power).	Verification: User / Technician Source Replacement: Technician
Cleaning	Every 90 days (minimum).	User / Technician
Micro-fan Cleaning	Every 90 days (minimum).	Technician
Routine Assessment/Gauging	6 months	Technician

7.6 - Disposal

Do not throw electronic devices or parts in the garbage when discarding the equipment or its parts.

To minimize pollution and ensure utmost protection of global environment, please recycle. For more information about "Waste Electrical and Electronic Equipment" visit related websites "Waste Electrical and Electronic Equipment – WEEE".

If there is a need to dispose the equipment, or parts, and they do not have a specific allocation set by the client, the specific item should be sent under the client's responsibility as of transport, manufacturer or its legal representative for appropriate action of disposal according to current national legislation.

Caution: The disposal of **batteries** must be carried out according to current legislation in the country.

Note: *The equipment and/or its parts must be sent on clean and aseptic conditions.*

Not fulfilling these conditions exempt the supplier of responsibility for possible impacts to the environment and/or people.

USER MANUAL
BILITRON™ SKY MODELO 5006
 THE SUPERLED PHOTOTHERAPY

7.7 Failure Diagnosis

Failure	Possible Causes	Correction
All LEDs are off and there is no specific message on the display.	a) Level of radiation at 0%; b) Bad connection of the Power supply cord; c) Performance of the safety thermostat due to overheating.	a) Adjust the radiation level; b) Connect to the mains power properly; c) Make sure that the holes for air circulation on the upper side of BILITRON™ are not blocked. If the problem persists, contact Authorized Technical Service.
One or more LEDs are off.	Burning of LED.	Contact Authorized Technical Service for LED changing.
Irradiance is very low.	a) BILITRON™ too distant from the patient; b) Dirty and/or dusty collimator; c) Very low mains power voltage; d) Module has already achieved limit; e) Error in the Radiometer reading; f) Intensity adjustment lower than 100%.	a) Check right height according to the manual; b) Clean Collimator, according to item 7 of the manual; c) Check network voltage; d) Change FANEM™ Power Supply Module for a new one - REF: 006.053.900; e) Check sensor placement of the Radiometer and its battery; f) Adjust a new value for the light intensity.
BILITRON™ does not stop at the right height of work.	a) Joints of the articulated arm with abrasion.	a) Check if the screws of the joints are tight; b) Change the arm.
Radiometer reading value does not change, even if positioning it at different places on the Mattress.	a) Disconnected optical probe; b) Uncalibrated optical probe or with broken wires.	a) Connect the Optical Probe; b) Contact Authorized Technical Service.


If it is not possible to solve the problem with these attitudes, call FANEM authorized technical service or its nearest local representative.

8- Warranty and Technical Assistance

- ♦ As all equipments from FANEM™, this also offers total guarantee of one (01) year against any possible defects in manufacturing parts (see warranty statement enclosed).
- ♦ For all types of maintenance, within or beyond the warranty, always look for Authorized Technical Service from FANEM™. Do not allow any third parties with improper technical qualification damage or change the original characteristics of your equipment.
- ♦ The SuperLED modules are guaranteed upon delivery.
- ♦ Always use original FANEM™ parts.
- ♦ ANVISA registration number 10.224.620.076
- ♦ Technician in Charge
Eng. Orlando Rossi Filho
CREA/SP 98.435/D

EC REP

Cinterqual LDA
 Travessa da Anunciada, 10 - 2 Esq-F
 2900-238
 Setúbal -- Portugal
 Phone: +351 265 238 237

Jose Oswaldo Caputto Flosi
 International Sales Manager
 FANEM LTDA

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