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Table of contents

	page
INTRODUCTION	3
1. PURPOSE	3
2. INDICATIONS	3
3. CONTRAINDICATIONS	3
4. SIDE EFFECTS	4
5. OPERATION CONDITIONS	4
6. PARTS OF THE MEDICAL DEVICE	6
7. TECHNICAL SPECIFICATIONS	11
8. DEVICE AND OPERATION	15
1. Principles of Operation	15
2. Application	15
3. Method of application	15
4. Device overview	17
5. Operation Plan for the TCM 300 BL – Morcellator for Surgical Procedures with Accessories	19
6. Before use	21
6.1 Connection to the power supply	21
6.2 Device preparation	21
6.2.1 Assembly of the morcellator set with protective tube (option 1)	22
6.2.2 Assembly of the Morcellator set with trocar cannula (option 2)	25
6.2.3 Assembly of the trocar cannula with myoma drill and mounting gripper (option 3)	28
6.3 Mains connection	29
7. Use	29
7.1 Using the control panel	30
7.2 Operation	31
7.2.1 Option 1: morcellation with protective tube	33
7.2.2 Option 2: morcellation with trocar cannula	34
7.2.3 Option 3: myomectomy	36
8. After use	37
9. SCOPE OF DELIVERY	39
10. CLEANING OF DEVICE COMPONENTS	43
General information about the disinfection, presterilization cleaning and	
Sterilization of MDs	44
General information in preparation MDs	44
Instrument preparation	46
10.1 Control unit and pedal	47
10.2 Preparation of the trocar cannula	47
10.3 Preparation of the electric motor	48
10.4 Preparation of the gear unit	50
10.5 Preparation of cutting tubes, protective tubes, obturators, grippers, seals and myoma drill	52

11.	MAINTENANCE AND REPAIR	53
11.1	General maintenance	53
11.2	Replacing the control unit fuses	54
11.3	Malfunction and troubleshooting	55
12.	SAFETY INFORMATION	56
12.1	EMC Manufacturer's Declaration of Conformity	56
12.2	Modification and misuse	56
12.3	Use precautions	56
12.4	List of operational hazards, risk analysis	57
12.5	Safety requirements, operation conditions and rules	57
12.6	Safety check	57
13.	MEDICAL DEVICE CLASSIFICATION	58
14.	COMPLIANCE WITH INTERNATIONAL STANDARDS	58
15.	PACKING AND MARKING	59
16.	TRANSPORTATION AND STORAGE CONDITIONS	70
17.	SERVICE LIFE	71
18.	INFORMATION ON DISPOSAL	71
19.	WARRANTY COVERAGE	71
	Appendix 1. Parts of the medical device	72
	Appendix 2. Table of spare parts and accessories to the medical device	74
	Appendix 3. Certificate of Analysis	87
	Appendix 4. EMC Report	88

1 General Information

INTRODUCTION

Nowadays newly developed medical devices make hysterectomy possible by a laparoscopic intervention. The advantage for the patient of the laparoscopic method is less pain and a much faster recovery. Patients appreciate much the minimal invasive character of this operation because it leaves almost no scars on the patient's skin, reduces the healing of the wounds and saves money and time for all involved parties. A morcellator being a special device for morcellation and further extraction of mass organs and tissues directly in the body cavity without additional incision of tissues to remove the organ is used for such procedures.

1. PURPOSE

The TCM 3000 BL - Morcellator for Surgical Procedures with Accessories is used for morcellation and extraction of tissue and for removing myomas of the uterus during laparoscopic procedures in gynecology.

2. INDICATIONS.

The TCM 3000 BL - Morcellator for Surgical Procedures with Accessories is used for surgical laparoscopic intervention during general surgical and gynecological procedures for morcellation and extraction of tissues in the following cases:

- myoma of the uterus;
- fibroma of the uterus;
- uterus neoplasms requiring partial hysterectomy;
- benign ovarian neoplasms.

3. CONTRAINDICATIONS

Relative or absolute contraindications may result from the general medical findings or in special cases in which the patient's risk for motor-driven tools is significantly increased. Ovaries, fallopian tubes, myomata and other structures must be devascularized and dissected before morcellation.

The laparoscopic deployment of morcellators is contraindicated at:

- Treatment of malignant tumors
- Treatment of vascularized tissue
- Preparation of tissue
- In gynecologic surgery in which tissue to be morcellated is known or suspected to contain malignancy
- For removal of uterine tissue containing fibroids in patients who are...
 - Peri or postmenopausal, or
 - Candidates for en bloc tissue removal, for example through the vagina or via a mini-laparotomy incision.

4. SIDE EFFECTS

Tissue morcellation with electric-powered surgical cutting devices is a well-established procedure. There have never been reported any unacceptable side effects provided the device is properly used as intended. Acceptable side effects include sensation of pain and anesthesia reactions.

5. OPERATION CONDITIONS

To ensure proper operation of the TCM 3000 BL - Morcellator for Surgical Procedures with Accessories one must meet the following conditions:

Ambient Conditions:

	Operation:
Relative humidity:	30% - 80%
Temperature:	10°C - 45°C
Atmospheric pressure:	800 hPa - 1060 hPa

	Before use read this manual and the manuals for all other equipment used in the procedure thoroughly.		The products described in this manual must be used by trained and qualified medical personnel only in accordance with the requirements and with reference to the manual.
	For protection of the patient, the user and third parties the products used must be prepared as specified in the manual.		Improper use or repair of the device, or failure to observe these instructions, will relieve us from any obligation arising from warranty provisions or other claims.
	The use of third-party products is the responsibility of the operator. Functionality and patient safety cannot be guaranteed with third-party accessories.		Prior to using the device, before startup, and before operation, the user must always ensure that the device and accessories are in good working order and are clean, sterile and operational.
	The device must be protected against direct sunrays, heat supplied by centrally batteries, excessive dust, vibration and mechanical shocks.		Repairs may only be performed by authorized NOUVAG service technicians.
	Uterine tissue may contain unsuspected cancer. The use of laparoscopic power morcellators during fibroid surgery may spread cancer, and decrease the long-term survival of patients. This information should be shared with patients when considering surgery with the use of these devices.		Use Nou-Clean spray (not included into the scope of delivery, ordered and delivered separately) to care for the electric motor and gear unit. Using other care products can result in malfunctions and/or cause the warranty to be revoked.

	The device is not sterile on delivery. All sterilizable parts must be sterilized before use (see Section Instrument Preparation).		Never operate the clamping mechanisms of the cutting tubes while the system is running. This could result in instrument damage.
	Operate the device only outside the danger zone of explosive, flammable mixtures or gases.		Keep hands clear of the cutting tubes when operating the device. Danger of serious injury.
	The electric motor ventilation slots must be kept clear in order to prevent the electric motor from overheating.		The TCM 3000 BL - Morcellator intended for Surgical Procedures with Accessories must be operated under the continuous supervision of medically qualified persons only.
	Ensure that the operating voltage setting corresponds to the local mains voltage.		Disregarding the guidelines on intermittent operation (1 min on, 1 min off for 8 cycles) may cause burns at contact with the gear unit.

Modification and misuse

-  Replacement of device parts with external ones or modification of the TCM 3000 BL - Morcellator and its accessories is prohibited. The manufacturer is not liable for any damages resulting from unauthorized modifications or manipulations. The warranty will be cancelled.
- Use of the TCM 3000 BL - Morcellator for Surgical Procedures with Accessories outside the indications described in Section **PURPOSE** is prohibited. The user or operator is solely responsible for any adverse consequences resulting from such use.

5

Requirements to the operator:

Education:

- a university-educated physician

Health:

- good physical and psychophysical condition

Language proficiency:


- English or
- the language of the Operation Manual.

Medical application:

The TCM 3000 BL - Morcellator for Surgical Procedures with Accessories is designed solely for use in special conditions of surgical emergency units in medical facilities by medically qualified persons. The medical device is stable during operation at the temperature of 10°C to 45°C and relative humidity of 30%

to 80%. The medical device is resistant to biological liquids and body tissues it comes in contact with during operation. The medical device is resistant to pre-sterilization cleaning, disinfection and sterilization. **The operator must wear medical gloves and a surgical mask while using this medical device.**

6. PARTS OF THE MEDICAL DEVICE

Picture 1. The TCM 3000 BL - Morcellator for Surgical Procedures with Accessories

Picture 1 demonstrates general view of the TCM 3000 BL - Morcellator for Surgical Procedures with Accessories. This medical device consists of a compact TCM 3000 BL - Morcellator control unit (1) with a simple and intuitive panel to control the speed of the attached instruments, a pedal with precise speed control (3), a gear unit (4), an ergonomic handle, attached to the gear unit (5) and an E-Motor 21 electric motor (6).

Picture 2 shows how to assemble the gear unit. Parts of the gear unit are shown and described in Picture 3.

Picture 2. How one can assemble the gear unit.

5136 Seal holder

5167 Cross seal for all cutting tube diameters,
10 units per a package

5166 Membrane seal for all cutting tube diameters,
10 units per a package

Seal unit with membrane screw

Seal unit

5163 Gear unit with seal unit incl. seals, seal holder and locking wrench (to be assembled as shown in picture 2). The electric motor speed is transmitted through the gear unit to the cutting tube at the appropriate ratio.


5183 The ergonomic handle that is slipped onto the gear unit protects the operator from excess heat of the electric motor and the gear unit.

2090 E-Motor 21 electric motor.
Autoclavable brushless electric motor with a 3 m long cable delivers adequate torque, while running smooth and consistent.

Picture 3. Gear unit with seals.

Picture 4 demonstrates how to assemble the Morcellator before its operation (the details of the assembling procedure are specified in Section 6 of Chapter 8. Device and Operation).

Picture 4. How one can assemble the Morcellator.

To suit the needs of every procedure the TCM 3000 BL - Morcellator for Surgical Procedures with Accessories has instruments with three different diameters. They are available in Ø12, Ø15 and Ø20 mm versions. By carefully composed and aligned instruments and accessories the system as a whole delivers excellent results at morcellation and guarantees for maximum safety at operation. Parts and accessories of the medical device are described in Appendix 2.

To ensure convenient and safe operation Nouvag AG offers the following design solutions:

Tissue Positioner

To prevent tissue pieces from rotating at morcellation, the protective tubes are equipped with a so called tissue positioner. The tissue positioner acts also as a guiding element at morcellation.

Cutting tube (Blade)

The high-tech cutting tubes manufactured by Nouvag AG, fulfill all requirements of Laparoscopic Hysterectomy as well as Myomenukleation. A high class refinement of the blade of the cutting tube enables multiple use with consistent cutting effect, without having to sharpen the blade.

Intended use of the parts and accessories of the TCM 3000 BL - Morcellator for Surgical Procedures with Accessories is listed in Table 1.

Table 1. Intended use of the parts and accessories of the TCM 3000 BL - Morcellator for Surgical Procedures with Accessories.

No	Part	Intended use
1.	Protective tube, Ø 12 mm Ø 15 mm Ø 20 mm	To insert the cutting tube and protect tissues during the procedure
2.	Obturator, Ø 12 mm Ø 15 mm Ø 20 mm	To forward the cutting tube to the tissues to be removed
3.	Cutting tube, Ø 12 mm Ø 15 mm Ø 20 mm	To resect the tissues to be removed
4.	Trocar cannula, Ø 12 mm Ø 15 mm	To insert the cutting tube and protect tissues during the procedure
5.	Myoma drill, Ø 10 mm	To conduct myomectomy (myoma removal)
6.	Grasping forceps, Claw gripper, Ø 10 mm	To grasp the resected tissue for its further removal
7.	Grasping forceps, Tenaculum gripper, Ø 10 mm Ø 5 mm Ø 13,5/15,5 mm	To remove small tissue pieces left

The control unit and the pedal don't come into contact with the patient, **the operator must wear medical gloves and a surgical mask while operating this medical device** (not included in the scope of delivery); the patient contacts only surgical instruments and indirectly seal units, fixing those instruments in the gear unit as well as Trocar sealing tubes. Materials these products are manufactured of are listed in Tables in Appendix 2. Mechanical performance of the materials is described in Table 2.

Table 2. Mechanical performance of the materials.

Material	Strength limit	Yield limit	Elongation	Hardness	Roughness	Density g/cm ³	Water absorption at 23°C, %	Permissible air operation temperature, °C
Stainless steel AISI 304	≥515 MPa	≥205 MPa	≥40 %	≤201 HB	Ra 0.8			
Stainless steel 1.4305	≥505 MPa	≥215 MPa	≥32%	≤240 HB	Ra 0.8			
Silicone 70 Shore A	≥12.0 MPa		32%	50+/-5 HRC				
Polytetrafluorethylene PTFE Turcon MF6	1.58 MPa	5 MPa	≥50 %	60 HRC		2.18		
Polyester-amide PEI MT blue		105 MPa		114 HRC	Ra 0.8	1.27	1.35	-50...170
Polyetheretherketone (PEEK)		< 110 MPa		105 HRC		1.31	0.45	-50...250

3. TECHNICAL SPECIFICATIONS

Main technical specifications of the TCM 3000 BL – Morcellator for Surgical Procedures with Accessories are listed below.

Picture 5. The TCM 3000 BL - Morcellator Control Unit

Table 3. Technical data: TCM 3000 BL Morcellator Control Unit

Voltage, variable	100 V / 115 V / 230 V, 50 – 60 Hz
Fuse	2 fuses, T 1 AL 250 V AC
Power consumption	110 VA
Maximum motor speed	40 000 rpm
Maximum torque, motor	6 Nm
Maximum torque at transmission outlet	80 Nm*
Number of motors	1
Number of pedals	1
Motor coupling	Intra Coupling ISO3964
Motor cable length	3 m
Pedal cable length	3 m
Morcellator speed	50 – 1000 rpm
Number of speed modes	3 (low speed (100 – 400 rpm) /medium speed (300-700 rpm)/high speed (500 – 1000 rpm)) (adjusted with the pedal)
Protection class	II
Method of protection against electric shock:	BF type
Maximum allowed time of operation mode selection	1 minute
Unit category	stationary
Lighting,	6 ultra-bright white leds
Video output	USB only
User interface	LCD display. 2 parameter setting lines, size 6x1.5 cm
Electronic control	Microprocessor: supply voltage – 230/115/110 V, programming voltage - 5 V
Software	Embedded, version 1.63/11.07.08, safety class A
Time for the request processing	2 sec
Software security level	B
Dimension of the program	8 573 lines of source code
Process frequency	8 MHz
Pedal	IPX68
Dimensions (W x D x H):	No more than 120 x 180 x 107 mm
Weight:	No more than 1,8 kg
Length of the power cable	3 m

Picture 6. E-Motor 21 electric motor

Table 4 . Technical data: E-Motor 21 electric motor

Rotation speed	2'000 to 40'000 rpm
Motor type	BLDC motor with Hall sensors
Weight, no more than	280 g without cable, 300 g with cable
Dimensions, no more than	Length: 8 cm, diameter: 2 cm
Output capacity	120 W
Connection	E-Type (ISO 3964)
Noise power level	50 dB
Cable length	3 m
Sterilization	Steam sterilization at max. temperature 135°C

13

Picture 7. Vario Pedal.

Table 5. Technical data: Vario Pedal

Type:	Vario type foot pedal
Number of speed modes	3 (low speed (100 – 400 rpm) /medium speed (300-700 rpm)/high speed (500 – 1000 rpm))
Protection class	IPX8 as per IEC 60529
Connection	E-Type (ISO 3964)
Weight, no more than	1 kg
Dimensions (W x D x H), : no more than	28 x 9 x 4.5 cm
Pedal pressure force	At least 10 N and up to 50 N maximum
Cable length	3m

Picture 8. Gear unit

Table 6. Technical data: gear unit

Reduction ratio	40:1
Maximum torque	100 Nm
Maximum input rotation speed	40'000 rpm
Maximum rotation speed at cutting tube (blade) outlet	1'000 rpm
Weight, no more than	450 g
Dimensions, no more than	handle: length: 7 cm, diameter: 1.5 cm. Body: total length: 11 cm, diameter: 4 cm
Cable length	3m
Connection	E-Type (ISO 3964)
Sterilization	Steam sterilization at max. temperature 135°C

Picture 9 .Gripper

Table 7. Technical data: gripper.

Full opening/closing effort (moment of force)	0,2-0.68 N, with two fingers
Angle of free gripper mouth movement	Up to 25-35 °
Head projection above the surface	No more than 0.5 mm, caulked
Radius of working part root face	No more than 0.1 mm

4. DEVICE AND OPERATION

1. Principles of Operation

The TCM 3000 BL – Morcellator for Surgical Procedures with Accessories consists of an electric motor, which drives a cutting tube (blade). You can adjust speed and direction (to the right or to the left) of cutting tube rotation. Cutting tubes with blades at the distal end are inserted through an abdominal wall that has been enlarged by an Insufflator (user-supplied), by the means of a protective tube with an obturator or through a trocar cannula. The tissue that is to be removed is guided to the cutting tube, cut off and removed through the tube. The cutting tube is installed in the gear unit attached to the electric motor which is located in an ergonomic handle. The cable connects the electric motor with the control panel displaying various morcellator operation modes for a physician. The physician adjusts morcellator speed using the control panel or a foot switch.

The set of morcellation instruments includes cutting tubes of various diameters, protective tubes, trocar cannulas, laparoscopic instruments (gripper, a myoma drill).

2. Application

The TCM 3000 BL – Morcellator for Surgical Procedures with Accessories is applied in surgical gynecology.

3. Method of application

Surgical procedure of morcellation is conducted under local or general anesthesia. The patient body shall be placed in Trendelenburg position.

The cutting tube of the morcellator is inserted through an abdominal wall with an obturator placed into the protective tube or through a trocar cannula. By rotating rapidly (rotation speed is adjusted by the foot pedal) the cutting tube cuts off the tissue that is to be removed. The cut off tissue is grasped with the

claw gripper. The tip of the cutting tube is always located inside the trocar cannula or the protective tube in the lower abdominal wall to avoid deflection or wandering of the cutting tube inside the abdominal cavity. The cut off tissues fall onto the cutting edge of the cutting tube without sucking inside. Large segments of the cut off tissue (10 mm long and larger) are removed with the gripper through an outlet of the protective tube or trocar cannula and placed into a sterile tray.

After large tissue segments are removed with Tenaculum gripper small pieces left are to be removed. The patient shall be moved out of Trendelenburg position so that intra-abdominal liquid starts coming into the lower pelvis. There intra-abdominal liquid is aspirated and replaced with heparinized Ringer's normal saline solution.

Due to effective cutting off the structure of removed tissues is just slightly damaged ensuring high reliability of histological analysis.

This manual describes three different options for removing tissue by a laparoscopic procedure using the TCM 3000 BL – Morcellator for Surgical Procedures with Accessories.

Option 1: Tissue removal by morcellation with a protective tube

Instruments (cutting tubes, gripper) are inserted through a protective tube sealed by a valve system. The valve can be covered by a thumb for removing the gripper from the sealed protective tube to prevent gas loss.

Option 2: Tissue removal by morcellation with a trocar cannula


Instruments (cutting tubes, gripper) are inserted through the trocar cannula. The trocar cannula contains a valve that prevents gas loss when withdrawing the instruments.

Option 3: Tissue removal (myomectomy) with a myoma drill through the trocar cannula

The myoma drill is inserted through the trocar cannula. The trocar cannula can be used with or without a cutting tube. No gas is lost after removal of instruments with the trocar cannula.

4. Device overview

Front view


1. Display
2. Motor LED
3. Pedal connection
4. Motor connection

Control panel


1. Vario button (variable speed range)
2. Low button (low speed range)
3. Medium button (medium speed range)
4. High button (high speed range)
5. Change rotational direction button
6. Up button (increase value)
7. Down button (reduce value)

Rear view

1. Mains connection
2. Power switch
3. Fuse module
4. Window with national voltage display
5. Type plate with type designation, reference number, serial number, information on power supply and device fuse

Gear unit


1. Flange for protective tube connect
2. Thread for seal holder
3. Cutting tube holder
4. Seal unit with membrane screw
5. Thread for membrane screw
6. Membrane screw
7. Seal unit thread
8. Motor connection coupling

Ergonomic handle

1. Release lever
2. Handle

E-Motor 21

1. Control unit connection plug
2. Electric motor cable
3. Electric motor cap
4. Electric motor
5. Electric motor clutch
6. Electric motor connection plug

Picture 10. Device overview


Picture 6 demonstrates configuration of the accessories to the TCM 3000 BL – Morcellator for Surgical Procedures with Accessories.

Protective tube with Obturator


1. Protective tube
2. Lock nut
3. Nose for tissue positioning (tissue positioner)
4. Obturator

Cutting tube

1. Cutting edge
2. Locking pins


Trocar Cannula

1. Protective tube
2. Seal cap holder
3. Valve actuator
4. Valve
5. Main body
6. Nose for tissue positioning (tissue positioner)

18

Gripper

1. Mouth insert
2. Gripper mouth
3. Cannula shaft
4. Fixed section
5. Flexible leg
6. Locking device
7. Spring component

} Handle


Picture. 11. Configuration of medical device accessories

the cutting tube to the tissue that has to be removed, the protective tube has to be unlocked by holding the gear unit with one hand and performing a counter movement of the protective tube with the other hand.

Option 2. Morcellation with a trocar cannula


When morcellating with a trocar cannula (Picture 12, option 2) the abdominal wall is kept in safe distance from the operation area by inflating the abdominal cavity. This gives also safe clearance for the surgical operation. The trocar cannula stays flexible while operating, which means it can be rotated against the organ or tissue to be protected. Only then the obturator can be removed and the cutting tube can be inserted.

The configuration of instruments needed for this procedure is shown on the System chart (Pic .13)

Picture 13. Image-guided morcellation (endoscope).

Option 3

Picture 14. The instruments used for option 3.

Option 3. Tissue removal (myomectomy) with a myoma drill through a trocar cannula

The myoma drill is inserted through the trocar cannula. Large segments of the cut off tissue are grasped and removed with the gripper through an outlet of trocar cannula and placed into a sterile tray. The Ø10 mm claw gripper are inserted through the trocar cannula to grasp a myomatous nodule and remove it.

6. Before use

The TCM 3000 BL - Morcellator control unit and all required accessories and instruments must be placed on a level nonslip surface. For reliable operation of the device, the latter requires enough space to avoid undesirable interference so remove all the limitations and everything affecting manipulations. The display, keypad and indicator lights must be fully visible at all times. It is very important to ensure that no objects can fall on the pedal. The power plug at the rear of the device must be accessible at all times.

6.1 Connection to the power supply

- Before plugging the power cable into the power socket for the first time, you must check the supply voltage setting next to the power switch.
- In order to prevent the risk of electric shock, the device may only be connected to a power network with a PE protective earth conductor.

If the voltage shown does not correspond to the local mains voltage, the grey fuse holder must be set to the correct voltage:

1. Switch the device off.
2. Unplug the power cable.
3. Use a screwdriver to open the fuse slot.
4. Remove the fuse holder.
5. Remove the grey fuse holder and reinsert it so that the local mains voltage setting is shown in the small window.
6. Slide the grey fuse holder back in and close the fuse slot.
7. Check the mains voltage shown on the fuse slot.
8. Plug the power cable back into the device..

21

6.2 Device preparation

The device is not sterile on delivery. All sterilizable parts must be sterilized before use (see Section Instrument Preparation).

If the accessories have already been sterilized: when removing them from the sterile packaging, ensure that the sterile packaging is not damaged and that the sterility indicator confirms sterility and shelf life of the sterile item has not expired.

The Morcellator (as per option 1 to option 3) must be assembled under sterile conditions (wear gloves and a mask, place Morcellator and accessories on a sterile surface).

6.2.1 Assembly of the morcellator set with protective tube (option 1)

There is a risk of injury to tissue by the sharp edge at the distal end of the cutting tube.

- Handle the cutting tube carefully.
- If the cutting tube is not used, rotate the protective tube to the NO CUT position.
- If the cutting tube is not used, it must be placed on a flat, sterile surface so it cannot fall and cause injury


- Check the cutting tube.
- Make sure that the cutting edge of the blade is sharp and undamaged (e.g. no cracks, deformation or burrs).
- Use only cutting tubes in perfect condition.

1. Inserting the cutting tube through the gear unit

Remove the seal unit from the gear unit.


Insert the cutting tube through the gear unit. The two guide pins on the cutting tube must fit into the slot on the flange of the gear unit.

2. Fixing the cutting tube to the gear unit

Hold the cutting tube with one hand.

Connect the sealing unit to the gear unit with the bayonet closure with the other hand.

1. Mounting protective tube

Select a suitable protective tube.

Push the protective tube over the distal end of the cutting tube. The guide pin of the protective tube must be in the distal slot of the gear unit.


Lock the protective tube in position with a click.

2. Protecting the cutting tube

The distal end of the cutting tube is very sharp.

To prevent injury, the cutting tube must be covered at all times when it is not in use.


To cover the cutting tube, rotate the protective tube until the words "NO CUT" are visible in the distal slot of the gear unit.

The cutting tube is completely covered by the protective tube.

3. Replacing seal holder

Remove the seal holder from the gear unit by rotating it slightly.


Unscrew the membrane screw from the seal unit.

Check the seal holder before use. Do not use damaged seals (frayed, discolored, yellowed, porous, non-elastic).

Position the edge of the cross seal (1) to point to the inside stop (2) of the seal holder (3) and press the cross seal in.


Position the edge of the membrane seal (4) to point to the slot (5) of the membrane screw (6) and press the membrane seal into the slot.

Screw the membrane screw (6) into the thread of the seal holder (3) to the stop and connect the seal unit to the gear unit.

4. Inserting obturator

Insert the obturator through the seal holder, gear unit and protective tube to the stop.

5. Attaching ergonomic handle

Lift and hold the release lever (1) of the ergonomic handle (2). Push the handle (2) over the motor connection of the gear unit. Release the handle release lever. The handle release lever must lock in position.

Do not open the handle release lever to remove the ergonomic handle during the procedure.

6. Connecting electric motor

Insert the electric motor coupling into the bottom end of the handle.

Continue to insert the electric motor until it locks into the gear unit's motor connection.

Align the red dot on the device connection plug of the electric motor with the red dot on the electric motor connection of the control unit.

Plug the electric motor device connection plug into the electric motor connection of the control unit. The plug must lock in position.

7. Connecting pedal

Align the red dot on the device connection plug of the pedal with the red dot on the pedal connection of the control unit.

Plug the pedal device connection plug into the pedal connection of the control unit.

The plug must lock in position.

Functional check

- Check Morcellator for loose parts
- Check that Morcellator is correctly assembled
- Are all labels easily readable?
- Check cutting tube: cutting tube evenly ground, no broken edges or bends?
- Switch on control unit
- Switch on electromotor (press foot switch).
- Run electric motor with Morcellator for about 20 seconds; it must reach the set speed.
- If it does not run evenly, switch off electric motor with Morcellator.
- Switch off control unit.

 The control unit will only operate if the foot switch has been connected to the control unit.

- Check that Morcellator is correctly assembled.
- Check that the gear unit runs smoothly; if necessary treat it with Nou-Clean care spray (shall be user-supplied, as it is not included into the scope of delivery).
- Check that the electric motor is clean and runs smoothly; if necessary clean it and treat with Nou-Clean care spray (shall be user-supplied, as it is not included into the scope of delivery).
- Repeat the functional check.

Risk of injury by electric motor failure

Problem: the electric motor becomes very hot.

Cause: the electric motor was damaged during preparation.

Remedy: do not continue to use the electric motor. **Always have a spare electric motor available.**

6.2.2 Assembly of the Morcellator set with trocar cannula (option 2)

Using the trocar cannula

- ✎ The trocar cannula is used instead of the protective tube. The integrated valve flap prevents the escape of gas when instruments are moved in and out.

The nose of the trocar cannula partially extends beyond the cutting tube to protect organs in the abdominal cavity from the sharp blade. This requires the Morcellator to be rotated under visual supervision (endoscopy) so the nose of the trocar cannula covers the cutting tube blade to protect the organ.

1. Visual check

- Check the trocar cannula for sealing, corrosion and sterility.
- Ensure that it is free from rust, dents and scratches.
- The distal end of the obturator must not be damaged.
- The sealing flap and the seal must not be damaged.


Risk of injury

Select blade. Make sure that the cutting edge of the blade is sharp and undamaged (e.g. no cracks or burrs).

2. Inserting the cutting tube through the gear unit


Remove the seal unit from the gear unit.
Insert the cutting tube through the gear unit. The two guide pins on the cutting tube must fit into the slot on the flange of the gear unit.

3. Fixing the cutting tube to the gear unit

Hold the cutting tube with one hand.


Connect the sealing unit to the gear unit with the bayonet closure with the other hand.

4. Mounting trocar cannula

Select a suitable trocar cannula. Open the trocar cannula valve by actuating the lever to prevent damage to the cutting tube by pushing the trocar cannula too far.

Push the trocar cannula over the distal end of the cutting tube. The trocar cannula must enter

and be in the distal slot of the gear unit.
Lock the trocar cannula in position with a click.


5. Replacing seal holder

Remove the seal holder from the gear unit by rotating it slightly.


Unscrew the membrane screw from the seal unit. Select a suitable seal holder.
Keep a cross seal and a membrane seal ready.

6. Position the edge of the cross seal (1) to point to the inside stop (2) of the seal holder (3) and press the cross seal in.

Position the edge of the membrane seal (4) to point to the slot (5) of the membrane screw (6) and press the membrane seal into the slot.

7. Screw the membrane screw (6) into the thread of the seal holder (3) to the stop and connect the seal unit to the gear unit.


8. Inserting obturator

Insert the obturator through the seal holder, gear unit and the trocar cannula.

9. Attaching handle

Lift and hold the release lever (1) of the ergonomic handle (2). Push the handle (2) over the motor connection of the gear unit. Release the handle release lever. The handle release lever must lock in position.


10. Connecting electric motor

Insert the electric motor coupling into the bottom end of the handle.

Continue to insert the electric motor until it locks into the gear unit's motor connection.

Align the red dot on the device connection plug of the electric motor with the red dot on the electric motor connection of the control unit.


Plug the electric motor device connection plug into the electric motor connection of the control unit. The plug must lock in position.

11. Connecting pedal

Align the red dot on the device connection plug of the pedal with the red dot on the pedal connection of the control unit.

Plug the pedal device connection plug into the pedal connection of the control unit. The plug must lock in position.

Functional check

- Check Morcellator for loose parts
- Check that Morcellator is correctly assembled
- Are all labels easily readable?
- Check cutting tube: cutting tube evenly ground, no broken edges or bends?
- Switch on control unit
- Switch on eletromotor (press foot switch).
- Run electric motor with Morcellator for about 20 seconds; it must reach the set speed.
- If it does not run evenly, switch off electric motor with Morcellator.
- Switch off control unit.

 The control unit will only operate if the foot switch has been connected to the control unit.

- Check that Morcellator is correctly assembled.
- Check that the gear unit runs smoothly; if necessary treat it with Nou-Clean care spray (shall be user-supplied, as it is not included into the scope of delivery).
- Check that the electric motor is clean and runs smoothly, if necessary clean it and treat with Nou-Clean care spray (shall be user-supplied, as it is not included into the scope of delivery).
- Repeat the functional check.

Risk of injury by electric motor failure

Problem: the electric motor becomes very hot.

Cause: the electric motor was damaged during preparation.

Remedy: do not continue to use the electric motor. **Always have a spare electric motor available.**

28

6.2.3 Assembly of the trocar cannula with myoma drill and mounting gripper (option 3)

Assembly of the trocar cannula with myoma drill is similar to that of the Morcellator set with the protective tube (option 1) and with the trocar cannula (option 2), except that myoma drill and gripper are used instead of cutting tubes.

6.3 Mains connection

Risk of electric shock

Connect the device to an earthed AC mains system that conforms to IEC guidelines only. Connect the device to a power supply that corresponds to the connection values on the type plate on the back of the device.

Risk of damage

Do not switch on the control unit and the other equipment until all cables have been connected. Otherwise the equipment may be damaged.

The device is not water-proof.

Keep the device away from liquids. Keep the mains connection socket dry. Do not use the device if liquid penetrates it.

Connecting the power cable

Switch off the power switch on the back of the device ("O" position) (Picture 10).

Plug the included power cable into the socket on the back of the device.

Plug the other end of the power cable into the power socket.

7. Use

Before use

Make sure that the product has been correctly prepared and inspected.

Risk of injury

The TCM 3000 BL - Morcellator control panel has been designed for operation solely with the instruments manufactured by Nuvag AG. The user is solely responsible if different products are used with the set. Incompatibility may risk injury to patients.

Switching on and off

Switch on the control panel. The display lights up. A short signal tone sounds. The control panel displays the value "200 rpm".

Operating mode (electric motor)

The motor is designed for intermittent operation (1 minute on, 1 minute off for 8 cycles).

7.1 Using the control panel

Setting the speed

The speed can be set manually in stages (within a range of 50 to 100 rpm).

Press ▲ to increase the speed.

Press ▼ to reduce the speed.

The set speed is shown on the display.

Starting electric motor

Press foot switch.

The blade rotates at the set speed.

The electric motor LED lights up.

To stop the electric motor remove foot from foot switch.

Risk of injury and damage due to heat!

The Morcellator may generate heat if operated for a long period.

Remedy:

Use Morcellator only in intermittent mode of "1 minute on, 1 minute off for 8 cycles".

Stop the Morcellator immediately if one or more components of the Morcellator set become very warm.

Precaution when changing parameters. Unusual behavior of instruments while operating may provoke false reactions and jeopardize the patient.

Every new setting and the new behavior of the instrument has to be verified and tested.

Available rotary modes

 (100 – 400 rpm)

 (300-700 rpm)

 (500 – 1000 rpm)


Activating mode setting

Press  to activate mode setting.

The Vario LED lights up.

The Low LED lights up (Default setting)

The display shows «VARIO LOW, 100-400 rpm».

Selecting mode

Press ,  or  to select the corresponding mode. The display shows the mode and the corresponding speed range that can be selected with the pedal. The corresponding LED lights up.

Starting electric motor

Press foot switch. The electric motor LED lights up. The blade rotates at a speed within the range set by how far down the foot switch is pressed. To stop the electric motor remove foot from foot switch.

Exit mode setting

Press .

Changing direction of rotation

To change the direction of rotation, press . The display shows the mode and the corresponding speed range that can be selected with the pedal. The corresponding LED lights up.

7.2 Operation

Risk of death or serious injury!

Unsupervised cutting may injure the abdominal wall and organs and cause the death of the patient.

Remedy:

- Before every procedure make sure that an endoscopic procedure is preferable to a conventional procedure.
- Use only the products described in this manual.
- During the procedure use the protective tube or trocar cannula as a precautionary measure. Follow the instructions for "Option 1: Morcellation with protective tube" and "Option 2: Morcellation with trocar cannula" in Section Device Preparation.
- Use Morcellator under visual control (endoscope) only.
- Use Morcellator only for morcellation of tissue that has been fully prepared and is clearly visible. Nouvag AG accepts no liability for any other applications.
- Nouvag also recommends the use of a second gripper or a similar retaining device to prevent uncontrolled movement of large morcellated pieces of tissue. This will require an additional percutaneous access.
- Fix the Morcellator in a ventrolateral position to prevent movements inside the uterus.
- Prevent movement of the lateral blood vessels, kidneys and retroperitoneum.

- Pull the tissue towards the rotating blade for morcellation with the gripper.

Risk of injury due to gas loss!

Gas may be lost during the procedure. This may result in the abdominal wall dropping.

- Select an Insufflator (not included into the scope of delivery) with adequate flow to prevent the abdominal wall from dropping.
- Use the Morcellator with protective tube as in option 1: when removing or replacing instruments, immediately close the proximal opening of the seal unit with the thumbs.
- Use the Morcellator with trocar cannula as in option 2: pull the Morcellator with cutting tube out of the trocar cannula. The trocar cannula valve closes to prevent gas loss.
- After removal of the morcellated tissue, insert the Morcellator into the trocar cannula and repeat the process.

Risk of injury!

The Morcellator system is designed for removing tissue with the greatest possible safety. For this purpose the protective tube and the trocar cannula are both fitted with a valve flap.

Morcellation without protective tube or trocar cannula is prohibited.

Tissue escapes the gripper.

Tissue may escape the gripper during the procedure and fall into the abdominal cavity.

Causes:

- The tissue is not gripped tightly enough.

Remedy:

- Grip the tissue at a different position.
- Press the jaws of the gripper together as tightly as possible.

Risk of injury due to heat generated by friction!

The user must ensure that as little friction heat as possible is generated during a procedure on a patient.

High speeds and high application pressure may cause thermal necrosis of tissue.

Risk of injury!

The distal end of the cutting tube is extremely sharp.

Precautions:

- Handle the cutting tube particularly cautiously and carefully.
- Where possible, have the Morcellator with protective tube in protected position (**NO CUT**), particularly if the cutting tube is not used for any period.
- Insert and remove the protective tube and cutting tube under visual control (endoscopy).

Risk of injury!

When inserting the obturator through the percutaneous access, there is risk of injury to blood vessels, intestinal loops or the bladder. This can result in significant complications, which will require a Laparotomy.

Risk of damage to the blade of the cutting tube!

If the jaws of the gripper are not completely closed when it is pulled from the protective tube, the blade will be damaged (see picture 15).

Picture 15. Use of gripper.

Risk of injury!

Damaged blades and cutting tubes must not be used.

Cutting tubes must not be sharpened.

7.2.1 Option 1: morcellation with protective tube.

1. Hysterectomy

Select the procedure for the hysterectomy in accordance with current medical standards.

2. Inserting cutting tube

Make sure that the obturator is inside the cutting tube.

Make sure that the protective tube covers the blade on the cutting tube (set position to NO CUT; "PROTECTED" is also shown on the black ring between the gear unit and protective tube).


Insert Morcellator (with obturator) into an existing percutaneous access under visual control (endoscopy).

3. Using cutting tube

To use the cutting tube, rotate the protective tube to uncover the blade of the cutting tube (unprotected CUT position is displayed). To cover the blade again, rotate the protective tube (protected NO CUT position is displayed; "PROTECTED" is also shown on the black ring between the gear unit and protective tube).

4. Inserting gripper

Pull the obturator out of the protective tube. Close the proximal opening of the gear unit with the thumbs immediately. Insert gripper carefully through the proximal opening of the gear unit.

5. Removing tissue


Grip the tissue with the gripper under good visual control (endoscopy).

The handle lock is set.

Start the electric motor of the Morcellator set. Pull the tissue towards the rotating blade for morcellation with the gripper.

As the cut progresses pull the gripper with the tissue further out of the Morcellator.

Carefully pull the complete gripper with the separated tissue out of the Morcellator and the patient.

The retracted mouth must not touch the blade of the cutting tube. Repeat the process until the tissue has been completely removed.

Make sure that the proximal opening of the seal unit is closed immediately with the thumb when there is no instrument in the Morcellator to prevent large quantities of gas from escaping.

7.2.2 Option 2: morcellation with trocar cannula


1. Hysterectomy

Select the procedure for the hysterectomy in accordance with current medical standards.

2. Inserting trocar cannula

Insert trocar cannula with obturator into an existing access under visual control (endoscopy).

Remove obturator.

3. Insert the cutting tube through the trocar cannula.

Actuate the valve flap lever of the trocar cannula and hold in this position.


Carefully insert prepared cutting tube with gear unit and seal unit into the proximal opening of the trocar cannula.

Cutting tube must be fully inserted. Release the valve flap lever of the trocar cannula.

The guide pin of the trocar cannula must be in the distal slot of the gear unit.


Lock the trocar cannula in position with a click.

Make sure that the cutting tube remains in the NO CUT position.

4. Inserting gripper

Insert gripper carefully through the proximal opening of the gear unit and trocar cannula.

5. Preparing cutting tube

Rotate the trocar cannula so the CUT position is visible from above.

If the cutting tube is completely inserted, the cutting edge is now only partially covered by the trocar cannula.


6. Protecting tissue

The nose on the trocar cannula, which partially covers the sharp edge of the blade, must protect tissue that must not be injured. Rotate the Morcellator to the desired position for this.

Protected area


Unprotected cutting area

7. Removing tissue

Grip the required tissue with the gripper.
The handle lock is set.
Start the Morcellator set motor.
Pull the tissue towards the rotating blade with the gripper.
Pull the gripper out of the patient through the gear unit. The retracted mouth of the gripper must not touch the blade.
Repeat the process until the morcellated tissue has been completely removed.

7.2.3 Option 3: myomectomy

The myoma drill has been designed exclusively for removing tissue from a myoma.

1. Inserting Morcellator

Insert Morcellator with protective tube through the percutaneous access.
Keep gripper ready to remove the tissue from the myoma drill.

2. Inserting myoma drill


Insert trocar cannula through an additional percutaneous access. Open valve flap of trocar cannula and insert myoma drill.

3. Uncovering myoma

Remove tissue from the myoma with the myoma drill.

4. Inserting gripper

Insert gripper through the second access and remove separated tissue from the myoma drill. Remove myoma drill.

5. Removing separated tissue

If the piece of tissue is larger than the diameter of the cutting tube, start the Morcellator set electric motor. Pull the tissue towards the rotating blade with the gripper.
(Smaller pieces of tissue can be removed without morcellation).
Carefully pull the complete gripper out of the Morcellator and the patient. The retracted mouth must not touch the blade of the cutting tube.

6. Removing additional tissue

If necessary, re-insert the myoma drill and repeat the above process.

8. After use

1. Switching off

Switch off the power switch on the back of the device ("O" position).

2. Removing protective tube (when using Morcellator with protective tube)

Rotate protective tube to NO CUT position.

Close the proximal opening of the gear unit with the thumb to prevent gas loss. Hold the protective tube and carefully disconnect gear unit with cutting tube and remove from patient. The protective tube can also be removed from the percutaneous access.

If applicable, remove any other instruments.

3. Removing cutting tube (when using Morcellator with trocar cannula)

Close the proximal opening of the gear unit with the thumb to prevent gas loss.

Hold the trocar cannula and carefully disconnect gear unit with cutting tube from the trocar cannula and remove from patient. The valve flap closes automatically as soon as the cutting tube is pulled out.

The trocar cannula can remain as access throughout the remaining procedure.


4. Disconnecting gear unit connections

Pull the electric motor coupling from the gear unit's motor connection.

Pull the electric motor device connection plug from the electric motor connection of the control unit.


Pull the electric motor out of the handle. Lift and hold the release lever of the handle.

Pull the handle from the gear unit's electric motor connection.

5. Disconnecting pedal connections

Pull the pedal device connection plug from the pedal connection of the control unit.

6. Preparing for preparation at the location of use

Take the product to the preparation area immediately. Prepare the product as specified in this manual.


7. Disassembling trocar cannula

Disassembling trocar cannula (four steps)

Disassemble trocar cannula before decontamination.

Disassemble trocar cannula and rinse in deionized water during disassembly.

1. Unscrew sealing cap holder from the proximal end of the trocar cannula.
2. Remove the internal seal from the sealing cap holder.
3. Unscrew the valve body from the main body.
4. Remove the red O-ring from the valve body.


8. Disassembly of electric motor

Disassembling electric motor (three steps).

Disassemble electric motor before contamination.

1. Unscrew handpiece holder from electric motor
2. Unscrew motor cap from electric motor
3. Pull motor plug from electric motor.


9. Disassembly of gear unit and Morcellator with protective tube

Disassembling seal unit

Remove membrane screw from seal unit and remove membrane seal.

Remove seal holder from gear unit and remove cross seal.

Risk of injury for the user!

The distal end of the cutting tube is very sharp, handle the cutting tube very carefully.


Removing protective tube and cutting tube

Make sure that the lock nut of the protective tube covers the distal end of the cutting tube (protected NO CUT position is visible).

Hold the gear unit with one hand so the distal end of the cutting tube does not point to the user or any other person.

With the other hand carefully disconnect the protective tube from the gear unit and pull it off.

Hold the cutting tube with the same hand and press it lightly against the gear unit until it releases itself from the handle and slides out smoothly. It can now be pulled out of the proximal opening of the gear unit.

9. SCOPE OF DELIVERY

39

The TCM 3000 BL – Morcellator for Surgical Procedures with Accessories

1.	3281	TCM 3000 BL Morcellator control unit	1 item,
2.	5163nou	Gear unit with seal unit,	1 item,
3.	5183nou	Ergonomic handle,	1 item,
4.	1507nou	Vario foot switch pedal with connecting cable 3 m,	1 item,
5.	22262	Power cable, 3 m long,	1 item,
6.	1928	Adapter for cleaning and lubricating gear unit,	1 item,
7.	1974	Adapter for cleaning and lubricating electric motor,	1 item,
8.	31946	Operation Manual on CD-ROM,	1 item,
9.	21606	Fuse T1AL 250WAC,	2 items
10.	5154nou	Cutting tube, Ø 12 mm (delivering on request),	up to 20 items
11.	5155nou	Cutting tube, Ø 15 mm (delivering on request),	up to 20 items
12.	5156nou	Cutting tube, Ø 20 mm (delivering on request) ,	up to 20 items

13.	5141nou	Trocar cannula for Ø 12 mm instruments (delivering on request) ,	up to 20 items
14.	5142nou	Trocar cannula for Ø 15 mm instruments (delivering on request),	up to 20 items
15.	5151nou	Obturator, Ø 12 mm (delivering on request),	up to 20 items
16.	5152nou	Obturator, Ø 15 mm (delivering on request),	up to 20 items
17.	5153nou	Obturator, Ø 20 mm (delivering on request),	up to 20 items
18.	5137nou	Protective tube, Ø 12 mm (delivering on request),	up to 20 items
19.	5138nou	Protective tube, Ø 15 mm (delivering on request),	up to 20 items
20.	5139nou	Protective tube, Ø 20 mm (delivering on request),	up to 20 items
21.	5193	Myoma drill, Ø 10 mm (delivering on request),	up to 20 items
22.	5194	Tenaculum gripper, Ø 5 mm (delivering on request),	up to 20 items
23.	5195	Claw gripper, Ø 10 mm (delivering on request),	up to 20 items
24.	5196	Tenaculum gripper, Ø 10 mm (delivering on request),	up to 20 items
25.	5192	Claw gripper, Ø 13.5/15 mm (delivering on request),	up to 20 items
26.	2090nou	E-Motor 21 electric motor with 3 m cable (delivering on request),	up to 20 items
27.	76052	Motor cable, 3 m long (delivering on request),	up to 20 items

On device receipt the following actions shall be completed:

- Ensure that the shipping package contains all parts listed in Section Scope of Delivery.
- If any parts are damaged or missing, contact the Nouvag agent from which the device was ordered or Nouvag AG.

To carry out immediate inspection of the TCM 3000 BL – Morcellator for Surgical Procedures with Accessories compare the received device with the below pictures of the set components:

- | | |
|---|---|
| 1. Control unit (3281) | 8. 2 fuses (21606) |
| 2. Power cable, 3 m long (22262) | 9. Vario pedal (1507nou) |
| 3. Gear unit (5163nou) | 10. Adapter for spray for cleaning and lubricating gear unit (1928) |
| 4. Seal holder (5136nou) | 11. Adapter for cleaning and lubricating electric motor (1974) |
| 5. Cross seal for all cutting tube diameters (5167nou) | 12. Ergonomic handle (5183nou) |
| 6. Membrane seal for all cutting tube diameters, 10 pcs per a package (5166nou) | 13. Operation Manual on CD-ROM (31946) |
| 7. E-Motor 21 electric motor (2090nou) | |

Cutting tube

Protective tube

Obturator

Trocar cannula

Myoma drill

Gripper


Picture 16. Visual check of the scope of delivery.

 The operation manual is not available in printed form with this product. The operation manual as PDF is on the CD-ROM that was included in the delivery. A CD-ROM drive and the Adobe Acrobat software is needed to read or print this operation manual. To order the operation manual in printed form please check on the type plate for the address or the website. If you send us the article number (31946 RUS) with the language code and the note "Printed copy", a printed manual will be sent to you free of charge.

 In line with regulations pertaining to hazardous materials, the item Nou-Clean care spray (Ref. 1984) is not delivered with the control unit but can be ordered separately from any official Nouvag service centre. www.nouvag.com

10. CLEANING OF DEVICE COMPONENTS

Schematic flow chart as per EN ISO 17664

General Information about the disinfection, presterilization cleaning and sterilization of MDs

- Multiple use MDs should be repeatedly sterilized and stored afterwards under conditions eliminating any potential for secondary contamination with microorganisms.
- Sterilization of such MDs should be performed immediately before the first use of a device, and after each patient to avoid cross-contamination.
- Disassembly, washing and rinsing of MDs which contacted patients blood and serum should be performed in rubber gloves, goggles, medical face mask, gown and wearing a waterproof apron over the gown. Personal hygiene rules should be strictly observed.
- Containers with working solutions of detergents and disinfectants should have lids, have clear labels with name of an agent, its concentration, intended use, date of preparation, best before date. For agents ready to use indicate name of an agent, and state date of beginning of its use.
- In case of disinfection and pre-sterilization cleaning with solutions of chemical agents MDs are submerged into a working solution of the agent to fill all channels and cavities. Dismountable devices should be submerged in a dismounted state, devices with locks should be soaked in open position, with making several operational moves with such devices.
- Disinfection by wiping may be performed for MDs which does not come into direct contact with a patient, or when their design prevents from using the submersion method.
- After disinfection multiple-use MDs should be cleaned from residues of a disinfectant according to recommendations provided in instructions for use of the specific disinfection agent.
- Degree of cleanliness should be confirmed by making tests: amidopyrinum test for presence of residual blood and phenolphthalein test for residues of a detergent. Control test must be put every day. Control should be performed on 1% of simultaneously treated devices of every name but not less then 3 pieces. Both external surfaces, and internal surfaces of tools should be checked! The document confirming the control test was fulfilled is a logbook for the registration the quality of the tests.
- This MD is sterilized by a physical method (autoclaving).
- Sterilization should be performed in a functionally detached premise in a "clean area".
- MD should be transported to the "clean area" in lockable sterilized containers.
- Staff performing operations of preliminary cleaning and disinfection of MDs should wear special garments, and a mask, goggles and rubber gloves, over the gown – a waterproof apron!
- Sterilized MDs should be stored in conditions preventing from secondary contamination with microorganisms, in a special cabinet protected from dust, moisture and temperature fluctuations.
- Shelf life should be stated on package and determined by type of a packing material according to instructions for use.
- It is not allowed to use sterilized MDs with an expired shelf life after sterilization.
- Sterilization of MDs should be recorded in a registration log.

44

General Recommendations in Preparation of MDs

- The preparation procedure includes cleaning, washing, disinfection, washing and drying.
- Devices should be cleaned only in a cleaning solution (e.g. Schülke perform or 3% DESOMEDAN ID), with observance of personal safety measures. When using solutions strictly follow instructions for use provided by the manufacturer!
- Use only deionized water (VE-water) for rinsing and washing.

- Use only fresh working solutions prepared every day. In cases of gross contamination it is recommended to change solution frequently.
- It is recommended to use soft fibreless tissues, plastic brushes or washing guns. It is forbidden to use metal items. Treatment should be performed in plastic or enamel-lined (without damaged enamel surface) containers with lockable lids.
- Manual cleaning of external surfaces should be performed by a soft disposable tissue, disinfected brushes. Brush cleaning should be performed over a solution surface, in order to avoid splashing of contaminated liquid
- For effective cleaning swivel tools (forceps) should be open, to reduce overlapping surfaces to the minimum. Dismountable devices should be dismantled.
- Hollow tools should be cleaned by a rinsing gun (rinsing pressure (3 bar – 4 bar).
- Do not clean morcellator with pressurized air.
- Do not clean electric motor with ultrasound or pressurized air.
- Disinfection of MD should be performed concurrently with cleaning.
- Machine cleaning and disinfection should be performed in a disinfecting washer (e.g. thermal disinfecter Miele G7835 CD) with performance of thermal disinfection step.
- For avoidance of damaging, parts of a MD should be properly fixed in the disinfecting washer. Loading devices should ensure proper rinsing, including internal cavities of devices, though the appropriate fittings
- Drying of a tool is a constituent part of the procedure of machine cleaning and disinfection.
- Upon completion devices are removed manually, using disposable sterile gloves, treated with an antiseptic agent.
- Control is performed visually. Critical parts, like handles, hinge joints, and cutting edges should be checked with special attention. All devices with internal channels should be checked for penetrability. If any devices fail to pass the check, such devices should be replaced.
- Insufficiently cleaned devices should be repeatedly cleaned, and then washed thoroughly.
- Tools with micro cracks near hinge joints and/or damaged, bent or otherwise worn devices should be replaced, as they cannot perform or may insufficiently perform their functions.
- If tools require any steps of care, the appropriate measures, generally, are performed before checking their functionality
- Before checking functionality it is necessary to grease hinge and threaded tools. It is necessary to ensure proper functioning of a tool by checking. For this purpose all dismantled tools should be re-assembled. If necessary, after the check tools may be dismantled again for sterilization. Dismounting/assembling should be performed according to instructions provided by the manufacturer.
- Damaged or faulty medical tools should be sent for repair or disposed.
- It is recommended to pack MDs into a sterilization bag. Follow the directions of the manufacture on the using and the time for keeping sterility!
- Package should provide for the possibility to label with the following information:
 - sterilization date,
 - packer,
 - shelf-life,
 - contents.
- If sterilization is performed without bag, shelf life of sterilized MDs is: 6 hours if sterilized in lattice boxes, and 20 days if sterilized in boxes with filters; kraft-bags can be kept in sterile conditions up to three days. Tools were sterilized without bags must be used immediately after sterilization.

- Sterilize in vacuum autoclave (class B or S by EN 13060) at 135°C. Sterilization parameters:

Sterilization mode					
Sterilization chamber steam pressure, MPa (kgf/cm ²)		Sterilization chamber operational temperature, °C		Sterilization hold time, min	
nominal value	limit	nominal value	limit	nominal value	maximum deviation
0.20	±0,02	135	±2	5-25	±2

- To control the sterilization the manufacture recommend to use the sterility indicator (e.g. indication bands Comply™ by the company 3M).
- The manufacturer confirms validation of the process of repeated sterilization by a Certificate (see Appendix 3)
- Returned tools are medical devices which, regardless of whether they were used or not, are sent back to the manufacturer.
- Possible causes of return of a tool may be the need to repair or maintain it. If a MD is returned to the manufacturer or authorized representative, they may be sent back to a user only if they:
 - were cleaned and disinfected according to instructions of a producer, are dry and declared "sanitary clean", and/or
 - clearly marked as not contaminated, and packed properly.
- Sanitation of returned goods, just like in general circulation, should be performed without delays, to avoid damaging of a tool. Sanitation should not be performed if it may modify or destruct a product, tamper results of a test or make such test impossible. In cases of doubt one should consult with the manufacturer.
- The following information should be included to accompanying documents to each individual returned medical device:
 - intended use of the medical device;
 - sanitation method;
 - date of treatment;
 - surname of a person who performed treatment;
 - reason for return.

Instrument preparation

Risk of infection!

Prepare the product before its first use and before every subsequent use as specified in this manual.

Inadequate and/or incomplete preparation of the product may cause an infection of the patient. Observe the following specific instructions:

- Do not clean Morcellator with compressed air.
- Use only cleaning and care spray recommended by Nuvag AG – other care sprays may destroy the Morcellator.
- At unpacked steam sterilization of single items it must be regarded that they have a stable

seat in the tray to prevent them from rolling around, causing damage.

- Allow the components of the 3000 BL – Morcellator for Surgical Procedures with Accessories to cool down after steam sterilization!

Immersing instruments in cleaning solution

- Never exceed the maximum concentration and exposure time specified by the manufacturer of the cleaning and disinfection solution.
- Air bubbles must not adhere to the instruments.
- All components of the instruments must be completely immersed in the cleaning solution.
- All lumina of the instruments must be completely filled with cleaning solution and free from bubbles.

- Do not use cleaning agents that contain solvents.
- If sterilized device components are not used immediately, a sterilization indicator must be placed on the packaging and it must be labeled with the sterilization date.

- Perform cleaning, disinfection and sterilization after every treatment.
- To autoclave single items, such as cutting tubes, Nouvag AG recommends the use of sterilization packages to prevent of damage.
- Make sure that sterilization packaging is no more than 80 % full.
- Autoclave material for at least five minutes at max. 135 °C.
- If sterilized device components are not used immediately, the sterilization packaging must be labeled with the sterilization date.
- Nouvag AG recommends including a sterility indicator.

47

10.1. Control unit and pedal

The control unit and pedal do not come into contact with the patient.

Wipe the outside using microbiologically tested surface disinfectant or a 70 % -isopropyl alcohol solution. The front plate of the control unit is sealed accordingly for this purpose and can be wiped clean.

Risk of damage!

Do not use incompatible preparation methods. The product may be damaged.
Before disinfection and cleaning using liquids unplug the device.

Risk of damage to the TCM 3000 BL - Morcellator control panel

- Wipe away dirt or dust with a soft cloth.
- Use a moist cloth for stubborn dirt.

10.2 Preparation of the trocar cannula

Essential requirements

Avoid long wait times until preparation due to risk of drying and corrosion. Maximum time between use and preparation: 2 hours. The trocar cannula is designed for 500 sterilization cycles.

Disassembly

Before decontamination disassemble the trocar cannula as described in section "Disassembling trocar cannula". Remove seal cap, internal seal and O-ring.

Manual cleaning

- Rinse and clean all surfaces and lumina thoroughly with cleaning solution and suitable brushes.
- Place trocar cannula in an ultrasonic bath with cleaning solution, if necessary.
- Rinse all components of the trocar cannula thoroughly with cleaning gun and deionized water.

Machine cleaning and disinfection

- Arrange all components of the trocar cannula in the cleaning machine.
- Use a fine-mesh basket for small parts, such as seals and O-rings.
- Run a standard machine cleaning and thermal disinfection cycle.

Autoclaving

Autoclave in vacuum autoclave (class B or S as per EN 13060) at 135°C for at least five minutes*. Do not exceed the maximum load of the sterilizer when sterilizing multiple instruments in one sterilization cycle. Autoclaving without vacuum must be followed by a drying phase. Allow the trocar cannula to dry in a bag for at least one hour at room temperature with the paper side up.

** The temperature holding times conform to the national guidelines and standards with the maximum holding time 25 minutes.*

48

Assembly

- Make sure that all parts of the trocar cannula are completely dry after preparation.
- The flap in the valve body must be open when screwing the spacer with internal seal to the valve body.
- Assembly the trocar cannula in reverse order of disassembly (see Section "Disassembling trocar cannula").

10.3. Preparation of E-Motor 21 Electric Motor

Frequent cleaning has little effect on the electric motor. The end of the product service life normally results from wear and damage caused by use. The electric motor is designed for 250 sterilization cycles.

Risk of damage

- Electric motor must not be cleaned with ultrasonic or compressed air.
- Do not kink the electric motor cable, otherwise it may be broken.
- Use Nouvag AG cleaning and care spray exclusively, other care sprays may destroy the electric motor.
- Steam sterilize the electric motor in a sterilization bag and then allow to cool..

At the site of use

Remove surface dirt with a napkin or paper towel.

Storage and transport

Avoid long wait times before preparation because of the risk of drying and corrosion; maximum period between use and preparation is eight hours.

Preparation for cleaning

Remove dirt on the electric motor with a napkin or paper towel. Unscrew motor cap and electric motor clutch and remove the cable. Do not place electric motor in ultrasonic bath.

Disassembly

Before decontamination, disassemble motor as described in Section "Disassembling electric motor".

Manual cleaning

Use a neutral cleaning agent, cleaning brushes and demineralized water (< 38°C).

- Rinse and brush off dirt on the electric motor, motor cap, cable and electric motor clutch.
- Apply cleaning agent to all surfaces and intermediate spaces with brush and brush thoroughly.
- Rinse electric motor, motor cap with cable and electric motor clutch thoroughly under running water.

Automatic cleaning and disinfection

Equipment: cleaning and disinfecting machine with special racks that connect the electric motor to the cleaning and disinfecting machine and the channels. Wash the electric motor from the front. Use only neutral cleaning agents.

- Position electric motor in rack (ensure that channels are clean and connected to the cleaning and disinfecting machine).
- Place the motor cap and cable, and electric motor clutch, in the basket.
- Set a cleaning cycle that offers sufficient cleaning and rinsing. Perform the final rinse with deionized water.
- Perform a 10-minute rinse cycle at 93°C to facilitate thermal disinfection.
- When removing the electric motor, check the motor cap and cable and electric motor clutch to verify whether soiling is still visible in the gaps and the grooves. If necessary, repeat the cycle or clean manually.

Manual disinfection

For manual disinfection, spray the electric motor, motor cap with the plug and cable underneath it, and the electric motor clutch with isopropyl-based disinfectant spray.

Drying

If a drying programme is not provided by the washer-disinfector, the electric motor must be dried in a drying cabinet at 60°C.

Inspection, maintenance and assembly

- Perform a visual inspection to check for damage, corrosion and wear and tear.
- Screw the electric motor to the cleaning and lubrication spray adapter (Ref. 1974) and spray with Nou-Clean spray (not included into the scope of delivery and is an option) for about three seconds and briefly wipe with a moist cloth (see instructions on spray can).
- After spraying, attach the electric motor clutch and screw the cable with motor cap onto the electric motor.

Packaging

- Storage of single unit: Pack the electric motor in individual packaging for sterile items. The bag must be large enough to ensure that the seal is not subject to strain. When using sterilization container, observe the instructions of the manufacturer.
- Storage of set: Sort electric motors on trays intended for this purpose or place them on all-purpose sterilization trays.

Sterilization

- Autoclave in the vacuum autoclave (Class B or S as per EN 13060) at 135°C for at least 5 minutes*.
- When sterilizing several instruments during one sterilization cycle, do not exceed the maximum sterilizer load.
- A drying cycle must be added in the case of autoclaves without a post-vacuum function. Allow the electric motor to dry in the bag for at least 1 hour at room temperature with the paper side facing upwards.

50

** Temperature exposure times are based on country-specific guidelines and standards. Maximum exposure time: 25 minutes.*

Storage

If the sterilized electric motor is not used immediately after sterilization, the material packaging must be labeled with the sterilization date. It is advisable to add a sterilization indicator.

10.4. Preparation of gear unit

The gear unit is designed for 250 sterilization cycles.

 Do not immerse gear unit in liquids.

Disassembly

Disassemble gear unit as described in section "Disassembly of gear unit and Morcellator with protective tube".

Manual cleaning

Use neutral cleaning agent, soft brushes and demineralized water.

- Rinse and brush surfaces and intermediate spaces thoroughly with neutral cleaning agent.
- Brush and rinse lumen thoroughly with neutral cleaning agent.
- Rinse all surfaces and lumina for at least 10 minutes. Make sure that all intermediate spaces are cleaned.

Machine cleaning and disinfection

- Equipment: cleaning and disinfecting machine with special racks that ensure that channels are cleaned. Use only neutral cleaning agents.
- Place the gear unit in the load carrier (ensure that the channels can be rinsed).
- Place seal holder and cross and membrane seals in basket.
- Set a cleaning cycle that offers sufficient cleaning and rinsing. Perform the final rinse with deionized water.
- Perform a 10-minute rinse cycle at 93°C to facilitate thermal disinfection.
- When removing the gear unit, the seal holder, the cross and membrane seals, check whether visible contamination still remains in the gaps and the grooves. If necessary, repeat the cycle or clean manually.

Maintenance and lubricating gear unit

- Perform a visual inspection to check for damage, corrosion and wear and tear.
- Use Nou-Clean spray (not included into the scope of delivery and is an option) exclusively.
- Attach Nou-Clean spray to the adapter for cleaning and lubricating the gear unit (Ref. 1928) to and spray gear unit for about three seconds and briefly wipe with a moist cloth (see instructions on spray can)).

51

Sterilizing gear unit

Autoclave in the vacuum autoclave (Class B or S as per EN 13060) at 135°C for at least 5 minutes*. When sterilizing several instruments during one sterilization cycle, do not exceed the maximum sterilizer load. A drying cycle must be added in the case of autoclaves without a post-vacuum function. Allow the gear unit to dry in the bag for at least 1 hour at room temperature with the paper side facing upwards.

** Temperature exposure times are based on country-specific guidelines and standards. Maximum exposure time: 25 minutes.*

Assembly

- Make sure that all parts of the gear unit are completely dry after preparation.
- Assemble gear unit in reverse order of disassembly (see Chapter Use).

10.5. Preparation of cutting tubes, protective tubes, obturators, grippers, seals and myoma drill

Clean all products with suitable cleaning brushes.

Reprocessing restrictions

The end of the product service life is normally determined by wear and tear and damage through use. The cutting tubes are designed for a service life of 20 cycles.

At the location of use

Remove surface soiling with a napkin/paper towel.

Storage and transport

No special requirements - avoid long periods before preparation to prevent the danger of drying-out and corrosion, maximum time between use and preparation is two hours. If sterilized instruments are not used immediately, the sterilization packaging must be labeled with the sterilization date. Nouvag AG recommends using sterility indicators.

Manual cleaning

- Use disinfectant solution, neutral cleaning agent, suitable cleaning brushes and demineralized water. Wipe all surfaces thoroughly with disinfectant (exposure time as specified by the disinfectant manufacturer).
- Immerse all parts in the neutral cleaning agent (exposure time as specified by the cleaning agent manufacturer). Use suitable brushes that can reach all parts of the lumen for cavities.
- Rinse all parts thoroughly under running water.

Automatic cleaning and disinfection

Equipment: Cleaning/disinfection device with a special load carrier which facilitates the connection of the instruments to the cleaning/disinfection device and the channels. Use a fine-mesh basket for small parts. Use neutral cleaning agent at the concentration recommended by the manufacturer.

- Place the cutting tubes, protective tubes, myoma drills, grippers and ergonomic handle in the load carrier (ensure that the channels can be rinsed).
- Place seal holder, membrane screw, membrane seal, cross seal, obturator in the fine-mesh basket and close the basket.
- Set a cleaning cycle that offers sufficient cleaning and rinsing. Perform the final rinse with deionized water.
- Perform a 10-minute rinse cycle at 93°C to facilitate thermal disinfection.
- When removing the products, check whether visible contamination still remains in the gaps and the grooves. If necessary, repeat the cycle or clean manually.

Drying

If a drying program is not provided by the washer-disinfector, all parts must be dried in a drying cabinet at 60°C.

Inspection and maintenance

- Perform a visual inspection to check for damage, corrosion and wear and tear.
- Check cutting edge for damage and rotary components for bending.
- Check that rotary components move smoothly.
- Check seals and seal holders for elasticity, cracking, deformation and porosity.

Packaging

- Storage of single unit: Pack the parts in individual packaging for sterile items. The bags must be large enough to ensure that the seal is not subject to strain. When using sterilization container, observe the instructions of the manufacturer.
- Storage of set: Sort the instruments on trays intended for this purpose or place them on all-purpose sterilization trays.

Risk of damage

The blade of the cutting tube may damage the sterile packaging.

Sterilization of ergonomic handle, gripper, cutting tubes, protective tubes, obturator, myoma drills, seal and membrane screw

Autoclave in the vacuum autoclave (Class B or S as per EN 13060) at 135°C for at least 5 minutes*. When sterilizing several instruments during one sterilization cycle, do not exceed the maximum sterilizer load. A drying cycle must be added in the case of autoclaves without a post-vacuum function. Allow the instruments to dry in the bag for at least 1 hour at room temperature with the paper side facing upwards.

** Temperature exposure times are based on country-specific guidelines and standards. Maximum exposure time: 25 minutes.*

53

11. Maintenance and repair

11.1 General maintenance

Before every use perform a visual inspection of the device to check for damage and wear and tear. Check cutting edges and rotary components for bending. Check rotary components for gaps.

Users can perform the following maintenance procedures themselves:

Maintenance, electric motor cable

- Replace a damaged electric motor cable with a new one (Ref. 76052).

Maintenance, protective tube

- Replace damaged protective tubes with new ones.
Protective tube, Ø 12 mm (Ref. 5137nou)

Protective tube, Ø 15 mm (Ref. 5138nou)

Protective tube, Ø 20 mm (Ref. 5139nou)

Maintenance, obturator

- Replace damaged obturators.
Obturator, Ø 12 mm (Ref. 5151nou)
Obturator, Ø 15 mm (Ref. 5152nou)
Obturator, Ø 20 mm (Ref. 5153nou)

Maintenance, trocar cannula seal caps

- Replace damaged seal caps with new ones.
Seal cap set for trocar cannula Ø 12/15 mm (Ref. 5177nou), 10 items per package
- Replace damaged seal caps with new ones.
O-ring for trocar cannula Ø 12/15 mm (Ref. 5180nou), 10 items per package
- Replace damaged seal cap holder with a new one.
Seal cap holder Ø 12 mm (Ref. 51484nou), 1 item per package
Seal cap holder Ø 15 mm (Ref. 51502nou), 1 item per package.

Maintenance, seal holder

- Replace damaged parts with new ones.
Seal holder (Ref. 5136nou), 1 item per package

Maintenance, seals

- Replace damaged seals with new ones.
Membrane seal for suitable for all cutting tube diameters (Ref. 5166nou), 10 items per package.
Cross seal for suitable for all cutting tube diameters (Ref. 5167nou), 10 items per package.

54

Maintenance, adapter for Nou-Clean spray (spray is not included into the scope of delivery and shall be ordered and supplied separately)

- Replace necessary adapters with new ones.
Adapter for cleaning and lubricating electric motor (Ref. 1974)
Adapter for cleaning and lubricating gear unit (Ref. 1928)

11.2 Replacing the control unit fuses

Users can replace faulty control unit fuses themselves. These are located at the rear of the device in the fuse slot beside the power switch:

1. Unplug the power plug.
2. Open the fuse slot using a screw driver.
3. Replace the faulty fuse T 1 AL 250 VAC.
4. Slide the fuse holder back in and close the fuse slot.
5. Check the mains voltage shown on the fuse slot.
6. Plug in the power plug again.

Picture 17. Fuses and their replacement.

Picture symbols:

1. Fuse slot locking mechanism
2. Display window for voltage setting
3. Fuse slot
4. Fuse 1
5. Fuse 2

11.3. Malfunction and troubleshooting

Table 8.

Malfunction	Cause	Solution	Refer to the Operation Manual
Device doesn't operate	Control panel not switched on	Set the power switch "I/O" to "I"	Section "Assembly of the Morcellator set with protective tube" (option 1) Section "Assembly of the Morcellator set with trocar cannula" (option 2)
	Power connection not established	Connect the control panel	Section "Assembly of the Morcellator set with protective tube" (option 1) Section "Assembly of the Morcellator set with trocar cannula" (option 2)
	Incorrect operating voltage	Check the mains voltage	Section "Connection to the power supply"
	Faulty fuse	Replace the fuse	Section "Replacing the control panel fuses"
Electric motor doesn't run	Electric motor not switched on	Switch on the electric motor using the treadplate	Section "Using the control panel"
	Electric motor not connected	Connect the motor cable to the control unit	Section "Assembly of the Morcellator set with protective tube" (option 1) Section "Assembly of the Morcellator set with trocar cannula"

			(option 2)
	Gear unit not correctly installed	Press the electric motor firmly onto the gear unit coupling until it clicks into position and check that it is secure by moving it slightly in the opposite direction.	Section "Assembly of the Morcellator set with protective tube" (option 1) Section "Assembly of the Morcellator set with trocar cannula" (option 2)
Foot switch doesn't operate	Foot switch not connected	Plug foot switch cable into the foot switch socket.	Section "Assembly of the Morcellator set with protective tube" (option 1) Section "Assembly of the Morcellator set with trocar cannula" (option 2)

12 Safety information

Your safety, the safety of your team, and of course that of your patients is very important to us. It is therefore essential to bear the following information in mind:

12.1 EMC Manufacturer's Declaration of Conformity

Please observe the information on electromagnetic compatibility provided with these Operation Manual (see Appendix 4).

12.2 Modification and misuse

- Replacement of device parts with external ones or modification of the TCM 3000 BL - Morcellator and its accessories is prohibited. The manufacturer is not liable for any damages resulting from unauthorized modifications or manipulations. The warranty will be cancelled.
- Use of the TCM 3000 BL - Morcellator for Surgical Procedures with Accessories outside the indications described in Section **PURPOSE** is prohibited. The user or operator is solely responsible for any adverse consequences resulting from such use.

12.3 Use precautions

In order to avoid risks and complications you must follow the below precautions:

- The cutting tube of the morcellator must be moved against the abdomen in one direction only;
- Cut off tissues to be removed must be extracted through a hollow tube;
- Uncontrolled movement of morcellator tubes and/or instruments towards lateral vessels, bowels or retroperitoneum must be excluded;

- The morcellator must be ventrolaterally positioned during the whole procedure in order to avoid bowels or vessels injuries;
- Sharp edge of the morcellator cutting tube must be visible during the whole procedure;
- The device must be used by qualified and experienced surgeons only;
- The device must be operated as per the Operation Manual and be functional at all times;
- Rotating tube must be set into motion under continuous visual control only.

12.4 List of operational hazards, risk analysis

The manufacturer runs the system of risk analysis and monitoring. Risk analysis is performed during post marketing monitoring of results of device use. Risk analysis is a part of the device technical documents. Residual risks of use of the TCM 3000 BL – Morcellator for Surgical Procedures with Accessories include injuries of vessels and other anatomical sites.

12.5 Safety requirements, operation conditions and rules

This product is designed solely for use by medically qualified, experienced and trained persons. Surgical procedure using the TCM 3000 BL – Morcellator for Surgical Procedures with Accessories must be completed under continuous supervision of medically qualified staff able to react to any possible adverse events and complications immediately.

To avoid injuries no device manipulations must be done while the electric motor is switched on.

When the device is not used, the protective tube must be pulled out to completely cover the blade of the cutting tube to avoid injuries.

57

Ambient conditions

	Operation:
Relative humidity:	30% - 80%
Temperature:	10°C - 45°C
Atmospheric pressure:	800 hPa - 1060 hPa

The medical device is always safe for the environment while being used..

12.6. Safety check

Safety check is a regular inspection being mandatory for valid medical devices. This procedure is intended for early detection of any device defects and risks for patients, users and third parties.

Qualified specialists inspect the TCM 3000 BL – Morcellator for Surgical Procedures with Accessories **every two years** and report inspection results in a special Device Maintenance Record.

If you have any questions regarding device maintenance and repair as well as after-sales servicing please contact our authorized representative in the Russian Federation: **ООО MedTechKomplekt, Russia, 125252, Moscow, 19 A, Kuusinena Str., office 118, phone/fax +7 (499) 198-81-09, E-mail: medkomplekt@bk.ru.**

Operation manual for TCM 3000 BL - Morcellator

13. MEDICAL DEVICE CLASSIFICATION

Electrical class	II
Method of protection against electric shock:	BF type
Purpose	Cutting instruments
Frequency of application	Multiple application
Type of perceived mechanical impact	Stationary
Type of human body contact	Short-term contact (less than 24 h) with milieu interieur

14. COMPLIANCE WITH INTERNATIONAL STANDARDS

The TCM 3000 BL – Morcellator for Surgical Procedures with Accessories complies with Council Directive 93/42/EEC on Medical Devices:

EN60601.1:2006	General electrical safety standards
EN60601-1-2/07.2007	General requirements for basic safety and essential performance. Collateral standard. Electromagnetic compatibility. Requirements and testing.
EN 62366-2008	Medical devices. Application of usability engineering to medical devices.
EN61000.4.2	(EMC). Part 4: Testing and measuring techniques.– Section 2: Electrostatic discharge immunity test. Basic EMC Publication.
EN61000.4.3	(EMC). Part 4: Testing and measuring techniques.– Section 3: Radiated, radio-frequency, electromagnetic field immunity test. Basic EMC Publication.
EN61000.4.4	(EMC). Part 4: Testing and measuring techniques.– Section 4: Electrical fast transient/burst immunity test. Basic EMC Publication.
EN61000.4.5	Electromagnetic compatibility (EMC). Part 4: Testing and measuring techniques. Section 5. Surge immunity test.
EN61000.4.6	Electromagnetic compatibility (EMC). Part 4: Testing and measuring techniques. Section 6. Immunity to conducted disturbances, induced by radio-frequency fields.
EN61000.4.8	(Electromagnetic compatibility (EMC). Part 4: Testing and measuring techniques. Section 8. Power frequency magnetic field immunity test.

EN61000.4.11	Electromagnetic compatibility (EMC). Part 4: Testing and measuring techniques. Section 11. Voltage dips, short interruption and voltage variations immunity test.
EN ISO 10993-1:2009	Biological evaluation of medical devices. Part 1: Evaluation and testing within a risk management
ISO 3964	Dental handpieces. Coupling dimensions.
EN ISO 17664:2004	Sterilization of medical devices - Information to be provided by the manufacturer for the processing of resterilizable medical devices
ISO 13485:2012	Design and development, production, distribution and servicing of medical devices
CE certificate	Compliance with Council Directive 93/42/EEC on Medical Devices

15. PACKING AND MARKING

Table 9. Explanation of symbols


	Important information		Autoclavable at max. 135° C
	Warning		Suitable for thermal disinfection
	Manufacturer	IPX8	Protection against continuous submersion in water.
	Applied part of BF type is the Morcellator.		Protection class II
	Certified by the Canadian Standards Association (CSA).		Electrical and electronic devices that have reached the end of their service life comprise hazardous waste and may not be disposed of together with household waste. Local disposal regulations apply.
	CE mark with Notified Body	REF	Symbol indicating the order number
	Intermittent operation, 1 min on, 1 min off for 8 cycles	SN	Symbol indicating the serial number and the date of manufacture (year/month)
	Refer to the Operation Manual		Motor
	Pedal		Biohazard

Picture 18. The plate on the rear panel of the device

60

Symbols:

Name and address of the manufacturer

Medical device name

Symbol referring to the Operation Manual

Information on the country of origin: Switzerland

Operation manual for TCM 3000 BL - Morcellator

REF 3281 Reference number of the medical device

SN 4859U1504R Serial number of the medical device


Indicator of hazardous waste which may not be disposed of together with household waste

Power : 100V~/115V~/230V~ 50-60Hz 110VA

Symbol indicating power consumption

Fuse : T1AL 250VAC Fuses used

Protection class II

Symbol indicating CE Certificate obtained

Symbol indicating that the medical device is certified by the Canadian Standards Association (CSA).

2D datamatrix (QR-code)

Picture 19. Symbol for applied part of BF type.

2. Individual package is marked with the following information:


Picture 20. Individual package of the trocar cannula

Picture 21. Individual package of the obturator

Picture 22. Individual package of the protective tube

63


Picture 23. Marking on individual package

Name and address of the manufacturer

Information on the country of origin: Switzerland

Name of the medical device: Protection sleeve 15 mm

REF 5138nou

Reference number of the medical device

LOT 1936C1602R

Number of medical device lot

CE CE mark with Notified Body

Non-sterile

Autoclavable at max. 135° C

Suitable for thermal disinfection

2D datamatrix (QR-code)

Qty.: 1

Number of medical devices: 1 item.

Stainless steel

Material the medical device is made of: stainless steel

Picture 24. . Marking of the tools for the medical device.

The marking is electro-eroded directly onto the non-working surface of the medical device. The marking includes the following information:

- information on the manufacturer: NOUVAGAG
- Reference number of the medical device: 5193
- material, the medical device is made of: GermanyStainless (stainless steel, produced in Germany)
- CE mark.

3. Manufacturer's carton is marked with the following data:


A shipping package (picture 17) protects the medical device against mechanical stress associated with transportation. The medical device packed for transportation is resistant to vibration and shock. The shipping package is a 395x395x245 mm (WxDxH) carton marked with the below information (pictures 25, 26, 27, 28).

Picture 25. Shipping package

Picture 26. Information on the shipping package.

Picture 27. Marking on the shipping package

Symbols:

Name and address of the manufacturer

Type : TCM 3000 BL Morcellator

Name of the medical device

REF 3281

Reference number

SN

Serial number

Power : 100V~/115V~/230V~ 50-60Hz 110VA

Symbol indicating power consumption

Fuse : T1AL 250VAC

Fuses used

Symbol indicating that CE Certificate has been obtained and that the medical device is certified by the Canadian Standards Association (CSA)

67

2D datamatrix (QR-code)


Indicator of hazardous waste which may not be disposed of together with household waste

Protection class II

Information on the country of origin

Picture 28. Informational block on the shipping package.

Informational block on the shipping package includes:

1. Corporate logo with manufacturer's address and phone numbers.
2. The following notice: «The operation manual is not available in printed form with this product. The operation manual as PDF is on the CD-ROM that was included in the delivery. A CD-ROM drive and the Adobe Acrobat software is needed to read or print this operation manual. To order the operation manual in printed form please check on the type plate for the address or the website. If you send us the article number (31946 RUS) with the language code and the note "Printed copy", a printed manual will be sent to you free of charge».

The shipping package is marked with the following handling symbols.

Picture 29. Handling symbols on the shipping package.

Table 10. Handling symbols on the package.


Keep dry

The device must not be located near moisture sources or in a humid facility


Handle with care.

Fragile goods. Goods shall be handled with care. Such a symbol is used to mark goods which may be damaged in the result of tossing or shock.

Transportation conditions

The symbol is used to mark the transportation requirements (humidity and atmospheric pressure).

Transportation conditions

The symbol is used to mark the transportation requirements (temperature).

The medical device is packed as follows inside the shipping carton:

Picture 30. Shipping package, inside view.

The device is packed into a 370x220x130 mm polyamide box, its accessories are packed into a 375x155x130 mm cardboard box. The box also contains a CD-ROM with the Operation Manual and a warranty card and both boxes are packed into the shipping package.

16. TRANSPORTATION AND STORAGE CONDITIONS.

There are no special storage requirements. The device may be stored at the below conditions. The device packed in the manufacturer's package should be stored in the warehouse and put on the racks. Avoid putting any other equipment or goods on the device.

There are no special transportation requirements as well. The device doesn't require disassembling for transportation. The device may be shipped by any covered vehicles or in multi-purpose containers as per the relevant shipping regulations at the below ambient conditions. Arrangement and attachment of shipping packages must ensure their stable position preventing packages from moving and bumping into each other.

Transportation and storage:

Relative humidity:	Max. 90 %
Temperature:	0°C - 60°C
Atmospheric pressure:	700 gPa - 1060 gPa

17. SERVICE LIFE

Shelf life of the medical device is 5 years. Average service life is 7 years from the date of manufacture provided the device is operated according to the Manual. After the service life is over this medical device must be disposed of. After exploitation resource is over the medical device must be disposed of as per national regulations of the importing country as well as recommendations of the manufacturer.

18. INFORMATION ON DISPOSAL

The medical device shall be disposed of and destroyed as per the classification and rules of assembling, operation, decontamination, handling, storage, transportation, recording and disposal enforced by the competent authority of the importing country (Safety class of medical wastes as per Sanitary Regulations SanPin 2.1.7.2790-10 according to Russian Classification).

Electrical and electronic devices that have reached the end of their service life comprise hazardous waste and may not be disposed of together with household waste. Prevailing national and local disposal regulations apply.

When disposing of the device, device parts and accessories, the requirements specified in legislation must be followed.

To ensure environmental protection, old devices can be returned to the dealer or manufacturer

Electric motors that have reached the end of their service life must not be disposed of with household waste. Electric motors must be sterilized before disposal. The national and local regulations for disposal of electronic waste must be observed.

If they are not sterilized before disposal, please observe the national and local regulations for disposal of infectious waste.

19. WARRANTY COVERAGE

The manufacturer Nouvag AG provides warranty coverage of the medical device. Manufacturer's address is St.Gallerstrasse, 23-25, CH 9403 Goldach, Switzerland. Purchasing the TCM 3000 BL - Morcellator entitles you to a 1-year warranty. If you return the warranty card for registration within four weeks of the date of purchase, warranty coverage will be extended for a further 6 months. Improper use or repair, or failure to observe these instructions, will relieve the manufacturer from any obligation arising from warranty provisions or other claims.

Any claims or questions regarding handling of this medical device please send to:

ООО МедТехКомплект, Russia, 125252, Moscow, 19 A, Kuusinena Str., office 118, phone/fax +7 (499) 198-81-09, E-mail: medkomplekt@bk.ru.

Date: 01.12.2016

Leszek Kazana /Area Sales Manager

Operation manual for TCM 3000 BL - Morcellator

L. Kazana

Appendix 1. Parts of the TCM 3000 BL – Morcellator for Surgical Procedures with Accessories:

I. Main parts

1.	3281	TCM 3000 BL Morcellator control unit	1 item;
2.	5163nou	Gear unit with seal unit,	1 item;
3.	5183nou	Ergonomic handle,	1 item;
4.	1507nou	Vario foot switch pedal with connecting cable 3 m,	1 item;
5.	22262	Power cable, 3 m long,	1 item;
6.	1928	Adapter for cleaning and lubricating gear unit,	1 item;
7.	1974	Adapter for cleaning and lubricating electric motor,	1 item;
8.	31946	Operation Manual on CD-ROM,	1 item;
9.	21606	Fuse T1AL 250WAC,	2 items
10.	5154nou	Cutting tube, Ø 12 mm (delivering on request),	up to 20 items
11.	5155nou	Cutting tube, Ø 15 mm (delivering on request),	up to 20 items
12.	5156nou	Cutting tube, Ø 20 mm (delivering on request),	up to 20 items
13.	5141nou	Trocar cannula for Ø 12 mm instruments (delivering on request),	up to 20 items
14.	5142nou	Trocar cannula for Ø 15 mm instruments (delivering on request),	up to 20 items
15.	5151nou	Obturator, Ø 12 mm (delivering on request),	up to 20 items
16.	5152nou	Obturator, Ø 15 mm (delivering on request),	up to 20 items
17.	5153nou	Obturator, Ø 20 mm (delivering on request),	up to 20 items
18.	5137nou	Protective tube, Ø 12 mm (delivering on request),	up to 20 items
19.	5138nou	Protective tube, Ø 15 mm (delivering on request),	up to 20 items
20.	5139nou	Protective tube, Ø 20 mm (delivering on request),	up to 20 items
21.	5193	Myoma drill, Ø 10 mm (delivering on request),	up to 20 items
22.	5194	Tenaculum gripper, Ø 5 mm (delivering on request),	up to 20 items

23.	5195	Claw gripper, Ø 10 mm (delivering on request),	up to 20 items
24.	5196	Tenaculum gripper, Ø 10 mm (delivering on request),	up to 20 items
25.	5192	Claw gripper, Ø 13.5/15 mm (delivering on request),	up to 20 items
26.	2090nou	E-Motor 21 electric motor with 3 m cable (delivering on request),	up to 20 items
27.	76052	Motor cable, 3 m long (delivering on request),	up to 20 items

II. Accessories

1.	1928	Adapter for cleaning and lubricating gear unit,	up to 20 items
2.	1974	Adapter for cleaning and lubricating electric motor,	up to 20 items
3.	5136nou	Seal holder, 1 item per package,	up to 20 packages
4.	5167nou	Cross seal for all cutting tube diameters, 10 items per package,	up to 20 packages
5.	5166nou	Membrane seal for all cutting tube diameters, 10 items per package,	up to 20 packages
6.	5177nou	Seal caps for trocar cannula for Ø 12 mm and Ø 15 mm instruments, 10 items per package,	up to 20 packages
7.	5180nou	O-seal for Ø 12 mm and Ø 15 mm trocar cannula, 10 items per package,	up to 20 packages
8.	51484nou	Seal cap holder for Ø 12 mm trocar cannula, 1 item per package,	up to 20 packages
9.	51502nou	Seal cap holder for Ø 15 mm trocar cannula, 1 item per package,	up to 20 packages

73

Date: 01.12.2016

Leszek Kazana /Area Sales Manager

L. Kazana

Appendix 2. Table of spare parts and accessories to the medical device.

Table 1. Dimensions and weight of the medical device parts

No	Part	Appearance	Weight (grams), no more than	Dimensions (mm), no more than
1	1928 Adapter for cleaning and lubricating gear unit		8	Diameter - 8 Length - 55 Width - 20
2	1974 Adapter for cleaning and lubricating electric motor		11	Diameter - 10 Length - 64 Width - 21 Thread M20


74

3	5183nou Ergonomic handle		50	Length - 120 Diameter - 28
3	21606 Fuse T1AL 250WAC		7	Length - 20 Diameter - 5

75

4	22262 Power cable, 3 m long		320	Length - 3000
5	31946 Operation Manual on CD-ROM		20	Diameter - 120

Table 2. Component parts and accessories of the medical device.



No	Part	Appearance	Weight (g), no more than	Dimensions (mm), no more than	Material brand
1.	5139nou Protective tube, Ø 20 mm		130	Length 185	Medical stainless steel 1.4305
2.	5138nou Protective tube, Ø 15 mm		100	Length 185	Medical stainless steel 1.4305


3.	5137nou Protective tube, Ø 12 mm		75	Length 185	Medical stainless steel 1.4305
4.	5153nou Obturator Ø 20 mm		83	Length 340	Polyetherimide PEI MT
5.	5152nou Obturator Ø 15 mm		82	Length 340	Polyetherimide PEI MT 78
6.	5151nou Obturator Ø 12 mm		81	Length 340	Polyetherimide PEI MT

7.	5156nou Cutting tube, Ø 20 mm		630	Length 260 Blade thickness 0,4	Medical stainless steel AISI 304
8.	5155nou Cutting tube, Ø 15 mm		610	Length 251 Blade thickness 0,4	Medical stainless steel AISI 304
9.	5154nou Cutting tube, Ø 12 mm		510	Length 251 Blade thickness 0,4	Medical stainless steel AISI 304
10	5142nou Trocac cannula for instruments Ø 15 mm		320	Length 185	Medical stainless steel 1.4305

11	5141nou Trocac cannula for instruments Ø 12 mm		312	Length 185	Medical stainless steel 1.4305
12	5193 Myoma drill, Ø 10 mm		100	Total length 330, Handle diameter 35, Handle length 100, Working length of the drill 30 Right-hand drill	Медицинская нержавеющая сталь AISI 304

80

13	5195 Claw gripper, Ø 10 mm		215	Total length 525 Working length 420 Mouth length 35	Medical stainless steel AISI 304
14	5196 Tenaculum gripper, Ø 10 mm		217	Total length 525 Working length 420 Mouth length 35	Medical stainless steel AISI 304

15	5194 Tenaculum gripper, Ø 5 mm		198	Total length 525 Working length 420 Mouth length 35	Medical stainless steel AISI 304
16	5192 Tenaculum gripper, Ø 13.5/ 15.0 mm		232	Total length 525 Working length 420 Mouth length 35	Medical stainless steel AISI 304

17	2090nou E-Motor 21 with 3 m long cable			300	E-Motor length - 80, diameter - 20, E-Motor clutch, length - 30, diameter - 8	Medical stainless steel 1.4305 Cable Polyetheretherketone (PEEK)
18	76052 Motor cable, 3 m long			20	Length 3000	Polyetheretherketone (PEEK)


21	5136nou Seal holder, 1 item per package		80	Outside diameter 30 Internal thread M30 Inside diameter - 20 Height- 33	Medical stainless steel 1.4305
23	5167nou Cross seal for all cutting tube diameters, 10 items per package		2	Diameter - 32,4 Height - 13 Edge- 25	Silicone 70 Shore A 84
24	5166nou Membrane seal for all cutting tube diameters, 10 items per package		1	Outside diameter - 28 Inside diameter - 5 Height - 4	Silicone 70 Shore A

25	5177nou Seal caps for trocar cannula Ø 12 mm and 15 mm, 10 items per package		2	Diameter - 20 Height - 8	Polytetrafluorethylene PTFE Turcon MF6
26	5180nou O-ring for trocar cannula Ø 12 mm and 15 mm, 10 items per package		1	Diameter - 32 Thickness - 1	Polytetrafluorethylene PTFE Turcon MF6
27	51484nou Seal cap holder for trocar cannula Ø 12 mm, 1 item per package		60	Outside diameter - 30 Inside diameter - 12 Height - 29 Thread M24	Medical stainless steel 1.4305

28	51502nou Seal cap holder for trocar cannula Ø 15 mm, 1 item per package		60	Outside diameter – 30 Inside diameter - 15 Height - 28 Thread M24 Medical stainless steel 1.4305
----	--	---	----	---

Date: 01.12.2016

Leszek Kazana /Area Sales Manager

CERTIFICATE OF ANALYSIS

MORCELLATOR

REVISION:2

STERILITY VALIDATION ACCORD. ISO 17664:2004

Product Codes: Morcellator
Manufacturer: NOUVAG AG

According to AAMI TIR12:2004 /ISO 17665-1:2006/ 11737-2:2009 Morcellators were tested for efficacy of the sponsors sterilization process (31932 41/10) for reuse (Zwisler report 1010.0717). The test specimen were analysed reassembled. For the sterilisation process #1 and #2, two test specimen each were contaminated with spore stripes of *Geobacillus stearothermophilus* (with worst case by placement of 10 biological indicators on each device). For the sterilization process #3, three test specimen were contaminated with a challenge suspension of *Geobacillus stearothermophilus* spores and Soybean-casein digest medium (with worst case by placement of 10 biological indicators on each device). This microbial challenge was let dry for approx. one hour. Thereafter, the test samples were sterilized in a partial cycle (2.5 min, 134°C) of the recommended sterilizing instructions, whereas reference test samples were kept unclean and the bacterial load was determined by bioburden.


Recommended Sterilization Process

- reassembled test specimen were double packed and sealed in foil for sterilization
- sterilization parameters (moisted heat, type B autoclave): two pre-vacuum steps, holding time 5 minutes at 134°C and drying step for ten minutes (validated after: partial cycle, holding time 2.5 minutes at 134°C and drying step for ten minutes)

After the entire procedure the test specimen were tested for sterility by direct immersion of the samples for 14 days in Soybean-casein digest medium. The test was replicated three times concerning sterilization efficiency (product code: Morcellator). In all cases the sterilization process applied resulted in a reduction of more than 10^6 bacteria per device.

Date: 01.12.2016

Leszek Kazana /Area Sales Manager

The sterilizing procedure applied is adequate for the tested medical device. 

Dr. Christian Draing, Study Director

Date: 05.11.2010

Christine Sauter, Quality Assurance Personnel

Date: 5.11.10

EMC - R e p o r t

Emission and Immunity Tests

Test Item :

Control unit for electronic motor

TCM 3000 BL Morcellator

Manufacturer:

NOUVAG AG

Date: 01.12.2016

Leszek Kazana /Area Sales Manager

L. Kazana

Without permission of the test centre this report is not permitted to be duplicated in extract.

Report No.: PLE130919

Test Item: Control unit for electronic motor

Type or Model: TCM 3000 BL Morcellator

Applicant: NOUVAG AG
St. Gallerstraße 23-25
CH-9403 Goldach

Arrival Date: 25. September 2013

Place of Testing: PRO EMV Labor Strausberg GmbH
Garzauer Chaussee
15344 Strausberg
e-mail: proemv@proemv.de

Phone: +493341 - 335255
Fax: +493341 - 335368
www.proemv.de

Standards: EN 60601-1-2 /07.2007

Date of Testing: 25. September to 09. October 2013

Procedure: Section 6 (Electromagnetic compatibility)

Test Result: The device meets the requirements for class A.

Tested by: A. Koch
Lab. Manager

2013-10-09
Date, Signature

Inspected : R. Erxleben
General Manager

2013-10-09
Date, Signature

Contents

1. Instruments And Equipment Used	4
2. Equipment under Test (EUT) Description	6
3. Test Set-up and Mode of Operation During the Tests	7
4. Test Program / Summary and Test Results	8
5. Performance Criteria	9
5.1. Description of EUT quality of performance	9
6. Measurements	10
6.1. Conducted Continuous Disturbance at Mains Port	10
6.1.1. Standards	10
6.1.2. Test Description	10
6.1.3. Limits	10
6.1.4. Measured Values	10
6.1.5. Test Result	10
6.2. Radiated Disturbance	11
6.2.1. Standards	11
6.2.2. Test Description	11
6.2.3. Limits	11
6.2.4. Set-up and Operating Conditions During the Measurements	11
6.2.5. Measured Values	11
6.2.6. Test Result	12
6.3. Immunity to Radiated Electromagnetic Fields	13
6.3.1. Standards	13
6.3.2. Set-up and Operating Conditions During the Measurements	13
6.3.3. Test Description	13
6.3.4. Test Result	13
6.4. Immunity to Conducted Disturbance by RF fields, injected Current	14
6.4.1. Standards	14
6.4.2. Test Description	14
6.4.3. Set-up and Operating Conditions During the Measurements	14
6.4.4. Test Result	14
6.5. Immunity to Fast Transient Disturbances (Bursts)	15
6.5.1. Standards	15
6.5.2. Test Description	15
6.5.3. Set-up and Operating Conditions During the Measurements	15
6.5.4. Test Result	15
6.6. Immunity to Electrostatic Discharge (ESD)	16
6.6.1. Standard	16
6.6.2. Test Description	16
6.6.3. Test Result	16
6.7. Immunity to Surge Pulses	17
6.7.1. Standards	17
6.7.2. Test Description	17
6.7.3. Measurements	17
6.7.4. Test Result	17
6.8. Voltage dips, interruptions and fluctuations immunity tests	18
6.8.1. Standards	18
6.8.2. Test levels	18
6.8.3. Test Result	18
6.9. Power frequency magnetic field immunity test	19
6.9.1. Standards	19
6.9.2. Test Description	19
6.9.3. Set-up and Operating Conditions During the Measurements	19
6.9.4. Test Result	19

6.10. Immunity to variations of power frequency	20
6.10.1. Standard	20
6.10.2. Test Description	20
6.10.3. Set-up and Operating Conditions During the Measurements	20
6.10.4. Test Result	20
7. Photo documentation	29

Other Aspects: Diagrams

Page 21 to 28

1. Instruments And Equipment Used

Conducted disturbance:

	Type	Make	Ser. / Inv. No.
EMI test receiver	ESIB 26	Rohde&Schwarz	SN 100172
EMI test software	ES-K1 V1.60	Rohde&Schwarz	
Network simulation	ESH2-Z5	Rohde&Schwarz	SN 832769/008

Radiated disturbance:

	Type	Make	Ser. / Inv. No.
EMI test receiver	ESIB 26	Rohde&Schwarz	SN 100172
EMI test software	ES-K1 V1.60	Rohde&Schwarz	
Bilog antenna	CBL 6111	Chase	SN 1134

Immunity test (RF field):

	Type	Make	Ser. / Inv. No.
Power amplifier	KAW3020	amplifier research	SN 10479-1
Power amplifier	30S1G3	amplifier research	SN 30778
Signal generator	SML03	Rohde&Schwarz	SN 100097
Test software (RF field)	OATS-sX V6.0.9	CONFORMITAS	
Power meter	NRVD	Rohde&Schwarz	SN 836519/011
Direction coupler	DC-6180	amplifier research	SN 14375
Thermal power head	NRV-Z 51	Rohde&Schwarz	SN 834519/023
Field strength meter	PMM 8051	PMM	SN 0106
Probe	BA 01	PMM	SN 059
Broadband antenna	BTA-L	Frankonia	SN 97061002
Horn antenna	BBHA 9120 E	Schwarzbeck	SN 0899

Immunity test (Conducted disturbance by RF field, current injection):

	Type	Make	Ser. / Inv. No.
Power amplifier	50W1000A	amplifier research	SN 14506
Signal generator	SMY01	Rohde&Schwarz	SN 833913/006
Test software (current injection)	OATS-sX V6.0.9	CONFORMITAS	
Power meter	NRVD	Rohde&Schwarz	SN 836519/011
HF-Voltage head	URV5-Z7	Rohde&Schwarz	SN 860620
Coupling networks	CDN801	MEB	
Coupling clamp	KEMZ 801	MEB	

Immunity test (Burst pulses):

	Type	Make	Ser. / Inv. No.
Burst generator	NSG 3025	Schaffner	SN 19431
Software Burst Test	NSG 3025 V02.01	Schaffner	

Immunity test (ESD):

	Type	Make	Ser. / Inv. No.
ESD simulator	NSG 435	Schaffner	SN 00000599

Immunity test (Surge pulses):

	Type	Make	Ser. / Inv. No.
Surge simulator	NSG 651	Schaffner	SN 00000221
Test software (Surge)	NSG 651 V02.05	Schaffner	
Coupling network	CDN 110	Schaffner	SN 00000309

Immunity test (Voltage dips, interruptions and fluctuations):

	Type	Make	Ser. / Inv. No.
Simulator	NSG 1003	Schaffner	SN 00000182
Software (Voltage fluctuations)	NSG1003 V1.22	Schaffner	

Immunity magnetic LF fields:

	Type	Make	Ser. / Inv. No.
AC-Source	ACS 500	EMTEST	SN 0804-01
Transformer	T 500	peco	SN 6018
Helmholtz coils	HHS1	PRO EMV GmbH	IN 007/96

Immunity to variations of power frequency:

	Type	Make	Ser. / Inv. No.
AC-Source	ACS 500	EMTEST	SN 0804-01

Measuring equipment:

	Type	Make	Ser. / Inv. No.
Absorbing chamber	8,4 x 7,2 x 5,3 [m]	Frankonia	
Antenna mast	MA 240	Deisel	
Rotary disk	DS 415	Deisel	
Controller	HD 100	Deisel	

2. Equipment under Test (EUT) Description

The EUT is a microprocessor-controlled motor system for aesthetic-plastic surgery and dermatology.

The test item is designed for use in medical areas (hospital, clinic, doctor's office).

Classification of the EUT into group and class based on CISPR 11:

The TCM 3000 BL Morcellator is an equipment of group1, class A.

Class A equipment is equipment suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

Class B equipment is equipment suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes

Technical Data:

Mains voltage (V)	: 100 / 115 / 230 V AC
Mains frequency (Hz)	: 50 - 60
Power consumption	: 60 VA
Protecting class	: II
Dimensions (cm)	: 115 x 120 x 180
Case design	: Plastic
Interfaces	: • Power supply "In" • Foot switch (Footcontrol Vario, NOUVAG AG, SN: 2924F1301R) • Electronic motor (NOUVAG AG, SN: 4330M1208R)
Hardware release	09/2013
Serial number	2906U1301R

For any further information see the applicant's system documentation.

3. Test Set-up and Mode of Operation During the Tests

The measurements were performed in a typical test configuration providing a maximum interference capability.

The test set up was done according to the named EMC base standards.

Mode of operation during the tests	Normal operating conditions
Mains voltage at tests	• 230 V AC / 50 Hz • 100 V AC / 60 (50) Hz
Adjustments at the EUT	• Prog. High 500 - 1000 RPM -> 1000 RPM
Connected lines, Length (meters)	• Power supply "In", unshielded (2,5) • Foot switch, shielded (2,8) • Electronic motor, shielded (2,8)
Indicators for the evaluation of the EUT behaviour	• Display • Motor driving

The test item was tested as a table top equipment.

Climatic conditions:

	Required	Actual
Ambient temperature	15 to 35°C	✓
Relative humidity	30 to 60 %	✓
Atmospheric pressure	86 to 106 kPa	✓

Unless otherwise noted these conditions are valid for all following measurements.

Details of device settings and test arrangements may also be seen in the photo documentation.

4. Test Program / Summary and Test Results

The device presented for testing was tested for compliance with the mentioned standards. The following tests (measurements) have been performed on the test item:

Test	Test level	Result
Conducted continuous disturbance at main ports	class B	passed
Radiated disturbance	class A*	passed
Harmonic currents in low-voltage supply systems	not applicable ! (see EN 61000-3-2, section 7, power consumption < 75 W)	
Voltage fluctuations and flicker in low-voltage supply systems	not applicable ! (see EN 61000-3-3, section 6; no significant flicker)	
Immunity to radiated electromagnetic fields	3 V/m	passed
Conducted immunity, injected currents	3 V	passed
Immunity to fast transient disturbances	1 / 2 kV	passed
Immunity to electrostatic discharge	6 / 8 kV	passed
Immunity to surge pulses	1 / 2 kV	passed
Voltage dips, interruptions and fluctuations		passed
Immunity magnetic LF fields	3 A/m	passed
Immunity to variations of power frequency	49 / 51 59 / 61	passed

*) Remark to emission tests:

When a medical equipment of the class A is provided for the use in the domestic establishment or to the connection to the public low-voltage power network, the instructions must contain the following one or an equally significant warning:

Warning!

This device/system is intended solely for use by medical specialists.
This device/system may cause radio disturbances or can disturb the operation of devices nearby. It may be necessary to take appropriate remedial measures, such as e.g. a new orientation, a new arrangement of the control unit or the shielding of the equipment.

5. Performance Criteria

Performance criteria corresponding with EN 60601-1-2:

Under the test the equipment (EUT) shall be able to provide the intended clinical benefit and remain safe. The EUT may exhibit degradation of performance (e.g. deviation from manufacturer's specification). However, the following degradations shall not be allowed.

- component failures
- changes in programmable parameters
- reset to factory defaults (manufacturer's presets)
- change of operating mode
- false alarms
- cessation of any intended operation, even if accompanied by an alarm
- initiation of any unintended operation, including unintended or uncontrolled motion, even if accompanied by an alarm
- error of a displayed numerical value sufficiently large to affect diagnosis or treatment
- noise on a waveform in which the noise is indistinguishable from physiologically - produced signals or the noise interferes with interpretation of physiologically - produced signals
- artefact or distortion in an image in which the artefact is indistinguishable from physiologically - produced signals or the distortion interferes with interpretation of physiologically - produced signals
- failure of automatic diagnosis or treatment EUT to diagnose or treat, even if accompanied by an alarm

5.1. Description of EUT quality of performance

The manufacturer has designated performance criterion as specified by EN 60601-1-2:2007.

Following reactions of the EUT are according to manufacturer's specifications permissible:

- Immunity test "Voltage dips, interruptions and fluctuations":
A shutdown of the test item is permissible.

6. Measurements

Note:

The test results shall apply exclusively to the device under test. They shall not represent a generally valid opinion on the properties of the respective products from the running production process.

6.1. Conducted Continuous Disturbance at Mains Port

6.1.1. Standards

EN 60601-1-2 /07.2007

CISPR 11 /2003 + A1 /2004, mod. + A2 /2006

group 1

class B

6.1.2. Test Description

The conducted disturbance at the mains port was tested. It was measured in the frequency range from 0.15 MHz to 30 MHz.

6.1.3. Limits

The limiting values are:

Frequency range	Limit class B (Detector)
0,15 - 0,5 MHz	66 - 56 dB μ V (QP) 56 - 46 dB μ V (AV)
0,5 - 5 MHz	56 dB μ V (QP) 46 dB μ V (AV)
5 - 30 MHz	60 dB μ V (QP) 50 dB μ V (AV)

Remark: for example 66 to 56 dB μ V means, the limiting value falls linearly with the logarithm of the frequency.

6.1.4. Measured Values

Input voltage (U_{in}): 230 V AC / 50 Hz

Measurement	Diagram	Phase	Detector	Remark
1	page 21, 23	L1/N	PK/QP	Test ok
2	page 22, 23	L1/N	AV/AV	Test ok

6.1.5. Test Result

The test item meets the requirements.

Note: No exceeding of the limits was observed.

6.2. Radiated Disturbance

6.2.1. Standards

EN 60601-1-2 /07.2007
CISPR 11 /2003 + A1 /2004, mod. + A2 /2006 group 1 class A

6.2.2. Test Description

The radiated disturbance was tested. It was measured in the frequency range from 30 MHz to 1000 MHz at a measuring distance of 10 m.

6.2.3. Limits

The QP limiting values are:

Frequency range	Limit class A
30 - 230 MHz	40 dB μ V/m
230 - 1000 MHz	47 dB μ V/m

6.2.4. Set-up and Operating Conditions During the Measurements

The disturbance levels indicated in the measurement diagrams (measured in an absorbing chamber) shall be reduced by approximately 10 dB before comparing them with the limit (conversion of the measuring distance from 3 m to 10 m). The correction values: test assembly attenuation, antenna factor, and cable attenuation, are taken into consideration in the measurement result.

The measurement diagrams with QP verification (decisive for the rating) represent the maximum which could be reached by turning the test item, by variations in the height of the antenna (1 - 4 m), and by changing the antenna polarisation.

Input voltage (U_{in}): 230 V AC / 50 Hz

6.2.5. Measured Values

6.2.5.1. Measured values of the pre-measurements in the absorbing chamber

Antenna height 2.0 m; horizontal and vertical polarization; antenna-to-EUT azimuth 0, 90, 180, 270 [dg]; MaxHold scan

Measurement	Diagram	Detector	Remark
1	page 24	PK	horizontal antenna
2	page 25	PK	vertical antenna

6.2.5.2. Measured values of the final measurement

Measurement	Diagram	Detector	Remark
3	page 26, 28	QP	horizontal antenna, Test ok
4	page 27, 28	QP	vertical antenna, Test ok

Remark:

Due to the radio interference of the environment a final measurement was not possible on the open area test site in 10 m and 3 m measuring distance. Therefore it was measured according to the current edition CISPR 11 /2009, mod. + A1 /2010 on a test setup according to the CISPR 16-1-4 with a measuring distance of 3 m and up to 10 dB increased limiting value.

6.2.6. Test Result

The test item meets the requirements for class A.

Note: No exceeding of the limits was observed.

6.3. Immunity to Radiated Electromagnetic Fields

6.3.1. Standards

EN 60601-1-2 /07.2007

IEC 61000-4-3

6.3.2. Set-up and Operating Conditions During the Measurements

The EUT was placed on a wooden rotation desk so that it was positioned in a 0.8 m distance from the chamber's floor and in a 3 m distance in front of the antenna.

The EUT was exposed to radiation from the front, from the right side, from the left side, and from behind.

The frequency range was passed through once in each polarisation.

Input voltage (U_{in}): 230 V AC / 50 Hz

6.3.3. Test Description

Basic standard	IEC 61000-4-3
Frequency range	80 - 2500 MHz
Step size	1 % of previous step
Modulation	AM 80%, 1 kHz sinus
Field strength	3 V/m
Polarisation	horizontal and vertical
Performance criteria	see section 5 / 5.1
Dwell time	2 s for each step

6.3.4. Test Result

Test item reactions: No unacceptable loss of performance or loss of data was observed.

During and after the influence of the disturbance the test item fulfilled its normal function properly.

6.4. Immunity to Conducted Disturbance by RF fields, injected Current

6.4.1. Standards

EN 60601-1-2 /07.2007

IEC 61000-4-6

6.4.2. Test Description

Basic standard	IEC 61000-4-6
Frequency range	0,15 - 80 MHz
Frequency steps	1%
Modulation	AM 80%, 1 kHz sinus
RF-voltage	3 V (unmodulated wave) *
Injected Signal lines	3 V (modulated wave)
Injected DC lines	3 V (modulated wave)
Injected AC-lines	3 V (modulated wave)
Rating criteria	see section 5 / 5.1
Dwell time	2 s for each frequency step

*) Deviating to IEC 61000-4-6 the test level was calibrated according to
EN 60601-1-2 within: -0% to +25%

6.4.3. Set-up and Operating Conditions During the Measurements

For the tests the following connected lines were led through coupling networks (CDN).

Connected Line	Coupling network (0,15 - 80 MHz)
230 V AC-power supply	CDN 801-M3
Footswitch	KEMZ 801
Electronic motor	KEMZ 801

The coupling network were separately connected to a RF generator for each test run.

6.4.4. Test Result

Test item reactions: No unacceptable loss of performance or loss of data was observed.

During and after the influence of the disturbance the test item fulfilled its normal function properly.

6.5. Immunity to Fast Transient Disturbances (Bursts)

6.5.1. Standards

EN 60601-1-2 /07.2007

IEC 61000-4-4

6.5.2. Test Description

Basic standard	IEC 61000-4-4
Pulse shape	5 / 50 ns
Repetition frequency	5 kHz
Test duration	> 1 min
Signal lines * (cap. coupling pliers)	1,0 kV (0,5 kV start, 0,5 kV steps to 1,0 kV)
DC lines (direct coupling)	2,0 kV (0,5 kV start, 0,5 kV steps to 2,0 kV)
AC-lines (direct coupling)	2,0 kV (0,5 kV start, 0,5 kV steps to 2,0 kV)
Rating criteria	see section 5 / 5.1

*) Applicable only to ports interfacing with cables whose total length according to the manufacturers functional specification may exceed 3 m.

Input voltage (U_{in}): 100 V AC / 50 Hz and 230 V AC / 50 Hz

6.5.3. Set-up and Operating Conditions During the Measurements

The burst pulses were coupled into the following lines:

Test voltage [kV]	Input coupling
+/- 2,0	AC-power supply (L, N, PE)

The test voltage was increased in steps from the lowest up to the selected test level.

6.5.4. Test Result

Test item reactions:

Test level [kV]	Remark
+/- 2,0	The content of the display moves towards to the right and to the left.

During and after the influence of the disturbance the test item fulfilled its normal function properly.

6.6. Immunity to Electrostatic Discharge (ESD)

6.6.1. Standard

EN 60601-1-2 /07.2007

IEC 61000-4-2

6.6.2. Test Description

Basic standard	IEC 61000-4-2
Air discharge	8 kV (2 kV start, 2 kV steps to 8 kV)
Contact discharge	6 kV (2 kV start, 2 kV steps to 6 kV)
Discharges per pole	> 10
Discharge-R	330 R
Discharge-C	150 pF
Performance criteria	see section 5 / 5.1

Test points:	C	A	Test level
horizontal / vertical coupling area	X		a.l.
Housing (Plastic parts)		X	a.l.
Primary power receptacle	X	X	a.l.
Power switch		X	a.l.
Connector housing of the Interfaces (Front)	X		a.l.
Cables of the Interfaces		X	a.l.
Display		X	a.l.
Control panel		X	a.l.
Footswitch (Metal parts)	X		a.l.
Footswitch (Plastic parts)		X	a.l.
Electronic motor housing	X		a.l.

C: Contact discharge

A: Air discharge

a.l.: All test levels

Input voltage (U_{in}): 230 V AC / 50 Hz

6.6.3. Test Result

Test item reactions: The display was influenced temporally through a discharge.

No unacceptable loss of performance or loss of data was observed.

During and after the influence of the disturbance the test item fulfilled his normal function properly.

6.7. Immunity to Surge Pulses

6.7.1. Standards

EN 60601-1-2 /07.2007

IEC 61000-4-5

6.7.2. Test Description

Basic standard	IEC 61000-4-5
Open circuit pulse shape	1,2 μ s / 50 μ s
Short circuit pulse shape	8 μ s / 20 μ s
Number of pulses	5 pos. and 5 neg.
Repetition pulse	all 60 seconds pulse
Angles	0°, 90°, 270°
AC lines	1 kV sym. 2 kV asym.
Performance criteria	see section 5 / 5.1

6.7.3. Measurements

Input voltage (U_{in}). 100 V AC / 50 Hz and 230 V AC / 50 Hz

AC power port:

sym. (S) = differential mode :L against N

$R_i = 2 \text{ Ohm}$ and $C = 18 \mu\text{F}$ (LZ)

asym. (A) = common mode :L, N against PE

$R_i = 12 \text{ Ohm}$ and $C = 9 \mu\text{F}$ (HZ)

Number of pulses	Form	Pol.	Up	Angle	Surge Tot.	Result
5 x (L against N)	Surge LZ	+	1000	S 0	5	passed
5 x (L against N)	Surge LZ	-	1000	S 0	10	passed
5 x (L against N)	Surge LZ	+	1000	S 90	15	passed
5 x (L against N)	Surge LZ	-	1000	S 90	20	passed
5 x (L against N)	Surge LZ	+	1000	S 270	25	passed
5 x (L against N)	Surge LZ	-	1000	S 270	30	passed
5 x (L against PE)	Surge HZ	+	2000	A ---	35	passed
5 x (L against PE)	Surge HZ	-	2000	A ---	40	passed
5 x (N against PE)	Surge HZ	+	2000	A ---	45	passed
5 x (N against PE)	Surge HZ	-	2000	A ---	50	passed

The test voltage was increased in steps from the lowest up to the selected test level.

6.7.4. Test Result

Test item reactions: No unacceptable loss of performance or loss of data was observed.

During and after the influence of the disturbance the test item fulfilled its normal function properly.

6.8. Voltage dips, interruptions and fluctuations immunity tests

6.8.1. Standards

EN 60601-1-2 /07.2007

IEC 61000-4-11

6.8.2. Test levels

Input voltage: $U_{in} = 230 \text{ V AC} / 50 \text{ Hz}$
 $U_{in} = 100 \text{ V AC} / 50 \text{ Hz}$

Basic standard	IEC 61000-4- 11			
	dips	dips	interruptions	interruptions
Input voltage U_{in}	U_{in}	U_{in}	U_{in}	U_{in}
Variable voltage U_{var}	$U_{in} -30 \%$	$U_{in} -60 \%$	0 V (-100 %)	0 V (-100%)
Duration of the events t_{event}	500 ms	100 ms	10 ms	5000 ms
Repetition time t_{rep}	20 s	20 s	20 s	20 s
Phase angle t_{phase} [grd]	0, 180	0, 180	0, 180	0, 180
Test duration	3 x per phase position	3 x per phase position	3 x per phase position	3 x per phase position
Performance criteria	see section 5 / 5.1	see section 5 / 5.1	see section 5 / 5.1	see section 5 / 5.1

6.8.3. Test Result

The test item meets the stipulated requirements.

Test item reactions 230 V AC:

Test	Parameter U_{var} / t_{event}	Remarks
interruptions	0 V / 10 ms	During and after the test the EUT functioned as intended.
dips	160 V / 500 ms	
dips	90 V / 100 ms	
interruptions	0 V / 5000 ms	EUT shuts down. After the disturbance the normal operation can be recovered by off- and on-switching of the equipment.

Test item reactions 100 V AC:

Test	Parameter U_{var} / t_{event}	Remarks
interruptions	0 V / 10 ms	During and after the test the EUT functioned as intended.
dips	70 V / 500 ms	
dips	40 V / 100 ms	
interruptions	0 V / 5000 ms	EUT shuts down. After the disturbance the normal operation can be recovered by off- and on-switching of the equipment.

6.9. Power frequency magnetic field immunity test (immunity magnetic LF fields)

6.9.1. Standards

EN 60601-1-2 /07.2007

IEC 61000-4-8

6.9.2. Test Description

Basic standard	IEC 61000-4-8
Mag. field strength	3 A/m
Frequency	50 / 60 Hz
Test duration	> 1 min (every direction)
Rating criteria	see section 5 / 5.1

6.9.3. Set-up and Operating Conditions During the Measurements

The EUT was arranged in the focal point of a Helmholtz coils and exposed to the magnetic field.

The EUT was rotated by 90° in order to expose it to the test field with different orientations

Input voltage (U_{in}): 230 V AC / 50 Hz

6.9.4. Test Result

Test item reactions: No unacceptable loss of performance or loss of data was observed.

During and after the influence of the disturbance the test item fulfilled its normal function properly.

6.10. Immunity to variations of power frequency

6.10.1. Standard

EN 60601-1-2 /07.2007

IEC 60601-1

6.10.2. Test Description

Basic standard	IEC 60601-1	IEC 60601-1
Nominal voltage	230 V	100 V
Nominal frequency	50 Hz	60 Hz
Test parameter	± 1 Hz of nominal frequency	± 1 Hz of nominal frequency
Test frequency	49 Hz (50 Hz - 1 Hz)	59 Hz (60 Hz - 1 Hz)
Test frequency	51 Hz (50 Hz + 1 Hz)	61 Hz (60 Hz + 1 Hz)
Duration of test	> 1 min	> 1 min
Rating criteria	see section 5 / 5.1	see section 5 / 5.1

6.10.3. Set-up and Operating Conditions During the Measurements

Nominal voltage: 100 V AC and 230 V AC

The test was carried out at the nominal voltage and $\pm 10\%$ of the nominal voltage.

6.10.4. Test Result

Test item reactions: No unacceptable loss of performance or loss of data was observed.

During and after the influence of the disturbance the test item fulfilled its normal function properly.

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


40.000000 42.51
120.000000 36.76
160.050000 34.12

7. Photo documentation


Figure 1 : Test Set-up "Conducted continuous disturbance at mains port"

Figure 2 : Test Set-up "Radiated electromagnetic fields and immunity to radiated electromagnetic fields"


Figure 3 : Test Set-up "Radiated electromagnetic fields and immunity to radiated electromagnetic fields"


Figure 4 : Test Set-up "Radiated electromagnetic fields and immunity to radiated electromagnetic fields"


Figure 5 : Test Set-up "Immunity to variations of power frequency"


Figure 6 : Test Set-up "Immunity to conducted disturbance by RF field, current injection"


Figure 7 : Test Set-up "Immunity to fast transient disturbances (Bursts)"


Figure 8 : Test Set-up "Immunity to electrostatic discharge (ESD)"


Figure 9 : Test Set-up "Immunity to surge pulses (surge)"


Figure 10 : Test Set-up "Immunity to voltage dips, interruptions and fluctuations"




Figure 11 : Test Set-up "Immunity magnetic LF fields"

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