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
Fig 1. SES 2000 Autoclave

INTRODUCTION

There are two versions of the SES2000 Autoclave, one with a printer and one without a printer as shown on the front cover of this publication.

DO'S & DON'TS

DO'S

- DO close door properly prior to selecting programme.
- DO leave chamber door slightly open after use.
- DO use distilled or deionized water, or water treated by reverse osmosis for filling the reservoir.
- DO ensure instruments are free from debris before placing them in sterilizer chamber.
- DO use the 'tray lifter' to remove sterilized load.
- DO use pouch rack for sterilization of pouches or bags and place only on the second from bottom shelf.
- DO keep chamber trays and chamber face clean.
- DO clean door seal and chamber face with a lint free cloth at weekly intervals.
- DO switch off unit and disconnect from mains supply before checking or renewing a fuse, or cleaning the unit.

DON'TS

- DON'T use trays or cassettes without perforations as this will result in a non-sterile load.
- DON'T process wrapped, sealed or porous goods or fluids in this autoclave.
- DON'T use the unit in zones of risk associated with flammable anaesthetics.
- DON'T cover the ventilation louvres.
- DON'T fill reservoir during a sterilizing cycle.
- DON'T use tap water in reservoir.
- DON'T use abrasive powders, chemicals or concentrated chloride solutions for cleaning purposes.

GENERAL

1. The SES 2000 is an autoclave designed for the sterilization of unwrapped instruments, utensils and other items. Alternatively, such items can be packed in L.M.G. SMITH BROTHERS "VIEW-PACK SELF-SEAL" POUCHES, but only using the special pouch rack which is available as an optional extra. The SES 2000 will also take cassettes using a cassette rack.

CAUTION:

- (1) Never use trays or cassettes without perforations as this will result in a non-sterile load.
- (2) DO NOT PROCESS WRAPPED, SEALED OR POROUS GOODS OR FLUIDS IN THIS AUTOCLAVE.
- (3) NOT TO BE USED IN ZONES OF RISK ASSOCIATED WITH FLAMMABLE ANAESTHETICS.

WARNING

Sterilization of equipment
not of solid metal construction.

The manufacturers of such equipment should be consulted about suitability for autoclaving and the maximum temperature which the equipment can withstand, as the drying phase could exceed sterilizing temperatures under certain conditions.

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programme selector touch button and has four programmes, as follows:


- 134°C Without drying
 - 121°C Without drying
 - 134°C With drying
 - 121°C With drying
3. The SES 2000 is capable of 4 cycles per hour when the '134°C Without drying' programme is selected. For sterilizing cycle times, refer to 'TECHNICAL DATA'.

Safety Features

4. The SES 2000 is provided with the following safety features:
- (1) Pressure Door Lock. The door is mechanically locked at pressures above 0.14 bar.
 - (2) Door Interlock Switch. This switch prevents programme starting if door is not properly closed.
 - (3) Pressure Relief Valve (safety valve). This valve safely releases excess pressure (i.e. above 2.6 bar).
 - (4) Microcomputer. The microcomputer constantly monitors all key functions. If an error arises, it will immediately stop the cycle, discharge pressure and cause the appropriate message to be displayed (see ERROR CORRECTION para. 31).
 - (5) Overheat control. The microcomputer operates in conjunction with an automatic reset thermostat to prevent the heating element from overheating. The autoclave is fitted with a fusible link which, should overheating occur, will cut all power to the heating element and computer, leaving only the red overheat warning lamp illuminated on the front panel.

Unpacking and Installation

5. Remove the SES 2000 from its packing case by lifting with the hands placed under the base at each side and with the chamber door facing towards you. Place the unit on a flat, level surface and ensure there is a working clearance of 150mm all round for adequate ventilation. **NEVER COVER THE VENTILATION LOUVRES.**
6. Open chamber door Fig.2. Remove sterilizing trays, tray lifter, reservoir lid, drip tray and feet adjustment pack. Discard all packing pieces.


7. Fit reservoir cover and place plastic drip tray at front of unit, below door. Fit feet adjustment pack components as follows:
- (1) Plastic cups to protect work surface from rubber feet marks.
 - (2) If work surface is not level, place spacers, as required, inside cup(s).
8. Leave chamber door slightly open when unit is not in use.
9. Fit electrical plug, if not already fitted, to mains supply cable. The UK colour code for plug wiring is:
- Live (line) - Brown
 - Neutral - Blue
 - Earth/ground - Green/Yellow
10. Check that the plug fuse is appropriate for voltage of equipment (see rating plate on back of cabinet).

WARNING
THIS EQUIPMENT MUST BE EARTHED


12. Remove reservoir cover and fill reservoir up to 'MAX' mark with distilled or deionized water, or water treated by reverse osmosis process. **DO NOT USE TAP WATER.** Refit cover.

PREPARATION FOR USE

13. Check that water level in reservoir is between 'MAX' and 'MIN' marks Fig. 3. **DO NOT FILL DURING CYCLE.**
14. Open chamber door, if closed.

15. Select mains on/off switch to 'I' (on). This will initiate a single audible signal followed by the display 'SES', then the number of cycles and finally 'ready' in the display window. If the door is closed, the display will alternate between 'door' and 'test' when power is switched on. In order to continue with cycle, the door must be opened, at which point the display will change to 'ready' together with four green indicators.

16. Load trays and place in sterilizing chamber Fig. 4. Bowls and dishes should be loaded on their sides to ensure total sterilization. Note: See information provided for sterilizing using special pouches/bags.

IMPORTANT

The SES 2000 is designed to sterilize instruments, **NOT TO WASH THEM.**

If instruments are not cleaned this may compromise the efficiency of the sterilization process.

17. Instruments should be pre-washed, preferably in an ultrasonic cleaner to remove amalgam, debris, etc., then rinsed to remove all traces of proprietary cleaners, chemicals or disinfectants, before being placed in sterilizer chamber. Residue from these products could result in a blockage of the water re-cycling system.
18. Ensure door handle is turned to fully open position (counter-clockwise), then close door by turning handle 180° clockwise Fig. 2.

STARTING THE STERILIZING CYCLE

19. Select required programme by pressing appropriate button on control panel Fig. 5. The indicator lamp for the selected programme remains illuminated. The sterilizing cycle will now proceed automatically. As the cycle progresses various displays will appear in the display window to indicate programme status, as follows:

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
DISPLAY SEQUENCE DURING CYCLE

20. Display - Interpretation
- door/test - Test door interlock switch (see para 15)
 - door - Sterilizing chamber door open
 - rEAdy - Programme can now be selected
 - StArt - Programme selected and cycle started
 - FILL - Chamber being filled
 - HEAt - Heater temperature below 92°C
 - 92 to 136 - Heating to sterilizing stage
 - S-135.5 - 'S' flashing, indicates sterilizing begun and timing started
 - cond - Steam being discharged and condensed
 - dry - Load being dried
 - End - Cycle complete
 - Ser - Service set to display at 1500 cycles. This will not inhibit the use of the autoclave. (Arrange for service).
 - door - For access to cycle counter press 121 with drying programme button.

DRYING

21. Commencement of the drying phase is indicated by the display 'dry' accompanied by a rapid intermittent audible signal. When this signal is heard, maximum drying will be achieved if the chamber door is opened halfway and then closed until it just touches the metal catch on the front fascia panel, leaving a small gap (approx. 1cm) for ventilation. Instruments can be taken out any time after this, but best results will be achieved by waiting until the automatic cycle is complete (see UNLOADING, para 22).

UNLOADING

22. On completion of a sterilizing cycle 'without drying', the display 'End' will appear, accompanied by a brief audible signal. When the door is opened to remove load the display will change to 'door'. On completion of a sterilizing cycle 'with drying' (where chamber door is opened slightly during the drying phase), the display 'End' will appear briefly accompanied by an audible signal and followed, after a few seconds, by 'door'. Open door fully and remove load USING THE TRAY LIFTER.

STERILIZING WITH POUCHES/BAGS

23. Sterilization of instruments or other metal items in pouches/bags should only be undertaken using the special pouch rack available as an optional extra. (Fig. 6 below) in conjunction with pouches as specified by the manufacturer (see para 1).


Fig. 6 Pouch rack

more than 150g or six items per pack, whichever is the least.

Following washing or ultrasonic cleaning, a suitable number of instruments should be placed in each pouch, ensuring that the pouch is not overloaded. The sealed pouches should be loaded horizontally between the coils of the rack, with no more than one pouch between adjacent coils. When autoclaving handpieces, they should be placed in an upright position in the pouch.

Note: Paper/film pouches must be loaded such that the film sides of adjacent pouches are face to face.

25. Pouches must only be sterilized using the '134°C with drying' programme, at the end of which the pouches should be dry and can be stored for use at a later date. The pouch rack should only be placed on the second from bottom shelf. Note: See 'Warning' in para 1.

CLEANING AND CARE

Cleaning

26. Keep chamber trays and chamber face clean. Door seal and chamber face should be cleaned using a lint-free cloth only.

CAUTION

Do not use abrasive powders, chemicals or concentrated chloride solutions

27. The reservoir must be drained at two-weekly intervals (Fig. 7) and recharged with distilled or deionized water, or water treated by reverse osmosis process. **DO NOT USE TAP WATER.** Regular cleaning of reservoir will reduce the effects of excessive handpiece lubricant which could be detrimental to the function of the unit.


Fig. 7 Drain tap

Care

28. Attention to the following will increase the life of your SES 2000:
- (1) After use leave chamber door slightly open.
 - (2) At weekly intervals lightly clean door seal and chamber face with a lint-free cloth.

FUSE RENEWAL

WARNING

Switch off unit and disconnect mains power supply before renewing a fuse.

29. Fuses are to be found at the back of the unit (see Fig. 1 detail). For fuse ratings refer to **TECHNICAL DATA**. To extract a fuse, insert a screwdriver or small coin in the slot of the fuse holder and twist it counter-clockwise. Having inspected/renewed a fuse, reverse the above procedure to re-secure fuse holder.

PLANNED PREVENTIVE MAINTENANCE AGREEMENT

30. We would strongly recommend that a 'P.P.M. Agreement' is taken out on your SES 2000. Although the unit requires minimal maintenance, it is good practice to have the autoclave checked and calibrated at regular intervals. This ensures that the exacting conditions necessary for sterilization are maintained throughout the unit's working life.

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31. If after switching on power there is no visual or audible display, first check local power supply connections. Also check all 3 fuses at rear of unit (see fuse renewal). If any error occurs during a cycle (i.e. any time after selecting programme), the cycle will abort and provided power supply

to the unit is maintained the error will be indicated by one of the visual displays shown below. If power fails during a cycle, check supply conditions, fuses and connections as already indicated above.

WARNING

If an error display appears during a cycle, do not switch-off power until any discharge of hot water or steam into reservoir has stopped. Also do not attempt to open chamber door until unit has cooled and internal pressure has fallen sufficiently to release pressure safety bolt.

DISPLAY

CAUSE

REMEDY

ELECx (+ audible signal)	Temporary failure in mains supply to unit.	(a) Check local supply conditions. Check supply plug wiring and power cable for possible breaks. (b) Re-start cycle (see NOTE below). NOTE: Because there can be no guarantee that the load has been sterilized, 'ELECx' will continue to be displayed after power has been restored. To re-start cycle: (a) Wait until discharge of water or steam into reservoir has ceased. (b) Switch off power and wait approximately 10 seconds. (c) Push and hold-in any one of the four programme selector buttons and, at the same time, switch on power. (d) Release programme selector button as soon as 'SES' appears. (e) Select required programme to restart cycle.
LoH20 (+ audible signal)	Water failed to enter chamber from reservoir.	DO NOT OPEN DOOR, for safety reasons, switch off power and check water level in reservoir is above 'MIN' mark. DO NOT FILL TO 'MAX' MARK during a cycle, if reservoir already contains water above the 'MIN' mark. Switch on power and re-start cycle.
H20 (+ audible signal)	Water level in chamber has dropped slightly.	(a) Wait for any discharge of water or steam into reservoir to cease, then switch-off power and carry out water level check, as for 'LoH20'. Re-start cycle. (b) If 'H20' display persists, call a Service Engineer.
ERR 2	Door opened before 100°C is reached or Door not fully closed at commencement of cycle or Door switch faulty	Switch off power, wait 5 seconds, switch power on and restart cycle. WARNING: Should door be opened beware of possible hot water leaking from chamber. If ERR 2 persists call Service Engineer.
Door/test	Door closed with power on.	See para 15.
Error	(a) Heater not working (b) Excessive or too low temperatures. (c) Fault with microcomputer system. (d) Low temperature during sterilizing.	(a) Check fuses (see FUSE RENEWAL). Re-start cycle. (b) Ensure unit is level. Allow unit to cool, then re-start cycle. (c) Allow unit to cool and re-start cycle. (d) Re-start cycle. If display 'Error' persists, call Service Engineer.

NOTE: To clear error displays 'LoH20', 'H20' or 'Error', switch-off unit, wait for approximately 10 seconds, then switch-on power again.

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