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# Merit Maestro

## MICROCATHETER

INSTRUCTIONS FOR USE

*[Signature]*  
DAN CLARK  
AUGUST 18, 2016

MERIT MEDICAL SYSTEMS, INC  
1600 WEST MERIT PARKWAY  
SOUTH JORDAN, UTAH 84095  
PHONE 801-253-1600  
FAX 801-253-1688  
www.merit.com


Instructions for Use  
Merit Maestro Microcatheter  
Consisting of 6 pages  
Dated August 18, 2016


### ACKNOWLEDGMENT

State of Utah            )  
                                  ) ss.  
County of Salt Lake    )

On this 18th day of August, 2016, before me, Jennifer Wojciechowski, a notary public, personally appeared Dan Clark, Director of Regulatory Affairs, proved on the basis of satisfactory evidence to be the person whose name is subscribed to the attached instrument, and acknowledge he executed the same.

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English

## INSTRUCTIONS FOR USE

**DESCRIPTION:** The Merit Maestro Microcatheter is available in three French size configurations. The first configuration incorporates a change in its outside diameter along its length from a 2.8F (0.93mm) proximal region to a 2.4F (0.80mm) flexible distal region of 20cm in length. This configuration employs a nominal proximal inner diameter of 0.025" with a nominal tip inner diameter of 0.020".

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**INDICATIONS FOR USE:** The Microcatheter is intended for general intravascular use, including peripheral and coronary vasculature. Once the subselective region has been accessed, the Microcatheter can be used for the controlled and selective infusion of diagnostic, embolic, or therapeutic materials into vessels.

The catheter should not be used in the cerebral vessels.

**CONTRAINDICATIONS:** None known

### WARNINGS:

- Due to contractual agreements, this Microcatheter is not for neurovascular use at or above the common carotid artery or at or above the vertebral artery.
- This device is intended to be used only by physicians trained in percutaneous intravascular techniques and procedures.
- Contents supplied sterile using an ethylene oxide (EtO) process. Do not use if sterile barrier is damaged.
- For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. After use, dispose of product and packaging in accordance with hospital, administrative, and/or local government policy.
- The infusion dynamic pressure with this Microcatheter should not exceed 800psi (5515 kPa). Infusion pressure in excess of this maximum may result in Microcatheter rupture, possibly resulting in patient injury. If flow through the Microcatheter becomes restricted, do not attempt to clear the Microcatheter lumen by infusion. The static pressure with this Microcatheter should not exceed 2068 kPa/300psi. Static pressure in excess of this maximum may result in Microcatheter rupture, possibly resulting in patient injury. Identify and resolve the cause of the blockage or replace the Microcatheter with a new Microcatheter before resuming infusion.
- Make sure that the guiding catheter does not slip out of the vessel. If the guiding catheter should leave the vessel when the Microcatheter and/or the guide wire is moved, this may result in the damage of the Microcatheter system.
- Microcatheter advancement beyond the end of the guide wire may result in vessel trauma.
- Do not use a power injector to infuse agents other than contrast media, as the Microcatheter may become blocked. The safety setting of injection pressure must not exceed the maximum injection pressure of 800psi (5515 kPa). Exceeding injection pressure beyond the maximum injection pressure may cause Microcatheter rupture. (See Instructions For Using a Power Injector)

### PRECAUTIONS:

- Only:** Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.
- Ensure embolic material compatibility with Microcatheter prior to use.
- Always monitor infusion rates when using the Microcatheter.
- When injecting contrast for angiography, ensure that the Microcatheter is not kinked or occluded.
- The Microcatheter has a lubricious hydrophilic coating on the outside of the catheter. It must be kept hydrated prior to removal from its carrier and during the actual procedure in order to be lubricious. This can be accomplished by attaching the Y-connector to a continuous saline drip.
- Prior to a procedure, all equipment to be used for the procedure should be carefully examined to verify proper function and integrity.
- Inspect the Microcatheter prior to use for any bends or kinks. Any Microcatheter damage may decrease the desired performance characteristics.
- Exercise care in handling of the Microcatheter during a procedure to reduce the possibility of accidental breakage, bending or kinking.
- When the Microcatheter is in the body, it should be manipulated only under fluoroscopy. Do not attempt to move the Microcatheter without observing the resultant tip response.
- Exchange Microcatheters frequently during lengthy procedures that require extensive manipulation or multiple guide wire exchanges.
- Never advance or withdraw an intravascular device against resistance until the cause of the resistance is determined by fluoroscopy. Movement of the Microcatheter or guide wire against resistance may result in separation of the Microcatheter or guide wire tip, damage to the Microcatheter, or vessel perforation.
- Because the Microcatheter may be advanced into narrow subselective vasculature, repeatedly assure that the Microcatheter has not been advanced so far as to interfere with its removal.
- Excessive tightening of a hemostatic valve onto the Microcatheter shaft may result in damage to the catheter.
- Read and follow the manufacturer's IFU for diagnostic, embolic, or therapeutic agents to be used with this Microcatheter.
- Do not use opened or damaged packages. Use prior to the "use before" date. Store at controlled room temperature.

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Access site complications  
Vessel perforation  
Vessel Spasm  
Hemorrhage  
Pain and tenderness  
Embolism  
Allergic reaction

Vascular thrombosis  
Thrombosis  
Ischemia  
Infection  
Vessel Dissection  
Distal embolization  
Death

TABLE 1: COMPATIBILITY INFORMATION

| Microcatheter Q.D. | Nominal Microcatheter I.D. | Maximum Guide wire Q.D. | Recommended guiding catheter       |
|--------------------|----------------------------|-------------------------|------------------------------------|
| 2.8F / 2.4F        | 0.020" (0.52mm)            | 0.018" (0.46mm)         | 0.040" (1.02mm) to 0.041" (1.04mm) |
| 2.8F / 2.8F        | 0.024" (0.62mm)            | 0.021" (0.53mm)         | 0.040" (1.02mm) to 0.041" (1.04mm) |
| 2.9F / 2.9F        | 0.027" (0.68mm)            | 0.021" (0.53mm)         | 0.042" (1.07mm) to 0.043" (1.09mm) |
| EMBOLICS           |                            |                         |                                    |
|                    | Particles                  | Spherical               | Coils                              |
| 2.8F / 2.4F        | ≤ 700 µm Emboli            | ≤ 700 µm Microspheres   | 0.46mm / 0.018"                    |
| 2.8F / 2.8F        | ≤ 700 µm Emboli            | ≤ 700 µm Microspheres   | 0.46mm / 0.018"                    |
| 2.9F / 2.9F        | ≤ 900 µm Emboli            | ≤ 900 µm Microspheres   |                                    |

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- Place the appropriate guiding catheter using standard technique. A rotating hemostasis valve may be connected to the guiding catheter luer adapter to continuously flush the guiding catheter with saline.
- Remove the spiral holder housing the Microcatheter from its sealed packaging.
- Attach a syringe filled with heparinized saline solution or sterile water to the luer lock fitting of the Microcatheter holder.
- Inject enough solution to wet the Microcatheter surface entirely. This will activate the hydrophilic coating on the Microcatheter surface. Note: The surface of the Microcatheter may become dry after removal from the holder. Additional wetting with heparinized saline or sterile water will renew the hydrophilic effect.
- Upon removal of the Microcatheter from the spiral holder, inspect Microcatheter to verify there is no damage prior to insertion.
- Attach a second hemostasis valve with side-arm adapter to the Microcatheter, purge any air and flush with heparinized saline or sterile water.
- Carefully insert guide wire into the Microcatheter and completely close the valve around the guide wire.
- Introduce the Microcatheter and guide wire assembly into the guiding catheter via the hemostasis valve (if used). If rotating hemostatic valve is used, tighten the valve around the Microcatheter to prevent backflow, but allowing some movement through the valve by the Microcatheter.
- Using fluoroscopy, introduce the Microcatheter and guide wire assembly into the vascular system, making sure the guide wire is always ahead of the Microcatheter. Advance the guide wire and Microcatheter to a selected vascular site by alternatively advancing the guide wire and then tracking the Microcatheter over the guide wire. Note: To facilitate Microcatheter handling, the proximal portion of the Microcatheter is uncoated to ensure a non-slip grip.
- Final positioning is accomplished by short advances of the guide wire and Microcatheter until the desired position is achieved and then confirmed by fluoroscopic visualization.
- Monitor Microcatheter placement and position during use.
- To infuse, completely remove the guide wire from the Microcatheter. Connect a syringe with infusate to the Microcatheter manifold luer, and infuse as required.

### INSTRUCTION FOR USING A POWER INJECTOR WITH THE MICROCATHETER:

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
## ANNEX TO INSTRUCTIONS FOR USE

### DEVICE SPECIFICATIONS

The Maestro Microcatheter is available in working lengths of 110, 130 and 150 cm with the distal tip offered in pre-shaped straight, 45 degree and swan neck configurations.

### MATERIAL OF CONSTRUCTION

The materials used in the manufacture of the components of the Maestro Microcatheter are identified below.

| COMPONENT          |                                | MATERIAL DESCRIPTION                                                                                                                                                                                                                                                                                                                                                                                        |
|--------------------|--------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Hub [M1]           |                                | polyamide                                                                                                                                                                                                                                                                                                                                                                                                   |
| Strain Relief [M2] |                                | Nylon 12 (polyamide elastomer)<br>Titanium oxide<br>Phtalocyanine blue                                                                                                                                                                                                                                                                                                                                      |
| Catheter Shaft     | Inner Shaft Sleeve Tubing [M6] | Fluon Polytetrafluoroethylene (PTFE)                                                                                                                                                                                                                                                                                                                                                                        |
|                    | Braid filling polymer [M5]     | polyetherblockamide                                                                                                                                                                                                                                                                                                                                                                                         |
|                    | Braid Ribbon Fibers [M4]       | Nylon 66                                                                                                                                                                                                                                                                                                                                                                                                    |
|                    | Outer Layer [M3, M7]           | <u>Proximal Section of Tubing</u> <ul style="list-style-type: none"> <li><input type="checkbox"/> Nylon 12</li> <li><input type="checkbox"/> Phtalocyanine blue</li> <li><input type="checkbox"/> Titanium oxide</li> </ul> <u>Distal Section of Tubing</u> <ul style="list-style-type: none"> <li><input type="checkbox"/> polyetherblockamide</li> <li><input type="checkbox"/> Barium sulfate</li> </ul> |

|  |                                          |                                                                                |
|--|------------------------------------------|--------------------------------------------------------------------------------|
|  |                                          | <input type="checkbox"/> Titanium oxide<br><input type="checkbox"/> Stearamide |
|  | <b>Coating [M9]</b>                      | Hydrophilic Coating                                                            |
|  | <b>Marker Band [M8]</b>                  | Pt                                                                             |
|  | <b>and Marker Band Coating Solution</b>  | polyetherblockamide<br>Stearamide                                              |
|  | <b>Catheter curve straightener [M10]</b> | Nylon 12<br>Titanium oxide                                                     |

#### **ACCESSORIES:**

The Maestro Microcatheter is packaged with two (2) 3 ml syringes, currently CE marked under Certificate CE 542368.

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#### **PACKAGING**

The Merit Microcatheter is sold in a single packaging configuration. The catheter is contained within a circular polyethylene dispenser hoop (holder). The two (2) 3ml syringes that are provided with the catheter are

packaged inside a 4" (10.2cm) x 8" (20.3cm) pouch Tyvek pouch which is nestled on top of the microcatheter holder. The catheter hoop (with the catheter inside) and the pouched syringes are placed inside a 10.5" (26.7cm) x 14.25 (36.2cm) Tyvek/polyethylene pouch that is sealed using a validated heat seal process. One pouch is placed inside an SBS carton and one carton is placed inside a Kraft E-flute shipper carton.

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##### **Storage**

The Maestro Microcatheters must be stored in an air-conditioned area (temperature: 0 – 38C degree, humidity: no condensing).

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### List of Product Label Part Numbers

|                                                                                     |                                                                                             |
|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|
|    | This symbol means that the device fully complies with European Directive 93/42/EEC.         |
|    | Manufacturer                                                                                |
|    | Authorized Representative                                                                   |
|    | Radiopaque marker                                                                           |
|    | Maximum guide wire                                                                          |
|    | Maximum pressure                                                                            |
|    | Single use                                                                                  |
|    | Reference Number                                                                            |
|  | Lot Number                                                                                  |
|  | Sterilization: Ethylene-Oxide Gas                                                           |
|  | Caution: U. S. federal law restricts this device to sale by or on the order of a physician. |
|  | Use by                                                                                      |
|  | Caution: Consult Accompanying Documents                                                     |

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
Instructions for Use  
Merit Maestro Microcatheter  
Dated 13 December 2016  
Consisting of 6 pages


### ACKNOWLEDGEMENT

State of Utah                    )  
                                          ) ss.  
County of Salt Lake        )

On this 13th day of December, 2016, before me, Jennifer Wojciechowski, a notary public, personally appeared Dan Clark, Director, Regulatory Affairs, proved on the basis of satisfactory evidence to be the person whose name is subscribed to the attached instrument, and acknowledge he executed the same.

Witness my hand an official seal.

  
\_\_\_\_\_  
Notary Public  
My Commission Expires: March 6, 2019

# Merit Maestro

## MICROCATHETER

INSTRUCTIONS FOR USE

*Don Clark*  
Don Clark  
12 DECEMBER 2016

# Merit Maestro MICROCATHETER

English

## INSTRUCTIONS FOR USE

**DESCRIPTION:** The Merit Maestro Microcatheter is available in three French size configurations. The first configuration incorporates a change in its outside diameter along its length from a 2.8F (0.93mm) proximal region to a 2.4F (0.80mm) distal region of 20cm in length. This configuration employs a nominal proximal inner diameter of 0.025" with a minimal tip inner diameter of 0.020".

The second configuration maintains a 2.8F (0.93mm) diameter throughout its length, but also incorporates a 20cm flexible distal region. The I.D. of this configuration has a nominal 0.025" inner diameter in the proximal region and a distal region inner diameter of 0.024".

The third configuration maintains a 2.9F (0.96mm) diameter throughout its length, but also incorporates a 20cm flexible distal region. The I.D. of this configuration has a nominal 0.027" inner diameter in the proximal and distal regions.

The Microcatheter lumen is able to accommodate steerable guide wires. A lubricious, hydrophilic coating is applied to the distal 80cm outer surface of all catheter configurations. The Microcatheter has a radio-opaque marker at the distal tip, located approximately 1.3mm proximal of the distal limit of the tip, to facilitate fluoroscopic visualization. The proximal end of the Microcatheter incorporates a standard Luer adapter to facilitate the attachment of accessories.

**INDICATIONS FOR USE:** The Microcatheter is intended for general intravascular use, including peripheral and coronary vasculature. Once the subselective region has been accessed, the Microcatheter can be used for the controlled and selective infusion of diagnostic, embolic, or therapeutic materials into vessels. The catheter should not be used in the cerebral vessels.

**CONTRAINDICATIONS:** None known

### WARNINGS:

- Due to contractual agreements, this Microcatheter is not for neurovascular use at or above the common carotid artery or at or above the vertebral artery.
- This device is intended to be used only by physicians trained in percutaneous intravascular techniques and procedures.
- Contents supplied sterile using an ethylene oxide (ETO) process. Do not use if sterile barrier is damaged.
- For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. After use, dispose of product and packaging in accordance with hospital, administrative, and/or local government policy.
- The infusion dynamic pressure with this Microcatheter should not exceed 800psi (55.15 kPa). Infusion pressure in excess of this maximum may result in Microcatheter rupture, possibly resulting in patient injury. If flow through the Microcatheter becomes restricted, do not attempt to clear the Microcatheter lumen by infusion. The static pressure with this Microcatheter should not exceed 2068 kPa/300psi. Static pressure in excess of this maximum may result in Microcatheter rupture, possibly resulting in patient injury. Identify and resolve the cause of the blockage or replace the Microcatheter with a new Microcatheter before resuming infusion.
- Make sure that the guiding catheter does not slip out of the vessel. If the guiding catheter should leave the vessel when the Microcatheter and/or the guide wire is moved, this may result in the damage of the Microcatheter system.
- Microcatheter advancement beyond the end of the guide wire may result in vessel trauma.
- Do not use a power injector to infuse agents other than contrast media, as the Microcatheter may become blocked. The safety setting of injection pressure must not exceed the maximum injection pressure of 800psi (55.15 kPa). Exceeding injection pressure beyond the maximum injection pressure may cause Microcatheter rupture. (See Instructions For Using a Power Injector)

### PRECAUTIONS:

- Only Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.
- Ensure embolic material compatibility with Microcatheter prior to use.
- Always monitor infusion rates when using the Microcatheter.
- When injecting contrast for angiography, ensure that the Microcatheter is not kinked or occluded.
- The Microcatheter has a lubricious hydrophilic coating on the outside of the catheter. It must be kept hydrated prior to removal from its carrier and during the actual procedure in order to be lubricious. This can be accomplished by attaching the Y-connector to a continuous saline drip.
- Prior to a procedure, all equipment to be used for the procedure should be carefully examined to verify proper function and integrity.
- Inspect the Microcatheter prior to use for any bends or kinks. Any Microcatheter damage may decrease the desired performance characteristics.
- Exercise care in handling of the Microcatheter during a procedure to reduce the possibility of accidental breakage, bending or kinking.
- When the Microcatheter is in the body, it should be manipulated only under fluoroscopy. Do not attempt to move the Microcatheter without observing the resultant tip response.
- Exchange Microcatheters frequently during lengthy procedures that require extensive manipulation or multiple guide wire exchanges.
- Never advance or withdraw an intravascular device against resistance until the cause of the resistance is determined by fluoroscopy. Movement of the Microcatheter or guide wire against resistance may result in separation of the Microcatheter or guide wire tip, damage to the Microcatheter, or vessel perforation.
- Because the Microcatheter may be advanced into narrow subselective vasculature, repeatedly assure that the Microcatheter has not been advanced so far as to interfere with its removal.
- Excessive tightening of a hemostatic valve onto the Microcatheter shaft may result in damage to the catheter.
- Read and follow the manufacturer's IFU for diagnostic, embolic, or therapeutic agents to be used with this Microcatheter.
- Do not use opened or damaged packages. Use prior to the "use before" date. Store at controlled room temperature.

### POTENTIAL COMPLICATIONS:

Possible complications include, but are not limited to:

|                           |                     |
|---------------------------|---------------------|
| Access site complications | Vascular thrombosis |
| Vessel perforation        | Thrombosis          |
| Vessel Spasm              | Ischemia            |
| Hemorrhage                | Infection           |
| Pain and tenderness       | Vessel Dissection   |
| Embolism                  | Distal embolization |
| Allergic reaction         | Death               |

TABLE 1: COMPATABILITY INFORMATION

| Microcatheter O.D. | Nominal Microcatheter I.D. | Maximum Guide wire O.D. | Recommended guiding catheter                          |
|--------------------|----------------------------|-------------------------|-------------------------------------------------------|
| 2.8F / 2.4F        | 0.020" (0.52mm)            | 0.018" (0.46mm)         | 0.040" (1.02mm) to 0.041" (1.04mm)                    |
| 2.8F / 2.8F        | 0.024" (0.62mm)            | 0.021" (0.53mm)         | 0.040" (1.02mm) to 0.041" (1.04mm)                    |
| 2.9F / 2.9F        | 0.027" (0.68mm)            | 0.021" (0.53mm)         | 0.041" (1.04mm) to 0.042" (1.07mm) to 0.043" (1.09mm) |

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**INSTRUCTIONS FOR USE: NOTE:** It is recommended that the Microcatheter be used with a guiding catheter and a sheath introducer.

- Place the appropriate guiding catheter using standard technique. A rotating hemostasis valve may be connected to the guiding catheter luer adapter to continuously flush the guiding catheter with saline.
- Remove the spiral holder housing the Microcatheter from its sealed packaging.
- Attach a syringe filled with heparinized saline solution or sterile water to the luer lock fitting of the Microcatheter holder.
- Inject enough solution to wet the Microcatheter surface entirely. This will activate the hydrophilic coating on the Microcatheter surface. Note: The surface of the Microcatheter may become dry after removal from the holder. Additional wetting with heparinized saline or sterile water will renew the hydrophilic effect.
- Upon removal of the Microcatheter from the spiral holder, inspect Microcatheter to verify there is no damage prior to insertion.
- Attach a second hemostasis valve with side-arm adapter to the Microcatheter, purge any air and flush with heparinized saline or sterile water.
- Carefully insert guide wire into the Microcatheter and completely close the valve around the guide wire.
- Introduce the Microcatheter and guide wire assembly into the guiding catheter via the hemostasis valve (if used). If rotating hemostatic valve is used, tighten the valve around the Microcatheter to prevent backflow, but allowing some movement through the valve by the Microcatheter.
- Using fluoroscopy, introduce the Microcatheter and guide wire assembly into the vascular system, making sure the guide wire is always ahead of the Microcatheter. Advance the guide wire and Microcatheter to a selected vascular site by alternatively advancing the guide wire and then tracking the Microcatheter over the guide wire. Note: To facilitate Microcatheter handling, the proximal portion of the Microcatheter is uncased to ensure a non-slip grip.
- Final positioning is accomplished by short advances of the guide wire and Microcatheter until the desired position is achieved and then confirmed by fluoroscopic visualization.
- Monitor Microcatheter placement and position during use.
- To infuse, completely remove the guide wire from the Microcatheter. Connect a syringe with infusate to the Microcatheter manifold luer, and infuse as required.

### INSTRUCTION FOR USING A POWER INJECTOR WITH THE MICROCATETER:

A power injector can be used to infuse a contrast media through the Microcatheter. Observe the warnings and cautions given below. The flow rate depends upon such factors as the viscosity of the contrast media, which varies with the type and temperature of the media, the model and setting of the power injector, and how the injector is connected to the Microcatheter. The observed flow rate values indicated below are for reference only.

### Flow Rate Tables

| Merit Maestro Catheter Size Shaft/Tip | Usable Length (cm) | Contrast Media      | Iodine Content (mg/ml) | Viscosity (cP) at 37°C | MEDRAD Flow Setting Conditions With Linear Rise @ 0.3sec | Flow Rate (ml/sec) | Volume (ml) | Actual Contrast Delivery ml/sec With Safety Pressure Setting of: 55.15 kPa (800psi) | Dead Space (Priming) Volume (ml) |
|---------------------------------------|--------------------|---------------------|------------------------|------------------------|----------------------------------------------------------|--------------------|-------------|-------------------------------------------------------------------------------------|----------------------------------|
| 2.8/2.4F                              | 110                | ISOVUE- (Iopamidol) | 300                    | 4.7                    | 6.0                                                      | 10                 | 5.55        | 0.63                                                                                |                                  |
|                                       |                    | Iopamidol           | 370                    | 9.4                    | 3.0                                                      | 10                 | 2.54        |                                                                                     |                                  |
|                                       | 130                | ISOVUE- (Iopamidol) | 300                    | 4.7                    | 6.0                                                      | 10                 | 5.14        | 0.70                                                                                |                                  |
|                                       |                    | Iopamidol           | 370                    | 9.4                    | 3.0                                                      | 10                 | 2.21        |                                                                                     |                                  |
| 2.8/2.8F                              | 110                | ISOVUE- (Iopamidol) | 300                    | 4.7                    | 6.0                                                      | 10                 | 4.08        | 0.76                                                                                |                                  |
|                                       |                    | Iopamidol           | 370                    | 9.4                    | 3.0                                                      | 10                 | 2.00        |                                                                                     |                                  |
|                                       | 130                | ISOVUE- (Iopamidol) | 300                    | 4.7                    | 6.0                                                      | 10                 | 5.57        | 0.63                                                                                |                                  |
|                                       |                    | Iopamidol           | 370                    | 9.4                    | 3.0                                                      | 10                 | 2.63        |                                                                                     |                                  |
| 2.9/2.9F                              | 110                | ISOVUE- (Iopamidol) | 300                    | 4.7                    | 6.0                                                      | 10                 | 5.07        | 0.70                                                                                |                                  |
|                                       |                    | Iopamidol           | 370                    | 9.4                    | 3.0                                                      | 10                 | 2.37        |                                                                                     |                                  |
|                                       | 130                | ISOVUE- (Iopamidol) | 300                    | 4.7                    | 6.0                                                      | 10                 | 4.70        | 0.77                                                                                |                                  |
|                                       |                    | Iopamidol           | 370                    | 9.4                    | 3.0                                                      | 10                 | 2.19        |                                                                                     |                                  |
| 2.9/2.9F                              | 110                | ISOVUE- (Iopamidol) | 300                    | 4.7                    | 7.0                                                      | 10                 | 6.82        | 0.69                                                                                |                                  |
|                                       |                    | Iopamidol           | 370                    | 9.4                    | 4.0                                                      | 10                 | 3.44        |                                                                                     |                                  |
|                                       | 130                | ISOVUE- (Iopamidol) | 300                    | 4.7                    | 7.0                                                      | 10                 | 6.26        | 0.77                                                                                |                                  |
|                                       |                    | Iopamidol           | 370                    | 9.4                    | 4.0                                                      | 10                 | 3.43        |                                                                                     |                                  |
| 2.9/2.9F                              | 150                | ISOVUE- (Iopamidol) | 300                    | 4.7                    | 7.0                                                      | 10                 | 5.59        | 0.85                                                                                |                                  |
|                                       |                    | Iopamidol           | 370                    | 9.4                    | 4.0                                                      | 10                 | 3.20        |                                                                                     |                                  |

### REFERENCE DATA


- Injector used: MEDRAD MARK V
- Injection pressure monitor/limit setting: 55.15 kPa (800psi)
- Linear rise seconds: 0.3 sec.

- Contrast Media temperature: 37°C
- Flow scale: ml/sec

## ANNEX TO INSTRUCTIONS FOR USE

### DEVICE SPECIFICATIONS

The Merit Maestro Microcatheter is available in working lengths of 110, 130 and 150 cm with the distal tip offered in pre-shaped straight, 45 degree and swan neck configurations.

### MATERIAL OF CONSTRUCTION

The materials used in the manufacture of the components of the Merit Maestro Microcatheter are identified below.

| COMPONENT          |                                | MATERIAL DESCRIPTION                                                                                                                                                                                                                                     |
|--------------------|--------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Hub [M1]           |                                | polyamide                                                                                                                                                                                                                                                |
| Strain Relief [M2] |                                | Nylon 12 (polyamide elastomer)<br>Titanium oxide<br>Phthalocyanine blue                                                                                                                                                                                  |
| Catheter Shaft     | Inner Shaft Sleeve Tubing [M6] | Fluon Polytetrafluoroethylene (PTFE)                                                                                                                                                                                                                     |
|                    | Braid filling polymer [M5]     | polyetherblockamide                                                                                                                                                                                                                                      |
|                    | Braid Ribbon Fibers [M4]       | Nylon 66                                                                                                                                                                                                                                                 |
|                    | Outer Layer [M3, M7]           | <u>Proximal Section of Tubing</u><br><input type="checkbox"/> Nylon 12<br><input type="checkbox"/> Phthalocyanine blue<br><input type="checkbox"/> Titanium oxide<br><br><u>Distal Section of Tubing</u><br><input type="checkbox"/> polyetherblockamide |

|  |                                          |                                                                                                                           |
|--|------------------------------------------|---------------------------------------------------------------------------------------------------------------------------|
|  |                                          | <input type="checkbox"/> Barium sulfate<br><input type="checkbox"/> Titanium oxide<br><input type="checkbox"/> Stearamide |
|  | <b>Coating [M9]</b>                      | Hydrophilic Coating                                                                                                       |
|  | <b>Marker Band [M8]</b>                  | Pt                                                                                                                        |
|  | <b>and Marker Band Coating Solution</b>  | polyetherblockamide<br>Stearamide                                                                                         |
|  | <b>Catheter curve straightener [M10]</b> | Nylon 12<br>Titanium oxide                                                                                                |

### **3ml SYRINGES:**

The Merit Maestro Microcatheter is packaged with two (2) 3 ml syringes, currently CE marked under Certificate CE 542368.

The Medallion Series Syringes are standard 3 piece piston syringes constructed with clear barrels, color plungers, and a black piston seal.

The syringe barrels are one-piece polycarbonate molded construction with flat finger grips and fixed male luer connectors. The syringe barrels are marked for accurate, clear and easy reading. The safety stop ring prevents inadvertent removal of the plunger from the syringe barrel. The stop ring withstands back pressure of up to 200 PSI in order to make mechanical filling of the syringe possible.

The syringe is intended to be used to infuse a contrast media through the Microcatheter, including the activation of the hydrophilic coating on the Microcatheter surface. (2) 3 ml syringes included for ease of mixing embolics/spheres.

### **Luer Options:**

Fixed Male Luer

Merit's Medallion syringes are categorized below.

| Device Materials |                      |
|------------------|----------------------|
| COMPONENT        | MATERIAL DESCRIPTION |
| Body             | Polycarbonate        |
| Plunger tip      | Silicone             |
| Plunger          | ABS                  |

Fill volume - 3ml

Mass (weight) - ≤8g

Size (general dimensions no more than) - 31.75mm x 12.7mm x 107.95mm

Graduation intervals (max) - 1/16 in. (1.6mm) per grad line

### **PACKAGING**

The Merit Microcatheter is sold in a single packaging configuration. The catheter is contained within a circular

polyethylene dispenser hoop (holder). The two (2) 3ml syringes that are provided with the catheter are packaged inside a 4" (10.2cm) x 8" (20.3cm) pouch Tyvek pouch which is nestled on top of the microcatheter holder. The catheter hoop (with the catheter inside) and the pouched syringes are placed inside a 10.5" (26.7cm) x 14.25 (36.2cm) Tyvek/polyethylene pouch that is sealed using a validated heat seal process. One pouch is placed inside an SBS carton and one carton is placed inside a Kraft E-flute shipper carton.

#### **PRODUCT SHELF LIFE**

The Merit Microcatheter is labeled for a three (3) year shelf life

#### **TRANSPORTATION, STORAGE AND USAGE**

##### **Transportation**

Transportation should be conducted in accordance with the contract which two sides agreed with (temperature: 0 – 38C degree, humidity: no condensing).

##### **Storage**

The Merit Maestro Microcatheter s must be stored in an air-conditioned area (temperature: 0 – 38C degree, humidity: no condensing).

**Operation temperature:** 25°C±10°C

Contact with the internal environment of the body at a temperature of 37°C±1° C

**Use Conditions** - in medical and preventive treatment facility.

#### **DISPOSAL**

Dispose the devices according to local environmental requirements.

Class of hazard for medical wastes by reference to specification of morphological composition (in accordance with Russian Sanitary Regulations and Standards 2.1.7.2790-10) Class B

#### **WARRANTY**



##### **Limited Warranty and Disclaimer**

Merit represents and warrants that the products shall materially conform to the specifications at the time of shipment to Customer and will remain materially free of defects in materials and workmanship for the period specified on the written limited warranty separately provided with each product. This warranty shall not apply to products that have been (i) modified, changed, repaired, refurbished, reprocessed, or altered by anyone other than Merit; (ii) subjected to misuse, mishandling, re-use, re-sterilization, accident, abuse, neglect or tampering; (iii) damaged by excessive physical, environmental or electrical stress; or that have had a serial number altered, defaced or removed; (iv) combined with any other product; (v) used outside the approved "Indications for Use" as cleared by the relevant competent authority, used contrary to the use outlined in the product specifications or in an application or environment for which such product was not designed or contemplated. Any technical product advice furnished by Merit is provided without duty or compensation and Merit assumes no obligation or liability, all such advice being given and accepted at Customer's risk. EXCEPT AS EXPRESSLY SET FORTH HEREIN MERIT MAKES NO REPRESENTATION OR WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, TO CUSTOMER, ITS END-USERS OR ANY THIRD PARTIES WITH RESPECT TO THE PRODUCTS. MERIT DISCLAIMS AND EXCLUDES ANY AND ALL OTHER EXPRESS, IMPLIED, AND STATUTORY WARRANTIES, INCLUDING WITHOUT LIMITATION, WARRANTIES AGAINST INFRINGEMENT, AND THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

##### **Limitation of Liability**

UNDER NO CIRCUMSTANCES SHALL MERIT BE LIABLE TO CUSTOMER OR ANY OTHER PERSON OR ENTITY FOR ANY PUNITIVE, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, WHETHER BASED UPON LOST GOODWILL, LOST RESALE PROFITS, WORK STOPPAGE, PRODUCT FAILURE, IMPAIRMENT OF OTHER PRODUCTS OR OTHERWISE AND WHETHER IN A CONTRACT, TORT OR OTHER ACTION FOR OR ARISING OUT OF BREACH OF WARRANTY, BREACH OF CONTRACT, DELAY, NEGLIGENCE, STRICT LIABILITY, OR OTHERWISE.

#### List of Product Label Part Numbers

|                                                                                     |                                                                                             |
|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|
|    | This symbol means that the device fully complies with European Directive 93/42/EEC.         |
|    | Manufacturer                                                                                |
|    | Authorized Representative                                                                   |
|    | Radiopaque marker                                                                           |
|    | Maximum guide wire                                                                          |
|    | Maximum pressure                                                                            |
|   | Single use                                                                                  |
|  | Reference Number                                                                            |
|  | Lot Number                                                                                  |
|  | Sterilization: Ethylene-Oxide Gas                                                           |
|  | Caution: U. S. federal law restricts this device to sale by or on the order of a physician. |
|  | Use by                                                                                      |
|  | Caution: Consult Accompanying Documents                                                     |

#### Customer Service

Merit Medical Systems, Inc.  
1600 West Merit Parkway  
South Jordan, Utah 84095 U.S.A.  
+ 1 801 253 1600

Merit Technologies LLC  
115054, Moscow  
Bakhrushina Str., 32, bldg. 1, Russia  
+ 7 495 221 89 02

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