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I, the undersigned,

Aryeh Zetler, Notary,

of G Center, Hatamar Street, Yokneam Elite, Israel

do hereby certify that the attached document marked "A" is to the best of my knowledge a true copy of the original document which is drawn up in the English language and has been produced to me.

In witness whereof I hereto set my signature and seal on 30.11.2014.

Fees in the amount of 446 Shekels plus Value Added Tax have been paid.

APOSTILLE

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1. STATE OF ISRAEL

This public document

2. Has been signed by

Advocate

ARYEH ZETLER

3. Acting in capacity of Notary

4. Bears the seal/stamp of

the above Notary

Certified

5. At the Magistrates Court of Nazareth

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
Syneron

"A"₂

Copyright © 2014 Syneron Medical Ltd. All rights reserved.

Syneron Medical Ltd. reserves the right to make changes to its products or specifications to improve performance, reliability, or manufacturability. Information furnished by Syneron Medical Ltd. is believed to be accurate and reliable. However, Syneron Medical Ltd. assumes no responsibility for its use. No license is granted by its implication or otherwise under any patent or patent rights of Syneron Medical Ltd.

No part of this document may be produced or transmitted in any form or by any means, electronic or mechanical, for any purpose, without the express written permission of Syneron Medical Ltd.

Data is subject to change without notification.

Syneron Medical Ltd. has patents and pending patent applications, trademarks, copyrights, or other intellectual property rights covering subject matter in this document. The furnishing of this document does not give you any license to these patents, trademarks, copyrights, or other intellectual property rights except as expressly provided in any written agreement from Syneron Medical Ltd.

Specifications are subject to change without notice.

Manual Catalog Part Number: PB81498EN

Revision Release Date: June 2014

- Before using the device, please check local regulations. If any local legislation is violated, use cannot be authorized.
- In the United States, federal law restricts prescription medical devices to be sold by or on the order of a physician, or properly licensed practitioner. Syneron makes no representations regarding federal, state, or local laws or regulations that might apply to the use and operation of this device.

Contact information for equipment manufacturer and technical services:

- **Syneron, Inc.**
3 Goodyear, Suite A
Irvine, CA 92618
USA
Tel: + 1.866.259.6661 (toll free)
- **CEpartner4u B.V.**
Esdoornlaan 13
3951 DB Maarn
The Netherlands
Tel: + 31.343.442.524
Fax: + 31.343.442.162
- **Syneron Medical (HK) Ltd.**
Room 3103, 31/F, Hopewell Centre,
183 Queen's Road East, Wan Chai, Hong Kong
Tel: +852-2543-4326
+852-9406-3090
Fax: +852-2543-4327
- **Syneron Medical Ltd.**
P.O. Box 550
Industrial Zone, Tavor Building
Yokneam Illit 20692
Israel
Tel: + 972.73.244.2200
Fax: + 972.73.244.2202
- Website: www.Syneron-Candela.com
- E-mail: info@syneron.com

TABLE OF CONTENTS

1. Before You Start.....7

Conventions Used in this Manual7

Explanation of the Symbols Used in the System.....8

2. Introduction.....9

System Overview10

Console10

Applicators.....11

Anatomy of Cellulite.....13

The VelaShape III Procedure13

Heating of Dermis and Hypodermis with elōs.....15

Mechanical Effects of Vacuum on the Skin15

3. Safety16

Introduction16

System Safety Features17

Patient and Personnel Safety17

Grounding the System18

Operating the System18

Infrared Energy Safety Issues18

Radio Frequency Energy Safety Issues18

Mechanical Rollers Safety Issues (VSmooth)19

Electrical and Mechanical Safety.....19

Maintaining the System19

Fire Hazards20

System Classifications.....20

Standard Compliances20

System Nameplate21

4. System Installation22

Unpacking the System.....22

Equipment List23

Consumable Supplies.....23

System Dimensions24

Installation Process.....25

Moving the System within the Facility.....27

Moving the System to another Facility.....27

Electrical Requirements.....27

Environmental Requirements28

5. Operating Instructions29

Activating the System29

Automatic Vacuum and RF Test29

Setting Treatment Parameters – Tablet Control Panel30

Default Settings32

Setting Treatment Parameters – Applicator Control Panel32

Initiating Treatment.....34

Cover Packs34

6. Clinical Guide: Cellulite and Circumference Treatment35

Realistic Expectations.....35

Indications.....35

4

PB81493EN

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✉

INFO@CIRMI.RU

🌐

WWW.GOSZDRAVNADZOR.RU

2014-10-05 05:25:56

Contraindications	35
Possible Side Effects	37
Pre-Treatment Preparation	37
Selecting the Appropriate Treatment Type	38
VelaShape III Measurement	39
Circumference	40
Cellulite	41
Treatment Procedure	42
Applicator Positioning	43
Test Spot	43
Treatment Technique: Identifying the Applicator and Cover	43
Treatment Using the VContour Applicator	43
Treatment Protocol Using the Circumference Mode	44
Treatment Protocol of Large Areas Using Circumference Mode (Medium Cover)	44
Treatment Protocol of Small Areas Using the Circumference Mode (Small Cover)	47
Treatment Protocol Using the Deep Mode	48
Treatment Protocol of Focal Fat Deposits Using Deep Mode (Large and Medium Covers)	48
Cellulite Mode Treatment Protocol Using the VSmooth Applicator	54
Treatment Protocol Using Cellulite Mode	54
Treatment Options with Cellulite Mode	59
Endpoints (When to Stop Treating an Area)	60
Post-Treatment Care	61
Treatment Conclusion	61
7. Maintenance	62
Cleaning the Console	62
Cleaning and Maintaining the VSmooth Applicator	63
Replacing the VSmooth Applicator's Filters	65
Cleaning and Maintaining the VContour Applicator	66
Cleaning the VContour Applicator's Filter	67
Cleaning the System Air Filter	68
8. Troubleshooting Guide	69
9. System Specifications	72
Output Parameters	72
Treatment Area Footprints	72
Control Panel	72
Electrical Requirements	72
Environmental Requirements	73
Treatment Room	73
Storage/Transportation	73
Degree of Protection Against Ingress of Water	73
System Console Dimensions [W x D x H]	73
Weight	73
10. Electromagnetic Regulation	74
Electromagnetic Emissions	74
Electromagnetic Immunity	75
Recommended Separation Distances	77

LIST OF FIGURES

Figure 1: VelaShape III System	10
Figure 2: VContour Applicator [Left: Operator Panel / Right: Applicator Elements].....	11
Figure 3: VSmooth Applicator Elements	12
Figure 4: Anatomy of Cellulite	13
Figure 5: elōs Heating and Mechanical Manipulation Action	13
Figure 6: Before Treatment	14
Figure 7: After VelaShape III Treatment Reduction in the Size of the Fat Cells and Fat Chambers	14
Figure 8: System Nameplate Location	21
Figure 9: System Dimensions	24
Figure 10: Applicator on Handle.....	25
Figure 11: Umbilical Connection Port.....	25
Figure 12: Service Panel	26
Figure 13: Main Menu Screen.....	30
Figure 14: Treatment Settings Ribbon	31
Figure 15: Default Settings.....	32
Figure 16: Applicator Control Panel with Initial Temperature Display.....	32
Figure 17: Applicator Control Panel Display	33
Figure 18: Control Panel Buttons	33
Figure 19: Parameter Icons and Level Displays	33
Figure 20: Illustration of Circumference Measurements: A) Fixed Distance between the Legs; B) Labeling the circumference; C) Measuring Horizontally; D) Labeling the Arm Against a Wall.	40
Figure 21: Illustration of cellulite stages (A- Stage Zero; B- Stage One; C- Stage Two; D- Stage Three) ..	41
Figure 22: Sub-Area Divided into Segments.....	44
Figure 23: Identifying Sub-Areas and Segments for Anterior and Posterior Abdomen Circumference Reduction	44
Figure 24: Cycle 1 "Horizontal Snake-Like"	45
Figure 25: Cycle 2 "Vertical Snake-Like".....	46
Figure 26: Sub-Area Divided into Segments.....	47
Figure 27: Cycle 1 and Cycle 2 in Small Sub-Areas	48
Figure 28: Deep Mode protocol Segmentation	48
Figure 29: Flanks (A) and Middle Back (B) Segmentation for the Treatment of Focal Fat Deposit.....	49
Figure 30: Cycle 1-Step #1 "Horizontal Snake-Like".....	51
Figure 31: Cycle 1 – Step #2 "Vertical Snake-Like"	52
Figure 32: Cellulite Mode Protocol Sub-Area Segmentation	54
Figure 33: Segment A Treatment Illustration	56
Figure 34: Segment B Treatment Illustration	56
Figure 35: Segment C Treatment Illustration	57
Figure 36: Segment D Treatment Illustration	57
Figure 37: Horizontal Snake-Like	58
Figure 38: Vertical Snake-Like	58
Figure 39: Removing the VSmooth Replaceable Cover	63
Figure 40: VSmooth Applicator's Internal Components	64
Figure 41: Replacing the VSmooth Applicator's Filters	65
Figure 42: Applicator's Internal Components	66
Figure 43: Replacing the VelaShape III Applicator's Filter.....	67
Figure 44: Cleaning the Air Filter	68

1. BEFORE YOU START

Warning

- Read this manual to become familiar with all safety requirements and operating procedures before attempting to operate the system.
- Any intense optical energy or radio frequency (RF) device can cause injury if used improperly.
- High voltage is present inside the system. Always be aware of the possible dangers and take proper safeguards as described in this manual.
- The VelaShape III system must be serviced only by Syneron Medical Ltd. qualified personnel.

Conventions Used in this Manual

The following conventions in the form of notes, cautions and warnings are used in this manual:

Note

The content of this **Note** offers general information that is important to keep in mind.

Caution

A **Caution** alerts the user to the possibility of a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or damage to the equipment.

Warning

A **Warning** alerts the user to the possibility of injury, death, or serious adverse reactions associated with the use or misuse of the system.

2014-10-05 05:25:56

Explanation of the Symbols Used in the System

Symbol	Description
	Refer to instruction manual/booklet
	Caution!
	CE Compliance Symbol
	Waste of Electrical and Electronic Equipment (WEEE) compliance symbol
	RF Interference may occur in the vicinity of equipment with the following symbol
	Product Shelf-Life: product may not be stored for longer than this period after being opened for the first time
	Fuse
	Type BF Equipment
	Manufacturer
	Authorized Representative in the European Union

2014-10-05 05:25:56

2. INTRODUCTION

The VelaShape III is a device that combines **elōs** technology (electro optical synergy) with mechanical manipulation of the skin using gentle **vacuum** suction.

The system is designed for use in medical and aesthetic practices, to be user-friendly, safe and reliable. It integrates the following technologies:

- Infrared (IR) optical energy.
- Bipolar conductive radio frequency (RF) current.
- Mechanical manipulation (vacuum with or without massage rollers).

These three technologies work together to provide dermal and sub-dermal heating, and mechanical modifications resulting in improved blood circulation, enhanced lymphatic drainage, as well as an increased metabolic rate of the adipose tissue.

While this system is considered safe for both the operator and the patient, it should be operated only by trained, qualified treating attendants.

System Overview

The system includes three core modules: the **Console**, **Applicator** and **Touch-Screen Control Panel** (see Figure 1).


Figure 1: VelaShape III System

Console

The VelaShape III console houses the following components:

- Power supply
- Vacuum pump
- RF generator
- Other power components
- IR LEDs power supply
- Main CPU
- Tablet

2014-10-05 05:25:56

Applicators

The **VSmooth** and **VContour** applicators are used for the physical implementation of the RF current, IR energy and mechanical manipulation treatment parameters.

The **VContour** applicator (see Figure 2) utilizes IR, RF and mechanical manipulation by vacuum. It is used for cellulite and circumference reduction of curved areas, for localized fatty deposits, and to reinforce results of **VSmooth** applicator, thus leading to contouring of smaller areas such as the lower abdomen, flanks and arms.

The **VSmooth** applicator (see Figure 3) utilizes IR, RF and mechanical manipulation by vacuum and massage rollers. It is used for cellulite reduction in large areas such as the thighs, buttocks, or abdomen and for circumference reduction.

The VelaShape III applicators are composed of a permanent handpiece and a replaceable cover. The replaceable covers of the VelaShape III applicators are usable for limited periods of time.

The replaceable cover of the VelaShape III applicator should be replaced according to the schedule printed on the cover's storage box, or sooner if one of the following is observed:

- After cleaning the replaceable cover and applicator spring connectors, an **Improve Vacuum Seal** message appears on the screen.
- Electrodes appear damaged.

The filter should be cleaned or replaced according to the schedule printed on the cover's storage box, or sooner if a reduction in vacuum effectiveness is observed.

Figure 2 shows the display and selector controls that are located on the operator panel of the applicator as well as on the touch-screen control panel (see Figure 2).


Figure 2: VContour Applicator [Left: Operator Panel / Right: Applicator Elements]

Figure 3 shows the VSmooth applicator and its operating elements:


Figure 3: VSmooth Applicator Elements

Note

The images in the figures are for illustration purposes only.

The applicator is applied to the treatment area so that all four edges of the aperture make a seal with the skin. This allows the vacuum to be effective and ensures that the RF electrodes are fully coupled to the skin. Failure to achieve full contact between the RF electrodes and skin, or failure to achieve proper vacuum will result in disabled IR and RF energies, indicated by an **RF out of range** message on the system's control panel and applicator screen.

2014-10-05 05:25:56

Anatomy of Cellulite

Cellulite is characterized by **Dimpled Skin Surface** – mattress-like appearance (see Figure 4). The fat cells are located in chambers which are separated by vertical and horizontal **Connective Tissue Septae**. The **Septae** are anchored rigid fibers that do not stretch with the skin as the volume within the walls of the chamber increases. When the fat chamber becomes too distended by enlarged fat cells and by retained intercellular fluids, it protrudes beyond the length of the vertical septae, pushing up the skin's surface and resulting in a dimpled skin appearance.


Figure 4: Anatomy of Cellulite

The VelaShape III Procedure

The application of **elōs** (IR and RF energies) to the dermis/hypodermis is intended to deploy heat to the tissues that contribute to cellulite. The **vacuum** contributes to increased blood circulation (see Figure 5). The blood, which is one of the targets of the hyperthermic effect of the IR and RF, responds in a higher diffusivity of oxygen. This in turn increases the metabolic breakdown of stored fat by the fat cells (lipolysis), thereby reducing the size of the fat cells. Consequently, the bumpiness of the skin, which is largely caused by distended fat cells protruding up from the anchored connective tissue septae, is reduced.


Figure 5: *elōs* Heating and Mechanical Manipulation Action

2014-10-05 05:25:56

Both actions, reducing the size of the fat cells and draining excess fluids to the lymphatic system, are accompanied by shrinkage of the connective tissue septa and facilitate the overall shrinkage of the size of the fat chambers, thus leading to a reduction in the appearance of cellulite and a reduction in the circumference of the treated area.


Figure 6: Before Treatment


Figure 7: After VelaShape III Treatment Reduction in the Size of the Fat Cells and Fat Chambers

2014-10-05 05:25:56

Heating of Dermis and Hypodermis with elōs

Light Energy (see Figure 5):

- IR light at wavelength of 850 nm and maximum power of 3.3 Watts.
- Maximum electrical power of the IR light is 9 Watts (9 J/sec).
- Depth of heating from light reaches the deep dermis.

RF Energy (see Figure 5):

- Bi-polar conductive RF energy is emitted from the RF electrodes at a frequency of ~1 MHz.
- Maximum electrical power of the RF energy is 150 Watts.
- The heat generated by the RF energy reaches the depth of the hypodermis.

Mechanical Effects of Vacuum on the Skin

- When all edges of the applicator are in full contact with the skin, a seal is made, allowing the vacuum pump to apply gentle suction.
- Negative pressure (vacuum) dilates the blood vessels, increasing the perfusion in the local area.
- Vacuum levels are up to 390 mbar of negative pressure.

3. SAFETY

This chapter describes the safety issues regarding the use and maintenance of the VelaShape III system, with special emphasis on optical and electrical safety.

Introduction

The system is designed for safe and reliable treatment, when used in accordance with proper operation and maintenance procedures. Only trained practitioners are qualified to operate the system. The user and all other personnel operating or maintaining the system should be familiar with the safety information provided in this chapter. Syneron Medical Ltd. assumes no liability whatsoever for any damage or injury as a result of an application of a product which is not in strict accordance with the instructions provided with the product. The primary consideration should be to maximize safety for both treating attendant and the patient.

Warning

- Read this chapter to be familiar with all of its safety requirements and operating procedures prior to operating the system.
- Any intense optical energy or RF device can cause injury if used improperly.
- High voltage is present inside the system.
- Always be aware of the possible dangers and take proper safeguards as described in this manual.
- **In the event of adverse effects or injury/illness, immediately contact:**
 - In North America: **Syneron, Inc.**
 - In Europe: **CEpartner4u B.V.**
 - Outside North America and Europe: **Syneron Medical Ltd.**
- For complete contact information please refer to page 3 of this manual.

2014-10-05 05:25:56

System Safety Features

The VelaShape III system incorporates the following safety features:

- The power electronics cannot be activated unless an applicator has been connected to the system.
- Upon system activation, IR, RF, and Vacuum are at Level **0** on the control panel.
- During activation, the system performs a self-test of the hardware to detect any faults, so that you can shut down the system if necessary.
- IR and RF energies are not delivered unless full contact with the tissue is conducted and the applicator's trigger or footswitch is pressed.

Patient and Personnel Safety

When handling the VelaShape III system, the primary safety concern should be the safety of both the treating attendant and the patient. Follow these personal safety guidelines to maintain a high level of personnel safety:

- Make sure that only well-trained attendants perform treatment.
- Before activating the system, make sure any long hair of the treating attendant or the patient is safely pulled back and tied.
- Ascertain that all jewelry is removed.
- Prior to scheduling treatments, obtain a detailed patient history to verify that the patient is eligible and is in fit condition to undergo treatment.
- Fully inform the patient regarding the treatment protocol, the likely results and any risks associated with treatment.
- Improper treatment may cause local burns.
- Opaque eye protection should be used when working within a 30 cm (12") proximity to the patient's eyes.

Warning

Prior to system activation, any personnel operating or maintaining the system should be familiar with the safety information provided in this chapter.

Grounding the System

When installing the system, ensure that the system is grounded through the grounding conductor in the power cable. This protective grounding is essential for safe operation.

Operating the System

- Keep all system panels closed. Do not remove the panels and do not touch any of the system's inner components.
- Keep your hands away from the applicator during system start-up.
- Keep the applicator dry and clean.
- Do not allow the replaceable applicator covers to come in contact with hard materials, excessive heat or humidity; these may damage the RF electrodes.
- Never use the system in the presence of explosive or flammable materials.
- Make sure that all treatment personnel are familiar with the system's controls and know how to shut down the system in the event of an emergency.

Infrared Energy Safety Issues

The VelaShape III emits gentle heating infrared light in a continuous wave mode. Light emission is contained within the applicator chamber and occurs only when the operator presses the trigger button on the applicator or footswitch, when an adequate vacuum seal is achieved and the RF electrodes are in full contact with the treated surface. Opaque eye protection should be used when working within a 30 cm (12") proximity to the patient's eyes.

Do not aim the applicator at anything but the intended target area; it should always be oriented toward the treatment area, and the trigger button or footswitch pressed only when full contact is made with the skin. In rare cases superficial crusting may be formed on the skin from the IR effect.

Radio Frequency Energy Safety Issues

The RF energy provides gentle heating that is distributed deep into the skin via bi-polar electrodes. Full coupling of both electrodes to the skin, after applying VelaSpray Ease or any other recommended product, is necessary before pressing the trigger button or footswitch, so that the RF energy is safely conducted to penetrate the skin. If coupling with the skin is broken due to incomplete contact or removing the applicator from the treatment site while the trigger or footswitch is pressed, IR and RF modes are disabled and an error message on the control panel indicates insufficient coupling. The trigger or footswitch should be released immediately and the applicator should be re-oriented on the skin. In rare cases superficial crusting may be formed on the skin from the RF effect.

2014-10-05 05:25:56

Because of the heating effect by RF energy, highly conductive objects such as superficial metal or conductive implants, rods, plates or pins should be avoided in the treatment area.

Mechanical Rollers Safety Issues (VSmooth)

Beware of any loose material in or near the immediate treatment area that may get caught in the rollers. There should be no loose cords, clothing, hair, etc. in the treatment area. Ensure integrity of rollers' coating before and during treatment.

Electrical and Mechanical Safety

- Keep all system panels closed. Removing the panels creates a safety hazard.
- Keep hands away from the applicator during system start-up.
- Perform maintenance procedures when the system is shut down and disconnected from its electrical power source.
- The system's cables are equipped with isolation coating and the electronic circuitry is also isolated from the mains power supply. In addition, the system's plastic covers serve as a means of isolation.
- Move the system slowly and carefully. The system weighs approximately 20 kg (44 lbs.) and may cause injury if proper care is not used when moving it.
- The system is grounded through the grounding conductor in the power cable. This protective grounding is essential for safe operation. The RF electrodes are the only non-grounded component of the system, due to the nature of the treatment.
- Portable and mobile RF communications equipment can affect medical electrical equipment, Interference may occur in the vicinity of the equipment.
- The use of accessories and cables other than those specified may result in increased emissions or decreased immunity of the device. It is forbidden to connect any third-party equipment to the system without written approval from Syneron Medical.

Maintaining the System

- Only authorized Syneron technicians are qualified to service the system.
- Shut down the system and disconnect the power cable before performing any maintenance procedures (see Chapter 7).

2014-10-05 05:25:56

Fire Hazards

- The absorption of optical energy or conduction of RF raises the temperature of the absorbing material. Do not use the system in the presence of explosive or flammable materials.
- Do not use flammable substances when preparing the skin for treatment.
- If alcohol is used for disinfecting, it must be allowed to dry thoroughly before the system is used.

System Classifications

- Electric shock protection: Class A, Group I
- Protection against ingress of liquids: Ordinary equipment
- Not suitable for use in presence of flammable substances.

Standard Compliances

- EN 60601-1-1:2006 Medical Electronic Equipment Part 1: General Requirements for Safety.
- EN 60601-1:2007 Medical Electronic Equipment Electromagnetic compatibility and FCC part 18.
- UL94 Flammability requirements.
- ISO 14791:2012 Application of Risk Management to Medical devices.
- IEC 62304:2006 Medical device software - Software life cycle processes (If applicable).
- ISO 10993-1:2009 Biological evaluation of medical devices.
- ISO 13485:2003 Medical devices - Quality management systems -- Requirements for regulatory purposes.

Warning

No modifications to this equipment are allowed!

2014-10-05 05:25:56

System Nameplate

The following illustration (see Figure 8) displays the location of the system's nameplate, which includes the system's identification, regulatory and electrical requirement information:


Figure 8: System Nameplate Location

4. SYSTEM INSTALLATION

The system is designed for installation in a clinical environment. The customer will carry out the installation as described in the following sections.

Unpacking the System

After unpacking the system, Syneron recommends that you retain the shipping box and the internal Styrofoam covers and braces in the event that the system will have to be relocated anytime in the future.

Unpack the system as follows:

1. Cut off and remove the three plastic straps that encase the cardboard box.
2. Lift the cardboard box up and off the packaged components and set it aside.
3. You will find a "**VelaShape III Unpacking Instructions**" sheet taped to the side of the top Styrofoam cover; follow the instructions on this sheet, and then continue with unpacking and installation from this point.
4. Carefully remove the packed accessories and small Styrofoam boxes from inside the Styrofoam covers.
5. Finish removing all of the Styrofoam covers and braces from around the body of the system.
6. Remove the plastic bag that is wrapped around the system.
7. Compare the accessories you unpack to the equipment list (next section) to ensure that all accessories are present. Report any missing accessories to your Syneron representative.
8. Lift the system up and out of the bottom Styrofoam brace and carefully set it on the floor.

Warning

The system weighs 20 kg (44 lbs.) and should only be lifted by two people. Lift correctly and protect your back!

9. Finish unpacking the accessories and continue with the installation process.

2014-10-05 05:25:56

Equipment List

The VelaShape III system is delivered with the following:

- VelaShape III system console.
- VContour Applicator (includes one medium-size replaceable cover).
- VSmooth Applicator (includes one replaceable cover).
- **Accessories kit, including:**
 - ⇒ Footswitch
 - ⇒ Applicator umbilical cable.
 - ⇒ System air filters (x 3).
 - ⇒ VelaSpray Ease lotion (2 bottles of 300 ml).
 - ⇒ VContour small replaceable cover.
 - ⇒ VContour large replaceable cover.
 - ⇒ One pair of protective opaque goggles.
 - ⇒ User Manual.
 - ⇒ Quick-Reference Guide.
 - ⇒ Filter extractor.

Consumable Supplies

The following consumable supplies may be ordered from your Syneron distributor:

- VelaShape III applicator's covers packs.
- VelaSpray Ease lotion.

System Dimensions

Space should be allocated with adequate ventilation and free airflow. The working area for the system should be prepared according to the system dimensions presented in Figure 9. In order to guarantee proper ventilation, always keep the sides of the system at least 0.5m (20") from the wall or from other obstructions to air flow. After positioning the system, lock the brakes on the rear wheels by pressing the pedals on top of each wheel.


Figure 9: System Dimensions

2014-10-05 05:25:56

Installation Process

After unpacking the VelaShape III system, install the applicator in the following manner:

1. Refer to Figure 10: ensure that the rubber grips are installed on both handles. If they came in the accessories box, slip them on to the handles.
2. Unpack the applicator and hang it on of the system's handles, on the rubber grip to prevent the applicator from slipping off the handle and falling.


Figure 10: Applicator on Handle

3. Insert the applicator's umbilical connector into the applicator's connection port, and press it in until it 'clicks' into place. The connector and connection port are polarized to prevent improper connection.


Figure 11: Umbilical Connection Port

Refer to Figure 12:

4. Unpack the power cable, ensure that it has a proper wall connector (for your country) and connect it to the power connection port on the system's service panel (6), and plug the other end into the appropriate electrical outlet.
5. Unpack the footswitch and connect it to the footswitch connection port on the system's service panel (2). The connector and connection port are polarized to prevent improper connection.

Note

If the system has been removed from cold storage, allow it time to warm up to room temperature (15-30°C/59-86°F) before first activation.


Figure 12: Service Panel

Service Panel Connections:

1. USB communications port (used only by Syneron service engineers).
2. Footswitch connection port.
3. Umbilical connection port.
4. Main On/Off switch.
5. Fuse housing (see Chapter 7 – **Maintenance**).
6. Power cable inlet.

2014-10-05 05:25:56

Moving the System within the Facility

1. Turn the system off and disconnect the power cable.
2. Disconnect the applicator and place it in a safe location.
3. Release the wheel brakes at the base of the system.
4. Slowly push or pull the system using the system's handle.

Moving the System to another Facility

1. Turn the system off and disconnect the power cable.
2. Disconnect the applicator and pack it in its storage case.
3. Release the wheel brakes at the base of the system.
4. Slowly push or pull the system using the system's handle and carefully lift it up onto the vehicle. Spread padding on the wall of the transport vehicle to avoid scratching the system's panels and strap the system to the wall.
5. Upon arrival at the new facility, position the system in its new location, using the system handle.
6. Reconnect the applicator, footswitch and the power cable.

Caution

Do not lay the system on its back or sides! Doing so may damage the vacuum pump.

Electrical Requirements

The system incorporates a universal power supply module that allows it to be connected to standard electrical outlets anywhere. The system's electrical requirements are:

- 100-230 VAC; up to 4A; 50-60 Hz; single phase

The electrical requirements are printed on the system's nameplate, located on the system's rear panel.

Input power lines should be free of transients, voltage and current spikes, sags and surges. Consequently, the system power line should not be shared with other heavy variable loads such as elevators, air conditioning systems, large motors, etc.

It is recommended that the system be connected to a separate power line with separate circuit breakers. Syneron cannot guarantee performance unless the system is connected to a dedicated circuit.

**Warning**

- For continued protection against fire, replace the fuses only with those of the same type and rating.
- Proper grounding is essential for safe operation.
- Do not replace any electric part of the system except fuses!
- Only authorized Syneron technicians are qualified to service the system.

Environmental Requirements

- Corrosive materials can damage electronic parts; therefore the system should operate in a non-corrosive atmosphere.
- Metallic dust can damage electrical equipment. Dust particles absorb light and heat up. Hot particles located on the aperture guide can damage it; therefore dust particles should be kept to a minimum.
- For optimal operation of the system, maintain room temperature between 15°-30°C (59°-86°F) and relative humidity of up to 80%.
- If the system has been stored in an environment cooler than the above mentioned temperatures, allow it time to warm up to room temperature before first activation.

2014-10-05 05:25:56

5. OPERATING INSTRUCTIONS

This chapter will assist you to operate the system for the first time, and will guide you through the initial setup procedures needed to perform in order to start operating the system.

Activating the System

1. Before turning the system on, verify that the power cable is connected to the power connection port on the system's rear panel and to the wall power outlet, connect the umbilical cable to its connection port and ensure that the applicator is connected to the umbilical cable.

Warning

High voltage is present inside the system. Always be aware of the possible dangers and take proper safeguards as described in this manual.

2. Toggle the main power switch (on the service panel) to the **On** position; the system initiates an automatic vacuum and RF test. When this test is completed the system is ready for operation.

Automatic Vacuum and RF Test

The VelaShape III is equipped with an automatic vacuum and RF self-test routine.

The test is performed after the applicator is connected.

During the test the VelaShape III system is tested for radiofrequency and vacuum performance. During this test, the trigger or footswitch of the applicator is disabled.

The test is comprised of three consecutive vacuum pulses, and can take up to 3 seconds.

It is important to verify that the applicator's chamber is not blocked or positioned against a sealing surface during the test (i.e., verify that the applicator is not applied to the patient during the test).

If a fault situation is detected, it will be displayed as an error message on the control panel (refer to Chapter 8 – **Troubleshooting**).

2014-10-05 05:25:56

When you 'pull up' the **Options** tab the **Treatment Settings** ribbon appears. The numbered tags in Figure 14 correlate to the following numbered explanations and definitions of the screen's elements.


Figure 14: Treatment Settings Ribbon

1. **IR Level** – using your finger to 'tap' on the **IR** (Infra-Red) button in order to toggle between these intensity levels: **0, 1, 2** or **3**.
2. **RF Level** – using your finger to 'tap' on the **RF** (Radio-Frequency) button in order to toggle between these intensity levels: **0, 1, 2** or **3**.
3. **VAC Level** – using your finger to 'tap' on the **VAC** (Vacuum) button in order to toggle between these intensity levels: **1, 2** or **3**.
4. **End Treatment** – press this button to end the patient's session, reset the clocks and be ready for the next patient.
5. **Utilities** – for use only by Syneron-trained system engineers, or under direct instructions from the Syneron Service Center.

2014-10-05 05:25:56

Default Settings

The intensity level of each parameter setting is displayed in the following illustrations:

Figure 15: Default Settings

Setting Treatment Parameters – Applicator Control Panel

After activation (main On/Off switch), the system starts the automatic vacuum and RF self-test routine. When complete, the display on the connected applicator initializes while displaying an icon of the connected type of cover for several seconds.

Press the central **OK** button (see Figure 18) to toggle between views of the current temperature sensor (see Figure 16), the parameter setting icons (see Figure 19) and information display screen (see Figure 17).


Figure 16: Applicator Control Panels

2014-10-05 05:25:56

**Note**

Both control panels – tablet and applicator – operate in parallel and display the same information.

Press the central **OK** button until the display changes to the information display screen from which you may adjust the system's operating parameters.

Refer to Figure 17: the following information is displayed:

- **IR** level currently set
- **RF** level currently set
- **Vacuum** level currently set
- **mBar**: current vacuum reading
- **T** (reading on left side): current temperature reading (reading on right side is only for use of Syneron technical personnel).


Figure 17: Applicator Control Panel Display

Refer to Figure 18: press the ► or ◀ buttons to toggle between the parameter icons (**IR**, **RF** or **VAC**).

Press the ▲ or ▼ button to select the desired intensity level for that parameter.

You may now proceed with the treatment procedure.


Figure 18: Control Panel Buttons

Figure 19: Parameter Icons and Level Displays

Initiating Treatment

Now that you have prepared the system and defined your treatment settings, you are ready to begin treatment. Apply the applicator with moderate pressure to the first treatment section and press the trigger on the applicator handle or the footswitch. To stop the treatment, release the trigger or the footswitch.

If desired you may pull up the **Options** tab to keep the parameter setting icons in view while you are performing the treatment.

Warning

You may treat patients only after passing appropriate training and have adhered to all local and state regulations regarding the use of this system.

Cover Packs

Four types of cover packs are available:

- **VContour Large Cover Pack** (25x50 mm) – contains five large covers and appropriate filters.
- **VContour Medium Cover Pack** (25x30 mm) – contains five medium covers and appropriate filters.
- **VContour Small Cover Pack** (10x13 mm) – contains five small covers and appropriate filters.
- **VSmooth Cover Pack** – contains five covers for the VSmooth applicator and appropriate filters.

Caution

Always ensure that the applicator filters are clean, dry and have not passed their maximal number of operating hours before starting a treatment.

6. CLINICAL GUIDE: CELLULITE AND CIRCUMFERENCE TREATMENT

This chapter provides you with the information required to treat cellulite and circumferences with the VelaShape III system. The entire treatment protocol is to the cleared indication, and unless otherwise specified, for each application specific procedure should be applied.

Realistic Expectations

- There will be improvement in cellulite appearance, but not necessarily complete elimination. There will be improvement in circumferences. However the improvement of both cellulite and circumference reduction is based on individual being treated.
- The degree of response to the VelaShape III treatment, and the number of treatment sessions required will vary among patients and will depend on the clinical and physiological condition at the start of the treatment regimen. Some patients respond better than others.
- The treatment results are temporary and one maintenance treatment session every 2-3 months is essential to sustain them.
- A healthy lifestyle (diet and exercise) may help to obtain better results, but is not essential. However, weight gain may have a negative effect on the results.

Indications

USA: The VelaShape III system is indicated for temporary reduction in circumferences of the abdomen and thighs, as well as temporary reduction in the appearance of cellulite. It is also indicated for the relief of minor aches and muscle spasms, as well as for the improvement of blood circulation.

EU/Canada/Singapore: The VelaShape III system is indicated for non-invasive treatment for body contouring via temporary cellulite and circumference reduction.

Contraindications

- Pacemaker or internal defibrillator.
- Superficial metal or other implants in the treatment area.
- Current or history of skin cancer, or current condition of any other type of cancer, or pre-malignant moles.

2014-10-05 05:25:56

- History of any kind of cancer.*
- Severe concurrent conditions, such as cardiac disorders.
- Pregnancy and nursing as well as 4-5 months post- natural delivery, or at least 6 months post-cesarean section, or until normal hormonal balance is regained.
- Impaired immune system due to immunosuppressive diseases such as AIDS and HIV, or use of immunosuppressive medications.*
- Diseases which may be stimulated by light at the wavelengths used.*
- Patients with a history of diseases stimulated by heat, such as recurrent Herpes Simplex in the treatment area, may be treated only following a prophylactic regimen.
- Poorly controlled endocrine disorders, such as diabetes.
- Any active condition in the treatment area such as sores, psoriasis, eczema and rash.
- History of skin disorders, keloids, abnormal wound healing, as well as very dry and fragile skin.
- History of bleeding coagulopathies, or use of anticoagulants.*
- Use of medications, herbs, food supplements & vitamins known to induce photosensitivity to light exposure at the wavelengths used, such as Isotretinoin (Accutane) within the last 6 months, tetracyclines, or St. John's Wort within the last 2 weeks.*
- Any surgical procedure in the treatment area within the last 3 months or before complete healing.
- Treating over tattoo or permanent makeup.
- Area with varicose veins.
- As per the practitioner's discretion, refrain from treating any condition which might make it unsafe for the patient.

() Although not recommended, these conditions may be treated at the discretion and under the full responsibility of the medical director/physician. In such cases, a small area should be treated and assessed a few days later to determine if the patient will tolerate the treatment without developing short-term adverse effects. Test spots cannot indicate long-term adverse events.*

**Note**

- In case of uncertainty regarding applicability or potential side effects, have the patient consult his/her physician and bring a written consent for treatment.
- Always test a small area with the selected parameters, and wait several minutes to detect results of the tested area before proceeding with the whole area treatment.
- Additionally for dark skin types (IV-VI), a small area should be treated and assessed a few days later to determine if the patient will tolerate the treatment without developing adverse effects.
- It is not recommended to treat the abdomen sooner than one hour after a meal.
- It is not recommended to treat patients with BMI >30.

Possible Side Effects

Certain side effects may be experienced during treatment or shortly afterwards, usually as a result of improper use of the system. Although these side effects are rare and temporary, they should be reported immediately to a physician for proper treatment.

These are the side effects that may appear in the treatment area:

- Pain
- Excessive skin redness (erythema)
- Hair pulling (by mechanical rollers)
- Damage to natural skin texture (scratching, crusting, blister, burn)
- Bruising

Pre-Treatment Preparation

- Complete or update the patient's medical and physical history.
- Exclude from treatment patients presenting any of the contraindications for treatment.
- Determine why the patient is seeking treatment and what the expectations are.
- Inform the patient about the treatment protocol, typical treatment results and possible adverse effects or discomfort.
- Instruct the patient about the safety issues.

2014-10-05 05:25:56

- Advise the patient to exfoliate the skin with appropriate creams before treatment. The patient should come to the treatment with clean, dry skin (no perfume, no lotion).
- Hair may impair the movement of the applicator and/or overheat the treatment area. Removing hair before treatment is advised.
- Advise the patient to avoid anticoagulants, such as aspirin throughout the treatment regimen, if medical condition permits and pertinent to physician approval. Anticoagulants increase the possibility of bruising.
- Before the first treatment session it is important to weigh your patient, photograph and measure the circumference of the treatment area in order to establish a baseline (see VelaShape III Measurement, page 39).
- Re-photograph and re-measure before the last treatment and before any maintenance treatment thereafter, and document the results achieved. Additionally you may do so every two sessions to document any progress. Re-weigh your patient before each subsequent session to ascertain no weight gain.
- Mark the treatment areas when the patient is standing (use white or yellow marker), as cellulite or contour may change characteristics when the patient is lying down.
- Make sure that all the applicators' parts are clean and intact. Make sure that the RF electrodes are clean and intact (not chipped). When using the VSmooth applicator, verify that the rollers function smoothly.
- Make sure that the cover filters are clean and correctly inserted (replace them if not).
- Ensure that the replaceable cover is firmly inserted and secured in place.
- Fit the patient with opaque eye protectors if the treatment area is within a 30 cm (12") proximity to his/her eyes.

Selecting the Appropriate Treatment Type

The VelaShape III combines two types of energy, IR and RF, and mechanical manipulations of vacuum and/or rollers (massage) in two different applicators: VContour applicator and VSmooth applicator. Each applicator combined and synchronized different actions to optimize the energy delivery. The VContour applicator combines the RF energy and mechanical manipulation of vacuum, while the VSmooth applicator combines the IR and RF energy and mechanical manipulation via vacuum and rollers (massage).

You should examine the patient's skin before treatment to determine the suitable type of treatment to apply, to ensure effective treatment with minimal risk of skin damage. Selecting the correct treatment type is essential for successful treatment.

- The VContour applicator is used for circumference reduction and focal fat deposit treatment in different body areas such as: abdomen, flanks, thigh, arms, saddlebags, submental etc. The applicator has three types of covers: **Large, Medium** and **Small**.
- The VSmooth applicator is used for cellulite and skin treatment in different zones such as: abdomen, thigh, buttocks, arms, etc. The applicator has one type of cover.

Table 1: Selecting the Appropriate Treatment Type

#	Treatment Type	Treatment Area	Applicator	Cover	Treatment Protocol
1	Circumferential Reduction	Abdomen ^a	VContour	Medium	Option 1: three bi-weekly sessions. Option 2: one long treatment (repeating twice in the same session).
2	Circumferential Reduction – Small Areas	Very Small Area ^b	VContour	Small	Three bi-weekly sessions
3	Localized Fatty Deposits	Abdomen ^c	VContour	Large followed by Medium	Three bi-weekly sessions
4	Cellulite	Abdomen, Thighs, Buttocks	VSmooth	VSmooth	Three bi-weekly sessions
5	Localized Fatty Deposits plus Cellulite	Abdomen, Thighs, Buttocks	VContour followed by VSmooth	Large followed by VSmooth	Three bi-weekly sessions


VelaShape III Measurement

In order to establish a baseline for the VelaShape III treatment and follow up on the results, it is advised to weigh the patient, photograph, assess the cellulite degree and measure the circumference of the treatment area. Circumference and cellulite measurements should be performed before the first treatment session and at least once more before the last treatment, and before any maintenance treatment thereafter, in order to document the results achieved. Additionally you may do so after the first two sessions (before starting the third treatment) to document any progress. Measurements after the treatment session may be inaccurate due to transient edema.

^a Circumference Reduction treatment may be performed on similarly size areas such as: thighs, arms, etc.

^b Circumference Reduction treatment may be performed on areas such as Submental.

^c Localized Fatty Deposits treatment may be performed on similarly size areas such as: flanks, thighs, middle back, etc.

2014-10-05 05:25:56

Circumference

The suggested circumference measurement method is:

1. The patient should stand in an upright position adjacent to a wall, barefoot and with both legs attached or at a fixed distance (see Figure 20 [A]).
2. A measuring ruler of hard material (wood or plastic) of at least one meter should be placed against the wall. Alternatively, a vertical metal stand with an attached marker may be used (see Figure 20 [B]).
3. Measurements should be taken for each anatomical area, usually the "thickest" part of each area:
 - Thighs and buttocks – upper third.
 - Abdomen – under and above the navel.
 - Arms – midway between shoulder and elbow.


Figure 20: Illustration of Circumference Measurements: A) Fixed Distance between the Legs; B) Labeling the circumference; C) Measuring Horizontally; D) Labeling the Arm Against a Wall.

4. During the entire measurement procedure, the distance between the patient's feet must remain constant.
5. At least four marks should be made with a white make-up pencil on each measured area, usually front, back and both sides, while the patient is rotating.
6. A measuring tape should be placed across all marks which were made in each anatomical area, creating a straight and horizontal continuous line between them (see Figure 20 [C]).
7. The measurements and the distance from the floor should be recorded in the patient folder, for identical positioning in subsequent measurements.
8. When measuring the arm circumference, indicate the distance from the wall and always measure at the same distance while the arm is perpendicular to a wall (see Figure 20 [D]).

2014-10-05 05:25:56

Cellulite

The suggested cellulite classification method is according to the **Nurnberger-Muller Scale** that classifies cellulite into four categories using the pinch test.

By gently pinching a few inches of skin between your thumb and forefinger, in the area intended for treatment, we can distinguish four (4) stages of cellulite appearance (*Hexsel et al.*):

1. **Stage Zero** – there is no alteration to the skin surface (no dimpling is seen on the skin when the patient is standing naked in front of the mirror). When the skin is pinched no orange peel effect is seen (see Figure 21-A).
2. **Stage One** – the skin of the affected area is smooth while the patient is standing or lying down (no dimpling is seen on the skin when the patient is standing), but undulations (bumps and lumps) on the skin surface can be seen on pinching the skin or during muscle contraction (see Figure 21-B).
3. **Stage Two** – the "orange peel" or "mattress" (dimpling) appearance is evident when standing, without the use of any manipulation (skin pinching or gluteus muscle contraction, see Figure 21-C).
4. **Stage Three** – presence of alteration described in stage two, plus presence of raised and depressed areas and nodules (dimpling that can be seen both when the patient is standing and lying down, see Figure 21-D).

Table 2: Cellulite stages summary according to Nurnberger-Muller Scale

	Standing	Supine	Pinch Test
Stage Zero	–	–	–
Stage One	–	–	Cellulite / Dimpling
Stage Two	Cellulite / Dimpling	–	Cellulite / Dimpling
Stage Three	Cellulite / Dimpling	Cellulite / Dimpling	Cellulite / Dimpling


Figure 21: Illustration of cellulite stages (A- Stage Zero; B- Stage One; C- Stage Two; D- Stage Three)

2014-10-05 05:25:56

Treatment Procedure

- Always inspect the applicator for cleanliness and integrity. Make sure that there are no cracks or other damage to any of the applicator's components.
- Select the parameters as suggested in the **Treatment Parameters** section of each treatment procedure. Always consider the patient's skin response and tolerance.
- Select and mark the treatment area while the patient is standing. The area dimensions will correlate to the cover aperture.
- It is essential to use VelaSpray Ease lotion, (not gel, and not any other product) for optimal coupling and RF conductivity. Moisten the treatment area with VelaSpray Ease lotion and rub it well into the skin using massage movements with gloved hands. Ensure that a thin layer of lotion (~1 mm) remains on the skin. **You may need to re-apply the lotion during treatment procedure.**
- Always test a small area with the selected parameters, and wait a few minutes to detect results of the tested area before proceeding with the treatment of the whole area.
- Orient the applicator perpendicular to the treatment area and apply minimal pressure, just enough for complete contact with the skin. When pressing the applicator's trigger or footswitch the skin will be suctioned into the applicator's chamber by the vacuum (excluding treatment where the vacuum is turned off) and energy will be delivered into the skin. If you hear air being suctioned into the applicator chamber, this indicates that the applicator is not positioned properly on the skin. If adequate coupling between the RF electrodes and the skin is not achieved, the RF and IR will be disabled and an **Improve Vacuum Seal** message will appear on the screen.
- You may manually assist the treatment by slightly stretching curved or loose skin with your free hand, and by slightly pushing or pulling the skin to or from the side towards the applicator. You may also improve contact by placing a rolled towel under the treatment area, which will push it up towards the applicator.
- The goal of each treatment session is to heat the deeper skin layers leading to a rise of the patient's skin temperature to the **endpoint**. The heat is usually accompanied by erythema (at various levels, depending on skin type and anatomical area) and sometimes slight edema starting already after four to five (4-5) minutes of treatment. The patient should feel a tolerable heat sensation. This typically occurs at a temperature of 43-45°C (110-113°F).
- The endpoint is reached when substantial and persistent radiant heat can be distinguished and the skin surface temperature measured by the IR thermometer is between 43-45°C (110-113°F). Once this endpoint is reached, this heat must be retained **for at least 10 minutes** for treatment to be effective. Therefore the minimal treatment duration of a single sub-area is 15 minutes.

**Note**

The most accurate skin temperature is read directly underneath the IR thermometer aperture.

**Applicator Positioning**

Position the applicator perpendicular to the skin and apply light but sufficient pressure to ensure good contact.

Test Spot

Before starting treatment, test a small area with the selected parameters and wait several minutes to observe the results.

Treatment Technique: Identifying the Applicator and Cover**Treatment Using the VContour Applicator**

Treatment can be performed with two different protocol modes:

- **Circumference Mode**
- **Deep Mode**

The **Circumference Mode** treatment protocol for circumference reduction supports two types of covers:

- When treating larger treated areas (such as abdomen, arms, thighs) use the **Medium Cover** as a standalone.
- When treating smaller treated areas (such as submental) use the **Small Cover** as a standalone.

The **Deep Mode** treatment protocol for focal fat deposits supports a combination of two covers. Treatment is performed using the **Large Cover** followed by the **Medium Cover** and it is used for different body areas (such as: flanks, saddlebag, middle back, etc.).

2014-10-05 05:25:56

Treatment Protocol Using the Circumference Mode

Treatment Protocol of Large Areas Using Circumference Mode (Medium Cover)

Zone Segmentation

Select the appropriate zone for treatment of circumference reduction using the **Medium Cover**; define the size of the targeted sub-area, large, medium or small, and divide it into segments:

- Large size sub-area (25 x 20 cm) will be divided into 20 segments.
- Medium size sub-area (15 x 20 cm) will be divided into 12 segments.
- Small size sub-area (10 x 20 cm) will be divided into 8 segments.

Each segment of the sub-area is the size of the **Medium Cover** (5 x 5 cm).


Figure 22: Sub-Area Divided into Segments


Figure 23: Identifying Sub-Areas and Segments for Anterior and Posterior Abdomen Circumference Reduction

**Note**

If the patient has a high BMI, additional sub-areas will be added and treated according to the same treatment protocol. Each sub-area will be treated individually in the same session.

Treatment Parameters

Applicator	Cover	Parameter	Level	Notes
VContour Applicator	Medium Cover	Vacuum	1-3 Start at level #1	May increase if needed for increased depth. Do not increase >1 for loose skin areas.
		RF	3	Maintain RF Level at 3 whenever possible to achieve an optimized result.

Cycle # 1:

Treat each segment from A1 to A20 in the **"Horizontal Snake-Like"** treatment technique (see Figure 24) with five consecutive stacking pulses, then four consecutive pulses, then three consecutive pulses and then two consecutive pulses. Continue treating each sub-area from A1 to A20 with two consecutive stacking pulses until the surface temperature reaches 43-45°C or the patient cannot tolerate further heat.


Figure 24: Cycle 1 **"Horizontal Snake-Like"**

2014-10-05 05:25:56

Cycle # 2:

In this cycle, the same sub-area set (A1-A20) will be started and treated with two consecutive stacked pulses using the "**Vertical Snake-Like**" treatment technique (see Figure 25) until the surface temperature reaches 43-45°C or the patient cannot tolerate further heat.


Figure 25: Cycle 2 "**Vertical Snake-Like**"

Cycle # 3:

In this cycle, the same sub-area set (start from sub-area A1 to A20) will be started and continued with two consecutive stacking pulses (same as in Cycle # 2) using the "**Horizontal Snake-Like**" treatment technique until the surface temperature reaches 43-45°C or the patient cannot tolerate further heat.

Cycle # 4 (Optional):

If the total time duration of Cycles # 1, 2 and 3 is less than 15 minutes, continue with the "**Vertical Snake-Like**" treatment technique (see Figure 25) and apply one pulse per sub-area until the 15-minute treatment duration for all four cycles is reached.

**Note**

When coupler agent is added over the treatment area, continue treatment when an additional stacking pulse is added to the current number of consecutive pulses used prior to adding the coupler; for example, if coupler agent is added at the stage when two consecutive pulses are being applied at each sub-area, the treatment should be continued with three consecutive pulses. When the whole area is completed with the three consecutive pulses, the treatment will continue with two consecutive pulses.

2014-10-05 05:25:56

For the medium and small size sub-areas set (including segments A1-A12 and A1-A8 respectively) treatment protocol is similar, and will be performed with the same three consecutive cycles.

Treatment Time

The large size sub-area (25 x 20 cm) set includes 20 segments (A1-A20), the medium size sub-area (15 x 20 cm) set includes 12 segments (A1-A12) and the small size sub-area (10 x 20 cm) set includes 8 segments (A1-A8). The 3-4 consecutive cycles should take at least 15 minutes.

Treatment Protocol of Small Areas Using the Circumference Mode (Small Cover)

Zone Segmentation

Select the appropriate zone for treatment of circumference reduction using the **Small Cover**, and divide it into segments.

The appropriate size of the targeted sub-area (2 x 24 cm) contains two (2) sets of eight (8) segments (1 x 3 cm).


Figure 26: Sub-Area Divided into Segments

Treatment Parameters[†]

Applicator	Cover	Parameter	Level
VContour Applicator	Small Cover	Vacuum	1
		RF	1 – 2*

(*) Start at level 1, if enabled; increase the RF to level 2 during the subsequent treatments.

Cycle #1:

Start with a sub-area set in the higher line (A1-A8) followed by the lower line of the sub-area (B1-B8). Treat each sub-area from A1 to A8 and B1 to B8 with 5 consecutive stacking pulses, then 4 consecutive pulses, then 3 consecutive pulses and then with 2 consecutive pulses. Continue treating each sub-area from A1 to A8 and B1 to B8 with 2 stacking pulse until the surface temperature reaches 43-45°C or the patient cannot tolerate further heat.

[†] Start at level 1; if possible, increase the RF level to 2 at subsequent treatments.

2014-10-05 05:25:56

Cycle #2:

A second cycle of same sets of sub-areas should be re-treated in the same manner (sub-area set A1 to A8, then sub-area set B1 to B8) with 2 stacking pulse until the surface temperature reaches 43-45°C, or the patient cannot tolerate further heat.


Figure 27: Cycle 1 and Cycle 2 in Small Sub-Areas

Treatment Time and Number of Treatment Sessions

The size sub-area (2 x 24cm) set includes 16 segments that should be treated in two cycles. Each subject will be treated with three bi-weekly treatments.

Treatment Protocol Using the Deep Mode

Treatment Protocol of Focal Fat Deposits Using Deep Mode (Large and Medium Covers)

Zone Segmentation

Select the appropriate zone for treatment of focal fat deposits using the Deep mode protocol, combining the **Large Cover** followed by the **Medium Cover**. Define the size of the targeted sub-area as 15cm x 15cm, and divide it into nine (9) segments.

Each segment of the sub-area is at the size of the **Large** and **Medium Cover** (5x5 cm).


Figure 28: Deep Mode protocol Segmentation

2014-10-05 05:25:56

"A" 49

In some cases, large areas can contain the 15cm x 15cm sub-area more than once. In this case each sub-area will be treated according to the treatment protocol.


Figure 29: Flanks (A) and Middle Back (B) Segmentation for the Treatment of Focal Fat Deposit

The **Deep Mode treatment protocol** includes two cycles, performed with the VContour applicator and the **Large Cover** without vacuum manipulation in cycle #1, while in cycle #2 the applicator is used with the **Medium Cover** and vacuum manipulation.

In order to achieve good clinical results, it is recommended to reach 43-45°C on the skin surface, maintaining it for at least 5 minutes in each step of the treatment (at least 15 minutes of treatment for the cycle).

2014-10-05 05:25:56

Treatment Parameters

Applicator	Cover	Cycle	Step	Parameter	Level
VContour	Large	Cycle #1	Step #1	Vacuum	Off
				RF	3 ^(e)
			Step #2	Vacuum	Off
				RF	3 / 2 ⁽²⁾
			Step #3	Vacuum	Off
				RF	3 / 2 ⁽²⁾
	Medium	Cycle #2	Step #1	Vacuum	1 – 3 ^(f)
				RF	3 / 2 ⁽²⁾
			Step #2	Vacuum	1 – 3 ⁽³⁾
				RF	3 / 2 ⁽²⁾
			Step #3	Vacuum	1 – 3 ⁽³⁾
				RF	3 / 2 ⁽²⁾

If the patient feels uncomfortable you may:

1. Re-apply VelaSpray Ease lotion at minimum intervals of one minute, in each step of each cycle.
2. Use the **"Quick Technique"** manner^g when stacking one pulse.
3. Only after trying both options 1 & 2 above and the patient still feels uncomfortable, reduce the RF level to 2.

^e Maintain RF Level at 3 whenever possible to achieve an optimized result.

^f Option to increase vacuum from 1 -> 2 -> 3 based on patient tolerance (in areas where there is minimal laxity).

^g **"Quick Technique"** manner: stacking one pulse in each segment, while moving the applicator faster in the **"Horizontal snake-like"** or **"Vertical snake-like"** manner according to the step.

2014-10-05 05:25:56

"A"

Cycle #1 – Step #1:

- Attach the **Large Cover** to the VContour applicator.
- Set the RF level to 3.
- Apply a thin layer of VelaSpray Ease lotion on the sub-area.
- Treat each segment from 1 to 9 in the “**Horizontal Snake-Like**” manner (see Figure 30) with four consecutive stacking pulses, then three consecutive stacking pulses and then two consecutive stacking pulses. Continue treating each segment of the sub-area (from 1 to 9) with two stacking pulses in the “**Horizontal Snake-Like**” manner until the treated area reaches two full minutes of treatment.
- After two full minutes of treatments, continue treating each segment of the sub-area (from 1 to 9) with one pulse in the “**Horizontal Snake-Like**” manner until the treated area reaches five full minutes of treatment.


Figure 30: Cycle 1-Step #1 “Horizontal Snake-Like”

2014-10-05 05:25:56

Cycle #1 – Step #2

- RF level should be maintained at level 3.
- Apply a thin layer of VelaSpray Ease lotion on the sub-area.
- Treat each segment of the sub-area (from 1 to 9) stacking two (2) pulses in the “**Vertical Snake-Like**” manner (see Figure 31). Continue treat each segment of the sub-area (from 1 to 9) stacking one pulse in the “**Vertical Snake-Like**” manner until the treated area reaches five full minutes of treatment.


Figure 31: Cycle 1 – Step #2 “**Vertical Snake-Like**”

Cycle #1 – Step #3

- RF level should be maintained at level 3.
- Apply a thin layer of VelaSpray Ease lotion on the sub-area.
- Treat each segment of the sub-area (from 1 to 9) stacking two (2) pulses in the “**Horizontal Snake-Like**” manner. Continue treat each segment of the sub-area (from 1 to 9) stacking one pulse in the “**Horizontal Snake-Like**” manner until the treated area reaches five full minutes of treatment.

Cycle #2 – Step #1

- Remove the **Large Cover** and attach the **Medium Cover** to the VContour applicator.
- Set the RF level to 3.
- Set the Vacuum level to 1.
- Apply a thin layer of VelaSpray Ease lotion on the sub-area.

NOTARY
NOTAR

- Treat each segment from 1 to 9 in the “**Horizontal Snake-Like**” manner with four consecutive stacking pulses, then three consecutive stacking pulses and then two consecutive stacking pulses. Continue treating each segment of the sub-area (from 1 to 9) stacking two pulses in the “**Horizontal Snake-Like**” manner until the treated area reaches two full minutes of treatment.
- After two full minutes of treatments, continue treating each segment of the sub-area (from 1 to 9) stacking one pulse in the “**Horizontal Snake-Like**” manner until the treated area reaches five full minutes of treatment.

Cycle #2 – Step #2

- RF level should be maintained at level 3.
- Apply a thin layer of VelaSpray Ease lotion on the sub-area.
- Treat each segment of the sub-area (from 1 to 9) stacking two (2) pulses in the “**Vertical Snake-Like**” manner. Continue treat each segment of the sub-area (from 1 to 9) stacking one pulse in the “**Vertical Snake-Like**” manner until the treated area reaches five full minutes of treatment.

Cycle #2 – Step #3

- RF level should be maintained at level 3.
- Apply a thin layer of VelaSpray Ease lotion on the sub-area.
- Treat each segment of the sub-area (from 1 to 9) stacking two (2) pulses in the “**Horizontal Snake-Like**” manner. Continue treat each segment of the sub-area (from 1 to 9) stacking one pulse in the “**Horizontal Snake-Like**” manner until the treated area reaches five full minutes of treatment.

Treatment Time

The size of the sub-area (15 x 15cm) set includes 9 segments that should be treated in two cycles of at least 15 minutes each.

Cellulite Mode Treatment Protocol Using the VSmooth Applicator

Cellulite mode is the treatment protocol of choice for cellulite and skin laxity. It is supported by the elōs technology, combining RF and IR energy with mechanical manipulation via vacuum and massage. The Treatment is performed using the **Large Electrode Cover** attached to the **VSmooth** applicator, and it can be performed on different body areas, such as: thighs, abdomen, buttocks, etc.

Treatment Protocol Using Cellulite Mode

Zone Segmentation

Select the appropriate zone and target the sub-area for treatment of cellulite and skin laxity; first define the area to be treated in **Cycle # 1** according to **Deep Mode - Zone Segmentation**. For Cycle #2 define the size of the targeted sub-area as approximately 20 x 25cm, each sub-area will be divided into four elliptical, overlapping segments as shown in Figure 32.


Figure 32: Cellulite Mode Protocol Sub-Area Segmentation

Figure 32 above demonstrates single sub-areas segmentation:

- Segment **A** – superior segment
- Segment **B** – inferior segment
- Segment **C** – left segment
- Segment **D** – right segment

The **Cellulite Mode** treatment protocol includes six steps performed using the VSmooth applicator and the attached cover.

Note

If the patient has a high BMI or the treated area is large, more than one sub-area will be treated according to the same treatment protocol. Each sub-area will be treated individually in the same session.

2014-10-05 05:25:56

- **Example # 1:** Thighs area can be divided up into three sub-areas:
 - ⇒ Anterior lateral
 - ⇒ Anterior inner
 - ⇒ Posterior
- **Example # 2:** Anterior abdomen area can be divided up into two sub-areas:
 - ⇒ Right
 - ⇒ Left

Treatment Parameters

Applicator	Cover	Step	Parameter	Level	Note
VSmooth	Large Electrode	All Steps (Step # 1 to Step # 6)	RF	3	Maintain RF Level at 3 whenever possible to achieve optimized results.
			IR	1-3 Start at 1	Start with IR level 1 and increase according to patient's tolerance from 1 -> 2 -> 3.
			Vacuum	1-3 Start at 1	May increase if needed for increased depth. Do not increase >1 for loose skin areas.
			Massage	On	

2014-10-05 05:25:56

Step #1:

- Attach the Large electrode cover to the VSmooth applicator.
- Set the RF level on 3; IR level on 1; and Vacuum level on 1 on the touch screen or on the applicator.
- Apply a thin layer of VelaSpray Ease lotion on the sub-area.
- Treat segment **A** in a circular horizontal **Clock-Wise** manner (see Figure 33) in a gliding technique, until the skin surface at the treated segment reaches 43-45°C or the patient cannot tolerate the heat.


Figure 33: Segment A Treatment Illustration

Step #2:

- Maintain the RF at level 3; IR level at 1; and Vacuum level at 1 on the touch screen or on the applicator.
- Apply a thin layer of VelaSpray Ease lotion on the sub-area.
- Treat segment **B** in a circular horizontal **Clock-Wise** manner (see Figure 34) in a gliding technique, until the skin surface at the treated segment reaches 43-45°C or the patient cannot tolerate the heat.


Figure 34: Segment B Treatment Illustration

2014-10-05 05:25:56

Step #3:

- Maintain the RF at level 3; IR level at 1; and Vacuum level at 1 on the touch screen or on the applicator.
- Apply a thin layer of VelaSpray Ease lotion on the sub-area.
- Treat segment **C** in a circular horizontal **Clock-Wise** manner (see Figure 35) in a gliding technique, until the skin surface at the treated segment reaches 43-45°C or the patient cannot tolerate the heat.


Figure 35: Segment C Treatment Illustration

Step #4:

- Maintain the RF at level 3; IR level at 1; and Vacuum level at 1 on the touch screen or on the applicator.
- Apply a thin layer of VelaSpray Ease lotion on the sub-area.
- Treat segment **D** in a circular horizontal **Clock-Wise** manner (see Figure 36) in a gliding technique, until the skin surface at the treated segment reaches 43-45°C or the patient cannot tolerate the heat.



Figure 36: Segment D Treatment Illustration


2014-10-05 05:25:56

Step #5:

- Maintain the RF at level 3; IR level at 1; and Vacuum level at 1 on the touch screen or on the applicator.
- Apply a thin layer of VelaSpray Ease lotion on the sub-area.
- Treat all treated sub-area in the **Horizontal Snake-Like** manner (see Figure 37) in a gliding technique, until the skin surface at the treated segment reaches 43-45°C or the patient cannot tolerate the heat.

**Figure 37: Horizontal Snake-Like****Step #6:**

- Maintain the RF at level 3; IR level at 1; and Vacuum level at 1 on the touch screen or on the applicator.
- Apply a thin layer of VelaSpray Ease lotion on the sub-area.
- Treat all treated sub-area in the **Vertical Snake-Like** manner (see Figure 38) in a gliding technique, until the skin surface at the treated segment reaches 43-45°C or the patient cannot tolerate the heat.

**Figure 38: Vertical Snake-Like****Treatment Time**

In order to achieve good clinical results, it is recommended that the total treatment duration for each sub-area should be:

- For the thighs - each sub-area should be treated for at least 10minutes.
- For the abdomen and buttocks - each sub-area should be treated for at least 15minutes.

2014-10-05 05:25:56

Treatment Options with Cellulite Mode

The **Cellulite Mode** treatment protocol can be used as:

- Standalone for the treatment of cellulite and skin laxity.
- Combined treatment with **Circumference Mode** for circumference reduction in addition to the treatment of cellulite and skin laxity.
- Combined treatment with **Deep Mode** for localized fatty deposits in addition to the treatment of cellulite and skin laxity.

Treatment should perform as follows (see Table 3):

Table 3: Treatment Options with Cellulite Mode

Treatment Protocol	Treated Area	# of Sub-Areas	Treatment Duration (Recommended)
Cellulite Mode	Thighs	Three (3) sub-areas: • Anterior lateral • Anterior inner • Posterior	10 minutes per sub-area.
	Abdomen	Normal Size - one (1) sub-area	15 minutes per sub-area, wait 5 minutes and repeat.
		Large Size - divide abdomen into half, two (2) sub-areas	15 minutes per sub-area, wait 5 minutes and repeat.
Circumference Mode plus Cellulite Mode	Thighs	Three (3) sub-areas: • Anterior lateral • Anterior inner • Posterior	15 minutes per sub-area with Circumference Mode. 10 minutes per sub-area with Cellulite Mode.
	Abdomen	Normal Size - one (1) sub-area	15 minutes with Circumference Mode. 15 minutes with Cellulite Mode.
		Large Size - divide abdomen into half, two (2) sub-areas	15 minutes per sub-area with Circumference Mode. 15 minutes per sub-area with Cellulite Mode.
Deep Mode plus Cellulite Mode	Thighs	Three (3) sub-areas: • Anterior lateral • Anterior inner • Posterior	15 minutes per sub-area with Deep Mode (Cycle #1) **. 10 minutes per sub-area with Cellulite Mode.
	Abdomen	Normal Size - one (1) sub-area	15 minutes with Deep Mode (Cycle #1). 15 minutes with Cellulite Mode.
		Large Size - divide abdomen into half, two (2) sub-areas	15 minutes per sub-area with Deep Mode (Cycle #1). 15 minutes per sub-area with Cellulite Mode.

(*) Perform Deep Mode Cycle #1 only in the zones that require treatment.

(**) When choosing to treat **Localized Fatty Deposits and Cellulite with Deep Mode followed by Cellulite Mode** on the thigh area: the sub-area, anterior inner thigh, should be treated with 15 minutes of Cellulite Mode only (Deep Mode - Cycle #1 will be skipped).

2014-10-05 05:25:56

Endpoints (When to Stop Treating an Area)

- Use the radiant thermometer to measure the temperature on the skin surface, skin temperature should be 43-45°C.
- If the patient feels an intolerable heat sensation before reaching 43-45°C, try to re-apply some VelaSpray Ease lotion, if patient still cannot tolerate the heat, lower the RF and continue treatment. If the patient tolerates this well, complete the indicated pulses and phases.
- Mild erythema is normal.
- Heat should be persistent.

Warning

In the event of excessive heat sensation or any side effects (i.e. burn), stop treatment and apply cooling.

Caution

Release the trigger or footswitch **before** you lift the applicator from the skin.

Note

- Compliance with the treatment protocol is essential to achieving the desired results.
- Always maintain complete contact between the applicator chamber and the skin.
- For small zones such as the navel area, treat twice, but with less stacks until the endpoint is reached, and ensure a 5-minute interval between treatments.
- When reinforcing results, treat the small area (such as localized fat deposits) only once.
- After treatment, the treated area may show temporary signs of erythema and slight edema.
- A **Bad Coupling** message on the control panel indicates inadequate coupling, which results in RF being disabled.
- VelaShape III treatment performed 3–5 weeks post-lipolysis can enhance smooth recovery from the lipolysis treatment.

2014-10-05 05:25:56

Post-Treatment Care

- If heat sensation is excessive, cool the area with conventional cooling methods such as cold (not frozen) packs, air cooling with SynerCool, etc.
- In case of side effects such as burns, stop treatment; apply cooling and appropriate creams. Resume treating the area only after the side effects have subsided, and with reduced parameter levels.
- Avoid post-treatment overheating (e.g., Jacuzzi, sauna, etc.), or risk of physical damage.

Treatment Conclusion

- Treatment should be concluded when the results are satisfactory to the patient or according to the operator's discretion.
- Improvement will take some time to be apparent, and is affected by clinical conditions at the start of the treatment regimen. Main impacts may be:
 - ⇒ Transient skin texture that may show improvement after each session.
 - ⇒ Circumference reduction may show improvement after 2–3 sessions.
 - ⇒ Cellulite may show improvement after 4 sessions.
- Take photos and measurements of the treatment areas after the last session.

7. MAINTENANCE

While the VelaShape III system should not be serviced by anyone other than an authorized Syneron representative, the system requires regular preventive maintenance and cleaning to ensure optimal system performance. This chapter describes the routine maintenance procedures that operators may perform:

- Cleaning the VelaShape III system.
- Cleaning the VelaShape III applicator and its replaceable covers.
- Cleaning the air filter.
- Replacing the applicator's replaceable covers and filters when required

Note

Each patient should be assigned with **his/her own set of covers** in order to provide safe and hygienic treatments. This also allows the patient to keep track of the number of treatment hours administered.

Warning

The  symbol located near the power cable's connection port indicates:

WARNING! DISCONNECT POWER SUPPLY CABLE BEFORE MAINTENANCE OR SERVICE PROCEDURES TO AVOID ELECTRIC SHOCK.

For continued protection against fire, replace the fuses only with ones of the same type and rating.

Cleaning the Console

Clean the system at least once a week. Turn the system off and wipe all surfaces with a soft, damp, non-abrasive cloth.

Clean the umbilical cable with pre-moistened wipes.

Warning

Be careful not to spill any liquids on the system.

Cleaning and Maintaining the VSmooth Applicator

Caring for the VSmooth applicator involves the following maintenance procedures before and after each treatment:

- Removing, cleaning and replacing the replaceable cover.
- Cleaning the interior of the applicator.
- Replacing the applicator's filters.

Caution

While gently cleaning the interior of the applicator, avoid moisture and dirt from getting into the applicator's filters.

The replaceable cover of the VSmooth applicator should be replaced according to the schedule printed on the new cover's box label, or sooner if one of the following is observed:

- Rollers appear damaged, discolored or their movement is impaired.
- Roller axis appears to wear out (color turns from silver to copper).

Proceed as follows:

1. Turn the system off and disconnect the replaceable applicator cover by pressing the two spring-loaded buttons, then pull the cover straight off the applicator body, as shown in Figure 39.


Figure 39: Removing the VSmooth Replaceable Cover

**Note**

The image in the figure is for illustration purposes only.

2. If it becomes necessary to clean the applicator cover during treatment:
 - While the system is on, clean the rollers of the VSmooth applicator while pressing the trigger, with all modes at level **0**, while rolling over a mild disinfectant-dampened tissue (that does not contain alcohol).
 - Turn the system off and remove the cover. Use hot tap water (no more than 122°F/50°C). The cover should be flushed and washed thoroughly in hot water and mild disinfectant soap or detergent that do not contain alcohol, until all residual VelaSpray Ease is removed. Wipe and dry the cover with a soft cloth. Reconnect the cover only when completely dry.
3. Clean the VSmooth applicator's interior components at the end of every day as follows (see Figure 40):
 - Replace the two filters (see Figure 41) according to the schedule printed on the cover's box. Never use a pair of filters for more than 10 accumulative hours.
 - Clean the RF contacts and ID contacts with a soft cloth.
 - Gently clean the IR lamp aperture using a gentle optical lens cleaning solution, such as Fisherbrand® Optical Lens Cleaner or equivalent; be careful not to scratch it.
 - Gently dry the IR lamp aperture with lint-free lens cleaning tissue.
 - Reconnect the clean and dry (or new) applicator's cover.


Figure 40: VSmooth Applicator's Internal Components

2014-10-05 05:25:56

**Caution**

Verify that the vacuum's aperture does not become clogged during cleaning. Avoid inserting spray residue into the aperture.

Replacing the VSmooth Applicator's Filters

Refer to Figure 41:

1. Use the supplied filter extractor to pull the two filters out of the body of the applicator, as shown in view **A**. Discard the used filters.
2. Unpack a pair of new filters and firmly insert them into the filter chambers, as shown in view **B**.


Figure 41: Replacing the VSmooth Applicator's Filters

**Caution**

Make sure that the interior of the applicator is completely dry before inserting the filters.

Caution

Make sure that the interior of the applicator and the cover are completely dry before inserting the filter and reconnecting the cover.

Cleaning the VContour Applicator's Filter

Remove the filter from the cover if it becomes necessary to clean or replace it.

1. Use the supplied filter extractor to pull the filter out of the cover, as shown in Figure 43.
2. Wash the filter and allow it to dry completely before reusing it.
3. Firmly insert the clean filter into the filter chamber.


Figure 43: Replacing the VelaShape III Applicator's Filter

Note

The image in the figure is for illustration purposes only.

Cleaning the System Air Filter

While the system air filter should be cleaned at least once a week, the filter cup needs to be checked daily and cleaned according to the amount of dust and debris that accumulates.

To clean the system air filter (see Figure 44):

1. Turn the system off.
2. Open the cover of the system's rear utility cabinet to access the filter assembly.
3. Release Screw # 1 on the bottom of the filter cup to open the cup and wipe off the dust and debris.
4. Release Screw # 2 to remove the filter, wash it with water and soap and set it aside to dry.
5. When the filter is completely dry, fit it back into place, or put a new filter and secure the filter cup using the screws.
6. The same filter may be re-used after washing. It should be discarded when it is no longer white after cleaning.

Caution

Be sure that the filter is completely dry before re-assembling.


Figure 44: Cleaning the Air Filter

2014-10-05 05:25:56

8. TROUBLESHOOTING GUIDE

The VelaShape III system is equipped with a self-testing software routine that continuously monitors system operation. If a system error is detected, the modules are disabled until the problem is solved.

The following troubleshooting guides do not attempt to list all possible system failures. Any error not listed should be referred to your Syneron representative.

Table 4 offers a list of possible system errors that will be displayed on the control panel's display, and how to correct them. If the corrective actions listed in the table do not solve the problem, contact your Syneron representative.

Table 5 offers several fault situations that may occur with no operator panel indication or control panel display, and their appropriate corrective actions. If these do not resolve the problem, contact your Syneron representative.

Most times a system error can be cleared by restarting the VelaShape III system. Always try this remedy before contacting your Syneron representative.

2014-10-05 05:25:56

Table 4: System Errors Displayed by Terminal Messages and/or Operator Panel Indicators

Message	Symptom	Corrective Action
Improve Vacuum Seal	Vacuum is not implemented correctly upon the skin.	- Remove and clean the applicator cover. - Replace the cover filter. - Put on a new cover. - Make sure the patient skin is cleaned and shaved.
RF Out of Range	The system cannot apply RF correctly.	- Make sure there is a good vacuum seal between the applicator cover and patient skin, - Make sure there is an adequate amount of gel. - Make sure the applicator is placed on the skin appropriately.
No applicator, No Tip	System freezes, no RF or IR pulses emitted to the patient's skin even if footswitch/trigger is pressed.	- Let go of footswitch/trigger and make sure that the applicator and tip are well connected and sealed. - Let go of footswitch/trigger and wait 10 minutes for the system to recognize the applicator and tip. - If after 10 minutes the system does not recognize the applicator and tip restart the system. - If after restart the system does not recognize the applicator and tip unplug other medical devices in the room and restart the system again.

2014-10-05 05:25:56

Table 5: Other Fault Situations with No Error Messages

Symptom	Corrective Action
No vacuum suction	1. Turn the system off. 2. Check that air filter is clean and the filter cover is properly fastened (the cover's sealing O-ring should be properly installed). 3. Check the applicator connection to the umbilical cable, and the umbilical cable's connection to the system. 4. Verify that the cover is correctly attached to the applicator. 5. Replace applicator filters or cover. 
No IR output	1. Turn the system off. 2. Call your Syneron representative for further instructions.

2014-10-05 05:25:56

9. SYSTEM SPECIFICATIONS

Output Parameters

<i>Console</i>	
Peak RF Power	150 Watts
RF Frequency	1 MHz
RF Mode	Bi-Polar
Light Wavelength	850 nm
Vacuum	(-200) to (-390) Millibars Pulsed
<i>VContour Applicator</i>	
IR Electrical Power	Up to 1.7 W
RF Power – Small Cover	Up to 10 W (average per duty cycle)
RF Power – Medium Cover	Up to 50 W (average per duty cycle)
RF Power – Large Cover	Up to 75 W (average per duty cycle)
<i>VSmooth Applicator</i>	
IR Electrical Power	Up to 3.2 W
RF Power	Up to 75 W (average per duty cycle)

Treatment Area Footprints

<i>VContour Applicator</i>	
Large Cover	25 mm x 50 mm
Medium Cover	25 mm x 30 mm
Small Cover	10 mm x 13 mm
<i>VSmooth Applicator</i>	
Cover with Rollers	40 mm x 40 mm

Control Panel

- Tablet computer with 10.1 inch, high-definition GUI display
- Android operating system

Electrical Requirements

- 100-230 VAC; up to 4A; 50-60 Hz; single phase

2014-10-05 05:25:56

Environmental Requirements

Treatment Room

- Temperature: 15° – 30°C [59° – 86°F]
- Relative humidity: Up to 80%
- Air pressure range: 70 – 106 kPa

Storage/Transportation

- Temperature: 10° – 55°C [50° – 131°F]
- Relative humidity: 10% – 80%
- Air pressure range: 70 – 106 kPa

Degree of Protection Against Ingress of Water

- IPX0 – Ordinary Equipment

System Console Dimensions [W x D x H]

- 38 x 49 x 132 cm / 15 x 19 x 51.8 inches

Weight

- VelaShape III system: 20 kg. / 44.1 lbs.
- VContour Applicator: 0.4 kg. / 0.9 lbs.
- VSmooth Applicator: 1.0 kg. / 2.2 lbs.

2014-10-05 05:25:56

10. ELECTROMAGNETIC REGULATION

Electromagnetic Emissions

The VelaShape III system is intended for use in the following electromagnetic environment. The owner of the system and/or the user must ensure that it is operated the correct environment.

<i>Emissions Test</i>	<i>Compliance</i>
CISPR 11: RF Emissions	Group 1

The VelaShape III system must emit electromagnetic energy in order to perform its intended function. Nearby electronic equipment may be affected.

<i>Emissions Test</i>	<i>Compliance</i>
CISPR 11: RF Emissions	Class A
IEC 61000-3-2: Harmonic Emissions	N/A
IEC 61000-3-3: Voltage fluctuations/flicker emissions	N/A

The VelaShape III system is suitable for use in all establishments other than domestic, and may be used in domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, provided the following warning statement is heeded:

Warning

This device is intended for use by health care professionals only. This device may cause interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigating measures, such as re-orienting or relocating the device or shielding the location.

2014-10-05 05:25:56

Electromagnetic Immunity

Immunity Test	IEC 60601 Test Level	Compliance
IEC 61000-4-2: Electrostatic Discharge (ESD)	± 6 kV contact ± 8 kV air	$\pm 2, 4, 6$ kV Contact $\pm 2, 4, 6, 8$ kV Air
IEC 61000-4-4: Electrical Fast Transient/Burst	± 2 kV for power supply lines	± 2 kV
	± 1 kV for input/output lines	N/A
IEC 61000-4-5: Surge	± 1 kV differential Mode	1.0 DM
	± 2 kV common mode	2.0 CM
IEC 61000-4-8: Power Frequency (50/60 Hz) Magnetic Field	3 A/m	3 A/m
IEC 61000-4-11: Voltage Dips, Short Interruptions and Voltage Variations on Power Supply Input Lines	$<5\% U_T$ ($>95\%$ dip in U_T) for 0.5 cycle	$<5\% U_T$ ($>95\%$ dip in U_T) for 0.5 cycle
	$40\% U_T$ (60% dip in U_T) for 5 cycles	$40\% U_T$ (60% dip in U_T) for 5 cycles
	$70\% U_T$ (30% dip in U_T) for 25 cycles	$70\% U_T$ (30% dip in U_T) for 25 cycles
	$5\% U_T$ (95% dip in U_T) for 5 cycles	$5\% U_T$ (95% dip in U_T) for 5 cycles

- Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
- Mains power quality should be that of a typical commercial or hospital environment.
- If the VelaShape III system is required to operate during power interruptions, the system should be powered from an uninterruptible power supply.
- Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

**Note**

U_T is the AC mains voltage prior to application of the test level.

2014-10-05 05:25:56

Electromagnetic Immunity (continued)

Immunity Test	IEC 60601 Test Level	Compliance	Guidelines
61000-4-3: Radiated RF	3 V/m 80 MHz to 2.5 GHz	3 V/m 80 MHz to 2.5 GHz	$d = 1.17 * \text{SQRT}(P)$ 80 MHz to 800 MHz $d = 2.33 * \text{SQRT}(P)$ 800 MHz to 2.5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ¹ , should be less than the compliance level in each frequency range ² . Interference may occur in the vicinity of equipment marked with this symbol: 
61000-4-6: Conducted RF	3 V _{RMS} 150 kHz to 80 MHz	3 V _{RMS} 150 kHz to 80 MHz	Portable and mobile RF communications equipment should be used no closer to any part of the VelaShape III system, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.17 * \text{SQRT}(P)$

(1) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the VelaShape III system is used exceeds the applicable RF compliance level above, the VelaShape III system should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the VelaShape III system.

(2) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

**Note**

- At 80 and 800 MHz, the higher frequency range applies.
- These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Recommended Separation Distances

The following are the recommended separation distances between portable and mobile RF communications equipment and the VelaShape III system.

The VelaShape III system is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The owner of the system and/or the user can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the VelaShape III system as recommended below, according to the maximum output power of the communications equipment.

For transmitters rated at a maximum output power not listed below, the recommended separation distance (*d*) in meters (M) can be estimated using the equation applicable to the frequency of the transmitter, where *P* is the maximum output power rating of the transmitter in watts (W) according to the transmitter's manufacturer.

Rated Maximum Output Power of Transmitter (W)	Separation Distance According to Frequency of Transmitter (M)		
	150 KHz to 80 MHz <i>d</i> = 1.17 SQRT(<i>P</i>)	80 MHz to 800 MHz <i>d</i> = 1.17 SQRT(<i>P</i>)	800 MHz to 2.5 GHz <i>d</i> = 2.33 SQRT(<i>P</i>)
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.69	3.69	7.38
100	11.67	11.67	23.33
150	14.32	14.32	28.53

Note

At 80 and 800 MHz, the higher frequency range applies.

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