Instructions for Use of Expanders

The physical properties of the POLYTECH Health & Aesthetics implants have been optimized as to closely resemble body tissue, thus allowing implants to be ideal for use in aesthetic, corrective or reconstructive surgery. Every type of implant is available in different sizes and shapes (see current price list). The present instructions for use applies to temporary implants (expanders): expanders equipped with a valve are delivered empty and are gradually filled intra- and postoperatively after their implantation. They have a smooth or textured surface. For more information on the shell and the surface, please refer to the chapter "Information on materials". Design as well as production, delivery and customer service are subject to our quality management control according to EN ISO 9001 / EN ISO 13485. The products are CE-certified and comply with the requirements of the European Medical Device Directive (MDD) and other national regulations.

All materials used are tested for biocompatibility and conform to the medical quality standards. Furthermore, using different standard tests, the materials' mechanical properties such as elongation, tear resistance and fatigue have been investigated.

INFORMATION ON MATERIALS

Implant shells

POLYTECH implants and expanders have a chemically and mechanically resistant, soft silicone elastomer shell which consists of several polydimethylsiloxane layers. The shell surface is smooth or textured. Wall thickness of the shell depends on the size of the product: for smooth shells it is approximately 0.3 – 0.9mm, for textured (POLYtxt®) shells it is roughly 0.45 – 1.2mm. The average surface roughness of textured shells is approx. 250µm. Mechanical stability: elongation $\geq 450\%$, tensile strength $> 5\text{MPa}$, tensile set $\leq 10\%$; Strength of seams and seals must not fail when elongated to 300% and held for a period of 10s.

Filler material: for the filling of expanders, we recommend to use a sterile, pyrogen-free, 0.9% physiological saline solution.

Allowable maximum filling of the expander for the safe use of it after implantation is not more than two-fold excess of the nominal volume.

Expanders, sterile

Packaging includes:

1. Instruction on use for expanders - 1 (ea);
2. Instruction on use for Delivery Sleeve - (1ea);
3. Labels – 8 (ea);
4. Expander passport – 1 (ea);
5. Delivery sleeve for POLYtxt implants , Reference ID 72001 -1(ea) (if necessary);
6. Magnet finder for expanders with integrated valve , Reference ID 20738-012 – 1 (ea) (if necessary);
7. Injection valve with silicone lengthened tube and connector for expanders with remote valve. Reference ID 1409 – 1 (ea) (if necessary);

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Materials and technical characteristics of components:

1. For textured expanders: 1 polyethylene (open-ended) insertion sleeve.

Low density polyethylene, CLEANFLEX®; Technical characteristics of the sleeve:

Length: 370 mm ± 10 mm, width: 165 mm ± 10 mm; thickness: 0.05 mm ± 0.01 mm,

weight : 4.1 g - 7.4 g, tensile strength : warp direction ≥ 15 MPa; cross direction ≥ 10 MPa; Manufacturer: Polytech Health & Aesthetics GmbH, Germany;

2. Magnet Finder for Expanders with integrated valve

Following implantation via a magnet finder (polycarbonate Makrolon 2405) the physician determines the location of the integrated valve, which is equipped with a permanent magnet. Magnet finder length is 122mm ± 2 mm, weight is not more than 6

grams. Integrated valve (inner diameter is $18\text{mm} \pm 2\text{ mm}$, outer diameter is $50\text{mm} \pm 2\text{ mm}$) fixed on the front surface side of the expander. Integrated valve consist of silicone elastomers MED- 4830 and MED-4870 manufactured using an injection moulding procedure. Integrated magnet valve (Nd2Fe14B) is placed into housing of stainless steel material number 1.4401;

Manufacturer: Polytech Health & Aesthetics GmbH, Germany;

3. Injection valve with silicone lengthened tube and connector for expanders with remote valve. Remote valve (inner diameter is $17\text{mm} \pm 2\text{ mm}$, outer diameter is $34\text{ mm} \pm 2\text{ mm}$) consist of the silicone elastomer MED-4830 manufactured using an injection moulding procedure. To avoid that the injection needle pierces the valve when filling the expander, a disk made of titanium (material numberb3.7035) has been integrated in the valve. The valve is connected to the expander by a silicone tube with connector made of PVDF HD1000 (polyvinylidene fluoride). Length of tube is $50 \pm 2\text{ mm}$, inner diameter is $1.58 \pm 2\text{ mm}$, outer diameter is $3.18 \pm 2\text{ mm}$, silicone MED-4755, Nusil. One tube end is glued to a silicone flange (MED-4830), which is fixed inside the shell. Injection valve weigh is not more than 9 g. Manufacturer: Polytech Health & Aesthetics GmbH, Germany;

INDICATIONS AND POSSIBLE CONTRAINDICATIONS

Indications: temporary extension of tissue and skin (e.g. for reconstruction after ablation mammae, burns, posttraumatic loss of tissue) and correction of defects (e.g. Poland Syndrome). We recommend to replace expanders by final implants after 6 months implantation at the latest.

Possible contra-indications

General clinical contra-indications In addition to general invasive measures can include: patients with mental instabilities, repeated attempts or failures of contour correction, clinically persistent infections or systemic diseases , expected allergies or extraordinary immune response to implants abscesses, cysts or tumours in the area of expanders, especially cancer and recurrent metastases, advanced fibrocystic diseases, severe radiation injuries in the prospective area of expanders, heavy burn scars in the prospective area of implants Insufficient tissue covering in the prospective area of implantation (e.g. after preceding breast reduction) or reduced vascularization, wound healing impairments, existing costal injuries, pregnancy, breast feeding (concerning implantation in the area of mammae), pace-makers in case of skin expanders with

integrated magnet valves, MRI (magnetic resonance imaging) examinations scheduled for the implantation period of the expander

Custom-made implants

Customer-specific expanders (normally mouldbased) are subject to the same quality standards as applied to the above mentioned products. Any information given in this brochure also refers to custom-made implants.

TRANSPORT AND PACKAGING

Products are delivered individually packed using sterile double packaging systems. Care should be taken that the inner package is neither opened nor damaged before the implant is to be definitely implanted. Any sterile product packaging is provided with self-adhesive labels showing manufacturer, volume or size, order and serial numbers. Labels should be affixed to the patient's records and to the implant passport.

Each package includes:

- Double sterile blister packaging with one expander
- Delivery sleeve (for POLYbxt implants) Reference ID 72001 ;
- Magnet finder for expanders with integrated valve , Reference ID 20738-012 ;
- Injection valve with silicone lengthened tube and connector for expanders with remote valve. Reference ID 1409 ;
- Instruction on use for expanders
- Instruction on use for Delivery Sleeve
- Label (8ea)
- Expander passport

Expander passport includes the following information and is completed by physician:

Symbols:

Order number

Serial number

„Use-by date“ (month/year)


Sterilisation with steam or dry heat

single use only

Consult instructions for use



Humidity limitation (0% - 70% rH)

Upper limit of temperature: 35°C maximum for storage

Dimensions A, B, C and if applicable dimension D; as shown in the schematic drawings.

CE marking

STERILIZATION, STORAGE, SHELF-LIFE

WARNING: The expanders are delivered in a sterile double-packaging. In the operating theatre, product sterility is ensured as long as the sterile inner package has neither been damaged nor opened. A replacement expander should be available during implantation surgery. Product sterility is not guaranteed when the sterile packaging is damaged. We do not recommend to use such an expander. When in doubt, please contact POLYTECH Health & Aesthetics or your local distributor.

Products are sterilized using dry heat. Sterility is guaranteed for 5 years from sterilization date as long as the sterile double packaging has not been compromised (see label for maximum storage period ☒ YYYY-MM).

Sterility of products is not guaranteed with damaged packaging. In the operating theatre, product sterility is ensured as long as the sterilized inner package has neither been damaged nor opened.

The components are placed in a standard cardboard packing box with expanders.

Components

Delivery sleeve: delivered sterile by gamma irradiation, disposable.

Magnet finder is delivered non sterile, reusable. For cleaning it can be used usual antiseptic spirituous wipe.

Injection valve with silicone lengthened tube and connector: is delivered sterile by dry heat, disposable.

Materials of packaging for components:

Delivery sleeve is packed inside a double sterile packaging (two sealed pouches) made from Perfecseal, supplier is Perfecseal Ltd, Germany.

Magnet finder is packed inside a polyethylene packaging Hekofolit, supplier is Hekofolit Heinz Koch GMBH, Germany.

Injection valve with silicone lengthened tube and connector is packed inside a peel off pouch (made from Stericlin, supplier is Stericlin, Germany).

Components are packed into standard cardboard box together with expanders.

Products must be stored under dry conditions at room temperature so as not to exceed a relative ambient humidity of 70% and a temperature of 35°C.

Transportation conditions: $-29^{\circ}\text{C} \pm 2^{\circ}\text{C}$ and $60^{\circ}\text{C} \pm 2^{\circ}$, relative humidity
 $30\% \pm 5\% - 85\% \pm 5\%$.

Since the product packaging is gas-permeable, do not allow exposure of the product package to steam or gas. Implants should be stored in horizontal position so that labels can be easily read. Products are made for single use only and must neither be reused nor resterilized. This is due to the following: the surface quality does not permit safe decontamination and no adequate decontamination procedure can be recommended., the products are low pressure / vacuum-sensitive and may not be subjected to the usual sterilization procedures.

- The product integrity especially of the shell is subject to a maximum of mechanical stress during implantation and explantation and cannot be sufficiently rated by the user.
- The mechanical safety limits as required by international standards for new products cannot be guaranteed after implantation / explantation of the product because of the mechanical forces applied during the procedure. Therefore, the required product safety cannot be guaranteed at reuse.

PATIENT INFORMATION

Every patient who has to undergo surgery, or any person authorized to act in his / her name on the occasion of surgery, must be informed by the attending physician about any operative or postoperative risks in general as well as those related to expanders and implants. Patients should be advised in a professional and very responsible way so that no unrealistic expectations in terms of possibilities and cosmetic results can arise. It is necessary for patients to know that implants are not intended for life and may have to be replaced after some time. Expanders are for temporary implantation only; we recommend a maximum implantation time of 6 months.

Patients should be invited to consult their physicians at regular intervals (every 6 to 12 months) to have the implants and expanders including surrounding tissue checked. Also, they should be informed that abnormal strain, excessive movement or traumatization of surrounding tissue can result in implant rupture under the skin and, therefore, in a premature reoperation and explantation.

The physician in charge should advise patients on how to distinguish the expander (and the valve in case of expanders with remote valve) from tissue after surgery and how to

examine for tumours themselves. Also it is advisable to inform patients that magnet valves integrated into products may interfere with data carriers, computer screens and watches.

Expander passport

Expander passports are attached to each POLYTECH Health & Aesthetics GmbH product. This passport should be completed and given to patients to ensure availability of implantation date and implant type in case of need.

Screening diagnostics

Silicone implants/expanders can interfere with screening techniques and have a falsifying effect on results. It is, therefore, necessary that patients inform the medical staff involved about the existence of implants or expanders.

Patients having expanders with integrated valves should not undergo magnetic resonance imaging (MRI) examinations as there may be an interference with magnetic fields which may impair patients' health and lead to false MRI results.

ATTENTION: All of your patients who have an expander with integrated valve and permanent magnet have to be advised without fail that before undergoing an MRI (magnetic resonance imaging) examination it is absolutely compulsory to inform the corresponding physician of this magnetic valve.

POSSIBLE COMPLICATIONS

Possible complications associated with expanders may be due to individual anamnesis, medication, operation technique, and postoperative care. Intra- and postoperative risks related to implants and expanders are specified below:

- **Skin necroses, sloughings and extrusions** may occur due to excessive tension of the skin overlaying the implant (e.g. oversizing the implant in relation to the implant space or accelerated expander-triggered tissue extension) or by interrupted or reduced skin vascularization (e.g. traumatization of skin flap during the surgical intervention). Treatment with steroids may also provoke necroses. Smokers are high-risk patients as blood circulation in the overlaying skin may be impaired. Revisions may be necessary.
- **Haematomas** may occur prior to an infection, increasing the risk of capsular fibrosis, and therefore shall be avoided or treated immediately. Large postoperative oedemas or haematomas may, if untreated, provoke extrusion or

displacement of implants. Non-resorbable haematomas may require an immediate explantation. Meticulous haemostasis during a surgery is the best prevention.

- **Seromas**, due to traumatization or excessive postoperative mobility, may cause swellings and pain. Possible therapies: immobilization, compression, drainage or explantation, if necessary.
- **Excessive swellings** may be caused by serous fluid (referred to "seromas").
- **Displacement or continuous pain** may be due to improperly sized or placed implants. Over-sized implants, capsular contracture as well as irritations due to excessive implant movement may provoke such symptoms. Pain may occur in the operated area as well as in pectoralis and shoulder-arm areas after mammary surgery.
- **Folding** (also referred to as "wrinkling" in medical literature): palpable implant edges or distortions may be due to a discrepancy between implant size and cavity, insufficient tissue covering.
- **Infections** in connection with breast implants are very scarce (0.114%, Brand, 1993). Infections with unclear etiology after implantation should be treated immediately. The use of antibiotics, drainage or explantation can be necessary. A culture of the wound secretion should be made. Typical germs are: *Staphylococcus aureus*, *Staphylococcus epidermis* or *Streptococcus*. Not all infections can be treated with the product remaining implanted. The "toxic shock" syndrome has been reported in single cases in connection with breast implants.
- **Reddening of the skin**, in the literature described as "rash", can be observed with specific indications in different frequencies, depending on the surface structure (approx. 1–5% for Microthane®-surfaced mammary implants). This reddening of the skin should not be confused with an infection. It differs from an infection by itching and the absence of systemic infection symptoms. It usually occurs 7 to 10 days after the implantation and can last 2 to 3 weeks. The use of steroids is recommended in literature for the case that a treatment is absolutely necessary.
- **Capsular contracture / fibrosis**: formation of a fibrous capsular around an implant is a normal host response to the implant. Form and position of the implant can change due to a capsular contracture. A clinically relevant (Baker III, IV) contracture may incur aesthetic impairment and pain. An explantation may become necessary. A closed capsulotomy is not recommended since this may

damage the product. If such a capsulo-tomy cannot be avoided, a screening procedure should be employed to investigate possible implant damage.

- **Desensibilization or hypersensibilization** may occur transiently or permanently (e.g. totally or partially of breast or scrotum). There may be changes in sensibility in the presence of mammary implants, especially when a peri- or transareolar incision was made. Sometimes, a transient change in sensibility is reported when using facial implants.
- **Permeation of saline solution:** silicone elastomers are permeable to gas and steam and hence to saline solution. Permeation of saline solution through the expander shell may cause a loss of liquid and a decrease of volume.
- **Implant rupture or implant leakage** are absolute indications for an explantation. Damaging of the expander (intraoperatively, e.g. by too short incisions, by surgical instruments, or postoperatively, e.g. when filling expanders, in case of punctures, biopsies, closed capsulotomies etc.), traumata or material fatigue are, in our opinion, the most frequent reasons. Deflation of expanders: if the recommended 6-month implantation period is exceeded, the expander shell may fold – as a possible side-effect of the expansion phase – which may weaken the shell and lead to a deflation of expander.
- **Overfilling** of expanders may lead to ischaemia and other complications (e.g. necrosis, extrusion). Filling of expanders should be adapted to individual conditions.
- **Cancer risk:** epidemiologic studies of women having implants reveal that breast cancer risk for them is not higher if compared with control groups. Degradation of micropolyurethane-foam may, under certain conditions, incur a release of 2,4 toluenediamine (2,4 TDA) concentrations which are, however, extremely low.
- **Benign calcification** around the implant is possible.
- **The ability of breast feed** may, irrespective of breast implants, be affected by breast surgery.
- There is no evidence that breast implants/expanders could be detrimental to the health of infants.
- **Anaplastic large cell lymphoma (ALCL):** Investigations by the American FDA have shown a possible link between breast implants and the development of ALCL. The FDA believes that women with breast implants may have a very small

increased risk to develop ALCL. It is recommending that the responsible physicians inform their patients accordingly. There are no data for expanders.

SURGICAL TECHNIQUES

Products, as supplied, are exclusively intended for implantation by skilled surgeons.

The physician in charge should be experienced in filling expanders in order to avoid hurting the patient or damaging the product.

In no case should expanders be modified prior to or during implantation. Before implantation, each product should be checked for possible damage. Before implanting an expander, the shell integrity has to be tested, if need be, by partially filling it with sterile saline solution and applying moderate hand pressure to the shell. Damaged or modified products must not be implanted!

A replacement expander should be available during the implantation!

A direct combination of the expander with a medicinal preparation has to be avoided. Implantation may be performed by a variety of surgical techniques. It is only up to the surgeon to decide on the technique which should be in line with the medical requirements as well as the individual needs of patients. Patients should be examined prior to surgery in order to guarantee sufficient, tension-free tissue covering above the implant and the incision. Pressure, traction or other stress factors, to which implant and incision area may re-act, must be avoided.

Implantation must be performed carefully. Implants can be damaged by excessive manual strain, manipulation or surgical instruments. Should the elastomer shell be damaged, filler material might leak from the implant.

We recommend to use a polyethylene sleeve (as part of the delivery) for implantation of textured expanders (refer to "Dolsky Technique").

Following recommendations should be observed:

- Incisions should be made so as to correspond with the implant size and surface.
- Haemostasis should be carried out meticulously. Cauterizers may damage the implant shell.
- Implant size should be selected carefully. It is recommended to store different sizes of implants in order to have the appropriate one available. We recommend using sizers to determine the appropriate size during surgery (see current price list).

- Implant cavity should correspond to implant size. Implants should be placed flatly and without wrinkles.
- Appropriate tissue covering must be ensured or, if necessary, preference be given to submuscular placement of implants.

Prevention of contamination

The implant surface must not be contaminated with talcum, dust or skin oils. Cotton fibres, talcum or other surface contaminants can be attracted by electrostatic charging and, thus, provoke an increased foreign body reaction. Extraordinary care has to be taken to prevent contamination. Medical staff involved should, therefore, clean sterile gloves with sterile saline solution before touching the implant.

Recommended passing procedure

1. Open outer package under aseptic conditions. Remove outer sterile packaging. Attach self-adhesive product label to patient's records.
2. Open outer sterile packaging. Offer inner sterile package to the nurse for putting it on the sterile instrument table under sterile conditions. The same handling procedure is followed for separately packed components, if provided.
3. Open inner sterile package under sterile conditions shortly before implantation. Prepare the sterile implant according to the requested surgical technique and hand it to the surgeon.

Filling of expanders

ATTENTION: Expanders are intended for temporary implantation only. We recommend a maximum implantation time of 6 months.

Expanders are filled by injecting through a valve sterile, pyrogen-free physiological (0.9%) saline solution using a 23-gauge needle (e.g. Sterican® article no. 4657667 and 4665635, B.Braun, Melsungen). Also, cannulas with a smaller diameter, i.e. gauge 24, 25 etc., can be used. Inside the self-sealing valve is a cannula stop serving as a puncture barrier.

Prior to the implantation, integrity of shell and functionality of valve should be checked by filling the expander partially and by slightly pressing upon the shell. Damaged or modified products must not be implanted. Air bubbles, if any, should be aspirated before the implantation. First or subsequent filling volumes are dependent on the individual anatomical needs. Successive filling intervals should be agreed between physician and patient. Care should be taken not to overextend tissue or expander.

Before finally closing a wound it has to be ensured that the expander fits smoothly in with the surrounding tissue. After surgery, successive filling may incur formation of new folds

in the expander shell which, however, must be considered a possible side-effect of the expansion phase. Exceeding the recommended 6-month implantation period, folds may weaken the shell and lead to deflation of expander.

Such an effect is more likely to occur if the time intervals between individual fillings are extended. Therefore, expanders should be filled as quickly as circumstances allow. Once the expansion phase has been completed, the expander should be removed.

Expanders with remote valve


Please follow instructions below when handling the valve and trimming the tube:

1. The remote valve is permanent-ly connected to a tube of 5cm length that must not be shortened.
2. If necessary, cut the tube attached to the expander down to the necessary length.
3. Fit the expander tube onto the connector at the end of the 5-cm valve tube. The two tube ends should be arranged so as to reach at least the shoulder of the raised middle section of the connector.
4. Check the correct position and the firm hold of the two tubes on the connector.
5. We recommend fixing the two ends of the tube by using non-absorbable suture material to intensify the hold. Care should be taken not to perforate the tube.
6. Tube and expander should be kept free of tension when positioning the remote valve, taking into account the patient's normal movements.

Expanders with integrated valve and permanent magnet

ATTENTION: expanders with integrated magnets are not eligible for patients having pace-makers or patients who need an MRI examination during the implantation period.

Before the expander pocket is closed, please ensure that access to the integrated valve is not blocked by folds of the expander shell pushed over the valve. The use of permanent magnets in the valve guarantees an easy localization of the valves after implantation. To this effect, please use the enclosed magnet detector. The car-danic mechanism of the indicator magnet allows to exactly locate the magnet in the valve and determine the direction of injection: move the magnet detector up and down, right and left over the skin in the area of the magnet valve. When the detector's indicator is in a vertical position above the valve, press the magnet detector tip lightly against the skin to provoke an inflammatory reaction serving as injection mark.

Post-operative care : The patient's condition should be controlled after surgery during an appropriate period of time in order to follow up as quickly as possible on possible host response to implants or filling material (refer to "complications").

WARRANTY

All products are produced and tested for quality and consistency in conformity with the Medical Devices Act. Individual components have been selected thoroughly and comply with the latest standards for medical excellence, as does manufacturing.

Any information laid down in the "Instructions for use" must be followed. Replacement of products by POLYTECH Health & Aesthetics GmbH is only possible if manufacturing faults are definitely proven. There is no warranty whatever, expressed or implied, including any warranty under commercial law and any warranty for modifications the user may apply to products.

No representative or distributor of the manufacturer is entitled to change the above-listed warnings or precautions. Neither is any natural person nor legal party allowed to assume any liability or responsibility except the manufacturer.

Return of products: on all matters relating to the circulation of this medical product in Russia please apply to our Authorized representative in LLC "Beauty Matrix"; Russian Federation 119002, Moscow Kaloshin Pereulok, 4, build. 1.; Tel. 7499 918-77-70; info@bmatrix.ru

Disposal


Products that have expired should be destroyed as an epidemiologically non hazardous waste, close in composition to the municipal solid waste; Used, damaged parts thereof shall be destroyed as an epidemiologically hazardous waste.

Disposal of contaminated or used expanders and expanders of end shelf-life must be made in accordance with standard practice, adopted in the medical institution, and procedure established by local, federal and state regulations and the rules of law.

Our business activities are based on the DIN EN ISO quality standards guide. Our production effort, manufacturing and quality management is governed by the Medical Devices Act and EN ISO 9001 / EN ISO 13485, ranging from the choice of raw material, product design, production, sterilization, storage, final distribution to customer service.

Herewith I declare, that the present Instruction for Use is applicable to expanders. It is the latest version of this document at time of this declaration.

Dieburg, 05.12.2017

Dr. Oliver Bögershausen
Senior Director Operations

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