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

INDICATIONS

ULTRAPRO ADVANCED Mesh is indicated for the repair of abdominal wall hernias and abdominal wall deficiencies that require the addition of a reinforcing material to obtain the desired surgical result.

CONTRAINDICATIONS

ULTRAPRO ADVANCED Mesh should not be placed intraperitoneally.

DESCRIPTION OF THE MEDICAL DEVICE

ULTRAPRO ADVANCED Mesh is provided sterile.

ULTRAPRO ADVANCED Mesh is comprised of approximately equal parts of absorbable poliglecaprone-25 monofilament fiber and non-absorbable polypropylene monofilament fiber.

The polymer of the undyed and dyed (phthalocyanine blue, Color Index No. 74160) polypropylene fiber is identical to the material used for dyed/undyed PROLENE Suture.

The poliglecaprone 25 copolymer is identical to the material used in MONOCRYL Suture.

After absorption of the poliglecaprone-25 component, only the polypropylene mesh remains.

The remaining polypropylene mesh stretches to a greater degree perpendicular to the blue stripes.

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ACTIONS/PERFORMANCE

An Animal study showed that implantation of ULTRAPRO ADVANCED Mesh elicits an inflammatory reaction that stimulates the deposition of a thin fibrous layer of tissue, which can grow through the interstices of the mesh, thus incorporating the mesh into adjacent tissue. The absorbable component is essentially absorbed within 91 days, whereas the polypropylene material is not absorbed.

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TECHNICAL CHARACTERISTICS OF
MEDICAL DEVICE

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THE MATERIALS OF WHICH MADE THE MEDICAL DEVICE

The ULTRAPRO ADVANCED Mesh is made out of Prolene (un-dyed and dyed) monofilament and MONOCRYL(un-dyed) monofilament. The Prolene monofilament is an uncoated, synthetic, nonabsorbable thread composed of an isotactic crystalline stereoisomer of polypropylene, a synthetic linear polyolefin. The ULTRAPRO ADVANCED Mesh has approximately 17 % (by weight) pigmented polypropylene filaments and approximately 83 % natural polypropylene filaments by weight. The clear portions of the meshes are extruded from natural polypropylene resin. The blue strands are made from pigmented resin with copper phthalocyanine, blue 15, color index No. 74160. The Monocryl monofilament is an undyed, uncoated, synthetic, absorbable thread composed of poliglecaprone-25.

THE METHODS AND CONDITIONS OF USE

This medical device is intended for professional use in the medical and health care institutions. ULTRAPRO ADVANCED Mesh is intended for use only by physicians who are trained in the surgical procedures and techniques required for the hernia repairs and the implantation of synthetic meshes. The selection of mesh for any given patient is a function of numerous factors including, but not limited to, the patient's past medical and surgical history, current medical condition (i.e., comorbidities), surgical technique, and size and location of the hernia. The physician is advised to consult the medical literature regarding techniques,

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complications, and adverse reactions before selecting a mesh.

APPLICATION/INSTRUCTIONS FOR USE

Sizing and Placement

ULTRAPRO ADVANCED Mesh must always be placed extraperitoneally. ULTRAPRO ADVANCED Mesh can be trimmed at the surgeon's discretion while providing the necessary overlap to reduce the likelihood of recurrence.

The use of thermal cutting devices is not recommended because it has not been tested.

Ventral/Incisional and Inguinal Repair

In ventral or incisional hernia repair, the ULTRAPRO ADVANCED Mesh should extend at least 3 to 5 cm beyond the margins of the hernia defect, unless the surgeon decides that additional overlap onto healthy tissue is needed. When used in inguinal hernia repair, the mesh should provide sufficient overlap of the fascial defect on all sides.

The blue stripes on the mesh can be used for orientation and alignment purposes.

Place the mesh so that it lies flat to the tissue.

Fixation

The method of securing the implant (e.g., non-absorbable sutures, absorbable/no absorbable tackers) to provide for adequate mesh fixation and to minimize the risk of recurrence should be determined at the surgeon's discretion, based on the individual patient's needs and the nature of the repair.

Spacing and distribution between fixation points and technique should be determined at

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the surgeon's discretion to provide adequate mesh fixation, to reduce the risk of mesh migration, and to optimize mesh-to-tissue contact to foster tissue ingrowth.

Preclinical data and reported experience suggest having the fixation points at least 1 cm from the edge of the mesh and approximately 1 and 2 cm apart from each other.

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WARNINGS

- Failure to properly follow instructions may result in improper functioning of the device and could lead to injury. Please read all information carefully.
- If this device is used in patients with the potential for growth or tissue expansion (such as infants, children, or women who may become pregnant), the surgeon should be aware that the device will not stretch significantly as the patient grows.
- This device is indicated for abdominal wall hernia repair and not for gynecologic procedures. Gynecologic procedures should be performed using devices indicated for gynecologic repairs.
- It is recommended that the device not be used in a contaminated field, as contamination of the device may lead to infection that may require removal of the device.
- As with any implant, an acute and permanent foreign body response will occur. In some patients, this response can result in one or more of the adverse reactions listed below.
- The device is a permanent implant that is designed to integrate into the tissue.

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In cases in which the device needs to be removed, in part or in whole, significant dissection may be required.

- Insufficient overlap on any side of the defect may increase the risk of postoperative complications, including recurrence. Consult Application/Instructions for Use section.
- Insufficient or improper fixation may increase the risk of postoperative complications, including recurrence. Consult Application/Instructions for Use section.
- Do not resterilize/reuse. Reuse of this device (or portions of this device) may create a risk of product degradation, which may result in device failure and/or cross-contamination, which may lead to infection or transmission of blood-borne pathogens to patients and anyone coming in contact with the device.
- Inspect the mesh carefully before implantation. Do not use the device if it is damaged.

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PRECAUTIONS

The safety and effectiveness of pretreating of ULTRAPRO ADVANCED Mesh in combination with solutions (e.g., saline,

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medications) prior to implantation have not been studied.

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ADVERSE REACTIONS

Potential adverse reactions are those typically associated with surgically implantable materials, including:

- infection,
- inflammation,
- seroma formation,
- acute or chronic pain,
- foreign body sensation,
- hematoma,
- nerve damage, soft tissue injury, adhesion formation,
- fistula formation,
- extrusion/erosion, excessive contraction or shrinkage of the tissue surrounding the mesh,
- mesh failure/hernia recurrence.

One or more revision surgeries may be necessary to treat the above-mentioned adverse reactions. Revision surgery may not resolve the adverse reactions and may pose a risk of additional adverse reactions.

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STERILIZATION

The device is sterilized by Ethylene Oxide

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gas.
Do not resterilize.
Do not use if package is opened or damaged.

PACKAGING

Products are packed to accommodate flat screens without folds
ULTRAPRO ADVANCED Mesh device is secured in a paper folder, and a printed adhesive label is placed on the folder. The folder is then sealed within a foil pouch, comprised of a foil top layer and a foil bottom layer. A printed adhesive label is placed on the foil pouch. The sealed foil pouch and IFU are then placed within the secondary packaging (carton). A barcode label is applied to the back of the sales carton.

PACKAGING MATERIAL

Primary package – foil: Novelis; Amcor

Secondary package – carton: MPS (Multi Packaging Solutions)

Tertiary package: cardboard shipping container

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

Labeling of the Tertiary package:

- Name of Medical Device;
 - Transportation conditions;
 - Manufacturers Bar-codes;
 - Batch number;
 - Used before;
 - Quantity;
 - Cautions;
 - Address of manufacturer

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**USAGE, STORAGE AND
 TRANSPORTATION CONDITIONS**

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Usage:

Temperature +32°C - +42°C;
 Internal environment of human body

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PRODUCT EXPIRATION

ULTRAPRO ADVANCED Mesh: 2 years
Do not use beyond the expiration date

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WARRANTIES

The medical device is intact and unopened package stores all parameters of quality and sterility during shelf life.

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DISPOSAL

Dispose the product and its residual elements or waste items in accordance with local laws and regulations.

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AUTHORIZED REPRESENTATIVE IN THE RUSSIAN FEDERATION

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