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Abbott

Abbott Ireland Diagnostics
Division,
Lisnamuck,
Longford
Co. Longford
Ireland

Tel. +353 43 3331000
Fax +353 43 3331001

I, Lourda Fitzpatrick, Project Manager Regulatory Affairs, Abbott Ireland Diagnostics Division, Longford, Ireland, hereby certify that the attached documents that keep the number of pages indicated below are the Insert for Use. These documents are in a full accordance with actual manuals of manufacturer.

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Signature

Date: 19 OCT 2016

Lourda Fitzpatrick
Project Manager Regulatory Affairs
Abbott Ireland Diagnostics Division
Lisnamuck, Longford
Co. Longford
Ireland

ABBOTT DIAGNOSTICS DIVISION LONGFORD

INTENDED USE

The ARCHITECT 25-OH Vitamin D Controls are for the estimation of test precision and the detection of systematic analytical deviations of the ARCHITECT iSystem when used for the quantitative determination of 25-hydroxyvitamin D (25-OH vitamin D) in human serum and plasma.

Refer to the reagent package insert and the ARCHITECT System Operations Manual for additional information.

CONTENTS

- 1 Bottle (8.0 mL) Low Control contains 25-OH vitamin D prepared in PBS buffer with human serum. Preservatives: ProClin 950, Sodium Azide.
- 1 Bottle (8.0 mL) Medium Control contains 25-OH vitamin D prepared in PBS buffer with human serum. Preservatives: ProClin 950, Sodium Azide.
- 1 Bottle (8.0 mL) High Control contains 25-OH vitamin D prepared in PBS buffer with human serum. Preservatives: ProClin 950, Sodium Azide.

The controls are at the following ranges and approximate concentrations:

Control	Target Concentration		Control Range	
	(ng/mL)	(nmol/L)	(ng/mL)	(nmol/L)
CONTROL L	20.0	50.0	14.0 – 26.0	35.0 – 65.0
CONTROL M	40.0	100.0	28.0 – 52.0	70.0 – 130.0
CONTROL H	75.0	187.5	52.5 – 97.5	131.3 – 243.8

PRECAUTIONS

- **IVD**
- For *In Vitro* Diagnostic Use

-  **CAUTION:** This product contains human-sourced and/or potentially infectious components. Refer to the **CONTENTS** section of this package insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human-sourced materials should be considered potentially infectious. It is recommended that these reagents and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.¹⁻⁴
- The human serum used in the controls is nonreactive for HBsAg, HIV-1 RNA or HIV-1 Ag, anti-HCV and anti-HIV-1/HIV-2.

The following warnings and precautions apply to: **CONTROL L** / **CONTROL M** / **CONTROL H**

WARNING	Contain methylisothiazolone and sodium azide.
H317	May cause an allergic skin reaction.
EUH032	Contact with acids liberates very toxic gas.
Prevention	
P261	Avoid breathing mist / vapors / spray.
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves / protective clothing / eye protection.
Response	
P302+P352	IF ON SKIN: Wash with plenty of water.
P333+P313	If skin irritation or rash occurs: Get medical advice / attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
Disposal	
P501	Dispose of contents / container in accordance with local regulations.

Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 8.

STORAGE

- Controls are stable until the expiration date when stored and handled as directed.
- Do not use past expiration date.
- Controls must be stored at 2-8°C in an upright position and may be used immediately after removal from 2-8°C storage.

PREPARATION FOR ANALYSIS

Controls may be used immediately after removal from 2-8°C storage.

Prior to each use, mix by gentle inversion.

After each use, tightly close the caps and return the controls to 2-8°C storage.

BIBLIOGRAPHY

- 1. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- 2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- 3. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- 4. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.

Key to Symbols

	Caution
	Consult instructions for use
	Manufacturer
	Temperature limitation
	Use by/Expiration date
	Concentration
	Contains Sodium Azide. Contact with acids liberates very toxic gas.
	Control Low, Medium, High (L,M,H)
	In Vitro Diagnostic Medical Device
	Lot Number
	Product of Ireland
	Range
	List Number
	Warning: May cause an allergic reaction.

ICHITECT is a trademark of Abbott Laboratories in various jurisdictions. All other trademarks are property of their respective owners.

 Abbott Ireland
Diagnostics Division
Lisnamuck, Longford
Co. Longford
Ireland
+353-43-3331000

Customer Service: Contact your local representative
to find country-specific contact information on
www.abbottdiagnostics.com

Revised July 2015.
© 2015 Abbott Laboratories

ADDENDUM

to the Instruction for Use (for Russia)

SHIPPING

Recommended shipping conditions for device correspond to storage conditions provided in the Instructions for Use (2-8°C). Device is stable until the expiration date indicated on the package. Device is no longer suitable for use if shipping temperature limits have been violated (GOST R 51088-2013). Ship in an enclosed carrier vehicle of any kind in accordance with shipping rules, applicable to the respective carrier.

WARRANTY

Manufacturer guarantees that device remains stable until the expiration date when stored and handled as directed. Manufacturer set the warranty shelf life of device which corresponds to period between the date of manufacture and expiration date that are indicated on the label and/or in accompanying documentation. Device must be discarded if stored with deviations in manufacturer's requirements (GOST R 51088-2013).

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1. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
3. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
4. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.

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