TO THE
**FEDERAL SERVICE ON
SURVEILLANCE IN HEALTHCARE
(ROSZDRAVNADZOR)**
Department of Registration of
Medical Devices

Covering letter

With this letter we

B. Braun Surgical, S.A.
Carretera de Terrassa
121, 08191 Rubí (Barcelona)
Spain

herewith state that in the application of registration of

Histoacryl Flexible Pack

a instructions for use (IFU) for Russia which is drafted by LLC B. Braun Medical Russia according to Russian's related regulatory will be submitted. The updated contents are listed in below:

1. add the information of shelf life;
2. add the information about the transportation conditions of the product;
3. add the information about utilization of product;
4. add the information about warranty and claims.

Besides provided IFU was amended according to the recommendation of Roszdravnadzor and includes the mentioned data:

1. Use conditions of medical device
2. Solubility of Histoacryl Flexible
3. Viscosity of Histoacryl Flexible
4. Age restrictions of medical device;
5. Correct name of medical device;
6. Clarification of used symbols of medical device.

Moreover IFU has been recently updated and is already approved by B. Braun Surgical, S.A. Main changes done are:

1. Histoacryl Tip optional
2. Long incisions treatment by pre-suturing
3. Re-organisation of the information in the IFU

For and on behalf of

B. Braun Surgical, S.A. Rubí, 2 of March, 2016

Silvia Orus
Global Regulatory Affairs Manager
Closure Technologies
B. Braun Surgical, S.A.

B. BRAUN
B. Braun Surgical SA
Carretera de Terrassa, 121
08191 Rubí (Barcelona)
Tel.: 93 586 62 00*

Dr. Pau Tullon
R&D, Regulatory Affairs & Quality Director
CoE, Closure Technologies
B. Braun Surgical, S.A.

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Instructions for Use

Synthetic tissue adhesive with an application tip Histoacryl® *Flexible Pack* in a set:

1. 0.5 ml polypropylene ampoule of adhesive Histoacryl® *Flexible* – 5 pcs.; application tip Histoacryl® *Tip* – 5 pcs., instruction for use – 1 pcs.
2. 0.5 ml polypropylene ampoule of adhesive Histoacryl® *Flexible* – 10 pcs.; application tip Histoacryl® *Tip* – 10 pcs., instruction for use – 1 pcs.

Device description

Synthetic tissue adhesive with an application tip Histoacryl® *Flexible Pack* in a set is topical skin adhesive. Histoacryl® *Flexible* is a sterile, liquid consisting of n-butyl-2-cyanoacrylate and a softener. For easy visual assessment of the thickness of the layer which has been applied, Histoacryl® *Flexible* is coloured blue with D&C violet No. 2 dye.

Histoacryl® *Flexible* is supplied in 0.5 mL single use plastic ampoules. One ampoule should be used only in one patient. Each ampoule is sealed within an aluminium pouch in order to protect the ampoule from external contamination keeping at the same time the external sterility of the ampoule. Histoacryl® *Flexible* remains liquid until exposed to water or water containing substances, including tissues. Histoacryl® *Flexible*, when in contact with skin tissue, polymerizes forming a film that adheres to the underlying surface.

Studies⁽¹⁾ have shown that Histoacryl® *Flexible* acts as a barrier to microbial penetration as long as the adhesive film remains intact.

⁽¹⁾In vitro test assessing Histoacryl® *Flexible* microbial barrier during 7 days for the following bacteria: *Staphylococcus aureus*, *Staphylococcus epidermidis*, *Escherichia coli*, *Pseudomonas aeruginosa*, *Enterococcus faecium*, *Brevundimonas diminuta* and *Candida albicans*.

Indications of use


Synthetic tissue adhesive with an application tip Histoacryl® *Flexible Pack* in a set is indicated for closure of minimum-tension skin wounds up to 25 cm from clean surgical incisions and simple, thoroughly cleansed, trauma-induced lacerations.

The maximum length of the incision, without exceeding 25 cm, which can be treated with Histoacryl® *Flexible Pack* depends on the wound site, patient condition and surgeon's criterion. Histoacryl® *Flexible Pack* must be used in combination with subcuticular sutures to release the tension along the wound line.


Mode of application and dosage

Histoacryl® *Flexible Pack* is to be used only for the closure of smooth and fresh skin wounds.


- a) Prior to the application of Histoacryl® *Flexible Pack*, provide local anaesthesia if necessary, ensure that the edges of the wound have been thoroughly cleaned and debrided and always assure proper hemostasis. Ensure that any remaining topical medications/anaesthetics is removed prior to application.
- b) Histoacryl® *Flexible Pack* must be used in conjunction with and not in substitution of subcuticular sutures. Release tension along the wound line and ensure wound edges approximation with subcuticular sutures before applying Histoacryl® *Flexible Pack* to the skin surface. Appose skin edges avoiding irregularities on the wound surface.
- c) Pull open the aluminium pouch and remove the ampoule containing the adhesive in sterile condition. The ampoule containing the adhesive should only be removed from the aluminium pouch immediately prior to use.
- d) Hold the ampoule with the cannula pointing upwards, and flick it sharply to ensure no adhesive remains trapped within the cannula. Open the ampoule by twisting off the ribbed tip of the cannula. Hold the yellow part of the ampoule between two fingers, keeping it vertical with the thin end pointing upwards while opening the ampoule. This will prevent any loss of Histoacryl® *Flexible* by spillage from the ampoule when opening. Histoacryl® *Flexible* should be applied immediately after opening the ampoule.
- e) Histoacryl® *Flexible* may be used with an optional applicator tip Histoacryl® *Tip*. The applicator tip Histoacryl® *Tip* is a small cylindrical device intended to provide easy distribution of the adhesive over the wound or incision area. The applicator tip Histoacryl® *Tip* should be attached to the open cannula of the Histoacryl® *Flexible* ampoule as shown in the following pictures.


1. Open blister and take out the application tip

2. Open the ampoule by twisting off the ribbed tip

3. Attach application tip to the open ampoule

4. Apply the glue in a thin layer to the approximated wound edges and hold in apposition for approximately 30 seconds

f) Apply Histoacryl® *Flexible Pack* to the skin surface by carefully pressing the side walls of the plastic ampoule.

g) Hold the edges together for approximately 30 seconds after application to allow Histoacryl® *Flexible* to cure and to prevent displacement of the wound edges. Using forceps to hold skin edges in apposition while applying Histoacryl® *Flexible Pack* remains at surgeon's criterion.

Avoid heavy application. Unless otherwise prescribed, as less Histoacryl® *Flexible* skin adhesive as possible should be applied; the amount applied is sufficient if slight colouring is visible. A sparing application of thin layers or spots is required for undisturbed wound healing.

h) Discard any remaining Histoacryl® *Flexible Pack* without contaminating the surroundings. Prior to discarding the ampoule, replace the ribbed tip that was removed on opening the ampoule, putting it back on the end of the cannula.

Mode of Action

In the presence of tissue moisture, Histoacryl® *Flexible* immediately polymerises into a solid substance that attaches firmly to the tissue.

Contraindications

Only use Histoacryl® *Flexible* topically. Histoacryl® *Flexible* is not indicated for internal use. Thus, is not to be applied to any internal organs, blood vessels, nerve tissue, mucosal surfaces or mucocutaneous junctions, areas with dense natural hair or within the conjunctival sac of the eye. **Do not use for sclerosation therapy. Do not use for mesh fixation.**

Histoacryl® *Flexible* is not intended to be applied to the surface of the eye. If Histoacryl® *Flexible* accidentally drops in the eye, adhesion can be avoided by rinsing with water immediately. If adhesion or substance depositions occur, especially at the eye, they will dissolve without further action in a few days. Do not force the eye open.

Histoacryl® *Flexible* is not to be applied to wounds that show evidence of infection or gangrene, or to wounds of decubitus aetiology.

Histoacryl® *Flexible* is not to be used on patients with known preoperative systemic infections, uncontrolled diabetes, or diseases or conditions that are known to interfere with the wound healing process.

Histoacryl® *Flexible* has not been evaluated in patients with a history of hypertrophic scarring or keloid formation.

Histoacryl® *Flexible* is not indicated in patients with known hypersensitivity to cyanoacrylate, formaldehyde, triacetin or D&C violet No. 2 dye.

Warnings

Histoacryl® *Flexible* must not be introduced into the wound, since this would interfere with wound healing.

Histoacryl® *Flexible* generates a small amount of heat during polymerisation and should not be applied to tissues that may be affected by such heat.

A heavy application may cause thermal damage to tissues, resulting in delayed healing. Connective tissue healing can be impeded when too much tissue adhesive is applied.

Histoacryl® *Flexible* should not be applied to wet or bleeding wounds. Excess moisture (such as water or alcohol) or blood presence, may accelerate polymerisation, resulting in the generation of excess heat.

Histoacryl® *Flexible* is not to be applied below the surface of the skin because the polymerised adhesive is not absorbed by tissue and may elicit a foreign body reaction.

Histoacryl® *Flexible* is not intended to be applied to wounds subject to high skin tension or on areas of increased skin tension such as the elbows, knees, or knuckles. Histoacryl® *Flexible* is not to be used in areas of skin excision. Do not reuse. Reuse might cause an infection hazard for patients and/or users, impaired product functionality, and risk of injury, illness or death due to contamination.

Precautions

Only product of liquid consistency can be used. The condition of the adhesive can be judged with the plastic ampoule unopened.

The areas to be joined should be as dry as possible. Moisture may accelerate Histoacryl® *Flexible* polymerisation and might affect wound closure results. Any remaining topical medications/anesthetics must be removed prior to application.

If too much adhesive is accidentally applied, it can be removed during the first few seconds using a dry swab. Care must be taken, the dry swab could adhere.

The skin surfaces should be held together for approximately 30 seconds after application of Histoacryl® *Flexible*. Corrections cannot be made once the adhesive has hardened.

Wounds should be kept dry following closure with Histoacryl® *Flexible*. Do not apply liquid or topical medications onto wounds closed with Histoacryl® *Flexible*. The polymerised film might be weakened, leading to skin edge separation.

Histoacryl® *Flexible* should not be used on wound sites that will be subjected to repeated or prolonged moisture, movement or friction.

Accidental bonding of unwanted skin may occur. Do not pull apart skin. If necessary, adhesions in less sensitive skin areas can be carefully dissolved with acetone or by soaking in warm water until the skin can be separated.

Care must be taken to avoid direct contact of instruments, cloths or gloves with the adhesive applied, as they might adhere to the tissue. Instruments stained with Histoacryl® *Flexible* can be cleaned with dimethyl formamide or acetone.

The user must be familiar with the particular properties of liquid tissue adhesive before applying Histoacryl® *Flexible* skin adhesive.

Histoacryl® *Flexible* is intended for single use. Opened unused ampoules should be discarded.

Side effects

The use of this product leads to an exothermic polymerisation reaction and the released heat may result in a sensation of warmth.

Cyanoacrylates can be associated with a limited period of local sensitization or irritation reaction in the application area; a transient foreign body reaction can occasionally take the form of an inflammatory reaction.

Histoacryl® *Flexible* migration below the epidermal surface may result in "tattooing" of the underlying tissue.

Sterility

The sterility of Histoacryl® *Flexible* tissue adhesive is ensured by membrane filtration and aseptic filling. Histoacryl® *Tip* is ETO-sterilized. Do not re-sterilise.

Storage

Histoacryl® *Flexible Pack* should be stored in its original sealed aluminum pouch in temperatures at or below +25 °C away from moisture and direct light sources.

Product range

0.5 ml polypropylene ampoule of adhesive Histoacryl® *Flexible* – 5 pcs.; REF 1051250P
application tip Histoacryl® *Tip* – 5 pcs., instruction for use – 1 pcs.

0.5 ml polypropylene ampoule of adhesive Histoacryl® *Flexible* – 10 pcs.; REF 1051260P
application tip Histoacryl® *Tip* – 10 pcs., instruction for use – 1 pcs.

Symbols used on label

DIM Volume of Histoacryl® *Flexible*

REF Cat. No.

LOT Batch Number

 Use until

STERILE A Method of sterilisation of Histoacryl® *Flexible*: aseptic filling

STERILE EO Method of sterilization of Histoacryl® *Tip*: ETO-sterilization

 +25°C Upper limit of temperature (+25°C)

 Do not reuse

 See Instructions for Use

CE₀₁₂₃ CE-mark and identification number of notified body. Product conforms to the essential requirements of the Medical Device Directive 93/42/EEC of Histoacryl® *Flexible* and Histoacryl® *Tip*

PC Conformity marking of the Russian Federation

08/201

Synthetic tissue adhesive with an application tip Histoacryl® *Flexible Pack* in a set should be used by the professional's appointment in healthcare facilities.

Synthetic tissue adhesive Histoacryl® *Flexible* practically insoluble in water, ethanol and easily soluble in dimethyl sulfoxide (DMSO), and acetone.

Viscosity of synthetic tissue adhesive Histoacryl® *Flexible* is 3.8 - 4.6 mPaxs.

There is no age restrictions of the use of Synthetic tissue adhesive with an application tip Histoacryl® *Flexible Pack* in a set.

The shelf life of synthetic tissue adhesive with an application tip Histoacryl® *Flexible Pack* in a set in undamaged packing at proper storage conditions is 2 years from the production date.

The medical device must not be used after the expiry date.

Transportation conditions for synthetic tissue adhesive with an application tip Histoacryl® *Flexible Pack* in a set are equal to that of storage.

Expired or used synthetic tissue adhesive with an application tip Histoacryl® *Flexible Pack* in a set must be utilized as G-class waste in accordance with SanPin 2.1.7.2790-10.

B. Braun Surgical, S.A. guarantees that synthetic tissue adhesive with an application tip Histoacryl® *Flexible Pack* in a set complies with the requirements of all regulatory documents provided the storage and transportation conditions are observed. B. Braun Surgical, S.A. shall not be liable for the inappropriate (non intended) use of the medical device.

Date of information 02/2017

Version 03

R&D, Regulatory Affairs and Quality Director CoE,
Closure Technologies B. Braun Surgical, S.A.

Global Regulatory Affairs Manager CoE,
Closure Technologies B. Braun Surgical, S.A.

 Pau Turon
6/3/17
 Silvia Orús
3/3/2017

DG8079633

08/2016

.../...TICATION: I, CARLOS MATEO MARTINEZ DE BARTOLOME, Civil Law Notary in the City of Barcelona, Spain, certify that the signatures written on the first and third sheet are the true signatures of **Mrs. SILVIA ORUS PUIGVERT**, holder of a valid Spanish D.N.I number 33.906.328-G and **Mr. PAU TURON DOLS**, holder of a valid Spanish D.N.I number 46.041.264-W for its coincidence with the ones they commonly use and known by me. This document consists of FOUR sheets of plain paper sealed and initialed by me, with the security seals numbers this and previous numbering in sequence. I also note that my statement does not apply to the content of the document, which does not have public effect. In Barcelona, the 06 March 2017.-----

It stays written down and filed in my register with the number: **282**.

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This Apostille only certifies the authenticity of the signature and the capacity of the person who has signed the public document, and, where appropriate, the identity of the seal or stamp which the public document bears.

This Apostille does not certify the content of the document for which it was issued.

[This Apostille is not valid for use anywhere within Spain]

[To verify the issuance of this Apostille, see <https://eregister.justicia.es>]

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