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
Medtronic

MEDTRONIC Regulatory Affairs

Medtronic AF Solutions
8200 Coral Sea Street
Mounds View, MN 55112

www.medtronic.com

12/20/2013

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www.medtronic.com

UL 60601-1:2003.

Conformité Européenne (European Conformity). This symbol means that the device fully complies with European Directive 93/42/EEC (NB 0050).

Do not dispose of this product in the unsorted municipal waste stream. Dispose of this product according to local regulations. See <http://recycling.medtronic.com> for instructions on proper disposal of this product.

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Introduction to the generator 1

Description and intended use

The Model 990018 Medtronic GENius Multi-Channel RF Ablation Generator delivers temperature-controlled radiofrequency (RF) energy to user-selectable electrodes on compatible Medtronic multi-electrode cardiac ablation catheters.

The generator automatically recognizes the attached cardiac ablation catheter and programs preset default temperature, ablation duration, and energy mode setting parameters. Ablation parameters, including ablation duration, energy mode, target temperature, and channels (electrode pairs), can also be manually selected. A display screen monitor is available for viewing real-time system and ablation parameters.

The generator features the following RF energy mode selections: bipolar only, unipolar only, and combination ratios of bipolar:unipolar energy mode selections of 1:1, 2:1, and 4:1.

An ECG interface box connects the catheter with the generator, patient electrodes, the external monitoring equipment, and the pacing equipment that can be used during the procedure.

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Medtronic AF Solutions
Sea Street
MN 55112

medtronic.com

ablation. For PVAC, PVAC GOLD, MAAC, and MASC, deselect the affected electrode pair or discontinue ablation, if the electrode does not reach 50°C and is receiving maximum power delivery.

Use of RF energy near implanted devices – Implantable devices, such as pacemakers and implantable cardioverter-defibrillators (ICDs), may be inhibited or otherwise affected by RF energy. Since RF ablation therapy may affect implantable devices:

- Have temporary pacing and defibrillation equipment available.
- Program pacemaker sensing parameters to VOO or DOO to ensure that RF energy is not sensed as an intrinsic event.
- Program ICD detection parameters to OFF or MONITOR to ensure that the RF energy is not sensed and inappropriately treated as a ventricular fibrillation episode.
- Monitor surface and/or intracardiac electrograms and/or vital signs during RF energy delivery to assess for device interaction. Take appropriate action if any interaction is detected.
- Reprogram and confirm all parameters of implantable devices after an RF ablation procedure.
- Exercise extreme caution when delivering ablation energy in close proximity to implanted leads.
- Refer to the appropriate implantable device technical manual for additional information.

Electrical isolation during ablation – Do not allow the patient to contact grounded equipment that might produce electrical current leakage during ablation or DCCV (Direct Current Cardioversion). Patient should not contact earthed metal parts or parts with appreciable capacitance to earth. Electrical current leakage may induce arrhythmias that may result in the patient's death.

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Medtronic AF Solutions
10000 Grand Sea Street
Eden Prairie, MN 55112

www.medtronic.com

Placement of monitoring electrodes – Place all physiological monitoring electrodes as far away as possible from the patient return electrodes and their leads to avoid RF interference affecting the ability to interpret patient electrocardiograms (EGG).

Generator settings or equipment failure – Do not increase target temperature or ablation duration settings before checking for obvious defects or misapplication, as patient injury may occur. Apparent low power output or failure of the equipment to function correctly at normal settings may indicate faulty application of the patient return electrodes or failure of an electrical lead.

Electromagnetic interference (EMI) radiated – The generator emits RF energy during ablation at a frequency level that may cause EMI with unshielded electronic equipment. To minimize EMI, the generator should be moved away from any other electronic device. If EMI is apparent during the application of RF energy, EMI may be reduced by repositioning the generator or other equipment.

Electromagnetic interference (EMI) susceptibility – The generator has been designed to minimize electromagnetic interference (EMI). If interference should occur, move the generator away from the device generating the interference or place the generator at a different angle.

Cable positioning – Position all cables associated with the generator so they do not come in contact with the patient or other leads. RF interference may affect the ability to interpret patient EGMs.

Leakage current from connected devices – Use only isolated amplifiers, pacing equipment, and ECG equipment (IEC 60601-1 Type CF equipment, or equivalent), as patient injury or death may occur. Do not allow leakage current from any connected devices to the patient to exceed

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10 microamperes (μA) under any circumstances.

Equipment modification – Do not modify this equipment. Modifications may reduce system effectiveness and impact safety.

Flammable materials – Do not allow flammable material in the area where RF ablation procedures are performed. The risk of igniting flammable gases, flammable agents used for cleaning or disinfecting, or other materials is inherent in the application of RF energy.

Long-term risk – The long-term risks of lesions created by RF ablation have not been established. In particular, the long-term effects of lesions in proximity to the specialized conduction system or coronary vasculature are unknown.

No direct skin contact during system testing – Do not test the operation of the ablation system through direct contact with the electrodes. Testing through direct skin contact may cause damage to healthy tissue.

EP recording systems – Use only high-frequency, current-limiting EP recording systems or patient injury may occur.

Equipment stability during use – Do not disturb or move the generator, ECG interface box, or cables while the catheter is positioned within the patient, as patient injury may occur.

System compatibility – Use only the following components with the Medtronic GENius Multi-Channel RF Ablation Generator: the supplied Medtronic ECG interface box and cable, GENius Jr. Remote Control and cable, compatible Medtronic multi-electrode cardiac ablation catheters, compatible Medtronic catheter interface cables, the supplied Medtronic ECG amplifier cable, and Valleylab Patient Return Electrodes (Model E7506, two required). The safety and use of the system with other catheters or components have not been tested.

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Medtronic AF Solutions
8700 Coral Sea Street
Minnetonka, MN 55112

www.medtronic.com

Proper use of patient return electrodes – Use only the designated Valleylab Patient Return Electrodes (Model E7506, two required) when operating the ablation system. Prepare the area where the patient return electrodes are attached in accordance with the patient return electrode manufacturer documentation. For best contact, avoid areas with adipose tissue, bony prominences, fluid invasion, scar tissue, and excess hair (shave the patient if necessary). Do not use a patient return electrode that is beyond its expiration date, in damaged packaging. Incorrect use of the patient return electrodes may reduce efficacy and patient safety during ablation.

Grounded wall outlet – To avoid the risk of electric shock, this equipment must only be connected to supply mains with protective earth (grounded receptacle).

Unintentional power output – Failure of the ablation equipment can result in an unintentional power output increase. An unintentional power output increase may impact patient safety.

Stimulation equipment – Certain stimulation equipment may alter the RF energy waveform, which may trigger arrhythmias on rare occasions. To eliminate the influence of stimulation equipment during ablation, disconnect the stimulation equipment while ablating.

Precautions

Related product literature – Do not attempt to operate this Medtronic cardiac ablation system prior to completely reading and understanding the

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Medtronic

MEDTRONIC Regulatory Affairs

Medtronic AF Solutions
40 Coral Sea Street
Mounds View, MN 55112

www.medtronic.com

GENius Multi-Channel RF Generator operator's manual and the technical manual of the compatible Medtronic cardiac ablation catheter.

Required use environment – Cardiac ablation procedures should be performed only in a fully equipped electrophysiology laboratory.

Qualified users – This equipment should be used only by or under the supervision of physicians trained in left atrial ablation procedures using this Medtronic multi-channel RF ablation generator and compatible catheters.

Regular inspection -- Regularly inspect reusable cables for visual evidence of damage.

Disposal of system components – Contact Medtronic for the return of non-disposable components. The following components are considered non-disposable: generator, ECG interface box, remote control, and cables.

Generator placement – Place the generator on a level surface. Do not place other equipment on top of the generator. Do not place the generator on top of other equipment.

Overheating the generator – Do not place the generator near heat-generating equipment. Do not block cooling inlets or outlets on the generator.

Fluid incursion – The generator may not function correctly if the electronic circuitry or the connectors are wet:

- Do not allow any fluid or moisture into the generator, the ECG interface box, or the cables.
- Do not hang fluids above the generator.
- Do not immerse the cables into fluids.

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www.medtronic.com

Indications

See the technical manual for the catheter being used.

1. PVAC, PVAC GOLD Pulmonary Vein Ablation Catheter

The Pulmonary Vein Ablation Catheter (GOLD) is designed to be used for the creation of endocardial lesions (focal and linear) during cardiac ablation procedures for the treatment of symptomatic Atrial Fibrillation (AF).

The PVAC/PVAC GOLD is also intended to be used for cardiac electrophysiological (EP) mapping of pulmonary vein potentials, delivery of diagnostic pacing stimuli and verifying electrical isolation of the pulmonary veins post-treatment.

2. MAAC Multi-Array Ablation Catheter

The Multi-Array Ablation Catheter is designed to be used for the creation of endocardial lesions (focal and linear) in arrhythmogenic tissue in the left atrial body during cardiac ablation procedures for the treatment of symptomatic Persistent and Long Standing

Persistent Atrial Fibrillation (AF). The MAAC is also intended to be used for cardiac electrophysiological (EP) mapping and delivery of diagnostic pacing stimuli.

3. MASC Multi-Array Septal Catheter

The Multi-Array Septal Catheter is designed to be used for the creation of endocardial lesions (focal and linear) in arrhythmogenic tissue along the left atrial septum during cardiac ablation procedures for the treatment of symptomatic Persistent and Long Standing Persistent Atrial Fibrillation (AF).

The MASC is also intended to be used for cardiac electrophysiological (EP) mapping and delivery of diagnostic pacing stimuli.

4. PV-TRACKER Guide wire

The PV-Tracker guide wire is intended to facilitate the placement of devices during diagnostic and interventional procedures

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Medtronic AF Solutions
8201 Coral Sea Street
New, MN 55112

www.medtronic.com

5. ARRIVE Braided transseptal sheath

The ARRIVE Braided Transseptal Sheath is indicated for the percutaneous introduction of various types of cardiovascular catheters to all heart chambers including the left atrium via transseptal puncture.

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- Active sepsis
- Left atrial thrombus or myxoma
- Atrial Septal Patch or atrial septal defect repair with percutaneous device
- Known sensitivity to heparin
- Blood clotting abnormalities
- Venous filtering device (Greenfield Filter)
- Known or suspected left atrial myxoma
- Previous intra-atrial septal patch or device placement
- Previous or systemic embolization from the left side of the heart

Potential adverse events

Potential adverse events associated with cardiac catheter ablation procedures include, but are not limited to, the following conditions:

- Allergic reaction to x-ray contrast media
- Atrial septal defect
- Catheter entrapment
- Anesthesia reactions
- Arrhythmias, proarrhythmia
- Arterial esophageal fistula
- Bleeding related to anticoagulation

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Medtronic AF Solutions
Central Sea Street
Minneapolis, MN 55112

medtronic.com

- Body temperature elevation
- Bradycardia
- Cardiac perforation of the heart or other organs during transseptal puncture or other procedures
- Cardiac tamponade
- Cardiac thromboembolism
- Cerebrovascular accident (CVA) or transient ischemic attack (TIA)
- Chest discomfort
- Chronic cough
- Component damage to ICD or implanted pacemaker
- Death
- Dislodgement of implantable cardioverter defibrillator (ICD) or permanent pacing leads
- Heart failure
- Hematoma
- Hemoptysis
- High creatinine phosphokinase or troponin level
- Hypotension/hypertension
- Infections
- Myocardial infarction or ischemia
- Nerve injury or nerve damage
- Obstruction, perforation, damage, or spasm of the vascular system including the coronary circulation system
- Pericarditis or endocarditis
- Unintended complete or incomplete atrioventricular node (AV-Node) or sinus node block or damage
- Pleural or pericardial effusion
- Pneumonia
- Pneumothorax
- Pulmonary embolism
- Pulmonary infiltrates
- Pulmonary vein narrowing or stenosis
- Pseudoaneurysm in groin
- Radiation injury or damage and late malignancy

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Generator controls and connections 2

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System overview

Figure 2-1 shows connections among the generator, its components, the patient return electrodes, and a compatible catheter.

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
Figure 2-1. The generator, its components, patient return electrodes, and catheter

- 1 GENius cardiac ablation generator
- 2 Monitor
- 3 VGA monitor cables (not included)

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available for operation. The panel also provides receptacles for the ECG interface box and the optional remote control.


Figure 2-2. Generator front panel and monitor handle

- 1 Standby button
- 2 Start button
- 3 STOP button
- 4 Channel buttons
- 5 All Channels button
- 6 Target temperature adjustment buttons
- 7 Ablation duration adjustment buttons
- 8 Remote control cable receptacle
- 9 ECG interface box cable receptacle
- 10 Energy mode buttons
- 11 Monitor handle

Standby button

After the catheter has been positioned for ablation and prior to initiating an ablation, the Standby button is used to display electrode temperature information for all selected channels. Standby is also used to reset the Cumulative RF usage display (see "Cumulative RF usage display" on page 29).

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Shaw, Minnesota, MN 55112

www.medtronic.com

Note: Press and hold the Standby button for 3 seconds to reset the Cumulative RF usage display.

Start button

The Start button begins energy delivery according to the front panel settings. During the ablation, a continuous audible tone is emitted, indicating that energy is being delivered. When complete, an end-of-ablation audible tone (beep) sounds.

STOP button

The STOP button can be pressed anytime during an ablation sequence to immediately stop energy delivery. When the STOP button is pressed, the end-of-ablation audible tone (beep) sounds and the generator enters Setup mode.

Channel buttons

While in Setup mode, with the catheter connected, channels (electrode pairs) can selectively be turned On or Off for ablation. Each channel button (channels 1–6) corresponds to a pair of electrodes and includes an indicator that is illuminated when selected. Channels 7 and 8 are non-functional.

All Channels button

The All Channels button is available to select or de-select all channels with a single press. Channels can be de-selected during ablation. The default setting for channel selection is all electrode channels off.

Target temperature adjustment buttons

The target temperature can be adjusted (set) during Setup mode or Ablation mode. Push the **Temp Up** button to increase the target temperature and **Temp Down** button to decrease target temperature. Adjustments can be made from 50 °C to 65 °C in increments of 1 °C.

Ablation duration adjustment buttons

The value (in seconds) for ablation duration can be selected during Setup mode. Adjustments can be made from 45 to 120 seconds in 5-second increments.

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Medtronic AF Solutions
Sea Street
MN 55112
www.medtronic.com

Remote control cable receptacle

This receptacle accepts the connector from the optional remote control cable.

ECG interface box cable receptacle

This receptacle accepts the connector from the ECG interface box cable.

Energy mode buttons

The generator features these RF energy mode selections: bipolar only, unipolar only, and combination energy mode selections of 1:1, 2:1, and 4:1. The first number of the combined energy mode selection refers to the bipolar component, and the second number refers to the unipolar component. For example, in the 1:1 energy mode, both bipolar and unipolar are equally distributed over time. The generator also employs a software power limit based on the energy mode selected.

Bipolar button

The Bipolar button is used to select bipolar energy delivery. When selected, its indicator is illuminated.

Unipolar button

The Unipolar button is used to select unipolar energy delivery. When selected, its indicator is illuminated.

1:1 button

The 1:1 button is used to select bipolar-to-unipolar energy at a 1:1 ratio. When selected, its indicator is illuminated.

2:1 button

The 2:1 button is used to select bipolar-to-unipolar energy at a 2:1 ratio. When selected, its indicator is illuminated.

4:1 button

The 4:1 button is used to select bipolar-to-unipolar energy at a 4:1 ratio. When selected, its indicator is illuminated.

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Medtronic AF Solutions
Sea Street
Minneapolis, MN 55112

www.medtronic.com

If used, the remote control provides the same ablation controls as the front panel (see “Front panel” on page 19). Figure 2-3 shows the layout of the remote control panel.


Figure 2-3. Remote control panel

- 1 Standby button
- 2 Start button
- 3 STOP button
- 4 Channel buttons
- 5 All Channels button
- 6 Target temperature adjustment buttons
- 7 Ablation duration adjustment buttons
- 8 Energy mode buttons

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Medtronic


MEDTRONIC Regulatory Affairs

Medtronic AF Solutions
Coral Sea Street
Mounds View, MN 55112

www.medtronic.com

Rear panel

The rear panel of the generator provides the following receptacles and connectors:

- 1 Power input receptacle
- 2 Equipotential (chassis ground) receptacle
- 3 RS-232
- 4 USB connector
- 5 VGA Out connectors

Power input receptacle


The power input receptacle accepts the connector from the hospital-grade power cord (supplied) and includes the main On/Off switch.

Equipotential (chassis ground) receptacle

This single contact receptacle provides a positive connection to the generator chassis. Connecting a ground wire cable to this receptacle is recommended.

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Figure 3-1. Language selection instruction display

After the user follows on-screen instructions to press STOP, the generator displays a key to available languages (Figure 3-2).

CH1 = Danish CH6 = Polish
CH2 = Czech CH6 = Portuguese
CH3 = Finnish
CH4 = Greek
Press Stop Button for Language Menu
Appuyez sur Stop pour ouvrir le menu Langue
Stopp-Schaltfläche für Sprachenmenü drücken
Pulse el botón Stop para acceder al menú de idioma
Premere il pulsante Stop (Arresto) per accedere al menu delle lingue
Druk op de Stopknop voor het menu Taal
Tryck på stoppknappen för språkmeny
Trykk på Stop-knappen for å finne språkmenyen
Tryk på knappen Stop for sprogmenu
Stisknète tlačítko Stop pro nabídku jazyků
Valitse kielivalikko painamalla Stop
Πιέστε το κουμπί Stop (Διακοπή) για το μενού γλώσσας
Aby wejść do Menu języka, naciśnij przycisk Stop
Prima o botão de Stop para menu de idioma

Figure 3-2. Language selection key display

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Medtronic AF Solutions
Central Sea Street
MN 55112


www.medtronic.com

Instructions direct the user to press the corresponding Channel button to select a language, or to press STOP for more languages.

These screens are displayed at startup only.

Patient return electrode placement image

Prior to ablation, the display indicates patient electrode placement (see Figure 3-3).


This screen is displayed:

- After the generator has successfully initialized and after the language selection screen has timed out.
- After the Cumulative RF usage display has been reset.
- After one hour with no generator use.

Ensure that the patient return electrodes are properly placed and

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Medtronic AF Solutions
8900 Coral Sea Street
Minnetonka, MN 55112

www.medtronic.com

connected to the ECG box, and then press STOP to clear message and continue to the setup screen.

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Setup Screen

In Setup mode, the Setup screen allows for the selection of ablation settings, as shown in Figure 3-4.

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
Medtronic AF Solutions
8200 Coral Sea Street
Minneapolis, MN 55112

www.medtronic.com


- 1 Catheter image (will vary depending on catheter type used)
- 2 Ablation Time; starts incrementing after the initiation of ablation and counts up to the selected value for ablation duration
- 3 Energy mode identifier
- 4 Target temperature in °C
- 5 Target temperature range
- 6 Cumulative RF usage display - data may vary
- 7 Selectable channels (pairs)

Figure 3-4. Setup Screen

This screen is displayed:

- Upon initialization, if no error conditions are detected and after the patient return electrode placement image has been dismissed
- Immediately after an ablation sequence
- When the STOP button is pressed
- When a catheter is detected

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
www.medtronic.com

Standby Screen

In Standby mode, the Standby Screen (Figure 3-5) displays the following information.

Figure 3-5. Standby Screen

- 1 Catheter image showing selected pairs registering in blue
- 2 Ablation Time
- 3 Energy mode identifier
- 4 Target temperature
- 5 Target temperature range
- 6 Temperature readings for each selectable electrode in a selected channel

This screen is displayed:


- After the Standby button is pressed on the front panel or remote control box, if used.

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
Ablation Screen


During ablation, the Ablation Screen (Figure 3-6) displays the following information.

- 1 Catheter image with color-coded temperature indicators for each electrode
- 2 Ablation Time (increments up to preset time during ablation)

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Medtronic

MEDTRONIC
Regulatory Affairs

Medtronic AF Solutions
1000 1st Sea Street
Minneapolis, MN 55112

www.medtronic.com

3 Area for system status, user prompts, and system messages

4 Energy delivery mode identifier

5 Target temperature in °C

6 Target temperature range

7 Temperature readings for each electrode in a selected channel

8 Real-time power output data (watts) for each selected electrode within a channel

9 Real-time ablation duration at target temperature

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Figure 3-6. Sample Ablation screen

Note: The temperature readings for each electrode provide a 1-second running average that is calculated based on 8 temperature readings per second.

This screen is displayed:

■ Throughout ablation. It is displayed when the Start button has been pressed on the front panel or the remote control, if no error conditions are detected.

Temperature color-code key

The channel bars and catheter image use the following colors to provide temperature information:

Color	Significance
White	Temperature ramp-up
Green	Temperatures between 50°C and 65°C with power $\geq 3 \text{ W}$
Yellow	■ Temperatures between 66°C and 70°C at any power

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

Medtronic AF Solutions
12910 Coral Sea Street
Minnetonka, MN 55112

www.medtronic.com

- 1 Area for system status, user prompts, and system messages
- 2 Message popup
- 3 Cumulative RF usage display - data may vary

Figure 3-7. Sample system message display

A screen similar to this may be displayed when an error condition is detected.

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AF Solutions
Sea Street
MN 55112

www.medtronic.com

Prior to use, inspect all system components for visible damage. (See "Compatible components" on page 14). Replace any component if damage is found.

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Ensuring the catheter interface cable is sterile

The catheter interface cable is supplied sterile and ready for use. If repeat sterilization is required, it is recommended that the catheter interface cable be sterilized using the steam method with specific parameters. For steam sterilization, rely on standard hospital practice. Alternatively, a sterile probe cover can be used.

Steam method

Resterilization parameters

Method of sterilization	Steam
Steam cycle temperature	127 °C (260 °F)
Steam cycle time at temperature	35 ± 5 minutes
Cool-down time	30 minutes minimum

If the recommended sterilization parameters are followed, the catheter interface cable can be resterilized and reused up to 10 times. Continue to inspect for damage (bent pins, cracks in jacket, etc.) before each use.

Note: Following sterilization, allow the catheter interface cable to reach ambient temperature.

Sterile probe cover

If the steam sterilization parameters cannot be followed, it is recommended that a sterile probe cover be used to fully enclose the catheter interface cable before use. If a sterile probe cover is used in place of steam sterilization, there is no restriction for reuse of the catheter interface cable.

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Medtronic AF Solutions
4700 Cornl Sea Street
Mounds View, MN 55112

medtronic.com

Note: Inspect for damage (bent pins, cracks in jacket, etc.) before each use.

Connecting non-sterile cables

1. Connect the power cord.
 - a. Plug the power cord into the power entry receptacle on the back of the generator.
 - b. Plug the power cord into an AC outlet.
 - c. Position the generator and power cord so that the cord can be easily accessed and disconnected if necessary.
2. Connect the ECG interface box cable to the generator front panel receptacle labeled *ECG Interface Box*.
To do this, align the red connector key to the 12:00 position, and then push the connector into the receptacle firmly until it stops. To disconnect, pull the connector body until it separates from the receptacle.
3. Connect the other end of the ECG interface box cable to the receptacle labeled *Generator* on the ECG interface box.
4. Connect the ECG amplifier cable to the receptacle labeled *ECG Amplifier* on the ECG interface box. Connect the other end of the ECG amplifier cable to an external recording system.

Warning: Use only isolated amplifiers, pacing equipment, and ECG equipment (IEC 60601-1 Type CF equipment, or equivalent), as patient injury or death may occur. Do not allow leakage current from any connected devices to the patient to exceed 10 microamperes (μA) under any circumstances.

Connecting the remote control (optional)

1. If the remote control is used, connect one end of the remote control cable to the generator front panel receptacle labeled Remote Control.
2. If the remote control is to be used outside the sterile field, connect the other end of the remote control cable to the rear panel receptacle on the

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www.medtronic.com

remote control labeled Generator.

If the remote control is to be used within the sterile field, use sterile technique to cover the remote control with a sterile cover and connect the remote control cable to the rear panel receptacle on the remote control labeled Generator.

Initializing the generator

1. Pull out the handle on the monitor until it stops, and then rotate the monitor up until it stops in the vertical position.
2. Turn the generator power switch On (located on the rear panel). The generator initiates a power-up/self-test routine.
3. During the self-test routine, which may take up to 2 minutes, the screen intermittently displays test information and language selection menus. Upon completion of a successful self-test, the generator enters Setup mode and all default settings are enabled.

Default settings (no catheter attached)

Mode	Setup
Channels	All off (de-selected)
Target temperature	60 °C
Ablation duration	60 seconds
Energy delivery	Bipolar:Unipolar, 4:1

In the event of a power-up/self-test failure, the generator displays a system message, emits an audible tone (beep), and is locked from further use. For more information, see "System messages" on page 37.

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Medtronic

MEDTRONIC Regulatory Affairs

Medtronic Solutions
10000 Central Expressway
Minneapolis, MN 55112

medtronic.com

according to the manufacturer documentation. Refer to Figure 4-1 for proper placement. For best contact, avoid areas with adipose tissue, bony prominences, fluid invasion, scar tissue, and excess hair (shave the patient if necessary).

Note: Be sure to use only Valleylab Patient Return Electrodes (Model E7506).


Figure 4-1. Patient return electrode placement

For proper system operation, position the patient return electrodes according to the patient return electrode placement image (see Figure 4-1).
2. Connect the patient return electrodes to the receptacles labeled Patient Return Electrode on the ECG interface box.

Performing an ablation procedure 5

Attaching the catheter

Note: For catheter placement and use instructions, consult the technical manual of the compatible Medtronic cardiac ablation catheter being used.

Use sterile technique to attach the catheter to the sterile catheter interface

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Medtronic

MEDTRONIC Regulatory Affairs

Medtronic AP Solutions
1 Sea Street
Minneapolis, MN 55112

www.medtronic.com

cable.

When a compatible Medtronic multi-electrode cardiac ablation catheter is connected, the generator identifies the catheter (an hourglass is displayed during this process), then displays the appropriate catheter image and default system parameters.

Important: After identification, the catheter can be used for up to 6 hours.

Following the attachment of a catheter, the generator begins to monitor system conditions. For information about conditions it may detect, see "System messages" on page 37.

Setting generator parameters

During Setup mode, prior to delivery of RF energy, ensure that the following parameters are set according to physician preference:

- Channels (system default is Off)
- Duration (system default is 60 seconds)
- Target temperature (system default is 60 °C)
- Energy mode (system default exists for each catheter) - also see Table 5-1

For information about selecting ablation options, see Chapter 2, "Generator controls and connections" on page 15.

The generator software limits power to the electrodes depending on the selected energy mode.

Selected energy mode	Power limit
Unipolar	10 W
Bipolar	6 W
1:1	10 W
2:1	9 W

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Medtronic

MEDTRONIC Regulatory Affairs

AT Solutions
10000 Sea Street
Minnetonka, MN 55112
medtronic.com

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Performing ablation

Note: Consult the technical manual of the compatible Medtronic multi-electrode cardiac ablation catheter being used.

1. Verify that setup steps have been completed. (See Chapter 4, "Preparing for ablation" on page 30.)
2. Verify that the catheter has been inserted into the patient and attached to the system (see "Attaching the catheter" on page 34) and that ablation settings have been selected (see "Setting generator parameters" on page 34).
3. If necessary, press the STOP button on the generator to dismiss the Verify Return Electrodes Are Attached message and patient return electrode placement image.
4. Press the Standby button to automatically calibrate all thermocouples to the patient's body temperature. No energy is delivered to the catheter while in Standby mode.

Note: It is possible to begin ablation from Setup mode. To do so, press the Start button instead of Standby. Calibration will still occur.

5. Press the Start button to deliver RF energy and start the ablation sequence according to the front panel settings. The ablation feedback tone is enabled and sounds throughout the ablation. The Ablation Time display begins to increment.
6. Monitor the ablation readings provided on the Ablation Screen during the sequence. (See "Ablation Screen" on page 27 for more information about the display.)

Note: Target temperature can be adjusted and channel(s) de-selected during the ablation. The STOP button can be used at any time to end the ablation.

Following completion of the ablation sequence, the end-of-ablation beep sounds and the generator enters Setup mode. (See "Setup Screen" on page

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Medtronic

MEDTRONIC Regulatory Affairs

Medtronic RF Solutions
10000 Central Expressway
Minneapolis, MN 55112

www.medtronic.com

25 for more information about the display.)
7. Change the settings prior to the next ablation, if required. If no settings are changed, the generator maintains the last settings used during ablation.

8. Repeat steps 4 and 5 to perform additional ablation sequences, if required. (See "Cumulative RF usage display" on page 29 for information about total ablation time in the current session.)

Notes:

■ To prevent inaccurate temperature acquisition, the generator does not permit ablation until all selected electrode temperatures are below 45 °C.

■ See "System messages" on page 37 for a list of possible exceptions to normal procedure, and how the generator responds.

Shutting down and disconnecting the system

1. Following the final delivery of RF energy, turn off the generator at the main's On/Off switch located on the rear panel of the generator. On the power switch represents Off.
2. Disconnect all cables.
3. Ensure that all connections are clean and dry. If the catheter interface cable is heavily soiled, it is recommended that the cable be discarded.
4. Clean and store the generator and components as indicated in "Storage and handling" on page 43 and in Appendix A: "GENius generator components" on page 47

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The generator may display the following system messages. Follow the steps in the Troubleshooting column in the sequence specified until the system message clears. (Some messages may require another ablation to determine if the error condition has been corrected.). If the steps are unsuccessful, contact your Medtronic representative.

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Medtronic

MEDTRONIC Regulatory Affairs

Medtronic AT Solutions
1000 Central Expressway, Suite 200
St. Paul, MN 55112

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Medtronic

MEDTRONIC Regulatory Affairs

AF Solutions
1000 1st Sea Street
Bloomington, MN 55112

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Medtronic AF Solutions
1000 Coral Sea Street
Bridges View, MN 55112

www.medtronic.com

- Storage temperature limitations: 0 °C to 60 °C (32 °F to 140 °F)
- Storage relative humidity limitations: 15% to 95%
- Handle with care.
- Keep the generator dry.

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Cleaning

Wipe off the generator with a damp cloth. If necessary, use a mild detergent solution. Dry thoroughly.
Do not immerse in water.

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Technical specifications of the generator

Information about catheters is not included here. See the technical manual of the compatible Medtronic multi-electrode cardiac ablation catheter being used.

See Appendix A: GENius generator components on page 47 for specifications for generator components.

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GENius Multi-Channel RF Ablation Generator

Physical characteristics	
Weight	38 lbs (17.6 kg)
External dimensions	Height: 8.6 inches (21.8 cm) Width: 17.5 inches (44.4 cm) Depth: 17.1 inches (43.4 cm)
Features	Integral carrying handle Flat front panel Integral rubber feet Pull-out, 15-in display screen

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Medtronic AF Solutions
Sea Street
MN 55112

medtronic.com

	LCD resolution: 1024 x 768
Environmental characteristics	
Operating temperature	10 °C to 40 °C (50 °F to 104 °F)
Operating relative humidity	30% to 75% non-condensing
Storage temperature	0 °C to 60 °C (32 °F to 140 °F)
Storage humidity	15% to 95% non-condensing
Safety (Type CF) characteristics	
Safety class	Class I
System electrical leakage	≤10 µA
Grounding system	Floating (applied parts); integral Equipotential grounding connection
Defibrillator protection	Ability to withstand 5000 VDC defibrillator pulse
Operating parameters	
Rated output power	145 W
Rated load range	80 Ω to 250 Ω

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
Medtronic AF Solutions
Sea Street
Minneapolis, MN 55112

medtronic.com

The remote control is connected to either of two remote control cables that are available:

- Remote control cable Model 990041—15 ft (4.57 m)
- Remote control cable Model 990042—25 ft (7.62 m)

Instructions for Use

1. Ensure that the generator is turned Off (recommended).
2. If the remote control is used, connect one end of the remote control cable to the generator front panel receptacle labeled Remote Control.
3. If the remote control is to be used outside the sterile field, connect the other end of the remote control cable to the rear panel receptacle on the remote control labeled Generator.


If the remote control is to be used within the sterile field, use sterile technique to cover the remote control and connect the remote control cable to the rear panel receptacle on the remote control labeled Generator.

Cleaning and storage

Wipe off the remote control with a damp cloth. If necessary, use a mild detergent solution. Dry thoroughly.

Do not immerse in water.

Store in accordance with temperature described in technical specifications (page 49).

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Technical Specifications

Physical characteristics	
Weight	1.2 lbs (0.55 kg)
External dimensions	Height: 1.75 inches (4.4 cm) Width: 8.6 inches (21.8 cm) Depth: 4.3 inches (10.9 cm)
Features	Integral rubber feet Plastic case
Environmental characteristics	
Storage temperature	0 °C to 60 °C (32 °F to 140 °F)
Storage humidity	15% to 95% non-condensing
Operating temperature	10 °C to 40 °C (50 °F to 104 °F)
Operating relative humidity	30% to 75% non-condensing
Operating parameters	
Degree of protection against ingress of liquid	IPX0 (not protected)
Flammability	Not suitable for use in presence of flammable materials

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
www.medtronic.com

Intended Use

The remote control cable Model 990041 and Model 990042 are designed for use with the GENius Multi-Channel RF Ablation Generator Model 990018 and the GENius Jr. Remote Control Model 990029. It is **required only if the optional remote control is used**.

Description

The remote control cable connects the remote control and the generator.

Instructions for Use

1. Ensure that the generator is turned Off (recommended).
2. Plug the remote control cable into the remote control port.
3. Plug the cable into the Remote Control Cable port on the front of the generator.

Cleaning and Storage

Wipe off the remote control cable with a damp cloth.

If necessary, use a mild detergent solution. Dry thoroughly.

Do not immerse in water.


Store in accordance with conditions described in technical specifications (page 51).

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Technical Specifications

Physical characteristics	
Length (two lengths available)	15 ft (4.57 m), 25 ft (7.62 m)
Environmental characteristics	
Storage temperature	0 °C to 60 °C (32 °F to 140 °F)
Storage humidity	15% to 95% non-condensing

ECG interface box

Intended Use

The ECG interface box is designed for use with the GENius Multi-Channel RF Ablation Generator Model 990018. The ECG interface box is required for proper system operation.

Description

The ECG interface box connects the catheter with the generator, patient return electrodes, the external monitoring equipment, and the pacing equipment that can be used during the procedure.

The ECG interface box has the following receptacles:

■ Generator receptacle– This receptacle accepts the connector from the ECG interface box cable.

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
Medtronic AF Solutions
1000 Central Sea Street
St. Paul, MN 55112

medtronic.com

- **ECG Amplifier receptacle**— This receptacle accepts the connector from the ECG amplifier cable.
- **Catheter receptacle**— This receptacle accepts the connector from the catheter interface cable.
- **Patient Return Electrode receptacles**— These two receptacles accept the connector from the Valleylab Patient Return Electrodes (Model E7506) (delivered separately).


The following cables and electrodes are used with the ECG interface box:

- **ECG Interface Box Cable Model 990020**— This cable connects the ECG Interface Box to the generator.
- **ECG Amplifier Cable Model 990027**— This cable connects the ECG Interface Box to an external recording system.
- **Catheter Interface Cable Model 990004**— This cable connects the ECG Interface Box to the ablation catheter.
- **Patient Return Electrodes**—Refer to Valleylab Patient Return

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MEDTRONIC
Regulatory Affairs

Medtronic AF Solutions
Sea Street
Minneapolis, MN 55112
www.medtronic.com

Electrodes (Model E7506) manufacturer documentation for more information.

Instructions for Use

1. Ensure that the generator is turned Off (recommended).
2. Insert the connector of each interconnecting cable to the appropriate receptacle on the ECG interface box.

Cleaning and Storage

Wipe off the ECG interface box with a damp cloth. If necessary, use a mild detergent solution. Dry thoroughly.
Do not immerse in water.
Store in accordance with temperature described in technical specifications

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Technical Specifications

Physical characteristics	
Weight	1.2 lbs (0.55 kg)
External dimensions	Height: 2.5 inches (6.3 cm) Width: 4.75 inches (12 cm) Depth: 7.75 inches (19.6 cm)
Environmental characteristics	
Operating temperature	10 °C to 40 °C (50 °F to 104 °F)
Operating relative humidity	30% to 75% non-condensing

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
ECG Interface Box Model 990028.

Description

The ECG interface box cable connects the ECG interface box with the generator.

Instructions for Use

1. Ensure that the generator is turned Off (recommended).
2. Connect the ECG interface box cable to the receptacle on the ECG interface box labeled Generator.
3. Connect the other end to the receptacle labeled ECG Interface Box on the front panel of the generator.

Cleaning and Storage

Wipe off the ECG interface box cable with a damp cloth. If necessary, use a mild detergent solution. Dry thoroughly.

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
ECG amplifier cable

Intended Use

The ECG Amplifier Cable Model 990027 is designed for use with the GENius Multi-Channel RF Ablation Generator Model 990018 and the ECG Interface Box Model 990028.

Description

The ECG amplifier cable connects the ECG interface box to an external recording system.

Instructions for Use

1. Ensure that the generator is turned Off (recommended).

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
Medtronic AF Solutions
Coral Sea Street
MN 551

www.medtronic.com

2. Connect the ECG amplifier cable to the receptacle labeled ECG Amplifier on the ECG interface box.
3. Connect the other end of the ECG amplifier cable to an external recording system.

Cleaning and Storage

Wipe off the ECG amplifier cable with a damp cloth. If necessary, use a mild detergent solution. Dry thoroughly.
Do not immerse in water.
Store in accordance with temperature described in technical specifications

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Technical Specifications

Physical characteristics	
Length	6 ft (1.8 m)
Environmental characteristics	
Storage temperature	0 °C to 60 °C (32 °F to 140 °F)
Storage humidity	15% to 95% non-condensing

Catheter interface cable

Intended Use

The Catheter Interface Cable Model 990004 is designed for use with the GENius Multi-Channel RF Ablation Generator and with compatible

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
Medtronic AF Solutions
8000 Coral Sea Street
Minneapolis, MN 55112

www.medtronic.com

Medtronic multi-electrode cardiac ablation catheters.

Description

The catheter interface cable connects the ECG interface box and the catheter. It is shipped sterile and must be resterilized or covered with a sterile cover before reuse.

Instructions for Use

1. Connect the cable to the catheter.
2. Connect the catheter interface cable to the ECG interface box, at the receptacle labeled Catheter.

Cleaning and Storage

The catheter interface cable is shipped sterile. After use and prior to subsequent use, the catheter interface cable must be resterilized. (Alternatively, it may be covered in a sterile fashion during reuse.) Do not immerse in water.

Store in accordance with temperature described in technical specifications

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MEDTRONIC

Regulatory Affairs

Medtronic AF Solutions
R200 Sea Street
Mound View, MN 55112

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If the recommended sterilization parameters are followed, the catheter interface cable can be resterilized and reused up to 10 times. Continue to inspect for damage (for example, bent connection pins, damage to jacket, etc.) before each use.

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Note: Inspect for damage (for example, bent connection pins, damage to jacket, etc.) before each use.

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Power Cord

Intended Use

The Medtronic Model 990023, 990024, 990025, 990035, 990082, and 990085 power cords are designed for use with the GENius Multi-Channel RF Ablation Generator Model 990018.

Description

The power cord supplies electricity to the generator. It is required for operation of the generator.

The power cord connects to:

- The console – The cord plugs into the specified inlet on the back of the generator.
- A wall outlet – The cord connects to a source of line power.

Instructions for Use

1. Plug the power cord into the power entry receptacle on the back of the generator.
2. Plug the power cord into an AC outlet.
3. Position the generator and power cord so that the cord can be easily accessed and disconnected if necessary.

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www.medtronic.com

Cleaning and Storage

Wipe off the power cord with a damp cloth. If necessary, use a mild detergent solution or isopropyl alcohol. Dry thoroughly.

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Disposal

Do not dispose of this product in the unsorted municipal waste stream.

Dispose of this product according to local regulations. See

<http://recycling.medtronic.com> for instructions on proper disposal of this product.

Technical specifications

Physical characteristics	
Length	990023, 990024, 990035, 990082, 990085: 2.5 m (8.2 ft) 990025: 3.6 m (12 ft)
Plug Type	990023: CEE 7/7 990024: BS 1363/A 990025: NEMA-5-15 990035: SEV 1011 990082: AS/NZS 3112:2000 990085: Danish Medical DK-2-8A
Connector Type	IEC 60320 C13
Current Rating	10 A
Voltage Rating	250 VAC
Conductor Size	1.00 mm ²
Environmental characteristics	
Storage temperature	0 °C to 60 °C (32 °F to 140 °F)
Storage humidity	15% to 95% non-condensing

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www.medtronic.com

CH-1131 Tolochenaz
Switzerland
Tel. +41-21-802-7000
Fax +41-21-802-7900

Australia
Medtronic Australasia Pty Ltd
97 Waterloo Road
North Ryde, NSW 2113
Australia

World Headquarters
Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432-5604
USA
Tel. +1-763-514-4000

Manufacturer
Medtronic Ablation Frontiers LLC
2210 Faraday Ave., Suite 100
Carlsbad, CA 92008
USA
Tel. +1-760-827-0001
Fax +1-760-827-0020
Technical Support: +1-763-514-4000

Medtronic E.C. Authorized Representative
Medtronic B.V.
Earl Bakkenstraat 10
6422 PJ Heerlen
The Netherlands
Tel. +31-45-566-8000

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