ACKNOWLEDGED
as true, reliable, and correct

President of Spident Co., Ltd.

Mr. Jemo Ahn

SPIDENT CO., LTD

Spident Co., Ltd.

**203&312, Korea Industrial Complex, 722,
Gojan-Dong, Namdong-Gu,
Incheon, Korea 405-821**

www.spident.co.kr

Tel.: +82(32)821-0071

Fax: +82(32)821-0074

info@spident.co.kr

EXPLOITATIVE DOCUMENT
of the MEDICAL DEVICE

Light curing liquid dam FineDam
Accessories:

- 1. Prefilled syringe, 1,2 ml – 4 ea.;**
- 2. Tips – not more than 10 ea.;**
- 3. PE bag for tips with zipper**

2015

Contents

1. Name of the Medical Device.....	3
2. Technical Specifications	3
3. Recommended Use of the Medical Device	3
4. Area of Use	3
5. Indications for use	3
6. Contraindication for Use and Potential Adverse Effects (Precautions)	3
7. Composition of the Medical Device	4
8. Description of the Operating Principles	5
9. Directions of Use.....	5
11. Storage conditions	6
12. Transportation conditions.....	6
13. Packaging and package content.....	6
14. Labeling information	7
15. Safe Disposal and Utilization Instructions	7
16. Manufacturer's Warranty.....	7
17. Shelf life of the Medical Device.....	8
18. Manufacturer	8
19. Manufacturer's Representative and Customers Support Service in Russian Federation	8

1. Name of the Medical Device

Light curing liquid dam FineDam with accessories:

- Prefilled syringe, 1,2 ml – 4 ea.;
- Tips – not more than 10 ea.;
- PE bag for tips with zipper.

2. Technical Specifications

Light curing liquid dam FineDam has the following features and benefits:

- Easy to use;
- High elasticity;
- Has low adhesiveness and easy to remove as one piece from slots and holes;
- Due to the blue colour creates visible contract with gingiva, tooth and whitening material.

Main technical characteristics of FineDam:

Characteristic	Units	Specification
Appearance		Paste
Colour		Blue
Odour		Odourless
Weight (in a syringe)	g	2,00±10%
Curing depth	mm	4.24 (standard: ≥1,5)
Water sorbtion	µg/mm ³	12,74 (standard: ≤ 40)
Solubility	µg/mm ³	1,99 (standard: ≤ 7,5)
Curing time	sec	Not more than 30

3. Recommended Use of the Medical Device

FineDam is a light-cured resin barrier for isolating tissue adjacent to teeth being whitened. For single tooth bleaching, FineDam may be used to protect adjacent teeth.

FineDam has the following functions: 1) working field isolation from the rest oral tissues; 2) protection of soft tissues from dangerous liquids; 3) decrease the risk of infecting.

For single use only.

4. Area of Use

Medical and health care institutions. Dentistry.

5. Indications for use

Working field isolation from the rest oral tissues during treatment (for example – whitening).

6. Contraindication for Use and Potential Adverse Effects (Precautions)

Contraindication:

Allergic to any ingredients in FineDam.

Special precautions:

- In case of a sensitizing reaction to any of the FineDam component, stop using FineDam and consult a doctor. In these cases FineDam must no longer be used.
- FineDam is only to be used for oral tissues. Only in dentistry.
- Avoid contact with eyes. In case of accidental contact with eyes wash eyes with running water for 15 minutes.
- Avoid contact with skin. In case of accidental skin contact wash skin with running soap water then consult a medical specialist if necessary.
- Not for internal use in normal conditions.
- Please read the instruction before use.
- Close the cap immediately after use.
- Individual measures: safety glasses, gloves.
- Store only in original packaging.
- Do not use expired FineDam.

7. Composition of the Medical Device

Light curing liquid dam FineDam has the following composition:

Chemical Name	CAS No	Per 100 g
Urethane dimethacrylate mixture	72869-86-4	74.98
2-Hydroxyethyl methacrylate	868-77-9	19.99
Silicon Dioxide, Amorphous	7631-86-9	4.00
Camphorquinone	10373-78-1	0.20
Ethyl 4-dimethylaminobenzoate	10287-53-3	0.40
2,6-Di-tert-butyl-p-cresol	128-37-0	0.01
Benzotriazol	2440-22-4	0.01
Titanium dioxide	13463-67-7	0.40
Blue pigment	147-14-8	0.01

The list of accessories:

1. Syringe, 1,2 ml – 4 ea.;
2. Tips – not more than 10 ea.;
3. PE bag for tips with zipper

FineDam subsidiary materials content:

No	Accessories	Material, Type	Manufacturer
1	Syringe	Polypropylene Yuplene R370Y	Shina Corporation http://www.shinacor.co.kr/
2	Tip	Polypropylene Moplen HP640R Stainless steel AISI 304 (SUS 304) (needle)	Biodent Co., Ltd http://www.biodent.co.kr/
3	PE bag for tips with zipper	PE	

8. Description of the Operating Principles

Dental dam (Kofferdam – in German) – is a working field isolation system, which is usually a latex dam, that isolates the teeth to be treated from other oral tissues. Dental dam benefits:

- Protection against ingestion;
- Eliminating the possibilities to damage mucosal membranes;
- Protection of the staff from infection;
- Oral soft tissues retraction;
- Working field is kept dry and disinfected for a long time;
- Good visibility, the mouth is always opened;
- Discussion with the patient is stopped;
- No need in rinsing;
- No cotton rolls.

Using of latex dental dams has some shortcomings, the main ones are patient's sensitivity to latex and some difficulties in dam positioning, requiring certain skills for the dentist. Liquid dam is a light cured paste or jelly material containing methacrylate, prefilled in syringes. Time of curing (time of light polymerization) is 20 seconds. The liquid jelly or paste fills the interdental space, adheres to the enamel and gingival margins, creating good protection for the soft tissues. The liquid dam has high elasticity and easily removes as one piece.

9. Directions of Use

1. Remove the cap from the syringe.
2. Place the tip onto the syringe.
3. Turn tip 90 degree to lock in place.
4. Apply FineDam (see fig. 1) – 4-6 sm width and 1,5-2 sm depth and 5 mm put at the enamel.
5. Make sure that the material applied interproximally and embraces the last tooth to be whitened.
6. Light curing for 20 sec.


Fig. 1. Directions of use of FineDam.

10. Conditions of Use (Temperature, Relative Humidity, Atmospheric Pressure)

The medical device must be used in accordance with the good industrial hygiene and safety practice.

Conditions of use:

Temperature: $(+2^{\circ} - +27^{\circ}) \pm 1^{\circ}\text{C}$;

Relative humidity: (not more than 70%) $\pm 2\%$;

Atmospheric pressure: 760 ± 10 mm Hg.

11. Storage conditions

Temperature: $(+2 - +27^{\circ}\text{C}) \pm 1^{\circ}\text{C}$;

Relative humidity: (not more than 70%) $\pm 2\%$;

Atmospheric pressure: 101 ± 1 kPa (760 ± 10 mm Hg).

The device must only be kept in original packaging.

The devices must be stored in a dry, dark and cool place away from direct sunlight, atmospheric precipitations and soil moisture in stands or cases 30 cm above the floor, away from direct sunlight, heaters and other sources of heat (not less than 1 m).

The packed devices must be protected from dirt, dust, moist, and damage.

12. Transportation conditions

Temperature	Relative humidity	Atmospheric pressure	Period of time
$+2^{\circ}\text{C} - +27^{\circ}\text{C}$	not more than 70%	760 ± 10 mm Hg	Not more than 72 hours
Note: the device must be kept at the ambient temperature for at least 8 hours before use			

Safety of transportation:

The device is not recognized as dangerous and not regulated a dangerous good.

Avoid exposure to heat (above $+27 \pm 1^{\circ}\text{C}$) and light during transportation.

13. Packaging and package content

Package content:

Light curing liquid dam FineDam is supplied in prefilled syringes 1,2 ml. Each syringe is closed with a protective cap made from PP. 4 syringes along with tips (not more than 10 ea.), packed into a zipper PE bag, and IfU are packed into the customer package – a labeled PE-bag. The bag is sealed.

Labeling information:

The device package shows the following required information:

- manufacturing company name;
- brand name of the device;
- type of the devices;
- name and registered office address of the manufacturer;
- name and registered office address of the manufacturer's representatives;
- country of the manufacturer;
- manufacturing date;

- advertising slogan "You're the Artist".
- number and date of GOST-R Registration Certificate.
- information of customers' support service;
- barcode;
- batch code / LOT number;
- use before date;
- symbols, specifying the information necessary for the proper usage of the the medical device (specified in section 14. Labeling Information).

14. Labeling information

	Batch code		CE mark
	Manufacturer		Manufacturing date
	EC representative		Expiry date
	For single use		Do not use if the package id damaged
	Store between 2°C – 27°C		GMP mark
	Keep dry		Keep away from sunlight
	Non sterile		Rosstandard mark

15. Safe Disposal and Utilization Instructions

Medical device must be disposed after usage, after the expiry date is reached, or if the package was damaged and the product cannot be used.

Medical devices must be disposed in accordance with the medical practices, which are regular for the medical institutions, and with the local state and federal regulatory documents. The choice of methods of safe disposal depends on the medical institution capacity and profile, availability of the waste decontamination/deactivation installations, the local routine ways of the waste decontamination/deactivation (burning, landfilling, utilization).

In Russian Federation the medical device must be disposed in accordance with the actual Russian Federation Sanitary Regulations and Standards for medical wastes management.

16. Manufacturer's Warranty

Spident Co., Ltd. guarantees that this product has been manufactured and packaged with the due care prescribed in the EU Directive for Medical Devices 93/42/EEC. Each product is individually tested and inspected before it is released for packaging. Spident Co., Ltd. will replace any device if defects which can be proven to have originated during manufacture or packaging are immediately brought to the attention of Spident Co., Ltd. or its distributors. This warranty applies in lieu of and ruling out any other verbal or written, express, implicit or other warranty.

As a result of biological differences in individuals, no product is 100% effective under all circumstances. Due to this fact and since we have no control over the conditions, Spident Co.,

Ltd. and its distributors do not guarantee good effectiveness in application; unsatisfactory results also cannot be ruled out. Spident Co., Ltd. and its distributors shall not be liable for any, damage or losses of any nature arising directly or indirectly from the use of the device. The product liability expires if any damage caused by transport, storage or improper handling by the operator, or other factors outside the control of the manufacturer occurs. Spident Co., Ltd. and its distributors shall not be liable for damages arising from reuse of the product.

17. Shelf life of the Medical Device

Shelf life of the medical device – 2 years from the manufacturing date.

18. Manufacturer

Spident Co., Ltd.
203&312, Korea Industrial Complex, 722,
Gojan-Dong, Namdong-Gu,
Incheon, Korea 405-821
www.spident.co.kr
Tel.: +82(32)821-0071
Fax: +82(32)821-0074
info@spident.co.kr

19. Manufacturer's Representative and Customers Support Service in Russian Federation

LLC «Spident»
Per. Dmitrovskiy 13, room 7,
Saint-Petersburg, 191025, Russia
Tel. +7 (812) 670-80-80
Fax +7 (812) 670-91-39
info@spident.ru

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