Labo America Inc.

Address: 920 Auburn Court, Fremont, CA 94538, USA

Phone: +1 (510) 445-1257

Fax: +1 (510) 445-1317

E-mail: sales@laboamerica.com

Website: <http://www.laboamerica.com>

General Director

Gautam Aggarwal


February 10, 2017

A handwritten signature in black ink, appearing to read "Gautam Aggarwal", written over a horizontal line.

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Instruction manual

For registration of the product on the territory of the Russian Federation

To guarantee proper operation of the microscope, as well as to avoid injury in the process of work with the instrument, it is strongly recommended to read this manual in full and comprehend it before commencement of work.

Content

1. INTRODUCTION	3
2. SAFETY GUIDELINE	4
3. SPECIAL INSTRUCTIONS	5
4. DELIVERY SET	6
5. DESCRIPTION OF PARTS OF THE MICROSCOPE	11
6. ASSEMBLY PROCEDURE	15
7. ELECTRICAL INTERFACES	19
8. HANDLING	20
9. OPERATIONAL INSTRUCTIONS FOR MICROSCOPE	21
10. HOW TO FOCUS ON OBJECT	22
11. REPLACEMENT OF OBJECTIVE/EYEPIECE	22
12. REPLACEMENT OF LIGHT SOURCE	22
13. USE OF ACCESSORIES	23
14. FUSE REPLACEMENT	27
15. TENSION ADJUSTMENT IN THE PROCESS OF ACCESSORIES USE	27
16. DISINFECTION AND STERILIZATION	27
17. CLEANING AND MAINTENANCE	28
18. TROUBLESHOOTING	29
19. TECHNICAL SPECIFICATIONS	30
20. OPERATION, TRANSPORTATION AND STORAGE CONDITIONS	37
21. INSTALLATION AND USE OF Pixel Pro SOFTWARE	37
22. INSTALLATION AND MAINTENANCE OF MICROSCOPE	39
23. MARKING	40
24. INFORMATION ON ELECTROMAGNETIC COMPATIBILITY	40
25. MANUFACTURER WARRANTY	43
26. DISPOSAL CONSIDERATIONS	43
27. MANUFACTURER INFORMATION:	44

I. INTRODUCTION

This manual extends to operative microscopes Labomed PRIMA OPH, PRIMA ENT, PRIMA DNT, PRIMA GN (hereinafter referred to as the microscopes) designed for magnification of minute structures, during performance of delicate surgical procedures, which high-power magnification in transmitted light.

Area of use –referral medical and preventive treatment facilities.

Microscopes are manufactured in the following embodiments:

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Microscopes provide high quality of optical image, good depth of focus and wide field of view. Adjustment of light, built-in decline, adjustment of optical head helps to reduce fatigue and allow working in comfortable conditions for a long time.

Abilities of these microscopes:

1. Optical head can be easily positioned with the help of suspended arm.
2. Step five-fold magnifier gives optimal magnification at five different levels.
3. Lighting with cold light of LED lamp of 50 W is provided with the help of fiber optic cable. Lighting is adjusted to achieve the most suitable brightness with the help of adjusting handle, which is conveniently located on the suspended arm.
4. Electric power supply of the lamp runs down automatically, when suspended arm is raised in its highest position.
5. When microscope is not used, the arm can be bent towards the body of the microscope for the purpose of compact storage.
6. Mobile counter provides excellent steadiness along with the microscope mobility.

2. SAFETY GUIDELINE

1. This microscope is manufactured in compliance with CE safety regulations.
2. Microscope is designed for use only in accordance with this manual.
3. Service and repair are permissible only by authorized persons.
4. Replace blown-out fuses with similar fuses (the same voltage, current intensity and run-down characteristics).
5. Connect the equipment to power outlet and use plug with adjusted earthing.
6. Don't apply force in the process of connection of cables connectors. If the plug and outlet can't be connected, make sure they match each other. If any connector is damaged, please, contact the representative of the manufacturer.
7. Make sure that free passable ventilation air inlets/outlets of the body are not closed.
8. Microscopes are designed for use only in dry locations. Prevent microscope elements from water penetration. Don't put vessels with liquid on the microscope body.
9. Microscope has overheat protection by means of thermal switch.
10. Manufacturer is not responsible for damages caused by actions of unauthorized persons; it will also deprive you of any rights for warranty.
11. It is recommended to use microscope only with accessories that are included into delivery set. If you want to use other accessories, make sure they are certified by LABOMED with respect to safety operation together with the microscope.

Highly important: For identification Service Department workers must know series number of the microscope.

3. SPECIAL INSTRUCTIONS

3.1 Before every use and after reassembly of the microscope

1. Check out electric connectors (connections).
2. Install aseptic caps, panels or leads, which were removed or open.
3. Pay your attention to the stickers on the microscope, such as warning signs, warning triangles with exclamation mark or notes.
4. Don't close ventilation outlets.

3.2 When microscope is functioning

1. Don't look directly at the light source and objective lens of the microscope.
2. When light is on, both ends of fiber optic cable should be connected. Otherwise there is risk of fire or burns.
3. Any type of emission has a heavy toll on biological tissue. It is valid for illuminated field; in this regard adjust brightness and length of light of the field at the absolute bare minimum.
4. Adjust resistance of the suspended arm for your comfort.

3.3 After each use of the microscope

1. Shut down the microscope from stationary power source.
2. When microscope is not used, suspended arm can be bent towards the microscope for the purpose of compact storage.

3.4 Contact during operation

1. During operation of microscope to perform required functions microscope and its component parts are in direct contact with a patient.

4. DELIVERY SET

Microscope is delivered disassembled.

While unpacking, please, check delivery set according to Table 1.

Table 1

No.	Name	Quantity, pieces
I	Operative microscope Labomed PRIMA OPH, PRIMA ENT, PRIMA DNT, PRIMA GN*:	
1	Operative microscope PRIMA OPH, consisting of:	
1.1	Microscope stand	1
1.2	Rotary arm	1
1.3	Suspended arm	1
1.4	Inclining arm with the head	1
1.5	Wide-field eyepiece with adjustment 10×/Fn 18 mm.	2
1.6	Objective f=200 mm with fine focus knob and sterilized cap	1
1.7	Binocular tube 45° or 90° or 0°- 210°	1
1.8	Network cable	1
1.9	Pedal lever	1
1.10	Fuse 1A/250 V	2
1.11	H-shape mobile pillar with 4 wheels	1
1.12	Sterilized caps for knobs of arms blocking	6
1.13	Sterilized caps for magnifier knob	2
1.14	Sterilized cap for control knob	2
1.15	Sterilized cap for adjustment handle of lightning	1
1.16	Objective cover	1
1.17	Pointing device XYZ	1
1.18	Operational documentation:	
1.18.1	User manual	1
1.18.2	Warranty card	1
	Accessories:	
1	Adjusted wide-field eyepiece 10×/FN 16 mm with detachable cap	2
2	Adjusted wide-field eyepiece 12.5×/FN 18 mm with fixing cap	2
3	Bracket on the wall	1
4	Overhead bracket with column	1
5	Table-top bracket with short column	1
6	Key set for installation (5 in pack)	1
7	Supervision bridge for the assistant with inclined head, without eyepieces	No more than 2
8	Objective f=175 mm with focusing mechanism with sterilized cap	1
9	Objective f=200 mm with focusing mechanism with sterilized cap	1
10	Objective f=250 mm with focusing mechanism with sterilized cap	1
11	Objective f=300 mm with focusing mechanism with sterilized cap	1
12	Objective f=400 m with focusing mechanism with sterilized cap	1
13	Dust cover	1
14	Beam splitter	1
15	Double beam splitter 50:50	1
16	Double beam splitter 80:20	1
17	Inclined double beam splitter 50:50	1
18	Inclined double beam splitter 80:20	1

19	Adapter for digital camera DSLR	1
20	Adapter for CCD camera	1
21	Adapter for Sony Handy Cam	1
22	Adapter for Canon camera	1
23	Adapter for Nikon camera	1
24	Adapter for Sony camera	1
25	Camera with built-in beam splitter	1
26	Feeding cable for camera with built-in beam splitter	1
27	Camera	1
28	USB cable/microUSB.	1
29	HDMI cable/miniHDMI.	1
30	Extra LED 50 W	1
31	CD for Pixel Pro software	1
32	USB key for software Pixel Pro	1
33	Binocular tube 45°	1
34	Binocular tube 90°	1
35	Binocular tube 0°- 210°	1
2	Operative microscope PRIMA ENT, consisting of:	
2.1	Microscope counter	1
2.2	Rotary arm	1
2.3	Suspended arm	1
2.4	Inclining arm with the head	1
2.5	Wide-field eyepiece with adjustment 10×/Fn 18 mm.	2
2.6	Objective f=250 mm with fine focus knob and sterilized cap	1
2.7	Binocular tube 45° or 90° or 0°- 210°	1
2.8	Network cable	1
2.9	Pedal lever	1
2.10	Fuse 1A/250 V	2
2.11	H-shape mobile pillar with 4 wheels	1
2.12	Sterilized caps for knobs of arms blocking	6
2.13	Sterilized caps for magnifier knob	2
2.14	Sterilized cap for control knob	2
2.15	Sterilized cap for adjustment knob of lightning	1
2.16	Objective cover	1
2.17	Operational documentation:	
2.17.1	User manual	1
2.17.2	Warranty card	1
	Accessories:	
1	Adjusted wide-field eyepiece 10×/FN 16 mm with detachable cap	2
2	Adjusted wide-field eyepiece 12.5×/FN 18 mm with fixing cap	2
3	Bracket on the wall	1
4	Overhead bracket with column	1
5	Table-top bracket with short column	1
6	Key set for installation (5 in pack)	1
7	Supervision bridge for the assistant with inclined head, without eyepieces	No more than 2
8	Objective f=175 mm with focusing mechanism with sterilized cap	1

9	Objective f=200 mm with focusing mechanism with sterilized cap	1
10	Objective f=250 mm with focusing mechanism with sterilized cap	1
11	Objective f=300 mm with focusing mechanism with sterilized cap	1
12	Objective f=400 m with focusing mechanism with sterilized cap	1
13	Dust cover	1
14	Beam splitter	1
15	Double beam splitter 50:50	1
16	Double beam splitter 80:20	1
17	Inclined double beam splitter 50:50	1
18	Inclined double beam splitter 80:20	1
19	Adapter for digital camera DSLR	1
20	Adapter for CCD camera	1
21	Adapter for Sony Handy Cam	1
22	Adapter for Canon camera	1
23	Adapter for Nikon camera	1
24	Adapter for Sony camera	1
25	Camera with built-in beam splitter	1
26	Feeding cable for camera with built-in beam splitter	1
27	Camera	1
28	USB cable/microUSB.	1
29	HDMI cable/miniHDMI.	1
30	Extra LED 50 W	1
31	CD for Pixel Pro software	1
32	USB key for software Pixel Pro	1
33	Binocular tube 45°	1
34	Binocular tube 90°	1
35	Binocular tube 0°- 210°	1
3	Operative microscope PRIMA DNT, consisting of:	1
3.1	Microscope stand	1
3.2	Rotary arm	1
3.3	Suspended arm	1
3.4	Inclining arm with the head	1
3.5	Wide-field eyepiece with adjustment 10×/Fn 18 mm.	2
3.6	Objective f=250 mm with fine focus knob and sterilized cap	1
3.7	Binocular tube 45° or 90° or 0°- 210°	1
3.8	Network cable	1
3.9	Pedal lever	1
3.10	Fuse 1A/250 V	2
3.11	H-shape mobile pillar with 4 wheels	1
3.12	Sterilized caps for knobs of arms blocking	6
3.13	Sterilized caps for magnifier knob	2
3.14	Sterilized cap for control knob	2
3.15	Sterilized cap for adjustment knob of lightning	1
3.16	Protector of objective	1
3.17	Objective cover	1
3.18	Operational documentation:	
3.18.1	User manual	1

3.18.2	Warranty card	1
	Accessories:	
1	Adjusted wide-field eyepiece 10×/FN 16 mm with detachable cap	2
2	Adjusted wide-field eyepiece 12.5×/FN 18 mm with fixing cap	2
3	Bracket on the wall	1
4	Overhead bracket with column	1
5	Table-top bracket with short column	1
6	Key set for installation (5 in pack)	1
7	Supervision bridge for the assistant with inclined head, without eyepieces	No more than 2
8	Objective f=175 mm with focusing mechanism with sterilized cap	1
9	Objective f=200 mm with focusing mechanism with sterilized cap	1
10	Objective f=250 mm with focusing mechanism with sterilized cap	1
11	Objective f=300 mm with focusing mechanism with sterilized cap	1
12	Objective f=400 mm with focusing mechanism with sterilized cap	1
13	Dust cover	1
14	Beam splitter	1
15	Double beam splitter 50:50	1
16	Double beam splitter 80:20	1
17	Inclined double beam splitter 50:50	1
18	Inclined double beam splitter 80:20	1
19	Adapter for digital camera DSLR	1
20	Adapter for CCD camera	1
21	Adapter for Sony Handy Cam	1
22	Adapter for Canon camera	1
23	Adapter for Nikon camera	1
24	Adapter for Sony camera	1
25	Camera with built-in beam splitter	1
26	Feeding cable for camera with built-in beam splitter	1
27	Camera	1
28	USB cable/microUSB.	1
29	HDMI cable/miniHDMI.	1
30	Extra LED 50 W	1
31	CD for Pixel Pro software	1
32	USB key for software Pixel Pro	1
33	Binocular tube 45°	1
34	Binocular tube 90°	1
35	Binocular tube 0°- 210°	1
4	Operative microscope PRIMA GN, consisting of:	
4.1	Microscope stand	1
4.2	Rotary arm	1
4.3	Suspended arm	1
4.4	Inclining arm with the head	1
4.5	Wide-field eyepiece with adjustment 10×/Fn 18 mm.	2
4.6	Objective f=300 mm with fine focus knob and sterilized cap	1
4.7	Binocular tube 45° or 90° or 0°- 210°	1
4.8	Network cable	1

4.9	Fuse 1A/250 V	2
4.10	H-shape mobile pillar with 4 wheels	1
4.11	Sterilized caps for knobs of arms blocking	6
4.12	Sterilized caps for magnifier knob	2
4.13	Sterilized cap for control knob	2
4.14	Sterilized cap for adjustment knob of lightning	1
4.15	Objective cover	1
4.16	Operational documentation:	
4.16.1	User manual	1
4.16.2	Warranty card	1
	Accessories:	
1	Adjusted wide-field eyepiece 10×/FN 16 mm with detachable cap	2
2	Adjusted wide-field eyepiece 12.5×/FN 18 mm with fixing cap	2
3	Bracket on the wall	1
4	Overhead bracket with column	1
5	Table-top bracket with short column	1
6	Key set for installation (5 in pack)	1
7	Supervision bridge for the assistant with inclined head, without eyepieces	No more than 2
8	Objective f=175 mm with focusing mechanism with sterilized cap	1
9	Objective f=200 mm with focusing mechanism with sterilized cap	1
10	Objective f=250 mm with focusing mechanism with sterilized cap	1
11	Objective f=300 mm with focusing mechanism with sterilized cap	1
12	Objective f=400 m with focusing mechanism with sterilized cap	1
13	Dust cover	1
14	Beam splitter	1
15	Double beam splitter 50:50	1
16	Double beam splitter 80:20	1
17	Inclined double beam splitter 50:50	1
18	Inclined double beam splitter 80:20	1
19	Adapter for digital camera DSLR	1
20	Adapter for CCD camera	1
21	Adapter for Sony Handy Cam	1
22	Adapter for Canon camera	1
23	Adapter for Nikon camera	1
24	Adapter for Sony camera	1
25	Camera with built-in beam splitter	1
26	Feeding cable for camera with built-in beam splitter	1
27	Camera	1
28	USB cable/microUSB.	1
29	HDMI cable/miniHDMI.	1
30	Extra LED 50 W	1
31	CD for Pixel Pro software	1
32	USB key for software Pixel Pro	1
33	Binocular tube 45°	1
34	Binocular tube 90°	1
35	Binocular tube 0°- 210°	1

*- embodiment is chosen by the Customer from represented models: PRIMA OPH, PRIMA ENT, PRIMA DNT, PRIMA GN.

5. DESCRIPTION OF PARTS OF THE MICROSCOPE


Fig. 1 Microscope PRIMA DNT

1. Wheel with brake
2. Mobile H-shaped base
3. Microscope pillar
4. Rotary arm
5. Suspended arm
6. Knob of the rotary arm lock
7. Knob of the suspended arm lock.
8. Knob of the inclined arm lock
9. Inclined arm
10. Knob of light filter switching
11. Hydraulic lock of the suspended arm
12. Adjustment of the spring tension
13. Lead of the rotary arm
14. Binocular tube with eyepieces
15. Magnification switch
16. Handles
17. Objective
18. Knob of the head block
19. Knob of lightning adjustment
20. Optic fiber cable
21. Pedal lever


Fig. 2 Microscope PRIMA ENT

1. Wheel with brake
2. Mobile H-shaped base
3. Microscope pillar
4. Rotary arm
5. Suspended arm
6. Knob of the rotary arm lock
7. Knob of the suspended arm lock.
8. Knob of the inclined arm lock
9. Inclined arm
10. Knob of light filter switching
11. Hydraulic lock of the suspended arm
12. Adjustment of the spring tension
13. Lead of the rotary arm
14. Binocular tube with eyepieces
15. Magnification switch
16. Handles
17. Objective
18. Knob of the head block
19. Knob of lightning adjustment
20. Optic fiber cable
21. Pedal lever


Fig. 3 Microscope PRIMA GN

1. Wheel with brake
2. Mobile H-shaped base
3. Microscope pillar
4. Rotary arm
5. Suspended arm
6. Knob of the rotary arm lock
7. Knob of the suspended arm lock.
8. Knob of the inclined arm lock
9. Inclined arm
10. Knob of light filter switching
11. Hydraulic lock of the suspended arm
12. Adjustment of the spring tension
13. Lead of the rotary arm
14. Binocular tube with eyepieces
15. Magnification switch
16. Handles
17. Objective
18. Knob of the head block
19. Knob of lightning adjustment
20. Optic fiber cable


Fig. 4 Microscope PRIMA OPH

1. Wheel with brake
2. Mobile H-shaped base
3. Microscope pillar
4. Rotary arm
5. Suspended arm
6. Knob of the rotary arm lock
7. Knob of the suspended arm lock.
8. Knob of the inclined arm lock
9. Pointing device XYZ
10. Handle
11. Hydraulic lock of the suspended arm
12. Adjustment of the spring tension
13. Lead of the rotary arm
14. Binocular tube with eyepieces
15. Magnification switch
16. Handles
17. Objective
18. Knob of the head block
19. Knob of lightning adjustment
20. Optic fiber cable
21. Pedal lever

6. ASSEMBLY PROCEDURE

6.1 Mobile stand is installed by means of fixing of pillar to the base. Put the pillar on the lead screw of the base. Match screw hole with the screw. Hold the pillar and wind screws up at the bottom of the foundation with the help of wrenches included into delivery set.

See figure 5.

6.2 After installation of the column on the mobile base, make sure that all the elements are connected properly.

6.3 Take out rotary arm and suspended arm from the box.

Install rotary arm on lag hinge (1). See figure 6.

Fix rotary arm with the help of screw on top. See figure

● 1

Release retention screw of the rotary arm and screws of suspended arm to make rotation possible.


Fig. 5


Fig. 6


Fig. 6.1

6.4. Take out inclined arm with the head from the box and install it on the suspended arm, having inserted guiding axle (1) directly into suspended arm. Don't forget to release retention screw (2) before adjustment of guiding axle. Fix the locking (3). See figure 7, 7.1, 7.2.


Fig. 7 for PRIMA DNT and ENT


Fig. 7.1 for PRIMA GN


Fig. 7.2 for PRIMA OPH

6.5. Install binocular tube and eyepieces on the head. Fix binocular tube with the help of lock. See figure 12.


Fig. 8

6.6. Remove warning tape from suspended arm. Take protecting screw out of suspended arm with the help of wrench 5.0 mm. Replace protective screw with fixation knob, packing washer and cap (2). Tighten in full. See figure 9.


Fig. 9

6.7 Remove cable cover, having squeezing out two screws, located on the top. Pass optic fiber cable. See figure 10.


Fig. 10

6.8. Install cable cover on the suspended arm, having matched driving screw with the slot, situated on the lead, using wrench 3 mm.

6.9. Remove lock knobs A and B of the rotary arm. Remove self-tapping screws, fixed on the upper part and bottom of the lead of the rotary arm. Fix panels on the rotary arm and put back the screws that have been removed before. In the process of panel installation make sure that duct for optic fiber cable is made.

See fig. 11.


Fig. 11

6.10 Pass optic fiber cable into the hole, situated on the back wall of the head, pushing it through slightly as it's shown on figure 12. Audible click proves that the cable was installed correctly.


Fig. 12

6.11. Connect cable of the pedal lever to the microscope (for embodiments PRIMA DNT, PRIMA ENT, PRIMA OPH) See figure 13.


Fig. 13

Procedures, shown on figures 8 - 13 are identical for each embodiment.

7. ELECTRICAL INTERFACES

7.1 Connect feeding cable to the slot of alternating current (1) see fig. 14 and 14.1 situated on the rotary arm.

Turn on electric supply pressing power switch control (2).

Note: Electric line voltage on the place of microscope installation must correspond to the permissible voltage range. Otherwise, you mustn't work with the microscope.


Fig. 14 PRIMA OPH, DNT, ENT


Fig. 14.1 PRIMA GN

7.2 Additionally for PRIMA OPH:

- 7.2.1. Switch connector 7 pin, passing from XY junction to the terminal block of the electric power supply.
- 7.2.2. Switch connector 4 pin, passing from control junction Z of the axle to the terminal block of the electric power supply.
- 7.2.3. Switch connector 12 pin, passing from pedal control to the terminal block of the electric power supply.

Electrical interfaces are ready.

8. HANDLING

8.1 Power switch control ON/OFF (on Fig. 15) is situated under rotary arm. In 'ON' position, green LED switches on and cooling fan is launched. Before switching on the microscope, adjust light control device at the minimal level.

To increase life cycle of the lamp, switch off lantern, if the microscope is not used for a long time.

8.2 Light control device

Light control device is situated in front of the suspended arm. With the help of light control device vie field lightning can be adjusted. See figure 15.

8.3 Brakes

Prevent the pillar from uncontrolled movements by means of applying the brakes, located on the wheels. To release the brake, press the upper part of the brakes. See figure 16.

8.4 Knob of rotary arm lock

This knob helps to lock movements of the rotary arm, lock the position after initial focusing of the interest range by means of clockwise rotation of the knob. See fig. 15.

8.5 Knob of suspended arm lock

This knob locks suspended arm movement in any positional suitable for you by means of clockwise rotation of the knob. See fig. 15.

8.6 Knob of inclined arm lock

This knob is situated on suspended arm and helps to lock optic head in any position by means of clockwise rotation. See fig.17.

8.7 Hydraulic knob of suspended arm

Adjust tension of the suspended arm by means of adjusting screw in any position suitable for you, rotating knob clockwise or counter clockwise.


Fig. 15


Fig. 16


Fig. 17

9. OPERATIONAL INSTRUCTIONS FOR MICROSCOPE

Adjustment of magnification (see fig. 18)

1. Adjust the biggest magnification by means of one of the shift knobs of magnification (26a), located on the head.
2. Precise focusing is done with FOV knob (26b).
3. Absolute center adjusting of scope view in the view field can be done with the handle movement (26c).
4. Make sure that shift knob is in the index point in "stop" position.
5. Lightning is adjusted with the knob, figure 18.1. To adjust desirable level of lightning rotates the knob clockwise or counter clockwise.


Fig. 18


Fig. 18.1

10. HOW TO FOCUS ON OBJECT

1. Adjust eyepieces at '0' diopter.
2. Adjust IPD (interpupillary distance) of the optic head with the help of IPD scale.
3. Adjust highest factor of magnification in "stop" position with the help of one of the knobs situated on the head. Due to this scope view remains in focus at all levels of magnification.
4. Precise focusing is done with FOV knob (1, fig. 19) by means of its clockwise or counter clockwise rotation, or with the help of pedal putting the blades up and down.
(for embodiments OPH, ENT, DNT)


Fig. 19

11. REPLACEMENT OF OBJECTIVE/EYEPIECE

Fig 20.

1. Objectives can be disconnected, turning them counter clockwise.
Objectives can be installed back turning them clockwise.
Additionally for PRIMA DNT:
Use objective protector while working for its protection.
2. To install eyepieces, insert tube of the optic head.


Fig. 20

12. REPLACEMENT OF LIGHT SOURCE

Fig. 21.

Open leads of the rotary arm. Disconnect optic fiber cable and replace light source assembly with the new one. Put leads of the arm back.


Fig. 21

13. USE OF ACCESSORIES

13.1 Double beam splitter. (fig. 22)

Beam splitter splits main ray of light into several, one of them is directed to the binocular tube and two auxiliary ports, where adapters for photo/video devices and observation bridge for the assistant with inclined head. Correlation of light splitting can be 50:50 or 80:20. If the correlation is 80:20, 80% of light goes to the observed and 20% to the connected photo/video devices.


Fig. 22

13.2 Inclined beam splitter. (fig. 23)

Inclined beam splitter serves for the same purpose as beam splitter, but it provides inclined position of 60°.


Fig. 23

13.3 Adapter для Sony Handy cam. (fig. 24)

Adapter for Sony Handy cam permits to connect camcorder Sony Handy cam HDR-CX-130 to the beam splitter.


Fig. 24

13.4 Observation bridge for the assistant with inclined head. (fig.25)

Supervision bridge for the assistant with inclined head is used for observation of the object observed with a microscope by the assistant of the operator.


Fig. 25

13.5 Adapter for CCD camera. (fig. 26)

Adapter for CCD camera makes possible the connection of CCD camera to the beam splitter to work in real-time mode.


Fig. 26

13.6 Adapter for digital camera DSLR. (fig. 27)

Adapter for digital camera DSLR makes possible connection of DSLR camera to the microscope.


Fig. 27

Compatible types of cameras:

Canon DS126291, Nikon D60, Sony HDR-CX-130.

13.7 Adapters for cameras Nikon, Canon, Sony. (fig.28)

Allow connecting of DSLR camera to beam splitter with removed objective.


Fig. 28

13.8 Installation of beam splitter with observation bridge for the assistant with inclined head and DSLR camera (fig.29)

1. Release retention screw ① on the head.
2. Install beam splitter ③ on the head ② and tighten retention screw on the head.
3. To install observation bridge for the assistant ④ in beam splitter, insert in into the port and tighten locking ring.
4. To install adapter for digital camera DSLR ⑤ into beam splitter, insert it into the port of beam splitter and tighten locking ring ⑥.


Fig. 29

13.9 Installation of beam splitter with adapter for Sony Handy cam and adapter for CCD camera (fig.30):

1. Release retention screw ① on the head.
2. Install beam splitter ③ on the head ②.
3. Tighten retention screw ① on the head.
4. To install adapter for Sony Handycam Into beam splitter, insert adapter for Sony Handycam ④ into beam splitter and tighten locking ring.
5. To install adapter of CCD camera into beam splitter insert adapter ⑤ into beam splitter and tighten locking ring ⑥.


Fig. 30

13.10 Installation of Canon, Nikon and Sony DSLR cameras (fig.31):

1. Remove objective mechanism from the camera (as shown on the figure to the right).
2. Install adapter for camera (Canon, Nikon or Sony) on DSLR adapter ③, screwing the rung tightly.
- ④. Insert compatible camera into ② adapter for camera ①.
3. Insert adapter DSLR ③ into beam splitter and tighten locking ring.


Fig. 31

13.11 Installation of CCD camera and Sony Handy cam (fig.32):

1. Install Sony handy cam ① on the adapter ② and lock its rolling rings clockwise.
2. Install CCD camera ③ on adapter ④ and lock it clockwise.


Fig. 32

13.12 Operation manual of camera with built-in beam splitter.

Camera with built-in beam splitter and camera are digital cameras, developed for work with microscope, which make it possible to display image of the observed object on computer and take videos.

Camera control is performed with the help of software PixelPro (see clause 21).

CONNECTION TO MICROSCOPE.

1. Camera with built-in beam splitter is connected to microscope with the help of feeding cable. To connect camera to computer use attached cable USB/microUSB and cable HDMI/miniHDMI.

Installation of camera with built-in beam splitter is performed in the following way:

1.1 Release retention screw on the head.

1.2. Install camera in built-in beam splitter on the head and tighten retention screw on the head.

2. Feeding of camera and connection with computer is conducted with USB cable.

Installation of camera is carried out in the following manner:

2.1. Release retention screw on the head.

2.2. Install beam splitter on the head and tighten retention screw on the head.

2.3. Install adapter for camera in beam splitter and insert in into the port of beam splitter and tighten locking ring.


Fig. 33

- 1- Adjusting knob of image brightness
- 2- Button of video recording on SD- card
- 3 -Slot for SD card
- 4- HDMI connector
- 5- Power connector
- 6- USB connector (see fig.33)

MEASURES OF PRECAUTION:

- Never look into microscope to the sun, source of bright light and laser radiation – IT IS DANGEROUS FOR SIGHT!
- Don't disassemble camera without any professional help;
- Keep camera away from moisture, don't use under the rain;
- Don't drop camera, keep it away from excess load from other objects.

Don't apply undue exertion to lock and retention screws.

13.13 Built-in light filters.

Light filters are used to more in-depth and detailed study of the object of research. Thanks to its usage, it is possible to emphasize visually the more important object areas for the study.

Use of light filters:

- blue and green filters are used to enhance resolution;
- yellow are used to improve contrast.

14. FUSE REPLACEMENT

Fuses are located together with AC hole, situated close to power switch control. To open fuse block use flathead screwdriver. Two fuses are inside, one is functional and the other is spare accessory fuse. Replace burnt fuse with the functional one and put the block back.

15. TENSION ADJUSTMENT IN THE PROCESS OF ACCESSORIES USE


Figure 34.

After installation of additional accessories, new load on suspended arm, should be compensated by means of tension adjustment with the help of screw situated on suspended arm by means of it turning clockwise or counter clockwise.

16. DISINFECTION AND STERILIZATION

For diagnostics:

Moisten soft cotton cloth in aseptic solution (for example, Sagrotan – P); if it's necessary, wipe frequently contacted parts, such as control knobs, handles etc.

After operation:

Sterilized caps are contained in every part the surgeon contacts with in the process of operation. After use if necessary remove remains of blood or liquids from caps with the help of antiseptic solution. Sterilize them after each use of microscope.

Disinfection of camera and camera with built-in beam splitter.

In case of contamination of optical surfaces first it is necessary to blow off dust and small particles, then wipe with soft clean cloth, wetted with substance for lenses cleaning. Treatment of camera frame must be performed similar to cleaning of microscope surfaces.

Sterilization of caps:

ATTENTION: Don't use worn, damaged sterilized caps.

Sterilization of caps must be done by means of autoclaving in autoclave with specifications corresponding to sterilization.

The process of sterilization consists of some stages. Caps in sterilization caps are put into autoclave. Required parameters of sterilization are set (see manual for sterilizer)

and take the required time according to these specifications. At the end of exposure sterilization is considered finished.

Parameters of sterilization: Temperature: 134°C, time: 10 minutes; pressure of steam: 2.1±0.1 kgs/cm².

Shelf life of sterilized caps in sterilized boxes is 20 days.

Maximum level of contamination in steam condensate

Substance	Condensate, mg/kg
Evaporation residue	1,0
Silicon oxide SiO ₂	0,1
Iron	0,1
Cadmium	0,005
Lead	0,05
Residue of heavy metals (except iron, cadmium and lead)	0,1
Chlorine Cl	0,1
Phosphate P ₂ O ₅	0,1

17. CLEANING AND MAINTENANCE

17.1 Cleaning of optical surfaces

Coarse particles of dirt/dust should be removed from optic surfaces (objectives, eyepieces) with clean dry air.

Moisten soft cotton cloth with cleaning for lenses and wipe delicately the surface of lens, starting from center towards downside.

Install blanking cover on the objective and cover the head of the microscope with the dust cover while service interruptions.

17.2 Cleaning of mechanic surfaces

All mechanic surfaces of the equipment can be wiped with damp cloth. Don't use aggressive or abrasive agents for cleaning. To remove residual dirt it's possible to use any dishwashing detergent.

17.3 Maintenance

Maintenance of the product is allowed by skilled professional only.

Service maintenance is performed if necessary, contact aftersales service.

18. TROUBLESHOOTING

Problem	Possible reason	Solution
No light	Feeding cable is not connected	Connect cable
	Power switch control is turned off	Turn off the switch
	Fuse is blown	Replace the fuse
	Cable feeding defect	Replace feeding cable
	No voltage in electric line	Contact technical expert
	Disorders of electronic system	Contact service department
	Wrong connection of LED to the microscope or lamp	Adjust maximum intensity
Insufficient lightning	Low level of brightness	Adjust brightness with the help of brightness control
	Wrong connection of LED to the microscope or lamp	Adjust maximum intensity
	Optic fiber cable defect (no uniform lightning)	Replace LED
LED lamp of lightning system doesn't function, at the same time fan is working LED lamp constantly switches off and switches on in the process of work	Brightness control	
	Poor contact of the lamp mode	Install lamp module in the right way
	Faulty LED lamp	Use replacement lamp
	Vent slots are closed or dirty	Vent slots should be clean, clean them when necessary.
	Temperature fuse in the lamp housing is dirty	Clean fuse with dry brush; blow off if necessary.
	Fan defect. Malfunction of electronic system.	Contact service department
Movement of suspended arm up and down is too tough	Hydraulic lock of suspended arm is excessively tight.	Release hydraulic lock of suspended arm
Microscope is not stable	Brakes for wheels are not applied	Hit the brake pedal
Image is not seen in the image field	Magnifier is adjusted according to the index in the wrong way	Indicate magnification index in the right way

19. TECHNICAL SPECIFICATIONS

19.1 Zooming scales of microscopes (the way of changing for zooming scale for fixed objective is done by step-by-step method) are specified in Table 2.

Table 2

Eyepiece	Zooming switch position	Zooming scale, N-fold $\pm 5\%$
Eyepiece 10×	0,4 ^x	4
	0,63 ^x	6,3
	1,0 ^x	10
	1,6 ^x	16
	2,5 ^x	25
Eyepiece 12.5×	0,4 ^x	5
	0,63 ^x	7.875
	1,0 ^x	12.5
	1,6 ^x	20
	2,5 ^x	31,25

19.2 The diameter of illuminated area shall be not less than 85 mm.

19.3 Interpupillary distance should be regulated in the range from 55 mm p to and including 75 mm.

19.4 Angular inclination of eyepieces in binocular tubes should be in accordance with Table 3.

Table 3

Type	Angular inclination of eyepieces, $\pm 3^\circ$
Binocular tube 90	90
Binocular tube 45	45
Binocular tube 210	0-210

19.5 Interval for eyepieces dioptric adjustment shall be in the range from -5 to +5.

19.6 The focal length of objectives should be in accordance with Table 4.

Table 4

Type	Focal length, mm, $\pm 2\%$
Objective 175	175
Objective 200	200
Objective 250	250
Objective 300	300
Objective 400	400

19.7 Object-side linear field of microscopes should be in accordance with Table 5

Table 5

Zooming switch position	Object-side linear field of microscope, $\pm 0,1$ mm				
	Objective 175	Objective 200	Objective 250	Objective 300	Objective 400
0,4 ^x	46,39	52,94	66,18	79,41	105,88
0,63 ^x	31,00	35,30	44,10	53,00	70,60
1,0 ^x	18,58	21,21	26,51	31,81	42,42
1,6 ^x	11,60	13,25	16,70	19,90	26,50
2,5 ^x	7,42	8,47	10,59	12,71	16,94

19.8 Microscopes resolving power should be in accordance with Table 6.

Table 6

Eyepiece	Zooming switch position	Resolving power, lines/mm
Eyepiece 10×	0,4 ^x	16
	0,63 ^x	25
	1,0 ^x	40
	1,6 ^x	64
	2,5 ^x	100
Eyepiece 12.5×	0,4 ^x	20
	0,63 ^x	30
	1,0 ^x	50
	1,6 ^x	80
	2,5 ^x	125

19.9 Light source specification: L.E.D. (source voltage is 12V; capacity is 50 W).

19.10 Depends on type, built-in light filters should be in accordance with Table 7

Table 7

Type	Filter	Spectral transmission, nm
PRIMA OPH	Blue	460±5÷490±5
	Green	510±5÷560±5
PRIMA ENT	Green	510±5÷560±5
	Yellow	580±5÷600±5
PRIMA DNT	Green	510±5÷560±5
	Yellow	580±5÷600±5
PRIMA GN	Green	510±5÷560±5

19.11 Variations for illuminating intensity should range from 15000±500 to 80000±1000 lux.

19.12 Motorized XYZ motion:

- along Z-direction, it shall be equal to 40±0,5 mm (for OPH, ENT, DNT types);
- along X-direction it shall be equal to 50±0,5 mm (for OPH types);
- along Y-direction it shall be equal to 50±0,5 mm (for OPH types).

19.13 The difference in zooming of the right and the left optical systems should not exceed 1,5%.

19.14 The divergence of microscopes right and left optical systems shall be: no more than 15' in vertical direction; no more than 10' in horizontal direction.

19.15 Focal difference of microscopes right and left optical systems shall not exceed $1,5D_F$, Deviation of focal planes during zooming adjustment shall not exceed $3D_F$, where: $D_F = \lambda / 2NA^2 + 1/(7M_{TOTVIS}NA)$, where D_F – is depth of objective field, M_{TOTVIS} – is total visible zooming, λ – is wave-length, NA – is numerical aperture.

19.16 Images angling difference in microscopes right and left optical systems shall not exceed 2°.

19.17 Difference in height of exit eye pupil in right and left optical systems shall be equal to 1.5 mm for zero diopt. on diopter scale.

19.18 Calibration tolerance of diopter scale shall be in the range of: ±0,25 diopt. for zero diopt. on diopter scale.

19.19 Re-alignment of microscopes moving parts:

19.19.1 Vertical adjustment of microscopes should be equal to 550±5 mm.

19.19.2 Microscopes should be rotating for the entire circle ±360°.

19.19.3 Microscopes should be rotating left and right for not less than ±150° about the rotary arm.

19.19.4 Microscopes should be rotating for not less than ±90° about the inclining arm.

19.19.5 Microscopes should be rotating up and down about the suspended arm for not less than ±135°.

- 19.20 The force, required to move the microscope, should be no more than 80 n.
- 19.21 Wheels brakes of mobile base should be turned on when applying pressure is not exceeding 150 n.
- 19.22 The set up height of suspended arm should not be altered under load of 7 kg.
- 19.23 The degree of dust and moisture protection of the microscope must be of IP20 in accordance with GOST 14254-96.
- 19.24 The degree of dust and moisture protection of the pedal must be of IPX6 in accordance with GOST14254-96.
- 19.25 Fuse capacity for the fuse of 5 x 20 is equal to 1A/250V.
- 19.26 AC of 220V voltage and 50 Hz of frequency should power microscopes. Microscopes shall be functioning for power-supply variation in the range of 198 -242 V.
- 19.27 The maximum of microscopes power consumption should be no more than 130 Watts.
- 19.28 Beam splitter, dual beam splitter 50:50 and inclined double beam splitter 50:50 should ensure separation of luminous flux in a ratio of 50/50.
- 19.29 Dual beam splitter 80:20 and inclined double beam splitter 80:20 should ensure separation of luminous flux in a ratio of 80/20.
- 19.30 Measurements:
- 19.30.1 Overall dimensions of microscopes are:
- PRIMA OPH, PRIMA ENT, PRIMA DNT (L × W × H), the dimensions shall be: (1415x560x1756)±5 mm.
 - PRIMA GN (L × W × H), the dimensions shall be: (1415x560x934)±5 mm.
- 19.30.2 The height up to suspended arm of microscope is:
- for types of PRIMA OPH, PRIMA ENT, PRIMA DNT it shall be: 1466±3 mm;
 - for types of PRIMA GN it shall be: 650±3 mm.
- 19.30.3 The height up to the head end of microscope is:
- for types of PRIMA ENT, PRIMA DNT, PRIMA OPH it shall be 1146±3 mm
 - for types of PRIMA GN it shall be: 991±3 mm.
- 19.30.4 The length of rotary arm shall be: 525±2 mm.
- 19.30.5 The length of suspended arm shall be: 600±2 mm.
- 19.30.6 Dimensions of the wheeled mobile base (W × H) shall be: (580x560)±3 mm.
- 19.30.7 The diameter of the base shall be 60±2 mm.
- 19.30.8 The length of overhead bracket with column shall be: 1677±3 mm.
- 19.30.9 Dimensions of the wall bracket (W × H) shall be: (356x148)±3 mm.
- 19.30.10 The length of the desk-top bracket shall be: 768±2 mm.
- 19.30.11 Network cable length must be 2±0.1 m.
- 19.30.12 Cord length of the pedals must be 2±0.1 m.
- 19.30.13 Feeding cable for camera with built-in beam splitter should be 0,5±0.05 m.
- 19.30.14 USB/micro USB cable length must be 1±0.05 m.
- 19.30.15 HDMI /mini HDMI cable length should be 1±0.05 m.
- 19.30.16 Dimensional specifications of the camera (L × W × H), it shall be: (29x29x46) ±0,1 mm.
- 19.30.17 Dimensional specifications of camera with built-in beam splitter (L × W × H), it shall be: (150x25x100) ±0,1 mm.
- 19.31 The weight:
- 19.31.1 Microscope shall weight not more than 72 kg.
- 19.31.2 The pedal weight shall not exceed 800 g.
- 19.31.3 Weight of camera it shall be 200±10 g.
- 19.32.4 Weight of camera with built-in splitter it shall be 350±15 g.
- 19.32 Specification of cameras shall be:
- 19.32.1 Specification of camera with built-in beam splitter should be:
- Sensor device – is 1 / 2.5 "5,0 MП CMOS, active area – is 5.70 mm x 4,28 mm, resolution – is 2592 x 1944 pixels; numeralization – is 12 bits of RGB, ratio of signal/interference – is 68 dB (measured at 10 ms of exposure; video images frequency - is 30 fps (1920 x 1080 pixels), maximal exposure - is 3 sec, minimal exposure is 1 sec, interface – is SD card, it has USB 2.0, HDMI ports,

power supply – is 5V/2A, watt consumption – is 1.8 W; it is compatible with Microsoft Windows XP / Vista / 7/8/10.

19.32.2 Specification of camera shall be: max resolution – is 5 MP (2592x1932 pixels), it is USB connected; numeralization – is 12 bits of RGB, interface – is C-mount, 0.45 adapter, sensor device - is CMOS, it has pixel of 2.2 μm x 2.2 μm , ratio of signal/interference - is 40 dB, display space – is 5.70 mm x 4.28 mm , maximal exposure - is 76 sec; video images frequency - is 7 fps for full resolution; it is proportional to size of display, it is compatible with Microsoft Windows XP / Vista / 7/8/10.

19.33 The corrected sound power level, generated during use of microscopes should not exceed 55 dBA.

19.34 Software PixelPro is used for device functioning.

19.34.1 Specifications of the software:

- Software version: V11.0

- Production date: 2015

- Program size: 20 Mb

19.34.2 System requirements

- Operating system: Windows Vista, XP, 7, 8, 10;

- Processor speed: 1.7 GHz and higher;

- Core memory: 512 Mb or more;

- Drive mechanism: DVD-ROM;

- Hard disk drive: no less than 400 Mb.

19.35 The USB key is used to authenticate the software user.

USB key has the following characteristics:

- interface USB 2.0

- USB key size (LxWxH) should be: (70x21x10) $\pm$ 1 mm.

19.36 For recording images on camera with built-in beam splitter SD memory card slot is used.


19.37 The wrench set for assembling consists of 5 wrenches of the following sizes: 2.5; 3; 4; 5; 6 mm.

19.38 Setting time of operational mode of microscopes must be no more than 5 s.


19.39 Microscope must provide continuous operation no less than 8 hours per day.

19.40 Microscope must endure the impact of sinusoidal vibration in vertical range of frequency band from 10 to 55 Hz with displacement amplitude up to 0.15 mm.


Dimensions and movement of mobile parts PRIMA ENT, DNT

PRIMA OPH


PRIMA GN


Types of support


On the wall

Overhead

Desktop

20. OPERATION, TRANSPORTATION AND STORAGE CONDITIONS

For microscope

Operation	Temperature Humidity (without condensation) Pressure	+10°C....+40° C 30%.....90% 700gPa.....1,060gPa
Transportation and storage	Temperature Humidity (without condensation) Pressure	-50°C....+50° C 10%.....100% 500 gPa.....1,060 gPa

For camera and cameras with built-in beam splitter

Operation	Temperature Humidity (without condensation) Pressure	+10°C....+40° C 30%.....90% 700 gPa.....1,060 gPa
Transportation and storage	Temperature Humidity (without condensation) Pressure	-50C....+50° C 10%.....100% 500 gPa.....1,060 gPa

Wrap details and units of microscope into package paper and put them into stowage box. All stowage boxes must be sealed up in plastic bags with silica gel. Stowage boxed in plastic bags must be packed in transportation box. Pillar and pedal must be packed in a separate transportation box. When transported microscope and its accessories must be put into boxes in such a way not to mix them, when boxes are shaken. Transportation by all means of roofed transport is considered acceptable. After transportation (or storage) at low temperatures microscope in transportation pack must be left on the room with temperature from +15 to +25°C no less than 10 hours. Only after that it can be unpacked and used. Keep microscope under the cover. All replacement parts, accessories and instrument must be kept in boxes.

21. INSTALLATION AND USE OF Pixel Pro SOFTWARE

Software Pixel Pro is used for connection of the camera with built-in beam splitter and the camera without built-in beam splitter to the microscope and makes it possible to watch, record and record video in real-time mode, to do pictures of objects observed in the microscope.

1) Installation of software.

- Insert installation disk into CD-DVD drive of your personal computer.
- Click "setup.exe" file two times.
- Select the language for installation, the program is available in the following languages (English, Русский, Español, Français, Romana, Deutsch, Italiano)
- Touch the button "Accept", isnatll DotNet Framework 3.5. Please wait until installation is over and touch the button "Next".
- Choose "I agree" and touch the button "Next".
- Choose Folder for installation of default program "C:\program files\Labomed\PixelPro Software" and touch the button "Next".
- Touch the button "OK" and install FFDSHow CODEC package.
- Save FFDSHow default settings and then touch the button "Next" to continue.
- After successful installation, touch the button "Close" to complete the installation.

- Create a shortcut on the desktop to launch the program.
- Before open the program connect USB key to USB port of personal computer, which is used to authenticate the user of software.
- To launch the program click two times "PixelPro icon".

Image acquisition

- 1) Click button "Acquisition" to capture image. Decreased image will appear in the lower panel from the captured frame.
- 2) To view the image, go to the folder, in which image was saved.
- 3) Click close button (in the upper right corner) in the window of image acquisition to close the window.

Video capture

1. In the open window: click button "recording"  to record video of selected area.
2. Click "stop" button  to stop recording. If necessary you can click "start recording" button again and start new recording.
3. Video is stored in the selected folder, where you can view it.
4. Click close button (X in the upper right corner) to close the window.

Legal obligations of software manufacturer

- 1) Consumer claims for indemnification or indemnification of vainly incurred expenses, regardless legal nature of the claim, are considered in accordance with this provision. By virtue of Law on responsibility of manufacturer for the product, the responsibility remains in force.
- 2) Company Labo America Inc. bears unlimited liability for infliction of life harm, corporal hurt or personal injuries.
- 3) Company Labo America Inc. bears unlimited liability for losses, caused by deliberate or roughly negligent conduct of the employees of Labo America Inc. Company Labo America Inc. bears liability for losses cause by negligence, only if duties are violated, performance of which is essential to achieve the purpose of the Contract (duties essential for the contract). The duties essential for the contract are such duties, which must guarantee proper performance of the contract or which must be regularly performed by the Consumer. In case of essential duties violation the liability of company Labo America is limited with the damage, which can occur in the course of usual transmission of Software or works on technical and service maintenance (clause 22).
- 4) Liability for data loss is limited with expenses for recovery, which are possible in case of regular creating backup files corresponding to risks.
- 5) Clauses from 1 to 4 are valid for personal responsibility of employees, members, representatives and authorized companies of Labo America Inc.
- 6) Defects of supplied Software, including documentation, are corrected by company Labo America Inc. during guarantee period of one year from, the moment of first use of the Software after written notification of the User. Correction of defect is done by means of free correction or replacement at the choice of company Labo America Inc.

Consumer rights

Right for copying and access protection and program change.

- 1) The user can copy this Software, if it is necessary for the use of program. Installation of the program from original data carrier on mass tape storage of used equipment hardware, as well as program loading into frame memory belong to necessary types of copying.
- 2) The user has the right for copying of Software within the net used in clinic.

- 3) The consumer has the right for copying of Software for safety consideration. Creating and storage of one backup file is admitted.
- 4) Copying, drawing out of provisions given above, which printout of programming code is related to, is prohibited.
- 5) It is prohibited to delete or change copyright notice, series number and other program identifiers.
- 6) Conversion of program code into other forms of code (build), as well as other kinds of reproduction of other stages of Software manufacturing (reverse engineering), including program change is not admitted.

22. INSTALLATION AND MAINTENANCE OF MICROSCOPE

Assembling

Assembling of microscope is performed by personnel with corresponding education and skills authorized by LLC "HEALTHY WORLD" for microscope assembling.

Daily checks

To provide a good image avoid dirt contamination of ocular lenses and/or don't touch them. When the device is not used, cover it with dust sheet. If objective or ocular lenses are dirty, clean them following the instructions "CLEANING OF OPTIC SURFACES" on page 29. Check if air holes are not closed. Check electric connectors. When the device is not used for a long time, always switch it off (Power switch control ON/OFF).

User service

Name	Timing	Content
Cleaning	After every use of the device	• After usage of the device, remove rests of blood or liquids from its external parts.
Sterilization	After every use of the device	• After usage of the device, remove rests of blood or liquids from its sterilized caps. Then clean and sterilize them.
Cleaning of lenses	If necessary	• Objective lenses • Ocular lenses
Functional check	Before every use of the device	• Check if all knobs, screws and fixing mechanisms are tightened • Check balance of suspended arm • Before use of the device, after long-term storage, make sure it works safely.

Maintenance service

Name	Timing	Content
Cleaning of all units	Every 12 months	• Clean external parts of the device • Clean optical system
Check of safe operation	Every 12 months	• Suspended arm: up-down motion, balance • Functioning of knob of rotary arm lock • Functioning of knob of suspended arm lock • Functioning of focusing assembly • Mechanism light source change • Light system • Function of light adjustment • Functioning of pedal lever

23. MARKING

- Serial No.	S/N
- " Manufacturing date "	 (year/month)
- MANUFACTURER	
- sign "Current type"	
- "Type of working part"	
- Degree of dust and moisture protection	IP20
- «Do not throw out! Special disposal is required.»	
- CE MARKING	
- Safety sign, refer to user manual	

24 INFORMATION ON ELECTROMAGNETIC COMPATIBILITY

Information presented in this clause of extracts from the standard for electric and medical devices. It must be considered during assembling and combining of microscopes with other devices. In case of uncertainty refer to the full text of standard.

Table 1: Radiation

Recommendations and data of the manufacturer – electromagnetic emission		
Surgery microscope Labomed is intended for use in electromagnetic environment described below. User of surgery microscope Labomed must provide its use in such environment.		
Harmful radiation measurement	Level	Electromagnetic environment – recommendations
High-frequency radiation CISPr 11	group 1	Surgery microscope Labomed uses high-frequency energy for its functioning only. For this reason its high-frequency radiation is unperceivable, and probability of jamming for close electromagnetic devices is negligible.
High-frequency radiation CISPr 11	Class B	Surgery microscope Labomed is intended for use in institutions (non-residential area), as well as in buildings, directly connected to residential electric network, feeding dwelling houses.
Harmonic vibrations IEC 61000-3-2	Class A	
Voltage fluctuations/ pulse radiation IEC 61000-3-3	respectively	

Table 2: Interference immunity (all devices)

Recommendations and data of the manufacturer – electromagnetic emission			
Surgery microscope Labomed is intended for use in electromagnetic environment described below. User of surgery microscope Labomed must provide its use in such environment.			
High frequency disturbance test (standard)	IEC 60601 Test level	Level of correspondence	Electromagnetic environment – recommendations
Discharge of static electricity (ESD) IEC 61000-4-2	Contact discharge ±6 kV Air discharge ±8 kV	Contact discharge ±6 kV Air discharge ±8 kV	Floor – wooden, concrete or covered with ceramic tile. If the floor has synthetic coating, atmospheric relative humidity must be no less than 30%.
Instantaneous values of excess voltage IEC 61000-4-4	±2 kV for power networks ±1 kV for input and output wires	±2 kV for power networks ±1 kV for input and output wires	Power network quality must correspond to typical commercial or hospital environment.
Impulse voltage IEC 61000-4-5	±1 kV Differential-mode voltage ±2 kV Single-phase voltage	±1 kV Differential-mode voltage ±2 kV Single-phase voltage	Power network quality must correspond to typical commercial or hospital environment.
Voltage spikes, short breaks and electric power supply fluctuations IEC 61000-4-11	<5 % U_T (0.5 cycle) 40 % U_T (5 cycles) 70 % U_T (25 cycles) <5 % U_T for 5 s	5 % U_T (0.5 cycle) 40 % U_T (5 cycles) 70 % U_T (25 cycles) <5 % U_T for 5 s	Power network quality must correspond to typical commercial or hospital environment.
Power frequency (50/60 Hz) Magnetic field IEC 61000-4-8	3 A/m	Not applicable	

Note: U_T is alternating current system voltage before use of level check

Table 4: Interference immunity (except life sustaining devices)

Recommendations and data of the manufacturer – electromagnetic emission			
Surgery microscope Labomed is intended for use in electromagnetic environment described below. User of surgery microscope Labomed must provide its use in such environment.			
Electromagnetic environment – recommendations			
Portable and mobile electronic devices, including their wires, can be used at the distance from surgery microscope Labomed, no less than that calculated by equation for this carrier frequency of the transmitter.			
High frequency disturbance test (standard)	IEC60601 Test level	Level of correspondence	Recommended safety distance
Conductive perturbation influences HF IEC61000-4-6	$3B_{rms}$ 150kHz–80MHz	$3B_{rms}$	$d=2.4\sqrt{P}$ 150kHz–80MHz
Emitting perturbation influences HF IEC61000-4-3	3V/m 80MHz–2.5GHz	3V/m	$d=2.4\sqrt{P}$ 80MHz–2.5GHz
Where P is nominal rating power of the transmitter in watts (W) according to the data of manufacturer, and d is recommended safety distance in meters (m). Field intensity from fixed radiofrequency transmitters on all frequencies must be lower than applied level of radio frequency radiation according to electromagnetic mapping at the place of installation. ^b disturbances are possible close to devices with the following sign 			

- a Field intensity of fixed transmitters, for example, base stations, cellular networks and mobile ground radio stations, amateur radios, radio and television transmitters AM and FM theoretically can't be predicted with accuracy. To determine electromagnetic environment in regard to fixed transmitters, it is necessary to perform electromagnetic exploration of the place of installation. If measured field intensity in the place of operation of surgery microscope Labomed exceeds applied level of radiofrequency radiation (see above), microscope Labomed must be under observation to make sure it functions in the right way. If abnormalities are observed, additional measures must be required, for example, change of orientation or transportation of microscope Labomed to another place.
- b In frequency range 150 kHz – 80 MHz field intensity must be lower than 3 V/m.

Note:

1: for frequencies 80 MHz and 800 MHz value for higher frequency range is valid.

2: indicated data can't be applied in all situations. Extension of electromagnetic oscillation changes because of absorption and reflection properties of structures, people and subjects.

Table 6: Recommended safety distances (except life sustaining devices)**Recommended safety distances between portable and mobile high frequency communication devices and surgery microscope Labomed**

Surgery microscope Labomed is intended for use in electromagnetic environment, where perturbation influence of high frequency is under control. The owner of the user of surgery microscope Labomed can help to avert electromagnetic perturbation influence, keeping minimal distance between portable and mobile high frequency communication devices (transmitters) and surgery microscope Labomed, depending on power output of communication device (see below).

Maximum power output of a transmitter in W	Safety distance depending on frequency of a transmitter in meters
	150 kHz – 2.5 Ghz $d = 2.4 \sqrt{P}$ in m
0.01	0.24
0.1	0.76
1	2.4
10	7.6
100	24.0

For transmitters, maximum power of which is not presented in the table, it is possible to calculate recommended safety distance (d) in meters (m) by the equation for corresponding line, where P is maximum power of a transmitter in watts (W) in accordance with data of transmitter manufacturer.

Note: 1: indicated data can't be applied in all situations. Extension of electromagnetic oscillation changes because of absorption and reflection properties of structures, people and subjects.

25. MANUFACTURER WARRANTY

The manufacturer guarantees the compliance of the microscope with all the requirements of the present regulatory document while observing working, transportation and storage conditions by the customer.

Guarantee period of the microscope is 12 months from the date of sale.

Guarantee period of software is 12 months from the date of the first use.

For warranty or routine repair refer to address: Limited Liability Company "HEALTHY WORLD" 191014, Saint Petersburg, Ligovsky prospekt, 3/9, apt.1, phone: . 8(812) 579-74-97

26. DISPOSAL CONSIDERATIONS

Microscopes should be prepared and delivered to special licensed institutions dealing with disposal of electronic devices. Disposal should be performed in accordance with effective local regulations. According to Sanitary Rules and Regulations 2.1.7.2790 microscope belongs to class A – epidemiological non-hazardous wastes. Microscope must be sanitized in accordance with Methodological guidelines MU-287-113 dated 30.12.1998.

27. MANUFACTURER INFORMATION:

Name and address of the manufacturer:

Labo **America** Inc., 920 Auburn Court, Fremont, CA 94538, USA

Medical device's production site:

Labotech **Microscopes** India Private Ltd, 19, HSIDC Industrial Area, Ambala Cantt,
Haryana – 133006 India

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