Appendix to the IFU

Gyda L 14-Dec-2015

Gloria Nguyen
 Regulatory Affairs Specialist
 Cardiac Rhythm Heart Failure

Gloria Nguyen Dec-14-2015

Implantable Pulse Generator Vitatron A20SR, model A20A1

1. Shelf life	18 months
2. Indications	These Vitatron A-series implantable pulse generators (IPGs) are indicated for use to improve cardiac output, prevent symptoms, or protect against arrhythmias related to cardiac impulse formation or conduction disorders. These devices are indicated for use in patients who are experiencing exercise intolerance or exercise restrictions related to an arrhythmia. Using rate response modes may restore heart rate variability and improve cardiac output.
3. Requirements for the protection of the environment in the application of medical devices. Procedure for utilization/liquidation of MD	Explant the implantable device postmortem. In some countries, explanting battery-operated implantable devices is mandatory because of environmental concerns; please check the local regulations. In addition, if subjected to incineration or cremation temperatures, the device may explode. Please return explanted devices to Vitatron for analysis and disposal.
4. Classification	1) Class depending on the potential risk of a medical product, in accordance with the nomenclature classification - CLASS III / 3 2) Class of electrical safety, laser safety, and other types of hazards (if any) –Not applicable for implantable devices

Vitatron Holding B.V.
 Mailing address : P.O. Box 1220
 6201 MP Maastricht
 Visiting Address : Endepolsdomein 5
 6229 GW Maastricht


Tel : +31 43 356 6551
 Fax : +31 43 356 6542
 K.v.K. Zuid Limburg : 090.30.178
 V.A.T. : NL 0013.07.915.B.01

5. Package content

- Implantable pulse generator
- Torque wrench
- Product documentation*

* Product documentation : **Manual, Warranty, Individual implantation card**, Appendix to the Manual (IFU) (not in package content, shipped upon request), Reference Manual (not in package content, visit <http://manuals.medtronic.com/>)

6. Technical parameters and Leads compatibility

Physical characteristic	Vitatron A20A1
Dimensions (mm) [Height × Length × Thickness]	40.2 × 42.9 × 7.5
Mass (g)	21.5
Volume (cc)	9.7
Connector	IS-1 single chamber
Case material	Titanium
Connector materials	Polyurethane, silicone rubber
Power source	Lithium iodine
Average available capacity (Ah)	0.89
Initial voltage (V)	2.8

1) Battery and power consumption for Model A20

Vitatron Holding B.V.

Mailing address : P.O. Box 1220
6201 MP Maastricht
Visiting Address : Endepoelsdomein 5
6229 GW Maastricht

Tel : +31 43 356 6551

Fax : +31 43 356 6542

K.v.K. Zuid Limburg : 090.30.178

V.A.T. : NL 0013.07.915.B.01

Parameter	Capability	Notes
Model	Sigma 213	
Type	Lithium-Iodine	
Manufacturer	Medtronic Energy and Component Center	
Voltage	2.8 V	
Usable capacity	0.92 Ah	Capacity available from BOS to EOS
Minimum capacity at RRT/ERI	0.06 Ah	Minimum capacity is measured at nominal settings, 500 Ω load in uni-polar or bipolar polarity, 100% pacing
Device input impedance	> 150 k Ω	Device input impedance is measured to CENELEC standard (EN 45502-2-1)

3) Lead and connector compatibility

Warning: Verify lead and connector compatibility before using a lead with this device. Using an incompatible lead may damage the connector, result in electrical current leakage, or result in an intermittent electrical connection. Note: If you are using a lead that requires an adaptor for this device, please contact your Vitatron representative for information about compatible lead adaptors. Select a compatible lead. Refer to Table .

Table Lead and connector compatibility

Model	Connector	Polarity	Primary leads
A20A1A1		Bipolar/Unipolar	IS-1a BI

Vitatron Leads for use with Vitatron A-Series IPGs

Device Name	Device Model Number(s)
Impulse II	IHP09B, IHP09JB
Crystalline	ICL08, ICL08B, ICL08JB, ICM09, ICM09B, ICM09JB
Crystalline ActFix	ICF09B, ICQ09B
Excellence PS+	IMK49B, IMK49JB
Brilliant S+	IMW14Q, IMW15Q, IMW16Q, IMW17Q, IMW18Q

IS-1 refers to the International Connector Standard (see Document No. ISO 5841-3) whereby pulse generators and leads so designated are assured of meeting the electrical and mechanical parameters specified in the IS-1 International Standard.

Requirement ISO 5841	pass/ not pass
1. Connection lead:	
A) location and number of insulating rings	pass
B) size	pass

Vitatron Holding B.V.

Mailing address : P.O. Box 1220
6201 MP Maastricht
Visiting Address : Endepolsdomein 5
6229 GW Maastricht

Tel : +31 43 356 6551
Fax : +31 43 356 6542
K.v.K. Zuid Limburg : 090.30.178
V.A.T. : NL 0013.07.915.B.01

B) continuity of the electrical circuit	continuous - pass
2. strength of the introduction / extraction flow gauge	<14H - pass
3. electrical impedance between the conductive parts	≥50 kOhm - pass
4. Deformation and its influence on the connection terminal on	<14H - pass
5. Impact of bipolar clamping screw connector sockets on Unipolar connection terminal	pass
6. Socket outlet A) size B) max force introduction	pass ≤9 H - pass
7. Labeling sign («IS-1» clear marking)	pass

Magnet operation and Recommended replacement time

Magnet operation		Indicators of RRT/ERI status	
Without magnet	With magnet	Without magnet	With magnet
VV/AAI	VOO/AOO at 85 min ⁻¹	VVI at 65 min ⁻¹	VOO at 65 min

Note: The device does not respond to the application of a magnet for one hour after the use of a pro unless the session is ended with the command option to immediately clear data collected in the device. command for ending a session allows the device to retain collected data for one hour.

Magnet operation

Magnet operation may be initiated by performing one of the following actions:

- Placing a transtelephonic magnet over the IPG
- Performing a Magnet Test from the programmer

Note: The device does not respond to the application of a magnet for one hour after the use of a programmer unless the session is ended with a command option to immediately clear data collected in the device. The default command for ending a session allows the device to retain collected data for one hour.

The following responses occur automatically during magnet mode operation:

- The device performs a Threshold Margin Test (TMT) to allow a check for loss of capture. A TMT is performed at 100 min⁻¹ with amplitude reduced by 20% on the third pulse.

Warning: Loss of capture during a TMT at a 20% reduction in amplitude indicates an inadequate safety margin. As soon as possible, perform a pacing threshold test and reprogram outputs to establish a 2:1 stimulation safety margin.

Note: When Transtelephonic Monitor is set to On, TMT is delayed for five seconds to enhance communication with Transtelephonic equipment.

- The magnet rate becomes 65 min⁻¹ when the Recommended Replacement Time (RRT/ERI) is reached or when a full electrical reset occurs.
- The magnet rate for normal operation switches to a rate of 85 min⁻¹ and the appropriate asynchronous mode, see the following table for mode information.

Programmed mode	Magnet mode
Dual chamber	DOO
VDD	VOO
Single chamber	VOO/AOO

Recommended Replacement Time (RRT/ERI)

The Recommended Replacement Time (RRT/ERI) warns when the device battery is nearing depletion. When the battery voltage has dropped below a defined limit,

Vitatron Holding B.V.

Mailing address : P.O. Box 1220
6201 MP Maastricht
Visiting Address : Endepotsdomein 5
6229 GW Maastricht

Tel : +31 43 356 6551
Fax : +31 43 356 6542
K.v.K. Zuid Limburg : 090.30.178
V.A.T. : NL 0013.07.915.B.01

the device sets the RRT/ERI status and reverts to VVI operation at a rate of 65 min⁻¹. At most programmed settings, the majority of devices will function for a minimum of three months after RRT/ERI. The device's indicators of RRT/ERI status are listed in the following table:



Operation	Indicator
Nonmagnet mode	VVI mode at 65 min ⁻¹
Magnet mode	VOO mode at 65 min ⁻¹
Telemetry	Replacement message on the programmer
Battery/Lead information	Replacement message and displayed battery voltage on programmer

Projected service life Vitatron A20SR:

Pacing	A Amplitude, V Amplitude	Rate, Pulse Width	Lead impedance	
			500 Ω	1000 Ω
			Longevity (years)	
SSIR or SSI, 0%	2.0 V	60 min ⁻¹ , 0.4 ms	11.3	11.3
	2.5 V		10.4	10.4
	3.5 V		11.1	11.1
SSIR or SSI, 50%	2.0 V	60 min ⁻¹ , 0.4 ms	10.4	10.8
	2.5 V		9.5	9.9
	3.5 V		9.0	9.9
SSIR or SSI, 100%	2.0 V	60 min ⁻¹ , 0.4 ms	9.6	10.4
	2.5 V		8.7	9.5
	3.5 V		7.5	8.9
SSIR or SSI, 0%	2.5 V	70 min ⁻¹ , 0.5 ms	10.3	---
	5.0 V		10.2	---
SSIR or SSI, 100%	2.5 V	70 min ⁻¹ , 0.5 ms	8.0	---
	5.0 V		4.6	---
SSIR or SSI, 100%	5.0 V	70 min ⁻¹ , 1.0 ms	3.1	---
SSIR or SSI, 100%	5.0 V	100 min ⁻¹ , 1.0 ms	2.3	---

4)Wrench -Torque

Mass: 8g

The torque is limited to 12 to 17 in-oz. during first clockwise revolution; it will click when torque value in this range is achieved.

Vitatron Holding B.V.

Mailing address : P.O. Box 1220
6201 MP Maastricht
Visiting Address : Endepolsdomein 5
6229 GW Maastricht

Tel : +31 43 356 6551
Fax : +31 43 356 6542
K.v.K. Zuid Limburg : 090.30.178
V.A.T. : NL 0013.07.915.B.01

7. Device Materials	Material (longterm contact with internal environment of the body)	Device Component (IPG Vitatron A20)
	Titanium, grade 1 Fe to 0,12 C to 0,05 Si to 0,08 N to 0,04 Ti 99,5 - 99,99 O to 0,1 H to 0,008	Left-hand shield Right-hand shield
	Silicone Rubber MED 2000	Adhesive
	Polyurethane, 75 durometer Elasthane	Connector module
	ETR silicone rubber, 35 durometer	Split grommet
	Wrench -Torque (not in contact with the human body) Materials: 1)Haft - Polyetherimide (G.E. Ultem 1100F-8104) 2)Hex Shaft - 455 Corrosion resistant steel	
8. Software	Programmers and software – Use only Medtronic programmers and Vitatron application software to communicate with the device. Programmers and software from other manufacturers are not compatible with Vitatron devices. Version - VSF20 Date of issue – 2013 Safety class – C (regarding IEC 62304) User Characteristics - The primary audience for the programmer software is intended to be clinicians and nurses for the purpose of programming and interrogating the IPG. A secondary audience are the in-house Medtronic software and hardware development engineers and verification test engineers. The environment that the IPG application software and the 2090 programmer will be used in is a follow- up setting. Session Description This feature provides the capability of defining a patient session and controlling activities that occur at the start of a patient session. A patient session is generally defined as the period of use of the programmer for one patient involving one implanted device during which data collection and/or programming of the device is performed. This is termed an “active” or “programming” session. Common to all NG IPG sessions is the need to identify the device with which the programmer is communicating either through an auto-id or manual model selection. The device may also be identified during the process in which demographics information is provided to fill in registration form information. The information gathered during a session will be stored in a patient session record. This	

Vitatron Holding B.V.

Mailing address : P.O. Box 1220

6201 MP Maastricht

Visiting Address : Endepolsdomein 5

6229 GW Maastricht

Tel : +31 43 356 6551

Fax : +31 43 356 6542

K.v.K. Zuid Limburg : 090.30.178

V.A.T. : NL 0013.07.915.B.01

	data will include the patient demographic information. At the start of a NG IPG session, after the device has been identified and an initial interrogate has occurred, the state of the device (i.e. the state of the device parameters and its components relative to patient history) will be analyzed through the quick scan diagnostic procedure (TACHY). The state of the device also will be saved at this time within a diagnostics. The patient session record will become the repository of other data collected during the patient session (i.e. test results, episodes, etc.) If the session is communicating with an implantable device, model identification of the device will occur prior to the start of a session. The model information will be used by the session to make available only the model applicable aspects of things like protocol and nominal selection lists. Connection type Programmer - Device Device uses the Telemetry (175kHz) for connection with programmer through the programming head Range of coverage – 6cm.
9. Methods and conditions of disinfection and sterilization	Sterile. Single use only. Products cannot be resterilized
10. Operating and Storage conditions	<u>Checking and opening the package</u> – Before opening the sterile package tray, visually check for any signs of damage that might invalidate the sterility of the package contents. <u>If the package is damaged</u> – The device packaging consists of an outer tray and inner tray. Do not use the device or accessories if the outer packaging tray is wet, punctured, opened, or damaged. Return the device to Vitatron because the integrity of the sterile packaging or the device functionality may be compromised. This device is not intended to be resterilized. <u>Dropped device</u> – Do not implant the device if it has been dropped on a hard surface from a height of 30 cm or more after it is removed from its packaging. <u>Use by date</u> – Do not implant the device after the “Use by” date on the package label. Battery longevity may be reduced. <u>For single use only</u> – Do not resterilize and reimplant an explanted device that has been contaminated by contact with body fluids. <u>Sterilization</u> – Vitatron has sterilized the package contents with ethylene oxide before shipment. This device is for single use only and is not intended to be resterilized. <u>Avoid magnets</u> – To avoid damaging the device, store the device in a clean area away from magnets, kits containing magnets, and any sources of electromagnetic interference.

Vitatron Holding B.V.

Mailing address : P.O. Box 1220
 6201 MP Maastricht
 Visiting Address : Endepolsdomein 5
 6229 GW Maastricht

Tel : +31 43 356 6551
 Fax : +31 43 356 6542
 K.v.K. Zuid Limburg : 090.30.178
 V.A.T. : NL 0013.07.915.B.01

	Temperature limits – Store and transport the package between -18°C and $+55^{\circ}\text{C}$. Electrical reset may occur at temperatures below -18°C. Device longevity may decrease and performance may be affected at temperatures above $+55^{\circ}\text{C}$. Device temperature – Allow the device to reach room temperature before it is programmed or implanted. Device temperature above or below room temperature may affect initial device function.
11. Operational condition	External defibrillation equipment – Keep external defibrillation equipment nearby for immediate use whenever arrhythmias are possible or intentionally induced during device testing, implant procedures, or post-implant testing. External pacing instrument – Keep an external pacing instrument available for immediate use. When the lead is disconnected, pacemaker-dependent patients are without pacing support. Allow the device to reach room temperature before it is programmed or implanted. Device temperature above or below room temperature may affect initial device function. Cellular phones (Vitatron Models A20A1) This device contains a filter that prevents most cellular telephone transmissions from interacting with device operation. To further minimize the possibility of interaction, take the following precautions: • Maintain a minimum separation of 15 cm between the device and the cellular telephone, even if the cellular telephone is not on. • Maintain a minimum separation of 30 cm between the device and any antenna transmitting above 3 W. • Hold the cellular telephone to the ear farthest from the device. This device has been tested to the frequency ranges used by common cellular phone transmission technologies. Based on this testing, the device should not be affected by the normal operation of cellular phones using such technologies. Electromagnetic interference (EMI) (Vitatron Models A20A1) Instruct patients to avoid devices that generate strong EMI. Electromagnetic interference may cause device malfunction or damage such as prevention of programming, detection, or therapy delivery. The patient should move away from the EMI source or turn off the source because this usually allows the device to return to its normal mode of operation. EMI may be emitted from these sources: • high-voltage power lines • communication equipment such as microwave transmitters, linear power amplifiers, or high-powered amateur transmitters • commercial electrical equipment such as arc welders, induction furnaces, or resistance welders Home appliances that are in good working order and properly grounded do not usually produce enough EMI to interfere with device operation. There are reports of temporary disturbances caused by electric hand tools or electric razors used directly over the implant site. Electronic article surveillance (EAS) (Vitatron Models A20A1)

Vitatron Holding B.V.

Mailing address : P.O. Box 1220

6201 MP Maastricht

Visiting Address : Endepolsdomein 5

6229 GW Maastricht

Tel : +31 43 356 6551

Fax : +31 43 356 6542

K.v.K. Zuid Limburg : 090.30.178

V.A.T. : NL 0013.07.915.B.01

	Electronic article surveillance equipment, such as retail theft prevention systems, may interact with devices and result in inappropriate therapy delivery. Advise patients to walk directly through an EAS system and not remain near an EAS system longer than necessary. Static magnetic fields(Vitatron Models A20A1) Patients should avoid equipment or situations where they would be exposed to static magnetic fields greater than 10 gauss or 1 mT. Static magnetic fields may cause the device to operate a synchronously. Sources of static magnetic fields include, but are not limited to, stereo speakers, bingo wands, extractor wands, magnetic badges, or magnetic therapy products. Device temperature – Allow the device to reach room temperature before it is programmed or implanted. Device temperature above or below room temperature may affect initial device function. Mode of operation, and the implant technique used by the physician. Basic rate, test pulse rate, pulse duration, and pulse amplitude remain within expected tolerances when the device temperature is between 20 °C to 43 °C. Protection of the IPG from damage caused by atmospheric pressure changes: Implantable parts of the IPG withstand changes of pressure. all devices passed the requirements in accordance with ISO 14708 (for detailed information please see <u>nEw3 Mechanical DAU Design Verification Test Report</u> (section 9.6) During the procedure the product is inside of the human body and is operated at a body temperature from 32 to 42 °C, pressure 50-150 kPa; Humidity up to 100%.
12. Transportation conditions	Packed by the manufacturer products are transported by all modes of transport in closed vehicles in accordance with the rules of carriage of goods, operating in each mode of transport. Do not store and transport this product above + 55 °C. Transport in a manner that protects the integrity of the package and the sterile barrier. Keep in a dry location away from heat and sunlight.
13. Utilization	Explant the implantable device postmortem. In some countries, explanting battery-operated implantable devices is mandatory because of environmental concerns; please check the local regulations. In addition, if subjected to incineration or cremation temperatures, the device may explode. Please return explanted devices to Vitatron for analysis and disposal. Waste class – B (in accordance with Russian legislation)
14. Warranty	

Vitatron Holding B.V.

Mailing address : P.O. Box 1220

6201 MP Maastricht

Visiting Address : Endepolsdomein 5

6229 GW Maastricht

Tel : +31 43 356 6551

Fax : +31 43 356 6542

K.v.K. Zuid Limburg : 090.30.178

V.A.T. : NL 0013.07.915.B.01

Terms of Limited Warranty

- A. This Limited Warranty provides the following assurance to a patient who receives this Vitatron pulse generator:
- (1) Should the pulse generator fail to function within normal tolerances at any time due to a defect in materials or workmanship, Vitatron will issue a credit as defined in Subsection A.3. Battery depletion will occur with time and, as used in Subsection A.2., is not considered a defect in material or workmanship.
 - (2) A credit as defined in Subsection A.3. will be issued if the pulse generator is explanted within 60 months due to battery depletion, in accordance with the pulse generator's Recommended Replacement Time Indicator. Each warranty period commences on the implant date. The warranty on battery depletion is made void if the pulse generator has been programmed to output values higher than the standard factory settings for an accumulated period of longer than 3 months.
 - (3) The credit is equal to the purchase price and is issued against the purchase of any Vitatron pulse generator requested as a replacement, or, at the discretion of Vitatron, is a functionally comparable replacement Vitatron pulse generator at no charge. Normally the credit will be issued only to the hospital or physician purchasing the replacement pulse generator.
- B. To qualify for the Limited Warranty, the following conditions must be met:
- (1) The pulse generator must be implanted prior to its "Use By" date.
 - (2) The pulse generator for which credit is sought must be returned to Vitatron in a shockproof package within 30 days of its removal and shall then become the property of Vitatron. Furthermore, a "Request for analysis" form should be completed and returned with the pulse generator.
- C. This Limited Warranty is limited to its express terms. In particular:
- (1) Vitatron shall not be liable to any person for any medical expenses or any direct or consequential damages resulting from implant or removal of any pulse generator or caused by any defect, failure or malfunction of any pulse generator, whether a claim for such damage is based upon warranty, contract, tort or otherwise. The exclusive remedy for any defect, failure or malfunction of any pulse generator shall be as stated.
 - (2) This Limited Warranty is made only to the patient in whom the pulse generator was originally implanted and is expressly in lieu of any implied warranty of merchantability or fitness for a particular purpose, whether arising from statute, common law, custom or otherwise. The remedies set forth in this Limited Warranty shall be the exclusive remedies available to any person. No person has any authority to bind Vitatron to any representation, condition or warranty except as contained in this Limited Warranty.
 - (3) This Limited Warranty is only valid as to the pulse generator and is not applicable, in whole or part, to any other Vitatron product including, but not limited to, Vitatron pacing leads.
 - (4) Vitatron shall inspect the returned pulse generator and determine whether credit is due under the terms of this warranty. No credit will be issued if Vitatron finds evidence of improper handling or of material alteration to the pulse generator.
 - (5) This Limited Warranty is only applicable to pulse generators sold through a sales organization authorized by Vitatron. Vitatron as used in this warranty shall mean Vitatron B.V. in The Netherlands.
 - (6) The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. Users may benefit from statutory warranty rights under the legislation governing the sales of consumer goods. If any part or term of this Limited Warranty is held by any court of competent jurisdiction to be illegal, unenforceable, or in conflict with applicable law, the validity of the remaining portion of the Limited Warranty shall not be affected, and all rights and obligations shall be construed and enforced as of this Limited Warranty did not contain the particular part or term held to be invalid.

15. Labeling

Vitatron Holding B.V.


Mailing address : P.O. Box 1220
6201 MP Maastricht
Visiting Address : Endepotsdomein 5
6229 GW Maastricht

Tel : +31 43 356 6551
Fax : +31 43 356 6542
K.v.K. Zuid Limburg : 090.30.178
V.A.T. : NL 0013.07.915.B.01

Shelf Box label:

A20 SR
vitatron

SR

IS-1 BI
IS-1 UNI

IPG device

Model A20A1

NBG AAIR / VVIR

Adaptive

VVIR

2.8 mV
Sensitivity

3.5 V $\uparrow$
0.4 ms $\downarrow$
Amplitude / pulse width

AUTO
Pacing polarity

60 min⁻¹
Lower rate

330 ms
Refractory period

85 min⁻¹
Magnet rate

AUTO
Sensing polarity

Package contents 1x

IPG Device

Product documentation

Torque wrench

SN 360005148G

2017-11-28
YYYY-MM-DD

2016-06
YYYY-MM

CE0344

STERILE EO

-18 °C / 0 °F

56 °C / 133 °F

Vitatron Holding B.V.
Endepolsdomein 5
6229 GW Maastricht
The Netherlands

Manufactured in:
Singapore
Product of Singapore

MPA3R / A001 A MP54363A001 A

Sterile package label:

Vitatron Holding B.V.

Mailing address : P.O. Box 1220

6201 MP Maastricht

Visiting Address : Endepolsdomein 5


6229 GW Maastricht

Tel : +31 43 356 6551

Fax : +31 43 356 6542

K.v.K. Zuid Limburg : 090.30.178

V.A.T. : NL 0013.07.915.B.01

Symbol

Explanation

CE0344

Conformité Européenne (European Conformity). This symbol means that the device fully complies with AIMD Directive 90/385/EEC (NB 0344) and R&TTEDirective 1999/5/EC.

Do not use if package is damaged


Do not reuse

Sterilized using ethylene oxide

Consult instructions for use

Date of manufacture

Vitatron Holding B.V.

Mailing address : P.O. Box 1220

6201 MP Maastricht

Visiting Address : Endepolsdomein 5

6229 GW Maastricht

Tel : +31 43 356 6551

Fax : +31 43 356 6542

K.v.K. Zuid Limburg : 090.30.178

V.A.T. : NL 0013.07.915.B.01


Manufacturer

Use by

Reorder number


Serial number

Model number

NBG pacing code that describes device functionality by indicating chamberpaced, chamber sensed, response to sensed event, and rate response capability

Temperature limitation

Adaptive

Package contents

IPG device

Product documentation

Torque wrench

Amplitude / pulse width

Sensitivity

Vitatron Holding B.V.

Mailing address : P.O. Box 1220

6201 MP Maastricht

Visiting Address : Endepolsdomein 5

6229 GW Maastricht

Tel : +31 43 356 6551

Fax : +31 43 356 6542

K.v.K. Zuid Limburg : 090.30.178

V.A.T. : NL 0013.07.915.B.01

Pacing polarity

Sensing polarity

Refractory period

Refractory period

Magnet Rate

Primary Packaging – Sterile Package

Package part	Material
Inner Tray	Coextruded copolyester material consisting of a cap/core/cap construction. The core is virgin polyethylene terephthalate glycol (PETG) and the cap layers are virgin PETG with an antiblocking additive.
Outer Tray	Thermoformed polyethylene terephthalate glycol (PETG)
Inner Lid	Tyvek 1073B lid coated with Oliver 18B adhesive 74.32 cm ² surface area penetrable by EO gas
Outer Lid	Tyvek 1073B lid coated with Oliver 18B adhesive 144.77 cm ² surface area penetrable by EO gas

Sterile package dimensions:

Vitatron Holding B.V.

Mailing address : P.O. Box 1220

6201 MP Maastricht

Visiting Address : Endepotsdomein 5


6229 GW Maastricht

Tel : +31 43 356 6551

Fax : +31 43 356 6542

K.v.K. Zuid Limburg : 090.30.178

V.A.T. : NL 0013.07.915.B.01

Shelf box dimensions:

17,74 cm

4,31 cm

32,24 cm

Secondary Packaging –Shelf Box

Package part	Material
Shelf Box	White vat non-bending chip Box wrap – white litho Adhesive – Davis NW300C or Wisdom 6567, 6578 or equivalent Shrink-wrap of a polyethylene (PE) shrink film

Sticker for final package

Sticker will be added locally in Russia after approval with information below (information can be changed in accordance with Russian legislation):

	Product name (in Russian)	
	Manufacturer (name)	UPN (unique product number)
	Country of origin	CFN n(customer facing number)
	See IFU RC # _____ date of issue _____	

Vitatron Holding B.V.

Mailing address : P.O. Box 1220
6201 MP Maastricht
Visiting Address : Endepoisdomein 5
6229 GW Maastricht

Tel : +31 43 356 6551
Fax : +31 43 356 6542
K.v.K. Zuid Limburg : 090.30.178
V.A.T. : NL 0013.07.915.B.01

Mark of conformity with the declaration of conformity

16. Customer Service

Manufacturer

Vitatron Holding B.V.
Endepolsdomein 5
6229 GW Maastricht
The Netherlands
Tel. +31 43 356 6551

Representative company of Vitatron B.V. in Moscow Russia

123317 Moscow,
Presnenskaya Naberezhnaya , 10
Tel.: +7 495 580 7377; +7 495 580 7377 FREE
Fax: + 7 495 580 7378
E-mail: info.russia@medtronic.ru
Internet: www.medtronic.eu, www.medtronic.ru

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